Chemicals



Gündüz
Chemical Solutions

DASITGROUP | CARLO ERBA REAGENTS

CARLO ERBA Reagents has been a part of the DASIT Group since July 2013 and more than ever I am pleased and honored to lead a multinational Team that combines the virtuosity of its pharmacist founder with more than one century tradition of excellence, and the spirit of service typical of our Group.

CARLO ERBA Reagents is driven by passion, an entrepreneurial spirit, eagerness to innovate and, above all, the desire to consistently connect with and serve our Customers.

The 2020 edition of the CARLO ERBA Reagents catalog is a proof of how the pride of the past can be combined with the foresight for the future.

Special thanks to you, dear Customers, who continue to demonstrate your loyalty and affection in Europe and around the world.

Angelo Fracassi, President of DASIT Group



CONTACT US

DASITGROUP | CARLO ERBA REAGENTS

ITALIA

CARLO ERBA Reagents S.r.l.
Via Raffaele Merendi 22
20010 Cornaredo (MI)

Servizio Clienti
servizioclienticer@dgroup.it
Tel.: +39 02 93 99 190
Fax: +39 02 93 99 10 01

FRANCE

CARLO ERBA Reagents SAS
Chaussée du Vexin,
Parc d'affaire des Portes
27106 Val de Reuil

Service Client
serviceclient@cer.dgroup.it
Tél.: +33 2 32 09 20 00
Fax: +33 2 32 59 11 89

DEUTSCHLAND

CARLO ERBA Reagents GmbH
Denzlinger Str. 27
79312 Emmendingen

Kundendienst
info.de@cer.dgroup.it
Tel.: +49 07641 468 81 90
Fax: +49 07641 468 81 919

ESPAÑA

CARLO ERBA Reagents S.A.
Calle Filadors 35, 6ª
Planta Puerta 5
08208 Sabadell (BCN)

Servicio Cliente
serviciocliente@cer.dgroup.it
Tel.: +34 93 693 37 35
Fax: +34 93 724 31 68

ALL OTHER COUNTRIES

Customer Service
export@cer.dgroup.it
Ph.: +33 2 32 09 20 00
Fax: +33 2 32 59 11 89

www.carloerbareagents.com

- *Online catalogue (search by CAS, product name, product code) with updated specifications;*
- *Updated MSDS in multi-language and certificates of analysis available to download;*
- *Catalogues and brochures available to download;*
- *Contacts for technical information and requests for quotations.*

Ask us for an e-commerce account:

- *For a full management of your online order, with automatic data transfer in our database;*
- *To check our stock availability in real time;*
- *To check your quotations, orders, delivery notes and invoices online;*
- *To view our price list, your specific commercial conditions and all our promotional activities.*

PRODUCTION

CARLO ERBA Reagents ensures for its customers both quality and service thanks to its flexible production plants, its modern quality control laboratories and its efficient logistics organization.

Our 2 production sites in Val de Reuil and Peypin manufacture:

- Solvents, salts, organic and inorganic molecules, mixtures
- From technical to ultrapure grade
- From raw material to pharma excipient grade following the IPEC 2015 guidelines evaluated and authorized by French Ministry of Health (ANSM Agency), decret n°2011/62/UE

To meet your needs, our production tools include:

- Fully automated lines
- Distillation columns
- Purification columns
- Mixer
- Grinder
- Packaging under inert atmosphere
- Storage tanks
- ISO 8 Clean rooms








CUSTOM SOLUTIONS

CARLO ERBA Reagents is the “partner in your choice” for your specific needs.

Customization can be proposed under request. Our expertise allows us to give you the most suitable solution.
Our experienced Quality, Production and Logistics Teams are providing you with:

- Specific purity
- Outsourcing services
- Safety stocks
- Custom mixtures
- Purifications
- Specific QC tests
- Certifications
- Custom packaging
- Shuttle service
- Deliveries in large volumes







CARLO ERBA Reagents products are classified according to their level of purity and suitable for being used in laboratory analysis or industrial preparations.

A color system assigned to each grade helps the customer to find the suitable product for his application.

It is under the customer responsibility to make sure the purchased product is suitable for his use and /or application. CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of inappropriate use.

RS – Specific Grade Reagents

RS

High purity Solvents, Acids and Reagents specifically designed to be used in the following instrumental techniques:

- **Trace Analysis**
  Superpure Acids and Ultrapure Acids for sample mineralization,
  Standards Solutions for AAS, ICP-OES and ICP-MS for instrumental calibration
  Standards Solutions for Pesticides, PAHs, PCB analysis
- **Liquid Chromatography**
  High Purity Solvents suitable for UHPLC-MS, LC-MS, HPLC GOLD UltraGradient Grade, HPLC PLUS Gradient Grade, HPLC Isocratic Grade, HPLC Preparative,
  Ion Chromatography Standard Solutions,
  Ion Pair Chromatography Reagents,
  Silica Gel & Filter Aids available in different mesh sizes.
- **Gas Chromatography**
  High Purity Solvents suitable for GC-HEADSPACE, for GC-MS, for organic trace analysis (ATRASOL®) and for residue analysis (PESTIPUR®)
- **Spectroscopy UV-vis/IR**
  SPECTROSOL® Solvents
- **Organic Synthesis**
  Anhydrous Solvents
  Deuterated Solvents
- **Karl Fischer**
  KF Reagents (ERBAqua®) for volumetric & coulometric titration pyridine-free
- **Food Analysis:**
  Specific Reagents, Acids, Bases and Kjeldahl Catalysts specific for food matrix analysis
- **Histology, Hematology and Cytodiagnosis:**
  Fixatives, Solvents, Embedding Media, Staining Solutions, Dyes, Mounting Media and Immersion Media
- **Conductivity**
  Standard Solutions
- **Electronic Application**
  MOS, RSE and VLSI quality according to semiconductors production processes requirements
- **Physical Properties Standards**
  Density, Melting point, Osmolality, Brix, Viscosity & Colour

RPE – Analytical Grade Reagents

RPE

Standard grade in the laboratory complying with ACS and/or ISO. Many products comply with various pharmacopoeias, including BP, EP and USP for analytical purposes

- **Solvents, Acids, Salts & Reagents**
  A wide range of products guaranteed for analytical purposes
- **Volumetry**
  Titrated Solutions ready-to-use or concentrated in NORMEX vials
- **pHmetry**
  pH Buffer Solutions with and without color, available in different sizes and packaging
- **Indicators**
  Pure Substances and ready to use indicating solutions all provided with their color index

RE – Technical Grade Reagents

RE

Purified commercial products used in many chemicals laboratories and industries

- **Acids, Bases, Solvents, Salts, Organic and Inorganic substances**
  Provided with guaranteed purity for the basic applications in the industry and laboratory
- **Green Solvents**
- **Solvents for Octane Number Determination**

Xcipharm™ – Excipient Grade Products

Xcipharm™

Excipients used in the formulation of pharmaceutical products. Manufactured under IPEC guide 2015

- **Acids, Bases, Solvents, Salts, Organic and Inorganic substances, mixtures and dilutions**
  According to the main Pharmacopoeias: EP, USP, NF, JP, BP, DAB, FU, FP, to be used as Excipients in drug manufacturing.
  A range of more than 100 products with packaging from 1 mL to 30.000 L for liquids and from 100 g to 1 Ton for solids.

ERBApharm® – Pharmaceutical grade Products

ERBApharm®

Raw materials for the pharmaceutical production, complying with the application: EP, USP, NF, JP, BP, DAB, FU, FP.

- **Acids, Bases, Solvents, Salts, Organic and Inorganic substances**
  According to the main Pharmacopoeias: EP, USP, NF, JP, BP, DAB, FU, FP or prepared with raw materials according to them, to be used as raw materials, synthetic intermediates, buffering agents, solvents in the production of active ingredients, and not as excipients.
  A range of more than 500 products with packaging from 1 L to 30.000 L for liquids and from 100 g to 1 Ton for solids.

Xcipharm™ is the new product range dedicated to excipients. We manufacture these products according to the IPEC guidelines, assuring a high level of quality for your excipients.

We can adapt our offer to your needs with packaging sizes starting at 1 L up to bulk quantity.

We provide all of the tests, certifications and statements to make the registration of your products easier to the health authorities.

QUALITY ASSURANCE

- French Ministry of Health (ANSM Agency) registration
- Traceability of raw material and packaging
- Flow chart
- Change control

PROCESS

- Validated cleaning procedures and dedicated equipment
- Labels consolidation
- Certificate of analysis with the name of the producer and the date of production of the raw material
- Stability studies
- Sample library of the raw material (1 year) and of the finished product (shelf life + 1 year)

DOCUMENTATION

- BSE/TSE statement
- OGM statement
- Residual solvents statement
- ICH Q3D
- Risk assessment (2015/C95/02)



Product list (non exhaustive) of excipients. Custom mixtures and dilution are also available

- Acetic acid glacial
- Acetone
- Benzoic acid
- Benzyl alcohol
- Calcium chloride $2H_2O$
- Citric acid
- Dietanolamine
- Diethyl ether
- Diethyl phthalate
- D(+)-glucose
- EDTA
- Ethanol
- Ethyl acetate
- Glycerine
- Glycine
- Hydrochloric acid
- L(+) Lactic acid
- Lactose
- Magnesium chloride $6H_2O$
- Magnesium stearate
- Magnesium sulfate $7H_2O$
- D-Mannitol
- Paraffin oil
- 2-Phenylethanol
- Phosphoric acid
- Potassium chloride
- Potassium hydroxide
- Potassium permanganate
- Potassium Phosphate Monobasic
- Propanol-2
- Propylene glycol
- Ricin oil/Castor oil
- Saccharose
- Sodium acetate trihydrate
- Sodium alginate
- Sodium bicarbonate
- Sodium carbonate
- Sodium chloride
- Sodium citrate dibasic
- Sodium citrate tribasic
- Sodium hydroxide
- Sodium metabisulfite
- Sodium phosphate dibasic
- Sodium sulfite
- Sodium thiosulfate $5H_2O$
- Stearic acid
- Sulfuric acid
- Tannic acid
- Tartaric acid
- Titanium dioxide
- Triethanolamine

Do not hesitate to contact us if your product is not in the list.



iPEC EUROPE

QUALITY

CARLO ERBA Reagents is ISO 9001:2015 certified.

Our quality system is based on risk management in each step, from the selection of raw material to the finished product. Each plant includes a QC lab with experienced chemists controlling raw materials, semi-finished and finished products.



CARLO ERBA Reagents has defined in this catalogue different grades for its products: RE, RPE, RS, ERBApharm®, Xcipharm™, Tailor made. All products have defined documentation depending on their grade in order to answer customers requirements.

Document	RE	RPE	RS	ERBApharm®	Xcipharm™	Tailor made
Certificate of Analysis	Yes	Yes	Yes	Yes	Yes	Yes
Material Safety Data Sheet	Yes	Yes	Yes	Yes	Yes	Yes
Self-questionnaire	Standard	Standard	Standard	ERBApharm®	Xcipharm™	Depending on the quality defined in the Quality Agreement
BSE/TSE	-	-	-	Yes if available	Yes	Depending on the quality defined in the Quality Agreement
Residual Solvents	-	-	-	Yes if available	Yes	Depending on the quality defined in the Quality Agreement
OGM	-	-	-	Yes if available	Yes	Depending on the quality defined in the Quality Agreement
Supply Chain and Risk Assessment (including Flow Chart, ICHQ3D, endotoxins, allergens, latex, melamine, phthalates, glucose and lactose…)	-	-	-	-	Yes	Depending on the quality defined in the Quality Agreement
Non-Disclosure Agreement	-	-	-	-	Yes	Depending on the quality defined in the Quality Agreement


These documents do not exonerate the customer of his responsibility to check that the product supplied by CARLO ERBA Reagents is suitable for his use.



CERTIFICATES OF ANALYSIS

Certificates of Analysis **for CARLO ERBA Reagents products are available at www.carloerbareagents.com by entering the lot number and the product code printed on the label.**

The expiration date for each batch is printed on the label and the certificate of analysis. This information applies to products stored in their original, unopened packaging, away from heat and light, as specified in the safety data sheet.



Carlo Erba Reagents S.A.S.
Parc d'Activités des Portes
Chaussée du Vexin - BP 616
27106 Val de Reuil Cedex
Tél.:02.32.09.20.00
Fax 02.32.09.20.20

Z.I DE VALDONNE - BP 4
13124 PEYPIN
Tél.:04.42.32.41.41
Fax 04.42.72.41.62

Carlo Erba Reagents S.r.l.
Via R. Merendi, 22
20010 Cornaredo (Mi)
Tel. 02 93 991 90 Fax 02 93 991 001

ISO 9001: 2008

Certificate of Analysis

PRODUCT : Sodium chloride ERBApharm-According to pharmacopoeia: Ph.Eur.-FU-Ph.Franc.-BP-DAB-USP-JP

CODE : 368253

LOT N° : V0A589090A ***METHOD*** : 6572

EXPIRING DATE : 2023/01 ***EDITION*** : 6

TEST	U.M.	SPECIFICATION	RESULT
Description	-	White crystalline powder	Conform
Identification	-	Positive	Positive
Appearance of solution	Ph.Eur.	Conform	Conform
Acidity or alkalinity	Ph.Eur.	Conform	Conform
Residue solvents	USP	Conform	Conform
Barium	Ph.Eur.	Conform	Conform
Iodide	Ph.Eur.	Conform	Conform
Ferrocyanide	Ph.Eur.	Conform	Conform
Nitrite	Ph.Eur.	Conform	Conform
Loss on drying	%	<= 0.5	0.04
Mg,alkal.earth met.(Ca)	ppm	<= 100	<100
Bromide	ppm	<= 100	<100
	ppm	<= 25	<25
	ppm	<= 3	<3
	ppm	<= 200	<200
	ppm	<= 0.2	<0.2
	ppm	<= 1	<1
	ppm	<= 2	<2
	ppm	<= 500	<500
	% s.s.	99.0 ÷ 100.5	100.3

ompliance with the current legislation, raw material for pharmaceutical uses included.

Approve Date :01/01/2020

Not signed electronically issued document
QUALITY CONTROL RESPONSIBLE
B. COULANGE (VDR)


1. PRODUCT
2. CODE
3. LOT NUMBER
4. EXPIRATION DATE
5. CODE OF THE METHOD USED FOR THE PRODUCT QUALIFICATION
6. CERTIFICATE EDITION NUMBER WHICH CHANGES WHEN GUARANTEED SPECIFICATIONS ARE UPDATED
7. DESCRIPTION OF THE QUALITY CONTROL TESTS THE PRODUCT UNDERGOES
8. UNIT OF MEASURE
9. GUARANTEED PRODUCT SPECIFICATIONS
10. TEST RESULTS OBTAINED FOR THE SPECIFIC LOT, AND FOR EACH SINGLE TEST, IN COMPLIANCE WITH INTERNATIONAL STANDARDS, IF APPLICABLE
11. APPROVAL DATE FOR TEST RESULTS

GHS is an international system designed to harmonize the classification, labeling and packaging of chemicals, transposed into the Regulation No. 1272/2008 and known as CLP (Classification, Labeling and Packaging). Reg. CLP is, since June 1st 2015, the only legislation enforced in the EU for the classification and labeling of substances and mixtures and is periodically updated through “Adaptations to Technical Progress (ATP)”.

All CARLO ERBA Reagents-branded substances and mixtures are in compliance with the current regulations.



1. *PRODUCT NAME*
2. *PRODUCT CODE*
3. *SIZE*
4. *GRADE*
5. *APPLICATION*
6. *GHS HAZARD SYMBOLS*
7. *GHS REFERENCES OF WARNING AND PRECAUTIONARY*
8. *LOT NUMBER AND EXPIRATION DATE*
9. *INTERNATIONAL IDENTIFICATION*
10. *MAIN TECHNICAL SPECIFICATION*
11. *BARCODE TO READ PRODUCT CODE AND BATCH NUMBER*




SAFETY - DATA SHEET


Based on the new CLP classification system (Reg. 1272/2008), CARLO ERBA Reagents creates a safety data sheet for each product in compliance with Article 31 of the REACH Regulation No. 1907/2006 and its amendments.

CARLO ERBA Reagents safety data sheets are regularly updated and available at:


www.carloerbareagents.com




Download our CLP poster on
www.carloerbareagents.com



HOW TO READ


PRODUCT NAME
TRANSPORT CLASSIFICATION
ACE
Acetone
• Acetone • Acétone • Acetona • Aceton
INTERNATIONAL CLASSIFICATION
CH_3COCH_3
Molecular Weight: 58,01
CAS: 67-64-1
EEC-N: 200-662-2
Classification transport
ONU: 1090
Transport Hazard class: 3
Packing group II
Acetone > RS - For HPLC - Isocratic Grade
Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.790 ÷ 0.792
Refractive index at 20°C 1.3581 ÷ 1.3601
Boiling point ..... 55.8 ÷ 56.3
Water (K.F.) ..... ≤ 500
Residue on evaporation ..... ≤5
Acidity ..... ≤0.0005 me
APPLICATION
Code | Size | Packaging
412501 | 1 l | Glass bottle
412502 | 2.5 l | Glass bottle
Acetone > RS - For GC-MS
Appearance ..... Clear colourless liquid
Refractive index at 20°C ..... 1.357 - 1.361
Water (K.F.) ..... ≤ 500 ppm
Colour ..... ≤ 10 APHA
Residue on evaporation ..... ≤ 2 pp
Acidity (acetic acid) ..... ≤ 20 ppr
Assay (GC) ..... ≥ 99.95
Ethyl alcohol ..... ≤ 100 ppr
PRODUCT CODE
Code | Size | Packaging
400952 | 1 l | Glass bottle
Acetone > RS - ATRASOL - For traces analysis
Appearance ..... Clear colourless liquid
Free acid (as CH3COOH) ..... ≤ 20
UNIT SIZE
PACKAGING



Chemical specifications are subject to change.

Please see the updated specifications on www.carloerbareagents.com

PREMIUM CHOICE



Protection, Safety, Reliability
for any Laboratory

- Microbiological & Cytotoxic Safety Cabinets (Class I-II-III)
- Laminar Airflow Cabinets
- Pharma Industry Laminar Airflow systems
- Customized solutions for special applications
- Chemical Fume Cupboards
- FORLAB *4you* modular, safe, ergonomic, eco friendly furniture
- Ultra-freezers

Manufactured by DASITGROUP FASTER *in Ferrara (IT) in a brand new 10.000 sqm plant.*

www.carloerbareagents.com

APPLICATIONS

LIQUID CHROMATOGRAPHY

GAS CHROMATOGRAPHY

ION PAIR CHROMATOGRAPHY

ION CHROMATOGRAPHY

TRACE ANALYSIS: METALS

TRACE ANALYSIS: ORGANICS

COD ANALYSIS

PHARMACEUTICAL PRODUCTION

PHARMACEUTICAL QUALITY CONTROL

PHYSICALS AND CHEMICAL CHARACTERISTICS

pHMETRY

VOLUMETRY

FOOD ANALYSIS

ELECTRONICS

HISTOLOGY, HEMATOLOGY AND CYTODIAGNOSTIC

ORGANIC CHEMISTRY

NMR SPECTROSCOPY

UV SPECTROSCOPY

GREEN CHEMISTRY

PETROCHEMICAL ANALYSIS

KARL FISCHER TITRATION

INDICATORS

CONDUCTIMETRY

BEST CHOICE

Tailor made

Chemicals such as you imagine

CARLO ERBA Reagents thanks to its experience and production flexibility manufacture tailor made products for industry and laboratories.



- Custom mixture
- Purification
- Custom packaging
- Deliveries in bulk tanks and isotanks
- Shuttle service
- Analytical control
- Batch management
- Quality assurance

www.carloerbareagents.com

LIQUID CHROMATOGRAPHY

Leader in the market of solvents for chromatography and trace analysis, CARLO ERBA Reagents extended its range of solvents for HPLC in order to satisfy the ever increasing demand requirements in terms of equipment and detection methods, particularly for impurities which, by interaction, can affect the result's reliability.
Our solvents for HPLC meet perfectly the requirements of this analytical technique by guaranteeing the optimal specifications on the following elements: purity, non-volatile residue content and UV transmission.

UHPLC-MS Solvents

The UHPLC-MS is certainly the chromatographic technique for users who, besides being on the lookout for the best analytical performances, works at very high pressure, with minimum solvent consumption and guarantee of resolution and reproducibility of results. CARLO ERBA Reagents, recognized as a leading company in the manufacture of solvents for chromatography, is always attuned to the market needs and has developed a specific range of solvents dedicated to solvents for UHPLC-MS in order to meet the quality requirements of this refined analytical technique.
Our solvents are characterized by:
- Purity higher than 99.95 %
- High UV transmittance
- Excellent baseline quality in gradient, tested specifically for UHPLC
- Test with reserpine (30 < ppb), specific for LCMS
- Low content in inorganic and metallic ions
- Non volatile residue content less than 1 ppm
- Packaged in amber glass bottles pretreated with 1.1-difluoroethane to significantly reduce the potential formation of metals adduct

Description	Notes	Size	Code	Page
Acetonitrile		1 l	412041	142
Acetonitrile		2.5 l	412042	142
Methanol		1 l	414941	488
Methanol		2.5 l	414942	488
Water		1 l	412091	786
Water		2.5 l	412092	786

LC-MS Solvents

Liquid chromatography coupled with mass spectrometry is a recently developed technique which is suitable for the analysis of fairly polar, non-volatile and thermally stable compounds. It allows the separation and characterization of many more compounds than GC-MS chromatography and also provides information on the molecular weight and structure of the HPLC peaks. The significant advantages of this combined technique has led to its widespread use in many analytical fields, particularly in the pharmaceutical, environmental and biotechnological sectors.
CARLO ERBA Reagents guarantees specifications of a high-performance LC-MS solvents:
- High purity, low acidity, alkalinity and residue, ideal fluorescence, absorbance/transmittance and gradient test
- Low metal content, in the order of ppb, in order to prevent interactions with ionized species in the mass analyzer
- LC-MS suitability: no signal is greater than the molecular peak of reserpine (609 amu) at the concentration of 30 ppb, in a range from 50 to 2000 amu
- Packaged in amber glass bottles pretreated with 1.1-difluoroethane to reduce significantly the potential formation of metals adduct

Description	Notes	Size	Code	Page
Acetonitrile		1 l	412341	142
Acetonitrile		2.5 l	412342	142
Ethyl acetate		1 l	448383	340
Ethyl acetate		2.5 l	448384	340
Methanol		1 l	414831	488
Methanol		2.5 l	414832	488
Propan-2-ol		1 l	415183	623
Propan-2-ol		2.5 l	415184	623
Water		1 l	412111	786
Water		2.5 l	412112	786

LC-MS Acids, Salts and Blends

CARLO ERBA Reagents proposes a wide range of acids and salts, specifically tested for LC-MS coupling, in addition to ready-to-use blends solutions for mobile phases.

Description	Notes	Size	Code	Page
Acetic acid glacial		10 x 1 ml	401411	130
Acetic acid glacial		10 x 2.5 ml	401412	130
Acetic acid glacial		50 ml	401413	130
Acetic acid glacial		1 l	401414	130
Acetonitrile + 0.1% v/v formic acid		1 l	412331	145
Acetonitrile + 0.1% v/v formic acid		2.5 l	412332	145
Acetonitrile + 0.1% v/v trifluoroacetic acid		1 l	412321	145
Acetonitrile + 0.1% v/v trifluoroacetic acid		2.5 l	412322	145
Ammonium acetate		50 g	418781	166
Ammonium formate		50 g	419741	171
Formic acid 99%		10 x 1 ml	405821	362
Formic acid 99%		10 x 2.5 ml	405822	362
Formic acid 99%		50 ml	405823	362
Formic acid 99%		1 l	405824	362
Methanol + 0.1% v/v formic acid		1 l	414861	492
Methanol + 0.1% v/v formic acid		2.5 l	414862	492
Methanol + 0.1% v/v trifluoroacetic acid		1 l	414871	493
Methanol + 0.1% v/v trifluoroacetic acid		2.5 l	414872	493
Trifluoroacetic acid		10 x 1 ml	411541	773
Trifluoroacetic acid		10 x 2.5 ml	411542	773
Trifluoroacetic acid		50 ml	411543	773
Water + 0.1% v/v formic acid		1 l	412121	789
Water + 0.1% v/v formic acid		2.5 l	412122	789
Water + 0.1% v/v trifluoroacetic acid		1 l	412031	789

HPLC Ultragradient Grade Solvents

The gradient control of elution and drift at critical wavelengths of our HPLC solvents, Gold UltraGradient and Plus Gradient Grade, guarantee a peak free baseline. Their optimal sensitivity allows you to evaluate in the best possible way the impurities of your samples. To make sure that no particles in the mobile phase will hinder your analyses, we carry out a microfiltration of our GOLD solvents.

Description	Notes	Size	Code	Page
Acetonitrile		1 l	412371000	142
Acetonitrile		2.5 l	412372000	142
Acetonitrile		4 l	412374	142
Methanol		1 l	412721	488
Methanol		2.5 l	412722	488
Methanol		4 l	412724	488
Methanol		5 l	412725	488

HPLC Gradient Grade Solvents

The solvents of this product line guarantee excellent short-wavelength performance and limited drift, which makes them ideal for gradient and trace analysis.
To make sure that no particle in the mobile phase will hindler your analyses, we carry out a microfiltration of our HPLC gradient plus.

Description	Notes	Size	Code	Page
Acetonitrile		1 l	412393	142
Acetonitrile		1 l	412391000	142

Product specifications are subject to changes. Please visit our website for updates.

Acetonitrile		2.5 l	412392000	142
Acetonitrile		5 l	412395	142
Ethanol absolute anhydrous		500 ml	412704	331
Ethanol absolute anhydrous	Only for Italian market	1 l	412701	331
Ethanol absolute anhydrous	Only for Italian market	1 l	412703	331
Ethanol absolute anhydrous		1 l	4127012	331
Ethanol absolute anhydrous		1 l	4127032	331
Ethanol absolute anhydrous	Only for Italian market	2.5 l	412702	331
Ethanol absolute anhydrous		2.5 l	4127022	331
Methanol		1 l	412381	488
Methanol		2.5 l	412383	488
Propan-2-ol		1 l	412711000	623
Propan-2-ol		2.5 l	412712000	623
Water		1 l	412141	786
Water		2.5 l	412142	786

HPLC Isocratic Grade Solvents

These products, due to their high purity and strictly controlled chemical-physical parameters, adequately meet the needs of modern analytical HPLC.
Available in glass bottles (1l and 2.5l) or stainless steel shuttle drums (5 to 1000l), their characteristics satisfy the requirements of the most advanced HPLC techniques.

Description	Notes	Size	Code	Page
Acetic acid glacial		1 l	401431	130
Acetic acid glacial		2.5 l	401432	130
Acetone		1 l	412501	138
Acetone		2.5 l	412502	138
Acetonitrile		1 l	412411000	143
Acetonitrile		2.5 l	412412000	143
Acetonitrile		4 l	412413000	143
Butanol-1		1 l	412511000	226
Butanol-1		2.5 l	412512000	226
n-Butyl chloride		1 l	431821	229
tert-Butylmethylether		1 l	432031	230
tert-Butylmethylether		2.5 l	432032	230
tert-Butylmethylether		4 l	432034	230
Chloroform, stab. with Amylene		1 l	412571	259
Chloroform, stab. with Amylene		2.5 l	412572	259
Chloroform, stab. with Ethanol		1 l	412652	259
Chloroform, stab. with Ethanol		2.5 l	412653	259
Cyclohexane		1 l	412431000	286
Cyclohexane		2.5 l	412432000	286
1,2-Dichloroethane		1 l	447191	299
1,2-Dichloroethane		2.5 l	447192	299
Dichloromethane, stab. with Amylene		1 l	412621000	300
Dichloromethane, stab. with Amylene		2.5 l	412622000	300
Dichloromethane, stab. with Ethanol		1 l	412662	300
Dichloromethane, stab. with Ethanol		2.5 l	412661	300
Diethyl ether		1 l	412671	308
Diethyl ether		2.5 l	412672	308
N,N-Dimethylformamide		1 l	444981	314

N,N-Dimethylformamide		2.5 l	444982	314
Dimethylsulphoxide		1 l	445141	317
Dimethylsulphoxide		2.5 l	445142	317
1,4-Dioxane		1 l	443231	321
Ethanol absolute anhydrous	Only for Italian market	1 l	412521	331
Ethanol absolute anhydrous		1 l	4125212	331
Ethanol absolute anhydrous	Only for Italian market	2.5 l	412522	331
Ethanol absolute anhydrous		2.5 l	4125222	331
Ethanol 96°	Only for Italian market	1 l	414541	334
Ethanol 96°		1 l	4145412	334
Ethanol 96°	Only for Italian market	2.5 l	414542	334
Ethanol 96°		2.5 l	4145422	334
Ethyl acetate		1 l	412611000	340
Ethyl acetate		2.5 l	412612000	340
n-Heptane 99%		1 l	412591000	380
n-Heptane 99%		2.5 l	412592000	380
n-Heptane		1 l	446831	382
n-Heptane		2.5 l	446832	382
n-Hexane 99%		1 l	412691	385
n-Hexane 99%		2.5 l	412692	385
n-Hexane		1 l	412601000	386
n-Hexane		2.5 l	412602000	386
Hexane mixture of isomers		1 l	412632	388
Hexane mixture of isomers		2.5 l	412631	388
Isohexane		1 l	445152	438
Isohexane		2.5 l	445151	438
Isooctane		1 l	412441000	439
Isooctane		2.5 l	412442000	439
Methanol		1 l	412531	489
Methanol		1 l	412533	489
Methanol		2.5 l	412532	489
Methanol		2.5 l	412535	489
Methyl acetate		2.5 l	P0043721	495
2-Methyltetrahydrofuran		1 l	412681	505
2-Methyltetrahydrofuran		2.5 l	412682	505
n-Pentane		1 l	P0643716	558
n-Pentane		2.5 l	P0643721	558
Propan-1-ol		1 l	412541000	622
Propan-1-ol		2.5 l	412542000	622
Propan-2-ol		1 l	412821	623
Propan-2-ol		1 l	412421000	624
Propan-2-ol		2.5 l	412422000	624
Tetrahydrofuran		1 l	412451000	751
Tetrahydrofuran		1 l	412453000	751
Tetrahydrofuran		2.5 l	412452000	751
Tetrahydrofuran, stab. with BHT		1 l	412471	751
Tetrahydrofuran, stab. with BHT		2.5 l	412472	751
Toluene		1 l	412641000	765
Toluene		2.5 l	412642000	765
Triethylamine		1 l	489631	772

Triethylamine		2.5 l	489633	772

HPLC Preparative Solvents

Our range of solvents for HPLC preparative have been designed to satisfy the requirements of separations and purifications. Their low non-volatile residue content (from 5 to 10 ppm maximum) allows to optimize the operation conditions and to make impurity-free preparations. All these solvents are available in 2.5 l bottles and in stainless steel shuttle drums from 5 to 1 000 l.

Description	Notes	Size	Code	Page
Acetonitrile		2.5 l	412409	143
tert-Butylmethylether		2.5 l	432022000	230
Chloroform		2.5 l	438641	260
Dichloromethane		2.5 l	463281	300
Dichloromethane		2.5 l	463291	300
Ethyl acetate		2.5 l	448211	341
Propan-2-ol		2.5 l	415112	624
Tetrahydrofuran		2.5 l	487352	752
Toluene		2.5 l	488531	766

Silica Gel and Filter aids

Besides the widely used silica gel, other products with particular characteristics are also available and offer a series of valid alternatives for resolving numerous separation problems.

Description	Notes	Size	Code	Page
Aluminum oxide (acid)		250 g	417185	155
Aluminum oxide (acid)		1 kg	417182	155
Aluminum oxide (basic)		100 g	417214	155
Aluminum oxide (basic)		1 kg	417217	155
Aluminum oxide (neutral)		250 g	417245	156
Aluminum oxide (neutral)		1 kg	417241	156
Aluminum oxide activated		1 kg	312261	156
Calcium carbonate		250 g	433245	237
Cellulose, powder		250 g	436061	250
Charcoal activated		250 g	434455	255
Charcoal activated		1 kg	434454	255
Charcoal activated	Granular	1 kg	P4610017	255
Dicalite 4158		500 g	P8880014	297
Dicalite 4158		1 kg	P8880017	297
Dicalite 4158		5 kg	P8880027	297
Florisil 60-100 mesh		100 g	452331	356
Florisil 60-100 mesh		500 g	452333	356
Florisil 60-100 mesh		1 kg	452332	356
Florisil 60-100 mesh		100 g	452271	356
Florisil 60-100 mesh		500 g	452273	356
Florisil 100-200 mesh		100 g	452351	357
Florisil 100-200 mesh		500 g	452353	357
Kieselguhr composed		250 g	449895	449
Kieselguhr composed		1 kg	449897	449
Magnesium oxide		1 kg	459617	473
Sand purified		1 kg	477153	644
Silica gel 60A 6 - 35µ		1 kg	P2010017	648
Silica gel 60A 6 - 35µ		5 kg	P2010027	648
Silica gel 60A 6 - 35µ		25 kg	P2010044	648

Silica gel 60A 20 - 45µ	1 kg	P2200017	648
Silica gel 60A 20 - 45µ	5 kg	P2200027	648
Silica gel 60A 35 - 70µ	1 kg	P2000017	649
Silica gel 60A 35 - 70µ	2 kg	P2000026	649
Silica gel 60A 35 - 70µ	5 kg	P2000027	649
Silica gel 60A 35 - 70µ	25 kg	P2000044	649
Silica gel 60A 35 - 70µ	100 g	453351	649
Silica gel 60A 35 - 70µ	500 g	453352	649
Silica gel 60A 35 - 70µ	1 kg	453353	649
Silica gel 60A 35 - 70µ	5 kg	453355	649
Silica gel 60A 40 - 63µ	1 kg	P2050017	649
Silica gel 60A 40 - 63µ	5 kg	P2050027	649
Silica gel 60A 40 - 63µ	25 kg	P2050044	649
Silica gel 60A 70 - 200µ	1 kg	P2100017	650
Silica gel 60A 70 - 200µ	2 kg	P2100026	650
Silica gel 60A 70 - 200µ	5 kg	P2100027	650
Silica gel 60A 70 - 200µ	25 kg	P2100044	650
Silica gel 60A 0,06÷0,20 mm	500 g	453336	650
Silica gel 60A 0,06÷0,20 mm	1 kg	453337	650
Silica gel 60A 0,06÷0,20 mm	5 kg	453332	650
Silica gel 60A 0,06÷0,20 mm	20 kg	453331	650

GAS CHROMATOGRAPHY

Broad spectrum chemical analysis of trace level components is a continuing challenge for any analytical chemist. This challenge is further confounded when chemical impurities may be present in common organic solvents or when chemical artifacts may be formed, produced and introduced during an analytical procedure. Minimizing and understanding these chemical artifacts is critical for trace level detection and is crucial for accurate analytical conclusions.
CARLO ERBA Reagents GC Solvents are the best choice for your complex mixture challenges.

GC-MS Solvents

The recent technological advances of GC-MS, GC-MS/MS and 2D GC-MS have opened new analytical horizons, in terms of selectivity of the result, and allowed a reduction of detection limits, reducing the need for cleaning the sample and the introduction of faster methods for sample preparation.
The role and the choice of the quality of the solvent is consequently crucial for the production of a precise and accurate analytical data.
That is why we are introducing a new product range dedicated to the most demanding need for GC-MS. These products were specifically tested for GC/MS test for individual signals, with a retention range of C11 to C40 with a scanning area of 30-600 amu with a guarantee of less than 2µg/l of impurities.
The CARLO ERBA Reagents GC-MS solvents guarantee excellent performance, even for the analysis of the most complex mixtures.

Description	Notes	Size	Code	Page
Acetone		1 l	400952	138
Chloroform		1 l	438732	260
Dichloromethane		1 l	463332	301
Dichloromethane		1 l	463342	301
Ethyl acetate		1 l	448342	341
n-Hexane 99%		1 l	447212	385
Methanol		1 l	414952	489
n-Pentane 99%		1 l	468172	557
n-Pentane		1 l	468182	558

Head Space Solvents

The operating principle of this technique is based on the chromatographic analysis of the vapor phase in thermodynamic equilibrium above the sample enclosed in a sealed container.
Analysis of residue solvents using GC Heaspace techniques, has become a major control procedure in pharmaceutical and food related industries. CARLO ERBA Reagents has recently developed solvents specifically tested for GC-HS applications. Their purity and handling specifications meet the requirements of the latest Pharmacopoeia guidelines for the residual solvent content in pharmaceutical products.

Description	Notes	Size	Code	Page
N,N-Dimethylacetamide		1 l	444311	313
N,N-Dimethylformamide		1 l	444991	315
Dimethylsulphoxide		1 l	445121	318
N-Methyl-2-pyrrolidone		1 l	462881	502
Water		1 l	412011	786

PESTIPUR® Solvents for pesticides residue analysis

The control of pesticide residues in the food and environmental sectors is remarkably important today, as these substances represent a potential public health hazard.
The purity of the solvent is a determinant factor in obtaining reliable results. Thus it is essential to have products available with suitable parameters for this type of application.
To meet these needs, CARLO ERBA Reagents offers its PESTIPUR® line of solvents, specific for the extraction of pesticides and the analysis of chlorinated and nitrogenous residues, even at trace levels.
Our products are prepared according to the most advanced distillation techniques and strictly controlled in order to guarantee the highest level of quality.
Various functionality tests ensure a stable base line in gas chromatography.
For the entire PESTIPUR® line, the absence of critical impurities is ensured by means of precise functionality tests in GC-ECD and GC-NPD.

Description	Notes	Size	Code	Page
Acetone		1 l	400991	138
Acetone		2.5 l	400992000	138
Acetone		4 l	400994	138

Acetone	2.5 l	400932	138
Acetonitrile	1 l	401241	143
Acetonitrile	2.5 l	401242	143
Acetonitrile	4 l	401243	143
tert-Butylmethylether	1 l	432061	230
tert-Butylmethylether	2.5 l	432062	230
Chloroform, stab. with Amylene	1 l	438681	260
Chloroform, stab. with Amylene	2.5 l	438682	260
Chloroform, stab. with Ethanol	1 l	438651	260
Chloroform, stab. with Ethanol	2.5 l	438652	260
Cyclohexane	1 l	436931	286
Cyclohexane	2.5 l	436932	286
Dichloromethane, stab. with Amylene	1 l	442291	301
Dichloromethane, stab. with Amylene	2.5 l	442292000	301
Dichloromethane, stab. with Amylene	4 l	442294	301
Dichloromethane, stab. with Ethanol	1 l	442261	302
Dichloromethane, stab. with Ethanol	2.5 l	442262	302
Diethyl ether	1 l	447651	308
Diethyl ether	2.5 l	447652	308
N,N-Dimethylformamide	1 l	444941	315
N,N-Dimethylformamide	2.5 l	444942	315
Ethyl acetate	1 l	448351	341
Ethyl acetate	2.5 l	448352000	341
n-Heptane 99%	1 l	446951	380
n-Heptane 99%	2.5 l	446952	380
Heptane mixture of isomers	1 l	446841	382
Heptane mixture of isomers	2.5 l	446842	382
n-Hexane 99%	1 l	447111	385
n-Hexane 99%	2.5 l	447112000	385
n-Hexane	1 l	447011	386
n-Hexane	2.5 l	447012	386
n-Hexane	4 l	447013	386
Hexane mixture of isomers	1 l	447181	388
Hexane mixture of isomers	2.5 l	447182	388
Isohexane	1 l	447131	439
Isohexane	2.5 l	447132	439
Isooctane	1 l	456791	439
Isooctane	2.5 l	456792	439
Methanol	1 l	414930	489
Methanol	2.5 l	414932	489
n-Pentane	1 l	468161	558
n-Pentane	2.5 l	468162	558
Petroleum ether 40 - 65°C	1 l	447851	566
Petroleum ether 40 - 65°C	2.5 l	447852	566
Petroleum ether 35 - 60°C	1 l	447862	567
Petroleum ether 35 - 60°C	2.5 l	447861	567
Propan-2-ol	1 l	415281	624
Toluene	1 l	488591	766
Toluene	2.5 l	488592	766
Toluene	4 l	488594	766

Product specifications are subject to changes. Please visit our website for updates.

ATRASOL® Solvents for the detection of traces in organic compounds and hydrocarbons

Specific solvents for gas chromatographic analysis of trace pollutants.
High purity, guaranteed absence of extraneous peaks in gas chromatographic determinations and guarantee of reproducibility and repeatability of the result are the main feature of this line.
Furthermore, for all the ATRASOL® solvents, the absence of critical impurities is ensured by means of precise functionality tests in GC-ECD and GC-FID.

Description	Notes	Size	Code	Page
Acetone		1 l	P0053216	138
Acetone		2.5 l	P0053221	138
Acetone		4 l	P0053282	138
Chloroform		1 l	P02432E16	260
Chloroform		2.5 l	P02432E21	260
Dichloromethane, stab. with Amylene		1 l	P02932A16	301
Dichloromethane, stab. with Amylene		2.5 l	P02932A21	301
Dichloromethane, stab. with Amylene		4 l	P02932A82	301
Dichloromethane, stab. with Ethanol		1 l	P02932E16	301
Dichloromethane, stab. with Ethanol		2.5 l	P02932E21	301
N,N-Dimethylformamide		1 l	P0343216	315
N,N-Dimethylformamide		2.5 l	P0343221	315
Dimethylsulphoxide		1 l	P0353216	318
Dimethylsulphoxide		2.5 l	P0353221	318
Ethyl acetate		1 l	P0023216	341
Ethyl acetate		2.5 l	P0023221	341
n-Hexane 99%		1 l	P052323016	385
n-Hexane 99%		2.5 l	P052323021	385
Methanol		1 l	P0933216	489
Methanol		2.5 l	P0933221	489
n-Pentane 99%		1 l	P064323016	557
n-Pentane 99%		2.5 l	P064323021	557
Toluene		1 l	P0713216	766
Toluene		2.5 l	P0713221	766
Toluene		4 l	P0713282	766

ATRASOL® Solvents for Hydrocarbon index determination according to EN ISO 9377-2

The European regulation UNI ISO 9377-2 "Determination of hydrocarbon oil index - Method using solvent extraction and gas chromatography", established the criteria for the evaluation of the hydrocarbon index in water using gas chromatography. This procedure is suitable for surface water, wastewater and water from sewage treatment plants.
CARLO ERBA Reagents offer suitable extraction solvents, with their boiling range between 56 and 69°C.
Each production lot is specifically analyzed so that the hydrocarbon index is less than or equal to 0.1 mg/l, in the retention time window between n-decane and n-tetracontane.

Description	Notes	Size	Code	Page
n-Hexane		1 l	P0523216	386
n-Hexane		2.5 l	P0523221	386
Isohexane		1 l	P6263216	438
Isohexane		2.5 l	P6263221	438
n-Pentane		1 l	P0643216	558
n-Pentane		2.5 l	P0643221	558
Petroleum ether 35 - 60°C		1 l	P0883216	567
Petroleum ether 35 - 60°C		2.5 l	P0883221	567

ION PAIR CHROMATOGRAPHY

Ion Pair Chromatography was developed to allow the separation of complex mixtures of polar and ionic molecules, which often are not well separated by ion exchange chromatography. The selectivity is determined by the mobile phase: the organic eluent is supplemented with a specific ion-pairing reagent. The IPC reagents are large ionic molecules having a charge opposite to the targeted analyte, as well as a hydrophobic region to interact with the stationary phase. The counter-ion combines with the ions of the eluent, becoming ion pairs in the stationary phase. Ion pairs are then separated on Reverse-phase HPLC columns.

Reagents for Ion Pair Chromatography

Derivatives can be used to solve common analytical problems related to ionic or polarized products in chromatography. The following reagents are additives for the mobile phase that allow the separation in reversed-phase HPLC of ionic or highly polar substances (counter-ion tetraalkylammonium for anionic electrolytes, alkyl or aryl sulfonate for cationic electrolytes).
Benefits of CARLO ERBA ion-pair reagents are:

- The purity of the mobile phase and therefore the accuracy of the results depend on the quality of the additive
- The specifications of our ion pair reagents are in line with the requirements of Reverse-phase HPLC: high purity ≥ 99%, minimum UV absorption in the far UV, controlled pH and minimum loss on drying

Description	Notes	Size	Code	Page
1-Butanesulfonic acid sodium salt		25 g	405631	225
1-Butanesulfonic acid sodium salt		100 g	405632	225
1-Decanesulfonic acid sodium salt		25 g	405871	290
1-Decanesulfonic acid sodium salt		100 g	405872	290
1-Dodecanesulfonic acid sodium salt		25 g	405881	324
1-Dodecanesulfonic acid sodium salt		100 g	405882	324
Dodecyltrimethylammonium bromide		25 g	405941	325
1-Heptanesulphonic acid sodium salt		25 g	405851	383
1-Heptanesulphonic acid sodium salt		100 g	405852	383
1-Hexanesulphonic acid sodium salt		25 g	405621	389
1-Hexanesulphonic acid sodium salt		100 g	405622	389
1-Hexanesulphonic acid sodium salt monohydrate		25 g	405921	389
1-Hexanesulphonic acid sodium salt monohydrate		100 g	405922	389
1-Octanesulphonic acid sodium salt		25 g	405861	537
1-Octanesulphonic acid sodium salt		100 g	405862	537
1-Octanesulphonic acid sodium salt		1 kg	405863	537
1-Octanesulfonic acid sodium salt monohydrate		25 g	405931	538
1-Octanesulfonic acid sodium salt monohydrate		100 g	405932	538
1-Pentanesulphonic acid sodium salt		25 g	405841	559
1-Pentanesulphonic acid sodium salt		100 g	405842	559
1-Pentanesulphonic acid sodium salt monohydrate		25 g	405891	560
1-Pentanesulphonic acid sodium salt monohydrate		100 g	405892	560
1-Propanesulfonic acid sodium salt		25 g	405901	627
1-Propanesulfonic acid sodium salt		100 g	405902	627
Tetrabutylammonium bisulfate		25 g	405971	748
Tetrabutylammonium bisulfate		100 g	405972	748

Product specifications are subject to changes. Please visit our website for updates.

ION CHROMATOGRAPHY

Ion chromatography is a widely used technique that separates ions and polar molecules based on their affinity to the ion exchanger. It is often used in protein purification and water analysis. it works on almost any kind of charged molecule—including large proteins, small nucleotides, and amino acids.

Standard Solutions for Ion Chromatography

Our standard solutions for ion chromatography are obtained by dissolution of a high-purity salt (+99.9%) in water.
They are characterized by:
- Produced under a Quality Management System that is Accredited according to ISO/IEC 17025 and ISO 17034;
- Concentrations equal to 1000 ppm;
- Guaranteed titer with its uncertainty;
- Raw materials selected and verified against N.I.S.T. Standard Reference Materials;
- Available in HDPE bottles;
- Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval;

Description	Notes	Size	Code	Page
Ammonium standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503311	165
Ammonium standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503313	165
Bromate standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503171	210
Bromate standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503173	210
Bromide standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503211	210
Bromide standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503213	210
Calcium standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503221	235
Calcium standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503223	235
Calcium standard solution	conc. 1.000 ppm Matrix: Water and nitric acid	500 ml	503389	235
Chlorate standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503181	257
Chlorate standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503183	257
Chloride standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503231	257
Chloride standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503233	257
Chlorite standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503191	257
Chlorite standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503193	257
Chromate standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503241	264
Chromate standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503243	264
Cyanide standard solution	conc. 1.000 ppm Matrix: Water and nitric acid	100 ml	503358	285
Fluoride standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503251	358
Fluoride standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503253	358
Iodide standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503261	424
Iodide standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503263	424
Lithium standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503281	462
Lithium standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503283	462
Magnesium standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503291	469
Magnesium standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503293	469
Magnesium standard solution	conc. 1.000 ppm Matrix: Water and nitric acid	500 ml	503390	469
Nitrate standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503331	527
Nitrate standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503333	527
Nitrite standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503321	534
Nitrite standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503323	534
Phosphate standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503341	574
Phosphate standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503343	574
Potassium standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503271	583
Potassium standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503273	583
Sodium standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503301	661
Sodium standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503303	661

Strontium standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503361	723
Sulfate standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503351	729
Sulfate standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503353	729

Multistandard Solutions for Ion Chromatography

For instrument calibration, the following multi-ion standard solutions are available for ion chromatography, complete with certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval. CARLO ERBA Reagents offers to customers the possibility of requesting quotes and ordering custom-made multi-ion solutions.
Just send us your request, specifying the ions of interest, their respective concentrations and the volume requested.

Description	Notes	Size	Code	Page
Multianions standard for ion chromatography	7 elements: Br-, Cl-, NO3-, NO2-, PO43-, SO42-, F- 1g/l each - Matrix: Water	100 ml	504526	512
Multianions standard for ion chromatography	7 elements: Br- 100ppm, SO42- 150ppm, PO43- 50ppm, Cl- 30ppm, NO2- 30ppm, NO3- 20ppm, F- 20 ppm - Matrix: Water	100 ml	504527	512
Multianions standard for ion chromatography	7 elements: F- 20mg/l; Cl- 100mg/l; NO2- 100mg/l; Br- 100mg/l; NO3- 100mg/l; PO43- 200mg/l; SO42- 100mg/l; Matrix: Water	500 ml	504677	512

Concentrated mobile phases for Ion Chromatography

The following eluents are filtered at 0.2µm and prepared from ultra-pure salts and 18-megaohm deionized water.
These are concentrated solutions that should be diluted by a factor of 100.

Description	Notes	Size	Code	Page
Eluent sodium bicarbonate	0.5 M Sodium bicarbonate	100 ml	504534	327
Eluent sodium bicarbonate	0.5 M Sodium bicarbonate	1 l	507578	327
Eluent sodium carbonate	0.5 M Sodium carbonate	100 ml	504533	327
Eluent sodium carbonate	0.5 M Sodium carbonate	1 l	507577	327
Eluent sodium carbonate/sodium bicarbonate	0.18 M Sodium carbonate / 0.17 M Sodium bicarbonate	100 ml	504530	328
Eluent sodium carbonate/sodium bicarbonate	0.22 M Sodium carbonate / 0.28 M Sodium bicarbonate	100 ml	504531	328
Eluent sodium carbonate/sodium bicarbonate	0.35 M Sodium carbonate / 0.1 M Sodium bicarbonate	100 ml	504532	328

Product specifications are subject to changes. Please visit our website for updates.

TRACE ANALYSIS: METALS

In choosing the most appropriate analytical method to determine metals, each laboratory must consider the sample type and concentration levels, the number of elements to be determined and the costs the choice implies.
As a result, flame and graphite furnace atomic absorption spectrophotometry (AA) and inductively coupled plasma (ICP and ICP-MS) emission spectrometry are the most widely used analytical methods for determining trace elements.
Instrumental analysis, using ICP or AA, generally involves a preliminary treatment of the sample. This operation, known as acid mineralization, consists in a digestion process with hot concentrated acid in order to extract the elements of interest. CARLO ERBA Reagents offers two specific complete range of products (acids, bases and water) for sample and blank preparation. The purity of these products guarantees maximum reliability of the result.

Superpure Acids for trace metal analysis at ppb level

SUPERPURE range is characterized by blank values generally between 0.5 and 1 ppb, for the 60-plus declared impurities.
They are produced using the most advanced sub-boiling distillation techniques, in special equipment made of quartz or Teflon and packaged in a controlled environment.
In order to minimize the possibility for contamination of the resultant distillate, the packaging is performed in a clean room. They are available in a wide variety of molecules and sizes.

Description	Notes	Size	Code	Page
Acetic acid glacial		500 ml	401405	131
Acetic acid glacial		1 l	401406	131
Acetic acid glacial		2.5 l	401407	131
Ammonia solution 20 - 22%		500 ml	420175	163
Hydrochloric acid 34-37%		500 ml	403915	396
Hydrochloric acid 34-37%		1 l	403916	396
Hydrochloric acid 34-37%		2.5 l	403917	396
Hydrochloric acid 29-31%		1 l	403921	398
Hydrofluoric acid 47-51%		500 ml	405716	411
Nitric acid 67-70%		500 ml	408115	529
Nitric acid 67-70%		1 l	408116	529
Nitric acid 67-70%		2.5 l	408117	529
Perchloric acid 65-71%		1 l	409193	560
Sulfuric acid 93-98%		500 ml	410405	734
Sulfuric acid 93-98%		1 l	410406	734
Sulfuric acid 93-98%		2.5 l	410407	734

Ultrapure Acids for trace metal analysis at ppt level

ULTRAPURE range is characterized by blank values generally between 50 and 1 ppt, for the 60-plus declared impurities. This extreme level of purity is obtained using the double sub-boiling distillation process and preserved in Teflon packaging, preconditioned with hot acid for at least one week.

Description	Notes	Size	Code	Page
Acetic acid glacial		500 ml	401361	131
Ammonia solution 20 - 22%		500 ml	420161	163
Hydrochloric acid 32-35%		500 ml	403891	396
Hydrofluoric acid 47-51%		500 ml	405611	411
Hydrogen peroxide solution 30-32%		500 ml	412051	413
Nitric acid 67-69%		250 ml	408052	530
Nitric acid 67-69%		500 ml	408051	530
Sulfuric acid 93-98%		500 ml	410351	733
Water		500 ml	412185	787

Standard solutions 1.000 ppm for AAS

Atomic absorption is the most sensitive technique available to analysts for the determination of metal impurities. It is a technique based on a sequential system which is generally slow but achieves sensitivity limits unattainable with other instrumental techniques.
The following ready-to-use standard single-element solutions are obtained by dissolution of the metal, at a purity level of 99.9%, in hydrochloric acid.
They are characterized by:
- Concentration of the metal equal to 1.000 ppm
- Available in 100ml and 500ml bottles in polyethylene or glass depending on compatibility
- Certificate of analysis with references on the N.I.S.T. Standard Reference Materials and uncertainty

Description	Notes	Size	Code	Page
Aluminum standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	504190	152
Aluminum standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	504186	152
Aluminum standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	E497405	152
Aluminum standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497401	152
Antimony standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	507525	183
Antimony standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	500 ml	507479	183
Antimony standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	E497415	183
Antimony standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497411	183
Arsenic standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	504439	185
Arsenic standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507496	185
Barium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507527	191
Barium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507481	191
Barium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	E497445	191
Barium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497441	191
Beryllium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	506941	201
Beryllium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507497	201
Bismuth standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507528	202
Bismuth standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507482	202
Bismuth standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	E497455	202
Bismuth standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497451	202
Boron standard solution	conc. 1.000 ppm Matrix: Water	100 ml	E497465	208
Boron standard solution	conc. 1.000 ppm Matrix: Water	500 ml	E497461	208
Cadmium standard solution	conc. 1000 ppm Matrix: Nitric acid	100 ml	507529	232
Cadmium standard solution	conc. 1000 ppm Matrix: Nitric acid	500 ml	507483	232
Calcium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507530	235
Calcium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507476	235
Calcium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	E497485	235
Calcium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497481	235
Cerium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507531	250
Cerium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507498	250
Cesium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507532	253
Cesium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507499	253
Chromium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	504195	265
Chromium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507485	265
Chromium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497501	265
Cobalt standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507533	271
Cobalt standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507484	271
Cobalt standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	E497495	271
Cobalt standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497491	271
Copper standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	504545	276

Copper standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507478	276
Copper standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	E497615	276
Copper standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497611	276
Dysprosium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507734	326
Dysprosium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507500	326
Erbium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507735	329
Erbium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507501	329
Europium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507736	352
Europium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507502	352
Gadolinium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507737	366
Gadolinium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507504	366
Gallium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507739	367
Gallium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507503	367
Germanium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507740	369
Germanium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507505	369
Gold standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	E497585	375
Gold standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497581	375
Hafnium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507741	378
Hafnium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507506	378
Holmium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507742	391
Holmium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507507	391
Indium standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	507743	423
Indium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507508	423
Iron standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	504194	429
Iron standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507393	429
Iron standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	E497515	429
Iron standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497511	429
Lanthanum standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507744	454
Lanthanum standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507509	454
Lead standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507752	456
Lead standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507490	456
Lead standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	E497595	456
Lead standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	E497591	456
Lithium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507745	462
Lithium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507486	462
Lithium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	E497525	462
Lithium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497521	462
Magnesium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503718	469
Magnesium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503719	469
Magnesium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	E497535	469
Magnesium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497531	469
Magnesium standard solution	conc. 1.000 ppm Matrix: Water	500 ml	507039	469
Manganese standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507746	479
Manganese standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507488	479
Manganese standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	E497545	479
Manganese standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497541	479
Mercury standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	497555	484
Mercury standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	497551	484
Mercury standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503640	484

Mercury standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507489	484
Molybdenum standard solution	conc. 1.000 ppm Matrix: Hydrofluoric acid and nitric acid	100 ml	507747	511
Molybdenum standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	500 ml	507494	511
Molybdenum standard solution	conc. 1.000 ppm Matrix: Ammonium hydroxide	100 ml	E497565	511
Molybdenum standard solution	conc. 1.000 ppm Matrix: Ammonium hydroxide	500 ml	E497561	511
Neodymium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507748	520
Neodymium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507510	520
Nickel standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507749	522
Nickel standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507487	522
Nickel standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	E497575	522
Nickel standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497571	522
Niobium standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	507750	527
Niobium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507511	527
Palladium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507751	553
Palladium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507512	553
Potassium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507753	582
Potassium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	506960	582
Potassium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	E497605	582
Potassium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497601	582
Rhenium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507754	638
Rhenium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507513	638
Rubidium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507755	641
Rubidium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507514	641
Samarium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507756	643
Samarium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507515	643
Scandium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507757	645
Scandium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507516	645
Selenium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507758	647
Selenium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507491	647
Selenium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	E497625	647
Selenium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497621	647
Silicon standard solution	conc. 1.000 ppm Matrix: Water	100 ml	E497635	651
Silicon standard solution	conc. 1.000 ppm Matrix: Water	500 ml	E497631	651
Silver standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507526	653
Silver standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507480	653
Sodium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507759	661
Sodium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503749	661
Sodium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	E497645	661
Sodium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497641	661
Strontium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507760	723
Strontium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507493	723
Strontium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	E497665	723
Strontium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497661	723
Tantalum standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507761	745
Tantalum standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507517	745
Tellurium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507762	747
Tellurium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507518	747
Thulium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507763	757

Thulium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507519	757
Tin standard solution	conc. 1.000 ppm Matrix: Hydrofluoric acid and nitric acid	100 ml	503949	760
Tin standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	500 ml	507492	760
Tin standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	E497655	760
Tin standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497651	760
Titanium standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	507764	763
Titanium standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	500 ml	507520	763
Tungsten standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	507765	779
Tungsten standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	500 ml	507521	779
Vanadium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507766	783
Vanadium standard solution	conc. 1000 ppm Matrix: Nitric acid	500 ml	504187	783
Vanadium standard solution	conc. 1.000 ppm Matrix: Sulfuric acid	100 ml	E497675	783
Vanadium standard solution	conc. 1.000 ppm Matrix: Sulfuric acid	500 ml	E497671	783
Ytterbium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507768	794
Ytterbium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507523	794
Yttrium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507767	795
Yttrium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507522	795
Zinc standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	507769	798
Zinc standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	E497685	798
Zinc standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	507477	798
Zinc standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	E497681	798
Zirconium standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	507770	804
Zirconium standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	500 ml	507524	804

Auxiliary products for AAS: Blank, Ionisation standards and matrix modifiers

Ionization buffers, prepared from high-purity salts (99,999%), are used to eliminate ionization phenomena and chemical interference in AAS. Matrix modifiers for graphite furnace AAS can be used diluted or mixed depending on the user's needs, are useful in eliminating the matrix effect in AAS.

Description	Notes	Size	Code	Page
Ammonium di-hydrogen phosphate 25 mg/L solution	Matrix: 1% Nitric acid	50 ml	503194	172
Ammonium nitrate 200 mg/l solution	Matrix: Water	50 ml	503195	174
Cadmium standard solution	conc. 5 ppb - Matrix: 2% Nitric acid	50 ml	504360	232
Cesium chloride 25 g/l solution	Matrix: Water	500 ml	504536	254
Copper standard solution	conc. 10 ppb - Matrix: 2% Nitric acid	50 ml	504361	276
Lanthanum chloride 25 g/l solution	Matrix: Hydrochloric acid	500 ml	504537	454
Lead standard solution	conc. 10 ppb - Matrix: 1% Nitric acid	50 ml	504364	456
Magnesium nitrate 10 g/l solution	Matrix: Water	50 ml	503196	473
Manganese standard solution	conc. 20 ppb - Matrix: 2% Nitric acid	50 ml	504362	480
Mercury standard solution	conc. 0.5 ppm - Matrix: 2% Nitric acid	100 ml	504370	485
Nickel standard solution	conc. 10 ppb - Matrix: 2% Nitric acid	50 ml	504363	522
Nickel (II) nitrate 10g/l	Matrix: 1% Nitric acid	50 ml	503197	524
Palladium nitrate 2 g/l solution	Matrix: 1% Nitric acid	50 ml	503198	554
Palladium nitrate 2 g/l solution	Matrix: 15% Nitric acid	50 ml	503202	554
Potassium chloride 25g/l in HCl	Matrix: 2% Hydrochloric acid	500 ml	504538	589

Water deionized and acidified	Matrix: 2 % Nitric acid	1 l	504550	789
Water deionized and acidified	Matrix: 5 % Nitric acid	1 l	504551	789
Water deionized and acidified	Matrix: 10 % Nitric acid	1 l	504552	789
Water deionized and acidified	Matrix: 2 % Hydrochloric acid	1 l	504553	789
Water deionized and acidified	Matrix: 5 % Hydrochloric acid	1 l	504554	789
Water deionized and acidified	Matrix: 10 % Hydrochloric acid	1 l	504557	789

Monoelement Standard solutions 10.000 ppm for ICP

These standard solutions are obtained by dissolution of the metal, at a purity level of 99.99%, in an acid (usually nitric acid).
They are characterized by:
- Produced under a Quality Management System that is Accredited according to ISO/IEC 17025 and ISO 17034
- Concentrations of 10.000 ppm
- Guaranteed titer with its uncertainty
- Raw materials selected and verified against N.I.S.T. Standard Reference Materials;
- Available in 100ml and 500ml polyethylene bottles
- Certificate of analysis with references to the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Description	Notes	Size	Code	Page
Aluminum standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	100 ml	503415	152
Aluminum standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	500 ml	503417	152
Antimony standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	100 ml	503895	182
Antimony standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	500 ml	503897	182
Arsenic standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503425	185
Arsenic standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503427	185
Barium standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	100 ml	503455	191
Barium standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	500 ml	503457	191
Bismuth standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503475	202
Bismuth standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503477	202
Boron standard solution	conc. 10.000 ppm Matrix: Water	100 ml	503445	208
Boron standard solution	conc. 10.000 ppm Matrix: Water	500 ml	503447	208
Cadmium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503495	232
Cadmium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503497	232
Calcium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503485	235
Calcium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503487	235
Cerium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503505	250
Cerium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503507	250
Cesium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503535	253
Cesium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503537	253
Chromium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503525	265
Chromium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503527	265
Cobalt standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503515	271
Cobalt standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503517	271
Copper standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503545	276
Copper standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503547	276
Dysprosium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	504235	325
Dysprosium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	504237	325
Erbium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	504245	329
Erbium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	504247	329
Europium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503575	351
Europium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503577	351
Gadolinium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503605	366
Gadolinium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503607	366

Germanium standard solution	conc. 10.000 ppm Matrix: Hydrofluoric acid and nitric acid	100 ml	504255	369
Germanium standard solution	conc. 10.000 ppm Matrix: Hydrofluoric acid and nitric acid	500 ml	504257	369
Gold standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	100 ml	503435	375
Gold standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	500 ml	503437	375
Hafnium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	504225	378
Hafnium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	504227	378
Holmium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	504265	390
Holmium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	504267	390
Indium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503655	423
Indium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503657	423
Iron standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503585	428
Iron standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503587	428
Lanthanum standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503685	453
Lanthanum standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503687	453
Lead standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503805	456
Lead standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503807	456
Lithium standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	100 ml	503695	462
Lithium standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	500 ml	503697	462
Lutetium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503615	467
Lutetium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503617	467
Magnesium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503715	469
Magnesium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503717	469
Manganese standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503725	479
Manganese standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503727	479
Mercury standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503635	484
Mercury standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503637	484
Molybdenum standard solution	conc. 10.000 ppm Matrix: Ammonium hydroxide	100 ml	503735	510
Molybdenum standard solution	conc. 10.000 ppm Matrix: Ammonium hydroxide	500 ml	503737	510
Neodymium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503765	520
Neodymium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503767	520
Nickel standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503775	522
Nickel standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503777	522
Niobium standard solution	conc. 10.000 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	503755	527
Niobium standard solution	conc. 10.000 ppm Matrix: Nitric acid and hydrofluoric acid	500 ml	503757	527
Palladium standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	100 ml	503815	553
Palladium standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	500 ml	503817	553
Phosphorus standard solution	conc. 10.000 ppm Matrix: Water	100 ml	503795	577
Phosphorus standard solution	conc. 10.000 ppm Matrix: Water	500 ml	503797	577
Platinum standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	100 ml	503835	581
Platinum standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	500 ml	503837	581
Potassium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503675	582
Potassium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503677	582
Praseodymium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503825	621
Praseodymium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503827	621
Rhodium standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	100 ml	503865	639
Rubidium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503845	640

Rubidium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503847	640
Ruthenium standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	100 ml	503875	641
Ruthenium standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	500 ml	503877	641
Samarium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503935	643
Samarium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503937	643
Scandium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503905	645
Scandium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503907	645
Selenium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503915	647
Selenium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503917	647
Silicon standard solution	conc. 10.000 ppm Matrix: Water	100 ml	503925	651
Silicon standard solution	conc. 10.000 ppm Matrix: Water	500 ml	503927	651
Silicon standard solution	conc. 10.000 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	504275	651
Silicon standard solution	conc. 10.000 ppm Matrix: Nitric acid and hydrofluoric acid	500 ml	504277	651
Silver standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503405	653
Silver standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503407	653
Sodium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503745	661
Sodium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503747	661
Strontium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503955	722
Strontium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503957	722
Sulfur standard solution	conc. 10.000 ppm Matrix: Water	100 ml	504295	730
Sulfur standard solution	conc. 10.000 ppm Matrix: Water	500 ml	504297	730
Tantalum standard solution	conc. 10.000 ppm Matrix: Hydrofluoric acid and nitric acid	100 ml	503965	745
Tantalum standard solution	conc. 10.000 ppm Matrix: Hydrofluoric acid and nitric acid	500 ml	503967	745
Tellurium standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	100 ml	503985	747
Tellurium standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	500 ml	503987	747
Terbium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	503975	748
Terbium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	503977	748
Thallium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	504015	754
Thallium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	504017	754
Tin standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	100 ml	503945	760
Tin standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	500 ml	503947	760
Titanium standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid and hydrofluoric acid	100 ml	504005	763
Titanium standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid and hydrofluoric acid	500 ml	504007	763
Tungsten standard solution	conc. 10.000 ppm Matrix: 10% ammonia	100 ml	504055	779
Tungsten standard solution	conc. 10.000 ppm Matrix: 10% ammonia	500 ml	504057	779
Uranium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	504035	782
Uranium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	504037	782
Vanadium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	504045	783
Vanadium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	504047	783
Ytterbium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	504075	794
Ytterbium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	504077	794
Yttrium standard solution	conc. 10.000 ppm Matrix: Nitric acid	100 ml	504065	795
Yttrium standard solution	conc. 10.000 ppm Matrix: Nitric acid	500 ml	504067	795
Zinc standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	100 ml	504085	798
Zinc standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid	500 ml	504087	798

Product specifications are subject to changes. Please visit our website for updates.

Zirconium standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid and hydrofluoric acid	100 ml	504095	804
Zirconium standard solution	conc. 10.000 ppm Matrix: Hydrochloric acid and hydrofluoric acid	500 ml	504097	804

Monoelement Standard solutions 1.000 ppm for ICP

ICP is a widely used analytical technique for trace metal analysis. It is based on a simultaneous system which allows quick and convenient analyses for a large number of determinable elements. One of the latest technological advances in the area of coupled analytical techniques involves the optimization of the ICP-MS technique, a versatile and vital instrument for the quick and reliable analysis of trace and ultra-trace metals.
These standard solutions are obtained by dissolution of the metal, at a purity level of 99.99%, in an acid (usually nitric acid).
They are characterized by:
- Produced under a Quality Management System that is Accredited according to ISO/IEC 17025 and ISO 17034
- Concentrations of 1.000 ppm
- Guaranteed titer with its uncertainty
- Raw materials selected and verified against N.I.S.T. Standard Reference Materials
- Available in 100ml and 500ml polyethylene bottles
- Certificate of analysis with references to the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Description	Notes	Size	Code	Page
Aluminum standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	503411	152
Aluminum standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	503413	152
Antimony standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	503891	182
Antimony standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	503893	182
Antimony standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	503899	182
Antimony standard solution	conc. 1.000 ppm Matrix: Hydrofluoric acid and nitric acid	500 ml	503898	182
Arsenic standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503421	185
Arsenic standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503423	185
Barium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	503451	191
Barium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	503453	191
Beryllium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	503461	201
Beryllium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	503463	201
Bismuth standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503471	202
Bismuth standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503473	202
Boron standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503441	208
Boron standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503443	208
Cadmium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503491	232
Cadmium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503493	232
Calcium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503481	235
Calcium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503483	235
Cerium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503501	250
Cerium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503503	250
Cesium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503531	253
Cesium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503533	253
Chromium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503521	265
Chromium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503523	265
Cobalt standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503511	271
Cobalt standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503513	271
Copper standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503541	276
Copper standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503543	276
Dysprosium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	504231	325
Dysprosium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	504233	325
Erbium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	504241	329

Erbium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	504243	329
Europium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503571	351
Europium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503573	351
Gadolinium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503601	366
Gadolinium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503603	366
Germanium standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	504251	369
Germanium standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	500 ml	504253	369
Gold standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	503431	375
Gold standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	503433	375
Hafnium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	504221	378
Hafnium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	504223	378
Holmium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	504261	390
Holmium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	504263	390
Indium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503651	423
Indium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503653	423
Iron standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503581	428
Iron standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503583	428
Lanthanum standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503681	453
Lanthanum standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503683	453
Lead standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503801	456
Lead standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503803	456
Lithium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	503691	462
Lithium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	503693	462
Lutetium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503611	467
Lutetium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503613	467
Magnesium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503711	469
Magnesium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503713	469
Manganese standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503721	479
Manganese standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503723	479
Mercury standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503631	484
Mercury standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503633	484
Molybdenum standard solution	conc. 1.000 ppm Matrix: Ammonium hydroxide	100 ml	503731	510
Molybdenum standard solution	conc. 1.000 ppm Matrix: Ammonium hydroxide	500 ml	503733	510
Neodymium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503761	520
Neodymium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503763	520
Nickel standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503771	522
Nickel standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503773	522
Niobium standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	503751	527
Niobium standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	500 ml	503753	527
Palladium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	503811	553
Palladium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	503813	553
Phosphorus standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503791	577
Phosphorus standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503793	577
Platinum standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	503831	581
Platinum standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	503833	581
Potassium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503671	582
Potassium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503673	582

Praseodymium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503821	621
Praseodymium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503823	621
Rhodium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	503861	639
Rhodium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	503863	639
Rubidium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503841	640
Rubidium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503843	640
Ruthenium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	503871	641
Ruthenium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	503873	641
Samarium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503931	643
Samarium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503933	643
Scandium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503901	645
Scandium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503903	645
Selenium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503911	647
Selenium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503913	647
Silicon standard solution	conc. 1.000 ppm Matrix: Water	100 ml	503921	651
Silicon standard solution	conc. 1.000 ppm Matrix: Water	500 ml	503923	651
Silicon standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	504271	651
Silicon standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	500 ml	504273	651
Silver standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503401	653
Silver standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503403	653
Sodium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503741	661
Sodium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503743	661
Strontium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503951	722
Strontium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503953	722
Sulfur standard solution	conc. 1.000 ppm Matrix: Water	100 ml	504291	730
Sulfur standard solution	conc. 1.000 ppm Matrix: Water	500 ml	504293	730
Tantalum standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	503961	745
Tantalum standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	500 ml	503963	745
Tellurium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	503981	747
Tellurium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	503983	747
Terbium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	503971	748
Terbium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	503973	748
Thallium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	504011	754
Thallium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	504013	754
Thorium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	504281	756
Thorium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	504283	756
Tin standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	503941	760
Tin standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	503943	760
Titanium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid and hydrofluoric acid	100 ml	504001	763
Titanium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid and hydrofluoric acid	500 ml	504003	763
Tungsten standard solution	conc. 1.000 ppm Matrix: Water	100 ml	504058	779
Tungsten standard solution	conc. 1.000 ppm Matrix: 4% ammonia	100 ml	504051	779
Tungsten standard solution	conc. 1.000 ppm Matrix: 4% ammonia	500 ml	504053	779
Uranium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	504031	782
Uranium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	504033	782
Vanadium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	504041	783

Vanadium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	504043	783
Ytterbium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	504071	794
Ytterbium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	504073	794
Yttrium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	504061	795
Yttrium standard solution	conc. 1.000 ppm Matrix: Nitric acid	500 ml	504063	795
Zinc standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	504081	798
Zinc standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	500 ml	504083	798
Zirconium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid and hydrofluoric acid	100 ml	504091	804
Zirconium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid and hydrofluoric acid	500 ml	504093	804

Monoelement Standard solutions 1.000 ppm for ICP-MS

These standard solutions are obtained by dissolution of the metal, at a purity level of 99.999%, in an acid (usually nitric acid).
They are characterized by:
- Produced under a Quality Management System that is Accredited according to ISO/IEC 17025 and ISO 17034
- Concentrations of 1.000 ppm
- Guaranteed titer with its uncertainty
- Raw materials selected and verified against N.I.S.T. Standard Reference Materials and packaged in a cleanroom
- Available in 100ml LDPE bottles
- Certificate of analysis with references to the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Description	Notes	Size	Code	Page
Aluminum standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505309	151
Antimony standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505833	182
Arsenic standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505313	185
Barium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505329	191
Beryllium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505333	201
Bismuth standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505339	202
Boron standard solution	conc. 1.000 ppm Matrix: Water	100 ml	505323	208
Cadmium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505549	232
Calcium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505543	235
Chromium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505569	264
Cobalt standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505563	270
Copper standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505579	275
Gallium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505619	367
Germanium standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505633	369
Gold standard solution	conc. 1000 ppm Matrix: Hydrochloric acid	100 ml	505319	375
Indium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505663	422
Iron standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505613	428
Lanthanum standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505693	453
Lead standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505769	456
Lithium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505703	462
Lutetium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505709	467
Magnesium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505713	469
Manganese standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505719	479
Mercury standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505654	484
Molybdenum standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505723	510
Nickel standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505753	522
Phosphorus standard solution	conc. 1.000 ppm Matrix: Water	100 ml	505763	577
Platinum standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	505789	580
Potassium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505683	582

Rhodium standard solution	conc. 1.000 ppm Matrix: Hydrochloric acid	100 ml	505809	638
Scandium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505839	645
Selenium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505843	647
Silver standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505303	653
Sodium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505733	661
Strontium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505869	722
Sulfur standard solution	conc. 1.000 ppm Matrix: Water	100 ml	505823	730
Terbium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505883	748
Thallium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505913	754
Tin standard solution	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505863	760
Titanium standard solution	conc. 1.000 ppm Matrix: Hydrofluoric acid and nitric acid	100 ml	505909	763
Vanadium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505929	783
Yttrium standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505943	794
Zinc standard solution	conc. 1.000 ppm Matrix: Nitric acid	100 ml	505953	797

Monoelement Standard solutions 100 ppm for ICP-MS

These standard solutions are obtained by dissolution of the metal, at a purity level of 99.999%, in an acid (usually nitric acid).
They are characterized by:
- Produced under a Quality Management System that is Accredited according to ISO/IEC 17025 and ISO 17034
- Concentrations of 100 ppm
- Guaranteed titer with its uncertainty
- Raw materials selected and verified against N.I.S.T. Standard Reference Materials and packaged in a cleanroom
- Available in 100ml LDPE bottles
- Certificate of analysis with references to the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Description	Notes	Size	Code	Page
Aluminum standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505308	151
Antimony standard solution	conc. 100 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505835	182
Arsenic standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505315	185
Barium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505328	191
Beryllium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505335	201
Bismuth standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505338	202
Boron standard solution	conc. 100 ppm Matrix: Water	100 ml	505325	208
Cadmium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505548	232
Calcium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505545	235
Cerium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505555	250
Cesium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505575	253
Chromium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505568	264
Cobalt standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505565	270
Copper standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505578	275
Dysprosium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505585	325
Erbium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505595	329
Europium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505605	351
Gadolinium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505625	366
Gallium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505618	367
Germanium standard solution	conc. 100 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505635	369
Gold standard solution	conc. 100 ppm Matrix: Hydrochloric acid	100 ml	505318	375
Hafnium standard solution	conc. 100 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505645	378
Holmium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505658	390

Indium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505665	422
Iridium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505675	427
Iron standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505615	428
Lanthanum standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505695	453
Lead standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505768	456
Lithium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505705	462
Lutetium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505708	467
Magnesium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505715	469
Manganese standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505718	479
Mercury standard solution	conc. 100 ppm Matrix: Hydrochloric acid	100 ml	505655	484
Mercury standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	506918	484
Molybdenum standard solution	conc. 100 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505725	510
Neodymium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505745	519
Nickel standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505755	522
Niobium standard solution	conc. 100 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505738	526
Osmium standard solution	conc. 100 ppm Matrix: Hydrochloric acid	100 ml	505758	549
Palladium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505775	553
Phosphorus standard solution	conc. 100 ppm Matrix: Water	100 ml	505765	577
Platinum standard solution	conc. 100 ppm Matrix: Hydrochloric acid	100 ml	505788	580
Potassium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505685	582
Praseodymium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505785	621
Rhenium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505805	638
Rhodium standard solution	conc. 100 ppm Matrix: Hydrochloric acid	100 ml	505808	638
Rubidium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505795	640
Ruthenium standard solution	conc. 100 ppm Matrix: Hydrochloric acid	100 ml	505815	641
Samarium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505855	643
Scandium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505838	645
Selenium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505845	647
Silicon standard solution	conc. 100 ppm Matrix: Water	100 ml	505848	651
Silver standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505305	653
Sodium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505735	661
Strontium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505868	722
Sulfur standard solution	conc. 100 ppm Matrix: Water	100 ml	505825	730
Tantalum standard solution	conc. 100 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505875	744
Tellurium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505888	747
Terbium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505885	748
Thallium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505915	754
Thulium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505918	756
Tin standard solution	conc. 100 ppm Matrix: Hydrofluoric acid and nitric acid	100 ml	505865	760
Titanium standard solution	conc. 100 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505908	763
Tungsten standard solution	conc. 100 ppm. Matrix: Ammonium hydroxyde	100 ml	505935	779
Uranium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505923	782
Vanadium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505928	783
Ytterbium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505948	794
Yttrium standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505945	794
Zinc standard solution	conc. 100 ppm Matrix: Nitric acid	100 ml	505955	797

Zirconium standard solution	conc. 100 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505958	804

Monoelement Standard solutions 10 ppm for ICP-MS

These standard solutions are obtained by dissolution of the metal, at a purity level of 99.999%, in an acid (usually nitric acid).
They are characterized by:
- Produced under a Quality Management System that is Accredited according to ISO/IEC 17025 and ISO 17034
- Concentrations of 10 ppm
- Guaranteed titer with its uncertainty
- Raw materials selected and verified against N.I.S.T. Standard Reference Materials and packaged in a cleanroom
- Available in 100ml LDPE bottles
- Certificate of analysis with references to the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Description	Notes	Size	Code	Page
Aluminum standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505307	151
Antimony standard solution	conc. 10 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505832	182
Arsenic standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505312	185
Barium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505327	191
Beryllium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505332	201
Bismuth standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505337	202
Boron standard solution	conc. 10 ppm Matrix: Water	100 ml	505322	208
Cadmium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505547	232
Calcium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505542	235
Cerium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505552	250
Cesium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505572	253
Chromium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505567	264
Cobalt standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505562	270
Copper standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505577	275
Dysprosium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505582	325
Erbium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505592	329
Europium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505602	351
Gadolinium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505622	366
Gallium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505617	367
Germanium standard solution	conc. 10 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505632	369
Gold standard solution	conc. 10 ppm Matrix: Hydrochloric acid	100 ml	505317	375
Hafnium standard solution	conc. 10 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505642	378
Holmium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505657	390
Indium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505662	422
Iron standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505612	428
Lanthanum standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505692	453
Lead standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505767	456
Lithium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505702	462
Lutetium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505707	467
Magnesium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505712	469
Manganese standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505717	479
Mercury standard solution	conc. 10 ppm Matrix: Nitric acid	100 ml	505652	484
Molybdenum standard solution	conc. 10 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505722	510
Neodymium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505742	519
Nickel standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505752	522
Niobium standard solution	conc. 10 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505737	526

Palladium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505772	553
Phosphorus standard solution	conc. 10 ppm Matrix: Water	100 ml	505762	577
Platinum standard solution	conc. 10 ppm Matrix: Hydrochloric acid	100 ml	505787	580
Potassium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505682	582
Praseodymium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505782	621
Rhenium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505802	638
Rhodium standard solution	conc. 10 ppm Matrix: Hydrochloric acid	100 ml	505807	638
Rubidium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505792	640
Ruthenium standard solution	conc. 10 ppm Matrix: Hydrochloric acid	100 ml	505812	641
Samarium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505852	643
Scandium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505837	645
Selenium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505842	647
Silicon standard solution	conc. 10 ppm Matrix: Water	100 ml	505847	651
Silver standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505302	653
Sodium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505732	661
Strontium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505867	722
Sulfur standard solution	conc. 10 ppm Matrix: Water	100 ml	505822	730
Tantalum standard solution	conc. 10 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505872	744
Tellurium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505887	747
Terbium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505882	748
Thallium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505912	754
Thulium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505917	756
Tin standard solution	conc. 10 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505862	760
Titanium standard solution	conc. 10 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505907	763
Tungsten standard solution	conc. 10 ppm. Matrix: Ammonium hydroxyde	100 ml	505932	779
Uranium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505922	782
Vanadium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505927	783
Ytterbium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505947	794
Yttrium standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505942	794
Zinc standard solution	conc. 10 ppm. Matrix: Nitric acid	100 ml	505952	797
Zirconium standard solution	conc. 10 ppm Matrix: Nitric acid and hydrofluoric acid	100 ml	505957	804

Multielement Standard Solutions for ICP

These standard solutions are intended for calibration of the instrument and are ideal for checking the reproducibility of analytical results. They are obtained by dissolution of various metals in an acid, usually hydrochloric or nitric acid.
They are characterized by:
- Produced under a Quality Management System that is Accredited according to ISO/IEC 17025 and ISO 17034
- Guaranteed titer with a tolerance of 0.2% at the 95% confidence level
- Raw materials selected and verified against N.I.S.T. Standard Reference Materials
- Available in 100ml and 500ml polyethylene bottles

Contact us with your specific mixture (CAS, concentration, solvent, volume, packaging) and you will receive our offer according to your needs.

Description	Notes	Size	Code	Page
Multielement standard for ICP	13 elements: Al, Mg, Cr, Mn, Cu, Rh, In, Cd, Ce, Pb, Th, B, Ba 0,01mg/ml each - Matrix: Nitric acid	500 ml	504396	513
Multielement standard for ICP	22 elements: As, Be, Bi, Ca, Cd, Co, Cr, Cu, Fe, Li, Mg, Mn, Mo, Ni, Pb, Sb, Se, Sr, Ti, Tl, V, Zn 1ppm each - Matrix: Nitric acid	100 ml	504350	514

Multielement standard for ICP	33 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Cd, Co, Cs, Cr, Cu, Fe, In, K, Li, Mg, Mn, Mo, Na, Nb, Ni, Pb, Rb, Sb, Se, Sr, Ti, Tl, U, V, Zn 1ppm each - Matrix: Nitric acid	100 ml	504354	514
Multielement standard for ICP	40 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Cd, Co, Cs, Cr, Cu, Fe, Ga, Ge, In, K, Li, Mg, Mn, Mo, Na, Ni, Nb, Pb, Pd, Rb, Sb, Se, Sn, Sr, Ta, Ti, Tl, U, V, W, Zn, Zr 1ppm each - Matrix: Nitric acid	100 ml	504356	514
Multielement standard for ICP	22 elements: As, Be, Bi, Ca, Cd, Co, Cr, Cu, Fe, Li, Mg, Mn, Mo, Ni, Pb, Sb, Se, Sr, Ti, Tl, V, Zn 1ppm each - Matrix: Nitric acid	500 ml	504351	514
Multielement standard for ICP	28 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Ce, Co, Cr, Cu, Fe, K, Li, Mg, Mn, Mo, Na, Ni, Pb, Sb, Se, Sr, Ti, Tl, V, Zn 1ppm each - Matrix: Nitric acid	500 ml	504353	514
Multielement standard for ICP	33 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Cd, Co, Cs, Cr, Cu, Fe, In, K, Li, Mg, Mn, Mo, Na, Nb, Ni, Pb, Rb, Sb, Se, Sr, Ti, Tl, U, V, Zn 1ppm each - Matrix: Nitric acid	500 ml	504355	514
Multielement standard for ICP	40 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Cd, Co, Cs, Cr, Cu, Fe, Ga, Ge, In, K, Li, Mg, Mn, Mo, Na, Ni, Nb, Pb, Pd, Rb, Sb, Se, Sn, Sr, Ta, Ti, Tl, U, V, W, Zn, Zr 1ppm each - Matrix: Nitric acid	500 ml	504357	514
Multielement standard for ICP	6 elements: Au, Ir, Pb, Pt, Rh, Ru 100ppm each - Matrix: Nitric acid	100 ml	504301	514
Multielement standard for ICP	16 elements: Al, As, Ba, Be, Bi, Ca, Cs, Ga, In, K, Li, Mg, Na, Rb, Se, Sr 100ppm each - Matrix: Nitric acid	100 ml	504303	514
Multielement standard for ICP	13 elements: Ag, Cd, Co, Cr, Cu, Fe, Hg, Mn, Ni, Pb, Tl, V, Zn 100ppm each - Matrix: Nitric acid	100 ml	504305	514
Multielement standard for ICP	24 components: Ag 1mg/kg; Sb 1mg/kg; As 1mg/kg; Mn 1mg/kg; Cd 1mg/kg; Cr 1mg/kg; Ti 1mg/kg; Pb 1mg/kg; Co 1mg/kg; Ni 1mg/kg; Se 1mg/kg; V 1mg/kg; Mo 1mg/kg; Sn 1mg/kg; Ba 1mg/kg; Be 1mg/kg; Li 1mg/kg; Tl 1mg/kg; Bi 1mg/kg; Al 10mg/kg; Cu 10mg/kg; Fe 10mg/kg; B 10mg/kg; Zn 10mg/kg - Matrix: nitric acid	100 ml	504480	514
Multielement standard for ICP	22 elements: As, Be, Bi, Ca, Cd, Co, Cr, Cu, Fe, Li, Mg, Mn, Mo, Ni, Pb, Sb, Se, Sr, Ti, Tl, V, Zn 100ppm each - Matrix: Nitric acid	100 ml	504306	515
Multielement standard for ICP	28 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Cd, Co, Cr, Cu, Fe, K, Li, Mg, Mn, Mo, Na, Ni, Pb, Sb, Se, Sr, Ti, Tl, V, Zn 100ppm each - Matrix: Nitric acid	100 ml	504308	515
Multielement standard for ICP	33 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Cd, Co, Cs, Cr, Cu, Fe, In, K, Li, Mg, Mn, Mo, Na, Ni, Nb, Pb, Rb, Sb, Se, Sr, Ti, Tl, U, V, Zn 100ppm - Matrix: Nitric acid	100 ml	504310	515
Multielement standard for ICP	9 elements: Au, Ir, Os, Pb, Pt, Rh, Ru, Sn, Te 100ppm each - Matrix: Hydrochloric acid	100 ml	504312	515
Multielement standard for ICP	22 elements: As, Be, Bi, Ca, Cd, Co, Cr, Cu, Fe, Li, Mg, Mn, Mo, Ni, Pb, Sb, Se, Sr, Ti, Tl, V, Zn 100ppm each - Matrix: Nitric acid	500 ml	504307	515
Multielement standard for ICP	28 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Cd, Co, Cr, Cu, Fe, K, Li, Mg, Mn, Mo, Na, Ni, Pb, Sb, Se, Sr, Ti, Tl, V, Zn 100ppm each - Matrix: Nitric acid	500 ml	504309	515
Multielement standard for ICP	33 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Cd, Co, Cs, Cr, Cu, Fe, In, K, Li, Mg, Mn, Mo, Na, Ni, Nb, Pb, Rb, Sb, Se, Sr, Ti, Tl, U, V, Zn 100ppm each - Matrix: Nitric acid	500 ml	504311	515
Multielement standard for ICP	9 elements: Au, Ir, Os, Pb, Pt, Rh, Ru, Sn, Te 100ppm each - Matrix: Hydrochloric acid	500 ml	504313	515

Multielement standard for ICP and ICP-MS	28 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Ce, Co, Cr, Cu, Fe, K, Li, Mg, Mn, Mo, Na, Ni, Pb, Sb, Se, Sr, Ti, Tl, V, Zn 1ppm each - Matrix: Nitric acid	100 ml	504352	515

Multielement Standard Solutions for ICP-MS

These standard solutions are commonly used for tuning ICP-MS instruments and for quality control. They are characterized by:
• Produced under a Quality Management System that is Accredited according to ISO/IEC 17025 and ISO 17034
• Concentrations of 10 ppm
• Guaranteed titer with its uncertainty
• Raw materials selected and verified against N.I.S.T. Standard Reference Materials and packaged in a cleanroom
• Available in 100ml LDPE bottles
• Certificate of analysis with references to the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval
Contact us with your specific mixture (Elements, concentration, matrix, volume, quantity) and you will receive our offer according to your needs.

Description	Notes	Size	Code	Page
Multielement standard for ICP and ICP-MS	9 elements: Be, Mg, Co, In, Rh, Ce, Ba, Pb, U 10ppm each - Matrix: Nitric acid	100 ml	504392	515
Multielement standard for ICP and ICP-MS	13 elements: Ba, Be, Bi, Ce, Cu, Ho, In, Li, Mg, Pb, Tl, U, Y 10ppm each - Matrix: Nitric acid	100 ml	504393	515

TRACE ANALYSIS: ORGANICS

CARLO ERBA Reagents offers tailored formulation of organic substances according to the ISO 17025 accredited Quality Management System and ISO Guide 34.
These organic mixtures are specific for instrument calibration and traceable to NIST.
Specific mixtures could be prepared if you send us CAS number, concentration, solvent, volume and packaging.

Standard Solutions for Pesticides analysis

Description	Notes	Size	Code	Page
Organic standard: Pesticide mixture	45 components 10mg/ml each in cyclohexane/ acetone: Acetochlor [CAS:34256-82-1]10mg/l; Aclonifen [CAS:74070-46-5] 10mg/l; Alachlor [CAS:15972-60-8] 10mg/l; Bifenthrin [CAS:82657-04-3] 10mg/l; Cadusafos [CAS:95465-99-9] 10mg/l; Captan [CAS:133-06-2] 10mg/l; Carbofuran [CAS:1563-66-2] 10mg/l; Chlorfenvinphos [CAS:470-90-6]10mg/l; Chlormephos [CAS:24934-91-6] 10mg/l; Chlorothalonil [CAS:1897-45-6] 10mg/l Chlorpyrifos [CAS:2921-88-2] 10mg/l; Chlorpyrifos methyl [CAS:5598-13-0] 10mg/l; lambda-Cyhalothrin [CAS:91465-08-6] 10mg/l; Cypermethrin [CAS:52315-07-8]10mg/l; Delta-HCH [CAS:319-86-8] 10mg/l; Diazinon [CAS:333-41-5] 10mg/l; Dichlobenil [CAS:1194-65-6] 10mg/l; Dinoterb [CAS:1420-07-1] 10mg/l; Endosulfan-total (sulfate) [CAS:1031-07-8] 10mg/l; Fipronil [CAS:120068-37-3] 10mg/l; Folpet [CAS:133-07-3] 10mg/l; Heptachlor-endo-epoxide [CAS:28044-83-9] 10mg/l; Hexachloro-1,3-butadiene [CAS:87-68-3]10mg/l; Iprodione [CAS:36734-19-7] 10mg/l; Isofenphos [CAS:25311-71-1]10mg/l; Malathion [CAS:121-75-5] 10mg/l; Metazachlor [CAS:67129-08-2] 10mg/l; Oxadiazon [CAS:19666-30-9] 10mg/l; Oxyfluorfen [CAS:42874-03-3] 10mg/l; Parathion (Parathion-ethyl) [CAS:56-38-2] 10mg/l; Parathion-methyl [CAS:298-00-0] 10mg/l; Pendimethalin [CAS:40487-42-1] 10mg/l; Pentachlorobenzene [CAS:608-93-5] 10mg/l; Procymidone [CAS:32809-16-8] 10mg/l; Propachlor [CAS:1918-16-7] 10mg/l; Tebutam [CAS:35256-85-0] 10mg/l; Tefluthrin [CAS:79538-32-2] 10mg/l; Terbufos [CAS:13071-79-9] 10mg/l; Tolylfluanid [CAS:731-27-1] 10mg/l; Triazophos [CAS:24017-47-8] 10mg/l; Trifluralin [CAS:1582-09-8] 10mg/l; Vinclozolin [CAS:50471-44-8] 10mg/l; Piperonyl butoxide [CAS:51-03-6] 10mg/l; Metolachlor [CAS:51218-45-2] 5mg/l; S-Metolachlor [CAS:87392-12-9] 5mg/l	1 ml	506897	546

Organic standard: Pesticide mixture	79 components in acetone: Bifenthrin [CAS:82657-04-3] 120µg/ml; lambda-Cyhalothrin [CAS:91465-08-6] 100µg/ml; Cypermethrin [CAS:5231 5-07-8] 130µg/ml; Deltamethrin [CAS:52918-63-5] 130µg/ml; Fenvalerate [CAS:51630-58-1] 105µg/ml; Permethrin [CAS:52645-53-1] 100µg/ml; tau-Fluvalinate [CAS:102851-06-9] 100µg/ml; Tetramethrin [CAS:7696-12-0] 100µg/ml; Aldrin [CAS:309-00-2] 20µg/ml; cis-Chlordane [CAS:5103-71-9] 20µg/ml; trans-Chlordane [CAS:5103-74-2] 20µg/ml; 2,4'-DDD [CAS:53-19-0] 20µg/ml; 4,4'-DDD (TDE) [CAS:72-54-8] 20µg/ml; 2,4'-DDE [CAS:3424-82-6] 20µg/ml; 4,4'-DDE [CAS:72-55-9] 20µg/ml; 2,4'-DDT [CAS:789-02-6] 20µg/ml; 4,4'-DDT [CAS:50-29-3] 20µg/ml; Dieldrin [CAS:60-57-1] 20µg/ml; Endosulfan-alpha [CAS:959-98-8] 20µg/ml; Endosulfan-beta [CAS:33213-65-9] 20µg/ml; Endosulfan-total (sulfate) [CAS:1031-07-8] 20µg/ml; Endrin [CAS:72-20-8] 20µg/ml; Endrin aldehyde [CAS:7421-93-4] 20µg/ml; Alp ha-HCH [CAS:319-84-6] 20µg/ml; Beta-HCH [CAS:319-85-7] 20µg/ml; Delta-HCH [CAS:319-86-8] 20µg/ml; Gamma-HCH (Lindane) [CAS:58-89-9] 20µg/ml; Heptachlor [CAS:76-44-8] 20µg/ml; Heptachlor-endo-epoxide [CAS:28044-83-9] 20µg/ml; Heptachlor-exo-epoxide [CAS:1024-57-3] 20µg/ml;Hexachlorobenzene [CAS:118-74-1] 20µg/ml; PCB 209 [CAS:2051-24-3] 20µg/ml; PCB 29 [CAS:15862-07-4] 20µg/ml; Vinclozolin [CAS:50471-44-8] 20µg/ml; Alachlor [CAS:15972-60-8] 100µg/ml; Bromopropylate [CAS:18181-80-1] 50µg/ml; Chlorothalonil [CAS:1897-45-6] 25µg/ml; Dicofol [CAS:115-32-2] 75µg/ml; Iprodione [CAS:36734-19-7] 200µg/ml; Nitrofen [CAS:1836-75-5] 20µg/ml; oxy-Chlordane [CAS:27304-13-8] 20µg/ml; Phosalone [CAS:2310-17-0] 20µg/ml; Procymidone [CAS:32809-16-8] 150µg/ml; Tetradifon [CAS:116-29-0] 20µg/ml; Bromophos-ethyl [CAS:4824-78-6] 100µg/ml; Bromophos-methyl [CAS:2104-96-3] 100µg/ml; Chlorfenvinphos [CAS:470-90-6] 100µg/ml; Chlorpyrifos (Chlorpyrifos-ethyl) [CAS:2921-88-2] 100µg/ml; Chlorpyrifos methyl [CAS:5598-13-0] 100µg/ml; Diazinon [CAS:333-41-5] 100µg/ml; Dichlorvos [CAS:62-73-7] 100µg/ml; Dimethoate [CAS:60-51-5] 100µg/ml; Disulfoton [CAS:298-04-4] 50µg/ml; Fenchlorphos [CAS:299-84-3] 100µg/ml; Fenthion [CAS:55-38-9] 100µg/ml; Malathion [CAS:121-75-5] 100µg/ml; Parathion (Parathion-ethyl) [CAS:56-38-2] 100µg/ml; Parathion-methyl [CAS:298-00-0] 100µg/ml; Pirimiphos-methyl [CAS:29232-93-7] 100µg/ml; Terbufos [CAS:13071-79-9] 100µg/ml; Acephate [CAS:30560-19-1] 100µg/ml; Azinphos-ethyl [CAS:2642-71-9] 400µg/ml; Azinphos-methyl [CAS:86-50-0] 400µg/ml; Demeton-S-methyl [CAS:919-86-8] 100µg/ml; Ethion [CAS:563-12-2] 20µg/ml; Fenamiphos [CAS:22224-92-6] 50µg/ml; Fenitrothion [CAS:122-14-5] 50µg/ml; Fonofos [CAS:944-22-9] 40µg/ml; Metalaxyl [CAS:57837-19-1] 600µg/ml; Methamidophos [CAS:10265-92-6] 100µg/ml; Methidathion [CAS:950-37-8] 100µg/ml; Mevinphos [CAS:7786-34-7] 100µg/ml; Monocrotophos [CAS:6923-22-4] 100µg/ml; Oxa dixyl [CAS:77732-09-3] 400µg/ml; Phorate [CAS:298-02-2] 50µg/ml; Phosphamidon [CAS:13171-21-6] 100µg/ml; Pirimiphos-ethyl [CAS:23505-41-1] 50µg/ml; Triazophos [CAS:24017-47-8] 100µg/ml; Tefluthrin [CAS:79538-32-2] 10µg/ml	1 ml	506905	546

Product specifications are subject to changes. Please visit our website for updates.

Organic standard: Pesticide mixture	29 component 20µg/ml each in toluene/acetone: Aldrin [CAS:309-00-2]; Alpha-HCH [CAS:319-84-6]; Beta-HCH [CAS:319-85-7]; Delta-HCH [CAS:319-86-8]; cis-Chlordane (alpha-Chlordane) [CAS:5103-71-9]; Dieldrin [CAS:60-57-1]; Endosulfan-alpha [CAS:959-98-8]; Endos ulfan-beta [CAS:33213-65-9]; Endosulfan-total (sulfate) [CAS:1031-07-8]; Endrin [CAS:72-20-8]; Endrin aldehyde [CAS:7421-93-4]; Endrin ketone [CAS:53494-70-5]; Gamma-HCH (Lindane) [CAS:58-89-9]; trans-Chlordane (Gamma-Chlordane) [CAS:5103-74-2]; Heptachlor [CAS:76-44-8]; Heptachlor-exo-epoxide (cis-Heptachlorepoxide (cis-, exo-,) [CAS:1024-57-3]; Methoxychlor (DMTD) [CAS:72-43-5]; 4,4'-DDD (TDE) [CAS:72-54-8]; 4,4'-DDE [CAS:72-55-9]; 4,4'-DDT [CAS:50-29-3]; Dicofol [CAS:115-32-2]; Nitrofen [CAS:1836-75-5]; Isodrin [CAS:465-73-6]; Alachlor [CAS:15972-60-8]; Hexachlorobenzene (HCB) [CAS:118-74-1]; 2,4'-DDE [CAS:3424-82-6]; 2,4'-DDD [CAS:53-19-0]; 2,4'-DDT [CAS:789-02-6]; oxy-Chlordane [CAS:27304-13-8]; trans-Nonachlor [CAS:39765-80-5]	1 ml	506948	546
Organic standard: Pesticide mixture	17 components 20µg/ml each in toluene/acetone: Cyfluthrin [CAS:68359-37-5]; Cypermethrin [CAS:52315-07-8]; Fenvalerate [CAS:51630-58-1]; Permethrin [CAS:52645-53-1]; Phenothrin [CAS:26002-80-2]; Tetramethrin [CAS:7696-12-0]; lambda-Cyhalothrin [CAS:91465-08-6]; Piperonyl butoxide [CAS:51-03-6]; Bifenthrin [CAS:82657-04-3]; Chlorothalonil [CAS:1897-45-6]; Quintozene [CAS:82-68-8]; Tecnazene [CAS:117-18-0]; Chlorobenzilate [CAS:510-15-6]; Vinclozolin [CAS:50471-44-8]; Chlordecone hydrate [CAS:143-50-0]; Captan [CAS:133-06-2]	1 ml	506950	546
Organic standard: Pesticide mixture	12 components 10µg/ml each in acetonitrile: Azoxystrobin [CAS:131860-33-8]; Boscalid [CAS:188425-85-6]; Carbendazim [CAS:10605-21-7]; Chlorpyrifos [CAS:2921-88-2]; Cyprodinil [CAS:121552-61-2]; Linuron [CAS:330-55-2]; Metalaxyl [CAS:57837-19-1]; Methomyl [CAS:16752-77-5]; Myclobutanil [CAS:88671-89-0]; Pyrimethanil [CAS:53112-28-0]; Pirimicarb [CAS:23103-98-2]; Thiabendazole [CAS:148-79-8]	10 ml	506803	546

Standard Solutions for PCBs analysis

Description	Notes	Size	Code	Page
Organic standard PCB	PCB 29 [15862-07-4] 10µg/ml in isooctane	1 ml	507127	540
Organic standard PCB	PCB 30 [35693-92-6] 10µg/ml in isooctane	1 ml	507128	540
Organic standard PCB	PCB 73 [74338-23-1] 10µg/ml in isooctane	1 ml	507129	540
Organic standard PCB	PCB 89 [73575-57-2] 10µg/ml in isooctane	1 ml	507131	540
Organic standard PCB	PCB 90 [68194-07-0] 10µg/ml in isooctane	1 ml	507132	540
Organic standard PCB	PCB 106 [70424-69-0] 10µg/ml in isooctane	1 ml	507133	540
Organic standard PCB	PCB 164 [74472-45-0] 10µg/ml in isooctane	1 ml	507134	540
Organic standard PCB	PCB 143 [68194-15-0] 10µg/ml in isooctane	1 ml	507135	540
Organic standard PCB	PCB 155[33979-03-2] 10µg/ml in isooctane	1 ml	507136	540
Organic standard PCB	PCB 198 [68194-17-2] 10µg/ml in isooctane	1 ml	507137	540
Organic standard PCB	PCB 207 [52663-79-3] 10µg/ml in isooctane	1 ml	507138	540
Organic standard PCB	PCB 209 [2051-24-3] 10µg/ml in isooctane	1 ml	507139	540
Organic standard PCB	PCB 209 [2051-24-3] 100µg/ml in isooctane	1 ml	507154	540
Organic Standard: PCB multielement mixture	7 components 10 µg/ml each in isooctane: PCB 28 [CAS:7012-37-5]; PCB 52 [CAS:35693-99-3]; PCB 101 [CAS:37680-73-2]; PCB 118 [CAS:31508-00-6]; PCB 138 [CAS:35065-28-2]; PCB 153 [CAS:35065-27-1]; PCB 180 [CAS:35065-29-3]	5 x 1 ml	507103	543

Organic Standard: PCB multielement mixture	2 components 100 µg/ml each in isooctane: PCB 30 [CAS:35693-92-6]; PCB 155 [CAS:33979-03-2]	1 ml	507609	543
Organic Standard: PCB multielement mixture	8 components 100 µg/ml each in hexane: PCB 28 [CAS:7012-37-5]; PCB 52 [CAS:35693-99-3]; PCB 101 [CAS:37680-73-2]; PCB 118 [CAS:31508-00-6]; PCB 138 [CAS:35065-28-2]; PCB 153 [CAS:35065-27-1]; PCB 180 [CAS:35065-29-3]; PCB 194 [CAS:35694-08-7]	1 ml	507679	543
Organic Standard: PCB multielement mixture	14 components 10 mg/l each in hexane: PCB 28 [CAS:7012-37-5]; PCB 31 [CAS:16606-02-3]; PCB 52 [CAS:35693-99-3]; PCB 101 [CAS:37680-73-2]; PCB 105 [CAS:32598-14-4]; PCB 118 [CAS:31508-00-6]; PCB 132 [CAS:38380-05-1]; PCB 138 [CAS:35065-28-2]; PCB 149 [CAS:38380-04-0]; PCB 153 [CAS:35065-27-1]; PCB 160 [CAS:41411-62-5]; PCB 163 [CAS:74472-44-9]; PCB 180 [CAS:35065-29-3]; PCB 193 [CAS:69782-91-8]	1.2 ml	507062	543
Organic Standard: PCB multielement mixture	19 components in ethyle acetate: PCB 18 [CAS:37680-65-2] 0.34mg/l; PCB 28 [CAS:7012-37-5] 0.6mg/l; PCB 52 [CAS:35693-99-3] 0.9mg/l; PCB 77 [CAS:32598-13-3] 2.2mg/l; PCB 81 [CAS:70362-50-4] 3.3mg/l; PCB 101 [CAS:37680-73-2] 2.28mg/l; PCB 105 [CAS:32598-14-4] 2.6mg/l; PCB 114 [CAS:74472-37-0] 9.6mg/l; PCB 118 [CAS:31508-00-6] 2.6mg/l; PCB 123 [CAS:65510-44-3] 2.7mg/l; PCB 126 [CAS:57465-28-8] 3mg/l; PCB 138 [CAS:35065-28-2] 6mg/l; PCB 153 [CAS:35065-27-1] 5mg/l; PCB 156 [CAS:38380-08-4] 5mg/l; PCB 157 [CAS:69782-90-7] 7mg/l; PCB 167 [CAS:52663-72-6] 8mg/l; PCB 169 [CAS:32774-16-6] 10mg/l; PCB 180 [CAS:35065-29-3] 10mg/l; PCB 189 [CAS:39635-31-9] 7mg/l	5 ml	506732	543
Organic Standard: PCB multielement mixture	14 components 1 mg/l each in hexane: PCB 28 [CAS:7012-37-5]; PCB 52 [CAS:35693-99-3]; PCB 101 [CAS:37680-73-2]; PCB 153 [CAS:35065-27-1]; PCB 138 [CAS:35065-28-2]; PCB 180 [CAS:35065-29-3] PCB 194 [CAS:35694-08-7]; 1,2,4-Trichlorobenzene [CAS:120-82-1]; Hexachloro-1, 3-butadiene [CAS:87-68-3]; Hexachlorobenzene [CAS:118-74-1]; Alpha-HCH [CAS:319-84-6]; Beta-HCH [CAS:319-85-7]; Gamma-HCH (Lindane) [CAS:58-89-9]; Delta-HCH [CAS:319-86-8]	5 ml	507889	543
Organic Standard: PCB multielement mixture	14 components 10 µg/ml each in isooctane according to EN 61619: PCB 18 [CAS:37680-65-2]; PCB 28 [CAS:7012-37-5]; PCB 31 [CAS:16606-02-3]; PCB 44 [CAS:41464-39-5]; PCB 52 [CAS:35693-99-3]; PCB 101 [CAS:37680-73-2]; PCB 118 [CAS:31508-00-6]; PCB 138 [CAS:35065-28-2]; PCB 149 [CAS:38380-04-0]; PCB 153 [CAS:35065-27-1]; PCB 170 [CAS:35065-30-6]; PCB 180 [CAS:35065-29-3]; PCB 194 [CAS:35694-08-7]; PCB 209 [CAS:2051-24-3]	10 ml	507115	543
Organic Standard: PCB multielement mixture	PCB 138 [35065-28-2] 50µg/ml in isooctane	10 ml	509144	543
Organic Standard: PCB multielement mixture	PCB 153 [35065-27-1] 50µg/ml in isooctane	10 ml	509145	543
Organic Standard: PCB multielement mixture	PCB 18 [37680-65-2] 50µg/ml in isooctane	10 ml	509146	543
Organic Standard: PCB multielement mixture	PCB 28 [7012-37-5] 50µg/ml in isooctane	10 ml	509147	543
Organic Standard: PCB multielement mixture	PCB 52 [35693-99-3] 50µg/ml in isooctane	10 ml	509148	543

Product specifications are subject to changes. Please visit our website for updates.

Standard Solutions for PAHs analysis

Description	Notes	Size	Code	Page
Organic Standard: PAH multielement mixture	15 components 50ug/ml each in acetonitrile: Acenaphthene [CAS:83-32-9]; Anthracene [CAS:120-12-7]; Benzo(a)anthracene [CAS:56-55-3]; Benzo(a)pyrene [CAS:50-32-8]; Benzo(b) fluoranthene [CAS:205-99-2]; Benzo(g,h,i) perylene [CAS:191-24-2]; Benzo(k)fluoranthene [CAS:207-08-9]; Chrysene [CAS:218-01-9]; Dibenzo(a,h)anthracene [CAS:53-70-3]; Fluoranthene [CAS:206-44-0]; Fluorene [CAS:86-73-7]; Indeno(1,2,3-c,d)pyrene [CAS:193-39-5]; Naphthalene [CAS:91-20-3]; Phenanthrene [CAS:85-01-8]; Pyrene [CAS:129-00-0]	5 x 1 ml	507063	542
Organic Standard: PAH multielement mixture	13 components in acetonitrile/ acetone 95/5: Phenanthrene [CAS:85-01-8] 600µg/ml; Anthracene [CAS:120-12-7] 40µg/ml; Fluoranthene [CAS:206-44-0] 160µg/ml; Pyrene [CAS:129-00-0] 160µg/ml; Benzo(a) anthracene [CAS:56-55-3] 20µg/ml; Chrysene [CAS:218-01-9] 80µg/ml; (95/5) Benzo(b) fluoranthene [CAS:205-99-2] 20µg/ml; Benzo(k)fluoranthene [CAS:207-08-9] 10µg/ml; Benzo(j)fluoranthene [CAS:205-82-3] 20µg/ml; Benzo(a)pyrene [CAS:50-32-8] 20µg/ml; Dibenzo(a,h)anthracene [CAS:53-70-3] 10µg/ml; Indeno(1,2,3-c,d)pyrene [CAS:193-39-5] 12µg/ml; Benzo(g,h,i)perylene [CAS:191-24-2] 20µg/ml	10 x 1 ml	506835	542
Organic Standard: PAH multielement mixture	19 components 100 µg/ml each in acetonitrile: Naphthalene [CAS:91-20-3]; Acenaphthylene [CAS:208-96-8]; 1-Methylnaphthalene [CAS:90-12-0]; 2-Methylnaphthalene [CAS:91-57-6]; Acenaphthene [CAS:83-32-9]; Fluorene [CAS:86-73-7]; Phenanthrene [CAS:85-01-8]; Anthracene [CAS:120-12-7]; Fluoranthene [CAS:206-44-0]; Pyrene [CAS:129-00-0]; 2-Methyl-Fluoranthene [CAS:33543-31-6]; Benzo(a)anthracene [CAS:56-55-3;] Chrysene [CAS:218-01-9]; Benzo(b)fluoranthene [CAS:205-99-2]; Benzo(k)fluoranthene [CAS:207-08-9]; Benzo(a)pyrene [CAS:50-32-8]; Benzo(g,h,i)perylene [CAS:191-24-2]; Indeno(1,2,3-c,d)pyrene [CAS:193-39-5]; Dibenzo(a,h)anthracene [CAS:53-70-3]	1 ml	506878	542
Organic Standard: PAH multielement mixture	23 components 1 µg/ml each in methanol: Anthracene [CAS:120-12-7]; Fluoranthene [CAS:206-44-0]; Naphthalene [CAS:91-20-3]; Benzo(a)pyrene [CAS:50-32-8]; Benzo(b) fluoranthene [CAS:205-99-2]; Benzo(g,h,i) perylene [CAS:191-24-2]; Benzo(k)fluoranthene [CAS:207-08-9]; Indeno(1,2,3-c,d)pyrene [CAS:193-39-5]; Biphenyl [CAS:92-52-4]; Acenaphthene [CAS:83-32-9]; Acenaphthylene [CAS:208-96-8]; Ben zo(a)anthracene [CAS:56-55-3]; Chrysene [CAS:218-01-9]; Dibenzo(a,h) anthracene [CAS:53-70-3]; Fluorene [CAS:86-73-7]; 2-Methylnaphthalene [CAS:91-57-6]; 2-Methyl-Fluoranthene [CAS:33543-31-6]; Phenanthrene [CAS:85-01-8]; Pyrene [CAS:129-00-0]; 1-Benzothiophene [CAS:95-15-8]; Dibenzothiophene [CAS:132-65-0]; Benzo(e)pyrene [CAS:192-97-2]; Perylene [CAS:198-55-0]	1 ml	506938	542

Organic Standard: PAH multielement mixture	6 components in acetonitrile: Fluoranthene [CAS:206-44-0] 2mg/l; Benzo(b)fluoranthene [CAS:205-99-2] 2mg/l; Benzo(k)fluoranthene [CAS:207-08-9] 2mg/l; Benzo(a)pyrene [CAS:50-32-8] 2mg/l; Benzo(g,h,i)perylene [CAS:191-24-2] 10mg/l; Indeno(1,2,3-c,d)pyrene [CAS:193-39-5] 10mg/l	1 ml	506979	542
Organic Standard: PAH multielement mixture	15 components in acetonitrile: Acenaphthene [CAS:83-32-9] 5mg/l; Fluorene [CAS:86-73-7] 5mg/l; Fluoranthene [CAS:206-44-0] 5mg/l; Benzo(a)anthracene [CAS:56-55-3] 5mg/l; Chrysene [CAS:218-01-9] 5mg/l; Benzo(b)fluoranthene [CAS:205-99-2] 5mg/l; Benzo(a)pyrene [CAS:50-32-8] 5mg/l; Dibenzo(a,h)anthracene [CAS:53-70-3] 5mg/l; Benzo(g,h,i)perylene [CAS:191-24-2] 10mg/l; Indeno(1,2,3-c,d)pyrene [CAS:193-39-5]10mg/l; 2-Methylnaphthalene [CAS:91-57-6] 10mg/l; 2-Methyl-Fluoranthene [CAS:33543-31-6] 10mg/l; Anthracene [CAS:120-12-7] 2mg/l; Benzo(k)fluoranthene [CAS:207-08-9] 2mg/l; Pyrene [CAS:129-00-0] 20mg/l	1 ml	506980	542
Organic Standard: PAH multielement mixture	19 components in Methanol: 2-Methylnaphthalene [CAS:91-57-6] 40mg/l; Anthracene [CAS:120-12-7] 20mg/l; Fluoranthene [CAS:206-44-0] 20mg/l; 2-Methyl-Fluoranthene [CAS:33543316] 20mg/l Benzo(a) anthracene [CAS:56-55-3] 20mg/l; Benzo(b) fluoranthene [CAS:205-99-2] 20mg/l Benzo(k) fluoranthene [CAS:207-08-9] 20mg/l; Benzo(a) pyrene [CAS:50-32-8] 20mg/l; Dibenzo(a,h) anthracene [CAS:53-70-3] 20mg/l; Benzo(g,h,i) perylene [CAS:191-24-2] 20mg/l; Indeno(1,2,3-c,d)pyrene [CAS:193-39-5] 20mg/l; Benzo(b) chrysene [CAS:214-17-5] 2mg/l; Naphthalene [CAS:91-20-3] 40mg/l; Acenaphthene [CAS:83-32-9] 40mg/l; Fluorene [CAS:86-73-7] 20mg/l; Phenanthrene [CAS:85-01-8] 20mg/l; Pyrene [CAS:129-00-0] 20mg/l; Chrysene [CAS:218-01-9] 20mg/l; Acenaphthylene [CAS:208-96-8] 400mg/l	1 ml	507094	542
Organic Standard: PAH multielement mixture	15 components 1 µg/ml each in dichloromethane: Naphthalene [CAS:91-20-3]; Acenaphthene [CAS:83-32-9]; Fluorene [CAS:86-73-7]; Phenanthrene [CAS:85-01-8]; Anthracene [CAS:120-12-7]; Fluoranthene [CAS:206-44-0]; Pyrene [CAS:129-00-0]; Benzo(a)anthracene [CAS:56-55-3]; Chrysene [CAS:218-01-9] Benzo(b)fluoranthene [CAS:205-99-2]; Benzo(k)fluoranthene [CAS:207-08-9]; Benzo(a)pyrene [CAS:50-32-8]; Dibenzo(a,h)anthracene [CAS:53-70-3]; Benzo(g,h,i)perylene [CAS:191-24-2]; Indeno(1,2,3-c,d)pyrene [CAS:193-39-5]	1.5 ml	507859	542
Organic Standard: PAH multielement mixture	15 components 1 µg/ml each in dichloromethane: Naphthalene [CAS:91-20-3]; Acenaphthene [CAS:83-32-9]; Fluorene [CAS:86-73-7]; Phenanthrene [CAS:85-01-8]; Anthracene [CAS:120-12-7]; Fluoranthene [CAS:206-44-0]; Pyrene [CAS:129-00-0]; Benzo(a)anthracene [CAS:56-55-3]; Chrysene [CAS:218-01-9] Benzo(b)fluoranthene [CAS:205-99-2]; Benzo(k)fluoranthene [CAS:207-08-9]; Benzo(a)pyrene [CAS:50-32-8]; Dibenzo(a,h)anthracene [CAS:53-70-3]; Benzo(g,h,i)perylene [CAS:191-24-2]; Indeno(1,2,3-c,d)pyrene [CAS:193-39-5]	1.5 ml	507899	542

Organic Standard: PAH multielement mixture	19 components 10mg/l each in methanol: Naphthalene [CAS:91-20-3]; Acenaphthylene [CAS:208-96-8]; Acenaphthene [CAS:83-32-9]; Fluorene [CAS:86-73-7]; Phenanthrene [CAS:85-01-8]; Anthracene [CAS:120-12-7]; Fluoranthene [CAS:206-44-0]; Pyrene [CAS:129-00-0]; Benzo(a)anthracene [CAS:56-55-3]; Chrysene [CAS:218-01-9]; Benzo(b) fluoranthene [CAS:205-99-2]; Benzo(k) fluoranthene [CAS:207-08-9]; Benzo(a) pyrene [CAS:50-32-8]; Dibenzo(a,h)anthracene [CAS:53-70-3]; Benzo(g,h,i)perylene [CAS:191-24-2]; Indeno(1,2,3-c,d)pyrene [CAS:193-39-5]; Benzo(e)pyrene [CAS:192-97-2]; 2-Methylnaphthalene [CAS:91-57-6]; 2-Methyl-Fluoranthene [C AS:33543-31-6]	10 ml	506821	542

Standard Solutions for Hydrocarbons analysis

Description	Notes	Size	Code	Page
Organic standard: Mixture for hydrocarbon analysis	5 components 5000µg/ml each in methanol: Benzene; Toluene; o-Xylene; m-Xylene; p-Xylene	1 ml	506736	546
Organic standard: Mixture for hydrocarbon analysis	14 components in methanol: 1,1-Dichloroethene 1000µg/ml; Dichloromethane 5000µg/ml; trans-1,2-Dichloroethene 5000µg/ml; 1,1-Dichloroethane 5000µg/ml; cis-1,2-Dichloroethene 5000µg/ml; 1,2-Dichloroethane 5000µg/ml; Chloroform 500µg/ml; 1,1,1-Trichloroethane 500µg/ml; Trichloroethene 500µg/ml; Bromodichloromethane 500µg/ml; Dibromochloromethane 500µg/ml;Tribromomethane 500µg/ml; Tetrachloromethane 100µg/ml; Tetrachloroethene 100µg/ml; Hydrocarbons Mixture Benzene; 1-Bromo-2-chloroethane; Chlorobenzene; 1,1-Dichloroethane; 1,2-dichloroethane	1 ml	506742	546
Organic standard: Mixture for hydrocarbon analysis	6 components 1000 µg/ml each in methanol: Benzene; Toluene; o-Xylene; m-Xylene; p-Xylene; Ethylbenzene	1 ml	507189	546
Organic standard: Mixture for hydrocarbon analysis	6 components 1000 µg/ml each in methanol: Trichloroethene; Tetrachloroethene; 1,2-Dichloroethane; Tetrachloromethane; 1,1,1-Trich loroethane; 1,1,2-Trichloroethane	1 ml	507190	546
Organic standard: Mixture for hydrocarbon analysis	4 components 1000 µg/ml each in methanol: Tribromomethane; Chloroform; Bromodichloromethane; Dibromochloromethane	1 ml	507191	546
Organic standard: Mixture for hydrocarbon analysis	6 components 1000 µg/ml each in carbon disulfide: Benzene; Ethylbenzene; Toluene; m-Xylene; o-Xylene; p-Xylene	2 ml	507474	546
Organic standard: Mixture for hydrocarbon analysis	22 components 2500mg/Kg each 1-Bromo-2-chloroethane; Chlorobenzene; 1,1-Dichloroethane; 1,1-Dichloroethene; cis-1,2-Dichloroethene; trans-1,2-Dichloroethene; Dichloromethane; Pentachloroethane; 1,1,2,2-Tetrachloroethane; Tetrachloroethene; Tetrachloromethane; 1,1,2-Trichloroethane; Trichloroethene: Chloroprene; Chloromethane; Vinylchloride; 1,3-Butadiene; Chloroethane; 1,2-Dichlorobutane; Ethylene; Chloroform; Matrix: Benzene	100 ml	506614	546

Standard Solutions for Method ISO 9377-2

For the determination of mineral oils according to UNI ISO 9377-2 "Determination of hydrocarbon oil index - Method using solvent extraction and gas chromatography", the current regulations impose specific mixtures of standard solutions. CARLO ERBA Reagents offers a complete range of such products, each supplied with an analysis certificate containing complete information on the composition and gravimetric verification performed in reference to NIST standards.

Description	Notes	Size	Code	Page
Organic standard: Mixture for hydrocarbon analysis	Standard quality control of two mineral oils in acetone	1 ml	506002	547
Organic standard: Mixture for hydrocarbon analysis	Mixture of mineral oil without additive 2 to 5 mg / ml each in hexane	1 ml	506010	547
Organic standard: Mixture for hydrocarbon analysis	Standard mixture of n-alkanes (C10 to C40 in pairs) of 50 mcg / ml each in hexane	1 ml	506020	547
Organic standard: Mixture for hydrocarbon analysis	2 component: 5g/l each of Mineral Oil [CAS:8042-47-5]; Diesel Oil in n-Hexane	5 ml	506012	547
Organic standard: Mixture for hydrocarbon analysis	N-tetracontane mixture (20 mg / l) and n-decane (20 mg / l) in hexane	5 ml	506040	547
Organic standard: Mixture for hydrocarbon analysis	Mixture of mineral oil without additive 2 to 1 mg / ml each in hexane	10 ml	506011	547
Organic standard: Mixture for hydrocarbon analysis	2 components: 5g/l each of Mineral Oil [CAS:8042-47-5]; Diesel Oil in n-Hexane	10 ml	506013	547
Organic standard: Mixture for hydrocarbon analysis	Standard mixture of n-alkanes (C10 to C40 in pairs) of 50 µg / ml each in hexane	10 ml	506021	547
Organic standard: Mixture for hydrocarbon analysis	Mother solution stearyl stearate 2 g / l in hexane	10 ml	506030	547

COD ANALYSIS

The COD (chemical oxygen demand) is the quantity of oxygen needed for the complete chemical oxidation of organic and inorganic compounds present in a water samples, especially for the control of pollutions levels. CARLO ERBA Reagents offers the strong oxidants and acids used for COD determination.

Reagents for COD Analysis

Description	Notes	Size	Code	Page
Ferroin 0.025 mol/l solution		100 ml	526751	354
Iron (II) ammonium sulfate 0.12N		1 l	526761	429
o-Phenantroline-Iron (II) sulphate solution in sulphuric acid		100 ml	E450043	569
Potassium dichromate 0.0417 mol/l (0.25N)		1 l	470451	593
Potassium dichromate 0.04 mol/l (0.24N) in 80 g/l HgSO4		1 l	526711	593
Potassium dichromate 0.04 mol/l (0.24N) in 80 g/l HgSO4		2.5 l	526712	593
Silver sulfate solution 0.7% in sulfuric acid		1 l	424191	659
Silver sulfate solution 0.7% in sulfuric acid		2.5 l	424192	659
Sulfuric acid with 10 g/l Silver sulfate		1 l	526605	742
Sulfuric acid with 10 g/l Silver sulfate		2.5 l	526606	742
Sulfuric acid with 6.6 g/l Silver sulfate		2.5 l	526602	743

PHARMACEUTICAL PRODUCTION

CARLO ERBA Reagents has always carefully monitored developments in the quality of raw materials and reagents used in the pharmaceutical industry. The experience gained in over 160 years as a supplier of raw materials and excipients and the countless approvals issued by the most important pharmaceutical multinationals have allowed CARLO ERBA Reagents to become a benchmark in the pharmaceutical industry worldwide, distinguishing itself for its premium quality and reliability.

Xcipharm™: excipient use

Xcipharm™ is the new product range dedicated to excipients. We manufacture these products according to IPEC guidelines assuring a high level of quality for your excipients. We can adapt our offer to your needs with SKU sizes starting at 100 ml up to bulk quantity. We provide all of the necessary tests, certification and statements to make the registration of your products to the health authority easier.

QUALITY ASSURANCE
- Registration to the French Minister of Health (ANSM)
- Traceability of raw material and packaging
- Flow chart
- Change control

PROCESS
- Validated cleaning procedures and or dedicated equipment
- Labels consolidation
- Certificate of analysis with the name of the producer and the date of production of the raw material
- Stability studies
- Sample library of the raw material (1 year) and of the finish product (shelf life + 1 year)

DOCUMENTATION
- BSE/TSE statement
- GMO statement
- Residual solvents statement ICH Q3C
- Elemental Impurities ICH Q3D
- Risk assessment (2015/C95/02)

ERBApharm®: raw material use

CARLO ERBA Reagents has developed the ERBApharm® product line on the specific requests of the pharmaceutical market. Their specifications comply with the effective requirements of Pharmacopoeia or - in the absence of those requirements - with strict sales specifications. These products are designed to be used as raw materials and reagents for the production of active principles.
The ERBApharm® line includes every kind of product: solvents, acids and basis, titrated solutions, organic and inorganic powders.
The documents available for these products comply with the need of information related to their use (residual solvents, BSE/TSE declaration and GMO).

Description	Notes	Size	Code	Page
Acetic acid glacial		1 l	302016	133
Acetic acid glacial		2.5 l	302011	133
Acetic acid glacial		5 l	302014	133
Acetic acid glacial		30 kg	302015	133
Acetic acid glacial		200 kg	302013	133
Acetone		1 l	301505	140
Acetone		2.5 l	301506	140
Acetone		5 l	301502	140
Acetone		5 l	301503	140
Acetone		16 kg	301501	140
Acetone		22 kg	301504	140
Acetone		160 kg	301507	140
Aluminum chloride hexahydrate		1 kg	311257	153
Aluminum chloride hexahydrate		5 kg	311252	153
Aluminum chloride hexahydrate		25 kg	311256	153
Aluminum chloride hexahydrate		50 kg	311254	153
Aluminum potassium sulfate dodecahydrate		1 kg	312401	156
Aluminum potassium sulfate dodecahydrate		5 kg	312404	156
Aluminum potassium sulfate dodecahydrate		10 kg	312402	156

p-Aminobenzoic acid	100 g	391804	158
p-Aminobenzoic acid	1 kg	391805	158
Ammonia solution 28%	1 l	314861	161
Ammonia solution 28%	2 l	314863	161
Ammonia solution 28%	25 kg	314866	161
Ammonium chloride	1 kg	313957	169
Ammonium chloride	2.5 kg	313952	169
Ammonium chloride	5 kg	313956	169
Ammonium chloride	25 kg	313951	169
Ammonium chloride	50 kg	313954	169
Benzalkonium chloride	1 kg	322737	197
Benzalkonium chloride	5 kg	322738	197
Benzoic acid	1 kg	302087	199
Benzoic acid	5 kg	302089	199
Benzoic acid	25 kg	302082	199
Benzyl alcohol	1 l	308131	200
Benzyl alcohol	2.5 l	308132	200
Benzyl alcohol	23 kg	308138	200
Benzyl alcohol	200 l	308137	200
Benzyl benzoate	1 l	323101	201
Benzyl benzoate	2.5 l	323102	201
Boric acid	25 kg	302185	205
Boric acid	1 kg	302177	206
Boric acid	5 kg	302179	206
Boric acid	25 kg	302178	206
Caffeine anhydrous	500 g	326356	234
Caffeine anhydrous	1 kg	326357	234
Caffeine anhydrous	25 kg	326358	234
Calcium acetate anhydrous	1 kg	326511	236
Calcium acetate anhydrous	5 kg	326512	236
Calcium acetate anhydrous	25 kg	326513	236
Calcium carbonate	1 kg	327101	237
Calcium carbonate	25 kg	327105	237
Calcium chloride dihydrate	1 kg	327607	239
Calcium chloride dihydrate	5 kg	327609	239
Calcium chloride dihydrate	25 kg	327603	239
Calcium chloride hexahydrate	1 kg	327507	239
Calcium chloride hexahydrate	5 kg	327509	239
Calcium gluconate	1 kg	330608	240
Calcium gluconate	5 kg	330609	240
Calcium gluconate	25 kg	330601	240
Calcium hydroxide	1 kg	331007	241
Calcium hydroxide	5 kg	331008	241
Calcium hydroxide	25 kg	331003	241
Calcium lactate	1 kg	331407	241
Calcium lactate	5 kg	331408	241
Calcium pantothenate	100 g	331602	242
Calcium phosphate dibasic dihydrate	1 kg	330307	243
Calcium phosphate dibasic dihydrate	25 kg	330303	243
Calcium phosphate tribasic	1 kg	330407	244

Calcium phosphate tribasic	5 kg	330409	244
Calcium phosphate tribasic	25 kg	330403	244
Calcium stearate	2.5 kg	332262	244
Calcium stearate	10 kg	332261	244
Calcium stearate	25 kg	332265	244
Calcium sulfate dihydrate	5 kg	331752	244
Calcium sulfate dihydrate	25 kg	331751	244
Camphor natural	500 g	332356	246
Camphor synthetic	500 g	332406	246
Camphor synthetic	5 kg	332401	246
Camphor synthetic	25 kg	332402	246
Castor oil	1 l	356351	249
Castor oil	5 l	356352	249
Castor oil	28 kg	356353	249
Cetyl alcohol	1 kg	308357	255
Cetyl alcohol	5 kg	308358	255
Cetyl alcohol	25 kg	308359	255
Chlorobutanol hemihydrate	1 kg	301357	259
Chlorobutanol hemihydrate	5 kg	301356	259
Chloroform	1 l	334351	262
Chloroform	2.5 l	334353	262
Chloroform	25 kg	334356	262
Chloroform	200 l	529301	262
Chloroform	250 kg	334354	262
Citric acid anhydrous	500 g	302486	269
Citric acid anhydrous	1 kg	302487	269
Citric acid anhydrous	5 kg	302485	269
Citric acid anhydrous	25 kg	302488	269
Citric acid anhydrous	50 kg	302484	269
Citric acid monohydrate	1 kg	302557	269
Citric acid monohydrate	5 kg	302559	269
Citric acid monohydrate	25 kg	302551	269
Citric acid monohydrate	50 kg	302554	269
Citric acid monohydrate	1 kg	302507	270
Citric acid monohydrate	5 kg	302509	270
Citric acid monohydrate	25 kg	302501	270
Citric acid monohydrate	50 kg	302504	270
Copper (II) sulfate pentahydrate	1 kg	364757	281
Copper (II) sulfate pentahydrate	5 kg	364759	281
Copper (II) sulfate pentahydrate	25 kg	364752	281
Di-n-butylphthalate	26 kg	325701	297
Dichloromethane, stab. with Ethanol	1 l	354501	303
Dichloromethane, stab. with Amylene	1 l	337331	304
Dichloromethane, stab. with Amylene	2.5 l	337333	304
Dichloromethane, stab. with Amylene	25 l	337335	304
Dichloromethane, stab. with Amylene	200 l	337337	304
Dichloromethane, stab. with Ethanol	2.5 l	525320	304
Dichloromethane, stab. with Ethanol	200 l	525321	304
Diethanolamine	215 kg	337801	306
Diethyl ether	1 l	340751	309

Product specifications are subject to changes. Please visit our website for updates.

Diethyl ether		40 x 100 g	340731	309
Diethyl ether		20 kg	340752	309
Diethyl ether		140 kg	340759	309
Diethyl phthalate		1 l	338112	310
Diethyl phthalate		2.5 l	338115	310
Diethyl phthalate		30 kg	338113	310
Diethyl phthalate		200 l	338114	310
Ethanol absolute anhydrous		1 l	529121	333
Ethanol absolute anhydrous	Only for Italian market	2.5 l	308661	333
Ethanol absolute anhydrous	Only for Italian market	2.5 l	308662	333
Ethanol absolute anhydrous		2.5 l	3086612	333
Ethanol absolute anhydrous		2.5 l	3086622	333
Ethanol absolute anhydrous		5 l	529122	333
Ethanol absolute anhydrous	Untaxed, for Italian license holders only	5 l	5291222	333
Ethanol absolute anhydrous		10 l	529124	333
Ethanol absolute anhydrous		25 l	308664	333
Ethanol absolute anhydrous		25 l	308667	333
Ethanol absolute anhydrous	Untaxed, for italian license holders only	25 l	3086642	333
Ethanol absolute anhydrous		200 l	308663	333
Ethanol absolute anhydrous		200 l	529125	333
Ethanol 96°		5 l	524135	335
Ethanol 96°		25 l	524132	335
Ethanol 96°	Only for Italian market	1 l	308644	336
Ethanol 96°	Only for Italian market	1 l	308647	336
Ethanol 96°		1 l	3086442	336
Ethanol 96°		1 l	3086472	336
Ethanol 96°	Only for Italian market	2.5 l	308641	336
Ethanol 96°	Only for Italian market	2.5 l	308649	336
Ethanol 96°		2.5 l	3086412	336
Ethanol 96°		2.5 l	3086492	336
Ethanol 96°		5 l	529141	336
Ethanol 96°	Untaxed, for Italian license holders only	5 l	5291412	336
Ethanol 96°		10 l	308646	336
Ethanol 96°	Untaxed, for italian license holders only	10 l	3086462	336
Ethanol 96°		25 l	308645	336
Ethanol 96°	Untaxed, for italian license holders only	25 l	3086452	336
Ethanol 96°		27 l	308648	336
Ethanol 96°	Untaxed, for italian license holders only	27 l	3086482	336
Ethanol 96°		200 l	308643	336
Ethanol 70% v/v	Sold by box: 20 bottles + 4 sprayers	20 x 500 ml	529187000	336
Ethanol 70% v/v		5 l	529189	336
Ethanol 70% v/v		200 l	529183	336
Ethanol 70% v/v	Sold by box 6 bottles + 3 sprayers	1 l	529184	337
Ethanol 70% v/v		1 l	529184000	337
Ethyl acetate		1 l	341506	342
Ethyl acetate		2.5 l	341503	342
Ethyl acetate		24 kg	341502	342
Ethyl acetate		25 l	529221	342
Ethyl acetate		200 l	529222	342
Ethylenediaminetetraacetic acid		5 kg	303251	344

Ethylenediaminetetraacetic acid	25 kg	303252	344
Ethylenediaminetetraacetic acid disodium salt	1 kg	303201	345
Ethylenediaminetetraacetic acid disodium salt	5 kg	303203	345
Ethylenediaminetetraacetic acid disodium salt	25 kg	303202	345
Ethylenediaminetetraacetic acid disodium salt	1 kg	303227	346
Ethylenediaminetetraacetic acid disodium salt	5 kg	303226	346
Ethylenediaminetetraacetic acid disodium salt	25 kg	303225	346
Formaldehyde 35% w/w	1 l	310351	359
Formaldehyde 35% w/w	2.5 l	310356	359
Formaldehyde 35% w/w	5 l	310358	359
Formaldehyde 35% w/w	10 kg	310349	359
Formaldehyde 35% w/w	30 kg	310348	359
Formaldehyde 35% w/w	55 kg	310355	359
Formic acid 99%	1 l	303911	363
Formic acid 99%	2.5 l	303912	363
Formic acid 99%	30 kg	303913	363
D(+)Glucose anhydrous	1 kg	346987	370
D(+)Glucose anhydrous	5 kg	346989	370
D(+)Glucose anhydrous	25 kg	346983	370
D(+)Glucose monohydrate	1 kg	346971	371
D(+)Glucose monohydrate	5 kg	346972	371
D(+)Glucose monohydrate	25 kg	346973	371
Glycerol (30°Bé)	1 l	346161	373
Glycerol (30°Bé)	2.5 l	346165	373
Glycerol (30°Bé)	5 l	346162	373
Glycerol (30°Bé)	35 kg	346164	373
Glycerol (30°Bé)	250 kg	346167	373
Glycine	1 kg	346207	374
Glycine	5 kg	346205	374
Glycine	25 kg	346208	374
Gum arabic	1 kg	347107	377
Hydrochloric acid 37%	1 l	302621	395
Hydrochloric acid 37%	2.5 l	302626	395
Hydrochloric acid 37%	5 l	302643	395
Hydrochloric acid 37%	10 l	302624	395
Hydrochloric acid 37%	25 kg	302623	395
Hydrochloric acid 37%	40 kg	302622	395
Hydrochloric acid 37%	55 kg	302627	395
Hydrochloric acid 37%	220 kg	302625	395
Hydrochloric acid 10%	10 kg	302591	400
Hydrogen peroxide solution 35%	2.5 l	307742	413
Hydrogen peroxide solution 30%	25 kg	307685	415
Hydrogen peroxide solution 3%	1 l	307671	416
Hydrogen peroxide solution 3%	50 kg	307678	416
Iodine	100 g	348454	424
Iodine	250 g	348455	424
Iodine	1 kg	348457	424
Iodine	5 kg	348451	424
Iodine	20 kg	348452	424
Iron (II) sulfate heptahydrate	1 kg	344957	431

Product specifications are subject to changes. Please visit our website for updates.

Iron (II) sulfate heptahydrate	5 kg	344959	431
L(+)-Lactic acid	1 l	304652	452
L(+)-Lactic acid	2.5 l	304651	452
L(+)-Lactic acid	25 kg	304653	452
Lactose monohydrate	1 kg	348707	453
Lactose monohydrate	5 kg	348708	453
Lactose monohydrate	10 kg	348702	453
Lactose monohydrate	25 kg	348703	453
Lanolin anhydrous	1 kg	347357	453
Lanolin anhydrous	5 kg	347359	453
Magnesium carbonate basic	1 kg	349257	470
Magnesium carbonate basic	5 kg	349279	470
Magnesium carbonate basic	25 kg	349272	470
Magnesium chloride hexahydrate	1 kg	349357	471
Magnesium chloride hexahydrate	5 kg	349359	471
Magnesium chloride hexahydrate	25 kg	349355	471
Magnesium hydroxide	1 kg	349455	472
Magnesium oxide heavy	1 kg	349655	473
Magnesium oxide heavy	5 kg	349656	473
Magnesium oxide heavy	25 kg	349653	473
Magnesium stearate	2.5 kg	350032	474
Magnesium stearate	20 kg	350033	474
Magnesium stearate	25 kg	350035	474
Magnesium sulfate heptahydrate	1 kg	349852	475
Magnesium sulfate heptahydrate	5 kg	349859	475
Magnesium sulfate heptahydrate	25 kg	349851	475
Maize starch	1 kg	313071	476
Maize starch	5 kg	313072	476
Maize starch	25 kg	313073	476
Maleic acid	500 g	407266	477
Maleic acid	5 kg	407261	477
Maleic acid	25 kg	407263	477
D-Mannitol	1 kg	352051	482
D-Mannitol	5 kg	352052	482
D-Mannitol	25 kg	352053	482
L-Menthol	50 g	352103	483
L-Menthol	500 g	352106	483
Methanol	1 l	309204	492
Methanol	2.5 l	309203	492
Methanol	25 l	309201	492
Methanol	200 l	529100	492
Methyl 4-hydroxybenzoate	1 kg	354007	499
Methyl 4-hydroxybenzoate	5 kg	354008	499
Methyl salicylate	1 l	354152	505
Methyl salicylate	25 kg	354155	505
Nicotinamide	100 g	392304	525
Nicotinamide	1 kg	392307	525
Oil refined of almonds	1 l	356251	539
Orthophosphoric acid 85%	1 l	304061	548
Orthophosphoric acid 85%	2.5 l	304062	548

Orthophosphoric acid 85%	40 kg	304063	548
Paraffin oil	1 l	356601	555
Paraffin oil	5 l	356608	555
Paraffin oil	23 kg	356603	555
Paraffin oil	185 kg	356607	555
Paraffin white soft	1 kg	388407	556
Paraffin white soft	5 kg	388409	556
Phenol	1 kg	343407	570
2-Phenylethanol	1 l	529022	573
2-Phenylethanol	2.4 l	529021	573
Potassium acetate	1 kg	358907	583
Potassium acetate	5 kg	358908	583
Potassium acetate	25 kg	358903	583
Potassium bromide	1 kg	359707	586
Potassium bromide	5 kg	359702	586
Potassium chloride	1 kg	360107	587
Potassium chloride	5 kg	360109	587
Potassium chloride	25 kg	360106	587
Potassium chloride	50 kg	360104	587
Potassium citrate tribasic monohydrate	500 g	359956	591
Potassium citrate tribasic monohydrate	1 kg	359957	591
Potassium citrate tribasic monohydrate	2.5 kg	359958	591
Potassium citrate tribasic monohydrate	5 kg	359959	591
Potassium hydroxide, flakes	25 kg	362201	599
Potassium hydroxide, flakes	25 kg	362202	599
Potassium hydroxide, pellets	1 kg	362237	600
Potassium hydroxide, pellets	5 kg	362239	600
Potassium hydroxide, pellets	25 kg	362235	600
Potassium iodide	250 g	362405	607
Potassium iodide	1 kg	362407	607
Potassium iodide	5 kg	362409	607
Potassium iodide	10 kg	362403	607
Potassium iodide	25 kg	362402	607
Potassium metabisulfite	1 kg	362627	609
Potassium metabisulfite	5 kg	362629	609
Potassium metabisulfite	10 kg	362622	609
Potassium metabisulfite	25 kg	362623	609
Potassium nitrate	1 kg	363007	610
Potassium nitrate	5 kg	363009	610
Potassium nitrate	25 kg	363002	610
Potassium permanganate	1 kg	363107	611
Potassium permanganate	5 kg	363109	611
Potassium permanganate	25 kg	363101	611
Potassium phosphate monobasic	1 kg	361507	615
Potassium phosphate monobasic	5 kg	361509	615
Potassium phosphate monobasic	25 kg	361503	615
Potassium sodium tartrate tetrahydrate	1 kg	363457	617
Potassium sodium tartrate tetrahydrate	5 kg	363459	617
Propan-2-ol	1 l	309501	626
Propan-2-ol	2.5 l	309505	626

Product specifications are subject to changes. Please visit our website for updates.

Propan-2-ol	5 l	529165	626
Propan-2-ol	10 l	309506	626
Propan-2-ol	25 l	309504	626
Propan-2-ol	25 l	309507	626
Propan-2-ol	200 l	309500	626
Propan-2-ol	200 l	309509	626
Propan-2-ol 70%	5 l	524195	626
Propyl p-hydroxybenzoate	50 g	363953	629
Propyl p-hydroxybenzoate	500 g	363956	629
Propylene glycol	1 l	346701	630
Propylene glycol	2.5 l	346703	630
Propylene glycol	60 kg	346705	630
Propylene glycol	200 kg	346708	630
Rice starch	1 kg	313107	640
Rice starch	2.5 kg	313108	640
Rice starch	5 kg	313109	640
Rice starch	25 kg	313102	640
Salicylic acid	1 kg	306381	643
Salicylic acid	1 kg	306377	643
Silver nitrate	100 g	320904	655
Silver nitrate	1 kg	320907	655
Sodium acetate anhydrous	1 kg	366377	662
Sodium acetate anhydrous	5 kg	366372	662
Sodium acetate anhydrous	25 kg	366371	662
Sodium acetate trihydrate	1 kg	366207	663
Sodium acetate trihydrate	5 kg	366209	663
Sodium acetate trihydrate	25 kg	366205	663
Sodium alginate	100 g	366551	663
Sodium alginate	1 kg	366552	663
Sodium alginate	5 kg	366553	663
Sodium benzoate	1 kg	366757	665
Sodium benzoate	5 kg	366759	665
Sodium benzoate	25 kg	366754	665
Sodium bicarbonate	1 kg	366908	666
Sodium bicarbonate	5 kg	366909	666
Sodium bicarbonate	25 kg	366902	666
Sodium bicarbonate	50 kg	366904	666
Sodium bromide	1 kg	367357	667
Sodium bromide	5 kg	367359	667
Sodium carbonate anhydrous	1 kg	367707	668
Sodium carbonate anhydrous	5 kg	367703	668
Sodium carbonate anhydrous	25 kg	367705	668
Sodium carbonate anhydrous	50 kg	367704	668
Sodium carbonate decahydrate	1 kg	367608	668
Sodium carbonate decahydrate	5 kg	367609	668
Sodium carbonate decahydrate	25 kg	367601	668
Sodium carbonate monohydrate	1 kg	367691	669
Sodium carbonate monohydrate	5 kg	367692	669
Sodium carbonate monohydrate	25 kg	367693	669
Sodium carbonate monohydrate	50 kg	367694	669

Sodium chloride	10 kg	368281	671
Sodium chloride	1 kg	368257	671
Sodium chloride	5 kg	368259	671
Sodium chloride	25 kg	368253	671
Sodium citrate dibasic sesquihydrate	1 kg	367951	672
Sodium citrate tribasic anhydrous	1 kg	368107	672
Sodium citrate tribasic anhydrous	20 kg	368102	672
Sodium citrate tribasic dihydrate	1 kg	368057	673
Sodium citrate tribasic dihydrate	5 kg	368058	673
Sodium citrate tribasic dihydrate	10 kg	368052	673
Sodium citrate tribasic dihydrate	25 kg	368051	673
Sodium citrate tribasic dihydrate	50 kg	368054	673
Sodium glycerophosphate pentahydrate	1 kg	369447	676
Sodium glycerophosphate pentahydrate	5 kg	369449	676
Sodium hydroxide, pearls	1 kg	369743	678
Sodium hydroxide, pearls	5 kg	369741	678
Sodium hydroxide, pearls	25 kg	369742	678
Sodium hydroxide, pearls	25 kg	369744	678
Sodium hydroxide, pellets	1 kg	369777	679
Sodium hydroxide, pellets	5 kg	369772	679
Sodium hydroxide, pellets	20 kg	369771	679
Sodium hydroxide, pellets	25 kg	369774	679
Sodium iodide	250 g	370305	693
Sodium iodide	1 kg	370307	693
Sodium iodide	5 kg	370309	693
Sodium metabisulfite	1 kg	370751	694
Sodium metabisulfite	2.5 kg	370752	694
Sodium metabisulfite	25 kg	370753	694
Sodium nitrite	1 kg	371901	696
Sodium nitrite	5 kg	371902	696
Sodium nitrite	25 kg	371903	696
Sodium phosphate dibasic anhydrous	1 kg	369212	700
Sodium phosphate dibasic anhydrous	5 kg	369213	700
Sodium phosphate dibasic anhydrous	25 kg	369211	700
Sodium phosphate dibasic anhydrous	25 kg	369275	700
Sodium phosphate dibasic dihydrate	5 kg	369185	701
Sodium phosphate dibasic dodecahydrate	1 kg	369158	702
Sodium phosphate dibasic dodecahydrate	5 kg	369159	702
Sodium phosphate dibasic dodecahydrate	25 kg	369152	702
Sodium phosphate dibasic dodecahydrate	50 Kg	369154	702
Sodium phosphate monobasic dihydrate	1 kg	369138	703
Sodium phosphate monobasic dihydrate	5 kg	369139	703
Sodium phosphate monobasic dihydrate	25 kg	369132	703
Sodium phosphate monobasic monohydrate	1 kg	369143	703
Sodium phosphate monobasic monohydrate	5 kg	369141	703
Sodium phosphate monobasic monohydrate	25 kg	369142	703
Sodium salicylate	1 kg	373607	705
Sodium salicylate	5 kg	373608	705
Sodium salicylate	25 kg	373603	705
Sodium sulfate anhydrous	1 kg	375713	707

Product specifications are subject to changes. Please visit our website for updates.

Sodium sulfate anhydrous	25 kg	375716	707
Sodium sulfite anhydrous	1 kg	376006	708
Sodium sulfite anhydrous	2.5 kg	376008	708
Sodium sulfite anhydrous	5 kg	376009	708
Sodium sulfite anhydrous	10 kg	376002	708
Sodium sulfite anhydrous	25 kg	376003	708
Sodium tetraborate decahydrate	1 kg	367207	710
Sodium tetraborate decahydrate	5 kg	367209	710
Sodium tetraborate decahydrate	25 kg	367201	710
Sodium thiosulfate pentahydrate	1 kg	377907	711
Sodium thiosulfate pentahydrate	5 kg	377909	711
Sodium thiosulfate pentahydrate	25 kg	377901	711
Sorbitol (no crystallizable) solution 70%	1 l	379021	715
Sorbitol (no crystallizable) solution 70%	5 l	379022	715
Stearic acid	2.5 kg	307112	722
Stearic acid	25 kg	307115	722
D(+)-Sucrose	1 kg	365157	726
D(+)-Sucrose	5 kg	365158	726
D(+)-Sucrose	25 kg	365152	726
Sulfuric acid 96%	1 l	306651	733
Sulfuric acid 96%	2.5 l	306657	733
Sulfuric acid 96%	50 kg	306653	733
Talc	1 kg	382107	744
Talc	5 kg	382109	744
Talc	25 kg	382105	744
Tannic acid	1 kg	307157	744
Tannic acid	5 kg	307152	744
Tannic acid	25 kg	307153	744
L(+) Tartaric Acid	1 kg	307357	745
L(+) Tartaric Acid	5 kg	307359	745
L(+) Tartaric Acid	1 kg	307307	745
L(+) Tartaric Acid	5 kg	307309	745
L(+) Tartaric Acid	25 kg	307301	745
Thymol	250 g	384205	757
Thymol	1 kg	384201	757
Thymol	2.5 kg	384202	757
Titanium dioxide	1 kg	385751	764
Titanium dioxide	5 kg	385752	764
Titanium dioxide	25 kg	385753	764
Triethanolamine	1 l	386301	772
Triethanolamine	2.5 l	386303	772
Triethanolamine	30 kg	386304	772
Triethanolamine	220 kg	386305	772
Tris (hydroxymethyl)-aminomethane	25 kg	313441	777
Vanillin	100 g	388104	784
Vanillin	1 kg	388107	784
Vanillin	5 kg	388108	784
Vanillin	10 kg	388102	784
Water purified	1 l	307606	788
Water purified	5 l	307601	788

Water purified	10 l	307602	788
Water purified	25 kg	307603	788
Water purified	200 l	307604	788
Zinc oxide	1 kg	393507	801
Zinc oxide	5 kg	393509	801
Zinc oxide	25 kg	393503	801
Zinc stearate	1 kg	395451	801
Zinc stearate	10 kg	395452	801
Zinc sulfate heptahydrate	1 kg	394007	802
Zinc sulfate heptahydrate	5 kg	394009	802
Zinc sulfate heptahydrate	25 kg	394001	802

ERBApharm®: Titrated and Diluted solutions manufactured from raw material according to Ph.Eur

CARLO ERBA Reagents provides a range of titrated and diluted solutions described (or not) in the monographs of European Pharmacopoeia. If there are not listed in it, CARLO ERBA Reagents is able to propose you the solution you need manufactured with raw materials conform to the required Pharmacopoeia, and with same guarantee of quality in terms of traceability, documentation and analysis. A feasibility study is carried out to meet your requirements.

Description	Notes	Size	Code	Page
Ethanol 50% v/v		5 l	529261	337
Hydrochloric acid 5%		10 l	PS0864/41	400
Hydrochloric acid 6 mol/l (6N)		25 l	528651	401
Hydrochloric acid 4 mol/l (4N)		1 l	528681	402
Hydrochloric acid 2 mol/l (2N)		1 l	528691	403
Hydrochloric acid 1 mol/l (1N)		1 l	528583	404
Hydrochloric acid 1 mol/l (1N)		5 l	528584	404
Hydrochloric acid 0.1 mol/l (0.1N)		1 l	528661	406
Hydrochloric acid 0.1 mol/l (0.1N)		5 l	528662	406
Propan-2-ol 70%		1 l	524182	627
Propan-2-ol 70%	6 units / box	1 l	524183	627
Propan-2-ol 70%		2.5 l	524184	627
Propan-2-ol 70%		5 l	524181	627
Sodium hydroxide solution 30%		1 l	369704	682
Sodium hydroxide solution 30%		20 l	369702	682
Sodium hydroxide solution 30%		10 kg	369701000	682
Sodium hydroxide solution 30%		200 l	369706	682
Sodium hydroxide solution 20% w/w		1 l	480631	683
Sodium hydroxide 6 mol/l (6N)		2 l	524651	684
Sodium hydroxide 3 mol/l (3N)		500 ml	524732	685
Sodium hydroxide 2 mol/l (2N)		1 l	524671	686
Sodium hydroxide 1 mol/l (1N)		1 l	524761	687
Sodium hydroxide 1 mol/l (1N)		1 l	524621	687
Sodium hydroxide 0.25 mol/l (0.25N)		5 l	369812	689
Sodium hydroxide 0.1 mol/l (0.1N)		1 l	524631	690

PHARMACEUTICAL QUALITY CONTROL

The European Pharmacopoeia defines, in chapters "2. Methods of analysis" and "4. Reagents", the procedures and products suitable for the analysis of pharmaceutical raw materials or finished products.
CARLO ERBA Reagents offers a wide range of products, prepared according to the instructions set forth in the European Pharmacopoeia currently in place, which meet the specific needs of analysis and quality control laboratories in the pharmaceutical industry.
Our solutions are ready to use or to be diluted immediately prior to use, allowing significant time savings with guaranteed conformity with the European Pharmacopoeia (Ph.Eur.) requirements. To facilitate their identification, our reagents are named as indicated in the Pharmacopoeia and coded with their specific reference number preceded by the number 61.
The following is guaranteed for all our pharmacopoeia reagents and solutions:
- Shelf life ranging from 2 to 24 months
- No particular storage or transport conditions
- Certificate of analysis with the lot number, expiration date and declared conformity with the European Pharmacopoeia requirements

These reagents and solutions are to be used during the analysis of pharmaceutical raw materials or finished products, and therefore should not be confused with the pharmacopoeia products.

USP Reagents

These reagents are specially intended for use in testing USP and monographs. These solutions are produced according USP specifications using USP reagent quality components.

Description	Notes	Size	Code	Page
Ammonia solution 6N		1 l	617000151	164
Arsenic trioxide solution	Arsenic trioxyde stock solution	100 ml	617000001	185
Barium chloride 30 g/l	Barium chloride TS	100 ml	617000161	194
Bromine solution	Brome TS	100 ml	617000141	210
Cobalt (II) chloride in solution	Cobaltous Chloride CS	100 ml	616001028	272
Cobalt (II) chloride in solution	Cobaltous Chloride CS	500 ml	616001057	272
Copper (II) sulphate solution	Cupric sulfate CS	100 ml	616001038	281
Copper (II) sulphate solution	Cupric sulfate CS	500 ml	616001037	281
Degree of coloration of liquids: primary solutions	Yellow primary solution	100 ml	612002100	291
Degree of coloration of liquids: primary solutions	Red primary solution	100 ml	612002200	291
Degree of coloration of liquids: primary solutions	Blue primary solution	100 ml	612002300	291
Hydrochloric acid 1 mol/l (1N)		1 l	617000191	404
Iron chloride in solution	Ferric Chloride CS	100 ml	616001048	434
Iron chloride in solution	Ferric Chloride CS	500 ml	616001047	434
Lead (II) acetate cotton	Lead acetate cotton	10 g	617000301	457
Lead (II) nitrate	Stock Solution TS	100 ml	617000321	458
Methyl red solution	Methyl red solution TS	100 ml	617000111	504
Mixtures for residual solvents analysis	5 elements (Class 1): Benzene 10mg/ml; Tetrachloromethane (Carbon tetrachloride) 20mg/ml; 1,2-Dichloroethane 25mg/ml; 1,1-Dichloroethene 40mg/ml; 1,1,1-Trichloroethane 50mg/ml; Matrix: Dimethylsulphoxide	1 ml	507692	509
Mixtures for residual solvents analysis	16 elements (Class 2): Acetonitrile 2.05mg/ml; Chlorobenzene 1.8mg/ml; Cumene 0.34mg/ml; Cyclohexane 19.4mg/ml; cis-1,2-Dichloroethene 4.7mg/ml; trans-1,2-Dichloroethene 4.7mg/ml; 1,4-Dioxan 1.9mg/ml; Ethylbenzene 1.84mg/ml; Methanol 15mg/ml; Methylcyclohexane 5.9mg/ml; Dichloromethane 3mg/ml; Tetrahydrofuran 3.6mg/ml; Toluene 4.45mg/ml; m-Xylene 6.51mg/ml; o-Xylene 0.98mg/ml; p-Xylene 1.52mg/ml; Matrix: Dimethylsulphoxide	1 ml	507693	509
Mixtures for residual solvents analysis	8 elements (Class 2): Chloroform 60µg/ml; 1,2-Dimethoxyethane 100µg/ml; n-Hexane 290µg/ml; 2-Hexanone 50µg/ml; Nitromethane 50µg/ml; Pyridine 200µg/ml; 1,2,3,4-Tetrahydronaphthalene (Tetralin) 100µg/ml; Trichloroethene 80µg/ml; Matrix: Dimethylsulphoxide	1 ml	507694	509

Phenolphthalein solution 1% in ethanol	Phenolphtalein TS	100 ml	617000131	571
Silver nitrate solution	Silver nitrate TS	1 l	617000201	659
Sodium hydroxide 1 mol/l (1N)		500 ml	617000121	687
Thioacetamide	Thioacetamide TS	100 ml	617000211	754
Titanium trichloride-sulfuric acid reagent	Titanium Trichloride-Sulfuric Acid TS	100 ml	617000221	765

Elemental Impurities - Multielement Standard for ICP analysis according to USP

Description	Notes	Size	Code	Page
Multielement standard for ICP	4 elements: Cd 5mg/l; Pb 10mg/l; As 15mg/l; Hg 15mg/l. Matrix: 7% HNO3	100 ml	506120	513
Multielement standard for ICP	8 elements: Cd 25mg/l; Pb 5mg/l; As 1.5mg/l; Hg 15mg/l; Mo 100mg/l; Ni 500mg/l; V 100mg/l; Cu 1000mg/l. Matrix: 7% HNO3	100 ml	506110	513
Multielement standard for ICP	Precious metals - 6 elements: Ir 100 mg/l, Pt 100 mg/l; Os 100 mg/l; Rh 100 mg/l; Pd 100 mg/l; Ru 100 mg/l. Matrix: 15% HCl	100 ml	506150	513
Multielement standard for ICP	Precious metals - 6 elements: Ir 10 mg/l, Pt 10 mg/l; Os 10 mg/l; Rh 10 mg/l; Pd 10 mg/l; Ru 10 mg/l. Matrix: 15% HCl	100 ml	506130	513
Multielement standard for ICP	8 elements: Cd 2.5mg/l; Pb 5mg/l; As 1.5mg/l; Hg 1.5mg/l; Mo 10mg/l; Ni 50mg/l; V 10mg/l; Cu 100mg/l. Matrix: 7% HNO3	100 ml	506140	513

Ph. Eur. Reagents, Chapter 2.2.1: Clarity and degree of opalescence of liquids

Description	Notes	Size	Code	Page
Primary opalescent suspension	Formazin suspension	100 ml	612201100	621
Primary opalescent suspension	Formazin suspension	1 l	612201101	621

Ph. Eur. Reagents, Chapter 2.2.2: Degree of coloration of liquids

Description	Notes	Size	Code	Page
Degree of coloration of liquids: primary solutions	Yellow primary solution	100 ml	612202100	290
Degree of coloration of liquids: primary solutions	Red primary solution	100 ml	612202200	290
Degree of coloration of liquids: primary solutions	Blue primary solution	100 ml	612202300	290
Degree of coloration of liquids: standard solutions	Standard solution B (brown)	125 ml	612202510	291
Degree of coloration of liquids: standard solutions	Standard solution BY (brownish-yellow)	125 ml	612202520	291
Degree of coloration of liquids: standard solutions	Standard solution Y (yellow)	125 ml	612202530	291
Degree of coloration of liquids: standard solutions	Standard solution GY (greenish-yellow)	125 ml	612202540	291
Degree of coloration of liquids: standard solutions	Standard solution R (red)	125 ml	612202550	291
Hydrochloric acid, dilute	Dilution matrix HCl 10g/L	1 l	612202400	409

Product specifications are subject to changes. Please visit our website for updates.

Ph. Eur. Reagents, Chapter 2.2.25: Absorption spectrophotometry, ultraviolet and visible

The European Pharmacopoeia defines in its “Analytical Methods”, chapter 2.2.25, the calibration of spectrophotometers absorption in the ultraviolet and visible.
For each parameter, the method of analysis recommends to use diluted solutions.
Among the multiple benefits, here are the most importants:
- Reduces associated risk of handling some hazardous reagents
- Traceability according to NIST for raw material and instrumentation used to prepare and control the finish product
- Available in 100 ml amber bottles or permanently sealed cuvettes (possibility of re-calibration on request)
- Can be used with all UV-VIS Spectrophotometers

Description	Notes	Size	Code	Page
Holmium perchlorate in solution		100 ml	506473	391
Potassium chloride 12g/l		10 ml	506432	589
Potassium chloride 12g/l		100 ml	506433	589
Potassium dichromate - Sulfuric acid solution	conc. 60 mg/l	2 x 10 ml	506442	592
Potassium dichromate - Sulfuric acid solution	conc. 600 mg/l	2 x 10 ml	506452	592
Potassium dichromate - Sulfuric acid solution	conc. 60 mg/l	100 ml	506443	592
Potassium dichromate - Sulfuric acid solution	conc. 600 mg/l	100 ml	506453	592
Toluene in solution in hexane		10 ml	506462	768
Toluene in solution in hexane		100 ml	506463	768
Water	Spectrophotometry Stray Light Blank	100 ml	506411	787

Ph. Eur. Reagents, Chapter 2.2.3: Potentiometric determination of pH

These buffer solutions are intended for:
- Calibration of pH-meters
- Validation of analytical methods
- Preparation of "working reference samples"
- Detection limit and linearity studies

These Certified Reference Materials are produced under a quality management system that is:
- Registered to ISO 9001 Quality Management System
- Accredited according to ISO/IEC 17025
- Accredited according to ISO 17034

Description	Notes	Size	Code	Page
Buffer pH 1.68	Certified Reference Material	500 ml	612203168	216
Buffer pH 4	pH 4.01 at 25°C - Certified Reference Material	500 ml	612203401	218
Buffer pH 6.88	pH 6.87 at 25°C - Certified Reference Material	500 ml	612203687	220
Buffer pH 7.4	pH 7.41 at 25°C - Certified Reference Material	500 ml	612203741	222
Buffer pH 9.22	pH 9.18 at 25°C - Certified Reference Material	500 ml	612203918	223

Ph. Eur. Reagents, Chapter 2.4.24: Identification and control of residual solvents

The International Conference on Harmonization (ICH) adopted the “Impurities: guidelines for residual solvents”, which sets the pharmaceutical limits for residual solvents in the manufacturing of drug substances, excipients, and drug products. The methods used to determine these limits are described in chapter 467 of the USP and 2.4.24 of the European Pharmacopeia.
Three classes of solvents are categorized:
Class 1: Solvents to be avoided (Known human carcinogens, strongly suspected human carcinogens, environmental hazards)
Class 2: Residual solvents (Solvents to be limited, nongenotoxic animal carcinogens or possible causative agents of other irreversible toxicity, such as neurotoxicity or teratogenicity. Solvents suspected of other significant but reversible toxicities.)
Class 3: Residual Solvents (Solvents with Low Toxic Potential. Solvents with low toxic potential to humans; no health-based exposure limit is needed. Class 3 residual solvents may have PDEs of up to 50 mg or more per day.)
CARLO ERBA Reagents developed a product line of mixture of class 1 and 2 solvents, offered at concentrations within the acceptable limit mentioned in chapter “5.4 - Residual solvents”. This allows you to control the amount of residual solvents in your starting material used in the manufacturing of drug products.
Advantages: delivered with a certificate of analysis guaranteeing total traceability:
- Lot number
- Expiration date
- CAS number of each component
- Molecular formula of each component
- Lot number of each starting material
- Concentration of each starting material
- Concentration variation limits

Description	Notes	Size	Code	Page
Mixtures for residual solvents analysis	5 elements (Class 1): Benzene 4mg/ml; 1,2-Dichloroethane 10mg/ml; Tetrachloromethane 8mg/ml; 1,1-Dichloroethene 16mg/ml; 1,1,1-Tric hloroethane 20mg/ml; Matrix: Dimethylsulphoxide	1 ml	507688	509
Mixtures for residual solvents analysis	14 elements (Class 2): Chlorobenzene 360µg/ml; Cyclohexane 3880µg/ml; cis-1,2-Dichloroethene 1870µg/ml; Dichloromethane 600µg/ml; Ethylbenzene 369µg/ml; n-Hexane 290µg/ml; Methylcyclohexane 1180µg/ml; n,n-Dimethylformamide 880µg/ml; Toluene 890µg/ml; 1,1,2-Trichloroethene 80µg/ml; m-Xylene 1302µg/ml; o-Xylene 195µg/ml; p-Xylene 304µg/ml; Tetrahydrofuran 720µg/ml; Matrix: Dimethylsulfoxide	1 ml	507689	509
Mixtures for residual solvents analysis	11 elements (Class 2): Acetonitrile 410mg/l; Chloroform 60mg/l; 1,2-Dimethoxyethane 100mg/l; n,n-Dimethylacetamide 1090mg/l; Dioxan 380mg/l; 2-Hexanone 50mg/l; Methanol 3000mg/l; Nitromethane 50mg/l; Pyridine 200mg/l; 1,2,3,4-Tetrahydronaphthalene 100mg/l; Isop ropylbenzene (Cumene) 70mg/l; Matrix: Dimethylsulphoxide/ Water	1 ml	507690	509
Mixtures for residual solvents analysis	6 elements (Class 2): Ethyleneglycol-monoethyl ether (2-Ethoxyethanol) 160mg/l; Ethyleneglycol 620mg/l; Formamide 220mg/l; Ethylene glycol-monomethyl ether (2-Methoxyethanol) 50mg/l; 1-Methyl-2-pyrrolidon 4840mg/l; Sulfolan 160mg/l; Matrix: Water	1 ml	507691	509

Ph. Eur. Reagents, Chapter 4.1.1: Reagents

Description	Notes	Size	Code	Page
Acetic acid 30%	Ref Ph.Eur 1000401	1 l	611000401	134
Acetic acid 12%	Ref Ph.Eur 1000402	1 l	611000402	135
Acetic anhydride	Acetic anhydride solution R1 Ref Ph.Eur 1000501	1 l	611000501	137
Aminohippuric acid reagent	Ref Ph.Eur 1003701	100 ml	611003701	158
Aminopyrazolone solution		100 ml	611004601	159
Ammonia solution 17%	Ref Ph.Eur 1004701	250 ml	611004701	164
Ammonia solution diluted	Ammonia, dilute R1 Ref Ph.Eur 1004702	1 l	611004702	164

Product specifications are subject to changes. Please visit our website for updates.

Ammonia solution diluted	Ammonia, dilute R2 Ref Ph.Eur 1004703	1 l	611004703	164
Ammoniacal solution of copper tetrammine	Ref Ph.Eur 1022600	100 ml	611022600	165
Ammonium carbonate solution 158 g/l	Ref Ph.Eur 1005201	1 l	611005201	168
Anisaldehyde solution	Ref Ph.Eur 1007301	100 ml	611007301	181
Anisaldehyde solution	Anisaldehyde solution R1 Ref Ph.Eur 1007302	100 ml	611007302	181
Antimony trichloride	Ref Ph.Eur 1007701	100 ml	611007701	183
Barium chloride dihydrate	Ref Ph.Eur 1009300	100 g	611009300	192
Barium chloride solution 61 g/l	Barium chloride solution R1 Ref Ph.Eur 1009301	100 ml	611009303	193
Barium chloride solution 61 g/l	Barium chloride solution R1 Ref Ph.Eur 1009301	250 ml	611009309	193
Barium chloride solution 61 g/l	Barium chloride solution R1 Ref Ph.Eur 1009301	1 l	611009301	193
Barium hydroxide solution 47.3 g/l	Ref Ph.Eur 1009401	250 ml	611009409	195
Barium hydroxide solution 47.3 g/l	Ref Ph.Eur 1009401	1 l	611009401	195
Biuret reagent	Biuret reagent Ref Ph.Eur 1011601	1 l	611011601	205
Bromine solution	Ref Ph.Eur 1012401	100 ml	611012401	210
Bromine water	Ref Ph.Eur 1012402	50 ml	611012409	211
Bromine water	Ref Ph.Eur 1012402	100 ml	611012402	211
Bromocresol green solution	Ref Ph.Eur 1012601 / Colur change: pH 3.6 (yellow) to pH 5.2 (blue)	100 ml	611012601	211
Bromocresol green - Methyl red solution	Ref Ph.Eur 1012602	100 ml	611012602	212
Bromocresol purple solution	Ref Ph.Eur 1012701	100 ml	611012701	213
Bromophenol blue solution	Ref Ph.Eur 1012801/Color change: pH 2.8 (yellow) to pH 4.4 (blue-violet)	100 ml	611012801	214
Bromophenol blue solution	Bromophenol blue solution R1 Ref Ph.Eur 1012802	100 ml	611012802	214
Bromophenol blue solution	Bromophenol blue solution R2 Ref Ph.Eur 1012803	100 ml	611012803	214
Bromothymol blue solution	Bromothymol blue solution R1 Ref Ph.Eur 1012901/ Color change: pH 5.8 (yellow) to pH 7.4 (blue)	100 ml	611012901	215
Bromothymol blue solution	Bromothymol blue solution R3 Ref Ph.Eur 1012903	100 ml	611012903	215
Calcium sulfate hemihydrate solution	Ref Ph.Eur 1015201	100 ml	611015202	245
Calcium sulfate hemihydrate solution	Ref Ph.Eur 1015201	1 l	611015201	245
Chloral hydrate	Ref Ph.Eur 1017901	100 ml	611017901	256
Citric acid monohydrate	Ref Ph.Eur 1021000	100 g	611021000	269
Congo red solution	Ref Ph.Eur 1022001	100 ml	611022001	274
Copper (II) sulfate solution 12.5%	Ref Ph.Eur 1022500	100 ml	611022501	281
Cresol red R solution in Sodium hydroxide 0.1N/ Ethanol 96% / Water	Ref Ph.Eur 1022801	100 ml	611022801	283
Crystal violet solution 0.5% in anhydrous acetic acid	Ref Ph.Eur 1022901	100 ml	611022901	284
Cupri-citric solution	Ref Ph.Eur 1023100	1 l	611023100	284
Cupri-tartaric solution	Ref Ph.Eur 1023300	2 x 500 ml	611023300	285
Dichloromethane acidified with 1% hydrochloric acid	Ref Ph.Eur 1055901	100 ml	611055901	305
1,4-Dioxane	Dioxane stock solution 1.0 mg/ml Ref Ph.Eur 1032001	50 ml	611032001	321
1,4-Dioxane	Dioxane solution R1 0.1 mg/ml Ref Ph.Eur 1032003	50 ml	611032003	321
1,4-Dioxane	Dioxane solution 0.5 mg/ml Ref Ph.Eur 1032002	100 ml	611032002	321
Diphenylamine solution 1% in sulfuric acid	Ref Ph.Eur 1032101	100 ml	611032109	322

Diphenylamine solution 1% in sulfuric acid	Ref Ph.Eur 1032101	1 l	611032101	322
Diphenylamine solution 1% in sulfuric acid	Diphenylamine solution R1 Ref Ph.Eur 1032102	1 l	611032102	322
Eriochrome black T	Ref Ph.Eur 1056801	100 g	611056801	329
Ethanol 96°	Ref Ph.Eur 1002501	1 l	611002501	334
Ethylene oxide solution	Ethylene oxide solution R2 Ref Ph.Eur 1036408	1 ml	611036408	349
Ethylene oxide solution	Ethylene oxide stock solution Ref Ph.Eur 1036401	10 ml	611036401	349
Ferroin 0.025 mol/l solution	Ref Ph.Eur 1038100	100 ml	611038100	354
Formaldehyde 35% w/w	Ref Ph.Eur 1039101	100 ml	611039101	359
Fuchsin solution decolorised	Ref Ph.Eur 1039401	100 ml	611039401	365
Fuchsin solution decolorised	Fuchsin solution, decolorised R1 Ref Ph.Eur 1039402	100 ml	611039402	365
Holmium perchlorate in solution	Ref Ph.Eur 1043101	1 l	611043101	391
Hydrochloric acid 25% w/v	Hydrochloric acid R1 Ref Ph.Eur 1043501	1 l	611043501	398
Hydrochloric acid, dilute	Ref Ph.Eur 1043503	1 l	611043503	409
Hydrochloric acid, dilute	Hydrochloric acid, dilute R1 Ref Ph.Eur 1043504	1 l	611043504	409
Hydrochloric acid, brominated	Ref Ph.Eur 1043507	1 l	611043507	409
Hydroxylamine solution, alcoholic	Ref Ph.Eur 1044301	100 ml	611044301	417
Indigo carmine solution	Ref Ph.Eur 1045601	1 l	611045601	422
Iodine bromide solution	Ref Ph.Eur 1045901	1 l	611045901	426
Iodoplatinate reagent	Ref Ph.Eur 1046300	200 ml	611046309	427
Iodoplatinate reagent	Ref Ph.Eur 1046300	1 l	611046300	427
Iron (III) ammonium sulfate solution 100 g/l	Ferric ammonium sulfate solution R2 Ref Ph.Eur 1037702	100 ml	611037703	433
Iron (III) ammonium sulfate solution 100 g/l	Ferric ammonium sulfate solution R2 Ref Ph.Eur 1037702	1 l	611037702	433
Lanthanum nitrate solution 50 g/l	Ref Ph.Eur 1048001	1 l	611048001	454
Lead (II) acetate basic solution	Ref Ph.Eur 1048400	100 ml	611048400	457
Lead (II) acetate cotton	Ref Ph.Eur 1048101	10 g	611048101	457
Lead (II) acetate paper	Ref Ph.Eur 1048102	50 pc	611048102	457
Lead (II) acetate solution 95 g/l	Ref Ph.Eur 1048103	100 ml	611048103	458
Lead (II) nitrate solution 33 g/l	Ref Ph.Eur 1048301	1 l	611048301	459
Malachite green solution 0.5% in acetic acid anhydrous	Ref Ph.Eur 1050501	1 l	611050501	476
Mercuric bromide paper	Ref Ph.Eur 1052101	50 pc	611052101	484
Mercury chloride solution 54 g/l	Ref Ph.Eur 1052201	100 ml	611052201	485
Mercury (II) sulfate solution	Ref Ph.Eur 1052600	100 ml	611052600	486
Methanol, hydrochloric	Ref Ph.Eur 1053203	100 ml	611053203	493
Methyl Orange solution 0.1%	Solution in ethanol Ref Ph.Eur 1054802	100 ml	611054802	501
Methyl orange mixed solution	Ref Ph.Eur 1054801	100 ml	611054801	502
Methyl red solution	Ref Ph.Eur 1055102	100 ml	611055102	504
Methyl red mixed solution	Ref Ph.Eur 1055101	100 ml	611055101	504
Molybdovanadic reagent	Ref Ph.Eur 1056700	100 ml	611056700	511
alpha-Naphtholbenzein solution 0.2% in acetic acid	Ref Ph.Eur 1057601	100 ml	611057601	518
Ninhydrin solution	Ninhydrin solution R1 Ref Ph.Eur 1058304	100 ml	611058304	526
Ninhydrin solution	Ninhydrin solution R2 Ref Ph.Eur 1058305	100 ml	611058305	526
Ninhydrin and Tin (II) chloride reagent	Ref Ph.Eur 1058301	100 ml	611058301	526
Nitric acid, dilute	Ref Ph.Eur 1058402	100 ml	611058402	534
Nitric acid, dilute	Ref Ph.Eur 1058402	250 ml	611058409	534
Pararosaniline solution, decolorised	Ref Ph.Eur 1062201	100 ml	611062201	556

Product specifications are subject to changes. Please visit our website for updates.

Perchloric acid solution	Ref Ph.Eur 1062901	100 ml	611062901	562
Phenol red solution	Ref Ph.Eur 1063601 / Colour change: pH 6.8 (yellow) to pH 8.4 (reddish-violet)	100 ml	611063601	571
Phenol red solution	Phenol red solution R2 Ref Ph.Eur 1063603	500 ml	611063603	571
Phenolphthalein solution 1% in ethanol	Phenolphtalein solution R1 Ref Ph.Eur 1063703	100 ml	611063703	571
Phenolphthalein solution 0.1%	Ref Ph.Eur 1063702	100 ml	611063702	572
Phenolphtalein solution 0.1% in ethanol	Ref Ph.Eur 1063702	1 l	611063709	572
Phenylhydrazine hydrochloride solution	Ref Ph.Eur 1064501	100 ml	611064501	573
Phosphomolybdotungstic reagent	Ref Ph.Eur 1065000	100 ml	611065000	576
Phosphotungstic acid solution	Ref Ph.Eur 1065200	100 ml	611065200	578
Picric acid solution	Ref Ph.Eur 1065801	100 ml	611065801	579
Picric acid solution	Picric acid solution R1 Ref Ph.Eur 1065802	100 ml	611065802	579
Potassium carbonate	Ref Ph.Eur 1068900	100 g	611068900	586
Potassium chloride 0.1 mol/l (0.1N)	Ref Ph.Eur 1069101	1 l	611069101	589
Potassium chromate 5% solution	Ref Ph.Eur 1069201	1 l	611069201	591
Potassium dichromate solution 0.5%	Potassium dichromate solution R1 Ref Ph.Eur 1069502	100 ml	611069509	592
Potassium dichromate solution 0.5%	Potassium dichromate solution R1 Ref Ph.Eur 1069502	1 l	611069502	592
Potassium dichromate solution 106 g/l	Ref Ph.Eur 1069501	1 l	611069501	594
Potassium ferricyanide	Ref Ph.Eur 1069800	500 g	611069700	594
Potassium ferrocyanide trihydrate	Ref Ph.Eur 1069800	500 g	611069800	595
Potassium ferrocyanide solution 53 g/l	Ref Ph.Eur 1069801	100 ml	611069801	596
Potassium hydrogen phthalate 0.2 mol/l (0.2N)	Ref Ph.Eur 1070001	1 l	611070001	598
Potassium hydroxide solution 3% in ethanol	Ref Ph.Eur 1070303	100 ml	611070303	601
Potassium hydroxide 2 mol/l (2N) in ethanol	Ref Ph.Eur 1070301	100 ml	611070301	602
Potassium hydroxide 0.5 mol/l (0.5N) in ethanol	Ref Ph.Eur 1070302	1 l	611070302	603
Potassium hydroxide in solution	Ref Ph.Eur 3005001	500 ml	613005001	606
Potassium iodide solution	Solution saturated Ref Ph.Eur 1070504	100 ml	611070504	608
Potassium iodide solution	Solution iodinated R1 Ref Ph.Eur 1070505	100 ml	611070505	608
Potassium iodide solution	Potassium iodide solution 166 g/l Ref Ph.Eur 1070502	1 l	611070502	608
Potassium iodobismuthate solution	Ref Ph.Eur 1070600	100 ml	611070600	609
Potassium iodobismuthate solution	Potassium iodobismuthate solution R2 Ref Ph.Eur 1070602	100 ml	611070602	609
Potassium permanganate solution 3%	Ref Ph.Eur 1070902	1 l	611070902	612
Potassium permanganate and phosphoric acid solution	Ref Ph.Eur 1070901	100 ml	611070901	613
Potassium phosphate monobasic 0.2 mol/l (0.2N)	Ref Ph.Eur 1069601	1 l	611069601	615
Potassium pyroantimonate solution	Ref Ph.Eur 1071300	100 ml	611071303	616
Potassium pyroantimonate solution	Ref Ph.Eur 1071300	250 ml	611071309	616
Potassium pyroantimonate solution	Ref Ph.Eur 1071300	1 l	611071302	616
Potassium tetraiodomercurate solution, alkaline	Ref Ph.Eur 1071600	200 ml	611071600	619
Potassium thiocyanate solution	A 97 g/l solution Ref Ph.Eur 1071801	1 l	611071801	621
Salicylaldehyde azine	Ref Ph.Eur 1075500	100 ml	611075500	642
Silver manganese paper	Ref Ph.Eur 1078200	50 stripes	611078200	655
Silver nitrate solution	Silver nitrate R1 Ref Ph.Eur 1078301	100 ml	611078307	658
Silver nitrate solution	Silver nitrate R2 Ref Ph.Eur 1078302	100 ml	611078306	658
Silver nitrate solution	Silver nitrate R1 Ref Ph.Eur 1078301	1 l	611078301	658
Silver nitrate solution	Silver nitrate R2 Ref Ph.Eur 1078302	1 l	611078302	658
Sodium carbonate solution	A 106 g/l solution ref Ph.Eur 1079301	1 l	611079301	670
Sodium hydroxide solution 20% w/v	Ref Ph.Eur 1081401	1 l	611081401	683

Sodium hydroxide solution	Ref Ph.Eur 1081402	1 l	611081402	691
Sodium hydroxide solution	Sodium hydroxide solution, strong Ref Ph.Eur 1081404	1 l	611081404	691
Sodium hydroxide solution, methanolic	Sodium hydroxide solution, methanolic R1 Ref Ph.Eur 1081405	100 ml	611081405	691
Sodium hypochlorite solution in water	Sodium hypochlorite solution, strong Ref Ph.Eur 1081600	250 ml	611081609	692
Sodium hypochlorite solution in water	Sodium hypochlorite solution, strong Ref Ph.Eur 1081600	1 l	611081600	692
Sodium sulfide nonahydrate solution	Ref Ph.Eur 1083901	100 ml	611083901	707
Sodium sulfide nonahydrate solution	Sodium sulfide solution R1 Ref Ph.Eur 1083902	100 ml	611083902	707
Starch soluble solution	Ref Ph.Eur 1085103	100 ml	611085103	721
Starch soluble solution	Ref Ph.Eur 1085103	1 l	611085104	721
Sulfomolybdic reagent	Sulfomolybdic reagent R3 Ref Ph.Eur 1086500	100 ml	611086500	729
Sulfuric acid, dilute	A 98 g/l solution Ref Ph.Eur 1086804	1 l	611086804	743
Thioacetamide solution 40 g/l	Ref Ph.Eur 1089602	100 ml	611089603	755
Thioacetamide solution 40 g/l	Ref Ph.Eur 1089602	1 l	611089602	755
Thioglycolic acid 80%	Ref Ph.Eur 1089700	10 ml	611089700	755
Thymol blue solution	Ref Ph.Eur 1090601	100 ml	611090601	758
Thymolphthalein solution 0.1% in ethanol	Ref Ph.Eur 1090701	100 ml	611090701	759
Tin (II) chloride solution	Ref Ph.Eur 1085001	100 ml	611085001	761
Titanium trichloride-sulfuric acid reagent	Ref Ph.Eur 1091202	100 ml	611091202	765
o-Tolidine solution	Ref Ph.Eur 1123001	500 ml	611123001	765
Tris(hydroxymethyl)aminomethane solution	Ref Ph.Eur 1094201	100 ml	611094201	777
Vanillin solution, phosphoric	Ref Ph.Eur 1095302	100 ml	611095302	784
Water	Water ammonium-free Ref Ph.Eur 1095501	1 l	611095501	787
Water	Water nitrate-free Ref Ph.Eur 1095506	1 l	611095506	787
Zinc, activated	Ref Ph.Eur 1096501	100 g	611096501	797
Zinc chloride solution, iodinated	Ref Ph.Eur 1096602	500 ml	611096603	800
Zinc chloride-formic acid solution	Ref Ph.Eur 1096601	1 l	611096601	800

Ph. Eur. Reagents Chapter 4.1.2: Standard solutions for limit tests

To be diluted immediately prior to use as described in the pharmacopoeia.

Description	Notes	Size	Code	Page
Aluminum standard solution	A 200 ppm solution Ref Ph.Eur 5000200	100 ml	615000200	151
Aluminum standard solution	A 10 ppm solution: to dilute according to Ref Ph.Eur 5000201	100 ml	615000201	151
Aluminum standard solution	A 2 ppm solution: to dilute according to Ref Ph.Eur 5000202	100 ml	615000202	151
Aluminum standard solution	A 100 ppm solution: to dilute according to Ref Ph.Eur 5000203	100 ml	615000203	151
Ammonium standard solution	A 2.5 ppm solution: to dilute according to Ref Ph.Eur 5000301	100 ml	615000301	165
Ammonium standard solution	A 1 ppm solution: to dilute according to Ref Ph.Eur 5000302	100 ml	615000302	165
Ammonium standard solution	A 100 ppm solution: to dilute according to Ref Ph.Eur 5000300	100 ml	615000309	165
Antimony standard solution	A 100 ppm solution Ref Ph.Eur 5000400	100 ml	615000400	182
Arsenic standard solution	A 1 ppm solution: to dilute according to Ref Ph.Eur 5000501	100 ml	615000501	185
Arsenic standard solution	A 10 ppm solution: to dilute according to Ref Ph.Eur 5000500	100 ml	615000509	185
Barium standard solution	A 0.1 % solution Ref Ph.Eur 5000601	100 ml	615000601	191

Product specifications are subject to changes. Please visit our website for updates.

Barium standard solution	A 50 ppm solution: to dilute according to Ref Ph.Eur 5000600	100 ml	615000609	191
Bismuth standard solution	A 100 ppm solution: to dilute according to Ref Ph.Eur 5005300	100 ml	615005300	202
Cadmium standard solution	A 0.1 % solution Ref Ph.Eur 5000700	100 ml	615000700	232
Cadmium standard solution	A 10 ppm solution: to dilute according to Ref Ph.Eur 5000701	100 ml	615000709	232
Calcium standard solution	A 100 ppm solution: to dilute according to Ref Ph.Eur 5000801	100 ml	615000801	234
Calcium standard solution	A 100 ppm alcoholic solution: to dilute according to Ph.Eur 5000802	100 ml	615000802	234
Calcium standard solution	A 10 ppm solution: to dilute according to Ref Ph.Eur 5000803	100 ml	615000803	234
Calcium standard solution	A 100 ppm solution R1: to dilute according to Ref Ph.Eur 5000804	100 ml	615000804	234
Calcium standard solution	A 400 ppm solution: to dilute according to Ref Ph.Eur 5000800	100 ml	615000809	234
Chloride standard solution	A 5 ppm solution: to dilute according to Ref Ph.Eur 5000901	100 ml	615000901	257
Chloride standard solution	A 8 ppm solution: to dilute according to Ref Ph.Eur 5000900	100 ml	615000909	257
Chloride standard solution	A 50 ppm solution: to dilute according to Ref Ph.Eur 5004100	100 ml	615004100	257
Chromium standard solution	A 0.1 % solution Ref Ph.Eur 5001002	100 ml	615001002	264
Chromium standard solution	A 100 ppm solution Ref Ph.Eur 5001000	1 l	615001000	264
Cobalt standard solution	A 100 ppm solution Ref Ph.Eur 5004300	1 l	615004300	270
Copper standard solution	A 0.1 % solution Ref Ph.Eur 5001100	100 ml	615001100	275
Ferricyanide standard solution	A 50 ppm solution: to dilute according to Ph.Eur 5001300	100 ml	615001300	353
Ferrocyanide standard solution	A 100 ppm solution: to dilute according to Ph.Eur 5001200	100 ml	615001209	354
Fluoride standard solution	A 1 ppm solution: to dilute according to Ref Ph.Eur 5001401	100 ml	615001401	358
Fluoride standard solution	A 10 ppm solution: to dilute according to Ref Ph.Eur 5001400	100 ml	615001409	358
Germanium standard solution	A 100 ppm solution Ref Ph.Eur 5004400	1 l	615004400	369
Glyoxal standard solution	A 20 ppm solution: to dilute according to Ph.Eur 5003700	100 ml	615003700	374
Iodine 10 ppm	Concentrated solution: to dilute according to Ref Ph.Eur 5003800	100 ml	615003809	426
Iron standard solution	A 10 ppm solution: to dilute according to Ref Ph.Eur 5001601	100 ml	615001601	428
Iron standard solution	A 8 ppm solution: to dilute according to Ref Ph.Eur 5001602	100 ml	615001602	428
Iron standard solution	A 2 ppm solution: to dilute according to Ref Ph.Eur 5001603	100 ml	615001603	428
Iron standard solution	A 0.1 % solution Ref Ph.Eur 5001605	100 ml	615001605	428
Iron standard solution	A 250 ppm solution: to dilute according to Ref Ph.Eur 5001606	100 ml	615001606	428
Iron standard solution	A 20 ppm solution: to dilute according to Ref Ph.Eur 5001600	100 ml	615001609	428
Lead standard solution	A 0.1 % solution Ref Ph.Eur 5001700	100 ml	615001700	455
Lead standard solution	A 10 ppm solution: to dilute according to Ref Ph.Eur 5001702	100 ml	615001702	455
Lead standard solution	A 2 ppm solution: to dilute according to Ref Ph.Eur 5001703	100 ml	615001703	455
Lead standard solution	A 1 ppm solution: to dilute according to Ref Ph.Eur 5001704	100 ml	615001704	455

Lead standard solution	A 0.1 ppm solution: to dilute according to Ref Ph.Eur 5001705	100 ml	615001705	455
Lead standard solution	A 10 ppm solution R1: to dilute according to Ref Ph.Eur 5001706	100 ml	615001706	455
Lead standard solution	A 100 ppm solution: to dilute according to Ref Ph.Eur 5001701	100 ml	615001709	455
Lead standard solution	A 0.1 % solution Ref Ph.Eur 5001700	500 ml	615001701	455
Magnesium standard solution	A 10 ppm solution: to dilute according to Ref Ph.Eur 5001801	100 ml	615001801	468
Magnesium standard solution	A 10 ppm solution: to dilute according to Ref Ph.Eur 5001802	100 ml	615001802	468
Magnesium standard solution	A 0.1 % solution Ref Ph.Eur 5001803	100 ml	615001803	468
Magnesium standard solution	A 100 ppm solution: to dilute according to Ref Ph.Eur 5001800	100 ml	615001809	468
Manganese standard solution	A 100 ppm solution Ref Ph.Eur 5004500	1 l	615004500	479
Manganese standard solution	A 1.000 ppm solution Ref Ph.Eur 5005800	1 l	615005800	479
Mercury standard solution	A 10 ppm solution: to dilute according to Ref Ph.Eur 5001901	100 ml	615001901	484
Mercury standard solution	A 1000 ppm solution Ref Ph.Eur 5001900	1 l	615001900	484
Nickel standard solution	A 0.1 ppm solution: to dilute according to Ref Ph.Eur 5002001	100 ml	615002001	521
Nickel standard solution	A 0.2 ppm solution: to dilute according to Ref Ph.Eur 5002002	100 ml	615002002	521
Nickel standard solution	A 10 ppm solution: to dilute according to Ref Ph.Eur 5002000	100 ml	615002009	521
Nitrate standard solution	A 10 ppm solution: to dilute according to Ref Ph.Eur 5002101	100 ml	615002101	527
Nitrate standard solution	A 2 ppm solution: to dilute according to Ref Ph.Eur 5002102	100 ml	615002102	527
Nitrate standard solution	A 100 ppm solution: to dilute according to Ref Ph.Eur 5002100	100 ml	615002109	527
Palladium standard solution	A 500 ppm solution Ref Ph.Eur 5003600	100 ml	615003600	553
Phosphate standard solution	A 5 ppm solution: to dilute according to Ref Ph.Eur 5002200	100 ml	615002200	574
Phosphate standard solution	A 200 ppm solution Ref Ph.Eur 5004200	1 l	615004200	574
Potassium standard solution	A 20 ppm solution: to dilute according to Ref Ph.Eur 5002401	100 ml	615002401	582
Potassium standard solution	A 0.1 % solution Ref Ph.Eur 5002402	100 ml	615002402	582
Potassium standard solution	A 100 ppm solution: to dilute according to Ref Ph.Eur 5002400	100 ml	615002409	582
Potassium standard solution	A 600 ppm solution: to dilute according to Ref Ph.Eur 5005100	100 ml	615005100	582
Selenium standard solution	A 1 ppm solution Ref Ph.Eur 5002501	100 ml	615002501	647
Selenium standard solution	A 100 ppm solution Ref Ph.Eur 5002500	1 l	615002500	647
Silver standard solution	A 5 ppm solution: to dilute according to Ref Ph.Eur 5002600	100 ml	615002609	653
Sodium standard solution	A 50 ppm solution: to dilute according to Ref Ph.Eur 5002701	100 ml	615002701	661
Sodium standard solution	A 200 ppm solution: to dilute according to Ref Ph.Eur 5002700	100 ml	615002709	661
Sodium standard solution	A 1000 ppm solution Ref Ph.Eur 5005700	1 l	615005700	661
Sulfate standard solution	A 10 ppm solution R1: to dilute according to Ref Ph.Eur 5002801	100 ml	615002801	728
Sulfate standard solution	A 100 ppm solution: to dilute according to Ref Ph.Eur 5002802	100 ml	615002802	728
Sulfate standard solution	A 10 ppm solution: to dilute according to Ref Ph.Eur 5002800	100 ml	615002809	728
Sulfite standard solution	A 1,5 ppm solution Ref Ph.Eur 5002900	100 ml	615002900	729

Thallium standard solution	A 10 ppm solution Ref Ph.Eur 5003000	100 ml	615003000	754
Tin standard solution	A 0.1 ppm solution: to dilute according to Ref Ph.Eur 5003101	100 ml	615003101	760
Tin standard solution	A 5 ppm solution: to dilute according to Ref Ph.Eur 5003100	100 ml	615003109	760
Titanium standard solution	A 100 ppm solution Ref Ph.Eur 5003200	1 l	615003200	763
Vanadium standard solution	A 1 g/l solution Ref Ph.Eur 5003300	100 ml	615003300	783
Zinc standard solution	A 10 ppm solution: to dilute according to Ref Ph.Eur 5003402	100 ml	615003402	797
Zinc standard solution	A 5 ppm solution: to dilute according to Ref Ph.Eur 5003403	100 ml	615003403	797
Zinc standard solution	A 100 ppm solution: to dilute according to Ref Ph.Eur 5003401	100 ml	615003409	797
Zirconium standard solution	A 1 g/l solution Ref Ph.Eur 5003500	100 ml	615003500	804

Ph. Eur. Reagents Chapter 4.1.3: Buffers Solutions

Description	Notes	Size	Code	Page
Acetate buffer pH 6.0	Ref Ph.Eur 4002200	1 l	614002200	130
Acetate buffer pH 4.6	Ref Ph.Eur 4001400	1 l	614001400	130
Acetone	Buffered acetone solution Ref Ph.Eur 4000100	1 l	614000100	140
Ammonium chloride buffer solution pH 10.7	Ref Ph.Eur 4013400	1 l	614013400	169
Ammonium chloride buffer solution pH 10.0	Ref Ph.Eur 4007300	100 ml	614007301	170
Ammonium chloride buffer solution pH 10.0	Ref Ph.Eur 4007300	1 l	614007300	170
Ammonium chloride buffer solution pH 9.5	Ref Ph.Eur 4007200	1 l	614007200	170
Buffer pH 2	Ref Ph.Eur 4000200	1 l	614000200	217
Buffer pH 3.5	Ref Ph.Eur 4000600	250 ml	614000601	217
Buffer pH 3.5	Ref Ph.Eur 4000600	1 l	614000600	217
Buffer pH 3.7	Ref Ph.Eur 4000900	1 l	614000900	218
Buffer pH 5.2	Ref Ph.Eur 4001700	1 l	614001700	219
Buffer pH 7	Ref Ph.Eur 4003500	1 l	614003500	220
Buffer pH 7.4	Ref Ph.Eur 4004600	1 l	614004600	222
Buffer pH 9	Ref Ph.Eur 4000700	1 l	614007000	222
Phosphate buffer pH 9.0	Ref Ph.Eur 4008300	1 l	614008300	574
Phosphate buffer pH 7.4	Ref Ph.Eur 4004800	1 l	614004800	575
Phosphate buffer pH 6.8	Ref Ph.Eur 4003400	1 l	614003400	575
Phosphate buffer pH 6.0	Ref Ph.Eur 4002400	1 l	614002400	575
Phosphate buffer pH 5.5	Ref Ph.Eur 4002000	1 l	614002000	575
Phosphate buffer pH 3.0	Ref Ph.Eur 4000500	100 ml	614000501	576
Phosphate buffer pH 3.0	Ref Ph.Eur 4000500	1 l	614000500	576
Phosphate buffer pH 2.0	Ref Ph.Eur 4007900	1 l	614007900	576
Total-ionic-strength-adjustment buffer	Ref Ph.Eur 4007700	1 l	614007700	769
Total-ionic-strength-adjustment buffer	Ref Ph.Eur 4008800	1 l	614008800	769
Tris(hydroxymethyl)aminomethane buffer solution pH 8.1	Ref Ph.Eur 4006200	1 l	614006200	777
Tris(hydroxymethyl)aminomethane-EDTA buffer solution pH 8.4	Ref Ph.Eur 4006600	1 l	614006600	777
Tris(hydroxymethyl)aminomethane sodium chloride buffer solution pH 7.4	Ref Ph.Eur 4004900	1 l	614004900	778

Ph. Eur. Reagents Chapter 4.2.1: Primary standards for volumetric solutions

Description	Notes	Size	Code	Page
Benzoic acid	Ref Ph.Eur 2000200	100 g	612000200	199
Potassium bromate	Ref Ph.Eur 2000300	50 g	612000300	584
Potassium hydrogen phthalate	Ref Ph.Eur 2000400	50 g	612000400	597
Sodium carbonate monohydrate	Ref Ph.Eur 2000500	50 g	612000500	668
Sodium chloride	Ref Ph.Eur 2000600	250 g	612000600	670
Sulfanilic acid	Ref Ph.Eur 2000700	50 g	612000700	728
Zinc standard solution	Ref Ph.Eur 2000800	100 g	612000800	797

Ph. Eur. Reagents Chapter 4.2.2: Volumetric solutions

Description	Notes	Size	Code	Page
Ammonium thiocyanate 0.1 mol/l (0.1N)	Ref Ph.Eur 3000500	500 ml	613000501	179
Ammonium thiocyanate 0.1 mol/l (0.1N)	Ref Ph.Eur 3000500	1 l	613000500	179
Barium chloride 0.1 mol/l (0.2N)	Ref Ph.Eur 3000600	1 l	613000600	193
Barium perchlorate 0.05 mol/l	Ref Ph.Eur 3000700	1 l	613000700	196
Barium perchlorate 0.025 mol/l	Ref Ph.Eur 3009600	500 ml	613009601	196
Bromide - bromate 0.0167 mol/l	Ref Ph.Eur 3001000	1 l	613001000	210
Cerium (IV) ammonium nitrate 0.1 mol/l	Ref Ph.Eur 3000100	1 l	613000100	251
Cerium (IV) ammonium nitrate 0.01 mol/l	Ref Ph.Eur 3000200	1 l	613000200	251
Cerium (IV) ammonium sulfate 0.1 mol/l	Ref Ph.Eur 3000300	250 ml	613000301	252
Cerium (IV) ammonium sulfate 0.1 mol/l	Ref Ph.Eur 3000300	1 l	613000300	252
Cerium (IV) ammonium sulfate 0.01 mol/l	Ref Ph.Eur 3000400	1 l	613000400	252
Cerium (IV) sulfate 0.1 mol/l	Ref Ph.Eur 3001100	500 ml	613001101	253
Cerium (IV) sulfate 0.1 mol/l	Ref Ph.Eur 3001100	1 l	613001100	253
Cupriethylenediamine hydroxide 1 mol/l	Ref Ph.Eur 3008700	1 l	613008700	285
Ethylenediaminetetraacetic acid disodium salt 0.1 mol/l (0.2N)	Ref Ph.Eur 3005900	500 ml	613005901	346
Ethylenediaminetetraacetic acid disodium salt 0.1 mol/l (0.2N)	Ref Ph.Eur 3005900	1 l	613005900	346
Hyamine 1622 solution 0.004M	Ref Ph.Eur 3000900	100 ml	613000901	392
Hyamine 1622 solution 0.004M	Ref Ph.Eur 3000900	1 l	613000900	392
Hydrochloric acid 6 mol/l (6N)	Ref Ph.Eur 3001500	1 l	613001500	401
Hydrochloric acid 2 mol/l (2N)	Ref Ph.Eur 3001700	1 l	613001700	403
Hydrochloric acid 1 mol/l (1N)	Ref Ph.Eur 3001800	500 ml	613001801	403
Hydrochloric acid 1 mol/l (1N)	Ref Ph.Eur 3001800	1 l	613001800	403
Hydrochloric acid 0.1 mol/l (0.1N)	Ref Ph.Eur 3002100	500 ml	613002101	406
Hydrochloric acid 0.1 mol/l (0.1N)	Ref Ph.Eur 3002100	1 l	613002100	406
Hydrochloric acid 0.1 mol/l (0.1N) in ethanol	Ref Ph.Eur 3008800	1 l	613008800	407
Iodine 0.5 mol/l (1N)	Ref Ph.Eur 3009400	1 l	613009400	425
Iodine 0.05 mol/l (0.1N)	Ref Ph.Eur 3002700	1 l	613002700	425
Iodine 0.01 mol/l (0.02N)	Ref Ph.Eur 3002900	1 l	613002900	426
Iron (II) sulfate 0.1mol/l	Ref Ph.Eur 3001400	1 l	613001400	431
Iron (III) ammonium sulfate 0.1 mol/l	Ref Ph.Eur 3001300	1 l	613001300	433
Lanthanum nitrate 0.1 mol/l	Ref Ph.Eur 3010100	1 l	613010100	454
Lead (II) nitrate 0.1 mol/l	Ref Ph.Eur 3003100	500 ml	613003101	459
Lead (II) nitrate 0.1 mol/l	Ref Ph.Eur 3003100	1 l	613003100	459
Lead (II) nitrate 0.05 mol/l	Ref Ph.Eur 3009700	100 ml	613009700	459
Perchloric acid 0.1 mol/l (0.1N)	Ref Ph.Eur 3003900	1 l	613003900	561
Perchloric acid 0.05 mol/l (0.05N)	Ref Ph.Eur 3004000	1 l	613004000	562

Product specifications are subject to changes. Please visit our website for updates.

Potassium bromate 0.033 mol/l (0.198N)	Ref Ph.Eur 3004200	1 l	613004200	585
Potassium bromate 0.02 mol/l (0.12N)	Ref Ph.Eur 3004300	1 l	613004300	585
Potassium dichromate 0.0167 mol/l (0.1N)	Ref Ph.Eur 3004600	1 l	613004600	593
Potassium hydroxide 1 mol/l (1N)	Ref Ph.Eur 3009100	1 l	613009100	602
Potassium hydroxide 0.5 mol/l (0.5N) in ethanol	Ref Ph.Eur 3004900	1 l	613004900	603
Potassium hydroxide 0.1 mol/l (0.1N)	Ref Ph.Eur 3004800	1 l	613004800	604
Potassium hydroxide 0.1 mol/l (0.1N) in ethanol	Ref Ph.Eur 3005100	1 l	613005100	605
Potassium iodate 0.05 mol/l (0.3N)	Ref Ph.Eur 3005200	1 l	613005200	606
Potassium permanganate 0.02 mol/l (0.1N)	Ref Ph.Eur 3005300	100 ml	613005301	612
Potassium permanganate 0.02 mol/l (0.1N)	Ref Ph.Eur 3005300	250 ml	613005309	612
Potassium permanganate 0.02 mol/l (0.1N)	Ref Ph.Eur 3005300	1 l	613005300	612
Silver nitrate 0.1 mol/l (0.1N)	Ref Ph.Eur 3005600	1 l	613005600	657
Sodium arsenite 0.1 mol/l (0.2N)	Ref Ph.Eur 3005800	1 l	613005800	664
Sodium hydroxide 1 mol/l (1N)	Ref Ph.Eur 3006300	500 ml	613006301	686
Sodium hydroxide 1 mol/l (1N)	Ref Ph.Eur 3006300	1 l	613006300	686
Sodium hydroxide 0.1 mol/l (0.1N)	Ref Ph.Eur 3006600	500 ml	613006601	689
Sodium hydroxide 0.1 mol/l (0.1N)	Ref Ph.Eur 3006600	1 l	613006600	689
Sodium hydroxide 0.1 mol/l (0.1N) in ethanol	Ref Ph.Eur 3007000	100 ml	613007001	690
Sodium hydroxide 0.1 mol/l (0.1N) in ethanol	Ref Ph.Eur 3007000	1 l	613007000	690
Sodium methoxide 0.1 mol/l	Ref Ph.Eur 3007100	100 ml	613007101	695
Sodium methoxide 0.1 mol/l	Ref Ph.Eur 3007100	1 l	613007100	695
Sodium nitrite 0.1 mol/l (0.1N)	Ref Ph.Eur 3007200	1 l	613007200	697
Sodium thiosulfate 0.1 mol/l (0.1N)	Ref Ph.Eur 3007300	500 ml	613007301	712
Sodium thiosulfate 0.1 mol/l (0.1N)	Ref Ph.Eur 3007300	1 l	613007300	712
Sulfuric acid 0.5 mol/l (1N)	Ref Ph.Eur 3007800	1 l	613007800	738
Sulfuric acid 0.05 mol/l (0.1N)	Ref Ph.Eur 3008000	500 ml	613008001	740
Sulfuric acid 0.05 mol/l (0.1N)	Ref Ph.Eur 3008000	1 l	613008000	740
Tetrabutylammonium hydroxide solution 0.1 mol/l (0.1N)	Ref Ph.Eur 3008300	1 l	613008300	749
Tetrabutylammonium hydroxide 0.1 mol/l (0.1N) in 2-propanol	Ref Ph.Eur 3008400	1 l	613008400	749
Zinc sulfate 0.1 mol/l (0.2N)	Ref Ph.Eur 3008600	500 ml	613008601	802
Zinc sulfate 0.1 mol/l (0.2N)	Ref Ph.Eur 3008600	1 l	613008600	802

PHYSICALS AND CHEMICAL CHARACTERISTICS

In order to help you save time, CARLO ERBA Reagents enlarges its range of standards for quality control laboratories: melting point standards, colour standard, osmolality standards, density standards, brix standards, refractive index standards. Traceability, reliability are the key word of these products. You will find the most commonly used items. Others are available, don't hesitate to contact us.

Melting point, Colour, Osmolality, Density, Brix & Refractive Index Standards

Description	Notes	Size	Code	Page
ASTM colour standards	Sample A1	100 ml	540601	187
ASTM colour standards	Sample A3	100 ml	540602	187
ASTM colour standards	Sample A5	100 ml	540603	187
ASTM colour standards	Sample A7	100 ml	540604	187
Brix Standards	Sucrose Stabilised 0%, Nominal Refractive Index @ 20°C: 1.332986	15 ml	540201	209
Brix Standards	Sucrose Stabilised 5%, Nominal Refractive Index @ 20°C: 1.340264	15 ml	540202	209
Brix Standards	Sucrose Stabilised 7%, Nominal Refractive Index @ 20°C: 1.343253	15 ml	540203	209
Brix Standards	Sucrose Stabilised 10%, Nominal Refractive Index @ 20°C: 1.347824	15 ml	540204	209
Brix Standards	Sucrose Stabilised 11.5%, Nominal Refractive Index @ 20°C: 1.350149	15 ml	540205	209
Brix Standards	Sucrose Stabilised 12%, Nominal Refractive Index @ 20°C: 1.350930	15 ml	540206	209
Brix Standards	Sucrose Stabilised 15%, Nominal Refractive Index @ 20°C: 1.355679	15 ml	540207	209
Brix Standards	Sucrose Stabilised 20%, Nominal Refractive Index @ 20°C: 1.363842	15 ml	540208	209
Brix Standards	Sucrose Stabilised 25%, Nominal Refractive Index @ 20°C: 1.372328	15 ml	540209	209
Brix Standards	Sucrose Stabilised 30%, Nominal Refractive Index @ 20°C: 1.381149	15 ml	540210	209
Brix Standards	Sucrose Stabilised 35%, Nominal Refractive Index @ 20°C: 1.390322	15 ml	540220	209
Brix Standards	Sucrose Stabilised 40%, Nominal Refractive Index @ 20°C: 1.39986	15 ml	540221	209
Brix Standards	Sucrose Stabilised 45%, Nominal Refractive Index @ 20°C: 1.409777	15 ml	540222	209
Brix Standards	Sucrose Stabilised 50%, Nominal Refractive Index @ 20°C: 1.420087	15 ml	540223	209
Brix Standards	Sucrose Stabilised 60%, Nominal Refractive Index @ 20°C: 1.441928	15 ml	540224	209
Brix Standards	Sucrose Stabilised 67.5%, Nominal Refractive Index @ 20°C: 1.459290	15 ml	540225	209
Density Standard	0.6960g/ml at 15°C	100 ml	540401	292
Density Standard	0.8715g/ml at 15°C	100 ml	540402	292
Density Standard	1.0040g/ml at 15°C	100 ml	540403	292
Density Standard	1.2498g/ml at 15°C	100 ml	540404	292
Density Standard	0.6919g/ml at 20°C	100 ml	540405	292
Density Standard	0.7033g/ml at 20°C	100 ml	540406	292
Density Standard	0.7488g/ml at 20°C	100 ml	540407	292
Density Standard	0.8668g/ml at 20°C	100 ml	540408	292
Density Standard	1.0005g/ml at 20°C	100 ml	540409	292
Density Standard	1.0301g/ml at 20°C	100 ml	540410	292
Density Standard	1.0792g/ml at 20°C	100 ml	540411	292
Density Standard	1.1651g/ml at 20°C	100 ml	540412	292

Product specifications are subject to changes. Please visit our website for updates.

Density Standard	1.2486g/ml at 20°C	100 ml	540413	292
Density Standard	1.3304g/ml at 20°C	100 ml	540414	292
Density Standard	1.5799g/ml at 20°C	100 ml	540415	292
Density Standard	1.7470g/ml at 20°C	100 ml	540416	292
Density Standard	1.9141g/ml at 20°C	100 ml	540417	292
Density Standard	0.8207g/ml at 40°C	100 ml	540418	292
Density Standard	0.9323g/ml at 40°C	100 ml	540420	292
Density Standard	1.2408g/ml at 40°C	100 ml	540421	292
Density Standard	0.9990g/ml at 60°C	100 ml	540422	292
Density Standard	0.7524g/ml at 15°C	100 ml	540451	293
Density Standard	0.7721g/ml at 15°C	100 ml	540452	293
Density Standard	0.7933g/ml at 15°C	100 ml	540453	293
Density Standard	0.8168g/ml at 15°C	100 ml	540454	293
Density Standard	0.8428g/ml at 15°C	100 ml	540455	293
Density Standard	0.8715g/ml at 15°C	100 ml	540456	293
Density Standard	0.6919g/ml at 20°C	100 ml	540457	293
Density Standard	0.7033g/ml at 20°C	100 ml	540458	293
Density Standard	0.7488g/ml at 20°C	100 ml	540459	293
Density Standard	0.7893g/ml at 20°C	100 ml	540460	293
Density Standard	0.8126g/ml at 20°C	100 ml	540461	293
Density Standard	0.8384g/ml at 20°C	100 ml	540462	293
Density Standard	0.8668g/ml at 20°C	100 ml	540463	293
Density Standard	0.9098g/ml at 20°C	100 ml	540464	293
Density Standard	0.9476g/ml at 20°C	100 ml	540465	293
Density Standard	1.0005g/ml at 20°C	100 ml	540566	293
Density Standard	1.0301g/ml at 20°C	100 ml	540567	293
Density Standard	0.8622g/ml at 25°C	100 ml	540568	293
Density Standard	0.9438g/ml at 25°C	100 ml	540569	293
Density Standard	0.9969g/ml at 25°C	100 ml	540470	293
Density Standard	0.9245g/ml at 50°C	100 ml	540471	293
Density Standard	0.9695g/ml at 60°C	100 ml	540472	293
Density Standard	0.9815g/ml at 80°C	100 ml	540473	293
Gardner Colour Standards	Colour 2	100 ml	540701	368
Gardner Colour Standards	Colour 4	100 ml	540702	368
Gardner Colour Standards	Colour 6	100 ml	540703	368
Gardner Colour Standards	Colour 8	100 ml	540704	368
Gardner Colour Standards	Colour 10	100 ml	540705	368
Gardner Colour Standards	Colour 12	100 ml	540706	368
Gardner Colour Standards	Colour 14	100 ml	540707	368
Gardner Colour Standards	Colour 16	100 ml	540708	368
Melting point standards	Benzophenone 47 to 49°C	1 g	540001	483
Melting point standards	p-Nitrotoluene 52 to 54°C	1 g	540002	483
Melting point standards	Vanillin 81 to 83°C	1 g	540003	483
Melting point standards	Acetanilide 113 to 116°C	1 g	540014	483
Melting point standards	Benzoic Acid 121 to 123°C	1 g	540004	483
Melting point standards	Phenacetin 133 to 135°C	1 g	540005	483
Melting point standards	Salicylic Acid 158 to 160°C	1 g	540006	483
Melting point standards	Sulfanilamide 164 to 166°C	1 g	540007	483
Melting point standards	Caffeine 235 to 238°C	1 g	540008	483
Melting point standards	Carbazole 243 to 247°C	1 g	540009	483

Melting point standards	Anthraquinone 283 to 286°C	1 g	540010	483
Melting point standards	Set Sulphanilamide Caffeine Vanillin	3 x 1 g	540011	483
Melting point standards	Set Benzophenone (4749°C) Benzoic Acid (121-123°C) Anthraquinone (283-286°C)	3 x 1 g	540012	483
Melting point standards	Set Vanillin (81-83°C) Phenacetin (134-136°C) Caffeine (235-237°C)	3 x 1 g	540013	483
Osmolality Standards	100mOsm/Kg H2O	12 x 5 ml	540301	550
Osmolality Standards	1500mOsm/Kg H2O	12 x 5 ml	540302	550
Osmolality Standards	200mOsm/Kg H2O	12 x 5 ml	540303	550
Osmolality Standards	2000mOsm/Kg H2O	12 x 5 ml	540304	550
Osmolality Standards	290mOsm/Kg H2O	12 x 5 ml	540305	550
Osmolality Standards	300mOsm/Kg H2O	12 x 5 ml	540306	550
Osmolality Standards	400mOsm/Kg H2O	12 x 5 ml	540307	550
Osmolality Standards	500mOsm/Kg H2O	12 x 5 ml	540308	550
Osmolality Standards	850mOsm/Kg H2O	12 x 5 ml	540309	550
Osmolality Standards	900mOsm/Kg H2O	12 x 5 ml	540310	550
Osmolality Standards Protein Based	240mOsm/Kg H2O	12 x 5 ml	540351	550
Osmolality Standards Protein Based	280mOsm/Kg H2O	12 x 5 ml	540352	550
Osmolality Standards Protein Based	320mOsm/Kg H2O	12 x 5 ml	540353	550
Osmolality Standards Urine Based	300mOsm/Kg H2O	12 x 5 ml	540354	550
Osmolality Standards Urine Based	800mOsm/Kg H2O	12 x 5 ml	540355	550
Refractive Index standards	1.34325 at 20°C	15 ml	540101	637
Refractive Index standards	1.34782 at 20°C	15 ml	540102	637
Refractive Index standards	1.35171 at 20°C	15 ml	540103	637
Refractive Index standards	1.37233 at 20°C	15 ml	540104	637
Refractive Index standards	1.38115 at 20°C	15 ml	540105	637
Refractive Index standards	1.40978 at 20°C	15 ml	540106	637
Refractive Index standards	1.42009 at 20°C	15 ml	540107	637
Refractive Index standards	1.44193 at 20°C	15 ml	540108	637
Saybolt Colour Standards	-15	100 ml	540709	644
Saybolt Colour Standards	+0	100 ml	540710	644
Saybolt Colour Standards	+12	100 ml	540711	644
Saybolt Colour Standards	+15	100 ml	540712	644
Saybolt Colour Standards	+19	100 ml	540713	644
Saybolt Colour Standards	+25	100 ml	540714	644
Saybolt Colour Standards	+30	100 ml	540715	644
Viscosity standards	6.7cSt@20°C, 5.8cSt@25°C, 4.2cSt@37.78°C, 4cSt@40°C, 3.2cSt@50°C	500 ml	540801	785
Viscosity standards	14cSt@20°C, 12cSt@25°C, 8cSt@37.78°C, 7.5cSt@40°C, 5.8cSt@50°C	500 ml	540802	785
Viscosity standards	20cSt@20°C, 16cSt@25°C, 11cSt@37.78°C, 10cSt@40°C, 7.5cSt@50°C	500 ml	540803	785
Viscosity standards	30cSt@20°C, 24cSt@25°C, 15cSt@37.78°C, 14cSt@40°C, 10cSt@50°C	500 ml	540804	785
Viscosity standards	88cSt@20°C, 66cSt@25°C, 35cSt@37.78°C, 32cSt@40°C, 21cSt@50°C	500 ml	540805	785
Viscosity standards	160cSt@20°C, 120cSt@25°C, 60cSt@37.78°C, 54cSt@40°C, 35cSt@50°C	500 ml	540806	785
Viscosity standards	790cSt@20°C, 580cSt@25°C, 280cSt@37.78°C, 250cSt@40°C, 160cSt@50°C	500 ml	540807	785
Viscosity standards	3300cSt@20°C, 2300cSt@25°C, 1100cSt@37.78°C, 940cSt@40°C, 560cSt@50°C	500 ml	540808	785

Product specifications are subject to changes. Please visit our website for updates.

Viscosity standards	19000cSt@20°C, 12000cSt@25°C, 4000cSt@37.78°C, 3400cSt@40°C, 1700cSt@50°C	500 ml	540809	785
Viscosity standards	28000cSt@20°C, 17000cSt@25°C, 6000cSt@37.78°C, 5100cSt@40°C, 2500cSt@50°C	500 ml	540810	785
Viscosity standards	41000cSt@20°C, 25000cSt@25°C, 8000cSt@37.78°C, 6700cSt@40°C, 3200cSt@50°C	500 ml	540811	785

pHMETRY

In industrial sectors as well as scientific analysis and research, the growing need for specific reagents for calibration of pH-meters, pH determination and buffering needs has led CARLO ERBA Reagents to offer the following product lines:

- Buffers in ready-to-use solution: colorless & colored
- Buffers in concentrated solution - Normex vial: colorless & colored
- pH indicators in solution
- pH indicator papers

Buffers in colorless solution, ready-to-use

These solutions are traceable to N.I.S.T and precisely standardized at 20°C, which makes them ideal for solving calibration problems with other solutions and preventing errors due to various factors such as the kind of salt or water used.

Description	Notes	Size	Code	Page
Ammonia buffer solution pH 10		5 l	PS0194/22	165
Ammonia buffer solution pH 10		5 l	PS0194/95	165
Buffer acetate pH 4.5		5 l	PS0784/95	216
Buffer pH 1		500 ml	486211	216
Buffer pH 1.68		500 ml	486751	216
Buffer pH 2		500 ml	486231	217
Buffer pH 3		500 ml	486251	217
Buffer pH 3		1 l	486252	217
Buffer pH 3.56		500 ml	486741	218
Buffer pH 4		500 ml	486271	218
Buffer pH 4		1 l	486273	218
Buffer pH 4		5 l	486274	218
Buffer pH 4		10 l	486276	218
Buffer pH 4.62		500 ml	486841	219
Buffer pH 5		500 ml	486311	219
Buffer pH 6		500 ml	486331	220
Buffer pH 6.88		500 ml	486871	220
Buffer pH 7		500 ml	486451	221
Buffer pH 7		1 l	486453	221
Buffer pH 7		5 l	486454	221
Buffer pH 7		10 l	486456	221
Buffer pH 7		25 l	486455	221
Buffer pH 7.20 Weise		500 ml	486411	221
Buffer pH 8		500 ml	486541	222
Buffer pH 8		1 l	486542	222
Buffer pH 9		500 ml	486591	223
Buffer pH 9		1 l	486593	223
Buffer pH 9		5 l	486594	223
Buffer pH 9.22		500 ml	486881	223
Buffer pH 10		500 ml	486611	224
Buffer pH 10		1 l	486613	224
Buffer pH 10		5 l	486614	224
Buffer pH 10		10 l	486615	224
Buffer pH 11		500 ml	486771	224
Buffer pH 11		1 l	486772	224
Buffer pH 12		500 ml	486691	225
Buffer pH 13		500 ml	486701	225
Buffer pH 13		1 l	486702	225
Phosphate buffer pH 7.4		5 l	524965	575

Product specifications are subject to changes. Please visit our website for updates.

Phosphate buffer pH 7.4	5 l	PS0740/95	575
Phosphate buffer pH 7.2	2.5 l	525925	575
Phosphate buffer pH 7.2	25 l	525921	575
Phosphate buffer pH 6.8	10 l	524952	575

Buffers in colored solution, ready-to-use

To enable immediate identification of the buffer, the following colored solutions, traceable to N.I.S.T, are available.

Description	Notes	Size	Code	Page
Buffer pH 4	Color: Red	500 ml	486761	218
Buffer pH 4	Color: Red	1 l	486762	218
Buffer pH 7	Color: Green	500 ml	486791	221
Buffer pH 7	Color: Green	1 l	486792	221
Buffer pH 9	Color: Blue	500 ml	PS0427/19	223
Buffer pH 10.06	Color: Blue	500 ml	486811	224
Buffer pH 10.06	Color: Blue	1 l	486812	224

Buffers in colorless solution, NORMEX

Buffers usually consist of a diluted solution of the buffering system.
Given their usually low concentration, long-term storage of these solutions may lead to the development of mold or other inconveniences.
To avoid these problems, concentrated buffer solutions packaged in polythene NORMEX vials are available.
With the aid of a 500ml volumetric flask, a funnel and a glass stirrer, each Normex buffer package can be used to prepare 500ml of buffer solution at the desired pH with a maximum error of ±0.05 pH units. Detailed instructions for use are printed on the package of each individual Normex buffer. The correct preparation procedure provides for the use of boiled distilled water with a recommended preparation temperature of 20°C.
The characteristics of this line are:
• vial specifically designed by CARLO ERBA Reagents to meet the user's needs
• quick and simple preparation of 500ml of solution at a known pH

Description	Notes	Size	Code	Page
Buffer pH 1	To dilute to 500 ml		486221	216
Buffer pH 2	To dilute to 500 ml		486241	217
Buffer pH 3	To dilute to 500 ml		486261	217
Buffer pH 4	To dilute to 500 ml		486281	218
Buffer pH 5	To dilute to 500 ml		486301	219
Buffer pH 6	To dilute to 500 ml		486321	220
Buffer pH 6.8	To dilute to 500 ml		486401	220
Buffer pH 7	To dilute to 500 ml		486421	221
Buffer pH 7.2	To dilute to 500 ml		486441	221
Buffer pH 7.4	To dilute to 500 ml		486461	222
Buffer pH 8	To dilute to 500 ml		486531	222
Buffer pH 9	To dilute to 500 ml		486571	223
Buffer pH 10	To dilute to 500 ml		486601	224
Buffer pH 11	To dilute to 500 ml		486631	224
Buffer pH 12	To dilute to 500 ml		486621	225
Buffer pH 13	To dilute to 500 ml		486641	225

Buffers in colored solution, NORMEX

To enable immediate identification of the buffer, the following colored concentrated buffer solutions packaged in polythene NORMEX vials are available.
The characteristics of this line are:
• vial specifically designed by CARLO ERBA Reagents to meet the user's needs
• quick and simple preparation of 500ml of solution at a known pH

Description	Notes	Size	Code	Page
Buffer pH 4	Color: red - To dilute to 500 mL		486291	219
Buffer pH 7	Color: yellow - To dilute to 500 mL		486431	221
Buffer pH 10.06	Color: blue - To dilute to 500 mL		486581	224

pH indicators in solution

The direct method, i.e., placing the indicator in the solution being analyzed, is the most efficient and thus the most widely adopted method for acid-base titration.
Since the indicator competes with the species being titrated, its use in significant quantities may alter the result of the titration; therefore it is critically important to choose the most appropriate indicator for the type of analysis being performed.
For the determination of pH values in aqueous solutions, universal indicators in solution are available for measurements in a variety of pH ranges. These products are supplied complete with a color scale and detailed instructions for use. The sensitivity of these solutions is higher than that of pH indicator papers.

Description	Notes	Size	Code	Page
Alizarin saturated solution in ethanol		250 ml	E415932	149
Alkali Blue 6B solution 2% in ethanol		250 ml	E428541	150
Bromocresol green 0.04% hydroalcoholic solution		250 ml	E491255	211
Bromocresol purple solution 0.4% in ethanol		250 ml	E470045	212
Bromophenol blue solution 0.4% in ethanol		250 ml	E428665	213
Bromophenol blue solution 0.02%		100 ml	428691	214
Bromothymol blue 0.4% in ethanol		250 ml	E428715	215
Bromothymol blue 0.02%		100 ml	428731	215
o-Cresol Red solution 0.2% in ethanol		250 ml	E476805	283
Crystal violet solution 0.5% in anhydrous acetic acid		500 ml	E491551	284
Indicator universal pH 0-5 hydroalcoholic solution		25 ml	E455661	421
Indicator universal pH 0-5 hydroalcoholic solution		500 ml	E455662	421
Indicator universal pH 1-11 hydroalcoholic solution		25 ml	E455702	421
Indicator universal pH 1-11 hydroalcoholic solution		500 ml	E455706	421
Indicator universal pH 1-11 water solution		25 ml	E455711	422
Indicator universal pH 1-11 water solution		500 ml	E455712	422
Methylene blue solution 1%		500 ml	E429011	499
Methyl Orange solution 0.1%		500 ml	E423562	501
Methyl red solution water/ethanol 0.2%		250 ml	E476915	503
Methyl red solution 0.1% in ethanol		250 ml	E476921	504
Phenol Red solution 0.2% in ethanol		250 ml	E476845	571
Phenolphthalein solution 1% in ethanol		250 ml	451191	572
Phenolphthalein solution 1% in ethanol	Only for italian market	250 ml	E451191	572
Phenolphthalein solution 1% in ethanol		1 l	451192	572
Phenolphthalein solution 1% in ethanol	Only for italian market	1 l	E451192	572
Thymol blue 0.4% in ethanol		250 ml	E429235	757
Thymolphthalein 0.1% hydroalcoholic solution		250 ml	E487755	759
o-Tolidine solution 0.1%		1 l	488461	765
Tropaeolin O solution 0.1%		500 ml	E490056	779

Product specifications are subject to changes. Please visit our website for updates.

pH paper, litmus and indicator paper in rolls for special uses

Indicator papers represent a convenient and particularly simple instrument for measuring pH.
Indicator papers are actually filter paper impregnated with indicator solutions.
A wide variety of indicator papers are available, allowing the user to select the product most suited to the type of sample being analyzed and the degree of certainty required on the final result.
This is the most basic type of paper and does not require a color scale for pH determination.
This type of test paper provides a relative indication and is useful in determining whether a solution is acid, neutral or basic.

Description	Notes	Size	Code	Page
Congo red paper	Congo red paper, Color change: red --> Blue, Change pH 5.0-->3.0	1 roll	435220000	273
Litmus paper	Blue litmus paper, Color change: blue --> red, change pH 8.0 - 5.0	1 roll	435260000	466
Litmus paper	Neutral litmus paper, Color change: red <-- purple --> Blue, Change pH 5.0 - 8.0	1 roll	435300000	466
Litmus paper	Red litmus paper, Color change: red --> blue, Change pH 5.0 - 8.0	1 roll	435340000	466
Potassium iodide starch paper	Paper starch iodide, Color change: White --> Blue-purple	1 roll	434980000	608

pH paper, high-sensitivity indicator paper in rolls

High-sensitivity test papers cover small pH ranges with a detection range of between 0.2 - 0.5 units, while the attached color scale provides the reference colors.
These test papers are suitable for measuring the pH of unbuffered or weakly buffered solutions.

Description	Notes	Size	Code	Page
Indicator papers	Monochromatic scale, pH range 1.0 - 11.0, Sensitivity 1,0 pH	1 roll	435140000	420
Indicator papers	Monochromatic scale, pH range 1.0 - 14.0, Sensitivity 1,0/2,0 pH	1 roll	435150000	420
Indicator papers	Monochromatic scale, pH range 3.8 - 5.8, Sensitivity 0.2/0.3 pH	1 roll	435161000	420
Indicator papers	Monochromatic scale. pH range 0.5 - 5.5. Sensitivity 0.5 pH	1 roll	435421000	420
Indicator papers	Monochromatic scale. pH range 4.0 - 7.0. Sensitivity 0.3 pH	1 roll	435431000	420
Indicator papers	Monochromatic scale. pH range 6.4 - 8.0. Sensitivity 0.2 pH	1 roll	435441000	420
Indicator papers	Monochromatic scale. pH range 7.2 - 9.7. Sensitivity 0.3 pH	1 roll	435451000	420
Indicator papers	Monochromatic scale. pH range 5.5 - 9.0. Sensitivity 0.5 pH	1 roll	435511000	420

pH paper, in roll three-colored

Three-color test papers have three parallel sections with different colors, formed by three different mixtures of indicators separated by a fine hydrophobic strip to prevent any color migration. The three indicator sections ensure high-precision results that are remarkably easy to read. The reference color chart provides three different color indications for each pH unit.

Description	Notes	Size	Code	Page
Indicator papers	Trichromic scale, pH range 1.0 - 11.0, Sensitivity 1,0 pH	1 roll	435131000	420

pH paper strips, with built-in color scale

These indicator papers are unique in that they have the color scale provided directly on each strip, usually divided into 8 color notches with the corresponding pH indicated for each.
This allows the user to perform quick and precise pH measurements without the need for an external color guide. These test papers are also suitable for colored solutions or suspensions, since the color of the solution to be analyzed acts in the same manner on the indicator scale as it does on the result of the solution being tested.
An invisible hydrophobic strip placed just above the last color notch on the color scale forms a barrier that prevents the test liquid, which may be corrosive, toxic or otherwise dangerous, from coming into contact with the user's fingers through capillary action.

Description	Notes	Size	Code	Page
Indicator papers	pH range 1.0 - 2.8. Sensitivity 0.2/0.3	100 stripes	435493000	421
Indicator papers	pH range 1.8 - 3.8. Sensitivity 0.2/0.3	100 stripes	435494000	421
Indicator papers	pH range 3.8 - 5.5. Sensitivity 0.2/0.3	100 stripes	435496000	421
Indicator papers	pH range 6.0 - 8.1. Sensitivity 0.2/0.3	100 stripes	435498000	421
Indicator papers	pH range 8.0 - 9.7. Sensitivity 0.2/0.3	100 stripes	435502000	421

pH paper strips, color-fixed

These indicator papers are indelible. This property depends on the chemical bond between the indicator and the cellulose fiber, and it allows the test papers to be used even in highly alkaline solutions.
The strip length allows a wide operating margin, thus preventing contact between the operator's fingers and the solution being analyzed.

Description	Notes	Size	Code	Page
Indicator papers	pH range 0.0 - 14.0, Sensitivity 1.0	100 stripes	435121000	421
Indicator papers	pH range 0.0 - 6.0. Sensitivity 0.5	100 stripes	435642000	421
Indicator papers	pH range 2.0 - 9.0. Sensitivity 0.5	100 stripes	435643000	421
Indicator papers	pH range 4.5 - 10.0. Sensitivity 0.5	100 stripes	435644000	421
Indicator papers	pH range 7.0 - 14.0. Sensitivity 0.3/0.4	100 stripes	435645000	421

Solutions for the electrode conservation

The lifetime of a pH electrode depends on the way it is being maintained and kept. A set of solutions are used to keep electrodes accurate and durable.

Description	Notes	Size	Code	Page
Hydrofluoric acid diluted		250 ml	405775	412
Potassium chloride 3.5 mol/l (3.5N)		250 ml	471225	587
Potassium chloride 3.5 mol/l (3.5N) + silver chloride		250 ml	471245	587
Potassium chloride 3 mol/l (3N)		250 ml	471215	588
Potassium chloride 3 mol/l (3N) water-glycerol solution		250 ml	471275	588
Potassium chloride 3 mol/l (3N) + silver chloride		250 ml	471235	588
Potassium chloride saturated solution		250 ml	471265	590
Potassium chloride solution		250 ml	471285	590
Potassium nitrate 1 mol/l (1N)		250 ml	473045	610

Product specifications are subject to changes. Please visit our website for updates.

VOLUMETRY

Volumetric titration is a standard analytical technique used specifically for quantitative determination of analytes in solution. Despite the numerous sophisticated instrumental innovations resulting from recent technological developments, the classical volumetric analysis technique has maintained its importance in the area of analytical chemistry.
CARLO ERBA Reagents offers different lines of volumetric solutions to meet the needs of a wide variety of users:
- Ready-to-use solutions
- Concentrated solutions in Normex vial

Volumetric solutions, ready-to-use

For routine analyses, our line of ready-to-use volumetric solutions is ideal for quick and convenient use with guaranteed quality. CARLO ERBA Reagents offers a wide range with exact title of the batch written on the label and certificate of analysis, provided with traceability to S.R.M. from N.I.S.T., a +/- 0.1% precision, where noted.
These solutions are available in both polythene or glass bottles, with ISO 45 bottle mouths which are adaptable to automatic titrators, as well as 10 and 5 liters Kubidos® packages supplied with tap and cap, ideal for excellent storage and dosing of the product.
Kubidos®. consists of a cubic box containing a HDPE, a tap and cap ensuring easy flow of the required volume.
With the low volume of product in contact with air, the risks of carbonation of alkaline solutions and microbial contamination are limited.

Description	Notes	Size	Code	Page
Acetic acid 1 mol/l (1N)		5 l	524605	136
Acetic acid 1 mol/l (1N)		1 l	502000	136
Acetic acid 0.1 mol/l (0.1N)		1 l	P3100015	136
Ammonium thiocyanate 1 mol/l (1N)		500 ml	420946	179
Ammonium thiocyanate 0.1 mol/l (0.1N)		1 l	420977	179
Ethylenediaminetetraacetic acid disodium salt 0.1 mol/l (0.2N)	Certified with NIST traceability	1 l	405511000	346
Ethylenediaminetetraacetic acid disodium salt 0.1 mol/l (0.2N)	Certified with NIST traceability	5 l	405513000	346
Ethylenediaminetetraacetic acid disodium salt 0.1 mol/l (0.2N)	Certified with NIST traceability	5 l	405514000	346
Ethylenediaminetetraacetic acid disodium salt 0.1 mol/l (0.2N)	Certified with NIST traceability	10 l	405512000	346
Ethylenediaminetetraacetic acid disodium salt 0.05 mol/l (0.1N)	Certified with NIST traceability	1 l	405501000	347
Ethylenediaminetetraacetic acid disodium salt 0.05 mol/l (0.1N)	Certified with NIST traceability	5 l	405502000	347
Ethylenediaminetetraacetic acid disodium salt 0.01 mol/l (0.02N)	Certified with NIST traceability	1 l	405442000	347
Ethylenediaminetetraacetic acid disodium salt 0.01 mol/l (0.02N)	Certified with NIST traceability	5 l	405443000	347
Hydrochloric acid 2 mol/l (2N)	Certified with NIST traceability	1 l	404067000	403
Hydrochloric acid 2 mol/l (2N)	Certified with NIST traceability	5 l	404062000	403
Hydrochloric acid 2 mol/l (2N)	Certified with NIST traceability	10 l	404061000	403
Hydrochloric acid 1 mol/l (1N)		5 l	528673	404
Hydrochloric acid 1 mol/l (1N)	Certified with NIST traceability	1 l	404097000	404
Hydrochloric acid 1 mol/l (1N)	Certified with NIST traceability	5 l	404092000	404
Hydrochloric acid 1 mol/l (1N)	Certified with NIST traceability	10 l	404091000	404
Hydrochloric acid 0.5 mol/l (0.5N)	Certified with NIST traceability	1 l	404147000	405
Hydrochloric acid 0.5 mol/l (0.5N)	Certified with NIST traceability	5 l	404142000	405
Hydrochloric acid 0.5 mol/l (0.5N)	Certified with NIST traceability	10 l	404141000	405
Hydrochloric acid 0.2 mol/l (0.2N) in propanol-2		1 l	526535	405
Hydrochloric acid 0.1 mol/l (0.1N)	Certified with NIST traceability	1 l	404197000	406
Hydrochloric acid 0.1 mol/l (0.1N)	Certified with NIST traceability	5 l	404192000	406
Hydrochloric acid 0.1 mol/l (0.1N)	Certified with NIST traceability	10 l	404191000	406
Hydrochloric acid 0.1 mol/l (0.1N)	Certified with NIST traceability	10 l	404195000	406
Hydrochloric acid 0.1 mol/l (0.1N) in propanol-2		1 l	526536	407
Hydrochloric acid 0.05 mol/l (0.05N)		1 l	PS0587/15	407

Hydrochloric acid 0.04 mol/l (0.04N)		10 l	PS0206/41	408
Hydrochloric acid 0.02 mol/l (0.02N)		1 l	PS0342/15	408
Hydrochloric acid 0.02 mol/l (0.02N)		5 l	526537	408
Hydrochloric acid 0.01 mol/l (0.01N)	Certified with NIST traceability	1 l	404267	408
Hydrofluoric acid 0.1 mol/l		1 l	507410	412
Iodine 0.5 mol/l (1N)	Certified with NIST traceability	500 ml	456135000	425
Iodine 0.5 mol/l (1N)	Certified with NIST traceability	1 l	456137000	425
Iodine 0.05 mol/l (0.1N)	Certified with NIST traceability	500 ml	456036000	425
Iodine 0.05 mol/l (0.1N)	Certified with NIST traceability	1 l	456037000	425
Lithium methoxide 0.1 mol/l (0.1N)		500 ml	E458321	464
Nitric acid 2 mol/l (2N)	Certified with NIST traceability	5 l	408185000	533
Nitric acid 1 mol/l (1N)	Certified with NIST traceability	500 ml	408176000	533
Nitric acid 1 mol/l (1N)	Certified with NIST traceability	1 l	408171000	533
Nitric acid 0.1 mol/l (0.1N)	Certified with NIST traceability	500 ml	408206000	533
Oxalic acid 0.5 mol/l (1N)		500 ml	408826	551
Oxalic acid 0.05 mol/l (0.1N)		500 ml	408856	551
Perchloric acid 0.1 mol/l (0.1N) in acetic acid	Certified with NIST traceability	500 ml	409136	562
Perchloric acid 0.1 mol/l (0.1N) in acetic acid	Certified with NIST traceability	1 l	409131	562
Perchloric acid 0.01 mol/l (0.01N)		500 ml	E409141	562
Potassium hydroxide 1 mol/l (1N)	Certified with NIST traceability	1 l	472287000	602
Potassium hydroxide 1 mol/l (1N)	Certified with NIST traceability	5 l	472282000	602
Potassium hydroxide 1 mol/l (1N)	Certified with NIST traceability	10 l	472281000	602
Potassium hydroxide 0.5 mol/l (0.5N)	Certified with NIST traceability	1 l	472337000	602
Potassium hydroxide 0.5 mol/l (0.5N)	Certified with NIST traceability	5 l	472332000	602
Potassium hydroxide 0.5 mol/l (0.5N)	Certified with NIST traceability	10 l	472331000	602
Potassium hydroxide 0.5 mol/l (0.5N) in ethanol	Certified with NIST traceability	1 l	472021000	603
Potassium hydroxide 0.5 mol/l (0.5N) in ethanol	Certified with NIST traceability	1 l	472022000	603
Potassium hydroxide 0.5 mol/l (0.5N) in methanol	Certified with NIST traceability	500 ml	472366000	603
Potassium hydroxide 0.5 mol/l (0.5N) in methanol	Certified with NIST traceability	1 l	472364000	603
Potassium hydroxide 0.5 mol/l (0.5N) in methanol	Certified with NIST traceability	5 l	472367000	603
Potassium hydroxide 0.25 mol/l (0.25N)	Certified with NIST traceability	1 l	472427000	604
Potassium hydroxide 0.25 mol/l (0.25N)	Certified with NIST traceability	5 l	472422000	604
Potassium hydroxide 0.25 mol/l (0.25N)	Certified with NIST traceability	10 l	472421000	604
Potassium hydroxide 0.1 mol/l (0.1N)	Certified with NIST traceability	1 l	472457000	605
Potassium hydroxide 0.1 mol/l (0.1N)	Certified with NIST traceability	5 l	472452000	605
Potassium hydroxide 0.1 mol/l (0.1N)	Certified with NIST traceability	10 l	472451000	605
Potassium hydroxide 0.1 mol/l (0.1N) in ethanol	Certified with NIST traceability	1 l	472041000	605
Potassium hydroxide 0.1 mol/l (0.1N) in ethanol	Certified with NIST traceability	1 l	472042000	605
Potassium hydroxide 0.1 mol/l (0.1N) in methanol	Certified with NIST traceability	500 ml	472486000	605
Potassium hydroxide 0.1 mol/l (0.1N) in methanol	Certified with NIST traceability	1 l	472484000	605
Potassium permanganate 0.2 mol/l (1N)	Certified with NIST traceability	1 l	473514000	612
Potassium permanganate 0.02 mol/l (0.1N)	Certified with NIST traceability	1 l	473567000	612
Potassium permanganate 0.02 mol/l (0.1N)	Certified with NIST traceability	5 l	473565000	612
Potassium thiocyanate 0.1 mol/l (0.1N)		1 l	E474417	620
Silver nitrate 1 mol/l (1N)	Certified with NIST traceability	500 ml	424036000	656
Silver nitrate 1 mol/l (1N)	Certified with NIST traceability	1 l	424035000	656
Silver nitrate 0.5 mol/l (0.5N)	Certified with NIST traceability	1 l	424051000	656
Silver nitrate 0.1 mol/l (0.1N)	Certified with NIST traceability	1 l	424067000	657
Silver nitrate 0.1 mol/l (0.1N)	Certified with NIST traceability	5 l	424062000	657
Silver nitrate 0.1 mol/l (0.1N)	Certified with NIST traceability	5 l	424063000	657

Silver nitrate 0.1 mol/l (0.1N)	Certified with NIST traceability	10 l	424061000	657
Silver nitrate 0.1 mol/l (0.1N) in 2-propanol		1 l	PS0250/16	657
Silver nitrate 0.05 mol/l (0.05N)	Certified with NIST traceability	1 l	424101000	657
Silver nitrate 0.01 mol/l (0.01N)		1 l	PS0030/15	658
Silver nitrate 0.01 mol/l (0.01N) in propanol-2		1 l	PS0252/16	658
Sodium carbonate 0.5 mol/l (1N)		500 ml	479186	669
Sodium hydroxide 2 mol/l (2N)	Certified with NIST traceability	500 ml	480686000	686
Sodium hydroxide 2 mol/l (2N)	Certified with NIST traceability	1 l	480687000	686
Sodium hydroxide 2 mol/l (2N)	Certified with NIST traceability	5 l	480682000	686
Sodium hydroxide 2 mol/l (2N)	Certified with NIST traceability	10 l	480681000	686
Sodium hydroxide 2 mol/l (2N)	Certified with NIST traceability	20 l	480684000	686
Sodium hydroxide 1 mol/l (1N)	Certified with NIST traceability	1 l	480717000	687
Sodium hydroxide 1 mol/l (1N)	Certified with NIST traceability	5 l	480711000	687
Sodium hydroxide 1 mol/l (1N)	Certified with NIST traceability	5 l	480714000	687
Sodium hydroxide 1 mol/l (1N)	Certified with NIST traceability	10 l	480713000	687
Sodium hydroxide 0.5 mol/l (0.5N)	Certified with NIST traceability	1 l	480777000	688
Sodium hydroxide 0.5 mol/l (0.5N)	Certified with NIST traceability	5 l	480771000	688
Sodium hydroxide 0.5 mol/l (0.5N)	Certified with NIST traceability	10 l	480772000	688
Sodium hydroxide 0.5 mol/l (0.5N)	Certified with NIST traceability	10 l	480773000	688
Sodium hydroxide 0.25 mol/l (0.25N)	Certified with NIST traceability	1 l	480867000	688
Sodium hydroxide 0.25 mol/l (0.25N)	Certified with NIST traceability	5 l	480861000	688
Sodium hydroxide 0.25 mol/l (0.25N)	Certified with NIST traceability	10 l	480862000	688
Sodium hydroxide 0.2 mol/l (0.2N)		1 l	P3440015	689
Sodium hydroxide 0.1 mol/l (0.1N)	Certified with NIST traceability	1 l	480897000	690
Sodium hydroxide 0.1 mol/l (0.1N)	Certified with NIST traceability	5 l	480891000	690
Sodium hydroxide 0.1 mol/l (0.1N)	Certified with NIST traceability	10 l	480892000	690
Sodium hydroxide 0.1 mol/l (0.1N)	Certified with NIST traceability	10 l	480893000	690
Sodium hydroxide 0.01 mol/l (0.01N)		1 l	PS0215/15	690
Sodium thiosulfate 0.5 mol/l (0.5N)		1 l	P3530015	712
Sodium thiosulfate 0.2 mol/l (0.2N)		5 l	P3520022	712
Sodium thiosulfate 0.1 mol/l (0.1N)	Certified with NIST traceability	1 l	484077000	712
Sodium thiosulfate 0.1 mol/l (0.1N)	Certified with NIST traceability	5 l	484072000	712
Sodium thiosulfate 0.1 mol/l (0.1N)	Certified with NIST traceability	10 l	484071000	712
Sulfuric acid 1 mol/l (2N)	Certified with NIST traceability	1 l	410547000	738
Sulfuric acid 1 mol/l (2N)	Certified with NIST traceability	10 l	410548000	738
Sulfuric acid 0.5 mol/l (1N)	Certified with NIST traceability	1 l	410577000	738
Sulfuric acid 0.5 mol/l (1N)	Certified with NIST traceability	5 l	410572000	738
Sulfuric acid 0.5 mol/l (1N)	Certified with NIST traceability	5 l	410575000	738
Sulfuric acid 0.5 mol/l (1N)	Certified with NIST traceability	10 l	410571000	738
Sulfuric acid 0.25 mol/l (0.5N)	Certified with NIST traceability	1 l	410667000	739
Sulfuric acid 0.25 mol/l (0.5N)	Certified with NIST traceability	5 l	410663000	739
Sulfuric acid 0.25 mol/l (0.5N)	Certified with NIST traceability	10 l	410662000	739
Sulfuric acid 0.125 mol/l (0.25N)		5 l	PS0445/22	740
Sulfuric acid 0.05 mol/l (0.1N)	Certified with NIST traceability	1 l	410717000	741
Sulfuric acid 0.05 mol/l (0.1N)	Certified with NIST traceability	5 l	410712000	741
Sulfuric acid 0.05 mol/l (0.1N)	Certified with NIST traceability	10 l	410711000	741
Sulfuric acid 0.05 mol/l (0.1N)	Certified with NIST traceability	10 l	410715000	741
Sulfuric acid 0.025 mol/l (0.05N)		10 l	PS0016/96	741
Sulfuric acid 0.02 mol/l (0.04N)		1 l	PS0219/15	741
Sulfuric acid 0.02 mol/l (0.04N)		5 l	PS0219/95	741

Sulfuric acid 0.02 mol/l (0.04N)	10 l	PS0219/96	741
Sulfuric acid 0.01 mol/l (0.02N)	1 l	PS0047/15	741
Sulfuric acid 0.005 mol/l (0.01N)	5 l	PS0026/95	742
Tetrabutylammonium hydroxide 0.1 mol/l (0.1N) in 2-propanol	500 ml	E487031	749
Tetrabutylammonium hydroxide 0.1 mol/l (0.1N) in methanol / propanol-2 (50/50)	1 l	P3840016	750
Zinc sulfate 0.1 mol/l (0.2N)	1 l	494921	803
Zinc sulfate 0.05 mol/l (0.05N)	1 l	494931	803

NORMEX, concentrated volumetric solutions

CARLO ERBA Reagents also offers a series of concentrated volumetric solutions packaged in convenient NORMEX vials.
These are ideal for users who would like to prepare solutions immediately before the analysis in a rapid and precise manner. The contents of each vial, brought to a volume of 1000ml with distilled water, allows the user to prepare volumetric solutions at a known concentration with a guaranteed titration factor equal to ± 0.005.

Description	Notes	Size	Code	Page
Acetic acid 0.1 mol/l (0.1N)	Volume: 55 ml		401561	136
Ammonium thiocyanate 0.1 mol/l (0.1N)	Volume: 55 ml		421001	179
Ammonium thiocyanate 0.01 mol/l (0.01N)	Volume: 55 ml		421061	179
Ethylenediaminetetraacetic acid disodium salt 0.1 mol/l (0.2N)	Volume: 165 ml		405421	346
Ethylenediaminetetraacetic acid disodium salt 0.01 mol/l (0.02N)	Volume: 55 ml		405431	347
Hydrochloric acid 1 mol/l (1N)	Volume: 165 ml		404111	404
Hydrochloric acid 0.5 mol/l (0.5N)	Volume: 165 ml		404161	405
Hydrochloric acid 0.1 mol/l (0.1N)	Volume: 55 ml		404211	406
Hydrochloric acid 0.01 mol/l (0.01N)	Volume: 55 ml		404251	408
Iodine 0.05 mol/l (0.1N)	Volume: 60 ml		456051	426
Iodine 0.005 mol/l (0.01N)	Volume: 60 ml		456121	426
Nitric acid 0.1 mol/l (0.1N)	Volume: 55 ml		408231	533
Oxalic acid 0.05 mol/l (0.1N)	Volume: 165 ml		408871	551
Oxalic acid 0.005 mol/l (0.01N)	Volume: 55 ml		408901	551
Potassium bromate 0.0167 mol/l (0.1N)	Volume: 60 ml		470681	585
Potassium dichromate 0.0167 mol/l (0.1N)	Volume: 60 ml		470501	593
Potassium hydroxide 1 mol/l (1N)	Volume: 165 ml		472311	602
Potassium hydroxide 0.5 mol/l (0.5N)	Volume: 55 ml		472391	603
Potassium hydroxide 0.1 mol/l (0.1N)	Volume: 55 ml		472511	605
Potassium iodate 0.0167 mol/l (0.1N)	Volume: 60 ml		472601	606
Potassium iodate 0.00167 mol/l (0.01N)	Volume: 60 ml		472631	607
Potassium permanganate 0.02 mol/l (0.1N)	Volume: 65 ml		473591	612
Potassium permanganate 0.002 mol/l (0.01N)	Volume: 60 ml		473661	613
Silver nitrate 0.1 mol/l (0.1N)	Volume: 60 ml		424081	657
Silver nitrate 0.01 mol/l (0.01N)	Volume: 60 ml		424161	658
Sodium arsenite 0.05 mol/l (0.1N)	Volume: 60 ml		402381	664
Sodium carbonate 0.05 mol/l (0.1N)	Volume: 55 ml		479211	669
Sodium chloride 0.1 mol/l (0,1N)	Volume: 55 ml		479781	672
Sodium hydroxide 1 mol/l (1N)	Volume: 165 ml		480741	687
Sodium hydroxide 0.5 mol/l (0.5N)	Volume: 55 ml		480801	688
Sodium hydroxide 0.1 mol/l (0.1N)	Volume: 55 ml		480921	690
Sodium hydroxide 0.01 mol/l (0.01N)	Volume: 55 ml		481001	690
Sodium thiosulfate 0.1 mol/l (0.1N)	Volume: 55 ml		484121	713
Sodium thiosulfate 0.01 mol/l (0.01N)	Volume: 55 ml		484161	713

Product specifications are subject to changes. Please visit our website for updates.

Sulfuric acid 0.5 mol/l (1N)	Volume: 165 ml	410591	738
Sulfuric acid 0.25 mol/l (0.5N)	Volume: 55 ml	410681	739
Sulfuric acid 0.05 mol/l (0.1N)	Volume: 55 ml	410731	741
Sulfuric acid 0.005 mol/l (0.01N)	Volume: 55 ml	410791	742

FOOD ANALYSIS

Rigorous control of all the products in the food chain (milk, cereals, meats, fruits and vegetables) has recently become even more important to prevent problems associated with the adulteration of foodstuffs as well as to guarantee that all the necessary elements for a healthy diet are present in the food we eat every day.

Reagents for food analysis: Milk, Olive Oil, Wine & more

In order to offer its customers ready-to-use reagents prepared in compliance with the specifications set forth by international standards and regulations for food analysis, CARLO ERBA Reagents has developed a wide range of products dedicated to the agriculture and food sector.

Description	Notes	Size	Code	Page
Acetic acid 27%		5 l	508645	135
Acetic acid 20%		10 l	PS0237/41	135
Acetone / water 98/2 (v/v) with bromophenol blue 0.02 g/l		5 l	PS0852/29	141
ADF Solution		2.5 l	526625	147
ADF Solution		10 l	526623	147
Amidoschwarz 10B solution		5 l	502050	157
Amidoschwarz 10B solution		10 l	502051	157
n-Amyl alcohol		1 l	413783	180
tert-Amyl alcohol		250 ml	413941	180
tert-Amyl alcohol		1 l	413944	180
Boric acid 4%		5 l	502002	206
Boric acid 4% with indicator		5 l	502601	206
Boric acid 3%		2.5 l	PS0563/21	206
Boric acid 1% with indicator		5 l	502611	207
Boric acid 1% with indicator		10 l	502612	207
Boric acid 20g/l		5 l	PS0703/22	207
Boric acid 20 g/l with indicator		5 l	PS0562/22	207
Boric buffer solution		10 l	PS0226/41	207
Carrez reagent potassium salt		1 l	502711	249
Carrez reagent zinc salt		1 l	502701	249
Fehling's A reagent		500 ml	449926	353
Fehling's A reagent		1 l	449927	353
Fehling's B reagent		500 ml	E449936	353
Fehling's B reagent		1 l	E449937	353
Folin-Ciocalteu's reagent		500 ml	E463562	358
Griess' reagent		1 l	454481	376
Griess' reagent A		500 ml	454452	376
Griess' reagent B		500 ml	454462	376
Hanus's reagent		1 l	E454872	379
Hydrochloric acid 50% v/v		1 l	504571	393
Hydrochloric acid 26%		2.5 l	PS0769/20	398
Hydrochloric acid 20%		5 l	PS0751/29	399
Hydrochloric acid 12%		5 l	PS0347/22	399
Hydrochloric acid 12%		25 l	PS0347/49	399
Hydrochloric acid 10%		10 l	PS0768/41	400
Hydrochloric acid 1.128% m/v		1 l	502761	401
Hydrochloric Acid 9 mol/l (9N)		2.5 l	PS0313/20	401
Hydrochloric acid 6 mol/l (6N)		1 l	502831	401
Hydrochloric acid 6 mol/l (6N)		18 l	502832	401
Hydrochloric acid 6 mol/l (6N)		5 l	528550000	401

Product specifications are subject to changes. Please visit our website for updates.

Hydrochloric acid 5 mol/l (5N)		1 l	P3160015	402
Hydrochloric acid 5 mol/l (5N)		5 l	P3160095	402
Hydrochloric acid 4 mol/l (4N)		1 l	502010	402
Hydrochloric acid 4 mol/l (4N)		1 l	PS0589/15	402
Hydrochloric acid 4 mol/l (4N)		5 l	PS0589/22	402
Hydrochloric acid 3 mol/l (3N)		1 l	502621	403
Hydrochloric acid 3 mol/l (3N)		2.5 l	502622	403
Hydrochloric acid 3 mol/l (3N)		25 l	502011	403
Hydrochloric acid 0.714 mol/l (N/1.4)		10 l	526531	404
Hydrochloric acid 0.2 mol/l (0.2N)		1 l	502631	405
Hydrochloric acid 0.0714 mol/l (N/14)		10 l	526533	407
Hydrogen peroxide solution 30%		5 l	502044	414
Isoamyl alcohol	With indicator	500 ml	E413903	436
Isoamyl alcohol	Without indicator	1 l	413892	436
Isoamyl alcohol		500 ml	413801	436
Mixture C.H.M.		2.5 l	524411	507
Mixture C.H.M.		5 l	524412	507
Mix Diethyl ether/Ethanol 70/30 w/w		1 l	463251	507
Mix Diethyl ether/Ethanol 70/30 w/w		2.5 l	463255	507
Mix Diethyl ether / Ethanol 95% 50/50(w/w) - 30 mg/l phenolphtaleine		5 l	529371	507
Mix Ethanol absolute/Diethyl ether 50/50 (w/w)		5 l	529311	507
Mix Ethanol absolute/Diethyl ether 50/50 (v/v)		5 l	529381	508
Mixture for checking solderings		5 l	502671	509
NDF Plus solution		25 l	526941	519
NDF Solution		2.5 l	526920	519
NDF Solution		25 l	526921	519
Nessler's reagent single solution		500 ml	464231	520
Nessler's reagent single solution		1 l	464232	520
Nessler's reagent solution A		500 ml	464422	520
Nessler's reagent solution B		500 ml	464432	521
Nitric acid 8 mol/l (8N)		2.5 l	PS0311/20	532
Potassium chromate 5% solution		1 l	502681	591
Potassium hydroxide solution 33%		5 l	PS0766/22	601
Potassium hydroxide 0.46 mol/l (0.46N)		5 l	502212	604
Potassium hydroxide 0.23 mol/l (0.23N)		5 l	502092	604
Reagent for lipolysis		2.5 l	524910	635
Sand of Fontainebleau		1 kg	502064	644
Sand of Fontainebleau		5 kg	502063	644
Sand of Fontainebleau		25 kg	502062	644
Selenic mixture		250 g	463421	646
Selenic mixture		1 kg	463422	646
Sodium chloride 5 mol/l (5N)		1 l	502131	672
Sodium hydroxide solution 50%		10 l	P4540041	680
Sodium hydroxide solution 50%		25 l	P4540049	680
Sodium hydroxide solution 40%		5 l	502721	680
Sodium hydroxide solution 40%		10 l	502722	680
Sodium hydroxide solution 35-37%		5 l	502112	681
Sodium hydroxide solution 35%		1 l	480591	681
Sodium hydroxide solution 35%		25 kg	480593	681

Sodium hydroxide solution 32%		25 kg	524510	681
Sodium hydroxide solution 32%		1 l	480561	681
Sodium hydroxide solution 32%		2.5 l	480566	681
Sodium hydroxide solution 32%		5 l	526521	681
Sodium hydroxide solution 32%		10 l	480564	681
Sodium hydroxide solution 32%		25 kg	480562	681
Sodium hydroxide solution 32%		30 kg	480563	681
Sodium hydroxide solution 30%		5 l	502741	682
Sodium hydroxide solution 30%		1 l	502731	682
Sodium hydroxide solution 20% w/v		10 l	524505	683
Sodium hydroxide solution 20% w/w		1 l	480621	683
Sodium hydroxide solution 20% w/w		30 kg	480622	683
Sodium hydroxide solution 10% w/v		5 l	508615	683
Sodium hydroxide solution 10% w/v		5 l	524506	684
Sodium hydroxide solution 10% w/v		10 l	524507	684
Sodium hydroxide solution 10% w/v		5 l	526642	684
Sodium hydroxide solution 10% w/v		10 l	526641	684
Sodium hydroxide solution 10% w/v		50 l	526644	684
Sodium hydroxide solution 5% w/v		5 l	524502	684
Sodium hydroxide solution 5% w/v		10 l	524501	684
Sodium hydroxide solution 5% w/v		5 l	526632	684
Sodium hydroxide solution 5% w/v		10 l	526634	684
Sodium hydroxide 5 mol/l (5N)		1 l	526513	685
Sodium hydroxide 5 mol/l (5N)		5 l	526512	685
Sodium hydroxide 4 mol/l (4N)		2.5 l	502662	685
Sodium hydroxide 4 mol/l (4N)		10 l	502664	685
Sodium hydroxide 0.7 mol/l (N/1.4)		10 l	526511	687
Sodium hydroxide 0.357 mol/l (0.357N)	Certified with NIST traceability	1 l	480837000	688
Sodium hydroxide 0.2 mol/l (0.2N)		500 ml	502782	689
Sodium hydroxide 0.2 mol/l (0.2N)		10 l	502781000	689
Sodium thiosulfate 0.0394 mol/l (0.0394N)		2.5 l	484141	713
Sodium thiosulfate 0.0197 mol/l (0.0197N)		2.5 l	484155	713
Sulfuric acid 98%		1 l	410421	730
Sulfuric acid 98%		2.5 l	502641	730
Sulfuric acid 90%		1 l	410391	734
Sulfuric acid 90%		2.5 l	410394	734
Sulfuric acid 85%		1 l	PS0433/15	734
Sulfuric acid 72%		2.5 l	502771	735
Sulfuric acid 69%		2.5 l	PS0893/21	735
Sulfuric acid 62%		2.5 l	PS0894/21	735
Sulfuric acid 50%		1 l	E306702	735
Sulfuric acid 50%		5 l	528541	735
Sulfuric acid 50%		35 kg	E306704	735
Sulfuric acid 30%		1 l	PS0009/15	736
Sulfuric acid 25%		1 l	504562	736
Sulfuric acid 25%		2.5 l	PS0212/21	736
Sulfuric acid 20%		1 l	410511000	736
Sulfuric acid 10% v/v		1 l	502591	737
Sulfuric acid 4 mol/l (8N)		1 l	526741	737
Sulfuric acid 2.5 mol/l (5N)		1 l	P3240015	737

Product specifications are subject to changes. Please visit our website for updates.

Sulfuric acid 0.33 mol/l (2N/3)	1 l	410634	738
Sulfuric acid 0.26 mol/l (0.52N)	5 l	502202	739
Sulfuric acid 0.166 mol/l (0.333N)	1 l	PS0217/15	739
Sulfuric acid 0.13 mol/l (0.26N)	5 l	502651	740
Sulfuric acid 0.1 mol/l (0.2N)	1 l	502100000	740
Sulfuric acid d=1.820	5 l	502020	743
Wijs' reagent	250 ml	E491901	790
Wijs' reagent	1 l	E491902	790

Kjeldahl, Nitrogen content

For the determination of the nitrogen content using the Kjeldahl method, the following mineralization catalysts are available.

Description	Notes	Size	Code	Page
Kjeldahl antifoam	Composition: Sodium sulfate 0.97 g/Silicone antifoam 0.03 g	1000 x 1 g	502811	449
Kjeldahl catalyst according to Wieninger	Composition: Sodium sulfate 4.88 g/ Copper sulfate 0.07 g/ Selenium 0.05 g	1000 x 5 g	502821	450
Kjeldahl catalyst for water analysis	Composition: Potassium sulfate 5.0 g/Selenium 5 mg	1000 x 5 g	502121	450
Kjeldahl catalyst for water analysis	Composition: Potassium sulfate 5.0 g/Selenium 50 mg	1000 x 5 g	502122	450
Kjeldahl catalyst without selenium and titanium	Composition: Potassium sulfate 3.50 g/Copper sulfate 0.40 g	1000 x 3.9 g	502791	450
Kjeldahl catalyst without selenium and titanium	Composition: Potassium sulfate 5.0 g/Copper sulfate 0.50 g	1000 x 5 g	502792	450
Kjeldahl selenium catalyst	Composition: Potassium sulfate 4.63 g/Copper sulfate 0.28 g/Selenium 0.09g	1000 x 5 g	502120	450
Kjeldahl titanium catalyst	Composition: Potassium sulfate 3.5 g/Copper sulfate 0.105 g/Titane dioxide 0.105 g	1000 x 3.5 g	502123	451
Kjeldahl titanium catalyst	Composition: Potassium sulfate 5.00 g/Copper sulfate 0.15 g/Titane dioxide 0.15 g	500 x 5 g	502802	451

ELECTRONICS

CARLO ERBA Reagents line of products for the electronics industry, is characterized by high chemical purity. As they are subjected to meticulous filtration processes, these solvents guarantee a particulate content less than 250 ppm, for particles with diameters larger than 0.5 micron.

RSE, Special Solvents & Reagents

RSE (Special Reagents for Electronics) for all applications which do not require control of particle content.

Description	Notes	Size	Code	Page
Acetic acid glacial		1 l	401463	131
Acetic acid glacial		2.5 l	401462	131
Acetone		1 l	401051	139
Acetone		2.5 l	401058	139
Acetone		5 l	401054	139
Acetone		5 l	401055	139
Acetone		22 kg	401052	139
Ammonia solution 30%		1 l	420071	160
Ammonia solution 30%		2 l	420073	160
Ammonia solution 30%		5 l	420077	160
Ammonia solution 30%		25 kg	420075	160
Ammonia solution 25%		5 l	420085	162
Ammonia solution 25%		25 kg	420084	162
Dichloromethane		1 l	463162	303
Dichloromethane		2.5 l	463161	303
Ethanol absolute anhydrous	Only for Italian market	1 l	414587	332
Ethanol absolute anhydrous		1 l	4145872	332
Ethanol absolute anhydrous	Only for Italian market	2.5 l	414583	332
Ethanol absolute anhydrous		2.5 l	4145832	332
Ethyl acetate		1 l	448307	342
Ethyl acetate		2.5 l	448308	342
Ethyl acetate		5 l	448306	342
Glycerol (30°Bé)		1 l	453771	372
Glycerol (30°Bé)		2.5 l	453772	372
Hydrochloric acid 37%		1 l	403977	394
Hydrochloric acid 37%		2.5 l	403971	394
Hydrofluoric acid 50%		1 l	405737	410
Hydrogen peroxide solution 30%		1 l	412161	414
Hydrogen peroxide solution 30%		5 l	412162	414
Hydrogen peroxide solution 30%		25 kg	412163	414
Isopar G		2.5 l	526151	440
Methanol		1 l	414917	490
Methanol		2.5 l	414914	490
Nitric acid 69.5%		1 l	408097	528
Nitric acid 69.5%		2.5 l	408098	528
Nitric acid 65%		1 l	408101	531
Nitric acid 65%		2.5 l	408102	531
Nitric acid 18%		1 l	408191	532
Orthophosphoric acid 85%		1 l	406022	548
Orthophosphoric acid 85%		2.5 l	406021	548
Potassium hydroxide, pellets		1 kg	472097	599

Product specifications are subject to changes. Please visit our website for updates.

Potassium hydroxide solution 45%	5 l	472103	600
Propan-2-ol	1 l	415237	625
Propan-2-ol	2.5 l	415235	625
Propan-2-ol	5 l	415231	625
Propan-2-ol	5 l	415238	625
Propan-2-ol	27 l	415236	625
Propan-2-ol	200 l	415233	625
Sodium hydroxide, pellets	1 kg	480527	679
Sodium hydroxide, pellets	5 kg	480522	679
Sodium hydroxide, pellets	25 kg	480525	679
Sulfuric acid 96%	1 l	410374	731
Sulfuric acid 96%	2.5 l	410371	731
Xylene, mix of isomers	1 l	492358	791
Xylene, mix of isomers	2.5 l	492359	791

MOS (Metal Oxide Semiconductor) Solvents & Reagents

MOS (Metal Oxide Semiconductor) for MOS circuit production processes.

Description	Notes	Size	Code	Page
Acetone		1 l	401042	140
Acetone		2.5 l	401041	140
Ammonia solution 25%		1 l	420051	162
Ammonia solution 25%		2.5 l	420052	162
Hydrochloric acid 37%		1 l	403942	394
Hydrochloric acid 37%		2.5 l	403941	394
Hydrofluoric acid 50%		1 l	405653	410
Hydrogen peroxide solution 30%		1 l	412081	414
Methanol		1 l	414822	491
Methanol		2.5 l	414821	491
Nitric acid 69.5%		1 l	408151	528
Nitric acid 69.5%		2.5 l	408152	528
Propan-2-ol		1 l	415162	625
Propan-2-ol		2.5 l	415161	625
Sulfuric acid 96%		1 l	410382	732
Sulfuric acid 96%		2.5 l	410381	732

VLSI (Very Large Scale Integration) Solvents & Reagents

VLSI (Very Large Scale Integration) for the production of microcircuits using VLSI technology.

Description	Notes	Size	Code	Page
Acetone		1 l	527651	139
Acetone		2.5 l	527650	139
Acetone		5 l	527655	139
Ethanol absolute anhydrous		1 l	527681	332
Ethanol absolute anhydrous		2.5 l	527680	332
Hydrochloric acid 37%		1 l	527601	394
Hydrochloric acid 37%		2.5 l	527600	394
Hydrogen peroxide solution 30%		1 l	527621	414
Hydrogen peroxide solution 30%		2.5 l	527620	414
Methanol		1 l	527641	490
Methanol		2.5 l	527640	490

Nitric acid 69.5%	1 l	527671	528
Nitric acid 69.5%	2.5 l	527670	528
Orthophosphoric acid 85%	1 l	527592	547
Orthophosphoric acid 85%	2.5 l	527591	547
Propan-2-ol	1 l	527696	624
Propan-2-ol	2.5 l	527690	624
Propan-2-ol	30 l	527691	624
Sulfuric acid 96%	1 l	527631	731
Sulfuric acid 96%	2.5 l	527630	731

Product specifications are subject to changes. Please visit our website for updates.

HISTOLOGY, HEMATOLOGY AND CYTODIAGNOSTIC

A wide range of products for the sample preparation in histology, hematology and cytodiagnosis.

Fixing media

CARLO ERBA Reagents offers a wide range of fixatives, in concentrated or ready-to-use diluted form.

Description	Notes	Size	Code	Page
Fixative AFA liquid	60 ml jars filled at 30 ml. Box of 500	480 x 30 ml	508840	354
Fixative AFA liquid	In Vitro Diagnostic Medical Device	1 l	526267	354
Fixative AFA liquid	In Vitro Diagnostic Medical Device	5 l	526263001	354
Fixative Bouin Hollande liquid		1 l	526268	355
Fixative Bouin liquid	In Vitro Diagnostic Medical Device	1 l	526270	355
Fixative Bouin liquid	In Vitro Diagnostic Medical Device	5 l	526261	355
Fixative Bouin liquid	In Vitro Diagnostic Medical Device	25 l	526311	355
Fixative Bouin liquid	In Vitro Diagnostic Medical Device	1 l	429751	355
Fixative Davidson liquid	60 ml jars filled at 30 ml. Box of 500	30 ml	508881	355
Fixative Davidson liquid		5 l	526277	355
Fixative FIXALL-HIS liquid		5 l	526274	356
Fixative liquid without acetic acid		10 l	526264	356
Formaldehyde 37% w/v		1 l	415661	358
Formaldehyde 37% w/v		2.5 l	415666	358
Formaldehyde 37% w/v		5 l	415667	358
Formaldehyde 37% w/v neutralized		1 l	415686	359
Formaldehyde 37% w/v neutralized		5 l	415682	359
Formaldehyde 37% w/v neutralized		10 kg	415683	359
Formaldehyde 37% w/v neutralized		30 kg	415684	359
Formaldehyde 37% w/v neutralized		55 kg	415685	359
Formaldehyde 10% v/v according to Lillie		5 l	526912	360
Formaldehyde 10% v/v according to Lillie		25 l	526911	360
Formaldehyde 5% w/v buffered at pH 6.9		5 l	415674	360
Formaldehyde 5% w/v buffered at pH 6.9		10 l	415672	360
Formaldehyde 4% w/v buffered at pH 6.9	In Vitro Diagnostic Medical Device	1 l	415634	361
Formaldehyde 4% w/v buffered at pH 6.9	In Vitro Diagnostic Medical Device	5 l	415631	361
Formaldehyde 4% w/v buffered at pH 6.9	In Vitro Diagnostic Medical Device	10 l	415633	361
Formaldehyde 4% w/v buffered at pH 6.9	In Vitro Diagnostic Medical Device	20 l	415636	361
Formaldehyde 4% w/v buffered at pH 6.9	60 ml jars filled at 30 ml. Box of 500	30 ml	508861	361
Formaldehyde 4% w/v buffered at pH 6.9	500 ml jars filled at 300 ml. Box of 32	300 ml	508863	361
Formaldehyde 4% w/v buffered at pH 6.9		500 ml	524920	361
Formaldehyde 4% w/v buffered at pH 6.9	1 l jar filled at 800 ml	800 ml	526937	361
Formaldehyde 4% w/v buffered at pH 6.9		1 l	415694	361
Formaldehyde 4% w/v buffered at pH 6.9	5 l bucket filled at 2.5 l	2.5 l	526931	361
Formaldehyde 4% w/v buffered at pH 6.9		5 l	415691	361
Formaldehyde 4% w/v buffered at pH 6.9		5 l	415695	361
Formaldehyde 4% w/v buffered at pH 6.9		5 l	526936	361
Formaldehyde 4% w/v buffered at pH 6.9		10 l	415693	361
Formaldehyde 4% w/v buffered at pH 6.9		10 l	526933	361
Formaldehyde 4% w/v buffered at pH 6.9		20 l	415696	361
Formaldehyde 4% w/v buffered at pH 6.9		30 kg	415692	361
Formaldehyde 4% w/v with sodium chloride		1 l	526934	361
Formaldehyde acetic	60 ml jars filled at 30 ml. Box of 500	480 x 30 ml	508871	362

Formaldehyde acetic	1 l	526231	362
Formaldehyde acetic	5 l	526273	362
Paraformaldehyde	1 kg	387507	556

Solvents for dehydration, de-waxing and diaphanization

Dehydration is obtained by bathing the tissue in an alcoholic series of increasing concentration.
Other solvents can be used instead of alcohols, as long as they mix with water, xylene and paraffins.
Solvent Plus, a mixture of Isoparaffins is widely and efficiently used for diaphanization and deparaffinization as a substitute of xylene.
Histolemon is a natural, non-toxic solvent for histology.
It is used as both a diaphanizing agent and a deparaffinizing agent, and it can replace xylene, toluene and all other solvents used for the same purposes. It is a product of natural origin, citrus-fruit scented, with the following characteristics: totally non-toxic, scarcely volatile, mixable with alcohols, compatible with all kinds of paraffins. It completely solves the problem of toxicity in the workplace.

Description	Notes	Size	Code	Page
Ethanol absolute anhydrous	Only for Italian market	1 l	414601	332
Ethanol absolute anhydrous	Only for Italian market	1 l	414607	332
Ethanol absolute anhydrous		1 l	4146012	332
Ethanol absolute anhydrous		1 l	4146072	332
Ethanol absolute anhydrous	Only for Italian market	2.5 l	414605	332
Ethanol absolute anhydrous	Only for Italian market	2.5 l	414608	332
Ethanol absolute anhydrous		2.5 l	4146052	332
Ethanol absolute anhydrous		2.5 l	4146082	332
Ethanol absolute anhydrous		5 l	414603	332
Ethanol absolute anhydrous		5 l	414606	332
Ethanol absolute anhydrous	Untaxed, for Italian license holders only	5 l	4146032	332
Ethanol absolute anhydrous	Untaxed, for Italian license holders only	5 l	4146062	332
Ethanol absolute anhydrous		10 l	414604	332
Ethanol 96°	Only for Italian market	1 l	414634	335
Ethanol 96°	Only for Italian market	1 l	414637	335
Ethanol 96°		1 l	4146342	335
Ethanol 96°		1 l	4146372	335
Ethanol 96°	Only for Italian market	2.5 l	414631	335
Ethanol 96°	Only for Italian market	2.5 l	414632	335
Ethanol 96°		2.5 l	4146312	335
Ethanol 96°		2.5 l	4146322	335
Ethanol 96°		5 l	414635	335
Ethanol 96°	Untaxed, for Italian license holders only	5 l	4146352	335
Ethanol 96°		10 l	414638	335
Ethanol 70% v/v		2.5 l	308771	337
Ethanol 70% v/v	Untaxed, for italian license holders only	25 l	3087752	337
Ethanol absolute denaturated		1 l	528761	338
Ethanol absolute denaturated		2.5 l	528765	338
Ethanol absolute denaturated		5 l	528763	338
Ethanol absolute denaturated		5 l	528764	338
Ethanol absolute denaturated		10 l	528766	338
Ethanol absolute denaturated		25 l	528762	338
Ethanol 95° denaturated		1 l	528771	338
Ethanol 95° denaturated		5 l	528775	338
Ethanol 95° denaturated		10 l	528772	338
Ethanol 95° denaturated		25 l	528773	338
Ethanol 95° denaturated		200 l	528774	338
Histolemon		1 l	454911	390

Product specifications are subject to changes. Please visit our website for updates.

Histolemon		2.5 l	454912	390
Histolemon		5 l	454915	390
Methanol		1 l	414814	491
Methanol		1 l	414819	491
Methanol		2.5 l	414815	491
Methanol		2.5 l	414816	491
Methanol		5 l	524102	491
Methanol		5 l	524103	491
Methanol		10 l	414818	491
Mixture Ethanol 95° / Isopropanol	In Vitro Diagnostic Medical Device	5 l	414551	508
Mixture Ethanol 99° / Isopropanol	In Vitro Diagnostic Medical Device	5 l	414511	508
Propan-1-ol		1 l	415104	622
Propan-1-ol		2.5 l	415102	622
Propan-1-ol		10 l	415108	622
Propan-2-ol		1 l	415154	625
Propan-2-ol		2.5 l	415156	625
Propan-2-ol		2.5 l	415158	625
Propan-2-ol		5 l	415173	625
Propan-2-ol		5 l	529174	625
Propan-2-ol		10 l	415153	625
Solvent Plus		2.5 l	446187	714
Solvent Plus		5 l	446181	714
Tetrahydrofuran		1 l	487308	752
Tetrahydrofuran		2.5 l	487303	752
Tetrahydrofuran		5 l	487305	752
Tetrahydrofuran		5 l	487307	752
Toluene		1 l	488551	767
Toluene		2.5 l	488555	767
Toluene		5 l	488552	767
Water		1 l	307592	788
Water		2.5 l	307593	788
Water		5 l	307582	788
Water		10 l	307586	788
Xylene, mix of isomers		1 l	492301	791
Xylene, mix of isomers		2.5 l	492306	791
Xylene, mix of isomers		5 l	492305	791

Embedding media

CARLO ERBA Reagents paraffins, in addition to possessing the typical properties required of an inclusion medium, offer the advantages of being highly purified and filtered, with a melting point between 56 and 58°C.
In order to optimize the infiltration times and guarantee better preparation of the sample, paraffin is available without DMSO additive.

Description	Notes	Size	Code	Page
Paraffin 56-58°C - Erbaplast (without DMSO)	In Vitro Diagnostic Medical Device	4 x 2 kg	467958	555

Reagents and staining solutions

Ready-to-use staining solutions produce bright and well-contrasted colors.
Ideal for easy microscope readings, they allow valid results to be obtained with shorter staining times and greater certainty.
Solutions for Papanicolaou staining are especially suited for laboratory use in cytological diagnostics, thanks to their quality and reproducibility characteristics: bright and contrasted colors, high staining capacity, short staining time,long solution life, reproducibility of the colors, and perfect long-term preservation of the stained sections.
Panoptic staining achieved with Giemsa and May Grünwald reagents provides a greater body of information with a single staining. Indeed, blood smear staining allows the diagnosis of all pathologies related to malformations of red and white blood cells as well as distribution

imbalances of the different kinds of white blood cells.

Description	Notes	Size	Code	Page
Alcian Blue 8GS 1%		250 ml	428551	148
Amman's lactophenol solution		100 ml	457531	160
Benedict's reagent		1 l	E425742	197
Carbolated Methylene Blue hydroalcoholic solution		100 ml	428991	248
Carbolated Toluidine Blue hydroalcoholic solution		100 ml	429291	248
Crystal violet oxalate for Gram-Hucker Kit	In Vitro Diagnostic Medical Device	250 ml	491561	284
Decalcifying agent	In Vitro Diagnostic Medical Device	1 l	441221	290
Differentiator for kit Gram-Hucker	In Vitro Diagnostic Medical Device	250 ml	444131	311
Ehrlich's reagent		500 ml	E446302	327
Eosin Y 1% solution acqueous	In Vitro Diagnostic Medical Device	1 l	446644	328
Eosin Y 0.5% solution alcoholic	In Vitro Diagnostic Medical Device	1 l	446664	329
Esbach's reagent		1 l	446981	330
Gentian violet carbolated solution	Dye for microscopy (bacteriology) according to GRAM	250 ml	E491651	368
Gentian violet carbolated solution	Dye for microscopy (bacteriology) according to NICOLLE	250 ml	E491661	368
Giemsa's reagent	In Vitro Diagnostic Medical Device	100 ml	453614	369
Giemsa's reagent	In Vitro Diagnostic Medical Device	6 x 100 ml	E453612	369
Giemsa's reagent	In Vitro Diagnostic Medical Device	500 ml	453616	369
Giemsa's reagent	In Vitro Diagnostic Medical Device	6 x 500 ml	E453613	369
Giemsa's reagent	In Vitro Diagnostic Medical Device	2.5 l	453611	369
Giemsa's reagent	In Vitro Diagnostic Medical Device	4 x 2.5 l	E453615	369
Gram - Hucker Kit	In Vitro Diagnostic Medical Device	4 x 250 ml	454441	376
Haemalum solution according to Carazzi	In Vitro Diagnostic Medical Device	250 ml	434351	378
Haemalum solution according to Mayer	In Vitro Diagnostic Medical Device	1 l	446372	378
Haemalum solution according to Mayer	In Vitro Diagnostic Medical Device	1 l	446377	378
Hematoxylin solution according to Mayer	In Vitro Diagnostic Medical Device	100 ml	460511	379
Hematoxylin solution according to Mayer	In Vitro Diagnostic Medical Device	6 x 100 ml	460512	379
Hematoxylin solution according to Mayer	In Vitro Diagnostic Medical Device	1 l	460513	379
Hematoxylin solution according to Mayer	In Vitro Diagnostic Medical Device	6 x 1 l	460515	379
Lactophenol blue solution	In Vitro Diagnostic Medical Device	100 ml	428901	452
Lugol concentrated solution		1 l	458741	466
Lugol solution for Gram-Hucker kit	In Vitro Diagnostic Medical Device	250 ml	458751	466
Lugol's Reagent Iodine-Iodide Solution	Medical Device	250 ml	458762	466
Lugol's Reagent Iodine-Iodide Solution	Medical Device	6 x 250 ml	E458761	466
Lugol's Reagent Iodine-Iodide Solution	Medical Device	1 l	458763	466
May Grünwald reagent	In Vitro Diagnostic Medical Device	100 ml	460584	482
May Grünwald reagent	In Vitro Diagnostic Medical Device	6 x 100 ml	E460582	482
May Grünwald reagent	In Vitro Diagnostic Medical Device	500 ml	460586	482
May Grünwald reagent	In Vitro Diagnostic Medical Device	6 x 500 ml	E460583	482
May Grünwald reagent	In Vitro Diagnostic Medical Device	2.5 l	460581	482
May Grünwald reagent	In Vitro Diagnostic Medical Device	4 x 2.5 l	E460585	482
Mayer's reagent		500 ml	460502	482
Methylene blue saturated solution		250 ml	E429031	498
Papanicolaou Haematoxylin solution according to Harris	In Vitro Diagnostic Medical Device	500 ml	446462	554
Papanicolaou Haematoxylin solution according to Harris	In Vitro Diagnostic Medical Device	6 x 500 ml	446464	554

Product specifications are subject to changes. Please visit our website for updates.

Papanicolaou Haematoxylin solution according to Harris	In Vitro Diagnostic Medical Device	1 l	446461	554
Papanicolaou Haematoxylin solution according to Harris	In Vitro Diagnostic Medical Device	6 x 1 l	446465	554
Papanicolaou Haematoxylin solution according to Harris	In Vitro Diagnostic Medical Device	2.5 l	446463	554
Papanicolaou Haematoxylin solution according to Harris	In Vitro Diagnostic Medical Device	4 x 2.5 l	446466	554
Papanicolaou solution EA 50	In Vitro Diagnostic Medical Device	500 ml	467782	554
Papanicolaou solution EA 50	In Vitro Diagnostic Medical Device	6 x 500 ml	E467784	554
Papanicolaou solution EA 50	In Vitro Diagnostic Medical Device	1 l	467781	554
Papanicolaou solution EA 50	In Vitro Diagnostic Medical Device	6 x 1 l	E467785	554
Papanicolaou solution EA 50	In Vitro Diagnostic Medical Device	2.5 l	467783	554
Papanicolaou solution EA 50	In Vitro Diagnostic Medical Device	4 x 2.5 l	E467786	554
Papanicolaou solution OG 6	In Vitro Diagnostic Medical Device	500 ml	467792	555
Papanicolaou solution OG 6	In Vitro Diagnostic Medical Device	6 x 500 ml	E467794	555
Papanicolaou solution OG 6	In Vitro Diagnostic Medical Device	1 l	467791	555
Papanicolaou solution OG 6	In Vitro Diagnostic Medical Device	6 x 1 l	E467795	555
Papanicolaou solution OG 6	In Vitro Diagnostic Medical Device	2.5 l	467793	555
Papanicolaou solution OG 6	In Vitro Diagnostic Medical Device	4 x 2.5 l	E467796	555
Safranine T hydroalcoholic solution for Gram-Hucker Kit	In Vitro Diagnostic Medical Device	250 ml	477241	642
Schiff's reagent for PAS coloration	In Vitro Diagnostic Medical Device	500 ml	477591	646
Schiff's reagent for PAS coloration	In Vitro Diagnostic Medical Device	6 x 500 ml	477592	646
Tauber reagent		500 ml	490422	746
Turk's reagent	In Vitro Diagnostic Medical Device	500 ml	E490451	780
Wright's stain solution in methanol	In Vitro Diagnostic Medical Device	100 ml	492011	790
Ziehl-Neelsen's reagent	In Vitro Diagnostic Medical Device	250 ml	493101	796
Ziehl-Neelsen's reagent	In Vitro Diagnostic Medical Device	1 l	493102	796

Dyes

To aid identification of the finest details of cellular structures, very pure products must be used. For this purpose the CARLO ERBA Reagents dyes in powder, in addition to possessing such characteristics, are also classified and certified through the combined use of TLC and UV-Visible spectrophotometry.

Description	Notes	Size	Code	Page
Alcian blue 8GX		25 g	428561	149
Alkali blue 6B		25 g	428532	150
Aniline blue soluble in water		25 g	428582	180
Azure II		5 g	424721	189
Azure II eosin		5 g	424731	190
Bismarck brown R		25 g	431252	201
Biuret 97%		25 g	428432	205
Brillant cresyl blue		10 g	428811	208
Brillant cresyl blue		25 g	428812	208
Brilliant green		25 g	491152	209
Chrysoidine Y		25 g	440572	268
Congo red		25 g	476762	273
Congo red		100 g	476764	273
Coomassie brilliant blue R 250		25 g	428642	274
Eosin B		25 g	446602	328
Eosin Y		25 g	446632	328
Eosin Y		100 g	446634	328

Erythrosin extra B	25 g	446972	330
Erythrosin extra B	100 g	446971	330
Fast green FCF	25 g	491391	353
Fuchsin acid	25 g	452812	364
Fuchsin acid	100 g	452814	364
Fuchsin basic	25 g	452842	365
Fuchsin basic	100 g	452844	365
Gentian violet	25 g	388703	368
Gentian violet	50 g	388701	368
Gold(III) chloride trihydrate	1 g	467007	375
Hematoxylin	25 g	446472	379
Hematoxylin	100 g	446473	379
Hematoxylin	1 kg	446475	379
Light green	10 g	491371	461
Light green	25 g	491372	461
Malachite green	25 g	491303	476
Malachite green	100 g	491304	476
Methyl blue	25 g	428932	497
Methyl green	10 g	491351	499
Methyl green	25 g	491352	499
Naphthol yellow S	25 g	453562	518
Neutral red	10 g	476951	521
Nigrosine	25 g	464852	525
Nigrosine	50 g	464853	525
Nuclear fast red	10 g	477011	536
Nuclear fast red	25 g	477012	536
Orange G	25 g	423432	539
Orange II	10 g	423341	539
Orcein	5 g	466858	540
Orcein	25 g	466859	540
Phloxin B	10 g	452051	574
Phloxin B	25 g	452052	574
Ponceau red BS	10 g	476941	581
Ponceau red S	5 g	476981	581
Ponceau red S	25 g	476982	581
Rosolic acid	25 g	409702	640
Safranine T	25 g	477232	642
Sudan black B	25 g	464241	726
Sudan III	25 g	485902	726
Sudan yellow	10 g	453581	727
Tartrazine	50 g	486903	746
Toluidine blue	25 g	429282	769
Victoria blue	10 g	429381	785
Victoria blue	25 g	429382	785

Mounting media

Mounting media of natural or synthetic origin in various solvents, depending on the specific needs of the users, with all the typical characteristics required of common mounting media.

Description	Notes	Size	Code	Page
Canada balsam		100 g	321553	247
Canada balsam		250 g	321554	247
Eukitt		100 ml	554194	351
Eukitt		250 ml	554193	351
Eukitt		500 ml	554192	351
Oil of cedar wood		100 ml	466753	538
Oil of cedar wood		1 l	466757	538

Immersion media

An especially useful range of products for optical microscopy.

Description	Notes	Size	Code	Page
Immersion oil		100 ml	466782	420
Immersion oil		1 l	466783	420

ORGANIC CHEMISTRY

Anhydrous solvents

In analytical and synthetic chemistry many operations must take place in an anhydrous environment in order to prevent alterations, the occurrence of side reactions and a reduced yield. For this purpose, high-quality solvents with minimum water content must be used. CARLO ERBA Reagents anhydrous solvents eliminate the time-consuming need for anhydrification of the solvents while ensuring results that are certain and reliable over time. These solvents are the result of specific production processes, optimized and controlled to obtain the highest degree of purity and minimum water content. The packaging materials used are ideal for preserving product quality, and the products are available in various sizes, in septum bottles or glass bottles.

Description	Notes	Size	Code	Page
Acetone		200 ml	P0051010	139
Acetone		1 l	P0051016	139
Acetonitrile		200 ml	P0061010	144
Acetonitrile	Water content < 50 ppm	200 ml	P00610S10	144
Acetonitrile		1 l	P0061016	144
Acetonitrile	Water content < 50 ppm	1 l	P00610S16	144
Acetonitrile		2.5 l	P0061021	144
Acetonitrile	Water content < 50 ppm	2.5 l	P00610S21	144
Butanol-1		200 ml	P0171010	226
Butanol-1		1 l	P0171016	226
tert-Butanol		1 l	P0191016	227
n-Butyl acetate		1 l	P0011016	228
n-Butyl acetate		2.5 l	P0011021	228
tert-Butylmethylether		1 l	P0921016	231
Chloroform		200 ml	P02410A10	261
Chloroform		1 l	P02410A16	261
Chloroform		2.5 l	P02410A21	261
Chloroform		200 ml	P02410E10	261
Chloroform		2.5 l	P02410E21	261
Cyclohexane		200 ml	P0251010	286
Cyclohexane		1 l	P0251016	286
1,2-Dichloroethane		200 ml	P0281010	299
1,2-Dichloroethane		1 l	P0281016	299
1,2-Dichloroethane		2.5 l	P0281021	299
Dichloromethane, stab. with Amylene		200 ml	P02910A10	302
Dichloromethane, stab. with Amylene		1 l	P02910A16	302
Dichloromethane, stab. with Amylene		2.5 l	P02910A21	302
Dichloromethane, stab. with Ethanol		200 ml	P02910E10	302
Dichloromethane, stab. with Ethanol		1 l	P02910E16	302
Dichloromethane, stab. with Ethanol		2.5 l	P02910E21	302
Diethyl ether		200 ml	P0441010	309
Diethyl ether		1 l	P0441008	309
Diethyl ether		1 l	P0441016	309
Diethyl ether		2.5 l	P0441021	309
Diisopropylether		1 l	P0431016	312
N,N-Dimethylformamide		200 ml	P0341010	315
N,N-Dimethylformamide		1 l	P0341016	315
N,N-Dimethylformamide		2.5 l	P0341021	315
1,4-Dioxane		200 ml	P0361010	321
1,4-Dioxane		1 l	P0361016	321

Product specifications are subject to changes. Please visit our website for updates.

1,4-Dioxane	2.5 l	P0361021	321
Ethanol absolute anhydrous	200 ml	P013A1010	331
Ethanol absolute anhydrous	1 l	P013A1016	331
Ethanol absolute anhydrous	2.5 l	P013A1021	331
Ethyl acetate	200 ml	P0021010	341
Ethyl acetate	1 l	P0021016	341
Ethyl acetate	2.5 l	P0021021	341
Ethyl methyl ketone	1 l	P0201016	350
Formamide	200 ml	P6151010	362
n-Heptane 99%	1 l	P0501016	381
n-Heptane 99%	2.5 l	P0501021	381
n-Hexane	1 l	P0521016	387
Isobutanol	1 l	P0531016	437
Isopentane	1 l	P0651016	441
Methanol	200 ml	P0931010	490
Methanol	1 l	P0931016	490
Methanol	2.5 l	P0931021	490
Methylcyclohexane	1 l	P0581016	497
Methyl isobutyl ketone	1 l	P0601016	500
N-Methyl-2-pyrrolidone	200 ml	P0871010	502
n-Pentane	1 l	P0641016	559
Propan-1-ol	1 l	P0941016	622
Propan-1-ol	2.5 l	P0941021	622
Propan-2-ol	200 ml	P0951010	624
Propan-2-ol	1 l	P0951016	624
Pyridine	200 ml	P0671010	630
Pyridine	1 l	P0671016	630
Pyridine	2.5 l	P0671021	630
Sulfolane	1 l	P932SP16	729
Tetrahydrofuran	200 ml	P0701010	752
Tetrahydrofuran	1 l	P0701016	752
Tetrahydrofuran	2.5 l	P0701021	752
Toluene	200 ml	P0711010	766
Toluene	1 l	P0711016	766
Toluene	2.5 l	P0711021	766

Anhydrous solvents, with molecular sieves

To improve prevention of contamination from external humidity, these solvents are supplied with molecular sieves.

Description	Notes	Size	Code	Page
Acetonitrile	On molecular sieves 3A	200 ml	P00610T10	144
Chloroform	On molecular sieves 4A	2.5 l	P02410AT21	261
Dichloromethane	On molecular sieves 4A	200 ml	P02910AT10	302
Dichloromethane	On molecular sieves 4A	1 l	P02910AT16	302
Diethyl ether	On molecular sieves 4A, Water content < 20ppm	200 ml	P04410T10	309
Diethyl ether	On molecular sieves 4A, Water content < 20ppm	1 l	P04410T16	309
N,N-Dimethylformamide	On molecular sieves 4A	200 ml	P03410T10	315
N,N-Dimethylformamide	On molecular sieves 4A	1 l	P03410T16	315

Ethyl acetate	On molecular sieves 4A, Water content < 20ppm	2.5 l	P00210T21	341
Tetrahydrofuran	On molecular sieves 4A	200 ml	P07010T10	752
Tetrahydrofuran	On molecular sieves 4A	1 l	P07010T16	752
Toluene	On molecular sieves 4A	200 ml	P07110T10	766

NMR SPECTROSCOPY

To meet the various needs of NMR spectroscopy, CARLO ERBA Reagents offers a wide range of deuterated solvents, characterized by high production standards resulting from the use of spectroscopically pure raw materials and accurate analyses and controls in line with our long tradition of quality.
Our solvents are available with various degrees of isotopic purity and a vast assortment of packages, with an attention to detail that complements the overall reliability and safety guaranteed by our products.

NMR solvents

The choice of deuterated solvents is very important in NMR spectroscopy.
These solvents contribute to shielding the substance they dissolve, as well as making a significant contribution to expanding the electric reaction field generated by the solvent as a result of the polarization produced by the molecules in the solute.
To meet the varied needs of this sophisticated analytical technique, CARLO ERBA Reagents offers a wide range of deuterated solvents, characterized by high production standards resulting from the use of spectroscopically pure raw materials and accurate analyses and controls in line with our long tradition of quality.
Our solvents are available in various degrees of isotopic purity and a vast assortment of packages, with an attention to detail that complements the overall reliability and safety guaranteed by the products.

Description	Notes	Size	Code	Page
Acetic acid-d4		10 x 0.75 ml	P5039	137
Acetone-d6		10 x 0.6 ml	P5060	141
Acetone-d6		10 x 0.75 ml	P5049	141
Acetone-d6		10 ml	P5044A	141
Acetone-d6		25 ml	P5045	141
Acetone-d6		100 ml	P5046	141
Acetonitrile-d3		2 x 0.6 ml	P5070	145
Acetonitrile-d3		10 x 0.75 ml	P5079	145
Acetonitrile-d3		5 ml	P5073A	145
Benzene-d6		10 x 0.75 ml	P5089	198
Benzene-d6		25 ml	P5085	198
Benzene-d6		100 ml	P5086	198
Chloroform-d		10 x 0.6 ml	P5130	263
Chloroform-d		100 ml	P5505	263
Chloroform-d		10 x 0.75 ml	P5119	263
Chloroform-d		25 ml	P5115	263
Chloroform-d		100 ml	P5116	263
Chloroform-d		500 ml	P5117	263
Chloroform-d		1 l	P5118	263
Chloroform-d		100 ml	P5325	263
Chloroform-d + 0.03% TMS		100 ml	P5006	263
Cyclohexane-d12		2 x 0.5 ml	P5151A	287
Deuterium oxide-d2		10 x 0.75 ml	P5179	293
Deuterium oxide-d2		25 ml	P5175	293
Deuterium oxide-d2		10 x 0.75 ml	P5169	293
Deuterium oxide-d2		5 x 10 ml	P5164	293
Deuterium oxide-d2		25 ml	P5165	293
Deuterium oxide-d2		25 ml	P5165S	293
Deuterium oxide-d2		100 ml	P5166	293
Deuterium oxide-d2		1 l	P5168	293
Deuterium oxide-d2 + 0.01% DMSO		10 x 0.6 ml	P5170D	294
Deuterium oxide-d2 + 0.5% TSP d4		10 x 0.6 ml	P5161T	294
Deuterium oxide-d2 + 0.03% TSP d4		10 x 0.6 ml	P5160T	294
1,2-Dichlorobenzene-d4		5 ml	P5533A	299
Dichloromethane-d2		10 x 0.6 ml	P5330	305

Dichloromethane-d2	10 x 0.75 ml	P5339	305
Dichloromethane-d2	25 ml	P5335	305
N,N-Dimethylformamide-d7	2 x 0.75 ml	P5189A	316
Dimethylsulphoxide-d6	10 x 0.6 ml	P5220	319
Dimethylsulphoxide-d6	10 x 0.75 ml	P5229	319
Dimethylsulphoxide-d6	10 x 0.6 ml	P5200	319
Dimethylsulphoxide-d6	10 x 0.75 ml	P5209	319
Dimethylsulphoxide-d6	10 ml	P5204A	319
Dimethylsulphoxide-d6	5 x 10 ml	P5204S	319
Dimethylsulphoxide-d6	25 ml	P5205	319
Dimethylsulphoxide-d6	100 ml	P5206	319
Dimethylsulphoxide-d6 + 0.03% TMS	10 x 0.75 ml	P5541	319
Dimethylsulphoxide-d6 + 0.03% TMS	25 ml	P5545	319
Dimethylsulphoxide-d6 + 0.03% TMS	10 x 0.6 ml	P5602	320
Dimethylsulphoxide-d6 + 0.03% TMS	25 ml	P5605	320
Ethanol-d6 anhydrous	2 x 1 ml	P5262A	340
Formic acid-d	5 ml	P5733	364
Hydrochloric acid-d 20%	25 ml	P5685	409
Hydrochloric acid-d 1 mol/l	25 ml	P5695	409
Methanol-d4	10 x 0.6 ml	P5310	493
Methanol-d4	10 x 0.75 ml	P5319	493
Methanol-d4	10 x 0.6 ml	P5280	493
Methanol-d4	10 x 0.75 ml	P5289	493
Methanol-d4	5 ml	P5283A	493
Methanol-d4	5 x 10 ml	P5284	493
Methanol-d4	5 x 10 ml	P5284S	493
Methanol-d4	25 ml	P5285	493
Methanol-d4 + 0.03% TMS	10 x 0.6 ml	P5140	494
Methanol-d3	10 x 0.75 ml	P5309	494
Methanol-d1	25 ml	P5275	494
Orthophosphoric acid-d3 85% in D2O	25 ml	P5055	549
Pyridine-d5	2 x 0.6 ml	P5370	631
Pyridine-d5	2 x 0.75 ml	P5369A	631
Pyridine-d5	10 ml	P5364A	631
Sodium hydroxide-d 1 30%	25 ml	P5675	691
Sodium hydroxide-d 1 mol/l	25 ml	P5665	691
Tetrachloroethane-d2	25 ml	P5435	751
Tetrahydrofuran-d8	2 x 0.6 ml	P5380	753
Tetrahydrofuran-d8	25 ml	P5385	753
Toluene-d8	2 x 0.75 ml	P5399A	768
Toluene-d8	5 ml	P5393A	768
Toluene-d8	25 ml	P5395	768
Trifluoroacetic acid-d	2 x 0.75 ml	P5419A	774
Trifluoroacetic acid-d	5 ml	P5413A	774

UV SPECTROSCOPY

In molecular structure research and equilibrium studies, as well as studies on kinetics and steric effects, it is important to use solvents with high UV transmittance which contain no interfering substances absorbent in the IR spectral band and are suited for fluorescence analysis.

SPECTROSOL® Solvents for optical spectroscopy

CARLO ERBA Reagents offers a line of special solvents for spectroscopy. These are high-purity products obtained through specific production processes, controlled and packaged in order to adequately meet the needs of modern ultraviolet, infrared and fluorescence analytical techniques.

Description	Notes	Size	Code	Page
Acetone		1 l	401034	139
Acetone		2.5 l	401032	139
Acetonitrile		1 l	401216	143
Acetonitrile		2.5 l	401212	143
tert-Butylmethylether		1 l	432001	230
tert-Butylmethylether		2.5 l	432002	230
Chloroform, stab. with Amylene		1 l	438591	260
Chloroform, stab. with Amylene		2.5 l	438592	260
Chloroform, stab. with Ethanol		1 l	438664	261
Chloroform, stab. with Ethanol		2.5 l	438662	261
Cyclohexane		1 l	436967	286
Cyclohexane		2.5 l	436963	286
1,2-Dichloroethane		1 l	P0282716	299
Dichloromethane, stab. with Amylene		1 l	442371	302
Dichloromethane, stab. with Amylene		2.5 l	P02927A21	302
Dichloromethane, stab. with Ethanol		1 l	463025	302
Diethyl ether		1 l	447593	308
N,N-Dimethylformamide		1 l	444957	315
N,N-Dimethylformamide		2.5 l	444956	315
Dimethylsulphoxide		1 l	445112	318
Dimethylsulphoxide		2.5 l	445111	318
Ethanol absolute anhydrous	Only for Italian market	1 l	414677	331
Ethanol absolute anhydrous		1 l	4146772	331
Ethanol 96°	Only for Italian market	1 l	414667	334
Ethanol 96°		1 l	4146672	334
Ethyl acetate		1 l	448271	341
Ethyl acetate		2.5 l	448272	341
n-Heptane 99%		1 l	446824	380
n-Heptane 99%		2.5 l	P0502721	380
n-Hexane 99%		1 l	447051	385
n-Hexane 99%		2.5 l	447052	385
n-Hexane		1 l	446934	387
n-Hexane		2.5 l	446932	387
Isooctane		1 l	456754	439
Isooctane		2.5 l	456753	439
Methanol		1 l	414902	490
Methanol		2.5 l	414903	490
Methylcyclohexane		1 l	P0582716	497
Paraffin oil		100 ml	466792	555
n-Pentane 99%		1 l	468142	557
n-Pentane 99%		2.5 l	468141	557

Potassium bromide	100 g	470701	585
Propan-2-ol	1 l	415213	624
Propan-2-ol	2.5 l	P0952721	624
Tetrachloroethylene	1 l	P0682716	750
Tetrachloroethylene	2.5 l	P0682721	750
Tetrahydrofuran	1 l	487345	752
Tetrahydrofuran	2.5 l	487346	752
Toluene	1 l	488601	766
Toluene	2.5 l	488602	766
1,2,4-Trichlorobenzene	2.5 l	P0722721	771
Trifluoroacetic acid	1 l	P0082746	773
Trifluoroacetic acid	2.5 l	P0082747	773

Product specifications are subject to changes. Please visit our website for updates.

PEPTIDE SYNTHESIS

CARLO ERBA Reagents propose a specific range of the most commonly used solvents for peptide and/or DNA synthesis. These solvents are tested specifically to guarantee the absence of amines.

Solvents for Peptides Synthesis

Description	Notes	Size	Code	Page
Acetonitrile		200 ml	P0063510	144
Acetonitrile		1 l	P0063516	144
Acetonitrile		2.5 l	P0063521	144
N,N-Dimethylformamide		1 l	P0343516	316
N,N-Dimethylformamide		2.5 l	P0343521	316
N,N-Dimethylformamide		5 l	P0343522	316
N,N-Dimethylformamide		10 l	P0343541	316
N,N-Dimethylformamide		25 l	P0343549	316
N,N-Dimethylformamide		25 l	P0343550	316
N,N-Dimethylformamide		200 l	P0343567	316
N-Methyl-2-pyrrolidone		1 l	P0873516	502
N-Methyl-2-pyrrolidone		2.5 l	P0873521	502
N-Methyl-2-pyrrolidone		10 l	P0873541	502
N-Methyl-2-pyrrolidone		25 l	P0873549	502
N-Methyl-2-pyrrolidone		200 l	P0873566	502
Piperidine		500 ml	P0663518	580
Piperidine		1 l	P0663516	580
Piperidine		2.5 l	P0663521	580
Pyridine		1 l	P0673516	630
Pyridine		2.5 l	P0673521	630
Trifluoroacetic acid		100 ml	P0082103	773
Trifluoroacetic acid		2.5 l	P0082147	773

GREEN CHEMISTRY

Growing awareness of the environmental consequences of chemical products and the processes by which they are produced has led to the development of the concept of "Sustainable (Green) Chemistry" in the United States during the early nineties. The definition given by its founder, Paul T. Anastas, is the following:
"Green Chemistry is the utilization of a set of principles that reduces or eliminates the use and generation of hazardous substances in the design, manufacture and application of chemical products."
It is based on 12 principles, which take into consideration the environmental, economic and safety aspects of chemistry.

Green solvents

In the interest of Green Chemistry, CARLO ERBA Reagents offers the following Green Solvents.

Description	Notes	Size	Code	Page
Cyclopentyl methyl ether		1 l	P8010216	289
Cyclopentyl methyl ether		5 l	P8010229	289
N,N'-Dimethylpropylene urea		500 ml	P8020218	317
N,N'-Dimethylpropylene urea		1 l	P8020216	317
N,N'-Dimethylpropylene urea		5 l	P8020229	317
N,N'-Dimethylpropylene urea		25 l	P8020248	317
N,N'-Dimethylpropylene urea		200 l	P8020268	317
1,3-Dioxolane		1 l	P8030216	322
1,3-Dioxolane		5 l	P8030222	322
1,3-Dioxolane		25 l	P8030249	322
1,3-Dioxolane		200 l	P8030268	322
2-Methyltetrahydrofuran		1 l	P9960216	505
2-Methyltetrahydrofuran		2.5 l	P9960221	505
2-Methyltetrahydrofuran		5 l	P9960229	505
2-Methyltetrahydrofuran		25 l	P9960248	505
2-Methyltetrahydrofuran		200 l	P9960268	505
4-Methyltetrahydropyran		500 ml	P9990218	506
4-Methyltetrahydropyran		1 l	P9990216	506
4-Methyltetrahydropyran		2.5 l	P9990221	506
1,3-Propanediol		1 l	P8040216	627
1,3-Propanediol		5 l	P8040222	627
1,3-Propanediol		190 l	P8040268	627

Product specifications are subject to changes. Please visit our website for updates.

PETROCHEMICAL ANALYSIS

CARLO ERBA Reagent offers a new range of ASTM solvents for testing engines. These solvents can be used in RON and MON methods (D2699 and D2700).

Solvents for the Octane Number

The determination of the octane number of a gasoline is a measure of its quality and performence as a fuel. CARLO ERBA Reagents offers solvents for the measure of the octane rate.

Description	Notes	Size	Code	Page
n-Heptane 99%		5 l	524263	381
n-Heptane 99%		25 l	524265	381
n-Heptane 99%		200 l	524267	381
Isooctane		5 l	528960	440
Isooctane		25 l	528961	440
Octane 80 blend		5 l	525992	537
Octane 80 blend		25 l	525993	537
Octane 80 blend		140 kg	525994	537
Toluene		5 l	386102	767
Toluene		25 l	386104	767
Toluene		180 kg	386106	767

Reagents for TAN & TBN Analysis

The determination of TAN, the measure of acid concentration and of TBN, the measure of alkaline concentration present in a lubricant are essential tests in industrial machinery and engine applications. CARLO ERBA Reagents offers reagents for TAN and TBN analysis.

Description	Notes	Size	Code	Page
Reagent TAN		2.5 l	PS0327/21	635
Reagent TAN		5 l	PS0327/29	635
Reagent TAN		10 l	PS0327/39	635
Reagent TBN ASTM D2896		2.5 l	PS0423/21	636
Reagent TBN ASTM D2896		5 l	PS0423/29	636
Reagent TBN ASTM D2896		10 l	PS0423/39	636
Reagent TBN ASTM D4739		5 l	526615	636

KARL FISCHER TITRATION

ERBAqua® is the CARLO ERBA Reagents brand, for its complete range of pyridine-free reagents for the volumetric and coulometric Karl Fischer determination of water.
The main features of this range are: more safety due to py-free and one component reagents non-hazardous, fast and stable endpoints and long term titre stability.

ERBAqua® One component volumetric reagents

In one-component volumetric Karl Fischer Titration, the titrant contains all the reagents required by the reaction: iodine, sulfur dioxide, base and an alcohol.
Available in two different titer strengths, 5 mg/ml and 2 mg/ml, they are suitable for routine analysis, and thanks to their methanol-free formulation, they can be used also if the sample contains aldehydes and ketones.

Description	Notes	Size	Code	Page
Karl Fischer reagent 1 component 2 mg H2O/ml		1 l	570021	443
Karl Fischer reagent 1 component 5 mg H2O/ml		1 l	570011	443
Karl Fischer titrant 1 component 5 mg H2O/ml for aldehydes and ketones		1 l	570081	443

ERBAqua® Solvents for One component volumetric reagents

Used to dissolve the sample, when a one component reagent is used for titration. In addition to methanol and anhydrous chloroform, other solvents are available:
• 570031, for sample containing oils and fats
• 570041, for sample containing aldehydes and ketones. It is especially suitable for high molecular weight products and non polar constituents mixture

Description	Notes	Size	Code	Page
Chloroform		1 l	P02410E16	261
Karl Fischer solvent for aldehydes and ketones one component		1 l	570041	444
Karl Fischer solvent for oils one component		1 l	570031	444
Methanol	Water content max 50 ppm	1 l	414981	490

ERBAqua® Two components volumetric reagents

For users who do frequent Karl Fischer analysis and need a higher degree of accuracy than one-component volumetric titration can provide, we recommend the two components range.
They give more accuracy in the results and longer shelf life because the reagents required by the Karl Fischer reaction are separated between the titrant (iodine) and the working medium (sulfur dioxide and base), but this implies that they must be coupled used.

Description	Notes	Size	Code	Page
Karl Fischer titrant 2 component 2 mg H2O/ml		1 l	570061	444
Karl Fischer titrant 2 component 5 mg H2O/ml		1 l	570051	444

ERBAqua® Solvents for Two component volumetric reagents

Together with the two-component titrants, there are different types of solvents.
For oils or other non-polar compounds, 570101 is the suitable working medium. In case of need for an extra buffering capacity of 5 mmole of acid/ml, 570111 is the suitable working medium.

Description	Notes	Size	Code	Page
Karl Fischer solvent 2 component		1 l	570071	445
Karl Fischer solvent 2 component for aldehydes and ketones - Methanol free		1 l	570091	445
Karl Fischer solvent for oils 2 component		1 l	570101	445
Karl Fischer 2 component buffered solvent		500 ml	570111	445

Product specifications are subject to changes. Please visit our website for updates.

ERBAqua® Coulometric reagents - Anolite solutions for cells with diaphragm

Coulometric Karl Fischer method is highly effective and suitable when water content in the sample is less than 0.1% and high accurancy on the result is needed.
Designed for use with all titrators and titration cell types, a wide range is available in order to allow the necessary reagents both when a diaphragm or a diaphragmless cell is used.
For units with a diaphragm, two reagent solutions are required: an anolyte and a catolyte.

Description	Notes	Size	Code	Page
Karl Fischer anolyte solution		500 ml	570121	446
Karl Fischer anolyte solution for ketones and aldehydes - Methanol free		500 ml	570161	446
Karl Fischer anolyte solution - CFC free		500 ml	570141	446
Karl Fischer anolyte solution, oven		500 ml	570151	446
Karl Fischer anolyte solution for oils		500 ml	570171	447

ERBAqua® Coulometric reagents - Catolite solutions for cells with diaphragm

When coulometric titration is performed with a cell with a diaphragm there is a need for anolyte or catholyte solutions.

Description	Notes	Size	Code	Page
Karl Fischer catholyte solution		125 ml	570181	447
Karl Fischer catholyte solution for aldehydes and ketones		125 ml	570191	447

ERBAqua® Coulometric reagents - Anolite solutions for cells without diaphragm

When the coulometric titration is done with a cell without a diapragm, an anolyte solution is sufficient.

Description	Notes	Size	Code	Page
Karl Fischer anolyte solution for cells with and without diagrams		500 ml	570131	447

ERBAqua® Water standards - Gravimetric

Standardization of a Karl Fischer reagent is necessary in order to determine its water equivalency.
Gravimetric standards with content in water in mg/g.

Description	Notes	Size	Code	Page
Karl Fischer water standard 10.0 mg/g		10 x 5 ml	570221	448
Karl Fischer water standard 1.0 mg/g		10 x 5 ml	570211	448
Karl Fischer water standard 0.10 mg/g		10 x 5 ml	570201	449
Sodium tartrate dihydrate		100 g	483561	709

ERBAqua® Water standards - Volumetric

Volumetric standards with content in water in mg/ml.

Description	Notes	Size	Code	Page
Karl Fischer water standard 5.0 mg/ml		10 x 5 ml	570231	448

INDICATORS

Indicators represent a practical and important tool for monitoring the progress of a reaction of an aqueous solution, operations which are often essential for obtaining correct analytical data.

Indicators, for UV-fluorescence, Redox, Precipitation and Complexometry

- UV Fluorescence indicators
  The use of chromatic indicators may not be equally effective with turbid or colored solutions, or when the change in concentration is not fast enough. In these cases, instrumental methods or fluorescence indicators may be adopted to identify the end point.
- Oxidation-reduction indicators
  These indicators are substances that vary in color depending on whether they are in oxidized or reduced form. Their behavior is very similar to that of the indicators used in acid-base titration; however, while the latter are sensitive to changes in the solution's pH, oxidation-reduction indicators are sensitive to changes in the system's potential. The color changes are usually very clear and well-defined.
- Precipitation indicators
  Precipitation titration methods have very limited applications compared to other types of volumetric analysis, but the few that are still employed are very useful in practical terms.
  The titration process is based on the formation of an insoluble compound between the titrating agent and the substance being titrated, which gradually results from the reaction that occurs during titration.
  Precipitation indicators allow visual identification of the end point of titration thanks to a change in color, which corresponds to the variation of a key characteristic, such as a change in the precipitate's electric charge (isoelectric point).
- Complexometry indicators
  These are organic colorings, mainly of the azo group, which form stable complexes with metals and are characterized by different colors depending on whether they are in free form or complex form in the solution.

Description	Notes	Size	Code	Page
Acridine orange		25 g	423461	147
Alizarin		25 g	415892	149
Alizarin red		25 g	416002	149
Alizarin yellow R		10 g	453451	150
Anthrone		25 g	423282	182
Arsenazo III		1 g	424281	184
Arsenazo III		25 g	424282	184
Azomethine H		10 g	424691	189
Azomethine H		25 g	424692	189
Bromocresol green		1 g	491207	211
Bromocresol green		25 g	491208	211
Bromocresol purple		5 g	470038	212
Bromocresol purple		25 g	470039	212
Bromophenol blue		5 g	428658	213
Bromophenol blue		25 g	428659	213
Bromophenol blue		50 g	428653	213
Bromophenol blue		500 g	428655	213
Bromophenol blue indicator		1 l	PS0269/15	214
Bromophenol blue TAC indicator		1 l	PS0189/15	214
Bromophenol blue TAC indicator		1 l	PS0189/16	214
Bromothymol blue		5 g	428708	215
Bromothymol blue		25 g	428702	215
Bromothymol blue		50 g	428703	215
Calcon		25 g	434171	245
Calconcarbonic acid		5 g	403308	245
Calmagite		5 g	434181	246
Chloramine T sodium salt		25 g	437555	256
Chloranil		50 g	437601	256
Chromotropic acid disodium salt		25 g	404872	268
Clayton's yellow		5 g	453518	270

Product specifications are subject to changes. Please visit our website for updates.

Clayton's yellow	25 g	453519	270
m-Cresol purple	1 g	470067	282
m-Cresol purple	25 g	470068	282
o-Cresol red	5 g	476778	283
Crystal violet	25 g	491502	284
Diacetyldioxime	50 g	441553	295
Diacetyldioxime sodium salt	50 g	441623	295
Diacetyldioxime sodium salt	250 g	441625	295
2,6-Dichlorophenolindophenol sodium salt	5 g	442508	305
2,6-Dichloroquinone-4-chlorimide	5 g	442458	305
Diethylenetriaminepentacetic acid	250 g	405192	308
Dimedone	25 g	444252	313
p-Dimethylaminobenzaldehyde	100 g	444604	314
p-Dimethylaminobenzaldehyde	250 g	444603	314
p-Dimethylaminobenzylidenerhodanine	5 g	444678	314
Dimidium bromide	1 g	445232	320
Dimidium bromide	5 g	445231	320
Dimidium bromide	25 g	445233	320
sym-Diphenylcarbazide	25 g	443752	323
sym-Diphenylcarbazide	100 g	443754	323
sym-Diphenylcarbazone	10 g	443801	323
Diphenylthiocarbazone	50 g	444053	323
Dodecylbenzenesulphonic acid sodium salt	10 g	405351	325
Dodecylbenzenesulphonic acid sodium salt	25 g	405352	325
Eriochrome black T	10 g	464221	330
Eriochrome black T	25 g	464222	330
Eriochromocyanine R	10 g	446811	330
Eriochromocyanine R	25 g	446812	330
Fluorescein	25 g	452086	357
Fluorescein	50 g	452083	357
Fluorescein	500 g	452087	357
Fluorescein sodium salt	25 g	452112	357
Fluorescein sodium salt	50 g	452113	357
Fluorescein sodium salt	1 kg	452117	357
Fluorescein sodium salt	25 g	345356	357
Fluorescein sodium salt	1 kg	345357	357
Glyoxal-bis-(2-hydroxyanil)	10 g	454131	374
Idrimeter Erba Solution A	500 ml	E455256	419
Idrimeter Erba Solution A	1 l	E455257	419
Idrimeter Erba Solution B	500 ml	E455266	419
Idrimeter Erba Solution B	1 l	E455267	419
Idrimeter Erba Indicator C	10 g	E455271	419
Idrimeter Erba Indicator C	100 g	E455274	419
Indicator for ammoniacal nitrogen solution	250 ml	E455651	420
Indicator for iodometry	25 g	455622	420
Indicator for iodometry	250 g	455621	420
Indigo carmine dried	25 g	434932	422
Inulin	10 g	455901	424
Inulin	25 g	455902	424
Inulin	100 g	455903	424

Litmus	100 g	489054	465
Luminol	25 g	458772	467
Metanil yellow	25 g	453542	487
3-Methyl-2-benzothiazolinone hydrazone hydrochloride	5 g	462238	497
Methylene blue	100 g	428984	498
Methylene blue	25 g	429982	498
Methylene blue	500 g	429981	498
Methyl orange	25 g	423504	501
Methyl orange	50 g	423503	501
Methyl orange	250 g	423505	501
Methyl orange	500 g	423501	501
Methyl red	25 g	476882	503
Methyl red	50 g	476883	503
Methyl red	250 g	476881	503
Methylthymol blue sodium salt	1 g	429021	506
Methylthymol blue sodium salt	25 g	429022	506
Methyl yellow	25 g	444552	506
Murexide	5 g	463608	516
Murexide	25 g	463609	516
Neocuproine hydrochloride	1 g	444731	519
Phenol red	5 g	476838	570
Phenol red	25 g	476839	570
Phenolphthalein	100 g	451154	571
Phenolphthalein	500 g	451156	571
Pyrocatechol violet	1 g	491871	632
Pyrocatechol violet	25 g	491872	632
Quinaldine red	25 g	476688	634
Red for oils O	25 g	476961	636
Starch paste solution 1%	250 ml	E477301	721
Starch paste solution 1%	1 l	E477302	721
Sudan III hydroalcoholic saturated solution	250 ml	E485952	727
Sudan yellow	25 g	453582	727
Thymol blue	5 g	429228	757
Thymol blue	25 g	429222	757
Thymol blue	50 g	429223	757
Thymolphthalein	5 g	487728	758
Thymolphthalein	25 g	487729	758
Tropaeolin O	10 g	490001	778
Tropaeolin O	25 g	490002	778
Xylenecyanol	1 g	492211	792
Xylenecyanol	25 g	492212	792
Xylenol orange	1 g	423597	793
Xylenol orange	5 g	423598	793
Xylenol orange	25 g	423599	793

Product specifications are subject to changes. Please visit our website for updates.

CONDUCTIMETRY

The determination of the electrical conductivity or conductance is a key physical-chemical parameter for water analysis.


Standard Solutions

CARLO ERBA Reagents offers the following reference solutions, which are certified with NIST traceability.

Description	Notes	Size	Code	Page
Standard solution 1.30 µS/cm		250 ml	575231	715
Standard solution 5 µS/cm		250 ml	575001	716
Standard solution 10 µS/cm		250 ml	575011	716
Standard solution 20 µS/cm		500 ml	575021	716
Standard solution 50 µS/cm		500 ml	575031	716
Standard solution 84 µS/cm		500 ml	575041	717
Standard solution 100 µS/cm		500 ml	575051	717
Standard solution 147 µS/cm		500 ml	575061	717
Standard solution 200 µS/cm		500 ml	575071	717
Standard solution 500 µS/cm		500 ml	575081	717
Standard solution 1000 µS/cm		500 ml	575091	718
Standard solution 1413 µS/cm		500 ml	575101	718
Standard solution 5000 µS/cm		500 ml	575111	718
Standard solution 10000 µS/cm		500 ml	575121	718
Standard solution 12880 µS/cm		500 ml	575131	718
Standard solution 20000 µS/cm		500 ml	575141	719
Standard solution 50000 µS/cm		500 ml	575151	719
Standard solution 100000 µS/cm		500 ml	575161	719
Standard solution 150000 µS/cm		500 ml	575171	719
Standard solution 200000 µS/cm		500 ml	575181	720
Standard solution 300000 µS/cm		500 ml	575191	720
Standard solution 350000 µS/cm		500 ml	575201	720
Standard solution 450000 µS/cm		500 ml	575211	720
Standard solution 500000 µS/cm		500 ml	575221	721

PARTNER IN YOUR E-CHOICE

Web & e-commerce



easy to use

quick access to products

e-commerce & punch-out

www.carloerbareagents.com

Product specifications are subject to changes. Please visit our website for updates.

PRODUCTS



BEST CHOICE

Labware Consumable & Instruments

Lab Logistics Group



- Glassware
- Plasticware
- Liquid Handling
- Chromatography
- Analitycal Instruments

www.carloerbareagents.com

Absorbent for split liquids

• Assorbente per liquidi versati • Absorbant pour liquides répandus • Absorbente para líquidos derramados • Absorbens für verschüttete Flüssigkeiten

HEU210

Absorbent for split liquids > RE - Pure

RE

Description .....Polvere nocciola Identification ..... Positive

Code	Size	Packaging	Notes
300101	1 kg	Plastic bottle	
300102	5 kg	Plastic tank	

Acetal

• Acetale • Acétal • Acetal • Acetal

Synonym:
• *1,1-Diethoxyacetal*
• *1,1-Diethoxyethane*

$CH_3CH(OC_2H_5)_2$
Molecular Weight: 118,18
CAS: 105-57-7
EEC-N: 203-310-6

Classification transport
ONU: 1088
Transport Hazard class: 3
Packing group II

Danger
H225-H315-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P332+P313

Acetal > RE - Pure - Reagent Ph.Eur.

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.822 ÷ 0.829
Refractive index at 20°C ..... 1.3790 ÷ 1.3850
Assay (GLC) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
400155	100 ml	Glass bottle	

Acetamide

• Acetammide • Acétamide • Acetamida • Acetamid

Synonym:
Amide C2

CH_3CONH_2
Molecular Weight: 59,07
CAS: 60-35-5
EEC-N: 200-473-5

Warning
H351
P201-P202-P280-P308+P313-P405-P501a

Acetamide > RPE - For analysis

RPE

Description ..... Colourless crystals
Identification ..... Positive
Melting point ..... 78.5 ÷ 81.5 °C
Acetate ..... ≤0.2 %
Chloride ..... ≤20 ppm
Alcohol-benzene insol. ..... ≤100 ppm
Heavy metals (Pb) ..... ≤5 ppm
Residue on ignition ..... ≤100 ppm
Sulphate ..... ≤20 ppm
Fe ..... ≤5 ppm
Assay (ex nitrogen) ..... 99 ÷ 100 %

Code	Size	Packaging	Notes
400204	100 g	Plastic bottle	

Acetanilide

• Acetanilide • Acétanilide • Acetanilida • Acetanilid

Synonym:
N-Phenylacetamide

$CH_3CONHC_6H_5$
Molecular Weight: 135,17
CAS: 103-84-4
EEC-N: 203-150-7

Warning
H302
P264-P270-P301+P312a-P330-P501a

Acetanilide > RPE - For analysis

RPE

Description ....white to beige to grey to light brown .....powder or crystals or flakes
Melting point ..... 111 ÷ 115 ° C
Assay (GLC) ..... ≥ 99.0 %

Code	Size	Packaging	Notes
400255	250 g	Plastic bottle	

Acetate buffer pH 6.0

• Tampone acetato pH 6.0 • Tampon acétate pH 6.0 • Tampón acetato pH 6.0 • Puffer acetate pH 6.0

Acetate buffer pH 6.0 > RS - For analysis according to Ph. Eur. Chap. 4.1.3

RS

Code	Size	Packaging	Notes
614002200	1 l	Plastic bottle	Ref Ph.Eur 4002200

Acetate buffer pH 4.6

• Tampone acetato pH 4.6 • Tampon acétate pH 4.6 • Tampón acetato pH 4.6 • Puffer acetate pH 4.6

Acetate buffer pH 4.6 > RS - For analysis according to Ph. Eur. Chap. 4.1.3

RS

Code	Size	Packaging	Notes
614001400	1 l	Plastic bottle	Ref Ph.Eur 4001400

Acetic acid glacial

• Acido acetico glaciale • Acide acétique glacial • Acido acético glacial • Essigsäure

CH_3COOH
Molecular Weight: 60,05
CAS: 64-19-7
EEC-N: 200-580-7

Classification transport
ONU: 2789
Transport Hazard class: 8
Packing group II



Danger
H226-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Acetic acid glacial > RS - For LC/MS

RS

Description	Clear colourless liquid	Subst. reducing KMnO4	Conform	at 260 nm	≥ 80 %	Ca	≤ 0.1 ppm
Colour (APHA)	≤10	Subst.reducing dichromate	Conform	at 270 nm	≥ 95 %	Mg	≤ 0.1 ppm
Refractive index at 20°C	1.3711 ÷ 1.3731	Residue on evaporation	≤ 5 ppm	at 280 nm	≥ 97 %	K	≤ 0.1 ppm
Density at 20° C	1.0501 ÷ 1.0521	Assay (GC)	≥ 99.95 %	Al	≤ 0.05 ppm		
Boiling point	118.3 - 118.8. °C	U.V. Transmittance		Fe	≤ 0.2 ppm		
Water (K.F.)	≤ 0.1 %	at 254 nm	≥ 30 %	Na	≤ 0.5 ppm		

Code	Size	Packaging	Notes
401411	10 x 1 ml	Glass ampoule	
401412	10 x 2.5 ml	Bottle	
401413	50 ml	Plastic bottle	
401414	1 l	Glass bottle	

Additive for eluent phase for LC-MS. Store at temperature > 20 °C

Acetic acid glacial > RS - For HPLC - Isocratic Grade

RS

Appearance	Clear colourless liquid	Water content (K.F.)	≤ 0.05 % m/m	UV transmittance at 260 nm	≥ 80 %	Assay (GC)	≥ 99.8 %
Colour	≤ 10 Hazen	UV transmittance at 254 nm	≥ 25 %	Non volatile residue	≤ 10 mg/Kg		

Code	Size	Packaging	Notes
401431	1 l	Glass bottle	
401432	2.5 l	Glass bottle	

Store at temperature > 20 °C

Product specifications are subject to changes. Please visit our website for updates.

Acetic acid glacial > RS - Ultrapure - For trace analysis at ppt level

RS

Description	Clear colourless liquid	Mn	≤ 10 ppt	Cs	≤ 10 ppt	Re	≤ 10 ppt
Identification	Positive	Mo	≤ 10 ppt	Dy	≤ 1 ppt	Rh	≤ 50 ppt
Ag	≤ 50 ppt	Na	≤ 100 ppt	Er	≤ 1 ppt	Rb	≤ 10 ppt
Al	≤ 50 ppt	Ni	≤ 50 ppt	Eu	≤ 1 ppt	Ru	≤ 50 ppt
As	≤ 50 ppt	Pb	≤ 10 ppt	Gd	≤ 1 ppt	Sm	≤ 1 ppt
Ba	≤ 10 ppt	Sb	≤ 50 ppt	Ga	≤ 10 ppt	Sc	≤ 10 ppt
Be	≤ 10 ppt	Sn	≤ 50 ppt	Ge	≤ 10 ppt	Te	≤ 1 ppt
Bi	≤ 10 ppt	Sr	≤ 10 ppt	Hf	≤ 10 ppt	Tb	≤ 1 ppt
Ca	≤ 50 ppt	Ti	≤ 10 ppt	Ho	≤ 1 ppt	Tm	≤ 1 ppt
Cd	≤ 10 ppt	V	≤ 10 ppt	In	≤ 1 ppt	W	≤ 10 ppt
Co	≤ 10 ppt	Zn	≤ 50 ppt	La	≤ 1 ppt	Yb	≤ 1 ppt
Cr	≤ 10 ppt	Zr	≤ 10 ppt	Li	≤ 10 ppt	Y	≤ 1 ppt
Cu	≤ 10 ppt	Assay (acidimetric)	≥ 99 %	Lu	≤ 10 ppt	Tl	≤ 10 ppt
Fe	≤ 50 ppt	U	≤ 1 ppt	Nd	≤ 1 ppt		
K	≤ 50 ppt	Th	≤ 1 ppt	Pt	≤ 50 ppt		
Mg	≤ 50 ppt	Ce	≤ 10 ppt	Pr	≤ 1 ppt		

Code	Size	Packaging	Notes
401361	500 ml	Plastic bottle	

Store at temperature > 20 °C

Acetic acid glacial > RS - Superpure - For trace analysis at ppb level

RS

Description	Clear colourless liquid	Cr	≤ 1 ppb	Mg	≤ 0.5 ppb	Sr	≤ 0.5 ppb
Colour (APHA)	≤ 10	Cs	≤ 0.1 ppb	Mn	≤ 0.5 ppb	Tb	≤ 0.1 ppb
Identification	Positive	Cu	≤ 0.5 ppb	Mo	≤ 0.5 ppb	Te	≤ 0.5 ppb
Chloride	≤ 1 ppm	Dy	≤ 0.1 ppb	Na	≤ 1 ppb	Th	≤ 0.1 ppb
Phosphate	≤ 1 ppm	Er	≤ 0.1 ppb	Nd	≤ 0.1 ppb	Ti	≤ 0.5 ppb
Sulphate	≤ 0.5 ppm	Eu	≤ 0.1 ppb	Ni	≤ 0.5 ppb	Tl	≤ 0.1 ppb
Reducing dichromate	Conform	Fe	≤ 1 ppb	Pb	≤ 0.1 ppb	Tm	≤ 0.1 ppb
Subst. reducing KMn04	Conform	Ga	≤ 0.1 ppb	Pr	≤ 0.1 ppb	U	≤ 0.1 ppb
Al	≤ 1 ppb	Ge	≤ 0.5 ppb	Pt	≤ 0.5 ppb	V	≤ 0.5 ppb
Ag	≤ 1 ppb	Gd	≤ 0.1 ppb	Rb	≤ 0.1 ppb	W	≤ 0.5 ppb
As	≤ 0.5 ppb	Hf	≤ 0.1 ppb	Re	≤ 0.1 ppb	Y	≤ 0.1 ppb
Ba	≤ 0.5 ppb	Hg	≤ 1 ppb	Rh	≤ 0.5 ppb	Yb	≤ 0.1 ppb
Be	≤ 0.1 ppb	Ho	≤ 0.1 ppb	Ru	≤ 0.5 ppb	Zn	≤ 1 ppb
Bi	≤ 0.1 ppb	In	≤ 0.1 ppb	Sb	≤ 0.5 ppb	Zr	≤ 0.1 ppb
Ca	≤ 1 ppb	K	≤ 1 ppb	Sc	≤ 0.1 ppb	Assay (acidimetric)	≥ 99 %
Cd	≤ 0.5 ppb	La	≤ 0.1 ppb	Se	≤ 1 ppb		
Ce	≤ 0.1 ppb	Li	≤ 0.1 ppb	Sm	≤ 0.1 ppb		
Co	≤ 0.1 ppb	Lu	≤ 0.1 ppb	Sn	≤ 0.5 ppb		

Code	Size	Packaging	Notes
401405	500 ml	Plastic bottle	
401406	1 l	Plastic bottle	
401407	2.5 l	Plastic bottle	

Store at temperature > 20 °C

Acetic acid glacial > RS - RSE - For electronic use

RS

Description	Clear liquid	Phosphate	≤0.1 ppm	Ca	≤0.2 ppm	Na	≤0.2 ppm
Colour (APHA)	≤10	Heavy metals (Pb)	≤0.2 ppm	Cd	≤0.01 ppm	Ni	≤0.03 ppm
Identification	Positive	Reducing chromic acid	≤100 ppm	Co	≤0.01 ppm	Pb	≤0.01 ppm
Water miscibility	Conform	Subst. reducing KMn04	≤10 ppm	Cr	≤0.03 ppm	Pt	≤0.05 ppm
Freezing point	≥16.24 °C	Sulphate	≤0.5 ppm	Cu	≤0.01 ppm	Sb	≤0.005 ppm
Assay (acidimetric)	≥99.9 %	Ag	≤0.02 ppm	Fe	≤0.1 ppm	Sn	≤0.02 ppm
Density at 20° C	1.0495 ÷ 1.0503	Al	≤0.01 ppm	Ga	≤0.02 ppm	Sr	≤0.02 ppm
Boiling point	118.3 ÷ 118.8 °C	As	≤0.005 ppm	In	≤0.02 ppm	Ta	≤0.1 ppm
Residue on evaporation	≤5 ppm	Au	≤0.05 ppm	K	≤0.1 ppm	Ti	≤0.05 ppm
Formic acid	≤0.1 %	B	≤0.01 ppm	Li	≤0.02 ppm	Tl	≤0.05 ppm
Acetic anhydride	≤100 ppm	Ba	≤0.1 ppm	Mg	≤0.05 ppm	V	≤0.05 ppm
Chloride	≤1 ppm	Be	≤0.02 ppm	Mn	≤0.01 ppm	Zn	≤0.05 ppm
Carbonyl Compounds (CO)	≤2 ppm	Bi	≤0.02 ppm	Mo	≤0.05 ppm	Zr	≤0.05 ppm

Code	Size	Packaging	Notes
401463	1 l	Glass bottle	
401462	2.5 l	Glass bottle	

Store at temperature > 20 °C

Acetic acid glacial > RS - For potentiometry

RS

Refractive index at 20°C ..... 1.371 - 1.374	Reducing substances ..... Conform	Aluminium (Al) ..... ≤ 0.05 mg/Kg	Magnesium (Mg) ..... ≤ 0.1 mg/Kg
Water content (K.F.) ..... ≤ 1000 mg/Kg	Assay (GC) ..... ≥ 99.8 %	Barium (Ba) ..... ≤ 0.1 mg/Kg	Manganese (Mn) ..... ≤ 0.05 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg	Acetic anhydride ..... ≤ 0.025 %	Cadmium (Cd) ..... ≤ 0.05 mg/Kg	Lead (Pb) ..... ≤ 0.05 mg/Kg
Colour ..... ≤ 10 Hazen	Chloride (Cl-) ..... ≤ 1 mg/Kg	Cobalt (Co) ..... ≤ 0.05 mg/Kg	Zinc (Zn) ..... ≤ 0.1 mg/Kg
Titrable base ..... Conform	Sulphate (SO4--) ..... ≤ 1 mg/Kg	Iron (Fe) ..... ≤ 1 mg/Kg	Arsenic (As) ..... ≤ 0.02 mg/Kg

Code	Size	Packaging	Notes
P00725P15	1 l	Plastic bottle	
P00725P21	2.5 l	Glass bottle	

Store at temperature > 20°C

Acetic acid glacial > RS - For titration in non-aqueous medium

RS

Description ..... Clear colourless liquid	Identification ..... Positive	Density at 20° C ..... 1.049 ÷ 1.051	Assay ..... ≥ 99.8 %

Code	Size	Packaging	Notes
401453	1 l	Glass bottle	
401455	2.5 l	Glass bottle	

Store at temperature > 20 °C

Acetic acid glacial > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description . Clear liquid or crystalline mass	Formic acid ..... ≤0.05 %	Bi ..... ≤0.02 ppm	Na ..... ≤0.5 ppm
Colour (APHA) ..... ≤10	Acetic anhydride ..... ≤100 ppm	Ca ..... ≤0.2 ppm	Ni ..... ≤0.03 ppm
Identification ..... Positive	Chloride ..... ≤1 ppm	Cd ..... ≤0.05 ppm	Pb ..... ≤0.02 ppm
Water miscibility ..... Conform	Carbonyl Compounds (CO) ..... ≤50 ppm	Co ..... ≤0.01 ppm	Sr ..... ≤0.02 ppm
Titrable base ..... Conform	Phosphate ..... ≤0.5 ppm	Cr ..... ≤0.03 ppm	Ti ..... ≤0.1 ppm
Subst. reducing KMn04 ..... Conform	Heavy metals (Pb) ..... ≤0.5 ppm	Cu ..... ≤0.01 ppm	Tl ..... ≤0.05 ppm
Density at 20° C ..... 1.0501 ÷ 1.0521	Sulphate ..... ≤0.5 ppm	Fe ..... ≤0.2 ppm	V ..... ≤0.05 ppm
Refractive index at 20°C . 1.3711 ÷ 1.3731	Ag ..... ≤0.02 ppm	K ..... ≤0.1 ppm	Zn ..... ≤0.05 ppm
Boiling point ..... 118.3 ÷ 118.8 °C	Al ..... ≤0.05 ppm	Li ..... ≤0.02 ppm	Zr ..... ≤0.1 ppm
Freezing point ..... ≥16 °C	As ..... ≤0.01 ppm	Mg ..... ≤0.1 ppm	Reducing chromic acid ..... Conform ACS
Assay (acidimetric) ..... 99.5 ÷ 100.5 %	Ba ..... ≤0.1 ppm	Mn ..... ≤0.01 ppm	Water (K.F.) ..... ≤ 1500 ppm
Residue on evaporation ..... ≤ 10 ppm	Be ..... ≤0.02 ppm	Mo ..... ≤0.05 ppm	Acetaldehyde ..... ≤ 500 ppm

Code	Size	Packaging	Notes
401421	1 l	Glass bottle PVC coated	
401422	1 l	Glass bottle	
401424	2.5 l	Glass bottle	
401425	30 kg	Plastic drum	

Store at temperature > 20 °C

Acetic acid glacial > RPE - For analysis

RPE

Description . Clear liquid or crystalline mass	Acetic anhydride ..... ≤100 ppm	Ca ..... ≤0.2 ppm	Ni ..... ≤0.03 ppm
Colour (APHA) ..... ≤10	Chloride ..... ≤1 ppm	Cd ..... ≤0.05 ppm	Pb ..... ≤0.02 ppm
Identification ..... Positive	Carbonyl Compounds (CO) ..... ≤50 ppm	Co ..... ≤0.01 ppm	Sr ..... ≤0.02 ppm
Water miscibility ..... Conform	Phosphate ..... ≤0.5 ppm	Cr ..... ≤0.03 ppm	Ti ..... ≤0.1 ppm
Subst. reducing KMn04 ..... Conform	Heavy metals (Pb) ..... ≤0.5 ppm	Cu ..... ≤0.01 ppm	Tl ..... ≤0.05 ppm
Density at 20° C ..... 1.0501 ÷ 1.0521	Sulphate ..... ≤0.5 ppm	Fe ..... ≤0.5 ppm	V ..... ≤0.05 ppm
Refractive index at 20°C . 1.3711 ÷ 1.3731	Ag ..... ≤0.02 ppm	K ..... ≤0.1 ppm	Zn ..... ≤0.05 ppm
Boiling point ..... 118.3 ÷ 118.8 °C	Al ..... ≤0.05 ppm	Li ..... ≤0.02 ppm	Zr ..... ≤0.1 ppm
Freezing point ..... ≥16 °C	As ..... ≤0.01 ppm	Mg ..... ≤0.1 ppm	Water (K.F.) ..... ≤ 1500 ppm
Assay (acidimetric) ..... 99.5÷100.5 %	Ba ..... ≤0.1 ppm	Mn ..... ≤0.01 ppm	Assay (CPG) ..... ≥ 99.8 %
Residue on evaporation ..... ≤10 ppm	Be ..... ≤0.02 ppm	Mo ..... ≤0.05 ppm	Acetaldehyde ..... ≤ 500 ppm
Formic acid ..... ≤0.05 %	Bi ..... ≤0.02 ppm	Na ..... ≤0.5 ppm	

Code	Size	Packaging	Notes
401391	1 l	Glass bottle	
524520	1 l	Plastic bottle	
401392	2.5 l	Glass bottle	
524521	2.5 l	Plastic bottle	
401396	30 kg	Plastic drum	
401397	200 kg	Polythene-metal drum	

Store at temperature > 20°C

Product specifications are subject to changes. Please visit our website for updates.

Acetic acid glacial > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-NF-DAB-JP

ERBApharm

Description .....Clear, colourless liquid or crystalline mass
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Colour ..... ≤ 10 APHA
Reducing impurities.....Conform Ph.Eur.
Ready oxidizable substances.....Conform USP-NF
Freezing point ..... ≥ 15.6 °C
Refractive index at 20°C..... 1.370 - 1.374
Water (K.F.) .....≤ 3000 ppm
Residue on evaporation .....≤ 50 mg/l
Chloride.....≤ 2 mg/l
Sulphate .....≤ 2 mg/l
Sulfate..... Conform USP-NF
Heavy metals (Pb).....≤ 5 ppm
Fe .....≤ 5 ppm
Assay (acidimetric) .....99.5 ÷ 100.5 %
Origin (BSE/TSE)..... Synthesis

Code	Size	Packaging	Notes
302016	1 l	Glass bottle	
302011	2.5 l	Glass bottle	
302014	5 l	Plastic tank	
302015	30 kg	Plastic drum	
302013	200 kg	Polythene-metal drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade. Store at tempera ture > 20 °C

Acetic acid glacial > RE - Pure

RE

Description .....Clear liquid
Density at 20°C .....1.049 ÷ 1.051
Colour (APHA) ..... ≤ 10
Identification ..... Positive
Water miscibility.....Conform
Assay ..... ≥ 99 %
Subst. reducing KMnO4.....≤ 40 ppm
Density at 20° C .....1.050 ÷ 1.052
Refractive index at 20°C.....1.371 ÷ 1.373
Boiling point.....118.3 ÷ 118.8 °C
Freezing point ..... ≥ 16.24 °C
Residue on evaporation .....≤ 30 ppm
Water .....≤ 2500 ppm
Formic acid.....≤ 100 ppm
Chloride.....≤ 1 ppm
Carbonyl Compounds (CO) .....≤ 2 ppm
Heavy metals (Pb).....≤ 1 ppm
Reducing chromic acid .....≤ 30 ppm
Sulphate .....≤ 2 ppm
As .....≤ 0.5 ppm
Cu.....≤ 0.1 ppm
Fe .....≤ 1 ppm
Ni.....≤ 0.1 ppm
Assay (acidimetric) ..... ≥ 99.8 %

Code	Size	Packaging	Notes
302031	1 l	Plastic bottle	
302032	2.5 l	Glass bottle	
302034	5 l	Plastic tank	
302033	10 l	Plastic tank	

Store at temperature > 20 °C

Acetic acid glacial > RE - Pure - For fibers

RE

Description .....Clear liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Water miscibility.....Conform
Density at 20° C .....1.050 ÷ 1.052
Refractive index at 20°C.1.3711 ÷ 1.3731
Boiling point.....118.3 ÷ 118.8 °C
Subst. reducing KMnO4.....≤40 ppm
Residue on evaporation .....≤8 ppm
Formic acid..... ≤0.1 %
Acetic anhydride .....≤500 ppm
Chloride.....≤1 ppm
Carbonyl Compounds (CO).....≤50 ppm
Heavy metals (Pb).....≤1 ppm
Sulphate .....≤1 ppm
Al .....≤0.1 ppm
As .....≤0.5 ppm
Cu.....≤0.1 ppm
Fe .....≤0.5 ppm
Ni.....≤0.1 ppm
Assay ..... ≥ 99.8 %
Water (K.F.) .....≤ 2500 ppm
Assay (CPG) ..... ≥ 99.5 %
Freezing point .....16.2 ÷ 16.6 °C

Code	Size	Packaging	Notes
302021	200 kg	Polythene-metal drum	

Store at temperature > 20 °C

Acetic acid 96%

• Acido acetico 96% • Acide acétique 96% • Acido acético 96% • Essigsäure 96%

CH_3COOH
Molecular Weight: 60,05
CAS: 64-19-7
EEC-N: 200-580-7

Classification transport
ONU: 2789
Transport Hazard class: 8
Packing group II



Danger
H226-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Acetic acid 96% > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C .....1.049 ÷ 1.051
Residue on evaporation .....≤100 ppm
Heavy metals (Pb).....≤50 ppm
Sulphate .....≤1000 ppm
Iron .....≤50 ppm
Assay (acidimetric) ..... ≥96 %

Code	Size	Packaging	Notes
302002	1 l	Glass bottle	
302003	2.5 l	Glass bottle	
302005	25 kg	Plastic tank	
302007	50 kg	Plastic tank	

Store at temperature > 20 °C

Acetic acid 80%

• Acido acetico 80% • Acide acétique 80% • Acido acético 80% • Essigsäure 80%

CH_3COOH
Molecular Weight: 60,05
CAS: 64-19-7

Classification transport
ONU: 2789
Transport Hazard class: 8
Packing group II



Danger
H226-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Acetic acid 80% > RE - Pure

RE

Description ..... Clear liquid	Density at 20° C ..... 1.069 ÷ 1.071	Acetic anhydride ..... ≤500 ppm	As ..... ≤0.4 ppm
Acetic acid content ..... 80.5 - 82.5 %	Subst. reducing KMn04 ..... ≤40 ppm	Chloride ..... ≤1 ppm	Cu ..... ≤0.1 ppm
Colour (APHA) ..... ≤10	Residue on evaporation ..... ≤30 ppm	Heavy metals (Pb) ..... ≤2 ppm	Fe ..... ≤0.5 ppm
Identification ..... Positive	Formic acid ..... ≤200 ppm	Sulphate ..... ≤1 ppm	Ni ..... ≤0.1 ppm
Water miscibility ..... Conform	Alcoh acetone acetaldeh ..... ≤50 ppm	Al ..... ≤0.1 ppm	Assay (acidimetric) ..... 80.5 ÷ 82.5 %

Code	Size	Packaging	Notes
301855	5 l	Plastic tank	
301852	50 kg	Plastic tank	
301853	200 kg	Plastic drum	

Store at temperature > 20°C

Acetic acid 45%

• Acido acetico 45% • Acide acétique 45% • Acido acético 45 % • Essigsäure 45%

CH_3COOH
Molecular Weight: 60,05
CAS: 64-19-7

Classification transport
ONU: 2790
Transport Hazard class: 8
Packing group III



Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Acetic acid 45% > RE - Pure - For glassware washing

RE

Description ..... Clear colourless liquid Assay (acetic acid) ..... 44 ÷ 46 %

Code	Size	Packaging	Notes
526545	5 l	Plastic tank	
526546	10 l	Plastic tank	

Acetic acid 30%

• Acido acetico 30% • Acide acétique 30% • Acido acético 30% • Essigsäure 30%

CH_3COOH
Molecular Weight: 60,05
CAS: 64-19-7

Classification transport
ONU: 2790
Transport Hazard class: 8
Packing group III



Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Acetic acid 30% > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Assay ..... 290 - 310 g/L

Code	Size	Packaging	Notes
611000401	1 l	Glass bottle	Ref Ph.Eur 1000401

Product specifications are subject to changes. Please visit our website for updates.

Acetic acid 27%

• Acido acetico 27% • Acide acétique 27% • Acido acético 27 % • Essigsäure 27%

CH_3COOH
Molecular Weight: 60,05
CAS: 64-19-7

Classification transport
ONU: 2790
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Acetic acid 27% > RS - For glassware washing

RS

Assay (acidimetric) ..... 26 ÷ 28 %

Code	Size	Packaging	Notes
508645	5 l	Plastic tank	

Content is guaranteed for standardized volumes at 20 °C. Keep tightly closed

Acetic acid 25%

• Acido acetico 25% • Acide acétique 25% • Acido acético 25 % • Essigsäure 25%

CH_3COOH
Molecular Weight: 60,05
CAS: 64-19-7

Classification transport
ONU: 2790
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Acetic acid 25% > RE - Pure

RE

Density d20/4 ..... 1.03 - 1.032 Acetic acid content ..... 24 - 26 %

Code	Size	Packaging	Notes
PS0222/52	30 l	Plastic tank	

Acetic acid 20%

• Acido acetico 20% • Acide acétique 20% • Acido acético 20 % • Essigsäure 20%

CH_3COOH
Molecular Weight: 60,05
CAS: 64-19-7

Classification transport
ONU: 2790
Transport Hazard class: 8
Packing group III

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Acetic acid 20% > RPE - For analysis

RPE

Density d20/4 ..... 1.024 - 1.026 Acetic acid content ..... 19 - 21 %

Code	Size	Packaging	Notes
PS0237/41	10 l	Plastic tank	

Acetic acid 12%

• Acido acetico 12% • Acide acétique 12% • Acido acético 12 % • Essigsäure 12%

CH_3COOH
Molecular Weight: 60,05
CAS: 64-19-7

Classification transport
ONU: 2790
Transport Hazard class: 8
Packing group III

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Acetic acid 12% > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611000402	1 l	Plastic bottle	Ref Ph.Eur 1000402

Acetic acid 12% > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Density at 20° C ..... 1.014 ÷ 1.016 Assay (acidimetric) ..... 11.5 - 12.5 %

Code	Size	Packaging	Notes
401531	1 l	Plastic bottle	

Content is guaranteed for standardized volumes at 20 °C. Keep tightly closed

Acetic acid 12% > RE - Pure

RE

Density d20/4 ..... 1.014 - 1.015 Acetic acid content ..... 11.5 - 12.5 %

Code	Size	Packaging	Notes
PS0221/29	5 l	Plastic tank	

Acetic acid 1 mol/l (1N)

• Acido acetico 1 mol/l (1N) • Acide acétique 1 mol/l (1N) • Acido acético 1 mol/l (1N) • Essigsäure 1 mol/l (1N)

CH_3COOH
Molecular Weight: 60,05
CAS: 64-19-7

Acetic acid 1 mol/l (1N) > RPE - For analysis - Reag. Ph.Eur. - Reag. USP

RPE

Assay (potentiometric) ..... 0.998 ÷ 1.002 mol/L

Code	Size	Packaging	Notes
524605	5 l	Plastic tank	

Content is guaranteed for standardized volumes at 20°C. Keep tightly closed

Acetic acid 1 mol/l (1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.999 - 1.001 N

Code	Size	Packaging	Notes
502000	1 l	Plastic bottle	

Content is guaranteed for standardized volumes at 20 °C. Keep tightly closed

Acetic acid 0.1 mol/l (0.1N)

• Acido acetico 0.1 mol/l (0.1N) • Acide acétique 0.1 mol/l (0.1N) • Acido acético 0.1 mol/l (0.1N) • Essigsäure 0.1 mol/l (0.1N)

CH_3COOH
Molecular Weight: 60,05
CAS: 64-19-7

Acetic acid 0.1 mol/l (0.1N) > RPE - For analysis

RPE

Assay (potentiomety) ..... 0.0999 - 0.1001 N

Code	Size	Packaging	Notes
P3100015	1 l	Plastic bottle	

Acetic acid 0.1 mol/l (0.1N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
401561		Plastic ampoule	Volume: 55 ml

6,005 g of CH3COOH. Volumetric concentrated solution to prepare 1 L of solution 0,1 N

Acetic acid 0.03 mol/l (0.03N)

• Acido acetico 0.03 mol/l (0.03N) • Acide acétique 0.03 mol/l (0.03N) • Acido acético 0.03 mol/l (0.03N) • Essigsäure 0.03 mol/l (0.03N)

CH_3COOH
Molecular Weight: 60,05
CAS: 64-19-7

Acetic acid 0.03 mol/l (0.03N) > RPE - For analysis - Reag. Ph.Eur. - Reag. USP

RPE

Assay (potentiometric) ..... 0.02994 ÷ 0.03006 mol/L

Code	Size	Packaging	Notes
524611	10 l	Kubidos	

Product specifications are subject to changes. Please visit our website for updates.

Acetic acid-d4

• Acido acetico-d4 • Acide acétique-d4 • Acido acético-d4 • Essigsäure d-4

Synonym:
Tetradeuteroacetic acid

CD_3COOD
Molecular Weight: 64,08
CAS: 1186-52-3
EEC-N: 214-693-4

Classification transport
ONU: 2789
Transport Hazard class: 8
Packing group II

Danger
H226-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Acetic acid-d4 > RS - For NMR - min 99,9%

RS

Code	Size	Packaging	Notes
P5039	10 x 0.75 ml	Glass ampoule	

For specifications, contact our customer service for a certificate of analysis

Acetic acid isopropyl ester ▶ Isopropyl acetate

Acetic acid magnesium salt ▶ Magnesium acetate tetrahydrate

Acetic anhydride

• Anidride acetica • Anhydride acétique • Anhídrido acético • Essigsäureanhydrid

$(CH_3CO)_2O$
Molecular Weight: 102,09
CAS: 108-24-7

Classification transport
ONU: 2733
Transport Hazard class: 3
Packing group II

Danger
H225-H302-H312-H332-H314-H335
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P362+P364-P403+P233

Acetic anhydride > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611000501	1 l	Glass bottle	Acetic anhydride solution R1 Ref Ph.Eur 1000501

Storage: protected from light and air

Acetic anhydride > RPE - For analysis - ACS

RPE

Description	Clear liquid	Refractive index at 20°C	1.3881 ÷ 1.3931	Chloride	≤5 ppm	Fe	≤5 ppm
Colour (APHA)	≤20	Boiling point	136 ÷ 142 °C	Phosphate	≤10 ppm	Assay (GLC)	≥97.0 %
Identification	Positive	Residue on evaporation	≤30 ppm	Heavy metals (Pb)	≤2 ppm		
Density at 20° C	1.080 ÷ 1.084	Subst. reducing KMnO4	Conform ACS	Sulphate	≤5 ppm		

Code	Size	Packaging	Notes
421491	1 l	Glass bottle	
421496	2.5 l	Glass bottle	
421493	30 kg	Plastic drum	

Storage: protected from light and air

Acetic anhydride > RE - Pure

RE

Description	Clear liquid	Refractive index at 20°C	1.3856 ÷ 1.3956	Chloride	≤5 ppm	Assay (GLC)	≥99.5 %
Identification	Positive	Residue on evaporation	≤100 ppm	Heavy metals (Pb)	≤ 2 ppm	Assay (acidimetric)	≥99.5 % (m/m)
Colour	≤ 10 APHA	Subs. reducing KMnO4	≤ 0.02 %	Al	≤ 1 ppm		
Density at 20° C	1.079 ÷ 1.085	Acetic acid	≤ 0.5 %	Fe	≤ 1 ppm		

Code	Size	Packaging	Notes
316501	1 l	Glass bottle	
316503	30 kg	Plastic drum	
316502	210 kg	Metal drum	

Storage: protected from light and air

Acetoacetic ester ▶ Ethyl acetoacetate

Acetone

• Acetone • Acétone • Acetona • Aceton

Synonym:
2-Propanone

CH_3COCH_3
Molecular Weight: 58,01
CAS: 67-64-1
EEC-N: 200-662-2

Classification transport
ONU: 1090
Transport Hazard class: 3
Packing group II

Danger
H225-H319-H336-HEU066
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Acetone > RS - For HPLC - Isocratic Grade

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.790 ÷ 0.792
Refractive index at 20°C. 1.3581 ÷ 1.3601
Boiling point..... 55.8 ÷ 56.3 ° C
Water (K.F.) ..... ≤ 500 ppm
Residue on evaporation ..... ≤5 ppm
Acidity ..... ≤0.0005 meq/g
Alcalinity ..... ≤0.0002 meq/g
Assay (GLC) ..... ≥99.9 %
U.V. Transmittance
At 335 nm ..... ≥ 60 %
at 340 nm ..... ≥ 85 %
At 345 nm ..... ≥ 90 %
at 350 nm ..... ≥ 98 %
at 360 nm ..... ≥ 99 %

Code	Size	Packaging	Notes
412501	1 l	Glass bottle	
412502	2.5 l	Glass bottle	

Acetone > RS - For GC-MS

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C..... 1.357 - 1.361
Water (K.F.) ..... ≤ 500 ppm
Colour ..... ≤ 10 APHA
Residue on evaporation ..... ≤ 2 ppm
Acidity (acetic acid)..... ≤ 20 ppm
Assay (GC) ..... ≥ 99.95 %
Ethyl alcohol ..... ≤ 100 ppm
Methyl alcohol..... ≤ 500 ppm
Isopropyl alcohol ..... ≤ 500 ppm
GC-MS.Individual peak (n-hexadecane). ≤ 2 µg/L
Ret.range n-undecane to n-tetracontane

Code	Size	Packaging	Notes
400952	1 l	Glass bottle	

Acetone > RS - ATRASOL - For traces analysis

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C..... 1.357 - 1.361
Water content (K.F.) ..... ≤ 500 mg/Kg
Non volatile residue..... ≤ 2 mg/Kg
Colour ..... ≤ 10 Hazen
Free acid (as CH3COOH)..... ≤ 20 mg/Kg
Assay (GC) ..... ≥ 99.9 %
GC (FID) - NC Atrasol ..... Conform
Ethanol ..... ≤ 100 mg/Kg
Methanol ..... ≤ 500 mg/Kg
2-Propanol..... ≤ 500 mg/Kg
GC-ECD.Individual peak (CCl4)..... ≤ 1 µg/l
Ret.range dichloromethane to 1,2,4-trichlorobenzene
GC-ECD.Individual peak (Lindane) .≤ 2 ng/L
Ret.range 1,2,4-trichlorobenzene to decachlorobiphenyle
GC-FID.Individual peak (n-hexadecane). ≤ 2 µg/L
Ret.range n-undecane to n-tetracontane

Code	Size	Packaging	Notes
P0053216	1 l	Glass bottle	
P0053221	2.5 l	Glass bottle	
P0053282	4 l	Glass bottle	

Acetone > RS - PESTIPUR - For pesticide analysis

RS

Description ..... Clear liquid
Identification ..... Positive
Colour ..... ≤ 10 hazen
Assay (GLC) ..... ≥ 99.8 %
Water ..... ≤ 0.05 %
Not volatile residue..... ≤ 2 mg/kg
Free acid (as CH3COOH)..... ≤ 20 mg/kg
GC-ECD (Lindane standard) ..... ≤ 3 ng/l
GC-NPD (Ethylparathion standard) .. ≤ 3 ng/l
Refractive index at 20°C..... 1.357 ÷ 1.361

Code	Size	Packaging	Notes
400991	1 l	Glass bottle	
400992000	2.5 l	Glass bottle	
400994	4 l	Glass bottle	

Acetone > RS - PESTIPUR - For pesticide analysis - PAH guaranteed

RS

Description ..... Clear liquid
Identification ..... Positive
Colour ..... ≤ 10 hazen
Refractive index at 20°C..... 1.357 ÷ 1.361
Assay (GLC) ..... ≥ 99.8 %
Water ..... ≤ 0.05 %
Not volatile residue..... ≤ 2 mg/kg
Free acid (as CH3COOH)..... ≤ 20 mg/kg
GC-ECD (Lindane standard) ..... ≤ 3 ng/l
GC-NPD (Ethylparathion standard) .. ≤ 3 ng/l
PAH test (according to ISO 17993) .. Passes test
Test against standard (each 16 PAH) . ≤ 0.1 µg/l

Code	Size	Packaging	Notes
400932	2.5 l	Glass bottle	

16 selected PAHs tested according to ISO 17993:2002, each max 0,1µg/l

Product specifications are subject to changes. Please visit our website for updates.

Acetone > RS - SPECTROSOL - For optical spectroscopy

RS

Description ..... Clear liquid	Boiling point ..... 55.8 ÷ 56.3 ° C	Assay (GLC) ..... ≥99.9 %	at 335 nm ..... ≥60 %
Colour (APHA) ..... ≤10	Water (K.F.) ..... ≤ 500 ppm	Fluorescence	at 340 nm ..... ≥85 %
Identification ..... Positive	Residue on evaporation ..... ≤5 ppm	at 365 nm ..... ≤2 ppb	at 345 nm ..... ≥95 %
Density at 20° C ..... 0.790 ÷ 0.792	Acidity ..... ≤0.0005 meq/g	U.V. Transmittance	at 350 nm ..... ≥98 %
Refractive index at 20°C . 1.3581 ÷ 1.3601	Alcalinity ..... ≤0.0002 meq/g	at 330 nm ..... ≥16 %	

Code	Size	Packaging	Notes
401034	1 l	Glass bottle	
401032	2.5 l	Glass bottle	

Acetone > RS - Anhydrous - For analysis

RS

Clear,colourless liq.appearance ..... Conform	Non volatile residue ..... ≤ 10 mg/Kg	Diacetyl alcohol ..... ≤ 500 mg/Kg	2-Propanol ..... ≤ 500 mg/Kg
Refractive index at 20°C ..... 1.357 - 1.361	Assay (GC) ..... ≥ 99.8 %	Benzene ..... ≤ 2 mg/Kg	
Density d20/20 ..... 0.790 - 0.793	Colour ..... ≤ 10 Hazen	Methanol ..... ≤ 500 mg/Kg	
Water content (K.F.) ..... ≤ 100 mg/Kg	Free acid (as CH3COOH) ..... ≤ 20 mg/Kg	Ethanol ..... ≤ 100 mg/Kg	

Code	Size	Packaging	Notes
P0051010	200 ml	Bottle with septum	
P0051016	1 l	Glass bottle	

Acetone > RS - VLSI - For electronic use

RS

Code	Size	Packaging	Notes
527651	1 l	Plastic bottle	
527650	2.5 l	Plastic bottle	
527655	5 l	Plastic bottle	

Particles control < 250 particles 0.5 µm/ml
For specifications, contact our customer service for a certificate of analysis

Acetone > RS - RSE - For electronic use

RS

Description ..... Clear liquid	Methyl alcohol ..... ≤500 ppm	Bi ..... ≤0.02 ppm	Na ..... ≤0.2 ppm
Colour (APHA) ..... ≤10	Aldehyde ..... ≤10 ppm	Ca ..... ≤0.1 ppm	Ni ..... ≤0.01 ppm
Identification ..... Positive	Chloride ..... ≤0.1 ppm	Cd ..... ≤0.005 ppm	Pb ..... ≤0.01 ppm
Water miscibility ..... Conform	Phosphate ..... ≤0.1 ppm	Co ..... ≤0.005 ppm	Pt ..... ≤0.02 ppm
Assay (GLC) ..... ≥99.8 %	Heavy metals (Pb) ..... ≤0.2 ppm	Cr ..... ≤0.01 ppm	Sb ..... ≤0.01 ppm
Resistivity ..... ≥5 Mohm.cm	Sulphate ..... ≤0.5 ppm	Cu ..... ≤0.01 ppm	Sn ..... ≤0.02 ppm
Density at 20° C ..... 0.790 ÷ 0.792	Subst. reducing KMn04 ..... ≤2 ppm	Fe ..... ≤0.05 ppm	Sr ..... ≤0.02 ppm
Boiling point ..... 55.8 ÷ 56.3 °	Ag ..... ≤0.02 ppm	Ga ..... ≤0.02 ppm	Ti ..... ≤0.05 ppm
Water (K.F.) ..... ≤0.2 %	Al ..... ≤0.05 ppm	In ..... ≤0.02 ppm	Tl ..... ≤0.05 ppm
Residue on evaporation ..... ≤5 ppm	As ..... ≤0.01 ppm	K ..... ≤0.1 ppm	V ..... ≤0.05 ppm
Acidity (formic acid) ..... ≤15 ppm	Au ..... ≤0.05 ppm	Li ..... ≤0.02 ppm	Zn ..... ≤0.01 ppm
Alcalinity (NH3) ..... ≤2 ppm	B ..... ≤0.01 ppm	Mg ..... ≤0.1 ppm	Zr ..... ≤0.05 ppm
Ethyl alcohol ..... ≤100 ppm	Ba ..... ≤0.1 ppm	Mn ..... ≤0.01 ppm	
Isopropyl alcohol ..... ≤500 ppm	Be ..... ≤0.02 ppm	Mo ..... ≤0.05 ppm	

Code	Size	Packaging	Notes
401051	1 l	Glass bottle	
401058	2.5 l	Glass bottle	
401054	5 l	Plastic tank	
401055	5 l	Metal tank	
401052	22 kg	Metal drum	

Acetone > RS - MOS - For electronic use

RS

Description	Clear liquid	Methyl alcohol	≤500 ppm	Bi	≤0.02 ppm	Na	≤0.2 ppm
Colour (APHA)	≤10	Aldehyde	≤10 ppm	Ca	≤0.1 ppm	Ni	≤0.01 ppm
Identification	Positive	Chloride	≤0.1 ppm	Cd	≤0.005 ppm	Pb	≤0.01 ppm
Water miscibility	Conform	Phosphate	≤0.1 ppm	Co	≤0.005 ppm	Pt	≤0.02 ppm
Assay (GLC)	≥99.8 %	Heavy metals (Pb)	≤0.2 ppm	Cr	≤0.01 ppm	Sb	≤0.01 ppm
Resistivity	≥5 Mohm.cm	Sulphate	≤0.5 ppm	Cu	≤0.01 ppm	Sn	≤0.02 ppm
Density at 20° C	0.790 ÷ 0.792	Subst. reducing KMn04	≤2 ppm	Fe	≤0.05 ppm	Sr	≤0.02 ppm
Boiling point	55.8 ÷ 56.3 °	Ag	≤0.02 ppm	Ga	≤0.02 ppm	Ti	≤0.05 ppm
Water (K.F.)	≤0.2 %	Al	≤0.05 ppm	In	≤0.02 ppm	Tl	≤0.05 ppm
Residue on evaporation	≤5 ppm	As	≤0.01 ppm	K	≤0.1 ppm	V	≤0.05 ppm
Acidity (formic acid)	≤15 ppm	Au	≤0.05 ppm	Li	≤0.02 ppm	Zn	≤0.01 ppm
Alcalinity (NH3)	≤2 ppm	B	≤0.01 ppm	Mg	≤0.1 ppm	Zr	≤0.05 ppm
Ethyl alcohol	≤100 ppm	Ba	≤0.1 ppm	Mn	≤0.01 ppm		
Isopropyl alcohol	≤500 ppm	Be	≤0.02 ppm	Mo	≤0.05 ppm		

Code	Size	Packaging	Notes
401042	1 l	Glass bottle	
401041	2.5 l	Glass bottle	

Acetone > RS - For analysis according to Ph. Eur. Chap. 4.1.3

RS

Code	Size	Packaging	Notes
614000100	1 l	Glass bottle	Buffered acetone solution Ref Ph.Eur 4000100

Acetone > RPE - For analysis - ISO - ACS - Reag.Ph.Eur. - Reag.USP

RPE

Description	Clear liquid	Acidity	≤0.0003 meq/g	B	≤0.1 ppm	Mn	≤0.02 ppm
Colour (APHA)	≤10	Alcalinity	≤0.0006 meq/g	Ba	≤0.1 ppm	Ni	≤0.01 ppm
Identification (I.R.)	Conform	Ethyl alcohol	≤100 ppm	Ca	≤0.5 ppm	Pb	≤0.01 ppm
Water solubility	Conform	Isopropyl alcohol	≤500 ppm	Cd	≤0.05 ppm	Sn	≤0.1 ppm
Density at 20° C	0.790 ÷ 0.792	Methyl alcohol	≤500 ppm	Co	≤0.05 ppm	Zn	≤0.2 ppm
Refractive index at 20°C	1.3581 ÷ 1.3601	Aldehyde	≤10 ppm	Cr	≤0.02 ppm	Assay (GLC)	≥99.8 %
Boiling point	55.8 ÷ 56.3 ° C	Heavy metals (Pb)	≤0.2 ppm	Cu	≤0.01 ppm	Related substances (GLC)	Conform
Water (K.F.)	≤0.2 %	Subst. reducing KMn04	≤2 ppm	Fe	≤0.1 ppm	Benzene	≤ 2 ppm
Residue on evaporation	≤10 ppm	Al	≤0.5 ppm	Mg	≤0.02 ppm	Diacetyl alcohol	≤ 500 ppm

Code	Size	Packaging	Notes
400961	1 l	Plastic bottle	
400971	1 l	Glass bottle	
400965	2.5 l	Plastic bottle	
400974	2.5 l	Glass bottle	
400962	5 l	Plastic tank	
400963	10 l	Plastic tank	
400978	16 kg	Plastic tank	
400972	22 kg	Metal drum	
400979	160 kg	Metal drum	

Acetone > ERBApharm - According to pharmacopoeia: BP-NF-Ph.Eur.

ERBApharm

Description	Clear colourless liquid	Appearance of solution	Conform Ph. Eur.	Impurity C	≤ 2 ppm (v/v)	Assay (GLC)	≥ 99.0 %
Density at 20°C	0.790 ÷ 0.793	Acidity or alcalinity	Conform Ph.Eur.	Other impurities (GC)	≤ 0.05 %(v/v)	Density at 25°C	≤ 0.789
Identification B (Ph.Eur)	Positive	Water insoluble substances	Conform Ph.Eur.	Reducing substances	Conform Ph.Eur.	Water (GLC)	≤ 0.5 %
Identification C (Ph.Eur)	Positive	Related substances (GLC)	Conform Ph.Eur.	Water (K.F.)	≤ 3 g/l	Origin (BSE/TSE)	Synthesis
Identification (I.R.)	Positive	Impurities A,B	≤ 0.05 %v/v	Residue on evaporation	≤ 40 ppm(p/v)	Residual solvents (Current ICH)	Conform

Code	Size	Packaging	Notes
301505	1 l	Glass bottle	
301506	2.5 l	Glass bottle	
301502	5 l	Aluminium can	
301503	5 l	Plastic tank	
301501	16 kg	Plastic tank	
301504	22 kg	Metal drum	
301507	160 kg	Metal drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

DASIT GROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

Acetone > RE - Pure

RE

Description ..... Clear liquid	Density at 20°C ..... 0.788÷0.792	Water (K.F.) ..... ≤0.25 % m/m	Assay (GLC) ..... ≥99.8 % (GLC)
Colour ..... ≤10 APHA	Refractive index at 20°C... 1.3601÷1.3581	Residue on evaporation ..... ≤15 ppm	Acidity (acetic acid)..... ≤ 200 ppm
Identity (IR)..... Positive	Boiling point..... 55.7÷56.7 °C	Water miscibility..... Complete	Diacetyl alcohol ..... ≤ 500 ppm

Code	Size	Packaging	Notes
508200	1 l	Glass bottle	
508201	2.5 l	Glass bottle	
528203	5 l	Plastic tank	
528206	10 l	Plastic tank	
528201	25 l	Plastic tank	
528204	200 l	Metal drum	

Acetone / water 98/2 (v/v) with bromophenol blue 0.02 g/l

• Acetone / Acqua 98/2 (v/v) + Blu di Bromofenolo 0.02 g/l • Acétone / eau 98/2 (v/v) + Bleu de Bromophénol 0.02 g/l
• Acetona / Agua 98/2 (v/v) + Azul de bromofenol 0.02 g/l • Aceton / Wasser 98/2 (v/v) mit Bromphenolblau 0.02 g/l

Classification transport
ONU: 1090
Transport Hazard class: 3
Packing group II

Danger
H225-H319-H336
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Acetone / water 98/2 (v/v) with bromophenol blue 0.02 g/l > RS - For analysis

RS

Code	Size	Packaging	Notes
PS0852/29	5 l	Plastic tank	

Acetone-d6

• Acetone-d6 • Acétone-d6 • Acetona-d6 • Aceton-d6

Synonym:
Hexadeuteroacetone

CD_3COCD_3
Molecular Weight: 64,12
CAS: 666-52-4
EEC-N: 211-563-9

Classification transport
ONU: 1090
Transport Hazard class: 3
Packing group II

Danger
H225-H319-H336-HEU066
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Acetone-d6 > RS - For NMR - min 99.96%

RS

Code	Size	Packaging	Notes
P5060	10 x 0.6 ml	Glass ampoule	

For specifications, contact our customer service for a certificate of analysis

Acetone-d6 > RS - For NMR - min 99.8%

RS

Code	Size	Packaging	Notes
P5049	10 x 0.75 ml	Glass ampoule	
P5044A	10 ml	Glass ampoule	
P5045	25 ml	Glass bottle	
P5046	100 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Acetonitrile

• Acetonitrile • Acétonitrile • Acetonitrilo • Acetonitril

Synonym:
Methyl cyanide

CH_3CN
Molecular Weight: 41,05
CAS: 75-05-8
EEC-N: 200-835-2

Classification transport
ONU: 1648
Transport Hazard class: 3
Packing group II

Danger
H225-H302-H312-H332-H319
P210-P241-P280-P303+P361+P353-P304+P340-P305+P351+P338

Acetonitrile > RS - For UHPLC-MS

RS

Description ..... Clear colourless liquid
Colour ..... ≤ 5 APHA
Identification (I.R.) ..... Positive
Refractive index at 20°C ..... 1.342 - 1.346
Residue on evaporation ..... ≤ 1 ppm
Acidity ..... ≤ 0.0003 meq/g
Alkalinity ..... ≤ 0.0002 meq/g
Assay (CPG) ..... ≥ 99.99 %
Water (K.F.) ..... ≤ 100 ppm

Transmittance
At 191 nm ..... ≥ 40 %
At 195 nm ..... ≥ 80 %
At 200 nm ..... ≥ 95 %
At 215 nm ..... ≥ 97 %
≥ 230 nm ..... ≥ 99 %
Absorbance
At 220 nm ..... ≤ 0.01 AU
At 254 nm ..... ≤ 0.005 AU

Fluorescence (quinine)
At 365 nm ..... ≤ 0.5 ppb
At 450 nm ..... ≤ 0.5 ppb
UHPLC gradient peak
At 210 nm ..... ≤ 0.4 mAU
Drift at 210 nm ..... ≤ 6 mAU
Drift at 254 nm ..... ≤ 2 mAU
Test LC-MS TIC (50-2000m/z) ES I(+)
Sensitive Impurities (reserpine) ..... ≤ 30 ppb

Metals compounds
Al ..... ≤ 20 ppb
Fe ..... ≤ 20 ppb
Na ..... ≤ 50 ppb
Ca ..... ≤ 50 ppb
Mg ..... ≤ 20 ppb
K ..... ≤ 50 ppb

Code	Size	Packaging	Notes
412041	1 l	Glass bottle	
412042	2.5 l	Glass bottle	

Acetonitrile > RS - For LC/MS

RS

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Identification (I.R.) ..... Conform
Refractive index at 20°C ..... 1.342 ÷ 1.346
Residue on evaporation ..... ≤ 2 ppm
Acidity ..... ≤ 0.0005 meq/g
Alkalinity ..... ≤ 0.0002 meq/g

Assay (GLC) ..... ≥ 99.95 %
Transmission UV (1cm, ref water)
At 195 nm ..... ≥ 80 %
At 200 nm ..... ≥ 95 %
At 220 nm ..... ≥ 98 %
≥ 230 nm ..... ≥ 99 %
Fluorescence (quinine)

At 254 nm ..... ≤ 1 ppb
At 365 nm ..... ≤ 0.5 ppb
HPLC gradient
At 210 nm ..... ≤ 1 mAU
Test LC-MS TIC (50-2000m/z) ESI (+)
Sensitive Impurities (reserpine) ..... ≤ 50 ppb
Metals compounds

Al ..... ≤ 50 ppb
Fe ..... ≤ 50 ppb
Na ..... ≤ 50 ppb
Ca ..... ≤ 50 ppb
Mg ..... ≤ 50 ppb
K ..... ≤ 50 ppb
Water (K.F.) ..... ≤ 100 ppm

Code	Size	Packaging	Notes
412341	1 l	Glass bottle	
412342	2.5 l	Glass bottle	

Acetonitrile > RS - For HPLC - GOLD - Ultragradient grade

RS

Description ..... Clear liquid
Colour (APHA) ..... ≤ 10
Identification ..... Positive
Miscb. with Acetone ..... Conform
Water miscibility ..... Conform
Miscibility in ether ..... Conform
Miscibility in methanol ..... Conform
Density at 20°C ..... 0.781 ÷ 0.785

Refractive index at 20°C ..... 1.342 ÷ 1.344
Distillation range ..... 80.5 ÷ 82.5 °C
Water (KF) ..... ≤ 100 ppm
Residue on evaporation ..... ≤ 2 ppm
Acidity ..... ≤ 0.0003 meq/g
Alcalinity ..... ≤ 0.0002 meq/g
Assay (GLC) ..... ≥ 99.9 %
Fluorescence (quinine)

At 254 nm ..... ≤ 1 ppb
At 365 nm ..... ≤ 0.5 ppb
At 450 nm ..... ≤ 0.5 ppb
Absorbance
At 190 nm ..... ≤ 0.6 AU
At 200 nm ..... ≤ 0.03 AU
At 220 nm ..... ≤ 0.01 AU
Transmittance

At 195 nm ..... ≥ 80 %
At 228 nm ..... ≥ 99 %
From 230 to 420 nm ..... ≥ 99 %
Functionality for HPLC
At 210 nm ..... ≤ 1 mAU
drift at 210 nm ..... ≤ 12 mA.U.
HPLC Gradient ..... Passed test
UV cut off ..... ≤ 190 nm

Code	Size	Packaging	Notes
412371000	1 l	Glass bottle	
412372000	2.5 l	Glass bottle	
412374	4 l	Glass bottle	

Acetonitrile > RS - For HPLC PLUS Gradient grade - ACS - Reag.Ph.Eur. - Reag.USP

RS

Description ..... Clear colourless liquid
Colour (APHA) ..... ≤ 10
Identification ..... Positive
Density at 20° C ..... 0.781 ÷ 0.785
Refractive index at 20°C 1.3420 ÷ 1.3443
Distillation range 95% distils between 80 ÷ 82 °C
Titrable acid ..... ≤ 0.0008 meq/g

Titrable base ..... ≤ 0.0002 meq/g
Residue on evaporation ..... ≤ 0.0002 %
Water (K.F.) ..... ≤ 0.01 %
Litmus paper test ..... Conform
Assay (GLC) ..... ≥ 99.9 %
Fluorescence
at 254 nm ..... ≤ 1 ppb

at 365 nm ..... ≤ 0.5 ppb
at 450 nm ..... ≤ 0.5 ppb
U.V. Transmittance
at 195 nm ..... ≥ 79 %
at 200 nm ..... ≥ 90 %
at 210 nm ..... ≥ 95 %
at 220 nm ..... ≥ 98 %

From 240 to 420 nm ..... ≥ 98 %
Absorbance
At 190 nm ..... ≤ 1.00 AU
Absorbance ACS ..... Pass test
Gradient elution ACS ..... pass test
Functionality for HPLC

Code	Size	Packaging	Notes
412393	1 l	Glass bottle PVC coated	
412391000	1 l	Glass bottle	
412392000	2.5 l	Glass bottle	
412395	5 l	Aluminium can	

Product specifications are subject to changes. Please visit our website for updates.

Acetonitrile > RS - For HPLC - Isocratic Grade

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20 °C ..... 0.781 ÷ 0.785
Iodometric test ..... Conform
Refractive index at 20°C. 1.3420 ÷ 1.3440
Boiling point ..... 81.1 ÷ 82.1 °C
Acidity or alkalinity ..... ≤0.0008 meq/g
Water (K.F.) ..... ≤200 ppm
Residue on evaporation ..... ≤2 ppm
UV transmittance at 220 nm ..... ≥ 94 %
Assay (GLC) ..... ≥99.9 %
UV transmittance at 230 nm ..... ≥ 97 %
UV Transmittance
At 197 nm ..... ≥ 82 %
UV transmittance at 240 nm ..... ≥ 98 %
at 200 nm ..... ≥85 %
Free acid (as CH3COOH)... ≤ 0.002 % m/m
at 205 nm ..... ≥89 %
at 210 nm ..... ≥92 %
at 220 nm ..... ≥94 %
at 230 nm ..... ≥98 %
at 240 nm ..... ≥99 %
Fluorescence quinine 254 nm ..... ≤ 1 ppb

Code	Size	Packaging	Notes
412411000	1 l	Glass bottle	
412412000	2.5 l	Glass bottle	
412413000	4 l	Glass bottle	

Acetonitrile > RS - For HPLC 230

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C ..... 1.342 - 1.346
Density d20/20 ..... 0.780 - 0.785
Boiling point ..... 80.0 - 82.5 °C
Colour ..... ≤ 10 Hazen
Non volatile residue ..... ≤ 5 mg/Kg
Free acid (as CH3COOH) ..... ≤ 20 mg/Kg
Assay (GC) ..... ≥ 99.9 %
UV transmittance at 230 nm ..... ≥ 80 %
UV transmittance at 250 nm ..... ≥ 98 %

Code	Size	Packaging	Notes
P00637S16	1 l	Glass bottle	
P00637S21	2.5 l	Glass bottle	

Acetonitrile > RS - For preparative HPLC

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Colour ..... ≤10 APHA
Density at 20° C ..... 0.781 ÷ 0.785
Refractive index at 20°C. 1.3420 ÷ 1.3440
Boiling point ..... 81.1 ÷ 82.1 ° C
Acidity or alkalinity ..... ≤0.0008 meq/g
Water (K.F.) ..... ≤300 ppm
Residue on evaporation ..... ≤5 ppm
Assay (GLC) ..... ≥99.9 %
U.V. Transmittance
at 230 nm ..... ≥50 %
at 250 nm ..... ≥98 %

Code	Size	Packaging	Notes
412409	2.5 l	Glass bottle	
412407	23 l	Metal drum	

Acetonitrile > RS - PESTIPUR - For pesticide analysis

RS

Description ..... Clear liquid
Colour ..... ≤ 10 hazen
Identification ..... Positive
Water ..... ≤ 0.03 %
Acidity (acetic acid) ..... ≤ 20 ppm
Not volatile residue ..... ≤ 5 ppm
GC-ECD (Lindano) ..... ≤ 3 ng/l
Assay (GLC) ..... ≥ 99.9 %

Code	Size	Packaging	Notes
401241	1 l	Glass bottle	
401242	2.5 l	Glass bottle	
401243	4 l	Glass bottle	

Acetonitrile > RS - SPECTROSOL - For optical spectroscopy - Reag.Ph.Eur.

RS

Description ..... Clear colourless liquid
Density at 20° C ..... 0.781 ÷ 0.785
Refractive index at 20°C. 1.3410 ÷ 1.3450
Litmus paper test ..... Conform
Distillation range 95% distils between ..... 80 - 82 °C
Acidity or alkalinity ..... ≤0.0008 meq/g
Water (K.F.) ..... ≤300 ppm
Residue on evaporation ..... ≤5 ppm
Assay (GLC) ..... ≥99.9 %
Fluorescence
at 254 nm ..... ≤1 ppb
U.V. Transmittance
at 200 nm ..... ≥90 %
at 210 nm ..... ≥94 %
at 220 nm ..... ≥96 %
Colour ..... ≤ 10 APHA
at 230 nm ..... ≥98 %
Identification ..... Positive
UV Absorbance from 255 nm to 420 nm.. ≤ 0.01 AU

Code	Size	Packaging	Notes
401216	1 l	Glass bottle	
401212	2.5 l	Glass bottle	

Acetonitrile > RS - Anhydrous - For analysis

RS

Appearance ..... Clear colourless liquid	Water content (K.F.) .....≤ 100 mg/Kg	Free acid (as CH3COOH).....≤ 20 mg/Kg	Free alkali .....≤ 0.0003 meq/g
Refractive index at 20°C.....1.342 - 1.346	Non volatile residue.....≤ 10 mg/Kg	Density d20/20 .....0.780 - 0.785	Propionitrile .....≤ 300 mg/Kg
Identification (IR).....Conform	Colour ..... ≤ 10 Hazen	Iron (Fe).....≤ 0.5 mg/Kg	
Boiling point..... 80.0 - 82.5 °C	Assay (GC) ..... ≥ 99.9 %	Copper (Cu) .....≤ 0.5 mg/Kg	

Code	Size	Packaging	Notes
P0061010	200 ml	Bottle with septum	
P00610S10	200 ml	Bottle with septum	Water content < 50 ppm
P00610T10	200 ml	Bottle with septum	On molecular sieves 3A
P0061016	1 l	Glass bottle	
P00610S16	1 l	Glass bottle	Water content < 50 ppm
P0061021	2.5 l	Glass bottle	
P00610S21	2.5 l	Glass bottle	Water content < 50 ppm

Acetonitrile > RS - For peptide synthesis

RS

Appearance ..... Clear colourless liquid	Density d20/20 .....0.780 - 0.785	Water content (K.F.) .....≤ 30 mg/Kg	Assay (GC) ..... ≥ 99.9 %
Refractive index at 20°C.....1.342 - 1.346	Boiling point..... 80.0 - 82.5 °C	Colour ..... ≤ 10 Hazen	Non volatile residue .....≤ 5 mg/Kg

Code	Size	Packaging	Notes
P0063510	200 ml	Bottle with septum	
P0063516	1 l	Glass bottle	
P0063521	2.5 l	Glass bottle	

Acetonitrile > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Clear colourless liquid	Water (K.F.) .....≤ 500 ppm	K.....≤0.05 ppm	Assay (GLC) ..... ≥99.8 %
Identification ..... Positive	Residue on evaporation .....≤10 ppm	Mg .....≤0.05 ppm	Colour (APHA) ..... ≤10
Density at 20° C .....0.781 ÷ 0.785	Ca.....≤0.2 ppm	Na.....≤1 ppm	Acidity .....≤8 µeq/g
Refractive index at 20°C.1.3410 ÷ 1.3450	Cu.....≤0.05 ppm	Pb.....≤0.05 ppm	Alcalinity.....≤0.6 µeq/g
Distillation range 95% distils between ..... 80 - 82 °C	Fe .....≤0.2 ppm	Zn .....≤0.5 ppn	Litmus paper test .....Conform

Code	Size	Packaging	Notes
401183000	1 l	Glass bottle	
401185000	2.5 l	Glass bottle	

Acetonitrile > RE - Pure

RE

Appearance ..... Clear colourless liquid	Colour ..... Colourless	Water content (K.F.) .....≤ 500 mg/Kg	Assay (GC) ..... ≥ 99.9 %
Refractive index at 20°C.....1.342 - 1.346	Density d20/20 .....0.780 - 0.785	Free acid (as CH3COOH).....≤ 20 mg/Kg	
Identification (IR).....Conform	Boiling point..... 80.0 - 82.5 °C	Non volatile residue .....≤ 20 mg/Kg	

Code	Size	Packaging	Notes
P0060228	5 l	Plastic tank	
P0060248	25 l	Metal drum	
P0060268	200 l	Metal drum	

Product specifications are subject to changes. Please visit our website for updates.

Acetonitrile + 0.1% v/v formic acid

• Acetonitrile + 0.1% v/v acido formico • Acétonitrile + 0.1% v/v d'acide formique • Acetonitrilo+ 0.1% v/v acido formico
• Acetonitril + 0,1% v/v Ameisensäure

CH_3CN
Molecular Weight: 41,05
CAS: 75-05-8

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II

Danger
H225-H302-H332-H319
P210-P241-P280-P303+P361+P353-P304+P340-P305+P351+P338

Acetonitrile + 0.1% v/v formic acid > RS - For LC/MS

RS

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Acidity (formic acid) ..... 0.095 ÷ 0.105 %
HPLC Gradient
At 254 nm ..... ≤ 50 mAU
Transmittance
At 210 nm ..... ≥ 5 %
At 230 nm ..... ≥ 15 %
Assay (CPG) ..... ≥ 99.5 %
Test LC-MS TIC (100-2000m/z)
Sensitive Impurities (reserpine) ..... ≤ 50 ppb
Metals content
Al ..... ≤ 0.5 ppm
Fe ..... ≤ 0.5 ppm
Ca ..... ≤ 0.5 ppm
Mg ..... ≤ 0.5 ppm
Na ..... ≤ 2 ppm
K ..... ≤ 0.5 ppm
Raw material used
Acetonitrile LC-MS(code 412342) ..... Batch number
Formic acid 98-99% (code 405820) . Batch number

Code	Size	Packaging	Notes
412331	1 l	Glass bottle	
412332	2.5 l	Glass bottle	

Acetonitrile + 0.1% v/v trifluoroacetic acid

• Acetonitrile + 0.1% v/v acido trifluoroacetico • Acétonitrile + 0.1% v/v d'acide trifluoroacétique • Acetonitrilo+ 0.1% v/v acido trifluoroacético
• Acetonitril + 0.1% v/v trifluoressigsäure

CH_3CN
Molecular Weight: 41,05
CAS: 75-05-8

Classification transport
ONU: 1648
Transport Hazard class: 3
Packing group II

Danger
H225-H302-H332-H319
P210-P241-P280-P303+P361+P353-P304+P340-P305+P351+P338

Acetonitrile + 0.1% v/v trifluoroacetic acid > RS - For LC/MS

RS

Description ..... Clear colourless liquid
Assay (GC) (without TFA) ..... ≥ 99.9 %
Trifluoroacetic acid content (V/V) 0.095 - 0.105 %
Water (K.F.) ..... ≤ 150 ppm
Residue on evaporation ..... ≤ 2 ppm
UV transmittance (1 cm, ref. water)
At 195 nm ..... ≥ 20 %
At 230 nm ..... ≥ 50 %
At 254 nm ..... ≥ 90 %
At 260 nm ..... ≥ 95 %
Fluorescence (quinine)
At 365 nm ..... ≤ 0.5 ppb
HPLC gradient
Drift at 254 nm ..... ≤ 30 mAU
Test LC-MS TIC (50-2000m/z) ES I(+)
Sensitive Impurities (reserpine) ..... ≤ 50 ppb
Metals compounds
Al ..... ≤ 30 ppb
Fe ..... ≤ 50 ppb
Na ..... ≤ 50 ppb
Ca ..... ≤ 50 ppb
Mg ..... ≤ 30 ppb
K ..... ≤ 50 ppb

Code	Size	Packaging	Notes
412321	1 l	Glass bottle	
412322	2.5 l	Glass bottle	

Acetonitrile-d3

• Acetonitrile-d3 • Acétonitrile-d3 • Acetonitrilo-d3 • Acetonitril-d3

Synonym:
Trideuteroacetonitrile

CD_3CN
Molecular Weight: 44,07
CAS: 2206-26-0
EEC-N: 218-616-5

Classification transport
ONU: 1648
Transport Hazard class: 3
Packing group II

Danger
H225-H302-H312-H332-H319
P210-P241-P280-P303+P361+P353-P304+P340-P305+P351+P338

Acetonitrile-d3 > RS - For NMR - min 99.8%

RS

Code	Size	Packaging	Notes
P5070	2 x 0.6 ml	Glass ampoule	
P5079	10 x 0.75 ml	Glass ampoule	
P5073A	5 ml	Glass ampoule	

For specifications, contact our customer service for a certificate of analysis

Acetophenone

• Acetofenone • Acétophénone • Acetofenona • Acetophenon

Synonym:
Methyl phenyl ketone

$C_6H_5COCH_3$
Molecular Weight: 120,15
CAS: 98-86-2
EEC-N: 202-708-7

Warning
H302-H319
P264-P280i-P301+P312a-P305+P351+P338-P337+P313-P501a

Acetophenone > RE - Pure

RE

Description ..... Clear, colourless to pale yellow
Identification ..... Positive
Density at 20° C ..... 1.018 ÷ 1.038
Refractive index at 20°C ..... 1.5283 ÷ 1.5383
Boiling point ..... 200.5 ÷ 203.5 ° C
Assay (GLC) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
301251	1 l	Glass bottle	

Acetylacetone

• Acetilacetone • Acetylacétone • Acetilacetona • Acetylaceton

Synonym:
• *2,4-Pentanedione*
• *Diacetylmethane*

$CH_3COCH_2COCH_3$
Molecular Weight: 100,12
CAS: 123-54-6
EEC-N: 204-634-0

Classification transport
ONU: 2310
Transport Hazard class: 3
Packing group III

Danger
H226-H302-H311-H331
P210-P280-P303+P361+P353-P304+P340-P311a-P330-P361+P364-P403+P233

Acetylacetone > RPE - For analysis

RPE

Description ..... Clear liquid
Identification (I.R.) ..... Conform
Water miscibility ..... Conform
Density at 20° C ..... 0.971 ÷ 0.981
Refractive index at 20°C ..... 1.4510 ÷ 1.4540
Water (K.F.) ..... ≤0.1 %
Residue on evaporation ..... ≤100 ppm
Assay (GLC) ..... ≥99.5 %

Code	Size	Packaging	Notes
400305	100 ml	Glass bottle	
400307	1 l	Glass bottle	

n-Acetyl-L-cysteine

• n-Acetil-L-cisteina • n-Acétyl-L-cystéine • n-Acetil-L-cisteína • n-Acetyl-L-cystein

Synonym:
LNAC

$HSCH_2CH(NHCOCH_3)COOH$
Molecular Weight: 163,19
CAS: 616-91-1
EEC-N: 210-498-3

n-Acetyl-L-cysteine > RPE - For analysis

RPE

Description ..... White crystals
Identification ..... Positive
Ash ..... < 0.5 %
Assay (acidimetric) ..... > 97.5 %

Code	Size	Packaging	Notes
400522	100 g	Glass bottle	

p-Acetylphenetidine

• p-Acetilfenetidina • p-Acétylphénétidine • p-Acetilfenetidina • Phenacetin

Synonym:
• *Phenacetin*
• *Acetophenelide*

$C_2H_5OC_6H_4NHCOCH_3$
Molecular Weight: 179,22
CAS: 62-44-2
EEC-N: 200-533-0

Warning
H302-H312-H332
P261-P264-P271-P280h-P301+P312a-P304+P340

p-Acetylphenetidine > RE - Pure

RE

Description ..... White crystalline powder
Identification ..... Positive
Melting point ..... 133 ÷ 138 ° C
Assay (HPLC) ..... > 96.0 %

Code	Size	Packaging	Notes
300857	1 kg	Plastic bottle	

Acid Red 51 ▶ Erythrosin extra B

Product specifications are subject to changes. Please visit our website for updates.

a b c d e f g h i j k l m n o p q r s t u v w x y z

Acid Red 87 ▶ Eosin Y

Acid Red 91 ▶ Eosin B

Acid Violet 19 ▶ Fuchsin acid

ACN ▶ Acetonitrile

Acridine orange
• Arancio acridina • Orange acridine • Anaranjado de acridina • Acridinorange

Synonym:
3,6-Bis(dimethylamino)acridine hydrochloride

$C_{17}H_{20}ClN_3$
Molecular Weight: 301,82
CAS: 65-61-2
EEC-N: 200-614-0

Acridine orange > RS - For microscopy - C.I. 46005 — RS

Description ..... Orange powder Identification ..... Positive

Code	Size	Packaging	Notes
423461	25 g	Plastic bottle	

ADF Solution
• Soluzione ADF • Solution ADF • Solución ADF • ADF Lösung

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group III

Warning
H315-H319-H412
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

ADF Solution > RPE - For analysis — RPE

Density at 20°C ..... 1.020 ÷ 1.040

Code	Size	Packaging	Notes
526625	2.5 l	Plastic bottle	
526623	10 l	Plastic tank	

Composition: Trimethylcethylammonium bromure: 20 g Sulfuric acid 1N: QSP 1 L according to NF V18-122

Adipic acid
• Acido adipico • Acide adipique • Acido adípico • Adipinsäure

Synonym:
Hexanedioic acid

$HOOC(CH_2)_4COOH$
Molecular Weight: 146,14
CAS: 124-04-9
EEC-N: 204-673-3

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Adipic acid > RPE - For analysis — RPE

Description ..... White crystalline powder
Identification ..... Positive
Solution colour ..... ≤ 5 APHA
Melting point ..... 151 ÷ 153 °C
Water ..... ≤ 0.2 %
Residue on ignition ..... ≤ 0.0002 %
HNO3 ..... ≤ 4 ppm
Assay (acidimetric) ..... ≥ 99 %

Code	Size	Packaging	Notes
401785	250 g	Plastic bottle	

ß-Alanine

• ß-Alanina • ß-Alanine • ß-Alanina • ß-Alanin

Synonym:
3-Aminopropionic acid

$NH_2CH_2CH_2COOH$
Molecular Weight: 89,09
CAS: 107-95-9
EEC-N: 203-536-5

ß-Alanine > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Loss on drying ..... ≤ 0.5 %
Ammonium ..... ≤ 1000 ppm
Chloride ..... ≤ 400 ppm
Heavy metals (Pb) ..... ≤ 10 ppm
Residue on ignition ..... ≤ 0.2 %
Sulphate ..... ≤ 480 ppm
Fe ..... ≤ 30 ppm
Assay (non-aqueos medium) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
413603	50 g	Glass bottle	

Albumin from eggs powder

• Albumina d'uovo polvere • Albumine d'oeuf poudre • Albúmina de huevo polvo • Albumin aus Eipulver

CAS: 9006-59-1
EEC-N: 232-692-7

Albumin from eggs powder > RE - Pure

RE

Description ..... White powder yellowish
Identification ..... Positive
pH ..... 6.0 ÷ 8.0
Water (K.F.) ..... ≤ 8.0 %

Code	Size	Packaging	Notes
413671	1 kg	Plastic bottle	
413672	5 kg	Plastic bottle	

Albumin from eggs, dried

• Albumina d'uovo secca • Albumine d'oeuf sèche • Albúmina de huevo desecada • Eiweiß aus Eiern, getrocknet

CAS: 9006-59-1
EEC-N: 232-692-7

Albumin from eggs, dried > RS - For biochemistry

RS

Description ..... White powder yellowish
Identification ..... Positive

Code	Size	Packaging	Notes
413654	100 g	Plastic bottle	
413656	500 g	Plastic bottle	

Alcian Blue 8GS 1%

• Blu alcian 8GS 1% • Bleu alcian 8GS 1% • Azul de alcian 8GS 1% • Alcianblau 8GS 1%

Synonym:
Ingrain Blue 1

$C_{56}H_{68}N_{16}CuS_4Cl_4$
Molecular Weight: 1298,88
CAS: 75881-23-1

HEU210

Alcian Blue 8GS 1% > RS - For microscopy

RS

Description ..... Blue clear liquid
Maximum lambda absorption max (DMSO). 672 - 680 nm
Absorbance(Lambda max;DMSO;dil.0. 005g/L) ..... 0.35 - 0.55
pH at 20°C ..... 2.4 - 2.6

Code	Size	Packaging	Notes
428551	250 ml	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Alcian blue 8GX

• Blu alcian 8GX • Bleu alcian 8GX • Azul de alcian 8GX • Alcianblau 8GX

Synonym:
Ingrain Blue 1

$C_{56}H_{68}N_{16}CuS_4Cl_4$
Molecular Weight: 1298,88
CAS: 33864-99-2
EEC-N: 251-705-7

Alcian blue 8GX > RS - For microscopy - C.I. 74240

RS

Description ..... Violet crystalline powder Identification ..... Positive

Code	Size	Packaging	Notes
428561	25 g	Glass bottle	

Dye for histochemistry

Alizarin

• Alizarina • Alizarine • Alizarina • Alizarin

Synonym:
• *1,2-Dihydroxyanthraquinone*
• *Mordant Red 11*

$C_{14}H_8O_4$
Molecular Weight: 240,21
CAS: 72-48-0
EEC-N: 200-782-5

Warning
H302
P264-P270-P301+P312a-P330-P501a

Alizarin > RPE - For analysis - C.I. 58000

RPE

Description ..... Orange red powder Loss on drying ..... ≤2 % Assay ..... ≥96.0 %
Identification ..... Positive Alluminium sensitivity ..... ≥0.1 µg/ml

Code	Size	Packaging	Notes
415892	25 g	Glass bottle	

Dye for microscopy. Indicator acid - base (pH 5.8 to 7.2 - 11.0 to 13.0)

Alizarin red

• Rosso alizarina • Rouge Alizarine • Rojo Alizarina • Alizarinrot S

Synonym:
• *Alizarin sulfonic acid sodium salt*
• *3,4-Dihydroxy-9,10-dioxo-2-anthracenesulfonic acid sodium salt*

$C_{14}H_7NaO_7S$
Molecular Weight: 342,26
CAS: 130-22-3
EEC-N: 204-981-8

Alizarin red > RPE - For analysis - C.I. 58005

RPE

Description ..... Brown orange powder Identification ..... Positive Colour change ..... yellow violet pH range ..... 5.0 - 6.6

Code	Size	Packaging	Notes
416002	25 g	Glass bottle	

Alizarin saturated solution in ethanol

• Alizarina soluzione satura in alcole etilico • Alizarine solution saturée dans l'éthanol
• Alizarina solución saturada en alcohol etilico • Alizarin gesättigte Lösung in Ethanol

Synonym:
• *1,2-Dihydroxyanthraquinone*
• *Mordant Red 11*

$C_{14}H_8O_4$
Molecular Weight: 240,21
CAS: 72-48-0

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II

Danger
H225-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

Alizarin saturated solution in ethanol > RPE - For analysis

RPE

Description ..... Orange brown liquid Identification ..... Positive

Code	Size	Packaging	Notes
E415932	250 ml	Bottle	

Indicator acid - base. Indicator for absorption and complexometry. Saturated alcoholic solution

Alizarin yellow R

• Giallo alizarina R • Jaune d'alizarine R • Amarillo de alizarina R • Alizaringelb R

Synonym:
• *Mordant Orange 1*
• *5-(4-Nitrophenylazo)salicylic acid*

$C_{13}H_9N_3O_5$
Molecular Weight: 287,23
CAS: 2243-76-7
EEC-N: 218-818-3

Warning
H302-H319
P264-P280i-P301+P312a-P305+P351+P338-P337+P313-P501a

Alizarin yellow R > RPE - For analysis - C.I. 14030

RPE

Description ..... Brown crystalline powder
Identification ..... Positive
Loss on drying ..... ≤15 %
Residue on ignition ..... 18.0 ÷ 28.0 %
Colour change ..... yellow-orange
E (1% ö 1 cm) a 492 nm ..... 700 ÷ 1000

Code	Size	Packaging	Notes
453451	10 g	Glass bottle	

Alkali blue 6B

• Blu alcali 6B • Bleu alcalin 6B • Azul de alcali 6B • Alkaliblau 6B

Synonym:
Acid blue 10

$C_{37}H_{29}N_3O_3S$
Molecular Weight: 595,72
CAS: 1324-76-1
EEC-N: 215-385-2

Alkali blue 6B > RS - For microscopy - C.I. 42765

RS

Description ..... Brown violet powder
Identification ..... Positive

Code	Size	Packaging	Notes
428532	25 g	Glass bottle	

Dye for cytology

Alkali Blue 6B solution 2% in ethanol

• Blu alcali 6B soluzione 2% in alcole etilico • Bleu alcalin 6B solution 2% dans l'éthanol
• Azul de alcali 6B solución 2% en alcohol etilico • Alkaliblaue 6B 2% ige Lösung in Ethanol

Synonym:
Acid blue 10

$C_{32}H_{28}N_3NaO_4S$
Molecular Weight: 613,72
CAS: 1324-76-1

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II

Danger
H225-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

Alkali Blue 6B solution 2% in ethanol > RPE - For analysis

RPE

Description ..... Blue liquid
Identification ..... Positive

Code	Size	Packaging	Notes
E428541	250 ml	Bottle	

Alkylbenzyldimethylammonium chloride ▶ Benzalkonium chloride

Allylthiourea

• Alliltiourea • Allylthiourée • Aliltiourea • Allylthioharnstoff

Synonym:
• *1-Allyl-2-thiourea*
• *Thiosinamine*

CH_2:$CHCH_2NHCSNH_2$
Molecular Weight: 116,19
CAS: 109-57-9
EEC-N: 203-683-5

Classification transport
ONU: 2811
Transport Hazard class: 6.1
Packing group III

Danger
H301
P264-P270-P301+P310a-P330-P405-P501a

Allylthiourea > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Melting point ..... 70 ÷ 78 ° C
Perdita essiccamento (50°C) ..... ≤ 1.5 %
Assay (ex nitrogen) ..... 97.5 ÷ 102.5 %

Code	Size	Packaging	Notes
416281	25 g	Glass bottle	
416283	100 g	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Almond oil ▶ Oil refined of almonds

Alumina ▶ Aluminum hydroxide

Alumina white

• Allumina bianca • Alumine blanche • Alúmina blanca • Weißes Aluminiumoxid

Synonym:
Alumina

Al_2O_3
Molecular Weight: 101,96
CAS: 1344-28-1
EEC-N: 215-691-6

Alumina white > RS - For metallography (Medium metals) — RS

Description ..... White suspension
Identification ..... Positive

Code	Size	Packaging	Notes
416531	120 g	Plastic bottle	

Aluminum, powder

• Alluminio, polvere • Aluminium, poudre • Aluminio, polvo • Aluminium, pulver

Al
Molecular Weight: 26,98
CAS: 7429-90-5
EEC-N: 231-072-3

Classification transport
ONU: 1396
Transport Hazard class: 4.3
Packing group II

Danger
H250-H261
P210-P222-P223-P231a+P232-P280-P402+P404

Aluminum, powder > RPE - For analysis — RPE

Description ..... Grey powder
Identification ..... Positive
Cu ..... ≤ 0.03 %
Fe ..... ≤ 0.6 %
Zn ..... ≤ 0.08 %
Assay (complexometric) ..... ≥ 95 %

Code	Size	Packaging	Notes
416817	1 kg	Metallic can	
416815	25 kg	Plastic bucket	

Aluminum standard solution

• Alluminio standard soluzione • Aluminium standard solution • Aluminum, solución patrón • Aluminium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Aluminum standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2 — RS

Code	Size	Packaging	Notes
615000200	100 ml	Plastic bottle	A 200 ppm solution Ref Ph.Eur 5000200
615000201	100 ml	Plastic bottle	A 10 ppm solution: to dilute according to Ref Ph.Eur 5000201
615000202	100 ml	Plastic bottle	A 2 ppm solution: to dilute according to Ref Ph.Eur 5000202
615000203	100 ml	Plastic bottle	A 100 ppm solution: to dilute according to Ref Ph.Eur 5000203

Aluminum standard solution > RS - Standard solution for ICP-MS — RS

Code	Size	Packaging	Notes
505307	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505308	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505309	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Aluminum standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503411	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503413	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503415	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid
503417	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Aluminum standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
504190	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504186	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
E497405	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
E497401	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Aluminum standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
416581		Plastic ampoule	conc. 1.000 ppm Matrix: Nitric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Aluminum ammonium sulfate dodecahydrate

• Alluminio ammonio solfato dodecaidrato • Aluminum ammonium sulfate dodécahydraté
• Aluminio y amonio sulfato dodecahidrato • Aluminiumammoniumsulfat-Dodecahydrat

Synonym:
• *Ammonium alum*
• *Ammonium aluminum sulfate dodecahydrate*

$AlNH_4(SO_4)_2.12H_2O$
Molecular Weight: 453,34
CAS: 7784-26-1

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Aluminum ammonium sulfate dodecahydrate > RPE - For analysis

RPE

Description .....White crystals Identification ..... Positive Assay (complexometric) ..... ≥ 98.5 % Loss on drying (300°C)..... 45 ÷ 48 %

Code	Size	Packaging	Notes
416895	500 g	Plastic bottle	
416897	1 kg	Plastic bottle	
416892	25 kg	Plastic bucket	

Aluminum ammonium sulfate dodecahydrate > RE - Pure

RE

Description . White semitransparent crystals Chloride..... ≤50 ppm Heavy metals (Pb)..... ≤100 ppm Assay (complexometric) ..... ≥97 %
Identification ..... Positive Water-insoluble matter ..... ≤500 ppm Fe ..... ≤30 ppm

Code	Size	Packaging	Notes
311009	5 kg	Plastic tank	
311002	25 kg	Plastic bucket	

Product specifications are subject to changes. Please visit our website for updates.

Aluminum chloride anhydrous

• Alluminio cloruro anidro • Aluminium chlorure anhydre • Aluminio cloruro anhidro • Aluminiumchlorid. wasserfrei

$AlCl_3$
Molecular Weight: 133,34
CAS: 7446-70-0
EEC-N: 231-208-1

Classification transport
ONU: 1726
Transport Hazard class: 8
Packing group II

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Aluminum chloride anhydrous > RE - Pure

RE

Description ... Yellow grey cristalline powder
Identification ..... Positive
Fe ..... ≤ 100 ppm
Assay (complexometric) ..... ≥ 99.0 %

Code	Size	Packaging	Notes
416996	500 g	Glass bottle	

Aluminum chloride hexahydrate

• Alluminio cloruro esaidrato • Aluminium chlorure hexahydraté • Aluminio cloruro hexahidrato • Aluminiumchlorid-Hexahydrat

$AlCl_3.6H_2O$
Molecular Weight: 241,44
CAS: 7784-13-6
EEC-N: 231-208-1

Warning
H302-H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Aluminum chloride hexahydrate > RPE - For analysis - Reag. Ph. Eur.

RPE

Description ..... white or slightly yellow, crystalline powder
..... or colourless crystals, deliquescent.
Identification ..... Positive
pH sol. 5% at 25° C ..... 2.5 ÷ 3.5
Sulphate ..... ≤ 50 ppm
Cd ..... ≤ 50 ppm
Cu ..... ≤ 50 ppm
Fe ..... ≤ 50 ppm
Zn ..... ≤ 50 ppm
Assay (complexometric) ..... ≥ 99 %
Ca ..... ≤ 100 ppm
Co ..... ≤ 50 ppm
K ..... ≤ 100 ppm
Na ..... ≤ 0.05 %
Ni ..... ≤ 50 ppm
Pb ..... ≤ 50 ppm

Code	Size	Packaging	Notes
416942	100 g	Plastic bottle	
416943	500 g	Plastic bottle	
416947	1 kg	Plastic bottle	
416949	5 kg	Plastic tank	
416945	25 kg	Plastic bucket	

Aluminum chloride hexahydrate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP

ERBApharm

Description ..... white or slightly yellow, crystalline powder
Identification A ..... Positive
Identification B ..... Positive
Appearance of solution ..... Conform Ph. Eur.
Water (K.F.) ..... 42.0 ÷ 48.0 %
Alkali and alkaline-earth metals ..... ≤ 0.5 %
Sulfate ..... ≤ 100 ppm
Heavy metals (Pb) ..... ≤ 20 ppm
Fe ..... ≤ 10 ppm
Assay (complexometry) ..... 95.0 - 101.0 %

Code	Size	Packaging	Notes
311257	1 kg	Plastic bottle	
311252	5 kg	Plastic tank	
311256	25 kg	Plastic bucket	
311254	50 kg	Fibre drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Aluminum hydroxide

• Alluminio idrossido • Aluminium hydroxyde • Aluminio hidróxido • Aluminiumhydroxid

$Al(OH)_3$
Molecular Weight: 78
CAS: 21645-51-2
EEC-N: 244-492-7

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Aluminum hydroxide > RPE - For analysis

RPE

Description ..... White powder
Identification (I.R.)..... Positive
Loss on ignition..... ≤ 36.5 %
pH sol. 5%..... 8.5 - 10.0
Soluble matter in water..... ≤ 0.2 %
Assay (Al_2O_3) ..... ≥ 63.5 %

Code	Size	Packaging	Notes
417046	500 g	Plastic bottle	
417047	1 kg	Plastic bottle	

Aluminum hydroxide > RE - Pure

RE

Description ..... White powder
Identification (I.R.)..... Positive
Assay (Al_2O_3) ..... ≥ 60 %

Code	Size	Packaging	Notes
311734	1 kg	Plastic bottle	

Aluminum nitrate nonahydrate

• Alluminio nitrato nonaidrato • Aluminium nitrate nonahydraté • Aluminio nitrato nonahydrate • Aluminiumnitrat-Nonahydrat

$Al(NO_3)_3.9H_2O$
Molecular Weight: 375,13
CAS: 7784-27-2
EEC-N: 236-751-8

Classification transport
ONU: 1438
Transport Hazard class: 5.1
Packing group III

Danger
H272-H315-H319
P210-P280-P305+P351+P338-P332+P313-P362+P364-P337+P313

Aluminum nitrate nonahydrate > RPE - For analysis

RPE

Description ..... Semitransparent crystals
Identification ..... Positive
pH sol. 5% at 20°C ..... 2.5 ÷ 3.5
Chloride..... ≤ 20 ppm
Sulphate..... ≤ 50 ppm
Fe ..... ≤ 50 ppm
Assay ..... ≥ 99 %

Code	Size	Packaging	Notes
417095	100 g	Plastic bottle	
417096	500 g	Plastic bottle	
417097	1 kg	Plastic bottle	

Aluminum nitrate nonahydrate > RE - Pure

RE

Description ..... Cristalli trasparenti
Identification ..... Positive
Assay (complexometric) ..... ≥ 98 %

Code	Size	Packaging	Notes
312007	1 kg	Plastic bottle	
312008	5 kg	Plastic tank	
312001	25 kg	Plastic bucket	

Product specifications are subject to changes. Please visit our website for updates.

Aluminum oxide

• Alluminio ossido • Aluminium oxyde • Aluminio óxido • Aluminiumoxid

Synonym: *Alumina*

Al_2O_3
Molecular Weight: 101,96
CAS: 1344-28-1
EEC-N: 215-691-6

Aluminum oxide > RPE - For analysis - Reag. Ph. Eur.

RPE

Description ..... White powder
Identification ..... Positive
Chloride ..... ≤ 0.005 %
Heavy metals (as Pb) ..... ≤ 0.005 %
Sulfate ..... ≤ 0.05 %
Fe ..... ≤ 0.03 %
Assay (complexometry) ..... ≥ 98 %
Loss on ignition ..... ≤ 0.5 %
Alkalis (as Na20) ..... ≤ 0.2 %
Silicate (as SiO2) ..... ≤ 0.1 %

Code	Size	Packaging	Notes
417144	100 g	Plastic bottle	
417145	250 g	Plastic bottle	
417146	500 g	Plastic bottle	
417147	1 kg	Plastic bottle	

Aluminum oxide > RE - Pure

RE

Description ..... White powder
Identification ..... Positive
Fe2O3 ..... ≤ 250 ppm
Na2O ..... ≤ 0.25 %
Si ..... ≤ 400 ppm
Assay (complexometric) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
312258	2.5 kg	Plastic bucket	
312259	5 kg	Plastic tank	
312252	25 kg	Plastic bucket	

Aluminum oxide (acid)

• Alluminio ossido (acido) • Aluminium oxyde (acide) • Aluminio óxido (ácido) • Aluminiumoxid (Säure)

Synonym: *Alumina*

Al_2O_3
Molecular Weight: 101,96
CAS: 1344-28-1
EEC-N: 215-691-6

Aluminum oxide (acid) > RS - For chromatography according to Brockmann

RS

Description ..... White crystalline powder
Identification ..... Positive
Activity grade 1 ..... Conform

Code	Size	Packaging	Notes
417185	250 g	Plastic bottle	
417182	1 kg	Plastic bottle	

Aluminum oxide (basic)

• Alluminio ossido (basico) • Aluminium oxyde (basique) • Aluminio óxido (básico) • Aluminiumoxid (basisch)

Synonym: *Alumina*

Al_2O_3
Molecular Weight: 101,96
CAS: 1344-28-1
EEC-N: 215-691-6

Aluminum oxide (basic) > RS - For chromatography according to Brockmann

RS

Description ..... White granular powder
Identification ..... Positive
Activity grade 1 ..... Conform

Code	Size	Packaging	Notes
417214	100 g	Plastic bottle	
417217	1 kg	Plastic bottle	

Aluminum oxide (neutral)

• Alluminio ossido (neutro) • Aluminium oxyde (neutre) • Aluminio óxido (neutro) • Aluminiumoxid (neutral)

Synonym: *Alumina*

Al_2O_3
Molecular Weight: 101,96
CAS: 1344-28-1
EEC-N: 215-691-6

Aluminum oxide (neutral) > RS - For chromatography according to Brockmann

RS

Description	White granular powder	Activity grade 1	Conform	Fe2O3	≤ 0.03 %
Identification	Positive	pH sospension 10% H2O	6.5 ÷ 7.5	SiO2	≤ 0.03 %
Na2O	≤ 0.4 %				

Code	Size	Packaging	Notes
417245	250 g	Plastic bottle	
417241	1 kg	Plastic bottle	

Aluminum oxide activated

• Alluminio ossido attivo • Aluminium oxyde actif • Aluminio óxido activo • Aluminiumoxid aktiviert

Synonym: *Alumina*

Al_2O_3
Molecular Weight: 101,96
CAS: 1344-28-1
EEC-N: 215-691-6

Aluminum oxide activated > RE - Pure

RE

Description	Whitish granules	Identification	Positive	Diameter	0.1 ÷ 0.5 mm

Code	Size	Packaging	Notes
312261	1 kg	Plastic bottle	

Aluminum potassium sulfate dodecahydrate

• Alluminio potassio solfato dodecaidrato • Aluminium potassium sulfate dodécahydraté
• Aluminio y potasio sulfato dodecahidrato • Aluminiumkaliumsulfat-Dodecahydrat

Synonym: *Potassium aluminum sulfate dodecahydrate*

$AlK(SO_4)_2.12H_2O$
Molecular Weight: 474,39
CAS: 7784-24-9

Aluminum potassium sulfate dodecahydrate > RPE - For analysis - ACS

RPE

Description	White crystals	Ammonium	≤50 ppm	Fe	≤10 ppm
Identification	Positive	Chloride	≤5 ppm	Na	≤200 ppm
Water-insoluble matter	≤50 ppm	Heavy metals (Pb)	≤10 ppm	Assay (complexometric)	98.0 ÷ 102.0 %

Code	Size	Packaging	Notes
417295	100 g	Plastic bottle	
417296	500 g	Plastic bottle	
417297	1 kg	Plastic bottle	

Aluminum potassium sulfate dodecahydrate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP

ERBApharm

Description	White crystalline powder	pH 10% at 25° C	3.0 ÷ 3.5	Heavy metals (Pb)	≤ 20 ppm
Identification	Positive	Loss on drying	43.0 ÷ 46.0 %	Fe	≤ 100 ppm
Appearance of solution	Conform Ph.Eur.	Ammonium	≤ 0.2 %	Assay (complexometric)	99.0 ÷ 100.5 %

Code	Size	Packaging	Notes
312401	1 kg	Plastic bottle	
312404	5 kg	Plastic tank	
312402	10 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Product specifications are subject to changes. Please visit our website for updates.

Aluminum potassium sulfate dodecahydrate > RE - Pure

RE

Description ..... White crystalline powder
Identification ..... Positive
Water-insoluble matter ..... ≤ 500 ppm
Heavy metals (Pb) ..... ≤ 50 ppm
Fe ..... ≤ 300 ppm
Assay (complexometric) ..... ≥ 97 %

Code	Size	Packaging	Notes
312508	2.5 kg	Plastic bottle	

Aluminum sulfate

• Alluminio solfato • Aluminium sulfate • Aluminio sulfato • Aluminiumsulfat

$Al_2(SO_4)_3.18H_2O$
Molecular Weight: 666,43
CAS: 7784-31-8
EEC-N: 233-135-0

Danger
H318
P280i-P305+P351+P338-P310a

Aluminum sulfate > RPE - For analysis - ACS

RPE

Description ..... White crystals
Identification ..... Positive
Ca ..... ≤ 100 ppm
Water-insoluble matter ..... ≤ 100 ppm
Heavy metals (Pb) ..... ≤ 10 ppm
Chloride ..... ≤ 50 ppm
Fe ..... ≤ 20 ppm
Assay (complexometric) .... 98.0 ÷ 102.0 %
K ..... ≤ 50 ppm
Mg ..... ≤ 20 ppm
Na ..... ≤ 0.02 %

Code	Size	Packaging	Notes
417424	100 g	Plastic bottle	
417425	500 g	Plastic bottle	
417427	1 kg	Plastic bottle	

Aluminum sulfate > RE - Pure

RE

Description ..... White crystals
Identification ..... Positive
Heavy metals (Pb) ..... ≤50 ppm
Fe ..... ≤100 ppm
Assay (complexometric) ..... 16 ÷ 18 %(Al2O3)

Code	Size	Packaging	Notes
312753	5 kg	Plastic tank	
312752	10 kg	Carton box	
312751	25 kg	Plastic bucket	

Amidoschwarz 10B solution

• Amido nero 10B in soluzione • Noir amido 10B solution • Negro Amido 10B solucion • Amidoschwarz 10B losung

CAS: 1064-48-8

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Amidoschwarz 10B solution > RS - For agroalimentary analysis

RS

Density at 20°C ..... 1.006 ÷ 1.012
pH at 20°C ..... 2.30 ÷ 2.50
Absorbance (sol 1/100) at 620 nm ..... 0.700 ÷ 0.730 AU

Code	Size	Packaging	Notes
502050	5 l	Plastic tank	
502051	10 l	Plastic tank	

4-Aminoantipyrine ▶ 4-Aminophenazone

p-Aminobenzenesulfonamide ▶ Sulfanilamide

4-Aminobenzenesulfonic acid ▶ Sulfanilic acid

p-Aminobenzoic acid

• Acido p-amminobenzoico • Acide p-aminobenzoique • Acido p-aminobenzoico • p-Aminobenzoesäure

Synonym:
• *4-Aminobenzoic acid*
• *PABA*

$H_2NC_6H_4COOH$
Molecular Weight: 137,14
CAS: 150-13-0
EEC-N: 205-753-0

H412
P273-P501a

p-Aminobenzoic acid > ERBApharm - According to pharmacopoeia: USP

ERBApharm

Description	Yellowish powder	Organic impurities	Conform USP	Loss on drying	≤0.2 %	Assay (acidimetric)	98.0 ÷ 102.0 % s.s.
Identification A (USP)	Positive	Aniline	≤ 10 ppm	Sulphated ash	≤0.1 %		
Identification B (USP)	Positive	p-Toluidine	≤ 10 ppm	Heavy metals (Pb)	≤20 ppm		

Code	Size	Packaging	Notes
391804	100 g	Plastic bottle	
391805	1 kg	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Aminocyclohexane ▶ Cyclohexylamine

2-Aminoethanol ▶ Ethanolamine

Aminohippuric acid reagent

• Acido amminoippurico reattivo • Réactif à l'acide aminohippurique • Reactivo de ácido aminohipúrico • Reagenz mit Aminohippursäure

Classification transport
ONU: 1993

Aminohippuric acid reagent > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611003701	100 ml	Glass bottle	Ref Ph.Eur 1003701

4-Amino-3-hydroxy-1-naphthalenesulfonic acid ▶ 1-Amino-2-naphthol-4-sulfonic acid

1-Amino-2-naphthol-4-sulfonic acid

• Acido 1-ammino-2-naftolo-4-solfonico • Acide 1-amino-2-naphtol-4-sulfonique
• Acido 1-amino-2-naftol-4-sulfonico • 1-Amino-2-naphthol-4-sulfonsäure

Synonym:
4-Amino-3-hydroxy-1-naphthalenesulfonic acid

$NH_2C_{10}H_5(OH)SO_3H$
Molecular Weight: 239,25
CAS: 116-63-2
EEC-N: 204-147-3

1-Amino-2-naphthol-4-sulfonic acid > RPE - For analysis

RPE

Description	Pink granular powder	Water (K.F.)	≤ 5 %	Assay	≥ 94 %
Identification	Positive	Sulphated ash	≤ 0.5 %		

Code	Size	Packaging	Notes
402032	25 g	Glass bottle	

For the determination of phosphates

Product specifications are subject to changes. Please visit our website for updates.

4-Aminophenazone

• 4-Amminofenazone • 4-Aminophénazone • 4-Aminofenazona • 4-Aminophenazon

Synonym:
• *4-Aminoantipyrine*
• *4-Amino-2,3-dimethyl-1-phenyl-3-pyrazolin-5-one*

$C_6H_5NN(CH_3)C(CH_3):C(NH_2)CO$
Molecular Weight: 203,25
CAS: 83-07-8
EEC-N: 201-452-3

Warning
H302
P264-P270-P301+P312a-P330-P501a

4-Aminophenazone > RS - For phenols detection — RS

Description ..... Yellowish crystalline powder
Identification ..... Positive
Melting point ..... 105.5 ÷ 110 ° C
Loss on drying ..... ≤ 1.5 %
Residue on ignition ..... ≤ 0.1 %
Assay (non-aqueos medium) ..... ≥ 97.5 % (s.s.)

Code	Size	Packaging	Notes
418381	25 g	Glass bottle	

m-Aminophenol

• m-Amminofenolo • m-Aminophenol • m-Aminofenol • m-Aminophenol

Synonym:
3-Aminophenol

$NH_2C_6H_4OH$
Molecular Weight: 109,13
CAS: 591-27-5
EEC-N: 209-711-2

Classification transport
ONU: 2512
Transport Hazard class: 6.1
Packing group III

Warning
H302-H332-H411
P261-P264-P271-P301+P312a-P304+P340-P501a

m-Aminophenol > RE - Pure — RE

Description ..... Yellowish powder
Identification ..... Positive
Melting point ..... 122 ÷ 126 ° C
Assay (ex nitrogen) ..... ≥98 %

Code	Size	Packaging	Notes
418564	100 g	Glass bottle	

p-Aminophenol

• p-Amminofenolo • p-Aminophénol • p-Aminofenol • p-Aminophenol

Synonym:
• *4-Aminophenol*
• *4-Hydroxyaniline*

$NH_2C_6H_4OH$
Molecular Weight: 109,13
CAS: 123-30-8
EEC-N: 204-616-2

Classification transport
ONU: 2512
Transport Hazard class: 6.1
Packing group III

Warning
H302-H332-H341-H410
P261-P264-P271-P280-P304+P340-P308+P313

p-Aminophenol > RE - Pure — RE

Description ..... White powder or yellow
Identification ..... Positive
Melting point ..... 185 ÷ 195 °C
Assay (GLC) ..... ≥ 96.0 %

Code	Size	Packaging	Notes
418594	100 g	Glass bottle	

2-Aminopropane ▶ Isopropylamine

Aminopyrazolone solution

• Amminopirazolone soluzione • Aminopyrazolone en solution • Aminopirazolone en solución • Aminopyrazol in Lösung

$C_{13}H_{17}N_3O$
Molecular Weight: 231,29
CAS: 58-15-1

Aminopyrazolone solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1 — RS

Code	Size	Packaging	Notes
611004601	100 ml	Glass bottle	

Amman's lactophenol solution

• Lattofenolo d'Amman soluzione • Lactophénol d'Amman en solution • Lactofenol de Amman solución • Lactophenol von Amman in Lösung

Classification transport
ONU: 2810
Transport Hazard class: 6.1
Packing group II



Danger
H302-H331-H314-H341-H373
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P308+P313-P403+P233

Amman's lactophenol solution > RS - For microscopy

RS

Description ..... Amber liquid
Identification ..... Positive
Density at 20°C ..... 1.155 ÷ 1.159

Code	Size	Packaging	Notes
457531	100 ml	Glass bottle	

Dye for bacteriology. Contains phenol

Ammonia solution 32%

• Ammoniaca soluzione 32% • Ammoniaque solution 32% • Amonio hidróxido 32% • Ammoniaklösung 32%

Synonym:
Ammonium hydroxide solution

NH_4OH
Molecular Weight: 35,05
CAS: 1336-21-6
EEC-N: 215-647-6

Classification transport
ONU: 2672
Transport Hazard class: 8
Packing group III



Danger
H302-H314-H335-H400
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Ammonia solution 32% > RE - Pure

RE

Description ..... Clear liquid
Color ..... ≤ 10 APHA
Identification ..... Positive
Density at 20°C ..... 0.880 ÷ 0.898
Assay (alcalimetric) ..... 28 ÷ 34 %

Code	Size	Packaging	Notes
528503	5 l	Plastic bottle	
528501	18 kg	Plastic tank	

Ammonia solution 30%

• Ammoniaca soluzione 30% • Ammoniaque solution 30% • Amonio hidróxido 30% • Ammoniaklösung 30%

Synonym:
Ammonium hydroxide solution

NH_4OH
Molecular Weight: 35,05
CAS: 1336-21-6
EEC-N: 215-647-6

Classification transport
ONU: 2672
Transport Hazard class: 8
Packing group III



Danger
H302-H314-H335-H400
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Ammonia solution 30% > RS - RSE - For electronic use

RS

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Density at 20° C ..... 0.892 ÷ 0.910
Carbonate ..... ≤5 ppm
Chloride ..... ≤0.25 ppm
Phosphate ..... ≤0.2 ppm
Heavy metals (Pb) ..... ≤0.2 ppm
Residue on ignition ..... ≤3 ppm
Subst. reducing KMn04 ..... ≤5 ppm
Sulphate ..... ≤1 ppm
Assay (alkalimetric) ..... 28 ÷ 32 %

Ag ..... ≤0.02 ppm
Al ..... ≤0.05 ppm
As ..... ≤0.025 ppm
Au ..... ≤0.05 ppm
B ..... ≤0.01 ppm
Ba ..... ≤0.1 ppm
Be ..... ≤0.02 ppm
Bi ..... ≤0.02 ppm
Ca ..... ≤0.2 ppm
Cd ..... ≤0.01 ppm
Co ..... ≤0.01 ppm
Cr ..... ≤0.01 ppm

Cu ..... ≤0.01 ppm
Fe ..... ≤0.03 ppm
Ga ..... ≤0.02 ppm
In ..... ≤0.02 ppm
K ..... ≤0.2 ppm
Li ..... ≤0.02 ppm
Mg ..... ≤0.1 ppm
Mn ..... ≤0.01 ppm
Mo ..... ≤0.05 ppm
Na ..... ≤0.5 ppm
Ni ..... ≤0.01 ppm
Pb ..... ≤0.01 ppm

Pt ..... ≤0.05 ppm
Sb ..... ≤0.05 ppm
Sn ..... ≤0.02 ppm
Sr ..... ≤0.02 ppm
Ta ..... ≤0.1 ppm
Ti ..... ≤0.05 ppm
Tl ..... ≤0.05 ppm
V ..... ≤0.05 ppm
Zn ..... ≤0.05 ppm
Zr ..... ≤0.05 ppm

Code	Size	Packaging	Notes
420071	1 l	Plastic bottle	
420073	2 l	Glass bottle	
420077	5 l	Plastic bottle	
420075	25 kg	Combined drum	

Product specifications are subject to changes. Please visit our website for updates.

Ammonia solution 30% > RPE - For analysis - ACS

RPE

Description	Clear colourless liquid	Heavy metals (Pb)	≤0.4 ppm	Cd	≤0.01 ppm	Na	≤1 ppm
Colour (APHA)	≤10	Subst. reducing KMnO4	≤8 ppm(5m)	Co	≤0.01 ppm	Ni	≤0.02 ppm
Identification	Positive	Silicate	≤10 ppm	Cr	≤0.02 ppm	Pb	≤0.02 ppm
Assay (alkalimetric)	28.0 ÷ 30.0 %	Sulphide	≤0.1 ppm	Cu	≤0.02 ppm	Zn	≤0.05 ppm
Density at 20° C	0.892 ÷ 0.910	Sulphate	≤2 ppm	Fe	≤0.05 ppm	Nitrate	≤2 ppm
Carbonate	≤10 ppm	Ag	≤0.02 ppm	K	≤0.2 ppm	Residue on ignition	≤20 ppm
Chloride	≤0.5 ppm	As	≤0.02 ppm	Mg	≤0.1 ppm		
Phosphate	≤0.3 ppm	Ca	≤0.5 ppm	Mn	≤0.01 ppm		

Code	Size	Packaging	Notes
419941	1 l	Glass bottle	
419943	2 l	Glass bottle	
419948	2 l	Plastic bottle	
419945	5 l	Plastic bottle	
419946	25 kg	Combined drum	

Ammonia solution 30% > RE - Pure

RE

Description	Clear liquid	Residue on evaporation	≤0.1 %	Sulphate	≤500 ppm
Identification	Positive	Chloride	≤300 ppm	Fe	≤20 ppm
Density at 15° C	0.89 ÷ 0.91	Heavy metals (Pb)	≤50 ppm	Assay (alkalimetric)	28 ÷ 32 %

Code	Size	Packaging	Notes
314873	2 l	Glass bottle	
314871	25 kg	Combined drum	

Ammonia solution 28%

• Ammoniaca soluzione 28% • Ammoniaque solution 28% • Amonio hidróxido 28% • Ammoniaklösung 28%

Synonym: *Ammonium hydroxide solution*

NH_4OH
Molecular Weight: 35,05
CAS: 1336-21-6
EEC-N: 215-647-6

Classification transport
ONU: 2672
Transport Hazard class: 8
Packing group III



Danger
H302-H314-H335-H400
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Ammonia solution 28% > ERBApharm - According to pharmacopoeia: NF-FU-Ph.Eur.

ERBApharm

Description	Clear colourless liquid	Carbonate	≤ 60 ppm	Fe	≤ 0.25 ppm	Ph.Eur.	
Identification	Positive	Chloride	≤ 1 ppm	Non volat.substances	≤ 0.002 % m/v	Pyridine and homologues	Conform Ph.Eur.
Appearance of solution	Conform Ph.Eur.	Sulphate	≤ 5 ppm	Assay (alkalimetric)	27.0 ÷ 30.0 % NH3	Origin (BSE/TSE)	Synthesis
Density at 20° C	0.892 ÷ 0.910	Heavy metals (Pb)	≤ 1 ppm	Ready oxidizable substances	Conform		

Code	Size	Packaging	Notes
314861	1 l	Glass bottle	
314863	2 l	Glass bottle	
314866	25 kg	Combined drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Ammonia solution 25%

• Ammoniaca soluzione 25% • Ammoniaque solution 25% • Amonio hidróxido 25% • Ammoniaklösung 25%

Synonym:
Ammonium hydroxide solution

NH_4OH
Molecular Weight: 35,05
CAS: 1336-21-6
EEC-N: 215-647-6

Classification transport
ONU: 2672
Transport Hazard class: 8
Packing group III

Danger
H302-H314-H335-H400
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Ammonia solution 25% > RS - RSE - For electronic use

RS

Description	Clear liquid	Assay (alkalimetric)	24.0 ÷ 26.0 %	Bi	≤0.02 ppm	Pb	≤0.01 ppm
Density at 20° C	0.901 ÷ 0907 g/ml	Heavy metals (Pb)	≤0.2 ppm	Ca	≤0.2 ppm	Pt	≤0.05 ppm
Colour (APHA)	≤10	Residue on ignition	≤3 ppm	Cd	≤0.01 ppm	Sb	≤0.05 ppm
Sulphated ash	≤ 5 ppm	Subst. reducing KMn04	≤5 ppm	Co	≤0.01 ppm	Sn	≤0.02 ppm
Identification	Positive	Total sulphur	≤1 ppm	Ga	≤0.02 ppm	Sr	≤0.02 ppm
Fe	≤ 50 ppb	Ag	≤0.02 ppm	In	≤0.02 ppm	Ta	≤0.1 ppm
Carbonate	≤5 ppm	Al	≤0.05 ppm	K	≤0.2 ppm	Ti	≤0.05 ppm
Cu	≤ 10 ppb	As	≤0.025 ppm	Li	≤0.02 ppm	Tl	≤0.05 ppm
Chloride	≤0.25 ppm	Au	≤0.05 ppm	Mg	≤0.1 ppm	V	≤0.05 ppm
Cr	≤ 10 ppb	B	≤0.01 ppm	Mn	≤0.01 ppm	Zn	≤0.05 ppm
Ni	≤ 10 ppb	Ba	≤0.1 ppm	Mo	≤0.05 ppm	Zr	≤0.05 ppm
Phosphate	≤0.2 ppm	Be	≤0.02 ppm	Na	≤0.5 ppm		

Code	Size	Packaging	Notes
420085	5 l	Plastic bottle	
420084	25 kg	Combined drum	
420031	170 kg	Plastic drum	

Ammonia solution 25% > RS - MOS - For electronic use

RS

Description	Clear liquid	Ag	≤0.02 ppm	Cu	≤0.01 ppm	Pt	≤0.05 ppm
Colour (APHA)	≤10	Al	≤0.05 ppm	Fe	≤0.03 ppm	Sb	≤0.05 ppm
Identification	Positive	As	≤0.025 ppm	Ga	≤0.02 ppm	Sn	≤0.02 ppm
Density at 20° C	0.901 ÷ 0.907	Au	≤0.05 ppm	In	≤0.02 ppm	Sr	≤0.02 ppm
Carbonate	≤5 ppm	B	≤0.01 ppm	K	≤0.2 ppm	Ta	≤0.1 ppm
Chloride	≤0.25 ppm	Ba	≤0.1 ppm	Li	≤0.02 ppm	Ti	≤0.05 ppm
Phosphate	≤0.2 ppm	Be	≤0.02 ppm	Mg	≤0.1 ppm	Tl	≤0.05 ppm
Heavy metals (Pb)	≤0.2 ppm	Bi	≤0.02 ppm	Mn	≤0.01 ppm	V	≤0.05 ppm
Residue on ignition	≤3 ppm	Ca	≤0.2 ppm	Mo	≤0.05 ppm	Zn	≤0.05 ppm
Subst. reducing KMn04	≤5 ppm	Cd	≤0.01 ppm	Na	≤0.5 ppm	Zr	≤0.05 ppm
Total sulphur	≤1 ppm	Co	≤0.01 ppm	Ni	≤0.01 ppm		
Assay (alkalimetric)	24.0 ÷ 26.0 %	Cr	≤0.01 ppm	Pb	≤0.01 ppm		

Code	Size	Packaging	Notes
420051	1 l	Plastic bottle	
420052	2.5 l	Plastic bottle	

Ammonia solution 25% > RPE - For analysis

RPE

Description	Clear liquid	Phosphate	≤0.3 ppm	As	≤0.02 ppm	K	≤0.2 ppm
Colour (APHA)	≤10	Heavy metals (Pb)	≤0.4 ppm	Ca	≤0.5 ppm	Mg	≤0.1 ppm
Identification	Positive	Residue on ignition	≤3 ppm	Cd	≤0.01 ppm	Mn	≤0.01 ppm
Density at 20° C	0.904 ÷ 0.910	Subst. reducing KMn04	≤8 ppm(5m)	Co	≤0.01 ppm	Na	≤1 ppm
Assay (alkalimetric)	24.0 ÷ 26.0 %	Sulphide	≤0.1 ppm	Cr	≤0.02 ppm	Ni	≤0.02 ppm
Carbonate	≤10 ppm	Total sulphur	≤1 ppm	Cu	≤0.02 ppm	Pb	≤0.02 ppm
Chloride	≤0.5 ppm	Ag	≤0.02 ppm	Fe	≤0.05 ppm	Zn	≤0.05 ppm

Code	Size	Packaging	Notes
419993	2 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Ammonia solution 20 - 22%

• Ammoniaca soluzione 20 - 22% • Ammoniaque solution 20 - 22% • Amonio hidróxido 20 - 22% • Ammoniaklösung 20 - 22 %

Synonym:
Ammonium hydroxide solution

NH_4OH
Molecular Weight: 35,05
CAS: 1336-21-6

Classification transport
ONU: 2672
Transport Hazard class: 8
Packing group III

Danger
H302-H314-H335-H400
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Ammonia solution 20 - 22% > RS - Ultrapure - For trace analysis at ppt level

RS

Description	Clear colourless liquid	Mg	≤ 20 ppt	Er	≤ 10 ppt	Rh	≤ 10 ppt
Identification	Positive	Mn	≤ 10 ppt	Eu	≤ 10 ppt	Rb	≤ 10 ppt
Ag	≤ 10 ppt	Na	≤ 20 ppt	Gd	≤ 10 ppt	Sm	≤ 10 ppt
Al	≤ 20 ppt	Ni	≤ 10 ppt	Ga	≤ 10 ppt	Sc	≤ 10 ppt
As	≤ 10 ppt	Pb	≤ 10 ppt	Ge	≤ 10 ppt	Te	≤ 10 ppt
Ba	≤ 10 ppt	Sn	≤ 10 ppt	Au	≤ 10 ppt	Tb	≤ 10 ppt
Be	≤ 10 ppt	Sr	≤ 10 ppt	Ho	≤ 10 ppt	Th	≤ 10 ppt
Bi	≤ 10 ppt	Tl	≤ 10 ppt	In	≤ 10 ppt	Tm	≤ 10 ppt
Ca	≤ 20 ppt	Ti	≤ 10 ppt	La	≤ 10 ppt	W	≤ 10 ppt
Cd	≤ 10 ppt	Zn	≤ 10 ppt	Li	≤ 10 ppt	U	≤ 10 ppt
Co	≤ 10 ppt	Assay (alkalimetric)	20 ÷ 22 %	Lu	≤ 10 ppt	V	≤ 10 ppt
Cr	≤ 10 ppt	Sb	≤ 10 ppt	Mo	≤ 10 ppt	Yb	≤ 10 ppt
Cu	≤ 20 ppt	Ce	≤ 10 ppt	Nd	≤ 10 ppt	Y	≤ 10 ppt
Fe	≤ 20 ppt	Cs	≤ 10 ppt	Nb	≤ 10 ppt	Zr	≤ 10 ppt
K	≤ 20 ppt	Dy	≤ 10 ppt	Pr	≤ 10 ppt		

Code	Size	Packaging	Notes
420161	500 ml	Plastic bottle	

Ammonia solution 20 - 22% > RS - Superpure - For trace analysis at ppb level

RS

Description	Clear liquid	Mg	≤ 1 ppb	Colour (APHA)	≤ 10	Rh	≤ 0.5 ppb
Identification	Positive	Mn	≤ 0.5 ppb	Ce	≤ 0.1 ppb	Rb	≤ 0.1 ppb
Ag	≤ 0.5 ppb	Mo	≤ 0.5 ppb	Cs	≤ 0.1 ppb	Sm	≤ 0.1 ppb
Al	≤ 1 ppb	Na	≤ 1 ppb	Dy	≤ 0.1 ppb	Sc	≤ 0.1 ppb
As	≤ 1 ppb	Ni	≤ 0.5 ppb	Er	≤ 0.1 ppb	Te	≤ 0.1 ppb
Ba	≤ 0.1 ppb	Pb	≤ 0.1 ppb	Eu	≤ 0.1 ppb	Tb	≤ 0.1 ppb
Be	≤ 0.1 ppb	Sb	≤ 0.5 ppb	Gd	≤ 0.1 ppb	Tl	≤ 0.1 ppb
Bi	≤ 0.1 ppb	Se	≤ 1 ppb	Ga	≤ 0.1 ppb	Tm	≤ 0.1 ppb
Ca	≤ 1 ppb	Sn	≤ 0.5 ppb	Ge	≤ 0.1 ppb	W	≤ 0.1 ppb
Cd	≤ 0.5 ppb	Sr	≤ 0.1 ppb	Au	≤ 0.5 ppb	Yb	≤ 0.1 ppb
Co	≤ 0.5 ppb	Ti	≤ 0.5 ppb	Ho	≤ 0.1 ppb	Y	≤ 0.1 ppb
Cr	≤ 0.5 ppb	V	≤ 0.5 ppb	In	≤ 0.1 ppb	Chloride	≤ 0.5 ppm
Cu	≤ 0.5 ppb	Zn	≤ 0.5 ppb	La	≤ 0.1 ppb	Phosphate	≤ 0.01 ppm
Fe	≤ 1 ppb	Zr	≤ 0.1 ppb	Lu	≤ 0.1 ppb	Sulphate	≤ 1 ppm
Hg	≤ 0.2 ppb	Assay (alkalimetric)	20 ÷ 22 %	Nd	≤ 0.1 ppb		
K	≤ 1 ppb	U	≤ 0.1 ppb	Nb	≤ 0.1 ppb		
Li	≤ 0.1 ppb	Th	≤ 0.1 ppb	Pr	≤ 0.1 ppb		

Code	Size	Packaging	Notes
420175	500 ml	Plastic bottle	

Ammonia solution 20 - 22% > RPE - For analysis

RPE

Description	Clear liquid	Phosphate	≤0.3 ppm	As	≤0.02 ppm	K	≤0.2 ppm
Colour (APHA)	≤10	Heavy metals (Pb)	≤0.4 ppm	Ca	≤0.5 ppm	Mg	≤0.1 ppm
Identification	Positive	Residue on calcination	≤3 ppm	Cd	≤0.01 ppm	Mn	≤0.01 ppm
Density at 20° C	0.917 ÷ 0.923	Subst. reducing KMn04	≤8 ppm(5m)	Co	≤0.01 ppm	Na	≤1 ppm
Assay (alkalimetric)	20 ÷ 22 %	Sulphide	≤0.1 ppm	Cr	≤0.02 ppm	Ni	≤0.02 ppm
Carbonate	≤10 ppm	Total sulphur	≤1 ppm	Cu	≤0.02 ppm	Pb	≤0.02 ppm
Chloride	≤0.5 ppm	Ag	≤0.02 ppm	Fe	≤0.05 ppm	Zn	≤0.05 ppm

Code	Size	Packaging	Notes
419981	1 l	Glass bottle	
419983	2 l	Glass bottle	
419984	25 kg	Combined drum	

Ammonia solution 17%

• Ammoniaca soluzione 17% • Ammoniaque solution 17% • Amonio hidróxido 17% • Ammoniaklösung 17%

Synonym: *Ammonium hydroxide solution*

NH_4OH
Molecular Weight: 35,05
CAS: 1336-21-6

Classification transport
ONU: 2672
Transport Hazard class: 8
Packing group III



Danger
H302-H314-H335-H400
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Ammonia solution 17% > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611004701	250 ml	Plastic bottle	Ref Ph.Eur 1004701

Storage: protected from atmospheric carbon dioxide, at a temperature below 20 °C

Ammonia solution 10%

• Ammoniaca soluzione 10% • Ammoniaque solution 10% • Amonio hidróxido 10% • Ammoniaklösung 10%

Synonym: *Ammonium hydroxide solution*

NH_4OH
Molecular Weight: 35,05
CAS: 1336-21-6

Classification transport
ONU: 2672
Transport Hazard class: 8
Packing group III



Danger
H302-H314-H335-H400
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Ammonia solution 10% > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Density at 15° C .....0.956 ÷ 0.962 Assay (alkalimetric).....9.5 ÷ 10.5 %

Code	Size	Packaging	Notes
E420001	1 l	Plastic bottle	
E420002	5 l	Plastic tank	

Ammonia solution 6N

• Ammonio soluzione 6N • Ammoniaque solution 6N • Amonio solución 6N • Ammoniaklösung 6N

Synonym: *Ammonium hydroxide solution*

NH_4OH
Molecular Weight: 35,05
CAS: 1336-21-6
EEC-N: 215-647-6

Classification transport
ONU: 2672
Transport Hazard class: 8
Packing group III



Danger
H302-H314-H335-H400
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Ammonia solution 6N > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
617000151	1 l	Plastic bottle	

Ammonia solution diluted

• Ammoniaca soluzione diluita • Ammoniaque solution diluée • Amonio hidróxido solución diluido
• Ammoniaklösung verdünnt

Synonym: *Ammonium hydroxide solution*

NH_4OH
Molecular Weight: 35,05
CAS: 1336-21-6



Danger
H315-H318
P264-P280a-P305+P351+P338-P310a-P362+P364-P332+P313

Ammonia solution diluted > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611004702	1 l	Plastic bottle	Ammonia, dilute R1 Ref Ph.Eur 1004702
611004703	1 l	Plastic bottle	Ammonia, dilute R2 Ref Ph.Eur 1004703

Product specifications are subject to changes. Please visit our website for updates.

Ammonia buffer solution pH 10

• Tampone ammoniacale pH 10 • Tampon ammoniacal pH10 • Tampón de amoníaco pH 10 • Ammoniakpufferlösung pH 10

Classification transport
ONU: 1719
Transport Hazard class: 8
Packing group III



Danger
H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Ammonia buffer solution pH 10 > RS - For analysis

RS

pH..... 10.5 - 10.9 unite pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
PS0194/22	5 l	Plastic tank	
PS0194/95	5 l	Kubidos	

Ammoniacal solution of copper tetrammine

• Rame tetrammina soluzione ammoniacale • Solution amoniacale de tétramminecuivre • Cobre tetrammina solución amoniacale • Ammoniaklösung von Tetramminecuivre

Classification transport
ONU: 3266
Transport Hazard class: 8
Packing group III



Danger
H302-H314-H335-H410
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Ammoniacal solution of copper tetrammine > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611022600	100 ml	Glass bottle	Ref Ph.Eur 1022600

Ammonium standard solution

• Ammonio standard soluzione • Ammonium standard solution • Amonio, solución patrón • Ammonium-Standardlösung



Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Ammonium standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615000301	100 ml	Plastic bottle	A 2.5 ppm solution: to dilute according to Ref Ph.Eur 5000301
615000302	100 ml	Plastic bottle	A 1 ppm solution: to dilute according to Ref Ph.Eur 5000302
615000309	100 ml	Plastic bottle	A 100 ppm solution: to dilute according to Ref Ph.Eur 5000300

Ammonium standard solution > RS - Standard solution for ion chromatography

RS

Code	Size	Packaging	Notes
503311	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503313	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Ammonium acetate

• Ammonio acetato • Ammonium acétate • Amonio acetato • Ammoniumacetat

CH_3COONH_4
Molecular Weight: 77,08
CAS: 631-61-8
EEC-N: 211-162-9

Ammonium acetate > RS - For LC/MS

RS

Assay	≥ 98 %	Al	≤ 1 ppm	Fe	≤ 2 ppm	Sr	≤ 1 ppm
Impurities	≤ 50 ppm	As	≤ 0.1 ppm	K	≤ 5 ppm	Zn	≤ 1 ppm
Water	≤ 2 %	Ba	≤ 1 ppm	Li	≤ 1 ppm	Grad. Elution H.Peak at 254 nm	< 0.001 AU
Residue on ignition	≤ 100 ppm	Bi	≤ 1 ppm	Mg	≤ 1 ppm	Grad. Elution drift at 254 nm	< 0.005 AU
pH	6.5 ÷ 7.4	Ca	≤ 5 ppm	Mn	≤ 1 ppm	T260nm (1M)	≥ 98 %
Melting point	112 ÷ 116 °C	Cd	≤ 1 ppm	Mo	≤ 1 ppm	T280nm (1M)	≥ 99 %
Chloride	≤ 5 ppm	Co	≤ 1 ppm	Na	≤ 5 ppm	Preparation	Filtered through 0.1
Nitrate	≤ 10 ppm	Cr	≤ 1 ppm	Ni	≤ 1 ppm		
Sulfate	≤ 10 ppm	Cu	≤ 1 ppm	Pb	≤ 1 ppm		

Code	Size	Packaging	Notes
418781	50 g	Glass bottle	

Additive for eluent phase for LC-MS

Ammonium acetate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description	Crist. bianchi	Water-insoluble matter	≤50 ppm	Nitrate	≤10 ppm	Fe	≤5 ppm
Identification	Positive	Residue on ignition	≤100 ppm	Sulphate	≤10 ppm	Assay (acc.to Sörensen)	≥98 %
pH sol. 5% at 25° C	6.7 ÷ 7.3	Chloride	≤5 ppm	Heavy metals (Pb)	≤5 ppm		

Code	Size	Packaging	Notes
418775	100 g	Plastic bottle	
418776	500 g	Plastic bottle	
418777	1 kg	Plastic bottle	
418772	5 kg	Plastic jar	
418773	5 kg	Plastic bucket	
418771	25 kg	Plastic bucket	

Ammonium acetate > RE - Pure

RE

Description	White semitransparent crystals	Identification	Positive	Water	≤ 2 %	Assay (non-aqueos medium)	≥ 97.5 % (s.s.)

Code	Size	Packaging	Notes
313507	1 kg	Plastic bottle	
313508	5 kg	Plastic tank	
313502	25 kg	Plastic bucket	
313504	50 kg	Fibre drum	

Ammonium aluminum sulfate dodecahydrate ▶ Aluminum ammonium sulfate dodecahydrate

Product specifications are subject to changes. Please visit our website for updates.

Ammonium bicarbonate

• Ammonio bicarbonato • Ammonium bicarbonate • Amonio bicarbonato • Ammoniumbicarbonat

Synonym:
Ammonium hydrogen carbonate

NH_4HCO_3
Molecular Weight: 79,06
CAS: 1066-33-7
EEC-N: 213-911-5

Warning
H302
P264-P270-P301+P312a-P330-P501a

Ammonium bicarbonate > RPE - For analysis

RPE

Description	White crystalline powder	Water-insoluble matter	≤30 ppm	Cu	≤10 ppm	Ni	≤10 ppm
Identification	Positive	Sulphate	≤50 ppm	Fe	≤10 ppm	Pb	≤10 ppm
pH sol. 5% at 25° C	7.0 ÷ 8.0	Residue on ignition	≤500 ppm	K	≤10 ppm	Assay (acidimetric)	≥98.5 %
Chloride	≤5 ppm	As	≤1 ppm	Mg	≤10 ppm		
Phosphate	≤5 ppm	Ca	≤100 ppm	Na	≤20 ppm		

Code	Size	Packaging	Notes
418925	100 g	Plastic bottle	
418926	500 g	Plastic bottle	
418927	1 kg	Plastic bottle	
418929	5 kg	Plastic tank	

Ammonium bicarbonate > RE - Pure

RE

Description	White crystalline powder	Sulphated ash	≤500 ppm	Sulphate	≤150 ppm
Identification	Positive	Chloride	≤50 ppm	Assay (acidimetric)	≥99 %

Code	Size	Packaging	Notes
313601	5 kg	Plastic tank	

Ammonium bifluoride ▶ Ammonium hydrogen difluoride

Ammonium bromide

• Ammonio bromuro • Ammonium bromure • Amonio bromuro • Ammoniumbromid

NH_4Br
Molecular Weight: 97,94
CAS: 12124-97-9
EEC-N: 235-183-8

Ammonium bromide > RPE - For analysis - ACS

RPE

Description	White crystalline powder	Residue on ignition	≤100 ppm	Sulphate	≤50 ppm	Assay (argentimetric)	≥99.0 %
Identification	Positive	Bromate	≤20 ppm	Heavy metals (Pb)	≤5 ppm		
pH sol. 5% at 25° C	4.5 ÷ 6.0	Chloride	≤0.2 %	Ba	≤20 ppm		
Water-insoluble matter	≤50 ppm	Iodide	Conform	Fe	≤5 ppm		

Code	Size	Packaging	Notes
419174	100 g	Plastic bottle	
419175	250 g	Plastic bottle	
419176	500 g	Plastic bottle	
419177	1 kg	Plastic bottle	

Ammonium carbamate

• Ammonio carbammato • Ammonium carbamate • Amonio carbamato • Ammoniumcarbamat

Synonym:
Carbamic acid ammonium salt

NH_2COONH_4
Molecular Weight: 78,07
CAS: 1111-78-0
EEC-N: 214-185-2

Warning
H302
P264-P270-P301+P312a-P330-P501a

Ammonium carbamate > RPE - For analysis

RPE

Description ..... Pezzi irregolari bianchi
Identification ..... Positive
Chloride ..... ≤ 5 ppm
Phosphate ..... ≤ 5 ppm
Nitrate ..... ≤ 10 ppm
Ca ..... ≤ 10 ppm
Cu ..... ≤ 5 ppm
Fe ..... ≤ 5 ppm
K ..... ≤ 50 ppm
Mg ..... ≤ 5 ppm
Na ..... ≤ 50 ppm
Ni ..... ≤ 5 ppm
Pb ..... ≤ 5 ppm
Zn ..... ≤ 5 ppm
Assay (alkalimetric) ..... ≥ 99.5 %
Sulphate ..... ≤ 10 ppm
Cd ..... ≤ 5 ppm
Co ..... ≤ 5 ppm
Cr ..... ≤ 5 ppm
Mn ..... ≤ 5 ppm
Residue on calcination ..... ≤ 20 ppm

Code	Size	Packaging	Notes
419201	100 g	Plastic bottle	
419204	500 g	Plastic bottle	
419202	1 kg	Plastic bottle	

Ammonium carbonate

• Ammonio carbonato • Ammonium carbonate • Amonio carbonato • Ammoniumcarbonat

Synonym:
Hartshorn salt

$(NH_4)_2CO_3$
Molecular Weight: 96,09
CAS: 10361-29-2
EEC-N: 233-786-0

Warning
H302
P264-P270-P301+P312a-P330-P501a

Ammonium carbonate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White crystals
Identification ..... Positive
Water-insoluble matter ..... ≤50 ppm
Non volat.substances ..... ≤100 ppm
Chloride ..... ≤5 ppm
Total sulphur ..... ≤20 ppm
Heavy metals (Pb) ..... ≤5 ppm
Fe ..... ≤5 ppm
Assay (alkalimetric) ..... ≥30.0 %

Code	Size	Packaging	Notes
419235	100 g	Plastic bottle	
419236	500 g	Plastic bottle	
419237	1 kg	Plastic bottle	
419239	5 kg	Plastic tank	
419232	25 kg	Plastic bucket	

Ammonium carbonate solution 158 g/l

• Ammonio carbonato soluzione 158 g/l • Ammonium carbonate solution 158 g/l • Amonio carbonato solución 158 g/l • Ammoniumcarbonatlösung 158 g/l

Molecular Weight: 96,09

HEU210

Ammonium carbonate solution 158 g/l > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611005201	1 l	Plastic bottle	Ref Ph.Eur 1005201

Product specifications are subject to changes. Please visit our website for updates.

Ammonium chloride

• Ammonio cloruro • Ammonium chlorure • Amonio cloruro • Ammoniumchlorid

Synonym:
Salmiac

NH_4Cl
Molecular Weight: 53,49
CAS: 12125-02-9
EEC-N: 235-186-4

Warning
H302-H319
P264-P280i-P301+P312a-P305+P351+P338-P337+P313-P501a

Ammonium chloride > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description	White crystalline powder	Water-insoluble matter	≤50 ppm	Sulphate	≤20 ppm	Fe	≤2 ppm
Identification	Positive	Residue on ignition	≤ 100 ppm	Heavy metals (Pb)	≤5 ppm	Mg	≤5 ppm
pH sol. 5% at 25° C	4.5 ÷ 5.5	Phosphate	≤2 ppm	Ca	≤10 ppm	Assay (argentimetric)	≥99.5 %

Code	Size	Packaging	Notes
419415	100 g	Plastic bottle	
419416	500 g	Plastic bottle	
419417	1 kg	Plastic bottle	
419419	5 kg	Plastic tank	
419412	25 kg	Drum	

Ammonium chloride > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP-DAB

ERBApharm

Description	White crystalline powder	Bromide and Iodide	Conform Ph.Eur.	Sulphated ash	≤0.1 %	Fe	≤20 ppm
Identification	Positive	Thiocyanate	Conform USP-NF	Heavy metals (Pb)	≤10 ppm	Assay (argentimetric)	99.5 ÷ 100.5 % s.s.
Appearance of solution	Conform Ph.Eur.	pH (1:20)	4.6 ÷ 6.0	Sulphate	≤150 ppm	Origin (BSE/TSE)	Synthesis
Acidity or alkalinity	Conform Ph.Eur.	Loss on drying	≤0.5 %	Ca	≤200 ppm	Residual solvents (Current ICH)	Conform

Code	Size	Packaging	Notes
313957	1 kg	Plastic bottle	
313952	2.5 kg	Plastic bottle	
313956	5 kg	Plastic tank	
313951	25 kg	Plastic bucket	
313954	50 kg	Fibre drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Ammonium chloride > RE - Pure

RE

Description	Polvere crist. bianca parz. ammassata	Identification	Positive	Not soluble matter	≤ 0.02 %	Fe	≤ 10 ppm
		Water (K.F.)	≤ 0.1 %	Sulphated ash	≤ 0.3 %	Assay (argentimetric)	≥ 99 %

Code	Size	Packaging	Notes
314002	2.5 kg	Plastic bottle	
314001	25 kg	Plastic bucket	

Ammonium chloride buffer solution pH 10.7

• Tampone cloruro di ammonio pH 10.7 • Tampon chlorure d'ammonium pH 10.7 • Tampón cloruro di amonio pH 10.7 • Puffer Ammoniumchlorid pH 10.7

Classification transport
ONU: 3266
Transport Hazard class: 8
Packing group III

Danger
H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Ammonium chloride buffer solution pH 10.7 > RS - For analysis according to Ph. Eur. Chap. 4.1.3

RS

Code	Size	Packaging	Notes
614013400	1 l	Plastic bottle	Ref Ph.Eur 4013400

Ammonium chloride buffer solution pH 10.0

• Tampone cloruro di ammonio pH 10.0 • Tampon chlorure d'ammonium pH 10.0 • Tampón cloruro di amonio pH 10.0 • Puffer Ammoniumchlorid pH 10.0

Classification transport
ONU: 3266
Transport Hazard class: 8
Packing group III

Danger
H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Ammonium chloride buffer solution pH 10.0 > RS - For analysis according to Ph. Eur. Chap. 4.1.3 — RS

Code	Size	Packaging	Notes
614007301	100 ml	Plastic bottle	Ref Ph.Eur 4007300
614007300	1 l	Plastic bottle	Ref Ph.Eur 4007300

Ammonium chloride buffer solution pH 9.5

• Tampone cloruro di ammonio pH 9.5 • Tampon chlorure d'ammonium pH 9.5 • Tampón cloruro di amonio pH 9.5 • Puffer Ammoniumchlorid pH 9.5

Ammonium chloride buffer solution pH 9.5 > RS - For analysis according to Ph. Eur. Chap. 4.1.3 — RS

Code	Size	Packaging	Notes
614007200	1 l	Plastic bottle	Ref Ph.Eur 4007200

tri-Ammonium citrate

• Ammonio citrato tribasico • tri-Ammonium citrate • Tri-amonio citrato • tri-Ammoniumcitrat

Synonym:
• *Ammonium citrate tribasic*
• *Citric acid triammonium salt*

$HOC(CO_2NH_4)(CH_2CO_2NH_4)_2$
Molecular Weight: 243,22
CAS: 3458-72-8
EEC-N: 222-394-5

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

tri-Ammonium citrate > RE - Pure — RE

Description ..... White crystalline powder
Identification ..... Positive
Chloride ..... ≤ 30 ppm
Sulphate ..... ≤ 150 ppm
Fe ..... ≤ 20 ppm
Assay ..... ≥ 97 %

Code	Size	Packaging	Notes
313895	25 kg	Plastic bucket	

Ammonium citrate dibasic

• Ammonio citrato bibasico • Ammonium citrate dibasique • Amonio citrato dibásico
• Ammoniumcitrat zweibasisch

Synonym:
• *Ammonium hydrogencitrate*
• *Citric acid ammonium salt*

$HOCCOOH(CH_2COONH_4)_2$
Molecular Weight: 226,18
CAS: 3012-65-5
EEC-N: 221-146-3

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Ammonium citrate dibasic > RPE - For analysis - ACS — RPE

Description ..... White crystals
Identification ..... Positive
Water-insoluble matter ..... ≤ 50 ppm
Residue on ignition ..... ≤ 50 ppm
Chloride ..... ≤ 10 ppm
Phosphate ..... ≤ 5 ppm
Oxalate ..... ≤ 500 ppm
Total sulphur ..... ≤ 50 ppm
Heavy metals (Pb) ..... ≤ 5 ppm
Fe ..... ≤ 10 ppm
Assay (acc.to Sörensen) ..... 98.0 ÷ 103.0 %

Code	Size	Packaging	Notes
419315	250 g	Plastic bottle	
419313	500 g	Plastic bottle	
419317	1 kg	Plastic bottle	
419312	25 kg	Plastic bucket	
419316	50 kg	Fibre drum	

Product specifications are subject to changes. Please visit our website for updates.

Ammonium citrate solution 20%

• Ammonio citrato soluzione 20% • Ammonium citrate solution 20% • Amonio citrato solución 20%
• Ammoniumcitratlösung 20%

Synonym:
• *Ammonium hydrogencitrate*
• *Citric acid ammonium salt*

$HOCCOOH(CH_2COONH_4)_2$
Molecular Weight: 226,18
CAS: 3012-65-5

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Ammonium citrate solution 20% > RPE - For analysis

RPE

Description ..... Clear colourless liquid
pH of the substance ..... 7 - 7.3
Assay ..... 19 - 21 %

Code	Size	Packaging	Notes
E419361	1 l	Bottle	

For the determination of phosphates

Ammonium citrate tribasic ▶ tri-Ammonium citrate

Ammonium cobalt(II) sulfate hexahydrate ▶ Cobalt (II) ammonium sulfate hexahydrate

Ammonium dihydrogenphosphate ▶ Ammonium phosphate monobasic

Ammonium fluoride

• Ammonio fluoruro • Ammonium fluorure • Amonio fluoruro • Ammoniumfluorid

NH_4F
Molecular Weight: 37,04
CAS: 12125-01-8
EEC-N: 235-185-9

Classification transport
ONU: 2505
Transport Hazard class: 6.1
Packing group III

Danger
H301-H311-H331
P261-P304+P340-P311a-P330-P361+P364-P403+P233

Ammonium fluoride > RPE - For analysis - ACS

RPE

Description ..... White crystals
Identification ..... Positive
Water-insoluble matter ..... ≤50 ppm
Residue on ignition ..... ≤100 ppm
Chloride ..... ≤10 ppm
Sulphate ..... ≤50 ppm
Heavy metals (Pb) ..... ≤5 ppm
Fe ..... ≤5 ppm
Assay (acc.to Sörensen) ..... ≥98.0 %

Code	Size	Packaging	Notes
419634	100 g	Plastic bottle	
419638	250 g	Plastic bottle	
419635	500 g	Plastic bottle	
419637	1 kg	Plastic bottle	

Ammonium formate

• Ammonio formiato • Ammonium formiate • Amonio formiato • Ammoniumformiat

Synonym:
Formic acid ammonium salt

$HCOONH_4$
Molecular Weight: 63,06
CAS: 540-69-2
EEC-N: 208-753-9

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Ammonium formate > RS - For LC/MS

RS

Assay ..... ≥ 98 %
Impurities ..... ≤ 50 ppm
Water ..... ≤ 0.5 %
pH ..... 5.5 ÷ 7.6
Melting point ..... 119 ÷ 121 °C
Chloride ..... ≤ 5 ppm
Sulfate ..... ≤ 10 ppm
Al ..... ≤ 1 ppm
As ..... ≤ 0.1 ppm
Ba ..... ≤ 1 ppm
Bi ..... ≤ 1 ppm
Ca ..... ≤ 5 ppm
Cd ..... ≤ 1 ppm
Co ..... ≤ 1 ppm
Cr ..... ≤ 1 ppm
Cu ..... ≤ 1 ppm
Fe ..... ≤ 3 ppm
K ..... ≤ 5 ppm
Li ..... ≤ 1 ppm
Mg ..... ≤ 1 ppm
Mn ..... ≤ 1 ppm
Mo ..... ≤ 1 ppm
Na ..... ≤ 5 ppm
Ni ..... ≤ 1 ppm
Pb ..... ≤ 1 ppm
Sr ..... ≤ 1 ppm
Zn ..... ≤ 1 ppm
Grad. Elution H.Peak at 254 nm ≤ 0.001 AU
Grad. Elution drift at 254 nm .... ≤ 0.005 AU
T260nm (1M) ..... ≥ 97 %
Preparation ..... Filtered through 0.1

Code	Size	Packaging	Notes
419741	50 g	Plastic bottle	

Additive for eluent phase for LC-MS

Ammonium formate > RPE - For analysis - Reag. Ph. Eur.

RPE

Description ..... Colourless crystals	Water ..... ≤ 0.5 %	Water-insoluble matter ..... ≤ 50 ppm	Assay (acc.to Sörensen)..... ≥ 98.5 %
Identification ..... Positive	Chloride..... ≤ 5 ppm	Heavy metals (Pb)..... ≤ 5 ppm	Melting point..... 119 ÷ 121 °C
pH sol. 5% at 25° C ..... 5.5 ÷ 7.0	Sulphate..... ≤ 10 ppm	Fe ..... ≤ 5 ppm	

Code	Size	Packaging	Notes
419734	100 g	Plastic bottle	
419735	250 g	Plastic bottle	
419736	500 g	Plastic bottle	
419737	1 kg	Plastic bottle	
419733	25 kg	Plastic bucket	

Ammonium heptamolybdate tetrahydrate ▶ Ammonium molybdate tetrahydrate

Ammonium hydrogen carbonate ▶ Ammonium bicarbonate

di-Ammonium hydrogen citrate ▶ Ammonium citrate dibasic

Ammonium hydrogen difluoride

• Ammonio bifluoruro • Ammonium bifluorure • Amonio bifluoruro • Ammoniumhydrogendifluorid

Synonym:
Ammonium bifluoride

$NH_4F.HF$
Molecular Weight: 57,04
CAS: 1341-49-7
EEC-N: 215-676-4

Classification transport
ONU: 1727
Transport Hazard class: 8
Packing group II

Danger
H301-H314
P280-P301+P310a-P301+P330+P331-
P303+P361+P353-P304+P340-P305+P351+P338

Ammonium hydrogen difluoride > RE - Pure

RE

Description ..... White flakes	Water (K.F.) ..... ≤ 0.5 %	Sulphite..... ≤ 100 ppm	Fe ..... ≤ 500 ppm
Identification ..... Positive	Sulphate..... ≤ 2000 ppm	Heavy metals (Pb)..... ≤ 200 ppm	Assay (acidimetric) ..... ≥ 94 %

Code	Size	Packaging	Notes
314261	1 kg	Plastic bottle	

di-Ammonium hydrogenphosphate ▶ Ammonium phosphate dibasic

Ammonium di-hydrogen phosphate 25 mg/L solution

• Ammonio diidrogeno fosfato 25 mg/l soluzione • Ammonium dihydrogénophosphate 25mg/l
• Amonio di-Hidrógeno Fosfato solución 25 mg/l • Ammoniumdihydrogenphosphat 25 mg/l

Synonym:
• *di-Ammonium hydrogenphosphate*
• *Monoammonium phosphate*

$NH_4H_2PO_4$
Molecular Weight: 115,03
CAS: 7722-76-1

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-
P362+P364-P337+P313

Ammonium di-hydrogen phosphate 25 mg/L solution > RS - Matrix modifiers for AAS-GTA

RS

Code	Size	Packaging	Notes
503194	50 ml	Plastic bottle	Matrix: 1% Nitric acid

Ammonium hydroxide solution ▶ Ammonia solution 30%

Product specifications are subject to changes. Please visit our website for updates.

Ammonium iodide

• Ammonio ioduro • Ammonium iodure • Amonio ioduro • Ammoniumiodid

NH_4I
Molecular Weight: 144,94
CAS: 12027-06-4
EEC-N: 234-717-7

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Ammonium iodide > RPE - For analysis - ACS

RPE

Description ..... White granules	Chloride + bromide (Cl) .....≤ 50 ppm	Heavy metals (Pb).....≤ 10 ppm	Assay (oxidimetric) ..... ≥ 99.0 %
Identification ..... Positive	Water-insoluble matter .....≤ 50 ppm	Ba.....≤ 20 ppm	Sulphate ..... ≤ 0.05 %
Residue on ignition..... ≤ 0.05 %	Phosphate .....≤ 10 ppm	Fe .....≤ 5 ppm	

Code	Size	Packaging	Notes
420133	50 g	Glass bottle	
420135	250 g	Glass bottle	

Stabilized with ~1,5% of NH4H2PO2

Ammonium iron(III) citrate green ▶ Iron (III) ammonium citrate green

Ammonium iron(III) citrate red ▶ Iron (III) ammonium citrate red

Ammonium iron(III) sulfate dodecahydrate ▶ Iron (III) ammonium sulfate dodecahydrate

Ammonium iron(II) sulfate hexahydrate ▶ Iron (II) ammonium sulfate hexahydrate

Ammonium molybdate tetrahydrate

• Ammonio molibdato tetraidrato • Ammonium molybdate tétrahydraté • Amonio molibdato tetrahidratado
• Ammoniummolybdat tetrahydrat

Synonym:
Ammonium heptamolybdate tetrahydrate

$(NH_4)_6Mo_7O_{24}.4H_2O$
Molecular Weight: 1235,86
CAS: 12054-85-2
EEC-N: 234-320-9

Warning
H302-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Ammonium molybdate tetrahydrate > RS - For microanalysis

RS

Description ..... White crystalline powder	AsO4, PO4, SiO4 (SiO2)..... ≤ 0.001 %	Heavy metals (Pb).....≤ 10 ppm	Assay (complexometric) .....81.0 ÷ 83.0 %
Identification ..... Positive	Chloride.....≤ 20 ppm	K.....≤ 100 ppm	
Water-insoluble matter ..... ≤ 0.005 %	Sulphate.....≤ 200 ppm	Mg .....≤ 50 ppm	
Nitrate ..... ≤ 0.003 %	Phosphate .....≤ 5 ppm	Na.....≤ 100 ppm	

Code	Size	Packaging	Notes
420391	100 g	Glass bottle	

Ammonium molybdate tetrahydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Green crystals	ppm	Sulphate.....≤ 200 ppm	Assay (oxidimetric) . 81.0 ÷ 83.0 % (MoO3)
Identification ..... Positive	Chloride.....≤ 20 ppm	Heavy metals (Pb).....≤ 10 ppm	Mg .....≤ 50 ppm
Water-insoluble matter .....≤ 50 ppm	Phosphate .....≤ 5 ppm	K.....≤ 100 ppm	
Arsenate, phosphate and silicate (SiO2)≤ 10	Nitrate .....≤ 30 ppm	Na.....≤ 100 ppm	

Code	Size	Packaging	Notes
420234	100 g	Glass bottle	
420236	500 g	Plastic bottle	
420238	2.5 kg	Plastic bottle	

Ammonium molybdate solution 2.5% in nitric acid

• Ammonio molibdato soluzione 2.5% in acido nitrico • Ammonium molybdate solution 2.5% dans l'acide nitrique
• Amonio molibdato solución 2.5% en acido nítrico • Ammonium molybdat Losung 2.5 %

$(NH_4)_6Mo_7O_{24}.4H_2O$
Molecular Weight: 1235,86
CAS: 12054-85-2

Classification transport
ONU: 2031
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Ammonium molybdate solution 2.5% in nitric acid > RPE - For analysis

RPE

Description ..... Clear liquid
Density at 20° C ..... 1.0 ÷ 1.2

Code	Size	Packaging	Notes
E420371	1 l	Bottle	

Ammonium nitrate

• Ammonio nitrato • Ammonium nitrate • Amonio nitrato • Ammonsalpeter

NH_4NO_3
Molecular Weight: 80,04
CAS: 6484-52-2
EEC-N: 229-347-8

Classification transport
ONU: 1942
Transport Hazard class: 5.1
Packing group III

Warning
H272-H319
P210-P220-P264-P280-P305+P351+P338-P337+P313

Ammonium nitrate > RPE - For analysis - ISO - ACS - Reag.Ph.Eur. - Reag.USP

RPE

Description ..... White crystals
Identification ..... Positive
pH sol. 5% at 25° C ..... 4.5 ÷ 6.0
Water-insoluble matter ..... ≤ 50 ppm
Residue on ignition ..... ≤ 100 ppm
Chloride ..... ≤ 5 ppm
Phosphate ..... ≤ 5 ppm
Nitrite ..... ≤ 5 ppm
Sulphate ..... ≤ 20 ppm
Heavy metals (Pb) ..... ≤ 5 ppm
Fe ..... ≤ 2 ppm
Assay (alkalimetric) ..... ≥ 95 %
Acidity ..... Conform

Code	Size	Packaging	Notes
420425	100 g	Plastic bottle	
420426	500 g	Plastic bottle	
420427	1 kg	Plastic bottle	
420429	5 kg	Plastic tank	
420422	25 kg	Plastic bucket	
420424	50 kg	Fibre drum	

Ammonium nitrate > RE - Pure

RE

Description ..... White pearls
Identification ..... Positive
pH sol. 5% at 25° C ..... 4.6 ÷ 7
Residue on ignition ..... ≤ 0.5 %
Total nitrogen ..... ≥ 34.2 %

Code	Size	Packaging	Notes
315509	5 kg	Plastic tank	
315502	25 kg	Plastic bucket	

Ammonium nitrate 200 mg/l solution

• Ammonio nitrato 200 mg/l soluzione • Ammonium nitrate 200 mg/l • Amonio nitrato 200 mg/l solución • Ammonsalpeter 200 mg/l

NH_4NO_3
Molecular Weight: 80,04
CAS: 6484-52-2

HEU210

Ammonium nitrate 200 mg/l solution > RS - Matrix modifiers for AAS-GTA

RS

Code	Size	Packaging	Notes
503195	50 ml	Plastic bottle	Matrix: Water

Product specifications are subject to changes. Please visit our website for updates.

Ammonium oxalate monohydrate

• Ammonio ossalato monoidrato • Ammonium oxalate monohydraté • Amonio oxalato monohidrato • Ammoniumoxalatmonohydrat

Synonym:
Oxalic acid diammonium salt

$(NH_4)_2C_2O_4.H_2O$
Molecular Weight: 142,11
CAS: 6009-70-7
EEC-N: 238-135-4

Warning
H302-H312
P264-P270-P280h-P301+P312a-P330-P501a

Ammonium oxalate monohydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description	White crystals	Residue on ignition	≤200 ppm	Heavy metals (Pb)	≤5 ppm
Identification	Positive	Chloride	≤20 ppm	Fe	≤2 ppm
Water-insoluble matter	≤50 ppm	Sulphate	≤20 ppm	Assay (oxidimetric)	99.0 ÷ 101.0 %

Code	Size	Packaging	Notes
420475	250 g	Plastic bottle	
420476	500 g	Plastic bottle	
420477	1 kg	Plastic bottle	
420478	2.5 kg	Plastic bottle	
420473	25 kg	Plastic bucket	

Ammonium oxalate solution 4%

• Ammonio ossalato soluzione 4% • Ammonium oxalate solution 4% • Amonio oxalato solución 4% • Ammoniumoxalat 4%

$(NH_4)_2C_2O_4.H_2O$
Molecular Weight: 142,11
CAS: 6009-70-7

HEU210

Ammonium oxalate solution 4% > RPE - For analysis

RPE

Description	Clear colourless liquid	Identification	Positive	Density at 20° C	1.00 ÷ 1.02	Assay (oxidimetric)	3.8 ÷ 4.2 %

Code	Size	Packaging	Notes
E420521	1 l	Bottle	

Ammonium persulfate

• Ammonio persolfato • Ammonium persulfate • Amonio persulfato • Ammoniumpersulfat

Synonym:
• *APS*
• *Ammonium peroxydisulfate*

$(NH_4)_2S_2O_8$
Molecular Weight: 228,2
CAS: 7727-54-0
EEC-N: 231-786-5

Classification transport
ONU: 1444
Transport Hazard class: 5.1
Packing group III

Danger
H272-H302-H315-H319-H334-H317-H335
P210-P280-P284-P304+P340-P305+P351+P338-P342+P311a-P403+P233

Ammonium persulfate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description	Yellowish crystals	Water-insoluble matter	≤50 ppm	Heavy metals (Pb)	≤50 ppm	Assay (oxidimetric)	≥98.0 %
Identification	Positive	Residue on ignition	≤500 ppm	Fe	≤10 ppm		
Acidity (H2SO4)	≤0.04 meq/g	Chloride & Chlorate(Cl)	≤10 ppm	Mn	≤0.5 ppm		

Code	Size	Packaging	Notes
420625	100 g	Plastic bottle	
420626	500 g	Plastic bottle	
420627	1 kg	Plastic bottle	
420629	5 kg	Plastic jar	
420623	25 kg	Plastic bucket	

Ammonium persulfate > RE - Pure

RE

Description	White crystals or yellowish	Heavy metals (Pb)	≤ 50 ppm	Assay (oxidimetric)	≥ 97.5 %
Identification	Positive	Fe	≤ 10 ppm		

Code	Size	Packaging	Notes
316008	500 g	Plastic bottle	
316002	25 kg	Plastic bucket	

Ammonium phosphate dibasic

• Ammonio fosfato bibasico • Ammonium phosphate dibasique • Amonio fosfato dibásico
• Ammoniumphosphat dibasisch

Synonym:
• *Ammonium hydrogen phosphate*
• *Diammonium hydrogenphosphate*

$(NH_4)_2HPO_4$
Molecular Weight: 132,06
CAS: 7783-28-0
EEC-N: 231-987-8

Warning
H312-H332
P261-P271-P280h-P304+P340-P312a-P501a

Ammonium phosphate dibasic > RPE - For analysis - ACS - Reag. USP

RPE

Description ..... White crystals	Sulphate ..... ≤ 100 ppm	Fe ..... ≤ 10 ppm	Ca ..... ≤ 10 ppm
Identification ..... Positive	Chloride ..... ≤ 10 ppm	K ..... ≤ 50 ppm	Mg ..... ≤ 5 ppm
pH sol. 5% at 25° C ..... 7.7 ÷ 8.1	Nitrate ..... ≤ 30 ppm	Na ..... ≤ 50 ppm	
Water-insoluble matter ..... ≤ 50 ppm	Heavy metals (Pb) ..... ≤ 10 ppm	Assay (alkalimetric) ..... ≥ 98.0 %	

Code	Size	Packaging	Notes
419835	100 g	Plastic bottle	
419836	500 g	Plastic bottle	
419837	1 kg	Plastic bottle	
419831	5 kg	Plastic tank	
419832	25 kg	Plastic bucket	
419834	50 kg	Fibre drum	

Ammonium phosphate dibasic > RE - Pure

RE

Description ..... white crystals	pH sol. 1M ..... 7.6 - 8.2	Assay (P2O5) ..... ≥ 53.4 %	Assay (alcalimetric) ..... 98 - 102 %
Identification ..... Positive	Insoluble in water ..... ≤ 0.2 %	Assay (nitrogen) ..... ≥ 21 %	

Code	Size	Packaging	Notes
314757	1 kg	Plastic bottle	
314758	2.5 kg	Plastic bottle	

Ammonium phosphate monobasic

• Ammonio fosfato monobasico • Ammonium phosphate monobasique • Amonio fosfato monobásico
• Ammoniumphosphat einbasig

Synonym:
• *Ammonium dihydrogenphosphate*
• *Monoammonium phosphate*

$NH_4H_2PO_4$
Molecular Weight: 115,03
CAS: 7722-76-1
EEC-N: 231-764-5

Ammonium phosphate monobasic > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White crystals	Water-insoluble matter ..... ≤ 50 ppm	Ca ..... ≤ 10 ppm	Na ..... ≤ 50 ppm
Identification ..... Positive	Heavy metals (Pb) ..... ≤ 5 ppm	Mg ..... ≤ 5 ppm	Assay (acidimetric) ..... ≥ 98.0 %
pH sol. 5% at 25° C ..... 3.8 ÷ 4.4	Nitrate ..... ≤ 0.001 %	Fe ..... ≤ 10 ppm	
Chloride ..... ≤ 5 ppm	Sulphate ..... ≤ 100 ppm	K ..... ≤ 50 ppm	

Code	Size	Packaging	Notes
419785	100 g	Plastic bottle	
419786	500 g	Plastic bottle	
419787	1 kg	Plastic bottle	

Ammonium phosphate monobasic > RE - Pure

RE

Appearance ..... White crystalline	(w/w)	Fluoride ..... ≤ 10 ppm	Hg ..... ≤ 1 ppm
P2O5 ..... ≥ 61.3 % (w/w)	pH (1% solution) ..... 4.3 - 4.7	As ..... ≤ 1 ppm	Cd ..... ≤ 1 ppm
Assay (as(NH4)H2PO4) ..... 98.0 - 102.2 %	Water insoluble Matter ..... ≤ 0.2 % (w/w)	Pb ..... ≤ 1 ppm	

Code	Size	Packaging	Notes
314507	1 kg	Plastic bottle	
314505	5 kg	Plastic tank	
314506	10 kg	Plastic tank	
314504	25 kg	Sack	

Product specifications are subject to changes. Please visit our website for updates.

a b c d e f g h i j k l m n o p q r s t u v w x y z

Ammonium rhodanide ▶ Ammonium thiocyanate

Ammonium sodium hydrogen phosphate ▶ Sodium ammonium hydrogen phosphate

Ammonium sulfamate

• Ammonio solfammato • Ammonium sulfamate • Amonio sulfamato • Ammoniumsulfamat

$NH_4SO_3NH_2$
Molecular Weight: 114,12
CAS: 7773-06-0
EEC-N: 231-871-7

Warning
H302
P264-P270-P301+P312a-P330-P501a

Ammonium sulfamate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White crystals	Melting point..... 131.0 ÷ 135.0 °C	Water-insoluble matter ..... ≤200 ppm	Assay (oxidimetric) ..... ≥98.0 %
Identification ..... Positive	Residue on ignition..... ≤0.1 %	Heavy metals (Pb)..... ≤5 ppm	

Code	Size	Packaging	Notes
420724	100 g	Glass bottle	
420725	500 g	Glass bottle	

Ammonium sulfate

• Ammonio solfato • Ammonium sulfate • Amonio sulfato • Ammoniumsulfat

Synonym:
Mascagnite

$(NH_4)_2SO_4$
Molecular Weight: 132,14
CAS: 7783-20-2
EEC-N: 231-984-1

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Ammonium sulfate > RPE - For analysis - ISO - Reag.Ph.Eur.

RPE

Description ..... White crystals	Water-insoluble matter ..... ≤50 ppm	Phosphate ..... ≤5 ppm	Fe ..... ≤5 ppm
Identification ..... Positive	Residue on ignition..... ≤100 ppm	Nitrate ..... ≤10 ppm	Assay (acc.to Sörensen)..... ≥99.0 %
pH sol. 5% at 25° C ..... 5.0 ÷ 6.0	Chloride..... ≤5 ppm	Heavy metals (Pb)..... ≤5 ppm	

Code	Size	Packaging	Notes
420775	100 g	Plastic bottle	
420776	500 g	Plastic bottle	
420777	1 kg	Plastic bottle	
420772	5 kg	Plastic tank	
420771	25 kg	Plastic bucket	
420774	50 kg	Fibre drum	

Ammonium sulfate > RE - Pure

RE

Description ..... White crystals	Chloride..... ≤ 3 ppm	Ca..... ≤ 10 ppm	Residue on calcination..... ≤ 0.01 % (SO4)
Identification ..... Positive	Nitrate ..... ≤ 10 ppm	Fe ..... ≤ 5 ppm	Assay ..... ≥ 99.0 %
pH sol. 5% at 20°C ..... 5 ÷ 6	Phosphate ..... ≤ 5 ppm	Heavy metals (Pb)..... ≤ 5 ppm	

Code	Size	Packaging	Notes
316257	1 kg	Plastic bottle	
316251	5 kg	Plastic tank	
316252	25 kg	Plastic bucket	

Ammonium sulfide solution 20%

• Ammonio solfuro soluzione 20% • Ammonium sulfure en solution à 20% • Amonio sulfuro solución 20% • Ammoniumsulfid in 20% iger Lösung

$(NH_4)_2S$
Molecular Weight: 68,141
CAS: 12135-76-1
EEC-N: 235-223-4

Classification transport
ONU: 2683
Transport Hazard class: 8
Packing group II



Danger
H225-H301-H314-HEU031
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338

Ammonium sulfide solution 20% > RPE - For analysis

RPE

Description	Yellow clear liquid	Residue on ignition	≤50 ppm	Fe	≤5 ppm	Sb	≤5 ppm
Identification	Positive	Sulphate	≤100 ppm	K	≤10 ppm	Sn	≤5 ppm
Density at 15° C	0.960 ÷ 1.000	As	≤1 ppm	Mg	≤5 ppm	Zn	≤5 ppm
Carbonate	≤50 ppm	Ca	≤30 ppm	Na	≤50 ppm	Assay (ex ammonium)	≥20 %
Chloride	≤50 ppm	Cd	≤5 ppm	Ni	≤5 ppm	Assay (argentimetric)	≥20 %
Heavy metals (Pb)	≤10 ppm	Cu	≤5 ppm	Pb	≤5 ppm		

Code	Size	Packaging	Notes
421101	1 l	Glass bottle	

Ammonium L(+)tartrate

• Ammonio L(+)tartrato • Ammonium L(+)tartrate • Amonio L(+)tartrato • Ammoniumtartrat

Synonym:
• *L-(+)-Tartaric acid diammonium salt*
• *Diammonium tartrate*

$(CHOHCOONH_4)_2$
Molecular Weight: 184,15
CAS: 3164-29-2
EEC-N: 221-618-9

Ammonium L(+)tartrate > RPE - For analysis

RPE

Description	White crystalline powder	Sulfated ashes	≤ 0.02 %	Water insoluble Matter	≤ 0.005 %	Heavy metals (as Pb)	≤ 0.0005 %
Identification	Positive	Phosphate	≤ 0.001 %	Chloride	≤ 0.001 %		
Assay (acidimetric)	≥ 99.0 %	Fe	≤ 0.0005 %	Sulfate	≤ 0.005 %		

Code	Size	Packaging	Notes
421206	500 g	Plastic bottle	

Ammonium thiocyanate

• Ammonio solfocianuro • Ammonium thiocyanate • Amonio sulfocianuro • Ammoniumthiocyanat

Synonym:
Ammonium rhodanide

NH_4SCN
Molecular Weight: 76,12
CAS: 1762-95-4
EEC-N: 217-175-6



Warning
H302-H312-H332-HEU032
P261-P264-P271-P280h-P301+P312a-P304+P340

Ammonium thiocyanate > RPE - For analysis - ISO - ACS - Reag.Ph.Eur. - Reag.USP

RPE

Description	White crystals	Water-insoluble matter	≤50 ppm	Chloride	≤50 ppm	Fe	≤3 ppm
Identification	Positive	Residue on ignition	≤250 ppm	Sulphate	≤50 ppm	Assay (argentimetric)	≥97.5 %
pH sol. 5% at 25° C	4.5 ÷ 6.0	Reducing iodine	≤0.004 meq/g	Heavy metals (Pb)	≤5 ppm		

Code	Size	Packaging	Notes
420885	250 g	Plastic bottle	
420886	500 g	Plastic bottle	
420887	1 kg	Plastic bottle	

Ammonium thiocyanate > RE - Pure

RE

Description	White crystals	Sulphated ash	≤ 0.03 %	Sulphate	≤ 50 ppm	Assay (argentimetric)	≥ 99 %
Identification	Positive	S	≤ 20 ppm	Fe	≤ 2 ppm		

Code	Size	Packaging	Notes
316307	1 kg	Plastic bottle	
316308	5 kg	Plastic tank	
316303	25 kg	Plastic bucket	

Product specifications are subject to changes. Please visit our website for updates.

Ammonium thiocyanate 1 mol/l (1N)

• Ammonio solfocianuro 1 mol/l (1N) • Ammonium thiocyanate 1 mol/l (1N) • Amonio sulfocianuro 1 mol/l (1N) • Ammoniumthiocyanat 1 mol/l (1N)

NH_4SCN
CAS: 1762-95-4

Ammonium thiocyanate 1 mol/l (1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.998 - 1.002 N

Code	Size	Packaging	Notes
420946	500 ml	Plastic bottle	

76,12 g of NH4SCN. Volumetric solution ready-to-use: 1 N. Stabilized with p-oxybenzoate

Ammonium thiocyanate 0.1 mol/l (0.1N)

• Ammonio solfocianuro 0.1 mol/l (0.1N) • Ammonium thiocyanate 0.1 mol/l (0.1N) • Amonio sulfocianuro 0.1 mol/l (0.1N)
• Ammoniumthiocyanat 0.1 mol/l (0.1N)

NH_4SCN
Molecular Weight: 76,12
CAS: 1762-95-4

Ammonium thiocyanate 0.1 mol/l (0.1N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613000501	500 ml	Plastic bottle	Ref Ph.Eur 3000500
613000500	1 l	Plastic bottle	Ref Ph.Eur 3000500

Ammonium thiocyanate 0.1 mol/l (0.1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.0998 - 0.1002 N

Code	Size	Packaging	Notes
420977	1 l	Plastic bottle	

7,612 g of NH4SCN. Volumetric solution ready-to-use: 0,1 N. Stabilized with p-oxybenzoate

Ammonium thiocyanate 0.1 mol/l (0.1N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
421001		Plastic ampoule	Volume: 55 ml

7,612 g of NH4SCN. Volumetric concentrated solution to prepare 1 L of solution 0,1 N

Ammonium thiocyanate 0.01 mol/l (0.01N)

• Ammonio solfocianuro 0.01 mol/l (0.01N) • Ammonium thiocyanate 0.01 mol/l (0.01N) • Amonio sulfocianuro 0.01 mol/l (0.01N)
• Ammoniumthiocyanat 0.01 mol/l (0.01N)

NH_4SCN
Molecular Weight: 76,12
CAS: 1762-95-4

Ammonium thiocyanate 0.01 mol/l (0.01N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
421061		Plastic ampoule	Volume: 55 ml

0,7612 g of NH4SCN. Volumetric concentrated solution to prepare 1 L of solution 0,01 N

n-Amyl alcohol

• Alcole n-amilico • Alcool n-amylique • Alcohol n-amílico • n-Amylalkohol

Synonym:
• *1-Pentanol*
• *Pentyl alcohol*

$CH_3(CH_2)_3CH_2OH$
Molecular Weight: 88,15
CAS: 71-41-0
EEC-N: 200-752-1

Classification transport
ONU: 1105
Transport Hazard class: 3
Packing group III

Warning
H226-H332-H315-H335
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

n-Amyl alcohol > RPE - For analysis

RPE

Description	Clear colourless liquid	Boiling point	136 - 138 °C	Cu	≤0.05 ppm	Na	≤1 ppm
Identification	Positive	Water (K.F.)	≤0.2 %	Fe	≤0.5 ppm	Pb	≤0.1 ppm
Density at 25°C	~ 0.811	Residue on evaporation	≤10 ppm	K	≤0.5 ppm	Zn	≤0.5 ppm
Refractive index at 20°C	1.4081 ÷ 1.4121	Ca	≤0.5 ppm	Mg	≤0.5 ppm	Assay (GLC)	≥99 %

Code	Size	Packaging	Notes
413783	1 l	Glass bottle	

n-Amyl alcohol > RE - Pure

RE

Description	Clear colourless liquid	Density at 20° C	0.815 ÷ 0.819	Water (K.F.)	≤0.2 %
Identification	Positive	Refractive index at 20°C	1.4061 ÷ 1.4141	Assay (GLC)	≥98.5 %

Code	Size	Packaging	Notes
307901	1 l	Glass bottle	

tert-Amyl alcohol

• Alcole ter-amilico • Alcool tert-amylique • Alcohol ter-amílico • 2-Methylbutanol-2

Synonym:
• *Tert-pentyl alcohol*
• *2-Methyl-2-butanol*

$(CH_3)_2C(OH)CH_2CH_3$
Molecular Weight: 88.151
CAS: 75-85-4
EEC-N: 200-908-9

Classification transport
ONU: 1105
Transport Hazard class: 3
Packing group II

Danger
H225-H332-H315-H335
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

tert-Amyl alcohol > RPE - For analysis

RPE

Description	Clear colourless liquid	Boiling point	101.3 ÷ 102.8 °C	Alcalinity (NH3)	≤5 ppm	Subst. reducing KMnO4	≤200 ppm
Identification	Positive	Residue on evaporation	≤10 ppm	Indole base	≤0.1 ppm	Cu	≤0.5 ppm
Ready carbonizable substances	Conform	Acids and esthers	≤500 ppm	Organic base (N)	≤7 ppm	Fe	≤0.5 ppm
Density at 20° C	0.808 ÷ 0.810	Acidity(valerianic.ac)	≤50 ppm	Carbonyl Compounds (CO)	≤80 ppm	Assay (GLC)	≥99.5 %
Refractive index at 20°C	1.4027 ÷ 1.4077	Water (K.F.)	≤0.2 %	Pyridine and homologues	≤30 ppm		

Code	Size	Packaging	Notes
413941	250 ml	Glass bottle	
413944	1 l	Glass bottle	
413945	25 l	Plastic tank	

Aniline blue soluble in water

• Blu anilina solubile in acqua • Bleu d'aniline soluble dans l'eau • Azul de anilina soluble en agua • Blaues Anilin. wasserlöslich

$C_{32}H_{25}N_3Na_2O_9S_3$
Molecular Weight: 737,74
CAS: 28631-66-5
EEC-N: 249-113-9

Warning
H302-H312-H332
P261-P264-P271-P280h-P301+P312a-P304+P340

Aniline blue soluble in water > RPE - For analysis - C.I. 42755

RPE

Description	Violet powder	Identification	Positive

Code	Size	Packaging	Notes
428582	25 g	Glass bottle	

Dye for microscopy (botanical-cytology-histology). Indicator acid - base (pH 9.4 ÷ 14.0)

Product specifications are subject to changes. Please visit our website for updates.

Aniline hydrochloride

• Anilina cloridrato • Aniline chlorhydrate • Anilina clorhidrato • Aniliniumchlorid

$C_6H_5NH_2.HCl$
Molecular Weight: 129,59
CAS: 142-04-1
EEC-N: 205-519-8

Classification transport
ONU: 1548
Transport Hazard class: 6.1
Packing group III



Danger
H301-H311-H331-H318-H317-H341-H351-H372-H400
P280-P304+P340-P305+P351+P338-P308+P313-P330-P361+P364-P403+P233

Aniline hydrochloride > RPE - For analysis

RPE

Description ..... White crystalline powder | Identification ..... Positive | Melting point..... 196 ÷ 199 ° C | Assay (acidimetric) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
422376	500 g	Plastic bottle	

Anisaldehyde

• Aldeide anisica • Aldéhyde anisique • Aldehído anísico • Anisaldehyd

Synonym:
• *p-Anisaldehyde*
• *4-Methoxybenzaldehyde*

4-$CH_3OC_6H_4CHO$
Molecular Weight: 136,15
CAS: 123-11-5
EEC-N: 204-602-6



Warning
H302
P264-P270-P301+P312a-P330-P501a

Anisaldehyde > RE - Pure

RE

Description ..... Clear yellow liquid | Identification ..... Positive | Density at 20° C ..... 1.121 ÷ 1.125 | Refractive index at 20°C..... 1.5710 ÷ 1.5750 | Assay (GLC) .....≥99 %

Code	Size	Packaging	Notes
415312	100 ml	Glass bottle	

Anisaldehyde solution

• Aldeide anisica soluzione • Aldéhyde anisique solution • Aldehído anísico solución • Anisaldehydlösung

Classification transport
ONU: 1993

Anisaldehyde solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611007301	100 ml	Glass bottle	Ref Ph.Eur 1007301
611007302	100 ml	Glass bottle	Anisaldehyde solution R1 Ref Ph.Eur 1007302

Anisic acid

• Acido anisico • Acide anisique • Acido anísico • Anissäure

Synonym:
• *4-Methoxybenzoic acid*
• *Draconic acid*

$CH_3OC_6H_4COOH$
Molecular Weight: 152,15
CAS: 100-09-4
EEC-N: 202-818-5



Warning
H302-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Anisic acid > RPE - For analysis

RPE

Description .... White to slightly grey - beige powder | Identification ..... Positive | Melting point..... 181 ÷ 186 °C | Assay (GLC) ..... ≥ 97.5 %

Code	Size	Packaging	Notes
402133	100 g	Glass bottle	

Anthrone

• Antrone • Anthrone • Antrona • Anthron

Synonym:
9(10H)-Anthracenone

$C_6H_4COC_6H_4CH_2$
Molecular Weight: 194,23
CAS: 90-44-8
EEC-N: 201-994-0

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Anthrone > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Yellow crystals
Identification ..... Positive
Melting point ..... 151 - 159 °C
Absorbance of reag.sol. ..... Conform
Solub.ty in Dieth.Ether ..... Conform
Sens. to carbohydrates ..... Conform

Code	Size	Packaging	Notes
423281	10 g	Glass bottle	
423282	25 g	Glass bottle	

For the determination of carbohydrates

Antimony standard solution

• Antimonio standard soluzione • Antimoine standard solution • Antimonio, solución patrón • Antimonstandardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group II

Danger
H290-H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Antimony standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615000400	100 ml	Plastic bottle	A 100 ppm solution Ref Ph.Eur 5000400

Antimony standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505832	100 ml	Plastic bottle	conc. 10 ppm Matrix: Nitric acid and hydrofluoric acid
505835	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid and hydrofluoric acid
505833	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Antimony standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503891	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503893	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503895	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid
503897	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid
503899	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid
503898	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrofluoric acid and nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Antimony standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
507525	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid
507479	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid
E497415	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
E497411	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Antimony standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
422731		Plastic ampoule	conc. 1.000 ppm Matrix: Hydrochloric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Antimony potassium tartrate

• Antimonio potassio tartrato • Antimoine potassium tartrate • Antimonio y potasio tartrato
• Antimon-Kalium-Tartrat

Synonym:
Potassium antimonyl tartrate

$C_4H_2KO_6Sb_6.1,5H_2O$
Molecular Weight: 333,93
CAS: 28300-74-5
EEC-N: 234-293-3

Classification transport
ONU: 1551
Transport Hazard class: 6.1
Packing group III



Warning
H302-H332-H411
P261-P264-P271-P301+P312a-P304+P340-P501a

Antimony potassium tartrate > RPE - For analysis - Reag. Ph. Eur.

RPE

Description ..... White crystalline powder
Identification ..... Positive
pH sol. 2% at 20°C ..... 3.8 - 4.5
Chloride ..... ≤ 0.005 %
Fe ..... ≤ 0.001 %
Water insoluble substances ..... ≤ 0.05 %
Loss on drying (105°C) ..... ≤ 2.7 %
As ..... ≤ 0.015 %
Assay (iodometric) ..... ≥ 99.0 %

Code	Size	Packaging	Notes
423035	250 g	Plastic bottle	
423036	500 g	Plastic bottle	
423037	1 kg	Plastic bottle	

Antimony trichloride

• Antimonio tricloruro • Antimoine trichlorure • Antimonio tricloruro • Antimontrichlorid

Synonym:
Antimony(III) chloride

$SbCl_3$
Molecular Weight: 228,11
CAS: 10025-91-9

Classification transport
ONU: 2922
Transport Hazard class: 8
Packing group II



Danger
H302-H332-H314-H351-H361d-H335-H372-H412-HEU301
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Antimony trichloride > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611007701	100 ml	Glass bottle	Ref Ph.Eur 1007701

Antimony trichloride > RPE - For analysis - ACS

RPE

Description ..... White crystals
Identification ..... Positive
Chloroform insoluble ..... ≤ 500 ppm
Sulphate ..... ≤ 50 ppm
As ..... ≤ 200 ppm
Ca ..... ≤ 50 ppm
Cu ..... ≤ 10 ppm
Fe ..... ≤ 20 ppm
K ..... ≤ 100 ppm
Na ..... ≤ 200 ppm
Pb ..... ≤ 50 ppm
Assay (iodometric) ..... ≥ 99.0 %

Code	Size	Packaging	Notes
422834	100 g	Glass bottle	
422835	250 g	Glass bottle	

Aqueous calcium hydroxide

• Acqua di calce • Eau de chaux • Calcio hidróxido acuoso • Kalkwasser

Aqueous calcium hydroxide > RPE - For analysis

RPE

Description ..... Clear liquid Identification ..... Positive Assay ..... 0.140 ÷ 0.169 % (p/p)

Code	Size	Packaging	Notes
E411921	1 l	Plastic bottle	

L(+)Arginine

• L(+)Arginina • L(+)Arginine • L(+)Arginina • L(+)Arginin

Synonym:
(S)-2-Amino-5-guanidinopentanoic acid

$C_6H_{14}N_4O_2$
Molecular Weight: 174,2
CAS: 74-79-3
EEC-N: 200-811-1

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

L(+)Arginine > RPE - For analysis

RPE

Description ..... White crystalline powder Melting point ..... ≥230 ° C Assay (non-aqueos medium) ..... ≥98 %
Identification ..... Positive Specific optical rotation ..... +26.3 ÷ +27.7 °

Code	Size	Packaging	Notes
424271	100 g	Glass bottle	

L(+)Arginine monohydrochloride

• L(+)Arginina monocloridrato • L(+)Arginine monochlorhydratée • L(+)Arginina monoclorhidrato
• L-Argininmonohydrochlorid

Synonym:
(S)-(+)-2-Amino-5-[(aminoiminomethyl)amino] pentanoic acid monohydrochloride

$C_6H_{14}N_4O_2 \cdot HCl$
Molecular Weight: 210,66
CAS: 1119-34-2
EEC-N: 214-275-1

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

L(+)Arginine monohydrochloride > RPE - For analysis

RPE

Description ..... White powder Loss on drying ..... ≤ 0.2 % Heavy metals (Pb) ..... ≤ 10 ppm Fe ..... ≤ 10 ppm
Identification ..... Positive Ammonium ..... ≤ 200 ppm Residue on ignition ..... ≤ 0.1 % Assay (non-aqueos medium) 98.5 ÷ 101.0 % s.s.
Potere rotat. spec. (c=8; HCl 6N) .. +21.4 ÷ +23.5 ° s.s. Water-insoluble matter ..... ≤ 100 ppm Sulphate ..... ≤ 300 ppm

Code	Size	Packaging	Notes
424268	5 g	Glass bottle	

Arsenazo III

• Arsenazo III • Arsenazo III • Arsenazo III • Arsenazo III

Synonym:
2,7-Bis(2-arsonophenylazo)chromotropic acid

$C_{22}H_{16}As_2N_4O_{14}S_2Na_2$
Molecular Weight: 774,36
CAS: 1668-00-4
EEC-N: 216-788-6

Classification transport
ONU: 3465
Transport Hazard class: 6.1
Packing group II

Danger
H301-H331-H410
P261-P271-P304+P340-P311a-P330-P403+P233

Arsenazo III > RPE - For analysis

RPE

Description ..... Brown granular powder Identification ..... Positive Sens.as complex.indicat ..... Conform

Code	Size	Packaging	Notes
424281	1 g	Glass bottle	
424282	25 g	Glass bottle	

Suitable for determination of Th, Zr, U, Cd, Zn, Ca

Product specifications are subject to changes. Please visit our website for updates.

Arsenic standard solution

• Arsenico standard soluzione • Arsenic standard solution • Arsenic, solución patrón • Arsen-Standardlösung

Danger
H350-HA26
P201-P202-P280-P308+P313-P405-P501a

Arsenic standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615000501	100 ml	Plastic bottle	A 1 ppm solution: to dilute according to Ref Ph.Eur 5000501
615000509	100 ml	Plastic bottle	A 10 ppm solution: to dilute according to Ref Ph.Eur 5000500

Arsenic standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505312	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505315	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505313	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Arsenic standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503421	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503423	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503425	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503427	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Arsenic standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
504439	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507496	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Arsenic trioxide solution

• Arsenico triossido soluzione • Arsenic trioxyde solution • Arsénico trióxido solución • Arsentrioxid-Lösung

Synonym:
• *Arsenic(III) oxide*
• *Arsenic trioxide*

As_2O_3
Molecular Weight: 197,84
CAS: 1327-53-3

Arsenic trioxide solution > RS - For analysis according to JP

RS

Code	Size	Packaging	Notes
616000001	100 ml	Plastic bottle	Arsenic (III) stock solution

Arsenic trioxide solution > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
617000001	100 ml	Plastic bottle	Arsenic trioxyde stock solution

L(+)Ascorbic acid

• Acido L(+)ascorbico • Acide L(+)ascorbique • Acido L(+)ascórbico • L(+)-Ascorbinsäure

Synonym:
• *Vitamin C*
• *L-Threoascorbic acid*

$COCOH:COHCHCHOHCH_2OH$
Molecular Weight: 176,13
CAS: 50-81-7
EEC-N: 200-066-2

L(+)Ascorbic acid > RPE - For analysis - ISO

RPE

Description .....White crystal. powder
Identification ..... Positive
Melting point.....190.5 ÷ 192.0 °C
Specific optical rotation... +20.5 ÷ +21.5 °
Loss on drying ..... ≤ 0.1 %
Chloride.....≤ 50 ppm
Water-insoluble matter .....≤ 30 ppm
Heavy metals (Pb).....≤10 ppm
Residue on ignition.....≤300 ppm
Sulphate.....≤ 20 ppm
Cu.....≤0.3 ppm
Fe .....≤2 ppm
Pb.....≤0.5 ppm
Assay (oxidimetric) ..... ≥99.0 %

Code	Size	Packaging	Notes
402404	100 g	Plastic bottle	
402406	500 g	Plastic bottle	
402407	1 kg	Plastic bottle	

L(+)Asparagine

• L(+)Asparagina • L(+)Asparagine • L(+)Asparagina • L(+)Asparagin

Synonym:
• *(S)-(+)-2-Aminosuccinamic acid*
• *L-Aspartic acid 4-amide*

$C_4H_8N_2O_3.H_2O$
Molecular Weight: 150,14
CAS: 5794-13-8

L(+)Asparagine > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification (I.R.)..... Positive
Specific optical rotation... +33.5 ÷ +36.5 °
Loss on drying .....11.5 ÷ 12.5 %
Residue on ignition..... ≤ 0.1 %
Heavy metals (Pb).....≤ 10 ppm
As .....≤ 1 ppm
Fe .....≤ 10 ppm
Ammonium ..... ≤ 0.1 %
Chloride..... ≤ 0.02 %
Sulphate..... ≤ 0.03 %
Assay (non-aqueos medium)98.5÷101.0 %
(s.s.)

Code	Size	Packaging	Notes
424544	100 g	Glass bottle	
424547	1 kg	Plastic bottle	

L(+)Aspartic acid

• Acido L(+)aspartico • Acide L(+)aspartique • Acido L(+)aspártico • L(+)Asparaginsäure

Synonym:
• *(S)-(+)-Aminosuccinic acid*
• *(S)-Aminobutanedioic acid*

$HOOCCHNH_2CH_2COOH$
Molecular Weight: 133,1
CAS: 56-84-8
EEC-N: 200-291-6

L(+)Aspartic acid > RE - Pure

RE

Description ..... White crystalline powder
Identification ..... Positive
Potere rotat. spec. (c=8; HCl 6N)..... +24 ÷ +26 °
Loss on drying ..... ≤ 0.2 %
Ammonium.....≤ 500 ppm
Chloride.....≤ 200 ppm
Sulphate.....≤ 300 ppm
As .....≤ 1 ppm
Fe .....≤ 10 ppm
Heavy metals (Pb).....≤ 10 ppm
Residue on ignition..... ≤ 0.1 %
Assay (non-aqueos medium)≥ 98.0 % (s.s.)

Code	Size	Packaging	Notes
402442	25 g	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

ASTM colour standards

• Standard del colore ASTM • Etalons couleurs ASTM • Patrones de color ASTM • ASTM-Farbstandards



Danger
H315-H304
P264-P280g-P301+P310a-P331-P362+P364-P332+P313

ASTM colour standards > RS - For calibration

RS

Code	Size	Packaging	Notes
540601	100 ml	Glass bottle	Sample A1
540602	100 ml	Glass bottle	Sample A3
540603	100 ml	Glass bottle	Sample A5
540604	100 ml	Glass bottle	Sample A7

Ausilab 110

• Ausilab 110 • Ausilab 110 • Ausilab 110 • Ausilab 110



Warning
H319-HEU208
P264-P280i-P305+P351+P338-P337+P313

Ausilab 110 > RE - Pure - For glassware manual washing

RE

Code	Size	Packaging	Notes
FG201149	4 x 5 kg	Plastic tank	

Multi-purpose detergent concentrate liquid, neutral. Contains: Anionic surfactants (15-30%), less than 5% emollients. Use: diluted in water, very foaming

Ausilab 140

• Ausilab 140 • Ausilab 140 • Ausilab 140 • Ausilab 140



Danger
H315-H318
P264-P280a-P305+P351+P338-P310a-P362+P364-P332+P313

Ausilab 140 > RE - Pure - For glassware manual washing

RE

Code	Size	Packaging	Notes
FG201220	4 x 5 kg	Plastic tank	

Degreaser concentrate, liquid, alkaline. Contains: EDTA and salts (5-15%), non-ionic surfactants (5-15%). Use: diluted in water

Ausilab 210

• Ausilab 210 • Ausilab 210 • Ausilab 210 • Ausilab 210



Danger
H315-H318
P264-P280a-P305+P351+P338-P310a-P362+P364-P332+P313

Ausilab 210 > RE - Pure - For glassware washers

RE

Code	Size	Packaging	Notes
FG201138	8 kg	Plastic bucket	

Detergent for glassware washers, powder, alkaline (suitable for use in the food industry). Contains: over than 30% of phosphates, less than 5% non-ionic surfactants. Use: concentration as dirt and water hardness

Ausilab 250
• Ausilab 250 • Ausilab 250 • Ausilab 250 • Ausilab 250



Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Ausilab 250 > RE - Pure - For glassware washers — RE

Code	Size	Packaging	Notes
FG201139	4 x 5 kg	Plastic tank	

Neutralizing detergent for glassware washers, liquid, acid. Contains: over than 30% of organics acids. Use: for washers with automatic metering

Ausilab 260
• Ausilab 260 • Ausilab 260 • Ausilab 260 • Ausilab 260

NaCl
Molecular Weight: 58,44
CAS: 7647-14-5

Ausilab 260 > RE - Pure - For glassware washers — RE

Code	Size	Packaging	Notes
FG20C1156C5	5 kg	Plastic bucket	

Regenerating salt for dishwashers

Ausilab 280
• Ausilab 280 • Ausilab 280 • Ausilab 280 • Ausilab 280

Classification transport
ONU: 1719
Transport Hazard class: 8
Packing group III



Danger
H314
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Ausilab 280 > RE - Pure - For glassware washers — RE

Code	Size	Packaging	Notes
FG201151	4 x 5 kg	Plastic tank	

Dishwashing detergent, liquid, alkaline (suitable for use in the food industry). Contains: over than 30% E.D.T.A., NaOH (5-15%). Use: For dishwashers with automatic metering

Ausilab 290
• Ausilab 290 • Ausilab 290 • Ausilab 290 • Ausilab 290

Classification transport
ONU: 1719
Transport Hazard class: 8
Packing group III



Danger
H314
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Ausilab 290 > RE - Pure - For glassware washers — RE

Code	Size	Packaging	Notes
FG201200	4 x 6 kg	Plastic tank	

Dishwashing detergent, liquid, alkaline, manufactured according to Ecocert standard. Contains: over than 30% Citrates and gluconates, NaOH (5-15%). Use: For dishwashers with automatic metering

Ausilab 400
• Ausilab 400 • Ausilab 400 • Ausilab 400 • Ausilab 400

Classification transport
ONU: 1987
Transport Hazard class: 3
Packing group III



Warning
H226-H319-H412
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

Ausilab 400 > RE - Pure - For hand washing cream — RE

Code	Size	Packaging	Notes
FG201146	6 x 800 ml	Box	

Alcoholic sanitizing gel, neutral. Contains: over 50% than alcohols, moisturizing agents (5-15%), benzalkonium chloride (0-1%). Use: Pure

Ausilab 500
• Ausilab 500 • Ausilab 500 • Ausilab 500 • Ausilab 500



Danger
H315-H318-HEU208
P264-P280a-P305+P351+P338-P310a-P362+P364-P332+P313

Ausilab 500 > RE - Pure - For washing — RE

Code	Size	Packaging	Notes
FG201140	15 x 750 ml	Box	

Purpose cleaner. Contains between 5 and 15% of alcohols, less than 5% non-ionic surfactants. Use: Direct spray of pure product

Azomethine H
• Azometina H • Azométhine H • Azometino H • Azomethin H

Synonym:
4-Hydroxy-5-(2-hydroxybenzylideneamino)-naphthalene-2,7-disulfonic acid monosodium salt hydrate

$C_{17}H_{12}NNaO_8S_2 \cdot xH_2O$
Molecular Weight: 445,4
CAS: 206752-32-1



Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Azomethine H > RPE - For analysis — RPE

Description ..... Yellow orange powder Identification ..... Positive

Code	Size	Packaging	Notes
424691	10 g	Glass bottle	
424692	25 g	Glass bottle	

Suitable for the determination of boron

Azure II
• Azzurro II • Azur II • Azur II • Azur II

CAS: 37247-10-2

Classification transport
ONU: 2811
Transport Hazard class: 6.1
Packing group III



Danger
H318
P280i-P305+P351+P338-P310a

Azure II > RS - For microscopy - C.I. 52010/52015 — RS

Description ..... Polvere nera con riflessi verdi Maximum absorption ..... 645 ÷ 650 nm Loss on drying at 110°C ..... ≤ 15 %
Identification ..... Positive E 1% / 1 cm on dried substance ..... 1850 ÷ 2100

Code	Size	Packaging	Notes
424721	5 g	Glass bottle	

Dye for bacteriology and hematology. Mix Azur - methylene blue

Azure II eosin

• Azzurro II eosina • Azur II Eosine • Azur II Eosina • Azur II Eosin

Synonym:
Azure II eosinate

CAS: 53092-85-6

Classification transport
ONU: 2811
Transport Hazard class: 6.1
Packing group III



Danger
H332-H318
P261-P271-P280i-P304+P340-P310a-P305+P351+P338

Azure II eosin > RS - For microscopy - C.I. 52010/52015/45380

RS

Description ..... Dark blue powder
Identification ..... Positive
Loss on drying (110°C)..... ≤ 10 %
Maximum absorption lambda max1645.0 ÷ 655.0 nm
Spec. absorption (E1%/1cm) lambda max1≥ 1000
Maximum absorption lambda max2515.0 ÷ 525.0 nm
Spec. absorption (E1%/1cm) lambda max2≥ 400

Code	Size	Packaging	Notes
424731	5 g	Glass bottle	

Dye for bacteriology, hematology histopathology. Mix Azur - Methylene blue - eosin

Product specifications are subject to changes. Please visit our website for updates.

Barbituric acid

• Acido barbiturico • Acide barbiturique • Acido barbitúrico • Barbitursäure

Synonym:
2,4,6-Trihydroxypyrimidine

$NHCONHCOCH_2CO$
Molecular Weight: 128,12
CAS: 67-52-7
EEC-N: 200-658-0

Barbituric acid > RPE - For analysis - Reag. Ph. Eur.

RPE

Description ..... White crystalline powder Identification ..... Positive Melting point..... 250 ÷ 256 ° C Assay (acidimetric) ..... 99.0 ÷ 101.0 %

Code	Size	Packaging	Notes
402532	25 g	Glass bottle	
402535	250 g	Glass bottle	

Barium standard solution

• Bario standard soluzione • Baryum standard solution • Bario, solución patrón • Barium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Warning
H290
P234-P390-P406

Barium standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615000601	100 ml	Plastic bottle	A 0.1 % solution Ref Ph.Eur 5000601
615000609	100 ml	Plastic bottle	A 50 ppm solution: to dilute according to Ref Ph.Eur 5000600

Barium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505327	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505328	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505329	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Barium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503451	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503453	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503455	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid
503457	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Barium standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
507527	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507481	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
E497445	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
E497441	500 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Barium standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
424861		Plastic ampoule	conc. 1.000 ppm Matrix: Nitric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Barium acetate

• Bario acetato • Baryum acétate • Bario acetato • Bariumacetat

$Ba(CH_3COO)_2$
Molecular Weight: 255,42
CAS: 543-80-6
EEC-N: 208-849-0

Warning
H302-H332
P261-P264-P271-P301+P312a-P304+P340-P501a

Barium acetate > RPE - For analysis - ACS

RPE

K.....≤ 30 ppm
Description ..... White crystalline powder
Identification ..... Positive
Water-insoluble matter .....≤ 100 ppm
Oxidizing subst.(NO3) .....≤ 50 ppm
Chloride.....≤ 10 ppm
Heavy metals (Pb).....≤ 5 ppm
Ca.....≤ 500 ppm
Na.....≤ 50 ppm
Fe .....≤ 10 ppm
Sr.....≤ 0.2 %
Assay (complexometric) ....99.0 ÷ 102.0 %

Code	Size	Packaging	Notes
424895	100 g	Plastic bottle	
424896	500 g	Plastic bottle	
424897	1 kg	Plastic bottle	

Barium carbonate

• Bario carbonato • Baryum carbonate • Bario carbonato • Bariumcarbonat

$BaCO_3$
Molecular Weight: 197,34
CAS: 513-77-9
EEC-N: 208-167-3

Warning
H302
P264-P270-P301+P312a-P330-P501a

Barium carbonate > RPE - For analysis - ACS - Reag. Ph.Eur.

RPE

Description ..... White powder
Identification ..... Positive
Water-soluble titrable base .≤ 0.002 meq/g
HCl-insoluble matter.....≤ 150 ppm
Oxidizing subst.(NO3) .....≤ 50 ppm
Chloride.....≤ 20 ppm
Sulphide .....≤ 10 ppm
Heavy metals (Pb).....≤ 10 ppm
Ca.....≤ 500 ppm
Fe .....≤ 20 ppm
Sr..... ≤ 0.7 %
Assay (alkalimetric).....99.0 ÷ 101.0 %
K.....≤ 50 ppm
Na.....≤ 200 ppm

Code	Size	Packaging	Notes
424945	250 g	Plastic bottle	
424943	25 kg	Plastic bucket	

Barium carbonate > RE - Pure

RE

Description ..... Hazel-white powder Identification ..... Positive Fe .....≤ 50 ppm Assay (alkalimetric)..... ≥ 99.0 %

Code	Size	Packaging	Notes
321507	1 kg	Plastic bottle	
321502	25 kg	Plastic bucket	

Barium chloride dihydrate

• Bario cloruro diidrato • Baryum chlorure dihydraté • Bario cloruro dihidratado • Bariumchlorid-Dihydrat

$BaCl_2$
Molecular Weight: 208,23
CAS: 10361-37-2

Classification transport
ONU: 1564
Transport Hazard class: 6.1
Packing group III

Danger
H301-H332-H319
P261-P271-P301+P310a-P304+P340-P305+P351+P338-P337+P313

Barium chloride dihydrate > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611009300	100 g	Glass bottle	Ref Ph.Eur 1009300

Product specifications are subject to changes. Please visit our website for updates.

Barium chloride dihydrate > RPE - For analysis - ACS

RPE

Appearance	White crystalline powder	Ca	≤ 0.05 %	Na	≤ 0.005 %	Loss on drying	14.0 - 16.0 %
pH	5.2 - 8.2	Fe	≤ 2 ppm	Heavy metals (as Pb)	≤ 5 ppm	Insoluble matter	≤ 0.005 %
BaCl2.2H20	≥ 99.0 %	Sr	≤ 0.1 %	K	≤ 0.0025 %	Oxidizing substances (as NO3)	≤ 0.005 %

Code	Size	Packaging	Notes
425025	100 g	Plastic bottle	
425026	500 g	Plastic bottle	
425027	1 kg	Plastic bottle	
425029	5 kg	Plastic jar	
425022	25 kg	Plastic bucket	

Barium chloride dihydrate > RE - Pure

RE

Description ..... White crystals Identification ..... Positive Assay (argentimetric) ..... ≥98 %

Code	Size	Packaging	Notes
321757	1 kg	Plastic bottle	
321758	5 kg	Plastic tank	
321752	25 kg	Plastic bucket	

Barium chloride solution 10%

• Bario cloruro soluzione 10% • Baryum chlorure solution 10% • Bario cloruro solución 10% • Bariumchlorid 10%

$BaCl_2.2H_2O$
Molecular Weight: 244,27
CAS: 10326-27-9

Warning
H302-H319
P264-P280i-P301+P312a-P305+P351+P338-P337+P313-P501a

Barium chloride solution 10% > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Density at 15° C ..... 1.08 ÷ 1.10 Assay (argentimetric) ..... 9.5 ÷ 10.5 %

Code	Size	Packaging	Notes
E425101	1 l	Plastic bottle	

Suitable for water analysis

Barium chloride 0.1 mol/l (0.2N)

• Bario cloruro 0.1 mol/l (0.2N) • Baryum chlorure 0.1 mol/l (0.2N) • Bario cloruro 0.1 mol/l (0.2N) • Bariumchlorid 0.1 mol / l (0.2 N)

$BaCl_2.2H_2O$
CAS: 10326-27-9

HEU210

Barium chloride 0.1 mol/l (0.2N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613000600	1 l	Plastic bottle	Ref Ph.Eur 3000600

Barium chloride solution 61 g/l

• Bario cloruro soluzione 61 g/l • Baryum chlorure 61g/l • Bario cloruro solución 61 g/l • Bariumchlorid 61 g/l

$BaCl_2$
Molecular Weight: 208,25
CAS: 10361-37-2

Barium chloride solution 61 g/l > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611009303	100 ml	Plastic bottle	Barium chloride solution R1 Ref Ph.Eur 1009301
611009309	250 ml	Plastic bottle	Barium chloride solution R1 Ref Ph.Eur 1009301
611009301	1 l	Plastic bottle	Barium chloride solution R1 Ref Ph.Eur 1009301

Barium chloride 30 g/l

• Bario cloruro 30 g/l • Baryum chlorure 30 g/l • Bario cloruro 30 g/l • Bariumchlorid 30 g/l

$BaCl_2.2H_2O$
Molecular Weight: 244,27
CAS: 10326-27-9

Barium chloride 30 g/l > RS - For analysis according to JP

RS

Code	Size	Packaging	Notes
616001018	100 ml	Plastic bottle	Barium chloride TS

Barium chloride 30 g/l > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
617000161	100 ml	Plastic bottle	Barium chloride TS

Barium chromate

• Bario cromato • Baryum chromate • Bario cromato • Bariumchromat

$BaCrO_4$
Molecular Weight: 253,33
CAS: 10294-40-3
EEC-N: 233-660-5

Classification transport
ONU: 1479
Transport Hazard class: 5.1
Packing group II



Warning
H272-H302-H332
P210-P220-P261-P271-P280-P304+P340

Barium chromate > RE - Pure

RE

Description ..... Yellow powder
Identification ..... Positive
HCl-insoluble matter ..... ≤ 0.1 %
Chloride ..... ≤ 100 ppm
Fe ..... ≤ 20 ppm
Assay (oxidimetric) ..... ≥ 99 %

Code	Size	Packaging	Notes
425245	250 g	Glass bottle	

Barium hydroxide octahydrate

• Bario idrossido ottaidrato • Baryum hydroxyde octahydraté • Bario hidróxido octahidrato • Bariumhydroxid-Octahydrat

$Ba(OH)_2.8H_2O$
Molecular Weight: 315,48
CAS: 12230-71-6
EEC-N: 241-234-5

Classification transport
ONU: 1564
Transport Hazard class: 6.1
Packing group III



Danger
H302-H314-HEU071
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Barium hydroxide octahydrate > RPE - For analysis

RPE

Description ..... White crystals
Identification ..... Positive
Carbonate ..... ≤2.0 %
Chloride ..... ≤10 ppm
HCl-insoluble matter ..... ≤50 ppm
Heavy metals (Pb) ..... ≤5 ppm
Sulphide ..... ≤5 ppm
Ca ..... ≤50 ppm
Fe ..... ≤5 ppm
Sr ..... ≤1.5 %
Assay (alkalimetric) ..... ≥98.0 %

Code	Size	Packaging	Notes
425295	100 g	Plastic bottle	
425296	500 g	Plastic bottle	
425297	1 kg	Plastic bottle	
425292	25 kg	Plastic bucket	

Barium hydroxide octahydrate > RE - Pure

RE

Description ..... White crystals
Identification ..... Positive
Chloride ..... ≤500 ppm
HCl-insoluble matter ..... ≤500 ppm
Fe ..... ≤50 ppm
Assay (alkalimetric) ..... ≥95 %

Code	Size	Packaging	Notes
322007	1 kg	Plastic bottle	
322009	5 kg	Plastic tank	
322001	25 kg	Plastic bucket	
322004	50 kg	Plastic bucket	

Product specifications are subject to changes. Please visit our website for updates.

Barium hydroxide solution 5%

• Bario idrossido soluzione 5% • Baryum hydroxyde solution 5% • Bario hidróxido solución 5% • Bariumhydroxidlösung 5%

$Ba(OH)_2.8H_2O$
Molecular Weight: 315,48
CAS: 12230-71-6

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group II

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Barium hydroxide solution 5% > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... ≥1.03
Assay ..... 4.7 ÷ 5.3 % (p/p)

Code	Size	Packaging	Notes
E425301	1 l	Bottle	

Barium hydroxide solution 47.3 g/l

• Bario idrossido soluzione 47.3 g/l • Baryum hydroxyde solution 47.3 g/l • Bario hidróxido solución 47.3 g/l • Bariumhydroxidlösung 47.3 g / l

$Ba(OH)_2.8H_2O$
Molecular Weight: 315,48
CAS: 12230-71-6

Barium hydroxide solution 47.3 g/l > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611009409	250 ml	Plastic bottle	Ref Ph.Eur 1009401
611009401	1 l	Plastic bottle	Ref Ph.Eur 1009401

Barium nitrate

• Bario nitrato • Baryum nitrate • Bario nitrato • Bariumnitrat

$Ba(NO_3)_2$
Molecular Weight: 261,34
CAS: 10022-31-8
EEC-N: 233-020-5

Classification transport
ONU: 1446
Transport Hazard class: 5.1
Packing group II

Danger
H272-H301-H332-H319
P210-P261-P271-P280-P304+P340-P305+P351+P338

Barium nitrate > RPE - For analysis - ACS

RPE

Description ..... White crystals
Identification ..... Positive
pH sol. 5% at 25° C ..... 5.0 ÷ 8.0
Chloride ..... ≤ 10 ppm
Fe ..... ≤ 5 ppm
Ca ..... ≤ 0.005 %
Assay ..... ≥ 99.0 %
Heavy metals (Pb) ..... ≤ 10 ppm

Code	Size	Packaging	Notes
425342	500 g	Plastic bottle	
425347	1 kg	Plastic bottle	
425341	25 kg	Plastic bucket	

Barium perchlorate trihydrate

• Bario perclorato triidrato • Baryum perchlorate trihydraté • Bario perclorato trihidrato • Bariumperchlorat-trihydrat

$Ba(ClO_4)_2.3H_2O$
Molecular Weight: 390,29
CAS: 10294-39-0
EEC-N: 236-710-4

Classification transport
ONU: 1447
Transport Hazard class: 5.1
Packing group II

Danger
H271-H302-H332
P210-P261-P271-P280-P283-P304+P340

Barium perchlorate trihydrate > RPE - For analysis

RPE

Description ..... White crystals
Identification ..... Positive
Total nitrogen ..... ≤20 ppm
Chlorate ..... ≤300 ppm
Chloride ..... ≤10 ppm
Water-insoluble matter ..... ≤100 ppm
Methyl alcohol insolub ..... ≤100 ppm
Heavy metals (Pb) ..... ≤5 ppm
Substances not ppt. H2S ..... ≤0.1 %
Ca ..... ≤100 ppm
Cu ..... ≤25 ppm
Fe ..... ≤3 ppm
K ..... ≤100 ppm
Na ..... ≤100 ppm
Ni ..... ≤25 ppm
Pb ..... ≤25 ppm
Sr ..... ≤0.6 %
Zn ..... ≤20 ppm
Assay (complexometric) ..... ≥99.8 %

Code	Size	Packaging	Notes
425411	50 g	Glass bottle	

Barium perchlorate 0.05 mol/l

• Bario perclorato 0.05 mol/l • Baryum perchlorate 0.05 mol/l • Bario perclorato 0.05 mol/l • Bariumperchlorat 0.05 mol/l

$Ba(ClO_4)_2$
Molecular Weight: 336,23
CAS: 13465-95-7

Classification transport
ONU: 2924
Transport Hazard class: 3
Packing group II

Danger
H225-H315-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P332+P313

Barium perchlorate 0.05 mol/l > RS - For analysis according to Ph. Eur. Chap. 4.2.2 — RS

Code	Size	Packaging	Notes
613000700	1 l	Plastic bottle	Ref Ph.Eur 3000700

Barium perchlorate 0.025 mol/l

• Bario perclorato 0.025 mol/l • Baryum perchlorate 0.025 mol/l • Bario perclorato 0.025 mol/l • Bariumperchlorat 0.025 mol/l

Synonym:
Barium diperchlorate

$Ba(ClO_4)_2$
Molecular Weight: 336.23
CAS: 13465-95-7

Classification transport
ONU: 2924
Transport Hazard class: 3
Packing group II

Danger
H225-H315-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P332+P313

Barium perchlorate 0.025 mol/l > RS - For analysis according to Ph. Eur. Chap. 4.2.2 — RS

Code	Size	Packaging	Notes
613009601	500 ml	Plastic bottle	Ref Ph.Eur 3009600

Barium sulfate

• Bario solfato • Baryum sulfate • Bario sulfato • Bariumsulfat

Synonym:
Baryte

$BaSO_4$
Molecular Weight: 233,4
CAS: 7727-43-7
EEC-N: 231-784-4

Barium sulfate > RPE - For analysis — RPE

Description	White powder	Total nitrogen	≤30 ppm	Silicate	≤20 ppm	Na	≤500 ppm
Identification	Positive	Chloride	≤300 ppm	As	≤1 ppm	Ni	≤5 ppm
Organic substances	Conform	Phosphate	≤10 ppm	Cd	≤5 ppm	Pb	≤5 ppm
Loss on ignition	≤1.5 %	Heavy metals (Pb)	≤10 ppm	Cu	≤5 ppm	Zn	≤10 ppm
Acidity (H2SO4)	≤100 ppm	Soluble barium salts	≤50 ppm	Fe	≤20 ppm	Assay (complexometric)	≥97 %
Alkalinity(Ba idroxide)	≤40 ppm	Soluble salts	≤0.2 %	K	≤100 ppm		

Code	Size	Packaging	Notes
425497	1 kg	Plastic bottle	

Barium sulfate > RE - Pure — RE

Description	White powder	Chloride	≤0.1 %	Subst. reducing KMnO4	≤40 ppm(10m)
Identification	Positive	Phosphate	≤200 ppm	Fe	≤20 ppm
Loss on ignition	≤2 %	Nitrate	≤50 ppm	Zn	≤20 ppm

Code	Size	Packaging	Notes
322607	1 kg	Plastic bottle	

Baryte ▶ Barium sulfate

Basic violet 3 ▶ Crystal violet

Basic Violet 14 ▶ Fuchsin basic

Product specifications are subject to changes. Please visit our website for updates.

Benedict's reagent

• Benedict reattivo • Réactif de Bénédict • Reactivo de Benedict • Reagenz von Benedict

H412
P273-P501a

Benedict's reagent > RS - For microscopy

RS

Description ..... Clear blue liquid
Identification ..... Positive

Code	Size	Packaging	Notes
E425742	1 l	Glass bottle	

Suitable for the determination of glucose

Benzaldehyde

• Aldeide benzoica • Aldéhyde benzoïque • Aldehído benzoico • Benzaldehyd

Synonym:
Bitter almond

C_6H_5CHO
Molecular Weight: 106,12
CAS: 100-52-7
EEC-N: 202-860-4

Classification transport
ONU: 1990
Transport Hazard class: 9
Packing group III

Warning
H302
P264-P270-P301+P312a-P330-P501a

Benzaldehyde > RPE - For analysis - Reag. Ph. Eur.

RPE

Description ..... Clear colorless to pale yellow
Identification ..... Positive
Refractive index at 20°C ..... 1.5440 ÷ 1.5460
Acidity (benzoic acid) ..... ≤ 0.5 %
Assay (GLC) ..... ≥ 99.0 %

Code	Size	Packaging	Notes
415362	500 ml	Glass bottle	

Benzalkonium chloride

• Benzalconio cloruro • Benzalkonium chlorure • Benzalconio cloruro • Benzalkoniumchlorid

Synonym:
Alkylbenzyldimethylammonium chloride

$C_6H_5CH_2N(CH_3)_2RCl$
Molecular Weight: 365
CAS: 63449-41-2
EEC-N: 264-151-6

Classification transport
ONU: 3259
Transport Hazard class: 8
Packing group II

Danger
H302-H312-H314-H400
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P362+P364

Benzalkonium chloride > ERBApharm - According to pharmacopoeia: NF

ERBApharm

Description ..... Polvere o polvere gelatinosa bianco-giallognola
Identification ..... Positive
Water not sol. matter ..... Conform NF
Ammine estranee ..... Conform NF
Rapporto dei comp. alchilici (HPLC) Conform NF
Water ..... ≤ 15.0 %
Residue on calcination ..... ≤ 2.0 %
Assay ..... 97.0 ÷ 103.0 % s.s.

Code	Size	Packaging	Notes
322737	1 kg	Plastic bottle	
322738	5 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Benzalkonium chloride > RE - Pure

RE

Description ..... Yellow liqui-solid
Identification ..... Positive
Aspetto soluzione 1% ..... Conform
Water (K.F.) ..... ≤ 10 %
Sulphated ash ..... ≤ 0.1 %
Assay (non-aqueos medium) ..... 98.5 ÷ 103.0 % (s.s.)

Code	Size	Packaging	Notes
322721	250 g	Plastic bottle	

Benzene
• Benzene • Benzène • Benceno • Benzol

C_6H_6
Molecular Weight: 78,11
CAS: 71-43-2
EEC-N: 200-753-7

Classification transport
ONU: 1114
Transport Hazard class: 3
Packing group II



Danger
H225-H310-H315-H319-H340-H350-H372-H304-HA26
P210-P241-P280-P301+P310a-P303+P361+P353-P305+P351+P338

Benzene > RPE - For analysis - ACS — RPE

Description ..... Clear liquid	Ready carbonizable substances ..... Conform	Residue on evaporation ..... ≤10 ppm	Total sulphur ..... ≤5 ppm
Colour (APHA) ..... ≤10	Water (K.F.) ..... ≤500 ppm	Tiophene ..... Conform	Assay (GLC) ..... ≥99.0 %

Code	Size	Packaging	Notes
426113	2.5 l	Glass bottle	

Benzene-d6
• Benzene-d6 • Benzène-d6 • Benceno-d6 • Benzol-d6

Synonym:
Hexadeuterobenzene

C_6D_6
Molecular Weight: 84,07
CAS: 1076-43-3
EEC-N: 214-061-8

Classification transport
ONU: 1114
Transport Hazard class: 3
Packing group II



Danger
H225-H301-H311-H330-H350-HA26
P210-P280-P284-P303+P361+P353-P304+P340-P403+P233

Benzene-d6 > RS - For NMR - min 99.5% — RS

Code	Size	Packaging	Notes
P5089	10 x 0.75 ml	Glass ampoule	
P5085	25 ml	Glass bottle	
P5086	100 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Benzenemethanol ▶ Benzyl alcohol

Benzenesulfonyl chloride
• Benzenesolfonile cloruro • Benzène sulfonyle chlorure • Bencenosulfonilo cloruro • Benzolsulfonylchlorid

$C_6H_5SO_2Cl$
Molecular Weight: 176,62
CAS: 98-09-9
EEC-N: 202-636-6

Classification transport
ONU: 2225
Transport Hazard class: 8
Packing group III



Danger
H302-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Benzenesulfonyl chloride > RPE - For analysis — RPE

Description ..... Clear colourless liquid	Density at 20° C ..... ~ 1.377	Melting point ..... 15 ÷ 17 ° C
Identification ..... Positive	Refractive index at 20°C ..... ~ 1.552	Assay ..... ≥99 %

Code	Size	Packaging	Notes
426231	10 ml	Glass bottle	

For derivatization

Benzethonium chloride ▶ Hyamine 1622

Product specifications are subject to changes. Please visit our website for updates.

Benzoic acid

• Acido benzoico • Acide benzoique • Acido benzóico • Benzoesäure

C_6H_5COOH
Molecular Weight: 122,12
CAS: 65-85-0
EEC-N: 200-618-2



Danger
H315-H318-H372
P260-P264-P305+P351+P338-P310a-P362+P364-P332+P313

Benzoic acid > RS - For analysis according to Ph. Eur. Chap. 4.2.1 — RS

Code	Size	Packaging	Notes
612000200	100 g	Plastic bottle	Ref Ph.Eur 2000200

Benzoic acid > RS - Standard for volumetry — RS

Description .....White crystalline powder
Identification ..... Positive
Assay (acidimetric) .....≥99.8 %

Code	Size	Packaging	Notes
402621	10 g	Glass bottle	

Benzoic acid > RPE - For analysis - ACS — RPE

Description .....White crystalline powder
Identification (I.R.)..... Positive
Insol. in alcol metilico .....≤ 0.005 %
Subst. reducing KMnO4..... Conform
Heavy metals (Pb)..... ≤ 5 ppm
Chlorinated compounds.....≤ 0.005 %
Sulphur compounds .....≤ 0.002 %
Residue on calcination.....≤ 0.005 %
Assay (acidimetric) .....≥ 99.5 %

Code	Size	Packaging	Notes
402636	100 g	Plastic bottle	
402635	250 g	Plastic bottle	
402637	1 kg	Plastic bottle	

Benzoic acid > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP — ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Carbonizable substances ...Conform Ph.Eur.
Oxidizing substances .....Conform Ph.Eur.
Freezing point ..... 121 ÷ 123 ° C
Water (K.F.) ..... ≤0.7 %
Sulphated ash ..... ≤0.05 %
Halogenated and halides.....≤300 ppm
Heavy metals (Pb).....≤10 ppm
Assay (acidimetric) .....99.5 ÷ 100.5 %s.s.
Origin (BSE/TSE)..... Synthesis
Residual solvents (Current ICH).....Conform

Code	Size	Packaging	Notes
302087	1 kg	Plastic bottle	
302089	5 kg	Plastic tank	
302082	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Benzoic acid benzyl ester ▶ Benzyl benzoate

p-Benzoquinone

• p-Benzochinone • p-Benzoquinone • p-benzoquinona • p-Benzochinon

Synonym:
Quinone

$C_6H_4O_2$
Molecular Weight: 108,09
CAS: 106-51-4
EEC-N: 203-405-2

Classification transport
ONU: 2587
Transport Hazard class: 6.1
Packing group II



Danger
H301-H331-H315-H319-H335-H400
P304+P340-P311a-P305+P351+P338-P330-P362+P364-P403+P233

p-Benzoquinone > RPE - For analysis — RPE

Description .....Yellow-green powder
Identification ..... Positive
Melting point..... 112.0 ÷ 116.0 ° C
Water .....≤ 0.5 %
Assay (iodometric).....≥ 98.5 %

Code	Size	Packaging	Notes
436853	50 g	Glass bottle	
436854	100 g	Plastic bottle	

For spectrophotometric microdetermination of amines

Benzyl alcohol

• Alcole benzilico • Alcool benzylique • Alcohol bencílico • Benzylalkohol

Synonym:
Benzenemethanol

$C_6H_5CH_2OH$
Molecular Weight: 108,14
CAS: 100-51-6
EEC-N: 202-859-9

Warning
H302-H312-H332-H319
P261-P264-P271-P304+P340-P305+P351+P338-P337+P313

Benzyl alcohol > RPE - For analysis - Stabilized with 0,02% BHA

RPE

Description	Clear colourless liquid	Boiling point	204.9 ÷ 205.9 °C	Benzaldheide(GLC)	≤0.1 %	Peroxides (H2O2)	≤10 ppm
Identification (I.R.)	Conform	Acidity (benzoic acid)	≤200 ppm	Total chlorine	≤50 ppm	Residue on ignition	≤20 ppm
Water miscibility	Conform	Water (K.F.)	≤0.1 %	Carbonyl Compounds (CO)	≤100 ppm	Fe	≤10 ppm
Refractive index at 20°C	1.5376 ÷ 1.5416	Alcalinity (NaOH)	≤34 ppm	Heavy metals (Pb)	≤5 ppm	Assay (GLC)	≥99.5 %

Code	Size	Packaging	Notes
414052	1 l	Glass bottle	

Store between 2-8 °C

Benzyl alcohol > RPE - For analysis

RPE

Description	Clear liquid	Density at 20° C	~ 1.04	Water (K.F.)	≤0.1 %	Assay (GLC)	≥99.8 %
Colour (APHA)	≤10	Refractive index at 20° C	~ 1.54	Total chlorine	≤100 ppm		
Identification	Positive	Boiling point	204.5 ÷ 205.5 °C	Benzaldeide	≤0.1 %		

Code	Size	Packaging	Notes
414022	1 l	Glass bottle	
414024	2.5 l	Glass bottle	

Store between 2-8 °C

Benzyl alcohol > ERBApharm - According to pharmacopoeia: BP-FU-NF-Ph.Eur.-Ph.Franc.

ERBApharm

Description	Clear colourless liquid	Refractive index at 20°C	1.538 ÷ 1.541	Solubility	Conform Ph. Eur.	0.3 %	
Identification	Positive	Residue on evaporation	≤500 ppm	Cyclohexylmethanol (GLC)	≤0.10 %	Assay (acidimetric)	98.0 ÷ 100.5 %
Appearance of solution	Conform Ph.Eur.	Analogous subst. GLC	Conform Ph.Eur.	Peaks sum rel.ret. less than C6H5CH2OH ≤ 0.04 %			
Acidity	Conform Ph.Eur.	Benzaldehyde (GLC)	≤0.15 %	Peaks sum rel.ret. great than C6H5CH2OH≤			
Density at 20° C	1.043 ÷ 1.049	Peroxide value	≤5				

Code	Size	Packaging	Notes
308131	1 l	Glass bottle	
308132	2.5 l	Glass bottle	
308138	23 kg	Drum	
308137	200 l	Drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Benzyl benzoate

• Benzile benzoato • Benzyle benzoate • Bencilo benzoato • Benzylbenzoat

Synonym:
Benzoic acid benzyl ester

$C_6H_5COOCH_2C_6H_5$
Molecular Weight: 212,25
CAS: 120-51-4
EEC-N: 204-402-9

Warning
H302-H411
P264-P270-P301+P312a-P330-P391-P501a

Benzyl benzoate > RPE - For analysis

RPE

Description	Clear liquid	Refractive index at 20°C	1.5651 ÷ 1.5711	Acidity (benzoic acid)	≤ 0.12 %
Identification	Positive	Melting point	18.6 ÷ 20.2 °C	Carbonyl Compounds (CO)	≤ 0.04 %
Density at 20°C	1.118 ÷ 1.122	Water (K.F.)	≤ 0.05 %	Assay (GLC)	≥ 99 %

Code	Size	Packaging	Notes
426761	250 ml	Glass bottle	
426763	1 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Benzyl benzoate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP

ERBApharm

Description ..... Clear colourless liquid
Identification ..... Positive
Acidity ..... Conform Ph.Eur.
Organic volatile impurities Conform USP-NF
Density at 20° C ..... 1.118 ÷ 1.122
Density at 25° C ..... 1.116 ÷ 1.120
Refractive index at 20°C ..... 1.568 ÷ 1.570
Freezing point ..... ≥ 17.0 ° C
Aldehyde ..... ≤ 0.05 %
Sulphated ash ..... ≤ 0.1 %
Assay (alkalimetric) ..... 99.0 ÷ 100.5 %

Code	Size	Packaging	Notes
323101	1 l	Glass bottle	
323102	2.5 l	Glass bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Beryllium standard solution

• Berillio standard soluzione • Béryllium standard solution • Berilio, solución patrón • Beryllium-Standardlösung

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group II

Beryllium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505332	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505335	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505333	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Beryllium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503461	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503463	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Beryllium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
506941	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507497	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Bismarck brown R

• Bruno Bismarck R • Brun Bismarck R • Pardo de Bismarck R • Bismarck brown R

$C_{21}H_{24}N_8.2HCl$
Molecular Weight: 461,39
CAS: 8005-78-5
EEC-N: 232-341-8

Warning
H302-H312-H332
P261-P264-P271-P280h-P301+P312a-P304+P340

Bismarck brown R > RS - For microscopy - C.I. 21010

RS

Description ..... Brown powder
Identification ..... Positive
E 1% / 1 cm a 460 nm ..... ≥250
Loss on drying (110°C) ..... ≤10 %

Code	Size	Packaging	Notes
431252	25 g	Glass bottle	

Dye for bacteriology, histology

Bismuth, granules

• Bismuto, granelli • Bismuth, granules • Bismuto, gránulos • Wismut, Granulat

Bi
Molecular Weight: 208,98
CAS: 7440-69-9
EEC-N: 231-177-4

Bismuth, granules > RPE - For analysis

RPE

Description ..... Metallic granules Identification ..... Positive Melting point ..... ~ 272 °C

Code	Size	Packaging	Notes
428064	100 g	Glass bottle	

Bismuth standard solution

• Bismuto standard soluzione • Bismuth standard solution • Bismuto, solución patrón • Wismut-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group II

Danger
H290-H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Bismuth standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615005300	100 ml	Plastic bottle	A 100 ppm solution: to dilute according to Ref Ph.Eur 5005300

Bismuth standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505337	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505338	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505339	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Bismuth standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503471	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503473	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503475	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503477	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Bismuth standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
507528	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507482	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
E497455	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
E497451	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Bismuth standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
428071		Plastic ampoule	conc. 1.000 ppm Matrix: Nitric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Bismuth(III) carbonate basic

• Bismuto carbonato basico • Bismuth (III) carbonate basique • Bismuto (III) carbonato basico
• Wismut (III) basisches Carbonat

Synonym:
Bismuth subcarbonate

$(BiO)_2CO_3$
Molecular Weight: 509,97
CAS: 5892-10-4
EEC-N: 227-567-9

Bismuth(III) carbonate basic > RPE - For analysis

RPE

Description ..... Yellowish powder
Identification ..... Positive
Loss on drying ..... ≤ 1.0 %
Ag ..... ≤ 25 ppm
As ..... ≤ 5 ppm
Cu ..... ≤ 50 ppm
Pb ..... ≤ 20 ppm
Assay (complexometric) .... 97.6 ÷ 100.7 %
Alkaly-alkaline earth ..... ≤ 1.0 %
Chloride ..... ≤ 0.05 %
Nitrate ..... ≤ 0.4 %

Code	Size	Packaging	Notes
428103	50 g	Glass bottle	
428105	250 g	Plastic bottle	

Bismuth(III) nitrate basic

• Bismuto nitrato basico • Bismuth (III) nitrate basique • Bismuto (III) nitrato basico
• Wismut (III) nitrat basisch

Synonym:
Bismuth subnitrate

$Bi_5O(OH)_9(NO_3)_4$
Molecular Weight: 1461,99
CAS: 1304-85-4
EEC-N: 215-136-8

Classification transport
ONU: 1477
Transport Hazard class: 5.1
Packing group II

Danger
H272-H315-H319-H335
P210-P261-P280-P304+P340-P305+P351+P338-P403+P233

Bismuth(III) nitrate basic > RPE - For analysis - Reag. Ph. Eur.

RPE

Description ..... White powder Identification ..... Positive Assay (Bi)..... 71.0 - 74.0 %

Code	Size	Packaging	Notes
428294	100 g	Glass bottle	

Bismuth(III) nitrate basic > RPE - For analysis

RPE

Description ..... White powder
Identification ..... Positive
Chloride ..... ≤ 200 ppm
Assay (Bi)..... 71.0 - 74.0 %
Ag ..... ≤ 25 ppm
Cu ..... ≤ 50 ppm
Pb ..... ≤30 ppm

Code	Size	Packaging	Notes
428284	100 g	Glass bottle	
428286	500 g	Glass bottle	

For Dragendorff reagent

Bismuth(III) nitrate pentahydrate

• Bismuto (III) nitrato pentaidrato • Bismuth (III) nitrate pentahydraté • Bismuto (III) nitrato pentahidratado • Wismut (III) nitratpentahydrat

$Bi(NO_3)_3.5H_2O$
Molecular Weight: 485,07
CAS: 10035-06-0
EEC-N: 233-791-8

Classification transport
ONU: 1477
Transport Hazard class: 5.1
Packing group II

Danger
H272-H315-H319-H335
P210-P261-P280-P304+P340-P305+P351+P338-P403+P233

Bismuth(III) nitrate pentahydrate > RPE - For analysis

RPE

Description ..... White crystals or crystalline powder
Identification ..... Positive
Chloride ..... ≤ 10 ppm
Insoluble matter (in Nitric Acid) .. ≤ 0.005 %
Sulphate ..... ≤ 50 ppm
Ag ..... ≤ 10 ppm
As ..... ≤ 10 ppm
Cu ..... ≤ 20 ppm
Fe ..... ≤ 10 ppm
Ni ..... ≤ 10 ppm
Pb ..... ≤ 10 ppm
Zn ..... ≤ 10 ppm
Ca ..... ≤ 20 ppm
Mg ..... ≤ 10 ppm
Assay ..... ≥ 98.0 %

Code	Size	Packaging	Notes
428234	100 g	Glass bottle	
428236	500 g	Glass bottle	

Bismuth(III) nitrate pentahydrate > RE - Pure

RE

Description ..... White crystals
Identification ..... Positive
Chloride ..... ≤350 ppm
Sulphate ..... ≤600 ppm
As ..... ≤10 ppm
Assay (complexometric) ..... ≥98 %

Code	Size	Packaging	Notes
324185	250 g	Glass bottle	

N,O-Bis(trimethylsilyl)acetamide

• N,O-Bis (trimetilsilil)acetammide • N,O-Bis(triméthylsilyl)acétamide • N,O-Bis (trimetilsilil) acetamida
• N,O-Bis (trimethylsilyl) acetamid

Synonym:
BSA

$CH_3C[:NSi(CH_3)_3]OSi(CH_3)_3$
Molecular Weight: 203,43
CAS: 10416-59-8
EEC-N: 233-892-7

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group III

Warning
H226-H302
P210-P241-P264-P280-P303+P361+P353-P403+P235

N,O-Bis(trimethylsilyl)acetamide > RPE - For analysis

RPE

Description ..... Clear yellow liquid
Identification ..... Positive
Density at 20° C ..... ~ 0.835
Assay (GC) ..... ≥ 97 %

Code	Size	Packaging	Notes
489934	25 ml	Glass bottle	

For derivatization

N,O-Bis(trimethylsilyl)-trifluoroacetamide

• N,O-Bis (trimetilsilil)trifluoroacetammide • N,O-Bis(triméthylsilyl)trifluoroacétamide
• N,O-Bis (trimetilsilil) trifluoroacetamida • N,O-Bis-(trimethylsilyl) trifluoracetamid

Synonym:
BSTFA

$C_8H_{18}F_3NOSi_2$
Molecular Weight: 257,39
CAS: 25561-30-2
EEC-N: 247-103-9

Classification transport
ONU: 2920
Transport Hazard class: 8
Packing group II

Danger
H226-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

N,O-Bis(trimethylsilyl)-trifluoroacetamide > RPE - For analysis

RPE

Description ..... Yellow clear liquid
Identification ..... Positive
Density at 20° C ..... ≤ 0.985
Assay (GLC) ..... ≥ 98 %

Code	Size	Packaging	Notes
489561	25 ml	Glass bottle	

For derivatization

Product specifications are subject to changes. Please visit our website for updates.

Biuret 97%

• Biureto 97% • Biuret 97% • Biuret 97% • Biuret 97%

Synonym:
• *Allophanic acid amide*
• *Carbamoyl urea*

$NH_2CONHCONH_2$
Molecular Weight: 103,08
CAS: 108-19-0
EEC-N: 203-559-0

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Biuret 97% > RS - For microscopy

RS

Description ..... White crystalline powder Identification ..... Positive Water (K.F.) ..... ≤ 2 %

Code	Size	Packaging	Notes
428432	25 g	Glass bottle	

Biuret reagent

• Biureto reattivo • Réactif au biuret • Biuret reactivo • Biuret-Reagenz

Synonym:
• *Allophanic acid amide*
• *Carbamoyl urea*

$NH_2CONHCONH_2$
Molecular Weight: 103,08
CAS: 108-19-0
EEC-N: 203-559-0

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Biuret reagent > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611011601	1 l	Plastic bottle	Biuret reagent Ref Ph.Eur 1011601

Boric acid

• Acido borico • Acide borique • Acido bórico • Borsäure

H_3BO_3
Molecular Weight: 61,83
CAS: 10043-35-3
EEC-N: 233-139-2

Danger
H360FD-HA26
P201-P202-P280-P308+P313-P405-P501a

Boric acid > RPE - For analysis - ISO - ACS - Reag.Ph.Eur. - Reag.USP

RPE

Description ..... White crystalline powder
Identification ..... Positive
Chloride ..... ≤10 ppm
Phosphate ..... ≤10 ppm
Sulphate ..... ≤100 ppm
Heavy metals (Pb) ..... ≤10 ppm
Ca ..... ≤50 ppm
Fe ..... ≤10 ppm
Insoluble in methanol ..... ≤50 ppm
Nonvolatile with methanol ..... ≤500 ppm
Appearance of solution ..... Conform
Alcohol solubility ..... Conform
Organic substances ..... Conform
pH solution 3.3% ..... 3.8 ÷ 4.8
Loss on drying ..... ≤0.5 %
Assay (acidimetric) ..... 99.5 ÷ 100.5 % s.s.

Code	Size	Packaging	Notes
402763	100 g	Plastic bottle	
402766	500 g	Plastic bottle	
402767	1 kg	Plastic bottle	
402765	10 kg	Plastic tank	
402762	25 kg	Plastic bucket	
402764	50 kg	Fibre drum	

Boric acid > ERBApharm - According to pharmacopoeia: NF

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution ..... Conform Ph.Eur.
Alcohol solubility ..... Conform Ph.Eur.
pH (sol. 3.3.%) ..... 3.8 ÷ 4.8
Organic substances ..... Conform Ph.Eur.
Loss on drying ..... ≤ 0.5 %
Sulphate ..... ≤ 450 ppm
Fe ..... ≤ 1 ppm
Heavy metals (Pb) ..... ≤ 15 ppm
Assay (acidimetric) ..... 99.5 ÷ 100.5 % s.s.
Origin (BSE/TSE) ..... Synthesis
Residual solvents (Current ICH) ..... Conform

Code	Size	Packaging	Notes
302185	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Boric acid > ERBApharm - According to pharmacopoeia: Ph.Eur.-FU-Ph.Franc.-DAB-USP

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Alcohol solubility .....Conform Ph.Eur.
Organic substances.....Conform Ph.Eur.
pH (sol. 3.3.%).....3.8 ÷ 4.8
Loss (silica gel) ..... ≤0.5 %
Sulphate.....≤450 ppm
Fe .....≤10 ppm
Heavy metals (Pb).....≤15 ppm
Assay (acidimetric) .....99.5 ÷ 100.5 % s.s.
Origin (BSE/TSE)..... Synthesis
Residual solvents (Current ICH).....Conform

Code	Size	Packaging	Notes
302177	1 kg	Plastic bottle	
302179	5 kg	Plastic tank	
302178	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Boric acid 4%

• Acido borico 4% • Acide borique 4% • Acido bórico 4% • Borsäure 4%

H_3BO_3 HEU210
Molecular Weight: 61,83
CAS: 10043-35-3

Boric acid 4% > RS - For agroalimentary analysis

RS

Description ..... Clear colourless liquid
Assay ..... 3.975 ÷ 4.025 %

Code	Size	Packaging	Notes
502002	5 l	Plastic tank	

According to NF V04-387: H3BO3: 40 g/L water QSP 1 L

Boric acid 4% with indicator

• Acido borico 4% con indicatore • Acide borique 4% avec indicateur • Acido bórico 4% con indicador • Borsäure 4% mit Indikator

H_3BO_3 HEU210
Molecular Weight: 61,83
CAS: 10043-35-3

Boric acid 4% with indicator > RS - For agroalimentary analysis

RS

Description ..... Clear mauve liquid
Assay ..... 3.975 ÷ 4.025 %

Code	Size	Packaging	Notes
502601	5 l	Plastic tank	

According to NF V04-211: H3BO3: 40 g/L water QSP 1L. Mixed indicator (Methyl red+Bromocresol green): 10 ml

Boric acid 3%

• Acido borico 3% • Acide borique 3% • Acido bórico 3% • Borsäure 3%

H_3BO_3 HEU210
Molecular Weight: 61,83
CAS: 10043-35-3

Boric acid 3% > RPE - For analysis

RPE

Boric acid content ..... 2.9 - 3.1 %

Code	Size	Packaging	Notes
PS0563/21	2.5 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Boric acid 1% with indicator

• Acido borico 1% con indicatore • Acide borique 1% avec indicateur • Acido bórico 1% con indicador • Borsäure 1%

H_3BO_3 HEU210
Molecular Weight: 61,83
CAS: 10043-35-3

Boric acid 1% with indicator > RS - For agroalimentary analysis

RS

Description ..... Clear grey/green liquid Assay ..... 0.993 ÷ 1.007 %

Code	Size	Packaging	Notes
502611	5 l	Plastic tank	
502612	10 l	Plastic tank	

According to NF V04-211: H3BO3: 10 g/L water QSP 1L. Mixed indicator (Methyl red+Bromocresol green): 10 ml

Boric acid 20g/l

• Acido borico 20g/l • Acide borique 20g/l • Acido bórico 20g/l • Borsäure 20 g/l

H_3BO_3 HEU210
Molecular Weight: 61,83
CAS: 10043-35-3

Boric acid 20g/l > RS - For analysis

RS

Concentration .....19 - 21 g/l

Code	Size	Packaging	Notes
PS0703/22	5 l	Plastic tank	

Boric acid 20 g/l with indicator

• Acido borico 20 g/l con indicatore • Acide borique 20 g/l avec indicateur • Acido bórico 20 g/l con indicador • Borsäure 20 g/l mit Indikator

H_3BO_3 HEU210
Molecular Weight: 61,83
CAS: 10043-35-3

Boric acid 20 g/l with indicator > RS - For nitrogen detection according to Kjeldahl

RS

Refractive index at 20°C..... 1.331 - 1.335 Density d20/4 ..... 1.003 - 1.007

Code	Size	Packaging	Notes
PS0562/22	5 l	Plastic tank	

Composition: Boric acid 20g, Red methyl solution 0.5 g/l: 10ml, Methylene blue solution 1.5 g/l: 2ml, water: QSP 1l according to ISO 5663-1984

Boric buffer solution

• Tampone borico soluzione • Solution tampon borique • Tampón bórico solución • Boric Pufferlösung

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Boric buffer solution > RS - For analysis

RS

Temperature of measurement ..... 19 - 21 °C pH..... 9.9 - 10.9 unite pH

Code	Size	Packaging	Notes
PS0226/41	10 l	Plastic tank	

Composition: boric acid 33g/l, potassium chloride 39.7g/l, sodium hydroxide 18.75g/l

2-Bornanone ▶ Camphor natural

Boron standard solution
• Boro standard soluzione • Bore standard solution • Boro, solución patrón • Bor Standardlösung

HEU210

Boron standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505322	100 ml	Plastic bottle	conc. 10 ppm Matrix: Water
505325	100 ml	Plastic bottle	conc. 100 ppm Matrix: Water
505323	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Boron standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503441	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503443	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503445	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Water
503447	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Boron standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
E497465	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Water
E497461	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Boron standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
429641		Plastic ampoule	conc. 1.000 ppm Matrix: Water - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Brilliant Blue R ▶ Coomassie brilliant blue R 250

Brillant cresyl blue
• Blu cresile brillante • Bleu de crésyl brillant • Azul de cresol brillante • Brillantes Kresylblau

$C_{17}H_{21}N_4O$. 0.5 $ZnCl_4$
Molecular Weight: 400,97
CAS: 81029-05-2
EEC-N: 279-675-0

Brillant cresyl blue > RS - For microscopy

RS

Description ..... Polvere verde-bluastra scura Identification ..... Positive E (1%/1cm) a 622 nm (in ETOH 50%) ..... 1500 ÷ 1700

Code	Size	Packaging	Notes
428811	10 g	Glass bottle	
428812	25 g	Glass bottle	

Dye for cytology and hematology

Product specifications are subject to changes. Please visit our website for updates.

Brilliant green

• Verde brillante • Vert brillant • Verde brillante • Brillantgrün

Synonym:
• *Basic Green 1*
• *Diamond green*

$C_{27}H_{34}N_2O_4S$
Molecular Weight: 482,64
CAS: 633-03-4
EEC-N: 211-190-1

Warning
H302
P264-P270-P301+P312a-P330-P501a

Brilliant green > RPE - For analysis - C.I. 42040

RPE

Description ..... Green crystline powder Identification ..... Positive Loss on drying ..... ≤ 10 % Colour change..... yellow green

Code	Size	Packaging	Notes
491152	25 g	Glass bottle	

Dye for microscopy (bacteriology-Botanical-histology). Indicator acid - base (pH 0.1 ÷ 2.6)

Brix Standards

• Brix Standards • Etalons en degré Brix • Brix patrones • Hengste im Brix-Grad

Brix Standards > RS - For calibration

RS

Code	Size	Packaging	Notes
540201	15 ml	Ampoule	Sucrose Stabilised 0%, Nominal Refractive Index @ 20°C: 1.332986
540202	15 ml	Ampoule	Sucrose Stabilised 5%, Nominal Refractive Index @ 20°C: 1.340264
540203	15 ml	Ampoule	Sucrose Stabilised 7%, Nominal Refractive Index @ 20°C: 1.343253
540204	15 ml	Ampoule	Sucrose Stabilised 10%, Nominal Refractive Index @ 20°C: 1.347824
540205	15 ml	Ampoule	Sucrose Stabilised 11.5%, Nominal Refractive Index @ 20°C: 1.350149
540206	15 ml	Ampoule	Sucrose Stabilised 12%, Nominal Refractive Index @ 20°C: 1.350930
540207	15 ml	Ampoule	Sucrose Stabilised 15%, Nominal Refractive Index @ 20°C: 1.355679
540208	15 ml	Ampoule	Sucrose Stabilised 20%, Nominal Refractive Index @ 20°C: 1.363842
540209	15 ml	Ampoule	Sucrose Stabilised 25%, Nominal Refractive Index @ 20°C: 1.372328
540210	15 ml	Ampoule	Sucrose Stabilised 30%, Nominal Refractive Index @ 20°C: 1.381149
540220	15 ml	Ampoule	Sucrose Stabilised 35%, Nominal Refractive Index @ 20°C: 1.390322
540221	15 ml	Ampoule	Sucrose Stabilised 40%, Nominal Refractive Index @ 20°C: 1.39986
540222	15 ml	Ampoule	Sucrose Stabilised 45%, Nominal Refractive Index @ 20°C: 1.409777
540223	15 ml	Ampoule	Sucrose Stabilised 50%, Nominal Refractive Index @ 20°C: 1.420087
540224	15 ml	Ampoule	Sucrose Stabilised 60%, Nominal Refractive Index @ 20°C: 1.441928
540225	15 ml	Ampoule	Sucrose Stabilised 67.5%, Nominal Refractive Index @ 20°C: 1.459290

Store between 2 - 8 °C

Bromate standard solution

• Bromati standard soluzione • Bromate standard solution • Bromato, solución patrón • Bromat-Standardlösung



Danger
H350-HA26
P201-P202-P280-P308+P313-P405-P501a

Bromate standard solution > RS - Standard solution for ion chromatography

RS

Code	Size	Packaging	Notes
503171	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503173	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Bromide - bromate 0.0167 mol/l

• Bromuri - bromati 0.0167 mol/l • Bromure - bromate 0.0167 mol/l • Bromuro - bromato 0.0167 mol/l • Bromidbromat 0.0167 mol/l



Danger
H350-HA26
P201-P202-P280-P308+P313-P405-P501a

Bromide - bromate 0.0167 mol/l > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613001000	1 l	Plastic bottle	Ref Ph.Eur 3001000

Bromide standard solution

• Bromuri standard soluzione • Bromure standard solution • Bromuro, solución patrón • Bromid-Standardlösung

Bromide standard solution > RS - Standard solution for ion chromatography

RS

Code	Size	Packaging	Notes
503211	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503213	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Bromine solution

• Bromo soluzione • Brome solution • Bromo, solución patrón • Brome lösung

Br_2
Molecular Weight: 159,82
CAS: 7726-95-6

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group I



Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Bromine solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611012401	100 ml	Glass bottle	Ref Ph.Eur 1012401

Bromine solution > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
617000141	100 ml	Glass bottle	Brome TS

Product specifications are subject to changes. Please visit our website for updates.

Bromine water

• Acqua di bromo • Eau de brome • Agua de bromo • Bromwasser

Br_2
Molecular Weight: 159,82
CAS: 7726-95-6

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group I

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Bromine water > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611012409	50 ml	Glass bottle	Ref Ph.Eur 1012402
611012402	100 ml	Glass bottle	Ref Ph.Eur 1012402

Storage: protected from light

Bromocresol green

• Verde bromocresolo • Vert de bromocrésol • Verde de bromocresol • Bromkresolgrün

$C_{21}H_{14}Br_4O_5S$
Molecular Weight: 698,05
CAS: 76-60-8
EEC-N: 200-972-8

Bromocresol green > RPE - For analysis - ACS

RPE

Description ..... Tan powder
Identification ..... Positive
Appearance of solution ..... Conform
Colour change ..... yellow-blue
pH range ..... 3.8 ÷ 5.4

Code	Size	Packaging	Notes
491207	1 g	Glass bottle	
491208	25 g	Glass bottle	

Clark indicator series. Complexometric indicator. Dye for microscopy (botanical-histology)

Bromocresol green 0.04% hydroalcoholic solution

• Verde bromocresolo soluzione 0.04% idroalcolica • Vert de bromocrésol 0.04% solution hydroalcoolique
• Verde de bromocresol 0.04% solución hidroalcohólica • Bromkresolgrün 0.04% ige hydroalkoholische Lösung

$C_{21}H_{14}Br_4O_5S$
Molecular Weight: 698,05
CAS: 76-60-8

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Bromocresol green 0.04% hydroalcoholic solution > RPE - For analysis

RPE

Description ..... Clear green liquid
Identification ..... Positive
Sensitivity(pH 4.0-5.4) ..... Conform
Colour change ..... yellow blue

Code	Size	Packaging	Notes
E491255	250 ml	Glass bottle	

Clark indicator series. Acid-base indicator (pH 3.8 ÷ 5.4)

Bromocresol green solution

• Verde bromocresolo soluzione • Vert de bromocrésol solution • Verde de bromocresol solución • Bromkresolgrünlösung

$C_{21}H_{14}Br_4O_5S$
Molecular Weight: 698,05
CAS: 76-60-8

Bromocresol green solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611012601	100 ml	Plastic bottle	Ref Ph.Eur 1012601 / Colur change: pH 3.6 (yellow) to pH 5.2 (blue)

Bromocresol green solution > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
617000101	100 ml	Plastic bottle	Bromocresol green TS

Bromocresol green - Methyl red solution

• Verde bromocresolo - Rosso metile soluzione • Vert de bromocrésol - Rouge de méthyle solution • Verde de bromocresol - Rojo de metilo solución • Bromkresolgrün - Methylrotlösung

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II

Danger
H225-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

Bromocresol green - Methyl red solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611012602	100 ml	Plastic bottle	Ref Ph.Eur 1012602

Bromocresol purple

• Porpora bromocresolo • Pourpre de bromocrésol • Púrpura de bromocresol • Bromkresolpurpur

Synonym:
5,5'-Dibromo-o-cresolsulfonphthalein

$C_{21}H_{16}Br_2O_5S$
Molecular Weight: 540,24
CAS: 115-40-2
EEC-N: 204-087-8

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Bromocresol purple > RPE - For analysis

RPE

Description ..... Pink powder Identification ..... Positive Colour change..... Purple - yellow pH range ..... 5.2 ÷ 6.8

Code	Size	Packaging	Notes
470038	5 g	Glass bottle	
470039	25 g	Glass bottle	

Clark indicator series. Acid-base indicator (pH 5.2 ÷ 6.8) Dye for microscopy (histology)

Bromocresol purple solution 0.4% in ethanol

• Porpora bromocresolo soluzione 0.4% in etanolo • Pourpre de bromocrésol solution 0.4% dans l'éthanol • Púrpura de bromocresol solución 0.4% en alcohol etilico • Bromkresol-Purpur-Lösung 0.4% in Ethanol

Synonym:
5,5'-Dibromo-o-cresolsulfonphthalein

$C_{21}H_{16}Br_2O_5S$
Molecular Weight: 540,24
CAS: 115-40-2

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Bromocresol purple solution 0.4% in ethanol > RPE - For analysis

RPE

Description ..... Red clear liquid Identification ..... Positive pH range ..... 5.2 - 6.8

Code	Size	Packaging	Notes
E470045	250 ml	Glass bottle	

Clark indicator series. Acid-base indicator (pH 5.2 ÷ 6.8) Dye for microscopy (histology)

Product specifications are subject to changes. Please visit our website for updates.

Bromocresol purple solution

• Porpora bromocresolo soluzione • Pourpre de bromocrésol solution • Púrpura de bromocresol solución • Bromkresol-Purpur-Lösung

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group III

Warning
H226-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

Bromocresol purple solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611012701	100 ml	Plastic bottle	Ref Ph.Eur 1012701

Colour change: pH 5.2 (yellow) to pH 6.8 (bluish-violet)

alpha-Bromonaphthalene

• alfa-Bromonaftalene • alpha-Bromonaphtalène • alpha-Bromonaftaleno • alpha-Bromnaphthalin

$C_{10}H_7Br$
Molecular Weight: 207,07
CAS: 90-11-9
EEC-N: 201-965-2

Warning
H302-H319
P264-P280i-P301+P312a-P305+P351+P338-P337+P313-P501a

alpha-Bromonaphthalene > RPE - For analysis

RPE

Description ..... Yellow or brown clear liquid
Identification ..... Positive
Density at 20° C ..... 1.479 ÷ 1.485
Refractive index at 20°C ..... 1.6562 ÷ 1.6602
Melting point ..... 2.0 ÷ 4.0 ° C
Assay (GLC) ..... ≥99 %

Code	Size	Packaging	Notes
430652	500 ml	Glass bottle	

Bromophenol blue

• Blu bromofenolo • Bleu de bromophénol • Azul de bromofenol • Bromphenolblau

Synonym:
3',3'',5',5''-Tetrabromophenolsulfonephthalein

$C_{19}H_{10}Br_4O_5S$
Molecular Weight: 669,96
CAS: 115-39-9
EEC-N: 204-086-2

Bromophenol blue > RPE - For analysis - ACS

RPE

Description ..... Pink powder
Identification ..... Positive
Appearance of solution ..... Conform
Colour change ..... yellow blue
pH range ..... 3.0 - 4.6

Code	Size	Packaging	Notes
428658	5 g	Glass bottle	
428659	25 g	Glass bottle	
428653	50 g	Glass bottle	
428655	500 g	Plastic bottle	

Clark indicator series. Dye for microscopy (histology)

Bromophenol blue solution 0.4% in ethanol

• Blu bromofenolo soluzione 0.4% in alcole etilico • Bleu de bromophénol solution 0.4% dans l'éthanol • Azul de bromofenol solución 0.4% en alcohol etilico • Bromphenolblau Lösung 0.4% in Ethanol

$C_{19}H_{10}Br_4O_5S$
Molecular Weight: 669,96
CAS: 115-39-9

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Bromophenol blue solution 0.4% in ethanol > RPE - For analysis

RPE

Description ..... Red liquid
Identification ..... Positive
Colour change ..... yellow blue
pH range ..... 3.0 - 4.6

Code	Size	Packaging	Notes
E428665	250 ml	Glass bottle	

Clark indicator series. Acid-base indicator (pH 3.0 to 4.6) indicator absorbance

Bromophenol blue solution 0.02%

• Blu bromofenolo soluzione 0.02% • Bleu de bromophénol solution 0.02% • Azul de bromofenol solución 0.02% • Bromphenolblau 0.02%

$C_{19}H_{10}Br_4O_5S$
Molecular Weight: 669,96
CAS: 115-39-9

Bromophenol blue solution 0.02% > RPE - For analysis

RPE

Description ..... Clear purple liquid Identification ..... Positive

Code	Size	Packaging	Notes
428691	100 ml	Glass bottle	

Dye for microscopy

Bromophenol blue solution

• Blu bromofenolo soluzione • Bleu de bromophénol solution • Azul de bromofenol solución • Bromphenolblau-Lösung

$C_{19}H_{10}Br_4O_5S$
Molecular Weight: 669,96
CAS: 115-39-9

Classification transport
ONU: 2733
Transport Hazard class: 3
Packing group II



Danger
H225-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

Bromophenol blue solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611012801	100 ml	Plastic bottle	Ref Ph.Eur 1012801/Color change: pH 2.8 (yellow) to pH 4.4 (blue-violet)
611012802	100 ml	Plastic bottle	Bromophenol blue solution R1 Ref Ph.Eur 1012802
611012803	100 ml	Plastic bottle	Bromophenol blue solution R2 Ref Ph.Eur 1012803

Bromophenol blue indicator

• Blu bromofenolo indicatore • Indicateur bleu de bromophénol • Indicador azul de bromofenol • Bromphenolblau-Indikator

Bromophenol blue indicator > RS - For analysis

RS

Temperature of measurement ..... 15 - 25 °C pH ..... 2.1 - 2.3 unite pH

Code	Size	Packaging	Notes
PS0269/15	1 l	Plastic bottle	

Composition: 97,3% NaOH 0.1N 2,7% H2SO4 96%, 0.8% bromophenol blue

Bromophenol blue TAC indicator

• Blu bromofenolo indicatore TAC • Indicateur TAC au bleu de bromophénol • TAC indicador azul de bromofenol • Bromphenolblau-TAC-Indikator

Bromophenol blue TAC indicator > RS - For analysis

RS

Temperature of measurement ..... 15 - 25 °C pH ..... 2 - 2.4 unite pH Colour ..... coca-cola Hazen

Code	Size	Packaging	Notes
PS0189/15	1 l	Plastic bottle	
PS0189/16	1 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Bromothymol blue

• Blu bromotimolo • Bleu de bromothymol • Azul de bromotimol • Bromthymolblau

Synonym:
3',3'-Dibromothymolsulfonphthalein

$C_{27}H_{28}Br_2O_5S$
Molecular Weight: 624,39
CAS: 76-59-5
EEC-N: 200-971-2

Warning
H302-H312-H332
P261-P264-P271-P280h-P301+P312a-P304+P340

Bromothymol blue > RPE - For analysis - ACS

RPE

Description ..... Brown powder
Identification ..... Positive
Appearance of solution ..... Conform
Colour change..... yellow-blue
pH range ..... 6.0 - 7.6

Code	Size	Packaging	Notes
428708	5 g	Glass bottle	
428702	25 g	Glass bottle	
428703	50 g	Glass bottle	

Clark indicator series

Bromothymol blue 0.4% in ethanol

• Blu bromotimolo soluzione 0.4% in alcole etilico • Bleu de bromothymol solution 0.4% dans l'éthanol • Azul de bromotimol solución 0.4% en alcohol etilico • Bromthymolblau 0.4% in Ethanol

$C_{27}H_{28}Br_2O_5S$
Molecular Weight: 624,39
CAS: 76-59-5

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Bromothymol blue 0.4% in ethanol > RPE - For analysis

RPE

Description ..... Dark green liquid
Identification ..... Positive
Colour change..... yellow blue
pH (Hydralcoholic sol.) ..... Conform 6.0 7.6

Code	Size	Packaging	Notes
E428715	250 ml	Glass bottle	

Clark indicator series

Bromothymol blue 0.02%

• Blu bromotimolo 0.02% • Bleu de bromothymol solution 0.02% • Azul de bromotimol 0.02% • Bromthymolblau 0.02%

Synonym:
3',3'-Dibromothymolsulfonphthalein

$C_{27}H_{28}Br_2O_5S$
Molecular Weight: 624,39
CAS: 76-59-5

Bromothymol blue 0.02% > RPE - For analysis

RPE

Description ..... Dark green liquid
Identification ..... Positive
pH at 20° C ..... 6.7 ÷ 6.9

Code	Size	Packaging	Notes
428731	100 ml	Glass bottle	

Dye for microscopy

Bromothymol blue solution

• Blu bromotimolo soluzione • Bleu de bromothymol solution • Azul de bromotimol solución • Bromthymolblau-Lösung

Synonym:
3',3'-Dibromothymolsulfonphthalein

$C_{27}H_{28}Br_2O_5S$
Molecular Weight: 624,39
CAS: 76-59-5

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II

Danger
H225-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

Bromothymol blue solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611012901	100 ml	Plastic bottle	Bromothymol blue solution R1 Ref Ph.Eur 1012901/ Color change: pH 5.8 (yellow) to pH 7.4 (blue)
611012903	100 ml	Plastic bottle	Bromothymol blue solution R3 Ref Ph.Eur 1012903

BSA ▶ N,O-Bis(trimethylsilyl)acetamide

BSTFA ▶ N,O-Bis(trimethylsilyl)-trifluoroacetamide

Buffer acetate pH 4.5
• Tampone acetato pH 4.5 • Tampon acétate pH 4.5 • Tampón acetato pH 4,5 • Puffer acetate pH 4.5

Buffer acetate pH 4.5 > RS - For analysis — RS

Temperature of measurement ..... 19 - 21 °C pH..... 4.3 - 4.7 unite pH

Code	Size	Packaging	Notes
PS0784/95	5 l	Kubidos	

Composition: sodium acetate anhydrous 164g/l, acetic acid 168g/l, deionized water 763g/l

Buffer pH 1
• Tampone pH 1 • Tampon pH 1 • Tampón pH 1 • Puffer pH 1

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Buffer pH 1 > RPE - For analysis — RPE

pH..... 0.98 - 1.02 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
486211	500 ml	Plastic bottle	

Composition: Glycocolle/ Sodium Chloride/Hydrochloric acid. Standardized against NIST reference materials

Buffer pH 1 > RPE - NORMEX - For analysis — RPE

Code	Size	Packaging	Notes
486221		Plastic ampoule	To dilute to 500 ml

Composition: Glycocolle/ Sodium Chloride/Hydrochloric acid. Volumetric concentrated solution to be diluted accordingly with the inst ructions for use printed on the package

Buffer pH 1.68
• Tampone pH 1.68 • Tampon pH 1.68 • Tampón pH 1.68 • Puffer pH 1.68

HEU210

Buffer pH 1.68 > RS - For analysis according to Ph. Eur. Chap. 2.2.3 — RS

Code	Size	Packaging	Notes
612203168	500 ml	Plastic bottle	Certified Reference Material

Composition: Potassium tetraoxalate 0.05 M

Buffer pH 1.68 > RPE - For analysis — RPE

Clear,colourless solution ..... Conform pH..... 1.66 - 1.70 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
486751	500 ml	Plastic bottle	

Composition: Potassium oxalate tetra-acid. Standardized against NIST reference materials

Product specifications are subject to changes. Please visit our website for updates.

Buffer pH 2
• Tampone pH 2 • Tampon pH 2 • Tampón pH 2 • Puffer pH 2

Buffer pH 2 > RS - For analysis according to Ph. Eur. Chap. 4.1.3 — RS

Code	Size	Packaging	Notes
614000200	1 l	Plastic bottle	Ref Ph.Eur 4000200

Buffer pH 2 > RPE - For analysis — RPE

pH ..... 1.98 - 2.02 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
486231	500 ml	Plastic bottle	

Composition: Citric acid/Hydrochloric acid/Sodium hydroxide. Standardized against NIST reference materials

Buffer pH 2 > RPE - NORMEX - For analysis — RPE

Description ..... Clear colourless liquid Identification ..... Positive pH at 20° C ..... 1.95 ÷ 2.05

Code	Size	Packaging	Notes
486241		Plastic ampoule	To dilute to 500 ml

Composition: Citric acid/Hydrochloric acid/Sodium hydroxide. Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Buffer pH 3
• Tampone pH 3 • Tampon pH 3 • Tampón pH 3 • Puffer pH 3



Warning
H290-H319
P234-P264-P280i-P305+P351+P338-P337+P313-P406

Buffer pH 3 > RPE - For analysis — RPE

pH ..... 2.98 - 3.02 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
486251	500 ml	Plastic bottle	
486252	1 l	Plastic bottle	

Composition: Citric acid/Hydrochloric acid/Sodium hydroxide. Standardized against NIST reference materials

Buffer pH 3 > RPE - NORMEX - For analysis — RPE

Description ..... Clear colourless liquid Identification ..... Positive pH at 20° C ..... 2.95 ÷ 3.05

Code	Size	Packaging	Notes
486261		Plastic ampoule	To dilute to 500 ml

Composition: Citric acid/Hydrochloric acid/Sodium hydroxide. Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Buffer pH 3.5
• Tampone pH 3.5 • Tampon pH 3.5 • Tampón pH 3.5 • Puffer pH 3.5

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group II



Warning
H290-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Buffer pH 3.5 > RS - For analysis according to Ph. Eur. Chap. 4.1.3 — RS

Code	Size	Packaging	Notes
614000601	250 ml	Plastic bottle	Ref Ph.Eur 4000600
614000600	1 l	Plastic bottle	Ref Ph.Eur 4000600

Buffer pH 3.56
• Tampone pH 3.56 • Tampon pH 3.56 • Tampón pH 3.56 • Puffer pH 3.56

Buffer pH 3.56 > RPE - For analysis — RPE

pH..... 3.54 - 3.58 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
486741	500 ml	Plastic bottle	

Composition: Potassium tartrate acide. Standardized against NIST reference materials

Buffer pH 3.7
• Tampone pH 3.7 • Tampon pH 3.7 • Tampón pH 3.7 • Puffer pH 3.7

Buffer pH 3.7 > RS - For analysis according to Ph. Eur. Chap. 4.1.3 — RS

Code	Size	Packaging	Notes
614000900	1 l	Plastic bottle	Ref Ph.Eur 4000900

Buffer pH 4
• Tampone pH 4 • Tampon pH 4 • Tampón pH 4 • Puffer pH 4

Buffer pH 4 > RS - For analysis according to Ph. Eur. Chap. 2.2.3 — RS

Code	Size	Packaging	Notes
612203401	500 ml	Plastic bottle	pH 4.01 at 25°C - Certified Reference Material

Composition: Potassium hydrogen phthalate 0.05 M

Buffer pH 4 > RPE - For analysis — RPE

pH..... 3.98 - 4.02 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
486271	500 ml	Plastic bottle	
486273	1 l	Plastic bottle	
486274	5 l	Kubidos	
486276	10 l	Kubidos	

Composition: Citric acid/Hydrochloric acid/Sodium hydroxide. Standardized against NIST reference materials

Buffer pH 4 > RPE - For analysis - Colored solution — RPE

pH..... 3.98 - 4.02 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
486761	500 ml	Plastic bottle	Color: Red
486762	1 l	Plastic bottle	Color: Red

Composition: Potassium phtalate acid. Standardized against NIST reference materials

Buffer pH 4 > RPE - NORMEX - For analysis — RPE

Description ..... Clear colourless liquid Identification ..... Positive pH at 20° C ..... 3.95 ÷ 4.05

Code	Size	Packaging	Notes
486281		Plastic ampoule	To dilute to 500 ml

Composition: Citric acid/Hydrochloric acid/Sodium hydroxide. Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Product specifications are subject to changes. Please visit our website for updates.

Buffer pH 4 > RPE - NORMEX - For analysis - Colored solution

RPE

Description ..... Red clear liquid Identification ..... Positive pH at 20° C ..... 3.95 ÷ 4.05

Code	Size	Packaging	Notes
486291		Plastic ampoule	Color: red - To dilute to 500 mL

Composition: Potassium phtalate acid. Traceable to NIST. Volumetric concentrated solution to be diluted accordingly with the instruc tions for use printed on the package

Buffer pH 4.5

• Tampone pH 4.5 • Tampon pH 4.5 • Tampón pH 4.5 • Puffer pH 4.5

Buffer pH 4.5 > RS - For HPLC

RS

pH ..... 4.45 - 4.55 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
524971	10 l	Plastic tank	

Buffer pH 4.62

• Tampone pH 4.62 • Tampon pH 4.62 • Tampón pH 4.62 • Puffer pH 4.62

HEU210

Buffer pH 4.62 > RPE - For analysis

RPE

pH ..... 4.60 - 4.64 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
486841	500 ml	Plastic bottle	

Composition: Sodium acetate / Acetic acid. Standardized against NIST reference materials

Buffer pH 5

• Tampone pH 5 • Tampon pH 5 • Tampón pH 5 • Puffer pH 5

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Buffer pH 5 > RPE - For analysis

RPE

pH ..... 4.98 - 5.02 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
486311	500 ml	Plastic bottle	

Composition: Citric acid/Sodium hydroxide. Standardized against NIST reference materials

Buffer pH 5 > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive pH at 20° C ..... 4.95 ÷ 5.05

Code	Size	Packaging	Notes
486301		Plastic ampoule	To dilute to 500 ml

Composition: Citric acid/Sodium hydroxide. Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Buffer pH 5.2

• Tampone pH 5.2 • Tampon pH 5.2 • Tampón pH 5.2 • Puffer pH 5.2

Buffer pH 5.2 > RS - For analysis according to Ph. Eur. Chap. 4.1.3

RS

Code	Size	Packaging	Notes
614001700	1 l	Plastic bottle	Ref Ph.Eur 4001700

Buffer pH 6
• Tampone pH 6 • Tampon pH 6 • Tampón pH 6 • Puffer pH 6

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Buffer pH 6 > RPE - For analysis — RPE

pH ..... 5.98 - 6.02 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
486331	500 ml	Plastic bottle	

Composition: Citric acid/Sodium hydroxide. Standardized against NIST reference materials

Buffer pH 6 > RPE - NORMEX - For analysis — RPE

Description ..... Clear colourless liquid Identification ..... Positive pH at 20° C ..... 5.95 ÷ 6.05

Code	Size	Packaging	Notes
486321		Plastic ampoule	To dilute to 500 ml

Composition: Citric acid/Sodium hydroxide. Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Buffer pH 6.8
• Tampone pH 6.8 • Tampon pH 6.8 • Tampón pH 6.8 • Puffer pH 6.8

Buffer pH 6.8 > RPE - NORMEX - For analysis — RPE

Description ..... Clear colourless liquid Identification ..... Positive pH at 20° C ..... 6.75 ÷ 6.85

Code	Size	Packaging	Notes
486401		Plastic ampoule	To dilute to 500 ml

Composition: Potassium phosphate monobasic/ Sodium hydroxide. Volumetric concentrated solution to be diluted accordingly with the in structions for use printed on the package

Buffer pH 6.88
• Tampone pH 6.88 • Tampon pH 6.88 • Tampón pH 6.88 • Puffer pH 6.88

Buffer pH 6.88 > RS - For analysis according to Ph. Eur. Chap. 2.2.3 — RS

Code	Size	Packaging	Notes
612203687	500 ml	Plastic bottle	pH 6.87 at 25°C - Certified Reference Material

Composition: Potassium dihydrogen phosphate 0.025 M + disodium hydrogen phosphate 0.025 M

Buffer pH 6.88 > RPE - For analysis — RPE

pH ..... 6.86 - 6.90 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
486871	500 ml	Plastic bottle	

Composition: Potassium phosphate monobasic / sodium phosphate dibasic. Standardized against NIST reference materials

Buffer pH 7
• Tampone pH 7 • Tampon pH 7 • Tampón pH 7 • Puffer pH 7

Buffer pH 7 > RS - For analysis according to Ph. Eur. Chap. 4.1.3 — RS

Code	Size	Packaging	Notes
614003500	1 l	Plastic bottle	Ref Ph.Eur 4003500

DASIT GROUP | CARLO ERBA REAGENTS
Product specifications are subject to changes. Please visit our website for updates.

Buffer pH 7 > RPE - For analysis

RPE

pH ..... 6.98 - 7.02 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
486451	500 ml	Plastic bottle	
486453	1 l	Plastic bottle	
486454	5 l	Kubidos	
486456	10 l	Kubidos	
486455	25 l	Plastic tank	

Composition: Potassium phosphate monobasic / sodium phosphate dibasic. Standardized against NIST reference materials

Buffer pH 7 > RPE - For analysis - Colored solution

RPE

Appearance ..... Green clear solution pH ..... 6.98 - 7.02 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
486791	500 ml	Plastic bottle	Color: Green
486792	1 l	Plastic bottle	Color: Green

Composition: Potassium phosphate monobasic/Sodium phosphate dibasic / Color: Green Standardized against NIST reference materials

Buffer pH 7 > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive pH at 20° C ..... 6.95 ÷ 7.05

Code	Size	Packaging	Notes
486421		Plastic ampoule	To dilute to 500 ml

Composition: Potassium phosphate monobasic / sodium phosphate dibasic. Volumetric concentrated solution to be diluted accordingly wi th the instructions for use printed on the package

Buffer pH 7 > RPE - NORMEX - For analysis - Colored solution

RPE

Description ..... Yellow clear liquid Identification ..... Positive pH at 20° C ..... 6.95 ÷ 7.05

Code	Size	Packaging	Notes
486431		Plastic ampoule	Color: yellow - To dilute to 500 mL

Composition: Potassium phosphate monobasic / sodium phosphate dibasic. Volumetric concentrated solution to be diluted accordingly wi th the instructions for use printed on the package

Buffer pH 7.2

• Tampone pH 7.2 • Tampon pH 7.2 • Tampón pH 7.2 • Puffer pH 7.2

Buffer pH 7.2 > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive pH at 20° C ..... 7.15 ÷ 7.25

Code	Size	Packaging	Notes
486441		Plastic ampoule	To dilute to 500 ml

Composition: Potassium phosphate monobasic / sodium hydroxide. Volumetric concentrated solution to be diluted accordingly with the i nstructions for use printed on the package

Buffer pH 7.20 Weise

• Tampone pH 7.20 secondo Weise • Tampon pH 7.20 selon Weise • Tampón pH 7.20 según Weise • Puffer pH 7.2 Weise

Buffer pH 7.20 Weise > RS - For analysis

RS

Description ..... Clear colourless liquid Identification ..... Conform pH at 20° C ..... 7.2 ± 0.05

Code	Size	Packaging	Notes
486411	500 ml	Plastic bottle	

Buffer pH 7.4
• Tampone pH 7.4 • Tampon pH 7.4 • Tampón pH 7.4 • Puffer pH 7.4

Buffer pH 7.4 > RS - For analysis according to Ph. Eur. Chap. 2.2.3 — RS

Code	Size	Packaging	Notes
612203741	500 ml	Plastic bottle	pH 7.41 at 25°C - Certified Reference Material

Composition: Potassium dihydrogen phosphate 0.0087 M + disodium hydrogen phosphate 0.0303 M

Buffer pH 7.4 > RS - For analysis according to Ph. Eur. Chap. 4.1.3 — RS

Code	Size	Packaging	Notes
614004600	1 l	Plastic bottle	Ref Ph.Eur 4004600

Buffer pH 7.4 > RPE - NORMEX - For analysis — RPE

Description ..... Clear colourless liquid Identification ..... Positive pH at 20° C ..... 7.35 ÷ 7.45

Code	Size	Packaging	Notes
486461		Plastic ampoule	To dilute to 500 ml

Composition: Potassium phosphate monobasic / sodium hydroxide. Volumetric concentrated solution to be diluted accordingly with the i nstructions for use printed on the package

Buffer pH 8
• Tampone pH 8 • Tampon pH 8 • Tampón pH 8 • Puffer pH 8



Warning
H290
P234-P390-P406

Buffer pH 8 > RPE - For analysis — RPE

pH ..... 7.98 - 8.02 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
486541	500 ml	Plastic bottle	
486542	1 l	Plastic bottle	

Composition: Boric acid/Sodium hydroxide/Hydrochloric acid. Standardized against NIST reference materials

Buffer pH 8 > RPE - NORMEX - For analysis — RPE

Description ..... Clear colourless liquid Identification ..... Positive pH at 20° C ..... 7.95 ÷ 8.05

Code	Size	Packaging	Notes
486531		Plastic ampoule	To dilute to 500 ml

Composition: Boric acid/Sodium hydroxide/Hydrochloric acid. Volumetric concentrated solution to be diluted accordingly with the inst ructions for use printed on the package

Buffer pH 9
• Tampone pH 9 • Tampon pH 9 • Tampón pH 9 • Puffer pH 9

HEU210

Buffer pH 9 > RS - For analysis according to Ph. Eur. Chap. 4.1.3 — RS

Code	Size	Packaging	Notes
614007000	1 l	Plastic bottle	Ref Ph.Eur 4000700

Product specifications are subject to changes. Please visit our website for updates.

Buffer pH 9 > RPE - For analysis

RPE

pH..... 8.98 - 9.02 unité pH Temperature of measurement..... 19 - 21 °C

Code	Size	Packaging	Notes
486591	500 ml	Plastic bottle	
486593	1 l	Plastic bottle	
486594	5 l	Kubidos	

Composition: Boric acid/Sodium hydroxide/Potassium chloride. Standardized against NIST reference materials

Buffer pH 9 > RPE - For analysis - Colored solution

RPE

pH..... 8.98 - 9.02 unité pH Temperature of measurement..... 19 - 21 °C

Code	Size	Packaging	Notes
PS0427/19	500 ml	Plastic bottle	Color: Blue

Composition: Boric acid/Sodium hydroxide/Potassium chloride/Methylene blue. Traceable to NIST

Buffer pH 9 > RPE - NORMEX - For analysis

RPE

Description..... Clear colourless liquid Identification..... Positive pH at 20° C..... 8.95 ÷ 9.05

Code	Size	Packaging	Notes
486571		Plastic ampoule	To dilute to 500 ml

Composition: Boric acid/Sodium hydroxide/Potassium chloride. Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Buffer pH 9.22

• Tampone pH 9.22 • Tampon pH 9.22 • Tampón pH 9.22 • Puffer pH 9.22

HEU210

Buffer pH 9.22 > RS - For analysis according to Ph. Eur. Chap. 2.2.3

RS

Code	Size	Packaging	Notes
612203918	500 ml	Plastic bottle	pH 9.18 at 25°C - Certified Reference Material

Composition: Disodium tetraborate 0.01 M

Buffer pH 9.22 > RPE - For analysis

RPE

pH..... 9.20 - 9.24 unité pH Temperature of measurement..... 19 - 21 °C

Code	Size	Packaging	Notes
486881	500 ml	Plastic bottle	

Composition: Sodium tetraborate. Standardized against NIST reference materials

Buffer pH 10

• Tampone pH 10 • Tampon pH 10 • Tampón pH 10 • Puffer pH 10



Danger
H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Buffer pH 10 > RS - For analysis according to AFNOR T90-003 normative

RS

pH..... 9 - 11 unite pH Temperature of measurement..... 19 - 21 °C

Code	Size	Packaging	Notes
PS0200/15	1 l	Plastic bottle	

Determination of the total concentration of Ca and Mg. Composition: ammonium chloride 64.5 g / l, ammonia 28% 440g / l EDTA-Mg 4.8 g / l deionized water 461.5 g / l

Buffer pH 10 > RPE - For analysis

RPE

Clear,colourless solution ..... Conform pH ..... 9.95 - 10.05 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
486611	500 ml	Plastic bottle	
486613	1 l	Plastic bottle	
486614	5 l	Kubidos	
486615	10 l	Kubidos	

Composition: Boric acid/Sodium hydroxide/Potassium chloride. Standardized against NIST reference materials

Buffer pH 10 > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive pH at 20° C ..... 9.95 ÷ 10.05

Code	Size	Packaging	Notes
486601		Plastic ampoule	To dilute to 500 ml

Composition: Boric acid/Sodium hydroxide/Potassium chloride. Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Buffer pH 10.06

• Tampone pH 10.06 • Tampon pH 10.06 • Tampón pH 10.06 • Puffer pH 10.06

HEU210

Buffer pH 10.06 > RPE - For analysis - Colored solution

RPE

Appearance ..... Blue clear solution pH ..... 10.01 - 10.11 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
486811	500 ml	Plastic bottle	Color: Blue
486812	1 l	Plastic bottle	Color: Blue

Composition: Boric acid/Sodium hydroxide/Potassium chloride. Standardized against NIST reference materials

Buffer pH 10.06 > RPE - NORMEX - For analysis - Colored solution

RPE

Description ..... Blue clear liquid Identification ..... Positive pH at 20° C ..... 10.01 ÷ 10.11

Code	Size	Packaging	Notes
486581		Plastic ampoule	Color: blue - To dilute to 500 mL

Composition: Sodium carbonate/Sodium bicarbonate. Volumetric concentrated solution to be diluted accordingly with the instructions f or use printed on the package

Buffer pH 11

• Tampone pH 11 • Tampon pH 11 • Tampón pH 11 • Puffer pH 11

Buffer pH 11 > RPE - For analysis

RPE

Appearance ..... Clear colourless solution pH ..... 10.98 - 11.02 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
486771	500 ml	Plastic bottle	
486772	1 l	Plastic bottle	

Composition: Sodium phosphate dibasic / Sodium hydroxide. Standardized against NIST reference materials

Buffer pH 11 > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive pH at 20° C ..... 10.95 ÷ 11.05

Code	Size	Packaging	Notes
486631		Plastic ampoule	To dilute to 500 ml

Composition: Boric acid/Sodium hydroxide/Potassium chloride. Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

DASITGROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

Buffer pH 12
• Tampone pH 12 • Tampon pH 12 • Tampón pH 12 • Puffer pH 12

Buffer pH 12 > RPE - For analysis
RPE

pH ..... 11.95 - 12.05 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
486691	500 ml	Plastic bottle	

Composition: Sodium phosphate/Sodium hydroxide. Standardized against NIST reference materials

Buffer pH 12 > RPE - NORMEX - For analysis
RPE

Description ..... Clear colourless liquid Identification ..... Positive pH at 20° C ..... 11.95 ÷ 12.05

Code	Size	Packaging	Notes
486621		Plastic ampoule	To dilute to 500 ml

Composition: Sodium phosphate/Sodium hydroxide. Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Buffer pH 13
• Tampone pH 13 • Tampon pH 13 • Tampón pH 13 • Puffer pH 13

Buffer pH 13 > RPE - For analysis
RPE

Clear,colourless solution ..... Conform pH ..... 12.95 - 13.05 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
486701	500 ml	Plastic bottle	
486702	1 l	Plastic bottle	

Composition: Potassium Chloride/Sodium hydroxide. Standardized against NIST reference materials

Buffer pH 13 > RPE - NORMEX - For analysis
RPE

Description ..... Clear colourless liquid Identification ..... Positive pH at 20° C ..... 12.95 ÷ 13.05

Code	Size	Packaging	Notes
486641		Plastic ampoule	To dilute to 500 ml

Composition: Potassium Chloride/Sodium hydroxide. Volumetric concentrated solution to be diluted accordingly with the instructions f or use printed on the package

Butanedioic acid ▶ Succinic acid

1-Butanesulfonic acid sodium salt
• Acido 1-butansolfonico sale sodico • Acide 1-butanesulfonique sel sodique
• Acido 1-butanosulfónico sal sódica • 1-Butansulfonsäure-Natriumsalz

Synonym:
• *Sodium 1-butanesulfonate*

$CH_3(CH_2)_3SO_3Na$
Molecular Weight: 160,17
CAS: 2386-54-1

1-Butanesulfonic acid sodium salt > RS - For ion pair chromatography
RS

Description ..... White powder Assay ..... ≥ 99.0 % At 210 nm ..... ≤ 1.0 AU At 250 nm ..... ≤ 0.04 AU
Water (K.F.) ..... ≤ 1.0 % Absorbance (5.5% in water) At 220 nm ..... ≤ 0.1 AU At 280 nm ..... ≤ 0.04 AU

Code	Size	Packaging	Notes
405631	25 g	Glass bottle	
405632	100 g	Plastic bottle	

Butanol-1

• 1-Butanolo • Butanol-1 • Butanol-1 • Butanol-1

Synonym:
Butyl alcohol

$CH_3(CH_2)_2CH_2OH$
Molecular Weight: 74,12
CAS: 71-36-3
EEC-N: 200-751-6

Classification transport
ONU: 1120
Transport Hazard class: 3
Packing group III

Danger
H226-H302-H315-H318-H335-H336
P210-P280-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P330-P362+P364-P403+P233

Butanol-1 > RS - For HPLC - Isocratic Grade

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.809 ÷ 0.811
Refractive index at 20°C . 1.3972 ÷ 1.4012
Boiling point ..... 117.0 ÷ 118.0 °C
Residue on evaporation ..... ≤5 ppm
Water (K.F.) ..... ≤0.1 %
Acidity or alkalinity ..... ≤0.0005 meq/g
Assay (GLC) ..... ≥99.8 %
U.V. Transmittance
at 210 nm ..... ≥10 %
at 235 nm ..... ≥80 %
at 310 nm ..... ≥98 %

Code	Size	Packaging	Notes
412511000	1 l	Glass bottle	
412512000	2.5 l	Glass bottle	

Filtered through 0.2 µm membrane

Butanol-1 > RS - Anhydrous - For analysis

RS

Refractive index at 20°C ..... 1.397 - 1.401
Water content (K.F.) ..... ≤ 400 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
Colour ..... ≤ 10 Hazen
Free acid (as CH3COOH) ..... ≤ 20 mg/Kg
Assay (GC) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
P0171010	200 ml	Bottle with septum	
P0171016	1 l	Glass bottle	

Butanol-1 > RPE - For analysis - ISO

RPE

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification (I.R.) ..... Conform
Water solubility ..... Conform
Density at 20° C ..... 0.808 ÷ 0.810
Boiling point ..... 116 ÷ 119 ° C
Water (K.F.) ..... ≤ 0.1 %
Residue on evaporation ..... ≤ 10 ppm
Acidity ..... ≤ 0.0008 meq/g
Alcalinity (NaOH) ..... ≤ 5 ppm
Carbonyl Compounds (CO) ..... ≤ 200 ppm
Acetone ..... ≤ 0.01 %
Isobutanolo ..... ≤ 0.15 %
2-Butanolo ..... ≤ 0.05 %
di-Butiletere ..... ≤ 0.1 %
Aldeide butirrica ..... ≤ 0.02 %
Ready carbonizable substances ..... Conform
Al ..... ≤ 0.5 ppm
Ba ..... ≤ 0.1 ppm
Ca ..... ≤ 0.5 ppm
Cd ..... ≤ 0.05 ppm
Co ..... ≤ 0.02 ppm
Cr ..... ≤ 0.02 ppm
Cu ..... ≤ 0.02 ppm
Fe ..... ≤ 0.1 ppm
Mg ..... ≤ 0.1 ppm
Mn ..... ≤ 0.02 ppm
Ni ..... ≤ 0.02 ppm
Pb ..... ≤ 0.1 ppm
Zn ..... ≤ 0.1 ppm
Assay (GLC) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
414131	1 l	Glass bottle	
414133	2.5 l	Glass bottle	
414132	22 kg	Metal drum	

Butanol-1 > RE - Pure

RE

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Identification ..... Positive
Density at 20° C ..... 0.808 ÷ 0.812
Refractive index at 20°C ..... 1.397 - 1.401
Boiling point ..... 116.8 ÷ 118.3 °C
Water (K.F.) ..... ≤0.1 %
Acidity (butirric acid) ..... ≤ 50 ppm
Residue on evaporation ..... ≤ 50 ppm
Acidity ..... ≤ 0.0008 meq/g
Assay (GLC) ..... ≥99 %

Code	Size	Packaging	Notes
308251	1 l	Glass bottle	
528300	5 l	Plastic tank	
308257	22 kg	Metal drum	
528301	25 l	Metal drum	
308259	160 kg	Metal drum	

Product specifications are subject to changes. Please visit our website for updates.

Butanol-2

• 2-Butanolo • Butanol-2 • Butan-2-olo • Butanol-2

Synonym:
sec-Butyl alcohol

$CH_3CHOHCH_2CH_3$
Molecular Weight: 74,12
CAS: 78-92-2
EEC-N: 201-158-5

Classification transport
ONU: 1120
Transport Hazard class: 3
Packing group III



Warning
H226-H319-H335-H336
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Butanol-2 > RPE - For analysis - Reag. Ph. Eur.

RPE

Description	Clear liquid	Density at 20° C	0.801 ÷ 0.811	Residue on evaporation	≤20 ppm	Assay (GLC)	≥99.0 %
Colour (APHA)	≤10	Refractive index at 20°C	1.3944 ÷ 1.3984	Acidity (butirric acid)	≤20 ppm		
Identification	Positive	Boiling point	99 - 100 °C	Alcalinity (NaOH)	≤10 ppm		
Water solubility	Conform	Water (K.F.)	≤0.2 %	Indole base	≤0.1 ppm		

Code	Size	Packaging	Notes
414264	1 l	Glass bottle	
414266	15 kg	Metal drum	
414261	160 kg	Metal drum	

Butanol-2 > RE - Pure - For synthesis

RE

Refractive index at 20°C	1.395 - 1.399	Non volatile residue	≤ 50 mg/Kg	Assay (GC)	≥ 99 %
Water content (K.F.)	≤ 2000 mg/Kg	Colour	≤ 10 Hazen	Free acid (as CH3COOH)	≤ 20 mg/Kg

Code	Size	Packaging	Notes
P0180241	10 l	Plastic tank	

tert-Butanol

• Alcole ter-butilico • tert-Butanol • Alcohol ter-butílico • tert-Butanol

Synonym:
• *2-Methyl-2-propanol*
• *tert-Butyl alcohol*

$(CH_3)_3COH$
Molecular Weight: 74,12
CAS: 75-65-0
EEC-N: 200-889-7

Classification transport
ONU: 1120
Transport Hazard class: 3
Packing group II



Danger
H225-H332-H319-H335
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

tert-Butanol > RS - Anhydrous - For analysis

RS

Refractive index at 25°C	1.383 - 1.387	Water content (K.F.)	≤ 200 mg/Kg	Assay (GC)	≥ 99.7 %
Colour	≤ 10 Hazen	Non volatile residue	≤ 10 mg/Kg	2- Propanol	≤ 0.30 %

Code	Size	Packaging	Notes
P0191016	1 l	Glass bottle	

tert-Butanol > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description	Clear colourless liq. or solid	Water solubility	Conform	Melting point	25 ÷ 26 ° C	Acidity	≤ 0.001 meq/g
Colour (APHA)	≤10	Alcohol miscibility	Conform	Water (K.F.)	≤0.1 %	Assay (GLC)	≥99.5 %
Identification (I.R.)	Conform	Boiling point	81.7 ÷ 82.7 ° C	Residue on evaporation	≤30 ppm	Carbonyl compounds (as HCHO)	≤ 0.01 %

Code	Size	Packaging	Notes
414343	500 ml	Plastic bottle	
414341	1 l	Glass bottle	
414346	2.5 l	Plastic bottle	
414342	25 l	Plastic tank	

tert-Butanol > RE - Pure

RE

Refractive index at 25°C	1.383 - 1.387	Water content (K.F.)	≤ 800 mg/Kg	Assay (GC)	≥ 99.5 %
Colour	≤ 10 Hazen	Non volatile residue	≤ 20 mg/Kg		

Code	Size	Packaging	Notes
P0190222	5 l	Plastic tank	
P0190268	200 l	Metal drum	

2-Butanone ▶ Ethyl methyl ketone

2-Butoxy ethanol

• 2-Butossietanolo • 2-Butoxyéthanol • 2-Butoxietanol • 2-Butoxyethanol

Synonym:
• *Ethylene glycol butyl ether*
• *Butyl glycol*

$CH_3(CH_2)_3OCH_2CH_2OH$
Molecular Weight: 118,18
CAS: 111-76-2
EEC-N: 203-905-0

Warning
H302-H312-H332-H315-H319
P261-P271-P304+P340-P305+P351+P338-P332+P313-P337+P313

2-Butoxy ethanol > RPE - For analysis

RPE

Description ..... Clear colourless liquid	Density at 20° C .....0.897 ÷ 0.905	Boiling point.....167 ÷ 172 ° C	Residue on evaporation .....≤50 ppm
Identification .....Positive	Refractive index at 20°C.1.4167 ÷ 1.4207	Water (K.F.) ..... ≤ 0.1 %	Assay (GLC) ..... ≥ 99 %

Code	Size	Packaging	Notes
453941	1 l	Glass bottle	

2-(2-Butoxyethoxy)ethanol

• 2-(2-Butossietossi)etanolo • 2-(2-Butoxyéthoxy)éthanol • 2-(2-Butoxietoxi) etanol
• 2-(2-Butoxyethoxy)ethanol

Synonym:
• *Diethylene glycol butyl ether*
• *BDG*

$C_8H_{18}O_3$
Molecular Weight: 162,23
CAS: 112-34-5
EEC-N: 203-961-6

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

2-(2-Butoxyethoxy)ethanol > RPE - For analysis

RPE

Description ..... Clear colourless liquid	Refractive index at 20°C.1.4296 ÷ 1.4346	Alcalinity (NH3).....≤0.85 ppm	Residue on ignition.....≤20 ppm
Identification .....Positive	Boiling point.....230.0 ÷ 232.0 °C	Carbonyl Compounds (CO) .....≤500 ppm	Fe .....≤2 ppm
Water miscibility.....Conform	Water (K.F.) ..... ≤0.2 %	Heavy metals (Pb).....≤2 ppm	Assay (GLC) ..... ≥99 %
Density at 20° C .....0.951 ÷ 0.959	Acidity (acetic acid).....≤300 ppm	Peroxides (H2O2) .....≤250 ppm	

Code	Size	Packaging	Notes
453881	1 l	Glass bottle	
453883	25 kg	Metal drum	

n-Butyl acetate

• n-Butile acetato • n-Butyle acétate • n-Butilo acetate • n-Butylacetat

$CH_3COO(CH_2)_3CH_3$
Molecular Weight: 116,16
CAS: 123-86-4
EEC-N: 204-658-1

Classification transport
ONU: 1123
Transport Hazard class: 3
Packing group III

Warning
H226-H336-HEU066
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

n-Butyl acetate > RS - Anhydrous - For analysis

RS

Refractive index at 20°C..... 1.392 - 1.396	Non volatile residue .....≤ 10 mg/Kg	Assay (GC).....≥ 99.5 %
Water content (K.F.).....≤ 100 mg/Kg	Colour .....≤ 10 Hazen	Free acid (as CH3COOH).....≤ 50 mg/Kg

Code	Size	Packaging	Notes
P0011016	1 l	Glass bottle	
P0011021	2.5 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

n-Butyl acetate > RPE - For analysis

RPE

Description	Clear colourless liquid	Water (K.F.)	≤0.1 %	Total silicon	≤0.02 ppm	Na	≤0.2 ppm
Identification (I.R.)	Conform	Butan-1-ol	≤0.5 %	Total sulphur	≤0.5 ppm	Pb	≤0.05 ppm
Ready carbonizable substances	Conform	n-Butyl formate	≤0.1 %	Ca	≤0.2 ppm	Zn	≤0.1 ppm
Density at 20° C	0.878 ÷ 0.884	Butyl propanoate	≤0.1 %	Cu	≤0.2 ppm	Assay (GLC)	≥99 %
Refractive index at 20°C	1.3926 ÷ 1.3976	iso-Butyl acetate	≤0.5 %	Fe	≤0.1 ppm		
Boiling point	126.0 ÷ 127.0 °C	Residue on evaporation	≤10 ppm	K	≤0.2 ppm		
Acidity or alkalinity	≤0.001 meq/g	Total phosphorus	≤0.2 ppm	Mg	≤0.02 ppm		

Code	Size	Packaging	Notes
431601000	1 l	Glass bottle	
431602000	2.5 l	Glass bottle	
431604	200 l	Metal drum	

n-Butyl acetate > RE - Pure

RE

Description	Clear colourless liquid	Refractive index at 20°C	1.3911 ÷ 1.3991	Residue on evaporation	≤100 ppm
Identification	Positive	Boiling point	123 ÷ 128 °C	Acidity (acetic acid)	≤300 ppm
Density at 20° C	0.876 ÷ 0.886	Water (K.F.)	≤0.1 %	Assay (GLC)	≥99 %

Code	Size	Packaging	Notes
325602	1 l	Glass bottle	
325601	2.5 l	Glass bottle	
325604	5 l	Plastic tank	
325603	24 kg	Metal drum	

n-Butyl alcohol ▶ Butanol-1

sec-Butyl alcohol ▶ Butanol-2

n-Butyl chloride

• n-Butile cloruro • n-Butyle chlorure • n-Butilo cloruro • n-Butylchlorid

Synonym:
1-Chlorobutane

$CH_3(CH_2)_3Cl$
Molecular Weight: 92,57
CAS: 109-69-3
EEC-N: 203-696-6

Classification transport
ONU: 1127
Transport Hazard class: 3
Packing group II

Danger
H225-H304-H412
P210-P241-P280-P301+P310a-P303+P361+P353-P403+P235

n-Butyl chloride > RS - For HPLC - Isocratic Grade

RS

Clear,colourless liq.appearance	Conform	Refractive index at 20°C	1.400 - 1.404	Non volatile residue	≤ 10 mg/Kg	UV transmittance at 250 nm	≥ 98 %
Identification	Conform	Water content (K.F.)	≤ 100 mg/Kg	UV transmittance at 230 nm	≥ 65 %	UV transmittance at 240 nm	≥ 95 %
Colour	≤ 10 Apha	Free acid (as HCl)	≤ 10 mg/Kg	UV transmittance at 235 nm	≥ 90 %	Assay (GC)	≥ 99.8 %

Code	Size	Packaging	Notes
431821	1 l	Glass bottle	

n-Butyl chloride > RPE - For analysis

RPE

Description	Clear colourless liquid	Refractive index at 20°C	1.3981 ÷ 1.4061	Colour	≤ 10 APHA	Water (K.F.)	≤ 150 ppm
Identification	Positive	Assay (GLC)	≥ 99.5 %	1-Butanol	≤ 0.05 %	Residue on evaporation	≤ 10 ppm

Code	Size	Packaging	Notes
431811	100 ml	Glass bottle	
431817	1 l	Glass bottle	

Butylhydroxytoluene

• Butile idrossitoluene • Butylhydroxytoluène • Butil hidroxitolueno • Butylhydroxytoluol

Synonym:
• *2,6-Di-tert-butyl-p-cresol*
• *BHT*

$[(CH_3)_3C]_2C_6H_2(CH_3)OH$
Molecular Weight: 220,36
CAS: 128-37-0
EEC-N: 204-881-4

Warning
H410
P273-P391-P501a

Butylhydroxytoluene > RPE - For analysis

RPE

Description ..... White semitransparent crystals
Identification ..... Positive
Melting point ..... 69.0 ÷ 71.0 ° C
Water (K.F.) ..... ≤ 0.12 %
Assay (GLC) ..... 99 ÷ 100 %

Code	Size	Packaging	Notes
432121	25 g	Glass bottle	

tert-Butylmethylether

• ter-Butilmetiletere • tert-Butylméthyléther • ter-Butilmetiléter • tert-Butylmethylether

Synonym:
• *MTBE*
• *Methyl tert-butyl ether*

$CH_3OC_4H_9$
Molecular Weight: 88,15
CAS: 1634-04-4
EEC-N: 216-653-1

Classification transport
ONU: 2398
Transport Hazard class: 3
Packing group II

Danger
H225-H315
P210-P241-P280-P303+P361+P353-P332+P313-P403+P235

tert-Butylmethylether > RS - For HPLC - Isocratic Grade

RS

Refractive index at 20°C ..... 1.367 - 1.371
Clear,colourless liq.appearance ..... Conform
Identification ..... Conform
Colour ..... ≤ 10 Apha
Water content (K.F.) ..... ≤ 100 mg/Kg
UV transmittance at 210 nm ..... ≥ 10 %
UV transmittance at 230 nm ..... ≥ 40 %
UV transmittance at 250 nm ..... ≥ 75 %
UV transmittance at 280 nm ..... ≥ 92 %
UV transmittance at 300 nm ..... ≥ 98 %
Assay (GC) ..... ≥ 99.8 %
Methanol + tert-butanol ..... ≤ 0.05 %
Non volatile residue ..... ≤ 10 mg/Kg
Hydrocarbons up to C8 ..... ≤ 0.05 %

Code	Size	Packaging	Notes
432031	1 l	Glass bottle	
432032	2.5 l	Glass bottle	
432034	4 l	Glass bottle	

tert-Butylmethylether > RS - For preparative HPLC

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.730 ÷ 0.750
Boiling point ..... 54.8 ÷ 55.8 ° C
Acidity or alkalinity ..... ≤0.0002 meq/g
Water (K.F.) ..... ≤200 ppm
Residue on evaporation ..... ≤2 ppm
Assay (GLC) ..... ≥99.5 %
U.V. Transmittance
at 240 nm ..... ≥60 %
at 270 nm ..... ≥95 %

Code	Size	Packaging	Notes
432022000	2.5 l	Glass bottle	

tert-Butylmethylether > RS - PESTIPUR - For pesticide analysis

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C ..... 1.367 - 1.371
Water content (K.F.) ..... ≤ 100 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 99.8 %
Non volatile residue ..... ≤ 5 mg/Kg
GC-ECD.Individual peak (Lindane) .. ≤ 3 ng/l
Retention time trichlorobenzene to mirex
Ret.range 1,2,4-trichlorobenzene
GC-NPD.Individual peak (Ethylparathion) ≤ 3 ng/l
to decachlorobiphenyle
Retention time Atrazin to Coumaphos

Code	Size	Packaging	Notes
432061	1 l	Glass bottle	
432062	2.5 l	Glass bottle	

For chlorinated and nitrogenous compounds analysis

tert-Butylmethylether > RS - SPECTROSOL - For optical spectroscopy

RS

Clear,colourless liq.appearance ..... Conform
Identification ..... Conform
Colour ..... ≤ 10 Apha
Refractive index at 20°C ..... 1.367 - 1.371
Water content (K.F.) ..... ≤ 100 mg/Kg
UV transmittance at 210 nm ..... ≥ 10 %
UV transmittance at 230 nm ..... ≥ 40 %
UV transmittance at 250 nm ..... ≥ 75 %
UV transmittance at 280 nm ..... ≥ 92 %
UV transmittance at 300 nm ..... ≥ 98 %
Methanol + tert-butanol ..... ≤ 0.05 %
Assay (GC) ..... ≥ 99.8 %
Non volatile residue ..... ≤ 10 mg/Kg
Hydrocarbons up to C8 ..... ≤ 0.05 %

Code	Size	Packaging	Notes
432001	1 l	Glass bottle	
432002	2.5 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

tert-Butylmethylether > RS - Anhydrous - For analysis

RS

Refractive index at 20°C ..... 1.367 - 1.371
Water content (K.F.) ..... ≤ 100 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 99.8 %
Methanol + tert-butanol ..... ≤ 0.05 %
Hydrocarbons up to C8 ..... ≤ 0.05 %

Code	Size	Packaging	Notes
P0921016	1 l	Glass bottle	

tert-Butylmethylether > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.739 ÷ 0.742
Boiling point ..... 53 ÷ 56 ° C
Water (K.F.) ..... ≤100 ppm
Peroxides (H2O2) ..... ≤10 ppm
Methyl alcohol ..... ≤0.1 %
Refractive index at 20°C ..... 1.368 ÷ 1.370
Assay (GLC) ..... ≥99.5 %

Code	Size	Packaging	Notes
432011	500 ml	Glass bottle	
432013	2.5 l	Glass bottle	
432015	20 kg	Aluminium can	

tert-Butylmethylether > RE - Pure

RE

Description ..... Clear liquid
Density at 20°C ..... 0.739 ÷ 0.741
Refractive index at 20°C. 1.3635 ÷ 1.3735
Boiling point ..... 54.8 ÷ 55.8 °C
Water (K.F.) ..... ≤ 300 ppm
Residue on evaporation ..... ≤ 30 ppm
Assay (GLC) ..... ≥ 99.9 %
Methanol + tert-butanol ..... ≤ 0.1 %

Code	Size	Packaging	Notes
528974	1 l	Glass bottle	
528970	5 l	Plastic tank	
528971	25 l	Metal drum	
528979	200 l	Metal drum	

n-Butyric acid

• Acido n-butirrico • Acide n-butyrique • Acido n-butírico • Buttersäure

$CH_3CH_2CH_2COOH$
Molecular Weight: 88,11
CAS: 107-92-6
EEC-N: 203-532-3

Classification transport
ONU: 2820
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

n-Butyric acid > RPE - For analysis - Reag. Ph. Eur.

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.953 ÷ 0.957
Refractive index at 20°C. 1.3961 ÷ 1.4021
Boiling point ..... 163.0 ÷ 164.0 °C
As ..... ≤ 3 ppm
Cd ..... ≤ 1 ppm
Hg ..... ≤ 1 ppm
Pb ..... ≤ 10 ppm
Assay (GLC) ..... ≥99 %

Code	Size	Packaging	Notes
403236	250 ml	Glass bottle	

Cadmium standard solution

• Cadmio standard soluzione • Cadmium standard solution • Cadmio, solución patrón • Cadmium-Standardlösung



Danger
H290-H340-H350-H373-HA26-HEU207
P260-P280-P308+P313-P314-P406-P501a

Cadmium standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615000700	100 ml	Plastic bottle	A 0.1 % solution Ref Ph.Eur 5000700
615000709	100 ml	Plastic bottle	A 10 ppm solution: to dilute according to Ref Ph.Eur 5000701

Cadmium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505547	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505548	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505549	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Cadmium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503491	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503493	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503495	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503497	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Cadmium standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
507529	100 ml	Plastic bottle	conc. 1000 ppm Matrix: Nitric acid
507483	500 ml	Plastic bottle	conc. 1000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Cadmium standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
432311		Plastic ampoule	conc. 1.000 ppm Matrix: Nitric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Cadmium standard solution > RS - Quality control standard solution for AAS (graphite furnace)

RS

Code	Size	Packaging	Notes
504360	50 ml	Plastic bottle	conc. 5 ppb - Matrix: 2% Nitric acid

Cadmium acetate dihydrate

• Cadmio acetato diidrato • Cadmium acetate dihydrate • Cadmio acetato dihidrato • Cadmiumacetat-Dihydrat

$Cd(CH_3COO)_2.2H_2O$
Molecular Weight: 266,52
CAS: 5743-04-4
EEC-N: 208-853-2

Classification transport
ONU: 2570
Transport Hazard class: 6.1
Packing group III



Danger
H301-H312-H330-H340-H350-H361d-H372-H410-HA26
P271-P280-P284-P301+P310a-P304+P340-P403+P233

Cadmium acetate dihydrate > RPE - For analysis

RPE

Description ..... White crystals	Nitrate ..... ≤ 30 ppm	Cu ..... ≤ 10 ppm	Zn ..... ≤ 50 ppm
Identification ..... Positive	Substances not ppt. H2S ..... ≤ 0.3 %	Fe ..... ≤ 10 ppm	Assay (complexometric) ..... ≥ 97.5 %
Water-insoluble matter ..... ≤ 0.005 %	Sulphate ..... ≤ 50 ppm	Pb ..... ≤ 50 ppm	

Code	Size	Packaging	Notes
432344	100 g	Glass bottle	
432345	250 g	Plastic bottle	
432347	1 kg	Plastic bottle	

Cadmium carbonate

• Cadmio carbonato • Cadmium carbonate • Cadmio carbonato • Cadmiumcarbonat

Synonym:
Otavite

$CdCO_3$
Molecular Weight: 172,42
CAS: 513-78-0
EEC-N: 208-168-9

Classification transport
ONU: 2570
Transport Hazard class: 6.1
Packing group II



Danger
H302-H312-H330-H340-H350-H372-H410-HA26
P280-P284-P304+P340-P310a-P320-P330-P362+P364-P403+P233

Cadmium carbonate > RPE - For analysis

RPE

Description ..... White powder	HCl-insoluble matter ..... ≤100 ppm	Pb ..... ≥100 ppm
Identification ..... Positive	Cu ..... ≤100 ppm	Assay (complexometric) ..... 97 ÷ 100 %

Code	Size	Packaging	Notes
432444	100 g	Glass bottle	
432446	500 g	Plastic bottle	

Cadmium chloride monohydrate

• Cadmio cloruro monoidrato • Cadmium chlorure monohydraté • Cadmio cloruro monohidrato • Cadmiumchlorid monohydrat

$CdCl_2.H_2O$
Molecular Weight: 201,32
CAS: 35658-65-2
EEC-N: 233-296-7

Classification transport
ONU: 2570
Transport Hazard class: 6.1
Packing group III



Danger
H301-H332-H350-H372-H410-HA26
P271-P280-P301+P310a-P330-P304+P340-P308+P313

Cadmium chloride monohydrate > RE - Pure

RE

Description ..... White crystals	Sulphate ..... ≤200 ppm	Pb ..... ≤50 ppm	
Identification ..... Positive	Cu ..... ≤20 ppm	Zn ..... ≤200 ppm	
Substances not ppt. H2S ..... ≤0.2 %	Fe ..... ≤10 ppm	Assay ..... ≥98 %	

Code	Size	Packaging	Notes
325741	100 g	Glass bottle	

Cadmium nitrate tetrahydrate

• Cadmio nitrato tetraidrato • Cadmium nitrate tétrahydraté • Cadmio nitrato tetrahidratado • Cadmiumnitrat-Tetrahydrat

Synonym: *Nitric acid,cadmium salt tetrahydrate*

$Cd(NO_3)_2.4H_2O$
Molecular Weight: 308,47
CAS: 10022-68-1
EEC-N: 233-710-6

Classification transport
ONU: 2570
Transport Hazard class: 6.1
Packing group III

Danger
H302-H312-H332-H340-H350-H372-H410-HA26
P271-P280-P304+P340-P308+P313-P330-P362+P364

Cadmium nitrate tetrahydrate > RPE - For analysis

RPE

Description	White crystals	Chloride	≤ 50 ppm	Pb	≤ 50 ppm
Identification	Positive	Cu	≤ 30 ppm	Sulfate	≤ 0.003 %
Water-insoluble matter	≤ 100 ppm	Fe	≤ 10 ppm	Assay (complexometric)	≥ 99 %

Code	Size	Packaging	Notes
432644	100 g	Glass bottle	
432645	500 g	Glass bottle	

Caffeine anhydrous

• Caffeina anidra • Caféine anhydre • Cafeína anhidra • Wasserfreies Koffein

Synonym: *1,3,7-Trimethylxanthine*

$C_8H_{10}N_4O_2$
Molecular Weight: 194,19
CAS: 58-08-2
EEC-N: 200-362-1

Warning
H302
P264-P270-P301+P312a-P330-P501a

Caffeine anhydrous > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP-DAB

ERBApharm

Description	White crystalline powder	Related substances	Conform Ph.Eur.	Melting point	235 ÷ 239 °C s.s.	Sulphate	≤ 500 ppm
Identification	Positive	Organic volatile impurities	Conform USP-NF	Loss on drying	≤ 0.5 %	Impurity (HPLC)	≤ 0.1 %
Appearance of solution	Conform Ph.Eur.	Ready carbonizable substances	Conform USP-NF	Sulphated ash	≤ 0.1 %	Assay (non-aqueos medium)	98.5 ÷ 101.0 % s.s.
Acidity	Conform Ph.Eur.	Other alkaloids	Conform USP-NF	Heavy metals (Pb)	≤ 10 ppm	Assay (HPLC)	98.5 ÷ 101.0 % s.s.

Code	Size	Packaging	Notes
326356	500 g	Plastic bottle	
326357	1 kg	Plastic bottle	
326358	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Calcium standard solution

• Calcio standard soluzione • Calcium standard solution • Calcio, solución patrón • Calcium-Standardlösung

Ca
CAS: 7440-70-2

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Calcium standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615000801	100 ml	Plastic bottle	A 100 ppm solution: to dilute according to Ref Ph.Eur 5000801
615000802	100 ml	Plastic bottle	A 100 ppm alcoholic solution: to dilute according to Ph.Eur 5000802
615000803	100 ml	Plastic bottle	A 10 ppm solution: to dilute according to Ref Ph.Eur 5000803
615000804	100 ml	Plastic bottle	A 100 ppm solution R1: to dilute according to Ref Ph.Eur 5000804
615000809	100 ml	Plastic bottle	A 400 ppm solution: to dilute according to Ref Ph.Eur 5000800

Product specifications are subject to changes. Please visit our website for updates.

Calcium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505542	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505545	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505543	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Calcium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503481	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503483	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503485	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503487	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Calcium standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
507530	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507476	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
E497485	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
E497481	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Calcium standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
432941		Plastic ampoule	conc. 1.000 ppm Matrix: Nitric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Calcium standard solution > RS - Standard solution for ion chromatography

RS

Code	Size	Packaging	Notes
503221	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503223	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503389	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water and nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Calcium acetate anhydrous

• Calcio acetato anidro • Calcium acétate anhydre • Calcio acetato anhidro • Calciumacetat wasserfrei

$Ca(CH_3COO)_2$
Molecular Weight: 158,17
CAS: 62-54-4
EEC-N: 200-540-9

Calcium acetate anhydrous > ERBApharm - According to pharmacopoeia: BP

ERBApharm

Description	White powder	Water (K.F.)	≤7 %	Sulphate	≤600 ppm	K	≤0.1 %
Identification	Positive	Alcalinity	Conform	Al	≤1 ppm	Na	≤0.5 %
Nitrate	Conform BP	Chloride	≤330 ppm	As	≤2 ppm	Mg	≤500 ppm
Ready oxidizable substances	Conform BP	Heavy metals (Pb)	≤20 ppm	Ba	≤50 ppm	Assay (complexometric)	98.0 ÷ 100.5 % s.s.

Code	Size	Packaging	Notes
326511	1 kg	Plastic bottle	
326512	5 kg	Plastic tank	
326513	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Calcium acetate anhydrous > RE - Pure

RE

Description	White powder	Chloride	≤ 500 ppm	Heavy metals (Pb)	≤ 5 ppm	Fe	≤ 10 ppm
Identification	Positive	Water-insoluble matter	≤ 500 ppm	Sulphate	≤ 500 ppm	Assay	≥ 99 %
pH sol. 1%	7 ÷ 8	Loss on drying	≤ 6 %	As	≤ 2 ppm		

Code	Size	Packaging	Notes
326507	1 kg	Plastic bottle	
326503	25 kg	Plastic bucket	

Calcium acetate monohydrate

• Calcio acetato monoidrato • Calcium acétate monohydrate • Calcio acetato monohidrato • Calciumacetat monohydrat

$Ca(CH_3COO)_2.H_2O$
Molecular Weight: 158,17
CAS: 5743-26-0
EEC-N: 611-528-1

Calcium acetate monohydrate > RPE - For analysis

RPE

Description	White crystalline powder	Chloride	≤ 20 ppm	Sulphate	≤ 50 ppm	Na	≤ 500 ppm
Identification	Positive	Water insoluble substances	≤ 0.01 %	Ba	≤ 30 ppm	Sr	≤ 0.1 %
pH sol. 5% at 25° C	7.5 ÷ 8.5	Heavy metals (Pb)	≤ 10 ppm	Fe	≤ 10 ppm	Assay (complexometric)	≥ 99 %
Loss on drying	≤ 10 %	Nitrate	≤ 20 ppm	K	≤ 100 ppm		

Code	Size	Packaging	Notes
432985	250 g	Plastic bottle	
432987	1 kg	Plastic bottle	
432982	25 kg	Plastic bucket	

Calcium carbonate

• Calcio carbonato • Calcium carbonate • Calcio carbonato • Calciumcarbonat

$CaCO_3$
Molecular Weight: 100,09
CAS: 471-34-1

Calcium carbonate > RS - For enviromental analysis

RS

Description	White powder	Sulphate	< 2000 ppm	As	< 4 ppm
Identification	Positive	Heavy metals (Pb)	< 20 ppm	Acid insoluble	< 0.2 %
Chloride	< 200 ppm	Perdita essiccamento (200°C)	< 0.5 %	Assay (complexometric)	> 98.5 % s.s.

Code	Size	Packaging	Notes
433216	500 g	Plastic bottle	

Low content in alkali

Product specifications are subject to changes. Please visit our website for updates.

Calcium carbonate > RS - For chromatography

RS

Description ..... White powder Identification ..... Positive

Code	Size	Packaging	Notes
433245	250 g	Plastic bottle	

Calcium carbonate > RS - TOC standard

RS

Code	Size	Packaging	Notes
505008	100 ml	Plastic bottle	conc. 50 mg/l
505009	100 ml	Plastic bottle	conc. 100 mg/l

Calcium carbonate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White powder
Identification ..... Positive
Diluted hydrochloric acid insoluble.... ≤ 100 ppm
Ammonium ..... ≤ 30 ppm
Chloride ..... ≤ 10 ppm
Fluoride ..... ≤ 15 ppm
Sulphate ..... ≤ 100 ppm
Heavy metals (Pb) ..... ≤ 10 ppm
Ba ..... ≤ 100 ppm
Fe ..... ≤ 30 ppm
K ..... ≤ 100 ppm
Mg ..... ≤ 200 ppm
Na ..... ≤ 0.1 %
Sr ..... ≤ 0.1 %
Assay (complexometric) 99.0 ÷ 100.5 % s.s.
As ..... ≤ 4 ppm

Code	Size	Packaging	Notes
433185	250 g	Plastic bottle	
433187	1 kg	Plastic bottle	
433183	25 kg	Sack	

Calcium carbonate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP

ERBApharm

Description ..... White powder
Identification ..... Positive
Barium ..... Conform Ph.Eur.
Organic volatile impurities Conform USP-NF
Loss on drying ..... ≤ 2.0 %
Magnesium and alkali metals..... ≤ 1.0 %
Substances insoluble in acetic acid ≤ 0.2 %
Sustances insoluble in hydrochloric acid .. ≤ 0.2 %
Chloride ..... ≤ 330 ppm
Fluoride ..... ≤ 50 ppm
Sulphate ..... ≤ 0.25 %
Heavy metals (Pb) ..... ≤ 20 ppm
As ..... ≤ 3 ppm
Fe ..... ≤ 200 ppm
Hg ..... ≤ 0.5 ppm
Pb ..... ≤ 3 ppm
Assay (complexometric) 98.5 ÷ 100.5 % s.s.
Origin (BSE/TSE) ..... Synthesis
Residual solvents (Current ICH) ..... Conform

Code	Size	Packaging	Notes
327101	1 kg	Plastic bottle	
327105	25 kg	Sack	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Calcium carbonate > RE - Pure

RE

Description ..... White powder
Identification ..... Positive
Loss on drying ..... ≤2 %
Assay (complexometric) ..... ≥ 98.0 %
Soluble alcalies(NaCO3) ..... ≤0.26 %
Chloride ..... ≤350 ppm
Apparent density ..... 400 ÷ 600 g/l
Heavy metals (Pb) ..... ≤30 ppm
Sulphate ..... ≤0.8 %
As ..... ≤5 ppm
Fe ..... ≤200 ppm
Assay (alkalimetric) ..... 98.5 ÷ 100.5 % s.s.

Code	Size	Packaging	Notes
327002	2.5 kg	Plastic bottle	Light powder
327059	5 kg	Plastic tank	Heavy powder
327003	25 kg	Plastic bucket	Light powder

Calcium chloride anhydrous

• Calcio cloruro anidro • Calcium chlorure anhydre • Calcio cloruro anhidro • Calciumchlorid wasserfrei

$CaCl_2$
Molecular Weight: 110,99
CAS: 10043-52-4
EEC-N: 233-140-8

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Calcium chloride anhydrous > RS - For microanalysis

RS

Description ..... White granules Identification ..... Positive Diameter ..... 1 - 6 mm Assay (complexometric) ..... ≥96 %

Code	Size	Packaging	Notes
433535	250 g	Plastic bottle	Granular

Calcium chloride anhydrous > RPE - For analysis

RPE

Description ..... White powder
Identification ..... Positive
Acidity ..... < 50 ppm
Alcalinity ..... < 0.3 %
Water-insoluble matter ..... < 0.1 %
Sulphate ..... < 500 ppm
Heavy metals (Pb) ..... < 50 ppm
Ba ..... < 100 ppm
Fe ..... < 50 ppm
Mg ..... < 0.1 %
Assay (complexometric) ..... ≥ 96.0 %

Code	Size	Packaging	Notes
433403	100 g	Plastic bottle	
433406	500 g	Plastic bottle	
433407	1 kg	Plastic bottle	
433405	25 kg	Drum	

Calcium chloride anhydrous > RE - Pure - Powder

RE

Description ..... White powder
Identification ..... Positive
Water insoluble Matter ..... ≤ 0.1 %
Heavy metals (Pb) ..... ≤ 50 ppm
Assay (complexometric) ..... ≥ 96 %

Code	Size	Packaging	Notes
328257	1 kg	Plastic bottle	
328258	5 kg	Plastic tank	
328252	25 kg	Plastic bucket	

Calcium chloride anhydrous > RE - Pure - Granular

RE

Description ..... white pellets
Identification ..... Positive
Diameter ..... 15 - 30 mm
Granulo 1,25-5 mm ..... ≥ 75 %
Assay (complexometric) ..... ≥ 90 %

Code	Size	Packaging	Notes
328757	1 kg	Plastic bottle	Ø 1.25 - 4 mm
328807	1 kg	Plastic bottle	Ø 15-30 mm
328759	5 kg	Plastic tank	Ø 1.25 - 4 mm
328809	5 kg	Plastic tank	Ø 15-30 mm

Calcium chloride dihydrate

• Calcio cloruro biidrato • Calcium chlorure dihydraté • Calcio cloruro dihidrato • Calciumchloriddihydrat

$CaCl_2.2H_2O$
Molecular Weight: 147,02
CAS: 10035-04-8
EEC-N: 233-140-8

Warning
H302-H319
P264-P280i-P301+P312a-P305+P351+P338-P337+P313-P501a

Calcium chloride dihydrate > RPE - For analysis - ACS

RPE

Description ..... White crystalline powder
Identification ..... Positive
pH sol. 5% at 25° C ..... 4.5 ÷ 8.5
Water-insoluble matter ..... ≤ 100 ppm
Heavy metals (Pb) ..... ≤ 5 ppm
Sulphate ..... ≤ 100 ppm
Ba ..... ≤ 50 ppm
Fe ..... ≤ 10 ppm
K ..... ≤ 100 ppm
Mg ..... ≤ 50 ppm
Assay (complexometric) ..... 99.0 ÷ 105.0 %
Na ..... ≤ 0.02 %
Sr ..... ≤ 0.1 %
Ammonium ..... ≤ 50 ppm
Oxidizing substances ..... ≤ 30 ppm

Code	Size	Packaging	Notes
433386	100 g	Plastic bottle	
433387	500 g	Plastic bottle	
433381	1 kg	Plastic bottle	
433382	5 kg	Plastic tank	
433384	25 kg	Plastic bucket	

Product specifications are subject to changes. Please visit our website for updates.

Calcium chloride dihydrate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP-DAB

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Acidity or alkalinity.....Conform Ph.Eur.
Ba .....Conform Ph.Eur.
Fe,Al and Phosphate ..... Conform USP-NF
pH (1:20) .....4.5 ÷ 9.2
Heavy metals (Pb).....≤ 10 ppm
Sulphate .....≤ 300 ppm
Mg and alkaline metals..... ≤ 0.5 %
Al .....≤ 1 ppm
Fe .....≤ 10 ppm
Assay (complexometric) ....99.0 ÷ 103.0 %
Origin (BSE/TSE)..... Synthesis
Residual solvents (Current ICH)..... Conform

Code	Size	Packaging	Notes
327607	1 kg	Plastic bottle	
327609	5 kg	Plastic tank	
327603	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Calcium chloride hexahydrate

• Calcio cloruro esaidrato • Calcium chlorure hexahydraté • Calcio cloruro hexahidrato • Calciumchloridhexahydrat

$CaCl_2.6H_2O$
Molecular Weight: 219,08
CAS: 7774-34-7
EEC-N: 233-140-8

Warning
H302-H319
P264-P280i-P301+P312a-P305+P351+P338-P337+P313-P501a

Calcium chloride hexahydrate > RPE - For analysis

RPE

Description .White semitransparent crystals
Identification ..... Positive
Oxidizing substances (NO3).....≤ 30 ppm
pH sol. 5% at 25° C .....4.5 ÷ 8.5
Ammonium .....≤ 50 ppm
Phosphate .....≤ 10 ppm
Water ins.ble/ppt NH4OH.....≤ 100 ppm
Heavy metals (Pb).....≤ 5 ppm
Sulphate .....≤ 100 ppm
As .....≤ 1 ppm
Ba .....≤ 50 ppm
Cu .....≤ 5 ppm
Fe .....≤ 25 ppm
Mg .....≤ 50 ppm
Ni.....≤ 5 ppm
Pb .....≤ 5 ppm
Zn .....≤ 10 ppm
Assay (complexometric) ..... ≥ 98.0 %
Mn .....≤ 5 ppm
Sr.....≤ 100 ppm

Code	Size	Packaging	Notes
433377	1 kg	Plastic bottle	
433371	5 kg	Plastic tank	
433375	25 kg	Drum	

Calcium chloride hexahydrate > ERBApharm - According to pharmacopoeia: Ph.Eur.

ERBApharm

Description ..... White crystalline mass
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Acidity or alkalinity.....Conform Ph.Eur.
Sulphate .....≤ 200 ppm
Al .....Conform Ph.Eur.
Ba .....Conform Ph.Eur.
Heavy metals (Pb).....≤ 15 ppm
Fe .....≤ 7 ppm
Mg and alkaline metals..... ≤ 0.3 %
Assay (complexometric) ....97.0 ÷ 103.0 %

Code	Size	Packaging	Notes
327507	1 kg	Plastic bottle	
327509	5 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Calcium chloride solution 0.025%

• Calcio cloruro soluzione 0.025% • Calcium chlorure solution 0.025% • Calcio cloruro solución 0.025% • Calciumchlorid 0.025%

Calcium chloride solution 0.025% > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Assay (complexometric) ..... 0.022 ÷ 0.028 %

Code	Size	Packaging	Notes
E433427	1 l	Plastic bottle	

Calcium dihydrogenphosphate monohydrate ▶ Calcium phosphate monobasic monohydrate

Calcium fluoride

• Calcio fluoruro • Calcium fluorure • Calcio fluoruro • Calciumfluorid

Synonym:
Fluorite

CaF_2
Molecular Weight: 78,08
CAS: 7789-75-5
EEC-N: 232-188-7

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Calcium fluoride > RPE - For analysis

RPE

Description	White powder	Chloride	≤ 0.1 %	Cr	≤ 5 ppm	Mn	≤ 1 ppm
Identification	Positive	Sulphate	≤ 50 ppm	Cu	≤ 1 ppm	Ni	≤ 5 ppm
Loss on ignition	≤ 0.5 %	Co	≤ 5 ppm	Fe	≤ 10 ppm	Assay (complexometric)	≥ 98 %

Code	Size	Packaging	Notes
433585	250 g	Plastic bottle	
433587	1 kg	Plastic bottle	

Calcium gluconate

• Calcio gluconato • Calcium gluconate • Calcio gluconato • Calciumgluconat

$C_{12}H_{22}CaO_{14}.H_2O$
Molecular Weight: 448,39
CAS: 18016-24-5

Calcium gluconate > ERBApharm - According to pharmacopoeia: FU-Ph.Eur.

ERBApharm

Description	White crystalline powder	Saccarose and red. sug.	Conform Ph.Eur.	Sulphate	≤ 100 ppm	TYMC	≤ 100 CFU/g
Identification	Positive	Mg and alkaline metals	≤ 0.4 %	Assay (complexometric)	98.5 ÷ 102.0 %		
Organic imp., boric ac.	Conform Ph.Eur.	Heavy metals (Pb)	≤ 10 ppm	TAMC	≤ 1000 CFU/g		

Code	Size	Packaging	Notes
330608	1 kg	Plastic bottle	
330609	5 kg	Plastic tank	
330601	25 kg	Sack	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Calcium hydroxide

• Calcio idrossido • Calcium hydroxyde • Calcio hidróxido • Calciumhydroxid

$Ca(OH)_2$
Molecular Weight: 74,09
CAS: 1305-62-0
EEC-N: 215-137-3

Danger
H315-H318
P264-P280a-P305+P351+P338-P310a-P362+P364-P332+P313

Calcium hydroxide > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description	White powder	Total sulphur	≤ 0.1 %	Assay (alkalimetric)	≥ 95.0 %	Mg	≤ 0.5 %
Identification	Positive	HCl-insoluble matter	≤ 300 ppm	K	≤ 500 ppm		
Carbonate	≤ 3.0 %	Heavy metals (Pb)	≤ 30 ppm	Na	≤ 500 ppm		
Chloride	≤ 300 ppm	Fe	≤ 500 ppm	Sr	≤ 500 ppm		

Code	Size	Packaging	Notes
433875	250 g	Plastic bottle	
433877	1 kg	Plastic bottle	
433874	25 kg	Plastic bucket	
433873	50 kg	Plastic drum	

Product specifications are subject to changes. Please visit our website for updates.

Calcium hydroxide > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP

ERBApharm

Description ..... White powder
Identification ..... Positive
Carbonate..... Conform USP-NF
Organic volatile impurities Conform USP-NF
HCl-insoluble matter..... ≤ 0.5 %
Heavy metals (Pb)..... ≤ 20 ppm
Mg and alkaline metals..... ≤ 4.0 %
Assay (complexometric) ....95.0 ÷ 100.5 %
Chloride..... ≤ 330 ppm
Sulphate ..... ≤ 0.4 %
As ..... ≤ 4 ppm
Assay (alkalimetric).....95.0 ÷ 100.5 %

Code	Size	Packaging	Notes
331007	1 kg	Plastic bottle	
331008	5 kg	Plastic tank	
331003	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Calcium hydroxide > RE - Pure

RE

Description ..... White powder
Identification ..... Positive
As ..... ≤ 3 ppm
Heavy metals (Pb)..... ≤ 20 ppm
Assay (complexometric) ..... ≥ 95 %

Code	Size	Packaging	Notes
326454	1 kg	Plastic bottle	
326458	25 kg	Plastic bucket	

Calcium lactate

• Calcio lattato • Calcium lactate • Calcio lactato • Calciumlactat

Synonym:
L-Lactic acid calcium salt

$(CH_3CHOHCOO)_2Ca.5H_2O$
Molecular Weight: 308,29
CAS: 5743-47-5
EEC-N: 248-953-3

Calcium lactate > ERBApharm - According to pharmacopoeia: BP-FU-Ph.Eur.-Ph.Franc.

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Acidity or alkalinity.....Conform Ph.Eur.
Ba .....Conform Ph.Eur.
Loss on drying ..... 22.0 ÷ 27.0 %
Chloride..... ≤ 200 ppm
Heavy metals (Pb)..... ≤ 10 ppm
Mg and alkaline metals..... ≤ 1 %
Sulphate ..... ≤ 400 ppm
Fe ..... ≤ 50 ppm
Assay (complexometric)98.0 ÷ 102.0 % s.s.

Code	Size	Packaging	Notes
331407	1 kg	Plastic bottle	
331408	5 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Calcium nitrate tetrahydrate

• Calcio nitrato tetraidrato • Calcium nitrate tétrahydraté • Calcio nitrato tetrahidrato • Calcium nitrate tetrahydrat

$Ca(NO_3)_2.4H_2O$
Molecular Weight: 236,15
CAS: 13477-34-4
EEC-N: 233-332-1

Classification transport
ONU: 1454
Transport Hazard class: 5.1
Packing group III



Danger
H272-H315-H319
P210-P280-P305+P351+P338-P332+P313-P362+P364-P337+P313

Calcium nitrate tetrahydrate > RPE - For analysis - ACS

RPE

Description .....White crystals
Identification ..... Positive
pH sol. 5% at 25° C .....5.0 ÷ 7.0
Water-insoluble matter .....≤ 50 ppm
Chloride..... ≤ 50 ppm
Nitrite ..... ≤ 10 ppm
Sulphate..... ≤ 20 ppm
Heavy metals (Pb)..... ≤ 5 ppm
Ba ..... ≤ 50 ppm
Fe ..... ≤ 5 ppm
Assay (complexometric) ....99.0 ÷ 103.0 %
Mg ..... ≤ 0.05 %
K..... ≤ 50 ppm
Na..... ≤ 100 ppm
Sr..... ≤ 0.05 %

Code	Size	Packaging	Notes
433955	100 g	Plastic bottle	
433956	500 g	Plastic bottle	
433957	1 kg	Plastic bottle	
433951	5 kg	Plastic tank	

Calcium nitrate tetrahydrate > RE - Pure

RE

Description ..... White crystals	Ammonium ..... ≤0.5 %	Heavy metals (Pb) ..... ≤50 ppm	Assay (complexometric) ..... 97 ÷ 100 %
Identification ..... Positive	Chloride ..... ≤500 ppm	Sulphate ..... ≤0.1 %	
pH sol. 5% at 25° C ..... 4.0 ÷ 7.0	Water-insoluble matter ..... ≤100 ppm	Fe ..... ≤100 ppm	

Code	Size	Packaging	Notes
331509	5 kg	Plastic tank	
331501	25 kg	Drum	

Calcium oxide, lumps

• Calcio ossido, pezzi • Calcium oxyde, morceaux • Calcio óxido, trozos • Calciumoxid, Klumpen

CaO
Molecular Weight: 56,08
CAS: 1305-78-8
EEC-N: 215-138-9



Danger
H315-H318-H335
P261-P304+P340-P310a-P305+P351+P338-P362+P364-P403+P233

Calcium oxide, lumps > RE - Pure

RE

Description ..... Whitish lumps	Loss on ignition ..... ≤5 %	HCl-insoluble matter ..... ≤0.5 %	Fe ..... ≤0.2 %
Identification ..... Positive	Carbonate ..... ≤5.0 %	As ..... ≤10 ppm	Assay (alkalimetric) ..... ≥95 %

Code	Size	Packaging	Notes
331557	1 kg	Plastic bottle	
331555	25 kg	Plastic bucket	

Calcium oxide, powder

• Calcio ossido, polvere • Calcium oxyde, poudre • Calcio óxido, polvo • Calciumoxid, Pulver

CaO
Molecular Weight: 56,08
CAS: 1305-78-8
EEC-N: 215-138-9



Danger
H315-H318-H335
P261-P304+P340-P310a-P305+P351+P338-P362+P364-P403+P233

Calcium oxide, powder > RE - Pure

RE

Description ..... White powder	Identification ..... Positive	Calcium carbonate ..... ≤10 %	Loss on ignition ..... ≤5 %

Code	Size	Packaging	Notes
331567	1 kg	Plastic bottle	
331564	25 kg	Plastic bucket	

Calcium pantothenate

• Calcio pantotenato • Calcium pantothénate • Calcio pantotenato • Calciumpantothenat

$C_{18}H_{32}O_{10}N_2Ca$
Molecular Weight: 476,6
CAS: 137-08-6
EEC-N: 205-278-9

Calcium pantothenate > ERBAPharm - According to pharmacopoeia: Ph.Eur.-FU-Ph.Franc.-DAB-USP

ERBApharm

Description ..... White powder	pH sol. 5% at 25° C ..... 6.8 ÷ 8.0	3-Aminopropionic acid ..... ≤ 0.5 %	Ca ..... 8.2 ÷ 8.6 % s.s.
Identification ..... Positive	Specific optical rotation on dry ..... +25.5 ÷ +27.5 °	Chloride ..... ≤ 200 ppm	Assay (protonometric) 98.0 ÷ 101.0 % s.s.
Appearance of solution ..... Conform Ph.Eur.	Loss on drying ..... ≤ 3.0 %	Heavy metals (Pb) ..... ≤ 20 ppm	
Organic volatile impurities Conform USP-NF		N ..... 5.7 ÷ 6.0 % s.s.	

Code	Size	Packaging	Notes
331602	100 g	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Product specifications are subject to changes. Please visit our website for updates.

Calcium phosphate dibasic dihydrate

• Calcio fosfato bibasico diidrato • Calcium phosphate dibasique dihydraté • Calcio fosfato dibásico dihidrato • Calciumphosphat zweiwertiges Dihydrat

$CaHPO_4.2H_2O$
Molecular Weight: 172,09
CAS: 7789-77-7
EEC-N: 231-826-1

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Calcium phosphate dibasic dihydrate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph. Franc.

ERBApharm

Description	White crystalline powder	Loss on ignition	24.5 ÷ 26.5 %	Heavy metals (Pb)	≤30 ppm	Assay (complexometric)	98.0 ÷ 105.0 %
Identification	Positive	Chloride	≤0.25 %	Sulphate	≤0.5 %		
Carbonate	Conform Ph.Eur.	Fluoride	≤50 ppm	As	≤3 ppm		
Ba	Conform Ph.Eur.	HCl-insoluble matter	≤0.2 %	Fe	≤400 ppm		

Code	Size	Packaging	Notes
330307	1 kg	Plastic bottle	
330303	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Calcium phosphate monobasic monohydrate

• Calcio fosfato monobasico monoidrato • Calcium phosphate monobasique monohydraté
• Calcio fosfato monobásico monohidrato • Einbasiges Calciumphosphat-Monohydrat

Synonym:
Calcium dihydrogenphosphate monohydrate

$Ca(H_2PO_4)_2.H_2O$
Molecular Weight: 252,07
CAS: 7758-23-8
EEC-N: 231-837-1

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Calcium phosphate monobasic monohydrate > RPE - For analysis

RPE

Description	White crystalline powder	As	≤ 1 ppm	Pb	≤ 1 ppm
Identification (I.R.)	Positive	Cd	≤ 1 ppm	Fluoride	≤ 30 ppm
Al	≤ 200 ppm	Hg	≤ 1 ppm	Assay (on dry)	≥ 95.0 %

Code	Size	Packaging	Notes
433685	250 g	Plastic bottle	

Calcium phosphate tribasic

• Calcio fosfato tribasico • Calcium phosphate tribasique • Calcio fosfato tribásico
• Calciumphosphat tribasisch

Synonym:
tri-Calcium (ortho)phosphate

$Ca_3(PO_4)_2$
Molecular Weight: 310,18
CAS: 7758-87-4
EEC-N: 231-840-8

Calcium phosphate tribasic > RPE - For analysis

RPE

Description	White powder	Loss on ignition	≤ 4 %	Pb	≤ 1 ppm	Zn	≤ 25 ppm
Identification	Positive	Heavy metals (Pb)	≤ 10 ppm	Assay (acidimetric)	≥ 90 %	Assay (P_2O_5)	40.5 ÷ 42 %
pH sol 10%	6.5 ÷ 7.5	As	≤ 1 ppm	Cu + Zn	≤ 50 ppm		

Code	Size	Packaging	Notes
433774	100 g	Plastic bottle	
433776	500 g	Plastic bottle	

Calcium phosphate tribasic > ERBApharm - According to pharmacopoeia: Ph.Eur.

ERBApharm

Description ..... White powder	Loss on ignition ..... ≤ 8.0 %	Heavy metals (Pb) ..... ≤ 30 ppm	Fe ..... ≤ 400 ppm
Identification ..... Positive	Chloride ..... ≤ 0.15 %	Sulphate ..... ≤ 0.5 %	Assay (complexometric) . 35.0 ÷ 40.0 % Ca
HCl-insoluble matter ..... ≤ 0.2 %	Fluoride ..... ≤ 75 ppm	As ..... ≤ 4 ppm	

Code	Size	Packaging	Notes
330407	1 kg	Plastic bottle	
330409	5 kg	Plastic bucket	
330403	25 kg	Fibre drum	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Calcium stearate

• Calcio stearato • Calcium stéarate • Calcio estearato • Calciumstearat

Synonym:
- *Octadecanoic acid calcium salt*
- *Stearic acid calcium salt*

$C_{36}H_{70}CaO_4$
Molecular Weight: 607,04
CAS: 1592-23-0
EEC-N: 216-472-8

Calcium stearate > ERBApharm - Vegetal origin - According to pharmacopoeia: USP-NF

ERBApharm

Description ..... White powder	Heavy metals (Pb) ..... ≤10 ppm	Stearic acid ..... ≥ 40.0 %
Identification ..... Positive	Origin (BSE/TSE) ..... Vegetable	Stearic + palmitic acid ..... ≥ 90.0 %
Loss on drying ..... ≤4.0 %	Assay (complexometric) (d.s) ..... 6.4 ÷ 7.4 % (Ca)	Residual solvents (Current ICH) ..... Conform

Code	Size	Packaging	Notes
332262	2.5 kg	Plastic bottle	
332261	10 kg	Plastic bucket	
332265	25 kg	Fibre drum	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Calcium sulfate dihydrate

• Calcio solfato biidrato • Calcium sulfate dihydraté • Calcio sulfato dihidrato • Calciumsulfat-Dihydrat

$CaSO_4.2H_2O$
Molecular Weight: 172,17
CAS: 10101-41-4
EEC-N: 231-900-3

Calcium sulfate dihydrate > RPE - For analysis - ACS

RPE

Description ..... White powder	HCl-insoluble matter ..... ≤ 0.02 %	K ..... ≤ 50 ppm	Assay (complexometric) ..... 98 ÷ 102 %
Identification ..... Positive	Heavy metals (Pb) ..... ≤ 20 ppm	Mg ..... ≤ 0.02 %	
Carbonate ..... Conform ACS	Nitrate ..... Conform ACS	Na ..... ≤ 0.02 %	
Chloride ..... ≤ 50 ppm	Fe ..... ≤ 10 ppm	Sr ..... ≤ 0.05 %	

Code	Size	Packaging	Notes
434155	100 g	Plastic bottle	
434156	500 g	Plastic bottle	
434151	25 kg	Plastic bucket	

Calcium sulfate dihydrate > ERBApharm - According to pharmacopoeia: NF

ERBApharm

Description ..... White crystalline powder	Fe ..... ≤100 ppm	Loss on drying ..... 19.0 ÷ 23.0 %
Identification ..... Positive	Heavy metals (Pb) ..... ≤10 ppm	Assay ..... 98.0 ÷ 101.0 %s.s.

Code	Size	Packaging	Notes
331752	5 kg	Plastic tank	
331751	25 kg	Sack	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Product specifications are subject to changes. Please visit our website for updates.

Calcium sulfate hemihydrate

• Calcio solfato emidrato • Calcium sulfate hemihydraté • Calcio sulfato hemihidrato • Calciumsulfat-Halbhydrat

$CaSO_4.1/2H_2O$
Molecular Weight: 145,15
CAS: 10034-76-1

Calcium sulfate hemihydrate > RE - Pure

RE

Description ..... White powder Identification ..... Positive Assay ..... ≥ 97 %

Code	Size	Packaging	Notes
331761	1 kg	Plastic bottle	
331762	5 kg	Plastic tank	
331763	25 kg	Plastic bucket	

Calcium sulfate hemihydrate solution

• Calcio solfato emidrato soluzione • Calcium sulfate hemihydraté solution • Calcio sulfato hemihidrato solución • Calciumsulfat-Halbhydratlösung

Calcium sulfate hemihydrate solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611015202	100 ml	Plastic bottle	Ref Ph.Eur 1015201
611015201	1 l	Plastic bottle	Ref Ph.Eur 1015201

Calcon

• Calcone • Calcon • Calcón • Calcon

Synonym:
1-(2-Hydroxy-1-naphthylazo)-2-naphthol-4-sulfonic acid sodium salt

$C_{20}H_{13}N_2NaO_5S$
Molecular Weight: 416,39
CAS: 2538-85-4
EEC-N: 219-810-2

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Calcon > RPE - For analysis - C.I. 15705

RPE

Description ..... Deep purple powder Identification ..... Positive

Code	Size	Packaging	Notes
434171	25 g	Glass bottle	

Complexometric indicator for Al, Fe, Zr

Calconcarbonic acid

• Acido calconcarbonico • Acide calconecarbonique • Acido calconcarbónico • Calconcarbonsäure

Synonym:
• *Calconcarboxylic acid*
• *3-Hydroxy-4-(2-hydroxy-4-sulfo-1-naphthylazo) naphthalene-2-carboxylic acid*

$C_{21}H_{14}O_7N_2S$
Molecular Weight: 438,42
CAS: 3737-95-9
EEC-N: 223-117-0

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Calconcarbonic acid > RPE - For analysis - Reag. Ph. Eur.

RPE

Description ..... Polvere viola scura Identification ..... Positive Water ..... ≤ 10 %

Code	Size	Packaging	Notes
403308	5 g	Glass bottle	

Calmagite

• Calmagite • Calmagite • Calmagita • Calmagite

Synonym:
3-Hydroxy-4-(2-hydroxy-5methylphenylazo)-naphtalene-1-sulfonic aicd

$C_{17}H_{14}N_2O_5S$
Molecular Weight: 358,37
CAS: 3147-14-6
EEC-N: 221-563-0

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Calmagite > RPE - For analysis

RPE

Description .....Black powder Identification ..... Positive

Code	Size	Packaging	Notes
434181	5 g	Glass bottle	

Suitable for the spectrophotometric determination of lanthanide

Camphor natural

• Canfora naturale • Camphre naturel • Alcanfor natural • Kampfer, Natürlicher

Synonym:
• 2-Bornanone
• 2-Camphanone

$C_{10}H_{16}O$
Molecular Weight: 152,23
CAS: 464-49-3
EEC-N: 207-355-2

Classification transport
ONU: 2717
Transport Hazard class: 4.1
Packing group III

Warning
H228-H315-H319-H335
P210-P241-P280-P304+P340-P305+P351+P338-P403+P233

Camphor natural > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-BP

ERBApharm

Description .....White crystals
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Water .....Conform Ph.Eur.
Acidity or alkalinity.....Conform Ph.Eur.
Sostanze analoghe (GLC) ...Conform Ph.Eur.
Melting point..... 175 ÷ 179 ° C
Specific optical rotation... +41.0 ÷ +43.0 °
Halogenated compounds .....≤100 ppm
Non volat.substances .....≤500 ppm

Code	Size	Packaging	Notes
332356	500 g	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Camphor synthetic

• Canfora sintetica • Camphre synthétique • Alcanfor sintetico • Kampfer, Synthetischer

$C_{10}H_{16}O$
Molecular Weight: 152,24
CAS: 76-22-2
EEC-N: 200-945-0

Classification transport
ONU: 2717
Transport Hazard class: 4.1
Packing group III

Warning
H228-H302-H371
P210-P241-P260-P264-P280-P301+P312a

Camphor synthetic > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP-DAB

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Water .....Conform Ph.Eur.
Acidity or alkalinity.....Conform Ph.Eur.
Related substances .....Conform Ph.Eur.
Melting point..... 172 ÷ 180 °C
Specific optical rotation.... -0.15 ÷ +0.15 °
Halogenated compounds .....≤100 ppm
Non volat.substances .....≤500 ppm
Origin (BSE/TSE)..... Synthesis
Residual solvents (Current ICH).....Conform

Code	Size	Packaging	Notes
332406	500 g	Plastic bottle	
332401	5 kg	Carton box	
332402	25 kg	Fibre drum	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Product specifications are subject to changes. Please visit our website for updates.

Canada balsam

• Balsamo del Canada • Baume du Canada • Balsamo del Canada • Kanadabalsam

CAS: 8007-47-4
EEC-N: 232-362-2

Warning
H226
P210-P241-P280-P303+P361+P353-P403+P235-P501a

Canada balsam > RS - Mounting medium for microscopy

RS

Description .....Pale yellow dense liquid Identification ..... Positive

Code	Size	Packaging	Notes
321553	100 g	Glass bottle	
321554	250 g	Glass bottle	

n-Caproic acid

• Acido n-caproico • Acide n-caproïque • Acido n-caproico • n-Capronsäure

Synonym:
Hexanoic acid

$CH_3(CH_2)_4COOH$
Molecular Weight: 116,16
CAS: 142-62-1
EEC-N: 205-550-7

Classification transport
ONU: 2829
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

n-Caproic acid > RPE - For analysis

RPE

Description .....Yellow colourless liquid Identification ..... Positive Refractive index at 20°C.1.4150 ÷ 1.4180 Assay (GLC) ..... ≥99 %

Code	Size	Packaging	Notes
403473	100 ml	Glass bottle	

Capryl alcohol ▶ Octanol-1

n-Caprylic acid

• Acido n-caprilico • Acide n-caprylique • Acido n-caprílico • n-Caprylsäure

Synonym:
Octanoic acid

$CH_3(CH_2)_6COOH$
Molecular Weight: 144,21
CAS: 124-07-2
EEC-N: 204-677-5

Classification transport
ONU: 3265
Transport Hazard class: 8
Packing group III

Danger
H314-H412
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

n-Caprylic acid > RE - Pure

RE

Description ..... Yellow clear liquid Identification ..... Positive Refractive index at 20°C.1.4268 ÷ 1.4288 Assay (GLC) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
403421	250 ml	Glass bottle	

Carbamic acid ammonium salt ▶ Ammonium carbamate

Carbolated Methylene Blue hydroalcoholic solution

• Blu metilene fenicato soluzione idroalcolica • Bleu de méthylène phéniqué solution hydroalcoolique • Azul de metileno fenicado solución hidroalcohólica • Methylenblau-phenolische hydroalkoholische Lösung

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group III



Warning
H226-H315-H319-H341
P210-P241-P280-P303+P361+P353-P305+P351+P338-P308+P313

Carbolated Methylene Blue hydroalcoholic solution > RS - For microscopy

RS

Description ..... Blue clear liquid Identification ..... Positive

Code	Size	Packaging	Notes
428991	100 ml	Bottle	

Dye for bacteriology. Water / ethanol mixture (70:30). Contains phenol

Carbolated Toluidine Blue hydroalcoholic solution

• Blu toluidina fenicato soluzione idroalcolica • Bleu de toluidine phéniqué solution hydroalcoolique • Azul de metileno fenicado solución hidroalcohólica • Phenolische Toluidinblau-Hydroalkohol-Lösung



Warning
H315-H319-H341
P280-P305+P351+P338-P308+P313-P362+P364-P332+P313-P337+P313

Carbolated Toluidine Blue hydroalcoholic solution > RS - For microscopy

RS

Description ..... Blue clear liquid Identification ..... Positive

Code	Size	Packaging	Notes
429291	100 ml	Bottle	

Dye for histology. Ethanol - water (10:90). Contains phenol

Carbon ▶ Charcoal activated

Carborundum, granules

• Carborundo, granulare • Carborundum, granulés • Carborundo, gránulos • Carborundum-Granulat

Synonym:
Silicon carbide

CSi
Molecular Weight: 40,1
CAS: 409-21-2
EEC-N: 206-991-8



Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Carborundum, granules > RPE - For analysis

RPE

Description ..... Black granules
Identification ..... Positive
Loss on ignition ..... ≤0.5 %
Chloride ..... ≤100 ppm
Heavy metals (Pb) ..... ≤50 ppm
Sulphate ..... ≤200 ppm
HCl solubility ..... ≤0.5 %
Fe ..... ≤200 ppm

Code	Size	Packaging	Notes
434766	500 g	Plastic bottle	

Carborundum, powder

• Carborundo, polvere • Carborundum, poudre • Carborundo, polvo • Carborundum-Pulver

Synonym:
Silicon carbide

CSi
Molecular Weight: 40,1
CAS: 409-21-2
EEC-N: 206-991-8



Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Carborundum, powder > RPE - For analysis

RPE

Description ..... Black greyish powder
Identification ..... Positive
Loss on ignition ..... ≤0.5 %
HCl solubility ..... ≤1 %
Fe ..... ≤0.5 %

Code	Size	Packaging	Notes
434786	500 g	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Carrez reagent potassium salt

• Carrez reattivo sale di potassio • Réactif de Carrez sel de potassium • Carrez reactivo sal de potasio • Carrez Reagenz Kaliumsalz

Synonym: *Carrez II*

HEU032

Carrez reagent potassium salt > RS - For agroalimentary analysis

RS

Density at 20°C ..... 1.056 ÷ 1.062

Code	Size	Packaging	Notes
502711	1 l	Plastic bottle	

Composition: according to NF V04-233: K4Fe(CN)6 3 H2O 106g water QSP 1 L

Carrez reagent zinc salt

• Carrez reattivo sale di zinco • Réactif de Carrez sel de zinc • Carrez reactivo sal de zinc • Carrez Reagenz Zinksalz

Synonym: *Carrez I*

Warning
H319-H412
P264-P273-P280i-P305+P351+P338-P337+P313-P501a

Carrez reagent zinc salt > RS - For agroalimentary analysis

RS

Density at 20°C ..... 1.108 ÷ 1.114 pH at 20°C ..... 4.40 ÷ 4.60

Code	Size	Packaging	Notes
502701	1 l	Glass bottle	

Composition: Zinc acetate dihydrate 219 g, acetic acid 30ml, water qsp 1 L

Casein

• Caseina • Caséine • Caseina • Kasein

CAS: 9000-71-9
EEC-N: 232-555-1

Casein > RS - For microscopy

RS

Description ..... Yellowish powder
Identification ..... Positive
Water ..... ≤ 13 %
Acidi liberi (ac. lattico) ..... ≤ 1 %
Ash ..... ≤ 3 % (s.s.)
Assay (ex nitrogen) ..... ≥ 92 % s.s.

Code	Size	Packaging	Notes
435963	50 g	Plastic bottle	

Castor oil

• Olio di ricino • Huile de ricin • Aceite de ricino • Rizinusöl

Molecular Weight: 932
CAS: 8001-79-4
EEC-N: 232-293-8

Castor oil > ERBApharm - According to pharmacopoeia: BP-FU-Ph.Eur.-Ph.Franc.

ERBApharm

Description ..... Slightly yellow, viscous liq.
Appearance ..... Clear at 40°C
Appearance of solution ..... Conform
Identification ..... Positive
Relative density at 20°C ..... about 0.958
Refractive index at 20° C ..... about 1.479
Refractive index at 20°C ..... about 1.479
Optical rotation ..... +3.5 ÷ +6.0 °
Specific absorbance at 270nm ..... ≤ 0.7 AU
Acid value ..... ≤ 1.5
Hydroxyl value ..... ≥ 160
Peroxide value ..... ≤ 10.0
Unsaponifiable matter ..... ≤ 0.8 %
Water (K.F.) ..... ≤ 0.3 %
Composition of fatty acids (GC) ..... Conform Ph.Eur.
Palmitic acid ..... ≤ 2.0 %
Stearic acid ..... ≤ 2.5 %
Ricinoleic acid ..... 85.0 - 92.0 %
Oleic acid and isomers
Linoleic acid ..... 2.5 - 7.0 %
Linolenic acid ..... ≤ 1.0 %
Eicosenoic acid ..... ≤ 1.0 %
Any other fatty acid ..... ≤ 1.0 %
Origin (BSE/TSE) ..... Vegetable
Residual solvents (Current ICH) ..... Conform
Residual solvents (CPMP/ICH/283/95) ..... Conform

Code	Size	Packaging	Notes
356351	1 l	Glass bottle	
356352	5 l	Aluminium can	
356353	28 kg	Metal drum	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Cedarwood oil ▶ Oil of cedar wood

Cellulose, powder

• Cellulosa, polvere • Cellulose, poudre • Celulosa, polvo • Zellulose, Pulver

Synonym:
• *α-Cellulose*
• *Cotton linters*

CAS: 9004-34-6
EEC-N: 232-674-9

Cellulose, powder > RS - For chromatography

RS

Description ..... Yellowish powder Identification ..... Positive

Code	Size	Packaging	Notes
436061	250 g	Plastic bottle	

Ceric ammonium nitrate ▶ Cerium (IV) ammonium nitrate

Ceric ammonium sulfate dihydrate ▶ Cerium (IV) ammonium sulfate dihydrate

Cerium standard solution

• Cerio standard soluzione • Cérium standard solution • Cerio, solución patrón • Cer-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Cerium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505552	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505555	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Cerium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503501	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503503	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503505	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503507	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Cerium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507531	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507498	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Product specifications are subject to changes. Please visit our website for updates.

Cerium (III) nitrate hexahydrate

• Cerio nitrato oso esaidrato • Cérium (III) nitrate hexahydraté • Cerio (III) nitrato hexahidratado
• Cer(III)-nitrat-Hexahydrat

Synonym:
• *Cerium trinitrate*
• *Nitric acid cerium salt*

$Ce(NO_3)_3.6H_2O$
Molecular Weight: 434,25
CAS: 10294-41-4
EEC-N: 233-297-2

Classification transport
ONU: 1477
Transport Hazard class: 5.1
Packing group III

Danger
H272
P210-P220-P280-P370+P378a-P501a

Cerium (III) nitrate hexahydrate > RE - Pure

RE

Description ..... Colourless crystals
Identification ..... Positive
Sulphate ..... ≤ 200 ppm
Fe2O3 ..... ≤ 35 ppm
CaO ..... ≤ 50 ppm
Assay ..... ≥ 99.4 %

Code	Size	Packaging	Notes
436203	50 g	Glass bottle	

Cerium (IV) ammonium nitrate

• Cerio ammonio nitrato ico • Cérium (IV) ammonium nitrate • Cerio (IV) amonio nitrato
• Cer (IV) ammoniumnitrat

Synonym:
Ceric ammonium nitrate

$(NH_4)_2Ce(NO_3)_6$
Molecular Weight: 548,23
CAS: 16774-21-3
EEC-N: 240-827-6

Classification transport
ONU: 1477
Transport Hazard class: 5.1
Packing group II

Danger
H272-H318
P210-P220-P280-P305+P351+P338-P310a-P501a

Cerium (IV) ammonium nitrate > RPE - For analysis - ACS - Reag. Ph.Eur.

RPE

Description ..... Orange crystalline powder
Identification ..... Positive
H2SO4-insoluble matter ..... ≤ 500 ppm
Chloride ..... ≤ 100 ppm
Phosphate ..... ≤ 200 ppm
Assay (oxidimetric) ..... ≥ 98.5 %
Fe ..... ≤ 50 ppm

Code	Size	Packaging	Notes
436081	50 g	Glass bottle	
436082	500 g	Glass bottle	

Cerium (IV) ammonium nitrate 0.1 mol/l

• Cerio ammonio nitrato ico 0.1 mol/l • Cérium (IV) ammonium nitrate 0.1 mol/l • Cerio (IV) amonio nitrato 0.1 mol/l • Cer (IV) ammoniumnitrat 0.1 mol/l

$Ce(NH_4)_2(NO_3)_6$
Molecular Weight: 548,22
CAS: 16774-21-3

Classification transport
ONU: 3093
Transport Hazard class: 8
Packing group II

Danger
H272-H290-H315-H318
P210-P280-P305+P351+P338-P310a-P362+P364-P332+P313

Cerium (IV) ammonium nitrate 0.1 mol/l > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613000100	1 l	Glass bottle	Ref Ph.Eur 3000100

Storage: protected from light

Cerium (IV) ammonium nitrate 0.01 mol/l

• Cerio ammonio nitrato ico 0.01 mol/l • Cérium (IV) ammonium nitrate 0.01 mol/l • Cerio (IV) amonio nitrato 0.01 mol/l • Cer (IV) ammoniumnitrat 0.01 mol/l

$Ce(NH_4)_2(NO_3)_6$
Molecular Weight: 548,22
CAS: 16774-21-3

Cerium (IV) ammonium nitrate 0.01 mol/l > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613000200	1 l	Plastic bottle	Ref Ph.Eur 3000200

Storage: protected from light

Cerium (IV) ammonium sulfate dihydrate

• Cerio ammonio solfato ico diidrato • Cérium (IV) ammonium sulfate dihydraté
• Cerio (IV) amonio sulfato dihidrato • Cer (IV) ammoniumsulfatdihydrat

Synonym:
Ceric ammonium sulfate

$Ce(NH_4)_4 (SO_4)_4.2H_2O$
Molecular Weight: 632,54
CAS: 10378-47-9

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Cerium (IV) ammonium sulfate dihydrate > RPE - For analysis - ACS

RPE

Description ..... Yellow-orange powder
Identification ..... Positive
Insoluble in diluted sulphuric acid ..... ≤ 0.05 %
Phosphate ..... ≤ 0.03 %
Fe ..... ≤ 100 ppm
Assay (oxidimetric) ..... ≥ 94 %

Code	Size	Packaging	Notes
436091	100 g	Glass bottle	
436092	500 g	Glass bottle	

Cerium (IV) ammonium sulfate 0.1 mol/l

• Cerio ammonio solfato ico 0.1 mol/l • Cérium (IV) ammonium sulfate 0.1 mol/l • Cerio (IV) amonio sulfato 0.1 mol/l • Cer (IV) ammoniumsulfat 0.1 mol/l

$Ce(NH_4)_4(SO_4)_4 \cdot 2H_2O$
Molecular Weight: 632,54
CAS: 10378-47-9

Classification transport
ONU: 2796
Transport Hazard class: 8
Packing group II

Warning
H290-H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Cerium (IV) ammonium sulfate 0.1 mol/l > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613000301	250 ml	Plastic bottle	Ref Ph.Eur 3000300
613000300	1 l	Plastic bottle	Ref Ph.Eur 3000300

Cerium (IV) ammonium sulfate 0.01 mol/l

• Cerio ammonio solfato ico 0.01 mol/l • Cérium (IV) ammonium sulfate 0.01 mol/l • Cerio (IV) amonio sulfato 0.01 mol/l • Cer (IV) ammoniumsulfat 0.01 mol/l

$Ce(NH_4)_4(SO_4)_4 \cdot 2H_2O$
Molecular Weight: 632,54
CAS: 10378-47-9

Cerium (IV) ammonium sulfate 0.01 mol/l > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613000400	1 l	Plastic bottle	Ref Ph.Eur 3000400

Cerium (IV) sulfate tetrahydrate

• Cerio solfato ico tetraidrato • Cérium (IV) sulfate tétrahydraté • Cerio (IV) sulfato tetrahidrato
• Cer (IV) sulfattetrahydrat

Synonym:
Ceric sulfate tetrahydrate

$Ce(SO_4)_2 \cdot 4H_2O$
Molecular Weight: 404,3
CAS: 10294-42-5
EEC-N: 237-029-5

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Cerium (IV) sulfate tetrahydrate > RE - Pure

RE

Description ..... Yellow powder
Identification ..... Positive
Chloride ..... ≤200 ppm
Heavy metals (Pb) ..... ≤50 ppm
Fe ..... ≤500 ppm
Assay (oxidimetric) ..... ≥98 % s.s.

Code	Size	Packaging	Notes
436402	25 g	Glass bottle	
436404	100 g	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Cerium (IV) sulfate 0.1 mol/l

• Cerio solfato ico 0.1 mol/l • Cérium (IV) sulfate 0.1 mol/l • Cerio (IV) sulfato 0.1 mol/l • Cer (IV) sulfat 0.1 mol/l

$Ce(SO_4)_2 \cdot 4H_2O$
Molecular Weight: 404,3
CAS: 10294-42-5

Classification transport
ONU: 2796
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Cerium (IV) sulfate 0.1 mol/l > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613001101	500 ml	Plastic bottle	Ref Ph.Eur 3001100
613001100	1 l	Plastic bottle	Ref Ph.Eur 3001100

Cerium (IV) sulfate 0.1 mol/l > RPE - For analysis

RPE

Code	Size	Packaging	Notes
436426	500 ml	Plastic bottle	

For specifications, contact our customer service for a certificate of analysis

Cesium standard solution

• Cesio standard soluzione • Césium standard solution • Cesio, solución patrón • Caesium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Cesium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505572	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505575	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Cesium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503531	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503533	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503535	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503537	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Cesium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507532	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507499	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Cesium chloride

• Cesio cloruro • Césium chlorure • Cesio cloruro • Cäsiumchlorid

CsCl
Molecular Weight: 168,36
CAS: 7647-17-8
EEC-N: 231-600-2

Warning
H312
P280h-P302+P352a-P312a-P362+P364-P501a

Cesium chloride > RPE - For analysis

RPE

Description	White crystalline powder	Ba	≤ 10 ppm	K	≤ 30 ppm	Rb	≤ 50 ppm
Identification	Positive	Ca	≤ 10 ppm	Mg	≤ 5 ppm	Sr	≤ 10 ppm
Sulphate	≤ 25 ppm	Cr	≤ 2 ppm	Na	≤ 50 ppm	Assay (argentimetric)	≥ 99 %
SiO2	≤ 2 ppm	Fe	≤ 5 ppm	Mn	≤ 5 ppm		
Al	≤ 1 ppm	Li	≤ 1 ppm	Pb	≤ 5 ppm		

Code	Size	Packaging	Notes
436502	25 g	Glass bottle	
436501	50 g	Glass bottle	

Cesium chloride 25 g/l solution

• Cesio cloruro 25 g/L soluzione • Césium chlorure 25 g/l • Cesio cloruro solución 25 g/l • Caesiumchlorid 25 g/l Lösung

CsCl
Molecular Weight: 168,36
CAS: 7647-17-8

HEU210

Cesium chloride 25 g/l solution > RS - Ionisation standard solution for AAS

RS

Code	Size	Packaging	Notes
504536	500 ml	Plastic bottle	Matrix: Water

Cesium sulfate

• Cesio solfato • Césium sulfate • Cesio sulfato • Cäsiumsulfat

Cs_2SO_4
Molecular Weight: 361,87
CAS: 10294-54-9
EEC-N: 233-662-6

Cesium sulfate > RPE - For analysis

RPE

Description	White crystalline powder	Heavy metals (Pb)	≤20 ppm	K	≤500 ppm	Zn	≤50 ppm
Identification	Positive	Al	≤5 ppm	Mg	≤5 ppm	Assay (acidimetric)	≥99 %
Total nitrogen	≤30 ppm	Ca	≤100 ppm	Na	≤200 ppm		
Chloride	≤40 ppm	Cu	≤5 ppm	Ni	≤5 ppm		
Water-insoluble matter	≤50 ppm	Fe	≤3 ppm	Pb	≤5 ppm		

Code	Size	Packaging	Notes
436534	25 g	Glass bottle	

Cetyl alcohol

• Alcole cetilico • Alcool cétylique • Alcohol cetílico • Cetylalkohol

Synonym:
1-Hexadecanol

$CH_3(CH_2)_{14}CH_2OH$
Molecular Weight: 242,44
CAS: 36653-82-4
EEC-N: 253-149-0

Warning
H315
P264-P280g-P302+P352a-P332+P313-P362+P364

Cetyl alcohol > RPE - For analysis

RPE

Description	Polvere cristallina bianca o incolore	Residue on ignition	≤0.1 %	Assay (GLC)	≥95.0% %
Identification	Positive	Melting point	47 ÷ 50 °C		

Code	Size	Packaging	Notes
414427	1 kg	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Cetyl alcohol > ERBApharm - According to pharmacopoeia: NF-Ph.Eur.-Ph.Franc.

ERBApharm

Description ..... White powder and flakes
Identification ..... Positive
Appearance of solution ..... Conform Ph.Eur.
Melting point ..... 46 ÷ 52 °C
Acid value ..... ≤ 1.0
Iodine value ..... ≤ 2.0
Hydroxyl value ..... 218 ÷ 238
Saponification value ..... ≤ 2.0
Water content ..... ≤ 0.5 %
Assay (GLC) ..... ≥ 95 %
Origin (BSE/TSE) ..... Vegetable

Code	Size	Packaging	Notes
308357	1 kg	Plastic bottle	
308358	5 kg	Plastic tank	
308359	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Charcoal activated

• Carbone attivo • Charbon actif • Carbón activado • Holzkohle aktiviert

Synonym: *Carbon*

C
Molecular Weight: 12,01
CAS: 7440-44-0
EEC-N: 931-328-0

Charcoal activated > RS - For chromatography

RS

Description ..... Black powder
Identification ..... Positive
Loss on drying ..... ≤ 10 %

Code	Size	Packaging	Notes
434455	250 g	Bag	
434454	1 kg	Bag	

Charcoal activated > RS - For microanalysis

RS

Description ..... Black fine powder
Identification ..... Positive

Code	Size	Packaging	Notes
434462	50 g	Glass bottle	

Charcoal activated > RE - Pure - For synthesis

RE

Code	Size	Packaging	Notes
P4610017	1 kg	Plastic bottle	Granular

For specifications, contact our customer service for a certificate of analysis

Charcoal decolorizing

• Carbone decolorante • Charbon décolorant • Carbón decolorante • Entfärbung von Holzkohle

Synonym: *Carbon*

C
Molecular Weight: 12,01
CAS: 7440-44-0
EEC-N: 931-328-0

Charcoal decolorizing > RPE - For analysis

RPE

Description ..... Black fine powder
Identification ..... Positive
pH ..... 10 - 11
Loss on drying ..... ≤ 10 %
Ash ..... ≤ 7 %

Code	Size	Packaging	Notes
434507	1 kg	Bag	
434501	20 kg	Fibre drum	

Charcoal vegetable

• Carbone vegetale • Charbon végétal • Carbón vegetal • Holzkohlegemüse

Synonym:
Carbon

C
Molecular Weight: 12,01
CAS: 7440-44-0
EEC-N: 931-328-0

Charcoal vegetable > RE - Pure

RE

Description ..... Black fine powder Identification ..... Positive

Code	Size	Packaging	Notes
332658	2.5 kg	Bag	
332659	5 kg	Plastic bucket	

Chloral hydrate

• Cloralio idrato • Chloral hydraté • Cloral hidrato • Hydriertes Chloral

$Cl_3CCH(OH)_2$
Molecular Weight: 165,4
CAS: 302-17-0

Classification transport
ONU: 2810
Transport Hazard class: 6.1
Packing group III



Danger
H301-H315-H319
P301+P310a-P330-P305+P351+P338-P362+P364-P332+P313-P337+P313

Chloral hydrate > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611017901	100 ml	Plastic bottle	Ref Ph.Eur 1017901

Chloramine T sodium salt

• Cloramina T sale sodico • Chloramine T sel sodique • Cloramina T sódica • Chloramins-T-Natriumsalz

Synonym:
N-Chloro-p-toluenesulfonamide sodium salt

$CH_3C_6H_4SO_2NClNa.3H_2O$
Molecular Weight: 227,65
CAS: 7080-50-4
EEC-N: 615-172-8

Classification transport
ONU: 3263
Transport Hazard class: 8
Packing group III



Danger
H302-H314-H334-HEU031
P280-P284-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P342+P311a

Chloramine T sodium salt > RPE - For analysis

RPE

Description ..... White crystalline powder pH solution 5% ..... 8.0 ÷ 10.3 Torbidità (0.5% in acqua) ..... ≤ 5 FTU
Identification ..... Positive Colore soluzione 5% ..... ≤ 25 APHA Assay (iodometric) ..... ≥ 97.5 %

Code	Size	Packaging	Notes
437555	25 g	Glass bottle	
437554	100 g	Plastic bottle	
437557	1 kg	Plastic bottle	
437551	25 kg	Fibre drum	

For determination of: Co, Cr, Fe, Hg, Mn, Ni, Sb

Chloranil

• Cloranile • Chloranile • Cloranilo • Chloranil

Synonym:
• *Tetrachloro-1,4-benzoquinone*
• *Tetrachloro-p-benzoquinone*

COCCl:CClCOCCl:CCl
Molecular Weight: 245,89
CAS: 118-75-2
EEC-N: 204-274-4

Classification transport
ONU: 3077
Transport Hazard class: 9
Packing group III



Warning
H315-H319-H410
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Chloranil > RPE - For analysis

RPE

Description Yellow to yellow-green powder or crystals Identification ..... Positive Assay (HPLC) ..... ≥ 98.5 % Loss on drying ..... ≤ 0.5 %

Code	Size	Packaging	Notes
437601	50 g	Glass bottle	

Reagent for the dehydrogenation of hydroaromatic compounds

Product specifications are subject to changes. Please visit our website for updates.

Chloranilic acid

• Acido cloranilico • Acide chloranilique • Acido cloranílico • Chloranilsäure

Synonym:
2,5-Dichloro-3,6-dihydroxy-2,5-cyclohexadiene-1,4-dione

$C_6H_2Cl_2O_4$
Molecular Weight: 208,99
CAS: 87-88-7
EEC-N: 201-780-7

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Chloranilic acid > RPE - For analysis

RPE

Description .....Red brick powder
Identification ..... Positive
Loss on drying ..... ≤1 %
Chloride.....≤50 ppm
Heavy metals (Pb).....≤20 ppm
Residue on ignition..... ≤0.1 %
Assay (argentimetric)..... ≥99 %

Code	Size	Packaging	Notes
403821	10 g	Glass bottle	
403822	100 g	Glass bottle	

Chlorate standard solution

• Clorati standard soluzione • Chlorate standard solution • Clorato, solución patrón • Chlorat-Standardlösung

Chlorate standard solution > RS - Standard solution for ion chromatography

RS

Code	Size	Packaging	Notes
503181	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503183	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Chloride standard solution

• Cloruri standard soluzione • Chlorure standard solution • Cloruro, solución patrón • Chlorid-Standardlösung

Chloride standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615000901	100 ml	Plastic bottle	A 5 ppm solution: to dilute according to Ref Ph.Eur 5000901
615000909	100 ml	Plastic bottle	A 8 ppm solution: to dilute according to Ref Ph.Eur 5000900
615004100	100 ml	Plastic bottle	A 50 ppm solution: to dilute according to Ref Ph.Eur 5004100

Chloride standard solution > RS - Standard solution for ion chromatography

RS

Code	Size	Packaging	Notes
503231	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503233	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Chlorite standard solution

• Cloriti standard soluzione • Chlorite standard solution • Clorita, solución patrón • Chlorit-Standardlösung

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group III

Chlorite standard solution > RS - Standard solution for ion chromatography

RS

Code	Size	Packaging	Notes
503191	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503193	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Chloroacetamide

• 2-Cloroacetammide • Chloroacétamide • Cloroacetamida • Chloracetamid

$CH_2ClCONH_2$
Molecular Weight: 93,51
CAS: 79-07-2
EEC-N: 201-174-2

Classification transport
ONU: 2811
Transport Hazard class: 6.1
Packing group III



Danger
H301-H317-H361f
P261-P280-P301+P310a-P330-P308+P313-P362+P364

Chloroacetamide > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Melting point ..... 116 ÷ 120 ° C
Water ..... ≤ 0.2 %
Assay (GLC) ..... ≥ 97.5 %

Code	Size	Packaging	Notes
437704	100 g	Glass bottle	

Chloroacetic acid

• Acido cloroacetico • Acide chloroacétique • Acido cloroacetico • Chloressigsäure

Synonym:
Monochloroacetic acid

$CH_2ClCOOH$
Molecular Weight: 94,5
CAS: 79-11-8
EEC-N: 201-178-4

Classification transport
ONU: 1751
Transport Hazard class: 6.1
Packing group II



Danger
H301-H311-H331-H314-H400
P280-P301+P330+P331-P303+P361+P353-P304+P340-P311a-P305+P351+P338-P361+P364-P403+P233

Chloroacetic acid > RPE - For analysis - Reag. Ph. Eur.

RPE

Description ..... White flakes
Identification ..... Positive
Melting point ..... 61 ÷ 63 °C
Assay (acidimetric) ..... ≥98.5 %

Code	Size	Packaging	Notes
404308	500 g	Plastic bottle	

p-Chlorobenzaldehyde

• p-Clorobenzaldeide • p-Chlorobenzaldéhyde • p-Clorobenzaldehído • p-Chlorbenzaldehyd

ClC_6H_4CHO
Molecular Weight: 140,57
CAS: 104-88-1
EEC-N: 203-247-4



Warning
H302-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

p-Chlorobenzaldehyde > RPE - For analysis

RPE

Description ..... White crystalline mass
Identification ..... Positive
Melting point ..... 46.0 ÷ 48.0 ° C
Residue on ignition ..... ≤100 ppm
Assay (GLC) ..... 99 ÷ 100 %

Code	Size	Packaging	Notes
438203	50 g	Glass bottle	

Chlorobenzene

• Clorobenzene • Chlorobenzène • Clorobenceno • Chlorbenzol

C_6H_5Cl
Molecular Weight: 113
CAS: 108-90-7
EEC-N: 203-628-5

Classification transport
ONU: 1134
Transport Hazard class: 3
Packing group III



Warning
H226-H332-H315-H411
P210-P241-P261-P280-P303+P361+P353-P304+P340

Chlorobenzene > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Alcohol miscibility ..... Complete
Diethyl ether miscib. ..... Complete
Density at 20° C ..... 1.103 ÷ 1.109
Refractive index at 20°C 1.5198 ÷ 1.5298
Boiling point ..... 131 ÷ 133 °C
Water (K.F.) ..... ≤200 ppm
Residue on evaporation ..... ≤30 ppm
Acidity (HCl) ..... ≤3 ppm
Benzene ..... ≤200 ppm
p- Dichlorobenzene ..... ≤0.02 %
o- Dichlorobenzene ..... ≤0.01 %
Free chlorine ..... ≤0.1 ppm
Assay (GLC) ..... ≥99.9 %

Code	Size	Packaging	Notes
438251	1 l	Glass bottle	
438255	2.5 l	Glass bottle	
438253	25 kg	Drum	

Product specifications are subject to changes. Please visit our website for updates.

Chlorobenzene > RE - Pure

RE

Description ..... Yellow clear liquid	Density at 20° C ..... 1.103 ÷ 1.109	Free acid (HCl) ..... ≤10 ppm	Residue on evaporation ..... ≤50 ppm
Identification ..... Positive	Refractive index at 20°C 1.5198 ÷ 1.5298	Benzene ..... ≤200 ppm	Assay (GLC) ..... ≥99.9 %
Titrable base ..... Conform	Boiling point ..... 131.5 ÷ 132.5 °C	Water (K.F.) ..... ≤500 ppm	

Code	Size	Packaging	Notes
334251	1 l	Glass bottle	
334255	2.5 l	Glass bottle	
334254	30 kg	Metal drum	

1-Chlorobutane ▶ n-Butyl chloride

Chlorobutanol hemihydrate

• Clorobutanolo emiidrato • Chlorobutanol hémihydrate • Clorobutanol hemihidrato
• Chlorobutanol Hemihydrat

Synonym:
β,β,β-*Trichloro-t-butanol*

$(CH_3)_2COHCCl_3$ù.1/2H_2O
Molecular Weight: 186,5
CAS: 6001-64-5
EEC-N: 200-317-6

Warning
H302-H312-H332-H315-H319
P261-P271-P304+P340-P305+P351+P338-P332+P313-P337+P313

Chlorobutanol hemihydrate > ERBApharm - According to pharmacopoeia: BP-FU-NF-Ph.Eur.-Ph.Franc.

ERBApharm

Description ... Colourless crystalline powder	Impurity B (acetone) (Ph.Eur.) ..... ≤ 0.10 %	Chloride ..... ≤100 ppm	Residual solvents (Current ICH) ..... Conform
Identification ..... Positive	Acidity ..... Conform Ph.Eur.	Assay (argentimetric) .. 98.0 ÷ 101.0 % s.s.	Reaction ..... Conform USP-NF
Appearance of solution ..... Conform Ph.Eur.	Water (K.F.) ..... 4.5 ÷ 5.5 %	Assay (USP) ..... 98.0 - 100.5 % (d.s.)	
Impurity A (chloroform) (Ph.Eur.) .. ≤ 60 ppm	Sulphated ash ..... ≤0.1 %	Origin (BSE/TSE) ..... Synthesis	

Code	Size	Packaging	Notes
301357	1 kg	Plastic bottle	
301356	5 kg	Plastic tank	

This product should be used in compliance with the current legislation.

1-Chloro-2,4-dinitrobenzene ▶ 2,4-Dinitrochlorobenzene

Chloroform

• Cloroformio • Chloroforme • Cloroformo • Chloroform

Synonym:
• *Methylidine trichloride*
• *Trichloromethane*

$CHCl_3$
Molecular Weight: 119,38
CAS: 67-66-3
EEC-N: 200-663-8

Classification transport
ONU: 1888
Transport Hazard class: 6.1
Packing group III

Danger
H302-H331-H315-H319-H351-H361d-H372-HEU301
P280-P304+P340-P305+P351+P338-P308+P313-P330-P362+P364-P403+P233

Chloroform > RS - For HPLC - Isocratic grade - Stabilized with amylene

RS

Description ..... Clear colourless liquid	Water (K.F.) ..... ≤ 100 ppm	Transmittance		Methylene chlorure ..... ≤ 50 ppm
Identification (I.R.) ..... Positive	Residue on evaporation ..... ≤ 5 ppm	At 250 nm ..... ≥ 50 %		Stabilized with amylene ..... ≤ 60 ppm
Density at 20°C ..... 1.479 - 1.483	Acidity ..... ≤ 0.0005 meq/g	At 260 nm ..... ≥ 90 %		
Refractive index at 20°C .. 1.4456 - 1.4496	Alkalinity ..... ≤ 0.0002 meq/g	At 275 nm ..... ≥ 98 %		
Boiling point ..... 61.0 - 61.5 °C	Assay (CPG) ..... ≥ 99.9 %	Carbon tetrachloride ..... ≤ 100 ppm		

Code	Size	Packaging	Notes
412571	1 l	Glass bottle	
412572	2.5 l	Glass bottle	

Chloroform > RS - For HPLC - Isocratic grade - Stabilized with ethanol

RS

Description ..... Clear colourless liquid	Boiling point ..... 61.0 ÷ 61.5 ° C	Alcalinity ..... ≤0.0002 meq/g	at 260 nm ..... ≥90 %	
Identification ..... Positive	Water (K.F.) ..... ≤ 100 ppm	Assay (GLC) ..... ≥99.9 %	at 275 nm ..... ≥98 %	
Density at 20° C ..... 1.479 ÷ 1.483	Residue on evaporation ..... ≤5 ppm	U.V. Transmittance	Carbon tetrachloride ..... ≤ 100 ppm	
Refractive index at 20°C ..... 1.444 ÷ 1.448	Acidity ..... ≤0.0005 meq/g	at 250 nm ..... ≥50 %	Methylene chlorure ..... ≤ 50 ppm	

Code	Size	Packaging	Notes
412652	1 l	Glass bottle	
412653	2.5 l	Glass bottle	

Chloroform > RS - For preparative HPLC - Stabilized with ethanol

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 1.479 ÷ 1.483
Refractive index at 20°C. 1.4456 ÷ 1.4496
Boiling point..... 61.0 ÷ 61.5 ° C
Water (K.F.) ..... ≤500 ppm
Residue on evaporation ..... ≤5 ppm
Alcalinity..... ≤0.0002 meq/g
Assay (GLC) ..... ≥99.0 %
U.V. Transmittance
at 250 nm ..... ≥50 %
at 275 nm ..... ≥98 %
Stabilized with ethyl alcohol .... 0.6 ÷ 1.0 %

Code	Size	Packaging	Notes
438641	2.5 l	Glass bottle	

Chloroform > RS - For GC-MS - Stabilized with ethanol

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C..... 1.444 - 1.448
Water (K.F.) ..... ≤ 100 ppm
Colour ..... ≤ 10 APHA
Residue on evaporation ..... ≤ 3 ppm
Acidity (HCl) ..... ≤ 5 ppm
Assay (GC) (without stabilizer) ... ≥ 99.95 %
Stabilizer (Ethanol)..... 0.4 - 1.0 % (w/w)
Carbon tetrachloride..... ≤ 100 ppm
Dichloromethane..... ≤ 100 ppm
GC-MS.Individual peak (n-hexadecane). ≤ 2 µg/L
Ret.range n-undecane to n-tetracontane

Code	Size	Packaging	Notes
438732	1 l	Glass bottle	

Chloroform > RS - ATRASOL - For analysis of volatile traces - Stabilized with ethanol

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C..... 1.444 - 1.448
Water content (K.F.) ..... ≤ 100 mg/Kg
Colour ..... ≤ 10 Hazen
Stabilizer (Ethanol)..... 0.4 - 1 % m/m
Dichloromethane..... ≤ 100 mg/Kg
Non volatile residue..... ≤ 3 mg/Kg
GC (FID) - NC Atrasol..... Conform
Carbon tetrachloride..... ≤ 100 mg/Kg
Assay (GC) (without stabilizer) ... ≥ 99.95 %
Free acid (as HCl)..... ≤ 5 mg/Kg
GC-ECD.Individual peak (Lindane) .≤ 2 ng/L
Ret.range 1,2,4-trichlorobenzene
to decachlorobiphenyle
GC-FID.Individual peak (n-hexadecane). ≤ 2 µg/L
Ret.range n-undecane to n-tetracontane

Code	Size	Packaging	Notes
P02432E16	1 l	Glass bottle	
P02432E21	2.5 l	Glass bottle	

Chloroform > RS - PESTIPUR - For pesticide analysis - Stabilized with amylene

RS

Refractive index at 20°C..... 1.444 - 1.448
Water content (K.F.) ..... ≤ 100 mg/Kg
Colour ..... ≤ 10 Hazen
Stabilizer (Amylene)..... 5 - 50 mg/Kg
Free acid (as HCl)..... ≤ 5 mg/Kg
Non volatile residue..... ≤ 5 mg/Kg
Assay (GC)..... ≥ 99.9 %
GC-ECD.Individual peak (Lindane) ..≤ 3 ng/l
Retention time trichlorobenzene to mirex
GC-NPD.Individual peak (Ethylparathion)≤ 3 ng/l
Retention time Atrazin to Coumaphos

Code	Size	Packaging	Notes
438681	1 l	Glass bottle	
438682	2.5 l	Glass bottle	

For chlorinated and nitrogenous compounds analysis

Chloroform > RS - PESTIPUR - For pesticide analysis - Stabilized with ethanol

RS

Description ..... Clear liquid
Colour ..... ≤ 10 hazen
Identification ..... Positive
Water ..... ≤ 0.01 %
Acidity (HCl) ..... ≤ 5 ppm
Not volatile residue..... ≤ 5 ppm
GC-ECD (Lindano) ..... ≤ 3 ng/l
GC-NPD (Ethylparation)..... ≤ 3 ng/l
Assay (GLC) ..... ≥ 99.9 %

Code	Size	Packaging	Notes
438651	1 l	Glass bottle	
438652	2.5 l	Glass bottle	

Chloroform > RS - SPECTROSOL - For optical spectroscopy - Stabilized with amylene

RS

Clear,colourless liq.appearance ..... Conform
Identification ..... Conform
Colour ..... ≤ 10 Apha
Refractive index at 20°C..... 1.444 - 1.448
Water content (K.F.) ..... ≤ 100 mg/Kg
Free acid (as HCl)..... ≤ 5 mg/Kg
Non volatile residue..... ≤ 5 mg/Kg
Stabilizer (Amylene)..... 5 - 50 mg/Kg
Carbon tetrachloride..... ≤ 100 mg/Kg
Dichloromethane..... ≤ 100 mg/Kg
UV transmittance at 250 nm ..... ≥ 50 %
UV transmittance at 260 nm ..... ≥ 85 %
UV transmittance at 280 nm ..... ≥ 95 %
Assay (GC)..... ≥ 99.9 %

Code	Size	Packaging	Notes
438591	1 l	Glass bottle	
438592	2.5 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Chloroform > RS - SPECTROSOL - For optical spectroscopy - Stabilized with ethanol

RS

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Density at 20° C ..... 1.479 ÷ 1.483
Refractive index at 20°C. 1.4461 ÷ 1.4491
Boiling point ..... 61.0 ÷ 61.5 ° C
Water (K.F.) ..... ≤100 ppm
Residue on evaporation ..... ≤5 ppm
Acidity (HCl) ..... ≤0.0005 meq/g
Alcalinity ..... ≤0.0002 meq/g
Assay (GLC) ..... ≥99.9 %
Fluorescence
at 254 nm ..... ≤2 ppb
at 365 nm ..... ≤2 ppb
U.V. Transmittance
at 245 nm ..... ≥15 %
at 250 nm ..... ≥50 %
at 260 nm ..... ≥90 %
at 275 nm ..... ≥98 %
Stabilized with ethyl alcohol .... 0.6 ÷ 1.0 %
Carbon tetrachloride ..... ≤ 100 ppm

Code	Size	Packaging	Notes
438664	1 l	Glass bottle	
438662	2.5 l	Glass bottle	

Chloroform > RS - Anhydrous - For analysis - Stabilized with amylene

RS

Refractive index at 20°C ..... 1.444 - 1.448
Water content (K.F.) ..... ≤ 50 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
Colour ..... ≤ 10 Hazen
Free acid (as HCl) ..... ≤ 5 mg/Kg
Assay (GC) ..... ≥ 99.95 %
Stabilizer (Amylene) ..... 5 - 50 mg/Kg
Carbon tetrachloride ..... ≤ 80 mg/Kg
Dichloromethane ..... ≤ 50 mg/Kg
1,2-dichloroethane ..... ≤ 10 mg/Kg

Code	Size	Packaging	Notes
P02410A10	200 ml	Bottle with septum	
P02410A16	1 l	Glass bottle	
P02410A21	2.5 l	Glass bottle	
P02410AT21	2.5 l	Glass bottle	On molecular sieves 4A

Chloroform > RS - Anhydrous - For analysis - Stabilized with ethanol

RS

Refractive index at 20°C ..... 1.444 - 1.448
Water content (K.F.) ..... ≤ 50 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) (without stabilizer) ... ≥ 99.95 %
Stabilizer (Ethanol) ..... 0.6 - 1 % m/m
Free acid (as HCl) ..... ≤ 5 mg/Kg
Carbon tetrachloride ..... ≤ 80 mg/Kg
Dichloromethane ..... ≤ 100 mg/Kg
1,2-dichloroethane ..... ≤ 10 mg/Kg

Code	Size	Packaging	Notes
P02410E10	200 ml	Bottle with septum	
P02410E16	1 l	Glass bottle	
P02410E21	2.5 l	Glass bottle	

Chloroform > RPE - For analysis - ISO - ACS - Reag.Ph.Eur. - Reag.USP - Stabilized ethanol

RPE

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification (I.R.) ..... Conform
Phosgene ..... Conform
Ready carbonizable substances ..... Conform
Acetone and aldeydes ..... Conform ACS
Acidità e cloruri ..... Conform ACS
Suit.for Dithizone col. ..... Conform ACS
Density at 20° C ..... 1.477 ÷ 1.481
Refractive index at 20°C. 1.4461 ÷ 1.4491
Boiling point ..... 60.5 ÷ 61.5 ° C
Water (K.F.) ..... ≤100 ppm
Residue on evaporation ..... ≤5 ppm
Acidity (HCl) ..... ≤5 ppm
Ethyl alcohol ..... 0.6 ÷ 1 %
Chloride ..... ≤0.4 ppm
Free chlorine ..... ≤0.1 ppm
Carbonyl Compounds (CO) ..... ≤5 ppm
Cu ..... ≤0.01 ppm
Fe ..... ≤0.1 ppm
Pb ..... ≤0.01 ppm
Zn ..... ≤0.05 ppm
Assay (GLC) ..... ≥99.9 %

Code	Size	Packaging	Notes
438613	1 l	Glass bottle	
438614	2.5 l	Glass bottle	
438612	5 l	Plastic tank	

Chloroform > RPE - For analysis - ISO - Stabilized with ethanol

RPE

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification (I.R.) ..... Conform
Free chlorine ..... Conform
Phosgene ..... Conform
Ready carbonizable substances ..... Conform
Alcohol miscibility ..... Complete
Benzene miscibility ..... Complete
Diethyl ether miscib. ..... Complete
Density at 20° C ..... 1.479 ÷ 1.483
Refractive index at 20°C. 1.4461 ÷ 1.4491
Boiling point ..... 61.0 ÷ 61.5 ° C
Water (K.F.) ..... ≤300 ppm
Residue on evaporation ..... ≤5 ppm
Acidity (HCl) ..... ≤5 ppm
Ethyl alcohol ..... 0.6 ÷ 1.0 %
Chloride ..... ≤0.4 ppm
Carbonyl Compounds (CO) ..... ≤5 ppm
Cu ..... ≤0.01 ppm
Fe ..... ≤0.1 ppm
Pb ..... ≤0.01 ppm
Zn ..... ≤0.05 ppm
Assay (GLC) ..... ≥99.9 %

Code	Size	Packaging	Notes
438601	1 l	Glass bottle	
438603	2.5 l	Glass bottle	
438607	35 kg	Aluminium can	
438606	250 kg	Metal drum	

Chloroform > RPE - For analysis - Stabilized with amylene

RPE

Clear,colourless liq.appearance ..... Conform	Refractive index at 20°C ..... 1.444 - 1.448	Non volatile residue ..... ≤ 10 mg/Kg	Dichloromethane ..... ≤ 50 mg/Kg
Identification ..... Conform	Water content (K.F.) ..... ≤ 100 mg/Kg	Stabilizer (Amylene) ..... 5 - 50 mg/Kg	1,2-dichloroethane ..... ≤ 10 mg/Kg
Colour ..... ≤ 10 Apha	Free acid (as HCl) ..... ≤ 5 mg/Kg	Carbon tetrachloride ..... ≤ 80 mg/Kg	Assay (GC) ..... ≥ 99.95 %

Code	Size	Packaging	Notes
438581	1 l	Glass bottle	
438582	2.5 l	Glass bottle	
P02405A28	5 l	Plastic tank	

Chloroform > ERBApharm - According to pharmacopoeia: BP - Stabilized with ethanol

ERBApharm

Description ..... Clear colourless liquid	Chloride ..... Conform BP	Related substances (CPG)	Total impurities ..... ≤ 1.0 %v/v
Identification ..... Positive	Aldehyde ..... Conform BP	Bromochloromethane ..... ≤ 0.5 %v/v	Ethanol ..... 1.0 ÷ 2.0 % (v/v)
Acidity or alkalinity ..... Conform BP	Density at 20° C ..... 1.474 ÷ 1.479	Methylene chlorure ..... ≤ 0.2 %v/v	Origin (BSE/TSE) ..... Synthesis
Foreign chlorin. comp. ..... Conform BP	Residue on evaporation ..... ≤40 ppm p/v	Carbon tetrachloride ..... ≤ 0.2 %v/v	Assay GC (without stabilizer) ..... ≥ 99.9 %
Free chlorine ..... Conform BP	Distillation range ..... Pass test**	Any single impurity ..... ≤ 0.2 %v/v	

Code	Size	Packaging	Notes
334351	1 l	Glass bottle	
334353	2.5 l	Glass bottle	
334356	25 kg	Plastic tank	
529301	200 l	Metal drum	
334354	250 kg	Metal drum	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Chloroform > RE - Pure - Stabilized with amylene

RE

Description ..... Clear colourless liquid	Density at 20°C ..... 1.478 ÷ 1.488	Acidity (HCl) ..... <10 ppm	Assay (GLC) ..... ≥ 99.9 %
Color ..... <10 APHA	Refractive index at 20°C. 1.4461 ÷ 1.4491	Residue on evaporation ..... ≤20 ppm	Amylene ..... ≤ 60 ppm
Identity (IR) ..... Positive	Boiling point ..... 61 ÷ 61.5 °C	Water (K.F.) ..... ≤ 300 ppm	

Code	Size	Packaging	Notes
528326	1 l	Glass bottle	
528328	2.5 l	Glass bottle	
528325	5 l	Plastic tank	
528329	25 l	Metal drum	
528327	200 l	Metal drum	

Chloroform > RE - Pure - Stabilized with ethanol

RE

Description ..... Clear colourless liquid	Density at 20°C ..... 1,478 ÷ 1,488	Acidity (HCl) ..... ≤ 50 ppm	Water (K.F.) ..... ≤ 0.03 %
Color ..... ≤ 10 APHA	Refractive index at 20°C. 1,4461 ÷ 1,4491	Ethyl alcohol ..... 0,6 ÷ 1 %	Assay (GLC) ..... ≥ 99 %
Identity (IR) ..... Positive	Boiling point ..... 61,0 ÷ 61,5 °C	Residue on evaporation ..... ≤ 20 ppm	Stab with 0,8 % ethanol

Code	Size	Packaging	Notes
508320	1 l	Glass bottle	
508321	5 l	Plastic tank	
P02402E40	10 l	Metal drum	
P02402E48	25 l	Metal drum	
508322	200 l	Metal drum	

Chloroform-d

• Cloroformio-d • Chloroforme-d • Cloroformo-d • Chloroform-d1

Synonym:
Deuterochloroform

$CDCl_3$
Molecular Weight: 120,37
CAS: 865-49-6
EEC-N: 200-663-8

Classification transport
ONU: 1888
Transport Hazard class: 6.1
Packing group III



Danger
H302-H331-H315-H319-H351-H361d-H372-HEU301
P280-P304+P340-P305+P351+P338-P308+P313-
P330-P362+P364-P403+P233

Chloroform-d > RS - For NMR - min 99.95%

RS

Code	Size	Packaging	Notes
P5130	10 x 0.6 ml	Glass ampoule	

For specifications, contact our customer service for a certificate of analysis

Chloroform-d > RS - For NMR - min 99.95% - Stabilized with 0.12% Ag

RS

Code	Size	Packaging	Notes
P5505	100 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Chloroform-d > RS - For NMR - min 99.8%

RS

Code	Size	Packaging	Notes
P5119	10 x 0.75 ml	Glass ampoule	
P5115	25 ml	Glass bottle	
P5116	100 ml	Glass bottle	
P5117	500 ml	Glass bottle	
P5118	1 l	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Chloroform-d > RS - For NMR - min 99.8% - Stabilized with 0.12% Ag

RS

Code	Size	Packaging	Notes
P5325	100 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Chloroform-d + 0.03% TMS

• Cloroformio-d + 0.03% TMS • Chloroforme-d + 0.03% TMS • Cloroformo-d + 0.03% TMS
• Chloroform-d + 0.03% TMS

Synonym:
Deuterochloroform

$CDCl_3$
Molecular Weight: 120,37
CAS: 865-49-6
EEC-N: 200-663-8

Classification transport
ONU: 1888
Transport Hazard class: 6.1
Packing group III



Danger
H302-H331-H315-H319-H351-H361d-H372-HEU301
P280-P304+P340-P305+P351+P338-P308+P313-
P330-P362+P364-P403+P233

Chloroform-d + 0.03% TMS > RS - For NMR - min 99.8%

RS

Code	Size	Packaging	Notes
P5006	100 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

4-Chlorophenol

• 4-Clorofenolo • p-Chlorophénol • p-Clorofenol • p-Chlorophenol

ClC_6H_4OH
Molecular Weight: 128,56
CAS: 106-48-9
EEC-N: 203-402-6

Classification transport
ONU: 2020
Transport Hazard class: 6.1
Packing group III

Warning
H302-H312-H332-H411
P261-P264-P271-P280h-P301+P312a-P304+P340

4-Chlorophenol > RPE - For analysis - Reag. Ph. Eur. — RPE

Description ..... Yellowish crystals Identification ..... Positive Water ..... ≤ 0.3 % Assay (GLC) ..... ≥ 99 %

Code	Size	Packaging	Notes
438504	100 g	Glass bottle	

Chloroplatinic acid hexahydrate ▶ Hexachloroplatinic acid hexahydrate

n-Chloro-p-toluenesulfonamide sodium salt ▶ Chloramine T sodium salt

Chromate standard solution

• Cromati standard soluzione • Chromate standard solution • Cromato, solución patrón • Chromatierte Standardlösung

Danger
H340-H350-H412-HEU203-HEU208-HA26
P201-P273-P280-P308+P313-P405-P501a

Chromate standard solution > RS - Standard solution for ion chromatography — RS

Code	Size	Packaging	Notes
503241	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503243	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Chromic anhydride ▶ Chromium (VI) oxide

Chromium standard solution

• Cromo standard soluzione • Chrome standard solution • Cromo, solución patrón • Chrome-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Warning
H290
P234-P390-P406

Chromium standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2 — RS

Code	Size	Packaging	Notes
615001002	100 ml	Plastic bottle	A 0.1 % solution Ref Ph.Eur 5001002
615001000	1 l	Plastic bottle	A 100 ppm solution Ref Ph.Eur 5001000

Chromium standard solution > RS - Standard solution for ICP-MS — RS

Code	Size	Packaging	Notes
505567	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505568	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505569	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Product specifications are subject to changes. Please visit our website for updates.

Chromium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503521	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503523	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503525	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503527	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Chromium standard solution > RS - Standard solution for AAS

RS

Description ..... Green-grey clear liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
504195	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
507485	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
E497501	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Chromium standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Orange clear liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
440641		Plastic ampoule	conc. 1.000 ppm Matrix: Hydrochloric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Chromium (III) chloride hexahydrate

• Cromo cloruro ico esaidrato • Chrome (III) chlorure hexahydrate • Cromo (III) cloruro hexahidrato
• Chrom (III) chloridhexahydrat

Synonym:
• *Chromium trichloride hexahydrate*
• *Hexaaquachromium (III) chloride*

$CrCl_3.6H_2O$
Molecular Weight: 266,45
CAS: 10060-12-5
EEC-N: 233-038-3

Warning
H302
P264-P270-P301+P312a-P330-P501a

Chromium (III) chloride hexahydrate > RPE - For analysis

RPE

Description ..... Polvere cristallina verde scuro
Identification ..... Positive
pH solution 5% ..... 2.0 ÷ 3.5
Non precipit. con NH4OH (SO4) ..... ≤ 0.2 %
Ammonium ..... ≤ 100 ppm
Sulphate ..... ≤ 100 ppm
Fe ..... ≤ 100 ppm
Pb ..... ≤ 50 ppm
Assay (oxidimetric) ..... ≥ 95.0 %

Code	Size	Packaging	Notes
440724	100 g	Glass bottle	
440727	1 kg	Plastic bottle	

Chromium (III) nitrate nonahydrate

• Cromo nitrato ico nonaidrato • Chrome (III) nitrate nonahydrate • Cromo (III) nitrato nonahidrato • Chrom (III) nitrat Nonahydrat

$Cr(NO_3)_3.9H_2O$
Molecular Weight: 400,15
CAS: 7789-02-8

Classification transport
ONU: 2720
Transport Hazard class: 5.1
Packing group III

Warning
H272-H332-H317-H411
P210-P220-P261-P271-P280-P304+P340

Chromium (III) nitrate nonahydrate > RPE - For analysis

RPE

Description	Violet crystals	Subst. not ppt NH4OH	≤ 0.2 %	Cu	≤ 10 ppm	Ni	≤ 50 ppm
Identification	Positive	Sulphate	≤ 50 ppm	Fe	≤ 200 ppm	Pb	≤ 20 ppm
pH sol. 5% at 25° C	2.0 ÷ 3.0	Ca	≤ 50 ppm	Mg	≤ 50 ppm	Zn	≤ 10 ppm
Ammonium	≤ 10 ppm	Cd	≤ 10 ppm	Mn	≤ 10 ppm	Assay (oxidimetric)	≥ 12.5 % Cr
Chloride	≤ 20 ppm	Co	≤ 10 ppm	Na	≤ 50 ppm		

Code	Size	Packaging	Notes
440775	250 g	Glass bottle	
440776	1 kg	Glass bottle	

Chromium (III) oxide

• Cromo sesquiossido • Chrome (III) oxyde • Cromo (III) óxido • Chrom (III) oxid

Synonym:
Chromia

Cr_2O_3
Molecular Weight: 151,99
CAS: 1308-38-9
EEC-N: 215-160-9

Chromium (III) oxide > RPE - For analysis

RPE

Description	Green powder	Identification	Positive	Water solubility	≤ 0.2 %	Assay (oxidimetric)	≥ 99 %

Code	Size	Packaging	Notes
440825	250 g	Plastic bottle	
440827	1 kg	Plastic bottle	
440823	25 kg	Plastic bucket	

Chromium (III) potassium sulfate dodecahydrate

• Cromo (III) di potassio solfato dodecaidrato • Chrome (III) potassium sulfate dodécahydraté
• Cromo (III) de potasio sulfato dodecahidrato • Chrom (III) kaliumsulfatdodecahydrat

Synonym:
• *Chrome alum*
• *Potassium chromium (III)sulfate dodecahydrate*

$CrK(SO_4)_2.12H_2O$
Molecular Weight: 499,39
CAS: 7788-99-0

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Chromium (III) potassium sulfate dodecahydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description	Purple crystals	Ammonium	≤100 ppm	Al	≤200 ppm
Identification	Positive	Chloride	≤20 ppm	Fe	≤100 ppm
Water-insoluble matter	≤100 ppm	Heavy metals (Pb)	≤100 ppm	Assay (oxidimetric)	98.0 ÷ 102.0 %

Code	Size	Packaging	Notes
440876	100 g	Plastic bottle	
440877	1 kg	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Chromium (III) sulfate

• Cromo solfato ico • Chrome (III) sulfate • Cromo (III) sulfato • Chrom (III) sulfat

Synonym:
• *Chromium(III) sulfate hydrate*
• *Chromium trisulfate*

$Cr_2(SO_4)_3.nH_2O$
Molecular Weight: 392,18
CAS: 15244-38-9

Classification transport
ONU: 3260
Transport Hazard class: 8
Packing group II



Danger
H302-H312-H332-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P362+P364

Chromium (III) sulfate > RPE - For analysis

RPE

Description	Green Crystal	Assay	≥ 98 %	Na	≤ 0.01 %
Identification	Positive	Chloride	≤ 0.01 %	Fe	≤ 0.005 %

Code	Size	Packaging	Notes
440955	250 g	Plastic bottle	
440957	1 kg	Plastic bottle	

Chromium (VI) oxide

• Cromo (VI) ossido • Chrome (VI) oxyde • Cromo (VI) óxido • Chrom (VI) oxid

Synonym:
Chromic anhydride

CrO_3
Molecular Weight: 99,99
CAS: 1333-82-0
EEC-N: 215-607-8

Classification transport
ONU: 1463
Transport Hazard class: 5.1
Packing group II



Danger
H271-H301-H311-H330-H314-H334-H317-H340-H350-H361f-H372-H410-HA26
P210-P280-P283-P284-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338-P342+P311a-P403+P233

Chromium (VI) oxide > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description	Purple crystals	Chloride	≤50 ppm	Fe,Al,Ba	≤300 ppm
Identification	Positive	Nitrate	≤500 ppm	Na	≤0.2 %
Water-insoluble matter	≤100 ppm	Sulphate	≤50 ppm	Assay (oxidimetric)	≥98.0 %

Code	Size	Packaging	Notes
421735	250 g	Plastic bottle	

Chromium (VI) oxide solution

• Cromo (VI) ossido soluzione • Chrome (VI) oxyde solution • Cromo (VI) óxido solución
• Chrom (VI) oxidlösung

Synonym:
Chromic anhydride

CrO_3
Molecular Weight: 99,99
CAS: 1333-82-0

Classification transport
ONU: 1755
Transport Hazard class: 8
Packing group II



Danger
H271-H301-H312-H330-H314-H334-H317-H340-H350-H361f-H335-H372-H410-HA26
P210-P280-P283-P284-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338-P342+P311a-P403+P233

Chromium (VI) oxide solution > RE - Pure

RE

Description	Orange clear liquid	Chloride	≤500 ppm	Assay (oxidimetric)	294.6 ÷ 300.6 g/l
Identification	Positive	Sulphate	≤0.1 %		

Code	Size	Packaging	Notes
317511	2.5 l	Glass bottle	

Chromotropic acid disodium salt

• Acido cromotropico sale bisodico • Acide chromotropique sel disodique • Acido cromotrópico sal disódica
• Chromotropsäure Dinatriumsalz-Dihydrat

Synonym:
• *4,5-Dihydroxynaphtalene-2,7-disulfonic acid disodium salt*
• *1,8-Dihydroxynaphtalene-3,6-disulfonic acid disodium salt*

$C_{10}H_6Na_2O_8S_2.2H_2O$
Molecular Weight: 400,29
CAS: 5808-22-0
EEC-N: 204-972-9

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Chromotropic acid disodium salt > RPE - For analysis

RPE

Description ..... Whitish powder
Identification ..... Positive
Water-insoluble matter ..... ≤ 0.02 %
Nitrate sinsitivity ..... Conform
Formaldehyde sensit. ..... Conform
Assay (acidimetric) ..... ≥ 97.5 %

Code	Size	Packaging	Notes
404872	25 g	Glass bottle	

Chrysoidine Y

• Crisoidina Y • Chrysoïdine Y • Crisoidina Y • Chrysoidin Y

Synonym:
Basic orange 2,4-Phenylazo-m-phenylene diamine monohydrochloride

$C_{12}H_{13}ClN_4$
Molecular Weight: 248,71
CAS: 532-82-1
EEC-N: 208-545-8

Classification transport
ONU: 3143
Transport Hazard class: 6.1
Packing group III

Danger
H302-H315-H318-H341-H410
P280-P305+P351+P338-P310a-P308+P313-P330-P362+P364

Chrysoidine Y > RS - For microscopy - C.I. 11270

RS

Description ..... Dark red powder
Identification ..... Positive

Code	Size	Packaging	Notes
440572	25 g	Glass bottle	

Dye for bacteriology and botanic

Cinchonine

• Cinconina • Cinchonine • Cinconina • Cinchonin

$C_{19}H_{22}N_2O$
Molecular Weight: 294,4
CAS: 118-10-5
EEC-N: 204-234-6

Warning
H302-H332
P261-P264-P271-P301+P312a-P304+P340-P501a

Cinchonine > RPE - For analysis

RPE

Description ..... White powder
Identification ..... Positive
Ready carbonizable substances ..... Conform
Melting point ..... 260 ÷ 265 °C
Specific optical rotation (c=5 in Ethano ..... +225 ÷ +230 °
Loss on drying ..... ≤ 1 %
Total chlorine ..... ≤ 300 ppm
Dil. H2SO4-ins. matter ..... ≤ 200 ppm
Heavy metals (Pb) ..... ≤ 10 ppm
Residue on ignition ..... ≤ 0.1 %
Total sulphur ..... ≤ 50 ppm
Fe ..... ≤ 10 ppm
Assay (non-aqueos medium) ..... ≥ 99 %

Code	Size	Packaging	Notes
437251	10 g	Glass bottle	

Citral

• Citrale • Citral • Citral • Citral

Synonym:
3,7-Dimethyl-2,6-octadienal

$C_{10}H_{16}O$
Molecular Weight: 152,24
CAS: 5392-40-5
EEC-N: 226-394-6

Warning
H315-H317
P261-P264-P280g-P362+P364-P333+P313-P501a

Citral > RE - Pure

RE

Description ..... Yellow liquid
Identification ..... Positive
Density at 20° C ..... 0.886 ÷ 0.890
Refractive index at 20°C ..... 1.4870 ÷ 1.4910
Residue on ignition ..... ≤500 ppm
Assay (GLC) ..... ≥97 %

Code	Size	Packaging	Notes
437401	25 ml	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Citric acid anhydrous

• Acido citrico anidro • Acide citrique anhydre • Acido cítrico anhidro • Zitronensäure wasserfrei

$C_6H_8O_7$
Molecular Weight: 192,13
CAS: 77-92-9
EEC-N: 201-069-1

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Citric acid anhydrous > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-BP-DAB-JP

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Ready carbonizable substances.....Conform
USP-NF
Organic volatile impurities Conform USP-NF
Water (K.F.) ..... ≤ 1.0 %
Oxalic acid.....≤ 360 ppm
Sulphated ash..... ≤ 0.1 %
Sulphate.....≤ 150 ppm
Heavy metals (Pb).....≤ 10 ppm
Assay (acidimetric) .....99.5 ÷ 100.5 % s.s.
Origin (BSE/TSE).....Vegetable
Residual solvents (Current ICH).....Conform

Code	Size	Packaging	Notes
302486	500 g	Plastic bottle	
302487	1 kg	Plastic bottle	
302485	5 kg	Plastic tank	
302488	25 kg	Plastic bucket	
302484	50 kg	Fibre drum	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Citric acid monohydrate

• Acido citrico monoidrato • Acide citrique monohydraté • Acido cítrico monohidrato • Zitronensäure monohydrat

$C_6H_8O_7$
Molecular Weight: 192,12
CAS: 5949-29-1

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Citric acid monohydrate > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611021000	100 g	Glass bottle	Ref Ph.Eur 1021000

Citric acid monohydrate > RPE - For analysis - ISO - ACS - Reag.Ph.Eur. - Reag.USP

RPE

Description .....White crystals
Identification ..... Positive
Water-insoluble matter .....≤50 ppm
Chloride.....≤10 ppm
Phosphate .....≤10 ppm
Sulphate.....≤20 ppm
Fe .....≤3 ppm
Pb .....≤2 ppm
Substances darkened by sulphuric acid ..... Conform
Appearance of solution .....Conform
Organic volatile impurities.....Conform
Water (K.F.) .....7.5 ÷ 8.8 %
Oxalic acid.....≤350 ppm
Sulphated ash.....≤200 ppm
Heavy metals (Pb).....≤10 ppm
As .....≤3 ppm
Assay (acidimetric) .....99.5 ÷ 100.5 % s.s.

Code	Size	Packaging	Notes
403725	250 g	Plastic bottle	
403727	1 kg	Plastic bottle	
403721	5 kg	Plastic tank	
403724	25 kg	Plastic bucket	
403722	50 kg	Fibre drum	

Citric acid monohydrate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-BP-DAB

ERBApharm

Description .....White crystals
Identification (I.R.)..... Positive
Clarity of solution.....Conform
Color of solution .....Conform
Appearance of solution .....Conform Ph.Eur.
Ready carbonizable substances.....Conform
USP-NF
Water (K.F.) .....7.5 ÷ 9.0 %
Oxalic acid.....≤ 360 ppm
Sulphated ash..... ≤ 0.1 %
Sulphate.....≤ 150 ppm
Heavy metals (Pb).....≤ 10 ppm
Assay (acidimetric) .....99.5 ÷ 100.5 % s.s.
Residual solvents (Current ICH).....Conform

Code	Size	Packaging	Notes
302557	1 kg	Plastic bottle	
302559	5 kg	Plastic tank	
302551	25 kg	Plastic bucket	
302554	50 kg	Fibre drum	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Citric acid monohydrate > ERBApharm - Powder - According to pharmacopoeia: Ph.Eur.-USP-FU-BP-DAB

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Ready carbonizable substances.....Conform
USP-NF
Water (K.F.) .....7.5 ÷ 9.0 %
Oxalic acid.....≤ 360 ppm
Sulphated ash ..... ≤ 0.1 %
Sulphate .....≤ 150 ppm
Heavy metals (Pb).....≤ 10 ppm
Assay (acidimetric) .....99.5 ÷ 100.5 % s.s.
Origin (BSE/TSE).....Vegetable
Residual solvents (Current ICH).....Conform

Code	Size	Packaging	Notes
302507	1 kg	Plastic bottle	
302509	5 kg	Plastic tank	
302501	25 kg	Plastic bucket	
302504	50 kg	Fibre drum	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Clayton's yellow

• Giallo Clayton • Jaune Clayton • Amarillo Clayton • Titangelb

Synonym:
• *Thiazole Yellow G*
• *Direct Yellow 9*

$C_{28}H_{19}N_5O_6S_4Na_2$
Molecular Weight: 695,73
CAS: 1829-00-1
EEC-N: 217-377-4

Clayton's yellow > RPE - For analysis - C.I. 19540

RPE

Description .....Yellow brown powder
Identification ..... Positive
pH range ..... 1.2 - 13.2
Loss on drying .....≤10 %
Mg sensitivity.....≤0.5 µg/ml
Residue on ignition.....18.4 ÷ 22.4 %

Code	Size	Packaging	Notes
453518	5 g	Glass bottle	
453519	25 g	Glass bottle	

Dye for microscopy. Indicator acid - base (pH 12.0 ÷ 13.0). Fluorescence indicator. For the determination of magnesium

Cobalt standard solution

• Cobalto standard soluzione • Cobalt standard solution • Cobalto, solución patrón • Kobalt-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III



Danger
H290-H314-H350i-H412-HEU208-HA26
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Cobalt standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615004300	1 l	Plastic bottle	A 100 ppm solution Ref Ph.Eur 5004300

Cobalt standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505562	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505565	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505563	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Cobalt standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503511	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503513	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503515	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503517	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Cobalt standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
507533	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507484	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
E497495	100 ml	Bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
E497491	500 ml	Bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Cobalt standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
439131		Plastic ampoule	conc. 1.000 ppm Matrix: Nitric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Cobalt (II) acetate tetrahydrate

• Cobalto acetato oso tetraidrato • Cobalt (II) acétate tétrahydraté • Cobalto (II) acetato tetrahidrato
• Kobalt (II) acetattetrahydrat

Synonym:
Cobaltous acetate tetrahydrate

$Co(CH_3COO)_2.4H_2O$
Molecular Weight: 249,08
CAS: 6147-53-1
EEC-N: 200-755-8



Danger
H302-H334-H317-H351
P261-P280-P284-P304+P340-P308+P313-P342+P311a

Cobalt (II) acetate tetrahydrate > RPE - For analysis

RPE

Description ..... Brown-red crystal Identification ..... Positive Chloride..... ≤ 0.005 % Sulfate..... ≤ 0.005 % Cu..... ≤ 0.0005 % Fe ..... ≤ 0.001 % Assay ..... 98 - 100 % Water insolubles..... ≤ 0.01 %

Code	Size	Packaging	Notes
439154	100 g	Glass bottle	
439155	250 g	Glass bottle	
439156	1 kg	Glass bottle	

Cobalt (II) ammonium sulfate hexahydrate

• Cobalto ammonio solfato oso esaidrato • Cobalt (II) ammonium sulfate hexahydraté
• Cobalto (II) y amonio sulfato hexahidrato • Kobalt (II) ammoniumsulfathexahydrat

Synonym:
Ammonium cobalt(II) sulfate hexahydrate

$Co(NH_4)_2(SO_4)_2.6H_2O$
Molecular Weight: 395,23
CAS: 13596-46-8
EEC-N: 237-043-1



Warning
H302-H332-H315-H319-H317-H335
P261-P271-P304+P340-P305+P351+P338-P337+P313-P403+P233

Cobalt (II) ammonium sulfate hexahydrate > RPE - For analysis

RPE

Description ..... Red crystal. powder Identification ..... Positive pH sol. 5% at 25° C ..... 4.0 ÷ 7.0 Assay (complexometric) ..... ≥98 %

Code	Size	Packaging	Notes
439204	100 g	Glass bottle	
439207	1 kg	Glass bottle	

Cobalt (II) chloride hexahydrate

• Cobalto cloruro oso esaidrato • Cobalt (II) chlorure hexahydrate • Cobalto (II) cloruro hexahidrato • Kobalt (II) chloridhexahydrat

Synonym: *Cobaltous chloride hexahydrate*

$CoCl_2.6H_2O$
Molecular Weight: 237,93
CAS: 7791-13-1
EEC-N: 231-589-4

Classification transport
ONU: 3288
Transport Hazard class: 6.1
Packing group III



Danger
H302-H334-H317-H350i-H410-HA26
P261-P280-P284-P304+P340-P308+P313-P342+P311a

Cobalt (II) chloride hexahydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Red-violet crystals
Identification ..... Positive
Water-insoluble matter .....≤ 100 ppm
Nitrate .....≤ 100 ppm
Sulphate .....≤ 100 ppm
Ca .....≤ 50 ppm
Cu .....≤ 20 ppm
Fe .....≤ 50 ppm
K .....≤ 100 ppm
Mg .....≤ 50 ppm
Na .....≤ 500 ppm
Ni ..... ≤ 0.1 %
Zn .....≤ 300 ppm
Assay (complexometric) ....98.0 ÷ 102.0 %

Code	Size	Packaging	Notes
439355	250 g	Plastic bottle	
439357	1 kg	Plastic bottle	
439353	25 kg	Sack	

Cobalt (II) chloride in solution

• Cobalto (II) cloruro in soluzione • Cobalt (II) chlorure en solution • Cobalto (II) cloruro en solución • Kobalt (II) chlorid in Lösung

$CoCl_2$
Molecular Weight: 129,84
CAS: 7646-79-9

Cobalt (II) chloride in solution > RS - For analysis according to JP

RS

Code	Size	Packaging	Notes
616001028	100 ml	Plastic bottle	Cobaltous Chloride CS

Cobalt (II) chloride in solution > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
616001057	500 ml	Plastic bottle	Cobaltous Chloride CS

Cobalt (II) nitrate hexahydrate

• Cobalto nitrato oso esaidrato • Cobalt (II) nitrate hexahydrate • Cobalto (II) nitrato hexahidrato • Kobalt (II) nitrathexahydrat

Synonym: *Cobaltous nitrate hexahydrate*

$Co(NO_3)_2.6H_2O$
Molecular Weight: 291,04
CAS: 10026-22-9
EEC-N: 233-402-1

Classification transport
ONU: 1477
Transport Hazard class: 5.1
Packing group II



Danger
H272-H302-H334-H317-H341-H350i-H360F-H410-HA26
P210-P261-P280-P284-P304+P340-P342+P311a

Cobalt (II) nitrate hexahydrate > RS - For enviromental analysis - ACS

RS

Description ..... Red crystals
Identification ..... Positive
Water-insoluble matter .....≤100 ppm
Ammonium ..... ≤0.2 %
Chloride .....≤20 ppm
Sulphate .....≤50 ppm
Ca .....≤50 ppm
Cu .....≤20 ppm
Fe .....≤10 ppm
K .....≤100 ppm
Mg .....≤50 ppm
Na .....≤500 ppm
Ni ..... ≤0.15 %
Pb .....≤20 ppm
Zn .....≤100 ppm
Assay (complexometric) ....98.0 ÷ 102.0 %

Code	Size	Packaging	Notes
439504	100 g	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Cobalt (II) nitrate hexahydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description	Red brick crystals	Sulphate	≤ 50 ppm	K	≤ 100 ppm	Pb	≤ 20 ppm
Identification	Positive	Ca	≤ 50 ppm	Mg	≤ 50 ppm	Zn	≤ 100 ppm
Water-insoluble matter	≤ 100 ppm	Cu	≤ 20 ppm	Na	≤ 500 ppm	Assay (complexometric)	98.0 ÷ 102.0 %
Chloride	≤ 20 ppm	Fe	≤ 10 ppm	Ni	≤ 0.15 %		

Code	Size	Packaging	Notes
439455	250 g	Plastic bottle	
439457	1 kg	Plastic bottle	

Cobalt (II) sulfate heptahydrate

• Cobalto solfato oso eptaidrato • Cobalt (II) sulfate heptahydrate • Cobalto (II) sulfato heptahidrato • Kobalt (II) sulfatheptahydrat

Synonym: *Cobaltous sulfate heptahydrate*

$CoSO_4.7H_2O$
Molecular Weight: 281.10
CAS: 10026-24-1
EEC-N: 233-334-2

Danger
H302-H334-H317-H350-H410-HA26
P261-P280-P284-P304+P340-P308+P313-P342+P311a

Cobalt (II) sulfate heptahydrate > RPE - For analysis

RPE

Description	Red crystals	Chloride	≤20 ppm	Cu	≤10 ppm	Zn	≤20 ppm
Identification	Positive	Water-insoluble matter	≤30 ppm	Fe	≤10 ppm	Assay (complexometric)	≥99 %
pH sol. 5% at 25° C	3.5 ÷ 4.5	Nitrate	≤50 ppm	Ni	≤500 ppm		
Ammonium	≤100 ppm	Subst. not ppt. (NH4)2S	≤0.1 %	Pb	≤10 ppm		

Code	Size	Packaging	Notes
439705	250 g	Plastic bottle	

Congo red

• Rosso Congo • Rouge Congo • Rojo Congo • Kongo-Rot

Synonym: *Direct Red 28*

$C_{32}H_{22}N_6Na_2O_6S_2$
Molecular Weight: 696,66
CAS: 573-58-0
EEC-N: 209-358-4

Danger
H350-H361d-HA26
P201-P202-P280-P308+P313-P405-P501a

Congo red > RPE - For analysis - C.I. 22120

RPE

Description	Red brown powder	Identification	Positive	Colour change	blue red	pH range	3.0 - 5.2

Code	Size	Packaging	Notes
476762	25 g	Glass bottle	
476764	100 g	Plastic bottle	

Dye for microscopy (histology). Indicator acid - base

Congo red paper

• Cartina rosso Congo • Papier rouge Congo • Rojo Congo papel • Kongorotes Papier

Danger
H350-HA26
P201-P202-P280-P308+P313-P405-P501a

Congo red paper > RS - For pHmetry

RS

Code	Size	Packaging	Notes
435220000	1 roll	Dispenser	Congo red paper, Color change: red --> Blue, Change pH 5.0-->3.0

Congo red solution

• Rosso Congo soluzione • Rouge Congo solution • Rojo Congo solución • Rote Kongo-Lösung

Synonym:
Direct Red 28

$C_{32}H_{22}N_6Na_2O_6S_2$
Molecular Weight: 696,66
CAS: 573-58-0

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group III



Danger
H226-H350-HA26
P210-P241-P280-P303+P361+P353-P308+P313-P403+P235

Congo red solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611022001	100 ml	Plastic bottle	Ref Ph.Eur 1022001

Change colour: pH 3.0 (blue) to pH 5.0 (pink)

Coomassie brilliant blue R 250

• Blu Coomassie R250 • Bleu de Coomassie brillant R 250 • Azul Coomassie R 250
• Helles Coomassie Blue R 250

Synonym:
• *Brilliant Blue R*
• *Acid Blue 83*

$C_{45}H_{44}N_3NaO_7S_2$
Molecular Weight: 825,99
CAS: 6104-59-2
EEC-N: 228-060-5



Warning
H302-H312-H332
P261-P264-P271-P280h-P301+P312a-P304+P340

Coomassie brilliant blue R 250 > RS - For microscopy - C.I. 42660

RS

Description ..... Dark violet powder Identification ..... Positive E (1%/1cm lambda max) .....≥ 300

Code	Size	Packaging	Notes
428642	25 g	Glass bottle	

Dye for histochemistry

Copper electrolytic rebaked, sheet

• Rame elettrolitico lastra ricotto • Cuivre électrolytique recuit, feuilles • Cobre electrolítico recocido, hojas • Kupferelektrolyt überbacken, Blatt

Cu
Molecular Weight: 63,55
CAS: 7440-50-8
EEC-N: 231-159-6

H413
P273-P501a

Copper electrolytic rebaked, sheet > RPE - For analysis

RPE

Description ..... Metallic sheet Identification ..... Positive Assay .....≥99.8 %

Code	Size	Packaging	Notes
475215	250 g	Box	

~ 0,15 mm thickness

Copper electrolytic, turnings

• Rame elettrolitico, tornitura • Cuivre électrolytique, tournures • Cobre electrolítico, virutas • Kupferelektrolyt, Späne

Cu
Molecular Weight: 63,55
CAS: 7440-50-8
EEC-N: 231-159-6

H413
P273-P501a

Copper electrolytic, turnings > RPE - For analysis

RPE

Description ..... Trucioli metallici Identification ..... Positive Assay .....≥99 %

Code	Size	Packaging	Notes
475305	250 g	Carton box	
475307	1 kg	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Copper electrolytic, wire

• Rame elettrolitico, filo • Cuivre électrolytique, fils • Cobre electrolítico, hilos • Kupferelektrolyt, Draht

Cu
Molecular Weight: 63,55
CAS: 7440-50-8
EEC-N: 231-159-6

H413
P273-P501a

Copper electrolytic, wire > RPE - For analysis

RPE

Description ..... Filo
Identification ..... Positive
HNO3-insoluble matter ..... ≤300 ppm
Sn + Sb ..... ≤200 ppm
Ag ..... ≤100 ppm
As ..... ≤5 ppm
Fe ..... ≤50 ppm
Mn ..... ≤10 ppm
P ..... ≤10 ppm
Pb ..... ≤200 ppm
Assay ..... ≥99.9 %

Code	Size	Packaging	Notes
475185	250 g	Box	
475187	1 kg	Box	

~ 1 mm diameter

Copper reduced, powder

• Rame ridotto, polvere • Cuivre réduit, poudre • Cobre reducido, polvo • Kupfer reduziert, Pulver

Cu
Molecular Weight: 63,55
CAS: 7440-50-8
EEC-N: 231-159-6

H413
P273-P501a

Copper reduced, powder > RPE - For analysis - Reag. Ph. Eur.

RPE

Description ..... Red-brown metallic powder
Identification ..... Positive
Assay ..... ≥ 98.5 % (Cu)

Code	Size	Packaging	Notes
475334	100 g	Glass bottle	
475337	1 kg	Glass bottle	

Copper standard solution

• Rame standard soluzione • Cuivre standard solution • Cobre, solución patrón • Kupfer-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Warning
H290
P234-P390-P406

Copper standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615001100	100 ml	Plastic bottle	A 0.1 % solution Ref Ph.Eur 5001100

Copper standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505577	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505578	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505579	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Copper standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503541	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503543	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503545	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503547	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Copper standard solution > RS - Standard solution for AAS

RS

Description ..... Clear blue liquid Identification ..... Positive Titration factor ..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
504545	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507478	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
E497615	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
E497611	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Copper standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Green clear liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
475151		Plastic ampoule	conc. 1.000 ppm Matrix: Nitric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Copper standard solution > RS - Quality control standard solution for AAS (graphite furnace)

RS

Code	Size	Packaging	Notes
504361	50 ml	Plastic bottle	conc. 10 ppb - Matrix: 2% Nitric acid

Copper (I) chloride

• Rame cloruro oso • Cuivre (I) chlorure • Cobre (I) cloruro • Kupfer (I) chlorid

Synonym:
- *Cuprous chloride*
- *Copper monochloride*

CuCl
Molecular Weight: 98,99
CAS: 7758-89-6
EEC-N: 231-842-9

Classification transport
ONU: 2802
Transport Hazard class: 8
Packing group III



Danger
H302-H312-H315-H318-H410
P264-P305+P351+P338-P310a-P330-P362+P364-P332+P313

Copper (I) chloride > RPE - For analysis

RPE

Description ..... green-dark powder Identification ..... Positive
Acid insoluble ..... ≤ 0.5 % Sulphate ..... ≤ 0.5 %
Fe ..... ≤ 0.02 % Titolo (manganometric) ..... ≥ 97 %

Code	Size	Packaging	Notes
475605	250 g	Plastic bottle	
475607	1 kg	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Copper (I) chloride solution 7% in ammonia

• Rame cloruro oso soluzione 7% in ammoniaca • Cuivre (I) chlorure solution 7% dans l'ammoniaque
• Cobre (I) cloruro solución 7% en amonio hidróxido • Kupfer (I) chloridlösung 7% in Ammoniak

Synonym:
• *Cuprous chloride*
• *Copper monochloride*

CuCl
Molecular Weight: 98,99
CAS: 7758-89-6



Danger
H314-H335-H411
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Copper (I) chloride solution 7% in ammonia > RS - For gas analysis according to Orsat — RS

Description ..... Blue liquid | Identification ..... Positive | Density at 20° C ..... ~ 1.08 | Assay ..... 6.5 ÷ 7.5 %

Code	Size	Packaging	Notes
E475632	1 l	Bottle	

Stabilized electrolytic copper

Copper (I) iodide

• Rame ioduro oso • Cuivre (I) iodure • Cobre (I) ioduro • Kupfer (I) iodid

Synonym:
Cuprous iodide

CuI
Molecular Weight: 190,45
CAS: 7681-65-4
EEC-N: 231-674-6



Danger
H302-H315-H318-H317-H372-H400-H411
P260-P264-P305+P351+P338-P310a-P330-P362+P364

Copper (I) iodide > RE - Pure — RE

Description ..... Polvere nocciola ammass. | Identification ..... Positive | Assay (oxidimetric) ..... ≥98 %

Code	Size	Packaging	Notes
364611	25 kg	Plastic bucket	

Copper (I) iodide > RE - Pure - For fibers — RE

Description ..... Polvere nocciola
Identification ..... Positive
Chloride ..... ≤ 0.05 %
Sulphate ..... ≤ 0.05 %
Cd ..... ≤ 50 ppm
Fe ..... ≤ 50 ppm
K ..... ≤ 0.05 %
Ni ..... ≤ 50 ppm
Pb ..... ≤ 100 ppm
Zn ..... ≤ 50 ppm
Assay (oxidimetric) ..... ≥ 98 %

Code	Size	Packaging	Notes
364637	1 kg	Plastic bottle	
364631	25 kg	Plastic bucket	

Copper (II) acetate hydrate

• Rame acetato ico idrato • Cuivre (II) acétate hydrate • Cobre (II) acetato hidrato • Kupfer (II) acetathydrat

Synonym:
Cupric acetate monohydrate

$Cu(CH_3COO)_2.H_2O$
Molecular Weight: 199,65
CAS: 6046-93-1
EEC-N: 205-553-3



Warning
H302
P264-P270-P301+P312a-P330-P501a

Copper (II) acetate hydrate > RPE - For analysis - ACS — RPE

Description Green-azure crystalline powder
Identification ..... Positive
Chloride ..... ≤ 30 ppm
Water-insoluble matter ..... ≤ 100 ppm
Sulphate ..... ≤ 100 ppm
Ca ..... ≤50 ppm
K ..... ≤100 ppm
Ni ..... ≤100 ppm
Na ..... ≤500 ppm
Assay (oxidimetric) ..... 98.0÷102.0 %

Code	Size	Packaging	Notes
475405	250 g	Plastic bottle	
475407	1 kg	Plastic bottle	

Copper (II) acetate hydrate > RE - Pure

RE

Description .. Dark green crystalline powder
Identification ..... Positive
Chloride ..... ≤ 200 ppm
pH (5%m/m solution) ..... 5.0 - 6.0
Sulphate ..... ≤ 0.1 %
Fe ..... ≤ 250 ppm
Assay (oxydimetric) ..... ≥ 98.0 %

Code	Size	Packaging	Notes
364007	1 kg	Plastic bottle	
364008	25 kg	Plastic bucket	

Copper (II) carbonate (basic)

• Rame carbonato basico ico • Cuivre (II) carbonate basique • Cobre (II) carbonato básico
• Kupfer(II)-hydroxidcarbonat

Synonym:
Cupric carbonate basic

$CuCO_3.Cu(OH)_2$
Molecular Weight: 221,12
CAS: 12069-69-1
EEC-N: 235-113-6

Classification transport
ONU: 3077
Transport Hazard class: 9
Packing group III



Warning
H302-H332-H319-H410
P261-P264-P271-P304+P340-P305+P351+P338-P337+P313

Copper (II) carbonate (basic) > RPE - For analysis

RPE

Description ..... Green azure powder
Identification ..... Positive
Fe ..... ≤ 300 ppm
Pb ..... ≤ 100 ppm
Assay (oxidimetric) ..... 53 ÷ 57 % (Cu)
Cd ..... ≤ 10 ppm
Chlorine ..... ≤ 250 ppm
Specific gravity ..... 0.8 ÷ 1.05 g/cm3
Acidity solubility ..... Complete

Code	Size	Packaging	Notes
475555	250 g	Plastic bottle	
475557	1 kg	Plastic bottle	
475553	25 kg	Plastic bucket	

Copper (II) chloride dihydrate

• Rame cloruro ico diidrato • Cuivre (II) chlorure dihydraté • Cobre (II) cloruro dihidrato
• Kupfer(II)-chlorid-Dihydrat

Synonym:
Cupric chloride dihydrate

$CuCl_2.2H_2O$
Molecular Weight: 170,47
CAS: 10125-13-0
EEC-N: 215-704-5

Classification transport
ONU: 2802
Transport Hazard class: 8
Packing group III



Danger
H290-H302-H312-H318-H400-H411
P264-P280-P305+P351+P338-P310a-P330-P362+P364

Copper (II) chloride dihydrate > RPE - For analysis

RPE

Description ..... Green - azure crystals
Identification ..... Positive
pH sol. 5% at 20°C ..... 3.0 ÷ 3.8
Total nitrogen ..... ≤ 40 ppm
Assay (complexometric) ..... ≥ 99 %
Ca ..... ≤ 20 ppm
K ..... ≤ 20 ppm
Mg ..... ≤ 20 ppm
Na ..... ≤ 20 ppm
Sulphate ..... ≤ 50 ppm
As ..... ≤ 1 ppm
Fe ..... ≤ 10 ppm
Ni ..... ≤ 50 ppm
Pb ..... ≤ 40 ppm

Code	Size	Packaging	Notes
475685	250 g	Plastic bottle	
475687	1 kg	Plastic bottle	

Copper (II) chloride dihydrate > RE - Pure

RE

Description ..... Blue crystals
Identification ..... Positive
pH sol. 5% at 20°C ..... 3 ÷ 3.8
Sulphate ..... ≤ 50 ppm
Nitrogen compounds (N) ..... ≤ 0.004 %
As ..... ≤ 1 ppm
Ca ..... ≤ 20 ppm
Nitrate ..... ≤ 0.5 %
K ..... ≤ 20 ppm
Zn ..... ≤ 0.05 %
Solubility ..... Conform
Pb ..... ≤ 40 ppm
Assay ..... ≥ 99 %

Code	Size	Packaging	Notes
364507	1 kg	Plastic bottle	
364508	5 kg	Plastic tank	

Product specifications are subject to changes. Please visit our website for updates.

Copper (II) nitrate trihydrate

• Rame nitrato ico triidrato • Cuivre (II) nitrate trihydraté • Cobre (II) nitrato trihidrato
• Kupfer(II) nitrat-trihydrat

Synonym:
Cupric nitrate trihydrate

$Cu(NO_3)_2.3H_2O$
Molecular Weight: 241,6
CAS: 10031-43-3
EEC-N: 221-838-5

Classification transport
ONU: 1477
Transport Hazard class: 5.1
Packing group II

Danger
H272-H302
P210-P220-P264-P280-P301+P312a-P501a

Copper (II) nitrate trihydrate > RPE - For analysis

RPE

Description	Blue crystals	Chloride	≤10 ppm	Sulphate	≤25 ppm	Ni	≤10 ppm
Identification	Positive	Water-insoluble matter	≤30 ppm	As	≤1 ppm	Pb	≤10 ppm
pH sol. 5% at 25° C	3.0 ÷ 4.0	Substances not ppt. H2S	≤500 ppm	Ba	≤50 ppm	Zn	≤10 ppm
Ammonium	≤10 ppm	Subst. ppt. by (NH4)2S	≤50 ppm	Fe	≤20 ppm	Assay (oxidimetric)	≥99.5 %

Code	Size	Packaging	Notes
475782	100 g	Plastic bottle	
475786	500 g	Plastic bottle	
475783	1 kg	Plastic bottle	
475784	2.5 kg	Plastic bottle	

Copper (II) oxide

• Rame ossido ico • Cuivre (II) oxyde • Cobre (II) óxido • Kupfer(II) oxid

Synonym:
Cupric oxide

CuO
Molecular Weight: 79,55
CAS: 1317-38-0
EEC-N: 215-269-1

Warning
H302-H410
P264-P270-P301+P312a-P330-P391-P501a

Copper (II) oxide > RPE - For analysis

RPE

Description	Black powder	HNO3-insoluble matter	≤ 0.02 %	Zolfo totale (SO4)	≤ 0.01 %
Identification	Positive	Substances not ppt. H2S	≤ 0.2 %	C	≤ 500 ppm
Free alkalis	Conform	Nitrogen compounds (N)	≤ 0.002 %	Assay (iodometric)	≥ 99 %

Code	Size	Packaging	Notes
475994	100 g	Glass bottle	
475997	1 kg	Glass bottle	

Copper (II) oxide wire

• Rame ossido filo • Cuivre (II) oxyde en fils • Cobre (II) óxido hilos • Kupfer (II) oxiddraht

Synonym:
Cupric oxide

CuO
Molecular Weight: 79,55
CAS: 1317-38-0
EEC-N: 215-269-1

Warning
H302-H410
P264-P270-P301+P312a-P330-P391-P501a

Copper (II) oxide wire > RS - For microanalysis - ACS

RS

Description	Grey wire	Nitrogen compounds (N)	≤20 ppm	Sulphur compounds	≤120 ppm
Identification	Positive	Carbon compounds	≤20 ppm		

Code	Size	Packaging	Notes
475966	500 g	Plastic bottle	

Copper (II) sulfate anhydrous

• Rame solfato ico anidro • Cuivre (II) sulfate anhydre • Cobre (II) sulfato anhidro • Kupfer(II) sulfat wasserfrei

Synonym: *Cupric sulfate*

$CuSO_4$
Molecular Weight: 159,6
CAS: 7758-98-7
EEC-N: 231-847-6

Classification transport
ONU: 3077
Transport Hazard class: 9
Packing group III



Warning
H302-H315-H319-H410
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Copper (II) sulfate anhydrous > RPE - For analysis

RPE

Description .....Greyish powder
Identification ..... Positive
Loss on drying ..... ≤1 %
Chloride.....≤100 ppm
Dil. H2SO4-ins. matter.....≤50 ppm
Ca.....≤50 ppm
Fe .....≤30 ppm
K.....≤100 ppm
Mg .....≤100 ppm
Na.....≤200 ppm
Ni.....≤50 ppm
Pb .....≤80 ppm
Assay (oxidimetric) .....99 ÷ 100.5 % (s.s.)

Code	Size	Packaging	Notes
476245	250 g	Plastic bottle	
476247	1 kg	Plastic bottle	
476243	25 kg	Plastic bucket	

Copper (II) sulfate anhydrous > RE - Pure

RE

Description ..... Grey powder
Identification ..... Positive
Fe .....≤0.1 %
Water-insoluble matter .....≤0.1 %
Assay (oxidimetric) .....≥97 %

Code	Size	Packaging	Notes
365006	500 g	Plastic bottle	
365007	1 kg	Plastic bottle	
365002	25 kg	Plastic bucket	

Copper (II) sulfate pentahydrate

• Rame solfato ico pentaidrato • Cuivre (II) sulfate pentahydraté • Cobre (II) sulfato pentahidrato • Kupfer(II)-sulfat-Pentahydrat

Synonym: *Blue Vitriol*

$CuSO_4.5H_2O$
Molecular Weight: 249,68
CAS: 7758-99-8
EEC-N: 231-847-6



Danger
H302-H318-H410
P264-P280i-P301+P312a-P305+P351+P338-P310a-P501a

Copper (II) sulfate pentahydrate > RS - For microanalysis

RS

Description ..... Blue crystals
Identification ..... Positive
Not soluble matter .....≤ 0.005 %
Nitrogen compounds (N) .....≤ 0.002 %
Assay (iodometric).....98.0 ÷ 102.0 %

Code	Size	Packaging	Notes
476154	100 g	Glass bottle	

Copper (II) sulfate pentahydrate > RPE - For analysis - ACS

RPE

Description .....Blue crystals
Identification ..... Positive
Total nitrogen.....≤20 ppm
Chloride.....≤10 ppm
H2SO4-insoluble matter .....≤50 ppm
Ca.....≤50 ppm
Fe .....≤30 ppm
K.....≤100 ppm
Na.....≤200 ppm
Ni.....≤50 ppm
Assay (oxidimetric) .....98.0 ÷ 102.0 %

Code	Size	Packaging	Notes
476096	100 g	Plastic bottle	
476097	1 kg	Plastic bottle	
476099	5 kg	Plastic jar	
476092	25 kg	Plastic bucket	

Product specifications are subject to changes. Please visit our website for updates.

Copper (II) sulfate pentahydrate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-BP

ERBApharm

Description ..... Blue crystalline powder
Identification ..... Positive
Appearance of solution ..... Conform Ph.Eur.
Loss on drying ..... 35.0 ÷ 36.5 %
Chloride ..... ≤100 ppm
Ca ..... ≤50 ppm
Fe ..... ≤30 ppm
K ..... ≤100 ppm
Na ..... ≤200 ppm
Ni ..... ≤50 ppm
Pb ..... ≤50.0 ppm
Assay (oxidimetric) ..... 99.0 ÷ 100.5 %

Code	Size	Packaging	Notes
364757	1 kg	Plastic bottle	
364759	5 kg	Plastic tank	
364752	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in Xcipharm™ grade

Copper (II) sulfate solution 12.5%

• Rame solfato ico soluzione 12.5% • Cuivre (II) sulfate solution 12.5% • Cobre (II) sulfato solución 12.5% • Kupfersulfat 12.5 %

Synonym: *Cupric sulfate*

$CuSO_4$
Molecular Weight: 159,6
CAS: 7758-98-7



Danger
H318-H410
P273-P280i-P305+P351+P338-P310a-P391-P501a

Copper (II) sulfate solution 12.5% > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611022501	100 ml	Plastic bottle	Ref Ph.Eur 1022500

Copper (II) sulfate solution 12.5% > RPE - For analysis

RPE

Description ..... Clear blue liquid
Identification ..... Positive
Density at 20° C ..... ~ 1.08
Assay (oxidimetric) ..... 12.0 ÷ 13.0 %

Code	Size	Packaging	Notes
E476211	1 l	Bottle	

Copper (II) sulphate solution

• Rame (II) solfato in soluzione • Cuivre (II) sulfate en solution • Cobre (II) sulfato de en solución • Kupfer(II)-sulfat in Lösung

$CuSO_4$
Molecular Weight: 159,61
CAS: 7758-98-7

Copper (II) sulphate solution > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
616001038	100 ml	Plastic bottle	Cupric sulfate CS
616001037	500 ml	Plastic bottle	Cupric sulfate CS

CPME ▶ Cyclopentyl methyl ether

m-Cresol

• m-Cresolo • m-Crésol • m-Cresol • m-Kresol

Synonym:
• *3-Hydroxytoluene*
• *3-Methylphenol*

$CH_3C_6H_4OH$
Molecular Weight: 108
CAS: 108-39-4
EEC-N: 203-577-9

Classification transport
ONU: 2076
Transport Hazard class: 6.1
Packing group II



Danger
H301-H311-H314
P280-P301+P310a-P301+P330+P331-
P303+P361+P353-P304+P340-P305+P351+P338

m-Cresol > RE - Pure

RE

Refractive index at 20°C ..... 1.539 - 1.543
Description ..... Reddish liquid
Water content (K.F.) ..... ≤ 500 mg/Kg
Identification ..... Positive
Colour ..... ≤ 120 Hazen
Water (K.F.) ..... ≤0.05 %
m-Cresol ..... ≥ 99 %
Assay (GLC) ..... ≥99 %
o-cresol ..... ≤ 0.3 %
Boiling point ..... ~ 202 ° C
Density d20/4 ..... None
Melting point ..... ~ 12 ° C

Code	Size	Packaging	Notes
440435	250 g	Glass bottle	
440437	1 kg	Glass bottle	

m-Cresol > RE - Pure - For synthesis

RE

Refractive index at 20°C ..... 1.539 - 1.543
Water content (K.F.) ..... ≤ 500 mg/Kg
Colour ..... ≤ 120 Hazen
m-Cresol ..... ≥ 99 %
o-cresol ..... ≤ 0.3 %
Density d20/4 ..... None

Code	Size	Packaging	Notes
PA160021	2.5 l	Glass bottle	

o-Cresol

• o-Cresolo • o-Crésol • o-Cresol • o-Kresol

Synonym:
2-Methylphenol

$CH_3C_6H_4OH$
Molecular Weight: 108,14
CAS: 95-48-7
EEC-N: 202-423-8

Classification transport
ONU: 3455
Transport Hazard class: 6.1
Packing group II



Danger
H301-H311-H314
P280-P301+P310a-P301+P330+P331-
P303+P361+P353-P304+P340-P305+P351+P338

o-Cresol > RE - Pure

RE

Description ..... Colourless to yellow brown liquid or solid
Identification ..... Positive
Boiling point ..... 189.5 ÷ 192.5 °C
Melting point ..... 29 ÷ 31 ° C
Water (K.F.) ..... ≤ 0.3 %
Assay (GLC) ..... ≥ 98 %

Code	Size	Packaging	Notes
440385	250 g	Glass bottle	
440387	1 kg	Glass bottle	
440381	2.5 Kg	Plastic bottle	

m-Cresol purple

• Porpora m-cresolo • Pourpre de m-cresol • Púrpura de m-cresol • m-Kresolpurpur

Synonym:
m-Cresolsulfonphthalein

$C_{21}H_{18}O_5S$
Molecular Weight: 382,44
CAS: 2303-01-7
EEC-N: 218-960-6

m-Cresol purple > RPE - For analysis

RPE

Description ..... Dark green powder
Identification ..... Positive
Solubilità (0.1% in ETOH 50%) ..... Complete
E (1%/1cm) a 578nm in tamp. ..... ≥ 900
Colour change ..... yellow - violet
pH range ..... 7.4 ÷ 9.0

Code	Size	Packaging	Notes
470067	1 g	Glass bottle	
470068	25 g	Glass bottle	

Clark indicator series. Acid-base indicator

Product specifications are subject to changes. Please visit our website for updates.

o-Cresol red

• Rosso o-cresolo • Rouge de o-crésol • Rojo de o-cresol • o-Kresolrot

Synonym:
o-Cresolsulfonphthalein

$C_{21}H_{18}O_5S$
Molecular Weight: 382,44
CAS: 1733-12-6
EEC-N: 217-064-2

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

o-Cresol red > RPE - For analysis

RPE

Description ..... Red brown powder
Identification ..... Positive
Sensitivity(pH 7.2-8.8) ..... Conform
Loss on drying (110°C) ..... ≤ 3.0 %
Colour change ..... yellow red

Code	Size	Packaging	Notes
476778	5 g	Glass bottle	

Clark indicator series

o-Cresol Red solution 0.2% in ethanol

• Rosso o-cresolo soluzione 0.2% in alcole etilico • Rouge de o-crésol solution 0.2% dans l'éthanol
• Rojo de o-cresol solución 0.2% en alcohol etilico • O-Kresol-rote Lösung 0.2% in Ethanol

Synonym:
o-Cresolsulfonphthalein

$C_{21}H_{18}O_5S$
Molecular Weight: 382,44
CAS: 1733-12-6

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group III

Warning
H226-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

o-Cresol Red solution 0.2% in ethanol > RPE - For analysis

RPE

Description ..... Red liquid
Identification ..... Positive
Sensitivity(pH 7.2-8.8) ..... Conform
Colour change ..... Yellow-red violet

Code	Size	Packaging	Notes
E476805	250 ml	Bottle	

Indicator series Clark. Indicator acid-base

Cresol red R solution in Sodium hydroxide 0.1N/ Ethanol 96% / Water

• Soluzione Rosso cresolo/Sodio idrossido 0.1N / Etanolo 96%/ Acqua • Rouge de crésol solution dans l'hydroxyde de sodium 0.1N / éthanol 96% / eau
• Rojo de o-cresol solución en sodio hidróxido 0.1N / alcohol etilico 96°/agua • Kresol-rote Lösung Natronlauge 0.1N/ ethanol 96% / wasser

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group III

Warning
H226
P210-P241-P280-P303+P361+P353-P403+P235-P501a

Cresol red R solution in Sodium hydroxide 0.1N/ Ethanol 96% / Water > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611022801	100 ml	Bottle	Ref Ph.Eur 1022801

o-Cresolsulfonphthalein ▶ o-Cresol red

Crystal violet

• Violetto cristalli • Violet cristallisé • Violeta cristal • Kristallviolett

Synonym:
• *Crystal violet solution*
• *Basic violet 3*

$C_{25}H_{30}ClN_3$
Molecular Weight: 407,99
CAS: 548-62-9
EEC-N: 208-953-6



Danger
H302-H318-H351-H410
P264-P280-P301+P312a-P305+P351+P338-P310a-P308+P313

Crystal violet > RPE - For analysis - C.I. 42555

RPE

Description ..... Dark green powder
Identification ..... Positive
Suitability for anhydrous titration ...Conform
Loss on drying ..... ≤9 %
Alcohol-insolub. matter..... ≤0.5 %
Residue on ignition..... ≤2.5 %
Colour change..... yellow blue
pH range ..... 0.1 - 2.0

Code	Size	Packaging	Notes
491502	25 g	Glass bottle	

Dye for microscopy (bacteriology-Botanical-histology). Indicator acid - base

Crystal violet oxalate for Gram-Hucker Kit

• Violetto cristalli ossalato soluzione per kit Gram-Hucker • Violet cristallisé oxalate pour kit de Gram-Hucker
• Violeta cristal oxalato solución para kit Gram-Hucker • Kristallviolettoxalat für Gram-Hucker Kit



Warning
H319-H351-H373-H412
P260-P264-P280-P305+P351+P338-P308+P313-P337+P313

Crystal violet oxalate for Gram-Hucker Kit > RS - For bacteriology

RS

Description ..... Violet clear liquid
Identification ..... Positive

Code	Size	Packaging	Notes
491561	250 ml	Bottle	In Vitro Diagnostic Medical Device

Crystal violet solution 0.5% in anhydrous acetic acid

• Violetto cristalli soluzione 0.5% in acido acetico anidro
• Violet cristallisé solution 0.5% dans l'acide acétique anhydre
• Violeta cristal solución 0.5% en acido acético anhidro • Kristallviolettlösung 0.5% in wasserfreier Essigsäure

Synonym:
• *Crystal violet solution*
• *Basic violet 3*

$C_{25}H_{30}ClN_3$
Molecular Weight: 407,99
CAS: 548-62-9

Classification transport
ONU: 2920
Transport Hazard class: 8
Packing group II



Danger
H226-H314-H412
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Crystal violet solution 0.5% in anhydrous acetic acid > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611022901	100 ml	Glass bottle	Ref Ph.Eur 1022901

Crystal violet solution 0.5% in anhydrous acetic acid > RPE - For analysis

RPE

Description ..... Violet clear liquid
Identification ..... Positive

Code	Size	Packaging	Notes
E491551	500 ml	Bottle	

Cupri-citric solution

• Cupri-citrico soluzione • Solution cupri-citrique • Cobre-cítrico solución • Cupri-Zitronensäure-Lösung



H411
P273-P391-P501a

Cupri-citric solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611023100	1 l	Plastic bottle	Ref Ph.Eur 1023100

Product specifications are subject to changes. Please visit our website for updates.

Cupriethylenediamine solution

• Cuprietilendiammina soluzione • Cuprié thylènediamine solution • Cobre etilendiamina solución
• Cupriethylendiamin-Lösung

Synonym:
• *Bis(ethylenediamine)copper(II) hydroxide solution*
• *Copper(II)-ethylenediame complex*

$Cu(NH_2CH_2CH_2NH_2)_2(OH)_2$
Molecular Weight: 217,76
CAS: 14552-35-3

Classification transport
ONU: 1761
Transport Hazard class: 8
Packing group II

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Cupriethylenediamine solution > RPE - For analysis

RPE

Description ..... Blue liquid
Identification ..... Positive
Assay (ex Cu) ..... 0.98 ÷ 1.02 M/l
Assay(Ethylendiamine) ..... 1.96 ÷ 2.04 M/l
Assay (cupriethylenediamine) ..... 0.98 ÷ 1.02 M/l

Code	Size	Packaging	Notes
E441071	1 l	Glass bottle	

Cupriethylenediamine hydroxide 1 mol/l

• Cuprietilendiammina idrossido 1 mol/l • Cuivre (II) éthylènediamine hydroxyde 1 mol/l • Cobre (II) dietilendiaminico hidroxido 1 mol/l
• Kupfer (II) ethylendiaminhydroxid 1mol/l

Classification transport
ONU: 1761
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Cupriethylenediamine hydroxide 1 mol/l > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613008700	1 l	Glass bottle	Ref Ph.Eur 3008700

Cupri-tartaric solution

• Cupri-tartarico soluzione • Solution cupri-tartrique • Cobre-tartárico solución • Cupri-Weinsäure-Lösung

Classification transport
ONU: 1719
Transport Hazard class: 8
Packing group II

Danger
H314-H411
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Cupri-tartaric solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611023300	2 x 500 ml	Plastic bottle	Ref Ph.Eur 1023300

Cyanide standard solution

• Cianuri standard soluzione • Cyanure solution standard • Cianuro, solución patrón • Cyanid-Standardlösung

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group III

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Cyanide standard solution > RS - Standard solution for ion chromatography

RS

Code	Size	Packaging	Notes
503358	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water and nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Cyclohexane

• Cicloesano • Cyclohexane • Ciclohexano • Cyclohexan

$CH_2(CH_2)_4CH_2$
Molecular Weight: 84,16
CAS: 110-82-7
EEC-N: 203-806-2

Classification transport
ONU: 1145
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H336-H304-H410
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Cyclohexane > RS - For HPLC - Isocratic Grade

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.776 ÷ 0.782
Refractive index at 20°C. 1.4229 ÷ 1.4299
Boiling point..... 80.2 ÷ 81.2 ° C
Water (K.F.) ..... ≤100 ppm
Residue on evaporation ..... ≤2 ppm
Acidity or alkalinity..... ≤0.0002 meq/g
Assay (GLC) ..... ≥99.9 %
U.V. Transmittance
at 200 nm ..... ≥4 %
at 210 nm ..... ≥15 %
at 220 nm ..... ≥50 %
at 230 nm ..... ≥80 %
at 240 nm ..... ≥92 %
At 250 nm ..... ≥ 98 %
at 260 nm ..... ≥98.5 %
Aromatic compounds..... ≤ 5 ppm

Code	Size	Packaging	Notes
412431000	1 l	Glass bottle	
412432000	2.5 l	Glass bottle	

Cyclohexane > RS - PESTIPUR - For pesticide analysis

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Colour ..... ≤ 10 hazen
Water ..... ≤ 100 ppm
Not volatile residue..... ≤ 2 ppm
GC-ECD (Lindano) ..... ≤ 3 ng/l
GC-NPD (Ethylparation)..... ≤ 3 ng/l
Assay (GLC) ..... ≥ 99.8 %
Refractive index at 20°C..... 1.424 ÷ 1.428

Code	Size	Packaging	Notes
436931	1 l	Glass bottle	
436932	2.5 l	Glass bottle	

Cyclohexane > RS - SPECTROSOL - For optical spectroscopy - Reag.Ph.Eur.

RS

Description ..... Clear colourless liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Density at 20° C ..... 0.776 ÷ 0.782
Refractive index at 20°C. 1.4229 ÷ 1.4299
Boiling point..... 80.2 ÷ 81.2 ° C
Melting point..... 5.5 ÷ 7.5 ° C
Acidity or alkalinity..... ≤0.0002 meq/g
Water (K.F.) ..... ≤100 ppm
Residue on evaporation ..... ≤5 ppm
Aromatic compounds..... ≤5 ppm
Assay (GLC) ..... ≥99.9 %
U.V. Transmittance
at 210 nm ..... ≥15 %
at 220 nm ..... ≥45 %
at 230 nm ..... ≥75 %
at 250 nm ..... ≥98 %
UV Absorbance at 220 nm ..... ≤ 0.35 AU
UV Absorbance at 235 nm ..... ≤ 0.16 AU
UV Absorbance at 240 nm ..... ≤ 0.05 AU
UV Absorbance at 250 nm ..... ≤ 0.01 AU

Code	Size	Packaging	Notes
436967	1 l	Glass bottle	
436963	2.5 l	Glass bottle	

Cyclohexane > RS - Anhydrous - For analysis

RS

Refractive index at 20°C..... 1.424 - 1.428
Water content (K.F.) ..... ≤ 50 mg/Kg
Non volatile residue..... ≤ 10 mg/Kg
Colour ..... ≤ 10 Hazen
Aromatic compounds..... ≤ 150 mg/Kg
Assay (GC) ..... ≥ 99.8 %
Methylcyclohexane ..... ≤ 1000 mg/Kg
Clear liquid appearance ..... Conform
Identification (IR)..... Conform
Density d20/4 ..... 0.775 - 0.782
Total sulphur (S) ..... ≤ 1 ppm

Code	Size	Packaging	Notes
P0251010	200 ml	Bottle with septum	
P0251016	1 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Cyclohexane > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Clear liquid	Alcohol miscibility ..... Complete	Melting point ..... 5.5 ÷ 7.5 ° C	Total sulphur ..... ≤1 ppm
Colour (APHA) ..... ≤10	Diethyl ether miscib. ..... Complete	Water (K.F.) ..... ≤100 ppm	Assay (GLC) ..... ≥99.8 %
Identification (I.R.) ..... Conform	Density at 20° C ..... 0.776 ÷ 0.782	Residue on evaporation ..... ≤10 ppm	Aromatic compounds ..... ≤ 150 ppm
Water solubility ..... Conform	Refractive index at 20°C 1.4229 ÷ 1.4299	Acidity (acetic acid) ..... ≤10 ppm	Methylcyclohexane ..... ≤ 0.1 %
Ready carbonizable substances ..... Conform	Boiling point ..... 80 ÷ 82 ° C	Subst. reducing KMnO4 ..... ≤20 ppm(5m)	

Code	Size	Packaging	Notes
436903	1 l	Glass bottle	
436905	2.5 l	Glass bottle	
436906	5 l	Plastic tank	
436901	10 l	Metal tank	
436902	21 kg	Metal drum	
436908	150 kg	Metal drum	

Cyclohexane > RE - Pure

RE

Description ..... Clear liquid	Density at 20° C ..... 0.774 ÷ 0.784	Boiling point ..... 80 ÷ 82 ° C	Residue on evaporation ..... ≤50 ppm
Identification ..... Positive	Refractive index at 20°C 1.4214 ÷ 1.4314	Water (K.F.) ..... ≤150 ppm	Assay (GLC) ..... ≥ 99.8 %

Code	Size	Packaging	Notes
333752	1 l	Glass bottle	
333751	2.5 l	Glass bottle	
528215	5 l	Plastic tank	
508235	10 l	Plastic tank	
333753	21 kg	Metal drum	
528216	25 l	Metal drum	
528217	200 l	Metal drum	

Cyclohexane-d12

• Cicloesano-d12 • Cyclohexane-d12 • Ciclohexano-d12 • Cyclohexan-d12

C_6D_{12}
Molecular Weight: 96,07
CAS: 1735-17-7
EEC-N: 217-077-3



Danger
H225-H315-H336-H304-H410
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Cyclohexane-d12 > RS - For NMR - min 99.5%

RS

Code	Size	Packaging	Notes
P5151A	2 x 0.5 ml	Glass ampoule	

For specifications, contact our customer service for a certificate of analysis

Cyclohexanol

• Cicloesanolo • Cyclohexanol • Ciclohexanol • Cyclohexanol

$CH_2(CH_2)_4CHOH$
Molecular Weight: 100
CAS: 108-93-0
EEC-N: 203-630-6



Warning
H302-H332-H315-H335
P261-P271-P304+P340-P332+P313-P362+P364-P403+P233

Cyclohexanol > RE - Pure

RE

Description ..... Clear colourless liquid	Identification ..... Positive	Residue on evaporation ..... ≤0.1 %	Assay (GLC) ..... ≥98 %

Code	Size	Packaging	Notes
333801	1 l	Glass bottle	

Cyclohexanol > RE - Pure - For synthesis

RE

Refractive index at 25°C..1.4620 - 1.4660 | Colour ..... ≤ 15 Hazen | Assay (GC) ..... ≥ 98.5 % | Cyclohexanone.....≤ 6000 mg/Kg

Code	Size	Packaging	Notes
P0260248	25 l	Metal drum	
P0260268	200 l	Metal drum	

Cyclohexanone

• Cicloesanone • Cyclohexanone • Ciclohexanona • Cyclohexanon

$CH_2(CH_2)_4CO$
Molecular Weight: 98,15
CAS: 108-94-1
EEC-N: 203-631-1

Classification transport
ONU: 1915
Transport Hazard class: 3
Packing group III



Warning
H226-H332
P210-P241-P261-P280-P303+P361+P353-P304+P340

Cyclohexanone > RPE - For analysis

RPE

Description ..... Clear colourless liquid | Identification ..... Positive | Density at 20° C .....0.941 ÷ 0.951 | Refractive index at 20°C.1.4477 ÷ 1.4537 | Boiling point.....152 ÷ 157 °C | Water (K.F.) ..... ≤0.1 % | Residue on evaporation .....≤500 ppm | Aldehydes(Formaldehyde) ..... ≤ 0.1 % | Cyclohexanol..... ≤ 0.1 % | Heavy metals (Pb).....≤2 ppm | Fe .....≤10 ppm | Assay (GLC) ..... ≥99.5 %

Code	Size	Packaging	Notes
437053	1 l	Glass bottle	
437052	2.5 l	Glass bottle	
437055	25 kg	Metal drum	

Cyclohexanone > RE - Pure

RE

Description ..... Yellow clear liquid | Identification ..... Positive | Density at 20° C .....0.941 ÷ 0.951 | Refractive index at 20°C.1.4457 ÷ 1.4557 | Boiling point.....152 ÷ 157 °C | Water (K.F.) ..... ≤0.1 % | Acidity(CycloexilcarAc) ..... ≤0.2 % | Cyclohexanol..... ≤ 0.1 % | Residue on evaporation .....≤500 ppm | Assay (GLC) ..... ≥99 %

Code	Size	Packaging	Notes
333901	1 l	Glass bottle	
333905	5 l	Aluminium can	
333902	26 kg	Metal drum	
528332	200 l	Metal drum	

Cyclohexylamine

• Cicloesilammina • Cyclohexylamine • Ciclohexilamina • Cyclohexylamin

Synonym:
Aminocyclohexane

$CH_2(CH_2)_4CHNH_2$
Molecular Weight: 99,18
CAS: 108-91-8
EEC-N: 203-629-0

Classification transport
ONU: 2357
Transport Hazard class: 8
Packing group II



Danger
H226-H302-H312-H314-H361f
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P362+P364

Cyclohexylamine > RPE - For analysis

RPE

Description ..... Clear colourless liquid | Identification ..... Positive | Density at 20° C .....0.864 ÷ 0.870 | Refractive index at 20°C.1.4572 ÷ 1.4612 | Boiling point.....134.0 ÷ 135.0 ° C | Water (K.F.) ..... ≤ 0.3 % | Assay (alkalimetric)..... ≥99 %

Code	Size	Packaging	Notes
437104	1 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Cyclopentyl methyl ether

• Ciclopentil-metil-etere • Cyclopentylméthyléther • Ciclopentil metil eter • Cyclopentylmethyl Ether

Synonym:
CPME

$C_6H_{12}O$
Molecular Weight: 100,16
CAS: 5614-37-9

Classification transport
ONU: 3271
Transport Hazard class: 3
Packing group II

Danger
H225-H302-H315-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P332+P313

Cyclopentyl methyl ether > RE - Pure

RE

Refractive index at 20°C..... 1.4199 - 1.4219
Water content (K.F.).....≤ 100 mg/Kg
Colour .....≤ 10 Hazen
Assay (GC).....≥ 99.9 %
Peroxides.....≤ 50 meq/Kg

Code	Size	Packaging	Notes
P8010216	1 l	Glass bottle	
P8010229	5 l	Plastic tank	
P8010248	25 l	Metal drum	
P8010268	200 l	Metal drum	

L-Cysteine

• L-Cisteina • L-Cystéine • L-Cisteína • L-Cystein

Synonym:
(R)-2-Amino-3-mercaptopropionic acid

$HSCH_2CH(NH_2)COOH$
Molecular Weight: 121,16
CAS: 52-90-4
EEC-N: 200-158-2

Warning
H302
P264-P270-P301+P312a-P330-P501a

L-Cysteine > RPE - For analysis

RPE

Description ..... White powder
Identification ..... Positive
Specific optical rotation.....+8 ÷ +9.5 °
Loss on drying .....≤ 0.5 %
Chloride.....≤ 500 ppm
Heavy metals (Pb).....≤ 10 ppm
Residue on ignition.....≤ 0.1 %
Sulphate.....≤ 300 ppm
Assay with HClO4..... ≥ 98.0 % (s.s)

Code	Size	Packaging	Notes
437308	5 g	Glass bottle	
437309	100 g	Glass bottle	

L-Cystine

• L-Cistina • L-Cystine • L-Cistina • L-Cystin

Synonym:
(R,R)-3,3'-Dithiobis(2-aminopropionic acid)

$C_6H_{12}N_2O_4S_2$
Molecular Weight: 240,3
CAS: 56-89-3
EEC-N: 200-296-3

Danger
H301-H315-H319-H335
P261-P271-P301+P310a-P304+P340-
P305+P351+P338-P403+P233

L-Cystine > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Loss on drying ..... ≤ 0.5 %
Potere rotat. spec. (C=2; HCl 1M)... -209 ÷ -224 ° (s.s.)
Ammonium.....≤ 200 ppm
Chloride.....≤ 200 ppm
Sulphate.....≤ 300 ppm
Heavy metals (Pb).....≤ 10 ppm
Residue on ignition..... ≤ 0.1 %
Fe .....≤ 10 ppm
Assay (acidimetric) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
437351	10 g	Glass bottle	
437355	100 g	Glass bottle	

Decalcifying agent
• Decalcificante • Décalcifiant • Descalcificante • Entkalkung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III



Warning
H290
P234-P390-P406

Decalcifying agent > RS - For histology

RS

Description ..... Clear colourless liquid | Identification ..... Positive | Titolo EDTA Na2 (mom. preparaz.) ≥ 0.26 % (p/p) | Titolo HCl (mom. preparaz.) ..≥ 2.21 % (p/p)

Code	Size	Packaging	Notes
441221	1 l	Plastic bottle	In Vitro Diagnostic Medical Device

Decane
• Decano • Décane • Decano • Decan

$CH_3(CH_2)_8CH_3$
Molecular Weight: 142,28
CAS: 124-18-5
EEC-N: 204-686-4

Classification transport
ONU: 2247
Transport Hazard class: 3
Packing group III



Danger
H226-H304
P210-P241-P280-P301+P310a-P303+P361+P353-P403+P235

Decane > RE - Pure

RE

Refractive index at 20°C ..... 1.408 - 1.412 | Colour ..... ≤ 10 Hazen | Assay (GC) ..... ≥ 99 %

Code	Size	Packaging	Notes
P0093016	1 l	Glass bottle	

Decanedioic acid ▶ Sebacic acid

1-Decanesulfonic acid sodium salt
• Acido 1-decansolfonico sale sodico • Acide 1-decanesulfonique sel sodique
• Acido 1-decanosulfónico sal sódica • 1-Decansulfonsäure-Natriumsalz

Synonym:
• *Sodium 1-decanesulfonate*

$CH_3(CH_2)_9SO_3Na$
Molecular Weight: 244,33
CAS: 13419-61-9



Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

1-Decanesulfonic acid sodium salt > RS - For ion pair chromatography

RS

Description ..... White to cream powder
Water (K.F.) ..... ≤ 1.0 %
Assay ..... ≥ 99.0 %
Absorbance (0,2M)
At 210 nm ..... ≤ 0.05 AU
At 220 nm ..... ≤ 0.03 AU
At 230 nm ..... ≤ 0.02 AU
At 260 nm ..... ≤ 0.02 AU
Solubility 5% in water ..... Clear and colourless solution

Code	Size	Packaging	Notes
405871	25 g	Glass bottle	
405872	100 g	Plastic bottle	

Degree of coloration of liquids: primary solutions
• Soluzioni primarie per il grado di colorazione dei liquidi • Solutions primaires pour le degré de coloration des liquides
• Soluciónes primarias para el grado de coloración de líquidos • Primärlösungen für den Einfärbungsgrad von Flüssigkeiten

Degree of coloration of liquids: primary solutions > RS - For analysis according to Ph. Eur. Chap. 2.2.2

RS

Code	Size	Packaging	Notes
612202100	100 ml	Glass bottle	Yellow primary solution
612202200	100 ml	Plastic bottle	Red primary solution
612202300	100 ml	Plastic bottle	Blue primary solution

Product specifications are subject to changes. Please visit our website for updates.

Degree of coloration of liquids: primary solutions > RS - For analysis according to Ph.Ch.

RS

Code	Size	Packaging	Notes
618000001	100 ml	Plastic bottle	Red primary solution
618000002	100 ml	Plastic bottle	Yellow primary Solution
618000003	100 ml	Plastic bottle	Blue primary solution

Degree of coloration of liquids: primary solutions > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
612002100	100 ml	Glass bottle	Yellow primary solution
612002200	100 ml	Plastic bottle	Red primary solution
612002300	100 ml	Plastic bottle	Blue primary solution

Degree of coloration of liquids: standard solutions

• Soluzioni standard per il grado di colorazione dei liquidi • Solutions étalon pour le degré de coloration des liquides
• Soluciónes estándar para el grado de coloración de líquidos • Standardlösungen für den Einfärbungsgrad von Flüssigkeiten

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group III



Danger
H290-H314-H334-H317-H350i-H412-HA26
P280-P284-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P362+P364-P342+P311a

Degree of coloration of liquids: standard solutions > RS - For analysis according to Ph. Eur. Chap. 2.2.2

RS

Code	Size	Packaging	Notes
612202510	125 ml	Plastic bottle	Standard solution B (brown)
612202520	125 ml	Plastic bottle	Standard solution BY (brownish-yellow)
612202530	125 ml	Plastic bottle	Standard solution Y (yellow)
612202540	125 ml	Plastic bottle	Standard solution GY (greenish-yellow)
612202550	125 ml	Plastic bottle	Standard solution R (red)

Density Standard

• Densità standard • Etalons de densité • Patrones de densidad • Dichtestandards

Density Standard > RS - PREMIUM - For calibration

RS

Code	Size	Packaging	Notes
540401	100 ml	Glass bottle	0.6960g/ml at 15°C
540402	100 ml	Glass bottle	0.8715g/ml at 15°C
540403	100 ml	Glass bottle	1.0040g/ml at 15°C
540404	100 ml	Glass bottle	1.2498g/ml at 15°C
540405	100 ml	Glass bottle	0.6919g/ml at 20°C
540406	100 ml	Glass bottle	0.7033g/ml at 20°C
540407	100 ml	Glass bottle	0.7488g/ml at 20°C
540408	100 ml	Glass bottle	0.8668g/ml at 20°C
540409	100 ml	Glass bottle	1.0005g/ml at 20°C
540410	100 ml	Glass bottle	1.0301g/ml at 20°C
540411	100 ml	Glass bottle	1.0792g/ml at 20°C
540412	100 ml	Glass bottle	1.1651g/ml at 20°C
540413	100 ml	Glass bottle	1.2486g/ml at 20°C
540414	100 ml	Glass bottle	1.3304g/ml at 20°C
540415	100 ml	Glass bottle	1.5799g/ml at 20°C
540416	100 ml	Glass bottle	1.7470g/ml at 20°C
540417	100 ml	Glass bottle	1.9141g/ml at 20°C
540418	100 ml	Glass bottle	0.8207g/ml at 40°C
540420	100 ml	Glass bottle	0.9323g/ml at 40°C
540421	100 ml	Glass bottle	1.2408g/ml at 40°C
540422	100 ml	Glass bottle	0.9990g/ml at 60°C

In accordance with ASTM D1480-12 for testing of Density or Relative Density (specific and API gravity) by Bingham Pycnometer. For density measurement by pycnometric techniques, vibrational techniques or hydrometer based techniques

Product specifications are subject to changes. Please visit our website for updates.

Density Standard > RS - QUALITY - For calibration

RS

Code	Size	Packaging	Notes
540451	100 ml	Glass bottle	0.7524g/ml at 15°C
540452	100 ml	Glass bottle	0.7721g/ml at 15°C
540453	100 ml	Glass bottle	0.7933g/ml at 15°C
540454	100 ml	Glass bottle	0.8168g/ml at 15°C
540455	100 ml	Glass bottle	0.8428g/ml at 15°C
540456	100 ml	Glass bottle	0.8715g/ml at 15°C
540457	100 ml	Glass bottle	0.6919g/ml at 20°C
540458	100 ml	Glass bottle	0.7033g/ml at 20°C
540459	100 ml	Glass bottle	0.7488g/ml at 20°C
540460	100 ml	Glass bottle	0.7893g/ml at 20°C
540461	100 ml	Glass bottle	0.8126g/ml at 20°C
540462	100 ml	Glass bottle	0.8384g/ml at 20°C
540463	100 ml	Glass bottle	0.8668g/ml at 20°C
540464	100 ml	Glass bottle	0.9098g/ml at 20°C
540465	100 ml	Glass bottle	0.9476g/ml at 20°C
540566	100 ml	Glass bottle	1.0005g/ml at 20°C
540567	100 ml	Glass bottle	1.0301g/ml at 20°C
540568	100 ml	Glass bottle	0.8622g/ml at 25°C
540569	100 ml	Glass bottle	0.9438g/ml at 25°C
540470	100 ml	Glass bottle	0.9969g/ml at 25°C
540471	100 ml	Glass bottle	0.9245g/ml at 50°C
540472	100 ml	Glass bottle	0.9695g/ml at 60°C
540473	100 ml	Glass bottle	0.9815g/ml at 80°C

In accordance with ASTM D4052 for testing of Density, Relative Density and API Gravity of Liquids by Digital Density Meter
For density measurement by vibrational techniques or hydrometer based techniques

Deuterium oxide-d2

• Deuterio ossido-d2 • Deutérium oxyde-d2 • Deuterio oxido-d2 • Deuteriumoxid-d

Synonym:
Heavy water

D_2O
Molecular Weight: 20,03
CAS: 7789-20-0
EEC-N: 232-148-9

Deuterium oxide-d2 > RS - For NMR - min 99.98%

RS

Code	Size	Packaging	Notes
P5179	10 x 0.75 ml	Glass ampoule	
P5175	25 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Deuterium oxide-d2 > RS - For NMR - min 99.8%

RS

Code	Size	Packaging	Notes
P5169	10 x 0.75 ml	Glass ampoule	
P5164	5 x 10 ml	Glass bottle	
P5165	25 ml	Glass bottle	
P5165S	25 ml	Bottle with septum	
P5166	100 ml	Glass bottle	
P5168	1 l	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Deuterium oxide-d2 + 0.01% DMSO

• Deuterio ossido-d2 + 0.01% DMSO • Deutérium oxyde-d2 + 0.01% DMSO
• Deuterio oxido-d2 + 0.01% DMSO • Deuteriumoxid-d + 0.01% DMSO

Synonym:
Heavy water

D_2O
Molecular Weight: 20,03
CAS: 7789-20-0
EEC-N: 232-148-9

Deuterium oxide-d2 + 0.01% DMSO > RS - For NMR - min 99.98%

RS

Code	Size	Packaging	Notes
P5170D	10 x 0.6 ml	Glass ampoule	

For specifications, contact our customer service for a certificate of analysis

Deuterium oxide-d2 + 0.5% TSP d4

• Deuterio ossido-d2 + 0.5% TSP d4 • Deutérium oxyde-d2 + 0.5% TSP d4
• Deuterio oxido-d2 + 0.5% TSP d4 • Deuteriumoxid-d + 0.5% TSP d4

Synonym:
Heavy water

D_2O
Molecular Weight: 20,03
CAS: 7789-20-0
EEC-N: 232-148-9

Deuterium oxide-d2 + 0.5% TSP d4 > RS - For NMR - min 99,9%

RS

Code	Size	Packaging	Notes
P5161T	10 x 0.6 ml	Glass ampoule	

For specifications, contact our customer service for a certificate of analysis

Deuterium oxide-d2 + 0.03% TSP d4

• Deuterio ossido-d2 + 0.03% TSP d4 • Deutérium oxyde-d2 + 0.03% TSP d4
• Deuterio oxido-d2 + 0.03% TSP d4 • Deuteriumoxid-d + 0.03% TSP d4

Synonym:
Heavy water

D_2O
Molecular Weight: 20,03
CAS: 7789-20-0
EEC-N: 232-148-9

Deuterium oxide-d2 + 0.03% TSP d4 > RS - For NMR - min 99.8%

RS

Code	Size	Packaging	Notes
P5160T	10 x 0.6 ml	Glass ampoule	

For specifications, contact our customer service for a certificate of analysis

Devarda's alloy

• Lega Devarda • Alliage de Devarda • Aleación según Devarda • Devarda-Legierung

CAS: 8049-11-4

Classification transport
ONU: 3132
Transport Hazard class: 4.3
Packing group II



Danger
H228-H261-H411
P210-P223-P231a+P232-P241-P280-P402+P404

Devarda's alloy > RPE - For analysis

RPE

Description ..... Greyish metallic powder
Identification ..... Positive

Concentration in element
N..... < 10 ppm

Al..... ~ 45 %
Cu..... ~ 50 %
Zn..... ~ 5 %

Code	Size	Packaging	Notes
457625	250 g	Glass bottle	
457627	1 kg	Plastic bottle	

Dextrose ▶ D(+)Glucose anhydrous

Product specifications are subject to changes. Please visit our website for updates.

Diacetone alcohol

• Diacetonalcole • Diacétone alcool • Diacetonalcohol • Diacetonalkohol

Synonym:
4-Hydroxy-4-methyl-2-pentanone

$(CH_3)_2C(OH)CH_2COCH_3$
Molecular Weight: 116
CAS: 123-42-2
EEC-N: 204-626-7

Warning
H319-H361fd-H335
P261-P271-P280-P304+P340-P305+P351+P338-P403+P233

Diacetone alcohol > RPE - For analysis

RPE

Description ..... Yellow clear liquid
Identification ..... Positive
Water miscibility ..... Conform
Alcohol miscibility ..... Conform
Density at 25° C ..... 0.927 ÷ 0.935
Refractive index at 20°C. 1.4212 ÷ 1.4272
Water (K.F.) ..... ≤ 0.1 %
Residue on evaporation ..... ≤ 50 ppm
Acidity (acetic acid) ..... ≤ 100 ppm
Alcalinity (NH3) ..... ≤ 10 ppm
Heavy metals (Pb) ..... ≤ 2 ppm
Fe ..... ≤ 1 ppm
Assay (GLC) ..... ≥ 99 %

Code	Size	Packaging	Notes
441771	250 ml	Glass bottle	
441774	1 l	Glass bottle	

Diacetone alcohol > RE - Pure

RE

Description ..... Yellow clear liquid
Identification ..... Positive
Density at 25° C ..... 0.926 ÷ 0.936
Refractive index at 20°C. 1.4192 ÷ 1.4292
Water (K.F.) ..... ≤ 0.3 %
Residue on evaporation ..... ≤ 100 ppm
Acidity (acetic acid) ..... ≤ 100 ppm
Assay (GLC) ..... ≥ 98 %

Code	Size	Packaging	Notes
337001	1 l	Glass bottle	
337002	18 kg	Metal drum	

Diacetyldioxime

• Diacetildiossima • Diacétyldioxime • Diacetildioxima • Diacetyldioxim

Synonym:
• *Dimethylglyoxime*
• *2,3-Butanedione dioxime*

$CH_3C{:}NOHC{:}NOHCH_3$
Molecular Weight: 116,12
CAS: 95-45-4
EEC-N: 202-420-1

Warning
H302
P264-P270-P301+P312a-P330-P501a

Diacetyldioxime > RPE - For analysis - ACS

RPE

Description ..... White crystalline powder
Identification ..... Positive
Suitability for Ni det. ..... Conform
Melting point ..... 239 ÷ 241 °C
Alcohol-insolub. matter ..... ≤ 0.05 %

Code	Size	Packaging	Notes
441553	50 g	Glass bottle	

Reagent for the spectrophotometric determination of: Co (III), Fe (II), Ni (II), Pd (II)

Diacetyldioxime sodium salt

• Diacetildiossima sale sodico • Diacétyldioxime sel sodique • Diacetildioxima sal sódica
• Diacetyldioximnatriumsalz

Synonym:
Dimethylglyoxime disodium salt octahydrate

$C_4H_6N_2Na_2O_2.8H_2O$
Molecular Weight: 304,09
CAS: 75006-64-3
EEC-N: 262-523-2

Warning
H302
P264-P270-P301+P312a-P330-P501a

Diacetyldioxime sodium salt > RPE - For analysis

RPE

Description ..... White powder
Identification ..... Positive
Suitability for Ni det. ..... Conform
Water (K.F.) ..... 43 - 48 %
Solubility (50mg/ml,50%ETOH) ..... Conform
Assay (gravimetric) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
441623	50 g	Glass bottle	
441625	250 g	Glass bottle	

For precipitation of metals

Diacetyldioxime solution 1% in ethanol

• Diacetildiossima soluzione 1% in alcole etilico • Diacétyldioxime solution 1% dans l'éthanol
• Diacetildioxima solución 1% en alcohol etilico • Diacetyldioximlösung 1% in Ethanol

Synonym:
• *Dimethylglyoxime*
• *2,3-Butanedione dioxime*

$CH_3C:NOHC:NOHCH_3$
Molecular Weight: 116,12
CAS: 95-45-4

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II



Danger
H225-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

Diacetyldioxime solution 1% in ethanol > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 15° C ..... 0.813 ÷ 0.819

Code	Size	Packaging	Notes
E441581	500 ml	Bottle	

Indicator for determining metals

1,2-Diaminoethane ▶ Ethylenediamine

Diatomaceous earth ▶ Kieselguhr composed

Dibenzoylmethane

• Dibenzoilmetano • Dibenzoylméthane • Dibenzoilmetano • Dibenzoylmethan

$(C_6H_5CO)_2CH_2$
Molecular Weight: 224,26
CAS: 120-46-7
EEC-N: 204-398-9

Dibenzoylmethane > RPE - For analysis

RPE

Description ..... Yellowish crystalline powder
Identification ..... Positive
Melting point ..... 77.0 ÷ 79.0 ° C
Water (K.F.) ..... ≤0.1 %
Residue on ignition ..... ≤500 ppm
Assay (GLC) ..... ≥ 97 %

Code	Size	Packaging	Notes
441873	25 g	Glass bottle	

Dibenzoyl peroxide

• Dibenzoile perossido • Benzoyle péroxyde • Dibenzoilo peróxido • Dibenzoylperoxid

Synonym:
Benzoyl peroxide

$(C_6H_5CO)_2O_2$
Molecular Weight: 242,22
CAS: 94-36-0
EEC-N: 202-327-6

Classification transport
ONU: 3104
Transport Hazard class: 5.2
Packing group -



Danger
H241-H319-H317-H410
P210-P261-P280-P305+P351+P338-
P337+P313-P410

Dibenzoyl peroxide > RE - Pure

RE

Description ..... White granular powder
Identification ..... Positive
Assay (oxidimetric) ..... ≥ 62.4 %

Code	Size	Packaging	Notes
427345	250 g	Plastic bottle	
427347	1 kg	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Di-n-butylphthalate

• n-Dibutilftalato • n-Dibutyle phtalate • n-Dibutiloftalato • Dibutylphthalat

Synonym:
Phtalic acid dibutyl ester

$C_6H_4[COOC_4H_9]_2$
Molecular Weight: 278,35
CAS: 84-74-2
EEC-N: 201-557-4



Danger
H360Df-H400-HA26
P273-P280-P308+P313-P391-P405-P501a

Di-n-butylphthalate > ERBAPharm - According to pharmacopoeia: Ph.Eur.

ERBApharm

Description .Clear colourless to very slightly yellow liq.
Appearance of solution .....Conform Ph.Eur.
Identification ..... Positive
Acidity .....Conform Ph.Eur.
Related substances (GLC) ..... ≤ 1.0 %
Water (K.F.) ..... ≤ 0.2 %
Sulphated ash ..... ≤ 0.1 %
Assay (alkalimetric).....99.0 ÷ 101.0 %
Density at 20°C .....1.043 ÷ 1.048
Refractive index at 20°C.....1.490 ÷ 1.495
Origin (BSE/TSE)..... Synthesis
Residual solvents (Current ICH).....Conform

Code	Size	Packaging	Notes
325701	26 kg	Metal drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Dicalite 4158

• Dicalite 4158 • Dicalite 4158 • Dicalite 4158 • Dicalite 4158

CAS: 93763-70-3

Dicalite 4158 > RE - For filtration

RE

Permeability (PRFv)..... 90 - 165
Wet cake density.....≤ 13.6 lbs/ft³
Float.....≤ 30 ml/20g

Code	Size	Packaging	Notes
P8880014	500 g	Plastic bottle	
P8880017	1 kg	Plastic bottle	
P8880027	5 kg	Plastic bucket	

Substitute to CELITE 545

Dichloroacetic acid

• Acido dicloroacetico • Acide dichloroacétique • Acido dicloroacético • Dichloressigsäure

$CHCl_2COOH$
Molecular Weight: 128,94
CAS: 79-43-6
EEC-N: 201-207-0

Classification transport
ONU: 1764
Transport Hazard class: 8
Packing group II



Danger
H311-H314-H400
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P361+P364

Dichloroacetic acid > RPE - For analysis

RPE

Description ..... Clear liquid
Identification ..... Positive
Colour ..... ≤ 50 APHA
Density at 20° C .....1.562 ÷ 1.572
Refractive index at 20°C.....1.4648 ÷ 1.4668
Assay (acidimetric) .....99.0 ÷ 101.0 %

Code	Size	Packaging	Notes
405101	250 ml	Glass bottle	
405103	1 l	Glass bottle	

Dichloroacetic acid > RE - Pure

RE

Description ..... Clear liquid
Identification ..... Positive
Freezing point .....≥ 12 °C
Water (KF).....≤ 0.3 %
Assay (acidimetric) .....≥ 99 %

Code	Size	Packaging	Notes
303151	1 l	Glass bottle	

m-Dichlorobenzene

• m-Diclorobenzolo • m-Dichlorobenzène • m-Diclorobenceno • m-Dichlorbenzol

Synonym: *1,3-Dichlorobenzene*

$C_6H_4Cl_2$
Molecular Weight: 147
CAS: 541-73-1
EEC-N: 208-792-1

Classification transport
ONU: 2810
Transport Hazard class: 6.1
Packing group III



Warning
H302-H411
P264-P270-P301+P312a-P330-P391-P501a

m-Dichlorobenzene > RPE - For analysis

RPE

Description .... Colourless or yellowish liquid
Identification ..... Positive
Density at 20° C .....1.283 ÷ 1.293
Refractive index at 20°C.1.5435 ÷ 1.5485
Boiling point.....171 ÷ 173 °C
Residue on ignition.....≤ 100 ppm
Assay (GLC) ..... ≥ 99 %

Code	Size	Packaging	Notes
442353	100 ml	Glass bottle	

o-Dichlorobenzene

• o-Diclorobenzolo • o-Dichlorobenzène • o-Diclorobenceno • o-Dichlorbenzol

Synonym: *1,2-Dichlorobenzene*

$C_6H_4Cl_2$
Molecular Weight: 147
CAS: 95-50-1
EEC-N: 202-425-9

Classification transport
ONU: 1591
Transport Hazard class: 6.1
Packing group III



Warning
H302-H332-H315-H319-H335-H410
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

o-Dichlorobenzene > RPE - For analysis

RPE

Description ..... Clear liquid
Identification ..... Positive
Density at 20° C ..... 1.299 ÷ 1.311
Refractive index at 20°C..... 1.5450 ÷ 1.5520
Water (K.F.) .....≤500 ppm
Residue on evaporation .....≤50 ppm
Acidity (HCl) .....≤50 ppm
Free chlorine.....≤1 ppm
Assay (GLC) .....≥99 %

Code	Size	Packaging	Notes
442301	500 ml	Glass bottle	

o-Dichlorobenzene > RE - Pure

RE

Description ..... Clear liquid
Identification ..... Positive
Assay (GLC) .....≥98 %

Code	Size	Packaging	Notes
337251	1 l	Glass bottle	

p-Dichlorobenzene

• p-Diclorobenzolo • p-Dichlorobenzène • p-Diclorobenceno • p-Dichlorbenzol

Synonym: *1,4-Dichlorobenzene*

$C_6H_4Cl_2$
Molecular Weight: 147
CAS: 106-46-7
EEC-N: 203-400-5



Warning
H319-H351-H410
P264-P280-P305+P351+P338-P308+P313-P337+P313-P501a

p-Dichlorobenzene > RPE - For analysis

RPE

Description ..... White crystals
Identification ..... Positive
Melting point..... 52.0 ÷ 54.0 ° C
Residue on ignition..... ≤500 ppm
Assay (GLC) .....≥99 %

Code	Size	Packaging	Notes
442406	500 g	Plastic bottle	
442407	1 kg	Plastic bottle	

p-Dichlorobenzene > RE - Pure

RE

Description .....White crystals
Identification ..... Positive
Melting point..... 52 ÷ 56 ° C
Assay (GLC) ..... ≥ 96.0 %

Code	Size	Packaging	Notes
337307	1 kg	Plastic bottle	
337303	25 kg	Plastic bucket	

Product specifications are subject to changes. Please visit our website for updates.

1,2-Dichlorobenzene-d4

• 1,2-Diclorobenzene-d4 • 1,2-Dichlorobenzène-d4 • 1,2-Diclorobenceno-d4 • 1,2-Dichlorbenzol-d4

Synonym:
Tetradeutero-1,2-dichlorobenzene

$C_6D_4Cl_2$
Molecular Weight: 151,03
CAS: 2199-69-1

Classification transport
ONU: 1591
Transport Hazard class: 6.1
Packing group III

Warning
H302-H315-H319-H335-H410
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

1,2-Dichlorobenzene-d4 > RS - For NMR - min 99%

RS

Code	Size	Packaging	Notes
P5533A	5 ml	Glass ampoule	

For specifications, contact our customer service for a certificate of analysis

1,2-Dichloroethane

• 1,2-Etano dicloro • 1,2-Dichloroéthane • 1,2-Dicloroetano • 1,2-Dichlorethan

Synonym:
• *Ethylene dichloride*
• *sym-Dichlorethane*

CH_2ClCH_2Cl
Molecular Weight: 98,97
CAS: 107-06-2
EEC-N: 203-458-1

Classification transport
ONU: 1184
Transport Hazard class: 3
Packing group II

Danger
H225-H302-H331-H315-H319-H350-H335-HA26
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P308+P313-P330-P362+P364-P403+P233

1,2-Dichloroethane > RS - For HPLC - Isocratic Grade

RS

Clear,colourless liq.appearance ..... Conform	Water content (K.F.) ..... ≤ 100 mg/Kg	UV transmittance at 255 nm ..... ≥ 90 %	Non volatile residue ..... ≤ 5 mg/Kg
Identification ..... Conform	UV transmittance at 240 nm ..... ≥ 20 %	UV transmittance at 260 nm ..... ≥ 98 %	
Colour ..... ≤ 10 Apha	UV transmittance at 245 nm ..... ≥ 55 %	Free acid (as HCl) ..... ≤ 5 mg/Kg	
Refractive index at 20°C ..... 1.443 - 1.447	UV transmittance at 250 nm ..... ≥ 78 %	Assay (GC) ..... ≥ 99.8 %	

Code	Size	Packaging	Notes
447191	1 l	Glass bottle	
447192	2.5 l	Glass bottle	

1,2-Dichloroethane > RS - SPECTROSOL - For optical spectroscopy

RS

Clear,colourless liq.appearance ..... Conform	Colour ..... ≤ 10 Hazen	UV transmittance at 240 nm ..... ≥ 90 %	Free acid (as HCl) ..... ≤ 5 mg/Kg
Refractive index at 20°C ..... 1.443 - 1.447	UV transmittance at 225 nm ..... ≥ 10 %	UV transmittance at 250 nm ..... ≥ 95 %	Assay (GC) ..... ≥ 99.8 %
Water content (K.F.) ..... ≤ 100 mg/Kg	UV transmittance at 230 nm ..... ≥ 50 %	UV transmittance at 260 nm ..... ≥ 99 %	Non volatile residue ..... ≤ 5 mg/Kg

Code	Size	Packaging	Notes
P0282716	1 l	Glass bottle	

1,2-Dichloroethane > RS - Anhydrous - For analysis

RS

Refractive index at 20°C ..... 1.443 - 1.447	Non volatile residue ..... ≤ 10 mg/Kg	Free acid (as HCl) ..... ≤ 5 mg/Kg
Water content (K.F.) ..... ≤ 100 mg/Kg	Colour ..... ≤ 10 Hazen	Assay (GC) ..... ≥ 99.8 %

Code	Size	Packaging	Notes
P0281010	200 ml	Bottle with septum	
P0281016	1 l	Glass bottle	
P0281021	2.5 l	Glass bottle	

1,2-Dichloroethane > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Clear liquid	Distillation range (≥ 95%) ..... 82 ÷ 84 ° C	Total phosphorus ..... ≤0.5 ppm	K ..... ≤0.2 ppm
Colour (APHA) ..... ≤10	Acidity ..... ≤0.0003 meq/g	Total silicon ..... ≤0.05 ppm	Mg ..... ≤0.1 ppm
Identification ..... Positive	Water (K.F.) ..... ≤200 ppm	Total sulphur ..... ≤0.5 ppm	Na ..... ≤0.5 ppm
Density at 20° C ..... 1.248 ÷ 1.264	Residue on evaporation ..... ≤10 ppm	Ca ..... ≤0.5 ppm	Pb ..... ≤0.02 ppm
Refractive index at 20°C ..... 1.4418 ÷ 1.4478	Free chlorine ..... ≤1 ppm	Cu ..... ≤0.05 ppm	Zn ..... ≤0.2 ppm
Boiling point ..... 83.0 ÷ 84.0 ° C	Subst. reducing KMnO4 ..... ≤10 ppm	Fe ..... ≤0.5 ppm	Assay (GLC) ..... ≥99.8 %

Code	Size	Packaging	Notes
447121000	1 l	Glass bottle	

1,2-Dichloroethane > RE - Pure

RE

Description ..... Clear liquid
Water content (K.F.) ..... ≤ 300 mg/Kg
Identification ..... Positive
Non volatile residue ..... ≤ 50 mg/Kg
Refractive index at 20°C ..... 1.4398 ÷ 1.4498
Colour ..... ≤ 10 Hazen
Colour (APHA) ..... ≤20
Free acid (as HCl) ..... ≤ 10 mg/Kg
Water (K.F.) ..... ≤300 ppm
Assay (GC) ..... ≥ 99.8 %
Assay (GLC) ..... ≥99.8 %
Residue on evaporation ..... ≤ 50 ppm
Acidity (HCl) ..... ≤ 10 ppm

Code	Size	Packaging	Notes
340151	1 l	Glass bottle	
P0280228	5 l	Plastic tank	
340155	34 kg	Metal drum	
P0280268	200 l	Metal drum	

Dichloromethane

• Diclorometano • Dichlorométhane • Diclorometano • Dichlormethan

Synonym:
Methylene chloride

CH_2Cl_2
Molecular Weight: 84,93
CAS: 75-09-2
EEC-N: 200-838-9

Classification transport
ONU: 1593
Transport Hazard class: 6.1
Packing group III

Warning
H315-H319-H351-H335-H336-H373
P271-P280-P304+P340-P305+P351+P338-P308+P313-P403+P233

Dichloromethane > RS - For HPLC - Isocratic grade - Stabilized with amylene

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 1.322 ÷ 1.328
Refractive index at 20°C ..... 1.4214 ÷ 1.4274
Boiling point ..... 39.6 ÷ 40.1 ° C
Acidity ..... ≤ 0.0001 meq/g
Water (K.F.) ..... ≤100 ppm
Residue on evaporation ..... ≤2 ppm
Assay (GLC) ..... ≥99.9 %
U.V. Transmittance
At 235 nm ..... ≥ 40 %
at 240 nm ..... ≥70 %
at 250 nm ..... ≥96 %
at 260 nm ..... ≥99 %
Amylene ..... ≤ 60 ppm
Alkalinity ..... ≤ 0.0002 meq/g

Code	Size	Packaging	Notes
412621000	1 l	Glass bottle	
412622000	2.5 l	Glass bottle	

Dichloromethane > RS - For HPLC - Isocratic grade - Stabilized with ethanol

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Boiling point ..... 39.55 ÷ 40.05 °C
Density at 20°C ..... 1.322 ÷ 1.328
Refractive index at 20°C ..... 1.4214 ÷ 1.4274
Acidity (Formic acid) ..... ≤ 5 ppm
Residue on evaporation ..... ≤ 5 ppm
Ethanol ..... 0.1 ÷ 0.4 %
Water (K.F.) ..... ≤ 100 ppm
Transmittance U.V.
at 235 nm ..... ≥ 40 %
at 240 nm ..... ≥ 75 %
at 245 nm ..... ≥ 90 %
at 255 nm ..... ≥ 96 %
at 260 nm ..... ≥ 99 %
Assay (GLC) ..... ≥ 99.9 %

Code	Size	Packaging	Notes
412662	1 l	Glass bottle	
412661	2.5 l	Glass bottle	

Dichloromethane > RS - For preparative HPLC - Stabilized with amylene

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 1.322 ÷ 1.328
Refractive index at 20°C ..... 1.4214 ÷ 1.4274
Boiling point ..... 39.6 ÷ 40.1 ° C
Water (K.F.) ..... ≤500 ppm
Residue on evaporation ..... ≤5 ppm
Assay (GLC) ..... ≥99.5 %
U.V. Transmittance
at 240 nm ..... ≥50 %
at 260 nm ..... ≥98 %
Amylene ..... 20 ÷ 60 ppm
Alkalinity ..... ≤ 0.0002 meq/g

Code	Size	Packaging	Notes
463281	2.5 l	Glass bottle	

Dichloromethane > RS - For preparative HPLC - Stabilized with ethanol

RS

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Identity (IR) ..... Positive
Density at 20°C ..... 1.322 ÷ 1.328
Refractive index at 20°C ..... 1.4214 ÷ 1.4274
Boiling point ..... 39.3 ÷ 40.3 °C
Alkalinity (NH3) ..... ≤ 5 ppm
Water (K.F.) ..... ≤ 500 ppm
Residue on evaporation ..... ≤ 5 ppm
Assay (GLC) ..... ≥ 99.5 %
Transmittance U.V.
at 240 nm ..... ≥ 50 %
at 260 nm ..... ≥ 98 %

Code	Size	Packaging	Notes
463291	2.5 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Dichloromethane > RS - For GC-MS - Stabilized with ethanol

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C ..... 1.422 - 1.426
Water (K.F.) ..... ≤ 100 ppm
Colour ..... ≤ 10 APHA
Residue on evaporation ..... ≤ 3 ppm
Acidity (HCl) ..... ≤ 5 ppm
Assay (GC) (without stabilizer) ... ≥ 99.95 %
Stabilizer (Ethanol) ..... 0.1 - 0.4 % (w/w)
Carbon tetrachloride ..... ≤ 20 ppm
Chloroform ..... ≤ 40 ppm
GC-MS.Individual peak (n-hexadecane). ≤ 2 µg/L
Ret.range n-undecane to n-tetracontane

Code	Size	Packaging	Notes
463332	1 l	Glass bottle	

Dichloromethane > RS - For GC-MS - Stabilized with amylene

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C ..... 1.422 - 1.426
Water (K.F.) ..... ≤ 100 ppm
Colour ..... ≤ 10 APHA
Residue on evaporation ..... ≤ 3 ppm
Acidity (HCl) ..... ≤ 5 ppm
Assay (GC) ..... ≥ 99.95 %
Amylene ..... 30 - 60 ppm
Carbon tetrachloride ..... ≤ 20 ppm
Chloroform ..... ≤ 50 ppm
GC-MS.Individual peak (n-hexadecane). ≤ 2 µg/L
Ret.range n-undecane to n-tetracontane

Code	Size	Packaging	Notes
463342	1 l	Glass bottle	

Dichloromethane > RS - ATRASOL - For analysis of volatile traces - Stabilized with amylene

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C ..... 1.422 - 1.426
Water content (K.F.) ..... ≤ 100 mg/Kg
Colour ..... ≤ 10 Hazen
Free acid (as HCl) ..... ≤ 5 mg/Kg
Chloroform ..... ≤ 50 mg/Kg
Non volatile residue ..... ≤ 3 mg/Kg
Carbon tetrachloride ..... ≤ 20 mg/Kg
Assay (GC) ..... ≥ 99.95 %
Stabilizer (Amylene) ..... 20 - 60 mg/Kg
GC-ECD.Individual peak (Lindane) .≤ 2 ng/L
Ret.range 1,2,4-trichlorobenzene
GC-FID.Individual peak (n-hexadecane). ≤ 2 µg/L
Ret.range n-undecane to n-tetracontane
PAH test (according to ISO 17993) .. Passes test
Test against standard (each 16 PAH) . ≤ 0.1 µg/l to decachlorobiphenyle

Code	Size	Packaging	Notes
P02932A16	1 l	Glass bottle	
P02932A21	2.5 l	Glass bottle	
P02932A82	4 l	Glass bottle	

16 selected PAHs tested according to ISO 17993:2002, each max 0,1µg/l

Dichloromethane > RS - ATRASOL - For analysis of volatile traces - Stabilized with ethanol

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C ..... 1.422 - 1.426
Water content (K.F.) ..... ≤ 100 mg/Kg
Colour ..... ≤ 10 Hazen
Free acid (as HCl) ..... ≤ 5 mg/Kg
Chloroform ..... ≤ 40 mg/Kg
Non volatile residue ..... ≤ 3 mg/Kg
Carbon tetrachloride ..... ≤ 20 mg/Kg
Assay (GC) (without stabilizer) ... ≥ 99.95 %
Stabilizer (Ethanol) ..... 0.1 - 0.4 % m/m
GC (FID) - NC Atrasol ..... Conform
GC-ECD.Individual peak (Lindane) .≤ 2 ng/L
Ret.range 1,2,4-trichlorobenzene to decachlorobiphenyle
GC-FID.Individual peak (n-hexadecane). ≤ 2 µg/L
Ret.range n-undecane to n-tetracontane
PAH test (according to ISO 17993) .. Passes test
Test against standard (each 16 PAH) . ≤ 0.1 µg/l

Code	Size	Packaging	Notes
P02932E16	1 l	Glass bottle	
P02932E21	2.5 l	Glass bottle	

16 selected PAHs tested according to ISO 17993:2002, each max 0,1µg/l

Dichloromethane > RS - PESTIPUR - For pesticide analysis - Stabilized with amylene

RS

Description ..... Clear liquid
Identification ..... Positive
Colour ..... ≤ 10 hazen
Assay (GLC) ..... ≥ 99.9 %
Water ..... ≤ 0.01 %
Not volatile residue ..... ≤ 5 mg/Kg
Free acid (as HCl) ..... ≤ 5 mg/Kg
GC-ECD (Lindane standard) ..... ≤ 3 ng/l
GC-NPD (Ethylparathion standard) .. ≤ 3 ng/l
Stabilizer (Amylene) ..... 20 - 50 mg/Kg
PAH test (according to ISO 17993) .. Passes test
Test against standard (each 16 PAH) .. ≤0.1 µg/l

Code	Size	Packaging	Notes
442291	1 l	Glass bottle	
442292000	2.5 l	Glass bottle	
442294	4 l	Glass bottle	

16 selected PAHs tested according to ISO 17993:2002, each max 0,1µg/l

Dichloromethane > RS - PESTIPUR - For pesticide analysis - Stabilized with ethanol

RS

Refractive index at 20°C ..... 1.422 - 1.426
Water content (K.F.) ..... ≤ 100 mg/Kg
Colour ..... ≤ 10 Hazen
Stabilizer (Ethanol) ..... 0.1 - 0.4 % m/m
Free acid (as HCl) ..... ≤ 5 mg/Kg
Non volatile residue ..... ≤ 5 mg/Kg
Assay (GC) (without stabilizer) ..... ≥ 99.9 %
GC-ECD.Individual peak (Lindane) .. ≤ 3 ng/l
Retention time trichlorobenzene to mirex GC-NPD.Individual peak (Ethylparathion) ≤ 3 ng/l
Retention time Atrazin to Coumaphos
PAH test (according to ISO 17993) .. Passes test
Test against standard (each 16 PAH) .. ≤0.1 µg/l

Code	Size	Packaging	Notes
442261	1 l	Glass bottle	
442262	2.5 l	Glass bottle	

16 selected PAHs tested according to ISO 17993:2002

Dichloromethane > RS - SPECTROSOL - For optical spectroscopy - Stabilized with amylene

RS

Description ..... Clear colourless liquid
Colour (APHA) ..... ≤10
Identification (I.R.) ..... Positive
Density at 20° C ..... 1.322 ÷ 1.328
Refractive index at 20°C. 1.4214 ÷ 1.4274
Boiling point ..... 39.6 ÷ 40.1 ° C
Water (K.F.) ..... ≤100 ppm
Residue on evaporation ..... ≤5 ppm
Acidity ..... ≤0.0005 meq/g
Assay (GLC) ..... ≥99.9 %
Fluorescence
at 254 nm ..... ≤2 ppb
at 365 nm ..... ≤2 ppb
U.V. Transmittance
at 235 nm ..... ≥35 %
at 240 nm ..... ≥70 %
at 245 nm ..... ≥90 %
at 250 nm ..... ≥95 %
at 260 nm ..... ≥99 %
Amylene ..... 20 ÷ 60 ppm
Alkalinity ..... ≤ 0.0002 meq/g

Code	Size	Packaging	Notes
442371	1 l	Glass bottle	
P02927A21	2.5 l	Glass bottle	

Dichloromethane > RS - SPECTROSOL - For optical spectroscopy - Stabilized with ethanol

RS

Description ..... Clear colourless liquid
Colour (APHA) ..... ≤10
Identification (I.R.) ..... Positive
Density at 20° C ..... 1.322 ÷ 1.328
Refractive index at 20°C. 1.4214 ÷ 1.4274
Boiling point ..... 39.6 ÷ 40.1 ° C
Water (K.F.) ..... ≤100 ppm
Residue on evaporation ..... ≤5 ppm
Acidity ..... ≤0.0005 meq/g
Alcalinity ..... ≤0.0002 meq/g
Assay (GLC) ..... ≥99.9 %
Fluorescence
at 254 nm ..... ≤2 ppb
at 365 nm ..... ≤2 ppb
U.V. Transmittance
at 235 nm ..... ≥35 %
at 240 nm ..... ≥70 %
at 245 nm ..... ≥90 %
at 250 nm ..... ≥95 %
at 260 nm ..... ≥99 %
Stabilizer (Ethanol) ..... 0.1 ÷ 0.4 %

Code	Size	Packaging	Notes
463025	1 l	Glass bottle	

Dichloromethane > RS - Anhydrous - For analysis - Stabilized with amylene

RS

Refractive index at 20°C ..... 1.422 - 1.426
Water content (K.F.) ..... ≤ 50 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 99.95 %
Stabilizer (Amylene) ..... 20 - 60 mg/Kg
Free acid (as HCl) ..... ≤ 3 mg/Kg
Carbon tetrachloride ..... ≤ 20 mg/Kg
Chloroform ..... ≤ 50 mg/Kg

Code	Size	Packaging	Notes
P02910A10	200 ml	Bottle with septum	
P02910AT10	200 ml	Bottle with septum	On molecular sieves 4A
P02910A16	1 l	Glass bottle	
P02910AT16	1 l	Glass bottle	On molecular sieves 4A
P02910A21	2.5 l	Glass bottle	
P02910A48	25 l	Metal drum	

Dichloromethane > RS - Anhydrous - For analysis - Stabilized with ethanol

RS

Refractive index at 20°C ..... 1.422 - 1.426
Water content (K.F.) ..... ≤ 50 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) (without stabilizer) ..... ≥ 99.95 %
Stabilizer (Ethanol) ..... 0.1 - 0.4 % m/m
Free acid (as HCl) ..... ≤ 3 mg/Kg
Carbon tetrachloride ..... ≤ 20 mg/Kg
Chloroform ..... ≤ 40 mg/Kg

Code	Size	Packaging	Notes
P02910E10	200 ml	Bottle with septum	
P02910E16	1 l	Glass bottle	
P02910E21	2.5 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Dichloromethane > RS - RSE - For electronic use - Stabilized with amylene

RS

Description	Clear liquid	Heavy metals (Pb)	≤0.1 ppm	Cr	≤0.01 ppm	Pb	≤0.02 ppm
Colour (APHA)	≤10	Ag	≤0.02 ppm	Cu	≤0.01 ppm	Pt	≤0.05 ppm
Identification	Positive	Al	≤0.05 ppm	Fe	≤0.1 ppm	Sb	≤0.01 ppm
Density at 20° C	1.322 ÷ 1.328	As	≤0.01 ppm	Ga	≤0.02 ppm	Sn	≤0.02 ppm
Boiling point	39.6 ÷ 40.1 ° C	Au	≤0.05 ppm	In	≤0.02 ppm	Sr	≤0.02 ppm
Resistivity	≥1 Mohm.cm	B	≤0.01 ppm	K	≤0.1 ppm	Ti	≤0.05 ppm
Assay (GLC)	≥99.5 %	Ba	≤0.1 ppm	Li	≤0.02 ppm	Tl	≤0.05 ppm
Water (K.F.)	≤100 ppm	Be	≤0.02 ppm	Mg	≤0.1 ppm	V	≤0.05 ppm
Residue on evaporation	≤5 ppm	Bi	≤0.02 ppm	Mn	≤0.01 ppm	Zn	≤0.05 ppm
Acidity (HCl)	≤10 ppm	Ca	≤0.2 ppm	Mo	≤0.05 ppm	Zr	≤0.05 ppm
Chloride	≤0.5 ppm	Cd	≤0.01 ppm	Na	≤0.2 ppm	Stabilized with amylene	20 ÷ 60 ppm
Free chlorine	≤0.2 ppm	Co	≤0.01 ppm	Ni	≤0.01 ppm	Alcalinity (NH3)	≤ 1 ppm

Code	Size	Packaging	Notes
463162	1 l	Glass bottle	
463161	2.5 l	Glass bottle	

Dichloromethane > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP - Stabilized with amylene

RPE

Description	Clear liquid	Boiling point	39.6 ÷ 40.1 ° C	Ba	≤0.1 ppm	Mn	≤0.02 ppm
Colour (APHA)	≤10	Water (K.F.)	≤100 ppm	Ca	≤0.5 ppm	Ni	≤0.02 ppm
Identification (I.R.)	Conform	Residue on evaporation	≤10 ppm	Cd	≤0.05 ppm	Pb	≤0.1 ppm
Ready carbonizable substances	Conform	Acidity	≤0.0003 meq/g	Co	≤0.02 ppm	Sn	≤0.1 ppm
Alcohol miscibility	Complete	Amylene	20 ÷ 60 ppm	Cr	≤0.02 ppm	Zn	≤0.1 ppm
Diethyl ether miscib.	Complete	Chloride	≤1 ppm	Cu	≤0.02 ppm	Assay (GLC)	≥99.9 %
Density at 20° C	1.322 ÷ 1.328	Al	≤0.5 ppm	Fe	≤0.1 ppm	Free halogens	Conform
Refractive index at 20°C	1.4214 ÷ 1.4274	B	≤0.02 ppm	Mg	≤0.1 ppm		

Code	Size	Packaging	Notes
463311	1 l	Glass bottle	
463314	2.5 l	Glass bottle	
524319	10 l	Plastic tank	
463318	20 kg	Drum	
524314	200 l	Metal drum	

Dichloromethane > RPE - For analysis - ACS - Stabilized with ethanol

RPE

Description	Clear colourless liquid	Free chlorine	≤ 0.1 ppm	Ba (Barium)	≤ 0.1 ppm	Mn (manganese)	≤0.02 ppm
Color	≤10 APHA	Identity (IR)	Positive	Ca (Calcium)	≤ 0.5 ppm	Ni (Nickel)	≤ 0.02 ppm
Density at 20°C	1.322 ÷ 1.328	Read. carboniz. subs.	Conform	Cd (Cadmium)	≤ 0.05 ppm	Pb (Lead)	≤ 0.1 ppm
Refractive index at 20°C	1.4214 ÷ 1.4274	Residue on evaporation	≤ 10 ppm	Co (Cobalt)	≤ 0.02 ppm	Sn (Tin)	≤0.1 ppm
Boiling point	39.3 ÷ 40.3 °C	Water (K.F.)	≤ 100 ppm	Cr (Chromium)	≤ 0.02 ppm	Zn (Zinc)	≤ 0.1 ppm
Acidity (Hydrochloric ac.)	≤ 3 ppm	Assay (GLC)	≥ 99.9 %	Cu (Copper)	≤ 0.02 ppm		
Chloride (Cl)	≤ 1 ppm	Al (Aluminium)	≤ 0.5 ppm	Fe (Iron)	≤ 0.1 ppm		
Ethyl alcohol	≤ 0.2 %	B (Boron)	≤ 0.02 ppm	Mg (magnesium)	≤0.1 ppm		

Code	Size	Packaging	Notes
463001	1 l	Glass bottle	
463003	2.5 l	Glass bottle	
463002	5 l	Plastic tank	
463008	250 kg	Metal drum	

Dichloromethane > ERBApharm - According to pharmacopoeia: NF - Stabilized with ethanol

ERBApharm

Description	Conform	Residue on evaporation	≤ 20 ppm	Free chlorine	Conform USP-NF	Ethanol	# 0,2 %
Identification (I.R.)	Positive	Acidity	≤ 10 ppm	Water (K.F.)	≤ 200 ppm		
Density d25/25	1,318 ÷ 1,322	Heavy metals (Pb)	≤ 1 ppm	Assay (GLC)	≥ 99.0 %		

Code	Size	Packaging	Notes
354501	1 l	Glass bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Dichloromethane > ERBApharm - According to pharmacopoeia: Ph.Eur.-NF - Stabilized with amylene

ERBApharm

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Identification ..... Positive
Appearance of solution ..... Clear colourless liquid Ph. Eur.
Density at 20°C ..... 1.320 - 1.332
Density at 25°C ..... 1.318 ÷ 1.322
Refractive index at 20°C ..... 1.423 - 1.425
Acidity ..... Pass test Ph.Eur.
Water (K.F.) ..... ≤ 200 ppm
Residue on evaporation ..... ≤ 20 ppm
Heavy metals (Pb) ..... ≤ 1 ppm
Free chlorine ..... Conform USP-NF
Related substances (CPG) ..... ≤ 0.1 %
Assay (GLC) ..... ≥ 99.9 %
Amylene ..... 20 ÷ 60 ppm
Carbon tetrachloride ..... ≤ 10 ppm (v/v)
Chloroforme ..... ≤ 50 ppm (v/v)
Origin (BSE/TSE) ..... Synthesis

Code	Size	Packaging	Notes
337331	1 l	Glass bottle	
337333	2.5 l	Glass bottle	
337335	25 l	Metal drum	
337337	200 l	Metal drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Dichloromethane > ERBApharm - According to pharmacopoeia: Ph.Eur. - Stabilized with ethanol

ERBApharm

Description ..... Conform Ph.Eur.
Colour ..... ≤ 10 APHA
Identification ..... Positive Ph.Eur.
Appearance of solution ..... Clear colourless liquid Ph.Eur.
Density at 20°C ..... 1.320 ÷ 1.332
Indice di rifraz. ..... 1.423 ÷ 1.425
Acidity ..... Conform Ph.Eur.
Residue on evaporation ..... ≤ 10 ppm
Water (K.F.) ..... ≤ 100 ppm
Free chlorine ..... Conform Ph.Eur.
Heavy metals (Pb) ..... ≤ 1 ppm
Related compounds ..... ≤ 0.1 %
Assay (GLC) ..... ≥ 99.95 %
Ethanol ..... ≤ 2.0 %
Amylene ..... ≤ 0.03 %
Carbon tetrachloride ..... ≤ 10 ppm (v/v)
Chloroforme ..... ≤ 50 ppm (v/v)

Code	Size	Packaging	Notes
525320	2.5 l	Glass bottle	
525321	200 l	Metal drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Dichloromethane > RE - Pure - stabilized with amylene

RE

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Density at 20°C ..... 1.320 ÷ 1.330
Refractive index at 20°C ..... 1.4194 ÷ 1.4294
Boiling point ..... 39.1 ÷ 40.6 °C
Acidity (HCl) ..... ≤ 10 ppm
Residue on evaporation ..... ≤ 20 ppm
Water (K.F.) ..... ≤ 200 ppm
Assay (GLC) ..... ≥ 99.8 %
Amylene ..... 20 ÷ 60 ppm

Code	Size	Packaging	Notes
528461	5 l	Plastic tank	
528464	10 l	Plastic tank	
528463	25 l	Plastic tank	
528462	200 l	Metal drum	

Dichloromethane > RE - Pure - Stabilized with ethanol

RE

Description ..... Clear colourless liquid
Color ..... <10 APHA
Density at 20°C ..... 1.32÷1.33
Refractive index at 20°C ... 1.4194÷1.4294
Boiling point ..... 39.3÷40.3 °C
Acidity (hydrochloric acid) ..... ≤ 5 ppm
Residue on evaporation ..... ≤20 ppm
Water (K.F.) ..... ≤ 200 ppm
Alcohol and ether miscibility ..... Complete
Assay (GLC) ..... ≥99.8 %
Ethanol (stab) ..... ~0,2 %

Code	Size	Packaging	Notes
528377	2.5 l	Glass bottle	
337315	5 l	Metal tank	
528372	5 l	Plastic tank	
528379	10 l	Plastic tank	
528370	25 l	Metal drum	
528371	200 l	Metal drum	

Dichloromethane > RE - Pure - Stabilized with cyclohexane

RE

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Assay (GLC) ..... ≥ 99.5 %
Density at 20°C ..... 1.320 ÷ 1.330
Refractive index at 20°C 1.4194 ÷ 1.4294
Boiling point ..... 39.0 ÷ 40.5 °C
Acidity ..... ≤ 1 %
Residue on evaporation ..... ≤ 5 ppm
Water (K.F.) ..... ≤ 500 ppm
Methyl alcohol ..... ≤ 500 ppm
Ethyl alcohol ..... ≤ 500 ppm
Chloroform ..... ≤ 1 %
Other impurities ..... ≤ 0.5 %
Cyclohexane ..... 0.0165 ÷ 0.0660 %(v/v)

Code	Size	Packaging	Notes
508370	1 l	Glass bottle	
508374	5 l	Plastic tank	
508375	200 l	Metal drum	

DASIT GROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

Dichloromethane acidified with 1% hydrochloric acid

• Diclorometano acidificato con acido cloridrico 1% • Dichlorométhane acidifié avec 1% d'acide chlorhydrique
• Diclorometano acidificado con ácido clorhídrico 1% • Dichlormethan mit 1% Salzsäure angesäuert

Synonym:
Methylene chloride

CH_2Cl_2
Molecular Weight: 84,93
CAS: 75-09-2
EEC-N: 200-838-9

Classification transport
ONU: 1593
Transport Hazard class: 6.1
Packing group III

Warning
H315-H319-H351-H335-H336-H373
P271-P280-P304+P340-P305+P351+P338-P308+P313-P403+P233

Dichloromethane acidified with 1% hydrochloric acid > RS - For analysis according to Ph. Eur. Chap. 4.1.1 — RS

Code	Size	Packaging	Notes
611055901	100 ml	Glass bottle	Ref Ph.Eur 1055901

Dichloromethane-d2

• Diclorometano-d2 • Dichlorométhane-d2 • Diclorometano-d2 • Dichlormethan-d2

Synonym:
Dideuteromethylenechloride

CD_2Cl_2
Molecular Weight: 86,95
CAS: 1665-00-5
EEC-N: 216-776-0

Classification transport
ONU: 1593
Transport Hazard class: 6.1
Packing group III

Warning
H315-H319-H351-H335-H336-H373
P271-P280-P304+P340-P305+P351+P338-P308+P313-P403+P233

Dichloromethane-d2 > RS - For NMR - min 99.5% — RS

Code	Size	Packaging	Notes
P5330	10 x 0.6 ml	Glass ampoule	
P5339	10 x 0.75 ml	Glass ampoule	
P5335	25 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

2,6-Dichlorophenolindophenol sodium salt

• 2,6-Diclorofenoloindofenolo sale sodico • 2,6-Dichlorophénolindophénol sel de sodium
• 2,6-Diclorofenol-indofenol sal sódica • 2,6-Dichlorphenolindophenol-Natriumsalz

Synonym:
2,6-Dichloro-N-(4-hydroxyphenyl)-1,4-benzoquinoneimine sodium salt

$C_{12}H_7Cl_2NO_2.Na$
Molecular Weight: 326,09
CAS: 620-45-1
EEC-N: 210-640-4

2,6-Dichlorophenolindophenol sodium salt > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP — RPE

Description ..... Dark green powder | Identification ..... Positive | Colouring interferences..... Conform | Loss on drying 120° C..... ≤12 %

Code	Size	Packaging	Notes
442508	5 g	Glass bottle	

For the determination of ascorbic acid

2,6-Dichloroquinone-4-chlorimide

• 2,6-Diclorochinone-4-clorimmide • 2,6-Dichloroquinone-4-chloroimide • 2,6-Dicloroquinona-4-clorimida
• 2,6-Dichlorchinon-4-chlorimid

Synonym:
• *Gibb's reagent*
• *N-2,6-trichloro-p-benzoquinoneimide*

CIN:CCH:CClCOCCl:CH
Molecular Weight: 210,45
CAS: 101-38-2
EEC-N: 202-937-2

Classification transport
ONU: 3224
Transport Hazard class: 4.1
Packing group II

Warning
H242-H315-H319-H335
P210-P280-P304+P340-P305+P351+P338-P403+P233-P410

2,6-Dichloroquinone-4-chlorimide > RPE - For analysis — RPE

Description ..... Yellow orange powder | Identification ..... Positive | Melting point..... 65.0 ÷ 67.0 ° C | Assay (argentimetric)..... ≥98 %

Code	Size	Packaging	Notes
442458	5 g	Glass bottle	

For spectrophotometric determination of vitamin B6

Diethanolamine

• Dietanolammina • Diéthanolamine • Dietanolamina • Diethanolamin

Synonym:
• *2,2'-Iminodiethanol*
• *Bis(2-hydroxyethyl)amine*

$(CH_2OHCH_2)_2NH$
Molecular Weight: 105,14
CAS: 111-42-2
EEC-N: 203-868-0



Danger
H302-H315-H318-H373
P260-P305+P351+P338-P310a-P330-P362+P364-P332+P313

Diethanolamine > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 30° C ..... 1.085 ÷ 1.091
Refract. index at 30° C ... 1.4723 ÷ 1.4783
Melting point ..... 27.8 ÷ 28.3 ° C
Water (K.F.) ..... ≤0.5 %
Monoethanolamine ..... ≤0.5 %
Residue on ignition ..... ≤50 ppm
Triethanolamine ..... ≤0.5 %
Assay (alkalimetric) ..... ≥99 %
Assay (GLC) ..... ≥99 %

Code	Size	Packaging	Notes
442554	100 g	Glass bottle	
442557	1 kg	Glass bottle	
442558	200 l	Metal drum	

Diethanolamine > ERBApharm - According to pharmacopoeia: USP-NF

ERBApharm

Description ..... Clear colourless liquid
Identification ..... Positive
Refract. index at 30° C ..... 1.473 ÷ 1.476
Water (K.F.) ..... ≤0,15 %
Triethanolamine ..... ≤1.0 %
Assay (alkalimetric) ..... 98,5÷101,0 %
Origin (BSE/TSE) ..... Synthesis
Residual solvents (Current ICH) ..... Conform

Code	Size	Packaging	Notes
337801	215 kg	Metal drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Diethylamine

• Dietilammina • Diéthylamine • Dietilamina • Diethylamin

$(C_2H_5)_2NH$
Molecular Weight: 73,14
CAS: 109-89-7
EEC-N: 203-716-3

Classification transport
ONU: 1154
Transport Hazard class: 3
Packing group II



Danger
H225-H302-H311-H332-H314-H335
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P361+P364-P403+P233

Diethylamine > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.705 ÷ 0.708
Refractive index at 20°C. 1.3840 ÷ 1.3900
Boiling point ..... 55 ÷ 56 ° C
Residue on evaporation ..... ≤50 ppm
Assay (GLC) ..... ≥99 %

Code	Size	Packaging	Notes
442756	1 l	Glass bottle	

Diethylamine > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Refractive index at 20°C. 1.3840 ÷ 1.3900
Assay (GLC) ..... ≥ 99 %
Boiling point ..... 55 ÷ 56 °C
Water (K.F.) ..... ≤ 0.2 %
Density at 20°C ..... 0.705 ÷ 0.708

Code	Size	Packaging	Notes
337501	1 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Diethyl carbonate

• Dietilcarbonato • Diéthylcarbonate • Dietilcarbonato • Diethylcarbonat

Synonym:
Carbonic acid diethyl ester

$(C_2H_5)_2CO_3$
Molecular Weight: 118,13
CAS: 105-58-8
EEC-N: 203-311-1

Classification transport
ONU: 2366
Transport Hazard class: 3
Packing group III

Warning
H226
P210-P241-P280-P303+P361+P353-P403+P235-P501a

Diethyl carbonate > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.971 ÷ 0.979
Refractive index at 20°C ..... 1.3820 ÷ 1.3870
Boiling point ..... 125.3 ÷ 126.8 ° C
Water (K.F.) ..... ≤500 ppm
Residue on evaporation ..... ≤20 ppm
Chloride ..... ≤5 ppm
Assay (GLC) ..... 99 ÷ 100 %

Code	Size	Packaging	Notes
443056	1 l	Glass bottle	

Diethylene glycol

• Glicol dietilenico • Glycol diéthylénique • Glicol dietilénico • Diethylenglykol

Synonym:
• *2,2'-Oxydiethanol*
• *2-Hydroxyethyl ether*

$(HOCH_2CH_2)_2O$
Molecular Weight: 106,12
CAS: 111-46-6
EEC-N: 203-872-2

Warning
H302
P264-P270-P301+P312a-P330-P501a

Diethylene glycol > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 1.115 ÷ 1.119
Refractive index at 20°C. 1.4450 ÷ 1.4500
Boiling point ..... 240 ÷ 252 ° C
Water (K.F.) ..... ≤0.3 %
Assay (GLC) ..... ≥98 %

Code	Size	Packaging	Notes
443255	1 l	Glass bottle	
443253	2.5 l	Glass bottle	
443252	25 kg	Drum	

Hygroscopic product. Store well sealed in a dry place

Diethylene glycol > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 1.112 ÷ 1.122
Refractive index at 20°C. 1.4425 ÷ 1.4525
Boiling point ..... 240 ÷ 252 ° C
Assay (GLC) ..... ≥99 %
Water ..... ≤0.3 %

Code	Size	Packaging	Notes
346301	1 l	Glass bottle	
346303	2.5 l	Glass bottle	
346304	30 kg	Metal drum	

Hygroscopic product. Store well sealed in a dry place

Diethylene glycol butyl ether ▶ 2-(2-Butoxyethoxy)ethanol

Diethylene glycol dimethyl ether

• Glicol dietilenico dimetiletere • Glycol diéthylénique diméthyléther • Glicol dietilénico dimetiléter
• Diethylenglycoldimethylether

Synonym:
• *Diglyme*
• *2-Methoxyethyl ester*

$(CH_3OCH_2CH_2)_2O$
Molecular Weight: 134,17
CAS: 111-96-6
EEC-N: 203-924-4

Classification transport
ONU: 3271
Transport Hazard class: 3
Packing group III

Danger
H226-H360FD-HEU019-HA26
P210-P241-P280-P303+P361+P353-P308+P313-P403+P235

Diethylene glycol dimethyl ether > RE - Pure

RE

Refractive index at 20°C ..... 1.406 - 1.41
Water content (K.F.) ..... ≤ 500 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 99.8 %

Code	Size	Packaging	Notes
P0410228	5 l	Metal tank	
P0410248	25 l	Metal drum	

Diethylene oxide ▶ 1,4-Dioxane

Diethylenetriaminepentacetic acid

• Acido dietilentriamminopentacetico • Acide diéthylènetriaminepentacétique
• Acido dietilentriaminopentaacético • Diethylentriaminpentaessigsäure

Synonym:
• *(Carboxymethylimino)bis(ethylenenitrilo) tetraacetic acid*
• *Pentetic acid*

$C_{14}H_{23}N_3O_{10}$
Molecular Weight: 393,35
CAS: 67-43-6
EEC-N: 200-652-8

Classification transport
ONU: -

Warning
H319-H412
P264-P273-P280i-P305+P351+P338-P337+P313-P501a

Diethylenetriaminepentacetic acid > RPE - For analysis

RPE

Description .....White crystalline powder Identification ..... Positive Assay (acidimetric) .....≥ 98.0 %

Code	Size	Packaging	Notes
405192	250 g	Glass bottle	

For the preparation of complexes: Na2Fe(DPTA), Na2[Cr(DPTA)] et H2[Gd(DPTA)]

Diethyl ether

• Dietiletere • Ether éthylique • Dietileter • Diethylether

Synonym:
Ethyl ether

$CH_3CH_2OCH_2CH_3$
Molecular Weight: 74,12
CAS: 60-29-7
EEC-N: 200-467-2

Classification transport
ONU: 1155
Transport Hazard class: 3
Packing group I

Danger
H224-H302-H336-HEU019-HEU066
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Diethyl ether > RS - For HPLC - Isocratic grade - Not stabilized

RS

Clear,colourless liq.appearance .....Conform
Identification .....Conform
Colour .....≤ 10 Apha
Refractive index at 20°C.....1.35 - 1.354
Water content (K.F.) .....≤ 100 mg/Kg
Free acid (as CH3COOH).....≤ 30 mg/Kg
Peroxides (as H2O2).....≤ 60 mg/Kg
Non volatile residue.....≤ 5 mg/Kg
UV transmittance at 220 nm ..... ≥ 10 %
UV transmittance at 230 nm ..... ≥ 45 %
UV transmittance at 250 nm ..... ≥ 75 %
UV transmittance at 280 nm ..... ≥ 92 %
UV transmittance at 300 nm ..... ≥ 98 %
Assay (GC) ..... ≥ 99.7 %

Code	Size	Packaging	Notes
412671	1 l	Glass bottle	
412672	2.5 l	Glass bottle	
412674	4 l	Glass bottle	

Diethyl ether > RS - PESTIPUR - For pesticide analysis - Not stabilized

RS

Description ..... Clear liquid
Colour .....≤ 10 hazen
Identification ..... Positive
Water .....≤ 0.02 %
Acidity (acetic acid).....≤ 30 ppm
Not volatile residue.....≤ 5 ppm
Peroxides (H2O2) .....≤ 60 ppm
GC-ECD (Lindano) .....≤ 3 ng/l
Assay (GLC) .....≥ 99.7 %

Code	Size	Packaging	Notes
447651	1 l	Glass bottle	
447652	2.5 l	Glass bottle	

Diethyl ether > RS - SPECTROSOL - For optical spectroscopy - Not stabilized

RS

Description .....Clear liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Colour ..... ≤ 10 Hazen
Density at 20° C .....0.712 ÷ 0.714
Peroxides (as H_2O_2).....≤ 60 mg/Kg
Boiling point..... 34.4 ÷ 34.9 ° C
Water (K.F.) .....≤100 ppm
Residue on evaporation .....≤5 ppm
Acidity .....≤0.0005 meq/g
Alcalinity.....≤0.0002 meq/g
Peroxides (H2O2) .....≤60 ppm
Assay (GLC) ..... ≥99.7 %
U.V. Transmittance
at 220 nm ..... ≥10 %
at 230 nm ..... ≥45 %
at 250 nm ..... ≥75 %
at 280 nm ..... ≥92 %
at 300 nm ..... ≥98 %

Code	Size	Packaging	Notes
447593	1 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Diethyl ether > RS - Anhydrous - For analysis - Stabilized with BHT

RS

Refractive index at 20°C ..... 1.35 - 1.354
Water content (K.F.) ..... ≤ 50 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 99.7 %
Peroxides (as H2O2) ..... ≤ 1 mg/Kg
Free acid (as CH3COOH) ..... ≤ 30 mg/Kg
Ketone and Aldehyde ..... ≤ 100 mg/kg
Methanol ..... ≤ 200 mg/Kg
Stabilizer (ionol) ..... 5 - 7 mg/Kg
Matter darkened by H2SO4 ..... ≤ 10 Hazen
Non volatile residue (without stab.) ≤ 10 mg/Kg

Code	Size	Packaging	Notes
P0441010	200 ml	Bottle with septum	
P04410T10	200 ml	Bottle with septum	On molecular sieves 4A, Water content < 20ppm
P0441008	1 l	Aluminium can	
P0441016	1 l	Glass bottle	
P04410T16	1 l	Glass bottle	On molecular sieves 4A, Water content < 20ppm
P0441021	2.5 l	Glass bottle	

Diethyl ether > RPE - For analysis - ACS - Not stabilized

RPE

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Water (K.F.) ..... ≤300 ppm
Residue on evaporation ..... ≤10 ppm
Acidity ..... ≤0.0002 meq/g
Ethyl alcohol ..... Conform
Carbonyl Compounds (CO) ..... ≤10 ppm
Peroxides (H2O2) ..... ≤1 ppm
Assay (GLC) ..... ≥99.0 %

Code	Size	Packaging	Notes
447534	1 l	Glass bottle	
447539	5 l	Aluminium can	
447532	20 kg	Aluminium can	
447531	140 kg	Metal drum	

Diethyl ether > RPE - For analysis - Stabilized with BHT

RPE

Refractive index at 20°C ..... 1.35 - 1.354
Water content (K.F.) ..... ≤ 200 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 99.7 %
Peroxides (as H2O2) ..... ≤ 1 mg/Kg
Free acid (as CH3COOH) ..... ≤ 30 mg/Kg
Ketone and Aldehyde ..... ≤ 100 mg/kg
Methanol ..... ≤ 200 mg/Kg
Stabilizer (ionol) ..... 5 - 7 mg/Kg
Matter darkened by H2SO4 ..... ≤ 10 Hazen
Non volatile residue (without stab.) ≤ 10 mg/Kg

Code	Size	Packaging	Notes
P0440508	1 l	Aluminium bottle	

Diethyl ether > RPE - For analysis - ACS - Stabilized with BHT

RPE

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification (I.R.) ..... Conform
Foreign odours ..... Conform
Ready carbonizable substances ..... Conform
Density at 20° C ..... 0.714 ÷ 0.716
Refractive index at 20°C ..... 1.35 - 1.354
Boiling point ..... 34.0 ÷ 35.0 ° C
Water (K.F.) ..... ≤200 ppm
Residue on evaporation ..... ≤10 ppm
Acetone ..... ≤50 ppm
Acidity (acetic acid) ..... ≤5 ppm
Alcalinity (NH3) ..... ≤1.4 ppm
Ethyl alcohol ..... ≤100 ppm
Methyl alcohol ..... ≤200 ppm
Carbonyl Compounds (CO) ..... ≤10 ppm
Heavy metals (Pb) ..... ≤1 ppm
Peroxides (H2O2) ..... ≤1 ppm
Al ..... ≤0.5 ppm
Ca ..... ≤0.5 ppm
Cd ..... ≤0.05 ppm
Co ..... ≤0.02 ppm
Cr ..... ≤0.02 ppm
Cu ..... ≤0.02 ppm
Fe ..... ≤0.1 ppm
Mg ..... ≤0.1 ppm
Mn ..... ≤0.02 ppm
Ni ..... ≤0.02 ppm
Pb ..... ≤0.05 ppm
Sn ..... ≤0.1 ppm
Zn ..... ≤0.1 ppm
Assay (GLC) ..... ≥99.8 %
Stabilized with obout 6 ppm BHT

Code	Size	Packaging	Notes
447521	1 l	Glass bottle	
447523	2.5 l	Glass bottle	
447522	5 l	Aluminium can	
447525	20 kg	Aluminium can	

Diethyl ether > ERBApharm - According to pharmacopoeia: BP-Ph.Eur. - Stabilized with BHT

ERBApharm

Description ..... Clear colourless liquid
Identification ..... Positive
Acidity ..... Conform Ph.Eur.
Aldehydes ..... Conform Ph.Eur.
Foreign odours ..... Conform Ph.Eur.
Peroxide ..... Conform Ph.Eur.
Density at 20° C ..... 0.714 ÷ 0.716
Boiling point ..... 34.0 ÷ 35.0 ° C
Water (K.F.) ..... ≤0.2 %
Non volat.substances ..... ≤20 ppm p/v
Assay (CPG) ..... ≥ 99.5 %
Origin (BSE/TSE) ..... Synthesis
Residual solvents (Current ICH) ..... Conform

Code	Size	Packaging	Notes
340751	1 l	Glass bottle	
340731	40 x 100 g	Glass bottle	
340752	20 kg	Aluminium can	
340759	140 kg	Metal drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Diethyl ether > RE - Pure - Stabilized with BHT

RE

Description ..... Clear colourless liquid
Density at 20°C ..... 0.710 ÷ 0.716
Identification ..... Positive
Density at 20° C ..... 0.710 ÷ 0.716
Boiling point..... 34.1 ÷ 35.1 ° C
Water (K.F.) ..... ≤300 ppm
Residue on evaporation ..... ≤20 ppm
Acidity (acetic acid)..... ≤20 ppm
Peroxyde (H2O2) ..... ≤ 1 ppm
Assay (GLC) ..... ≥99.5 %
Stabilized with BHT..... 5 ÷ 7 ppm

Code	Size	Packaging	Notes
528275	5 l	Aluminium can	
340762	20 kg	Aluminium can	
528276	25 l	Metal drum	
340765	140 kg	Metal drum	

n,n-Diethyl-p-phenylenediamine sulfate

• n,n-Dietil-p-fenilendiammina solfato • n,n-Diéthyl-p-phénylènediamine sulfate
• n,n-Dietil-p-fenilendiamonio sulfato • n,n-Diethyl-p-phenylendiamin-sulfat

Synonym:
4-Amino-N,N-diethylaniline sulfate salt

$NH_2C_6H_4N(C_2H_5)_2.H_2SO_4$
Molecular Weight: 262,33
CAS: 6283-63-2
EEC-N: 228-500-6

Classification transport
ONU: 2811
Transport Hazard class: 6.1
Packing group III

Danger
H301
P264-P270-P301+P310a-P330-P405-P501a

n,n-Diethyl-p-phenylenediamine sulfate > RPE - For analysis

RPE

DescriptionWhite to beige crystalline powder
Identification ..... Positive
Melting point..... 184 - 187 °C
Assay (oxidimetric) ..... ≥98.5 %

Code	Size	Packaging	Notes
443341	100 g	Plastic bottle	

n,n-Diethyl-p-phenylenediamine sulfate > RE - Pure

RE

Description ..... White to beige crystalline powder
Identification ..... Positive
Assay (oxidimetric) ..... ≥98 %

Code	Size	Packaging	Notes
338124	100 g	Plastic bottle	
338121	1 kg	Plastic bottle	

Diethyl phthalate

• Dietilftalato • Diéthyle phtalate • Dietil ftalato • Diethylphthalat

$C_6H_4(COOC_2H_5)_2$
Molecular Weight: 222,24
CAS: 84-66-2
EEC-N: 201-550-6

Diethyl phthalate > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 1.118 ÷ 1.121
Refractive index at 20° C.... 1.500 ÷ 1.505
Water (K.F.) ..... ≤0.1 %
Residue on ignition..... ≤100 ppm
Heavy metals (Pb)..... ≤5 ppm
Assay (GLC) ..... ≥99 %

Code	Size	Packaging	Notes
443404	1 l	Glass bottle	

Diethyl phthalate > ERBApharm - According to pharmacopoeia: BP-NF-Ph.Eur.-JPE

ERBApharm

Description .....Clear colourless,very slightly yellow liq.
Identification (IR)..... Positive
Refractive index at 20°C..... 1.500 ÷ 1.504
Appearance ..... Conform Ph.Eur.
Acidity (ml NaOH 0.1M) ..... ≤ 0.1 ml
Acid value..... ≤ 0.1 mg KOH/g
Related substances(GLC) ...Conform Ph.Eur.
Density (d20/20) ..... 1.118 - 1.121 Ph.Eur.
Density at 20°C ..... 1.118 - 1.122 NF
Water (K.F.) ..... ≤ 0.2 %
Sulphated ash..... ≤0.02 %
Assay (saponification)99.0 - 101.0 % (m/m)
Origin (BSE/TSE)..... Synthesis
Heavy metals (Pb)..... ≤ 20 ppm
As ..... ≤ 2 ppm
Residual solvents (Current ICH)..... Conform

Code	Size	Packaging	Notes
338112	1 l	Glass bottle	
338115	2.5 l	Glass bottle	
338113	30 kg	Metal drum	
338114	200 l	Metal drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Product specifications are subject to changes. Please visit our website for updates.

Differentiator for kit Gram-Hucker

• Differenziatore per kit Gram-Hucker • Différenciateur pour kit de Gram-Hucker • Diferenciador para kit Gram-Hucker
• Unterscheidungsmerkmal für Gram-Hucker Kit

Classification transport
ONU: 1987
Transport Hazard class: LQ



Danger
H225-H319-H371
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

Differentiator for kit Gram-Hucker > RS - For bacteriology

RS

Description ..... Liquido incolore
Identification ..... Positive

Code	Size	Packaging	Notes
444131	250 ml	Glass bottle	In Vitro Diagnostic Medical Device

Digitonin

• Digitonina • Digitonine • Digitonina • Digitonin

Synonym:
Digitin

$C_{56}H_{92}O_{29}$
Molecular Weight: 1229,34
CAS: 11024-24-1
EEC-N: 234-255-6

Classification transport
ONU: 1544
Transport Hazard class: 6.1
Packing group I



Danger
H301-H373
P260-P264-P301+P310a-P330-P314-P501a

Digitonin > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Water (K.F.) ..... ≤ 6 %
Residue on calcination ..... ≤ 0.3 %
Potere rotat. spec. (C=10; CH3COOH 75%) (s.s.) ..... -47 ÷ -49 °

Code	Size	Packaging	Notes
444207	1 g	Glass bottle	

Diglyme ▶ Diethylene glycol dimethyl ether

1,4-Dihydroxybenzene ▶ Hydroquinone

Diisopropylamine

• Diisopropilammina • Diisopropylamine • Diisopropilamina • Diisopropylamin

Synonym:
DIPA

$[(CH_3)_2CH]_2NH$
Molecular Weight: 101,19
CAS: 108-18-9
EEC-N: 203-558-5

Classification transport
ONU: 1158
Transport Hazard class: 3
Packing group II



Danger
H225-H302-H331-H314-H335
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Diisopropylamine > RPE - For analysis

RPE

Description ..... Clear liquid
Identification ..... Positive
Colour ..... ≤ 10 APHA
Refractive index at 20°C ..... 1.3910 ÷ 1.3930
Water ..... ≤ 0.2 %
Assay (GLC) ..... ≥ 99 %

Code	Size	Packaging	Notes
445981	250 ml	Glass bottle	

Diisopropylether

• Diisopropiletere • Ether isopropylique • Eter di-isopropílico • Di-isopropylether

Synonym:
Isopropyl ether

$[(CH_3)_2CH]_2O$
Molecular Weight: 102,18
CAS: 108-20-3
EEC-N: 203-560-6

Classification transport
ONU: 1159
Transport Hazard class: 3
Packing group II



Danger
H225-H336-HEU019-HEU066
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Diisopropylether > RS - Anhydrous - For analysis

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C ..... 1.366 - 1.370
Identification (IR) ..... Conform
Water content (K.F.) ..... ≤ 50 mg/Kg
Colour ..... ≤ 10 Hazen
Density d20/4 ..... 0.722 - 0.726
Peroxides (as H2O2) ..... ≤ 10 mg/Kg
Assay (GC) ..... ≥ 99.0 %
Stabilizer (ionol) ..... 2 - 15 mg/Kg
Free acid (as CH3COOH) ..... ≤ 10 mg/Kg
Non volatile residue (without stab.) ≤ 10 mg/Kg

Code	Size	Packaging	Notes
P0431016	1 l	Glass bottle	

Diisopropylether > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification (I.R.) ..... Conform
Alcohol miscibility ..... Complete
Chloroform miscibility ..... Complete
Diethyl ether miscib. ..... Complete
Density at 20° C ..... 0.719 ÷ 0.729
Boiling point ..... 66.5 ÷ 69.5 ° C
Water (K.F.) ..... ≤0.1 %
Peroxides (H2O2) ..... ≤5 ppm
Assay (GLC) ..... ≥98.5 %

Code	Size	Packaging	Notes
447932	1 l	Glass bottle	
447931	2.5 l	Glass bottle	
447933	5 l	Plastic tank	
447935	21 kg	Aluminium can	

Stabilized with ~10 ppm of BHT

Diisopropylether > RE - Pure

RE

Appearance ..... Clear colourless liquid
Identification (IR) ..... Conform
Water content (K.F.) ..... ≤ 1000 mg/Kg
Colour ..... ≤ 10 Hazen
Density d20/4 ..... 0.722 - 0.726
Peroxides (as H2O2) ..... ≤ 50 mg/Kg
Assay (GC) ..... ≥ 99.0 %
Stabilizer (ionol) ..... 2 - 15 mg/Kg
Free acid (as CH3COOH) ..... ≤ 10 mg/Kg
Non volatile residue (without stab.) ≤ 10 mg/Kg

Code	Size	Packaging	Notes
P0430228	5 l	Plastic tank	
P0430240	10 l	Metal tank	
P0430248	25 l	Metal drum	
P0430268	200 l	Metal drum	

N,N-Diisopropylethylamine

• N,N-Diisopropiletilammina • N,N-Diisopropyléthylamine • N,N-Diisopropiletilamina
• N,N-Diisopropylethylamin

Synonym:
• *DIPEA*
• *Ethyldiisopropylamine*

$C_8H_{19}N$
Molecular Weight: 129,25
CAS: 7087-68-5
EEC-N: 230-392-0

Classification transport
ONU: 3384
Transport Hazard class: 6.1
Packing group I



Danger
H225-H302-H331-H318-H335-H412
P210-P280-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P330-P403+P233

N,N-Diisopropylethylamine > RE - Pure

RE

Appearance ..... Clear, colourless to yellowish liquid
Identification ..... Conform
Refractive index at 20°C ..... 1.411 - 1.415
Water content (K.F.) ..... ≤ 3000 mg/Kg
Assay (GC) ..... ≥ 98 %

Code	Size	Packaging	Notes
P0400272	200 ml	Bottle with septum	
P0400252	30 l	Plastic tank	

Product specifications are subject to changes. Please visit our website for updates.

Dimedone

• Dimedone • Dimédon • Dimedona • Dimedon

Synonym:
Methone

$(CH_3)_2CCH_2COCH_2COCH_2$
Molecular Weight: 140,18
CAS: 126-81-8
EEC-N: 204-804-4

Dimedone > RPE - For analysis

RPE

Description ..... White crystalline powder | Identification ..... Positive | Melting point..... 146 ÷ 150 ° C | Assay (GLC) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
444252	25 g	Glass bottle	

For the determination of aldehydes

1,2-Dimethoxyethane

• 1,2-Dimetossietano • 1,2-Diméthoxyéthane • 1,2-Dimetoxietano • 1,2-Dimethoxyethan

Synonym:
• *Dimethylglycol*
• *Dimethyl ether ethylene glycol*

$C_4H_{10}O_2$
Molecular Weight: 90,12
CAS: 110-71-4
EEC-N: 203-794-9

Classification transport
ONU: 2252
Transport Hazard class: 3
Packing group II

Danger
H225-H332-H360FD-HEU019-HA26
P210-P241-P261-P280-P303+P361+P353-P304+P340

1,2-Dimethoxyethane > RS - Anhydrous - For analysis

RS

Refractive index at 20°C..... 1.377 - 1.381 | Water content (K.F.) ..... ≤ 200 mg/Kg | Colour ..... ≤ 10 Hazen | Assay (GC) ..... ≥ 99 %

Code	Size	Packaging	Notes
P0301010	200 ml	Bottle with septum	
P03010T16	1 l	Glass bottle	On molecular sieves 3A

1,2-Dimethoxyethane > RE - Pure

RE

Refractive index at 20°C..... 1.377 - 1.381
Water content (K.F.) ..... ≤ 500 mg/Kg
Non volatile residue ..... ≤ 50 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 99.5 %
Free acid (as CH3COOH)..... ≤ 150 mg/Kg

Code	Size	Packaging	Notes
P0300221	2.5 l	Glass bottle	
P0300268	200 l	Metal drum	

N,N-Dimethylacetamide

• N,N-Dimetilacetammide • N,N-Diméthylacétamide • N,N-Dimetilacetamida • N,N-Diethylacetamid

$CH_3CON(CH_3)_2$
Molecular Weight: 87,12
CAS: 127-19-5
EEC-N: 204-826-4

Danger
H312-H332-H360D-HA26
P261-P271-P280-P304+P340-P308+P313-P362+P364

N,N-Dimethylacetamide > RS - For Headspace chromatography

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.940 ÷ 0.946
Refractive index at 20°C. 1.4363 ÷ 1.4403
Boiling point..... 164.0 ÷ 166.0 °C
Water (K.F.) ..... ≤0.03 %

Acidity (acetic acid)..... ≤ 50 ppm
Assay (GLC) ..... ≥99.95 %
UV cut off..... ≤ 268 nm
U.V. Transmittance
at 268 nm ..... ≥ 10 %
at 275 nm ..... ≥ 55 %

At 300 nm ..... ≥ 85 %
at 350 nm ..... ≥ 98 %
at 400 nm ..... ≥ 99 %
GC/HS
Residual solvent of class 1(acc. to ICH) . ≤ 1 µg/g

Residual solvent of class 2(acc. to ICH)≤ 10 µg/g
Residual solvent of class 3(acc. to ICH)≤ 50 µg/g

Code	Size	Packaging	Notes
444311	1 l	Glass bottle	

N,N-Dimethylacetamide > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.940 ÷ 0.946
Refractive index at 20°C. 1.4343 ÷ 1.4403
Boiling point..... 164.0 ÷ 166.0 °C
Water (K.F.) ..... ≤0.05 %
Residue on evaporation ..... ≤50 ppm
Acidity (acetic acid)..... ≤150 ppm
Chloride..... ≤10 ppm
Heavy metals (Pb)..... ≤5 ppm
Sulphate..... ≤10 ppm
Fe ..... ≤5 ppm
Assay (GLC) ..... ≥99.8 %

Code	Size	Packaging	Notes
444307	1 l	Glass bottle	
444308	25 kg	Combined drum	
444309	190 kg	Combined drum	

p-Dimethylaminobenzaldehyde

• p-Dimetilamminobenzaldeide • p-Diméthylaminobenzaldéhyde • p-Dimetilaminobenzaldehído • p-Dimethylaminobenzaldehyd

Synonym: *Ehrlich's reagent*

$(CH_3)_2NC_6H_4CHO$
Molecular Weight: 149,19
CAS: 100-10-7
EEC-N: 202-819-0

Warning
H319-H317-H412
P261-P264-P280a-P305+P351+P338-P333+P313-P337+P313

p-Dimethylaminobenzaldehyde > RPE - For analysis

RPE

Description ..... Yellow or beige cristalline powder, chunk or flakes
Identification ..... Positive
Melting point..... 74 ± 2 ° C
Fee acidity ..... ≤ 1 %
Assay (non-aqueos medium)..... ≥ 97.5 %

Code	Size	Packaging	Notes
444604	100 g	Plastic bottle	
444603	250 g	Plastic bottle	

p-Dimethylaminobenzylidenerhodanine

• p-Dimetilamminobenzalrodanina • p-Diméthylaminobenzalrhodanine • p-Dimetilaminobencilidenrodanina • p-Diméthylaminobenzalrhodanine

$(CH_3)_2NC_6H_4CH:CSC:SNHCO$
Molecular Weight: 264,37
CAS: 536-17-4
EEC-N: 208-625-2

Warning
H302
P264-P270-P301+P312a-P330-P501a

p-Dimethylaminobenzylidenerhodanine > RPE - For analysis

RPE

Description ..... Red powder
Identification ..... Positive
Assay ($HClO_4$) ..... ≥ 98.5 %
Solubility (0,03% in acetone)... Clear orange solution

Code	Size	Packaging	Notes
444678	5 g	Glass bottle	

Acid-base indicator (pH 2,9÷4,0)

N,N-Dimethylformamide

• N,N-Dimetilformammide • N,N-Diméthylformamide • N,N-Dimetil formamida • N,N-Dimethylformamid

Synonym: *DMF*

$(CH_3)_2NOCH$
Molecular Weight: 73,09
CAS: 68-12-2
EEC-N: 200-679-5

Classification transport
ONU: 2265
Transport Hazard class: 3
Packing group III

Danger
H226-H312-H332-H319-H360D-HA26
P210-P241-P280-P303+P361+P353-P304+P340-P305+P351+P338

N,N-Dimethylformamide > RS - For HPLC - Isocratic Grade

RS

Appearance ..... Clear colourless liquid
Colour ..... ≤ 10 Apha
Water content (K.F.) ..... ≤ 0.03 % m/m
Non volatile residue ..... ≤ 0.0005 % m/m
Assay (GC) (on anhydrous) ..... ≥ 99.9 %
Free acid (as CH_3COOH)... ≤ 0.003 % m/m
UV transmittance at 270 nm ..... ≥ 30 %
UV transmittance at 275 nm ..... ≥ 60 %
UV transmittance at 300 nm ..... ≥ 90 %
UV transmittance at 320 nm ..... ≥ 97 %

Code	Size	Packaging	Notes
444981	1 l	Glass bottle	
444982	2.5 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

N,N-Dimethylformamide > RS - For Headspace chromatography

RS

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Density at 20° C ..... 0.945 ÷ 0.955
Refractive index at 20°C ..... 1.428 - 1.432
Boiling point ..... 152.0 ÷ 154.0 ° C
Water (K.F.) ..... ≤200 ppm
Residue on evaporation ..... ≤10 ppm
Assay (GLC) ..... ≥99.99 %
UV cut off ..... ≤ 269 nm
U.V. Transmittance
at 270 nm ..... ≥ 20 %
At 275 nm ..... ≥ 55 %
at 300 nm ..... ≥ 85 %
at 320 nm ..... ≥ 95 %
Residual solvent of class 2(acc. to ICH) ≤ 10 µg/g
Residual solvent of class 3(acc. to ICH) ≤ 50 µg/g
Residual solvent of class 1(acc. to ICH) . ≤ 1 µg/g

Code	Size	Packaging	Notes
444991	1 l	Glass bottle	

N,N-Dimethylformamide > RS - ATRASOL - For traces analysis

RS

Appearance ..... Clear liquid
Refractive index at 20°C ..... 1.428 - 1.432
Density d20/4 ..... 0.945 - 0.955
Colour ..... ≤ 10 Hazen
Water content (K.F.) ..... ≤ 200 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
Free acid (as HCOOH) ..... ≤ 20 mg/Kg
Free alkali as HN(CH3)2 ..... ≤ 10 mg/Kg
GC-FID.Individ. peak (hexane) ..... ≤ 3 mg/l
Retention time range before DMF

Code	Size	Packaging	Notes
P0343216	1 l	Glass bottle	
P0343221	2.5 l	Glass bottle	

N,N-Dimethylformamide > RS - PESTIPUR - For pesticide analysis

RS

Refractive index at 20°C ..... 1.428 - 1.432
Water content (K.F.) ..... ≤ 500 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 99.8 %
Non volatile residue ..... ≤ 10 mg/Kg
GC-ECD.Individual peak (Lindane) .. ≤ 3 ng/l
Retention time trichlorobenzene to mirex

Code	Size	Packaging	Notes
444941	1 l	Glass bottle	
444942	2.5 l	Glass bottle	

For chlorinated compounds analysis

N,N-Dimethylformamide > RS - SPECTROSOL - For optical spectroscopy

RS

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Density at 20° C ..... 0.945 ÷ 0.955
Refract. index at 25° C ... 1.4224 ÷ 1.4314
Boiling point ..... 152.0 ÷ 154.0 ° C
Water (K.F.) ..... ≤400 ppm
Residue on evaporation ..... ≤10 ppm
Acidity or alkalinity ..... ≤0.001 meq/g
Methyl alcohol ..... ≤100 ppm
Assay (GLC) ..... ≥99.9 %
U.V. Transmittance
at 270 nm ..... ≥28 %
at 280 nm ..... ≥72 %
at 300 nm ..... ≥90 %
at 320 nm ..... ≥96 %
at 330 nm ..... ≥98 %

Code	Size	Packaging	Notes
444957	1 l	Glass bottle	
444956	2.5 l	Glass bottle	

N,N-Dimethylformamide > RS - Anhydrous - For analysis

RS

Appearance ..... Clear liquid
Identification (IR) ..... Conform
Refractive index at 20°C ..... 1.428 - 1.432
Colour ..... ≤ 10 Hazen
Water content (K.F.) ..... ≤ 50 mg/Kg
Density d20/4 ..... 0.945 - 0.955
Non volatile residue ..... ≤ 20 mg/Kg
Free acid (as HCOOH) ..... ≤ 20 mg/Kg
Free alkali as HN(CH3)2 ..... ≤ 10 mg/Kg
Assay (GC) ..... ≥ 99.9 %
Methanol ..... ≤ 100 mg/Kg
Iron (Fe) ..... ≤ 5 mg/Kg

Code	Size	Packaging	Notes
P0341010	200 ml	Bottle with septum	
P03410T10	200 ml	Bottle with septum	On molecular sieves 4A
P0341016	1 l	Glass bottle	
P03410T16	1 l	Glass bottle	On molecular sieves 4A
P0341021	2.5 l	Glass bottle	

N,N-Dimethylformamide > RS - For peptide synthesis

RS

Appearance ..... Clear colourless liquid	Bromophenol blue test..... Conform	Non volatile residue ..... ≤ 15 mg/Kg	Iron (Fe) ..... ≤ 0.05 mg/Kg
Colour ..... ≤ 10 Hazen	Amines content..... ≤ 5 mg/Kg	Copper (Cu) ..... ≤ 0.05 mg/Kg	Nickel (Ni) ..... ≤ 0.05 mg/Kg
Refractive index at 20°C..... 1.428 - 1.432	Assay (GC)..... ≥ 99.9 %	Cadmium (Cd)..... ≤ 0.05 mg/Kg	Lead (Pb) ..... ≤ 0.1 mg/Kg
Water content (K.F.) ..... ≤ 300 mg/Kg	Methanol ..... ≤ 100 mg/Kg	Chromium (Cr) ..... ≤ 0.05 mg/Kg	

Code	Size	Packaging	Notes
P0343516	1 l	Glass bottle	
P0343521	2.5 l	Glass bottle	
P0343522	5 l	Plastic tank	
P0343541	10 l	Plastic tank	
P0343549	25 l	Plastic drum	
P0343550	25 l	Polythene-metal drum	
P0343567	200 l	Plastic drum	

N,N-Dimethylformamide > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Clear liquid	Refract. index at 25° C... 1.4224 ÷ 1.4314	Chloride..... ≤10 ppm	Fe ..... ≤0.05 ppm
Colour (APHA) ..... ≤15	Boiling point..... 152.0 ÷ 154.0 ° C	Heavy metals (Pb)..... ≤1 ppm	Ni..... ≤0.02 ppm
Identification (I.R.)..... Conform	Water (K.F.) ..... ≤0.03 %	Sulphate ..... ≤10 ppm	Pb ..... ≤0.1 ppm
Water miscibility..... Conform	Residue on evaporation ..... ≤20 ppm	Cd..... ≤0.05 ppm	Assay (GLC) ..... ≥99.9 %
Chloroform miscibility ..... Complete	Acidity (formic acid)..... ≤20 ppm	Cr..... ≤0.02 ppm	Methyl alcohol..... ≤ 100 ppm
Density at 20° C ..... 0.949 ÷ 0.952	Alcalinity (NH3)..... ≤20 ppm	Cu..... ≤0.02 ppm	

Code	Size	Packaging	Notes
444926	1 l	Glass bottle	
444923	2.5 l	Glass bottle	
444928	20 kg	Drum	
444925	190 kg	Metal drum	

N,N-Dimethylformamide > RE - Pure

RE

Description ..... Clear colourless liquid	Density at 20°C ..... 0.945 ÷ 0.950	Residue on evaporation ..... ≤ 50 ppm	Free acid (as HCOOH)..... ≤ 0.0020 %
Identification ..... Positive	Refractive index at 20°C. 1.4229 ÷ 1.4329	Water (K.F.) ..... ≤ 0.03 %	Free alkali as HN(CH3)2..... ≤ 0.0010 %
Colour ..... ≤ 10 APHA	Boiling point..... 152.0 ÷ 154.0 °C	Methyl alcohol..... ≤ 100 ppm	Assay (GLC) ..... ≥ 99.9 %

Code	Size	Packaging	Notes
508801	1 l	Glass bottle	
508804	1 l	Plastic bottle	
508802	2.5 l	Glass bottle	
528221	5 l	Plastic tank	
528220	25 l	Metal drum	
508803	200 l	Metal drum	

N,N-Dimethylformamide-d7

• N,N-Dimetilformammide-d7 • N,N-Diméthylformamide-d7 • N,N-Dimetilformamida-d7
• N,N-Dimethylformamid-d7

Synonym:
• *DMF-d7*
• *Heptadeutero-N,N-dimethylformamide*

$(CD_3)_2NOCD$
Molecular Weight: 80,14
CAS: 4472-41-7
EEC-N: 224-745-8

Classification transport
ONU: 2265
Transport Hazard class: 3
Packing group III



Danger
H226-H312-H332-H319-H360D-HA26
P210-P241-P280-P303+P361+P353-P304+P340-P305+P351+P338

N,N-Dimethylformamide-d7 > RS - For NMR - min 99.5%

RS

Code	Size	Packaging	Notes
P5189A	2 x 0.75 ml	Glass ampoule	

For specifications, contact our customer service for a certificate of analysis

Product specifications are subject to changes. Please visit our website for updates.

N,N-Dimethylformamide dimethylacetal

• N,N-Dimetilformammide dimetilacetale • N,N-Diméthylformamide diméthylacétal
• N,N-Dimetilformamida dimetilacetal • N,N-Dimethylformamiddimethylacetal

Synonym:
1,1-Dimethoxy-N,N-dimethylmethylamine

$(CH_3)_2NCH(OCH_3)_2$
Molecular Weight: 119,16
CAS: 4637-24-5
EEC-N: 225-063-3

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II

Danger
H225-H332-H315-H319
P210-P241-P280-P303+P361+P353-P304+P340-P305+P351+P338

N,N-Dimethylformamide dimethylacetal > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.890 - 0.910
Boiling point ..... 102 - 104 ° C
Assay (GC) ..... ≥ 97 %

Code	Size	Packaging	Notes
444901	10 ml	Glass bottle	

For derivatization

Dimethylglycol ▶ 1,2-Dimethoxyethane

Dimethylglyoxime disodium salt octahydrate ▶ Diacetyldioxime sodium salt

3,7-Dimethyl-2,6-octadienal ▶ Citral

N,N'-Dimethylpropylene uree

• N,N'-Dimetilpropilene urea • N,N'-Diméthylpropylène urée • N,N'-Dimetilpropileno urea
• N,N'-Dimethylpropylenharnstoff

Synonym:
DMPU

$C_6H_{12}N_2O$
Molecular Weight: 128,17
CAS: 7226-23-5
EEC-N: 230-625-6

Danger
H302-H318-H361f
P264-P280-P301+P312a-P305+P351+P338-P310a-P308+P313

N,N'-Dimethylpropylene uree > RE - Pure

RE

Clear,colourless to light yellow liq. ..Conform
Refractive index at 20°C..1.4883 - 1.4913
Water content (K.F.) ..... ≤ 1000 mg/Kg
Assay (GC) ..... ≥ 99.0 %

Code	Size	Packaging	Notes
P8020218	500 ml	Glass bottle	
P8020216	1 l	Glass bottle	
P8020229	5 l	Plastic tank	
P8020248	25 l	Metal drum	
P8020268	200 l	Metal drum	

Dimethylsulphoxide

• Dimetilsolfossido • Diméthylsulfoxyde • Dimetilsulfóxido • Dimethylsulfoxyd

Synonym:
• *Methyl sulfoxide*
• *DMSO*

CH_3SOCH_3
Molecular Weight: 78,13
CAS: 67-68-5
EEC-N: 200-664-3

Dimethylsulphoxide > RS - For HPLC - Isocratic Grade

RS

Appearance ..... Clear colourless liquid
Colour ..... ≤ 10 Apha
Water content (K.F.) ..... ≤ 0.04 % m/m
Non volatile residue ..... ≤ 0.0005 % m/m
Assay (GC) (on anhydrous) ..... ≥ 99.95 %
Free acid ..... ≤ 0.0005 meq/g
UV transmittance at 263 nm ..... ≥ 10 %
UV transmittance at 270 nm ..... ≥ 40 %
UV transmittance at 290 nm ..... ≥ 70 %
UV transmittance at 310 nm ..... ≥ 90 %
UV transmittance at 330 nm ..... ≥ 95 %
UV transmittance from 350 nm to 400 nm≥ 98 %

Code	Size	Packaging	Notes
445141	1 l	Glass bottle	
445142	2.5 l	Glass bottle	

Dimethylsulphoxide > RS - For Headspace chromatography

RS

Description ..... Clear liquid
Identification ..... Positive
Density at 20° C ..... 1.100 - 1.104
Water (K.F.) ..... ≤ 200 ppm
Refractive index at 20°C ..... 1.477 - 1.480
Residue on evaporation ..... ≤ 2 ppm
Assay (GLC) ..... ≥ 99.99 %
UV cut off ..... ≤ 265 nm
U.V. Transmittance
at 268 nm ..... ≥ 30 %
at 275 nm ..... ≥ 60 %
At 300 nm ..... ≥ 85 %
at 350 nm ..... ≥ 95 %
at 400 nm ..... ≥ 98 %
GC/HS
Residual solvent of class 1(acc. to ICH) . ≤ 1 µg/g
Residual solvent of class 2(acc. to ICH) ≤ 10 µg/g
Residual solvent of class 3(acc. to ICH) ≤ 50 µg/g

Code	Size	Packaging	Notes
445121	1 l	Glass bottle	

Dimethylsulphoxide > RS - ATRASOL - For traces analysis

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C ..... 1.477 - 1.481
Density d20/4 ..... 1.096 - 1.106
Water content (K.F.) ..... ≤ 200 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 99.98 %
Non volatile residue ..... ≤ 2 mg/Kg
GC-FID.Individ. peak (hexane) ..... ≤ 3 mg/l
Retention time range before DMSO

Code	Size	Packaging	Notes
P0353216	1 l	Glass bottle	
P0353221	2.5 l	Glass bottle	

Dimethylsulphoxide > RS - SPECTROSOL - For optical spectroscopy - Reag.Ph.Eur.

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 1.100 ÷ 1.104
Refractive index at 20°C ..... 1.478 ÷ 1.479
Melting point ..... ≥ 18.3 ° C
Water (K.F.) ..... ≤500 ppm
Residue on evaporation ..... ≤10 ppm
Acidity ..... ≤0.0005 meq/g
Alcalinity ..... ≤0.0002 meq/g
Assay (GLC) ..... ≥99.8 %
Fluorescence
at 254 nm ..... ≤2 ppb
at 365 nm ..... ≤2 ppb
Dimethylsulphone ..... ≤ 0.1 %
UV Absorbance at 262 nm ..... ≤ 1.00 AU
UV Absorbance at 270 nm ..... ≤ 0.46 AU
UV Absorbance at 290 nm ..... ≤ 0.16 AU
UV Absorbance from 340 nm ..... ≤ 0.01 AU
U.V. Transmittance
at 265 nm ..... ≥ 10 %
at 275 nm ..... ≥ 60 %
At 290 nm ..... ≥ 71 %
at 315 nm ..... ≥ 90 %
At 340 nm ..... ≥ 98 %

Code	Size	Packaging	Notes
445112	1 l	Glass bottle	
445111	2.5 l	Glass bottle	

Dimethylsulphoxide > RS - Anhydrous - For analysis

RS

Appearance ..... Clear liquid
Identification ..... Conform
Colour ..... ≤ 10 Hazen
Refractive index at 20°C ..... 1.477 - 1.481
Density d20/4 ..... 1.096 - 1.106
Water content (K.F.) ..... ≤ 50 mg/Kg
Assay (GC) ..... ≥ 99.8 %
Dimethylsulphone ..... ≤ 0.1 %

Code	Size	Packaging	Notes
508111	2.5 l	Glass bottle	

Dimethylsulphoxide > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Ready carbonizable substances ..... Conform
Density at 20° C ..... 1.096 ÷ 1.106
Refractive index at 20°C ..... 1.4765 ÷ 1.4825
Boiling point ..... 188.5 ÷ 189.5 °C
Melting point ..... 18.3 ÷ 18.7 °C
Water (K.F.) ..... ≤0.05 %
Residue on evaporation ..... ≤50 ppm
Acidity ..... ≤0.001 meq/g
Chloride ..... ≤5 ppm
Heavy metals (Pb) ..... ≤20 ppm
Sulphate ..... ≤50 ppm
Assay (GLC) ..... ≥99.9 %

Code	Size	Packaging	Notes
445103	1 l	Glass bottle	
445106	2.5 l	Glass bottle	
445107	15 kg	Plastic tank	
445101	25 kg	Plastic tank	

Dimethylsulphoxide > ERBApharm - According to pharmacopoeia: Ph.Eur.

ERBApharm

Appearance ..Colourless liquid or colourless crystals
Water solubility ..... Miscible
Ethanol solubility ..... Miscible
Refractive index at 20°C ..... 1.478 - 1.480
Identification C (IR) ..... Conform
Water content (K.F.) ..... ≤ 0.2 % m/m
Density d20/20 ..... 1.100 - 1.104
Free acid (ml NaOH 0.01N) ..... ≤ 5.0 ml
Solidification point ..... ≥ 18.3 °C
Assay (GC) ..... ≥ 99.9 %
Related substances (GC) ..... Conform
UV Absorbance at 275 nm ..... ≤ 0.3 AU
UV Absorbance at 285 nm ..... ≤ 0.2 AU
UV Absorbance at 295 nm ..... ≤ 0.2 AU
UV spectrum 270 - 350 nm ..... Smooth

Code	Size	Packaging	Notes
P0355016	1 l	Glass bottle	
P0355041	10 l	Plastic tank	
P0355049	25 l	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Product specifications are subject to changes. Please visit our website for updates.

Dimethylsulphoxide > RE - Pure

RE

Description ..... Clear colourless liquid	Assay (GC) ..... ≥ 99.8 %	Boiling point ..... 188 ÷ 190 °C	Dimethylsulphone ..... ≤ 0.1 %
Water content (K.F.) ..... ≤ 1000 mg/Kg	Density at 20°C ..... 1.096 ÷ 1.106	Acidity ..... ≤ 100 ppm	Assay (GLC) ..... ≥ 99.8 %
Identification ..... Positive	Refractive index at 20°C 1.4765 ÷ 1.4825	Water (K.F.) ..... ≤ 0.1 %	
Colour ..... ≤ 10 APHA	Non volatile residue ..... ≤ 50 mg/Kg	Residue on evaporation ..... ≤ 50 ppm	

Code	Size	Packaging	Notes
P03502T10	200 ml	Bottle with septum	On molecular sieves 4A
508001	1 l	Glass bottle	
P03502T16	1 l	Glass bottle	On molecular sieves 4A
508002	2.5 l	Glass bottle	
P03502T21	2.5 l	Glass bottle	On molecular sieves 4A
528335	5 l	Plastic tank	
P0350239	10 l	Plastic tank	
P0350242	20 l	Plastic tank	
P0350266	200 l	Combined drum	

Dimethylsulphoxide-d6

• Dimetilsolfossido-d6 • Diméthylsulfoxyde-d6 • Dimetilsulfóxido-d6 • Dimethylsulfoxid-d6

Synonym:
Hexadeuterodimethyl sulfoxide

CD_3SOCD_3
Molecular Weight: 84,18
CAS: 2206-27-1
EEC-N: 218-617-0

Dimethylsulphoxide-d6 > RS - For NMR - min 99.96%

RS

Code	Size	Packaging	Notes
P5220	10 x 0.6 ml	Glass ampoule	
P5229	10 x 0.75 ml	Glass ampoule	

Hygroscopic product. Store well sealed in a dry place. For specifications, contact our customer service for a certificate of analysi s

Dimethylsulphoxide-d6 > RS - For NMR - min 99.8%

RS

Code	Size	Packaging	Notes
P5200	10 x 0.6 ml	Glass ampoule	
P5209	10 x 0.75 ml	Glass ampoule	
P5204A	10 ml	Glass bottle	
P5204S	5 x 10 ml	Bottle with septum	
P5205	25 ml	Glass bottle	
P5206	100 ml	Glass bottle	

Hygroscopic product. Store well sealed in a dry place. For specifications, contact our customer service for a certificate of analysi s

Dimethylsulphoxide-d6 + 0.03% TMS

• Dimetilsolfossido-d6 + 0.03% TMS • Diméthylsulfoxyde-d6 + 0.03% TMS
• Dimetilsulfóxido-d6 + 0.03% TMS • Dimethylsulfoxid-d6 + 0.03% TMS

Synonym:
Hexadeuterodimethyl sulfoxide

CD_3SOCD_3
Molecular Weight: 84,18
CAS: 2206-27-1
EEC-N: 218-617-0

Dimethylsulphoxide-d6 + 0.03% TMS > RS - For NMR - min 99.95%

RS

Code	Size	Packaging	Notes
P5541	10 x 0.75 ml	Glass ampoule	
P5545	25 ml	Glass bottle	

Hygroscopic product. Store well sealed in a dry place. For specifications, contact our customer service for a certificate of analysi s

Dimethylsulphoxide-d6 + 0.03% TMS > RS - For NMR - min 99.8%

RS

Code	Size	Packaging	Notes
P5602	10 x 0.6 ml	Glass ampoule	
P5605	25 ml	Glass bottle	

Hygroscopic product. Store well sealed in a dry place. For specifications, contact our customer service for a certificate of analysi s

Dimidium bromide

• Dimidio bromuro • Dimidium bromure • Dimidio bromuro • Dimidiumbromid

Synonym:
- *3,8-Diamino-5-methyl-6-phenylphenanthridinium bromide*
- *Trypadine*

$C_{20}H_{18}BrN_3$
Molecular Weight: 380,29
CAS: 518-67-2
EEC-N: 208-256-7

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Dimidium bromide > RPE - For analysis

RPE

Description ..... Red-brown powder
Identification ..... Positive
Loss on drying ..... ≤ 5 %
Absorption Maximum ..... 523 - 528 nm
Absorption (A1%/1cm;1 max; methanol) ..... 168 - 177
Assay ..... ≥ 98 %

Code	Size	Packaging	Notes
445232	1 g	Glass bottle	
445231	5 g	Glass bottle	
445233	25 g	Glass bottle	

For the determination of surfactants

2,4-Dinitrochlorobenzene

• 2,4-Dinitroclorobenzene • 2,4-Dinitrochlorobenzène • 2,4-Dinitro-1-clorobenceno • 2,4-Dinitrochlorbenzol

Synonym:
- *DNCB*
- *1-Chloro-2,4-dinitrobenzene*

$(NO_2)_2C_6H_3Cl$
Molecular Weight: 202,55
CAS: 97-00-7
EEC-N: 202-551-4

Classification transport
ONU: 3441
Transport Hazard class: 6.1
Packing group II

Danger
H301-H311-H331-H373-H410
P271-P304+P340-P311a-P330-P361+P364-P403+P233

2,4-Dinitrochlorobenzene > RE - Pure

RE

Description ..... Yellow crystal. powder
Identification ..... Positive
Melting point ..... 48 ÷ 52 ° C
Residue on ignition ..... ≤0.1 %
Assay (argentimetric) ..... ≥98 %

Code	Size	Packaging	Notes
445421	250 g	Glass bottle	

2,4-Dinitrophenylhydrazine (with 30% of water)

• 2,4-Dinitrofenilidrazina (con 30% di acqua) • 2,4-Dinitrophénylhydrazine (avec 30% d'eau) • 2,4-Dinitrofenilhidracina (con 30% de agua)
• 2,4-Dinitrophenylhydrazin (mit 30% Wasser)

$(NO_2)_2C_6H_3NHNH_2$
Molecular Weight: 198,14
CAS: 119-26-6
EEC-N: 204-309-3

Classification transport
ONU: 1325
Transport Hazard class: 4.1
Packing group II

Danger
H228-H302-H315-H319-HEU001
P210-P241-P280-P305+P351+P338-P332+P313-P337+P313

2,4-Dinitrophenylhydrazine (with 30% of water) > RPE - For analysis

RPE

Description ..... Reddish powder
Identification ..... Positive
Melting point ..... 198 ÷ 201 °C
Water (K.F.) ..... 30 ÷ 35 %
Residue on ignition ..... ≤ 500 ppm
Assay (HPLC) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
445524	100 g	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

1,4-Dioxane

• 1,4-Diossano • 1,4-Dioxane • 1,4-Dioxano • 1,4-Dioxan

Synonym:
Diethylene oxide

$OCH_2CH_2OCH_2CH_2$
Molecular Weight: 88,11
CAS: 123-91-1
EEC-N: 204-661-8

Classification transport
ONU: 1165
Transport Hazard class: 3
Packing group II

Danger
H225-H319-H351-H335-HEU019-HEU066
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

1,4-Dioxane > RS - For HPLC - Isocratic grade - Stabilized with BHT

RS

Description ..... Clear liquid
Colour ..... ≤ 20 APHA
Identification (I.R.) ..... Positive
Water (K.F.) ..... ≤ 250 ppm
Residue on evaporation ..... ≤ 5 ppm
Acidity (acetic acid) ..... ≤ 50 ppm
Alcalinity (NH3) ..... ≤ 3 ppm
Assay (CPG) ..... ≥ 99.5 %
Absorbance
At 245 nm ..... ≤ 0.229 AU
Transmittance
At 270 nm ..... ≥ 80 %

Code	Size	Packaging	Notes
443231	1 l	Glass bottle	

1,4-Dioxane > RS - Anhydrous - For analysis - Stabilized with BHT

RS

Refractive index at 20°C ..... 1.42 - 1.424
Water content (K.F.) ..... ≤ 100 mg/Kg
Colour ..... ≤ 10 Hazen
Peroxides (as H2O2) ..... ≤ 50 mg/Kg
Assay (GC) ..... ≥ 99.8 %
Stabilizer (ionol) ..... 20 - 80 mg/Kg
Non volatile residue (without stab.) ≤ 10 mg/Kg
Free acid (as CH3COOH) ..... ≤ 50 mg/Kg

Code	Size	Packaging	Notes
P0361010	200 ml	Bottle with septum	
P0361016	1 l	Glass bottle	
P0361021	2.5 l	Glass bottle	

1,4-Dioxane > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611032001	50 ml	Glass bottle	Dioxane stock solution 1.0 mg/ml Ref Ph.Eur 1032001
611032003	50 ml	Glass bottle	Dioxane solution R1 0.1 mg/ml Ref Ph.Eur 1032003
611032002	100 ml	Glass bottle	Dioxane solution 0.5 mg/ml Ref Ph.Eur 1032002

1,4-Dioxane > RPE - For analysis - ACS - Reag.Ph.Eur. - Reag.USP - Stabilized with BHT

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 1.032 ÷ 1.036
Refractive index at 20°C 1.4194 ÷ 1.4254
Boiling point ..... 100.5 ÷ 101.5 °C
Freezing point ..... 11.5 ÷ 12.1 ° C
Water (K.F.) ..... ≤500 ppm
Residue on evaporation ..... ≤20 ppm
Acetal ..... ≤50 ppm
Acidity ..... ≤0.0016 meq/g
Total phosphorus ..... ≤0.1 ppm
Total silicon ..... ≤0.05 ppm
Total sulphur ..... ≤0.2 ppm
Ca ..... ≤0.5 ppm
Cu ..... ≤0.02 ppm
Fe ..... ≤0.2 ppm
K ..... ≤0.1 ppm
Mg ..... ≤0.05 ppm
Na ..... ≤0.5 ppm
Pb ..... ≤0.05 ppm
Zn ..... ≤0.2 ppm
Assay (GLC) ..... ≥99.8 %
Peroxides (H2O2) ..... ≤50 ppm
Carbonyl (as HCHO) ..... ≤100 ppm

Code	Size	Packaging	Notes
443202000	1 l	Glass bottle	
443206000	2.5 l	Glass bottle	
443204000	5 l	Plastic tank	
443201000	28 kg	Metal drum	
443205	200 kg	Metal drum	

1,4-Dioxane > RE - Pure - Stabilized with BHT

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Colour ..... ≤ 10 APHA
Density at 20° C ..... 1.031 ÷ 1.037
Refractive index at 20°C 1.4174 ÷ 1.4274
Boiling point ..... 100.3 ÷ 101.8 ° C
Water (K.F.) ..... ≤0.1 %
Residue on evaporation ..... ≤50 ppm
Acetal ..... ≤0.2 %
Acidity (acetic acid) ..... ≤ 50 ppm
Peroxides (H2O2) ..... ≤ 50 ppm
Assay (GLC) ..... ≥99.5 %
Stabilized with BHT ..... 20 ÷ 80 ppm

Code	Size	Packaging	Notes
338001	1 l	Glass bottle	
338003	2.5 l	Glass bottle	
338002	5 l	Plastic tank	
338005	28 kg	Metal drum	
338004	25 l	Plastic tank	

1,3-Dioxolane

• 1,3-Diossolano • 1,3-Dioxolane • 1,3-Dioxolano • 1,3-Dioxolan

Synonym:
- *Ethylene glycol methylene ether*
- *Formaldehyde ethylene acetal*

$C_3H_6O_2$
Molecular Weight: 74,08
CAS: 646-06-0
EEC-N: 211-463-5

Classification transport
ONU: 1166
Transport Hazard class: 3
Packing group II



Danger
H225-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

1,3-Dioxolane > RE - Pure

RE

Refractive index at 20°C ..... 1.3980 - 1.4020
Water content (K.F.) ..... ≤ 150 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 99.9 %
Peroxides (as H2O2) ..... ≤ 10 mg/Kg
Stabilizer (ionol) ..... ~75 mg/kg

Code	Size	Packaging	Notes
P8030216	1 l	Glass bottle	
P8030222	5 l	Plastic tank	
P8030249	25 l	Plastic tank	
P8030268	200 l	Metal drum	

Diphenylamine solution 1% in sulfuric acid

• Difenilammina soluzione 1% in acido solforico • Diphénylamine solution 1% dans l'acide sulfurique • Difenilamina solución 1% en ácido sulfúrico
• Diphenylaminlösung 1% in Schwefelsäure

$(C_6H_5)_2NH$
Molecular Weight: 169,23
CAS: 122-39-4

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group II



Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Diphenylamine solution 1% in sulfuric acid > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611032109	100 ml	Glass bottle	Ref Ph.Eur 1032101
611032101	1 l	Glass bottle	Ref Ph.Eur 1032101
611032102	1 l	Glass bottle	Diphenylamine solution R1 Ref Ph.Eur 1032102

Storage: protected from light

4-Diphenylaminesulfonic acid sodium salt

• 4-Difenilammina solfonato sodico • 4-Diphénylamine sulfonate sodique
• 4-Difenilamina sulfonato sal sódica • 4-Diphenylamin-Natriumsulfonat

Synonym:
- *4-(Phenylamino)benzene sulfonic acid sodium salt*
- *Sodium diphenylamine-4-sulfonate*

$C_6H_5NHC_6H_4SO_3Na$
Molecular Weight: 271,27
CAS: 6152-67-6
EEC-N: 228-165-6



Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-
P332+P313-P403+P233

4-Diphenylaminesulfonic acid sodium salt > RPE - For analysis - ACS

RPE

Description ..... Whitish powder
Identification ..... Positive
Sensitivity as indicat. ..... Conform

Code	Size	Packaging	Notes
443671	10 g	Glass bottle	

Redox indicator. Purple / Red - Clear

Product specifications are subject to changes. Please visit our website for updates.

sym-Diphenylcarbazide

• sim-Difenilcarbazide • Sym-Diphénylcarbazide • sim-Difenilcarbacida • sym-Diphenylcarbazid

Synonym:
• *1,5-Diphenylcarbazide*
• *1,5-Diphenylcarbohydrazide*

$C_6H_5NHNHCONHNHC_6H_5$
Molecular Weight: 242,28
CAS: 140-22-7
EEC-N: 205-403-7

sym-Diphenylcarbazide > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White powder
Identification ..... Positive
Melting point ..... 173 ÷ 176 ° C
Loss on drying ..... ≤ 1 %
Chromate sensitivity ..... Conform
Sulphated ash ..... ≤ 0.05 %
Assay (HPLC) ..... ≥ 97.5 %

Code	Size	Packaging	Notes
443752	25 g	Glass bottle	
443754	100 g	Glass bottle	

Redox indicator

sym-Diphenylcarbazone

• sim-Difenilcarbazone • Sym-Diphénylcarbazone • sim-Difenilcarbazona • sym-Diphenylcarbazon

Synonym:
Phenylazoformic acid 2-phenylhydrazide

$C_6H_5N{:}NCONHNHC_6H_5$
Molecular Weight: 240,27
CAS: 538-62-5
EEC-N: 208-698-0

sym-Diphenylcarbazone > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Orange crystalline powder
Identification ..... Positive
Acetone solubility ..... Conform
Mercury sensitivity ..... Conform
Residue on ignition ..... ≤0.1 %

Code	Size	Packaging	Notes
443801	10 g	Glass bottle	

Contains sym-Diphenylcarbazide. Redox indicator

Diphenylthiocarbazone

• Difeniltiocarbazone • Diphénylthiocarbazone • Difeniltiocarbazona • Diphenylthiocarbazon

Synonym:
Dithizone

$C_6H_5NHNHCSN{:}NC_6H_5$
Molecular Weight: 256,32
CAS: 60-10-6
EEC-N: 200-454-1

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Diphenylthiocarbazone > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Black-violet powder
Identification ..... Positive
E 620/E 450 nm ..... ≥ 1.55
Heavy metals (Pb) ..... ≤ 20 ppm
Residue on ignition ..... ≤ 0.3 %
Assay (spectrophotom.) ..... ≥ 85.0 %

Code	Size	Packaging	Notes
444053	50 g	Glass bottle	

2,2'-Dipyridyl

• 2,2'-Dipiridile • 2,2'-Dipyridyle • 2,2'-Dipiridilo • 2,2'-Dipyridyl

Synonym:
2,2'-Bipyridine

$C_{10}H_8N_2$
Molecular Weight: 156,19
CAS: 366-18-7
EEC-N: 206-674-4

Classification transport
ONU: 2811
Transport Hazard class: 6.1
Packing group III

Danger
H301-H311
P264-P280h-P301+P310a-P330-P361+P364-P501a

2,2'-Dipyridyl > RPE - For analysis

RPE

Description ..... Polvere crist. quasi bianca
Identification ..... Positive
Melting point ..... 69 ÷ 72 ° C
Iron sensitivity ..... Conform
Assay (GLC) ..... ≥ 99.0 %

Code	Size	Packaging	Notes
445958	5 g	Glass bottle	

Reactive iron and molybdenum

2,6 - Di-tert-butyl-p-cresol ▶ Butylhydroxytoluene

Direct red 28 ▶ Congo red

Disodium hydrogen phosphate ▶ Sodium phosphate dibasic anhydrous

Disodium hydrogen phosphate dodecahydrate ▶ Sodium phosphate dibasic dodecahydrate

Dithiooxamide

• Ditioossammide • Dithiooxamide • Ditiooxamida • Dithiooxamids

$NH_2CSCSNH_2$
Molecular Weight: 120,2
CAS: 79-40-3
EEC-N: 201-203-9

Warning
H302-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Dithiooxamide > RPE - For analysis

RPE

Description ..... Red orange crystals
Identification ..... Positive
Residue on ignition ..... ≤500 ppm
Copper sensitivity ..... ≤ 0.05 µg/ml
Assay (ex nitrogen) ..... ≥98 %

Code	Size	Packaging	Notes
446008	5 g	Plastic bottle	

Dithizone ▶ Diphenylthiocarbazone

DMF ▶ N,N-Dimethylformamide

DMSO ▶ Dimethylsulphoxide

1-Dodecanesulfonic acid sodium salt

• Acido 1-dodecansolfonico sale sodico • Acide 1-dodecanesulfonique sel sodique
• Acido 1-dodecanosulfónico sal sódica • 1-Dodecansulfonsäure-Natriumsalz

Synonym:
Sodium 1-dodecanesulfonate

$CH_3(CH_2)_{11}SO_3Na$
Molecular Weight: 272,39
CAS: 2386-53-0

1-Dodecanesulfonic acid sodium salt > RS - For ion pair chromatography

RS

Description ..... White crystalline powder
Water (K.F.) ..... ≤ 1.0 %
Assay ..... ≥ 99.0 %
Absorbance (5% in water)
At 210 nm ..... ≤ 0.05 AU
At 220 nm ..... ≤ 0.03 AU
At 230 nm ..... ≤ 0.02 AU

Code	Size	Packaging	Notes
405881	25 g	Glass bottle	
405882	100 g	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Dodecylbenzenesulphonic acid sodium salt

• Acido dodecilbenzensolfonico sale sodico • Acide dodécylbenzènesulfonique sel sodique
• Acido dodecilbencensulfónico sal sódica • Dodecylbenzolsulfonsäure-Natrium-Salz

Synonym:
Sodium dodecylbenzenesulfonate

$C_{12}H_{25}C_6H_4.SO_3Na$
Molecular Weight: 348,49
CAS: 25155-30-0
EEC-N: 246-680-4

Warning
H302
P264-P270-P301+P312a-P330-P501a

Dodecylbenzenesulphonic acid sodium salt > RS - For surfactants detection

RS

Description ..... White-yellowish crystalline powder | Identification ..... Positive | Water ..... < 5.0 % | Assay ..... > 83 %

Code	Size	Packaging	Notes
405351	10 g	Glass bottle	
405352	25 g	Glass bottle	

Minimum 90% biodegradability

Dodecyltrimethylammonium bromide

• Dodeciltrimetilammonio bromuro • Dodécyltriméthylammonium bromure • Dodeciltrimetilamonio bromuro
• Dodecyltrimethylammoniumbromid

Synonym:
Lauryltrimethylammonium bromide

$CH_3(CH_2)_{11}N(CH_3)_3Br$
Molecular Weight: 308,34
CAS: 1119-94-4
EEC-N: 214-290-3

Classification transport
ONU: 3077
Transport Hazard class: 9
Packing group III

Warning
H302-H315-H319-H335-H410
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Dodecyltrimethylammonium bromide > RS - For ion pair chromatography

RS

Absorbance UV curve (10%)
A240nm (1M) ..... ≤ 0.2 AU | A250nm (1M) ..... ≤ 0.03 AU | A260nm (1M) ..... ≤ 0.02 AU | A500nm (1M) ..... ≤ 0.02 AU

Code	Size	Packaging	Notes
405941	25 g	Plastic bottle	

Dysprosium standard solution

• Disprosio standard soluzione • Dysprosium solution standard • Disprosio, solución patrón • Dysprosium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Dysprosium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505582	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505585	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Dysprosium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
504231	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504233	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504235	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
504237	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Dysprosium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507734	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507500	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Product specifications are subject to changes. Please visit our website for updates.

EDTA ▶ Ethylenediaminetetraacetic acid

EDTA disodium salt ▶ Ethylenediaminetetraacetic acid disodium salt

EDTA dipotassium salt dihydrate ▶ Ethylenediaminetetraacetic acid dipotassium salt dihydrate

EDTA potassium and magnesium salt dihydrate ▶ Ethylenediaminetetraacetic acid potassium and magnesium salt dihydrate

EDTA tetrasodium salt tetrahydrate ▶ Ethylenediaminetetraacetic acid tetrasodium salt tetrahydrate

EGTA ▶ Ethylene glycol bis(2-aminoethyl ether)-N,N,N',N'-tetraacetic acid

Ehrlich's reagent

• Ehrlich reattivo • Réactif d'Ehrlich • Ehrlich reactivo • Ehrlich Reagenz

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Danger
H290-H314-H317
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P362+P364

Ehrlich's reagent > RS - For microscopy

RS

Description ..... Yellow clear liquid Identification ..... Positive

Code	Size	Packaging	Notes
E446302	500 ml	Glass bottle	

For urobilinogen detection

Eluent sodium bicarbonate

• Eluente sodio bicarbonato • Eluant sodium dicarbonate • Eluyente sodio dicarbonato • Elutionsmittel Natriumbicarbonat

$NaHCO_3$
Molecular Weight: 84,01
CAS: 144-55-8

Eluent sodium bicarbonate > RS - Eluent concentrates for ion chromatography

RS

Code	Size	Packaging	Notes
504534	100 ml	Plastic bottle	0.5 M Sodium bicarbonate
507578	1 l	Plastic bottle	0.5 M Sodium bicarbonate

Eluent sodium carbonate

• Eluente sodio carbonato • Eluant sodium carbonate • Eluyente sodio carbonato • Elutionsmittel Natriumcarbonat

Synonym:
• *Calcined soda*
• *Carbonic acid disodium salt*

Na_2CO_3
Molecular Weight: 84.01
CAS: 497-19-8

HEU210

Eluent sodium carbonate > RS - Eluent concentrates for ion chromatography

RS

Code	Size	Packaging	Notes
504533	100 ml	Plastic bottle	0.5 M Sodium carbonate
507577	1 l	Plastic bottle	0.5 M Sodium carbonate

Eluent sodium carbonate/sodium bicarbonate

• Eluente sodio carbonato/sodio bicarbonato • Eluant sodium carbonate/sodium dicarbonate • Eluyente sodio carbonato/sodio dicarbonato
• Elutionsmittel Natriumcarbonat / Natriumbicarbonat

HEU210

Eluent sodium carbonate/sodium bicarbonate > RS - Eluent concentrates for ion chromatography

RS

Code	Size	Packaging	Notes
504530	100 ml	Plastic bottle	0.18 M Sodium carbonate / 0.17 M Sodium bicarbonate
504531	100 ml	Plastic bottle	0.22 M Sodium carbonate / 0.28 M Sodium bicarbonate
504532	100 ml	Plastic bottle	0.35 M Sodium carbonate / 0.1 M Sodium bicarbonate

Eosin B

• Eosina B • Eosine B • Eosina B • Eosin B

Synonym:
• *Acid Red 91*
• *4',5'-Dibromo-2',7'-dinitrofluorescein, disodium salt*

$C_{20}H_6Br_2N_2Na_2O_9$
Molecular Weight: 624,09
CAS: 548-24-3
EEC-N: 208-943-1

Eosin B > RPE - For analysis - C.I. 45400

RPE

Description ..... Brown greyish powder Identification ..... Positive Absorbtion ind.sensit. ..... Conform

Code	Size	Packaging	Notes
446602	25 g	Glass bottle	

Dye for microscopy (histology). Absorbance and fluorescence indicator

Eosin Y

• Eosina Y • Eosine Y • Eosina Y • Eosin Y

Synonym:
• *Acid Red 87*
• *2',4',5',7'-Tetrabromofluorescein disodium salt*

$C_{20}H_6Br_4Na_2O_5$
Molecular Weight: 691,86
CAS: 17372-87-1
EEC-N: 241-409-6

Warning
H312-H332
P261-P271-P280h-P304+P340-P312a-P501a

Eosin Y > RS - For microscopy - C.I. 45380

RS

Description ..... Red-brown powder Loss on drying ..... ≤ 10 % Assay ..... 80.00 ÷ 88.00 %
Identification ..... Positive Sens. as absorption indicator ..... Conform

Code	Size	Packaging	Notes
446632	25 g	Glass bottle	
446634	100 g	Plastic bottle	

Dye for histology

Eosin Y 1% solution acqueous

• Eosina Y 1% soluzione acquosa • Eosine Y 1% solution aqueuse • Eosina Y 1% solución acuosa
• Eosin Y 1% ige wässrige Lösung

Synonym:
• *Acid Red 87*
• *2',4',5',7'-Tetrabromofluorescein disodium salt*

$C_{20}H_6Br_4Na_2O_5$
Molecular Weight: 691,86
CAS: 17372-87-1

HEU210

Eosin Y 1% solution acqueous > RS - For histology

RS

Description ..... Brown liquid Identification ..... Positive Maximum absorption ..... 515÷518 nm A 1%/1cm (0.005 g/l) ..... 1200÷1400 nm

Code	Size	Packaging	Notes
446644	1 l	Plastic bottle	In Vitro Diagnostic Medical Device

Product specifications are subject to changes. Please visit our website for updates.

Eosin Y 0.5% solution alcoholic

• Eosina Y 0.5% soluzione alcolica • Eosine Y 0.5% solution alcoolique • Eosina Y 0.5% solución alcohólica • Eosin Y 0.5% alkoholische Lösung

Synonym:
• *Acid Red 87*
• *2',4',5',7'-Tetrabromofluorescein disodium salt*

$C_{20}H_6Br_4Na_2O_5$
Molecular Weight: 691,86
CAS: 17372-87-1

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II



Danger
H225-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

Eosin Y 0.5% solution alcoholic > RS - For microscopy

RS

Code	Size	Packaging	Notes
446664	1 l	Plastic bottle	In Vitro Diagnostic Medical Device

Erbium standard solution

• Erbio standard soluzione • Erbium solution standard • Erbio, solución patrón • Erbium-Standardlösung

Classification transport
ONU: 2837
Transport Hazard class: 8
Packing group III



Erbium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505592	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505595	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Erbium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
504241	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504243	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504245	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
504247	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Erbium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507735	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507501	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Eriochrome black T

• Nero eriocromo T • Noir ériochrome T • Negro de eriocromo T • Schwarzes Eriochrom T

Synonym:
Mordant Black 11

$C_{20}H_{12}N_3NaO_7S$
Molecular Weight: 461,39
CAS: 1787-61-7

HEU210

Eriochrome black T > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611056801	100 g	Plastic bottle	Ref Ph.Eur 1056801

Storage: protected from light

Eriochrome black T > RPE - For analysis - C.I. 14645

RPE

Description ..... Blackish brown powder Identification ..... Positive Loss on drying ..... ≤ 10 % Sensitivity ..... Conform

Code	Size	Packaging	Notes
464221	10 g	Glass bottle	
464222	25 g	Glass bottle	

Complexometric indicator

Eriochromocyanine R

• Eriocromo cianina R • Eriochrome cyanine R • Eriocromocianina R • Eriochromcyanin R

Synonym:
Chromoxane cyanine R

$C_{23}H_{15}Na_3O_9S$
Molecular Weight: 536,4
CAS: 3564-18-9
EEC-N: 222-641-7

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Eriochromocyanine R > RPE - For analysis - C.I. 43820

RPE

Description ..... Red brick powder Identification ..... Positive Alluminium sensitivity ..... ≥1 µg/ml

Code	Size	Packaging	Notes
446811	10 g	Glass bottle	
446812	25 g	Glass bottle	

For the determination of Al. Complexometric indicator

Erythrosin extra B

• Eritrosina extra B • Erythrosine extra B • Eritrosina extra B • Erythrosin extra B

Synonym:
• *Acid Red 51*
• *2',4',5',7'-Tetraiodofluorescein disodium salt*

$C_{20}H_6I_4O_5Na_2$
Molecular Weight: 879,87
CAS: 16423-68-0
EEC-N: 240-474-8

Warning
H302
P264-P270-P301+P312a-P330-P501a

Erythrosin extra B > RS - For microscopy - C.I. 45430

RS

Description ..... Red brown powder Identification ..... Positive Assorbanza 524-527 nm ..... Conform Assorbività specifica 1%/1cm 930-1170

Code	Size	Packaging	Notes
446972	25 g	Glass bottle	
446971	100 g	Plastic bottle	

Dye for histology

Esbach's reagent

• Esbach reattivo • Réactif d'Esbach • Esbach reactivo • Esbach-Reagenz

HEU210

Esbach's reagent > RS - For microscopy

RS

Description ..... Yellow clear liquid Identification ..... Positive

Code	Size	Packaging	Notes
446981	1 l	Plastic bottle	

For the determination of albumin

1,2-Ethanediol ▶ Ethylene glycol

Product specifications are subject to changes. Please visit our website for updates.

Ethanol absolute anhydrous

• Etanolo assoluto anidro • Ethanol absolu anhydre • Etanol absoluto anhidro • Ethanol absolut wasserfrei

Synonym:
Ethyl alcohol absolute anhydrous

C_2H_5OH
Molecular Weight: 46,07
CAS: 64-17-5
EEC-N: 200-578-6

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II

Danger
H225-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

Ethanol absolute anhydrous > RS - For HPLC PLUS Gradient grade

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.7893 ÷ 0.7899
Refractive index at 20°C. 1.3602 ÷ 1.3622
Boiling point..... 78.3 ÷ 78.8 ° C
Acidity or alkalinity..... ≤0.0002 meq/g
Water (K.F.) ..... ≤0.05 %
Residue on evaporation ..... ≤5 ppm
Assay (GLC) ..... ≥99.9 %
Fluorescence
at 254 nm ..... ≤2 ppb
at 365 nm ..... ≤2 ppb
U.V. Transmittance
at 210 nm ..... ≥25 %
at 220 nm ..... ≥50 %
at 230 nm ..... ≥75 %
At 240 nm ..... ≥ 85 %
at 250 nm ..... ≥90 %
At 260 nm ..... ≥ 98 %
≥ 270 nm ..... ≥ 99 %
Acetal + acetaldehyde ..... ≤ 10 ppm(v/v)
HPLC Gradient
At 235 nm ..... ≤ 5 mAU

Code	Size	Packaging	Notes
412704	500 ml	Glass bottle	
412701	1 l	Glass bottle	Only for Italian market
412703	1 l	Glass bottle PVC coated	Only for Italian market
4127012	1 l	Glass bottle	
4127032	1 l	Glass bottle PVC coated	
412702	2.5 l	Glass bottle	Only for Italian market
4127022	2.5 l	Glass bottle	

Ethanol absolute anhydrous > RS - For HPLC - Isocratic Grade

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.7893 ÷ 0.7899
Refractive index at 20°C. 1.3602 ÷ 1.3622
Boiling point..... 78.3 ÷ 78.8 ° C
Acidity or alkalinity..... ≤0.0002 meq/g
Water (K.F.) ..... ≤0.05 %
Residue on evaporation ..... ≤5 ppm
Assay (GLC) ..... ≥99.9 %
U.V. Transmittance
at 210 nm ..... ≥25 %
at 220 nm ..... ≥50 %
at 230 nm ..... ≥75 %
At 240 nm ..... ≥ 80 %
at 250 nm ..... ≥90 %
at 270 nm ..... ≥94 %
at 290 nm ..... ≥97 %

Code	Size	Packaging	Notes
412521	1 l	Glass bottle	Only for Italian market
4125212	1 l	Glass bottle	
412522	2.5 l	Glass bottle	Only for Italian market
4125222	2.5 l	Glass bottle	

Ethanol absolute anhydrous > RS - SPECTROSOL - For optical spectroscopy

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.7893 ÷ 0.7899
Boiling point..... 78.3 ÷ 78.8 ° C
Acidity or alkalinity..... ≤0.0002 meq/g
Water (K.F.) ..... ≤0.2 %
Residue on evaporation ..... ≤10 ppm
Assay (GLC) ..... ≥99.8 %
Fluorescence
at 254 nm ..... ≤2 ppb
at 365 nm ..... ≤2 ppb
U.V. Transmittance
at 210 nm ..... ≥25 %
at 220 nm ..... ≥50 %
at 230 nm ..... ≥75 %
at 250 nm ..... ≥89 %
at 270 nm ..... ≥94 %
at 290 nm ..... ≥97 %

Code	Size	Packaging	Notes
414677	1 l	Glass bottle	Only for Italian market
4146772	1 l	Glass bottle	

Ethanol absolute anhydrous > RS - Anhydrous - For analysis

RS

Refractive index at 20°C ..... 1.358 - 1.362
Water content (K.F.) ..... ≤ 200 mg/Kg
Colour ..... ≤ 10 Hazen
Alcohol content (20°C) ..... ≥ 99.9 % V/V
Non volatile residue ..... ≤ 10 mg/Kg
Assay (GC) ..... ≥ 99.8 %
Free acid (as CH3COOH) ..... ≤ 10 mg/Kg
Aldehydes (as acetaldehyde) ..... ≤ 3 mg/Kg
Esters (as CH3COOC2H5) ..... ≤ 25 mg/Kg

Code	Size	Packaging	Notes
P013A1010	200 ml	Bottle with septum	
P013A1016	1 l	Glass bottle	
P013A1021	2.5 l	Glass bottle	

Ethanol absolute anhydrous > RS - VLSI - For electronic use

RS

Code	Size	Packaging	Notes
527681	1 l	Plastic bottle	
527680	2.5 l	Plastic bottle	

For specifications, contact our customer service for a certificate of analysis

Ethanol absolute anhydrous > RS - RSE - For electronic use

RS

Description ..... Clear colourless liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Water miscibility ..... Conform
Ready carbonizable substances ..... Conform
Density at 20° C ..... 0.790 ÷ 0.793
Boiling point ..... 78.3 ÷ 78.8 ° C
Resistivity ..... ≥0.5 Mohm cm
Assay(alcohol.) at 20°C ..... ≥99.9 %
Water (K.F.) ..... ≤0.1 %
Residue on evaporation ..... ≤10 ppm
Acidity (acetic acid) ..... ≤10 ppm
Alcalinity (NH3) ..... ≤1 ppm
Methyl alcohol ..... ≤100 ppm

Chloride ..... ≤0.2 ppm
Carbonyl Compounds (CO) ..... ≤10 ppm
Total phosphorus ..... ≤0.1 ppm
Heavy metals (Pb) ..... ≤0.2 ppm
Subst. reducing KMnO4 ..... ≤2.5 ppm
Total sulphur ..... ≤1 ppm
Ag ..... ≤0.02 ppm
Al ..... ≤0.05 ppm
As ..... ≤0.01 ppm
Au ..... ≤0.05 ppm
B ..... ≤0.01 ppm
Ba ..... ≤0.1 ppm
Be ..... ≤0.02 ppm
Bi ..... ≤0.02 ppm

Ca ..... ≤0.5 ppm
Cd ..... ≤0.01 ppm
Co ..... ≤0.01 ppm
Cr ..... ≤0.02 ppm
Cu ..... ≤0.02 ppm
Fe ..... ≤0.1 ppm
Ga ..... ≤0.02 ppm
In ..... ≤0.02 ppm
K ..... ≤0.1 ppm
Li ..... ≤0.02 ppm
Mg ..... ≤0.1 ppm
Mn ..... ≤0.01 ppm
Mo ..... ≤0.05 ppm
Na ..... ≤0.5 ppm

Ni ..... ≤0.01 ppm
Pb ..... ≤0.02 ppm
Pt ..... ≤0.05 ppm
Sb ..... ≤0.01 ppm
Sn ..... ≤0.02 ppm
Sr ..... ≤0.02 ppm
Ti ..... ≤0.05 ppm
Tl ..... ≤0.05 ppm
V ..... ≤0.05 ppm
Zn ..... ≤0.03 ppm
Zr ..... ≤0.05 ppm

Code	Size	Packaging	Notes
414587	1 l	Glass bottle	Only for Italian market
4145872	1 l	Glass bottle	
414583	2.5 l	Glass bottle	Only for Italian market
4145832	2.5 l	Glass bottle	

Ethanol absolute anhydrous > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Identification (I.R.) ..... Positive
Water miscibility ..... Complete
Density at 20°C ..... 0.7893 ÷ 0.7899
Boiling point ..... 78.3 ÷ 78.8 °C
Refractive index at 20°C ..... 1.3602 ÷ 1.3622
Water (K.F.) ..... ≤ 0.1 %
Residue on evaporation ..... ≤ 10 ppm
Substances darkened by H2SO4 ..... Conform

Acidity (acetic acid) ..... ≤ 10 ppm
Alcalinity (NH3) ..... ≤ 1 ppm
Isopropyl alcohol ..... ≤ 30 ppm
Methyl alcohol ..... ≤ 50 ppm
Benzene ..... ≤ 2 ppm(v/v)
Carbonyl compounds (CO) ..... ≤ 5 ppm
Subs. reducing KMnO4 ..... ≤ 3 ppm
Heavy metals (Pb) ..... ≤ 1 ppm
Absorbance UV (5cm, ref. water) ..... Conform
Assay (CPG) ..... ≥ 99.9 %

Assay (alcohometric) at 20°C ..... ≥ 99.9 %v/v
Volatil impurities ..... Conform
Al ..... ≤ 0.5 ppm
B ..... ≤ 0.02 ppm
Ba ..... ≤ 0.1 ppm
Ca ..... ≤ 0.5 ppm
Cd ..... ≤ 0.05 ppm
Co ..... ≤ 0.02 ppm
Cr ..... ≤ 0.02 ppm
Cu ..... ≤ 0.02 ppm

Fe ..... ≤ 0.1 ppm
Mg ..... ≤ 0.1 ppm
Mn ..... ≤ 0.02 ppm
Ni ..... ≤ 0.02 ppm
Pb ..... ≤ 0.1 ppm
Sn ..... ≤ 0.1 ppm
Acetal + acetaldehyde ..... ≤ 10 ppm(v/v)

Code	Size	Packaging	Notes
414601	1 l	Plastic bottle	Only for Italian market
414607	1 l	Glass bottle	Only for Italian market
4146012	1 l	Plastic bottle	
4146072	1 l	Glass bottle	
414605	2.5 l	Plastic bottle	Only for Italian market
414608	2.5 l	Glass bottle	Only for Italian market
4146052	2.5 l	Plastic bottle	
4146082	2.5 l	Glass bottle	
414603	5 l	Aluminium can	
414606	5 l	Plastic bottle	
524125	5 l	Plastic tank	
4146032	5 l	Aluminium bottle	Untaxed, for Italian license holders only
4146062	5 l	Plastic bottle	Untaxed, for Italian license holders only
5241252	5 l	Plastic tank	Untaxed, for Italian license holders only
414604	10 l	Plastic tank	
414609	25 l	Plastic tank	
4146092	25 l	Plastic tank	Untaxed, for Italian license holders only
414602	200 l	Plastic drum	
4146022	200 l	Plastic drum	Untaxed, for Italian license holders only

DASIT GROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

Ethanol absolute anhydrous > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-BP-JP

ERBApharm

Description ..... Clear colourless liquid
Identification (I.R.)..... Positive
Color of solution ..... Pass test
Clarity of solution..... Pass test
Density at 20°C ..... 0.790 - 0.793
Density at 15.56°C ..... ≤ 0.7962
Boiling point..... 78 - 79 °C

Residue on evaporation ..... ≤ 25 ppm(m/v)
Assay (alcohometric) at 15,56°C..... ≥ 99.5 %(v/v)
Assay (alcohometric) at 20°C.. ≥ 99.5 %v/v
Acidity or alcalinity..... ≤ 30 ppm
Volatil impurities..... Pass test
Methyl alcohol..... ≤ 75 ppm(v/v)

Acetal + acetaldehyde..... ≤ 10 ppm(v/v)
Benzene ..... ≤ 2 ppm(v/v)
Total other impurities..... ≤ 300 ppm(v/v)
Water (K.F.) ..... ≤ 0.1 %
Absorbance UV (5cm, ref. water) .. Pass test
At 240 nm ..... ≤ 0.40 AU
From 250 to 260 nm..... ≤ 0.30 AU

From 270 to 340 nm..... ≤ 0.10 AU
235 - 340 nm..... Smooth curve
Origin (BSE/TSE)..... Vegetable
Residual solvents (Current ICH)..... Conform

Code	Size	Packaging	Notes
529121	1 l	Glass bottle	
308661	2.5 l	Plastic bottle	Only for Italian market
308662	2.5 l	Glass bottle	Only for Italian market
3086612	2.5 l	Plastic bottle	
3086622	2.5 l	Glass bottle	
529122	5 l	Plastic tank	
5291222	5 l	Plastic tank	Untaxed, for Italian license holders only
529124	10 l	Plastic tank	
308664	25 l	Combined drum	
308667	25 l	Plastic tank	
3086642	25 l	Combined drum	Untaxed, for italian license holders only
308663	200 l	Plastic drum	
529125	200 l	Plastic drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Ethanol absolute anhydrous > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive

Density at 20° C ..... 0.7866 ÷ 0.7926
Boiling point..... 78.0 ÷ 79.0 ° C

Water (K.F.) ..... ≤0.1 %
Residue on evaporation ..... ≤30 ppm

Acidity (acetic acid)..... ≤50 ppm
Assay(alcohol.) at 20°C ..... ≥99.9 % (v/v)

Code	Size	Packaging	Notes
308602	1 l	Plastic bottle	Only for Italian market
308607	1 l	Glass bottle	Only for Italian market
308608	1 l	Plastic bottle	Origin: synthesis
3086022	1 l	Plastic bottle	
3086072	1 l	Glass bottle	
308603	2.5 l	Glass bottle	Only for Italian market
308605	2.5 l	Plastic bottle	Only for Italian market
3086032	2.5 l	Glass bottle	
3086052	2.5 l	Plastic bottle	
528131	5 l	Plastic tank	
308609	10 l	Plastic tank	
3086092	10 l	Plastic tank	Untaxed, for italian license holders only
308601	25 l	Metal drum	
308604	25 l	Plastic tank	
308600	200 l	Metal drum	

Ethanol 96°

• Etanolo 96° • Ethanol 96° • Etanol 96° • Ethanol 96°

Synonym:
Ethyl alcohol 96°

C_2H_5OH
Molecular Weight: 46,07
CAS: 64-17-5
EEC-N: 200-578-6

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II

Danger
H225-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

Ethanol 96° > RS - For HPLC - Isocratic Grade

RS

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Identification (I.R.)..... Positive
Density at 20°C .....0.8050 ÷ 0.8124
Boiling point..... 78.0 ÷ 79.0 °C
Acidity ..... ≤ 0.0002 meq/g
Residue on evaporation ..... ≤ 5 ppm
Assay (alcohometric) at 20°C..96.0 ÷ 96.3
%v/v
Transmittance
At 210 nm ..... ≥ 30 %
At 254 nm ..... ≥ 96 %
At 280 nm ..... ≥ 98 %
Assay (CPG)..... ≥ 99.8 %

Code	Size	Packaging	Notes
414541	1 l	Glass bottle	Only for Italian market
4145412	1 l	Glass bottle	
414542	2.5 l	Glass bottle	Only for Italian market
4145422	2.5 l	Glass bottle	

Ethanol 96° > RS - SPECTROSOL - For optical spectroscopy

RS

Description ..... Clear colourless liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Density at 20° C .....0.8050 ÷ 0.8124
Boiling point..... 78.0 ÷ 79.0 ° C
Acidity or alkalinity.....≤0.0005 meq/g
Residue on evaporation .....≤10 ppm
Assay(alcohol.) at 20°C ..... ≥96.0 %(v/v)
Fluorescence
at 254 nm .....≤2 ppb
at 365 nm .....≤2 ppb
U.V. Transmittance
at 210 nm ..... ≥35 %
at 220 nm ..... ≥55 %
at 230 nm ..... ≥72 %
at 250 nm ..... ≥90 %
at 270 nm ..... ≥98 %

Code	Size	Packaging	Notes
414667	1 l	Glass bottle	Only for Italian market
4146672	1 l	Glass bottle	

Ethanol 96° > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611002501	1 l	Glass bottle	Ref Ph.Eur 1002501

Product specifications are subject to changes. Please visit our website for updates.

Ethanol 96° > RPE - For analysis - ACS - Reag. Ph.Eur.

RPE

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification (I.R.)..... Positive
Appearance of solution ..... Conform USP
Absorbance ..... Conform Ph Eur
Volatile impurities..... Conform Ph Eur
Water miscibility..... Conform ACS
Substances darkened by sulphuric acid ..... Conform ACS

Density at 20° C .....0.8050 ÷ 0.8124
Boiling point..... 78.0 ÷ 79.0 ° C
Residue on evaporation .....≤10 ppm
Acidity (acetic acid)..... ≤30 ppm
Alcalinity (NH3)..... ≤0.0002 meq/g
Alcole isopropilico-acetone .....Conform
Methyl alcohol..... ≤0.1 %
Carbonyl Compounds (CO) .....≤5 ppm
Acetal + acetaldehyde..... ≤ 10 ppm(v/v)

Subst. reducing KMnO4.....≤3 ppm
Ba.....≤0.1 ppm
Ca.....≤0.5 ppm
Cd.....≤0.05 ppm
Co.....≤0.02 ppm
Cr.....≤0.02 ppm
Cu.....≤0.02 ppm
Fe .....≤0.1 ppm
Mg .....≤0.1 ppm

Mn .....≤0.02 ppm
Ni.....≤0.02 ppm
Pb .....≤0.1 ppm
Sn .....≤0.1 ppm
Zn .....≤0.1 ppm
Assay(alcohol.) at 20°C96.0 ÷ 96.9 % (v/v)
Assay (GLC) ..... ≥95.0 %

Code	Size	Packaging	Notes
414634	1 l	Plastic bottle	Only for Italian market
414637	1 l	Glass bottle	Only for Italian market
4146342	1 l	Plastic bottle	
4146372	1 l	Glass bottle	
414631	2.5 l	Glass bottle	Only for Italian market
414632	2.5 l	Plastic bottle	Only for Italian market
4146312	2.5 l	Glass bottle	
4146322	2.5 l	Plastic bottle	
414635	5 l	Plastic tank	
4146352	5 l	Plastic tank	Untaxed, for Italian license holders only
414638	10 l	Plastic tank	
414639	25 l	Plastic tank	
414633	200 l	Plastic drum	

Ethanol 96° > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Clear liquid
Colour (APHA) ..... ≤ 10
Identification A ..... Conform
Identification B ..... Conform
Identification C ..... Conform
Absorbance
at 240 nm ..... ≤ 0.40
at 250 - 260 nm ..... ≤ 0.30
at 270 - 340 nm ..... ≤ 0.10
Methyl alcohol..... ≤ 200 ppm
Acetalduhyde + acetal..... ≤ 10 ppm

Benzene ..... ≤ 2 ppm
Other impurities ..... ≤ 300 ppm
Furfural..... Conform
Water miscibility..... Conform ACS
Tail products ..... Conform
Substances darkened by sulphuric acid ..... Conform ACS
Density at 20° C .....0.8050 ÷ 08124
Boiling point..... 78.0 ÷ 79.0 ° C
Residue on evaporation ..... ≤ 10 ppm
Acidity (acetic acid)..... ≤ 10 ppm

Alcalinity (NH3)..... ≤ 1 ppm
Isopropyl alcohol ..... ≤ 30 ppm
Carbonyl Compounds (CO) ..... ≤ 5 ppm
Subst. reducing KMnO4 ..... ≤ 3 ppm
Al ..... ≤ 0.5 ppm
B..... ≤ 0.02 ppm
Ba..... ≤ 0.1 ppm
Ca..... ≤ 0.5 ppm
Cd..... ≤ 0.05 ppm
Co..... ≤ 0.02 ppm
Cr..... ≤ 0.02 ppm

Cu..... ≤ 0.02 ppm
Fe ..... ≤ 0.1 ppm
Mg ..... ≤ 0.1 ppm
Mn ..... ≤ 0.02 ppm
Ni..... ≤ 0.02 ppm
Pb ..... ≤ 0.1 ppm
Sn ..... ≤ 0.1 ppm
Zn ..... ≤ 0.1 ppm
Assay(alcohol.) at 20°C96.0 ÷ 96.9 % (v/v)
Assay(densim.) at 15.5°94.9 ÷ 96.0 % (v/v)

Code	Size	Packaging	Notes
414622	200 l	Plastic drum	

Ethanol 96° > ERBApharm - According to pharmacopoeia: Ph.Eur.- Microbiological tested

ERBApharm

Description ..... Clear colourless liquid
Identification (I.R.)..... Positive
Appearance of solution ..... Pass test
Acidity or alcalinity..... Pass test
Absorbance UV (5cm, ref. water) .. Pass test

Volatil impurities..... Pass test
Density at 20°C ..... 0.805 - 0.812
Boiling point..... 78 - 79 °C
Residue on evaporation ... ≤ 0.0025 %(m/v)
Assay (alcohometric) at 20°C... 95.1 - 96.9

%v/v
Origin (BSE/TSE)..... Vegetable
Total aerobic microbial count (TAMC) .... ≤ 5 CFU/100ml
Total yeasts/mould count (TYMC)..... ≤ 5

CFU/100ml
Test of specified micro-organisms
Enterobacteriaceae..... Absent/100 ml
Staphylococcus aureus ..... Absent/100 ml
Residual solvents (Current ICH)..... Conform

Code	Size	Packaging	Notes
524135	5 l	Plastic tank	
524132	25 l	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Ethanol 96° > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP

ERBApharm

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.805 ÷ 0.812
Densità a 15.5°C ..... 0.812 ÷ 0.816
Acidity or alkalinity ..... Conform Ph.Eur.
Residue on evaporation ..... ≤0.0025 % m/v
Assay(densim.) at 15.5° 94.9 ÷ 96.0 % v/v
Assay(alcohol.) at 20°C .95.1 ÷ 96.9 % v/v
Appearance of solution ..... Conform Ph.Eur.
Boiling point ..... ~ 78 ° C
Origin (BSE/TSE) ..... Vegetable
Residual solvents (Current ICH) ..... Conform
Color of solution ..... Conform USP-NF
Clarity of solution ..... Conform USP-NF
Volatil impurities ..... Pass test Ph.Eur.
Methyl alcohol ..... ≤ 200 ppm(v/v)
Acetal + acetaldehyde ..... ≤ 10 ppm(v/v)
Benzene ..... ≤ 2 ppm(v/v)
Total other impurities ..... ≤ 300 ppm(v/v)
Absorbance UV (5cm, ref. water) .. Pass test
At 240 nm ..... ≤ 0.40 AU
From 250 to 260 nm ..... ≤ 0.30 AU
From 270 to 340 nm ..... ≤ 0.10 AU
235 - 340 nm ..... Smooth curve

Code	Size	Packaging	Notes
308644	1 l	Plastic bottle	Only for Italian market
308647	1 l	Glass bottle	Only for Italian market
3086442	1 l	Plastic bottle	
3086472	1 l	Glass bottle	
308641	2.5 l	Glass bottle	Only for Italian market
308649	2.5 l	Plastic bottle	Only for Italian market
3086412	2.5 l	Glass bottle	
3086492	2.5 l	Plastic bottle	
529141	5 l	Plastic tank	
5291412	5 l	Plastic tank	Untaxed, for Italian license holders only
308646	10 l	Plastic tank	
3086462	10 l	Plastic tank	Untaxed, for italian license holders only
308645	25 l	Plastic tank	
3086452	25 l	Plastic tank	Untaxed, for italian license holders only
308648	27 l	Combined drum	
3086482	27 l	Combined drum	Untaxed, for italian license holders only
308643	200 l	Plastic drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Ethanol 96° > RE - Pure

RE

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Identification ..... Positive
Water miscibility ..... Complete
Acidity (acetic ac) ..... ≤ 0.005 %
Boiling point ..... 78 ÷ 79.5 °C
Density ..... 0.805 ÷ 0.8125
Assay GLC ..... ≥ 99 %
Residue on evaporation ..... ≤ 25 ppm
Assay (alcohol) ..... 95 ÷ 96.9 % (v/v)
Assay (GLC) ..... ≥ 99 %

Code	Size	Packaging	Notes
528151	5 l	Plastic tank	
528152	10 l	Plastic tank	
528154	20 l	Plastic tank	
529152	25 l	Plastic tank	
528153	200 l	Metal drum	

Ethanol 70% v/v

• Etanolo 70% v/v • Ethanol 70% v/v • Etanol 70% v/v • Ethanol 70% v/v

Synonym:
Ethyl alcohol 70%

C_2H_5OH
Molecular Weight: 46,07
CAS: 64-17-5

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II

Danger
H225-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

Ethanol 70% v/v > ERBApharm - According to pharmacopoeia: Ph.Eur.

ERBApharm

Description ..... Clear colourless liquid
Identification (I.R.) ..... Positif
Identity (IR) ..... Positive
Assay (alcohometric) at 20°C 69 ÷ 71 %v/v
Assay (alcohol) ..... 69 ÷ 71 % (v/v)
Ethyl alcohol (or.Tereos) ..... Batch number
Purified water ..... Batch number

Code	Size	Packaging	Notes
529187000	20 x 500 ml	Spray bottle	Sold by box: 20 bottles + 4 sprayers
529189	5 l	Plastic tank	
529183	200 l	Plastic drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Ethanol 70% v/v > ERBApharm - According to pharmacopoeia: BP-Ph.Eur.

ERBApharm

Description ..... Clear colourless liquid
Volatil impurities (GC)..... Passed test
Residue on evaporation .....≤ 25 ppm(m/v)
Acidity (acetic acid)..... ≤ 30 ppm
Apparent density .....882.2 ÷ 887.1 kg/m3
Assay (alcohometric) at 20°C..... 69.0 ÷ 71.0 %v/v

Code	Size	Packaging	Notes
529184	1 l	Spray bottle	Sold by box 6 bottles + 3 sprayers
529184000	1 l	Spray bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Ethanol 70% v/v > RE - Pure

RE

Description ..... Clear liquid
Identity (IR)..... Positive
Colour ..... Incolore
Assay (alcohol)..... 69 ÷ 71 % (v/v)
Identification ..... Positive
Acidity ..... ≤ 0.003 %
Alcalinity..... ≤ 3 ppm
Methanol ..... ≤ 100 ppm
Residue on evaporation .....≤ 30 ppm
Water miscibility..... Conform ACS
Substances darkened by sulphuric acid ..... Conform ACS
Substances reducing KMn04 (O)...≤ 3 ppm
Assay(alcohol.) at 20°C ..... ≥ 70 % (v/v)

Code	Size	Packaging	Notes
308771	2.5 l	Plastic bottle	
528170	5 l	Plastic tank	
529186	10 l	Plastic tank	
308775	25 l	Plastic tank	
3087752	25 l	Plastic tank	Untaxed, for italian license holders only

Ethanol 60% v/v

• Etanolo 60% v/v • Ethanol 60% v/v • Etanol 60% v/v • Ethanol 60°

Synonym:
Ethyl alcohol 60%

C_2H_5OH
Molecular Weight: 46,07
CAS: 64-17-5
EEC-N: 200-578-6

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II

Danger
H225-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

Ethanol 60% v/v > RE - Pure

RE

Description ..... Clear colourless liquid
Assay (alcohol)..... 59.0 ÷ 61.0 % (v/v)

Code	Size	Packaging	Notes
529180	5 l	Plastic tank	
529181	5 l	Plastic bottle	

Ethanol 50% v/v

• Etanolo 50% v/v • Ethanol 50% v/v • Etanol 50% v/v • Ethanol 50% v/v

Synonym:
Ethyl alcohol 50%

C_2H_5OH
Molecular Weight: 46,07
CAS: 64-17-5

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II

Danger
H225-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

Ethanol 50% v/v > ERBApharm - Prepared from raw material according Ph.Eur

ERBApharm

Description ..... Clear colourless liquid
Residue on evaporation ..... ≤ 25 ppm
Acidity (acetic acid)..... ≤ 30 ppm
Volatil impurities (GC)..... Pass test
Assay (alcohometric) at 20°C..... 49.0 - 51.0 %v/v

Code	Size	Packaging	Notes
529261	5 l	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Ethanol absolute denaturated

• Etanolo assoluto denaturato • Ethanol absolu dénaturé • Etanol absoluto desnaturalizado • Ethanol absolut vergällt

Synonym:
Ethyl alcohol absolute denaturated

C_2H_5OH
Molecular Weight: 46,07
CAS: 64-17-5

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II



Danger
H225-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

Ethanol absolute denaturated > RE - Pure - According to European regulation

RE

Description ..... Clear colourless liquid Colour ..... ≤ 10 APHA Density at 20°C .....0.787 ÷ 0.793 Assay (alcohol) ..... ≥ 99.2 % (v/v)

Code	Size	Packaging	Notes
528761	1 l	Glass bottle	
528765	2.5 l	Glass bottle	
528763	5 l	Metal tank	
528764	5 l	Plastic tank	
528766	10 l	Plastic tank	
528762	25 l	Plastic tank	

According to regulation 2016/1867 of 20th october 2016 deleting 3199/93. Eurodenaturant: 1l isopropanol IPA), 1l methylethylketone (MEK), 1g of denatonium benzoate per hectolitre of absolute alcohol

Ethanol absolute denaturated > RE - Pure - According to italian denaturating procedure

RE

Description ..... Clear pink liquid Density at 15° C ..... ~ 0.796 Assay(alcohol.) at 20°C ..... ≥ 99.5 % (V/V)
Identification ..... Positive Boiling point..... 77 ÷ 79 °C

Code	Size	Packaging	Notes
308651	1 l	Plastic bottle	
308653	5 l	Plastic bottle	
308656	10 l	Plastic tank	
308655	25 l	Metal drum	
308652	200 l	Metal drum	

Denaturing procedures authorized only on Italian market

Ethanol 95° denaturated

• Etanolo 95° denaturato • Ethanol 95° dénaturé • Etanol 95° desnaturalizado • Ethanol vergällt 95°

Synonym:
Ethyl alcohol 95° denaturated

C_2H_5OH
Molecular Weight: 46,07
CAS: 64-17-5

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II



Danger
H225-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

Ethanol 95° denaturated > RE - Pure - According to European regulation

RE

Description ..... Clear colourless liquid Colour ..... ≤ 10 APHA Density at 20°C .....0.804 ÷ 0.812 Assay (alcohol) .....94.8 ÷ 96.8 % (v/v)

Code	Size	Packaging	Notes
528771	1 l	Glass bottle	
528775	5 l	Plastic tank	
528772	10 l	Plastic tank	
528773	25 l	Plastic tank	
528774	200 l	Metal drum	

According to regulation 2016/1867 of 20th october 2016 deleting 3199/93. Eurodenaturant: 1l isopropanol IPA), 1l methylethylketone MEK), 1g of denatonium benzoate per hectolitre of absolute alcohol

Product specifications are subject to changes. Please visit our website for updates.

Ethanol 94° denaturated

• Etanolo 94° denaturato • Ethanol 94° dénaturé • Etanol 94° desnaturalizado • Ethanol vergällt 94°

Synonym:
Ethyl alcohol 94° denaturated

C_2H_5OH
Molecular Weight: 46,07
CAS: 64-17-5

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II



Danger
H225-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

Ethanol 94° denaturated > RE - Pure - According to italian denaturating procedure

RE

Description .....Clear pink liquid Identification .....Positive Density at 15° C .....0.815 ÷ 0.825 Assay(alcohol.) at 20°C .....92 ÷ 96 %

Code	Size	Packaging	Notes
308621	1 l	Plastic bottle	
308623	5 l	Plastic bottle	
308625	10 l	Plastic tank	
308624	25 l	Metal drum	

Denaturing procedures authorized only on Italian market

Ethanol 90° denaturated

• Etanolo 90° denaturato • Ethanol 90° dénaturé • Etanol 90° desnaturalizado • Ethanol vergällt 90°

Synonym:
Ethyl alcohol 90° denaturated

C_2H_5OH
Molecular Weight: 46,07
CAS: 64-17-5

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II



Danger
H225-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

Ethanol 90° denaturated > RE - Pure - According to italian denaturating procedure

RE

Description .....Clear pink liquid Identification .....Positive Density at 15° C .....0.830 ÷ 0.840 Assay(alcohol.) at 20°C .....88 ÷ 92 %

Code	Size	Packaging	Notes
308681	1 l	Plastic bottle	
308683	5 l	Plastic bottle	
308682	25 l	Metal drum	
308687	160 kg	Metal drum	

Denaturing procedures authorized only on Italian market

Ethanol 70°modified

• Etanolo 70° modificato • Ethanol 70° modifié • Etanol 70° modificada • Ethanol 70° geändert

Synonym:
Ethyl alcohol 70° denaturated

C_2H_5OH
Molecular Weight: 46,07
CAS: 64-17-5

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group III



Warning
H226-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

Ethanol 70°modified > RE - Pure

RE

Assay .....69.5 ÷ 71.5 % v/v Density at 20°C .....0.881 ÷ 0.886

Code	Size	Packaging	Notes
528191	5 l	Plastic tank	
528192	20 l	Plastic tank	

Colored Yellow. Modified with camphor and tartrazine

Ethanol-d6 anhydrous

• Etanolo anidro-d6 • Ethanol anhydre-d6 • Etanol anhidro-d6 • Ethanol-d6

Synonym:
• *Ethanol-d6*
• *Ethyl alcohol-d6*

C_2D_5OD
Molecular Weight: 52,11
CAS: 1516-08-1
EEC-N: 216-162-2

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II

Danger
H225
P210-P241-P280-P303+P361+P353-P403+P235-P501a

Ethanol-d6 anhydrous > RS - For NMR - min 99%

RS

Code	Size	Packaging	Notes
P5262A	2 x 1 ml	Glass ampoule	

For specifications, contact our customer service for a certificate of analysis

Ethanolamine

• Etanolammina • Ethanolamine • Etanolamina • Ethanolamin

Synonym:
2-Aminoethanol

$NH_2CH_2CH_2OH$
Molecular Weight: 61,08
CAS: 141-43-5
EEC-N: 205-483-3

Classification transport
ONU: 2491
Transport Hazard class: 8
Packing group III

Danger
H302-H312-H332-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P362+P364

Ethanolamine > RPE - For analysis

RPE

Description ..Clear colourless to pale yellow liquid
Identification ..... Positive
Water miscibility.....Conform
Alcohol miscibility..... Complete
Refractive index at 20°C.1.4491 ÷ 1.4591
Boiling point.....169.5 ÷ 170.5 °C
Melting point.....9.8 ÷ 10.8 ° C
Chloride.....≤10 ppm
Diethanolamine ..... ≤0.5 %
Heavy metals (Pb).....≤2 ppm
Sulphate.....≤20 ppm
Triethanolamine..... ≤0.5 %
Fe .....≤1 ppm
Assay (alkalimetric)..... ≥99 %

Code	Size	Packaging	Notes
447351	1 l	Glass bottle	
447352	30 kg	Aluminium can	

Hygroscopic product. Store well sealed in a dry place

Ethyl acetate

• Etile acetato • Ethyle acétate • Etilo acetato • Ethylacetat

$CH_3COOC_2H_5$
Molecular Weight: 88,11
CAS: 141-78-6
EEC-N: 205-500-4

Classification transport
ONU: 1173
Transport Hazard class: 3
Packing group II

Danger
H225-H319-H336-HEU066
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Ethyl acetate > RS - For LC/MS

RS

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Identification (I.R.)..... Positive
Refractive index at 20°C..... 1.370 - 1.374
Water (K.F.) ..... ≤ 200 ppm
Residue on evaporation ..... ≤ 2 ppm
Acidity (acetic acid)..... ≤ 0.0030 %
Alcalinity (NH3)..... ≤ 0.0005 %
Assay (CPG) ..... ≥ 99.95 %
Transmittance
At 260 nm ..... ≥ 75 %
At 275 nm ..... ≥ 97 %
At 300 nm ..... ≥ 98 %
Metals compounds
Al ..... ≤ 50 ppb
Fe ..... ≤ 50 ppb
Na..... ≤ 50 ppb
Ca..... ≤ 50 ppb
Mg ..... ≤ 50 ppb
K..... ≤ 50 ppb

Code	Size	Packaging	Notes
448383	1 l	Glass bottle	
448384	2.5 l	Glass bottle	

Ethyl acetate > RS - For HPLC - Isocratic Grade

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C .....0.898 ÷ 0.902
Refractive index at 20°C.1.3699 ÷ 1.3739
Boiling point.....76.9 ÷ 77.4 ° C
Water (K.F.) .....≤300 ppm
Residue on evaporation .....≤10 ppm
Acidity or alkalinity.....≤0.0015 meq/g
Assay (GLC) ..... ≥99.9 %
U.V. Transmittance
at 260 nm ..... ≥76 %
at 270 nm ..... ≥94 %
at 300 nm ..... ≥ 97 %
Methyl alcohol.....≤ 100 ppm
Ethyl alcohol .....≤ 400 ppm

Code	Size	Packaging	Notes
412611000	1 l	Glass bottle	
412612000	2.5 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Ethyl acetate > RS - For preparative HPLC - Reag.Ph.Eur.

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.901 ÷ 0.902
Refractive index at 20°C. 1.3699 ÷ 1.3739
Boiling point ..... 76.9 ÷ 77.4 ° C
Water (K.F.) ..... ≤ 200 ppm
Residue on evaporation ..... ≤ 5 ppm
Assay (GLC) ..... ≥ 99.9 %
U.V. Transmittance
at 265 nm ..... ≥ 80 %
at 300 nm ..... ≥ 95 %

Code	Size	Packaging	Notes
448211	2.5 l	Glass bottle	

Ethyl acetate > RS - For GC-MS

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C ..... 1.370 - 1.374
Water (K.F.) ..... ≤ 150 ppm
Residue on evaporation ..... ≤ 2 ppm
Colour ..... ≤ 10 APHA
Acidity (acetic acid) ..... ≤ 30 ppm
Assay (GC) ..... ≥ 99.9 %
Ethyl alcohol ..... ≤ 200 ppm
Methyl alcohol ..... ≤ 100 ppm
GC-MS.Individual peak (n-hexadecane). ≤ 2 µg/L
Ret.range n-undecane to n-tetracontane
(scanning area 30-600amu)

Code	Size	Packaging	Notes
448342	1 l	Glass bottle	

Ethyl acetate > RS - ATRASOL - For traces analysis

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C ..... 1.370 - 1.374
Water content (K.F.) ..... ≤ 150 mg/Kg
Colour ..... ≤ 10 Hazen
Methanol ..... ≤ 100 mg/Kg
Non volatile residue ..... ≤ 2 mg/Kg
Ethanol ..... ≤ 200 mg/Kg
Assay (GC) ..... ≥ 99.9 %
Free acid (as CH3COOH) ..... ≤ 30 mg/Kg
GC (FID) - NC Atrasol ..... Conform
GC-ECD.Individual peak (Lindane) .≤ 2 ng/L
Ret.range 1,2,4-trichlorobenzene
to decachlorobiphenyle
GC-FID.Individual peak (n-hexadecane). ≤ 2 µg/L
Ret.range n-undecane to n-tetracontane

Code	Size	Packaging	Notes
P0023216	1 l	Glass bottle	
P0023221	2.5 l	Glass bottle	

Ethyl acetate > RS - PESTIPUR - For pesticide analysis

RS

Description ..... Clear liquid
Identification ..... Positive
Colour ..... ≤ 10 hazen
Assay (GLC) ..... ≥ 99.8 %
Water ..... ≤ 0.03 %
Not volatile residue ..... ≤ 2 mg/kg
Free acid (as CH3COOH) ..... ≤ 30 mg/kg
GC-ECD (Lindane standard) ..... ≤ 3 ng/l
GC-NPD (Ethylparathion standard) ..... ≤ 3 ng/l

Code	Size	Packaging	Notes
448351	1 l	Glass bottle	
448352000	2.5 l	Glass bottle	

Ethyl acetate > RS - SPECTROSOL - For optical spectroscopy

RS

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Density at 20° C ..... 0.898 ÷ 0.902
Refractive index at 20°C. 1.3699 ÷ 1.3739
Boiling point ..... 76.9 ÷ 77.4 ° C
Water (K.F.) ..... ≤100 ppm
Residue on evaporation ..... ≤10 ppm
Acidity ..... ≤0.0005 meq/g
Alcalinity ..... ≤0.0002 meq/g
Ethyl alcohol ..... ≤ 100 ppm
Methyl alcohol ..... ≤100 ppm
Assay (GLC) ..... ≥99.8 %
U.V. Transmittance
at 255 nm ..... ≥15 %
at 260 nm ..... ≥75 %
at 270 nm ..... ≥95 %
at 280 nm ..... ≥98 %

Code	Size	Packaging	Notes
448271	1 l	Glass bottle	
448272	2.5 l	Glass bottle	

Ethyl acetate > RS - Anhydrous - For analysis

RS

Refractive index at 20°C ..... 1.37 - 1.374
Water content (K.F.) ..... ≤ 100 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
Colour ..... ≤ 10 Hazen
Free acid (as CH3COOH) ..... ≤ 30 mg/Kg
Assay (GC) ..... ≥ 99.8 %
Methanol ..... ≤ 100 mg/Kg
Ethanol ..... ≤ 400 mg/Kg
Density d20/4 ..... 0.898 - 0.902
Methyl acetate ..... ≤ 0.10 %

Code	Size	Packaging	Notes
P0021010	200 ml	Bottle with septum	
P0021016	1 l	Glass bottle	
P0021021	2.5 l	Glass bottle	
P00210T21	2.5 l	Glass bottle	On molecular sieves 4A, Water content < 20ppm

Ethyl acetate > RS - RSE - For electronic use

RS

Description ..... Clear liquid	Ag ..... ≤0.02 ppm	Cr ..... ≤0.01 ppm	Ni ..... ≤0.01 ppm
Colour (APHA) ..... ≤10	Al ..... ≤0.05 ppm	Cu ..... ≤0.01 ppm	Pb ..... ≤0.02 ppm
Identification ..... Positive	As ..... ≤0.01 ppm	Fe ..... ≤0.02 ppm	Pt ..... ≤0.05 ppm
Ready carbonizable substances ..... Conform	Au ..... ≤0.05 ppm	Ga ..... ≤0.02 ppm	Sb ..... ≤0.01 ppm
Assay (GLC) ..... ≥99.9 %	B ..... ≤0.01 ppm	In ..... ≤0.02 ppm	Sn ..... ≤0.02 ppm
Resistivity ..... ≥20 Mohm.cm	Ba ..... ≤0.1 ppm	K ..... ≤0.1 ppm	Sr ..... ≤0.02 ppm
Density at 20° C ..... 0.898 ÷ 0.902	Be ..... ≤0.02 ppm	Li ..... ≤0.02 ppm	Ti ..... ≤0.05 ppm
Boiling point ..... 76.6 ÷ 77.6 ° C	Bi ..... ≤0.02 ppm	Mg ..... ≤0.1 ppm	Tl ..... ≤0.05 ppm
Water (K.F.) ..... ≤500 ppm	Ca ..... ≤0.2 ppm	Mn ..... ≤0.01 ppm	V ..... ≤0.05 ppm
Residue on evaporation ..... ≤10 ppm	Cd ..... ≤0.01 ppm	Mo ..... ≤0.05 ppm	Zn ..... ≤0.02 ppm
Acidity (acetic acid) ..... ≤50 ppm	Co ..... ≤0.01 ppm	Na ..... ≤0.2 ppm	Zr ..... ≤0.05 ppm

Code	Size	Packaging	Notes
448307	1 l	Glass bottle	
448308	2.5 l	Glass bottle	
448306	5 l	Metal tank	

Ethyl acetate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Clear liquid	Refractive index at 20°C 1.3699 ÷ 1.3739	Al ..... ≤0.5 ppm	Fe ..... ≤0.02 ppm
Colour (APHA) ..... ≤10	Boiling point ..... 76.9 ÷ 77.4 ° C	B ..... ≤0.02 ppm	Mg ..... ≤0.1 ppm
Identification (I.R.) ..... Conform	Water (K.F.) ..... ≤300 ppm	Ba ..... ≤0.1 ppm	Mn ..... ≤0.02 ppm
Alcohol miscibility ..... Complete	Residue on evaporation ..... ≤10 ppm	Ca ..... ≤0.5 ppm	Ni ..... ≤0.02 ppm
Foreign esters ..... Conform	Acidity ..... ≤0.0009 meq/g	Cd ..... ≤0.05 ppm	Pb ..... ≤0.02 ppm
Water solubility ..... Conform	Ethyl alcohol ..... ≤0.05 %	Co ..... ≤0.02 ppm	Sn ..... ≤0.1 ppm
Substances darkened by sulphuric acid ..... Conform	Methyl alcohol ..... ≤0.01 %	Cr ..... ≤0.02 ppm	Zn ..... ≤0.02 ppm
Density at 20° C ..... 0.901 ÷ 0.902	Methyl acetate ..... ≤0.1 %	Cu ..... ≤0.02 ppm	Assay (GLC) ..... ≥99.9 %

Code	Size	Packaging	Notes
448251	1 l	Glass bottle	
448256	2.5 l	Glass bottle	
448254	5 l	Plastic tank	
448252	10 l	Metal tank	
448258	10 l	Plastic tank	
448253	24 kg	Metal drum	
448255	180 kg	Metal drum	

Ethyl acetate > ERBApharm - According to pharmacopoeia: NF

ERBApharm

Description ..... Clear colourless liquid	Ready carbonizable substances ..... Conform USP-NF	Organic volatile impurities Conform USP-NF	Water (K.F.) ..... ≤0.1 %
Identification ..... Positive	Methyl compounds ..... Conform USP-NF	Density at 25° C ..... 0.894 ÷ 0.898	Cromatographic purity ..... ≥99.5 %
Acidity ..... Conform USP-NF		Residue on evaporation ..... ≤0.02 %	Assay (saponification) ..... 99.0 ÷ 100.5 %

Code	Size	Packaging	Notes
341511	1 l	Glass bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Ethyl acetate > ERBApharm - According to pharmacopoeia: DAB-NF-Ph.Eur.

ERBApharm

Description ..... Clear colourless liquid	Refractive index at 20°C ..... 1.370 ÷ 1.373	Readily carbonizable substances ... Conform USP-NF	Ethyl alcohol ..... ≤ 0.05 %
Colour ..... ≤ 10 APHA	Boiling point ..... 76 ÷ 78 °C	React. with ac. sulfuric ..... Conform Ph.Eur.	Titolo (saponificazione) ..... 99.0 ÷ 100.5 %
Appearance ..... Conform Ph.Eur.	Acidity (acetic acid) ..... ≤ 0.005 %	Methyl compounds ..... Conform USP-NF	Assay (GLC) ..... ≥ 99.9 %
Identification ..... Positive	Acidity ..... Conform USP-NF	Organic volatile impurities Conform USP-NF	Origin (BSE/TSE) ..... Synthesis
Density at 20°C ..... 0.898 ÷ 0.902	Residue on evaporation ..... ≤ 0.003 %	Related compounds ..... ≤ 0.2 %	Residual solvents (Current ICH) ..... Conform
Density at 25° C ..... 0.894 ÷ 0.898	Water (K.F.) ..... ≤ 0.04 %		

Code	Size	Packaging	Notes
341506	1 l	Glass bottle	
341503	2.5 l	Glass bottle	
341502	24 kg	Metal drum	
529221	25 l	Aluminium can	
529222	200 l	Metal drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Product specifications are subject to changes. Please visit our website for updates.

Ethyl acetate > RE - Pure

RE

Description	Clear colourless liquid	Density at 20°C	0.898 ÷ 0.902	Acidity (acetic ac)	≤ 30 ppm	Assay (GLC)	≥ 99.8 %
Colour	≤ 10 APHA	Boiling point	76.0 ÷ 77.5 °C	Residue on evaporation	≤ 20 ppm	Ethyl alcohol	≤ 0.04 %
Identity (IR)	Positive	Refractive index at 20°C	1.3699 ÷ 1.3739	Water (K.F.)	≤ 500 ppm		

Code	Size	Packaging	Notes
508221	1 l	Glass bottle	
508222	2.5 l	Glass bottle	
528295	5 l	Plastic tank	
528299	10 l	Metal tank	
528294	25 l	Metal drum	
528296	25 l	Plastic tank	
528297	200 l	Metal drum	

Ethyl acetoacetate

• Etile acetoacetato • Ethyle acétoacétate • Etilo acetoacetato • Ethylacetoacetat

Synonym:
Acetoacetic ester

$CH_3CH_2OCOCH_2COCH_3$
Molecular Weight: 130,14
CAS: 141-97-9
EEC-N: 205-516-1

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Ethyl acetoacetate > RE - Pure

RE

Description	Clear colourless liquid	Density at 20° C	1.021 ÷ 1.029	Boiling point	179 ÷ 181 ° C	Assay (GLC)	≥ 98 %
Identification	Positive	Refractive index at 20°C	1.415 ÷ 1.424	Water (K.F.)	≤ 0.1 %		

Code	Size	Packaging	Notes
341751	1 l	Glass bottle	

Ethyl alcohol absolute anhydrous ▶ Ethanol absolute anhydrous

Ethyl alcohol 96° ▶ Ethanol 96°

Ethyl alcohol 70% ▶ Ethanol 70% v/v

Ethyl alcohol 60% ▶ Ethanol 60% v/v

Ethyl alcohol 50% ▶ Ethanol 50% v/v

Ethyl alcohol absolute denaturated ▶ Ethanol absolute denaturated

Ethyl alcohol 95° denaturated ▶ Ethanol 95° denaturated

Ethyl alcohol 94° denaturated ▶ Ethanol 94° denaturated

Ethyl alcohol 90° denaturated ▶ Ethanol 90° denaturated

Ethyl alcohol 70° denaturated ▶ Ethanol 70°modified

Ethyl diisopropylamine ▶ N,N-Diisopropylethylamine

Ethylenediamine

• Etilendiammina • Ethylènediamine • Etilendiamina • Ethylendiamin

Synonym:
1,2-Diaminoethane

$NH_2CH_2CH_2NH_2$
Molecular Weight: 60,1
CAS: 107-15-3
EEC-N: 203-468-6

Classification transport
ONU: 1604
Transport Hazard class: 8
Packing group II



Danger
H226-H302-H311-H332-H314-H334-H317
P210-P280-P284-P301+P330+P331-
P303+P361+P353-P304+P340-P310a-
P305+P351+P338-P361+P364-P342+P311a

Ethylenediamine > RS - RSE - For electronic use

RS

Description	Clear liquid	Freezing point	10.2 ÷ 10.8 ° C	Cu	≤ 50 ppb	Assay (GLC)	≥ 98.0 %
Colour (APHA)	≤ 10	Residue on ignition	≤ 100 ppm	Fe	≤ 50 ppb		
Density at 20° C	0.890 ÷ 0.906 g/ml	Cr	≤ 50 ppb	Ni	≤ 50 ppb		

Code	Size	Packaging	Notes
449451	180 kg	Plastic drum	

Ethylenediamine > RPE - For analysis

RPE

Description	Clear colourless liquid	Density at 20° C	0.890 ÷ 0.906	Melting point	10 ÷ 12 °C	Fe	≤ 5 ppm
Identification	Positive	Refractive index at 20°C	1.4470 ÷ 1.4570	Heavy metals (Pb)	≤ 5 ppm	Assay (GLC)	≥ 98 %
Alcohol miscibility	Complete	Boiling point	115.8 ÷ 117.3 °C	Residue on ignition	≤ 100 ppm		

Code	Size	Packaging	Notes
449425	1 l	Glass bottle	
449426	5 l	Plastic tank	
449424	25 kg	Combined drum	

Ethylenediaminetetraacetic acid

• Acido etilendiamminotetracetico • Acide éthylènediaminetétracétique • Acido etilendiaminotetracético
• Ethylendiamintetraessigsäure

Synonym:
• *EDTA*
• *(Ethylenedinitrilo)tetraacetic acid*

$[CH_2N(CH_2COOH)_2]_2$
Molecular Weight: 292,24
CAS: 60-00-4
EEC-N: 200-449-4



Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Ethylenediaminetetraacetic acid > RPE - For analysis

RPE

Description	White powder	Loss on drying	≤0.1 %	NH4OH-Insoluble subst.	≤50 ppm	Cu	≤2 ppm
Identification	Positive	Nitrilotriacetic acid	≤0.15 %	Heavy metals (Pb)	≤5 ppm	Fe	≤5 ppm
Chelation power	Conform	Chloride	≤40 ppm	Residue on ignition	≤0.1 %	Assay (complexometric)	≥99 %

Code	Size	Packaging	Notes
405465	250 g	Plastic bottle	
405463	1 kg	Plastic bottle	

Suitable for complexometry

Ethylenediaminetetraacetic acid > ERBApharm - According to pharmacopoeia: NF

ERBApharm

Description	White powder	Nitrilotriacetic acid	≤0.3 %	Fe	≤50 ppm
Identification	Positive	Residue on ignition	≤0.2 %	Assay (complexometric)	98.0 ÷ 100.5 %

Code	Size	Packaging	Notes
303251	5 kg	Plastic tank	
303252	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Product specifications are subject to changes. Please visit our website for updates.

Ethylenediaminetetraacetic acid dipotassium salt dihydrate

- Acido etilendiamminotetracetico sale bipotassico diidrato
- Acide éthylènediaminetétracétique sel dipotassique dihydraté
- Acido etilendiaminotetracético sal dipotasica dihidrato
- Ethylendiamintetraessigsäure-Dikaliumsalz-Dihydrat

Synonym:
- *Dipotassium ethylenediaminetetraacetate dihydrate*
- *EDTA dipotassium salt*

$C_{10}H_{14}N_2O_8K_2.2H_2O$
Molecular Weight: 404,46
CAS: 25102-12-9

Warning
H302-H312-H332
P261-P264-P271-P280h-P301+P312a-P304+P340

Ethylenediaminetetraacetic acid dipotassium salt dihydrate > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
pH sol. 5% at 25° C ..... 4.0 ÷ 5.0
Chloride ..... ≤50 ppm
Heavy metals (Pb) ..... ≤10 ppm
Fe ..... ≤10 ppm
Assay (complexometric) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
405531	50 g	Glass bottle	

Suitable for complexometry

Ethylenediaminetetraacetic acid dipotassium salt dihydrate > RE - Pure

RE

Description ..... White granular powder
Identification ..... Positive
pH sol. 5% in H2O ..... 4.0÷5.0
Water-insoluble matter ..... ≤0.003 %
Fe ..... ≤ 0.001 %
Heavy metals (Pb) ..... ≤ 0.0005 %
Chloride ..... ≤0.004 %
Cyanide ..... ≤ 0.001 %
Sulphate ..... ≤0.01 %
Ca ..... ≤ 0.0005 %
Cu ..... ≤0.0001 %
Nitrilotriacetic acid ..... ≤0.05 %
Loss on drying 150° C ..... 8.5÷9.5 %
Assay (complexometric) ..... ≥99.0 %

Code	Size	Packaging	Notes
405582	25 kg	Plastic bucket	

Suitable for complexometry

Ethylenediaminetetraacetic acid disodium salt

- Acido etilendiamminotetracetico sale bisodico
- Acide éthylènediaminetétracétique sel disodique
- Acido etilendiaminotetracético sal disódica
- Ethylendiamintetraessigsäure dinatriumsalz

Synonym:
- *Disodium ethylenediaminetetraacetate dihydrate*
- *EDTA disodium salt*

$C_{10}H_{14}N_2Na_2O_8.2H_2O$
Molecular Weight: 372,24
CAS: 6381-92-6
EEC-N: 205-358-3

Warning
H332-H373
P260-P271-P304+P340-P312a-P314-P501a

Ethylenediaminetetraacetic acid disodium salt > RPE - For analysis - ACS

RPE

Description ..... White powder
Identification ..... Positive
pH sol. 5% at 25° C ..... 4.0 ÷ 6.0
Nitrilotriacetic acid ..... ≤0.1 %
Water-insoluble matter ..... ≤50 ppm
Heavy metals (Pb) ..... ≤50 ppm
Fe ..... ≤100 ppm
Assay (complexometric) ..... 99.0 ÷ 101.0 %

Code	Size	Packaging	Notes
405494	100 g	Plastic bottle	
405491	250 g	Plastic bottle	
405497	1 kg	Plastic bottle	
405492	25 kg	Plastic bucket	

Suitable for complexometry

Ethylenediaminetetraacetic acid disodium salt > ERBApharm - According to pharmacopoeia: FU-Ph.Eur.

ERBApharm

Description ..... White powder
Identification ..... Positive
Appearance of solution ..... Conform Ph.Eur.
pH sol. 5% at 25° C ..... 4.0 ÷ 5.5
Heavy metals (Pb) ..... ≤20 ppm
Fe ..... ≤80 ppm
Assay (complexometric) ..... 98.5 ÷ 101.0 %
Origin (BSE/TSE) ..... Synthesis
Residual solvents (Current ICH) ..... Conform
Impurity A ..... ≤ 0.1 %

Code	Size	Packaging	Notes
303201	1 kg	Plastic bottle	
303203	5 kg	Plastic tank	
303202	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Ethylenediaminetetraacetic acid disodium salt > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP

ERBApharm

Description ..... White powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Ca ..... Conform USP-NF
pH sol. 5% at 25° C .....4.0 ÷ 5.5
Loss on drying .....8.7 ÷ 11.4 %
Ac. nitrilotriacetico (Impurezza A) ... ≤ 0.1 %
Heavy metals (Pb).....≤ 20 ppm
Fe .....≤ 80 ppm
Assay (complexometric)99.0 ÷ 101.0 % s.s.
Origin (BSE/TSE)..... Synthesis
Residual solvents (Current ICH)..... Conform

Code	Size	Packaging	Notes
303227	1 kg	Plastic bottle	
303226	5 kg	Plastic tank	
303225	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Ethylenediaminetetraacetic acid disodium salt 0.1 mol/l (0.2N)

• Acido etilendiamminotetracetico sale bisodico 0.1 mol/l (0.2N) • Acide éthylènediaminetétracétique sel disodique 0.1 mol/l (0.2N) • Acido etilendiaminotetracético sal disódica 0.1 mol/l (0.2N) • Ethylendiamintetraessigsäure dinatriumsalz 0.1 mol/l (0.2N)

Synonym:
• *Disodium ethylenediaminetetraacetate dihydrate*
• *EDTA disodium salt*

$C_{10}H_{14}N_2Na_2O_8.2H_2O$
Molecular Weight: 372,24
CAS: 6381-92-6

HEU210

Ethylenediaminetetraacetic acid disodium salt 0.1 mol/l (0.2N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613005901	500 ml	Plastic bottle	Ref Ph.Eur 3005900
613005900	1 l	Plastic bottle	Ref Ph.Eur 3005900

Ethylenediaminetetraacetic acid disodium salt 0.1 mol/l (0.2N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid NIST 682
Assay (colorimetry) ..... 0.1996 - 0.2004 N

Code	Size	Packaging	Notes
405511000	1 l	Plastic bottle	Certified with NIST traceability
405513000	5 l	Kubidos	Certified with NIST traceability
405514000	5 l	Plastic tank	Certified with NIST traceability
405512000	10 l	Kubidos	Certified with NIST traceability

37.22 g of EDTA disodium salt. Volumetric solution ready-to-use

Ethylenediaminetetraacetic acid disodium salt 0.1 mol/l (0.2N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Conform Titration factor..... 1.000 ± 0.005

Code	Size	Packaging	Notes
405421		Plastic ampoule	Volume: 165 ml

Volumetric concentrated solution to prepare 1 L of solution 0,1 M

Product specifications are subject to changes. Please visit our website for updates.

Ethylenediaminetetraacetic acid disodium salt 0.05 mol/l (0.1N)

- Acido etilendiamminotetracetico sale bisodico 0.05 mol/l (0.1N)
- Acide éthylènediaminetétracétique sel disodique 0.05 mol/l (0.1N)
- Acido etilendiaminotetracético sal disódica 0.05 mol/l (0.1N)
- Ethylendiamintetraessigsäure dinatriumsalz 0.05 mol/l (0.1N)

Synonym:
- *Disodium ethylenediaminetetraacetate dihydrate*
- *EDTA disodium salt*

$C_{10}H_{14}N_2Na_2O_8.2H_2O$ HEU210
Molecular Weight: 372,24
CAS: 6381-92-6

Ethylenediaminetetraacetic acid disodium salt 0.05 mol/l (0.1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (colorimetry) ..... 0.0998 - 0.1002 N NIST 682

Code	Size	Packaging	Notes
405501000	1 l	Plastic bottle	Certified with NIST traceability
405502000	5 l	Plastic tank	Certified with NIST traceability

18.61 g of EDTA disodium salt. Volumetric solution ready-to-use

Ethylenediaminetetraacetic acid disodium salt 0.01 mol/l (0.02N)

- Acido etilendiamminotetracetico sale bisodico 0.01 mol/l (0.02N)
- Acide éthylènediaminetétracétique sel disodique 0.01 mol/l (0.02N)
- Acido etilendiaminotetracético sal disódica 0.01 mol/l (0.02N)
- Ethylendiamintetraessigsäure dinatriumsalz 0.01 mol/l (0.02N)

Synonym:
- *Disodium ethylenediaminetetraacetate dihydrate*
- *EDTA disodium salt*

$C_{10}H_{14}N_2Na_2O_8.2H_2O$ HEU210
Molecular Weight: 372,24
CAS: 6381-92-6

Ethylenediaminetetraacetic acid disodium salt 0.01 mol/l (0.02N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid NIST 682 Assay (colorimetry) ..... 0.01996 - 0.02004 N

Code	Size	Packaging	Notes
405442000	1 l	Plastic bottle	Certified with NIST traceability
405443000	5 l	Kubidos	Certified with NIST traceability

3.722 g of EDTA disodium salt. Volumetric solution ready-to-use

Ethylenediaminetetraacetic acid disodium salt 0.01 mol/l (0.02N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
405431		Plastic ampoule	Volume: 55 ml

3,3621 g EDTA. Volumetric concentrated solution to prepare 1 L of solution 0,01 M

Ethylenediaminetetraacetic acid potassium and magnesium salt dihydrate

- Acido etilendiamminotetracetico sale di potassio e magnesio diidrato
- Acide éthylènediaminetétracétique sel de potassium et de magnésium dihydraté
- Acido etilendiaminotetracético sal dipotasica-magnésica dihidrato
- Ethylendiamintetraessigsäure-Kaliumsalz und Magnesiumdihydrat

Synonym:
- *EDTA-K2Mg*
- *Ethylenediaminetetraacetic acid dipotassium magnesium salt*

$C_{10}H_{12}K_2MgN_2O_8 \cdot 2H_2O$
Molecular Weight: 462,8
CAS: 15708-48-2
EEC-N: 239-803-8

Warning
H302-H312-H332
P261-P264-P271-P280h-P301+P312a-P304+P340

Ethylenediaminetetraacetic acid potassium and magnesium salt dihydrate > RPE - For analysis

RPE

Description ..... White powder Heavy metals (as Pb) ..... ≤ 10 ppm Free magnesium ..... ≤ 0.01 % Water (K.F.) ..... 7.5 - 9.5 %
Identification (I.R.) ..... Conform Fe ..... ≤ 10 ppm Free EDTA ..... ≤ 0.05 % Assay (complexometry) ..... ≥ 98 % s.s

Code	Size	Packaging	Notes
405541	100 g	Plastic bottle	

Suitable for complexometry

Ethylenediaminetetraacetic acid tetrasodium salt tetrahydrate

- Acido etilendiamminotetracetico sale tetrasodico tetraidrato
- Acide éthylènediaminetétracétique sel tétracétique tétrahydraté
- Acido etilendiaminotetracético sal tetrasódica tetrahidrato
- Ethylendiamin-tetraessigsäure Tetranatriumsalz Tetrahydrat

Synonym:
- *EDTA tetrasodium salt*
- *Tetrasodium ethylenediaminetetraacetate tetrahydrate*

$[CH_2N(CH_2COONa)_2]_2.4H_2O$
Molecular Weight: 452,24
CAS: 13235-36-4
EEC-N: 200-573-9

Warning
H302-H312-H332
P261-P264-P271-P280h-P301+P312a-P304+P340

Ethylenediaminetetraacetic acid tetrasodium salt tetrahydrate > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
pH sol. 1% ..... 10.0 ÷ 12.0
Water (K.F.) ..... 15.0 - 17.0 %
Assay (complexometric) ..... ≥ 98.5 % (s.s.)

Code	Size	Packaging	Notes
405482	250 g	Plastic bottle	
405486	25 kg	Fibreboard box	

Suitable for complexometry

Ethylenediaminetetraacetic acid tripotassic salt

- Acido etilendiamminotetracetico sale tripotassico • Acide éthylènediaminetétracétique sel tripotassique
- Acido etilenodiaminatetracético sal tripotasíca • Ethylendiamintetraessigsäure Tripotassium salt

Synonym:
EDTA Tripotassium salt | Tripotassium ehylenediaminetetraacetate hydrate

$C_{10}H_{13}K_3N_2O_8.2H_2O$
Molecular Weight: 442,57
CAS: 65501-24-8

Ethylenediaminetetraacetic acid tripotassic salt > RE - Pure

RE

Code	Size	Packaging	Notes
405424	25 kg	Fibre drum	

Ethylene dichloride ▶ 1,2-Dichloroethane

Ethylene glycol

- Glicol etilenico • Ethylène glycol • Etilenglicol • Ethylenglycol

Synonym:
1,2-Ethanediol

CH_2OHCH_2OH
Molecular Weight: 62,07
CAS: 107-21-1
EEC-N: 203-473-3

Warning
H302-H373
P260-P264-P301+P312a-P330-P314-P501a

Ethylene glycol > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Water miscibility ..... Conform
Miscb. with Acetone ..... Complete
Alcohol miscibility ..... Complete
Reac.with Ammonium hydr ..... Conform
Sub reducing AgNO3 amm ..... Conform
Density at 20° C ..... 1.108 ÷ 1.118
Refractive index at 20°C ..... 1.4274 ÷ 1.4354
Boiling point ..... 194 ÷ 200 ° C
Water (K.F.) ..... ≤ 0.1 %
Acidity (acetic acid) ..... ≤ 3 ppm
Chloride ..... ≤ 2 ppm
Carbonyl Compounds (CO) ..... ≤ 100 ppm
Heavy metals (Pb) ..... ≤ 2 ppm
Peroxides (H2O2) ..... ≤ 5 ppm
Residue on ignition ..... ≤ 50 ppm
Sulphate ..... ≤ 20 ppm
Fe ..... ≤ 1 ppm
Assay (GLC) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
453905	1 l	Glass bottle	
453904	2.5 l	Glass bottle	
453906	5 l	Plastic tank	
453902	30 kg	Plastic drum	

DASIT GROUP | CARLO ERBA REAGENTS
Product specifications are subject to changes. Please visit our website for updates.

Ethylene glycol > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 1.108 ÷ 1.118
Refractive index at 20°C. 1.4264 ÷ 1.4364
Boiling point..... 194 - 200 ° C
Water (K.F.) ..... ≤0.3 %
Residue on ignition..... ≤100 ppm
Assay (GLC) ..... ≥98 %

Code	Size	Packaging	Notes
346501	1 l	Glass bottle	
346503	2.5 l	Glass bottle	
346502	5 l	Plastic tank	
346504	25 l	Plastic tank	
346509	60 kg	Plastic tank	
346508	230 kg	Metal drum	

Ethylene glycol bis(2-aminoethyl ether)-N,N,N',N'-tetraacetic acid

- Acido etilenglicole bis-(2-amminoetiletere) tetracetico
- Acide éthylèneglycol bis-(2-aminoethylether) tétracétique (EGTA)
- Acido etilenglicol-bis (2-aminoetil-éter) tetracético
- Ethylenglykol-bis-(2-aminoethyl)-n,n,n',n'-tetraessigsäure

Synonym:
EGTA

$C_{14}H_{24}N_2O_{10}$
Molecular Weight: 380,35
CAS: 67-42-5
EEC-N: 200-651-2

Ethylene glycol bis(2-aminoethyl ether)-n,n,n',n'-tetraacetic acid > RPE - For analysis

RPE

Description ..... White cryst. powder
Identification ..... Positive
Loss on drying ..... ≤1 %
Residue on ignition..... ≤0.1 %
Assay (complexometric) ..... ≥97 %

Code	Size	Packaging	Notes
405521	10 g	Glass bottle	
405522	100 g	Plastic bottle	

Suitable for complexometry

Ethylene glycol butyl ether ▶ 2-Butoxy ethanol

Ethylene glycol dimethyl ether ▶ 1,2-Dimethoxyethane

Ethylene glycol monomethyl ether ▶ 2-Methoxy ethanol

Ethylene oxide solution

• Soluzione ossido di etilene • Oxyde d'éthylène solution • Etileno óxido solución • Ethylenoxidlösung

Ethylene oxide solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611036408	1 ml	Glass ampoule	Ethylene oxide solution R2 Ref Ph.Eur 1036408
611036401	10 ml	Glass ampoule	Ethylene oxide stock solution Ref Ph.Eur 1036401

Ethyl ether ▶ Diethyl ether

Ethyl formate

• Etile formiato • Ethyle formiate • Etilo formiato • Ethylformiat

Synonym:
Formic acid ethyl ester

$C_3H_6O_2$
Molecular Weight: 74,08
CAS: 109-94-4
EEC-N: 203-721-0

Classification transport
ONU: 1190
Transport Hazard class: 3
Packing group II

Danger
H225-H302-H332-H319-H335
P210-P280-P303+P361+P353-P304+P340-
P305+P351+P338-P403+P233

Ethyl formate > RE - Pure

RE

Description ..... Clear liquid
Identification ..... Positive
Density at 20° C ..... 0.907 ÷ 0.927
Refractive index at 20°C 1.3547 ÷ 1.3647
Assay (GLC) ..... ≥ 98.0 %
Colour (APHA) ..... ≤ 20
Boiling point ..... 53.3 ÷ 55.3 °C
Acidity ..... ≤ 0.002 meq/g

Code	Size	Packaging	Notes
342101	1 l	Glass bottle	

Ethyl methyl ketone

• Metiletilchetone • Méthyléthylcétone • Metiletilcetona • Ethylmethylketon

Synonym:
• *2-Butanone*
• *MEK*

$CH_3CH_2COCH_3$
Molecular Weight: 72,11
CAS: 78-93-3
EEC-N: 201-159-0

Classification transport
ONU: 1193
Transport Hazard class: 3
Packing group II

Danger
H225-H319-H336-HEU066
P210-P280-P303+P361+P353-P304+P340-
P305+P351+P338-P403+P233

Ethyl methyl ketone > RS - Anhydrous - For analysis

RS

Refractive index at 20°C ..... 1.377 - 1.381
Water content (K.F.) ..... ≤ 200 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 99.5 %
Free acid (as CH3COOH) ..... ≤ 30 mg/Kg

Code	Size	Packaging	Notes
P0201016	1 l	Glass bottle	

Ethyl methyl ketone > RPE - For analysis - Reag. Ph. Eur.

RPE

Description ..... Clear colourless liquid
Identification (I.R.) ..... Positive
Colour ..... ≤ 10 APHA
Alcohol miscibility ..... Complete
Diethyl ether miscib. ..... Complete
Water solubility ..... Conform
Density at 20° C ..... 0.802 ÷ 0.808
Refractive index at 20°C 1.3784 ÷ 1.3834
Boiling point ..... 79 ÷ 80 °C
Water (K.F.) ..... ≤ 500 ppm
Residue on evaporation ..... ≤ 10 ppm
Acidity (acetic acid) ..... ≤ 30 ppm
Aldehydes(Formaldehyde) ..... ≤ 20 ppm
Heavy metals (Pb) ..... ≤ 1 ppm
Subst. reducing KMnO4 ..... ≤ 2 ppm(15m)
Al ..... ≤ 0.5 ppm
Ba ..... ≤ 0.1 ppm
Ca ..... ≤ 0.5 ppm
Cd ..... ≤ 0.05 ppm
Co ..... ≤ 0.02 ppm
Cr ..... ≤ 0.02 ppm
Cu ..... ≤ 0.02 ppm
Fe ..... ≤ 0.1 ppm
Mg ..... ≤ 0.1 ppm
Mn ..... ≤ 0.02 ppm
Ni ..... ≤ 0.02 ppm
Pb ..... ≤ 0.1 ppm
Sn ..... ≤ 0.1 ppm
Zn ..... ≤ 0.1 ppm
Assay (GLC) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
462701	1 l	Glass bottle	
462703	2.5 l	Glass bottle	
462704	10 l	Plastic tank	
462702	22 kg	Metal drum	

Ethyl methyl ketone > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Colour ..... ≤ 10 APHA
Density at 20° C ..... 0.800 ÷ 0.810
Refractive index at 20°C 1.3784 ÷ 1.3844
Boiling point ..... 79 ÷ 80 °C
Water (K.F.) ..... ≤ 0.1 %
Residue on evaporation ..... ≤ 20 ppm
Acidity (acetic acid) ..... ≤ 30 ppm
Total alcohol ..... ≤ 0.5 %
Assay (GLC) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
354254	1 l	Glass bottle	
354253	2.5 l	Glass bottle	
528975	5 l	Plastic tank	
354251	22 kg	Metal drum	
528976	25 l	Metal drum	
528977	200 l	Metal drum	

Product specifications are subject to changes. Please visit our website for updates.

Eugenol

• Eugenolo • Eugenol • Eugenol • Eugenol

Synonym:
• *2-Methoxy-4-(2-propenyl)phenol*
• *4-Allyl-2-methoxyphenol*

$C_{10}H_{12}O_2$
Molecular Weight: 164,21
CAS: 97-53-0
EEC-N: 202-589-1

Warning
H302-H317
P261-P264-P280g-P301+P312a-P333+P313-P501a

Eugenol > RPE - For analysis

RPE

Description ..... Yellow clear liquid
Identification ..... Positive
Density at 20° C ..... 1.05 ÷ 1.07
Refractive index at 20°C ..... 1.53 ÷ 1.54
Assay (GLC) ..... ≥ 96 %

Code	Size	Packaging	Notes
449773	100 ml	Glass bottle	

Eukitt

• Eukitt • EUKITT • Eukitt • EUKITT

Classification transport
ONU: 1307
Transport Hazard class: 3
Packing group III

Danger
H226-H315-H319-H335-H304
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Eukitt > RS - For microscopy

RS

Description ..... Viscous liquid
Identification ..... Positive

Code	Size	Packaging	Notes
554194	100 ml	Aluminium bottle	
554193	250 ml	Aluminium bottle	
554192	500 ml	Aluminium bottle	

Balm rapid inclusion

Europium standard solution

• Europio standard soluzione • Europium solution standard • Europio, solución patrón • Europium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Europium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505602	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505605	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Europium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503571	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503575	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503573	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503577	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Europium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507736	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507502	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Fast green FCF

• Verde solido FCF • Vert solide FCF • Verde sólido FCF • Schneller grüner FCF

Synonym:
Food green 3

$C_{37}H_{34}N_2Na_2O_{10}S_3$
Molecular Weight: 808,86
CAS: 2353-45-9
EEC-N: 219-091-5



Warning
H302-H312-H332-H351
P261-P271-P280-P304+P340-P308+P313-P330

Fast green FCF > RS - For microscopy - C.I. 42053

RS

Description ..... Red brown powder Identification ..... Positive

Code	Size	Packaging	Notes
491391	25 g	Glass bottle	

Dye for histology

Fehling's A reagent

• Fehling reattivo soluzione A • Réactif de Fehling solution A • Fehling reactivo solución A • Fehling-Reagenzlösung A

Classification transport
ONU: 3082
Transport Hazard class: 9
Packing group III



Danger
H318-H410
P273-P280i-P305+P351+P338-P310a-P391-P501a

Fehling's A reagent > RS - For glucose detection

RS

CuSO4.5H2O content ..... 69.16 - 69.43 g/l Description ..... Clear blue liquid Density at 20° C ..... 1.037 - 1.043 [CuSO4.5H2O] ..... 69.12 - 69.40 g/l

Code	Size	Packaging	Notes
449926	500 ml	Plastic bottle	
449927	1 l	Plastic bottle	
PS0492/41	10 l	Plastic tank	

Composition: Copper sulfate and benzoic acid

Fehling's B reagent

• Fehling reattivo soluzione B • Réactif de Fehling solution B • Fehling reactivo solución B • Fehling-Reagenzlösung B

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II



Danger
H314
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Fehling's B reagent > RS - For glucose detection

RS

Concentration ..... 345.3 - 346.7 g/l Description ..... Clear colourless liquid Density at 20° C ..... 1.259 - 1.265 Alcalinity (NaOH) ..... 119 - 121 0/00

Code	Size	Packaging	Notes
E449936	500 ml	Plastic bottle	
E449937	1 l	Plastic bottle	
PS0493/41	10 l	Plastic tank	

Composition: Potassium sodium tartrate and sodium hydroxide

Ferricyanide standard solution

• Ferricianuro standard soluzione • Ferricyanures standard solution • Ferricianuro, solución patrón • Ferricyanid-Standardlösung

Ferricyanide standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615001300	100 ml	Plastic bottle	A 50 ppm solution: to dilute according to Ph.Eur 5001300

Ferrocyanide standard solution

• Ferrocianuro standard soluzione • Ferrocyanures standard solution • Ferrocianuro, solución patrón • Ferrocyanid-Standardlösung

Ferrocyanide standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615001209	100 ml	Plastic bottle	A 100 ppm solution: to dilute according to Ph.Eur 5001200

Ferroin 0.025 mol/l solution

• Ferroina 0.025 mol/l soluzione • Indicateur Ferroïne 0.025 mol/l • Ferroína solución 0.025 mol/l
• Indikator Ferroin 0.025 mol/l

Synonym:
• *1,10-Phenanthroline iron(II) sulfate complex*
• *o-Phenanthroline ferrous sulfate complex*

$[Fe(C_{12}H_8N_2)_3]SO_4$
Molecular Weight: 692,52
CAS: 14634-91-4

H412
P273-P501a

Ferroin 0.025 mol/l solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611038100	100 ml	Plastic bottle	Ref Ph.Eur 1038100

Ferroin 0.025 mol/l solution > RS - For environmental analysis (COD determination)

RS

Description ..... Dark red liquid

Code	Size	Packaging	Notes
526751	100 ml	Bottle	

Ferron

• Ferron • Ferron • Ferron • Ferron

$IC_6H(OH)SO_3HN{:}CHCH{:}CH$
Molecular Weight: 351,12
CAS: 547-91-1
EEC-N: 208-938-4

Classification transport
ONU: 2585
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Ferron > RPE - For analysis

RPE

Description .....Yellow crystalline powder
Identification ..... Positive
Loss on drying ..... ≤0.5 %
Heavy metals (Pb).....≤20 ppm
Residue on ignition..... ≤0.1 %
Iron sensitivity .....≥5 µg/ml
Fe .....≤20 ppm

Code	Size	Packaging	Notes
406918	5 g	Glass bottle	

For extraction and spectrophotometric determination of Mo (VI), Pd (II), U (VI), V (III)

Fixative AFA liquid

• Liquido fissatore AFA • Liquide fixateur AFA • Líquido fijador AFA • AFA-Fixiererflüssigkeit

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II

Danger
H225-H319-H317-H350-HA26
P210-P241-P261-P280-P303+P361+P353-P305+P351+P338

Fixative AFA liquid > RS - For histology

RS

Description ..... Clear liquid

Code	Size	Packaging	Notes
508840	480 x 30 ml	Plastic bottle	60 ml jars filled at 30 ml. Box of 500
526267	1 l	Plastic bottle	In Vitro Diagnostic Medical Device
526263001	5 l	Plastic tank	In Vitro Diagnostic Medical Device

Contains Ethanol, formaldehyde and acetic acid

Product specifications are subject to changes. Please visit our website for updates.

Fixative Bouin Hollande liquid

• Liquido fissatore di Bouin Hollande • Liquide de Bouin Hollande • Líquido fijador de Bouin Holland • Alkohol Bouin Holland



Danger
H317-H350-HA26
P261-P280-P308+P313-P362+P364-P333+P313-P501a

Fixative Bouin Hollande liquid > RS - For histology

RS

Appearance ..... Clear liquid

Code	Size	Packaging	Notes
526268	1 l	Plastic bottle	

Fixative Bouin liquid

• Liquido fissatore di Bouin • Liquide de Bouin • Líquido fijador de Bouin • Bouins Pikroformol Wasserlisch



Danger
H315-H319-H317-H341-H350-H335-HA26
P261-P271-P280-P304+P340-P305+P351+P338-P403+P233

Fixative Bouin liquid > RS - For histology

RS

Appearance ..... Yellow liquid

Code	Size	Packaging	Notes
526270	1 l	Plastic bottle	In Vitro Diagnostic Medical Device
526261	5 l	Plastic tank	In Vitro Diagnostic Medical Device
526311	25 l	Plastic tank	In Vitro Diagnostic Medical Device

Contains formaldehyde, acetic acid and 2.4.6-trinitrophenol

Fixative Bouin liquid > RS - For hematology

RS

Description ..... Yellow clear liquid Identification ..... Positive

Code	Size	Packaging	Notes
429751	1 l	Plastic bottle	In Vitro Diagnostic Medical Device

Fixative Davidson liquid

• Liquido fissatore di Davidson • Liquide fixateur de Davidson • Líquido fijador de Davidson • Davidson Fixiererflüssigkeit

Classification transport
ONU: 2924
Transport Hazard class: 3
Packing group II



Danger
H225-H312-H332-H314-H317-H341-H350-H335-HA26
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P362+P364-P403+P233

Fixative Davidson liquid > RS - For histology

RS

Description ..... Clear colourless liquid

Code	Size	Packaging	Notes
508881	30 ml	Plastic bottle	60 ml jars filled at 30 ml. Box of 500
526277	5 l	Plastic tank	

Fixative FIXALL-HIS liquid

• Liquido fissatore FIXALL-HIS • Liquide fixateur FIXALL-HIS • Líquido fijador FIXALL-HIS • Fixierflüssigkeit FIXALL-HIS

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II



Danger
H225-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

Fixative FIXALL-HIS liquid > RS - For histology

RS

Code	Size	Packaging	Notes
526274	5 l	Plastic tank	

Formaldehyde substitute. Ready-to-use solution

Fixative liquid without acetic acid

• Liquido fissatore senza acido acetico • Liquide fixateur sans acide acétique • Liquido fijador sin acido acético • Fixierflüssigkeit ohne Essigsäure

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group III



Danger
H226-H312-H332-H315-H319-H317-H341-H350-H335-HA26
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Fixative liquid without acetic acid > RS - For histology

RS

Description ..... Clear colourless liquid

Code	Size	Packaging	Notes
526264	10 l	Plastic tank	

Florisil 60-100 mesh

• Florisil 60-100 mesh • Florisil 60-100 mesh • Florisil 60-100 mesh • Florisil 60-100 mesh

Synonym:
Magnesium silicate

MgO_3Si
Molecular Weight: 100,39
CAS: 1343-88-0
EEC-N: 215-681-1

Florisil 60-100 mesh > RS - Adsorbent for chromatography

RS

Description ..... White powder Identification ..... Positive

Code	Size	Packaging	Notes
452331	100 g	Plastic bottle	
452333	500 g	Plastic bottle	
452332	1 kg	Plastic bottle	

Florisil 60-100 mesh > RS - For residual pesticides analysis

RS

Description ..... White powder
Identification ..... Positive

Granulometry
60 mesh..... ca 2.0 %
70 mesh..... ca 20.0 %
80 mesh..... ca 41.0 %
100 mesh..... ca 30.0 %
140 mesh..... ca 6.9 %

Code	Size	Packaging	Notes
452271	100 g	Plastic bottle	
452273	500 g	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Florisil 100-200 mesh

• Florisil 100-200 mesh • Florisil 100-200 mesh • Florisil 100-200 mesh • Florisil 100-200 mesh

Synonym:
Magnesium silicate

MgO_3Si
Molecular Weight: 100,39
CAS: 1343-88-0
EEC-N: 215-681-1

Florisil 100-200 mesh > RS - Adsorbent for chromatography

RS

Description ..... White powder Identification ..... Positive

Code	Size	Packaging	Notes
452351	100 g	Plastic bottle	
452353	500 g	Plastic bottle	

Fluorescein

• Fluoresceina • Fluorescéine • Fluoresceína • Fluorescein

Synonym:
Acid Yellow 73

$C_{20}H_{12}O_5$
Molecular Weight: 332,32
CAS: 2321-07-5
EEC-N: 219-031-8

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Fluorescein > RPE - For analysis

RPE

Description ..... Red brown powder Identification ..... Positive

Code	Size	Packaging	Notes
452086	25 g	Glass bottle	
452083	50 g	Glass bottle	
452087	500 g	Plastic bottle	

Fluorescein sodium salt

• Fluoresceina sodica • Fluorescéine sodique • Fluoresceína sódica • Fluorescein Natriumsalz

Synonym:
• *Acid Yellow 73*
• *Uranine*

$C_{20}H_{10}Na_2O_5$
Molecular Weight: 376,28
CAS: 518-47-8
EEC-N: 208-253-0

Fluorescein sodium salt > RPE - For analysis - C.I. 45350

RPE

Description ..... Red brick powder Identification ..... Positive Absorbtion ind.sensit. ..... Conform Loss on drying ..... ≤10 %

Code	Size	Packaging	Notes
452112	25 g	Glass bottle	
452113	50 g	Plastic bottle	
452117	1 kg	Plastic bottle	

Fluorescein sodium salt > RE - Pure - C.I. 45350

RE

Description ..... Red brick powder Identification ..... Positive

Code	Size	Packaging	Notes
345356	25 g	Glass bottle	
345357	1 kg	Plastic bottle	
345355	5 kg	Plastic tank	

Fluoride standard solution

• Fluoruri standard soluzione • Fluorure solution standard • Fluoruro, solución patrón • Fluorid-Standardlösung

Fluoride standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615001401	100 ml	Plastic bottle	A 1 ppm solution: to dilute according to Ref Ph.Eur 5001401
615001409	100 ml	Plastic bottle	A 10 ppm solution: to dilute according to Ref Ph.Eur 5001400

Fluoride standard solution > RS - Standard solution for ion chromatography

RS

Code	Size	Packaging	Notes
503251	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503253	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Fluorite ▶ Calcium fluoride

Folin-Ciocalteu's reagent

• Folin-Ciocalteu reattivo • Réactif de Folin-Ciocalteu • Folin-Ciocalteu reactivo • Folin-Ciocalteu-Reagenz

Classification transport
ONU: 3264
Transport Hazard class: LQ



Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Folin-Ciocalteu's reagent > RS - For microscopy

RS

Description ..... Yellow clear liquid
Identification ..... Positive
Sensibilità ai fenoli (A a 650nm) ..... ≥ 0.26
Titolo (equiv. di acido) .....1.9 ÷ 2.1 N

Code	Size	Packaging	Notes
E463562	500 ml	Glass bottle	

For the determination of phenols. Store at +4 °C

Formaldehyde 37% w/v

• Aldeide formica 37% m/v • Aldéhyde formique 37% m/v • Formaldehído 37% p/v • Formaldehyd 37% w/v

Synonym:
Formalin

HCHO
Molecular Weight: 30,03
CAS: 50-00-0

Classification transport
ONU: 2209
Transport Hazard class: 8
Packing group III



Danger
H301-H311-H331-H314-H317-H341-H350-H370-H335-HA26
P280-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Formaldehyde 37% w/v > RPE - For analysis

RPE

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Density at 20°C ..... 1.085 ÷ 1.092
Acidity (formic acid) ..... ≤350 ppm
Chloride ..... ≤5 ppm
Heavy metals (Pb) ..... ≤2 ppm
Residue on ignition ..... ≤500 ppm
Sulphate ..... ≤20 ppm
Fe ..... ≤1 ppm
Pb ..... ≤1 ppm
Assay (oxidimetric) ..... ≥37 % (m/v)

Code	Size	Packaging	Notes
415661	1 l	Plastic bottle	
415666	2.5 l	Plastic bottle	
415667	5 l	Plastic bottle	
415669	30 kg	Plastic drum	
415662	55 kg	Plastic tank	

Stabilized with ~10% of methanol

Product specifications are subject to changes. Please visit our website for updates.

Formaldehyde 37% w/v neutralized

• Aldeide formica 37% p/v neutralizzata • Aldéhyde formique 37% m/v neutre
• Formaldehído 37% p/v neutralizata • Formaldehyd 37% m / v neutral

Synonym:
Formalin

HCHO
Molecular Weight: 30,03
CAS: 50-00-0

Classification transport
ONU: 2209
Transport Hazard class: 8
Packing group III



Danger
H301-H311-H331-H314-H317-H341-H350-H370-H335-HA26
P280-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Formaldehyde 37% w/v neutralized > RS - For histology

RS

Description ..... Clear liquid | Colour (APHA) ..... ≤ 10 | Identification ..... Positive | Assay (oxidimetric) ..... ≥ 37 % (m/v)

Code	Size	Packaging	Notes
415686	1 l	Plastic bottle	
415682	5 l	Plastic bottle	
415683	10 kg	Plastic tank	
415684	30 kg	Plastic drum	
415685	55 kg	Plastic tank	

Stabilized with ~10% of methanol. Neutralized with dolomite

Formaldehyde 35% w/w

• Aldeide formica 35% m/m • Aldéhyde formique 35% m/m • Formaldehyde 35% p/p
• Formaldehyd 35% m/m

Synonym:
Formalin

HCHO
Molecular Weight: 30,03
CAS: 50-00-0

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group III



Danger
H302-H311-H331-H314-H317-H341-H350-H335-HA26
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P308+P313-P361+P364-P403+P233

Formaldehyde 35% w/w > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611039101	100 ml	Plastic bottle	Ref Ph.Eur 1039101

Formaldehyde 35% w/w > ERBApharm - According to pharmacopoeia: BP-FU-Ph.Eur.

ERBApharm

Description ..... Clear colourless liquid | Appearance of solution ..... Conform Ph.Eur. | Sulphated ash ..... ≤0.1 % | Assay (iodometric) ..... 34.5 ÷ 38.0 % m/m
Identification ..... Positive | Acidity ..... Conform Ph.Eur. | Methyl alcohol ..... 9.0 ÷ 15.0 % v/v | Origin (BSE/TSE) ..... Synthesis

Code	Size	Packaging	Notes
310351	1 l	Plastic bottle	
310356	2.5 l	Plastic bottle	
310358	5 l	Plastic bottle	
310349	10 kg	Plastic tank	
310348	30 kg	Plastic drum	
310355	55 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Formaldehyde 30% w/v

• Aldeide formica 30 % m/v • Aldéhyde formique 30% m/v • Formaldehído 30% p/v • Formaldehyd 30% m/v

Synonym: *Formalin*

HCHO
Molecular Weight: 30,03
CAS: 50-00-0

Classification transport
ONU: 2209
Transport Hazard class: 8
Packing group III



Danger
H302-H311-H331-H314-H317-H341-H350-H371-H335-HA26
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P361+P364-P403+P233

Formaldehyde 30% w/v > RE - Pure

RE

Description ..... Clear colourless liquid Density at 20°C ..... 1.075 ÷ 1.100 Assay ..... 27.5 ÷ 32.5 % (m/v)

Code	Size	Packaging	Notes
524930	5 l	Plastic tank	

Stabilized with ~10% of methanol

Formaldehyde 10% v/v according to Lillie

• Aldeide formica 10% v/v (Liquido di Lille) • Aldéhyde formique 10% v/v selon Lillie
• Formaldehído 10% v/v según Lillie • Formaldehyd 10% v/v Lillie

Synonym: *Formalin*

HCHO
Molecular Weight: 30,03
CAS: 50-00-0



Danger
H317-H341-H350-HA26
P261-P280-P308+P313-P362+P364-P333+P313-P501a

Formaldehyde 10% v/v according to Lillie > RS - For histology

RS

Description ..... Clear colourless liquid pH at 20°C ..... 6.8 ÷ 7.0

Code	Size	Packaging	Notes
526912	5 l	Plastic tank	
526911	25 l	Plastic tank	

Stabilized with ~10% of methanol. Buffered at pH 7 with phosphate salts

Formaldehyde 5% w/v buffered at pH 6.9

• Aldeide formica 5% m/v tamponata pH 6.9 • Aldéhyde formique 5% m/v tamponné à pH 6.9
• Formaldehído 5% p/v tamponata pH 6.9 • Formaldehyd 5%m/v pH 6.9 gepuffert

Synonym: *Formalin*

HCHO
Molecular Weight: 30,03
CAS: 50-00-0



Danger
H317-H341-H350-HA26
P261-P280-P308+P313-P362+P364-P333+P313-P501a

Formaldehyde 5% w/v buffered at pH 6.9 > RS - For histology

RS

Description ..... Clear colourless liquid pH at 20°C ..... 6.9 - 7.1

Code	Size	Packaging	Notes
415674	5 l	Plastic tank	
415672	10 l	Plastic tank	

Buffered at pH 7 with phosphate salts

Product specifications are subject to changes. Please visit our website for updates.

Formaldehyde 4% w/v buffered at pH 6.9

• Aldeide formica 4% m/v tamponata pH 6.9 • Aldéhyde formique 4% m/v tamponné à pH 6.9
• Formaldehído 4% p/v tamponata pH 6.9 • Formaldehyd 4%m/v pH 6.9 gepuffert

Synonym:
Formalin

HCHO
Molecular Weight: 30,03
CAS: 50-00-0



Danger
H317-H341-H350-HA26
P261-P280-P308+P313-P362+P364-P333+P313-P501a

Formaldehyde 4% w/v buffered at pH 6.9 > RS - For histology - CE - IVD

RS

Code	Size	Packaging	Notes
415634	1 l	Plastic bottle	In Vitro Diagnostic Medical Device
415631	5 l	Plastic bottle	In Vitro Diagnostic Medical Device
415633	10 l	Kubidos	In Vitro Diagnostic Medical Device
415636	20 l	Plastic tank	In Vitro Diagnostic Medical Device

For specifications, contact our customer service for a certificate of analysis

Formaldehyde 4% w/v buffered at pH 6.9 > RS - For histology

RS

Description ..... Clear colourless liquid
pH at 20°C ..... 6.8 ÷ 7.0
Identification ..... Positive
Density at 20°C ..... ≥ 1.00
pH at 20° C ..... 6.8 ÷ 7.0
Assay (oxydimetric) ..... ≥ 4.0 %
Density at 20° C ..... ≥ 1.00
Assay (oxidimetric) ..... 3.9 - 4.1 %

Code	Size	Packaging	Notes
508861	30 ml	Plastic jar	60 ml jars filled at 30 ml. Box of 500
508863	300 ml	Plastic jar	500 ml jars filled at 300 ml. Box of 32
524920	500 ml	Plastic bottle	
526937	800 ml	Plastic jar	1 l jar filled at 800 ml
415694	1 l	Plastic bottle	
526931	2.5 l	Metal bucket	5 l bucket filled at 2.5 l
415691	5 l	Plastic bottle	
415695	5 l	Kubidos	
526936	5 l	Plastic tank	
415693	10 l	Kubidos	
526933	10 l	Plastic tank	
415696	20 l	Kubidos	
415692	30 kg	Plastic drum	

10% solution (v/v) buffered at pH 6.9 with phosphate ions. Stabilized with methanol

Formaldehyde 4% w/v with sodium chloride

• Aldeide formica 4% m/v con sodio cloruro • Aldéhyde formique 4% m/v salé
• Formaldehído 4% p/v con sodio cloruro • Formaldehyd 4% m/v mit natriumchlorid

Synonym:
Formalin

HCHO
Molecular Weight: 30,03
CAS: 50-00-0



Danger
H317-H341-H350-HA26
P261-P280-P308+P313-P362+P364-P333+P313-P501a

Formaldehyde 4% w/v with sodium chloride > RS - For histology

RS

Description ..... Clear colourless liquid
pH at 20°C ..... 7.30 ÷ 7.40
Stabilized with 1% of methyl alcohol

Code	Size	Packaging	Notes
526934	1 l	Plastic bottle	

Formaldehyde acetic

• Aldeide formica e acido acetico • Formol acétique • Formol acético • Formalin Essigsäure

Synonym: *Formalin*

HCHO
Molecular Weight: 30,03
CAS: 50-00-0



Danger
H315-H319-H317-H341-H350-H335-HA26
P261-P271-P280-P304+P340-P305+P351+P338-P403+P233

Formaldehyde acetic > RS - For histology

RS

Description ..... Clear colourless liquid

Code	Size	Packaging	Notes
508871	480 x 30 ml	Plastic jar	60 ml jars filled at 30 ml. Box of 500
526231	1 l	Plastic bottle	
526273	5 l	Plastic tank	

Formamide

• Formammide • Formamide • Formamida • Formamid

Synonym: *Formic amide*

$HCONH_2$
Molecular Weight: 45,02
CAS: 75-12-7
EEC-N: 200-842-0



Danger
H351-H360D-H373-HA26
P260-P280-P308+P313-P314-P405-P501a

Formamide > RS - Anhydrous - For analysis

RS

Refractive index at 20°C ..... 1.445 - 1.449
Water content (K.F.) ..... ≤ 500 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 99 %
Methanol ..... ≤ 0.2 %

Code	Size	Packaging	Notes
P6151010	200 ml	Bottle with septum	

Formamide > RPE - For analysis - ACS

RPE

Description ..... Clear liquid
Colour (APHA) ..... ≤ 10
Identification ..... Positive
Freezing point ..... 2.0 ÷ 3.0 °C
Assay (GLC) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
452282	250 ml	Glass bottle	
452286	1 l	Glass bottle	

Formic acid 99%

• Acido formico 99% • Acide formique 99% • Acido fórmico 99% • Ameisensäure 99%

HCOOH
Molecular Weight: 46,03
CAS: 64-18-6
EEC-N: 200-579-1

Classification transport
ONU: 1779
Transport Hazard class: 8
Packing group II



Danger
H226-H302-H331-H314-HEU071
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Formic acid 99% > RS - For LC/MS

RS

Description ..... Clear colourless liquid
Colour (APHA) ..... ≤10
Refractive index at 20°C 1.3709 ÷ 1.3719
Density at 20° C ..... 1.218 ÷ 1.222
Boiling point ..... 100.2 - 101.2 °C
Residue on evaporation ..... ≤10 ppm
Assay (acidimetric) ..... ≥99.0 %
U.V. Transmittance
at 260 nm ..... ≥ 20 %
at 270 nm ..... ≥ 85 %
at 280 nm ..... ≥ 95 %
at 300 nm ..... ≥ 98 %
at 320 nm ..... ≥ 99 %
Al ..... ≤ 0.05 ppm
Fe ..... ≤ 0.2 ppm
Na ..... ≤ 0.5 ppm
Ca ..... ≤ 0.2 ppm
Mg ..... ≤ 0.1 ppm
K ..... ≤ 0.1 ppm

Code	Size	Packaging	Notes
405821	10 x 1 ml	Glass ampoule	
405822	10 x 2.5 ml	Glass ampoule	
405823	50 ml	Plastic bottle	
405824	1 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Formic acid 99% > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description	Clear liquid	Chloride	≤5 ppm	Cd	≤0.05 ppm	Na	≤0.5 ppm
Colour (APHA)	≤10	Heavy metals (Pb)	≤2 ppm	Co	≤0.02 ppm	Ni	≤0.2 ppm
Water miscibility	Conform	Oxalate	≤50 ppm	Cr	≤0.05 ppm	Pb	≤0.02 ppm
Density at 20° C	1.218 ÷ 1.222	Sulphate	≤5 ppm	Cu	≤0.02 ppm	Sr	≤0.02 ppm
Refractive index at 20°C	1.3709 ÷ 1.3719	Sulphite	≤10 ppm	Fe	≤2 ppm	V	≤0.05 ppm
Boiling point	100.2 ÷ 101.2 °C	Ag	≤0.02 ppm	K	≤0.1 ppm	Zn	≤0.05 ppm
Residue on evaporation	≤20 ppm	Al	≤0.05 ppm	Li	≤0.02 ppm	Assay (acidimetric)	≥98 %
Acetic acid	≤500 ppm	Ba	≤0.05 ppm	Mg	≤0.5 ppm		
Ammonium	≤10 ppm	Bi	≤0.1 ppm	Mn	≤0.05 ppm		
Total nitrogen	≤20 ppm	Ca	≤0.2 ppm	Mo	≤0.02 ppm		

Code	Size	Packaging	Notes
405792	1 l	Glass bottle	
405793	5 l	Plastic tank	
405794	30 kg	Plastic drum	

Formic acid 99% > ERBApharm - According to pharmacopoeia: DAB

ERBApharm

Description	Clear colourless liquid	Density at 20° C	1.218 ÷ 1.222	Origin (BSE/TSE)	Synthesis
Identification	Positive	Assay (acidimetric)	≥ 98 % m/m		

Code	Size	Packaging	Notes
303911	1 l	Glass bottle	
303912	2.5 l	Glass bottle	
303913	30 kg	Plastic drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Formic acid 85%

• Acido formico 85% • Acide formique 85% • Acido fórmico 85% • Ameisensäure 85%

HCOOH
Molecular Weight: 46,03
CAS: 64-18-6

Classification transport
ONU: 1779
Transport Hazard class: 8
Packing group II

Danger
H226-H302-H331-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Formic acid 85% > RPE - For analysis

RPE

Description	Clear liquid	Total nitrogen	≤20 ppm	Bi	≤0.1 ppm	Mg	≤0.5 ppm
Colour (APHA)	≤10	Chloride	≤5 ppm	Cd	≤0.05 ppm	Mn	≤0.05 ppm
Identification	Positive	Heavy metals (Pb)	≤2 ppm	Co	≤0.02 ppm	Mo	≤0.02 ppm
Water miscibility	Conform	Oxalate	≤50 ppm	Cr	≤0.05 ppm	Ni	≤0.2 ppm
Density at 20° C	1.196 ÷ 1.199	Sulphate	≤5 ppm	Cu	≤0.02 ppm	Sr	≤0.02 ppm
Residue on evaporation	≤20 ppm	Sulphite	≤10 ppm	Fe	≤2 ppm	V	≤0.05 ppm
Acidity (acetic acid)	≤500 ppm	Ag	≤0.02 ppm	K	≤0.1 ppm	Zn	≤0.05 ppm
Ammonium	≤10 ppm	Ba	≤0.05 ppm	Li	≤0.02 ppm	Assay (acidimetric)	85 ÷ 87 %

Code	Size	Packaging	Notes
405832	1 l	Glass bottle	
405833	2.5 l	Glass bottle	
405835	30 kg	Plastic drum	

Formic acid 85% > RE - Pure

RE

Description	Clear liquid	Residue on evaporation	≤0.5 %	Sulphate	≤300 ppm
Identification	Positive	Chloride	≤100 ppm	Fe	≤50 ppm
Density at 20° C	1.196 ÷ 1.199	Heavy metals (Pb)	≤50 ppm	Assay (acidimetric)	85 ÷ 87 %

Code	Size	Packaging	Notes
303905	1 l	Glass bottle	
303901	30 kg	Plastic drum	

Formic acid-d

• Acido formico-d • Acide formique-d • Acido fórmico-d • Ameisensäure-d

HCOOD
Molecular Weight: 47,03
CAS: 925-94-0

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group I

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Formic acid-d > RS - For NMR - min 97%

RS

Code	Size	Packaging	Notes
P5733	5 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Formic acid ammonium salt ▶ Ammonium formate

Formic acid ethyl ester ▶ Ethyl formate

D(-)Fructose

• D(-)Fruttosio • D(-)Fructose • D(-)Fructosa • Lävulose

$C_6H_{12}O_6$
Molecular Weight: 180,16
CAS: 57-48-7
EEC-N: 200-333-3

D(-)Fructose > RPE - For analysis

RPE

Description ..... White powder
Identification ..... Positive
Melting point..... 101.5 ÷ 104.5 ° C
Specific optical rotation..... -93.0 ÷ -91.0 °
Water (K.F.) ..... ≤ 0.5 %
Chloride..... ≤40 ppm
Heavy metals (Pb)..... ≤10 ppm
Sulphate ..... ≤50 ppm
As ..... ≤1 ppm
Assay ..... ≥ 99.5 %

Code	Size	Packaging	Notes
452665	100 g	Plastic bottle	
452666	500 g	Plastic bottle	

Fuchsin acid

• Fucsina acida • Fuchsine acide • Fucsina ácida • Säurefuchsin

Synonym:
Acid Violet 19

$C_{20}H_{17}N_3Na_2O_9S_3$
Molecular Weight: 585,6
CAS: 3244-88-0
EEC-N: 221-816-5

Fuchsin acid > RPE - For analysis - C.I. 42685

RPE

Description ..... Dark green crystals
Identification ..... Positive
Decolorization with SO2..... Conform

Code	Size	Packaging	Notes
452812	25 g	Glass bottle	
452814	100 g	Plastic bottle	

Dye for microscopy (botanical-histology). Indicator acid - base (pH 12.0 ÷ 14.0)

Product specifications are subject to changes. Please visit our website for updates.

Fuchsin basic

• Fucsina basica • Fuchsine basique • Fucsina básica • Fuchsin

Synonym:
Basic Violet 14

$C_{20}H_{20}ClN_3$
Molecular Weight: 337,85
CAS: 632-99-5
EEC-N: 211-189-6



Warning
H315-H319-H351-H335
P261-P271-P280-P304+P340-P305+P351+P338-P403+P233

Fuchsin basic > RPE - For analysis - C.I. 42510

RPE

Description ..... Green crystals Identification ..... Positive Decolorization with SO2 ..... Conform

Code	Size	Packaging	Notes
452842	25 g	Glass bottle	
452844	100 g	Plastic bottle	

Dye for microscopy (bacteriology-Botanical-histology). Indicator acid - base (pH 1.0 ÷ 3.1)

Fuchsin solution decolorised

• Fucsina decolorata soluzione • Fuschine décolorée solution • Fucsina solución decolorizada
• Fuschine verfärbte Lösung

Synonym:
Basic Violet 14

$C_{20}H_{20}ClN_3$
Molecular Weight: 337,85
CAS: 632-99-5



Warning
H290
P234-P390-P406

Fuchsin solution decolorised > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611039401	100 ml	Glass bottle	Ref Ph.Eur 1039401
611039402	100 ml	Glass bottle	Fuchsin solution, decolorised R1 Ref Ph.Eur 1039402

Fumaric acid

• Acido fumarico • Acide fumarique • Acido fumárico • Fumarsäure

HOOCCH:CHCOOH
Molecular Weight: 116,07
CAS: 110-17-8
EEC-N: 203-743-0



Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Fumaric acid > RPE - For analysis

RPE

Description ..... White crystal powder Water (K.F.) ..... ≤ 0.5 % Heavy metals (Pb) ..... ≤ 10 ppm Maleic acid ..... ≤ 0.1 %
Identification ..... Positive Assay (dried base) ..... 99.5 ÷ 100.5 % Residue on ignition ..... ≤ 100 ppm

Code	Size	Packaging	Notes
406284	100 g	Glass bottle	
406287	1 kg	Plastic bottle	

Gadolinium standard solution

• Gadolinio standard soluzione • Gadolinium solution standard • Gadolinio, solución patrón • Gadolinium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III



Gadolinium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505622	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505625	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Gadolinium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503601	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503603	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503605	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503607	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Gadolinium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507737	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507504	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

D(+)Galactose

• D(+)Galattosio • D(+)Galactose • D(+)Galactosa • D(+)Galactose

$CH_2OHCH(CHOH)_3CHOH$
Molecular Weight: 180,16
CAS: 59-23-4
EEC-N: 200-416-4

D(+)Galactose > RPE - For analysis

RPE

Description ..... White powder
Identification ..... Positive
Potere rotat. spec. a20°C (C=10;H2O;NH3) ..+78 ÷ +81.5 ° (s.s.)
Water (K.F.) .....≤ 0.3 %
Sulphated ash .....≤ 0.1 %

Code	Size	Packaging	Notes
453125	250 g	Plastic bottle	
453126	1 kg	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Gallic acid monohydrate

• Acido gallico monoidrato • Acide gallique monohydraté • Acido gálico monohidrato
• Gallussäure monohydrat

Synonym:
3,4,5-Trihydroxybenzoic acid monohydrate

$3,4,5\text{-}(OH)_3C_6H_2COOH.H_2O$
Molecular Weight: 188,14
CAS: 5995-86-8
EEC-N: 205-749-9

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-
P332+P313-P403+P233

Gallic acid monohydrate > RPE - For analysis

RPE

Description ..... White powder
Identification ..... Positive
Loss on drying ..... ≤ 10 %
Sulphate ..... ≤200 ppm
Sulphated ash ..... ≤ 0.1 %
Assay (acidimetric) ..... ≥ 99.0 % (s.s.)

Code	Size	Packaging	Notes
406335	250 g	Plastic bottle	
406336	1 kg	Plastic bottle	

Gallic acid monohydrate > RE - Pure

RE

Description ..... Yellow crystalline powder
Identification ..... Positive
Loss on drying 100° C ..... ≤10 %
Residue on ignition ..... ≤0.1 %
Sulphate ..... ≤500 ppm

Code	Size	Packaging	Notes
304201	1 kg	Plastic bottle	

Gallium standard solution

• Gallio standard soluzione • Gallium standard solution • Galio solución estándar • Gallium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Gallium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505617	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505618	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505619	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Gallium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507739	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507503	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Gardner Colour Standards

• Standard Gardner del colore • Etalons couleurs Gardner • Patrones de color Gardner • Gardner Farbstandards

Gardner Colour Standards > RS - For calibration

RS

Code	Size	Packaging	Notes
540701	100 ml	Glass bottle	Colour 2
540702	100 ml	Glass bottle	Colour 4
540703	100 ml	Glass bottle	Colour 6
540704	100 ml	Glass bottle	Colour 8
540705	100 ml	Glass bottle	Colour 10
540706	100 ml	Glass bottle	Colour 12
540707	100 ml	Glass bottle	Colour 14
540708	100 ml	Glass bottle	Colour 16

Gelatine

• Gelatina • Gélatine • Gelatina • Gelatine

CAS: 9000-70-8
EEC-N: 232-554-6

Gelatine > RS - For microbiology

RS

Description ..... Yellowish crystalline powder Identification ..... Positive Loss on drying ..... ≤ 13 % Sulphated ash ..... ≤ 2 %

Code	Size	Packaging	Notes
453226	500 g	Plastic bottle	

Gentian violet

• Violetto genziana • Violet de gentiane • Violeta de genciana • Enzianviolett

Synonym:
• *Crystal violet solution*
• *Basic violet 3*

$C_{25}H_{30}ClN_3$
Molecular Weight: 407,99
CAS: 548-62-9
EEC-N: 208-953-6



Danger
H302-H318-H351-H410
P264-P280-P301+P312a-P305+P351+P338-P310a-P308+P313

Gentian violet > RE - Pure

RE

Description ..... Dark green powder Identification ..... Positive Loss on drying 100° C. ..... ≤10 %

Code	Size	Packaging	Notes
388703	25 g	Glass bottle	
388701	50 g	Glass bottle	

Dye for microscopy (bacteriology)

Gentian violet carbolated solution

• Violetto genziana soluzione fenicata • Violet de gentiane solution phéniquée
• Violeta de genciana solución fenicada • Enzianviolettkarbollösung

Synonym:
• *Crystal violet solution*
• *Basic violet 3*

$C_{25}H_{30}ClN_3$
Molecular Weight: 407,99
CAS: 548-62-9



Warning
H319-H412
P264-P273-P280i-P305+P351+P338-P337+P313-P501a

Gentian violet carbolated solution > RS - For microscopy

RS

Description ..... Violet clear liquid Identification ..... Positive

Code	Size	Packaging	Notes
E491651	250 ml	Glass bottle	Dye for microscopy (bacteriology) according to GRAM
E491661	250 ml	Glass bottle	Dye for microscopy (bacteriology) according to NICOLLE

Product specifications are subject to changes. Please visit our website for updates.

Germanium standard solution

• Germanio standard soluzione • Germanium solution standard • Germanio, solución patrón • Germanium-Standardlösung

Germanium standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615004400	1 l	Plastic bottle	A 100 ppm solution Ref Ph.Eur 5004400

Germanium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505632	100 ml	Plastic bottle	conc. 10 ppm Matrix: Nitric acid and hydrofluoric acid
505635	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid and hydrofluoric acid
505633	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Germanium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
504251	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid
504253	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid
504255	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrofluoric acid and nitric acid
504257	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrofluoric acid and nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Germanium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507740	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507505	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Gibb's reagent ▶ 2,6-Dichloroquinone-4-chlorimide

Giemsa's reagent

• Giemsa reattivo • Réactif de Giemsa • Giemsa reactivo • Giemsa-Reagenz

Classification transport
ONU: 1992
Transport Hazard class: 3
Packing group II



Danger
H225-H301-H370
P210-P241-P280-P301+P310a-P303+P361+P353-P403+P235

Giemsa's reagent > RS - For hematology

RS

Description ..... Blue clear liquid Identification ..... Positive Functionality ..... Conform

Code	Size	Packaging	Notes
453614	100 ml	Glass bottle	In Vitro Diagnostic Medical Device
E453612	6 x 100 ml	Bottle	In Vitro Diagnostic Medical Device
453616	500 ml	Bottle	In Vitro Diagnostic Medical Device
E453613	6 x 500 ml	Bottle	In Vitro Diagnostic Medical Device
453611	2.5 l	Bottle	In Vitro Diagnostic Medical Device
E453615	4 x 2.5 l	Bottle	In Vitro Diagnostic Medical Device

Dye according hematology Romanowski

Glass wool

• Lana di vetro • Laine de verre • Lana de vidrio • Glaswolle

Synonym:
• *Silica*
• *Silicon dioxide*

Molecular Weight: 60,09
CAS: 65997-17-3

Glass wool > RPE - For analysis

RPE

Description ..... Fine threads of glass

Code	Size	Packaging	Notes
457521	250 g	Carton box	

D-Glucitol ▶ Sorbitol

D(+)Glucose anhydrous

• D(+)Glucosio anidro • D(+)Glucose anhydre • D(+)Glucosa anhidra • D(+)Glucose wasserfrei

Synonym:
Dextrose

$CH_2OHCH(CHOH)_3CHOH$
Molecular Weight: 180,16
CAS: 50-99-7
EEC-N: 200-075-1

D(+)Glucose anhydrous > RPE - For analysis - ACS

RPE

Description ..... White crystalline powder
Identification ..... Positive
Specific optical rotation... +52.5 ÷ +53.0 °
Loss on drying ..... ≤0.2 %
Acidity ..... ≤0.002 meq/g
Starch ..... Conform
Chloride ..... ≤100 ppm
Water-insoluble matter ..... ≤50 ppm
Heavy metals (Pb) ..... ≤5 ppm
Residue on ignition ..... ≤200 ppm
Sulphat + sulphit (SO4) ..... ≤50 ppm
Fe ..... ≤5 ppm

Code	Size	Packaging	Notes
454335	100 g	Plastic bottle	
454336	500 g	Plastic bottle	
454337	1 kg	Plastic bottle	
454338	2.5 kg	Plastic bottle	
454333	25 kg	Plastic bucket	

D(+)Glucose anhydrous > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-BP-DAB-JP

ERBApharm

Description ..... White crystalline powder
Identifications ..... Positives
Identification IR (USP) ..... Positive
Appearance of solution ..... Conform
Soluble starch, sulfite ..... ≤ 15 ppm
Dextrin ..... Conform
Specific optical rotation at 20°C (anh.) .....
+52.5 ÷ +53.3 °
Sulphated ash ..... ≤ 0.1 %
Water (K.F.) ..... ≤ 1.0 %
Acidity ..... Conform
Conductivity at 20°C ..... ≤ 20 µS/cm
Conductivity at 25°C ..... ≤ 20 µS/cm
Heavy metals (Pb) ..... ≤ 4 ppm
Chloride ..... ≤ 0.018 %
Sulphate ..... ≤ 0.024 %
Loss on drying at 105°C ..... ≤ 1.0 %
As ..... ≤ 1 ppm
Assay (anh.) ..... 97.5 ÷ 102.0 %
Assay (optical rotation dried sub.)(JP) ≥ 99.5 %
Related substances (HPLC) ..... Conform
Maltose and isomaltose ..... ≤ 0.4 %
Maltotriose ..... ≤ 0.2 %
Fructose ..... ≤ 0.15 %
Unspecified ..... ≤ 0.10 %
Total impurities ..... ≤ 0.5 %

Code	Size	Packaging	Notes
346987	1 kg	Plastic bottle	
346989	5 kg	Plastic tank	
346983	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Product specifications are subject to changes. Please visit our website for updates.

D(+)Glucose monohydrate

• D(+)Glucosio monoidrato • D(+)Glucose monohydraté • D(+)Glucosa monohidrato
• D(+)-Glucose-Monohydrat

Synonym:
Dextrose monohydrate

$CH_2OHCH(CHOH)_3CHOH.H_2O$
CAS: 5996-10-1
EEC-N: 200-075-1

D(+)Glucose monohydrate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-BP-DAB-JP

ERBApharm

Description .....white crystalline powder
Identifications..... Positives
Identification IR (USP) ..... Positive
Appearance of solution .....Conform
Soluble starch, sulfite .....≤ 15 ppm
Dextrin.....Conform
Specific optical rotation at 20°C (anh.).....
+52.5 ÷ +53.3 °
Sulphated ash..... ≤ 0.1 %
Water (K.F.) .....7.5 ÷ 9.5 %
Acidity .....Conform
Conductivity at 20°C .....≤ 20 µS/cm
Conductivity at 25°C .....≤ 20 µS/cm
Heavy metals (Pb).....≤ 4 ppm
Chloride..... ≤ 0.018 %
Sulphate..... ≤ 0.025 %
Loss on drying at 105°C..... ≤ 1.0 %
As .....≤ 1 ppm
Assay (anh.).....97.5 ÷ 102.0 %
Assay (optical rotation dried sub.)(JP)≥ 99.5 %
Related substances (HPLC) .....Conform
Maltose and isomaltose ..... ≤ 0.4 %
Maltotriose..... ≤ 0.2 %
Fructose ..... ≤ 0.15 %
Unspecified..... ≤ 0.10 %
Total impurities..... ≤ 0.5 %
Origin (BSE-TSE) .....Vegetable

Code	Size	Packaging	Notes
346971	1 kg	Plastic bottle	
346972	5 kg	Plastic tank	
346973	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

L(+)Glutamic acid

• Acido L(+)glutammico • Acide L(+)glutamique • Acido L(+)glutámico • L-Glutaminsäure

Synonym:
(S)-2-Aminopentanedioic acid

$HOOCCH(NH_2)CH_2CH_2COOH$
Molecular Weight: 147,13
CAS: 56-86-0
EEC-N: 200-293-7

L(+)Glutamic acid > RPE - For analysis

RPE

Description ..... White powder
Identification ..... Positive
Potere rotat. specif. a 20°C+30.5 ÷ 32.5 ° (s.s.)
Loss on drying ..... ≤ 0.5 %
Ammonium.....≤ 200 ppm
Chloride.....≤ 200 ppm
Sulphate.....≤ 300 ppm
Heavy metals (Pb).....≤ 10 ppm
Residue on ignition..... ≤ 0.1 %
Trasmittance at 430nm (C=10; HCl 2N)≥ 98 %
Other amino-acids..... Not detectables
As203.....≤ 1 ppm
Fe .....≤ 30 ppm
Assay (non-aqueos medium) .98.5 ÷ 100.5 % (s.s.)

Code	Size	Packaging	Notes
406485	250 g	Plastic bottle	

L(+)Glutamic acid > RE - Pure

RE

Description ..... White powder
Identification ..... Positive
Loss on drying ..... ≤ 0.2 %
Other amino-acids..... Not detectables
Specific optical rotation at 20°C (C=10; .... +31.5 ÷ +32.5 °
Trasmittance at 430nm (C=10; HCl 2N)≥ 98 %
Chloride.....≤ 210 ppm
Sulphate.....≤ 280 ppm
Ammonium.....≤ 200 ppm
Heavy metals (Pb).....≤ 10 ppm
Sulphated ash..... ≤ 0.1 %
Fe .....≤ 30 ppm
Assay (acidimetric) .....98.5 ÷ 100.5 %s.s.

Code	Size	Packaging	Notes
304505	250 g	Plastic bottle	
304507	1 kg	Plastic bottle	

Glutardialdehyde solution 50%

• Aldeide glutarica soluzione 50% • Aldéhyde glutarique solution 50% • Glutaraldehído solución 50% • Glutaraldehydlösung 50%

Synonym: *Pentane-1,5-dial*

$C_5H_8O_2$
Molecular Weight: 100,12
CAS: 111-30-8

Classification transport
ONU: 2922
Transport Hazard class: 8
Packing group II



Danger
H301-H314-H334-H317-H335-H400-H411
P280-P284-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338-P342+P311a-P403+P233

Glutardialdehyde solution 50% > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 1.127 ÷ 1.133
Assay ..... 50.0 ÷ 52.0 %

Code	Size	Packaging	Notes
415721	10 l	Kubidos	

Glycerol (30°Bé)

• Glicerina (30°Bé) • Glycérine (30°Bé) • Glicerina (30°Bé) • Glycerin (30°Bé)

Synonym: *1,2,3-Propanetriol*

$CH_2OHCHOHCH_2OH$
Molecular Weight: 92,09
CAS: 56-81-5
EEC-N: 200-289-5

Glycerol (30°Bé) > RS - RSE - For electronic use

RS

Description ..... Clear liquid
Colour (APHA) ..... ≤ 10
Identification ..... Positive
Water miscibility ..... Complete
Alcohol miscibility ..... Complete
Acroleine, sugars and ammonia compounds Conform ACS
Ready carbonizable substances ..... Conform ACS
Density at 25°/25° C ..... ≥ 1.2570
Organic chlorine (Cl) ..... ≤ 30 ppm
Chloride ..... ≤ 2 ppm
Fatty acid esters(glyceryl trybutyrrate) ≤ 500 ppm
Heavy metals (Pb) ..... ≤ 2 ppm
Oxalate ..... ≤ 8 ppm
Peroxides (H2O2) ..... ≤ 5 ppm
Residue on ignition ..... ≤ 50 ppm
Sulphate ..... ≤ 10 ppm
As ..... ≤ 0.4 ppm
Ca ..... ≤ 10 ppm
Cr ..... ≤ 0.1 ppm
Cu ..... ≤ 0.1 ppm
Fe ..... ≤ 0.5 ppm
Ni ..... ≤ 0.1 ppm
Assay (densimetric) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
453771	1 l	Glass bottle	
453772	2.5 l	Glass bottle	

Glycerol (30°Bé) > RPE - For analysis - ACS - Reag. USP

RPE

Description ..... Clear liquid
Colour (APHA) ..... ≤ 10
Identification ..... Positive
Neutrality ..... Conform ACS
Ready carbonizable substances ..... Conform ACS
Water (K.F.) ..... ≤ 0.5 %
Organic chlorine (Cl) ..... ≤ 30 ppm
Heavy metals (Pb) ..... ≤ 2 ppm
Residue on ignition ..... ≤ 50 ppm
Sulphate ..... ≤ 10 ppm
Acrolein and glucose ..... Conform ACS
Esters of fatty acids (butyric acid) ≤ 500 ppm
Ca ..... ≤ 10 ppm
Fe ..... ≤ 0.5 ppm
Assay (GLC) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
453751	500 ml	Plastic bottle	
453752	1 l	Glass bottle	
453755	2.5 l	Glass bottle	
453759	35 kg	Metal drum	
453758	260 kg	Metal drum	

Glycerol (30°Bé) > ERBAPharm - According to pharmacopoeia: Ph.Eur.-USP

ERBApharm

Clear,colourless solution ..... Conform
Refractive index at 20°C ..... 1.47 - 1.475
Colour ..... ≤ 10 Hazen
Identification (IR) ..... Conform
Identification B ..... Conform
Density d25/25 ..... ≥ 1.249
Water content (K.F.) ..... ≤ 2 % m/m
Heavy metals (as Pb) ..... ≤ 5 mg/Kg
Residue on ignition ..... ≤ 100 mg/Kg
Acidity / alkalinity ..... Conform
Aldehydes ..... ≤ 10 mg/Kg
Esters ..... Conform
Chloride (Cl-) ..... ≤ 10 mg/Kg
Halogenated compounds ..... ≤ 35 mg/Kg
Sugars ..... Conform
Sulfuric ashes ..... ≤ 0.01 %
Volatile org. Impurities ..... Conform
Sulphate (SO4--) ..... ≤ 20 mg/Kg
Glycerol content ..... 99 - 101 %
Fatty acids and esters ..... Conform
Limit of chlorinated compounds ≤ 0.0030 %
Impurity A and related substances. Conform
Diethyleneglycol and ethyleneglycol imp..... Conform

Code	Size	Packaging	Notes
P6170541	10 l	Plastic tank	
P6170567	200 l	Plastic drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Product specifications are subject to changes. Please visit our website for updates.

Glycerol (30°Bé) > ERBApharm - Vegetal origin - According to Ph.Eur.-USP-FU-Ph.Franc.-BP-DAB

ERBApharm

Appearance ..... Clear colourless liquid
Identification A .....Conform Ph.Eur.
Identification B .....Conform Ph.Eur.
Identification C . Same Rt to standard by GC USP-NF
Appearance of solution S ... Clear,colourless Ph.Eur.
Acidity or alkalinity.....Conform Ph.Eur.
Sugars.....Conform Ph.Eur.
Esters.....Conform Ph.Eur.
Fatty acids and esters..... Conform USP-NF
Colour .... Not darker than standard USP-NF
Density at 25°C ..... ≥ 1.249
Refractive index at 20°C.....1.470 ÷ 1.475
Water (K.F.) ..... ≤ 2.0 %
Aldehydes.....≤ 10 ppm
Sulfated ashes ..... ≤ 0.01 %
Chloride.....≤ 10 ppm
Limit of chlorinated compounds (as Cl) ≤ 30 ppm
Halogenated compounds .....≤ 35 ppm
Heavy metals (Pb).....≤ 5 ppm
Sulphate .....≤ 20 ppm
Ethylene glycol ..... ≤ 0.1 %
Diethylene glycol ..... ≤ 0.1 %
Imp.A and related substances .....Conform Ph.Eur
Impurity Ret.time < glycerol ..... ≤ 0.1 %
Total imp.Ret.time > glycerol ..... ≤ 0.5 %
Related compounds..... Conform USP-NF
Any impurity..... ≤ 0.1 %
Total impurities..... ≤ 1.0 %
Assay (acidimetric) .....99.0 ÷ 101.0 % s.s.

Code	Size	Packaging	Notes
346161	1 l	Glass bottle	
346165	2.5 l	Glass bottle	
346162	5 l	Plastic tank	
346164	35 kg	Metal drum	
346167	250 kg	Metal drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Glycerol (30°Bé) > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C .....1.258 ÷ 1.264
Refractive index at 20°C.1.4710 ÷ 1.4740
Chloride.....≤10 ppm
Heavy metals (Pb).....≤5 ppm
Residue on ignition.....≤100 ppm
Sulphate.....≤20 ppm
As .....≤0.5 ppm
Fe .....≤1 ppm

Code	Size	Packaging	Notes
346102	1 l	Glass bottle	
346106	35 kg	Metal drum	

Glycerol 90% (28° Bé)

• Glicerina 90% (28° Bé) • Glycérine 90% (28°Bé) • Glicerina 90% (28° Bé) • Glycerin 90% (28° Bé)

Synonym:
1,2,3-Propanetriol

$CH_2OHCHOHCH_2OH$
Molecular Weight: 92,09
CAS: 56-81-5
EEC-N: 200-289-5

Glycerol 90% (28° Bé) > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C .....1.220 ÷ 1.233
Chlorinated compounds.....≤ 30 ppm
Heavy metals (Pb).....≤ 5 ppm
Residue on ignition.....≤ 100 ppm
Appearance of solution .....Conform
Acidity or alkalinity.....Conform
Refractive index at 20°C.1.4490 ÷ 1.4550
Aldehydes.....≤ 10 ppm
Esters.....Conform
Sugar .....Conform
Chloride.....≤ 10 ppm
Water (K.F.) .....11.5 ÷ 16.5 %
Assay .....83.5 ÷ 88.5 %

Code	Size	Packaging	Notes
346131	1 l	Glass bottle	
346132	2.5 l	Glass bottle	
346134	35 kg	Plastic drum	

Glycine

• Glicocolla • Glycocolle • Glicina • Glycin

Synonym:
• *Aminoacetic acid*
• *Glycocoll*

CH_2NH_2COOH
Molecular Weight: 75,07
CAS: 56-40-6
EEC-N: 200-272-2

Glycine > RPE - For analysis

RPE

Description .....White crystalline powder
Identification (I.R.)..... Positive
Hydrolyzable matter..... Conform
Loss on drying .....≤ 0.2 %
Chloride..... ≤ 70 ppm
Sulphate..... ≤ 65 ppm
Heavy metals (Pb)..... ≤ 20 ppm
Sulphated ash.....≤ 0.1 %
Assay (non-aqueos medium)..... ≥ 98.5 % (s.s.)

Code	Size	Packaging	Notes
453804	100 g	Glass bottle	
453807	1 kg	Plastic bottle	

Glycine > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-Ph.Franc.

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Hydrolyzable matter..... Conform USP-NF
pH solution 5% .....5.9 ÷ 6.4
Loss on drying ..... ≤0.2 %
Sulphated ash ..... ≤0.1 %
Chloride.....≤70 ppm
Heavy metals (Pb).....≤10 ppm
Sulphate.....≤65 ppm
Assay (non-aqueos medium) .98.5 ÷ 101.0 % s.s.
Sostanze ninidrina-positive Conform Ph.Eur.
Related substances (HPLC) Conform Ph.Eur.
Ammonium (NH4)..... ≤ 0.02 %
Ninhydrin positive substances .....Conform Ph.Eur.

Code	Size	Packaging	Notes
346207	1 kg	Plastic bottle	
346205	5 kg	Plastic tank	
346208	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Glycocoll ▶ Glycine

Glycolic acid

• Acido glicolico • Acide glycolique • Acido glicólico • Glycolsäure

Synonym:
Hydroxyacetic acid

$CH_2OHCOOH$
Molecular Weight: 76,052
CAS: 79-14-1
EEC-N: 201-180-5

Classification transport
ONU: 3261
Transport Hazard class: 8
Packing group II

Danger
H302-H332-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Glycolic acid > RPE - For analysis

RPE

Description ..... White crystals
Identification ..... Positive
Melting point..... 72 ÷ 80 °C
Assay (acidimetric) ..... ≥ 98.5 % s.s.
Water .....≤ 1 %

Code	Size	Packaging	Notes
406434	100 g	Glass bottle	

Glyoxal standard solution

• Gliossale standard soluzione • Glyoxal standard solution • Glioxal, solución patrón • Glyoxal-Standardlösung

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II

Danger
H225-HEU208
P210-P241-P280-P303+P361+P353-P403+P235-P501a

Glyoxal standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615003700	100 ml	Plastic bottle	A 20 ppm solution: to dilute according to Ph.Eur 5003700

Glyoxal-bis-(2-hydroxyanil)

• Gliossale-bis(2-idrossianile) • Glyoxal-bis-(2-hydroxyanile) • Glioxal-bis-(2-hidroxianilo)
• Glyoxalbis(2-hydroxyanil)

Synonym:
• *2,2'-(Ethanediylidenedinitrilo)diphenol*
• *GBHA*

$C_{14}H_{12}N_2O_2$
Molecular Weight: 240,26
CAS: 1149-16-2
EEC-N: 214-560-0

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Glyoxal-bis-(2-hydroxyanil) > RPE - For analysis

RPE

Description .....Beige powder
Identification ..... Positive
Loss on drying .....≤ 0.5 %
Sulphated ash.....≤ 0.1 %
Assay (HPLC) .....≥ 96.0 %

Code	Size	Packaging	Notes
454131	10 g	Glass bottle	
454132	25 g	Glass bottle	

Indicator for the complexometric titration of calcium

Product specifications are subject to changes. Please visit our website for updates.

Gold standard solution

• Oro standard soluzione • Or solution standard • Oro, solución patrón • Goldstandardlösung

Danger
H314-HEU208
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Gold standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505317	100 ml	Plastic bottle	conc. 10 ppm Matrix: Hydrochloric acid
505318	100 ml	Plastic bottle	conc. 100 ppm Matrix: Hydrochloric acid
505319	100 ml	Plastic bottle	conc. 1000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Gold standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503431	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503433	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503435	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid
503437	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Gold standard solution > RS - Standard solution for AAS

RS

Description ..... Yellow clear liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
E497585	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
E497581	500 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Gold standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Yellow clear liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
466961		Plastic ampoule	conc. 1.000 ppm Matrix: Hydrochloric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Gold(III) chloride trihydrate

• Oro (III) cloruro triidrato • Or (III) chlorure acide trihydraté • Oro tricloruro acido
• Goldenes (III) säurechloridtrohydrat

Synonym:
• *Tetrachloroauric(III) acid*
• *Hydrogen tetrachloroaurate(III)*

$HAuCl_4.3H_2O$
Molecular Weight: 393,83
CAS: 16961-25-4
EEC-N: 240-948-4

Classification transport
ONU: 3260
Transport Hazard class: 8
Packing group II

Danger
H314-H317
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P362+P364

Gold(III) chloride trihydrate > RPE - For analysis - ACS

RPE

Description ..... Orange crystalline powder Insol.in Diethyl ether ..... ≤ 0.1 % Titolo (Au) ..... ≥ 49.0 %
Identification ..... Positive Metalli (SO4) ..... ≤ 0.2 %

Code	Size	Packaging	Notes
467007	1 g	Glass ampoule	

Gram-Hucker Kit

• Kit Gram-Hucker • Kit Gram-Hucker • Kit Gram-Hucker • Gram-Hucker Kit

Classification transport
ONU: 1987
Transport Hazard class: 3
Packing group II



Danger
H225-H272-H290-H300-H310-H330-H314-H351-H361d-H371-H373-H304-H410-HEU071-HEU210
P210-P280-P284-P301+P310a-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338-P320-P361+P364-P403+P233

Gram - Hucker Kit > RS - For bacteriology

RS

Description ..... Identification ..... Positive

Code	Size	Packaging	Notes
454441	4 x 250 ml	Carton box	In Vitro Diagnostic Medical Device

Dye. Contains ethanol. 4 bottles of 250 ml. 1x 477241 Safranin T, 1x 458751 Lugol, 1x 444131 Differentiator 1x 491561 Violet oxalate

Griess' reagent

• Griess reattivo soluzione in acido acetico • Réactif de Griess • Griess reactivo solución en acido acético • Griess Reagenz

Classification transport
ONU: 3265
Transport Hazard class: 8
Packing group II



Danger
H314-HEU208
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Griess' reagent > RS - For nitrite detection

RS

Description ..... Clear pinkish liquid Identification ..... Positive

Code	Size	Packaging	Notes
454481	1 l	Glass bottle	

Griess' reagent A

• Griess reagente A • Réactif A de Griess • Griess reactivo A • Griess Reagenz A

Synonym:
• *4-Aminobenzenesulfonic acid*
• *Aniline-4-sulfonic acid*

Molecular Weight: 173,19

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III



Danger
H314-HEU208
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Griess' reagent A > RS - For nitrite detection

RS

Description ..... Clear pinkish liquid Identification ..... Positive Nitrite sensitivity ..... ≥ 1 mcg/ml

Code	Size	Packaging	Notes
454452	500 ml	Glass bottle	

Griess' reagent B

• Griess reagente B • Réactif B de Griess • Griess reactivo B • Griess Reagenz B

Classification transport
ONU: 3265
Transport Hazard class: 8
Packing group II



Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Griess' reagent B > RS - For nitrite detection

RS

Description ..... Clear pinkish liquid Identification ..... Positive Nitrite sensitivity ..... ≥ 1 mcg/ml

Code	Size	Packaging	Notes
454462	500 ml	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Gum arabic

• Gomma arabica • Gomme arabique • Goma arábica • Gummi arabisch

Synonym:
Acacia gum

CAS: 9000-01-5
EEC-N: 232-519-5

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Gum arabic > ERBApharm - According to pharmacopoeia: BP-FU-Ph.Eur.-Ph.Franc.

ERBApharm

Description	White-yellowish granules	Starch and Dextrins	Conform Ph.Eur.	Total ash	≤ 4.0 %	Escherichia coli	Absent Ph.Eur.
Identification	Positive	Saccarose and fructose	Conform Ph.Eur.	Microbial tests		Salmonella	Absent Ph.Eur.
Agar and tragacanth	Conform Ph.Eur.	Tannin	Conform Ph.Eur.	TAMC	≤ 10000 CFU/g		
Agar and sterculia	Conform Ph.Eur.	Loss on drying	≤ 10.0 %	TYMC	≤ 100 CFU/g		

Code	Size	Packaging	Notes
347107	1 kg	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Haemalum solution according to Carazzi

• Emallume reattivo soluzione secondo Carazzi • Hemalun en solution selon Carazzi • Emalume reactivo solución según Carazzi • Hemalun in Lösung nach Carazzi

HEU210

Haemalum solution according to Carazzi > RS - For histology

RS

Description ..... Liquido rosso-bruno Identification ..... Positive Density at 20°C ..... 1.086 ÷ 1.090 pH of the substance ..... 2.1 ÷ 2.3

Code	Size	Packaging	Notes
434351	250 ml	Glass bottle	In Vitro Diagnostic Medical Device

Haemalum solution according to Mayer

• Emallume reattivo secondo Mayer • Hemalun en solution selon Mayer • Emalume reactivo solución según Mayer • Hemalun in Lösung nach Mayer



Warning
H302-H371
P260-P264-P270-P301+P312a-P330-P501a

Haemalum solution according to Mayer > RS - For histology

RS

Description ..... Clear purple liquid Identification ..... Positive Assorbanza a 560 nm ..... ≥ 0.80

Code	Size	Packaging	Notes
446372	1 l	Plastic bottle	In Vitro Diagnostic Medical Device
446377	1 l	Plastic bottle	In Vitro Diagnostic Medical Device

Hafnium standard solution

• Afnio standard soluzione • Hafnium solution standard • Hafnio, solución patrón • Hafnium-Standardlösung

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group III



Hafnium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505642	100 ml	Plastic bottle	conc. 10 ppm Matrix: Nitric acid and hydrofluoric acid
505645	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid and hydrofluoric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Hafnium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
504221	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504223	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504225	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
504227	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Hafnium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507741	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507506	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Product specifications are subject to changes. Please visit our website for updates.

Hanus's reagent

• Hanus reattivo • Réactif de Hanus • Hanus reactivo • Hanus-Reagenz

Classification transport
ONU: 2920
Transport Hazard class: 8
Packing group II



Danger
H226-H314-H373
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Hanus's reagent > RS - For detection of iodine index

RS

Description ..... Brown red liquid
Identification ..... Positive
Iodine value ..... ≥ 90

Code	Size	Packaging	Notes
E454872	1 l	Glass bottle	

Hartshorn salt ▶ Ammonium carbonate

Heavy water ▶ Deuterium oxide-d2

Hematoxylin

• Ematossilina • Hématoxyline • Hematoxilina • Hämatoxylin

Synonym:
Natural Black 1

$C_{16}H_{14}O_6$
Molecular Weight: 302,29
CAS: 517-28-2
EEC-N: 208-237-3



Warning
H302-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Hematoxylin > RS - For microscopy - C.I. 75290

RS

Description ..... Brown powder
Identification ..... Positive
Loss on drying ..... 4 - 8 %
Assay ..... ≥ 75 %
Alcohol solubility ..... Conform
Melting point ..... ~ 140 °C

Code	Size	Packaging	Notes
446472	25 g	Glass bottle	
446473	100 g	Glass bottle	
446475	1 kg	Plastic bottle	

Dye for cytology

Hematoxylin solution according to Mayer

• Ematossilina soluzione secondo Mayer • Hématoxyline en solution selon Mayer
• Hematoxilina solución según Mayer • Hämatoxylinlösung nach Mayer

Synonym:
Natural Black 1

$C_{16}H_{14}O_6$
Molecular Weight: 302,29
CAS: 517-28-2

HEU210

Hematoxylin solution according to Mayer > RS - For histology

RS

Description ..... Red-violet liquid
Identification ..... Positive

Code	Size	Packaging	Notes
460511	100 ml	Bottle	In Vitro Diagnostic Medical Device
460512	6 x 100 ml	Bottle	In Vitro Diagnostic Medical Device
460513	1 l	Bottle	In Vitro Diagnostic Medical Device
460515	6 x 1 l	Bottle	In Vitro Diagnostic Medical Device

Heptafluorobutyric acid

• Acido eptafluorobutirrico • Acide heptafluorobutyrique • Acido heptafluorobutírico • Perfluorbuttersäure

$CF_3(CF_2)_2COOH$
Molecular Weight: 214,04
CAS: 375-22-4
EEC-N: 206-786-3

Classification transport
ONU: 3265
Transport Hazard class: 8
Packing group II

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Heptafluorobutyric acid > RPE - For analysis

RPE

Description ..... Colourless liquid
Identification ..... Positive
Density at 20° C ..... 1.64 ÷ 1.65
Boiling point ..... 120 ÷ 121 ° C
Assay (GC) ..... ≥ 98 %

Code	Size	Packaging	Notes
405451	10 ml	Glass bottle	

For derivatization

n-Heptane 99%

• n-Eptano 99% • n-Heptane 99% • n-Heptano 99% • n-Heptan 99%

$CH_3(CH_2)_5CH_3$
Molecular Weight: 100,21
CAS: 142-82-5
EEC-N: 205-563-8

Classification transport
ONU: 1206
Transport Hazard class: 3
Packing group II

Danger
H225-H315-H336-H304-H410
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

n-Heptane 99% > RS - For HPLC - Isocratic Grade

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.681 ÷ 0.687
Refractive index at 20°C ..... 1.3836 ÷ 1.3916
Boiling point ..... 97.9 ÷ 98.9 ° C
Water (K.F.) ..... ≤100 ppm
Residue on evaporation ..... ≤5 ppm
Acidity or alkalinity ..... ≤0.00015 meq/g
Assay (GLC) ..... ≥99.2 %
U.V. Transmittance
at 200 nm ..... ≥20 %
at 210 nm ..... ≥55 %
at 220 nm ..... ≥80 %
at 230 nm ..... ≥92 %
at 240 nm ..... ≥96 %
At 250 nm ..... ≥ 98 %
at 260 nm ..... ≥99 %
Aromatic compounds ..... ≤ 5 ppm

Code	Size	Packaging	Notes
412591000	1 l	Glass bottle	
412592000	2.5 l	Glass bottle	

n-Heptane 99% > RS - PESTIPUR - For pesticide analysis

RS

Description ..... Clear liquid
Identification ..... Positive
Colour ..... ≤ 10 APHA
Water ..... ≤ 0.005 %
Free acids (HCOOH) ..... ≤ 10 ppm
Non volatile residue ..... ≤ 5 mg/Kg
GC-ECD (Lindano) ..... ≤ 3 ng/l
Assay (GLC) ..... ≥ 99 %
Total sulphur (S) ..... ≤ 5 ppm
Ret.range 1,2,4-trichlorobenzene to decachlorobiphenyle

Code	Size	Packaging	Notes
446951	1 l	Glass bottle	
446952	2.5 l	Glass bottle	

n-Heptane 99% > RS - SPECTROSOL - For optical spectroscopy

RS

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Colour ..... ≤ 10 Hazen
Density at 20° C ..... 0.681 ÷ 0.687
Refractive index at 20°C ..... 1.3836 ÷ 1.3916
Boiling point ..... 97.9 ÷ 98.9 ° C
Water (K.F.) ..... ≤100 ppm
Residue on evaporation ..... ≤5 ppm
Acidity ..... ≤0.0005 meq/g
Alcalinity ..... ≤0.0002 meq/g
Aromatic compounds ..... ≤5 ppm
Assay (GLC) ..... ≥99 %
U.V. Transmittance
at 200 nm ..... ≥20 %
at 210 nm ..... ≥50 %
at 220 nm ..... ≥80 %
at 230 nm ..... ≥92 %
at 250 nm ..... ≥98 %

Code	Size	Packaging	Notes
446824	1 l	Glass bottle	
P0502721	2.5 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

n-Heptane 99% > RS - Anhydrous - For analysis

RS

Refractive index at 20°C ..... 1.386 - 1.390
Water content (K.F.) ..... ≤ 50 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 99 %
Aromatic compounds ..... ≤ 50 mg/Kg
Free acid (as CH3COOH) ..... ≤ 10 mg/Kg
Identification (IR) ..... Conform
Clear,colourless liq.appearance ..... Conform
Density d20/4 ..... 0.681 - 0.687
Total sulphur (S) ..... ≤ 5 ppm

Code	Size	Packaging	Notes
P0501016	1 l	Glass bottle	
P0501021	2.5 l	Glass bottle	

n-Heptane 99% > RPE - For analysis

RPE

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification ..... Conform
Chloroform miscibility ..... Complete
Diethyl ether miscib. ..... Complete
Ready carbonizable substances ..... Conform
Density at 20° C ..... 0.681 ÷ 0.687
Refractive index at 20°C 1.3836 ÷ 1.3916
Boiling point ..... 97.4 ÷ 99.4 ° C
Water (K.F.) ..... ≤100 ppm
Residue on evaporation ..... ≤ 10 ppm
Subst. reducing KMnO4 ..... ≤20 ppm(5m)
Tiophene ..... ≤10 ppm
Total sulphur ..... ≤ 5 ppm
Assay (GLC) ..... ≥99.0 %
Aromatic compounds ..... ≤ 50 ppm
Acidity (acetic acid) ..... ≤ 10 ppm

Code	Size	Packaging	Notes
446787	1 l	Glass bottle	
446785	2.5 l	Glass bottle	
446781	5 l	Aluminium can	
446783	5 l	Plastic tank	
446782	18 kg	Metal drum	
446789	135 kg	Metal drum	
446788	200 l	Metal drum	

n-Heptane 99% > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Colour ..... ≤ 10 APHA
Density at 20° C ..... 0.679 ÷ 0.689
Refractive index at 20°C 1.3826 ÷ 1.3926
Boiling point ..... 97.4 ÷ 99.4 ° C
Water (K.F.) ..... ≤ 150 ppm
Residue on evaporation ..... ≤ 20 ppm
Acidità (ac.acetico) ..... ≤ 10 ppm
Total sulphur ..... ≤ 5 ppm
Assay (GC) ..... ≥ 99 %
Aromatic compounds ..... ≤ 50 ppm

Code	Size	Packaging	Notes
339381	1 l	Glass bottle	
339385	2.5 l	Glass bottle	
339382	5 l	Aluminium can	
528224	5 l	Plastic tank	
528228	10 l	Metal tank	
339386	18 kg	Metal drum	
528225	25 l	Metal drum	
528226	200 l	Metal drum	

n-Heptane 99% > RE - ASTM

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.679 ÷ 0.689
Refractive index at 20°C 1.3826 ÷ 1.3926
Boiling point ..... 97.4 ÷ 99.4 ° C
Water (K.F.) ..... ≤ 150 ppm
Residue on evaporation ..... ≤ 20 ppm
Lead ..... ≤ 0.002 g/gal
Isooctane (ASTM) ..... ≤ 0.10 % v/v
n-heptane (ASTM) ..... ≥ 99.75 % v/v
Assay (GC) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
524263	5 l	Plastic tank	
524265	25 l	Metal drum	
524267	200 l	Metal drum	

Suitable for ASTM methods D2700 and D2699

n-Heptane
• n-Eptano • n-Heptane • n-Heptano • n-Heptan

$CH_3(CH_2)_5CH_3$
Molecular Weight: 100,21
CAS: 142-82-5
EEC-N: 205-563-8

Classification transport
ONU: 1206
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H336-H304-H410
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

n-Heptane > RS - For HPLC - Isocratic Grade — RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.681 - 0.687
Refractive index at 20°C ..... 1.386 - 1.390
Water (K.F.) ..... ≤ 100 ppm
Residue on evaporation ..... ≤ 5 ppm
Acidity or alkalinity ..... ≤ 0.00015 meq/g
Assay (GLC) ..... ≥ 95 %
U.V. Transmittance
at 200 nm ..... ≥ 20 %
at 210 nm ..... ≥ 45 %
at 220 nm ..... ≥ 80 %
At 250 nm ..... ≥ 98 %
Aromatic compounds ..... ≤ 10 ppm

Code	Size	Packaging	Notes
446831	1 l	Glass bottle	
446832	2.5 l	Glass bottle	

n-Heptane > RE - Pure — RE

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Identification ..... Positive
Residue on evaporation ..... ≤ 10 ppm
Water (K.F.) ..... ≤ 100 ppm
Assay (GLC) ..... ≥ 95.0 %
Benzene ..... ≤ 10 ppm
Aromatics ..... ≤ 100 ppm

Code	Size	Packaging	Notes
508212	1 l	Glass bottle	
508215	5 l	Plastic tank	
508216	25 l	Metal drum	
508217	200 l	Metal drum	

Heptane mixture of isomers
• Eptano miscela di isomeri • Heptane mélange d´isomères • Heptano mezcla de isómero • Heptan Isomerengemisch

$CH_3(CH_2)_5CH_3$
Molecular Weight: 100,21
CAS: 142-82-5
EEC-N: 927-510-4

Classification transport
ONU: 1206
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H336-H304-H411
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Heptane mixture of isomers > RS - PESTIPUR - For pesticide analysis — RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C ..... 1.387 - 1.407
Water content (K.F.) ..... ≤ 50 mg/Kg
Density at 15°C ..... 690 - 725 g/l
Colour ..... ≤ 10 Hazen
Distillation range ..... 85 - 105 °C
Aromatic compounds ..... ≤ 100 mg/Kg
Non volatile residue ..... ≤ 5 mg/Kg
GC-ECD.Individual peak (Lindane) .. ≤ 3 ng/l
Retention time trichlorobenzene to mirex
Ret.range 1,2,4-trichlorobenzene
to decachlorobiphenyle

Code	Size	Packaging	Notes
446841	1 l	Glass bottle	
446842	2.5 l	Glass bottle	

For chlorinated compounds analysis

Heptane mixture of isomers > RPE - For analysis — RPE

Appearance ..... Clear colourless liquid
Identification ..... Conform
Colour ..... ≤ 10 Apha
Refractive index at 20°C ..... 1.387 - 1.407
Density at 15°C ..... 690 - 725 g/l
Distillation range ..... 85 - 105 °C
End/initial boiling points difference ≤ 10 °C
Water content (K.F.) ..... ≤ 100 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
Aromatic compounds ..... ≤ 100 mg/Kg
n-hexane ..... None %
Toluene ..... None mg/Kg

Code	Size	Packaging	Notes
524381	5 l	Plastic tank	

Product specifications are subject to changes. Please visit our website for updates.

Heptane mixture of isomers > RE - Pure

RE

Description ..... Clear liquid
Identification ..... Positive
Color ..... ≤ 10 APHA
Refractive index at 20°C ..... 1,3870 ÷ 1,4070
Density at 15°C ..... 0.690 ÷ 0.725
Boiling point ..... 85 ÷ 105 °C
Water (K.F.) ..... ≤ 150 ppm
Residue on evaporation ..... ≤ 100 ppm
Aromatics (Thiophene) ..... ≤ 200 ppm

Code	Size	Packaging	Notes
528245	5 l	Plastic tank	
528246	25 l	Metal drum	
528247	200 l	Metal drum	

1-Heptanesulphonic acid sodium salt

• Acido 1-eptansolfonico sale sodico • Acide 1-heptanesulfonique sel sodique
• Acido 1-heptanosulfónico sal sódica • 1-Heptansulfonsäure-Natriumsalz

Synonym:
Sodium 1-heptanesulfonate

$C_7H_{15}O_3SNa$
Molecular Weight: 202,25
CAS: 22767-50-6
EEC-N: 245-210-5

1-Heptanesulphonic acid sodium salt > RS - For ion pair chromatography

RS

Description ..... White crystalline powder
Identification ..... Positive
Loss on drying ..... ≤ 2 %
Assay ..... ≥ 99 %
Assorbanza (Sol. 0.25M)
at 200 nm ..... ≤ 0.1 AU
at 210 nm ..... ≤ 0.05 AU
at 220 nm ..... ≤ 0.04 AU
at 230 nm ..... ≤ 0.03 AU
at 240 nm ..... ≤ 0.01 AU
at 250 nm ..... ≤ 0.01 AU
at 260 nm ..... ≤ 0.01 AU

Code	Size	Packaging	Notes
405851	25 g	Glass bottle	
405852	100 g	Plastic bottle	

Hexachloroplatinic acid hexahydrate

• Acido esacloroplatinico esaidrato • Acide hexachloroplatinique hexahydraté
• Acido hexacloroplatínico hexahidrato • Hexachloroplatinsäure-Hexahydrat

Synonym:
Chloroplatinic acid hexahydrate

$H_2PtCl_6.6H_2O$
Molecular Weight: 517,92
CAS: 18497-13-7
EEC-N: 241-010-7

Danger
H315-H319-H334-H317
P261-P284-P304+P340-P305+P351+P338-P337+P313-P342+P311a

Hexachloroplatinic acid hexahydrate > RPE - For analysis

RPE

Description ..... Red-orange mass
Identification ..... Positive
Tracce tot. di metalli ..... ≤ 0.1 %
Assay ..... 38 ÷ 40 % Pt

Code	Size	Packaging	Notes
470017	1 g	Glass bottle	

Hexadecane

• Esadecano • Hexadécane • Hexadecano • Hexadecan

Synonym:
Cetane

$C_{16}H_{34}$
Molecular Weight: 226,44
CAS: 544-76-3
EEC-N: 208-878-9

Danger
H304
P301+P310a-P331-P405-P501a

Hexadecane > RE - Pure

RE

Appearance ..... Clear colourless liquid
Colour ..... ≤ 10 Hazen
Identification ..... Conform
Refractive index at 20°C ..... 1.433 - 1.437
Assay (GC) ..... ≥ 99 %

Code	Size	Packaging	Notes
P0853016	1 l	Glass bottle	

1-Hexadecanol ▶ Cetyl alcohol

Hexadecyltrimethylammonium bromide ▶ Trimethylcetylammonium bromide

Hexafluoro-2-propanol

• Esafluoro-2-propanolo • Hexafluoro-2-Propanol • Hexafluoro-2-Propanol • Hexafluor-2-propanol

Synonym:
• *Hexafluoroisopropanol*
• *HFP*

$C_3H_2F_6O$
Molecular Weight: 168,04
CAS: 920-66-1
EEC-N: 213-059-4

Classification transport
ONU: 3265
Transport Hazard class: 8
Packing group II



Danger
H302-H312-H332-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P362+P364

Hexafluoro-2-propanol > RPE - For analysis

RPE

Identification (IR)..... Conform Colour ..... ≤ 10 Hazen Assay (GC) ..... ≥ 99 %

Code	Size	Packaging	Notes
P6080503	100 ml	Glass bottle	
P6080518	500 ml	Glass bottle	

Hexahydrotoluene ▶ Methylcyclohexane

Hexamethyldisilazane

• Esametildisilazano • Hexaméthyldisilazane • Hexametildisilazano • Hexamethyldisilazan

Synonym:
HMDS

$(CH_3)_3SiNHSi(CH_3)_3$
Molecular Weight: 161,4
CAS: 999-97-3
EEC-N: 213-668-5

Classification transport
ONU: 3286
Transport Hazard class: LQ



Danger
H225-H302-H311-H315-H319-H335
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Hexamethyldisilazane > RPE - For analysis

RPE

Description ..... Clear liquid Identification ..... Positive Assay (GLC) ..... ≥ 97.5 %
Colour ..... ≤ 10 APHA Refractive index at 20°C..... 1.4060 ÷ 1.4090

Code	Size	Packaging	Notes
446731	25 ml	Glass bottle	

For derivatization

Hexamethylenetetramine

• Esametilentetrammina • Hexaméthylènetétramine • Hexametilentetramina • Hexamethylentetramin

Synonym:
• *1,3,5,7-Tetraazatricyclo[3.3.1.13,7]decane*
• *Hexamine*

$(CH_2)_6N_4$
Molecular Weight: 140,19
CAS: 100-97-0
EEC-N: 202-905-8

Classification transport
ONU: 1328
Transport Hazard class: 4.1
Packing group III



Warning
H228-H317
P210-P241-P261-P280-P333+P313-P501a

Hexamethylenetetramine > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description .....colourless/white crystalline powder Identification ..... Positive Ash ..... ≤ 0.02 % Loss on drying 130°C..... ≤ 2 %
Water (K.F.) ..... ≤ 0.5 % Assay (alkalimetric)..... ≥ 99.0 % Heavy metals (Pb)..... ≤ 10 ppm

Code	Size	Packaging	Notes
446875	250 g	Plastic bottle	
446876	1 kg	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

n-Hexane 99%

• n-Esano 99% • n-Hexane 99% • n-Hexano 99% • n-Hexan 99%

$CH_3(CH_2)_4CH_3$
Molecular Weight: 86,18
CAS: 110-54-3
EEC-N: 203-777-6

Classification transport
ONU: 1208
Transport Hazard class: 3
Packing group II

Danger
H225-H315-H361f-H336-H373-H304-H411
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

n-Hexane 99% > RS - For HPLC - Isocratic Grade

RS

Clear,colourless liq.appearance .....Conform
Identification .....Conform
Colour .....≤ 10 Apha
Refractive index at 20°C.....1.373 - 1.377
Water content (K.F.) .....≤ 150 mg/Kg
Aromatic compounds.....≤ 5 mg/Kg
UV transmittance at 210 nm ..... ≥ 50 %
UV transmittance at 220 nm ..... ≥ 80 %
UV transmittance at 230 nm ..... ≥ 90 %
UV transmittance at 245 nm ..... ≥ 98 %
Non volatile residue.....≤ 5 mg/Kg
Assay (GC)..... ≥ 99 %
Total sulphur (S) .....≤ 5 ppm

Code	Size	Packaging	Notes
412691	1 l	Glass bottle	
412692	2.5 l	Glass bottle	

n-Hexane 99% > RS - For GC-MS

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C..... 1.373 - 1.377
Water (K.F.) .....≤ 50 ppm
Residue on evaporation .....≤ 2 ppm
Colour .....≤ 10 APHA
Assay (GC).....≥ 99.0 %
GC-MS.Individual peak (n-hexadecane).....≤ 2 µg/L
Ret.range n-undecane to n-tetracontane (scanning area 30-600amu)

Code	Size	Packaging	Notes
447212	1 l	Glass bottle	

n-Hexane 99% > RS - ATRASOL - For traces analysis

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C..... 1.373 - 1.377
Water content (K.F.) .....≤ 50 mg/Kg
Colour .....≤ 5 Hazen
Assay (GC) ..... ≥ 99 %
Non volatile residue .....≤ 2 mg/Kg
GC (FID) - NC Atrasol .....Conform
GC-ECD.Individual peak (CCl4) .....≤ 1 µg/l
Ret.range dichloromethane to 1,2,4-trichlorobenzene
GC-ECD.Individual peak (Lindane) .≤ 2 ng/L
Ret.range 1,2,4-trichlorobenzene
to decachlorobiphenyle
GC-FID.Individual peak (n-hexadecane). ≤ 2 µg/L
Ret.range n-undecane to n-tetracontane

Code	Size	Packaging	Notes
P052323016	1 l	Glass bottle	
P052323021	2.5 l	Glass bottle	

n-Hexane 99% > RS - PESTIPUR - For pesticide analysis

RS

Description .....Clear liquid
Identification ..... Positive
Colour .....≤ 10 hazen
Assay (GLC) ..... ≥ 99 %
Not volatile residue.....≤ 2 mg/kg
Water .....≤ 100 mg/kg
GC-ECD (Lindane standard) .....≤ 3 ng/l
GC-NPD (Ethylparathion standard) ..≤ 3 ng/l

Code	Size	Packaging	Notes
447111	1 l	Glass bottle	
447112000	2.5 l	Glass bottle	

n-Hexane 99% > RS - SPECTROSOL - For optical spectroscopy

RS

Clear,colourless liq.appearance .....Conform
Identification .....Conform
Colour .....≤ 10 Apha
Refractive index at 20°C.....1.373 - 1.377
Water content (K.F.) .....≤ 150 mg/Kg
Total sulphur (S) .....≤ 5 ppm
Aromatic compounds.....≤ 5 mg/Kg
Non volatile residue.....≤ 5 mg/Kg
Assay (GC)..... ≥ 99 %
UV transmittance at 210 nm ..... ≥ 60 %
UV transmittance at 215 nm ..... ≥ 70 %
UV transmittance at 220 nm ..... ≥ 80 %
UV transmittance at 230 nm ..... ≥ 94 %
UV transmittance at 245 nm ..... ≥ 98 %

Code	Size	Packaging	Notes
447051	1 l	Glass bottle	
447052	2.5 l	Glass bottle	

n-Hexane 99% > RS - Anhydrous - For HPLC

RS

Refractive index at 20°C.....1.373 - 1.377
Water content (K.F.) .....≤ 50 mg/Kg
Colour .....≤ 10 Hazen
Aromatic compounds.....≤ 5 mg/Kg
UV transmittance at 220 nm ..... ≥ 80 %
UV transmittance at 230 nm ..... ≥ 90 %
UV transmittance at 210 nm ..... ≥ 50 %
UV transmittance at 245 nm ..... ≥ 98 %
Non volatile residue.....≤ 5 mg/Kg
Assay (GC)..... ≥ 99 %
Total sulphur (S) .....≤ 5 ppm

Code	Size	Packaging	Notes
P05230S01/16	1 l	Glass bottle	
P05230S01/21	2.5 l	Glass bottle	

n-Hexane 99% > RPE - For analysis

RPE

Description ..... Clear colourless liquid	Non volatile residue .....≤ 10 mg/Kg	Water (K.F.) .....≤ 100 ppm	Total sulphur (S) .....≤ 5 ppm
Water content (K.F.) .....≤ 100 mg/Kg	Colour .....≤ 10 APHA	Assay (GC) ..... ≥ 99 %	Residue on evaporation .....≤ 10 ppm
Identification ..... Positive	Refractive index at 20°C.....1.373 ÷ 1.377	Aromatic compounds.....≤ 10 ppm	Assay (CPG) ..... ≥ 99 %

Code	Size	Packaging	Notes
447041	1 l	Glass bottle	
447042	2.5 l	Glass bottle	
P052053068	200 l	Metal drum	

n-Hexane 99% > RE - Pure

RE

Description ..... Clear liquid	Density at 20°C .....0.658 ÷ 0.662	Residue on evaporation .....≤ 20 ppm	Total sulphur .....≤ 5 ppm
Colour .....≤ 10 APHA	Refractive index at 20°C.....1.373 ÷ 1.377	Water (K.F.) .....≤ 150 ppm	Assay (GLC) ..... ≥ 99 %
Identification ..... Positive	Boiling point.....68.2 ÷ 69.2 °C	Bromine rating ..... ≤ 1	Aromatic compounds.....≤ 50 ppm

Code	Size	Packaging	Notes
528950	5 l	Plastic tank	
528951	25 l	Metal drum	
528952	200 l	Metal drum	

n-Hexane

• n-Esano • n-Hexane • n-Hexano • n-Hexan

$CH_3(CH_2)_4CH_3$
Molecular Weight: 86,18
CAS: 110-54-3
EEC-N: 203-777-6

Classification transport
ONU: 1208
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H361f-H336-H373-H304-H411
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

n-Hexane > RS - For HPLC - Isocratic grade - Reag.Ph.Eur.

RS

Description ..... Clear colourless liquid	Residue on evaporation .....≤2 ppm	at 230 nm ..... ≥92 %	Density at 20°C .....0.659 - 0.663
Identification ..... Positive	Acidity or alkalinity.....≤0.00015 meq/g	at 240 nm ..... ≥95 %	UV Absorbance from 260 nm to 420 nm.. ≤ 0.01 AU
Density at 25° C ..... ≥0.662	Assay (GLC) ..... ≥96 %	At 245 nm ..... ≥ 98 %	
Refractive index at 20°C.1.3750 ÷ 1.3760	U.V. Transmittance	at 250 nm ..... ≥99 %	
Boiling point..... 67 - 69 ° C	at 210 nm ..... ≥50 %	Aromatic compounds.....≤ 10 ppm	
Water (K.F.) .....≤100 ppm	at 220 nm ..... ≥82 %	Total sulphur (S) .....≤ 5 ppm	

Code	Size	Packaging	Notes
412601000	1 l	Glass bottle	
412602000	2.5 l	Glass bottle	

n-Hexane > RS - ATRASOL - For trace analysis, Suitable for Hydrocarbon index determination

RS

Appearance ..... Clear colourless liquid	Refractive index at 20°C.....1.373 - 1.377	Assay (GC) ..... ≥ 96.5 %	GC-ECD.Individual peak (Lindane) .≤ 2 ng/L
Identification .....Conform	Colour ..... ≤ 5 Hazen	Non volatile residue .....≤ 2 mg/Kg	GC-FID.Individual peak (C10-C40)..≤ 2 µg/l
Density d20/4 .....0.655 - 0.665	Water content (K.F.) .....≤ 50 mg/Kg	GC-FID.Hydrocarbon oil index...≤ 0.05 mg/l	to decachlorobiphenyle

Code	Size	Packaging	Notes
P0523216	1 l	Glass bottle	
P0523221	2.5 l	Glass bottle	

According to NF-EN-ISO9377-2 for hydrocarbon index determination

n-Hexane > RS - PESTIPUR - For pesticide analysis

RS

Description ..... Clear liquid	Refractive index at 20° C..... 1.373 ÷ 1.377	GC-ECD (Lindano) .....≤ 3 ng/l
Identification ..... Positive	Water (K.F.) .....≤ 150 ppm	GC-NPD (Ethylparation).....≤ 3 ng/l
Colour .....≤ 10 Hazen	Not volatile residue .....≤ 5 ppm	Assay (GLC) .....≥ 95 %

Code	Size	Packaging	Notes
447011	1 l	Glass bottle	
447012	2.5 l	Glass bottle	
447013	4 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

n-Hexane > RS - SPECTROSOL - For optical spectroscopy - Reag.Ph.Eur.

RS

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Density at 20°C ..... 0.659 ÷ 0.663
Refractive index at 20°C ..... 1.375 ÷ 1.376
Boiling point ..... 67 ÷ 69 ° C
Acidity or alkalinity ..... ≤0.00015 meq/g
Water (K.F.) ..... ≤100 ppm
Residue on evaporation ..... ≤2 ppm
Aromatic compounds ..... ≤5 ppm
Total sulphur ..... ≤5 ppm
Assay (GLC) ..... ≥95 %
U.V. Transmittance
at 220 nm ..... ≥82 %
at 230 nm ..... ≥92 %
at 240 nm ..... ≥95 %
at 250 nm ..... ≥99 %
UV Absorbance from 260 nm to 420 nm.. ≤ 0.01 AU

Code	Size	Packaging	Notes
446934	1 l	Glass bottle	
446932	2.5 l	Glass bottle	

n-Hexane > RS - Anhydrous - For analysis

RS

Refractive index at 20°C ..... 1.373 - 1.377
Water content (K.F.) ..... ≤ 50 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 95 %
Aromatic compounds ..... ≤ 10 mg/Kg
Clear,colourless liq.appearance ..... Conform
Identification (IR) ..... Conform
Density d20/4 ..... 0.655 - 0.665
Total sulphur (S) ..... ≤ 5 ppm

Code	Size	Packaging	Notes
P0521016	1 l	Glass bottle	

n-Hexane > RPE - For analysis - ACS - Reag. Ph.Eur.

RPE

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Alcohol miscibility ..... Complete
Diethyl ether miscib. ..... Complete
Chloroform miscibility ..... Complete
Ready carbonizable substances ..... Conform
Density at 20°C ..... 0.659 ÷ 0.663
Refractive index at 20°C ..... 1.375 ÷ 1.376
Distillation range ..... 67 ÷ 69 ° C
Water (K.F.) ..... ≤100 ppm
Residue on evaporation ..... ≤10 ppm
Water-soluble titrable acid ..≤0.0003 meq/g
Subst. reducing KMnO4 ..... ≤20 ppm (5m)
Tiophene ..... Conform ACS
Total sulphur ..... ≤ 5 ppm
Cu ..... ≤0.01 ppm
Fe ..... ≤0.1 ppm
Ni ..... ≤0.01 ppm
Pb ..... ≤0.01 ppm
Zn ..... ≤0.01 ppm
Assay (GLC) ..... ≥95 %
Assay (hexan isomer+methylcyclopentane)≥ 98.5 %
Aromatic compounds ..... ≤ 10 ppm
Al ..... ≤ 0.5 ppm

Code	Size	Packaging	Notes
446907	1 l	Glass bottle	
446902	2.5 l	Glass bottle PVC coated	
446903	2.5 l	Glass bottle	
446901	5 l	Aluminium can	
446991	5 l	Plastic tank	
446905	18 kg	Metal drum	
446904	130 kg	Metal drum	

n-Hexane > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Tiophene ..... Conform
Density at 20° C ..... 0.655 ÷ 0.665
Refractive index at 20°C. 1.3699 ÷ 1.3799
Boiling point ..... 68.2 ÷ 69.2 °C
Water (K.F.) ..... ≤ 100 ppm
Residue on evaporation ..... ≤ 20 ppm
Acidity (Caproic acid) ..... ≤ 40 ppm
Assay (GLC) ..... ≥ 95 %
Colour ..... ≤ 10 APHA
Aromatic compounds ..... ≤ 50 ppm
Total sulphur (S) ..... ≤ 5 ppm

Code	Size	Packaging	Notes
339751	1 l	Glass bottle	
339755	2.5 l	Glass bottle	
339752	5 l	Plastic tank	
339756	18 kg	Metal drum	
339758	25 l	Metal drum	
339757	130 kg	Metal drum	
339759	200 l	Metal drum	

Hexane mixture of isomers

• Esano miscela di isomeri • Hexane mélange d´isomères • Hexano mezcla de isómeros • Hexan Isomerengemisch

$CH_3(CH_2)_4CH_3$
Molecular Weight: 86,18
CAS: 110-54-3
EEC-N: 925-292-5

Classification transport
ONU: 1208
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H361-H336-H373-H304-H411
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Hexane mixture of isomers > RS - For HPLC - Isocratic Grade

RS

Clear,colourless liq.appearance ..... Conform
Identification ..... Conform
Colour ..... ≤ 10 Hazen
Refractive index at 20°C ..... 1.375 - 1.383
Water content (K.F.) ..... ≤ 150 mg/Kg
Aromatic compounds ..... ≤ 5 mg/Kg
Total sulphur (S) ..... ≤ 5 ppm
Non volatile residue ..... ≤ 5 mg/Kg
UV transmittance at 240 nm ..... ≥ 94 %
UV transmittance at 250 nm ..... ≥ 98 %
UV transmittance at 260 nm ..... ≥ 98 %

Code	Size	Packaging	Notes
412632	1 l	Glass bottle	
412631	2.5 l	Glass bottle	

Hexane mixture of isomers > RS - PESTIPUR - For pesticide analysis

RS

Description ..... Clear liquid
Identification ..... Positive
Colour ..... ≤ 10 Hazen
Refractive index at 20° C ..... 1.375 ÷ 1.383
Water (K.F.) ..... ≤ 150 ppm
Not volatile residue ..... ≤ 5 ppm
GC-ECD (Lindano) ..... ≤ 3 ng/l
GC-NPD (Ethylparation) ..... ≤ 3 ng/l
GC-ECD (CCl4) ..... ≤ 1 µg/l

Code	Size	Packaging	Notes
447181	1 l	Glass bottle	
447182	2.5 l	Glass bottle	

Hexane mixture of isomers > RPE - For analysis

RPE

Clear,colourless liq.appearance ..... Conform
Identification ..... Conform
Colour ..... ≤ 10 Apha
Refractive index at 20°C ..... 1.375 - 1.383
Water content (K.F.) ..... ≤ 100 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
GC chromatogram ..... Conform
Aromatic compounds ..... ≤ 50 mg/Kg
Total sulphur (S) ..... ≤ 5 ppm

Code	Size	Packaging	Notes
446892	1 l	Glass bottle	
446891	2.5 l	Glass bottle	
446893	25 l	Metal drum	

Hexane mixture of isomers > RE - Pure

RE

Description ..... Clear colourless liquid
Water content (K.F.) ..... ≤ 200 mg/Kg
Identification ..... Positive
Non volatile residue ..... ≤ 50 mg/Kg
Density at 20° C ..... 0.665 ÷ 0.675
Colour ..... ≤ 10 Hazen
Refractive index at 20°C. 1.3750 ÷ 1.3850
GC chromatogram ..... Conform
Boiling point ..... 63 ÷ 70 ° C
Total sulphur (S) ..... ≤ 5 ppm
Water (K.F.) ..... ≤200 ppm
Aromatics ..... ≤30 ppm
Residue on evaporation ..... ≤50 ppm
Total sulphur ..... ≤5 ppm

Code	Size	Packaging	Notes
339851	1 l	Glass bottle	
339852	2.5 l	Glass bottle	
528940	5 l	Plastic tank	
P0520040	10 l	Metal drum	
339856	18 kg	Metal drum	
528941	25 l	Metal drum	
528942	200 l	Metal drum	

Product specifications are subject to changes. Please visit our website for updates.

1-Hexanesulphonic acid sodium salt

• Acido 1-esansolfonico sale sodico • Acide 1-hexanesulfonique sel sodique
• Acido 1-hexanosulfónico sal sódica • 1-Hexansulfonsäure-Natriumsalz

Synonym:
Sodium hexanesulfonate

$CH_3(CH_2)_5SO_3Na$
Molecular Weight: 188,22
CAS: 2832-45-3
EEC-N: 220-601-3

Danger
H315-H319-H334-H317
P261-P284-P304+P340-P305+P351+P338-P337+P313-P342+P311a

1-Hexanesulphonic acid sodium salt > RS - For ion pair chromatography

RS

Description ..... White crystalline powder
Water (K.F.) ..... ≤ 2.0 %
Assay ..... ≥ 99.0 %

Absorbance (0,25M)
At 200 nm ..... ≤ 0.1 AU
At 210 nm ..... ≤ 0.05 AU
At 220 nm ..... ≤ 0.04 AU
At 230 nm ..... ≤ 0.03 AU
At 240 nm ..... ≤ 0.01 AU
At 250 nm ..... ≤ 0.01 AU
At 260 nm ..... ≤ 0.01 AU

Code	Size	Packaging	Notes
405621	25 g	Glass bottle	
405622	100 g	Plastic bottle	

1-Hexanesulphonic acid sodium salt monohydrate

• Acido 1-esansolfonico sale sodico monoidrato • Acide 1-hexanesulfonique sel sodique monohydrate
• Acido 1-hexanosulfónico sal sódica monohidrato • 1-Hexansulfonsäure-Natriumsalz-Monohydrat

Synonym:
• *Sodium 1-hexanesulfonate monohydrate*

$C_6H_{13}NaO_3S.H_2O$
Molecular Weight: 206,24
CAS: 207300-91-2

1-Hexanesulphonic acid sodium salt monohydrate > RS - For ion pair chromatography

RS

Description ..... White crystalline powder
Loss on drying ..... 7.0 - 9.0 %
Assay ..... ≥ 99.0 %

Absorbance (0,25M)
At 200 nm ..... ≤ 0.1 AU
At 210 nm ..... ≤ 0.05 AU
At 220 nm ..... ≤ 0.03 AU
At 230 nm ..... ≤ 0.02 AU
At 240 nm ..... ≤ 0.01 AU
At 250 nm ..... ≤ 0.01 AU

Code	Size	Packaging	Notes
405921	25 g	Glass bottle	
405922	100 g	Plastic bottle	

Hexanoic acid ▶ n-Caproic acid

Histamine dihydrochloride

• Istamina bicloridrato • Histamine dichlorhydratée • Histamina diclorhidrato • Histamindihydrochlorid

Synonym:
2-(4-Imidazolyl)ethylamine dihydrochloride

$C_5H_9N_3 \cdot 2HCl$
Molecular Weight: 184,07
CAS: 56-92-8
EEC-N: 200-298-4

Danger
H315-H319-H334-H317-H335
P261-P284-P304+P340-P305+P351+P338-P342+P311a-P403+P233

Histamine dihydrochloride > RPE - For analysis

RPE

Description ..... Yellowish crystalline powder
Identification ..... Positive
Loss on drying ..... ≤ 0.5 %
Residue on ignition ..... ≤ 0.1 %
Assay (non-aqueos medium) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
456851	10 g	Glass bottle	
456852	100 g	Glass bottle	

L-Histidine hydrochloride monohydrate

• L-Istidina monocloridrato monoidrato • L-Histidine monochlorhydratée monohydraté
• L-Histidina monoclorhidrato monohidrato • L-Histidinmonohydrochloridmonohydrat

Synonym:
• *L-Histidine monohydrochloride monohydrate*
• *L-alpha-Amino-beta-(4-imidazolyl)propanoic acid monohydrochloride*

$C_6H_9N_3O_2.HCl.H_2O$
Molecular Weight: 209,63
CAS: 5934-29-2

L-Histidine hydrochloride monohydrate > RPE - For analysis

RPE

Description White or almost white crystalline powder
Identification ..... Positive
Sulfate ..... ≤ 300 ppm ppm
Ammonium ..... ≤ 200 ppm
Fe ..... ≤ 10 ppm
Heavy metals (Pb) ..... ≤ 10 ppm
Loss on drying ..... 7.0 ÷ 10.0 %
Sulfated ashes ..... ≤ 0.1 %
pH sol. 5% ..... 3.0 ÷ 5.0
Assay (dried base) ..... 98.5 ÷ 101.0 %

Code	Size	Packaging	Notes
456952	25 g	Glass bottle	
456951	500 g	Plastic bottle	

Histolemon

• Histolemon • Histolemon • Histolemon • Histolemon

Synonym:
• *Orange oil*
• *Citrus sinensis*

CAS: 8028-48-6

Classification transport
ONU: 2052
Transport Hazard class: 3
Packing group III



Danger
H226-H315-H317-H304-H410
P210-P241-P261-P280-P301+P310a-P303+P361+P353

Histolemon > RS - For histology

RS

Description ..... Clear liquid
Colour (APHA) ..... ≤ 25
Identification ..... Positive
Density at 20° C ..... 0.835 ÷ 0.845
Assay (GLC) ..... ≥ 95 %

Code	Size	Packaging	Notes
454911	1 l	Glass bottle	
454912	2.5 l	Glass bottle	
454915	5 l	Plastic tank	

Citrus based histological clearing agent

Holmium standard solution

• Olmio standard soluzione • Holmium solution standard • Holmio, solución patrón • Holmium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group II

Holmium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505657	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505658	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Holmium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
504261	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504263	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504265	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
504267	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Product specifications are subject to changes. Please visit our website for updates.

Holmium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507742	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507507	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Holmium oxide

• Olmio ossido • Holmium oxyde • Holmio oxido • Holmiumoxid

Synonym:
Holmium(III) oxide

Ho_2O_3
Molecular Weight: 377,88
CAS: 12055-62-8
EEC-N: 235-015-3

Holmium oxide > RPE - For analysis

RPE

Description ..... Yellowish powder Identification ..... Positive Assay ..... ≥ 99.85 %

Code	Size	Packaging	Notes
466831	1 g	Glass bottle	

Holmium perchlorate in solution

• Olmio perclorato soluzione • Holmium perchlorate solution • Holmio perclorato solución • Holmiumperchloratlösung

$Ho(ClO_4)_3$
Molecular Weight: 481,3
CAS: 14017-54-0

Holmium perchlorate in solution > RS - For analysis according to Ph. Eur. Chap. 2.2.25

RS

Code	Size	Packaging	Notes
506473	100 ml	Glass bottle	

Holmium perchlorate in solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611043101	1 l	Glass bottle	Ref Ph.Eur 1043101

Hyamine 1622

• Hyamine 1622 • Hyamine 1622 • Hyamina 1622 • Hyamine 1622

Synonym:
• *Benzethonium chloride*
• *(Diisobutylphenoxyethoxyethyl)dimethylbenzylammonium chloride*

$C_{27}H_{42}ClNO_2$
Molecular Weight: 448,18
CAS: 121-54-0
EEC-N: 204-479-9

Classification transport
ONU: 2923
Transport Hazard class: 8
Packing group III



Danger
H301-H314-H410
P280-P301+P310a-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338

Hyamine 1622 > RPE - For analysis - Reag. Ph. Eur.

RPE

Description ..... White powder Loss on drying ..... ≤ 5 % Assay (argentimetric) ..... ≥ 96.0 %
Identification ..... Positive Melting point ..... 158 ÷ 163 °C

Code	Size	Packaging	Notes
454921	100 g	Plastic bottle	

Hyamine 1622 solution 0.004M

• Hyamine 1622 soluzione 0.004M • Hyamine 1622 solution 0.004M • Hyamina 1622 solución 0.004M
• Hyamin 1622-Lösung 0.004M

Synonym:
• *Benzethonium chloride*
• *(Diisobutylphenoxyethoxyethyl)dimethylbenzylammonium chloride*

$C_{27}H_{42}ClNO_2$
Molecular Weight: 448,18
CAS: 121-54-0

Hyamine 1622 solution 0.004M > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613000901	100 ml	Plastic bottle	Ref Ph.Eur 3000900
613000900	1 l	Plastic bottle	Ref Ph.Eur 3000900

Hyamine 1622 solution 0.004M > RS - For anionic surfactants detection

RS

Description ..... Clear colourless liquid Identification ..... Positive Assay (at production) ..... 0.0035 ÷ 0.0045 Mol/l

Code	Size	Packaging	Notes
E454972	1 l	Glass bottle	

Hydrazine dihydrochloride

• Idrazina dicloridrato • Hydrazine dichlorhydratée • Hidrazina diclorhidrato • Hydraziniumchlorid

Synonym:
Hydrazine dihydrochloride

$NH_2NH_2.2HCl$
Molecular Weight: 104,97
CAS: 5341-61-7
EEC-N: 226-283-2

Classification transport
ONU: 3288
Transport Hazard class: 6.1
Packing group III



Danger
H301-H311-H331-H317-H350-H410-HA26
P261-P280-P304+P340-P308+P313-P330-P361+P364-P403+P233

Hydrazine dihydrochloride > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Melting point ..... 193.5 ÷ 202.5 ° C
Water-insoluble matter ..... ≤50 ppm
Heavy metals (Pb) ..... ≤5 ppm
Residue on ignition ..... ≤100 ppm
Sulphate ..... ≤10 ppm
Fe ..... ≤5 ppm
Assay (oxidimetric) ..... ≥99 %

Code	Size	Packaging	Notes
455054	100 g	Plastic bottle	
455056	500 g	Plastic bottle	
455057	1 kg	Plastic bottle	

Hydrindantin

• Idrindantina • Hydrindantine • Hidrindantina • Hydrindantine

Synonym:
2,2,'3,3',3'-Hexahydro-2,2'-biindan-1,1'-dione

$C_{18}H_{10}O_6.2H_2O$
Molecular Weight: 358,26
CAS: 5950-69-6
EEC-N: 225-823-4



Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Hydrindantin > RPE - For analysis

RPE

Description ..... Beige powder Identification ..... Positive Melting point ..... ≥ 240 °C Assay ..... ≥ 95 %

Code	Size	Packaging	Notes
455291	5 g	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Hydriodic acid 57%

• Acido iodidrico 57% • Acide iodhydrique 57% • Acido idríodico 57% • Jodwasserstoffsäure 57%

HI
Molecular Weight: 127,92
CAS: 10034-85-2

Classification transport
ONU: 1787
Transport Hazard class: 8
Packing group II



Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Hydriodic acid 57% > RPE - For analysis - ACS

RPE

Description ..... Dark red-brown solution
Identification ..... Positive
Chloride + bromide (Cl) ..... ≤ 500 ppm
Free iodine (I) ..... ≤ 7500 ppm
Heavy metals (Pb) ..... ≤ 10 ppm
Residue on ignition ..... ≤ 100 ppm
Sulfates ..... ≤ 50 ppm
Fe ..... ≤ 10 ppm
Phosphates ..... ≤ 10 ppm

Code	Size	Packaging	Notes
406831	100 ml	Glass bottle	

Stabilized with ~1,5% of hypophosphorous acid

Hydrobromic acid 48%

• Acido bromidrico 48% • Acide bromhydrique 48% • Acido bromhídrico 48% • Bromwasserstoffsäure 48%

HBr
Molecular Weight: 80,92
CAS: 10035-10-6

Classification transport
ONU: 1788
Transport Hazard class: 8
Packing group III



Danger
H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Hydrobromic acid 48% > RPE - For analysis - ACS - ISO

RPE

Description ..... Yellow clear liquid
Identification ..... Positive
Chloride ..... ≤ 200 ppm
Total phosphorus ..... ≤ 2 ppm
Iodide ..... ≤ 20 ppm
Heavy metals (Pb) ..... ≤ 5 ppm
Residue on ignition ..... ≤ 20 ppm
Sulphate ..... ≤ 30 ppm
As ..... ≤ 0.5 ppm
Fe ..... ≤ 1 ppm
Se ..... ≤ 0.01 ppm
Assay (acidimetric) ..... 47.0 - 49.0 %
Cd ..... ≤ 0.5 ppm
Cu ..... ≤ 0.5 ppm
Pb ..... ≤ 0.5 ppm
Zn ..... ≤ 0.5 ppm
Sulphated ash ..... ≤ 0.005 %

Code	Size	Packaging	Notes
402925	250 ml	Glass bottle	
402922	1 l	Bottle	

Hydrochloric acid 50% v/v

• Acido cloridrico 50% v/v • Acide chlorhydrique 50% v/v • Acido clorhídrico 50% v/v • Salzsäure 50% v/v

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group II

Hydrochloric acid 50% v/v > RPE - For analysis

RPE

Code	Size	Packaging	Notes
504571	1 l	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Hydrochloric acid 37%

• Acido cloridrico 37% • Acide chlorhydrique 37% • Acido clorhidrico 37% • Salzsäure 37%

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group II

Danger
H290-H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Hydrochloric acid 37% > RS - VLSI - For electronic use

RS

Code	Size	Packaging	Notes
527601	1 l	Plastic bottle	
527600	2.5 l	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Hydrochloric acid 37% > RS - RSE - For electronic use

RS

Description	Clear liquid	Sulphite	≤0.5 ppm	Cr	≤0.01 ppm	Ni	≤0.01 ppm
Colour (APHA)	≤10	Ag	≤0.02 ppm	Cu	≤0.005 ppm	Pb	≤0.05 ppm
Identification	Positive	Al	≤0.05 ppm	Fe	≤0.1 ppm	Pt	≤0.05 ppm
Assay (acidimetric)	≥36.5 %	As	≤0.005 ppm	Ga	≤0.02 ppm	Sb	≤0.005 ppm
Density at 20° C	1.183 ÷ 1.189	Au	≤0.05 ppm	Hg	≤0.1 ppm	Sn	≤0.02 ppm
Ammonium	≤1 ppm	B	≤0.02 ppm	In	≤0.02 ppm	Sr	≤0.02 ppm
Bromide	≤2 ppm	Ba	≤0.05 ppm	K	≤0.1 ppm	Ta	≤0.1 ppm
Free chlorine	≤0.5 ppm	Be	≤0.02 ppm	Li	≤0.02 ppm	Ti	≤0.05 ppm
Phosphate	≤0.1 ppm	Bi	≤0.02 ppm	Mg	≤0.1 ppm	Tl	≤0.05 ppm
Heavy metals (Pb)	≤0.1 ppm	Ca	≤0.5 ppm	Mn	≤0.01 ppm	V	≤0.05 ppm
Residue on ignition	≤2 ppm	Cd	≤0.005 ppm	Mo	≤0.05 ppm	Zn	≤0.02 ppm
Sulphate	≤0.5 ppm	Co	≤0.01 ppm	Na	≤0.5 ppm	Zr	≤0.05 ppm

Code	Size	Packaging	Notes
403977	1 l	Glass bottle	
403971	2.5 l	Glass bottle	
403974	25 kg	Drum	

Hydrochloric acid 37% > RS - MOS - For electronic use

RS

Description	Clear liquid	Sulphite	≤0.5 ppm	Cr	≤0.01 ppm	Ni	≤0.01 ppm
Colour (APHA)	≤10	Ag	≤0.02 ppm	Cu	≤0.005 ppm	Pb	≤0.05 ppm
Identification	Positive	Al	≤0.05 ppm	Fe	≤0.1 ppm	Pt	≤0.05 ppm
Assay (acidimetric)	≥36.5 %	As	≤0.005 ppm	Ga	≤0.02 ppm	Sb	≤0.005 ppm
Density at 20° C	1.183 ÷ 1.189	Au	≤0.05 ppm	Hg	≤0.1 ppm	Sn	≤0.02 ppm
Ammonium	≤1 ppm	B	≤0.02 ppm	In	≤0.02 ppm	Sr	≤0.02 ppm
Bromide	≤2 ppm	Ba	≤0.05 ppm	K	≤0.1 ppm	Ta	≤0.1 ppm
Free chlorine	≤0.5 ppm	Be	≤0.02 ppm	Li	≤0.02 ppm	Ti	≤0.05 ppm
Phosphate	≤0.1 ppm	Bi	≤0.02 ppm	Mg	≤0.1 ppm	Tl	≤0.05 ppm
Heavy metals (Pb)	≤0.1 ppm	Ca	≤0.5 ppm	Mn	≤0.01 ppm	V	≤0.05 ppm
Residue on ignition	≤2 ppm	Cd	≤0.005 ppm	Mo	≤0.05 ppm	Zn	≤0.02 ppm
Sulphate	≤0.5 ppm	Co	≤0.01 ppm	Na	≤0.5 ppm	Zr	≤0.05 ppm

Code	Size	Packaging	Notes
403942	1 l	Glass bottle	
403941	2.5 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Hydrochloric acid 37% > RPE - For analysis - ISO

RPE

Description ..... Clear liquid	Residue on ignition ..... ≤5 ppm	Co ..... ≤0.01 ppm	Ni ..... ≤0.02 ppm
Colour (APHA) ..... ≤10	Sulphate ..... ≤1 ppm	Cr ..... ≤0.02 ppm	Pb ..... ≤0.05 ppm
Identification ..... Positive	Sulphite ..... ≤0.5 ppm	Cu ..... ≤0.01 ppm	Sr ..... ≤0.02 ppm
Density at 20° C ..... 1.181 ÷ 1.189	Al ..... ≤0.2 ppm	Fe ..... ≤0.2 ppm	Ti ..... ≤1 ppm
Residue on evaporation ..... ≤100 ppm	As ..... ≤0.01 ppm	Hg ..... ≤0.1 ppm	Tl ..... ≤0.05 ppm
Ammonium ..... ≤1 ppm	Ba ..... ≤0.1 ppm	Li ..... ≤0.05 ppm	V ..... ≤ 0.5 ppm
Bromide ..... ≤50 ppm	Be ..... ≤0.02 ppm	Mg ..... ≤0.3 ppm	Zn ..... ≤0.5 ppm
Free chlorine ..... ≤ 4 ppm	Bi ..... ≤0.05 ppm	Mn ..... ≤0.01 ppm	Zr ..... ≤0.05 ppm
Phosphate ..... ≤0.5 ppm	Ca ..... ≤0.5 ppm	Mo ..... ≤0.05 ppm	Assay (acidimetric) ..... ≥36.5 %
Heavy metals (Pb) ..... ≤ 1 ppm	Cd ..... ≤0.005 ppm	Na ..... ≤5 ppm	

Code	Size	Packaging	Notes
403871	1 l	Glass bottle	
403876	1 l	Glass bottle PVC coated	
524525	1 l	Plastic bottle	
403872	2.5 l	Glass bottle	
524526	2.5 l	Plastic bottle	
403878	5 l	Plastic tank	
403874	25 kg	Plastic drum	

Hydrochloric acid 37% > ERBApharm - According to pharmacopoeia: Ph.Eur.-NF-FU-Ph.Franc.-BP-JP

ERBApharm

Description ..... Clear colourless liquid	Free chlorine or bromide .. Conform USP-NF	Residue on evaporation ..... ≤ 100 ppm	Assay (acidimetric) ..... 36.5 ÷ 38.0 %
Identification ..... Positive	Sulphite ..... Conform USP-NF	Residue on ignition ..... ≤ 80 ppm	As ..... ≤ 1 ppm
Appearance of solution ..... Conform Ph.Eur.	Sulphate ..... Conform USP-NF	Free chlorine ..... ≤ 4 ppm	Hg ..... ≤ 0.04 ppm
Bromide or iodide ..... Conform USP-NF	Density at 20° C ..... ~ 1.18	Heavy metals (Pb) ..... ≤ 2 ppm	

Code	Size	Packaging	Notes
302621	1 l	Glass bottle	
302626	2.5 l	Glass bottle	
302643	5 l	Plastic tank	
302624	10 l	Plastic tank	
302623	25 kg	Plastic drum	
302622	40 kg	Plastic tank	
302627	55 kg	Plastic tank	
302625	220 kg	Plastic drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Hydrochloric acid 37% > RE - Pure

RE

Assay ..... 36.0 ÷ 38.5 %	Residue on ignition ..... ≤ 20 ppm	Heavy metals ..... ≤ 5 ppm	Free chloride (Cl) ..... ≤ 20 ppm
Description ..... Clear colourless liquid or yellowish	Density at 20°C ..... 1.183 ÷ 1.189	Assay (acidimetric) ..... 36.5 - 37.5 %	Iron (Fe) ..... ≤ 2 ppm

Code	Size	Packaging	Notes
528525	5 l	Tank	
303871	25 kg	Plastic tank	

Hydrochloric acid 34-37%

• Acido cloridrico 34-37% • Acide chlorhydrique 34-37% • Acido clorhidrico 34-37% • Salzsäure 34-37%

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group II

Danger
H290-H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Hydrochloric acid 34-37% > RS - Superpure - For trace analysis at ppb level

RS

Description	Clear liquid	Fe	≤ 1 ppb	Assay (acidimetric)	34 ÷ 37 %	Nb	≤ 0.1 ppb
Identification	Positive	Hg	≤ 0.1 ppb	Th	≤ 0.1 ppb	Pr	≤ 0.1 ppb
Bromide	≤ 10 ppm	K	≤ 1 ppb	U	≤ 0.1 ppb	Re	≤ 0.1 ppb
Free chlorine	≤ 0.5 ppm	Li	≤ 0.1 ppb	Colour (APHA)	≤ 10	Rh	≤ 0.1 ppb
Total phosphorus	≤ 0.01 ppm	Mg	≤ 0.5 ppb	Ce	≤ 0.1 ppb	Rb	≤ 0.1 ppb
Total sulphur	≤ 0.3 ppm	Mn	≤ 0.1 ppb	Cs	≤ 0.1 ppb	Ru	≤ 0.1 ppb
Ag	≤ 1 ppb	Mo	≤ 0.1 ppb	Dy	≤ 0.1 ppb	Sm	≤ 0.1 ppb
Al	≤ 1 ppb	Na	≤ 1 ppb	Er	≤ 0.1 ppb	Sc	≤ 0.1 ppb
As	≤ 0.5 ppb	Ni	≤ 0.5 ppb	Eu	≤ 0.1 ppb	Te	≤ 0.1 ppb
B	≤ 1 ppb	Pb	≤ 0.1 ppb	Gd	≤ 0.1 ppb	Tb	≤ 0.1 ppb
Ba	≤ 0.1 ppb	Sb	≤ 0.5 ppb	Ga	≤ 0.1 ppb	Tl	≤ 0.1 ppb
Be	≤ 0.1 ppb	Se	≤ 1 ppb	Au	≤ 0.5 ppb	Tm	≤ 0.1 ppb
Bi	≤ 0.1 ppb	Sn	≤ 0.5 ppb	Hf	≤ 0.1 ppb	W	≤ 0.1 ppb
Ca	≤ 1 ppb	Sr	≤ 0.1 ppb	Ho	≤ 0.1 ppb	Yb	≤ 0.1 ppb
Cd	≤ 0.1 ppb	Ti	≤ 0.5 ppb	In	≤ 0.1 ppb	Y	≤ 0.1 ppb
Co	≤ 0.1 ppb	V	≤ 0.5 ppb	La	≤ 0.1 ppb		
Cr	≤ 0.5 ppb	Zn	≤ 1 ppb	Lu	≤ 0.1 ppb		
Cu	≤ 0.5 ppb	Zr	≤ 0.1 ppb	Nd	≤ 0.1 ppb		

Code	Size	Packaging	Notes
403915	500 ml	Plastic bottle	
403916	1 l	Plastic bottle	
403917	2.5 l	Plastic bottle	

Hydrochloric acid 32-35%

• Acido cloridrico 32-35% • Acide chlorhydrique 32-35% • Acido clorhidrico 32-35% • Salzsäure 32-35%

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group II

Danger
H290-H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Hydrochloric acid 32-35% > RS - Ultrapure - For trace analysis at ppt level

RS

Description	Clear colourless liquid	Cs	≤ 10 ppt	U	≤ 1 ppt	Rh	≤ 10 ppt
Identification	Positive	Cu	≤ 10 ppt	Th	≤ 1 ppt	Rb	≤ 10 ppt
Ag	≤ 10 ppt	Dy	≤ 1 ppt	Sb	≤ 20 ppt	Ru	≤ 10 ppt
Al	≤ 20 ppt	Mn	≤ 10 ppt	Gd	≤ 1 ppt	Sm	≤ 1 ppt
As	≤ 50 ppt	Eu	≤ 1 ppt	Ga	≤ 10 ppt	Sc	≤ 10 ppt
Au	≤ 50 ppt	Fe	≤ 10 ppt	Hf	≤ 10 ppt	Te	≤ 1 ppt
B	≤ 100 ppt	Ni	≤ 20 ppt	Ho	≤ 1 ppt	Tb	≤ 1 ppt
Ba	≤ 10 ppt	Pb	≤ 10 ppt	In	≤ 1 ppt	Tm	≤ 1 ppt
Be	≤ 10 ppt	Sn	≤ 20 ppt	La	≤ 1 ppt	W	≤ 10 ppt
Bi	≤ 10 ppt	Sr	≤ 10 ppt	Li	≤ 10 ppt	Yb	≤ 1 ppt
Ca	≤ 10 ppt	Tl	≤ 10 ppt	Lu	≤ 10 ppt	Y	≤ 1 ppt
Cd	≤ 10 ppt	Ti	≤ 10 ppt	Nd	≤ 1 ppt	Zr	≤ 10 ppt
Ce	≤ 10 ppt	V	≤ 10 ppt	Nb	≤ 1 ppt		
Co	≤ 10 ppt	Zn	≤ 10 ppt	Pr	≤ 1 ppt		
Cr	≤ 10 ppt	Assay (acidimetric)	32 ÷ 35 %	Re	≤ 10 ppt		

Code	Size	Packaging	Notes
403891	500 ml	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Hydrochloric acid 32% (20°Bé)

• Acido cloridrico 32% (20°Bé) • Acide chlorhydrique 32% (20°Bé) • Acido clorhidrico 32% (20°Bé) • Salzsäure 32% (20°Be)

HCl
Molecular Weight: 36,461
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group II

Danger
H290-H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Hydrochloric acid 32% (20°Bé) > RE - Pure

RE

Description ..... Yellow clear liquid Density at 20° C .....≥1.154 Assay (acidimetric) .....≥31.0 %

Code	Size	Packaging	Notes
302664	25 kg	Plastic drum	(20°Bè)

Hydrochloric acid 32%

• Acido cloridrico 32% • Acide chlorhydrique 32% • Acido clorhidrico 32% • Salzsäure 32%

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group II

Danger
H290-H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Hydrochloric acid 32% > RPE - For analysis - ISO

RPE

Description .....Clear liquid
Colour (APHA) ..... ≤10
Identification .....Positive
Density at 20° C .....1.159 - 1.164
Ammonium .....≤1 ppm
Bromide.....≤50 ppm
Free chlorine.....≤ 4 ppm
Phosphate .....≤0.5 ppm
Heavy metals (Pb).....≤1 ppm
Residue on ignition.....≤5 ppm

Sulphate.....≤0.8 ppm
Sulphite.....≤0.5 ppm
Al.....≤0.05 ppm
As.....≤0.01 ppm
Ba.....≤0.1 ppm
Be.....≤0.02 ppm
Bi.....≤0.05 ppm
Ca.....≤0.5 ppm
Cd.....≤0.005 ppm
Co.....≤0.01 ppm

Cr.....≤0.02 ppm
Cu.....≤0.01 ppm
Fe .....≤0.1 ppm
Hg.....≤0.1 ppm
Li .....≤0.05 ppm
Mg .....≤0.3 ppm
Mn .....≤0.01 ppm
Mo .....≤0.05 ppm
Na.....≤5 ppm
Ni.....≤0.02 ppm

Pb.....≤0.05 ppm
Sr.....≤0.02 ppm
Ti .....≤1 ppm
Tl .....≤0.05 ppm
V.....≤0.5 ppm
Zn .....≤0.5 ppm
Zr.....≤0.05 ppm
Assay (acidimetric) .....32.0 - 33.0 %

Code	Size	Packaging	Notes
403981	2.5 l	Glass bottle	
403984	2.5 l	Plastic bottle	
403982	25 kg	Plastic drum	
403986	55 kg	Plastic tank	
403988	220 kg	Plastic drum	

Hydrochloric acid 32% > RE - Pure

RE

Description ...Clear colouless or light yellow liquid
Density at 20°C .....1.157 - 1.171

Free chlorine.....≤100 ppm
Heavy metals (Pb).....≤50 ppm

Residue on ignition..... ≤0.1 %
Sulphate.....≤200 ppm

Fe .....≤2 ppm
Assay .....32 ÷ 34 %

Code	Size	Packaging	Notes
302601	1 l	Glass bottle	
302602	25 kg	Plastic drum	
302604	30 kg	Plastic drum	

Hydrochloric acid 29-31%

• Acido cloridrico 29-31 % • Acide chlorhydrique 29-31 % • Acido clorhidrico 29-31 % • Salzsäure 29-31 %

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group II



Danger
H290-H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Hydrochloric acid 29-31% > RS - Superpure - For trace analysis at ppb level — RS

Code	Size	Packaging	Notes
403921	1 l	Plastic bottle	

For specifications, contact our customer service for a certificate of analysis

Hydrochloric acid 26%

• Acido cloridrico 26% • Acide chlorhydrique 26% • Acido clorhidrico 26% • Salzsaure 26%

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group II



Danger
H290-H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Hydrochloric acid 26% > RE - Pure — RE

HCl content..... 25 - 27 %

Code	Size	Packaging	Notes
PS0769/20	2.5 l	Plastic bottle	

Hydrochloric acid 25% w/v

• Acido cloridrico 25% m/v • Acide chlorhydrique 25% m/v • Acido clorhidrico 25% p/v • Salzsäure 25% m/v

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III



Warning
H290-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Hydrochloric acid 25% w/v > RS - For analysis according to Ph. Eur. Chap. 4.1.1 — RS

Code	Size	Packaging	Notes
611043501	1 l	Plastic bottle	Hydrochloric acid R1 Ref Ph.Eur 1043501

Product specifications are subject to changes. Please visit our website for updates.

Hydrochloric acid 23%

• Acido cloridrico 23% • Acide chlorhydrique 23% • Acido clorhidrico 23% • Salzsäure 23%

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III



Warning
H290-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Hydrochloric acid 23% > RPE - For analysis - ISO

RPE

Description	Clear liquid	Sulphate	≤0.8 ppm	Cu	≤0.01 ppm	Sr	≤0.02 ppm
Colour (APHA)	≤10	Sulphite	≤0.5 ppm	Fe	≤0.1 ppm	Ti	≤1 ppm
Identification	Positive	Al	≤0.05 ppm	Hg	≤0.1 ppm	Tl	≤0.05 ppm
Density at 20° C	1.113 - 1.119	As	≤0.01 ppm	Li	≤0.05 ppm	V	≤0.5 ppm
Ammonium	≤1 ppm	Be	≤0.02 ppm	Mg	≤0.3 ppm	Zn	≤0.5 ppm
Bromide	≤50 ppm	Bi	≤0.05 ppm	Mn	≤0.03 ppm	Zr	≤0.05 ppm
Free chlorine	≤4 ppm	Ca	≤0.5 ppm	Mo	≤0.05 ppm	Assay (acidimetric)	23.0 - 24.0 %
Phosphate	≤0.5 ppm	Cd	≤0.005 ppm	Na	≤5 ppm		
Heavy metals (Pb)	≤1 ppm	Co	≤0.01 ppm	Ni	≤0.02 ppm		
Residue on ignition	≤5 ppm	Cr	≤0.02 ppm	Pb	≤0.05 ppm		

Code	Size	Packaging	Notes
403901	1 l	Glass bottle	
403905	2.5 l	Glass bottle	
403909	25 kg	Plastic drum	

Hydrochloric acid 20%

• Acido cloridrico 20% • Acide chlorhydrique 20% • Acido clorhidrico 20% • Salzsäure 20%

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III



Warning
H290-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Hydrochloric acid 20% > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Assay ..... 19.50 - 20.50 %

Code	Size	Packaging	Notes
524561	10 l	Plastic tank	

Hydrochloric acid 20% > RE - Pure

RE

Assay ..... 19 - 21 %

Code	Size	Packaging	Notes
PS0751/29	5 l	Plastic tank	

Hydrochloric acid 12%

• Acido cloridrico 12% • Acide chlorhydrique 12% • Acido clorhidrico 12% • Salzsäure 12%

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III



Warning
H290-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Hydrochloric acid 12% > RS - For analysis

RS

Assay ..... 11.0 - 13.0 %
Density d20/4 ..... 1.052 - 1.063
Clear liquid appearance ..... Conform
Identification ..... Conform
Colour ..... ≤ J7

Code	Size	Packaging	Notes
PS0347/22	5 l	Plastic tank	
PS0347/49	25 l	Plastic tank	
PS0347/66	200 l	Plastic drum	

Hydrochloric acid 10%

• Acido cloridrico 10% • Acide chlorhydrique 10% • Acido clorhidrico 10% • Salzsäure 10%

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III



Warning
H290-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Hydrochloric acid 10% > ERBApharm - According to pharmacopoeia: Ph.Eur.

ERBApharm

Description ..... Clear colourless liquid
Identification ..... Positive
Appearance ..... Conform Ph.Eur.
Free chlorine ..... ≤ 1 ppm
Sulphate ..... ≤ 5 ppm
Heavy metals (Pb) ..... ≤ 2 ppm
Residue on evaporation ..... ≤ 100 ppm
Assay ..... 9.5 - 10.5 % (m/m)
Residual solvents (Current ICH) ..... Conform

Code	Size	Packaging	Notes
302591	10 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Hydrochloric acid 10% > RE - Pure

RE

Hydrochloric acid content ..... 9.9 - 10.1 %

Code	Size	Packaging	Notes
PS0768/41	10 l	Plastic tank	

Hydrochloric acid 8%

• Acido cloridrico 8% • Acide chlorhydrique 8% • Acido clorhidrico 8% • Salzsäure 8%

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III



Warning
H290
P234-P390-P406

Hydrochloric acid 8% > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Density at 20° C ..... 1.037 - 1.043
Ammonium ..... ≤30 ppm
Free chlorine ..... ≤0.5 ppm
Heavy metals (Pb) ..... ≤1 ppm
Sulphate ..... ≤1 ppm
Sulphite ..... ≤1 ppm
As ..... ≤0.01 ppm
Fe ..... ≤0.3 ppm
Assay (acidimetric) ..... 8.0 - 9.0 %

Code	Size	Packaging	Notes
404033	10 kg	Plastic tank	
404036	55 kg	Plastic tank	

Hydrochloric acid 5%

• Acido cloridrico 5% • Acide chlorhydrique 5% • Acido cloridrico 5% • Salzsäure 5%

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III



Warning
H290
P234-P390-P406

Hydrochloric acid 5% > ERBApharm - Prepared from raw material according Ph.Eur

ERBApharm

Hydrochloric acid content ..... 4.5 - 5.5 %

Code	Size	Packaging	Notes
PS0864/41	10 l	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Product specifications are subject to changes. Please visit our website for updates.

Hydrochloric acid 1.128% m/v

• Acido cloridrico 1.128% m/v • Acide chlorhydrique 1.128% m/v • Acido clorhidrico 1.128% p/v • Salzsäure 1.128% m/v

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III



Warning
H290
P234-P390-P406

Hydrochloric acid 1.128% m/v > RS - For agroalimentary analysis — RS

Description ..... Clear colourless liquid Assay .....1.123 ÷ 1.133 % m/v

Code	Size	Packaging	Notes
502761	1 l	Plastic bottle	

Hydrochloric Acid 9 mol/l (9N)

• Acido clorhidrico 9 mol/l (9N) • Acide chlorhydrique 9 mol/l (9N) • Acido cloridrico 9 mol/l (9N) • Salzsäure 9 mol/l (9N)

HCl
Molecular Weight: 36,46
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group II



Danger
H290-H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Hydrochloric Acid 9 mol/l (9N) > RS - For analysis — RS

Assay (potentiometry) ..... 8.982 - 9.018 N

Code	Size	Packaging	Notes
PS0313/20	2.5 l	Plastic bottle	

Hydrochloric acid 6 mol/l (6N)

• Acido cloridrico 6 mol/l (6N) • Acide chlorhydrique 6 mol/l (6N) • Acido clorhidrico 6 mol/l (6N) • Salzsäure 6 mol/l (6N)

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III



Warning
H290-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Hydrochloric acid 6 mol/l (6N) > RS - For agroalimentary analysis — RS

Description ..... Clear colourless liquid Assay .....5.97 ÷ 6.03 N Colour .....≤ 10 Hazen

Code	Size	Packaging	Notes
502831	1 l	Plastic bottle	
502832	18 l	Plastic tank	

Store between 15-30 °C

Hydrochloric acid 6 mol/l (6N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2 — RS

Code	Size	Packaging	Notes
613001500	1 l	Plastic bottle	Ref Ph.Eur 3001500

Hydrochloric acid 6 mol/l (6N) > ERBApharm - Prepared from raw material according Ph.Eur — ERBApharm

Assay (potentiometric) .....5.9 ÷ 6.1 mol/L

Code	Size	Packaging	Notes
528651	25 l	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Hydrochloric acid 6 mol/l (6N) > RE - Pure — RE

Description ..... Clear colourless liquid Colour .....≤ 10 Hazen Assay .....5.95 ÷ 6.05 N

Code	Size	Packaging	Notes
528550000	5 l	Plastic tank	

Hydrochloric acid 5 mol/l (5N)

• Acido cloridrico 5 mol/l (5N) • Acide chlorhydrique 5 mol/l (5N) • Acido clorhidrico 5 mol/l (5N) • Salzsäure 5 mol/l (5N)

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III



Warning
H290-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Hydrochloric acid 5 mol/l (5N) > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Assay (potentiometric) ..... 4.99 - 5.01 mol/L

Code	Size	Packaging	Notes
528731	25 l	Plastic tank	

Hydrochloric acid 5 mol/l (5N) > RPE - For analysis

RPE

Assay (potentiomety) ..... 4.995 - 5.005 N

Code	Size	Packaging	Notes
P3160015	1 l	Plastic bottle	
P3160095	5 l	Kubidos	

Hydrochloric acid 4 mol/l (4N)

• Acido cloridrico 4 mol/l (4N) • Acide chlorhydrique 4 mol/l (4N) • Acido clorhidrico 4 mol/l (4N) • Salzsäure 4 mol/l (4N)

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III



Warning
H290-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Hydrochloric acid 4 mol/l (4N) > RS - For agroalimentary analysis

RS

Description ..... Clear colourless liquid Assay ..... 3.992 ÷ 4.008 N

Code	Size	Packaging	Notes
502010	1 l	Plastic bottle	

145,84 g of HCl. Ready-to-use. Content is guaranteed for standardized volumes at 20 °C. Keep tightly closed

Hydrochloric acid 4 mol/l (4N) > RPE - For analysis

RPE

Assay (potentiomety) ..... 3.992 - 4.008 N

Code	Size	Packaging	Notes
PS0589/15	1 l	Plastic bottle	
PS0589/22	5 l	Plastic tank	
PS0589/49	25 l	Plastic tank	

Hydrochloric acid 4 mol/l (4N) > ERBApharm - Prepared from raw material according Ph.Eur

ERBApharm

Identification (Ph.Eur)..... Conform Assay (Ph.Eur)..... 3.8 - 4.2 N Origine (BSE-TSE)..... Conform Residual solvents (Current ICH)..... Conform

Code	Size	Packaging	Notes
528681	1 l	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Product specifications are subject to changes. Please visit our website for updates.

Hydrochloric acid 3 mol/l (3N)

• Acido cloridrico 3 mol/l (3N) • Acide chlorhydrique 3 mol/l (3N) • Acido clorhidrico 3 mol/l (3N) • Salzsäure 3 mol/l (3N)

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III



Warning
H290-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Hydrochloric acid 3 mol/l (3N) > RS - For agroalimentary analysis — RS

Description ..... Clear colourless liquid Assay ..... 2.95 ÷ 3.05 N

Code	Size	Packaging	Notes
502621	1 l	Plastic bottle	
502622	2.5 l	Plastic bottle	
502011	25 l	Plastic tank	

Hydrochloric acid 2 mol/l (2N)

• Acido cloridrico 2 mol/l (2N) • Acide chlorhydrique 2 mol/l (2N) • Acido clorhidrico 2 mol/l (2N) • Salzsäure 2 mol/l (2N)

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III



Warning
H290
P234-P390-P406

Hydrochloric acid 2 mol/l (2N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2 — RS

Code	Size	Packaging	Notes
613001700	1 l	Plastic bottle	Ref Ph.Eur 3001700

Hydrochloric acid 2 mol/l (2N) > RPE - For analysis — RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 1.998 - 2.002 N NIST 723 ..... e

Code	Size	Packaging	Notes
404067000	1 l	Plastic bottle	Certified with NIST traceability
404062000	5 l	Kubidos	Certified with NIST traceability
404061000	10 l	Kubidos	Certified with NIST traceability

72.92 g of HCl. Volumetric solution ready-to-use

Hydrochloric acid 2 mol/l (2N) > ERBApharm - Prepared from raw material according Ph.Eur — ERBApharm

Identification (Ph.Eur) ..... Conform Assay (Ph.Eur) ..... 1.9 -2.1 N Origine (BSE-TSE) ..... Conform Residual solvents (Current ICH)

Code	Size	Packaging	Notes
528691	1 l	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Hydrochloric acid 1 mol/l (1N)

• Acido cloridrico 1 mol/l (1N) • Acide chlorhydrique 1 mol/l (1N) • Acido clorhidrico 1 mol/l (1N) • Salzsäure 1 mol/l (1N)

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III



Warning
H290
P234-P390-P406

Hydrochloric acid 1 mol/l (1N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2 — RS

Code	Size	Packaging	Notes
613001801	500 ml	Plastic bottle	Ref Ph.Eur 3001800
613001800	1 l	Plastic bottle	Ref Ph.Eur 3001800

Hydrochloric acid 1 mol/l (1N) > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
617000191	1 l	Plastic bottle	

Hydrochloric acid 1 mol/l (1N) > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.998 - 1.002 N

Code	Size	Packaging	Notes
528673	5 l	Kubidos	
528671	10 l	Kubidos	
528672	200 l	Polythene-metal drum	

Content is guaranteed for standardized volumes at 20°C. Keep tightly closed

Hydrochloric acid 1 mol/l (1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.999 - 1.001 N NIST 723 ..... e

Code	Size	Packaging	Notes
404097000	1 l	Plastic bottle	Certified with NIST traceability
404092000	5 l	Kubidos	Certified with NIST traceability
404091000	10 l	Kubidos	Certified with NIST traceability
404094000	20 l	Plastic tank	Certified with NIST traceability

36.46 g of HCl. Volumetric solution ready-to-use

Hydrochloric acid 1 mol/l (1N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
404111		Plastic ampoule	Volume: 165 ml

36,46 g of HCl. Volumetric concentrated solution to prepare 1 L of solution 1 N

Hydrochloric acid 1 mol/l (1N) > ERBApharm - Prepared from raw material according Ph.Eur

ERBApharm

Clear,colourless liquid ..... Conform Assay (Ph.Eur) ..... 0.95 - 1.05 N Residual solvents (Current ICH) ..... Conform
Identification (Ph.Eur) ..... Conform Origine (BSE-TSE) ..... Conform

Code	Size	Packaging	Notes
528583	1 l	Plastic bottle	
528584	5 l	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Hydrochloric acid 0.714 mol/l (N/1.4)

• Acido cloridrico 0.714 mol/l (N/1.4) • Acide chlorhydrique 0.714 mol/l (N/1.4) • Acido clorhídrico 0.714 mol/l (N/1.4) • Salzsäure 0.714 mol/l (N/1.4)

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III



Warning
H290
P234-P390-P406

Hydrochloric acid 0.714 mol/l (N/1.4) > RS - For agroalimentary analysis

RS

Description ..... Clear colourless liquid Assay ..... 0.710 ÷ 0.718 N

Code	Size	Packaging	Notes
526531	10 l	Kubidos	

Content is guaranteed for standardized volumes at 20°C. Keep tightly closed

Product specifications are subject to changes. Please visit our website for updates.

Hydrochloric acid 0.5 mol/l (0.5N)

• Acido cloridrico 0.5 mol/l (0.5N) • Acide chlorhydrique 0.5 mol/l (0.5N) • Acido clorhidrico 0.5 mol/l (0.5N) • Salzsäure 0.5 mol/l (0.5N)

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III

Warning
H290
P234-P390-P406

Hydrochloric acid 0.5 mol/l (0.5N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.4995 - 0.5005 N NIST 723 ..... e

Code	Size	Packaging	Notes
404147000	1 l	Plastic bottle	Certified with NIST traceability
404142000	5 l	Kubidos	Certified with NIST traceability
404141000	10 l	Kubidos	Certified with NIST traceability

18.23 g of HCl. Volumetric solution ready-to-use

Hydrochloric acid 0.5 mol/l (0.5N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
404161		Plastic ampoule	Volume: 165 ml

18,230 g HCl. Volumetric concentrated solution to prepare 1 L of solution 0,5 N

Hydrochloric acid 0.2 mol/l (0.2N)

• Acido cloridrico 0.2 mol/l (0.2N) • Acide chlorhydrique 0.2 mol/l (0.2N) • Acido clorhidrico 0.2 mol/l (0.2N) • Salzsäure 0.2 mol/l (0.2N)

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III

Warning
H290
P234-P390-P406

Hydrochloric acid 0.2 mol/l (0.2N) > RS - For agroalimentary analysis

RS

Description ..... Clear colourless liquid Assay ..... 0.195 ÷ 0.205 N

Code	Size	Packaging	Notes
502631	1 l	Plastic bottle	

7,292 g of HCl. Ready-to-use solution according to NF V04-242. Content is guaranteed for standardized volumes at 20 °C. Keep tightly closed

Hydrochloric acid 0.2 mol/l (0.2N) in propanol-2

• Acido cloridrico 0.2 mol/l (0.2N) in isopropanolo • Acide chlorhydrique 0.2 mol/l (0.2N) dans propanol-2 • Acido clorhidrico 0.2 mol/l (0.2N) en propan-2-ol • Salzsäure 0.2 mol / l (0.2 N) in Propanol-2

HCl
Molecular Weight: 36,46
CAS: 7647-01-0

Classification transport
ONU: 2924
Transport Hazard class: 3
Packing group II

Danger
H225-H290-H319-H336
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Hydrochloric acid 0.2 mol/l (0.2N) in propanol-2 > RPE - For analysis

RPE

Description ..... Clear colourless liquid Color ..... ≤ 10 APHA Assay ..... 0,195 ÷ 0,205 N

Code	Size	Packaging	Notes
526535	1 l	Glass bottle	

Hydrochloric acid 0.1 mol/l (0.1N)

• Acido cloridrico 0.1 mol/l (0.1N) • Acide chlorhydrique 0.1 mol/l (0.1N) • Acido clorhidrico 0.1 mol/l (0.1N) • Salzsäure 0.1 mol/l (0.1N)

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III

Warning
H290
P234-P390-P406

Hydrochloric acid 0.1 mol/l (0.1N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2 — RS

Code	Size	Packaging	Notes
613002101	500 ml	Plastic bottle	Ref Ph.Eur 3002100
613002100	1 l	Plastic bottle	Ref Ph.Eur 3002100

Hydrochloric acid 0.1 mol/l (0.1N) > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP — RPE

Description ..... Clear colourless liquid Assay (potentiometric) ..... 0.0998 - 0.1002 mol/L

Code	Size	Packaging	Notes
528573	5 l	Kubidos	
528571	10 l	Kubidos	

Content is guaranteed for standardized volumes at 20°C. Keep tightly closed

Hydrochloric acid 0.1 mol/l (0.1N) > RPE - For analysis — RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.0999 - 0.1001 N NIST 723 ..... e

Code	Size	Packaging	Notes
404197000	1 l	Plastic bottle	Certified with NIST traceability
404192000	5 l	Kubidos	Certified with NIST traceability
404191000	10 l	Kubidos	Certified with NIST traceability
404195000	10 l	Plastic tank	Certified with NIST traceability

3.646 g of HCl. Volumetric solution ready-to-use: 0.1N

Hydrochloric acid 0.1 mol/l (0.1N) > RPE - NORMEX - For analysis — RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
404211		Plastic ampoule	Volume: 55 ml

3,646 g of HCl. Volumetric concentrated solution to prepare 1 L of solution 0,1 N

Hydrochloric acid 0.1 mol/l (0.1N) > ERBApharm - Prepared from raw material according Ph.Eur. - USP — ERBApharm

Clear,colourless liquid ..... Conform Assay (Ph.Eur) ..... 0.095 - 0.105 N Residual solvents (Current ICH) ..... Conform
Identification (Ph.Eur) ..... Conform Origine (BSE-TSE) ..... Conform

Code	Size	Packaging	Notes
528661	1 l	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Hydrochloric acid 0.1 mol/l (0.1N) > ERBApharm - Prepared from raw material according Ph.Eur — ERBApharm

Clear,colourless liquid ..... Conform Assay (Ph.Eur) ..... 0.095 - 0.105 N Residual solvents (Current ICH) ..... Conform
Identification (Ph.Eur) ..... Conform Origine (BSE-TSE) ..... Conform

Code	Size	Packaging	Notes
528662	5 l	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Product specifications are subject to changes. Please visit our website for updates.

Hydrochloric acid 0.1 mol/l (0.1N) in ethanol

• Acido cloridrico 0.1 mol/l (0.1N) in etanolo • Acide chlorhydrique 0.1 mol/l (0.1N) dans l'éthanol • Acido clorhidrico 0.1 mol/l (0.1N) en etanol • Salzsäure 0.1 mol / l (0.1 N) in Ethanol

HCl
Molecular Weight: 36,46
CAS: 7647-01-0

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II

Danger
H225
P210-P241-P280-P303+P361+P353-P403+P235-P501a

Hydrochloric acid 0.1 mol/l (0.1N) in ethanol > RS - For analysis according to Ph. Eur. Chap. 4.2.2 — RS

Code	Size	Packaging	Notes
613008800	1 l	Glass bottle	Ref Ph.Eur 3008800

Hydrochloric acid 0.1 mol/l (0.1N) in propanol-2

• Acido cloridrico 0.1 mol/l (0.1N) in isopropanolo • Acide chlorhydrique 0.1 mol/l (0.1N) dans propanol-2 • Acido clorhidrico 0.1 mol/l (0.1N) en propan-2-ol • Salzsäure 0.1 mol / l (0.1 N) in Propanol-2

HCl
Molecular Weight: 36,46
CAS: 7647-01-0

Classification transport
ONU: 1219
Transport Hazard class: 3
Packing group II

Danger
H225-H319-H336
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Hydrochloric acid 0.1 mol/l (0.1N) in propanol-2 > RPE - For analysis — RPE

Description ..... Clear liquid Color ..... ≤ 10 APHA Assay ..... 0,095 ÷ 0.105 N

Code	Size	Packaging	Notes
526536	1 l	Plastic bottle	

Hydrochloric acid 0.0714 mol/l (N/14)

• Acido cloridrico 0.0714 mol/l (N/14) • Acide chlorhydrique 0.0714 mol/l (N/14) • Acido clorhídrico 0.0714 mol/l (N/14) • Salzsäure 0.0714 mol/l (N/14)

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III

Warning
H290
P234-P390-P406

Hydrochloric acid 0.0714 mol/l (N/14) > RS - For agroalimentary analysis — RS

Description ..... Clear colourless liquid Assay ..... 0.0710 ÷ 0.0718 N

Code	Size	Packaging	Notes
526533	10 l	Kubidos	

Content is guaranteed for standardized volumes at 20°C. Keep tightly closed

Hydrochloric acid 0.05 mol/l (0.05N)

• Acido cloridrico 0.05 mol/l (0.05N) • Acide chlorhydrique 0.05 mol/l (0.05N) • Acido clorhidrico 0.05 mol/l (0.05N) • Salzsäure 0.05 mol / l (0.05 N)

HCl
Molecular Weight: 36.46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III

Warning
H290
P234-P390-P406

Hydrochloric acid 0.05 mol/l (0.05N) > RPE - For analysis — RPE

Assay (potentiometry) ..... 0.0499 - 0.0501 N

Code	Size	Packaging	Notes
PS0587/15	1 l	Plastic bottle	

Hydrochloric acid 0.04 mol/l (0.04N)

• Acido cloridrico 0.04 mol/l (0.04N) • Acide chlorhydrique 0.04 mol/l (0.04N) • Acido clorhídrico 0.04 mol/l (0.04N) • Salzsäure 0.04 mol/l (0.04N)

HCl
Molecular Weight: 36.46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III

Warning
H290
P234-P390-P406

Hydrochloric acid 0.04 mol/l (0.04N) > RS - For analysis — RS

Assay (potentiomety) ..... 0.03992 - 0.04008 N

Code	Size	Packaging	Notes
PS0206/41	10 l	Plastic tank	

Hydrochloric acid 0.02 mol/l (0.02N)

• Acido cloridrico 0.02 mol/l (0.02N) • Acide chlorhydrique 0.02 mol/l (0.02N) • Acido clorhídrico 0.02 mol/l (0.02N) • Salzsäure 0.02 mol/l (0.02N)

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Hydrochloric acid 0.02 mol/l (0.02N) > RPE - For analysis — RPE

Description ..... Clear colourless liquid Assay (acidimetric) ..... 0.01996 ÷ 0.02004 mol/L

Code	Size	Packaging	Notes
PS0342/15	1 l	Plastic bottle	
526537	5 l	Plastic tank	

Hydrochloric acid 0.01 mol/l (0.01N)

• Acido cloridrico 0.01 mol/l (0.01N) • Acide chlorhydrique 0.01 mol/l (0.01N) • Acido clorhídrico 0.01 mol/l (0.01N) • Salzsäure 0.01 mol/l (0.01N)

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Hydrochloric acid 0.01 mol/l (0.01N) > RPE - For analysis — RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.00998 - 0.01002 N NIST 84 l

Code	Size	Packaging	Notes
404267	1 l	Plastic bottle	Certified with NIST traceability

Provided with a Certificate of analysis with references on the analytical method

Hydrochloric acid 0.01 mol/l (0.01N) > RPE - NORMEX - For analysis — RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
404251		Plastic ampoule	Volume: 55 ml

0,3646 g HCl. Volumetric concentrated solution to prepare 1 L of solution 0,01 N

Product specifications are subject to changes. Please visit our website for updates.

Hydrochloric acid, dilute

• Acido cloridrico diluito • Acide chlorhydrique dilué • Acido clorhídrico diluido • Salzsäure verdünnen

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

HEU210

Hydrochloric acid, dilute > RS - For analysis according to Ph. Eur. Chap. 2.2.2

RS

Code	Size	Packaging	Notes
612202400	1 l	Plastic bottle	Dilution matrix HCl 10g/L

Hydrochloric acid, dilute > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611043503	1 l	Plastic bottle	Ref Ph.Eur 1043503
611043504	1 l	Plastic bottle	Hydrochloric acid, dilute R1 Ref Ph.Eur 1043504

Hydrochloric acid, brominated

• Acido cloridrico bromurato • Acide chlorhydrique bromé • Acido clorhídrico brominado • Salzsäure, bromiert

HCl
Molecular Weight: 36,46
CAS: 7647-01-0
EEC-N: 231-595-7

Classification transport
ONU: 2922
Transport Hazard class: 8
Packing group I



Danger
H290-H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Hydrochloric acid, brominated > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611043507	1 l	Glass bottle	Ref Ph.Eur 1043507

Hydrochloric acid-d 20%

• Acido cloridrico-d 20% • Acide chlorhydrique-d 20% • Acido clorhidrico-d 20% • Deuteriumchlorid 20%

Synonym:
• *Deuterium chloride*
• *Deutero-hydrochloric acid*

DCl
Molecular Weight: 37,47
CAS: 7698-05-7

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group II

Hydrochloric acid-d 20% > RS - For NMR - min 99.95%

RS

Code	Size	Packaging	Notes
P5685	25 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Hydrochloric acid-d 1 mol/l

• Acido cloridrico-d 1 mol/l • Acide chlorhydrique-d 1 mol/l • Acido clorhidrico-d 1 mol/l • Deuteriumchlorid 1 mol/l

DCl
Molecular Weight: 37,47
CAS: 7698-05-7
EEC-N: 231-715-8

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III



Warning
H290-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Hydrochloric acid-d 1 mol/l > RS - For NMR - min 99.8%

RS

Code	Size	Packaging	Notes
P5695	25 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Hydrofluoric acid 50%

• Acido fluoridrico 50% • Acide fluorhydrique 50% • Acido fluorhídrico 50% • Fluorwasserstoff 50%

HF
Molecular Weight: 20
CAS: 7664-39-3

Classification transport
ONU: 1790
Transport Hazard class: LQ



Danger
H300-H310-H330-H314
P280-P284-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Hydrofluoric acid 50% > RS - RSE - For electronic use

RS

Description	Clear liquid	Sulphate	≤1 ppm	Co	≤0.01 ppm	Ni	≤0.05 ppm
Colour (APHA)	≤10	Sulphite	≤2 ppm	Cr	≤0.01 ppm	Pb	≤0.02 ppm
Identification	Positive	Ag	≤0.02 ppm	Cu	≤0.01 ppm	Pt	≤0.02 ppm
Assay (acidimetric)	49.0 ÷ 51.0 %	Al	≤0.05 ppm	Fe	≤0.1 ppm	Sb	≤0.01 ppm
Density at 20°C	~ 1.17	As	≤0.03 ppm	Ga	≤0.02 ppm	Sn	≤0.02 ppm
Hydrofluosilicic acid	≤20 ppm	Au	≤0.02 ppm	In	≤0.02 ppm	Sr	≤0.02 ppm
Chloride	≤1 ppm	B	≤0.02 ppm	K	≤0.1 ppm	Ta	≤0.1 ppm
Phosphate	≤0.5 ppm	Ba	≤0.1 ppm	Li	≤0.02 ppm	Ti	≤0.1 ppm
Heavy metals (Pb)	≤0.1 ppm	Be	≤0.01 ppm	Mg	≤0.1 ppm	Tl	≤0.02 ppm
Nitrate	≤3 ppm	Bi	≤0.02 ppm	Mn	≤0.01 ppm	V	≤0.01 ppm
Residue on ignition	≤5 ppm	Ca	≤0.1 ppm	Mo	≤0.01 ppm	Zn	≤0.05 ppm
Subst. reducing KMn04	≤4 ppm	Cd	≤0.01 ppm	Na	≤0.2 ppm	Zr	≤0.01 ppm

Code	Size	Packaging	Notes
405737	1 l	Plastic bottle	

Considered as toxic gas

Hydrofluoric acid 50% > RS - MOS - For electronic use

RS

Description	Clear colourless liquid	Sulphate	≤1 ppm	Co	≤0.01 ppm	Ni	≤0.01 ppm
Colour (APHA)	≤10	Sulphite	≤2 ppm	Cr	≤0.01 ppm	Pb	≤0.02 ppm
Identification	Positive	Ag	≤0.02 ppm	Cu	≤0.01 ppm	Pt	≤0.02 ppm
Assay (acidimetric)	49.0 ÷ 51.0 %	Al	≤0.05 ppm	Fe	≤0.1 ppm	Sb	≤0.01 ppm
Density at 20° C	1.152 ÷ 1.158	As	≤0.03 ppm	Ga	≤0.02 ppm	Sn	≤0.02 ppm
Hydrofluosilicic acid	≤20 ppm	Au	≤0.02 ppm	In	≤0.02 ppm	Sr	≤0.02 ppm
Chloride	≤1 ppm	B	≤0.02 ppm	K	≤0.1 ppm	Ta	≤0.1 ppm
Phosphate	≤0.5 ppm	Ba	≤0.1 ppm	Li	≤0.02 ppm	Ti	≤0.1 ppm
Heavy metals (Pb)	≤0.1 ppm	Be	≤0.01 ppm	Mg	≤0.1 ppm	Tl	≤0.02 ppm
Nitrate	≤3 ppm	Bi	≤0.02 ppm	Mn	≤0.01 ppm	V	≤0.01 ppm
Residue on ignition	≤5 ppm	Ca	≤0.1 ppm	Mo	≤0.01 ppm	Zn	≤0.05 ppm
Subst. reducing KMn04	≤4 ppm	Cd	≤0.01 ppm	Na	≤0.2 ppm	Zr	≤0.01 ppm

Code	Size	Packaging	Notes
405653	1 l	Plastic bottle	

Considered as toxic gas

Hydrofluoric acid 50% > RPE - For analysis - ACS - ISO

RPE

Description	Clear liquid	Sulphate	≤2 ppm	Co	≤0.02 ppm	Ni	≤0.05 ppm
Colour (APHA)	≤10	Sulphite	≤2 ppm	Cr	≤0.05 ppm	Pb	≤0.2 ppm
Identification	Positive	Ag	≤0.02 ppm	Cu	≤0.02 ppm	Sr	≤0.02 ppm
Density at 20° C	1.152 ÷ 1.158	Al	≤0.05 ppm	Fe	≤0.2 ppm	Ti	≤0.1 ppm
Hydrofluosilicic acid	≤20 ppm	As	≤0.05 ppm	K	≤0.1 ppm	Tl	≤0.05 ppm
Chloride	≤1 ppm	Ba	≤0.1 ppm	Li	≤0.02 ppm	V	≤0.05 ppm
Phosphate	≤0.5 ppm	Be	≤0.02 ppm	Mg	≤0.2 ppm	Zn	≤0.05 ppm
Heavy metals (Pb)	≤0.5 ppm	Bi	≤0.1 ppm	Mn	≤0.05 ppm	Zr	≤0.1 ppm
Residue on ignition	≤5 ppm	Ca	≤0.5 ppm	Mo	≤0.05 ppm	Assay (acidimetric)	49 ÷ 51 %
Subst. reducing KMn04	≤4 ppm	Cd	≤0.01 ppm	Na	≤0.5 ppm		

Code	Size	Packaging	Notes
405722	1 l	Plastic bottle	

Considered as toxic gas

Product specifications are subject to changes. Please visit our website for updates.

Hydrofluoric acid 47-51%

• Acido fluoridrico 47-51% • Acide fluorhydrique 47-51% • Acido fluorhídrico 47-51% • Fluorwasserstoff 47-51%

HF
Molecular Weight: 20
CAS: 7664-39-3

Classification transport
ONU: 1790
Transport Hazard class: 8
Packing group II



Danger
H300-H310-H330-H314
P280-P284-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Hydrofluoric acid 47-51% > RS - Ultrapure - For trace analysis at ppt level

RS

Description	Clear liquid	Er	≤ 1 ppt	Mo	≤ 10 ppt	Te	≤ 1 ppt
Identification	Positive	Eu	≤ 1 ppt	Nd	≤ 1 ppt	Tb	≤ 1 ppt
Al	≤ 20 ppt	Gd	≤ 1 ppt	Ni	≤ 20 ppt	Tl	≤ 10 ppt
Sb	≤ 20 ppt	Ga	≤ 10 ppt	Nb	≤ 10 ppt	Th	≤ 1 ppt
As	≤ 50 ppt	Ge	≤ 10 ppt	Pd	≤ 20 ppt	Tm	≤ 1 ppt
Ba	≤ 10 ppt	Au	≤ 20 ppt	Pt	≤ 20 ppt	Sn	≤ 20 ppt
Be	≤ 10 ppt	Hf	≤ 10 ppt	K	≤ 10 ppt	Ti	≤ 20 ppt
Bi	≤ 10 ppt	Ho	≤ 1 ppt	Pr	≤ 1 ppt	W	≤ 20 ppt
B	≤ 100 ppt	In	≤ 1 ppt	Re	≤ 10 ppt	U	≤ 1 ppt
Cd	≤ 10 ppt	Fe	≤ 10 ppt	Rh	≤ 20 ppt	V	≤ 10 ppt
Ca	≤ 10 ppt	La	≤ 10 ppt	Rb	≤ 20 ppt	Yb	≤ 1 ppt
Ce	≤ 10 ppt	Pb	≤ 10 ppt	Ru	≤ 20 ppt	Y	≤ 1 ppt
Cs	≤ 10 ppt	Li	≤ 10 ppt	Sm	≤ 1 ppt	Total sulphur	≤ 100 ppb
Cr	≤ 10 ppt	Lu	≤ 1 ppt	Sc	≤ 10 ppt	Zr	≤ 10 ppt
Co	≤ 10 ppt	Mg	≤ 10 ppt	Ag	≤ 10 ppt	Assay	47 ÷ 51 %
Cu	≤ 10 ppt	Mn	≤ 10 ppt	Na	≤ 10 ppt		
Dy	≤ 1 ppt	Hg	≤ 50 ppt	Sr	≤ 10 ppt		

Code	Size	Packaging	Notes
405611	500 ml	Plastic bottle	

Hydrofluoric acid 47-51% > RS - Superpure - For trace analysis at ppb level

RS

Description	Clear liquid	Cr	≤ 1 ppb	Mg	≤ 1 ppb	Ag	≤ 0.5 ppb
Identification	Positive	Co	≤ 0.1 ppb	Mn	≤ 0.1 ppb	Na	≤ 1 ppb
Colour (APHA)	≤ 10	Cu	≤ 0.5 ppb	Hg	≤ 1 ppb	Sr	≤ 0.1 ppb
Chloride	≤ 4 ppm	Dy	≤ 0.1 ppb	Mo	≤ 0.1 ppb	Te	≤ 0.1 ppb
Total phosphorus	≤ 0.05 ppm	Er	≤ 0.1 ppb	Nd	≤ 0.1 ppb	Tb	≤ 0.1 ppb
Total sulphur	≤ 0.1 ppm	Eu	≤ 0.1 ppb	Ni	≤ 0.5 ppb	Tl	≤ 0.1 ppb
Hydrofluosilicic acid	≤ 20 ppm	Gd	≤ 0.1 ppb	Nb	≤ 0.1 ppb	Th	≤ 0.1 ppb
Al	≤ 1 ppb	Ga	≤ 0.1 ppb	Pd	≤ 0.2 ppb	Tm	≤ 0.1 ppb
Sb	≤ 0.2 ppb	Ge	≤ 0.1 ppb	Pt	≤ 0.2 ppb	Sn	≤ 0.5 ppb
As	≤ 0.5 ppb	Au	≤ 0.2 ppb	K	≤ 1 ppb	Ti	≤ 1 ppb
Ba	≤ 0.1 ppb	Hf	≤ 0.1 ppb	Pr	≤ 0.1 ppb	W	≤ 0.5 ppb
Be	≤ 0.1 ppb	Ho	≤ 0.1 ppb	Re	≤ 0.1 ppb	U	≤ 0.1 ppb
Bi	≤ 0.1 ppb	In	≤ 0.1 ppb	Rh	≤ 0.1 ppb	V	≤ 0.1 ppb
B	≤ 1 ppb	Fe	≤ 1 ppb	Rb	≤ 0.1 ppb	Yb	≤ 0.1 ppb
Cd	≤ 0.1 ppb	La	≤ 0.1 ppb	Ru	≤ 0.1 ppb	Y	≤ 0.1 ppb
Ca	≤ 1 ppb	Pb	≤ 0.1 ppb	Sm	≤ 0.1 ppb	Zn	≤ 1 ppb
Ce	≤ 0.1 ppb	Li	≤ 0.1 ppb	Sc	≤ 0.1 ppb	Zr	≤ 0.1 ppb
Cs	≤ 0.1 ppb	Lu	≤ 0.1 ppb	Se	≤ 1 ppb	Assay	47 ÷ 51 %

Code	Size	Packaging	Notes
405716	500 ml	Plastic bottle	

Hydrofluoric acid 39.5%

• Acido fluoridrico 39.5% • Acide fluorhydrique 39.5% • Acido fluorhídrico 39.5% • Fluorwasserstoff 39.5%

HF
Molecular Weight: 20
CAS: 7664-39-3

Classification transport
ONU: 1790
Transport Hazard class: 8
Packing group II



Danger
H300-H310-H330-H314
P280-P284-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Hydrofluoric acid 39.5% > RPE - For analysis - ACS - ISO

RPE

Description	Clear colourless liquid	Sulphate	≤ 2 ppm	Co	≤ 0.02 ppm	Ni	≤ 0.02 ppm
Identification	Positive	Sulphite	≤ 2 ppm	Cr	≤ 0.05 ppm	Pb	≤ 0.05 ppm
Density at 20° C	1.127 ÷ 1.133	Ag	≤ 0.02 ppm	Cu	≤ 0.02 ppm	Sr	≤ 0.02 ppm
Assay (acidimetric)	39.1 ÷ 39.9 %	Al	≤ 0.05 ppm	Fe	≤ 0.2 ppm	Ti	≤ 0.1 ppm
Hydrofluosilicic acid	≤ 20 ppm	As	≤ 0.05 ppm	K	≤ 0.1 ppm	Tl	≤ 0.05 ppm
Chloride	≤ 1 ppm	Ba	≤ 0.1 ppm	Li	≤ 0.02 ppm	V	≤ 0.05 ppm
Phosphate	≤ 0.5 ppm	Be	≤ 0.02 ppm	Mg	≤ 0.2 ppm	Zn	≤ 0.05 ppm
Heavy metals (Pb)	≤ 0.5 ppm	Bi	≤ 0.1 ppm	Mn	≤ 0.05 ppm	Zr	≤ 0.1 ppm
Residue on ignition	≤ 5 ppm	Ca	≤ 0.5 ppm	Mo	≤ 0.05 ppm		
Subst. reducing KMn04	≤ 4 ppm	Cd	≤ 0.01 ppm	Na	≤ 0.5 ppm		

Code	Size	Packaging	Notes
405761	1 l	Plastic bottle	
405765	5 l	Plastic bottle	

Hydrofluoric acid 39.5% > RE - Pure

RE

Description	Clear colourless liquid	Residue on ignition	≤ 0.5 %	Fe	≤ 500 ppm
Identification	Positive	Sulphate	≤ 1 %	Assay (acidimetric)	38.0 ÷ 39.9 %

Code	Size	Packaging	Notes
303731	1 l	Plastic bottle	

Hydrofluoric acid 0.1 mol/l

• Acido fluoridrico 0.1M • Acide fluorhydrique 0.1M • Acido fluorhídrico 0.1M • Fluorwasserstoff 0.1 mol/l

Hydrofluoric acid 0.1 mol/l > RPE - For analysis

RPE

Code	Size	Packaging	Notes
507410	1 l	Bottle	

Hydrofluoric acid diluted

• Acido fluoridrico diluito • Acide fluorhydrique dilué • Acido fluorhídrico diluido • Fluorwasserstoff Verdünnte

HF
CAS: 7664-39-3

Classification transport
ONU: 2922
Transport Hazard class: 8
Packing group II



Danger
H290-H301-H311-H314
P280-P301+P310a-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338

Hydrofluoric acid diluted > RPE - For analysis

RPE

Code	Size	Packaging	Notes
405775	250 ml	Bottle	

Contains ~ 2% of hydrofluoric acid, ~ 5% hydrochloric acid

Product specifications are subject to changes. Please visit our website for updates.

Hydrogen peroxide solution 40% w/v

• Perossido di idrogeno soluzione 40% m/v • Eau oxygénée solution 40% m/v • Hidrógeno peróxido solución 40% p/v • Wasserstoffperoxyd 40% m/v

H_2O_2
Molecular Weight: 34,01
CAS: 7722-84-1

Classification transport
ONU: 2014
Transport Hazard class: 5.1
Packing group II



Danger
H302-H332-H315-H318-H335-H336
P304+P340-P310a-P305+P351+P338-P330-P362+P364-P403+P233

Hydrogen peroxide solution 40% w/v > RE - Pure - Stabilized — RE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 18° C ..... 1.127 ÷ 1.137
Assay (oxidimetric) ..... ≥39 % m/v

Code	Size	Packaging	Notes
307701	1 l	Plastic bottle	
307708	5 l	Plastic bottle	
307709	60 kg	Plastic tank	

154 volumes

Hydrogen peroxide solution 35%

• Perossido di idrogeno soluzione 35% • Eau oxygénée solution 35% • Hidrógeno peróxido solución 35% • Wasserstoffperoxyd 35%

H_2O_2
Molecular Weight: 34,01
CAS: 7722-84-1

Classification transport
ONU: 2014
Transport Hazard class: 5.1
Packing group II



Danger
H302-H332-H315-H318-H335-H336
P304+P340-P310a-P305+P351+P338-P330-P362+P364-P403+P233

Hydrogen peroxide solution 35% > ERBApharm - Prepared from raw material according Ph.Eur — ERBApharm

Description ..... Clear colourless liquid
Identification ..... Positive
Appearance of solution ..... Passes test
Residue on ignition ..... ≤ 0.05 %
Free acid (H2SO4) ..... ≤ 0.025 %
Heavy metals (Pb) ..... ≤ 0.0002 %
Chloride ..... ≤ 0.005 %
Residue on evaporation ..... ≤ 0.10 %
Residual solvents (Current ICH) ..... Conform
Assay ..... 34.5 - 36.0 %

Code	Size	Packaging	Notes
307742	2.5 l	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade
130 volumes

Hydrogen peroxide solution 30-32%

• Perossido di idrogeno soluzione 30-32% • Eau oxygénée solution 30-32% • Hidrógeno peróxido solución 30-32% • Wasserstoffperoxyd 30-32%

H_2O_2
Molecular Weight: 34,01
CAS: 7722-84-1

Classification transport
ONU: 2014
Transport Hazard class: LQ



Danger
H302-H332-H318
P261-P264-P271-P304+P340-P310a-P305+P351+P338

Hydrogen peroxide solution 30-32% > RS - Ultrapure - For trace analysis at ppt level — RS

Description ..... Clear colourless liquid
Identification ..... Positive
Al ..... ≤ 50 ppt
Sb ..... ≤ 10 ppt
As ..... ≤ 100 ppt
Ba ..... ≤ 10 ppt
Be ..... ≤ 10 ppt
Bi ..... ≤ 10 ppt
B ..... ≤ 100 ppt
Cd ..... ≤ 10 ppt
Ca ..... ≤ 100 ppt
Ce ..... ≤ 1 ppt
Cs ..... ≤ 1 ppt
Cr ..... ≤ 10 ppt
Co ..... ≤ 10 ppt
Cu ..... ≤ 10 ppt
Dy ..... ≤ 1 ppt
Er ..... ≤ 1 ppt
Eu ..... ≤ 1 ppt
Gd ..... ≤ 1 ppt
Ga ..... ≤ 10 ppt
Ge ..... ≤ 10 ppt
Au ..... ≤ 10 ppt
Hf ..... ≤ 1 ppt
Ho ..... ≤ 1 ppt
In ..... ≤ 1 ppt
Fe ..... ≤ 20 ppt
La ..... ≤ 1 ppt
Pb ..... ≤ 10 ppt
Li ..... ≤ 10 ppt
Lu ..... ≤ 1 ppt
Mg ..... ≤ 20 ppt
Mn ..... ≤ 10 ppt
Hg ..... ≤ 50 ppt
Mo ..... ≤ 10 ppt
Nd ..... ≤ 1 ppt
Ni ..... ≤ 20 ppt
Nb ..... ≤ 10 ppt
Pd ..... ≤ 10 ppt
K ..... ≤ 20 ppt
Pr ..... ≤ 1 ppt
Re ..... ≤ 10 ppt
Rh ..... ≤ 10 ppt
Rb ..... ≤ 10 ppt
Ru ..... ≤ 10 ppt
Sm ..... ≤ 1 ppt
Sc ..... ≤ 10 ppt
Se ..... ≤ 100 ppt
Ag ..... ≤ 10 ppt
Na ..... ≤ 50 ppt
Sr ..... ≤ 10 ppt
Ta ..... ≤ 10 ppt
Te ..... ≤ 1 ppt
Tb ..... ≤ 1 ppt
Tl ..... ≤ 1 ppt
Th ..... ≤ 1 ppt
Tm ..... ≤ 1 ppt
Sn ..... ≤ 50 ppt
Ti ..... ≤ 20 ppt
W ..... ≤ 20 ppt
U ..... ≤ 1 ppt
V ..... ≤ 10 ppt
Yb ..... ≤ 1 ppt
Y ..... ≤ 1 ppt
Zn ..... ≤ 50 ppt
Zr ..... ≤ 10 ppt
Assay ..... 30 ÷ 32 % (p/p)

Code	Size	Packaging	Notes
412051	500 ml	Plastic bottle	

Hydrogen peroxide solution 30%

• Perossido di idrogeno soluzione 30% • Eau oxygénée solution 30% • Hidrógeno peróxido solución 30% • Wasserstoffperoxyd 30%

H_2O_2
Molecular Weight: 34,01
CAS: 7722-84-1

Classification transport
ONU: 2014
Transport Hazard class: LQ

Danger
H302-H332-H318
P261-P264-P271-P304+P340-P310a-P305+P351+P338

Hydrogen peroxide solution 30% > RS - VLSI - For electronic use

RS

Code	Size	Packaging	Notes
527621	1 l	Bottle	
527620	2.5 l	Bottle	

Stabilized; 110 volumes. For specifications, contact our customer service for a certificate of analysis

Hydrogen peroxide solution 30% > RS - RSE - For electronic use - Stabilized

RS

Description	Clear liquid	Chloride	≤0.5 ppm	Ca	≤0.05 ppm	Mn	≤0.01 ppm
Colour (APHA)	≤10	Phosphate	≤1 ppm	Cd	≤0.01 ppm	Na	≤0.05 ppm
Identification	Positive	Heavy metals (Pb)	≤1 ppm	Co	≤0.01 ppm	Ni	≤0.01 ppm
Density at 20° C	1.120 ÷ 1.124	Nitrate	≤2 ppm	Cr	≤0.01 ppm	Pb	≤0.01 ppm
Assay (oxidimetric)	29 ÷ 31 % m/m	Sulphate	≤1 ppm	Cu	≤0.01 ppm	Zn	≤0.02 ppm
Residue on evaporation	≤5 ppm	Al	≤0.1 ppm	Fe	≤0.03 ppm		
Acidity (H2SO4)	≤20 ppm	As	≤0.01 ppm	K	≤0.02 ppm		
Ammonium	≤1 ppm	Ba	≤0.01 ppm	Mg	≤0.02 ppm		

Code	Size	Packaging	Notes
412161	1 l	Plastic bottle	
412162	5 l	Plastic bottle	
412163	25 kg	Drum	

110 volumes

Hydrogen peroxide solution 30% > RS - MOS - For electronic use - Stabilized

RS

Description	Clear liquid	Sulphate	≤2 ppm	Cr	≤0.01 ppm	Pb	≤0.02 ppm
Colour (APHA)	≤10	Ag	≤0.02 ppm	Cu	≤0.01 ppm	Pt	≤0.05 ppm
Identification	Positive	Al	≤0.1 ppm	Fe	≤0.1 ppm	Sb	≤0.01 ppm
Density at 18° C	1.120 ÷ 1.124	As	≤0.01 ppm	Ga	≤0.02 ppm	Sn	≤0.25 ppm
Assay (oxidimetric)	29.0 ÷ 31.0 %	Au	≤0.05 ppm	In	≤0.02 ppm	Sr	≤0.02 ppm
Residue on evaporation	≤10 ppm	B	≤0.02 ppm	K	≤0.1 ppm	Ta	≤0.1 ppm
Acidity (H2SO4)	≤20 ppm	Ba	≤0.1 ppm	Li	≤0.01 ppm	Ti	≤0.05 ppm
Ammonium	≤1 ppm	Be	≤0.02 ppm	Mg	≤0.05 ppm	Tl	≤0.05 ppm
Chloride	≤0.5 ppm	Bi	≤0.02 ppm	Mn	≤0.01 ppm	V	≤0.05 ppm
Phosphate	≤1 ppm	Ca	≤0.2 ppm	Mo	≤0.05 ppm	Zn	≤0.05 ppm
Heavy metals (Pb)	≤0.2 ppm	Cd	≤0.01 ppm	Na	≤0.2 ppm	Zr	≤0.05 ppm
Nitrate	≤2 ppm	Co	≤0.01 ppm	Ni	≤0.01 ppm		

Code	Size	Packaging	Notes
412081	1 l	Plastic bottle	

110 volumes

Hydrogen peroxide solution 30% > RS - For agroalimentary analysis

RS

Aspect ..... Conform Assay ..... 29.0 ÷ 31.0 %

Code	Size	Packaging	Notes
502044	5 l	Plastic tank	

Stabilized; 110 volumes.

Hydrogen peroxide solution 30% > RS - For microanalysis - Stabilized

RS

Description ..... Clear colourless liquid Identification ..... Positive Density at 18° C ..... 1.120 ÷ 1.124

Code	Size	Packaging	Notes
412102	250 ml	Plastic bottle	

110 volumes

Product specifications are subject to changes. Please visit our website for updates.

Hydrogen peroxide solution 30% > RPE - For analysis - ACS - Reag.Ph.Eur. - Reag.USP - Stabilized

RPE

Description	Clear liquid	Nitrate	≤2 ppm	Ca	≤0.2 ppm	Mn	≤0.01 ppm
Colour (APHA)	≤10	Sulphate	≤2 ppm	Cd	≤0.01 ppm	Na	≤0.5 ppm
Identification	Positive	Phosphate	≤2 ppm	Co	≤0.01 ppm	Ni	≤0.02 ppm
Density at 20° C	1.10 ÷ 1.13	Ammonium	≤5 ppm	Cr	≤0.02 ppm	Pb	≤0.02 ppm
Organic stabilizers	≤500 ppm	Heavy metals (Pb)	≤1 ppm	Cu	≤0.01 ppm	Zn	≤0.1 ppm
Residue on evaporation	≤20 ppm	Al	≤0.5 ppm	Fe	≤0.1 ppm	Assay (oxidimetric)	29.0 ÷ 31.0 % m/m
Acidity	≤0.0006 meq/g	As	≤0.01 ppm	K	≤0.1 ppm		
Chloride	≤0.5 ppm	Ba	≤0.05 ppm	Mg	≤0.1 ppm		

Code	Size	Packaging	Notes
412077	175 ml	Plastic bottle	
412071	250 ml	Plastic bottle	
412072	1 l	Plastic bottle	
412074	25 kg	Plastic drum	
412076	200 kg	Plastic drum	

Hydrogen peroxide solution 30% > ERBApharm - According to pharmacopoeia: Ph.Eur. - Stabilized

ERBApharm

Description	Clear colourless liquid	Acidity	Conform Ph.Eur.	Organic stabilizers	≤ 500 ppm	Origin (BSE/TSE)	Synthesis
Identification	Positive	Non volat.substances	≤ 2 g/l	Assay (oxidimetric)	29.0 ÷ 31.0 %	Residual solvents (Current ICH)	Conform

Code	Size	Packaging	Notes
307685	25 kg	Drum	

110 volumes. In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Hydrogen peroxide solution 6%

• Perossido di idrogeno soluzione 6% • Eau oxygénée en solution à 6% • Hidrógeno peróxido solución 6% • Wasserstoffperoxyd 6%

H_2O_2
Molecular Weight: 34,01
CAS: 7722-84-1

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Hydrogen peroxide solution 6% > RE - Pure

RE

Description	Clear colourless liquid	Identification	Positive	Assay (oxidimetric)	5.75 ÷ 6.25 %

Code	Size	Packaging	Notes
307641	10 l	Bottle	

Stabilized; 20 volumes.

Hydrogen peroxide solution 3.5% w/v

• Perossido di idrogeno soluzione 3.5% m/v • Eau oxygénée solution 3.5% m/v • Hidrógeno peróxido solución 3.5% p/v • Wasserstoffperoxyd 3.5% w/v

H_2O_2
Molecular Weight: 34,01
CAS: 7722-84-1

HEU210

Hydrogen peroxide solution 3.5% w/v > RE - Pure

RE

Description	Clear colourless liquid	Residue on evaporation	≤0.2 %	Heavy metals (Pb)	≤5 ppm
Identification	Positive	Acidity (HCl)	≤360 ppm	As	≤1 ppm
Density at 20° C	~ 1.015	Fixative	≤500 ppm	Assay (oxidimetric)	3.45 ÷ 3.75 % m/v

Code	Size	Packaging	Notes
E307661	1 l	Bottle	

Stabilized; 12 volumes.

Hydrogen peroxide solution 3%

• Perossido di idrogeno soluzione 3% • Eau oxygénée solution 3% • Hidrógeno peróxido solución 3% • Wasserstoffperoxyd 3%

H_2O_2
Molecular Weight: 34,01
CAS: 7722-84-1

HEU210

Hydrogen peroxide solution 3% > ERBApharm - According to pharmacopoeia: Ph.Eur.-FU - Stabilized

ERBApharm

Description ..... Clear colourless liquid
Identification ..... Positive
Assay (oxidimetric) ..... 2.5 ÷ 3.5 %
Acidity ..... Conform Ph.Eur.
Non volat.substances ..... Conform Ph.Eur.
Organic stabilizers ..... Conform Ph.Eur.
Origin (BSE/TSE) ..... Synthesis

Code	Size	Packaging	Notes
307671	1 l	Plastic bottle	
307678	50 kg	Plastic tank	

10 volumes. In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Hydroquinone

• Idrochinone • Hydroquinone • Hidroquinona • Hydrochinon

Synonym:
• *1,4-Benzenediol*
• *1,4-Dihydroxybenzene*

$1,4\text{-}(OH)_2C_6H_4$
Molecular Weight: 110,11
CAS: 123-31-9
EEC-N: 204-617-8



Danger
H302-H318-H317-H341-H351-H400
P261-P280-P305+P351+P338-P310a-P330-P362+P364

Hydroquinone > RPE - For analysis

RPE

Description ..... White crystals
Identification ..... Positive
Melting point ..... 171 ÷ 174 ° C
Loss on drying ..... ≤1 %
Heavy metals (Pb) ..... ≤5 ppm
Residue on ignition ..... ≤200 ppm
Fe ..... ≤5 ppm
Assay (oxidimetric) ..... ≥99 %

Code	Size	Packaging	Notes
455325	250 g	Plastic bottle	

Hydroquinone > RE - Pure

RE

Description ..... White crystals
Identification (I.R.) ..... Conform
Melting point ..... 170 ÷ 175 ° C
Heavy metals (Pb) ..... ≤ 10 ppm
Resorcinol ..... ≤ 0.1 %
Fe ..... ≤ 10 ppm
Assay (oxidimetric) ..... ≥ 99 %

Code	Size	Packaging	Notes
348126	500 g	Plastic bottle	
348128	1 kg	Plastic bottle	
348129	5 kg	Plastic tank	
348124	25 kg	Plastic bucket	

p-Hydroxybenzaldehyde

• p-IdrossibenzAldeide • p-Hydroxybenzaldéhyde • p-Hidroxibenzaldehído • p-Hydroxybenzaldehyd

HOC_6H_4CHO
Molecular Weight: 122,12
CAS: 123-08-0
EEC-N: 204-599-1

p-Hydroxybenzaldehyde > RE - Pure

RE

Description ..... White crystalline powder
Identification ..... Positive
Melting point ..... 114 ÷ 117 ° C
Assay (acidimetric) ..... ≥98 %

Code	Size	Packaging	Notes
467254	100 g	Glass bottle	

Hydroxybenzene ▶ Phenol

2-Hydroxybenzoic acid ▶ Salicylic acid

Product specifications are subject to changes. Please visit our website for updates.

Hydroxylamine solution, alcoholic

• Idrossilammina soluzione, alcolico • Hydroxylamine solution alcoolique • Hidroxilamina solución, alcohólico • Hydroxylaminlösung, alkoholisch

Classification transport
ONU: 2733
Transport Hazard class: 3
Packing group II



Danger
H225-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Hydroxylamine solution, alcoholic > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611044301	100 ml	Glass bottle	Ref Ph.Eur 1044301

Hydroxylamine sulphate

• Idrossilammina solfato • Hydroxylamine sulfate • Hidroxilamonio sulfato • Hydroxylaminsulfat

Synonym:
Hydroxylammonium sulfate

$(NH_2OH)_2.H_2SO_4$
Molecular Weight: 164,14
CAS: 10039-54-0
EEC-N: 233-118-8

Classification transport
ONU: 2865
Transport Hazard class: 8
Packing group III



Warning
H290-H302-H312-H315-H319-H317-H351-H373-H400
P280-P305+P351+P338-P308+P313-P330-P362+P364-P337+P313

Hydroxylamine sulphate > RPE - For analysis

RPE

Description ..... White crystals
Identification ..... Positive
Ammonium ..... ≤0.1 %
Chloride ..... ≤10 ppm
Heavy metals (Pb) ..... ≤10 ppm
Fe ..... ≤5 ppm
Assay (oxidimetric) ..... ≥99 %

Code	Size	Packaging	Notes
455525	250 g	Plastic bottle	
455527	1 kg	Plastic bottle	
455523	25 kg	Plastic bucket	

4-Hydroxy-4-methyl-2-pentanone ▶ Diacetone alcohol

1-Hydroxynaphthalene ▶ 1-Naphthol

2-Hydroxynaphthalene ▶ 2-Naphthol

8-Hydroxyquinoline

• 8-Ossichinolina • 8-Hydroxyquinoléine • 8-Oxiquinoleína • 8-Hydroxychinolin

$HOC_6H_3N:CHCH:CH$
Molecular Weight: 145,16
CAS: 148-24-3
EEC-N: 205-711-1

Classification transport
ONU: 2811
Transport Hazard class: 6.1
Packing group III



Danger
H301-H318-H317-H360D-H410-HA26
P261-P280-P301+P310a-P330-P305+P351+P338-P308+P313

8-Hydroxyquinoline > RPE - For analysis

RPE

Description ..... Off white to light brown crystal
Identification ..... Positive
Melting point ..... 72.5 ÷ 74.5 °C
Residue on ignition ..... ≤2000 ppm
Sulphate ..... ≤200 ppm
Assay ..... ≥ 99.5 %

Code	Size	Packaging	Notes
467353	50 g	Glass bottle	
467355	250 g	Plastic bottle	
467356	1 kg	Plastic bottle	

3-Hydroxytoluene ▶ m-Cresol

Hypophosphorous acid 50%

• Acido ipofosforoso 50% • Acide hypophosphoreux 50% • Acido hipofosforoso 50% • Hypophosphorige Säure 50%

H_3PO_2
Molecular Weight: 66,04
CAS: 6303-21-5
EEC-N: 228-601-5

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Hypophosphorous acid 50% > RPE - For analysis

RPE

Description ..... Yellow colourless liquid	Fe ..... ≤50 ppm	Sulphate (SO4--) ..... ≤ 500 mg/Kg
Arsenic (As) ..... ≤ 0.05 mg/Kg	Lead (Pb) ..... ≤ 0.5 mg/Kg	Assay (acidimetric) ..... ≥49.5 %
Ca ..... ≤ 30 ppm	Chloride (Cl-) ..... ≤ 200 mg/Kg	Density at 20° C ..... ~ 1.22

Code	Size	Packaging	Notes
406961	100 ml	Glass bottle	
406962	1 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Idrimeter Erba Solution A

• Idrimeter Erba Soluzione A • Idrimeter Erba Solution A • Idrimeter Erba Solución A • Idrimeter Erba Lösung A

Idrimeter Erba Solution A > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
E455256	500 ml	Bottle	
E455257	1 l	Bottle	

For water hardness. Titrant 1 ml = 1 mg CaCO3

Idrimeter Erba Solution B

• Idrimeter Erba Soluzione B • Idrimeter Erba Solution B • Idrimeter Erba Solución B • Idrimeter Erba Lösung B



Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-
P362+P364-P337+P313

Idrimeter Erba Solution B > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive pH at 20° C ..... 9.8 ÷ 10.2

Code	Size	Packaging	Notes
E455266	500 ml	Bottle	
E455267	1 l	Bottle	

For water hardness. Buffer pH 10

Idrimeter Erba Indicator C

• Idrimeter Erba Indicatore C • Idrimeter Erba Indicateur C • Idrimeter Erba Indicador C • Idrimeter Erba Anzeige C

Idrimeter Erba Indicator C > RPE - For analysis

RPE

Description ..... Violet granular powder Identification ..... Positive

Code	Size	Packaging	Notes
E455271	10 g	Bottle	
E455274	100 g	Bottle	

For water hardness

Imidazole

• Imidazolo • Imidazole • Imidazol • Imidazol

Synonym:
• *1,3-Diaza-2,4-cyclopentadiene*
• *Glyoxaline*

$C_3H_4N_2$
Molecular Weight: 68,08
CAS: 288-32-4
EEC-N: 206-019-2

Classification transport
ONU: 3263
Transport Hazard class: 8
Packing group III



Danger
H302-H314-H360D-HA26
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Imidazole > RS - For microanalysis

RS

Description ..... White crystals Identification ..... Positive Assay (GC) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
445551	2 g	Glass bottle	
445552	25 g	Glass bottle	

Immersion oil
• Olio d'immersione • Huile pour immersion • Aceite de inmersión • Immersionsöl



Warning
H302-H411
P264-P270-P301+P312a-P330-P391-P501a

Immersion oil > RS - For microscopy

RS

Description ..... Clear liquid Identification ..... Positive Density at 20° C ..... ≥1.0

Code	Size	Packaging	Notes
466782	100 ml	Glass bottle	
466783	1 l	Glass bottle	

contains benzylbenzoate

Indicator for ammoniacal nitrogen solution
• Indicatore per azoto ammoniacale soluzione • Indicateur pour l'azote ammoniacal • Indicador para nitrógeno amoniacal solución
• Indikator für ammoniakalischen Stickstoff

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group III



Warning
H226-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

Indicator for ammoniacal nitrogen solution > RPE - For analysis

RPE

Description ..... Dark green liquid Identification ..... Positive pH range ..... 4.4 - 6.0

Code	Size	Packaging	Notes
E455651	250 ml	Glass bottle	

Indicator for iodometry
• Indicatore per iodometria • Indicateur pour iodométrie • Indicador para iodometría • Indikator für die Jodometrie

EEC-N: 232-679-6

Indicator for iodometry > RPE - For analysis

RPE

Description ..... White powder Identification ..... Positive

Code	Size	Packaging	Notes
455622	25 g	Glass bottle	
455621	250 g	Plastic bottle	

Indicator papers
• Cartine indicatrici di pH • Papier indicateur de pH • Papel indicador de pH • pH-Indikatorpapier

Indicator papers > RS - For pHmetry

RS

Code	Size	Packaging	Notes
435131000	1 roll	Dispenser	Trichromic scale, pH range 1.0 - 11.0, Sensitivity 1,0 pH
435140000	1 roll	Dispenser	Monochromatic scale, pH range 1.0 - 11.0, Sensitivity 1,0 pH
435150000	1 roll	Dispenser	Monochromatic scale, pH range 1.0 - 14.0, Sensitivity 1,0/2,0 pH
435161000	1 roll	Dispenser	Monochromatic scale, pH range 3.8 - 5.8, Sensitivity 0.2/0.3 pH
435421000	1 roll	Dispenser	Monochromatic scale. pH range 0.5 - 5.5. Sensitivity 0.5 pH
435431000	1 roll	Dispenser	Monochromatic scale. pH range 4.0 - 7.0. Sensitivity 0.3 pH
435441000	1 roll	Dispenser	Monochromatic scale. pH range 6.4 - 8.0. Sensitivity 0.2 pH
435451000	1 roll	Dispenser	Monochromatic scale. pH range 7.2 - 9.7. Sensitivity 0.3 pH
435511000	1 roll	Dispenser	Monochromatic scale. pH range 5.5 - 9.0. Sensitivity 0.5 pH

Product specifications are subject to changes. Please visit our website for updates.

Indicator papers > RS - High sensitivity-Integrated chromatic scale

RS

Code	Size	Packaging	Notes
435493000	100 stripes	Tube	pH range 1.0 - 2.8. Sensitivity 0.2/0.3
435494000	100 stripes	Tube	pH range 1.8 - 3.8. Sensitivity 0.2/0.3
435496000	100 stripes	Tube	pH range 3.8 - 5.5. Sensitivity 0.2/0.3
435498000	100 stripes	Tube	pH range 6.0 - 8.1. Sensitivity 0.2/0.3
435502000	100 stripes	Tube	pH range 8.0 - 9.7. Sensitivity 0.2/0.3

Strip 11 x 100 mm

Indicator papers > RS - Indelibles - with colour scale

RS

Code	Size	Packaging	Notes
435121000	100 stripes	Tube	pH range 0.0 - 14.0, Sensitivity 1.0
435642000	100 stripes	Tube	pH range 0.0 - 6.0. Sensitivity 0.5
435643000	100 stripes	Tube	pH range 2.0 - 9.0. Sensitivity 0.5
435644000	100 stripes	Tube	pH range 4.5 - 10.0. Sensitivity 0.5
435645000	100 stripes	Tube	pH range 7.0 - 14.0. Sensitivity 0.3/0.4

Strip 6 x 85 mm

Indicator universal pH 0-5 hydroalcoholic solution

• Indicatore universale pH 0-5 soluzione idroalcolica • Indicateur universel pH 0-5 solution hydroalcoolique
• Indicador universal pH 0-5 solución hidroalcohólica • Universalanzeige pH-Wert 0-5 hydroalkoholische Lösung

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II

Danger
H225-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

Indicator universal pH 0-5 hydroalcoholic solution > RPE - For analysis

RPE

Description .....Dark green liquid Identification ..... Positive Functionality ..... Conform

Code	Size	Packaging	Notes
E455661	25 ml	Glass bottle	
E455662	500 ml	Glass bottle	

With chromatic scale

Indicator universal pH 1-11 hydroalcoholic solution

• Indicatore universale pH 1-11 soluzione idroalcolica • Indicateur universel pH 1-11 solution hydroalcoolique
• Indicador universal pH 1-11 solución hidroalcohólica • Hydroalkoholische Universalanzeige für pH 1-11

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Indicator universal pH 1-11 hydroalcoholic solution > RPE - For analysis

RPE

Description .....Dark green liquid Identification ..... Positive Functionality ..... Conform

Code	Size	Packaging	Notes
E455702	25 ml	Glass bottle	
E455706	500 ml	Plastic bottle	

With chromatic scale

Indicator universal pH 1-11 water solution

• Indicatore universale pH 1-11 soluzione in acqua • Indicateur universel pH 1-11 solution aqueuse • Indicador universal pH 1-11 solución en agua • Universalindikator pH 1-11 wässrige Lösung

Indicator universal pH 1-11 water solution > RPE - For analysis

RPE

Description .....Dark green liquid Identification ..... Positive Functionality ..... Conform

Code	Size	Packaging	Notes
E455711	25 ml	Glass bottle	
E455712	500 ml	Glass bottle	

With chromatic scale

Indigo carmine dried

• Carminio indaco secco • Carmin d'indigo sec • Indigo Carmin seco • Indigo Carmine getrocknet

Synonym:
Acid Blue 74

$C_{16}H_8N_2Na_2O_8S_2$
Molecular Weight: 466,36
CAS: 860-22-0
EEC-N: 212-728-8

Indigo carmine dried > RPE - For analysis - C.I. 73015

RPE

Description ..... Polvere blu viola Identification ..... Positive Colour change.....(blue - yellow) Assay ..... ≥ 85 %

Code	Size	Packaging	Notes
434932	25 g	Glass bottle	

Dye for microscopy (bacteriology histology). Indicator acid - base (pH 11.6 ÷ 14)

Indigo carmine solution

• Carminio indaco soluzione • Carmin d'indigo solution • Indigo carmín solución • Indigo-Karmin-Lösung

Synonym:
Acid Blue 74

$C_{16}H_8N_2Na_2O_8S_2$
Molecular Weight: 466,36
CAS: 860-22-0

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Indigo carmine solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611045601	1 l	Plastic bottle	Ref Ph.Eur 1045601

Indium standard solution

• Indio standard soluzione • Indium solution standard • Indio, solución patrón • Indium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Indium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505662	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505665	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505663	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Product specifications are subject to changes. Please visit our website for updates.

Indium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503651	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503653	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503655	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503657	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Indium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507743	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid
507508	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Indole

• Indolo • Indole • Indol • indol

Synonym:
1H-Benzo[b]pyrrole

$C_6H_4NHCH:CH$
Molecular Weight: 117,15
CAS: 120-72-9
EEC-N: 204-420-7

Classification transport
ONU: 2811
Transport Hazard class: 6.1
Packing group III

Danger
H302-H311
P264-P270-P280h-P301+P312a-P330-P501a

Indole > RPE - For analysis

RPE

Description .....White powder or flakes
Identification ..... Positive
Melting point..... 51 ÷ 53 ° C
Assay (spectrophotom.) ..... ≥99 %

Code	Size	Packaging	Notes
455801	10 g	Glass bottle	

Inositol

• Inositolo • Inositol • Inositol • Inosit

Synonym:
1,2,3,4,5,6-Hexahydroxycyclohexane

$CHOH(CHOH)_4CHOH$
Molecular Weight: 180,16
CAS: 87-89-8
EEC-N: 201-781-2

Inositol > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Melting point..... 224.0 ÷ 227.0 ° C
Loss on drying ..... ≤0.5 %
Chloride..... ≤50 ppm
Heavy metals (Pb)..... ≤25 ppm
Residue on ignition..... ≤0.1 %
Sulphate..... ≤60 ppm
Calcium..... Negative
Fe ..... ≤5 ppm
Assay (gravimetric) ..... ≥98.0 %
As ..... ≤ 0.0003 %
Cd..... ≤ 0.00082 %
Hg..... ≤ 0.00006 %
Pb..... ≤ 0.001 %

Code	Size	Packaging	Notes
455853	50 g	Glass bottle	

Inositol > RE - Pure

RE

Description ..... White powder
Identification ..... Positive
Loss on drying ..... ≤0.5 %
Sulphated ash..... ≤0.1 %
Chloride..... ≤100 ppm
Heavy metals (Pb)..... ≤40 ppm
Sulphate..... ≤150 ppm
Fe ..... ≤10 ppm

Code	Size	Packaging	Notes
348354	100 g	Plastic bottle	

Inulin

• Inulina • Inuline • Inulina • Inulin

$(C_6H_{10}O_5)n$
CAS: 9005-80-5
EEC-N: 232-684-3

Inulin > RPE - For analysis

RPE

Description	White crystalline powder	(c= 2 in water)		Chloride	≤50 ppm	Fe	≤10 ppm
Identification	Positive	Loss on drying	≤10 %	Sulphate	≤50 ppm		
Specific optical rotation on dry	-32 ÷ -40 °	Residue on ignition	≤0.1 %	Heavy metals (Pb)	≤10 ppm		

Code	Size	Packaging	Notes
455901	10 g	Glass bottle	
455902	25 g	Glass bottle	
455903	100 g	Glass bottle	

Iodide standard solution

• Ioduri standard soluzione • Iodure solution standard • Yoduro, solución patrón • Jodid-Standardlösung

Iodide standard solution > RS - Standard solution for ion chromatography

RS

Code	Size	Packaging	Notes
503261	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503263	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Iodine

• Iodio • Iode • Yodo • Jod

I_2
Molecular Weight: 253,8
CAS: 7553-56-2
EEC-N: 231-442-4

Classification transport
ONU: 3495
Transport Hazard class: 8
Packing group III



Danger
H302-H312-H332-H315-H319-H335-H372-H400
P271-P304+P340-P305+P351+P338-P332+P313-P337+P313-P403+P233

Iodine > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP-DAB

ERBApharm

Description	Grey-violet crystals	Chloride + bromide (Cl)	≤ 0.025 %	Assay (iodometric)	99.8 ÷ 100.5 %
Identification	Positive	Residue on ignition	≤ 0.05 %		

Code	Size	Packaging	Notes
348454	100 g	Glass bottle	
348455	250 g	Glass bottle	
348457	1 kg	Glass bottle	
348451	5 kg	Metallic can	
348452	20 kg	Metal drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Product specifications are subject to changes. Please visit our website for updates.

Iodine resublimed

• Iodio bisublimato • Iode bisublimé • Yodo bisublimado • Jod wieder sublimiert

I_2
Molecular Weight: 253,8
CAS: 7553-56-2
EEC-N: 231-442-4

Classification transport
ONU: 3495
Transport Hazard class: 8
Packing group III

Danger
H302-H312-H332-H315-H319-H335-H372-H400
P271-P304+P340-P305+P351+P338-P332+P313-P337+P313-P403+P233

Iodine resublimed > RPE - For analysis

RPE

Description ..... Grey-violet crystals
Identification ..... Positive
Residue on evaporation ..... ≤100 ppm
Chlorine-Bromine ..... ≤100 ppm
Assay (oxidimetric) ..... ≥99.8 %

Code	Size	Packaging	Notes
455959	100 g	Glass bottle	
455955	250 g	Glass bottle	
455957	1 kg	Glass bottle	
455954	25 kg	Metal drum	

Iodine 0.5 mol/l (1N)

• Iodio 0.5 mol/l (1N) • Iode 0.5 mol/l (1N) • Yodo 0.5 mol/l (1N) • Jod 0.5 mol/l (1N)

I_2
Molecular Weight: 253,8
CAS: 7553-56-2

Warning
H373
P260-P314-P501a

Iodine 0.5 mol/l (1N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613009400	1 l	Glass bottle	Ref Ph.Eur 3009400

Storage: protected from light

Iodine 0.5 mol/l (1N) > RPE - For analysis

RPE

Description ..... Brown red liquid
Assay (potentiometry) ..... 0.99 - 1.01 N
NIST 136 ..... e

Code	Size	Packaging	Notes
456135000	500 ml	Glass bottle	Certified with NIST traceability
456137000	1 l	Glass bottle	Certified with NIST traceability

126.9 g of I2. Volumetric solution ready-to-use

Iodine 0.05 mol/l (0.1N)

• Iodio 0.05 mol/l (0.1N) • Iode 0.05 mol/l (0.1N) • Yodo 0.05 mol/l (0.1N) • Jod 0.05 mol/l (0.1N)

I_2
Molecular Weight: 253,8
CAS: 7553-56-2

Warning
H373
P260-P314-P501a

Iodine 0.05 mol/l (0.1N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613002700	1 l	Glass bottle	Ref Ph.Eur 3002700

Storage: protected from light

Iodine 0.05 mol/l (0.1N) > RPE - For analysis

RPE

Description ..... Brown red liquid
Assay (potentiometry) ..... 0.0998 - 0.1002 N
NIST 136 ..... e

Code	Size	Packaging	Notes
456036000	500 ml	Glass bottle	Certified with NIST traceability
456037000	1 l	Glass bottle	Certified with NIST traceability

12.69 g of I2. Volumetric solution ready-to-use

Iodine 0.05 mol/l (0.1N) > RPE - NORMEX - For analysis

RPE

Description ..... Brown red liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
456051		Glass ampoule	Volume: 60 ml

12,69 g of I2. Volumetric concentrated solution to prepare 1 L of solution 0,1 N

Iodine 0.01 mol/l (0.02N)

• Iodio 0.01 mol/l (0.02N) • Iode 0.01 mol/l (0.02N) • Yodo 0.01 mol/l (0.02N) • Jod 0.01 mol / l (0.02 N)

I_2
Molecular Weight: 253,8
CAS: 7553-56-2

Warning
H373
P260-P314-P501a

Iodine 0.01 mol/l (0.02N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613002900	1 l	Glass bottle	Ref Ph.Eur 3002900

Storage: protected from light

Iodine 0.005 mol/l (0.01N)

• Iodio 0.005 mol/l (0.01N) • Iode 0.005 mol/l (0.01N) • Yodo 0.005 mol/l (0.01N) • Jod 0.005 mol/l (0.01N)

I_2
Molecular Weight: 253,8
CAS: 7553-56-2

Danger
H315-H319-H372
P260-P264-P305+P351+P338-P332+P313-P362+P364-P337+P313

Iodine 0.005 mol/l (0.01N) > RPE - NORMEX - For analysis

RPE

Description ..... Brown clear liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
456121		Glass ampoule	Volume: 60 ml

1,269 g of I2. Volumetric concentrated solution to prepare 1 L of solution 0,01 N

Iodine 10 ppm

• Iodio 10 ppm • Iode 10 ppm • Yodo 10 ppm • Jod 10 ppm

Iodine 10 ppm > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615003809	100 ml	Glass bottle	Concentrated solution: to dilute according to Ref Ph.Eur 5003800

Iodine bromide solution

• Iodo bromuro soluzione • Iode bromure solution • Yodo bromuro solución • Jodbromidlösung

Classification transport
ONU: 2734
Transport Hazard class: 8
Packing group II

Danger
H226-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Iodine bromide solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611045901	1 l	Glass bottle	Ref Ph.Eur 1045901

Storage: protected from light

Product specifications are subject to changes. Please visit our website for updates.

Iodoform
• Iodoformio • Iodoforme • Yodoformo • Iodoform

Synonym:
Triiodomethane

CHI_3
Molecular Weight: 393,73
CAS: 75-47-8
EEC-N: 200-874-5

Warning
H302-H312-H332
P261-P264-P271-P280h-P301+P312a-P304+P340

Iodoform > RE - Pure

RE

Description .....Yellow crystalline powder
Identification .....Positive
Loss on drying .....≤ 1.0 %
Sulphated ash .....≤ 0.2 %
Assay (argentimetric) .....≥ 99.0 %
Appearance of solution .....Conform
Acidity or alkalinity.....Conform
Chloride.....≤ 50 ppm

Code	Size	Packaging	Notes
348554	100 g	Glass bottle	
348558	25 kg	Drum	

Iodomethane ▶ Methyl iodide

Iodoplatinate reagent
• Reattivo iodoplatinato • Réactif à l'iodoplatinate • Yodoplatinato reactivo • Iodoplatinat-Reagenz

Iodoplatinate reagent > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611046309	200 ml	Glass bottle	Ref Ph.Eur 1046300
611046300	1 l	Glass bottle	Ref Ph.Eur 1046300

Storage: protected from light

IPA ▶ Propan-2-ol

Iridium standard solution
• Iridio standard soluzione • Iridium solution standard • Iridio, solución patrón • Iridium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Iridium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505675	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Iron, powder
• Ferro, polvere • Fer, poudre • Hierro, polvo • Eisenpulver

Fe
Molecular Weight: 55,85
CAS: 7439-89-6
EEC-N: 231-096-4

Classification transport
ONU: 3089
Transport Hazard class: 4.1
Packing group III

Warning
H228
P210-P240-P241-P280-P370+P378a

Iron, powder > RPE - For analysis

RPE

Description .....Grey powder
Identification .....Positive
Assay .....≥97 %

Code	Size	Packaging	Notes
451377	1 kg	Plastic bottle	
451373	25 kg	Plastic bucket	

Iron, reduced by hydrogen

• Ferro, ridotto dall'idrogeno • Fer, réduit par l'hydrogène • Hierro, reducido para el hidrógeno • Eisen, reduziert durch Wasserstoff

Fe
Molecular Weight: 55,85
CAS: 7439-89-6
EEC-N: 231-096-4

Classification transport
ONU: 3089
Transport Hazard class: 4.1
Packing group III

Warning
H228
P210-P240-P241-P280-P370+P378a

Iron, reduced by hydrogen > RPE - For analysis

RPE

Description ..... Grey powder
Identification ..... Positive
Ferric ion ..... Conform
Nitrogen compounds (N) ..... ≤50 ppn
H2SO4-insoluble matter ..... ≤0.5 %
Water solubility ..... ≤0.1 %
Assay ..... ≥95 %

Code	Size	Packaging	Notes
451395	250 g	Glass bottle	
451397	1 kg	Plastic bottle	

Iron standard solution

• Ferro standard soluzione • Fer solution standard • Hierro, solución patrón • Eisenstandardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Warning
H290
P234-P390-P406

Iron standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615001601	100 ml	Plastic bottle	A 10 ppm solution: to dilute according to Ref Ph.Eur 5001601
615001602	100 ml	Plastic bottle	A 8 ppm solution: to dilute according to Ref Ph.Eur 5001602
615001603	100 ml	Plastic bottle	A 2 ppm solution: to dilute according to Ref Ph.Eur 5001603
615001605	100 ml	Plastic bottle	A 0.1 % solution Ref Ph.Eur 5001605
615001606	100 ml	Plastic bottle	A 250 ppm solution: to dilute according to Ref Ph.Eur 5001606
615001609	100 ml	Plastic bottle	A 20 ppm solution: to dilute according to Ref Ph.Eur 5001600

Iron standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505612	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505615	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505613	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Iron standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503581	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503583	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503585	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503587	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Iron standard solution > RS - Standard solution for AAS

RS

Description ..... Yellow clear liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
504194	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507393	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
E497515	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
E497511	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Iron standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Yellow clear liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
451311		Plastic ampoule	conc. 1.000 ppm Matrix: Nitric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Iron (II) ammonium sulfate hexahydrate

• Ferro ammonio solfato oso esaidrato • Fer(II) ammonium sulfate hexahydraté
• Hierro (II) amonio sulfato hexahidratado • Eisen (II) ammoniumsulfathexahydrat

Synonym:
• *Ammonium iron(II) sulfate hexahydrate*
• *Mohr's salt*

$Fe(NH_4)_2(SO_4)_2.6H_2O$
Molecular Weight: 392,14
CAS: 7783-85-9
EEC-N: 233-151-8

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Iron (II) ammonium sulfate hexahydrate > RPE - For analysis - ACS

RPE

Description ..... Green - azure crystals
Identification ..... Positive
Phosphate ..... ≤ 30 ppm
Ca ..... ≤ 50 ppm
Ferric salt ..... ≤ 100 ppm
Cu ..... ≤ 30 ppm
Mn ..... ≤ 100 ppm
Zn ..... ≤ 30 ppm
Assay (oxidimetric) ..... 98.5 ÷ 101.5 %
Water-insoluble matter ..... ≤ 100 ppm
Mg ..... ≤ 20 ppm
K ..... ≤ 20 ppm
Na ..... ≤ 0.02 %

Code	Size	Packaging	Notes
451453	100 g	Plastic bottle	
451451	500 g	Plastic bottle	
451457	1 kg	Plastic bottle	
451452	25 kg	Plastic bucket	

Iron (II) ammonium sulfate hexahydrate > RE - Pure

RE

Description ..... Green - azure crystals Identification ..... Positive Ferric ion ..... ≤ 0.01 % Assay (oxidimetric) ..... ≥ 98 %

Code	Size	Packaging	Notes
344007	1 kg	Plastic bottle	
344008	5 kg	Plastic tank	
344003	25 kg	Plastic bucket	

Iron (II) ammonium sulfate 0.12N

• Ammonio ferrosolfato (II) 0.12N • Fer (II) ammonium sulfate 0.12N • Amonio Hierro (II) sulfato 0.12N • Eisen (II) ammoniumsulfat 0.12 N

$(NH_4)_2Fe(SO_4)_2.6H_2O$
Molecular Weight: 392,14
CAS: 7783-85-9

HEU210

Iron (II) ammonium sulfate 0.12N > RS - For environmental analysis (COD determination)

RS

Code	Size	Packaging	Notes
526761	1 l	Bottle	

Iron (II) ammonium sulfate 0.1N

• Ferro ammonio solfato oso 0.1N • Fer(II) ammonium sulfate 0.1N • Amonio (II) hierro II sulfato 0.1 N
• Eisen (II) ammoniumsulfat 0.1 N

Synonym:
• *Ammonium iron(II) sulfate hexahydrate*
• *Mohr's salt*

$Fe(NH_4)_2(SO_4)_2.6H_2O$
Molecular Weight: 392,14
CAS: 7783-85-9

Classification transport
ONU: 2796
Transport Hazard class: 8
Packing group II

Warning
H290-H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Iron (II) ammonium sulfate 0.1N > RPE - For analysis

RPE

Assay (potentiometry) ..... 0.099 - 0.101 N

Code	Size	Packaging	Notes
P3250016	1 l	Glass bottle	

Iron (II) chloride tetrahydrate

• Ferro cloruro oso tetraidrato • Fer (II) chlorure tétrahydraté • Hierro (II) cloruro tetrahidratado
• Eisen(II)-chlorid-Tetrahydrat

Synonym:
Ferrous chloride tetrahydrate

$FeCl_2.4H_2O$
Molecular Weight: 198,81
CAS: 13478-10-9
EEC-N: 231-843-4

Classification transport
ONU: 3260
Transport Hazard class: 8
Packing group III

Danger
H302-H315-H318
P264-P305+P351+P338-P310a-P330-P362+P364-P332+P313

Iron (II) chloride tetrahydrate > RPE - For analysis

RPE

Description ..... Yellow-green crystals
Identification ..... Positive
Total nitrogen ..... ≤20 ppm
Phosphate ..... ≤10 ppm
Subst. not ppt NH4OH ..... ≤500 ppm
Ferric salt ..... ≤0.2 %
Sulphate ..... ≤50 ppm
As ..... ≤1 ppm
Cr ..... ≤20 ppm
Cu ..... ≤20 ppm
Mn ..... ≤0.1 %
Ni ..... ≤50 ppm
Pb ..... ≤20 ppm
Zn ..... ≤20 ppm
Assay (oxidimetric) ..... ≥ 99 %

Code	Size	Packaging	Notes
451574	100 g	Plastic bottle	
451575	500 g	Plastic bottle	
451576	1 kg	Plastic bottle	
451573	25 kg	Drum	

Iron (II) sulfate heptahydrate

• Ferro solfato oso eptaidrato • Fer (II) sulfate heptahydraté • Hierro (II) sulfato heptahidratado
• Eisen(II) sulfat heptahydrat

Synonym:
Ferrous sulfate heptahydrate

$FeSO_4.7H_2O$
Molecular Weight: 278,05
CAS: 7782-63-0
EEC-N: 231-753-5

Warning
H302-H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Iron (II) sulfate heptahydrate > RPE - For analysis - ACS

RPE

Description ..... Green - azure crystals
Identification ..... Positive
Ferric salt ..... ≤ 0.1 %
Water-insoluble matter ..... ≤ 100 ppm
Ca ..... ≤ 50 ppm
Chloride ..... ≤ 10 ppm
Phosphate ..... ≤ 10 ppm
Cu ..... ≤ 50 ppm
Mn ..... ≤ 0.05 %
Zn ..... ≤ 50 ppm
Assay (oxidimetric) ..... ≥ 99.0 %
Mg ..... ≤ 20 ppm
K ..... ≤ 20 ppm
Na ..... ≤ 0.02 %

Code	Size	Packaging	Notes
451878	100 g	Plastic bottle	
451877	1 kg	Plastic bottle	
451879	5 kg	Plastic jar	

Product specifications are subject to changes. Please visit our website for updates.

Iron (II) sulfate heptahydrate > ERBApharm - According to pharmacopoeia: BP-DAB-FU-Ph.Eur.-Ph.Franc.

ERBApharm

Description ..... Green azure crystals
Identification ..... Positive
pH solution 5% ..... 3.0 ÷ 4.0
Chloride ..... ≤ 200 ppm
Ferric ion ..... ≤ 0.3 %
Cr ..... ≤ 50 ppm
Cu ..... ≤ 50 ppm
Mn ..... ≤ 0.1 %
Ni ..... ≤ 50 ppm
Zn ..... ≤ 50 ppm
Assay (oxidimetric) ..... 98.0 ÷ 105.0 %

Code	Size	Packaging	Notes
344957	1 kg	Plastic bottle	
344959	5 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Iron (II) sulfate 0.1mol/l

• Ferro solfato oso 0.1mol/l • Fer (II) sulfate 0.1 mol/l • Hierro (II) sulfato 0.1mol/l • Eisen (II) sulfat 0.1 mol / l

Synonym:
Ferrous sulfate 0.1 mol/l

HEU210

Iron (II) sulfate 0.1mol/l > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613001400	1 l	Plastic bottle	Ref Ph.Eur 3001400

Iron (III) ammonium citrate green

• Ferro ammonio citrato verde • Fer (III) ammonium citrate vert • Hierro (III) y amonio citrato verde
• Eisen (III) -ammoniumcitratgrün

Synonym:
• *Ammonium iron(III) citrate*
• *Ammonium ferric citrate*

$FeNH_4(C_6H_5O_7)$
Molecular Weight: 482,19
CAS: 1185-57-5
EEC-N: 214-686-6

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Iron (III) ammonium citrate green > RE - Pure

RE

Description ..... Green-yellowish crystalline powder
Identification ..... Positive
Sulphate ..... ≤ 0.5 %
As ..... ≤ 4 ppm
Assay (oxidimetric) ..... 14.0 ÷ 16.0 % Fe

Code	Size	Packaging	Notes
343605	250 g	Plastic bottle	
343607	1 kg	Plastic bottle	
343606	5 kg	Plastic tank	

Iron (III) ammonium citrate red

• Ferro ammonio citrato rosso • Fer (III) ammonium citrate rouge • Hierro (III) y amonio citrato rojo
• Eisen (III) ammoniumcitrat rot

Synonym:
• *Ammonium iron(III) citrate*
• *Ammonium ferric citrate*

$FeNH_4(C_6H_5O_7)$
Molecular Weight: 482,19
CAS: 1185-57-5
EEC-N: 214-686-6

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Iron (III) ammonium citrate red > RE - Pure

RE

Description ..... Red-brown powder
Identification ..... Positive
Sulphate ..... ≤ 0.6 %
As ..... ≤ 5 ppm
Assay (oxidimetric) ..... 20 ÷ 23 % Fe

Code	Size	Packaging	Notes
343441	250 g	Plastic bottle	
343442	1 kg	Plastic bottle	
343443	5 kg	Plastic tank	

Iron (III) ammonium oxalate

• Ferro ammonio ossalato ico • Fer (III) ammonium oxalate • Hierro (III) y amonio oxalato • Eisen (III) ammoniumoxalat

$(NH_4)_3Fe(C_2O_4)_3.3H_2O$
Molecular Weight: 428,08
CAS: 15187-32-3

Iron (III) ammonium oxalate > RE - Pure

RE

Description ..... Green crystals
Identification ..... Positive
pH 10% at 25° C ..... 4.0 ÷ 6.0
Chloride ..... ≤50 ppm
Sulphate ..... ≤50 ppm
Assay (oxidimetric) ..... ≥98 %

Code	Size	Packaging	Notes
343757	1 kg	Plastic bottle	

Iron (III) ammonium sulfate dodecahydrate

• Ferro ammonio solfato ico dodecaidrato • Fer (III) ammonium sulfate dodécahydraté
• Hierro (III) amonio sulfato dodecahidrato • Eisen (III) ammoniumsulfat dodecahydrat

Synonym:
• *Ammonium iron(III) sulfate dodecahydrate*
• *Ammonium ferric sulfate dodecahydrate*

$FeNH_4(SO_4)_2.12H_2O$
Molecular Weight: 482,19
CAS: 7783-83-7
EEC-N: 616-517-5

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-
P362+P364-P337+P313

Iron (III) ammonium sulfate dodecahydrate > RPE - For analysis

RPE

Description Purple semitransparent crystals
Identification ..... Positive
pH sol. 5% at 25° C ..... 1.75 ÷ 2.75
Chloride ..... ≤3 ppm
Phosphate ..... ≤20 ppm
Diluted HCl-ins. matter ..... ≤50 ppm
Nitrate ..... ≤50 ppm
Subst. not ppt NH4OH ..... ≤0.1 %
As ..... ≤2 ppm
Ca ..... ≤10 ppm
Cr ..... ≤100 ppm
Cu ..... ≤10 ppm
K ..... ≤300 ppm
Mg ..... ≤50 ppm
Mn ..... ≤100 ppm
Na ..... ≤150 ppm
Ni ..... ≤20 ppm
Pb ..... ≤30 ppm
Zn ..... ≤20 ppm
Assay (oxidimetric) ..... ≥98 %

Code	Size	Packaging	Notes
451503	100 g	Plastic bottle	
451505	500 g	Plastic bottle	
451507	1 kg	Plastic bottle	
451502	25 kg	Plastic bucket	
451504	50 kg	Plastic bucket	

Iron (III) ammonium sulfate dodecahydrate > RE - Pure

RE

Description Cristalli semitrasparenti violacei
Identification ..... Positive
Assay (oxidimetric) ..... ≥ 98.0 %
Chloride ..... ≤ 0.05 %
Subst. not ppt NH4OH ..... ≤ 0.1 %
Cu ..... ≤ 20 ppm
Zn ..... ≤ 50 ppm

Code	Size	Packaging	Notes
344107	1 kg	Plastic bottle	
344108	5 kg	Plastic tank	

Iron (III) ammonium sulfate solution 33% in nitric acid

• Ferro ammonio solfato ico soluzione 33% in acido nitrico • Fer (III) ammonium sulfate 33% dans l'acide nitrique
• Hierro (III) amonio sulfato solución 33% en acido nítrico • Eisen (III) ammoniumsulfatlösung 33% in Salpetersäure

Classification transport
ONU: 2031
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Iron (III) ammonium sulfate solution 33% in nitric acid > RPE - For analysis

RPE

Description ..... Yellow - brown liquid
Identification ..... Positive
Assay ..... 32 ÷ 34 %

Code	Size	Packaging	Notes
E451521	500 ml	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Iron (III) ammonium sulfate 0.1 mol/l

• Ferro ammonio solfato ico 0.1 mol/l • Fer (III) ammonium sulfate 0.1 mol/l • Hierro (III) amonio sulfato 0.1 mol/l • Eisen (III) ammoniumsulfat 0.1 mol / l

Iron (III) ammonium sulfate 0.1 mol/l > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613001300	1 l	Plastic bottle	Ref Ph.Eur 3001300

Iron (III) ammonium sulfate solution 100 g/l

• Ferro ammonio solfato ico soluzione 100 g/l • Fer (III) ammonium sulfate 100 g/l • Hierro (III) amonio sulfato solución 100 g/l
• Eisen (III) ammoniumsulfat 100 g / l

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Iron (III) ammonium sulfate solution 100 g/l > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611037703	100 ml	Plastic bottle	Ferric ammonium sulfate solution R2 Ref Ph.Eur 1037702
611037702	1 l	Plastic bottle	Ferric ammonium sulfate solution R2 Ref Ph.Eur 1037702

Iron (III) chloride anhydrous sublimed

• Ferro cloruro ico anidro sublimato • Fer (III) chlorure anhydre sublimé • Hierro (III) cloruro anhidro sublimado
• Eisen(III) chlorid wasserfreies sublimiertes

Synonym:
Ferric chloride

$FeCl_3$
Molecular Weight: 162,21
CAS: 7705-08-0
EEC-N: 231-729-4

Classification transport
ONU: 1773
Transport Hazard class: 8
Packing group III

Danger
H302-H315-H318
P264-P305+P351+P338-P310a-P330-P362+P364-P332+P313

Iron (III) chloride anhydrous sublimed > RPE - For analysis

RPE

Description ..... Black powder
Identification ..... Positive
Water-insoluble matter ..... ≤ 1 %
As ..... ≤ 20 ppm
Cu ..... ≤ 0.1 %
Mn ..... ≤ 0.3 %
Ni ..... ≤ 500 ppm
Pb ..... ≤ 200 ppm
Zn ..... ≤ 0.1 %
Assay (oxidimetric) ..... ≥ 98 %

Code	Size	Packaging	Notes
451695	250 g	Glass bottle	
451696	1 kg	Glass bottle	
451692	25 kg	Plastic bucket	

Iron (III) chloride hexahydrate

• Ferro cloruro ico esaidrato • Fer (III) chlorure hexahydraté • Hierro (III) cloruro hexahidratado
• Eisen (III) chloridhexahydrat

Synonym:
Ferric chloride hexahydrate

$FeCl_3.6H_2O$
Molecular Weight: 270,3
CAS: 10025-77-1
EEC-N: 231-729-4

Classification transport
ONU: 3260
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Iron (III) chloride hexahydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Brown pieces
Identification ..... Positive
Ca ..... ≤ 100 ppm
Water-insoluble matter ..... ≤ 100 ppm
Ferrous salt ..... ≤ 20 ppm
Total phosphorus ..... ≤ 100 ppm
Nitrate ..... ≤ 100 ppm
Sulphate ..... ≤ 100 ppm
Cu ..... ≤ 30 ppm
Zn ..... ≤ 30 ppm
Assay (oxidimetric) ..... 97.0 ÷ 102.0 %
Mg ..... ≤ 50 ppm
K ..... ≤ 50 ppm
Na ..... ≤ 0.05 %

Code	Size	Packaging	Notes
451625	100 g	Plastic bottle	
451626	500 g	Plastic bottle	
451627	1 kg	Plastic bottle	

Iron (III) chloride hexahydrate > RE - Pure

RE

Description ..... Grani o blocco giallo scuro
Identification ..... Positive
Acidity (HCl) ..... ≤ 1.2 %
Fe (+2) ..... ≤ 0.9 %
Heavy metals (Pb) ..... ≤ 0.05 %
Assay (oxidimetric) ..... 59 ÷ 61 %
Assay (FeCl3.6H2O) ..... ≥ 99 %

Code	Size	Packaging	Notes
344507	1 kg	Plastic bottle	
344508	2.5 kg	Plastic bottle	
344509	5 kg	Plastic tank	
344504	25 kg	Plastic bucket	

Iron (III) chloride solution 4.5%

• Ferro cloruro ico soluzione 4.5% • Fer (III) chlorure solution 4.5% • Hierro (III) cloruro solución 4.5% • Eisen (III) chlorid 4.5%

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Danger
H290-H318
P234-P280i-P305+P351+P338-P310a-P390-P406

Iron (III) chloride solution 4.5% > RPE - For analysis

RPE

Description ..... Yellow clear liquid
Identification ..... Positive
Density at 20° C ..... 1.019 ÷ 1.025
Assay ..... 4.3 ÷ 4.7 % p/p

Code	Size	Packaging	Notes
E451653	1 l	Bottle	

Iron chloride in solution

• Ferro cloruro in soluzione • Fer chlorure en solution • Hierro Cloruro en solución • Eisen chlorid in Lösung

$FeCl_3$
Molecular Weight: 162,2
CAS: 7705-08-0

Iron chloride in solution > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
616001048	100 ml	Plastic bottle	Ferric Chloride CS
616001047	500 ml	Plastic bottle	Ferric Chloride CS

Iron (III) citrate

• Ferro citrato ico • Fer (III) citrate • Hierro (III) citrato • Eisen (III) citrat

Synonym:
• *Iron(III) citrate tribasic*
• *Ferric citrate*

$C_6H_8O_7.xFe$
Molecular Weight: 262,96
CAS: 2338-05-8
EEC-N: 219-045-4

Iron (III) citrate > RE - Pure

RE

Description ..... Red-brown crystals
Identification ..... Positive
Assay (iodometric) ..... 18 ÷ 20 % Fe

Code	Size	Packaging	Notes
344201	1 kg	Plastic bottle	
344202	5 kg	Plastic tank	

Product specifications are subject to changes. Please visit our website for updates.

Iron (III) oxide

• Ferro ossido ico • fer (III) oxyde • Hierro (III) óxido • Eisen (III) oxid

Synonym:
Ferric oxide

Fe_2O_3
Molecular Weight: 159,7
CAS: 1309-37-1
EEC-N: 215-168-2

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Iron (III) oxide > RPE - For analysis

RPE

Description	Red brown powder	Chloride	≤100 ppm	Cr	≤10 ppm	Na	≤200 ppm
Identification	Positive	Sulphate	≤100 ppm	Cu	≤50 ppm	Ni	≤50 ppm
Loss on drying	≤0.2 %	Ca	≤100 ppm	K	≤50 ppm	Pb	≤50 ppm
Diluted HCl-ins. matter	≤100 ppm	Cd	≤10 ppm	Mg	≤50 ppm	Zn	≤50 ppm
Total nitrogen	≤50 ppm	Co	≤50 ppm	Mn	≤50 ppm	Assay (oxidimetric)	≥99 %

Code	Size	Packaging	Notes
451824	100 g	Plastic bottle	
451826	500 g	Plastic bottle	

Iron (III) nitrate nonahydrate

• Ferro nitrato ico nonaidrato • Fer(III) nitrate nonahydraté • Hierro (III) nitrato nonahidratado • Eisen (III) nitratnonahydrat

Synonym:
Ferric nitrate nonahydrate

$Fe(NO_3)_3.9H_2O$
Molecular Weight: 404,04
CAS: 7782-61-8
EEC-N: 233-899-5

Classification transport
ONU: 1466
Transport Hazard class: 5.1
Packing group III

Danger
H272-H315-H319-H335
P210-P261-P280-P304+P340-P305+P351+P338-P403+P233

Iron (III) nitrate nonahydrate > RPE - For analysis - ACS

RPE

Description	Purplish crystals deliquescent	Subst. not ppt NH4OH	≤0.1 %	Chloride	≤5 ppm	Assay (oxidimetric)	98.0 ÷ 101.0 %
Identification	Positive	Water-insoluble matter	≤50 ppm	Sulphate	≤100 ppm		

Code	Size	Packaging	Notes
451723	100 g	Plastic bottle	
451725	500 g	Plastic bottle	
451727	1 kg	Plastic bottle	
451722	25 kg	Plastic bucket	

Iron (III) sulfate

• Ferro solfato ico • Fer (III) sulfate • Hierro (III) sulfato • Eisen (III) sulfat

Synonym:
Ferric sulfate hydrate

$Fe_2(SO_4)_3.nH_2O$
Molecular Weight: 399,88
CAS: 15244-10-7
EEC-N: 233-072-9

Danger
H302-H315-H318
P264-P305+P351+P338-P310a-P330-P362+P364-P332+P313

Iron (III) sulfate > RPE - For analysis

RPE

Description	Yellow-green powder	Chloride	≤ 0.1 %	Pb	≤20 ppm
Identification	Positive	As	≤3 ppm	Assay (oxidimetric)	20 ÷ 23 % Fe

Code	Size	Packaging	Notes
451925	100 g	Plastic bottle	
451926	500 g	Plastic bottle	
451927	1 kg	Plastic bottle	

Isoamyl acetate

• Isoamilo acetato • Isoamyle acétate • Isoamilo Acetato • Isoamylacetat

Synonym:
• *Acetic acid 3-methylbutyl ester*
• *Isopentyl acetate*

$C_7H_{14}O_2$
Molecular Weight: 130
CAS: 123-92-2
EEC-N: 204-662-3

Classification transport
ONU: 1104
Transport Hazard class: 3
Packing group III

Warning
H226-HEU066
P210-P241-P280-P303+P361+P353-P403+P235-P501a

Isoamyl acetate > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Water miscibility ..... Conform
Alcohol miscibility ..... Complete
Density at 20° C ..... 0.870 ÷ 0.874
Refractive index at 20°C 1.3963 ÷ 1.4043
Boiling point ..... 140.5 ÷ 143.5 ° C
Water (K.F.) ..... ≤500 ppm
Residue on evaporation ..... ≤20 ppm
Acidity (acetic acid) ..... ≤50 ppm
Assay (GLC) ..... ≥99.0 %

Code	Size	Packaging	Notes
417781	250 ml	Glass bottle	
417782	1 l	Glass bottle	

Isoamyl acetate > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.870 ÷ 0.874
Refractive index at 20°C ..... 1.3943 ÷ 1.4063
Boiling point ..... 140.5 ÷ 143.5 ° C
Water (K.F.) ..... ≤0.5 %
Residue on evaporation ..... ≤100 ppm
Acidity (acetic acid) ..... ≤0.1 %
Assay (GLC) ..... ≥ 98.0 %

Code	Size	Packaging	Notes
313251	1 l	Glass bottle	
313252	15 kg	Metal drum	

Isoamyl acetate > RE - Pure - For synthesis

RE

Appearance ..... Clear liquid
Refractive index at 20°C ..... 1.399 - 1.403
Density d20/4 ..... 0.871 - 0.874
Colour ..... ≤ 15 Hazen
Free acid (as CH3COOH) ..... ≤ 0.02 %
Assay (GC) ..... ≥ 99.0 %
Isoamyl alcohol ..... ≤ 0.6 %

Code	Size	Packaging	Notes
P6120028	5 l	Plastic tank	

Isoamyl alcohol

• Alcole isoamilico • Alcool isoamylique • Alcohol isoamílico • Isoamylalkohol

Synonym:
3-Methyl-1-butanol

$(CH_3)_2CHCH_2.CH_2OH$
Molecular Weight: 88,15
CAS: 123-51-3
EEC-N: 204-633-5

Classification transport
ONU: 1105
Transport Hazard class: 3
Packing group III

Danger
H226-H332-H315-H318-H335
P210-P280-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P362+P364-P403+P233

Isoamyl alcohol > RS - For analysis according to Gerber

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.810 ÷ 0.814

Code	Size	Packaging	Notes
E413903	500 ml	Glass bottle	With indicator
413892	1 l	Glass bottle	Without indicator

Isoamyl alcohol > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.805 ÷ 0.813
Refractive index at 20°C 1.4023 ÷ 1.4083
Boiling point ..... 130.5 ÷ 132.5 °C
Water (K.F.) ..... ≤0.5 %
Acidity ..... ≤0.002 meq/g
Acids and esthers ..... ≤0.2 %
Residue on evaporation ..... ≤30 ppm
Carbonyl (HCOH) ..... ≤0.1 %
Assay (GLC) ..... ≥98.5 %

Code	Size	Packaging	Notes
413801	500 ml	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Isoamyl alcohol > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.805 ÷ 0.813
Refractive index at 20°C 1.4023 ÷ 1.4083
Boiling point ..... 130.5 ÷ 132.5 °C
Water (K.F.) ..... ≤0.5 %
Acidity ..... ≤0.002 meq/g
Acids and esthers ..... ≤0.2 %
Residue on evaporation ..... ≤30 ppm
Carbonyl (HCOH) ..... ≤0.1 %
Assay (GLC) ..... ≥97 %

Code	Size	Packaging	Notes
413832	1 l	Glass bottle	
413836	2.5 l	Glass bottle	
413833	22 kg	Metal drum	

Isoamyl alcohol > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.807 ÷ 0.817
Refractive index at 20°C 1.4038 ÷ 1.4098
Boiling point ..... 132 ± 1.5 ° C
Water (K.F.) ..... ≤0.2 %
Residue on evaporation ..... ≤0.1 %
Acids and esthers ..... ≤0.5 %

Code	Size	Packaging	Notes
308001	1 l	Glass bottle	
308003	22 kg	Metal drum	

Isobutanol

• Isobutanolo • Isobutanol • Isobutanol • Isobutanol

Synonym:
Isobutyl alcohol

$(CH_3)_2CHCH_2OH$
Molecular Weight: 74,12
CAS: 78-83-1
EEC-N: 201-148-0

Classification transport
ONU: 1212
Transport Hazard class: 3
Packing group III

Danger
H226-H315-H318-H335-H336
P210-P280-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P362+P364-P403+P233

Isobutanol > RS - Anhydrous - For analysis

RS

Refractive index at 20°C ..... 1.394 - 1.398
Colour ..... ≤ 10 Hazen
Water content (K.F.) ..... ≤ 200 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
Assay (GC) ..... ≥ 99.5 %
1-butanol ..... ≤ 0.4 %
Free acid (as CH3COOH) ..... ≤ 40 mg/Kg

Code	Size	Packaging	Notes
P0531016	1 l	Glass bottle	

Isobutanol > RPE - For analysis - ACS

RPE

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification (I.R.) ..... Positive
Water solubility ..... Conform
Density at 20° C ..... 0.801 ÷ 0.803
Refractive index at 20°C 1.3945 ÷ 1.3975
Boiling point ..... 105 ÷ 109 ° C
Water (K.F.) ..... ≤0.1 %
Residue on evaporation ..... ≤10 ppm
Carbonyl Compounds (CO) ..... ≤100 ppm
Acidity (acetic acid) ..... ≤100 ppm
Alcalinity (NH3) ..... ≤10 ppm
Indole base ..... ≤0.1 ppm
Assay (GLC) ..... ≥99 %

Code	Size	Packaging	Notes
414211	1 l	Glass bottle	
414213	22 kg	Metal drum	
414214	200 l	Metal drum	

Isobutanol > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.800 ÷ 0.804
Refractive index at 20°C 1.3930 ÷ 1.3990
Boiling point ..... 105.0 ÷ 109.0 °C
Water (K.F.) ..... ≤0.1 %
Residue on evaporation ..... ≤50 ppm
Acidity (acetic acid) ..... ≤ 40 ppm
Assay (GLC) ..... ≥99 %
1-Butanol ..... ≤ 0.4 %

Code	Size	Packaging	Notes
308301	1 l	Glass bottle	
308303	22 kg	Metal drum	

Isobutyl acetate

• Isobutile acetato • Isobutyle acétate • Isobutilo acetato • Isobutylacetat

Synonym:
Acetic acid isobutyl ester

$(CH_3)_2CHCH_2OCOCH_3$
Molecular Weight: 116,16
CAS: 110-19-0
EEC-N: 203-745-1

Classification transport
ONU: 1213
Transport Hazard class: 3
Packing group II

Danger
H225-HEU066
P210-P241-P280-P303+P361+P353-P403+P235-P501a

Isobutyl acetate > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.867 ÷ 0.877
Refractive index at 20°C ..... 1.3851 ÷ 1.3951
Boiling point ..... 114 ÷ 118.5 ° C
Water (K.F.) ..... ≤0.1 %
Residue on evaporation ..... ≤100 ppm
Acidity (acetic acid) ..... ≤500 ppm
Assay (GLC) ..... ≥98.0 %

Code	Size	Packaging	Notes
431721	500 ml	Glass bottle	

Isobutyl acetate > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.867 ÷ 0.877
Refractive index at 20°C ..... 1.3851 ÷ 1.3951
Boiling point ..... 114 ÷ 118.5 ° C
Water (K.F.) ..... ≤0.2 %
Residue on evaporation ..... ≤100 ppm
Acidity (acetic acid) ..... ≤500 ppm
Assay (GLC) ..... ≥97 %

Code	Size	Packaging	Notes
325631	1 l	Glass bottle	
325633	24 kg	Metal drum	

Isobutylacetone ▶ Methylisoamyl ketone

Isobutyl alcohol ▶ Isobutanol

Isohexane

• iso-Esano • Isohexane • Isohexano • Isohexan

Synonym:
2-Methylpentane

C_6H_{14}
Molecular Weight: 86,18
CAS: 107-83-5
EEC-N: 203-523-4

Classification transport
ONU: 1208
Transport Hazard class: 3
Packing group II

Danger
H225-H315-H336-H304-H411
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Isohexane > RS - For HPLC - Isocratic Grade

RS

Clear,colourless liq.appearance ..... Conform
Colour ..... ≤ 10 Apha
Identification (IR) ..... Conform
Water content (K.F.) ..... ≤ 150 mg/Kg
Non volatile residue ..... ≤ 5 mg/Kg
n-hexane ..... ≤ 6 %
Total C6H14 isomers (GC) ..... ≥ 80 %
UV transmittance at 200 nm ..... ≥ 25 %
UV transmittance at 210 nm ..... ≥ 55 %
UV transmittance at 220 nm ..... ≥ 80 %
UV transmittance at 230 nm ..... ≥ 90 %
UV transmittance at 250 nm ..... ≥ 98 %

Code	Size	Packaging	Notes
445152	1 l	Glass bottle	
445151	2.5 l	Glass bottle	

Isohexane > RS - ATRASOL - For trace analysis, Suitable for Hydrocarbon index determination

RS

Appearance ..... Clear colourless liquid
Identification (IR) ..... Conform
Water content (K.F.) ..... ≤ 50 mg/Kg
Colour ..... ≤ 5 Hazen
n-hexane ..... ≤ 3 %
Total C6H14 isomers (GC) ..... ≥ 95 %
Non volatile residue ..... ≤ 2 mg/Kg
GC-ECD.Individual peak (Lindane) .≤ 3 ng/L
Ret.range 1,2,4-trichlorobenzene to decachlorobiphenyle
GC-FID.Hydrocarbon oil index...≤ 0.05 mg/l
Retention time n-decane - n-tetracontane
GC-FID.Individual peak (C10-C40)..≤ 5 µg/l
Total sulphur (S) ..... ≤ 1 ppm

Code	Size	Packaging	Notes
P6263216	1 l	Glass bottle	
P6263221	2.5 l	Glass bottle	

According to NF-EN-ISO9377-2 for hydrocarbon index determination

Product specifications are subject to changes. Please visit our website for updates.

Isohexane > RS - PESTIPUR - For pesticide analysis

RS

Appearance ..... Clear colourless liquid
Water content (K.F.) ..... ≤ 100 mg/Kg
Colour ..... ≤ 10 Hazen
Identification (IR) ..... Conform
n-hexane ..... ≤ 5 %
Total C6H14 isomers (GC) ..... ≥ 80 %
Non volatile residue ..... ≤ 5 mg/Kg
GC-ECD.Individual peak (Lindane) .. ≤ 3 ng/l
Retention time trichlorobenzene to mirex
Ret.range 1,2,4-trichlorobenzene to decachlorobiphenyle

Code	Size	Packaging	Notes
447131	1 l	Glass bottle	
447132	2.5 l	Glass bottle	

For chlorinated compounds analysis

Isohexane > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Identification (I.R.) ..... Conform
Water (K.F.) ..... ≤ 150 ppm
Residue on evaporation ..... ≤ 10 ppm
Aromatic compounds ..... ≤ 50 ppm
Total C6H14 isomers (GC) ..... ≥ 80 %
n-Hexane ..... ≤ 5 %
Methylcyclopentane ..... ≤ 16 %
Total sulphur (S) ..... ≤ 1 ppm

Code	Size	Packaging	Notes
447311	1 l	Glass bottle	
447312	2.5 l	Glass bottle	

Isooctane

• Isoottano • Isooctane • Isooctano • Isooctan

Synonym:
2,2,4-Trimethylpentane

$(CH_3)_3CCH_2CH(CH_3)_2$
Molecular Weight: 114,23
CAS: 540-84-1
EEC-N: 208-759-1

Classification transport
ONU: 1262
Transport Hazard class: 3
Packing group II

Danger
H225-H315-H336-H304-H410
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Isooctane > RS - For HPLC - Isocratic Grade

RS

Description ..... Clear liquid
Colour ..... ≤ 10 hazen
Density at 20° C ..... 0.687 ÷ 0.697
Identification ..... Positive
Refractive index at 20°C. 1.3885 ÷ 1.3945
Water ..... ≤ 50 ppm
Boiling point ..... 98.2 ÷ 100.2 ° C
Not volatile residue ..... ≤ 5 ppm
Water (K.F.) ..... ≤50 ppm
Free acids (CH3COOH) ..... ≤ 20 ppm
Residue on evaporation ..... ≤2 ppm
Aromatic compounds ..... ≤ 5 ppm
Acidity or alkalinity ..... ≤0.0002 meq/g
Assay (GLC) ..... ≥ 99.5 %
U.V. Transmittance
at 220 nm ..... ≥ 70 %
at 230 nm ..... ≥ 85 %
at 205 nm ..... ≥10 %
at 260 nm ..... ≥ 95 %
at 215 nm ..... ≥60 %
at 270 nm ..... ≥ 97 %
at 240 nm ..... ≥95 %
Filtered at 0.2 µm

Code	Size	Packaging	Notes
412441000	1 l	Glass bottle	
412442000	2.5 l	Glass bottle	

Isooctane > RS - PESTIPUR - For pesticide analysis

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Colour ..... ≤ 10 hazen
Water ..... ≤ 50 ppm
Free acids (CH3COOH) ..... ≤ 20 ppm
Not volatile residue ..... ≤ 2 ppm
GC-ECD (Lindano) ..... ≤ 3 ng/l
GC-NPD (Ethylparation) ..... ≤ 3 ng/l
Assay (GLC) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
456791	1 l	Glass bottle	
456792	2.5 l	Glass bottle	

Isooctane > RS - SPECTROSOL - For optical spectroscopy - Reag.Ph.Eur.

RS

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Density at 20° C ..... 0.691 ÷ 0.696
Refractive index at 20°C. 1.3910 ÷ 1.3930
Boiling point ..... 98 ÷ 100 ° C
Acidity or alkalinity ..... ≤0.0002 meq/g
Water (K.F.) ..... ≤50 ppm
Residue on evaporation ..... ≤5 ppm
Assay (GLC) ..... ≥99.5 %
Aromatic compounds ..... ≤5 ppm
U.V. Transmittance
at 255 nm ..... ≥98 %
at 210 nm ..... ≥35 %
at 220 nm ..... ≥72 %
at 230 nm ..... ≥85 %
UV Absorbance from 250 nm to 420 nm.. ≤ 0.01 AU

Code	Size	Packaging	Notes
456754	1 l	Glass bottle	
456753	2.5 l	Glass bottle	

Isooctane > RS - Standard for refractometry

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.687 ÷ 0.697
Refractive index at 20°C 1.3905 ÷ 1.3925

Code	Size	Packaging	Notes
456641	100 ml	Glass bottle	

Isooctane > RPE - For analysis - ACS

RPE

Description ..... Clear colourless liquid
Identification (I.R.) ..... Positive
Colour ..... ≤ 10 APHA
Density at 20°C ..... 0.691 - 0.696
Refractive index at 20°C ..... 1.391 - 1.393
Water-soluble titrable acid ≤ 0.0003 meq/g
Residue on evaporation ..... ≤ 0.001 %
Sulfur compounds (S) ..... ≤ 0.0005 %
Assay (CPG) ..... ≥ 99.5 %
Distillation range 95% distils between ..... 98 - 100 °C
Water (K.F.) ..... ≤ 100 ppm
Aromatics ..... ≤ 10 ppm

Code	Size	Packaging	Notes
456734	1 l	Glass bottle	
456732	2.5 l	Glass bottle	

Isooctane > RE - Pure - For synthesis

RE

Refractive index at 20°C ..... 1.389 - 1.393
Water content (K.F.) ..... ≤ 150 mg/Kg
Non volatile residue ..... ≤ 20 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 99 %
Aromatic compounds ..... ≤ 50 mg/Kg
Identification (IR) ..... Conform
Total sulphur (S) ..... ≤ 10 ppm

Code	Size	Packaging	Notes
P0630240	10 l	Metal drum	
P0630268	200 l	Metal drum	

Isooctane > RE - ASTM

RE

Description ..... Clear liquid
Identification ..... Positive
Density at 20°C ..... 0.687 ÷ 0.697
Refractive index at 20°C 1.3885 ÷ 1.3945
Boiling point ..... 98.2 ÷ 100.2 °C
Residue on evaporation ..... ≤ 20 ppm
Water (K.F.) ..... ≤ 150 ppm
Aromatics ..... ≤ 50 ppm
Total sulphur ..... ≤ 10 ppm
Lead ..... ≤ 0.5 mg/l
n-heptane (ASTM) ..... ≤ 0.10 % v/v
Isooctane (ASTM) ..... ≥ 99.75 % v/v
Assay (GLC) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
528963	2.5 l	Glass bottle	
528960	5 l	Plastic tank	
528961	25 l	Metal drum	
528962	200 l	Metal drum	

Suitable for ASTM methods D2700 and D2699

Isopar G

• Isopar G • Isopar G • Isopar G • Isopar G

CAS: 90622-57-4
EEC-N: 292-459-0

Classification transport
ONU: 3295
Transport Hazard class: 3
Packing group III



Danger
H226-H304-H413
P210-P241-P280-P301+P310a-P303+P361+P353-P403+P235

Isopar G > RS - RSE - For electronic use

RS

Appearance ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Density at 15°C ..... 0.745 ÷ 0.756
Refractive index at 20°C ..... 1.416 ÷ 1.418
Distillation range ..... 159 ÷ 176 °C
Residue on evaporation ..... ≤ 10 ppm
Aromatics ..... ≤ 0.01 %

Code	Size	Packaging	Notes
526151	2.5 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Isopentane

• Isopentano • Isopentane • Isopentano • Isopentan

Synonym:
2-Methylbutane

$(CH_3)_2CHCH_2CH_3$
Molecular Weight: 72,15
CAS: 78-78-4
EEC-N: 201-142-8

Classification transport
ONU: 1265
Transport Hazard class: 3
Packing group I



Danger
H224-H336-H304-H411-HEU066
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Isopentane > RS - Anhydrous - For analysis

RS

Refractive index at 20°C ..... 1.352 - 1.356
Water content (K.F.) ..... ≤ 50 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
Aromatic compounds ..... ≤ 20 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 95 %
Identification (IR) ..... Conform
Total sulphur (S) ..... ≤ 2 ppm

Code	Size	Packaging	Notes
P0651016	1 l	Glass bottle	

Isopentane > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Identification (I.R.) ..... Conform
Refractive index at 20°C ..... 1.352 ÷ 1.356
Water (K.F.) ..... ≤ 150 ppm
Residue on evaporation ..... ≤ 10 ppm
Aromatic compounds ..... ≤ 20 ppm
Assay (CPG) ..... ≥ 95 %
Total sulphur (S) ..... ≤ 2 ppm

Code	Size	Packaging	Notes
524391	1 l	Glass bottle	

Isopentane > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Colour ..... ≤ 10 APHA
Density at 20°C ..... 0.610 ÷ 0.630
Refractive index at 20°C. 1.3507 ÷ 1.3607
Boiling point ..... 27 ÷ 28.5 °C
Water (K.F.) ..... ≤ 200 ppm
Residue on evaporation ..... ≤ 20 ppm
Aromatics ..... ≤ 50 ppm
Total sulphur (S) ..... ≤ 2 ppm
Assay ..... ≥ 95.0 %

Code	Size	Packaging	Notes
528492	1 l	Glass bottle	
528491	5 l	Metal tank	

Isopentyl acetate ▶ Isoamyl acetate

Isopropanol ▶ Propan-2-ol

Isopropyl acetate

• Isopropile acetato • Isopropyle acétate • Isopropilo acetato • Isopropylacetat

Synonym:
Acetic acid isopropyl ester

$CH_3COOCH(CH_3)_2$
Molecular Weight: 102,13
CAS: 108-21-4
EEC-N: 203-561-1

Classification transport
ONU: 1220
Transport Hazard class: 3
Packing group II



Danger
H225-H319-H336-HEU066
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Isopropyl acetate > RPE - For analysis

RPE

Description ..... Clear liquid
Water content (K.F.) ..... ≤ 100 mg/Kg
Identification ..... Positive
Non volatile residue ..... ≤ 10 mg/Kg
Colour ..... ≤10 APHA
Density at 20° C ..... 0.869 ÷ 0.875
Assay (GC) ..... ≥ 99 %
Refractive index at 20°C. 1.3760 ÷ 1.3780
Free acid (as CH3COOH) ..... ≤ 50 mg/Kg
Boiling point ..... 88 ÷89 ° C
2-Propanol ..... ≤ 1000 mg/Kg
Water (K.F.) ..... ≤0.05 %
Residue on evaporation ..... ≤10 ppm
Acidity (acetic acid) ..... ≤50 ppm
Assay (GLC) ..... ≥99.8 %

Code	Size	Packaging	Notes
474821	250 ml	Glass bottle	

Isopropyl acetate > RE - Pure

RE

Refractive index at 20°C.....1.375 - 1.379
Appearance .....Clear liquid
Water content (K.F.) .....≤ 800 mg/Kg
Non volatile residue.....≤ 20 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC)..... ≥ 99 %
Free acid (as CH3COOH).....≤ 50 mg/Kg
2-Propanol.....≤ 1000 mg/Kg

Code	Size	Packaging	Notes
P0890228	5 l	Plastic tank	
P0890240	10 l	Metal tank	
P0890268	200 l	Metal drum	

Isopropylamine

• Isopropilammina • Isopropylamine • Isopropilamina • Isopropylamin

Synonym:
2-Aminopropane

$(CH_3)_2CHNH_2$
Molecular Weight: 59,11
CAS: 75-31-0
EEC-N: 200-860-9

Classification transport
ONU: 1221
Transport Hazard class: 3
Packing group I



Danger
H224-H301-H311-H331-H314-H335
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P311a-P305+P351+P338-P361+P364-P403+P233

Isopropylamine > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.687 ÷ 0.693
Residue on evaporation .....≤50 ppm
Water (K.F.).....≤0.1 %
Assay (GLC) .....≥99.5 %

Code	Size	Packaging	Notes
474756	1 l	Glass bottle	

Isopropyl ether ▶ Diisopropylether

Isopropyl myristate

• Miristato di isopropile • Myristate d'isopropyle • Isopropilo Miristato • Isopropylmyristat

Synonym:
• *Isopropyl tetradecanoate*
• *Myristic acid isopropyl ester*

$C_{17}H_{34}O_2$
Molecular Weight: 2.704.507
CAS: 110-27-0
EEC-N: 203-751-4

Isopropyl myristate > RS - For synthesis

RS

Clear,colourless liq.appearance .....Conform
Clear,colourless solution appearance ..... Conform
Alkohol miscibility..... Miscible
Dichloromethane miscibility..... Miscible
Identification B .....Conform
Refractive index at 20°C.....1.434 - 1.437
Assay (GC) (C17H34O2) ..... ≥ 90 %
Viscosity at 20°C.....5 - 6 mPa.s
Acid number ..... ≤ 1 mg KOH /g
Iodine number.....≤ 1 g I2/100g
Saponification number202 - 212 mg KOH/g
Water content (K.F.) ..... ≤ 0.1 % m/m
Sulfuric ashes ..... ≤ 0.1 %

Code	Size	Packaging	Notes
P6070268	205 l	Metal drum	

Product specifications are subject to changes. Please visit our website for updates.

Kaolin washed

• Caolino lavato • Kaolin lavé • Caolín lavado • Gewaschenes Kaolin

Synonym:
• *Halloysite nanoclay*
• *Kaolin clay*

CAS: 1332-58-7
EEC-N: 310-194-1

Kaolin washed > RE - Pure

RE

Description ..... White hazel powder Identification ..... Positive Loss on drying ..... ≤ 1 % Loss on ignition ..... 11 - 13 %

Code	Size	Packaging	Notes
332573	5 kg	Plastic tank	
332574	25 kg	Plastic bucket	

Karl Fischer reagent 1 component 2 mg H_2O/ml

• Karl Fischer reattivo 2 mg H_2O/ml monocomponente • Réactif de Karl Fischer 1 composant 2 mgH^2O/ml • Karl Fischer reactivo 1 compuesto 2 mgH_2O/ml
• Karl Fischer Reagenz 1 Komponente 2 mg H 2 0 / ml

Danger
H315-H319-H360D-H373-HA26
P280-P305+P351+P338-P308+P313-P362+P364-P332+P313-P337+P313

Karl Fischer reagent 1 component 2 mg H_2O/ml > RS - ERBAqua - Volumetric titration

RS

Code	Size	Packaging	Notes
570021	1 l	Glass bottle	
570251	1 l	Glass bottle	New formulation. Not regulated for transport.

Pyridine free. Suitable for general purpose or aldehydes and ketones - To be used with anhydrous solvents (methanol, chloroform)

Karl Fischer reagent 1 component 5 mg H_2O/ml

• Karl Fischer reattivo 5 mg H_2O/ml monocomponente • Réactif de Karl Fischer 1 composant 5 mg H_2O/ml • Karl Fischer reactivo 1 compuesto 5 mgH2O/ml
• Karl-Fischer-Reagenz 1 Komponente 5 mg H2O / ml

Danger
H315-H319-H360D-H373-HA26
P280-P305+P351+P338-P308+P313-P362+P364-P332+P313-P337+P313

Karl Fischer reagent 1 component 5 mg H2O/ml > RS - ERBAqua - Volumetric titration

RS

Code	Size	Packaging	Notes
570011	1 l	Glass bottle	
570241	1 l	Glass bottle	New formulation. Not regulated for transport.

Pyridine free. Suitable for general purpose or aldehydes and ketones - To be used with anhydrous solvents (methanol, chloroform)

Karl Fischer titrant 1 component 5 mg H_2O/ml for aldehydes and ketones

• Karl Fischer reattivo 5 mg H_2O/ml monocomponente per aldeidi e chetoni • Karl Fischer titrant 1 composant 5 mg H_2O/ml pour aldehydes et cétones
• Karl Fischer reactivo de valoración 1 compuesto 5 mg H_2O/ml para aldehídos y cet
• Karl Fischer titriert 1 Komponente 5 mg H_2O / ml für Aldehyde und Ketone

Classification transport
ONU: 1188
Transport Hazard class: 3
Packing group III

Danger
H226-H312-H315-H319-H360FD-H373-HA26
P210-P241-P280-P303+P361+P353-P305+P351+P338-P308+P313

Karl Fischer titrant 1 component 5 mg H_2O/ml for aldehydes and ketones > RS - ERBAqua - Volumetric titration

RS

Code	Size	Packaging	Notes
570081	1 l	Glass bottle	

To be used with one component Karl Fischer solvent for aldehydes and ketones

Karl Fischer solvent for aldehydes and ketones one component

• Karl Fischer solvente monocomponente per aldeidi e chetoni • Karl Fischer solvant 1 composant pour aldéhydes et cétones
• Karl Fischer solvente 1 compuesto para aldehídos y cetonas • Karl Fischer Lösungsmittel 1 Komponente für Aldehyde und Ketone

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group III



Danger
H226-H302-H315-H319-H351-H360FD-H372-HEU301-HA26
P210-P241-P280-P303+P361+P353-P305+P351+P338-P308+P313

Karl Fischer solvent for aldehydes and ketones one component > RS - ERBAqua - Volumetric titration

RS

Code	Size	Packaging	Notes
570041	1 l	Glass bottle	

To be used with one component Karl Fischer reagents. For determination of water in aldehydes and ketones. Especially suitable for higher molecular weight products and mixtures containing non - polar constituents.

Karl Fischer solvent for oils one component

• Karl Fischer solvente monocomponente per olii • Karl Fischer solvant pour huiles 1 composant • Karl Fischer solvente para los aceites 1 compuesto
• Karl Fischer Lösungsmittel für 1 Komponentenöle

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group III



Danger
H226-H302-H315-H318-H351-H360D-H370-H372-HEU301-HA26
P210-P280-P303+P361+P353-P305+P351+P338-P310a-P330-P362+P364

Karl Fischer solvent for oils one component > RS - ERBAqua - Volumetric titration

RS

Code	Size	Packaging	Notes
570031	1 l	Glass bottle	

To be used with Karl Fischer reagents. For determination of water in oils or other non polar compounds.

Karl Fischer titrant 2 component 2 mg H_2O/ml

• Karl Fischer titolante bicomponente 2 mg H_2O/ml • Karl Fischer titrant 2 composants 2 mg H_2O/ml
• Karl Fischer reactivo de valoración 2 compuestos 2 mg H_2O/ml • Karl Fischer titriert 2 Komponenten 2 mg H_2O / ml - nicht hygroskopisch

Classification transport
ONU: 1307
Transport Hazard class: 3
Packing group III



Danger
H226-H315-H319-H335-H373-H304
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Karl Fischer titrant 2 component 2 mg H_2O/ml > RS - ERBAqua - Volumetric titration

RS

Code	Size	Packaging	Notes
570061	1 l	Glass bottle	

To be used with two component Karl Fischer solvent

Karl Fischer titrant 2 component 5 mg H_2O/ml

• Karl Fischer titolante bicomponente 5 mg H_2O/ml • Karl Fischer titrant 2 composants 5 mg H_2O/ml
• Karl Fischer reactivo de valoración 2 compuestos 5 mg H_2O/ml • Karl Fischer Titriermittel 2 Komponenten 5 mg H_2O / ml - nicht hygroskopisch

Classification transport
ONU: 1307
Transport Hazard class: 3
Packing group III



Danger
H226-H315-H319-H335-H373-H304
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Karl Fischer titrant 2 component 5 mg H_2O/ml > RS - ERBAqua - Volumetric titration

RS

Description ..... Brown red liquid Identification ..... Positive Water equivalent ..... ≥ 5.0 mg/ml

Code	Size	Packaging	Notes
570051	1 l	Glass bottle	

To be used with two component Karl Fischer solvent

Product specifications are subject to changes. Please visit our website for updates.

Karl Fischer solvent 2 component

• Karl Fischer solvente bicomponente • Karl Fischer solvent 2 component • Karl Fischer solvente 2 componentes • Karl Fischer Lösungsmittel 2 Komponente

Classification transport
ONU: 1230
Transport Hazard class: 3
Packing group II



Danger
H225-H301-H331-H315-H319-H370-H373
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Karl Fischer solvent 2 component > RS - ERBAqua - Volumetric titration

RS

Code	Size	Packaging	Notes
570071	1 l	Glass bottle	

To be used with a two component Karl Fischer titrant 570051 or 570061

Karl Fischer solvent 2 component for aldehydes and ketones - Methanol free

• Karl Fischer solvente bicomponente per aldeidi e chetoni - senza metanolo • Karl Fischer solvant 2 composants pour aldéhydes et cétones - Sans méthanol
• Karl Fischer solvente 2 compuestos para aldehídos y cetonas - sin metanol • Karl Fischer Solvent 2 Komponenten für Aldehyde und Ketone - Methanolfrei

Classification transport
ONU: 1188
Transport Hazard class: 3
Packing group III



Danger
H226-H312-H315-H319-H360FD-HA26
P210-P241-P280-P303+P361+P353-P305+P351+P338-P308+P313

Karl Fischer solvent 2 component for aldehydes and ketones - Methanol free > RS - ERBAqua - Volumetric titration

RS

Code	Size	Packaging	Notes
570091	1 l	Glass bottle	

To be used with a two component Karl Fischer titrant 570051 or 570061. For samples containing aldehydes and ketones

Karl Fischer solvent for oils 2 component

• Karl Fischer solvente bicomponente per olii • Karl Fischer solvant pour huiles 2 composant • Karl Fischer solvente para los aceites 2 compuestos
• Karl Fischer Lösungsmittel für 2 Komponentenöle

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II



Danger
H225-H302-H315-H318-H351-H361d-H370-H372-H373-HEU301
P210-P280-P303+P361+P353-P305+P351+P338-P310a-P330-P362+P364

Karl Fischer solvent for oils 2 component > RS - ERBAqua - Volumetric titration

RS

Code	Size	Packaging	Notes
570101	1 l	Glass bottle	

To be used with a two component Karl Fischer titrant 570051 or 570061. For determination of water in oils or other non-polar coumpounds

Karl Fischer 2 component buffered solvent

• Karl Fischer solvente bicomponente tamponato • Karl Fischer solvant 2 composant avec tampon • Karl Fischer solvente 2 compuestos con tampón
• Karl Fischer Pufferlösung

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H319-H360D-H370-H373-HA26
P210-P241-P280-P303+P361+P353-P305+P351+P338-P308+P313

Karl Fischer 2 component buffered solvent > RS - ERBAqua - Volumetric titration

RS

Code	Size	Packaging	Notes
570111	500 ml	Glass bottle	

Karl Fischer anolyte solution

• Karl Fischer soluzione anolita • Karl Fischer anolyte solution • Karl Fischer catolito solución • Karl-Fischer-Anolytlösung - ohne Pyridin

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II



Danger
H225-H302-H315-H319-H351-H360D-H370-H372-HA26-HEU301
P210-P241-P280-P303+P361+P353-P305+P351+P338-P308+P313

Karl Fischer anolyte solution > RS - ERBAqua - Coulometric titration

RS

Code	Size	Packaging	Notes
570121	500 ml	Glass bottle	

Karl Fischer anolyte solution for ketones and aldehydes - Methanol free

• Karl Fischer soluzione anolita per aldeidi e chetoni • Karl Fischer anolyte solution pour cétones et aldéhydes - Sans méthanol
• Karl Fischer anolito solución para aldehídos y cetonas - libre de metanol • Karl-Fischer-Anolytlösung für Ketone und Aldehyde - ohne Methanol

Classification transport
ONU: 1992
Transport Hazard class: 3
Packing group III



Danger
H226-H302-H331-H315-H319-H351-H360FD-H372-HEU301-HA26
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P308+P313-P330-P362+P364-P403+P233

Karl Fischer anolyte solution for ketones and aldehydes - Methanol free > RS - ERBAqua - Coulometric titration

RS

Code	Size	Packaging	Notes
570161	500 ml	Glass bottle	

For cells with diagrams

Karl Fischer anolyte solution - CFC free

• Karl Fischer soluzione anolita • Karl Fischer anolyte solution - sans CFC • Karl Fischer anolito solución - sin CFC • Karl Fischer Anolyt-Lösung - ohne FCKW

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group III



Danger
H226-H315-H319-H360D-H370-H335-H373-HA26
P210-P241-P280-P303+P361+P353-P305+P351+P338-P308+P313

Karl Fischer anolyte solution - CFC free > RS - ERBAqua - Coulometric titration

RS

Code	Size	Packaging	Notes
570141	500 ml	Glass bottle	

Pyridin and chlorocarbon free vessel solution. Suitable for nonpolar samples. For cells with diagrams

Karl Fischer anolyte solution, oven

• Karl Fischer soluzione anolita, con fornetto • Karl Fischer anolyte solution, pour four • Karl Fischer anolito solución para horno
• Karl-Fischer-Anolytlösung für den Ofen

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II



Danger
H225-H302-H315-H319-H360D-H370-H373-HA26
P210-P241-P280-P303+P361+P353-P305+P351+P338-P308+P313

Karl Fischer anolyte solution, oven > RS - ERBAqua - Coulometric titration

RS

Code	Size	Packaging	Notes
570151	500 ml	Glass bottle	

Fritless reagent for use in conjuction with Karl Fischer oven-Pydidin and chlorocarbon free - For cells with & without diagrams

Product specifications are subject to changes. Please visit our website for updates.

Karl Fischer anolyte solution for oils

• Karl Fischer soluzione anolita per olii • Karl Fischer anolyte solution pour les huiles • Karl Fischer anolito solución para los aceites
• Karl-Fischer-Anolytlösung für Öle

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H319-H370-H304
P210-P241-P280-P301+P310a-P303+P361+P353-P305+P351+P338

Karl Fischer anolyte solution for oils > RS - ERBAqua - Coulometric titration

RS

Code	Size	Packaging	Notes
570171	500 ml	Glass bottle	

Suitable for oils and petroleum products.For cells with diagrams

Karl Fischer catholyte solution

• Karl Fischer soluzione catolita • Karl Fischer catholyte solution • Karl Fischer catolito solución • Karl-Fischer-Katholyt-Lösung - Pyridinfrei

Classification transport
ONU: 1188
Transport Hazard class: 3
Packing group III



Danger
H226-H312-H315-H319-H360FD-H373-HA26
P210-P241-P280-P303+P361+P353-P305+P351+P338-P308+P313

Karl Fischer catholyte solution > RS - ERBAqua - Coulometric titration

RS

Code	Size	Packaging	Notes
570181	125 ml	Glass bottle	

Karl Fischer catholyte solution for aldehydes and ketones

• Karl Fischer soluzione catolita per aldeidi e chetoni • Karl Fischer catholyte solution pour aldehydes et cétones
• Karl Fischer catolito solución para aldehídos y cetonas • Karl-Fischer-Katholytlösung

Classification transport
ONU: 1188
Transport Hazard class: 3
Packing group III



Danger
H226-H312-H315-H318-H360FD-H373-HA26
P210-P280-P303+P361+P353-P305+P351+P338-P310a-P362+P364

Karl Fischer catholyte solution for aldehydes and ketones > RS - ERBAqua - Coulometric titration

RS

Code	Size	Packaging	Notes
570191	125 ml	Glass bottle	

For cells with diagrams

Karl Fischer anolyte solution for cells with and without diagrams

• Karl Fischer reagente per celle con e senza diaframma • Karl Fischer réactif pour la cellule avec et sans diaphragme
• Karl Fischer reactivo para la celda con y sin diafragma • Karl-Fischer-Reagenz für Zellen ohne Diaphragma

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H319-H360FD-H370-H373-HA26
P210-P241-P280-P303+P361+P353-P305+P351+P338-P308+P313

Karl Fischer anolyte solution for cells with and without diagrams > RS - ERBAqua - Coulometric titration

RS

Code	Size	Packaging	Notes
570131	500 ml	Glass bottle	

Can also be used with titrators that contain diaphragms or ceramic frits

Karl Fischer water standard 10.0 mg/g

• Karl Fischer standard 10.0 mg/g acqua • Karl Fischer standard eau 10.0 mg/g • Karl Fischer agua estándar 10.0 mg/g
• Karl Fischer Standardwasser 10.0 mg / g

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group III



Danger
H226-H315-H318-H351-H360D-H335-H336-H304-HA26
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Karl Fischer water standard 10.0 mg/g > RS - ERBAqua - For checking calibration of coulometric Karl Fischer equipment

RS

Code	Size	Packaging	Notes
570221	10 x 5 ml	Glass ampoule	

Gravimetric standard

Karl Fischer water standard 5.0 mg/ml

• Karl Fischer standard 5.0 mg/ml acqua • Karl Fischer standard eau 5.0 mg/ml • Karl Fischer agua estándar 5.0 mg/ml
• Karl Fischer Standardwasser 5.0 mg / ml

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group III



Danger
H226-H315-H318-H351-H360D-H335-H336-H304-HA26
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Karl Fischer water standard 5.0 mg/ml > RS - ERBAqua - For checking calibration of coulometric Karl Fischer equipment

RS

Code	Size	Packaging	Notes
570231	10 x 5 ml	Glass ampoule	

Volumetric standard

Karl Fischer water standard 1.0 mg/g

• Karl Fischer standard 1.0 mg/g acqua • Karl Fischer standard eau 1.0 mg/g • Karl Fischer agua estándar 1.0 mg/g
• Karl Fischer Standardwasser 1.0 mg / g

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group III



Danger
H226-H302-H315-H318-H335-H336-H373-H304
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Karl Fischer water standard 1.0 mg/g > RS - ERBAqua - For checking calibration of coulometric Karl Fischer equipment

RS

Code	Size	Packaging	Notes
570211	10 x 5 ml	Glass ampoule	

Gravimetric standard

Product specifications are subject to changes. Please visit our website for updates.

Karl Fischer water standard 0.10 mg/g

• Karl Fischer standard 0.10 mg/g acqua • Karl Fischer standard eau 0.10 mg/g • Karl Fischer agua estándar 0.10 mg/g
• Karl Fischer Standardwasser 0.10 mg / g

Classification transport
ONU: 1307
Transport Hazard class: 3
Packing group II

Danger
H225-H315-H319-H335-H373-H304
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Karl Fischer water standard 0.10 mg/g > RS - ERBAqua - For checking calibration of coulometric Karl Fischer equipment

RS

Code	Size	Packaging	Notes
570201	10 x 5 ml	Glass ampoule	

Gravimetric standard

Kerosene

• Cherosene • Kérosène • Queroseno • Kerosin

CAS: 64742-47-8

Classification transport
ONU: 1223
Transport Hazard class: 3
Packing group III

Danger
H304
P301+P310a-P331-P405-P501a

Kerosene > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Boiling point min. ..... ≥ 175 °C
Boiling point max..... ≤ 245 °C
Water (K.F.) ..... ≤ 0.01 %
Refractive index at 20°C..... 1.441 ÷ 1.451
Residue on ignition..... ≤ 0.001 %

Code	Size	Packaging	Notes
302575	5 l	Plastic tank	

Kieselguhr composed

• Farina fossile composta • Terre de silice composée • Harina fósil compuesta • Kieselgur besteht

Synonym:
Diatomaceous earth

Molecular Weight: 1495
CAS: 91053-39-3
EEC-N: 293-303-4

Warning
H373
P260-P314-P501a

Kieselguhr composed > RS - For thin layer chromatography according to Stahl

RS

Description ..... Grey powder
Identification ..... Positive

Code	Size	Packaging	Notes
449895	250 g	Plastic bottle	
449897	1 kg	Plastic bottle	

Kjeldahl antifoam

• Kjeldahl antischiuma • Kjeldahl Antimousse • Kjeldahl antiespumante • Kjeldahl entschäumer

Kjeldahl antifoam > RS - For Kjeldahl

RS

Description ..... White round flat tablets

Code	Size	Packaging	Notes
502811	1000 x 1 g	Metallic can	Composition: Sodium sulfate 0.97 g/Silicone antifoam 0.03 g

Kjeldahl catalyst according to Wieninger

• Kjeldahl catalizzatore sec. Wieninger • Catalyseur Kjeldahl selon Wieninger • Catalizador Kjeldahl según Wieninger • Kjeldahl-Katalysator nach Wieninger



Warning
H319-H411
P264-P280i-P305+P351+P338-P337+P313-P391-P501a

Kjeldahl catalyst according to Wieninger > RS - For Kjeldahl

RS

Description ..... Grey round flat tablets

Code	Size	Packaging	Notes
502821	1000 x 5 g	Metallic can	Composition: Sodium sulfate 4.88 g/ Copper sulfate 0.07 g/ Selenium 0.05 g

Kjeldahl catalyst for water analysis

• Kjedahl catalizzatore per analisi dell'acqua • Catalyseur Kjedahl pour analyse de l'eau • Catalizador Kjeldahl para análisis del agua
• Kjedahl-Katalysator für die Wasseranalyse

HEU210

Kjeldahl catalyst for water analysis > RS - For Kjeldahl

RS

Description .....Dark grey round flat tablets

Code	Size	Packaging	Notes
502121	1000 x 5 g	Metallic can	Composition: Potassium sulfate 5.0 g/Selenium 5 mg
502122	1000 x 5 g	Metallic can	Composition: Potassium sulfate 5.0 g/Selenium 50 mg

Kjeldahl catalyst without selenium and titanium

• Kjeldahl catalizzatore • Catalyseur Kjeldahl exempt de sélénium et titane • Catalizador Kjeldahl sin selenio y potassio
• Kjeldahl-Katalysator Selen- und Titanfrei

Classification transport
ONU: 3077
Transport Hazard class: 9
Packing group III



Warning
H319-H410
P264-P280i-P305+P351+P338-P337+P313-P391-P501a

Kjeldahl catalyst without selenium and titanium > RS - For Kjeldahl

RS

Description ..... Blue speckeld round flat tablets

Code	Size	Packaging	Notes
502791	1000 x 3.9 g	Metallic can	Composition: Potassium sulfate 3.50 g/Copper sulfate 0.40 g
502792	1000 x 5 g	Metallic can	Composition: Potassium sulfate 5.0 g/Copper sulfate 0.50 g

Kjeldahl selenium catalyst

• Kjeldahl catalizzatore al selenio • Catalyseur Kjeldahl au sélénium • Catalizador Kjeldahl al selenio • Kjeldahl-Selen-Katalysator

Classification transport
ONU: 3077
Transport Hazard class: 9
Packing group III



Warning
H319-H410
P264-P280i-P305+P351+P338-P337+P313-P391-P501a

Kjeldahl selenium catalyst > RS - For Kjeldahl

RS

Description .....Dark grey round flat tablets

Code	Size	Packaging	Notes
502120	1000 x 5 g	Metallic can	Composition: Potassium sulfate 4.63 g/Copper sulfate 0.28 g/ Selenium 0.09g

Product specifications are subject to changes. Please visit our website for updates.

Kjeldahl titanium catalyst

• Kjeldahl catalizzatore al titanio • Catalyseur Kjeldahl au titane • Catalizador Keldahl al titanio • Kjeldahl-Titankatalysator

Classification transport
ONU: 3077
Transport Hazard class: 9
Packing group III



Warning
H319-H410
P264-P280i-P305+P351+P338-P337+P313-P391-P501a

Kjeldahl titanium catalyst > RS - For Kjeldahl

RS

Description ..... Blue speckled round flat tablets

Code	Size	Packaging	Notes
502123	1000 x 3.5 g	Metallic can	Composition: Potassium sulfate 3.5 g/Copper sulfate 0.105 g/ Titane dioxide 0.105 g
502802	500 x 5 g	Metallic can	Composition: Potassium sulfate 5.00 g/Copper sulfate 0.15 g/ Titane dioxide 0.15 g

Kovac reagent

• Kovac reattivo per indolo • Réactif pour l'indole selon Kovac • Kovac reactivo • Kovac-Reagenz

Classification transport
ONU: 2924
Transport Hazard class: 3
Packing group III



Warning
H226-H290-H332-H315-H319-H335
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Kovac reagent > RS - Reagent for indole

RS

Description .....Clear yellow green liquid Identification ..... Positive

Code	Size	Packaging	Notes
435922	100 ml	Glass bottle	

L(+)-Lactic acid

• Acido L(+)-lattico • Acide L(+)-lactique • Acido L(+)-láctico • L-(+)-Milchsäure

Synonym:
Sarcolactic acid

$CH_3CHOHCOOH$
Molecular Weight: 90,08
CAS: 79-33-4
EEC-N: 201-196-2

Danger
H315-H318
P264-P280a-P305+P351+P338-P310a-P362+P364-P332+P313

L(+)-Lactic acid > ERBApharm - According to pharmacopoeia: BP-DAB-FU-Ph.Eur.-Ph.Franc.

ERBApharm

Description .....Syrupy colourless to slightly yellow liq.
Identification .....Positive
Appearance .....Conform Ph.Eur.
Reducing substances .....Conform Ph.Eur.
Ether ins. substances .....Conform Ph.Eur.
Citric, oxalic, phosphoric acids .....Conform Ph.Eur
Density at 20° C .....1.20 ÷ 1.21
Sulphated ash..... ≤ 0.1 %
Heavy metals (Pb).....≤10 ppm
Sulphate.....≤200 ppm
Ca.....≤200 ppm
Assay (acidimetric) .....88.0 ÷ 92.0 %
Origin (BSE/TSE)..... Synthesis
Residual solvents (Current ICH).....Conform

Code	Size	Packaging	Notes
304652	1 l	Glass bottle	
304651	2.5 l	Glass bottle	
304653	25 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

L-Lactic acid calcium salt ▶ Calcium lactate

Lactophenol blue solution

• Blu lattofenolo soluzione • Bleu de lactophénol solution • Azul lactofenol solución • Lactophenolblau-Lösung

Classification transport
ONU: 2927
Transport Hazard class: 6.1
Packing group II

Danger
H302-H331-H314-H341-H373
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P308+P313-P403+P233

Lactophenol blue solution > RS - For microscopy

RS

Description .....Dark blue liquid
Identification ..... Positive

Code	Size	Packaging	Notes
428901	100 ml	Glass bottle	In Vitro Diagnostic Medical Device

Dye for bacteriology. Contains phenol and lactic acid

Lactose monohydrate

• Lattosio monoidrato • Lactose monohydrate • Lactosa monohidrato • Lactose Monohydrat

Synonym:
4-O-β-D-Galactopyranosyl-D-glucose

$C_{12}H_{22}O_{11}.H_2O$
Molecular Weight: 360,32
CAS: 10039-26-6
EEC-N: 200-559-2

Lactose monohydrate > RPE - For analysis - ACS

RPE

Description .....White crystalline powder
Identification ..... Positive
Dextrose.....Conform
Sucrose.....Conform
Water (K.F.)..... 4.0 ÷ 6.0 %
Water-insoluble matter .....≤50 ppm
Residue on ignition.....≤300 ppm
Heavy metals (Pb).....≤5 ppm
Fe .....≤5 ppm

Code	Size	Packaging	Notes
457551	100 g	Plastic bottle	
457552	250 g	Plastic bottle	
457557	1 kg	Plastic bottle	
457553	25 kg	Plastic bucket	

Product specifications are subject to changes. Please visit our website for updates.

Lactose monohydrate > ERBApharm - According to pharmacopoeia: DAB-BP-FU-NF-Ph.Eur.

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Acidity or alkalinity.....Conform Ph.Eur.
Absorbance .....Conform Ph.Eur.
Specific optical rotation on dry ..... +54.4 ÷ +55.9 °
Water (K.F.) .....4.5 ÷ 5.5 %
Loss on drying ..... ≤ 0.5 %
Sulphated ash..... ≤ 0.1 %
Heavy metals (Pb).....≤ 5 ppm
TAMC ..... ≤ 100 CFU/g
TYMC ..... ≤ 50 CFU/g
Escherichia coli ..... Absent Ph.Eur

Code	Size	Packaging	Notes
348707	1 kg	Plastic bottle	
348708	5 kg	Plastic tank	
348702	10 kg	Carton box	
348703	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Lanolin anhydrous

• Grasso di lana anidro • Lanoline anhydre • Grasa de lana anhidra • Lanolin Wasserfreies

Synonym:
Wool fat

CAS: 8006-54-0
EEC-N: 232-348-6

Lanolin anhydrous > ERBApharm - According to pharmacopoeia: FU-Ph.Eur.

ERBApharm

Description .....Soft yellow mass
Melt prod. description.....Clear liquid
Identification ..... Positive
Acids,alkal.water-solub.....Conform Ph.Eur.
Oxydable hydrosol.matt.....Conform Ph.Eur.
Water absorption capac. ....Conform Ph.Eur.
Acid value..... ≤1.0
Peroxide value..... ≤20
Saponification value.....90 ÷ 105
Loss on drying ..... ≤0.5 %
Sulphated ash..... ≤0.15 %
Chloride.....≤150 ppm
Paraffins..... ≤1.0 %
Residui pesticidi .....≤ 1 ppm

Code	Size	Packaging	Notes
347357	1 kg	Metallic can	
347359	5 kg	Plastic tank	
347355	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Lanthanum standard solution

• Lantanio standard soluzione • Lanthane solution standard • Lantano, solución patrón • Lanthan-Standardlösung

Classification transport
ONU: 3267
Transport Hazard class: 8
Packing group III

Lanthanum standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505692	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505695	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505693	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Lanthanum standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503681	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503683	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503685	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503687	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Lanthanum standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507744	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507509	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Lanthanum chloride 25 g/l solution

• Lantanio cloruro 25 g/l soluzione • Lanthane chlorure 25g/l • Lantano cloruro solución 25 g/l
• Lanthanchlorid 25 g/l

Synonym:
• *Lanthanum(III) chloride*
• *Lanthanum trichloride*

$LaCl_3$
Molecular Weight: 245,27
CAS: 10099-58-8

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Lanthanum chloride 25 g/l solution > RS - Ionisation standard solution for AAS

RS

Code	Size	Packaging	Notes
504537	500 ml	Plastic bottle	Matrix: Hydrochloric acid

Lanthanum nitrate hexahydrate

• Lantanio nitrato esaidrato • Lanthane nitrate hexahydraté • Lantano nitrato hexahidratado
• Lanthannitrat-Hexahydrat

Synonym:
Nitric acid, lanthanum (III) salt, hexahydrate

$La(NO_3)_3.6H_2O$
Molecular Weight: 433,02
CAS: 10277-43-7
EEC-N: 233-238-0

Classification transport
ONU: 1477
Transport Hazard class: 5.1
Packing group II

Danger
H272-H315-H319-H335
P210-P261-P280-P304+P340-P305+P351+P338-P403+P233

Lanthanum nitrate hexahydrate > RPE - For analysis

RPE

Description ..... White crystalline powder Identification ..... Positive Fe ..... ≤ 10 ppm Assay ..... ≥ 98.0 %

Code	Size	Packaging	Notes
457502	25 g	Glass bottle	
457506	250 g	Plastic bottle	

Lanthanum nitrate 0.1 mol/l

• Lantanio nitrato 0.1 mol/l • Lanthane nitrate 0.1 mol/l • Lantano nitrato 0.1 mol/l • Lanthannitrat 0.1 mol/l

Lanthanum nitrate 0.1 mol/l > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613010100	1 l	Plastic bottle	Ref Ph.Eur 3010100

Lanthanum nitrate solution 50 g/l

• Lantanio nitrato soluzione 50 g/l • Lanthane nitrate 50 g/l • Lantano nitrato solución 50 g/l • Lanthannitratlösung 50 g/l

HEU210

Lanthanum nitrate solution 50 g/l > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611048001	1 l	Plastic bottle	Ref Ph.Eur 1048001

Product specifications are subject to changes. Please visit our website for updates.

Lanthanum oxide
• Lantanio ossido • Lanthane oxyde • Lantano oxido • Lanthanoxid

La_2O_3
Molecular Weight: 325,81
CAS: 1312-81-8
EEC-N: 215-200-5

Lanthanum oxide > RPE - For analysis

RPE

Description ..... White or pink powder
Identification ..... Positive
Loss on ignition ..... ≤ 3.0 %
Fe2O3 ..... ≤ 250 ppm
Assay (complexometric) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
457511	100 g	Glass bottle	

Lauryltrimethylammonium bromide ▶ Dodecyltrimethylammonium bromide

Lead
• Piombo • Plomb • Plomo • Bleiblech

Pb
Molecular Weight: 207,2
CAS: 7439-92-1
EEC-N: 231-100-4

Danger
H360FD-H362-H372-HA26
P260h-P263-P264-P280-P308+P313-P501a

Lead > RPE - For analysis

RPE

Description ..... Grey foil
Identification ..... Positive
Assay ..... 99 ÷ 100 %

Code	Size	Packaging	Notes
468866	500 g	Box	

~ 0,7 mm thickness

Lead standard solution
• Piombo standard soluzione • Plomb solution standard • Plomo, solución patrón • Blei-Standardlösung

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-
P362+P364-P337+P313

Lead standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615001700	100 ml	Plastic bottle	A 0.1 % solution Ref Ph.Eur 5001700
615001702	100 ml	Plastic bottle	A 10 ppm solution: to dilute according to Ref Ph.Eur 5001702
615001703	100 ml	Plastic bottle	A 2 ppm solution: to dilute according to Ref Ph.Eur 5001703
615001704	100 ml	Plastic bottle	A 1 ppm solution: to dilute according to Ref Ph.Eur 5001704
615001705	100 ml	Plastic bottle	A 0.1 ppm solution: to dilute according to Ref Ph.Eur 5001705
615001706	100 ml	Plastic bottle	A 10 ppm solution R1: to dilute according to Ref Ph.Eur 5001706
615001709	100 ml	Plastic bottle	A 100 ppm solution: to dilute according to Ref Ph.Eur 5001701
615001701	500 ml	Plastic bottle	A 0.1 % solution Ref Ph.Eur 5001700

Lead standard solution > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
504901	100 ml	Plastic bottle	conc. 100 ppm
504902	500 ml	Plastic bottle	conc. 100 ppm

Lead standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505767	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505768	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505769	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Lead standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503801	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503803	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503805	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503807	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Lead standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
507752	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507490	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
E497595	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Nitric acid
E497591	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Lead standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
468791		Plastic ampoule	conc. 1.000 ppm Matrix: Nitric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Lead standard solution > RS - Quality control standard solution for AAS (graphite furnace)

RS

Code	Size	Packaging	Notes
504364	50 ml	Plastic bottle	conc. 10 ppb - Matrix: 1% Nitric acid

Product specifications are subject to changes. Please visit our website for updates.

Lead (II) acetate basic

• Piombo acetato basico • Plomb (II) acétate basique • Plomo (II) acetato básico • Blei (II) acetat basisch

Synonym:
Lead subacetate

$C_4H_6O_4Pb.2Pb(OH)_2$
Molecular Weight: 566,48
CAS: 1335-32-6
EEC-N: 215-630-3

Classification transport
ONU: 1616
Transport Hazard class: 6.1
Packing group III



Danger
H351-H360Df-H373-H410-HA26
P260-P280-P308+P313-P314-P405-P501a

Lead (II) acetate basic > RPE - For analysis - ACS

RPE

Description ..... White powder
Identification ..... Positive
Loss on drying ..... ≤1.5 %
Water-insoluble matter ..... ≤1.0 %
Insol.in dil.acetic ac. ..... ≤200 ppm
Chloride ..... ≤30 ppm
Nitrate-nitrite (NO3) ..... ≤30 ppm
Ca ..... ≤ 100 ppm
Cu ..... ≤20 ppm
Fe ..... ≤20 ppm
K ..... ≤ 200 ppm
Na ..... ≤ 500 ppm
Assay (alkalimetric) ..... ≥33.0 % PbO

Code	Size	Packaging	Notes
468984	100 g	Plastic bottle	
468985	250 g	Plastic bottle	
468987	1 kg	Plastic bottle	

Lead (II) acetate basic solution

• Piombo acetato basico soluzione • Plomb (II) acétate basique soluiton • Plomo (II) acetato básico solución • Blei (II) acetat basische Lösung

Synonym:
Lead subacetate

Classification transport
ONU: 1616
Transport Hazard class: 6.1
Packing group III



Danger
H351-H360Df-H373-H410-HA26
P260-P280-P308+P313-P314-P405-P501a

Lead (II) acetate basic solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611048400	100 ml	Glass bottle	Ref Ph.Eur 1048400

Lead (II) acetate cotton

• Cotone piombo (II) acetato • Coton à l'acétate de plomb (II) • Plomo (II) acetato algodón • Blei (II) acetat Baumwolle

$Pb(H_3CCOO)_2$
Molecular Weight: 325,3
CAS: 301-04-2



Danger
H360-H411-HA26-HEU201
P273-P280-P308+P313-P391-P405-P501a

Lead (II) acetate cotton > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611048101	10 g	Glass bottle	Ref Ph.Eur 1048101

Lead (II) acetate cotton > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
617000301	10 g	Plastic bottle	Lead acetate cotton

Lead (II) acetate paper

• Cartina piombo (II) acetato • Papier (II) à l'acétate de plomb • Plomo (II) acetato papel • Blei (II) acetatpapier



Danger
H360-H411-HA26-HEU201
P273-P280-P308+P313-P391-P405-P501a

Lead (II) acetate paper > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611048102	50 pc	Tube	Ref Ph.Eur 1048102

Lead (II) acetate solution 95 g/l

• Piombo acetato soluzione 95 g/l • Plomb (II) acétate solution 95g/l • Plomo (II) acetato solución 95 g/l
• Blei (II) acetatlösung 95 g / l

Synonym:
Lead subacetate



Danger
H360Df-H373-H411-HA26-HEU201
P260-P280-P308+P313-P314-P405-P501a

Lead (II) acetate solution 95 g/l > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611048103	100 ml	Plastic bottle	Ref Ph.Eur 1048103

Lead (II) acetate trihydrate

• Piombo di-acetato triidrato • Plomb (II) acétate trihydraté • Plomo (II) acetato 3-hidratado • Blei (II) acetat Trihydrat

$C_4H_6O_4Pb.3H_2O$
Molecular Weight: 379,33
CAS: 6080-56-4

Classification transport
ONU: 1616
Transport Hazard class: 6.1
Packing group III



Danger
H302-H332-H360Df-H373-H410-HA26-HEU201
P260-P264-P271-P280-P304+P340-P308+P313

Lead (II) acetate trihydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White crystals
Identification ..... Positive
Diluted acetic acid insoluble matter ... ≤100 ppm
Chloride ..... ≤5 ppm
Nitrate-nitrite (NO3) ..... ≤50 ppm
Cu ..... ≤20 ppm
Fe ..... ≤10 ppm
Assay (complexometric) .... 99.0 ÷ 103.0 %
Ca ..... ≤50 ppm
K ..... ≤50 ppm
Na ..... ≤100 ppm

Code	Size	Packaging	Notes
468934	100 g	Plastic bottle	
468935	250 g	Plastic bottle	
468937	1 kg	Plastic bottle	
468932	25 kg	Drum	

Lead (II) acetate trihydrate > RE - Pure

RE

Description ..... White crystals
Identification ..... Positive
Fe ..... ≤ 20 ppm
Zn ..... ≤ 50 ppm
Assay (complexometric) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
357253	25 kg	Plastic bucket	

Lead (II) nitrate

• Piombo nitrato • Plomb (II) nitrate • Plomo (II) nitrato • Blei (II) nitrat

Synonym:
Lead dinitrate

$Pb(NO_3)_2$
Molecular Weight: 331,21
CAS: 10099-74-8

Lead (II) nitrate > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
617000321	100 ml	Plastic bottle	Stock Solution TS

Lead (II) nitrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White crystals
Identification ..... Positive
Dil. HNO3-insol. matter ..... ≤ 50 ppm
Chloride ..... ≤ 10 ppm
Ca ..... ≤ 50 ppm
Cu ..... ≤ 20 ppm
Fe ..... ≤ 10 ppm
Assay (complexometric) ..... ≥ 99.0 %
K ..... ≤ 50 ppm
Na ..... ≤ 0.02 %

Code	Size	Packaging	Notes
469355	100 g	Plastic bottle	
469356	500 g	Plastic bottle	
469357	1 kg	Plastic bottle	
469353	25 kg	Plastic bucket	

Product specifications are subject to changes. Please visit our website for updates.

Lead (II) nitrate > RE - Pure

RE

Description ..... White crystals
Identification ..... Positive
Dil. HNO3-insol. matter ..... ≤ 0.005 %
Chloride ..... ≤ 50 ppm
pH solution 7.5% ..... 3 ÷ 4
Cu ..... ≤ 20 ppm
Fe ..... ≤ 10 ppm
Assay (complexometric) ..... ≥ 98.0 %

Code	Size	Packaging	Notes
358007	1 kg	Plastic bottle	
358008	5 kg	Plastic tank	

Lead (II) nitrate 0.1 mol/l

• Piombo nitrato 0.1 mol/l • Plomb (II) nitrate 0.1 mol/l • Plomo (II) nitrato 0.1 mol/l • Blei (II) nitrat 0.1 mol/l

Synonym:
Lead dinitrate



Danger
H360Df-H373-H411-HA26-HEU201
P260-P280-P308+P313-P314-P405-P501a

Lead (II) nitrate 0.1 mol/l > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613003101	500 ml	Plastic bottle	Ref Ph.Eur 3003100
613003100	1 l	Plastic bottle	Ref Ph.Eur 3003100

Lead (II) nitrate 0.05 mol/l

• Piombo nitrato 0.05 mol/l • Plomb (II) nitrate 0.05 mol/l • Plomo (II) nitrato 0.05 mol/l
• Blei (II) nitrat 0.05 mol/l

Synonym:
Lead dinitrate



Danger
H360Df-H373-H412-HA26-HEU201
P260-P280-P308+P313-P314-P405-P501a

Lead (II) nitrate 0.05 mol/l > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613009700	100 ml	Plastic bottle	Ref Ph.Eur 3009700

Lead (II) nitrate solution 33 g/l

• Piombo nitrato soluzione 33 g/l • Plomb (II) nitrate solution 33g/l • Plomo (II) nitrato solución 33 g/l
• Blei (II) nitratlösung 33 g/l

Synonym:
Lead dinitrate

Classification transport
ONU: 1935
Transport Hazard class: 6.1
Packing group III



Danger
H360Df-H373-H411-HA26-HEU201
P260-P280-P308+P313-P314-P405-P501a

Lead (II) nitrate solution 33 g/l > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611048301	1 l	Plastic bottle	Ref Ph.Eur 1048301

Lead (II) oxide

• Piombo ossido • Plomb (II) oxyde • Plomo (II) óxido • Blei (II) oxid

PbO
Molecular Weight: 223,2
CAS: 1317-36-8
EEC-N: 215-267-0



Danger
H302-H332-H360Df-H373-H410-HA26-HEU201
P260-P264-P271-P280-P304+P340-P308+P313

Lead (II) oxide > RPE - For analysis

RPE

Description	Polvere gialla o arancio	Insoluble in Acetic ac.	≤ 200 ppm	Ca	≤ 50 ppm	K	≤ 50 ppm
Identification	Positive	Nitrate	≤ 100 ppm	Cu	≤ 50 ppm	Na	≤ 0.02 %
Chloride	≤ 20 ppm	Ag	≤ 5 ppm	Fe	≤ 20 ppm	Assay (complexometric)	≥ 99.0 %

Code	Size	Packaging	Notes
469404	100 g	Glass bottle	

Lead (II) oxide > RE - Pure

RE

Description	Yellow powder	Identification	Positive	Assay	≥ 99.8 %

Code	Size	Packaging	Notes
358257	1 kg	Plastic bottle	
358259	5 kg	Plastic tank	
358252	25 kg	Plastic bucket	

Lead (II) sulfate

• Piombo solfato • Plomb (II) sulfate • Plomo (II) sulfato • Blei (II) sulfat

Synonym:
Anglesite

$PbSO_4$
Molecular Weight: 303,26
CAS: 7446-14-2
EEC-N: 231-198-9



Danger
H302-H332-H360Df-H373-H410-HA26-HEU201
P260-P264-P271-P280-P304+P340-P308+P313

Lead (II) sulfate > RPE - For analysis

RPE

Description	White powder	Chloride	≤20 ppm	Fe	≤20 ppm	Assay (complexometric)	≥99.0 %
Identification	Positive	Insol.in Ammonium acet.	≤500 ppm	K	≤500 ppm		
Loss on ignition	≤0.5 %	Nitrate	Conform	Na	≤0.1 %		

Code	Size	Packaging	Notes
469505	250 g	Plastic bottle	
469506	1 kg	Plastic bottle	

Lead (IV) oxide

• Piombo biossido • Plomb (IV) oxyde • Plomo (IV) óxido • Blei (IV) oxid

Synonym:
Lead peroxide

PbO_2
Molecular Weight: 239,2
CAS: 1309-60-0
EEC-N: 215-174-5

Classification transport
ONU: 1872
Transport Hazard class: 5.1
Packing group III



Danger
H302-H332-H360Df-H373-H410-HA26-HEU201
P260-P264-P271-P280-P304+P340-P308+P313

Lead (IV) oxide > RPE - For analysis

RPE

Description	Blackish powder	Sulphate	≤ 0.1 %	Loss on drying	≤ 0.5 %
Identification	Positive	HNO3-insoluble matter	≤ 0.05 %	Mn	≤ 20 ppm
Chloride	≤ 200 ppm	Substances not ppt. H2S	≤ 0.05 %	Assay (oxidimetric)	≥ 97.0 %

Code	Size	Packaging	Notes
469055	250 g	Plastic bottle	
469057	1 kg	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Lecithin of soya

• Lecitina di soia • Lécithine de soja • Lecitina de soja • Sojalecithin

Synonym:
• *L-α-Phosphatidylcholine*
• *1,2-Diacyl-sn-glycero-3-phosphocholine*

Molecular Weight: 750
CAS: 8002-43-5
EEC-N: 232-307-2

Lecithin of soya > RE - Pure

RE

Description ..... Powder	Water ..... < 1 %	Peroxide value..... < 5 meq O2/Kg	Assay (phospholipids as acetone insolubl.> 97 %
Identification ..... Positive	Acidity index ..... < 35 mg KOH/g	Insolubles in toluene..... < 0.3 %	Sulphated ash..... ≤ 12 %

Code	Size	Packaging	Notes
348754	1 kg	Plastic bottle	

L(+)-Leucine

• L(+)Leucina • L(+)-Leucine • L(+)-Leucina • L (+) - Leucin

Synonym:
(S)-2-Amino-4-methylpentanoic acid

$(CH_3)_2CHCH_2CHNH_2COOH$
Molecular Weight: 131,18
CAS: 61-90-5
EEC-N: 200-522-0

L(+)-Leucine > RPE - For analysis

RPE

Description ..... White crystals	Ammonium ..... ≤50 ppm	Residue on ignition..... ≤500 ppm	Fe ..... ≤10 ppm
Identification ..... Positive	Chloride..... ≤250 ppm	Tyrosine..... ≤200 ppm	Assay (non-aqueos medium)..... ≥99 %
Specific optical rotation... +15.1 ÷ +16.1 °	Total phosphorus..... ≤10 ppm	Tryptophan..... ≤300 ppm	
Loss on drying ..... ≤0.1 %	Heavy metals (Pb)..... ≤10 ppm	Total sulphur ..... ≤250 ppm	

Code	Size	Packaging	Notes
457928	5 g	Glass bottle	

Light green

• Verde luce SF • Vert lumière SF • Verde claro SF • Hellgrün

Synonym:
Acid Green 5

$C_{37}H_{34}N_2Na_2O_9S_3$
Molecular Weight: 792,86
CAS: 5141-20-8
EEC-N: 225-906-5

Warning
H312-H332
P261-P271-P280h-P304+P340-P312a-P501a

Light green > RS - For microscopy - C.I. 42095

RS

Description ..... Dark violet powder	Identification ..... Positive	Water solubility..... Conform	Assay ..... ≥ 95 %

Code	Size	Packaging	Notes
491371	10 g	Glass bottle	
491372	25 g	Glass bottle	

Dye for botanical-bacteriology-cytology

Lithium standard solution

• Litio standard soluzione • Lithium solution standard • Litio, solución patrón • Lithium-Standardlösung



Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Lithium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505702	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505705	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505703	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Lithium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503691	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503693	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503695	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid
503697	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Lithium standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
507745	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507486	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
E497525	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
E497521	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Lithium standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
458211		Plastic ampoule	conc. 1.000 ppm Matrix: Nitric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Lithium standard solution > RS - Standard solution for ion chromatography

RS

Code	Size	Packaging	Notes
503281	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503283	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Product specifications are subject to changes. Please visit our website for updates.

Lithium carbonate

• Litio carbonato • Lithium carbonate • Litio carbonato • Lithiumcarbonat

Synonym:
• *Carbolithium*
• *Carbonic acid lithium salt*

Li_2CO_3
Molecular Weight: 73,89
CAS: 554-13-2
EEC-N: 209-062-5

Warning
H302-H319
P264-P280i-P301+P312a-P305+P351+P338-P337+P313-P501a

Lithium carbonate > RPE - For analysis - ACS

RPE

Description ..... White powder
Identification ..... Positive
Chloride ..... ≤ 50 ppm
Diluted HCl-ins. matter ..... ≤ 100 ppm
Heavy metals (Pb) ..... ≤ 20 ppm
Nitrate ..... ≤ 5 ppm
Total sulphur ..... ≤ 0.2 %
Ca ..... ≤ 100 ppm
Fe ..... ≤ 20 ppm
K ..... ≤ 100 ppm
Assay (alkalimetric) ..... ≥ 99.0 %
Na ..... ≤ 0.1 %

Code	Size	Packaging	Notes
458204	100 g	Plastic bottle	
458207	1 kg	Plastic bottle	

Lithium carbonate > RE - Pure

RE

Description ..... White powder
Identification ..... Positive
Chloride ..... ≤300 ppm
Diluted HCl-ins. matter ..... ≤500 ppm
Sulphate ..... ≤0.2 %
Fe ..... ≤30 ppm
Assay (non-aqueos medium) ..... ≥98 %

Code	Size	Packaging	Notes
348955	250 g	Plastic bottle	
348957	1 kg	Plastic bottle	

Lithium chloride

• Litio cloruro • Lithium chlorure • Litio cloruro • Lithiumchlorid

LiCl
Molecular Weight: 42,39
CAS: 7447-41-8
EEC-N: 231-212-3

Warning
H302-H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Lithium chloride > RPE - For analysis - ACS

RPE

Description ..... White crystalline powder
Identification ..... Positive
Alcalinity ..... ≤ 0.008 meq/g
Loss on drying ..... ≤ 1.0 %
Water-insoluble matter ..... ≤ 100 ppm
Heavy metals (Pb) ..... ≤ 20 ppm
Sulphate ..... ≤ 100 ppm
Ba ..... ≤ 30 ppm
Ca ..... ≤ 100 ppm
Fe ..... ≤ 10 ppm
K ..... ≤ 100 ppm
Assay (argentimetric) ..... ≥ 99 %
Na ..... ≤ 0.2 %

Code	Size	Packaging	Notes
458254	100 g	Glass bottle	
458256	500 g	Plastic bottle	
458257	1 kg	Plastic bottle	

Lithium chloride > RE - Pure

RE

Description ..... White crystalline powder
Identification ..... Positive
Umidità (H2O) ..... ≤ 0.8 %
Sulphate ..... ≤ 500 ppm
Fe ..... ≤ 20 ppm
Assay (argentimetric) ..... ≥ 99 %

Code	Size	Packaging	Notes
458271	250 g	Plastic bottle	
458272	1 kg	Plastic bottle	
458275	5 kg	Plastic jar	
458273	25 kg	Plastic bucket	

Lithium hydride
• Litio idruro • Lithium hydrure • Litio hidruro • Lithiumhydrid

LiH
Molecular Weight: 7,95
CAS: 7580-67-8
EEC-N: 231-484-3

Classification transport
ONU: 1414
Transport Hazard class: 4.3
Packing group I



Danger
H260-H314
P223-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Lithium hydride > RE - Pure — RE

Description ..... Polvere bianco grigiast
Identification ..... Positive
Assay (gas volumetric) ..... ≥95 %

Code	Size	Packaging	Notes
458303	50 g	Glass bottle	

Lithium hydroxide anhydrous
• Litio idrossido anidro • Lithium hydroxyde anhydre • Litio hidróxido anhidro • Lithiumhydroxid wasserfrei

LiOH
Molecular Weight: 23,95
CAS: 1310-65-2
EEC-N: 215-183-4

Classification transport
ONU: 2680
Transport Hazard class: 8
Packing group II



Danger
H331-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Lithium hydroxide anhydrous > RPE - For analysis — RPE

Description ..... White powder
Identification ..... Positive
Fe ..... ≤ 0.0001 %
Chloride ..... ≤ 0.04 %
Sulfate ..... ≤ 0.03 %
Assay (LiOH) ..... ≥ 99.0 %
Na ..... ≤ 0.02 %
Ca ..... ≤ 0.005 %
Li2CO3 ..... ≤ 0.5 %
Insolubles ..... ≤ 0.01 %

Code	Size	Packaging	Notes
458281	250 g	Glass bottle	
458282	1 kg	Plastic bottle	

Lithium hydroxide monohydrate
• Litio idrossido monoidrato • Lithium hydroxyde monohydraté • Litio hidróxido monohidrato • Lithiumhydroxidmonohydrat

$LiOH.H_2O$
Molecular Weight: 41,96
CAS: 1310-66-3
EEC-N: 215-183-4

Classification transport
ONU: 2680
Transport Hazard class: 8
Packing group II



Danger
H302-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Lithium hydroxide monohydrate > RPE - For analysis — RPE

Description ..... White powder
Identification ..... Positive
Carbonate ..... ≤ 0.70 %
Chloride ..... ≤ 100 ppm
Sulphate ..... ≤ 300 ppm
Ca ..... ≤ 250 ppm
Fe ..... ≤ 15 ppm
Assay (acidimetric) ..... ≥98 %
Insoluble in water ..... ≤ 0.01 %
K + Na ..... ≤ 0.05 %
Insoluble in HCl ..... ≤ 0.01 %

Code	Size	Packaging	Notes
458292	1 kg	Plastic bottle	

Lithium methoxide 0.1 mol/l (0.1N)
• Litio metilato 0.1 mol/l (0.1N) • Lithium méthylate 0.1 mol/l (0.1N) • Litio metilato 0.1 mol/l (0.1N)
• Lithiummethylat 0.1 mol/l (0.1 N)

Synonym:
Lithium methylate

CH_3LiO
Molecular Weight: 38,02
CAS: 865-34-9

Classification transport
ONU: 1992
Transport Hazard class: 3
Packing group II



Danger
H225-H302-H314-H361d-H370-H336-H373
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Lithium methoxide 0.1 mol/l (0.1N) > RPE - For analysis — RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Titration factor ..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
E458321	500 ml	Glass bottle	

Toluene solution ready-to-use

Product specifications are subject to changes. Please visit our website for updates.

Lithium nitrate

• Litio nitrato • Lithium nitrate • Litio nitrato • Lithiumnitrat

$LiNO_3$
Molecular Weight: 68,95
CAS: 7790-69-4
EEC-N: 232-218-9

Classification transport
ONU: 2722
Transport Hazard class: 5.1
Packing group III

Danger
H272
P210-P220-P280-P370+P378a-P501a

Lithium nitrate > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Water ..... ≤ 1 %
Alcalinity ..... ≤ 0.05 %
Chloride ..... ≤ 100 ppm
Sulphate ..... ≤ 500 ppm
Fe2O3 ..... ≤ 50 ppm
Assay ..... ≥ 99.0 %

Code	Size	Packaging	Notes
458355	250 g	Plastic bottle	
458356	1 kg	Plastic bottle	

Lithium sulfate monohydrate

• Litio solfato monoidrato • Lithium sulfate monohydraté • Litio sulfato monohidrato • Lithiumsulfat-Monohydrat

$Li_2SO_4.H_2O$
Molecular Weight: 127,95
CAS: 10102-25-7
EEC-N: 233-802-4

Warning
H302
P264-P270-P301+P312a-P330-P501a

Lithium sulfate monohydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White crystalline powder
Identification ..... Positive
Loss on drying 150° C ..... 13.0 ÷ 15.0
Water-insoluble matter ..... ≤100 ppm
Chloride ..... ≤20 ppm
Nitrate ..... ≤10 ppm
Heavy metals (Pb) ..... ≤10 ppm
Fe ..... ≤10 ppm
K ..... ≤500 ppm
Na ..... ≤500 ppm
Assay (acidimetric) ..... ≥99.0 % s.s.

Code	Size	Packaging	Notes
458404	100 g	Plastic bottle	
458405	1 kg	Plastic bottle	

Lithium tetraborate anhydrous

• Litio tetraborato anidro • Lithium tétraborate anhydre • Litio tetraborato anhidro • Lithiumtetraborat wasserfrei

$Li_2B_4O_7$
Molecular Weight: 169,12
CAS: 12007-60-2
EEC-N: 234-514-3

Danger
H302-H318-H361d
P264-P280-P301+P312a-P305+P351+P338-P310a-P308+P313

Lithium tetraborate anhydrous > RE - Pure

RE

Description ..... White powder
Identification ..... Positive
Assay ..... ≥ 98 %

Code	Size	Packaging	Notes
458163	1 kg	Plastic bottle	
458164	5 kg	Plastic jar	

Litmus

• Tornasole • Tournesol • Tornasol • Lackmus

Molecular Weight: 3300
CAS: 1393-92-6
EEC-N: 215-739-6

Litmus > RPE - For analysis

RPE

Description ..... Dark blue granules
Identification ..... Positive
Colour change ..... rosso - blu
pH range ..... 4.8 ÷ 8.3

Code	Size	Packaging	Notes
489054	100 g	Plastic bottle	

Litmus paper
• Cartine tornasole • Papier tournesol • Papel tornasol • Sonnenblumenpapier

Classification transport
ONU: 2025

Litmus paper > RS - For pHmetry

RS

Code	Size	Packaging	Notes
435260000	1 roll	Dispenser	Blue litmus paper, Color change: blue --> red, change pH 8.0 - 5.0
435300000	1 roll	Dispenser	Neutral litmus paper, Color change: red <-- purple --> Blue, Change pH 5.0 - 8.0
435340000	1 roll	Dispenser	Red litmus paper, Color change: red --> blue, Change pH 5.0 - 8.0

Roll dispenser 5 m by 7 mm

Lugol concentrated solution
• Lugol "forte" soluzione concentrata • Réactif de Lugol solution concentrée • Lugol reactivo solución concentrada • Lugol Reagenz konzentrierte Lösung



Warning
H373
P260-P314-P501a

Lugol concentrated solution > RS - For microscopy

RS

Description ..... Liquido limpido rosso bruno Identification ..... Positive Assay as iodine (oxidimetric) ..... 4.9 ÷ 5.1 % (p/v)

Code	Size	Packaging	Notes
458741	1 l	Glass bottle	

Lugol solution for Gram-Hucker kit
• Lugol soluzione in acqua per kit Gram-Hucker • Réactif de Lugol solution pour kit de Gram-Hucker • Lugol solución para kit Gram - Hucker
• Lugol-Reagenzlösung für das Gram-Hucker-Kit



Danger
H318-H373-H411
P260-P280i-P305+P351+P338-P310a-P314-P501a

Lugol solution for Gram-Hucker kit > RS - For bacteriology

RS

Description ..... Brown clear liquid Identification ..... Positive Assay ..... 0.85 ÷ 0.95 %p/p (I2)

Code	Size	Packaging	Notes
458751	250 ml	Glass bottle	In Vitro Diagnostic Medical Device

Lugol's Reagent Iodine-Iodide Solution
• Lugol reattivo soluzione iodo-iodurata • Réactif de Lugol solution • Lugol reactivo solución yodo-yodurata • Lugol-Reagenzlösung



Warning
H373
P260-P314-P501a

Lugol's Reagent Iodine-Iodide Solution > RS - For colposcopy

RS

Description ..... Brown clear liquid Identification ..... Positive Assay as iodine (oxidimetric) ..... 1.9 ÷ 2.1 % p/v

Code	Size	Packaging	Notes
458762	250 ml	Glass bottle	Medical Device
E458761	6 x 250 ml	Glass bottle	Medical Device
458763	1 l	Glass bottle	Medical Device

Product specifications are subject to changes. Please visit our website for updates.

Luminol

• Luminol • Luminol • Luminol • Luminol

Synonym:
5-Amino-2,3-dihydro-1,4-phthalazinedione

$C_8H_7N_3O_2$
Molecular Weight: 177,17
CAS: 521-31-3
EEC-N: 208-309-4

Warning
H302-H332
P261-P264-P271-P301+P312a-P304+P340-P501a

Luminol > RPE - For analysis

RPE

Description ..... Polvere verde-giallog. Identification ..... Positive Assay (acidimetric) ..... 97.5 ÷ 102.5 %

Code	Size	Packaging	Notes
458772	25 g	Glass bottle	

For chemiluminescence

Lutetium standard solution

• Lutezio standard soluzione • Lutéthium solution standard • Lutecio, solución patrón • Lutethium-Standardlösung

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group III

Lutetium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505707	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505708	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505709	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Lutetium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503611	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503613	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503615	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503617	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

L(+)Lysine monohydrochloride

• L(+)Lisina monocloridrato • L(+)Lysine monochlorhydrate • L(+)Lisina monoclorhidrato
• L (+) Lysinmonohydrochlorid

Synonym:
(S)-2,6-Diaminohexanoic acid monohydrochloride

$C_6H_{14}N_2O_2.HCl$
Molecular Weight: 182,65
CAS: 657-27-2
EEC-N: 211-519-9

L(+)Lysine monohydrochloride > RPE - For analysis

RPE

Description ..... White powder
Identification ..... Positive
Specific optical rotation... +19.5 ÷ +21.5 °
Loss on drying ..... ≤0.3 %
Ammonium ..... ≤10 ppm
Total phosphorus ..... ≤5 ppm
Water-insoluble matter ..... ≤50 ppm
Heavy metals (Pb) ..... ≤10 ppm
Residue on ignition ..... ≤500 ppm
Total sulphur ..... ≤30 ppm
Fe ..... ≤10 ppm
Assay (ex nitrogen) ..... ≥99 %
Assay (argentimetric) ..... ≥99 %

Code	Size	Packaging	Notes
458122	25 g	Glass bottle	
458124	100 g	Plastic bottle	
458121	5 kg	Plastic tank	

Magnesium, powder

• Magnesio, polvere • Magnésium, poudre • Magnesio, polvo • Magnesiumpulver

Mg
Molecular Weight: 24,31
CAS: 7439-95-4
EEC-N: 231-104-6

Classification transport
ONU: 1869
Transport Hazard class: 4.1
Packing group III

Danger
H228-H251-H261
P210-P223-P231a+P232-P241-P280-P402+P404

Magnesium, powder > RPE - For analysis

RPE

Description ..... Grey powder Identification ..... Positive Titolo (AAS/ICP) ..... ≥ 99 % (Mg)

Code	Size	Packaging	Notes
459066	500 g	Metallic can	

Magnesium, ribbon

• Magnesio, nastro • Magnesium, rubans • Magnesio, tiras • Magnesiumbänder

Mg
Molecular Weight: 24,31
CAS: 7439-95-4
EEC-N: 231-104-6

Classification transport
ONU: 1869
Transport Hazard class: 4.1
Packing group III

Danger
H228-H251-H261
P210-P223-P231a+P232-P241-P280-P402+P404

Magnesium, ribbon > RPE - For analysis

RPE

Description ..... Ribbon Identification ..... Positive Assay ..... ≥ 99 %

Code	Size	Packaging	Notes
459044	100 g	Carton box	

Size ~ 0.2 x 3 mm

Magnesium, turnings

• Magnesio, tornitura • Magnésium, tournures • Magnesio, virutas • Magnesium dreht sich um

Mg
Molecular Weight: 24,31
CAS: 7439-95-4
EEC-N: 231-104-6

Classification transport
ONU: 1869
Transport Hazard class: 4.1
Packing group III

Danger
H228-H251-H261
P210-P223-P231a+P232-P241-P280-P402+P404

Magnesium, turnings > RPE - For analysis

RPE

Description ..... Silvery turnings Identification ..... Positive Other metals ..... ≤ 0.10 % Assay ..... ≥ 99.80 % (Mg)

Code	Size	Packaging	Notes
459085	250 g	Metallic can	

According to Grignard

Magnesium standard solution

• Magnesio standard soluzione • Magnésium solution standard • Magnesio, solución patrón • Magnesiumstandardlösung

Warning
H290
P234-P390-P406

Magnesium standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615001801	100 ml	Plastic bottle	A 10 ppm solution: to dilute according to Ref Ph.Eur 5001801
615001802	100 ml	Plastic bottle	A 10 ppm solution: to dilute according to Ref Ph.Eur 5001802
615001803	100 ml	Plastic bottle	A 0.1 % solution Ref Ph.Eur 5001803
615001809	100 ml	Plastic bottle	A 100 ppm solution: to dilute according to Ref Ph.Eur 5001800

Product specifications are subject to changes. Please visit our website for updates.

Magnesium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505712	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505715	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505713	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Magnesium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503711	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503713	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503715	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503717	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Magnesium standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
503718	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503719	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
E497535	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
E497531	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
507039	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Magnesium standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
458891		Plastic ampoule	conc. 1.000 ppm Matrix: Nitric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Magnesium standard solution > RS - Standard solution for ion chromatography

RS

Code	Size	Packaging	Notes
503291	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503293	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503390	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water and nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Magnesium acetate tetrahydrate

• Magnesio acetato tetraidrato • Magnésium acétate tétrahydraté • Magnesio acetato tetrahidratado
• Magnesiumacetat-Tetrahydrat

Synonym:
Acetic acid magnesium salt

$Mg(CH_3COO)_2.4H_2O$
Molecular Weight: 214,46
CAS: 16674-78-5
EEC-N: 205-554-9

Magnesium acetate tetrahydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description	White crystalline powder	Sulphate	≤ 50 ppm	Fe	≤ 5 ppm	Sr	≤ 50 ppm
Identification	Positive	Heavy metals (Pb)	≤ 5 ppm	K	≤ 50 ppm	Assay (complexometric)	98.0 ÷ 102.0 %
Water-insoluble matter	≤ 50 ppm	Ba	≤ 10 ppm	Mn	≤ 10 ppm		
Chloride	≤ 10 ppm	Ca	≤ 100 ppm	Na	≤ 50 ppm		

Code	Size	Packaging	Notes
459135	250 g	Plastic bottle	
459137	1 kg	Plastic bottle	
459131	25 kg	Drum	

Magnesium carbonate basic

• Magnesio carbonato basico • Magnésium carbonate basique • Magnesio carbonato básico
• Magnesiumcarbonat basisch

Synonym:
Magnesium hydroxide carbonate

$(MgCO_3)_4.Mg(OH)_2.5H_2O$
Molecular Weight: 485,69
CAS: 39409-82-0
EEC-N: 235-192-7

Magnesium carbonate basic > RPE - For analysis

RPE

Description	Light white powder	Chloride	≤300 ppm	Ca	≤ 0.2 %	Pb	≤ 0.5 ppm
Identification	Positive	Heavy metals (Pb)	≤10 ppm	Cd	≤ 0.5 ppm	Loss on calcination	55.0 - 60.0 %
Water solubility	≤0.5 %	Sulphate	≤ 0.20 %	Fe	≤400 ppm	Assay (complexometric) (MgO)	40.0 ÷ 45.0 %
Insoluble in acetic acid	≤ 0.05 %	As	≤2 ppm	Hg	≤ 0.1 ppm		

Code	Size	Packaging	Notes
459285	250 g	Plastic bottle	
459287	1 kg	Plastic bottle	

Magnesium carbonate basic > ERBApharm - According to pharmacopoeia: BP-FU-Ph.Eur.-Ph.Franc.

ERBApharm

Description	White powder	Chloride	≤0.07 %	Acet.ac.not soluble ma.	≤0.05 %	Assay (Mg oxide)	40.0 ÷ 45.0 %
Identification	Positive	Heavy metals (Pb)	≤20 ppm	As	≤2 ppm		
Appearance of solution	Conform Ph.Eur.	Sulphate	≤0.3 %	Ca	≤0.75 %		
Apparent density	60 ÷ 150 g/l	Water-soluble subst.	≤1.0 %	Fe	≤0.04 %		

Code	Size	Packaging	Notes
349257	1 kg	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Magnesium carbonate basic > ERBApharm - According to pharmacopoeia: USP

ERBApharm

Description	White powder	Heavy metals (Pb)	≤30 ppm	As	≤4 ppm	Assay (Mg oxide)	40.0 ÷ 43.5 %
Identification	Positive	Soluble salts	≤1.0 %	Ca	≤0.45 %	Escherichia coli	Absent
Apparent density	400 ÷ 500 g/l	Acid not soluble matter	≤0.05 %	Fe	≤0.02 %		

Code	Size	Packaging	Notes
349279	5 kg	Plastic tank	
349272	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade. Heavy powder

Product specifications are subject to changes. Please visit our website for updates.

Magnesium chloride hexahydrate

• Magnesio cloruro esaidrato • Magnésium chlorure hexahydraté • Magnesio cloruro hexahidrato • Magnesiumchlorid-Hexahydrat

$MgCl_2.6H_2O$
Molecular Weight: 203,31
CAS: 7791-18-6
EEC-N: 232-094-6

Magnesium chloride hexahydrate > RPE - For analysis - ACS - ISO

RPE

Description	White crystals	Phosphate	≤5 ppm	Ba	≤50 ppm	Mn	≤5 ppm
Identification	Positive	Nitrate	≤10 ppm	Ca	≤100 ppm	Na	≤50 ppm
Water-insoluble matter	≤50 ppm	Sulphate	≤20 ppm	Fe	≤5 ppm	Sr	≤50 ppm
Ammonium	≤20 ppm	Heavy metals (Pb)	≤5 ppm	K	≤50 ppm	Assay (complexometric)	99.0 ÷ 102.0 %

Code	Size	Packaging	Notes
459336	100 g	Plastic bottle	
459337	1 kg	Plastic bottle	
459331	5 kg	Plastic tank	
459332	25 kg	Plastic bucket	
459334	50 kg	Fibre drum	

Magnesium chloride hexahydrate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP

ERBApharm

Description	Colourless hygroscopic crystals	Ba	Conform USP-NF	Sulphate	≤50 ppm	Assay (complexometric)	98.0 ÷ 101.0 %
Identification B	Positive	K	Conform USP-NF	Not soluble matter	≤50 ppm	Origin (BSE/TSE)	Synthesis
Identification C	Positive	pH solution 5%	4.5 ÷ 7.0	Al	≤1 ppm	Residual solvents (Current ICH)	Conform
Identification (I.R.)	Positive	Water (K.F.)	51.0 ÷ 55.0 %	As	≤2 ppm		
Appearance of solution	Conform Ph.Eur.	Heavy metals (Pb)	≤10 ppm	Ca	≤100 ppm		
Acidity or alkalinity	Conform Ph.Eur.	Bromide	≤500 ppm	Fe	≤10 ppm		

Code	Size	Packaging	Notes
349357	1 kg	Plastic bottle	
349359	5 kg	Plastic tank	
349355	25 kg	Drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Magnesium chloride hexahydrate > RE - Pure

RE

Description	White cristalline powder	As	≤ 2 ppm	Cu	≤ 10 ppm	Pb	≤ 10 ppm
Identification	Positive	Al	≤ 1 ppm	Fe	≤ 5 ppm	Zn	≤ 10 ppm
Ammonium	≤ 50 ppm	Ba	≤ 20 ppm	K	≤ 3000 ppm	Assay (complexometric)	98.5 ÷ 101.5 %
Sulphate	≤ 100 ppm	Ca	≤ 1000 ppm	Na	≤ 5000 ppm		

Code	Size	Packaging	Notes
349377	1 kg	Plastic bottle	
349372	25 kg	Plastic bucket	

Magnesium glycerophosphate

• Magnesio glicerofosfato • Magnesium glycerophosphate • Magnesio glicerofosfato
• Magnesiumglycerophosphat

Synonym:
DL-alpha-glycerol phosphate magnesium salt hydrate

$C_3H_7O_2PO_4Mg$
Molecular Weight: 194,36
CAS: 927-20-8
EEC-N: 231-149-3

Magnesium glycerophosphate > RE - Pure

RE

Description	White crystalline powder	Residue on calcination	47 ÷ 52 %	Heavy metals (Pb)	≤ 10 ppm	Assay (Mg)	≥ 10.6 %
Identification	Positive	Glycerol-Alc.sol.impur.	≤ 2 %	As	≤ 4 ppm		

Code	Size	Packaging	Notes
349407	1 kg	Plastic bottle	

Magnesium hydrogen phosphate trihydrate

• Magnesio fosfato bibasico triidrato • Magnesium phosphate dibasique trihydraté
• Magnesio fosfato dibásico trihidrato • Magnesium hydrogenphosphat Trihydrat

Synonym:
• *Magnesium phosphate dibasic trihydrate*
• *Newberyite*

$MgHPO_4.3H_2O$
Molecular Weight: 174,34
CAS: 7782-75-4

Magnesium hydrogen phosphate trihydrate > RPE - For analysis

RPE

Description	White powder	HCl-insoluble matter	≤500 ppm	Cu	≤50 ppm	Pb	≤50 ppm
Identification	Positive	Sulphate	≤ 60 ppm	Fe	≤50 ppm	Assay (complexometric)	97 ÷ 100 %
Chloride	≤100 ppm	As	≤1 ppm	Ni	≤ 50 ppm		

Code	Size	Packaging	Notes
459435	250 g	Plastic bottle	
459437	1 kg	Plastic bottle	

Magnesium hydroxide

• Magnesio idrossido • Magnésium hydroxyde • Magnesio hidróxido • Magnesiumhydroxid

$Mg(OH)_2$
Molecular Weight: 58,33
CAS: 1309-42-8
EEC-N: 215-170-3

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Magnesium hydroxide > ERBAPharm - According to pharmacopoeia: BP-FU-Ph.Eur.-Ph.Franc.

ERBApharm

Description	White powder	Chloride	≤ 0.1 %	Acet.ac.not soluble ma.	≤ 0.1 %	Assay (complexometric)	95.0 ÷ 100.5 %
Identification	Positive	Heavy metals (Pb)	≤ 30 ppm	As	≤ 4 ppm		
Appearance of solution	Conform Ph.Eur.	Sulphate	≤ 1.0 %	Ca	≤ 1.5 %		
Loss on ignition	29.0 ÷ 32.5 %	Soluble matter	≤ 2.0 %	Fe	≤ 0.07 %		

Code	Size	Packaging	Notes
349455	1 kg	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Magnesium hydroxide carbonate ▶ Magnesium carbonate basic

Magnesium nitrate hexahydrate

• Magnesio nitrato esaidrato • Magnésium nitrate hexahydraté • Magnesio nitrato hexahidratado • Magnesiumnitrat-Hexahydrat

$Mg(NO_3)_2.6H_2O$
Molecular Weight: 256,41
CAS: 13446-18-9
EEC-N: 233-826-7

Classification transport
ONU: 1474
Transport Hazard class: 5.1
Packing group III

Danger
H272-H319
P210-P220-P264-P280-P305+P351+P338-P337+P313

Magnesium nitrate hexahydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description	White crystals	Chloride	≤10 ppm	Ca	≤100 ppm	Sr	≤50 ppm
Identification	Positive	Phosphate	≤5 ppm	Fe	≤5 ppm	Assay (complexometric)	98.0 ÷ 102.0 %
pH sol. 5% at 25° C	5.0 ÷ 8.2	Sulphate	≤50 ppm	K	≤50 ppm		
Water-insoluble matter	≤50 ppm	Heavy metals (Pb)	≤5 ppm	Mn	≤5 ppm		
Ammonium	≤30 ppm	Ba	≤50 ppm	Na	≤50 ppm		

Code	Size	Packaging	Notes
459535	100 g	Plastic bottle	
459536	500 g	Plastic bottle	
459537	1 kg	Plastic bottle	

Magnesium nitrate hexahydrate > RE - Pure

RE

Description	White or yellowish pellets	Identification	Positive	Fe	≤ 10 ppm	Assay (complexometric)	≥ 96 %

Code	Size	Packaging	Notes
349557	1 kg	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Magnesium nitrate 10 g/l solution

• Magnesio nitrato 10 g/l soluzione • Magnésium nitrate 10 g/l • Magnesio nitrato solución 10 g/l • Magnesiumnitrat 10 g/l

HEU210

Magnesium nitrate 10 g/l solution > RS - Matrix modifiers for AAS-GTA

RS

Code	Size	Packaging	Notes
503196	50 ml	Plastic bottle	Matrix: Water

Magnesium oxide

• Magnesio ossido • Magnésium oxyde • Magnesio óxido • Magnesiumoxid

MgO
Molecular Weight: 40,31
CAS: 1309-48-4
EEC-N: 215-171-9

Magnesium oxide > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611049901	100 g	Bottle	

Magnesium oxide > RS - For chromatography

RS

Description ..... White powder Identification ..... Positive

Code	Size	Packaging	Notes
459617	1 kg	Plastic bottle	

Magnesium oxide > RPE - For analysis

RPE

Description ..... White powder
Identification ..... Positive
Loss on ignition ..... ≤ 2 %
Assay (complexometric) ..... ≥ 99 %
Chloride ..... ≤ 0.25 %
Silicate ..... ≤ 500 ppm
Sulphat + sulphit (SO4) ..... ≤ 0.25 %
Ca ..... ≤ 0.3 %
Fe ..... ≤ 400 ppm
Na ..... ≤ 0.1 %

Code	Size	Packaging	Notes
459584	100 g	Plastic bottle	
459586	500 g	Plastic bottle	
459587	1 kg	Plastic bottle	

Magnesium oxide heavy

• Magnesio ossido pesante • Magnésium oxyde lourd • Magnesio óxido pesado • Magnesiumoxid schwer

MgO
Molecular Weight: 40,31
CAS: 1309-48-4
EEC-N: 215-171-9

Magnesium oxide heavy > ERBApharm - According to pharmacopoeia: Ph.Eur.

ERBApharm

Description ..... White powder
Identification ..... Positive
Appearance of solution ..... Conform Ph.Eur.
Apparent density ..... ≥0.25 g/ml
Loss on ignition ..... ≤8.0 %
Chloride ..... ≤0.1 %
Heavy metals (Pb) ..... ≤30 ppm
Sulphate ..... ≤1.0 %
Soluble matter ..... ≤2.0 %
Acet.ac.not soluble ma. ..... ≤0.1 %
As ..... ≤4 ppm
Ca ..... ≤1.5 %
Fe ..... ≤0.07 %
Assay (complexometric) Calc. .... 98.0 ÷ 100.5 %

Code	Size	Packaging	Notes
349655	1 kg	Plastic bottle	
349656	5 kg	Plastic tank	
349653	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Magnesium perchlorate

• Magnesio perclorato • Magnsium perchlorate • Magnesio perclorato • Magnesiumperchlorat

$Mg(ClO_4)_2$
Molecular Weight: 223,21
CAS: 10034-81-8
EEC-N: 233-108-3

Classification transport
ONU: 1475
Transport Hazard class: 5.1
Packing group II

Danger
H272-H315-H319-H335
P210-P261-P280-P304+P340-P305+P351+P338-P403+P233

Magnesium perchlorate > RPE - For analysis

RPE

Description White crystalline powder or flakes and/or chunks
Identification ..... Positive
Loss on drying (190°C) ..... ≤ 8 %
Titrable free acid ..... ≤ 0.005 meq/g
Titrable base ..... ≤ 0.025 meq/g
Suitability (for moisture absorption) ..... Passes test

Code	Size	Packaging	Notes
422254	100 g	Glass bottle	
422251	250 g	Glass bottle	
422252	1 kg	Glass bottle	

Magnesium peroxide

• Magnesio perossido • Magnésium péroxyde • Magnesio peróxido • Magnesiumperoxid

MgO_2
Molecular Weight: 56,3
CAS: 14452-57-4
EEC-N: 238-438-1

Classification transport
ONU: 1476
Transport Hazard class: 5.1
Packing group II

Danger
H272
P210-P220-P280-P370+P378a-P501a

Magnesium peroxide > RE - Pure

RE

Description ..... White powder
Identification ..... Positive
Assay (oxidimetric) ..... ≥15 %

Code	Size	Packaging	Notes
349757	1 kg	Plastic bottle	
349753	25 kg	Drum	

Magnesium phosphate dibasic trihydrate ▶ Magnesium hydrogen phosphate trihydrate

Magnesium stearate

• Magnesio stearato • Magnésium stéarate • Magnesio estearato • Magnesiumstearat

Synonym:
Stearic acid magnesium salt

$[CH_3(CH_2)_{16}CO_2]_2Mg$
Molecular Weight: 591,27
CAS: 557-04-0
EEC-N: 209-150-3

Magnesium stearate > ERBApharm - Vegetal origin-According to pharmacopoeia: Ph.Eur.-BP-FU-NF

ERBApharm

Description ..... White powder
Identification ..... Positive
Acidity or alkalinity ..... Conform Ph.Eur.
Loss on drying ..... ≤ 6.0 %
Stearic acid ..... ≥ 40.0 %
Stearic+Palmitic acid ..... ≥ 90.0 %
Chloride ..... ≤ 0.1 %
Sulphate ..... ≤ 1.0 %
Cd ..... ≤ 3 ppm
Ni ..... ≤ 5 ppm
Pb ..... ≤ 10 ppm
Assay (magnesium) ..... 4.0 ÷ 5.0 % s.s.
Microbial tests
TAMC ..... ≤ 1000 CFU/g
TYMC ..... ≤ 100 CFU/g
Escherichia coli ..... Absent
Salmonella ..... Absent

Code	Size	Packaging	Notes
350032	2.5 kg	Plastic bottle	
350033	20 kg	Fibre drum	
350035	25 kg	Fibre drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Product specifications are subject to changes. Please visit our website for updates.

Magnesium sulfate anhydrous

• Magnesio solfato anidro • Magnésium sulfate anhydre • Magnesio sulfato anhidro • Magnesiumsulfat wasserfrei

$MgSO_4$
Molecular Weight: 120,36
CAS: 7487-88-9
EEC-N: 231-298-2

Magnesium sulfate anhydrous > RE - Pure - For anhydrification

RE

Appearance ..... White powder
Identification ..... Conform
Water content ..... ≤ 2 % m/m
Assay (on dry) ..... ≥ 99 %
Bulk density ..... ~ 760 G/L

Code	Size	Packaging	Notes
P1460012	1 kg	Plastic bottle	
P1460027	5 kg	Plastic tank	
P1460044	25 kg	Plastic bucket	
P1460057	50 kg	Fibre drum	

Magnesium sulfate heptahydrate

• Magnesio solfato eptaidrato • Magnésium sulfate heptahydraté • Magnesio sulfato heptahidratado • Magnesiumsulfat heptahydrat

Synonym:
Epsom salts

$MgSO_4.7H_2O$
Molecular Weight: 246,48
CAS: 10034-99-8
EEC-N: 231-298-2

Magnesium sulfate heptahydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White crystalline powder or flakes
Identification ..... Positive
pH sol. 5% at 25° C ..... 5.0 ÷ 8.2
Water-insoluble matter ..... ≤50 ppm
Ammonium ..... ≤20 ppm
Chloride ..... ≤5 ppm
Nitrate ..... ≤20 ppm
Heavy metals (Pb) ..... ≤5 ppm
Ca ..... ≤200 ppm
Fe ..... ≤5 ppm
K ..... ≤50 ppm
Mn ..... ≤5 ppm
Na ..... ≤50 ppm
Sr ..... ≤50 ppm
Assay (complexometric) 98.0 ÷ 102.0 % s.s.
Appearance of solution ..... Conform
Acidity or alkalinity ..... Conform
Loss on drying ..... 48.0 ÷ 52.0 %
As ..... ≤2 ppm

Code	Size	Packaging	Notes
459665	100 g	Plastic bottle	
459666	500 g	Plastic bottle	
459667	1 kg	Plastic bottle	
459669	5 kg	Plastic tank	
459662	25 kg	Plastic bucket	

Magnesium sulfate heptahydrate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-BP-DAB

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution ..... Conform Ph.Eur.
Acidity or alkalinity ..... Conform Ph.Eur.
pH solution 5% ..... 5.0 ÷ 9.2
Loss on drying ..... 48.0 ÷ 52.0 %
Loss on ignition ..... 40.0 ÷ 52.0 %
Chloride ..... ≤140 ppm
Heavy metals (Pb) ..... ≤10 ppm
As ..... ≤2 ppm
Fe ..... ≤20 ppm
Se ..... ≤30 ppm
Assay (complexometric) 99.0 ÷ 100.5 % s.s.

Code	Size	Packaging	Notes
349852	1 kg	Plastic bottle	
349859	5 kg	Plastic tank	
349851	25 kg	Drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Maize starch

• Amido di mais • Amidon de maïs • Almidón de maiz • Maisstärke

$(C_6H_{10}O_5)n$
CAS: 9005-84-9
EEC-N: 232-686-4

Maize starch > ERBApharm - According to pharmacopoeia: BP-FU-NF-Ph.Eur.-Ph.Franc.

ERBApharm

Description ..... White or yellowish powder
Identification ..... Positive
Zolfo diossido ..... ≤ 50 ppm
Foreingn cellular elem. ..... Conform Ph.Eur.
pH at 20° C ..... 4.0 ÷ 7.0
Loss on drying ..... ≤ 15.0 %
Sulphated ash ..... ≤ 0.6 %
Oxidizing substances ..... ≤ 20 ppm
Fe ..... ≤ 10 ppm
Microbial tests
TAMC ..... ≤ 1000 CFU/g
TYMC ..... ≤ 100 CFU/g
Escherichia coli ..... Absent Ph.Eur.
Salmonella ..... Absent Ph. Eur.

Code	Size	Packaging	Notes
313071	1 kg	Plastic bottle	
313072	5 kg	Plastic tank	
313073	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Malachite green

• Verde malachite • Vert de malachite • Verde de malaquita • Malachitgrün

Synonym:
Basic green 4

$C_{23}H_{25}ClN_2$
Molecular Weight: 364,92
CAS: 569-64-2
EEC-N: 209-322-8

Classification transport
ONU: 3143
Transport Hazard class: 6.1
Packing group III



Danger
H302-H318-H361d-H410
P264-P280-P301+P312a-P305+P351+P338-P310a-P308+P313

Malachite green > RS - For microscopy - C.I. 42000

RS

Description ..... Green shining crystals
Identification ..... Positive

Code	Size	Packaging	Notes
491303	25 g	Glass bottle	
491304	100 g	Plastic bottle	

Dye for cytology

Malachite green solution 0.5% in acetic acid anhydrous

• Verde malachite soluzione 0.5% in acido acetico anidro
• Vert malachite solution 0.5% dans l'acide acétique anhydre
• Verde de malaquita solución 0.5% en acido acético anhidro
• Grüne Malachitlösung 0.5% in wasserfreier Essigsäure

Synonym:
Basic green 4

Classification transport
ONU: 3264

Malachite green solution 0.5% in acetic acid anhydrous > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611050501	1 l	Glass bottle	Ref Ph.Eur 1050501

Product specifications are subject to changes. Please visit our website for updates.

Maleic acid

• Acido maleico • Acide maléique • Acido maleico • Maleinsäure

Synonym:
• *cis-Butenedioic acid*
• *Toxilic acid*

HOOCCH:CHCOOH
Molecular Weight: 116,07
CAS: 110-16-7
EEC-N: 203-742-5

Classification transport
ONU: 3261
Transport Hazard class: 8
Packing group III

Warning
H302-H315-H319-H317-H335
P261-P271-P304+P340-P305+P351+P338-P337+P313-P403+P233

Maleic acid > ERBApharm - According to pharmacopoeia: BP-Ph.Eur.-USP-NF

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Fumaric acid .....Conform Ph.Eur.
Melting point..... 132 ÷ 135 °C
Water (K.F.) ..... ≤2.0 %
Sulphated ash ..... ≤0.1 %
Heavy metals (Pb).....≤10 ppm
Fe .....≤5 ppm
Assay (acidimetric) .....99.0 ÷ 101.0 % s s
Origin (BSE/TSE)..... Synthesis
Residual solvents (Current ICH).....Conform

Code	Size	Packaging	Notes
407266	500 g	Plastic bottle	
407261	5 kg	Plastic tank	
407263	25 kg	Fibre drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Maleic anhydride

• Anidride maleica • Anhydride maléique • Anhídrido maleico • Maleinsäureanhydrid

Synonym:
2,5-Furandione

OCOCH:CHCO
Molecular Weight: 98,06
CAS: 108-31-6
EEC-N: 203-571-6

Classification transport
ONU: 2215
Transport Hazard class: 8
Packing group III

Danger
H302-H314-H334-H317
P280-P284-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P362+P364-P342+P311a

Maleic anhydride > RPE - For analysis

RPE

Description .....White flakes
Identification ..... Positive
Melting point..... 52 ÷ 55 ° C
Acido maleico libero .....≤ 0.5 %
Assay (acidimetric) .....≥ 98.5 %

Code	Size	Packaging	Notes
421955	250 g	Plastic bottle	

DL-Malic acid

• Acido DL-malico • Acide DL-malique • Acido DL-málico • DL-Apfelsäure

Synonym:
• *(±)-2-Hydroxysuccinic acid*
• *DL-Hydroxybutanedioic acid*

HOOCCH(OH)CH_2COOH
Molecular Weight: 134,09
CAS: 617-48-1
EEC-N: 210-514-9

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

DL-Malic acid > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Water-insoluble matter ..... ≤ 0.1 %
Melting point..... 130 ÷ 132 ° C
Water ..... ≤ 0.3 %
Chloride.....≤ 5 ppm
Heavy metals (Pb).....≤ 20 ppm
Sulphated ash..... ≤ 0.05 %
As .....≤ 3 ppm
Pb .....≤ 10 ppm
Assay (acidimetric) ..... ≥ 99.0 %

Code	Size	Packaging	Notes
407314	100 g	Plastic bottle	
407316	500 g	Plastic bottle	

Malonic acid

• Acido malonico • Acide malonique • Acido malónico • Malonsäure

Synonym:
Propanedioic acid

$HOOCCH_2COOH$
Molecular Weight: 104,06
CAS: 141-82-2
EEC-N: 205-503-0

Warning
H302-H319
P264-P280i-P301+P312a-P305+P351+P338-P337+P313-P501a

Malonic acid > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Melting point..... 133 ÷ 136 °C
Chloride..... ≤ 100 ppm
Assay (acidimetric) ..... ≥ 98.5 %
Residue on ignition..... ≤ 0.5 %
Sulphate..... ≤ 0.1 %

Code	Size	Packaging	Notes
407363	50 g	Glass bottle	

Maltose monohydrate

• Maltosio monoidrato • Maltose monohydraté • Maltosa monohidrato • Maltose-Monohydrat

Synonym:
4-O-α-D-Glucopyranosyl-D-glucose

$C_{12}H_{22}O_{11}.H_2O$
Molecular Weight: 360,32
CAS: 6363-53-7
EEC-N: 200-716-5

Maltose monohydrate > RPE - For analysis

RPE

Description ..... Whitish powder
Identification ..... Positive
Potere rotat. spec. (C=2; H20 24h) .....130 ± 3 °
Water .....≤ 7.5 %
Melting point..... ca. 120 °C (dec)
Assay .....≥ 94.0 %

Code	Size	Packaging	Notes
459863	50 g	Glass bottle	
459865	250 g	Plastic bottle	

Maltose monohydrate > RE - Pure

RE

Description .....White crystalline powder
Identification ..... Positive
Water .....≤ 7 %
Sulphated ash.....≤ 0.1 %
Assay .....≥ 92 %

Code	Size	Packaging	Notes
350401	25 kg	Plastic bucket	

D(-)Mandelic acid

• Acido D(-)mandelico • Acide D(-)mandélique • Acido D(-)mandélico • D (-) Mandelsäure

Synonym:
(R)-alpha-hydroxyphenylacetic acid

$C_6H_5CH(OH)COOH$
Molecular Weight: 152,15
CAS: 611-71-2
EEC-N: 210-276-6

D(-)Mandelic acid > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Melting point..... 133 ÷ 135 ° C
Potere rotator. specif.(C=2 in Acqua)-155.4
÷ -153.4 °

Code	Size	Packaging	Notes
407421	25 g	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Manganese electrolytic

• Manganese elettrolitico • Manganèse électrolytique • Manganeso electrolítico • Elektrolytisches Mangan

Mn
Molecular Weight: 54,94
CAS: 7439-96-5
EEC-N: 231-105-1

Manganese electrolytic > RPE - For analysis

RPE

Description ..... Brown irregular flakes Identification ..... Positive Assay ..... 99.9 ÷ 100.0 %

Code	Size	Packaging	Notes
459965	250 g	Glass bottle	

Manganese standard solution

• Manganese standard soluzione • Manganèse solution standard • Manganeso, solución patrón • Mangan-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Manganese standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615004500	1 l	Plastic bottle	A 100 ppm solution Ref Ph.Eur 5004500
615005800	1 l	Plastic bottle	A 1.000 ppm solution Ref Ph.Eur 5005800

Manganese standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505717	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505718	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505719	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Manganese standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503721	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503723	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503725	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503727	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Manganese standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
507746	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507488	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
E497545	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
E497541	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Manganese standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
459911		Plastic ampoule	conc. 1.000 ppm Matrix: Hydrochloric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Manganese standard solution > RS - Quality control standard solution for AAS (graphite furnace)

RS

Code	Size	Packaging	Notes
504362	50 ml	Glass bottle	conc. 20 ppb - Matrix: 2% Nitric acid

Manganese (II) acetate tetrahydrate

• Manganese acetato oso tetraidrato • Manganèse (II) acétate tetrahydrate
• Manganeso (II) acetato tetraidrato • Mangan (II) acetat Tetrahydrat

Synonym:
Manganous acetate

$Mn(CH_3COO)_2.4H_2O$
Molecular Weight: 245,09
CAS: 6156-78-1

Warning
H315-H319-H361-H335
P261-P271-P280-P304+P340-P305+P351+P338-P403+P233

Manganese (II) acetate tetrahydrate > RPE - For analysis

RPE

Description ..... Light pink crystals
pH sol. 5%..... 6 - 7
Assay (complexometry)..... ≥ 99.0 %
Insoluble matter in water ..... ≤ 0.01 %
Chloride..... ≤ 20 ppm
Sulfate..... ≤ 50 ppm
Fe ..... ≤ 50 ppm
Cu..... ≤ 5 ppm
Ni..... ≤ 20 ppm

Code	Size	Packaging	Notes
460005	250 g	Plastic bottle	
460007	1 kg	Plastic bottle	
460001	25 kg	Plastic bucket	

Manganese (II) chloride tetrahydrate

• Manganese cloruro oso tetraidrato • Manganèse (II) chlorure tetrahydrate • Manganeso (II) cloruro tetrahidratado • Mangan (II) chlorid Tetrahydrat

$MnCl_2.4H_2O$
Molecular Weight: 197,91
CAS: 13446-34-9

Warning
H302
P264-P270-P301+P312a-P330-P501a

Manganese (II) chloride tetrahydrate > RPE - For analysis

RPE

Description ..... Pink crystalline powder
Identification ..... Positive
Substances reducing KMnO4 (O). ≤ 0.0005
%
Sulphate..... ≤ 50 ppm
Ca..... ≤ 10 ppm
Fe ..... ≤ 5 ppm
Mg ..... ≤ 50 ppm
Pb..... ≤ 5 ppm
Zn..... ≤ 10 ppm
Assay ..... ≥ 99 %

Code	Size	Packaging	Notes
460158	100 g	Plastic bottle	
460156	500 g	Plastic bottle	
460159	1 kg	Plastic bottle	

Manganese (II) chloride tetrahydrate > RE - Pure

RE

Description ..... Pink crystals
Identification ..... Positive
Water-insoluble matter ..... ≤500 ppm
Sulphate..... ≤500 ppm
Fe ..... ≤50 ppm
Assay (complexometric) ..... ≥98 %

Code	Size	Packaging	Notes
351507	1 kg	Plastic bottle	
351508	5 kg	Plastic tank	
351502	25 kg	Plastic bucket	

Product specifications are subject to changes. Please visit our website for updates.

Manganese (II) sulfate monohydrate

• Manganese solfato oso monoidrato • Manganèse (II) sulfate monohydraté • Manganeso (II) sulfato monohidrato • Mangan(II) sulfat Monohydrat

$MnSO_4.H_2O$
Molecular Weight: 162,09
CAS: 10034-96-5
EEC-N: 232-089-9



Warning
H373-H411
P260-P273-P314-P391-P501a

Manganese (II) sulfate monohydrate > RPE - For analysis - ACS

RPE

Description ..... Pink powder
Identification ..... Positive
Loss on ignition ..... 10.0 ÷ 12.0 %
Water-insoluble matter ..... ≤100 ppm
Subst. reducing KMn04 ..... Conform
Chloride ..... ≤50 ppm
Heavy metals (Pb) ..... ≤20 ppm
Ca ..... ≤50 ppm
Fe ..... ≤20 ppm
K ..... ≤100 ppm
Mg ..... ≤50 ppm
Na ..... ≤500 ppm
Ni ..... ≤200 ppm
Zn ..... ≤50 ppm
Assay (complexometric) .... 98.0 ÷ 101.0 %
Mn ..... 31.8 - 32.8 %

Code	Size	Packaging	Notes
460305	250 g	Plastic bottle	
460307	1 kg	Plastic bottle	

Manganese (II) sulfate monohydrate > RE - Pure

RE

Description ..... Polvere rosata
Identification ..... Positive
Water-insoluble matter ..... ≤ 500 ppm
Pb ..... ≤ 15 ppm
As ..... ≤ 5 ppm
Assay (complexometric) ..... ≥ 98 %
Mn ..... ≥ 31.8 %

Code	Size	Packaging	Notes
352007	1 kg	Plastic bottle	
352008	5 kg	Plastic tank	
352002	25 kg	Plastic bucket	

Manganese (IV) oxide

• Manganese biossido • Manganèse (IV) dioxyde • Manganeso (IV) dióxido • Mangan (IV) dioxid

Synonym:
Manganese dioxide

MnO_2
Molecular Weight: 86,94
CAS: 1313-13-9
EEC-N: 215-202-6



Warning
H302-H332
P261-P264-P271-P301+P312a-P304+P340-P501a

Manganese (IV) oxide > RPE - For analysis

RPE

Description ..... Black powder
Identification ..... Positive
Loss on drying 120° C ..... ≤1.5 %
Assay (oxidimetric) ..... ≥90.0 %

Code	Size	Packaging	Notes
460055	250 g	Plastic bottle	
460056	1 kg	Plastic bottle	
460052	25 kg	Plastic bucket	

D-Mannitol

• D-Mannitolo • D-Mannitol • D-Manitol • D - Mannit

Synonym:
Mannite

$CH_2OH(CHOH)_4CH_2OH$
Molecular Weight: 182,17
CAS: 69-65-8
EEC-N: 200-711-8

D-Mannitol > RPE - For analysis

RPE

Description ..... White crystal. powder
Identification (I.R.) ..... Positive
Melting point ..... 165 ÷ 167 ° C
Power rotat . specif . at 25 ° C (dry) +23.3 ÷ +24.3 °
Loss on drying 105°C ..... ≤ 0.05 %
Acidity (acetic acid) ..... ≤50 ppm
Water-insoluble matter ..... ≤100 ppm
Heavy metals (Pb) ..... ≤5 ppm
Residue on ignition ..... ≤100 ppm
Red.ing sugars(Glucose) ..... ≤0.1 %
Assay (oxidimetric) ..... 98.5 ÷ 101.5 %

Code	Size	Packaging	Notes
460355	250 g	Plastic bottle	
460357	1 kg	Plastic bottle	
460352	5 kg	Plastic tank	
460353	25 kg	Plastic bucket	

D-Mannitol > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-JP

ERBApharm

Description ..... White crystalline powder
Appearance of solution .....Conform Ph.Eur.
Acidity ..... Conform USP-NF
Melting point..... 165 ÷ 167 °C
Heavy metals (Pb).....≤ 5 ppm
Ni.....≤ 1 ppm
Assay (anhydrous) (HPLC) .97.0 ÷ 102.0 %
Reducing sugar..... ≤ 0.1 %
Conductivity..... ≤ 20 µS.cm-1
Microbial tests
TAMC ..... ≤ 1000 CFU/g
TYMC ..... ≤ 100 CFU/g
Escherichia coli ..... Absent Ph. Eur.
Salmonella..... Absent Ph. Eur.
Residue on ignition.....≤ 100 ppm
Identification (I.R.)..... Positive
Identification D ..... Positive
Related substances
Impurity A (HPLC) ..... ≤ 2.0 %
Impurities B+C (HPLC)..... ≤ 2.0 %
Other impurities (HPLC) ..... ≤ 0.10 %
Total impurities (HPLC)..... ≤ 2.0 %
Loss on drying 105°C..... ≤ 0.3 %
Specific optical rotation at 20°C (anh.).+23 ÷ +24 °

Code	Size	Packaging	Notes
352051	1 kg	Plastic bottle	
352052	5 kg	Plastic tank	
352053	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Mascagnite ▶ Ammonium sulfate

May Grünwald reagent

• May Grünwald reattivo • Réactif de May Grünwald • May Grünwald reactivo • Reagenz von Mai Grünwald

Classification transport
ONU: 1992
Transport Hazard class: 3
Packing group II



Danger
H225-H301-H370
P210-P241-P280-P301+P310a-P303+P361+P353-P403+P235

May Grünwald reagent > RS - For hematology

RS

Description .....Blue clear liquid
Identification ..... Positive
pH (20% in water) .....5.4 ÷ 7.5
E 1% / 1 cm a 522 nm.....130 ÷ 180
E 1% / 1 cm a 552 nm.....30 ÷ 40
E 1% / 1 cm a 650 nm.....250 ÷ 350
Functionality .....Conform

Code	Size	Packaging	Notes
460584	100 ml	Plastic bottle	In Vitro Diagnostic Medical Device
E460582	6 x 100 ml	Plastic bottle	In Vitro Diagnostic Medical Device
460586	500 ml	Plastic bottle	In Vitro Diagnostic Medical Device
E460583	6 x 500 ml	Plastic bottle	In Vitro Diagnostic Medical Device
460581	2.5 l	Plastic bottle	In Vitro Diagnostic Medical Device
E460585	4 x 2.5 l	Plastic bottle	In Vitro Diagnostic Medical Device

Dye according to Pappenheim hematology

Mayer's reagent

• Mayer reattivo soluzione in acqua • Réactif de Mayer • Mayer reactivo solución en agua • Reagenz von Mayer

Classification transport
ONU: 2024
Transport Hazard class: 6.1
Packing group II



Danger
H300-H315-H319-H341-H373-H412
P280-P301+P310a-P305+P351+P338-P308+P313-P332+P313-P337+P313

Mayer's reagent > RS - For alkaloids detection

RS

Description ..... Yellow clear liquid Identification ..... Positive

Code	Size	Packaging	Notes
460502	500 ml	Plastic bottle	

MEK ▶ Ethyl methyl ketone

Product specifications are subject to changes. Please visit our website for updates.

Melting point standards

• Punto di fusione standard • Standards de point de fusion • Patrones de puntos de fusión • Schmelzpunktstandards

$CH_3CONHC_6H_5$
Molecular Weight: 135,17
CAS: 103-84-4
EEC-N: 203-150-7

Warning
H302
P264-P270-P301+P312a-P330-P501a

Melting point standards > RS - For calibration

RS

Code	Size	Packaging	Notes
540001	1 g	Bottle	Benzophenone 47 to 49°C
540002	1 g	Bottle	p-Nitrotoluene 52 to 54°C
540003	1 g	Bottle	Vanillin 81 to 83°C
540014	1 g	Bottle	Acetanilide 113 to 116°C
540004	1 g	Bottle	Benzoic Acid 121 to 123°C
540005	1 g	Bottle	Phenacetin 133 to 135°C
540006	1 g	Bottle	Salicylic Acid 158 to 160°C
540007	1 g	Bottle	Sulfanilamide 164 to 166°C
540008	1 g	Bottle	Caffeine 235 to 238°C
540009	1 g	Bottle	Carbazole 243 to 247°C
540010	1 g	Bottle	Anthraquinone 283 to 286°C
540011	3 x 1 g	Bottle	Set Sulphanilamide Caffeine Vanillin
540012	3 x 1 g	Bottle	Set Benzophenone (4749°C) Benzoic Acid (121-123°C) Anthraquinone (283-286°C)
540013	3 x 1 g	Bottle	Set Vanillin (81-83°C) Phenacetin (134-136°C) Caffeine (235-237°C)

L-Menthol

• L-Mentolo • L-Menthol • L-Mentol • L-Menthol

Synonym:
5-Methyl-2-(1-methylethyl)cyclohexanol

$C_{10}H_{19}OH$
Molecular Weight: 156,27
CAS: 2216-51-5
EEC-N: 218-690-9

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

L-Menthol > ERBApharm - According to pharmacopoeia: USP

ERBApharm

Description ..... Colourless crystals
Identification ..... Positive
Specific optical rotation..... -51 ÷ -45 °
Melting point..... 41 ÷ 44 °C
Nonvolatil residue..... ≤ 0.05 %
Related compounds..... Conform USP-NF
Total impurities (GC) ..... ≤ 2.0 %
Individual impurities (GC) ..... ≤ 0.3 %

Code	Size	Packaging	Notes
352103	50 g	Glass bottle	
352106	500 g	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

2-Mercaptoethanol

• 2-Mercaptoetanolo • 2-Mercaptoéthanol • 2-Mercaptoetanol • 2-Mercaptoethanol

Synonym:
• *Thioethylene glycol*
• *2-Hydroxyethylmercaptan*

$HSCH_2CH_2OH$
Molecular Weight: 78,13
CAS: 60-24-2
EEC-N: 200-464-6

Classification transport
ONU: 1750
Transport Hazard class: 6.1
Packing group II

Danger
H301-H310-H314-H411
P280-P301+P310a-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338

2-Mercaptoethanol > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Refractive index at 20° C..... 1.4990 ÷ 1.5020
Water ..... ≤ 0.5 %
Assay (GLC) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
460691	10 ml	Glass bottle	

Mercuric bromide paper

• Carta al bromuro mercurico • Papier bromure mercurique • Mercurio bromuro papel • Quecksilberbromidpapier

Classification transport
ONU: 3082
Transport Hazard class: 9
Packing group III



Warning
H302-H332-H373-H411
P260-P264-P271-P301+P312a-P304+P340-P501a

Mercuric bromide paper > RS - For analysis according to Ph. Eur. Chap. 4.1.1 — RS

Code	Size	Packaging	Notes
611052101	50 pc	Glass bottle	Ref Ph.Eur 1052101

Storage: in a glass-stoppered container wrapped with black paper

Mercury standard solution

• Mercurio standard soluzione • Mercure solution standard • Mercurio, solución patrón • Quecksilber-Standardlösung

Mercury standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2 — RS

Code	Size	Packaging	Notes
615001901	100 ml	Plastic bottle	A 10 ppm solution: to dilute according to Ref Ph.Eur 5001901
615001900	1 l	Plastic bottle	A 1000 ppm solution Ref Ph.Eur 5001900

Mercury standard solution > RS - Standard solution for ICP-MS — RS

Code	Size	Packaging	Notes
505652	100 ml	Plastic bottle	conc. 10 ppm Matrix: Nitric acid
505655	100 ml	Plastic bottle	conc. 100 ppm Matrix: Hydrochloric acid
506918	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505654	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Mercury standard solution > RS - Standard solution for ICP — RS

Code	Size	Packaging	Notes
503631	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503633	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503635	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503637	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Mercury standard solution > RS - Standard solution for AAS — RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
497555	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
497551	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503640	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507489	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Mercury standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
460741		Plastic ampoule	conc. 1.000 ppm Matrix: Nitric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Mercury standard solution > RS - Quality control standard solution for AAS (graphite furnace)

RS

Code	Size	Packaging	Notes
504370	100 ml	Plastic bottle	conc. 0.5 ppm - Matrix: 2% Nitric acid

Mercury (I) chloride

• di-Mercurio dicloruro • Mercure (I) chlorure • Mercure (I) chlorure • Quecksilber (I) chlorid

Synonym:
• *Calomel*
• *Mercurous chloride*

Hg_2Cl_2
Molecular Weight: 472,09
CAS: 10112-91-1
EEC-N: 233-307-5



Warning
H302-H315-H319-H335-H410
P261-P271-P304+P340-P305+P351+P338-
P332+P313-P403+P233

Mercury (I) chloride > RE - Pure

RE

Description ..... White crystalline powder Residue on ignition ..... ≤ 0.02 % Assay ..... ≥ 99.5 %
Identification ..... Positive Sulphate ..... ≤ 0.01 %

Code	Size	Packaging	Notes
352654	100 g	Glass bottle	
352657	1 kg	Plastic bottle	

Mercury chloride solution 54 g/l

• Mercurio cloruro soluzione 54 g/l • Mercure (II) chlorure solution • Mercurio (II) dicloruro solución 54 g/l
• Quecksilber (II) chloridlösung

Synonym:
Mercuric chloride

Classification transport
ONU: 2024
Transport Hazard class: 6.1
Packing group II



Danger
H300-H310-H314-H341-H361f-H372-H410
P280-P301+P310a-P301+P330+P331-
P303+P361+P353-P304+P340-P305+P351+P338

Mercury chloride solution 54 g/l > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611052201	100 ml	Plastic bottle	Ref Ph.Eur 1052201

Mercury (II) chloride

• Mercurio dicloruro • Mercure (II) chlorure • Mercurio (II) dicloruro • Quecksilber (II) chlorid

Synonym:
Mercuric chloride

$HgCl_2$
Molecular Weight: 271,5
CAS: 7487-94-7
EEC-N: 231-299-8

Classification transport
ONU: 1624
Transport Hazard class: 6.1
Packing group II



Danger
H300-H310-H314-H341-H361f-H372-H410
P280-P301+P310a-P301+P330+P331-
P303+P361+P353-P304+P340-P305+P351+P338

Mercury (II) chloride > RPE - For analysis - ISO - ACS - Reag.Ph.Eur. - Reag.USP

RPE

Description ..... Pezzi bianchi Residue after reduction..... ≤200 ppm Appearance of solution ..... Conform Loss on drying ..... ≤1.0 %
Identification ..... Positive Fe ..... ≤20 ppm Acidity or alkalinity..... Conform
Solution in ethyl ether..... Conform Assay (complexometric) 99.5 ÷ 100.5 % s.s. Mercurous chloride ..... Conform

Code	Size	Packaging	Notes
461003	50 g	Glass bottle	

Mercury (II) iodide

• Mercurio ioduro ico • Mercure (II) iodure • Mercurio (II) yoduro • Quecksilber (II) iodid

Synonym: *Mercuric iodide red*

HgI_2
Molecular Weight: 454,45
CAS: 7774-29-0
EEC-N: 231-873-8

Classification transport
ONU: 1638
Transport Hazard class: 6.1
Packing group II



Danger
H300-H310-H330-H373-H410
P271-P284-P301+P310a-P304+P340-P361+P364-P403+P233

Mercury (II) iodide > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Red - brick powder	Solub. in KI solution ..... Conform	Mercurous compounds ..... ≤0.1 %
Identification ..... Positive	Soluble mercury salts ..... ≤500 ppm	Assay (oxidimetric) ..... ≥99.0 % s.s.

Code	Size	Packaging	Notes
461105	250 g	Glass bottle	

Mercury (II) oxide red

• Mercurio ossido rosso • Mercure (II) oxyde rouge • Mercurio (II) óxido rojo • Quecksilber (II) rotes Oxid

Synonym: *Mercuric oxide*

HgO
Molecular Weight: 216,61
CAS: 21908-53-2
EEC-N: 244-654-7

Classification transport
ONU: 1641
Transport Hazard class: 6.1
Packing group II



Danger
H300-H310-H330-H373-H410
P271-P284-P301+P310a-P304+P340-P361+P364-P403+P233

Mercury (II) oxide red > RPE - For analysis - ACS

RPE

Description ..... Red powder	Chloride ..... ≤ 250 ppm	Residue after reduction ..... ≤ 0.025 %
Identification ..... Positive	Sulphate ..... ≤ 150 ppm	Fe ..... ≤50 ppm
Diluted HCl-ins. matter ..... ≤ 0.03 %	Nitrogen compounds (N) ..... ≤ 50 ppm	Assay (complexometric) ..... ≥ 99.0 %

Code	Size	Packaging	Notes
461325	250 g	Glass bottle	

Mercury (II) sulfate

• Mercurio solfato ico • Mercure (II) sulfate • Mercurio (II) sulfato • Quecksilber (II) sulfat

Synonym: *Mercuric sulfate*

$HgSO_4$
Molecular Weight: 296,65
CAS: 7783-35-9
EEC-N: 231-992-5

Classification transport
ONU: 1645
Transport Hazard class: 6.1
Packing group II



Danger
H300-H310-H330-H373-H410
P271-P284-P301+P310a-P304+P340-P361+P364-P403+P233

Mercury (II) sulfate > RPE - For analysis - ACS

RPE

Description ..... White powder or yellow	Chloride ..... ≤30 ppm	Mercurous compounds ..... ≤0.15 %	Fe ..... ≤50 ppm
Identification ..... Positive	Residue on ignition ..... ≤200 ppm	Nitrate ..... Conform	Assay (oxidimetric) ..... ≥98.0 %

Code	Size	Packaging	Notes
461405	250 g	Plastic bottle	

Mercury (II) sulfate solution

• Mercurio solfato soluzione • Mercure (II) sulfate solution • Mercurio (II) sulfato solución
• Quecksilber (II) sulfatlösung

Synonym: *Mercuric sulfate*

Classification transport
ONU: 2922
Transport Hazard class: 8
Packing group II



Danger
H290-H302-H331-H314-H373-H411
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Mercury (II) sulfate solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611052600	100 ml	Glass bottle	Ref Ph.Eur 1052600

Product specifications are subject to changes. Please visit our website for updates.

Metanil yellow

• Giallo metanile • Jaune de méthanyle • Amarillo de metanilo • Metanilgelb

Synonym:
- *3-(4-Anilinophenylazo)benzenesulfonic acid sodium salt*
- *Acid yellow 36*

$C_{18}H_{14}N_3NaO_3S$
Molecular Weight: 375,38
CAS: 587-98-4
EEC-N: 209-608-2

Warning
H312-H332
P261-P271-P280h-P304+P340-P312a-P501a

Metanil yellow > RPE - For analysis - C.I. 13065

RPE

Description ..... Golden to orange powder | Identification ..... Positive | pH range ..... 1.2 - 2.3 | Colour change..... red-yellow

Code	Size	Packaging	Notes
453542	25 g	Glass bottle	

Metaphosphoric acid

• Acido metafosforico • Acide métaphosphorique • Acido metafosfórico • Metaphosphorsäure

Synonym:
meta-Phosphoric acid

HPO_3
Molecular Weight: 79,98
CAS: 37267-86-0

Classification transport
ONU: 3261
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Metaphosphoric acid > RPE - For analysis

RPE

Description Whitish semitransparent pieces | Heavy metals (Pb)..... ≤ 50 ppm | As ..... ≤2 ppm | Stabilizer(NaPO3) ..... 50 ÷ 60 %
Identification ..... Positive | Nitrate ..... ≤10 ppm | Fe ..... ≤100 ppm
Chloride..... ≤20 ppm | Subst. reducing KMnO4..... ≤100 ppm (5m) | Assay (acidimetric) ..... 40 ÷ 50 %

Code	Size	Packaging	Notes
407465	250 g	Plastic bottle	
407467	1 kg	Plastic bottle	

Methanesulfonic acid

• Acido metansolfonico • Acide méthanesulfonique • Acido metanosulfónico • Methansulfonsäure

CH_3SO_3H
Molecular Weight: 96,11
CAS: 75-75-2
EEC-N: 200-898-6

Classification transport
ONU: 2586
Transport Hazard class: 8
Packing group III

Danger
H302-H312-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P362+P364

Methanesulfonic acid > RE - Pure

RE

Description ..... Colourless to yellow liquid | Identification ..... Positive | Density at 20° C ..... 1.47 ÷ 1.48 | Assay (acidimetric) ..... ≥99 %

Code	Size	Packaging	Notes
407481	250 ml	Glass bottle	
407483	1 l	Glass bottle	

Methanol

• Metanolo • Méthanol • Metanol • Methanol

Synonym:
Methyl alcohol

CH_3OH
Molecular Weight: 32,04
CAS: 67-56-1
EEC-N: 200-659-6

Classification transport
ONU: 1230
Transport Hazard class: 3
Packing group II

Danger
H225-H301-H311-H331-H370
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Methanol > RS - For UHPLC-MS

RS

Description ..... Clear colourless liquid
Colour ..... ≤ 5 APHA
Identification (I.R.) ..... Positive
Refractive index at 20°C .. 1.3270 - 1.3300
Residue on evaporation ..... ≤ 1 ppm
Acidity ..... ≤ 0.0003 meq/g
Alkalinity ..... ≤ 0.00006 meq/g
Assay (CPG) ..... ≥ 99.99 %
Water (K.F.) ..... ≤ 200 ppm
Transmittance
At 210 nm ..... ≥ 40 %
At 225 nm ..... ≥ 70 %
At 230 nm ..... ≥ 80 %
≥ 260 nm ..... ≥ 98 %
Fluorescence (quinine)
At 254 nm ..... ≤ 1 ppb
At 365 nm ..... ≤ 1 ppb
UHPLC gradient peak
At 220 nm ..... ≤ 4 mAU
At 235 nm ..... ≤ 2 mAU
Drift at 220 nm ..... ≤ 30 mAU
Drift at 235 nm ..... ≤ 10 mAU
Sensitive Impurities (reserpine) ..... ≤ 30 ppb
Metals compounds
Al ..... ≤ 20 ppb
Fe ..... ≤ 20 ppb
Na ..... ≤ 50 ppb
Ca ..... ≤ 50 ppb
Mg ..... ≤ 20 ppb
K ..... ≤ 50 ppb

Code	Size	Packaging	Notes
414941	1 l	Glass bottle	
414942	2.5 l	Glass bottle	

Methanol > RS - For LC/MS

RS

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Identification (I.R.) ..... Conform
Refractive index at 20°C ..... 1.327 ÷ 1.331
Water (K.F.) ..... ≤ 200 ppm
Residue on evaporation ..... ≤ 2 ppm
Acidity ..... ≤ 0.0003 meq/g
Alkalinity ..... ≤ 0.00006 meq/g
Assay (GLC) ..... ≥ 99.95 %
Transmission UV (1cm, ref water)
At 210 nm ..... ≥ 30 %
At 225 nm ..... ≥ 65 %
At 235 nm ..... ≥ 85 %
At 250 nm ..... ≥ 95 %
≥ 260 nm ..... ≥ 98 %
Fluorescence (quinine)
At 254 nm ..... ≤ 1 ppb
At 365 nm ..... ≤ 1 ppb
HPLC gradient
Test LC-MS TIC (50-2000m/z) ESI (+)
Sensitive Impurities (reserpine) ..... ≤ 50 ppb
Metals compounds
Al ..... ≤ 50 ppb
Fe ..... ≤ 50 ppb
Na ..... ≤ 50 ppb
Ca ..... ≤ 50 ppb
Mg ..... ≤ 50 ppb
K ..... ≤ 50 ppb

Code	Size	Packaging	Notes
414831	1 l	Glass bottle	
414832	2.5 l	Glass bottle	

Methanol > RS - For HPLC - GOLD - Ultragradient grade

RS

Description ..... Clear liquid
Identification ..... Positive
Colour (APHA) ..... ≤ 10
Density at 20°C ..... 0.7910 ÷ 0.7930
Refractive index at 20°C . 1.3270 ÷ 1.3300
Distillation range ..... 64.1 ÷ 65.1 °C
Water (K.F.) ..... ≤ 0.02 %
Residue on evaporation ..... ≤ 5 ppm
Carbonyle compounds (CH3COCH3) ... ≤ 20 ppm
Substances reducing KMnO4 (O) ... ≤ 2 ppm
Acidity ..... ≤ 0.0003 meq/g
Alcalinity ..... ≤ 0.00006 meq/g
Ethanol ..... ≤ 50 ppm
Assay (GLC) ..... ≥ 99.9 %
Fluorescence
at 254 nm ..... ≤ 1 ppb
at 365 nm ..... ≤ 1 ppb
Trasmittance
at 210 nm ..... ≥ 30 %
at 220 nm ..... ≥ 55 %
at 225 nm ..... ≥ 65 %
at 235 nm ..... ≥ 85 %
at 240 nm ..... ≥ 90 %
at 250 nm ..... ≥ 95 %
at 260 nm ..... ≥ 98 %
Functionality for HPLC
HPLC Gradient ..... Passed test

Code	Size	Packaging	Notes
412721	1 l	Glass bottle	
412722	2.5 l	Glass bottle	
412724	4 l	Glass bottle	
412725	5 l	Aluminium can	

Methanol > RS - For HPLC PLUS Gradient grade

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.7917 ÷ 0.7921
Refractive index at 20°C . 1.3278 ÷ 1.3298
Boiling point ..... 64.1 ÷ 65.1 ° C
Water (K.F.) ..... ≤0.02 %
Residue on evaporation ..... ≤5 ppm
Acidity ..... ≤0.0005 meq/g
Alcalinity ..... ≤0.0002 meq/g
Assay (GLC) ..... ≥99.9 %
Fluorescence
at 254 nm ..... ≤1 ppb
at 365 nm ..... ≤1 ppb
U.V. Transmittance
at 210 nm ..... ≥30 %
at 220 nm ..... ≥50 %
at 235 nm ..... ≥80 %
at 260 nm ..... ≥98 %
Carbonyl compounds (CO) ..... ≤ 20 ppm
Ethyl alcohol ..... ≤ 200 ppm
HPLC Gradient

Code	Size	Packaging	Notes
412381	1 l	Glass bottle	
412383	2.5 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Methanol > RS - For HPLC Isocratic grade - ACS - Reag.Ph.Eur R2 - Reag.USP

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.791 ÷ 0.793
Refractive index at 20°C. 1.3278 ÷ 1.3298
Boiling point ..... 64 ÷ 65 °C
Residue on evaporation ..... ≤ 5 ppm
Water (K.F.) ..... ≤ 0.05 %
Acidity ..... ≤ 0.0005 meq/g
Alcalinity ..... ≤ 0.0002 meq/g
Substances darkened by H2SO4 ... Conform
Subs. reducing KMnO4 ..... Conform
Solubility in water ..... Conform
Carbonyl compounds (CO) ..... ≤ 0.001 %
Absorbance ACS ..... Conform
HPLC Gradient (ACS) ..... Conform
Assay (GLC) ..... ≥ 99.9 %
U.V. Transmittance
at 210 nm ..... ≥ 20 %
at 235 nm ..... ≥ 80 %
at 260 nm ..... ≥ 98 %
Absorbance
At 225 nm ..... ≤ 0.17 AU

Code	Size	Packaging	Notes
412531	1 l	Glass bottle PVC coated	
412533	1 l	Glass bottle	
412532	2.5 l	Glass bottle	
412535	2.5 l	Glass bottle PVC coated	

Methanol > RS - For HPLC - Isocratic Grade

RS

Description ..... Clear colourless liquid
Density at 20°C ..... 0.791 ÷ 0.793
Refractive index at 20°C. 1.3278 ÷ 1.3298
Boiling point ..... 63.6 ÷ 65.6 °C
Free alkali (as NH3) ..... ≤ 1 mg/Kg
Assay GLC ..... ≥ 99.9 %
Free acid (as HCOOH) ..... ≤ 10 mg/Kg
Acidity (Formic Ac) ..... ≤20 ppm
Residue on evaporation ..... ≤ 10 ppm
Water (KF) ..... ≤ 500 ppm
UV transmittance at 210 nm ..... ≥ 20 %
Transmittance UV
Transmittance UV at 210 nm ..... ≥20 %
Transmittance UV at 235 nm ..... ≥ 80 %
Transmittance UV at 260 nm ..... ≥ 98 %

Code	Size	Packaging	Notes
525101	1 l	Glass bottle	
525102	2.5 l	Glass bottle	

Methanol > RS - For GC-MS

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C ..... 1.327 - 1.331
Water (K.F.) ..... ≤ 300 ppm
Residue on evaporation ..... ≤ 2 ppm
Colour ..... ≤ 5 APHA
Acidity (formic acid) ..... ≤ 10 ppm
Alcalinity (NH3) ..... ≤ 1 ppm
Assay (GC) ..... ≥ 99.98 %
GC-MS.Individual peak (n-hexadecane). ≤ 2 µg/L
Ret.range n-undecane to n-tetracontane (scanning area 30-600amu)

Code	Size	Packaging	Notes
414952	1 l	Glass bottle	

Methanol > RS - ATRASOL - For trace analysis, Suitable for Volatile chlorinated compounds analysis

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C ..... 1.327 - 1.331
Water content (K.F.) ..... ≤ 300 mg/Kg
Colour ..... ≤ 5 Hazen
Non volatile residue ..... ≤ 2 mg/Kg
Assay (GC) ..... ≥ 99.98 %
Free alkali (as NH3) ..... ≤ 1 mg/Kg
Free acid (as HCOOH) ..... ≤ 10 mg/Kg
GC-ECD.Individual peak (CCl4) ..... ≤ 1 µg/l
Ret.range dichloromethane to 1,2,4-trichlorobenzene
GC-ECD.Individual peak (Lindane) .≤ 2 ng/L
Ret.range 1,2,4-trichlorobenzene to decachlorobiphenyle
GC-FID.Individual peak (n-hexadecane). ≤ 2 µg/L
Ret.range n-undecane to n-tetracontane

Code	Size	Packaging	Notes
P0933216	1 l	Glass bottle	
P0933221	2.5 l	Glass bottle	

Methanol > RS - PESTIPUR - For pesticide analysis

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Colour ..... ≤ 10 hazen
Water ..... ≤ 0.05 %
Assay (GC) ..... ≥ 99.9 %
Free acids (HCOOH) ..... ≤ 10 ppm
Non volatile residue ..... ≤ 5 mg/Kg
Free alkalies (NH3) ..... ≤ 1 ppm
Not volatile residue ..... ≤ 5 ppm
GC-ECD (Lindano) ..... ≤ 3 ng/l
GC-NPD (Ethylparation) ..... ≤ 3 ng/l
Assay (GLC) ..... ≥ 99.9 %

Code	Size	Packaging	Notes
414930	1 l	Glass bottle	
414932	2.5 l	Glass bottle	

Methanol > RS - SPECTROSOL - For optical spectroscopy

RS

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Density at 20° C ..... 0.791 ÷ 0.793
Refractive index at 20°C . 1.3280 ÷ 1.3296
Boiling point ..... 64.1 ÷ 65.1 ° C
Water (K.F.) ..... ≤300 ppm
Residue on evaporation ..... ≤5 ppm
Acidity ..... ≤0.0005 meq/g
Alcalinity ..... ≤0.0002 meq/g
Free alkali (as NH3) ..... ≤ 1 mg/Kg
Assay (GLC) ..... ≥99.9 %
Fluorescence
at 254 nm ..... ≤2 ppb
at 365 nm ..... ≤2 ppb
U.V. Transmittance
Ethanol ..... ≤ 200 mg/Kg
at 205 nm ..... ≥10 %
Density d20/20 ..... 0.791 - 0.793
at 220 nm ..... ≥50 %
at 230 nm ..... ≥75 %
at 240 nm ..... ≥90 %
at 260 nm ..... ≥98 %
Nessler test ..... Conform
Ethyl alcohol ..... ≤ 200 ppm
Carbonyl compounds (CO) ..... ≤ 20 ppm

Code	Size	Packaging	Notes
414902	1 l	Glass bottle	
414903	2.5 l	Glass bottle	

Methanol > RS - Anhydrous - For analysis

RS

Description ..... Clear liquid
Identification (I.R.) ..... Positive
Colour (APHA) ..... ≤ 10
Refractive index at 20°C ..... 1.327 - 1.331
Acidity ..... ≤ 0.002 %
Water (K.F.) ..... ≤ 0.005 %
Residue on evaporation ..... ≤ 0.001 %
Assay (GLC) ..... ≥ 99.9 %

Code	Size	Packaging	Notes
P0931010	200 ml	Bottle with septum	
414981	1 l	Glass bottle	Water content max 50 ppm
P0931016	1 l	Glass bottle	
P0931021	2.5 l	Glass bottle	

Methanol > RS - VLSI - For electronic use

RS

Code	Size	Packaging	Notes
527641	1 l	Plastic bottle	
527640	2.5 l	Plastic bottle	

For specifications, contact our customer service for a certificate of analysis

Methanol > RS - RSE - For electronic use

RS

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Water miscibility ..... Conform
Ready carbonizable substances ..... Conform
Assay (GLC) ..... ≥99.9 %
Resistivity ..... ≥0.5 Mohm.cm
Density at 20° C ..... 0.791 ÷ 0.793
Boiling point ..... 64.1 ÷ 65.1 ° C
Water (K.F.) ..... ≤500 ppm
Residue on evaporation ..... ≤10 ppm
Acidity (formic acid) ..... ≤15 ppm
Alcalinity (NH3) ..... ≤1 ppm
Ethyl alcohol ..... ≤200 ppm
Chloride ..... ≤0.2 ppm
Carbonyl Compounds (CO) ..... ≤5 ppm
Phosphate ..... ≤0.5 ppm
Heavy metals (Pb) ..... ≤0.2 ppm
Subst. reducing KMnO4 ..... ≤2.5 ppm
Total sulphur ..... ≤1 ppm
Ag ..... ≤0.2 ppm
Al ..... ≤0.05 ppm
As ..... ≤0.01 ppm
Au ..... ≤0.05 ppm
B ..... ≤0.01 ppm
Ba ..... ≤0.1 ppm
Be ..... ≤0.02 ppm
Bi ..... ≤0.02 ppm
Ca ..... ≤0.5 ppm
Cd ..... ≤0.005 ppm
Co ..... ≤0.01 ppm
Cr ..... ≤0.01 ppm
Cu ..... ≤0.01 ppm
Fe ..... ≤0.1 ppm
Ga ..... ≤0.02 ppm
In ..... ≤0.02 ppm
K ..... ≤0.1 ppm
Li ..... ≤0.02 ppm
Mg ..... ≤0.1 ppm
Mn ..... ≤0.01 ppm
Mo ..... ≤0.05 ppm
Na ..... ≤0.5 ppm
Ni ..... ≤0.01 ppm
Pb ..... ≤0.01 ppm
Pt ..... ≤0.05 ppm
Sb ..... ≤0.01 ppm
Sn ..... ≤0.02 ppm
Sr ..... ≤0.02 ppm
Ti ..... ≤0.05 ppm
Tl ..... ≤0.05 ppm
V ..... ≤0.05 ppm
Zn ..... ≤0.01 ppm
Zr ..... ≤0.05 ppm

Code	Size	Packaging	Notes
414917	1 l	Glass bottle	
414914	2.5 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Methanol > RS - MOS - For electronic use

RS

Description ..... Clear liquid	Chloride ..... ≤0.2 ppm	Ca ..... ≤0.5 ppm	Ni ..... ≤0.01 ppm
Colour (APHA) ..... ≤10	Carbonyl Compounds (CO) ..... ≤5 ppm	Cd ..... ≤0.005 ppm	Pb ..... ≤0.01 ppm
Identification ..... Positive	Phosphate ..... ≤0.5 ppm	Co ..... ≤0.01 ppm	Pt ..... ≤0.05 ppm
Water miscibility ..... Conform	Heavy metals (Pb) ..... ≤0.2 ppm	Cr ..... ≤0.01 ppm	Sb ..... ≤0.01 ppm
Ready carbonizable substances ..... Conform	Subst. reducing KMnO4 ..... ≤2.5 ppm	Cu ..... ≤0.01 ppm	Sn ..... ≤0.02 ppm
Assay (GLC) ..... ≥99.9 %	Total sulphur ..... ≤1 ppm	Fe ..... ≤0.1 ppm	Sr ..... ≤0.02 ppm
Resistivity ..... ≥0.5 Mohm.cm	Ag ..... ≤0.2 ppm	Ga ..... ≤0.02 ppm	Ti ..... ≤0.05 ppm
Density at 20° C ..... 0.791 ÷ 0.793	Al ..... ≤0.05 ppm	In ..... ≤0.02 ppm	Tl ..... ≤0.05 ppm
Boiling point ..... 64.1 ÷ 65.1 ° C	As ..... ≤0.01 ppm	K ..... ≤0.1 ppm	V ..... ≤0.05 ppm
Water (K.F.) ..... ≤0.05 %	Au ..... ≤0.05 ppm	Li ..... ≤0.02 ppm	Zn ..... ≤0.01 ppm
Residue on evaporation ..... ≤10 ppm	B ..... ≤0.01 ppm	Mg ..... ≤0.1 ppm	Zr ..... ≤0.05 ppm
Acidity (formic acid) ..... ≤15 ppm	Ba ..... ≤0.1 ppm	Mn ..... ≤0.01 ppm	
Alcalinity (NH3) ..... ≤1 ppm	Be ..... ≤0.02 ppm	Mo ..... ≤0.05 ppm	
Ethyl alcohol ..... ≤200 ppm	Bi ..... ≤0.02 ppm	Na ..... ≤0.5 ppm	

Code	Size	Packaging	Notes
414822	1 l	Glass bottle	
414821	2.5 l	Glass bottle	

Methanol > RS - For titration according to Karl Fischer

RS

Description ..... Clear colourless liquid	Density at 20° C ..... 0.791 ÷ 0.793	Assay (GLC) ..... ≥99.9 %
Identification ..... Positive	Water (K.F.) ..... ≤0.03 %	

Code	Size	Packaging	Notes
414881	1 l	Glass bottle	
414883	2.5 l	Glass bottle	

Methanol > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Clear liquid	Residue on evaporation ..... ≤ 8 ppm	B ..... ≤0.02 ppm	Mn ..... ≤0.01 ppm
Colour (APHA) ..... ≤10	Acidity ..... ≤ 0.0003 meq/g	Ba ..... ≤0.1 ppm	Ni ..... ≤0.01 ppm
Identification (I.R.) ..... Conform	Alkalinity ..... ≤ 0.00006 meq/g	Ca ..... ≤0.5 ppm	Pb ..... ≤0.01 ppm
Water miscib. (15:40) ..... Complete	Ethyl alcohol ..... ≤ 50 ppm	Cd ..... ≤0.05 ppm	Sn ..... ≤0.1 ppm
Ready carbonizable substances ..... Conform	Chloride ..... ≤0.5 ppm	Co ..... ≤0.01 ppm	Zn ..... ≤0.1 ppm
Density at 20° C ..... 0.791 ÷ 0.793	Carbonyl Compounds (CO) ..... ≤10 ppm	Cr ..... ≤0.02 ppm	Assay (GLC) ..... ≥99.9 %
Refractive index at 20°C 1.3280 ÷ 1.3296	Heavy metals (Pb) ..... ≤0.5 ppm	Cu ..... ≤0.01 ppm	Nessler test ..... Conform
Boiling point ..... 64 ÷ 65 ° C	Subst. reducing KMnO4 ..... ≤2 ppm	Fe ..... ≤0.1 ppm	
Water (K.F.) ..... ≤300 ppm	Al ..... ≤0.5 ppm	Mg ..... ≤0.1 ppm	

Code	Size	Packaging	Notes
414814	1 l	Glass bottle	
414819	1 l	Plastic bottle	
414815	2.5 l	Plastic bottle	
414816	2.5 l	Glass bottle	
524102	5 l	Plastic tank	
524103	5 l	Metal tank	
414818	10 l	Plastic tank	
414813	25 l	Plastic tank	
414817	200 l	Plastic drum	

Methanol > RPE - Anhydrous - For analysis

RPE

Description ..... Clear liquid	Density at 20° C ..... 0.791 ÷ 0.793	Acidity (formic acid) ..... ≤15 ppm	Subst. reducing KMnO4 ..... ≤3 ppm
Colour (APHA) ..... ≤10	Refractive index at 20°C 1.3280 ÷ 1.3297	Alcalinity (NH3) ..... ≤1 ppm	Cu ..... ≤0.01 ppm
Identification ..... Positive	Boiling point ..... 64.1 ÷ 65.1 ° C	Ethyl alcohol ..... ≤200 ppm	Ni ..... ≤0.01 ppm
Water miscib. (15:40) ..... Complete	Water (K.F.) ..... ≤0.01 %	Carbonyl Compounds (CO) ..... ≤5 ppm	Assay (GLC) ..... ≥99.9 %
Ready carbonizable substances ..... Conform	Residue on evaporation ..... ≤10 ppm	Heavy metals (Pb) ..... ≤0.5 ppm	Nessler test ..... Conform

Code	Size	Packaging	Notes
414854	1 l	Glass bottle	
414855	2.5 l	Glass bottle	

Methanol > ERBApharm - According to pharmacopoeia: DAB-NF-Ph.Eur.

ERBApharm

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Identification ..... Positive
Density at 20° C ..... 0.791 ÷ 0.793
Refractive index at 20°C ..... 1.328 ÷ 1.330
Acidity ..... Conform DAB
Acidity or alcalinity ..... Pass test Ph.Eur.
Acidity (formic acid) ..... ≤ 10 ppm
Alcalinity (NH3) ..... ≤1 ppm

Distillation range ..... 64 ÷ 65 ° C
Water (K.F.) ..... ≤ 300 ppm
Residue on evaporation ..... ≤10 ppm
Ready oxidizable substances ..... Conform USP-NF
Ready carbonizable substances ..... Conform USP-NF
Acetone ..... ≤10 ppm
Acetone and aldeydes ..... ≤ 20 ppm

Ethyl alcohol ..... ≤ 100 ppm
Heavy metals and Zn ..... ≤2 ppm
Fe ..... ≤1 ppm
Absorbance UV (1cm, ref. water)
At 230 nm ..... ≤ 0.15 AU
At 250 nm ..... ≤ 0.05 AU
At 270 nm ..... ≤ 0.02 AU
At 290 nm ..... ≤ 0.01 AU
Absorbance UV curve ..... Smooth Ph.Eur.

Benzene ..... ≤ 2 ppm(v/v)
Organic volatile impurities Conform USP-NF
Related substances (CPG) ..... Pass test
Assay (GLC) ..... ≥99.9 %
Origin (BSE/TSE) ..... Synthesis
Residual solvents (Current ICH) ..... Conform

Code	Size	Packaging	Notes
309204	1 l	Glass bottle	
309203	2.5 l	Glass bottle	
309201	25 l	Plastic tank	
529100	200 l	Plastic drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Methanol > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Colour ..... ≤10 APHA

Density at 20° C ..... 0.791 ÷ 0.793
Refractive index at 20°C. 1.3278 ÷ 1.3298
Boiling point ..... 64.1 ÷ 65.1 ° C

Residue on evaporation ..... ≤50 ppm
Acidity (formic acid) ..... ≤20 ppm
Water (K.F.) ..... ≤0.05 %

Subst. reducing KMnO4 ..... ≤3 ppm
Assay (GLC) ..... ≥99.9 %

Code	Size	Packaging	Notes
309004	1 l	Glass bottle	
309001	2.5 l	Glass bottle	
528101	5 l	Plastic tank	
309008	10 l	Plastic tank	
309002	25 l	Plastic tank	
309009	160 kg	Plastic drum	
309000	200 l	Plastic drum	

Methanol + 0.1% v/v formic acid

• Metanolo + 0.1% v/v acido formico • Méthanol + 0.1% v/v acide formique • Metanol + 0.1% v/v acido formico • Methanol + 0.1% v/v Ameisensäure

CH_3OH
Molecular Weight: 32,04
CAS: 67-56-1

Classification transport
ONU: 1992
Transport Hazard class: 3
Packing group II



Danger
H225-H301-H370
P210-P241-P280-P301+P310a-P303+P361+P353-P403+P235

Methanol + 0.1% v/v formic acid > RS - For LC/MS

RS

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
HPLC Gradient
At 254 nm ..... ≤10 mAU

Transmittance
At 230 nm ..... ≥10 %
Assay (CPG) ..... ≥99.5 %
Test LC-MS TIC (100-2000m/z)

Sensitive Impurities (reserpine) ..... ≤100 ppb
Metals content
Na ..... ≤2 ppm
Ca ..... ≤0.5 ppm

Mg ..... ≤0.5 ppm
K ..... ≤0.5 ppm
Raw material used

Code	Size	Packaging	Notes
414861	1 l	Glass bottle	
414862	2.5 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Methanol + 0.1% v/v trifluoroacetic acid

• Metanolo + 0.1% v/v acido trifluoroacetico • Méthanol + 0.1% v/v acide trifluoroacétique • Metanol + 0.1% v/v acido trifluoroacético • Methanol + 0.1% v/v trifluoressigsäure

CH_3OH
Molecular Weight: 32,04
CAS: 67-56-1

Classification transport
ONU: 1230
Transport Hazard class: 3
Packing group II



Danger
H225-H301-H370
P210-P241-P280-P301+P310a-P303+P361+P353-P403+P235

Methanol + 0.1% v/v trifluoroacetic acid > RS - For LC/MS

RS

Description ..... Clear colourless liquid	At 225 nm ..... ≥ 20 %	At 365 nm ..... ≤ 1 ppb	Al ..... ≤ 30 ppb
Assay (GC) (without TFA) ..... ≥ 99.9 %	At 240 nm ..... ≥ 50 %	HPLC gradient	Fe ..... ≤ 100 ppb
Trifluoroacetic acid content (V/V) 0.095 - 0.105 %	At 250 nm ..... ≥ 80 %	Drift at 254 nm ..... ≤ 40 mAU	Na ..... ≤ 50 ppb
Water (K.F.) ..... ≤ 500 ppm	At 260 nm ..... ≥ 95 %	Test LC-MS TIC (50-2000m/z) ES I(+)	Ca ..... ≤ 50 ppb
Residue on evaporation ..... ≤ 2 ppm	Fluorescence (quinine)	Sensitive Impurities (reserpine) ..... ≤ 50 ppb	Mg ..... ≤ 30 ppb
UV transmittance (1 cm, ref. water)	At 254 nm ..... ≤ 1 ppb	Metals compounds	K ..... ≤ 50 ppb

Code	Size	Packaging	Notes
414871	1 l	Glass bottle	
414872	2.5 l	Glass bottle	

Methanol, hydrochloric

• Metanolo, cloridrico • Méthanol - chlorhydrique • Metanol, clorhídrico • Methanol - Salzsäure

CH_3OH
Molecular Weight: 32,04
CAS: 67-56-1

Classification transport
ONU: 2929

Methanol, hydrochloric > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611053203	100 ml	Glass bottle	Ref Ph.Eur 1053203

Methanol-d4

• Alcole metilico-d4 • Méthanol-d4 • Metanol-d4 • Methanol-d4

Synonym:
Tetradeuteromethanol

CD_3OD
Molecular Weight: 36,07
CAS: 811-98-3
EEC-N: 212-378-6

Classification transport
ONU: 1230
Transport Hazard class: 3
Packing group II



Danger
H225-H301-H311-H331-H370
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Methanol-d4 > RS - For NMR - min 99.96%

RS

Code	Size	Packaging	Notes
P5310	10 x 0.6 ml	Glass ampoule	
P5319	10 x 0.75 ml	Glass ampoule	

For specifications, contact our customer service for a certificate of analysis

Methanol-d4 > RS - For NMR - min 99.8%

RS

Code	Size	Packaging	Notes
P5280	10 x 0.6 ml	Glass ampoule	
P5289	10 x 0.75 ml	Glass ampoule	
P5283A	5 ml	Glass ampoule	
P5284	5 x 10 ml	Glass ampoule	
P5284S	5 x 10 ml	Bottle with septum	
P5285	25 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Methanol-d4 + 0.03% TMS

• Alcole metilico-d4 + 0.03% TMS • Méthanol-d4 + 0.03% TMS • Metanol-d4 + 0.03% TMS • Methanol-d4 + 0.03% TMS

Synonym:
Tetradeuteromethanol

CD_3OD
Molecular Weight: 36,07
CAS: 811-98-3
EEC-N: 212-378-6

Classification transport
ONU: 1230
Transport Hazard class: 3
Packing group II



Danger
H225-H301-H311-H331-H370
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Methanol-d4 + 0.03% TMS > RS - For NMR - min 99.8%

RS

Code	Size	Packaging	Notes
P5140	10 x 0.6 ml	Glass ampoule	

For specifications, contact our customer service for a certificate of analysis

Methanol-d3

• Alcole metilico-d3 • Méthanol-d3 • Metanol-d3 • Methanol-d3

Synonym:
1,1,1-Trideuteromethanol

CD_3OH
Molecular Weight: 35,02
CAS: 1849-29-2
EEC-N: 217-435-9

Classification transport
ONU: 1230
Transport Hazard class: 3
Packing group II



Danger
H225-H301-H311-H331-H370
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Methanol-d3 > RS - For NMR - min 99.5%

RS

Code	Size	Packaging	Notes
P5309	10 x 0.75 ml	Glass ampoule	

For specifications, contact our customer service for a certificate of analysis

Methanol-d1

• Alcole metilico-d1 • Méthanol-d1 • Metanol-d1 • Methanol-d1

Synonym:
• *Methan(ol-d)*
• *Methyl alcohol-OD*

CH_3OD
Molecular Weight: 33,05
CAS: 1455-13-6
EEC-N: 215-933-0

Classification transport
ONU: 1230
Transport Hazard class: 3
Packing group II



Danger
H225-H301-H311-H331-H370
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Methanol-d1 > RS - For NMR - min 99.5%

RS

Code	Size	Packaging	Notes
P5275	25 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

2-MeTHF ▶ 2-Methyltetrahydrofuran

DL-Methionine

• DL-Metionina • DL-Méthionine • DL-Metionina • DL-Methionin

Synonym:
• *(±)-2-Amino-4-(methylmercapto)butyric acid*
• *DL-2-Amino-4(methylthio)butanoic acid*

$CH_3S(CH_2)_2CHNH_2COOH$
Molecular Weight: 149,21
CAS: 59-51-8
EEC-N: 200-432-1

DL-Methionine > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
pH (sol 2% in H2O) ..... 5.4 - 6.1
Loss on drying ..... ≤ 0.5 %
Ammonium ..... ≤1000 ppm
Chloride ..... ≤200 ppm
Heavy metals (Pb) ..... ≤10 ppm
Residue on ignition ..... ≤0.1 %
Sulphate ..... ≤200 ppm
Assay (non-aqueos medium) 99.0 ÷ 101.0 % (s.s.)

Code	Size	Packaging	Notes
463126	250 g	Plastic bottle	

DASIT GROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

Methone ▶ Dimedone

4-Methoxybenzaldehyde ▶ Anisaldehyde

4-Methoxybenzoic acid ▶ Anisic acid

2-Methoxy ethanol

• 2-Metossietanolo • Glycol éthylénique monométhylether • 2-Metoxietanol • 2-Methoxy-ethanol

Synonym:
Methyl glycol

$CH_2OHCH_2OCH_3$
Molecular Weight: 76,1
CAS: 109-86-4
EEC-N: 203-713-7

Classification transport
ONU: 1188
Transport Hazard class: 3
Packing group III

Danger
H226-H302-H312-H332-H360FD-HA26
P210-P241-P261-P280-P303+P361+P353-P304+P340

2-Methoxy ethanol > RPE - For analysis

RPE

Description	Clear colourless liquid	Water (K.F.)	≤0.1 %	Ca	≤0.5 ppm	Ni	≤0.02 ppm
Identification (I.R.)	Conform	Residue on evaporation	≤20 ppm	Cd	≤0.05 ppm	Pb	≤0.1 ppm
Water miscibility	Conform	Acidity (acetic acid)	≤50 ppm	Co	≤0.02 ppm	Sn	≤0.1 ppm
Benzene miscibility	Complete	Alcalinity (NH3)	≤0.85 ppm	Cr	≤0.02 ppm	Zn	≤0.1 ppm
Diethyl ether miscib.	Complete	Heavy metals (Pb)	≤2 ppm	Cu	≤0.02 ppm	Assay (GLC)	≥99.5 %
Density at 20° C	0.962 ÷ 0.968	Peroxides (H2O2)	≤10 ppm	Fe	≤1 ppm		
Refractive index at 20°C	1.4004 ÷ 1.4044	Al	≤0.5 ppm	Mg	≤0.01 ppm		
Boiling point	123.5 ÷ 124.5 ° C	Ba	≤0.1 ppm	Mn	≤0.02 ppm		

Code	Size	Packaging	Notes
454021	1 l	Glass bottle	
454024	2.5 l	Glass bottle	
454023	25 kg	Metal drum	
454028	180 kg	Metal drum	

alpha-Methoxyphenylacetic acid

• Acido alfa-metossifenilacetico • Acide alpha-méthoxyphénylacétique • Acido alfa-metoxifenilacético
• alpha-Methoxy(phenyl)essigsäure

Synonym:
• *MOPA*
• *alpha-Methyl-DL-mandilic acid*

$C_6H_5CH(OCH_3)COOH$
Molecular Weight: 166,18
CAS: 7021-09-2
EEC-N: 230-300-9

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

alpha-Methoxyphenylacetic acid > RPE - For analysis

RPE

Description	White powder	Identification	Positive	Melting point	69 ÷ 71 ° C	Assay	≥99 %

Code	Size	Packaging	Notes
407441	5 g	Glass bottle	

Methyl acetate

• Metile acetato • Méthyle acétate • Metilo acetato • Methylacetat

CH_3COOCH_3
Molecular Weight: 74,08
CAS: 79-20-9
EEC-N: 201-185-2

Classification transport
ONU: 1231
Transport Hazard class: 3
Packing group II

Danger
H225-H319-H336-HEU066
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Methyl acetate > RS - For HPLC - Isocratic Grade

RS

Refractive index at 20°C	1.359 - 1.363	Colour	≤ 10 Hazen	UV Absorbance at 255 nm	≤ 1 AU
Water content (K.F.)	≤ 500 mg/Kg	Assay (GC)	≥ 99.5 %	UV Absorbance at 275 nm	≤ 0.1 AU
Non volatile residue	≤ 10 mg/Kg	Free acid (as CH3COOH)	≤ 25 mg/Kg	UV Absorbance at 300 nm	≤ 0.01 AU

Code	Size	Packaging	Notes
P0043721	2.5 l	Glass bottle	

Methyl acetate > RPE - For analysis

RPE

Description ..... Clear liquid
Colour (APHA) ..... ≤ 10
Identification ..... Positive
Density at 20° C ..... 0.930 ÷ 0.936
Water (K.F.) ..... ≤ 0.1 %
Assay (GLC) ..... ≥ 99 %

Code	Size	Packaging	Notes
462017	1 l	Glass bottle	

Methyl acetate > RE - Pure

RE

Refractive index at 20°C ..... 1.359 - 1.363
Water content (K.F.) ..... ≤ 300 mg/Kg
Non volatile residue ..... ≤ 50 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 99 %
Methanol ..... ≤ 0.1 %
Free acid (as CH3COOH) ..... ≤ 50 mg/Kg

Code	Size	Packaging	Notes
P0040228	5 l	Plastic tank	
P0040240	10 l	Metal tank	
P0040248	25 l	Metal drum	
P0040268	200 l	Metal drum	

Methyl alcohol ▶ Methanol

4-Methylaminophenol sulfate

• Bis (4-idrossi-N-metilanilinio) solfato • 4-Méthylaminophénol sulfate • 4-Metil aminofenol sulfato • 4-Methylaminophenolsulfat

Synonym: *Metol*

$(CH_3NHC_6H_4OH)_2.H_2SO_4$
Molecular Weight: 344,39
CAS: 55-55-0
EEC-N: 200-237-1

Classification transport
ONU: 3077
Transport Hazard class: 9
Packing group III



Warning
H302-H317-H373-H410
P260-P264-P280g-P301+P312a-P333+P313-P501a

4-Methylaminophenol sulfate > RPE - For analysis - ACS

RPE

Description ..... beige crystalline powder
Identification ..... Positive
Suit.for phosphate det. ..... Conform
Residue on ignition ..... ≤ 0.15 %
Assay (oxidimetric) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
461805	250 g	Plastic bottle	

Methyl benzoate

• Metile benzoato • Méthyle benzoate • Metilo benzoato • Methylbenzoat

$C_6H_5COOCH_3$
Molecular Weight: 136,15
CAS: 93-58-3
EEC-N: 202-259-7



Warning
H302
P264-P270-P301+P312a-P330-P501a

Methyl benzoate > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 1.086 ÷ 1.090
Refractive index at 20° C ..... 1.51 ÷ 1.52
Water (K.F.) ..... ≤0.1 %
Residue on evaporation ..... ≤100 ppm
Acidity (benzoic acid) ..... ≤0.12 %
Cu ..... ≤ 0.2 ppm
Fe ..... ≤ 0.5 ppm
Ni ..... ≤ 0.2 ppm
Pb ..... ≤ 0.2 ppm
Assay (GLC) ..... ≥99 %

Code	Size	Packaging	Notes
462207	1 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

3-Methyl-2-benzothiazolinone hydrazone hydrochloride

• 3-Metile-2-benzotiazolinone idrazone cloridrato • 3-Méthyl-2-benzothiazoline hydrazone chlorhydrate
• 3-Metil-2-benzotiazolinona hidrazona clorhidrato • 3-Methyl-2-benzothiazolinhydrazonhydrochlorid

Synonym:
2-Hydrazono-3-methylbenzothiazoline hydrochloride

$C_8H_9N_3S.HCl.H_2O$
Molecular Weight: 233,72
CAS: 38894-11-0

Classification transport
ONU: 2811
Transport Hazard class: 6.1
Packing group III

Danger
H301-H315-H319-H351-H335
P261-P271-P280-P304+P340-P305+P351+P338-P403+P233

3-Methyl-2-benzothiazolinone hydrazone hydrochloride > RPE - For analysis

RPE

Description ..... White powder Identification ..... Positive Loss on drying ..... ≤ 8 % Assay (HPLC) ..... ≥ 97.5 % (d.s.)

Code	Size	Packaging	Notes
462238	5 g	Glass bottle	

Methyl blue

• Blu metile • Bleu de méthyle • Azul de metilo • Methylblau

Synonym:
Acid blue 93

$C_{37}H_{27}N_3Na_2O_9S_3$
Molecular Weight: 799,8
CAS: 28983-56-4
EEC-N: 249-352-9

Methyl blue > RS - For microscopy - C.I. 42780

RS

Description ..... Red-violet crystals Identification ..... Positive Sensib.(pH 9.0-11.0) ..... Conform

Code	Size	Packaging	Notes
428932	25 g	Glass bottle	

Dye for histology and microbiology

2-Methyl-2-butanol ▶ tert-Amyl alcohol

3-Methyl-1-butanol ▶ Isoamyl alcohol

2-Methylbutane ▶ Isopentane

Methyl cyanide ▶ Acetonitrile

Methylcyclohexane

• Metilcicloesano • Méthylcyclohexane • Metilciclohexano • Methylcyclohexan

Synonym:
Hexahydrotoluene

$CH_3CH(CH_2)_4CH_2$
Molecular Weight: 98,19
CAS: 108-87-2
EEC-N: 203-624-3

Classification transport
ONU: 2296
Transport Hazard class: 3
Packing group II

Danger
H225-H315-H336-H304-H411
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Methylcyclohexane > RS - SPECTROSOL - For optical spectroscopy

RS

Code	Size	Packaging	Notes
P0582716	1 l	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Methylcyclohexane > RS - Anhydrous - For analysis

RS

Refractive index at 20°C ..... 1.421 - 1.425 Non volatile residue ..... ≤ 10 mg/Kg Aromatic compounds ..... ≤ 300 mg/Kg Total sulphur (S) ..... ≤ 1 ppm
Water content (K.F.) ..... ≤ 50 mg/Kg Colour ..... ≤ 10 Hazen Assay (GC) ..... ≥ 99 %

Code	Size	Packaging	Notes
P0581016	1 l	Glass bottle	

Methylcyclohexane > RPE - For analysis

RPE

Refractive index at 20°C ..... 1.421 - 1.425
Water content (K.F.) ..... ≤ 150 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
Colour ..... ≤ 10 Hazen
Aromatic compounds ..... ≤ 300 mg/Kg
Assay (GC) ..... ≥ 99 %
Total sulphur (S) ..... ≤ 1 ppm

Code	Size	Packaging	Notes
P0580516	1 l	Glass bottle	

Methylcyclohexane > RE - Pure

RE

Description ..... Clear liquid
Density at 20°C ..... 0.767 ÷ 0.770
Refractive index at 20°C ..... 1.419 ÷ 1.427
Boiling point ..... 100.15 ÷ 100.65 °C
Benzene ..... ≤ 200 ppm
Total sulphur (S) ..... ≤ 2 ppm
Water (K.F.) ..... ≤ 0.03 %
Residue on evaporation ..... ≤ 50 ppm
Assay (GLC) ..... ≥ 99.0 % (GLC)
Acidity or alcalinity ..... Passes test

Code	Size	Packaging	Notes
528264	1 l	Glass bottle	
528261	5 l	Plastic tank	
528260	25 l	Metal drum	
528262	200 l	Metal drum	

Methylene blue

• Blu metilene • Bleu de méthylène • Azul de metileno • Methylenblau

Synonym:
• *Tetramethylthionine chloride*
• *3,7-bis(Dimethylamino)phenazathonium chloride*

$C_{16}H_{18}N_3SCl.3H_2O$
Molecular Weight: 373,9
CAS: 7220-79-3
EEC-N: 200-515-2

Warning
H302-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Methylene blue > RS - For analysis - Reag.Ph.Eur. - C.I. 52015

RS

Description ..... Dark green powder
Identification ..... Positive
Assay (dried base) ..... 96 ÷ 101 %
Solubility ..... Clear Blue Solution
Fe ..... Passes
Pb ..... ≤ 20 ppm
Loss on drying at 110°C ..... 10 -20 %

Code	Size	Packaging	Notes
428984	100 g	Plastic bottle	

Methylene blue > RPE - For analysis - C.I. 52015

RPE

Description ..... Polvere verde blu
Identification ..... Positive
Functionality ..... Conform

Code	Size	Packaging	Notes
429982	25 g	Glass bottle	
429981	500 g	Plastic bottle	

Redox indicator purple - colorless

Methylene blue saturated solution

• Blu metilene soluzione satura in alcole etilico • Bleu de méthylène solution saturée
• Azul de metileno solución saturada en alcohol etilico • Methylenblau gesättigte Lösung

Synonym:
Tetramethylthionine chloride

$C_{16}H_{18}N_3SCl$
Molecular Weight: 319,85 (an,)
CAS: 61-73-4

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II

Danger
H225-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

Methylene blue saturated solution > RPE - For analysis

RPE

Description ..... Blue liquid
Identification ..... Positive
Density at 20° C ..... 0.809 ÷ 0.815 g/ml
Assay ..... 1.0 ÷ 1.2 % (p/v)

Code	Size	Packaging	Notes
E429031	250 ml	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Methylene blue solution 1%

• Blu metilene soluzione 1% • Bleu de méthylène solution 1% • Azul de metileno solución 1%
• Methylenblaulösung 1%

Synonym:
Tetramethylthionine chloride

$C_{16}H_{18}N_3SCl$
Molecular Weight: 319,85 (an,)
CAS: 61-73-4

HEU210

Methylene blue solution 1% > RPE - For analysis

RPE

Description ..... Blue liquid Identification ..... Positive Density at 20° C ..... ~ 1.00 g/ml Assay (oxidimetric) ..... 0.9 ÷ 1.1 %

Code	Size	Packaging	Notes
E429011	500 ml	Plastic bottle	

Methylene chloride ▶ Dichloromethane

Methyl glycol ▶ 2-Methoxy ethanol

Methyl green

• Verde metile • Vert de méthyle • Verde de metilo • Methylgrün

$C_{26}H_{33}Cl_2N_3$
Molecular Weight: 458,48
CAS: 22383-16-0

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Methyl green > RS - For microscopy - C.I. 42585

RS

Description ..... Dark red powder Identification ..... Positive

Code	Size	Packaging	Notes
491351	10 g	Glass bottle	
491352	25 g	Glass bottle	

Dye for cytology

Methyl 4-hydroxybenzoate

• Metile p-ossibenzoato • Méthyle p-oxybenzoate • Metilo p-oxibenzoato • Methyl-4-hydroxybenzoat

Synonym:
4-Hydroxybenzoic acid propyl ester

$HOC_6H_4COOCH_3$
Molecular Weight: 152,15
CAS: 99-76-3
EEC-N: 202-785-7

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Methyl 4-hydroxybenzoate > ERBApharm - According to pharmacopoeia: DAB-BP-FU-NF-Ph.Eur.-Ph.Franc.

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution ..... Conform Ph.Eur.
Acidity (HCl) ..... Conform Ph.Eur.
Related substances ..... Conform Ph.Eur.
Organic volatile impurities Conform USP-NF
Melting point ..... 125 ÷ 128 °C
Sulphated ash ..... ≤ 0.1 %
Assay (acidimetric) ..... 99.0 ÷ 100.5 %
Origin (BSE/TSE) ..... Synthesis

Code	Size	Packaging	Notes
354007	1 kg	Plastic bottle	
354008	5 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Methyl iodide

• Metile ioduro • Méthyle iodure • Metilo yoduro • Methyljodid

Synonym:
Iodomethane

CH_3I
Molecular Weight: 141,94
CAS: 74-88-4
EEC-N: 200-819-5

Classification transport
ONU: 2644
Transport Hazard class: 6.1
Packing group I



Danger
H301-H312-H331-H315-H351-H335
P261-P280-P304+P340-P308+P313-P330-P362+P364-P403+P233

Methyl iodide > RPE - For analysis

RPE

Description ..... Clear liquid	Density at 20° C ..... 2.276 ÷ 2.284	Water (K.F.) ..... ≤ 500 ppm
Identification ..... Positive	Refractive index at 20°C ..... 1.5273 ÷ 1.5313	Residue on evaporation ..... ≤50 ppm
Colour ..... ≤150 APHA	Boiling point ..... 42.0 ÷ 43.0 °C	Assay (GLC) ..... ≥99 %

Code	Size	Packaging	Notes
462601	50 ml	Glass bottle	
462604	250 ml	Glass bottle	

Stabilized with silver

Methylisoamyl ketone

• Metilisoamilchetone • Méthylisoamylcétone • Metilo isoamilcetona • Methylisoamylketon

Synonym:
• *5-Methyl-2-hexanone*
• *Isobutylacetone*

$C_7H_{14}O$
Molecular Weight: 114,19
CAS: 110-12-3

Classification transport
ONU: 2302
Transport Hazard class: 3
Packing group III



Warning
H226-H332
P210-P241-P261-P280-P303+P361+P353-P304+P340

Methylisoamyl ketone > RE - Pure

RE

Appearance ..... Clear colourless liquid	Density d20/20 ..... 0.811 - 0.815	Non volatile residue ..... ≤ 50 mg/Kg	Assay (GC) ..... ≥ 98.5 %
Refractive index at 20°C ..... 1.404 - 1.408	Water content (K.F.) ..... ≤ 500 mg/Kg	Free acid (as CH3COOH) ..... ≤ 100 mg/Kg	

Code	Size	Packaging	Notes
P0900221	2.5 l	Glass bottle	

Methyl isobutyl ketone

• Metile isobutilchetone • Méthylisobutylcétone • Metilo isobutilcetona • Isobutylmethylketon

Synonym:
• *4-Methyl-2-pentanone*
• *MIBK*

$CH_3COCH_2CH(CH_3)_2$
Molecular Weight: 100,16
CAS: 108-10-1
EEC-N: 203-550-1

Classification transport
ONU: 1245
Transport Hazard class: 3
Packing group II



Danger
H225-H332-H319-H335-HEU066
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Methyl isobutyl ketone > RS - Anhydrous - For analysis

RS

Refractive index at 20°C ..... 1.394 - 1.398	Non volatile residue ..... ≤ 10 mg/Kg	Acetone ..... ≤ 0.1 %	Free acid (as CH3COOH) ..... ≤ 50 mg/Kg
Density d20/20 ..... 0.797 - 0.805	Colour ..... ≤ 10 Hazen	Mesityl and isomesityl oxide ..... ≤ 0.1 %	
Water content (K.F.) ..... ≤ 200 mg/Kg	Assay (GC) ..... ≥ 99.5 %	4-methyl-2-pentanol ..... ≤ 0.1 %	

Code	Size	Packaging	Notes
P0601016	1 l	Glass bottle	

Keep in a well-ventilated place

Methyl isobutyl ketone > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Clear liquid	Density at 20° C ..... 0.797 ÷ 0.805	Alcalinity (NH3) ..... ≤10 ppm	Zn ..... ≤0.1 ppm
Colour (APHA) ..... ≤10	Refractive index at 20°C. 1.3930 ÷ 1.3990	Cd ..... ≤0.05 ppm	Assay (GLC) ..... ≥99.5 %
Identification (I.R.) ..... Conform	Boiling point ..... 115.7 ÷ 116.7 °C	Cu ..... ≤0.1 ppm	Acetone ..... ≤ 0.1 %
Alcohol miscibility ..... Complete	Water (K.F.) ..... ≤0.05 %	Fe ..... ≤0.1 ppm	Mesityl and isomesityl oxide ..... ≤ 0.1 %
Benzene miscibility ..... Complete	Residue on evaporation ..... ≤10 ppm	Ni ..... ≤0.1 ppm	4-methyl-2-pentanol ..... ≤ 0.1 %
Diethyl ether miscib. ..... Complete	Acidity ..... ≤0.002 meq/g	Pb ..... ≤0.1 ppm	

Code	Size	Packaging	Notes
461945	1 l	Glass bottle	
461943	21 kg	Metal drum	

Keep in a well-ventilated place

DASIT GROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

Methyl isobutyl ketone > RE - Pure

RE

Description ..... Clear liquid	Density at 20°C ..... 0,797 ÷ 0,805	Assay (GC) ..... ≥ 99.5 %	Free acid (as CH3COOH) ..... ≤ 50 mg/Kg
Density d20/20 ..... 0.797 - 0.805	Non volatile residue ..... ≤ 30 mg/Kg	Acidity (Acetic ac.) ..... ≤ 50 ppm	Assay (GLC) ..... ≥ 99.5 %
Identification ..... Positive	Refractive index at 20°C 1.3930 ÷ 1.3990	Residue on evaporation ..... ≤ 30 ppm	Acetone ..... ≤ 0.1 %
Water content (K.F.) ..... ≤ 1000 mg/Kg	Boiling point ..... 115,7 ÷ 116,7 °C	Water (K.F.) ..... ≤ 1000 ppm	Colour ..... ≤ 15 APHA

Code	Size	Packaging	Notes
528980	5 l	Plastic tank	
528981	25 l	Metal drum	
P0600268	200 l	Metal drum	

Keep in a well-ventilated place

Methyl orange

• Arancio metile • Méthylorange • Naranja de metilo • Orangenschnaps

Synonym:
- *4-[4-(Dimethylamino)phenylazo]benzenesulfonic acid sodium salt*
- *Acid Orange 52*

$C_{14}H_{14}N_3NaO_3S$
Molecular Weight: 327,34
CAS: 547-58-0
EEC-N: 208-925-3

Classification transport
ONU: 3143
Transport Hazard class: 6.1
Packing group II

Danger
H301
P264-P270-P301+P310a-P330-P405-P501a

Methyl orange > RPE - For analysis - C.I. 13025 - ACS

RPE

Description ..... Yellow-orange powder | Identification ..... Positive | Colour change ..... rosso-giallo | pH range ..... 3.2 ÷ 4.4 pH

Code	Size	Packaging	Notes
423504	25 g	Glass bottle	
423503	50 g	Plastic bottle	
423505	250 g	Plastic bottle	
423501	500 g	Plastic bottle	
423502	25 kg	Drum	

Dye for microscopy (histology)

Methyl Orange solution 0.1%

• Arancio metile soluzione 0.1% • Méthylorange solution 0.1% • Naranja de metilo solución 0.1%
• Methylorangenlösung 0.1%

Synonym:
- *4-[4-(Dimethylamino)phenylazo]benzenesulfonic acid sodium salt*
- *Acid Orange 52*

$C_{14}H_{14}N_3NaO_3S$
Molecular Weight: 327,34
CAS: 547-58-0

Methyl Orange solution 0.1% > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611054802	100 ml	Plastic bottle	Solution in ethanol Ref Ph.Eur 1054802

Color change: pH 3.0 (red) to pH 4.4 (yellow)

Methyl Orange solution 0.1% > RPE - For analysis

RPE

Description ..... Orange clear liquid | Identification ..... Positive | Sensitivity(pH 3.1-4.4) ..... Conform | Colour change ..... red-yellow

Code	Size	Packaging	Notes
E423562	500 ml	Plastic bottle	

Acid-base indicator

Methyl orange mixed solution

• Indicatore misto di metile arancione • Indicateur mixte au méthylorange • Naranja de metilo solución mixta • Mischindikator mit Methylorange

Synonym:
- *4-[4-(Dimethylamino)phenylazo]benzenesulfonic acid sodium salt*
- *Acid Orange 52*

Methyl orange mixed solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611054801	100 ml	Plastic bottle	Ref Ph.Eur 1054801

Colour change: pH 3.0 (orange) to pH 4.4 (olive-green)

2-Methylpentane ▶ Isohexane

4-Methyl-2-pentanone ▶ Methyl isobutyl ketone

2-Methylphenol ▶ o-Cresol

Methyl phenyl ketone ▶ Acetophenone

2-Methyl-2-propanol ▶ tert-Butanol

N-Methyl-2-pyrrolidone

• N-Metile-2-pirrolidone • N-Méthyle-2-pyrrolidone • N-Metil-2-pirrolidona • N-Methyl-2-pyrrolidon

Synonym:
- *1-Methyl-2-pyrrolidone*
- *NMP*

$CH_2.(CH_2)_2.CON.CH_3$
CAS: 872-50-4
EEC-N: 212-828-1

Danger
H315-H319-H360D-H335-HA26
P261-P271-P280-P304+P340-P305+P351+P338-P403+P233

N-Methyl-2-pyrrolidone > RS - For Headspace chromatography

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 1.026 - 1.032
Refractive index at 20°C ..... 1.469 - 1.471
Water (K.F.) ..... ≤ 0.1 %

Assay (GLC) ..... ≥ 99.8 %
UV cut off ..... ≤ 269 nm
U.V. Transmittance
At 285 nm ..... ≥ 30 %
At 300 nm ..... ≥ 55 %

At 320 nm ..... ≥ 78 %
≥ 350 nm ..... ≥ 97 %
GC/HS
Residual solvent of class 1(acc. to ICH) . ≤ 1 µg/g

Residual solvent of class 2(acc. to ICH)≤ 10 µg/g
Residual solvent of class 3(acc. to ICH)≤ 50 µg/g

Code	Size	Packaging	Notes
462881	1 l	Glass bottle	

N-Methyl-2-pyrrolidone > RS - Anhydrous - For analysis

RS

Appearance ..... Clear liquid
Identification ..... Conform
Density d20/4 ..... 1.026 - 1.032

Refractive index at 20°C ..... 1.466 - 1.471
Colour ..... ≤ 20 Hazen
Water content (K.F.) ..... ≤ 200 mg/Kg

Assay (GC) ..... ≥ 99.8 %
Butyrolactone ..... ≤ 500 mg/Kg
Monomethylamine ..... ≤ 50 mg/Kg

Code	Size	Packaging	Notes
P0871010	200 ml	Bottle with septum	

N-Methyl-2-pyrrolidone > RS - For peptide synthesis

RS

Clear liquid appearance ..... Conform
Identification (IR) ..... Conform
Refractive index at 20°C ..... 1.466 - 1.471

Water content (K.F.) ..... ≤ 400 mg/Kg
Colour ..... ≤ 15 Hazen
Bromophenol blue test ..... Conform

Amines content ..... ≤ 5 mg/Kg
Assay (GC) ..... ≥ 99.7 %
Non volatile residue ..... ≤ 10 mg/Kg

Code	Size	Packaging	Notes
P0873516	1 l	Glass bottle	
P0873521	2.5 l	Glass bottle	
P0873541	10 l	Plastic tank	
P0873549	25 l	Plastic tank	
P0873566	200 l	Polythene-metal drum	

Product specifications are subject to changes. Please visit our website for updates.

N-Methyl-2-pyrrolidone > RPE - For analysis - ACS

RPE

Description ..... Clear liquid
Identification ..... Positive
Colour ..... ≤ 50 APHA
Density at 20° C ..... 1.026 ÷ 1.032
Refractive index at 20°C. 1.4670 ÷ 1.4710
Water (K.F.) ..... ≤0.05 %
Chloride ..... ≤ 1 ppm
Free amines (as CH3NH2) ..... ≤ 0.01 %
Butirrolattone ..... ≤ 500 ppm
Monometilammina ..... ≤ 50 ppm
Assay (GLC) ..... ≥99.8 %

Code	Size	Packaging	Notes
462872	1 l	Glass bottle	
462875	2.5 l	Glass bottle	
462874	23 kg	Drum	
462873	210 kg	Metal drum	

N-Methyl-2-pyrrolidone > RE - Pure

RE

Description ..... Clear liquid
Identification ..... Positive
Colour ..... ≤ 50 APHA
Density at 20°C ..... 1.026 ÷ 1.032
Refractive index at 20°C. 1.4670 ÷ 1.4710
Boiling point ..... 203.0 ÷ 205.0 °C
Water (K.F.) ..... ≤ 0.1 %
Assay (GLC) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
528341	1 l	Glass bottle	
528343	2.5 l	Glass bottle	
528340	5 l	Plastic tank	
528346	25 l	Metal drum	

Methyl red

• Rosso metile • Rouge méthyle • Rojo de metilo • Methylrot

Synonym:
• *2-(4-Dimethylaminophenylazo)benzoic acid*
• *4-Dimethylaminobenzene-2'-carboxylic acid*

$C_{15}H_{15}N_3O_2$
Molecular Weight: 269,31
CAS: 493-52-7
EEC-N: 207-776-1

Methyl red > RPE - For analysis - C.I. 13020

RPE

Description ..... Purple powder
Identification ..... Positive
Appear of water sol. ..... Conform
Colour change ..... red-yellow
pH range ..... 4.2 - 6.2

Code	Size	Packaging	Notes
476882	25 g	Glass bottle	
476883	50 g	Plastic bottle	
476881	250 g	Plastic bottle	

Methyl red solution water/ethanol 0.2%

• Rosso metile soluzione 0.2% in alcole etilico • Rouge méthyle solution 0.2% dans l'éthanol
• Rojo de metilo solución 0.2% en alcohol etilico • Methylrotlösung Wasser / Ethanol 0.2%

Synonym:
• *2-(4-Dimethylaminophenylazo)benzoic acid*
• *4-Dimethylaminobenzene-2'-carboxylic acid*

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Methyl red solution water/ethanol 0.2% > RPE - For analysis

RPE

Description ..... Purple liquid
Identification ..... Positive
pH sensitivity ..... 4.2 ÷ 6.2
Colour change ..... yellow red

Code	Size	Packaging	Notes
E476915	250 ml	Glass bottle	

Indicator series Clark indicator acid-base

Methyl red solution 0.1% in ethanol

• Rosso metile soluzione 0,1% in alcole etilico • Rouge méthyle solution 0.1% dans l'éthanol
• Rojo de metilo solución 0.1% en alcohol etilico • Methylrotlösung 0.1% in Ethanol

Synonym:
• *2-(4-Dimethylaminophenylazo)benzoic acid*
• *4-Dimethylaminobenzene-2'-carboxylic acid*

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II



Danger
H225-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

Methyl red solution 0.1% in ethanol > RPE - For analysis

RPE

Description ..... Purple liquid Identification ..... Positive pH range ..... 4.4 - 6.2

Code	Size	Packaging	Notes
E476921	250 ml	Glass bottle	

Indicator series Clark indicator acid-base

Methyl red solution

• Rosso metile soluzione • Rouge de méthyle solution • Rojo de metilo solución • Methylrotlösung

Synonym:
• *2-(4-Dimethylaminophenylazo)benzoic acid*
• *4-Dimethylaminobenzene-2'-carboxylic acid*

$C_{15}H_{15}N_3O_2$
Molecular Weight: 269,31
CAS: 493-52-7

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II



Danger
H225-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

Methyl red solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611055102	100 ml	Plastic bottle	Ref Ph.Eur 1055102

Colour change: pH 4.4 (red) to pH 6.0 (yellow)

Methyl red solution > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
617000111	100 ml	Plastic bottle	Methyl red solution TS

Methyl red mixed solution

• Indicatore misto di rosso di metile • Indicateur mixte au rouge de méthyle • Rojo de metilo solución mixta • Mischindikator mit Methylrot

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II



Danger
H225-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

Methyl red mixed solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611055101	100 ml	Plastic bottle	Ref Ph.Eur 1055101

Colour change: pH 5.2 (red-violet) to pH 5.6 (green)

Product specifications are subject to changes. Please visit our website for updates.

Methyl salicylate

• Metile salicilato • Méthyle salicylate • Metilo salicilato • Methylsalicylat

Synonym:
2-Hydroxybenzoic acid methyl ester

$HOC_6H_4COOCH_3$
Molecular Weight: 152,15
CAS: 119-36-8
EEC-N: 204-317-7

Warning
H302-H319
P264-P280i-P301+P312a-P305+P351+P338-P337+P313-P501a

Methyl salicylate > ERBApharm - According to pharmacopoeia: BP-DAB-FU-Ph.Eur.-Ph.Franc.

ERBApharm

Description ..... Yellowish liquid
Identification ..... Positive
Refractive index at 20°C ..... 1.535 ÷ 1.538
Assay (saponification) ..... 99.0 ÷ 101.0 %
Appearance of solution ..... Conform Ph.Eur.
Acidity ..... Conform Ph.Eur.
Density d20/20 ..... 1.182 - 1.188

Code	Size	Packaging	Notes
354152	1 l	Glass bottle	
354155	25 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Methyl sulfoxide ▶ Dimethylsulphoxide

Methyl tert-butyl ether ▶ tert-Butylmethylether

2-Methyltetrahydrofuran

• 2-Metiltetraidrofurano • 2-Méthyltétrahydrofurane • 2-Metiltetrahidrofurano • 2-Methyltetrahydrofuran

Synonym:
• *Tetrahydro-2-methylfuran*
• *2-MeTHF*

$C_5H_{10}O$
Molecular Weight: 86,14
CAS: 96-47-9
EEC-N: 202-507-4

Classification transport
ONU: 2536
Transport Hazard class: 3
Packing group II

Danger
H225-H302-H318-H335-HEU019
P210-P280-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

2-Methyltetrahydrofuran > RS - For HPLC - Isocratic Grade

RS

Clear,colourless liq.appearance ..... Conform
Density d20/4 ..... 0.853 - 0.857
Identification ..... Conform
Colour ..... ≤ 10 Apha
Refractive index at 20°C ..... 1.404 - 1.408
Water content (K.F.) ..... ≤ 200 mg/Kg
Non volatile residue ..... ≤ 5 mg/Kg
Peroxides (as H2O2) ..... ≤ 300 mg/Kg
Assay (GC) ..... ≥ 99.5 %
UV transmittance at 240 nm ..... ≥ 30 %
UV transmittance at 250 nm ..... ≥ 50 %
UV transmittance at 260 nm ..... ≥ 70 %
UV transmittance at 280 nm ..... ≥ 90 %
UV transmittance from 310 nm ..... ≥ 98 %

Code	Size	Packaging	Notes
412681	1 l	Glass bottle	
412682	2.5 l	Glass bottle	

2-Methyltetrahydrofuran > RE - Pure

RE

Appearance ..... Clear colourless liquid
Refractive index at 20°C ..... 1.404 - 1.408
Density d20/4 ..... 0.853 - 0.857
Water content (K.F.) ..... ≤ 300 mg/Kg
Assay (GC) ..... ≥ 99.9 %
Stabilizer (ionol) ..... 150 - 400 mg/Kg
Peroxides (as H2O2) ..... ≤ 100 mg/Kg

Code	Size	Packaging	Notes
P9960216	1 l	Glass bottle	
P9960221	2.5 l	Glass bottle	
P9960229	5 l	Plastic tank	
P9960248	25 l	Metal drum	
P9960268	200 l	Metal drum	

4-Methyltetrahydropyran

• 4-Metiltetraidropirano • 4-Méthyltétrahydropyrane • 4-Metiltetrahidropyrano • 4-Methyltetrahydropyran

$C_6H_{12}O$
Molecular Weight: 100,16
CAS: 4717-96-8
EEC-N: 225-207-5

Classification transport
ONU: 2924
Transport Hazard class: 3
Packing group II



Danger
H225-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

4-Methyltetrahydropyran > RE - Pure

RE

Appearance ..... Clear colourless liquid Colour ..... ≤ 10 Hazen Density d20/4 ..... 0.855 - 0.865 Assay (GC) ..... ≥ 99.0 %

Code	Size	Packaging	Notes
P9990218	500 ml	Glass bottle	
P9990216	1 l	Glass bottle	
P9990221	2.5 l	Glass bottle	

Methylthymol blue sodium salt

• Blu metiltimolo sale sodico • Bleu de méthylthymol sel de sodium • Azul de metiltimol sal sódica
• Methylthymolblau Natriumsalz

Synonym:
3,3'-Bis[n,n-di(carboxymethyl)aminomethyl thymolsulfonephthalein sodium salt

$C_{37}H_{40}O_{13}N_2Na_4S$
Molecular Weight: 844,76
CAS: 1945-77-3
EEC-N: 217-743-3



Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Methylthymol blue sodium salt > RPE - For analysis

RPE

Description ..... Greenish brown powder Identification ..... Positive

Code	Size	Packaging	Notes
429021	1 g	Glass bottle	
429022	25 g	Glass bottle	

Complexometric indicator

Methyl yellow

• Giallo metile • Jaune de méthyle • Amarillo de metilo • Methylgelb

Synonym:
• *4-(Dimethylamin)azobenzene*
• *N,N-Dimethyl-4-(phenylazo)aniline*

$(CH_3)_2NC_6H_4N{:}NC_6H_5$
Molecular Weight: 225,29
CAS: 60-11-7
EEC-N: 200-455-7

Classification transport
ONU: 3143
Transport Hazard class: 6.1
Packing group III



Danger
H301-H317
P261-P264-P280g-P301+P310a-P330-P362+P364

Methyl yellow > RPE - For analysis - C.I. 11020

RPE

Description ..... Yellow orange powder pH range ..... 2.9 - 4.0 Colour change ..... Red - yellow
Identification ..... Positive Loss on drying ..... ≤3 %

Code	Size	Packaging	Notes
444552	25 g	Glass bottle	

MIBK ▶ Methyl isobutyl ketone

Mineral oil ▶ Paraffin oil

Product specifications are subject to changes. Please visit our website for updates.

Mixture C.H.M.

• Miscela C.H.M. • Mélange CHM • Mezcla C.H.M. • Mischung C.H.M.

Classification transport
ONU: 1992
Transport Hazard class: 3
Packing group II



Danger
H225-H302-H315-H319-H351-H361d-H336-H372-H410-HEU301
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Mixture C.H.M. > RPE - For analysis

RPE

Description ..... Clear colourless liquid Colour ..... ≤ 10 APHA Density at 20°C ..... 1.070 ÷ 1.080 Assay (GLC) ..... Conform

Code	Size	Packaging	Notes
524411	2.5 l	Glass bottle	
524412	5 l	Plastic tank	

Composition: Chloroforme stab. Amylene: 49%(v/v) n-Heptane: 49%(v/v) Methanol: 2%(v/v)

Mix Diethyl ether/Ethanol 70/30 w/w

• Miscela alcole-etere 30/70 m/m • Mélange alcool/ether 30/70 m/m • Mezcla alcohol-éter 30/70 p/p • Mische Diethylether / Ethanol 70/30 m/m

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II



Danger
H225-H319-H336-HEU019
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Mix Diethyl ether/Ethanol 70/30 w/w > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Density at 20° C ..... 0.740 ÷ 0.750

Code	Size	Packaging	Notes
463251	1 l	Glass bottle	
463255	2.5 l	Glass bottle	

Contains Phenolphtalein

Mix Diethyl ether / Ethanol 95% 50/50(w/w) - 30 mg/l phenolphtaleine

• Etanolo 95%-Etere etilico 50/50(m/m) - 30 mg/l fenolftaleina • Ethanol 95%-Ether ethylique 50/50(m/m) - 30 mg/l phénolphtaléine
• Etanol 95% - Etil éter 50/50 (m/m) - 30 mg/l fenolftaleína • Ethanol 95% -Diethylether 50/50 m/m -30 mg/l Phenolphthalin

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group I



Danger
H224-H319-H336-HEU019
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Mix Diethyl ether / Ethanol 95% 50/50(w/w) - 30 mg/l phenolphtaleine > RS - For analysis

RS

Description ..... Clear colourless liquid Colour ..... < 10 APHA Density at 20°C ..... 0.757 ÷ 0.767

Code	Size	Packaging	Notes
529371	5 l	Aluminium can	

Mix Ethanol absolute/Diethyl ether 50/50 (w/w)

• Miscela Etanolo assoluto/Etere Etilico 50/50 (p/p) • Mélange Ethanol absolu/Ether Ethylique 50/50 (m/m) • Mezcla Etanol absoluto/Eter Etilico 50/50 (p/p)
• Mischung Ethanol absolut/Diethylether 50/50 (m/m)

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group I



Danger
H224-H319-H336-HEU019
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Mix Ethanol absolute/Diethyl ether 50/50 (w/w) > RE - Pure

RE

Description ..... Clear colourless liquid Colour ..... ≤ 10 APHA

Code	Size	Packaging	Notes
529311	5 l	Aluminium can	

Mix Ethanol absolute/Diethyl ether 50/50 (v/v)

• Miscela Etanolo assoluto/Etere Etilico 50/50 (v/v) • Mélange Ethanol absolu/Ether Ethylique 50/50 (v/v) • Mezcla Etanol absoluto/Eter Etilico 50/50 (v/v) • Mischung Ethanol absolut/Diethylether 50/50 (v/v)

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group -



Danger
H224-H319-H336-HEU019
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Mix Ethanol absolute/Diethyl ether 50/50 (v/v) > RE - Pure

RE

Description ..... Clear colourless liquid Colour ..... ≤ 10 APHA Density at 20°C ..... ~ 0.757

Code	Size	Packaging	Notes
529381	5 l	Aluminium can	

Mixture Ethanol 95° / Isopropanol

• Miscela Etanolo 95° / Isopropanolo • Mélange éthanol 95° / 2-propanol • Mezcla etanol 95°/ 2-propanol • Mischung Ethanol 95 ° / Isopropanol

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II



Danger
H225-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

Mixture Ethanol 95° / Isopropanol > RPE - For analysis

RPE

Description ..... Clear colourless liquid Density d20/4 ..... 0.78 - 0.80 Propan-2-ol content ..... ≤ 15 % (w/w)
Water (K.F.) ..... ≤ 5 % Ethanol content ..... ≥ 80 % (w/w) Assay (alcohometric) Total at 20°C ..... 94 - 96 % (V/V)

Code	Size	Packaging	Notes
414551	5 l	Plastic tank	In Vitro Diagnostic Medical Device

Mixture Ethanol 99° / Isopropanol

• Miscela Etanolo 99° / Isopropanolo • Mélange éthanol 99° / 2-propanol • Mezcla etanol 99°/2-propanol • Mischung Ethanol 99 ° / Isopropanol

Classification transport
ONU: 1987
Transport Hazard class: 3
Packing group II



Danger
H225-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

Mixture Ethanol 99° / Isopropanol > RPE - For analysis

RPE

Description ..... Clear colourless liquid Density d20/4 ..... 0.78 - 0.80 Propan-2-ol content ..... 5 - 15 % (w/w)
Water (K.F.) ..... ≤ 0.5 % Ethanol content ..... 85 - 95 % (w/w) Assay (alcohometric) Total at 20°C ..... ≥ 99 % (V/V)

Code	Size	Packaging	Notes
414511	5 l	Plastic bottle	In Vitro Diagnostic Medical Device

Mixture for bromine index determination

• Miscela per la determinazione dell'indice di bromo • Mélange pour indice de brome • Mezcla para la determinación de índice de bromo • Mischung für Bromindex

Classification transport
ONU: 2920
Transport Hazard class: 8
Packing group II



Danger
H226-H302-H314-H351-H370-H373
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Mixture for bromine index determination > RS - For analysis

RS

Refractive index at 20°C ..... 1.375 - 1.379 Density d20/4 ..... 1.060 - 1.068

Code	Size	Packaging	Notes
PS0573/21	2.5 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Mixture for checking solderings

• Miscela per il controllo delle saldature degli imballaggi • Mélange pour le contrôle des soudures des emballages
• Mezcla para el control de soldaduras de envases • Mischung zur Kontrolle von Lötungen

Classification transport
ONU: 1219
Transport Hazard class: 3
Packing group II



Danger
H225-H319-H336
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Mixture for checking solderings > RS - For agroalimentary analysis

RS

Description ..... Clear purple liquid Density at 20°C ..... 0.788 ÷ 0.798

Code	Size	Packaging	Notes
502671	5 l	Plastic tank	

Mixtures for residual solvents analysis

• Soluzioni standard per la ricerca di solventi residui • Mélanges pour la recherche des solvants résiduels • Mezclas para disolventes residuales analisis
• Gemische zum Nachweis von Restlösemitteln

Mixtures for residual solvents analysis > RS - For analysis according to Ph. Eur. Chap. 2.4.24

RS

Code	Size	Packaging	Notes
507688	1 ml	Glass ampoule	5 elements (Class 1): Benzene 4mg/ml; 1,2-Dichloroethane 10mg/ml; Tetrachloromethane 8mg/ml; 1,1-Dichloroethene 16mg/ml; 1,1,1-Tric hloroethane 20mg/ml; Matrix: Dimethylsulphoxide
507689	1 ml	Glass ampoule	14 elements (Class 2): Chlorobenzene 360µg/ml; Cyclohexane 3880µg/ml; cis-1,2-Dichloroethene 1870µg/ml; Dichloromethane 600µg/ml; Ethylbenzene 369µg/ml; n-Hexane 290µg/ml; Methylcyclohexane 1180µg/ml; n,n-Dimethylformamide 880µg/ml; Toluene 890µg/ml; 1,1,2-Trichloroethene 80µg/ml; m-Xylene 1302µg/ml; o-Xylene 195µg/ml; p-Xylene 304µg/ml; Tetrahydrofuran 720µg/ml; Matrix: Dimethylsulfoxide
507690	1 ml	Glass ampoule	11 elements (Class 2): Acetonitrile 410mg/l; Chloroform 60mg/l; 1,2-Dimethoxyethane 100mg/l; n,n-Dimethylacetamide 1090mg/l; Dioxan 380mg/l; 2-Hexanone 50mg/l; Methanol 3000mg/l; Nitromethane 50mg/l; Pyridine 200mg/l; 1,2,3,4-Tetrahydronaphthalene 100mg/l; Isop ropylbenzene (Cumene) 70mg/l; Matrix: Dimethylsulphoxide/ Water
507691	1 ml	Glass ampoule	6 elements (Class 2): Ethyleneglycol-monoethyl ether (2-Ethoxyethanol) 160mg/l; Ethyleneglycol 620mg/l; Formamide 220mg/l; Ethylene glycol-monomethyl ether (2-Methoxyethanol) 50mg/l; 1-Methyl-2-pyrrolidon 4840mg/l; Sulfolan 160mg/l; Matrix: Water

Mixtures for residual solvents analysis > RS - For analysis according to USP – Residual solvents

RS

Code	Size	Packaging	Notes
507692	1 ml	Glass ampoule	5 elements (Class 1): Benzene 10mg/ml; Tetrachloromethane (Carbon tetrachloride) 20mg/ml; 1,2-Dichloroethane 25mg/ml; 1,1-Dichloroethene 40mg/ml; 1,1,1-Trichloroethane 50mg/ml; Matrix: Dimethylsulphoxide
507693	1 ml	Glass ampoule	16 elements (Class 2): Acetonitrile 2.05mg/ml; Chlorobenzene 1.8mg/ml; Cumene 0.34mg/ml; Cyclohexane 19.4mg/ml; cis-1,2-Dichloroethene 4.7mg/ml; trans-1,2-Dichloroethene 4.7mg/ml; 1,4-Dioxan 1.9mg/ml; Ethylbenzene 1.84mg/ml; Methanol 15mg/ml; Methylcyclohexane 5.9mg/ml; Dichloromethane 3mg/ml; Tetrahydrofuran 3.6mg/ml; Toluene 4.45mg/ml; m-Xylene 6.51mg/ml; o-Xylene 0.98mg/ml; p-Xylene 1.52mg/ml; Matrix: Dimethylsulphoxide
507694	1 ml	Glass ampoule	8 elements (Class 2): Chloroform 60µg/ml; 1,2-Dimethoxyethane 100µg/ml; n-Hexane 290µg/ml; 2-Hexanone 50µg/ml; Nitromethane 50µg/ml; Pyridine 200µg/ml; 1,2,3,4-Tetrahydronaphthalene (Tetralin) 100µg/ml; Trichloroethene 80µg/ml; Matrix: Dimethylsulphoxide

Molecular sieves 3 A

• Setacci molecolari 3 A • Tamis moléculaire 3 A • Tamices moleculares 3 A • Molekularsieb 3 A

HEU210

Molecular sieves 3 A > RS - Pellets 1/16'' — RS

Description ..... Beige granules Identification ..... Positive Apparent density ..... 700 ÷ 800 g/l

Code	Size	Packaging	Notes
477731	250 g	Glass bottle	

Molecular sieves 3 A > RS - Pellets 1/8'' — RS

Description ..... Sferett e 1/8" Identification ..... Positive

Code	Size	Packaging	Notes
477721	250 g	Glass bottle	

Molecular sieves 3 A > RS - Pure — RS

Code	Size	Packaging	Notes
P1810017	1 kg	Plastic bottle	

Pellets between 1/8" and 1/16"

Molecular sieves 4 A

• Setacci molecolari 4 A • Tamis moléculaire 4 A • Tamices moleculares 4 A • Molekularsieb 4A

HEU210

Molecular sieves 4 A > RS - Pellets — RS

Code	Size	Packaging	Notes
P1820017	1 kg	Plastic bottle	
P1820027	5 kg	Plastic bucket	
P1820047	25 kg	Plastic bucket	

Granulometry: 1.6-2.5 mm

Molybdenum standard solution

• Molibdeno standard soluzione • Molybdène solution standard • Molibdeno, solución patrón • Molybdän-Standardlösung

Molybdenum standard solution > RS - Standard solution for ICP-MS — RS

Code	Size	Packaging	Notes
505722	100 ml	Plastic bottle	conc. 10 ppm Matrix: Nitric acid and hydrofluoric acid
505725	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid and hydrofluoric acid
505723	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Molybdenum standard solution > RS - Standard solution for ICP — RS

Code	Size	Packaging	Notes
503731	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Ammonium hydroxide
503733	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Ammonium hydroxide
503735	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Ammonium hydroxide
503737	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Ammonium hydroxide

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Product specifications are subject to changes. Please visit our website for updates.

Molybdenum standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
507747	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrofluoric acid and nitric acid
507494	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid
E497565	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Ammonium hydroxide
E497561	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Ammonium hydroxide

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Molybdenum standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
463431		Plastic ampoule	conc. 1.000 ppm Matrix: Ammonium hydroxide - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Molybdenum (VI) oxide

• Anidride molibdica • Anhydride molybdique • Anhídrido molibdénico • Molybdän (VI) oxid

Synonym:
Molybdenum trioxide

MoO_3
Molecular Weight: 143,94
CAS: 1313-27-5
EEC-N: 215-204-7

Classification transport
ONU: 3288
Transport Hazard class: 6.1
Packing group III



Warning
H319-H351-H335
P261-P271-P280-P304+P340-P305+P351+P338-P403+P233

Molybdenum (VI) oxide > RPE - For analysis

RPE

Description ..... White,slightly yellow, green, or gray solid
Identification ..... Positive
Ammonium ..... ≤ 100 ppm
Chloride ..... ≤ 50 ppm
Phosphate ..... ≤ 5 ppm
Nitrate ..... ≤ 100 ppm
Sulphate ..... ≤ 200 ppm
Fe ..... ≤ 50 ppm
Assay (oxidimetric) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
422004	100 g	Glass bottle	
422005	250 g	Glass bottle	

Molybdophosphoric acid ▶ Phosphomolybdic acid

Molybdovanadic reagent

• Reattivo molibdovanadico • Réactif molybdovanadique • Reactivo molibdenovanádico • Molybdovanadisches Reagens

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group III



Danger
H302-H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Molybdovanadic reagent > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611056700	100 ml	Plastic bottle	Ref Ph.Eur 1056700

Mordant Black 11 ▶ Eriochrome black T

Mordant Red 11 ▶ Alizarin

Morpholine

• Morfolina • Morpholine • Morfolina • Morpholin

Synonym:
Tetrahydro-1,4-oxazine

$NH(CH_2)_2OCH_2CH_2$
Molecular Weight: 87,12
CAS: 110-91-8
EEC-N: 203-815-1

Classification transport
ONU: 2054
Transport Hazard class: 8
Packing group I



Danger
H226-H302-H312-H332-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P362+P364

Morpholine > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP — RPE

Description ..... Clear liquid
Colour (APHA) ..... ≤ 10
Identification ..... Positive
Assay (acidimetric) ..... ≥ 99.0 %
Boiling point ..... 126.0 ÷ 130.0 °C
Density at 20°C ..... ~ 1.01

Code	Size	Packaging	Notes
463453	1 l	Glass bottle	

MTBE ▶ tert-Butylmethylether

Mucicarmine hydroalcoholic solution

• Mucicarminio soluzione idroalcolica • Mucicarmin en solution hydroalcoolique • Mucicarmin solución hidroalcohólica
• Mucicarmin in hydroalkoholischer Lösung

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II



Danger
H225-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

Mucicarmine hydroalcoholic solution > RS - For microscopy — RS

Description ..... Red clear liquid
Identification ..... Positive

Code	Size	Packaging	Notes
463531	100 ml	Glass bottle	

Ethanol-water mixture (50:50)

Multianions standard for ion chromatography

• Standard multianione per cromatografia ionica • Etalon multiéléments pour chromatographie ionique • Patrón multielementos para cromatografía iónica
• Multianionen-Standard für die Ionenchromatographie

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group II

Multianions standard for ion chromatography > RS - Standard solution according to EPA method — RS

Code	Size	Packaging	Notes
504526	100 ml	Plastic bottle	7 elements: Br-, Cl-, NO3-, NO2-, PO43-, SO42-, F- 1g/l each - Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Multianions standard for ion chromatography > RS - Quality control Standard — RS

Code	Size	Packaging	Notes
504527	100 ml	Plastic bottle	7 elements: Br- 100ppm, SO42- 150ppm, PO43- 50ppm, Cl- 30ppm, NO2- 30ppm, NO3- 20ppm, F- 20 ppm - Matrix: Water
504677	500 ml	Plastic bottle	7 elements: F- 20mg/l; Cl- 100mg/l; NO2- 100mg/l; Br- 100mg/l; NO3- 100mg/l; PO43- 200mg/l; SO42- 100mg/l; Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

DASIT GROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

Multielement standard for ICP

• Standard multielemento per ICP • Etalon multiéléments pour ICP • Patrón multielementos para ICP • Multielement-Standard für ICP

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Multielement standard for ICP > RS - For analysis according to USP – WK Dietary supplement

RS

Code	Size	Packaging	Notes
506120	100 ml	Plastic bottle	4 elements: Cd 5mg/l; Pb 10mg/l; As 15mg/l; Hg 15mg/l. Matrix: 7% HNO3

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Multielement standard for ICP > RS - For analysis according to USP – Oral elemental impurities

RS

Code	Size	Packaging	Notes
506110	100 ml	Plastic bottle	8 elements: Cd 25mg/l; Pb 5mg/l; As 1.5mg/l; Hg 15mg/l; Mo 100mg/l; Ni 500mg/l; V 100mg/l; Cu 1000mg/l. Matrix: 7% HNO3
506150	100 ml	Plastic bottle	Precious metals - 6 elements: Ir 100 mg/l, Pt 100 mg/l; Os 100 mg/l; Rh 100 mg/l; Pd 100 mg/l; Ru 100 mg/l. Matrix: 15% HCl

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Multielement standard for ICP > RS - For analysis according to USP – Parenteral elemental impurities

RS

Code	Size	Packaging	Notes
506130	100 ml	Plastic bottle	Precious metals - 6 elements: Ir 10 mg/l, Pt 10 mg/l; Os 10 mg/l; Rh 10 mg/l; Pd 10 mg/l; Ru 10 mg/l. Matrix: 15% HCl
506140	100 ml	Plastic bottle	8 elements: Cd 2.5mg/l; Pb 5mg/l; As 1.5mg/l; Hg 1.5mg/l; Mo 10mg/l; Ni 50mg/l; V 10mg/l; Cu 100mg/l. Matrix: 7% HNO3

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Multielement standard for ICP > RS - Optimisation solution

RS

Code	Size	Packaging	Notes
504396	500 ml	Plastic bottle	13 elements: Al, Mg, Cr, Mn, Cu, Rh, In, Cd, Ce, Pb, Th, B, Ba 0,01mg/ml each - Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Multielement standard for ICP > RS - Quality control Standard

RS

Code	Size	Packaging	Notes
504350	100 ml	Plastic bottle	22 elements: As, Be, Bi, Ca, Cd, Co, Cr, Cu, Fe, Li, Mg, Mn, Mo, Ni, Pb, Sb, Se, Sr, Ti, Tl, V, Zn 1ppm each - Matrix: Nitric acid
504354	100 ml	Plastic bottle	33 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Cd, Co, Cs, Cr, Cu, Fe, In, K, Li, Mg, Mn, Mo, Na, Nb, Ni, Pb, Rb, Sb, Se, Sr, Ti, Tl, U, V, Zn 1ppm each - Matrix: Nitric acid
504356	100 ml	Plastic bottle	40 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Cd, Co, Cs, Cr, Cu, Fe, Ga, Ge, In, K, Li, Mg, Mn, Mo, Na, Ni, Nb, Pb, Pd, Rb, Sb, Se, Sn, Sr, Ta, Ti, Tl, U, V, W, Zn, Zr 1ppm each - Matrix: Nitric acid
504351	500 ml	Plastic bottle	22 elements: As, Be, Bi, Ca, Cd, Co, Cr, Cu, Fe, Li, Mg, Mn, Mo, Ni, Pb, Sb, Se, Sr, Ti, Tl, V, Zn 1ppm each - Matrix: Nitric acid
504353	500 ml	Plastic bottle	28 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Ce, Co, Cr, Cu, Fe, K, Li, Mg, Mn, Mo, Na, Ni, Pb, Sb, Se, Sr, Ti, Tl, V, Zn 1ppm each - Matrix: Nitric acid
504355	500 ml	Plastic bottle	33 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Cd, Co, Cs, Cr, Cu, Fe, In, K, Li, Mg, Mn, Mo, Na, Nb, Ni, Pb, Rb, Sb, Se, Sr, Ti, Tl, U, V, Zn 1ppm each - Matrix: Nitric acid
504357	500 ml	Plastic bottle	40 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Cd, Co, Cs, Cr, Cu, Fe, Ga, Ge, In, K, Li, Mg, Mn, Mo, Na, Ni, Nb, Pb, Pd, Rb, Sb, Se, Sn, Sr, Ta, Ti, Tl, U, V, W, Zn, Zr 1ppm each - Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Multielement standard for ICP > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
504301	100 ml	Plastic bottle	6 elements: Au, Ir, Pb, Pt, Rh, Ru 100ppm each - Matrix: Nitric acid
504303	100 ml	Plastic bottle	16 elements: Al, As, Ba, Be, Bi, Ca, Cs, Ga, In, K, Li, Mg, Na, Rb, Se, Sr 100ppm each - Matrix: Nitric acid
504305	100 ml	Plastic bottle	13 elements: Ag, Cd, Co, Cr, Cu, Fe, Hg, Mn, Ni, Pb, Tl, V, Zn 100ppm each - Matrix: Nitric acid
504480	100 ml	Plastic bottle	24 components: Ag 1mg/kg; Sb 1mg/kg; As 1mg/kg; Mn 1mg/kg; Cd 1mg/kg; Cr 1mg/kg; Ti 1mg/kg; Pb 1mg/kg; Co 1mg/kg; Ni 1mg/kg; Se 1mg/kg; V 1mg/kg; Mo 1mg/kg; Sn 1mg/kg; Ba 1mg/kg; Be 1mg/kg; Li 1mg/kg; Tl 1mg/kg; Bi 1mg/kg; Al 10mg/kg; Cu 10mg/kg; Fe 10mg/kg; B 10mg/kg; Zn 10mg/kg - Matrix: nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Product specifications are subject to changes. Please visit our website for updates.

Multielement standard for ICP > RS - Calibrating solution

RS

Code	Size	Packaging	Notes
504306	100 ml	Plastic bottle	22 elements: As, Be, Bi, Ca, Cd, Co, Cr, Cu, Fe, Li, Mg, Mn, Mo, Ni, Pb, Sb, Se, Sr, Ti, Tl, V, Zn 100ppm each - Matrix: Nitric acid
504308	100 ml	Plastic bottle	28 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Cd, Co, Cr, Cu, Fe, K, Li, Mg, Mn, Mo, Na, Ni, Pb, Sb, Se, Sr, Ti, Tl, V, Zn 100ppm each - Matrix: Nitric acid
504310	100 ml	Plastic bottle	33 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Cd, Co, Cs, Cr, Cu, Fe, In, K, Li, Mg, Mn, Mo, Na, Ni, Nb, Pb, Rb, Sb, Se, Sr, Ti, Tl, U, V, Zn 100ppm - Matrix: Nitric acid
504312	100 ml	Plastic bottle	9 elements: Au, Ir, Os, Pb, Pt, Rh, Ru, Sn, Te 100ppm each - Matrix: Hydrochloric acid
504307	500 ml	Plastic bottle	22 elements: As, Be, Bi, Ca, Cd, Co, Cr, Cu, Fe, Li, Mg, Mn, Mo, Ni, Pb, Sb, Se, Sr, Ti, Tl, V, Zn 100ppm each - Matrix: Nitric acid
504309	500 ml	Plastic bottle	28 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Cd, Co, Cr, Cu, Fe, K, Li, Mg, Mn, Mo, Na, Ni, Pb, Sb, Se, Sr, Ti, Tl, V, Zn 100ppm each - Matrix: Nitric acid
504311	500 ml	Plastic bottle	33 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Cd, Co, Cs, Cr, Cu, Fe, In, K, Li, Mg, Mn, Mo, Na, Ni, Nb, Pb, Rb, Sb, Se, Sr, Ti, Tl, U, V, Zn 100ppm each - Matrix: Nitric acid
504313	500 ml	Plastic bottle	9 elements: Au, Ir, Os, Pb, Pt, Rh, Ru, Sn, Te 100ppm each - Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Multielement standard for ICP and ICP-MS

• Standard multielemento per ICP e ICP-MS • Etalon multiéléments pour ICP et ICP-MS • Patrón multielementos para ICP e ICP-MS
• Multielement-Standard für ICP und ICP-MS

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group II

Multielement standard for ICP and ICP-MS > RS - Solution de Tuning pour ICP-MS

RS

Code	Size	Packaging	Notes
504392	100 ml	Plastic bottle	9 elements: Be, Mg, Co, In, Rh, Ce, Ba, Pb, U 10ppm each - Matrix: Nitric acid
504393	100 ml	Plastic bottle	13 elements: Ba, Be, Bi, Ce, Cu, Ho, In, Li, Mg, Pb, Tl, U, Y 10ppm each - Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Multielement standard for ICP and ICP-MS > RS - Quality control standard solution

RS

Code	Size	Packaging	Notes
504352	100 ml	Plastic bottle	28 elements: Al, Ag, As, B, Ba, Be, Bi, Ca, Ce, Co, Cr, Cu, Fe, K, Li, Mg, Mn, Mo, Na, Ni, Pb, Sb, Se, Sr, Ti, Tl, V, Zn 1ppm each - Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Murexide

• Muresside • Murexide • Murexida • Murexid

Synonym:
5,5'-Nitrilodibarbituric acid monoammonium salt

$C_8H_8N_6O_6$
Molecular Weight: 284,19
CAS: 3051-09-0
EEC-N: 221-266-6

Murexide > RPE - For analysis - C.I. 56085

RPE

Description ..... Red violet powder Identification ..... Positive Loss on drying ..... ≤10 % Residue on ignition ..... ≤0.5 %

Code	Size	Packaging	Notes
463608	5 g	Glass bottle	
463609	25 g	Glass bottle	

Complexometric indicator

Product specifications are subject to changes. Please visit our website for updates.

Naphthalene
• Naftalene • Naphtalène • Naftaleno • Naphthalin

$C_{10}H_8$
Molecular Weight: 128,17
CAS: 91-20-3
EEC-N: 202-049-5

Classification transport
ONU: 1334
Transport Hazard class: 4.1
Packing group III



Warning
H302-H351-H410
P264-P280-P301+P312a-P330-P308+P313-P501a

Naphthalene > RPE - For analysis — RPE

Description ..... White flakes Identification ..... Positive Melting point..... 79 ÷ 82 ° C Assay (GLC) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
463654	100 g	Plastic bottle	
463655	250 g	Plastic bottle	
463651	1 kg	Plastic bottle	

1-Naphthol
• 1-Naftolo • 1-Naphtol • 1-Naftol • 1-Naphthol

Synonym:
1-Hydroxynaphthalene

$C_{10}H_7OH$
Molecular Weight: 144,17
CAS: 90-15-3
EEC-N: 201-969-4



Danger
H302-H312-H315-H318-H335
P304+P340-P310a-P305+P351+P338-P330-
P362+P364-P403+P233

1-Naphthol > RPE - For analysis — RPE

Description ..... Grey-brown flakes Identification ..... Positive Melting point..... 94 ÷ 98 ° C Assay (GLC) ..... ≥ 97.5 %

Code	Size	Packaging	Notes
463935	250 g	Plastic bottle	

1-Naphthol > RE - Pure — RE

Description ..... White pinkish crystals Melting point..... 94 ÷ 98 ° C Assay (GLC) ..... ≥98 %
Identification ..... Positive á-Naphthole..... ≤1 %

Code	Size	Packaging	Notes
354751	250 g	Plastic bottle	

2-Naphthol
• 2-Naftolo • 2-Naphtol • 2-Naftol • 2-Naphthol

Synonym:
2-Hydroxynaphthalene

$C_{10}H_7OH$
Molecular Weight: 144,17
CAS: 135-19-3
EEC-N: 205-182-7

Classification transport
ONU: 3077
Transport Hazard class: 9
Packing group III



Warning
H302-H332-H400
P261-P264-P271-P301+P312a-P304+P340-P501a

2-Naphthol > RPE - For analysis - C.I. 37500 — RPE

Description ..... Flakes pinky Identification ..... Positive Melting point..... ≥ 120 °C Assay (GLC) ..... ≥ 99.0 %

Code	Size	Packaging	Notes
463984	100 g	Plastic bottle	
463986	500 g	Plastic bottle	

alpha-Naphtholbenzein

• alfa-Naftolbenzeina • alpha-Naphtolbenzéine • alpha-Naftolbenceína • alpha-Naphtolbenzol

Synonym:
4,4'-(alpha-Hydroxybenzylidene)di-1-naphtol

$C_{27}H_{18}O_2$
Molecular Weight: 374,44
CAS: 145-50-6
EEC-N: 205-656-3

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

alpha-Naphtholbenzein > RPE - For analysis

RPE

Description ..... Red brown powder Identification ..... Positive Suitability for anhydrous titration ..... Conform

Code	Size	Packaging	Notes
463891	5 g	Glass bottle	

Acid-base indicator

alpha-Naphtholbenzein solution 0.2% in acetic acid

• alfa-Naftolbenzeina soluzione 0.2% in acido acetico
• alpha-Naphtolbenzéine solution 0.2% dans acide acétique
• alpha-Naftolbenceína solución 0.2% en acido acético • alpha-Naphtolbenzollösung 0.2% in Essigsäure

Synonym:
4,4'-(alpha-Hydroxybenzylidene)di-1-naphtol

$C_{27}H_{18}O_2$
Molecular Weight: 374,44
CAS: 145-50-6

Classification transport
ONU: 2789
Transport Hazard class: 8
Packing group II

Danger
H226-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

alpha-Naphtholbenzein solution 0.2% in acetic acid > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611057601	100 ml	Glass bottle	Ref Ph.Eur 1057601

Naphthol yellow S

• Giallo naftolo S • Jaune naphtol S • Amarillo naftol S • Gelbes Naphthol S

Synonym:
Acid yellow 1|Flavinic acid sodium salt|2,4-Dinitro-1-naphtol-7-sulfonic acid sodium salt

$C_{10}H_4N_2Na_2O_8S$
Molecular Weight: 358,19
CAS: 846-70-8
EEC-N: 212-690-2

Classification transport
ONU: 3143
Transport Hazard class: 6.1
Packing group III

Warning
H302-H312-H332
P261-P264-P271-P280h-P301+P312a-P304+P340

Naphthol yellow S > RS - For microscopy - C.I. 10316

RS

Description ..... Pv xx.na giallo arancio Identification ..... Positive

Code	Size	Packaging	Notes
453562	25 g	Glass bottle	

Dye for histology

N-(1-Naphtyl)ethylenediamine dihydrochloride

• N-(1-Naftil)etilendiammina dicloridrato • N-(1-Naphtyl)éthylènediamine dichlorhydraté
• N-(1-Naftil)etilendiamina diclorhidrato • N-(1-Naphthyl) ethylendiamin-dihydrochlorid

Synonym:
2-(1-Naphthylamino)ethylamine dihydrochloride

$C_{10}H_7NHCH_2CH_2NH_2 \cdot 2HCl$
Molecular Weight: 259,18
CAS: 1465-25-4
EEC-N: 215-981-2

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

N-(1-Naphtyl)ethylenediamine dihydrochloride > RPE - For analysis

RPE

Description White or beige powder, crystal or chunks Identification ..... Positive Water ..... ≤ 5 % Titolo (TLC) ..... ≥ 98 %

Code	Size	Packaging	Notes
463831	10 g	Glass bottle	

Search for sulfonamides in the blood and spectrophotometric determination of nitrites and nitrates for FIA middle

Product specifications are subject to changes. Please visit our website for updates.

NDF Plus solution

• Soluzione NDF PLUS • Solution NDF PLUS • Solución NDF PLUS • Lösung NDF PLUS

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

NDF Plus solution > RPE - For agroalimentary analysis

RPE

Appearance ..... Clear liquid

Code	Size	Packaging	Notes
526941	25 l	Plastic tank	

Composition: Sodium laurylsulfate: 30 g; EDTA: 18.61 g; Sodium borate: 6.81 g; Hydrogenophosphate disodique: 4.6 g; Glycol monoethylether: 10 ml; Water: QSP 1 L according to NF V18-122

NDF Solution

• Soluzione NDF • Solution NDF • Solución NDF • Lösung NDF

Danger
H319-H360FD-HA26
P264-P280-P305+P351+P338-P308+P313-P337+P313-P501a

NDF Solution > RPE - For agroalimentary analysis

RPE

Description ..... Clear liquid Density at 20°C ..... 1.010 ÷ 1.025 pH at 20°C ..... 6.10 ÷ 7.10

Code	Size	Packaging	Notes
526920	2.5 l	Glass bottle	
526921	25 l	Plastic tank	

Composition: Sodium laurylsulfate: 30 g; EDTA: 18.61 g; Sodium borate: 6.81 g; Hydrogenophosphate disodique: 4.6 g; Glycol monoethylether: 10 ml; Water: QSP 1l according to NF V18-122

Neocuproine hydrochloride

• Neocuproina cloridrato • Néocuproïne chlorhydraté • Neocuproina clorhidrato • Neocupronhydrochlorid

Synonym:
2,9-Dimethyl-1,10-phenanthroline hydrochloride

$C_{14}H_{12}N_2.HCl.H_2O$
Molecular Weight: 244,72 (an,)
CAS: 7296-20-0
EEC-N: 230-732-8

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Neocuproine hydrochloride > RPE - For analysis

RPE

Description ..... Yellowish powder Identification ..... Positive Copper sensitivity ..... ≥0.3 µg/ml Assay (non-aqueos medium) ..... ≥99 % s s

Code	Size	Packaging	Notes
444731	1 g	Glass bottle	

Neodymium standard solution

• Neodimio standard soluzione • Néodyme solution standard • Neodimio, solución patrón • Neodym-Standardlösung

Classification transport
ONU: 3267
Transport Hazard class: 8
Packing group III

Neodymium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505742	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505745	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Neodymium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503761	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503763	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503765	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503767	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Neodymium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507748	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507510	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Nessler's reagent single solution

• Nessler reattivo soluzione unica • Réactif de Nessler solution • Nessler reactivo solución única • Nesslers Reagenzlösung

Classification transport
ONU: 2922
Transport Hazard class: 8
Packing group II



Danger
H301-H314-H341-H373-H412
P280-P301+P310a-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338

Nessler's reagent single solution > RPE - For analysis

RPE

Description ..... Yellow clear liquid Identification ..... Positive Sensibilità all'azoto ..... Conform

Code	Size	Packaging	Notes
464231	500 ml	Plastic bottle	
464232	1 l	Plastic bottle	

For the determination of ammonia and ammonium salt

Nessler's reagent solution A

• Nessler reattivo soluzione A • Réactif de Nessler solution A • Nessler reactivo solución A • Nessler Reagenzlösung A

Classification transport
ONU: 3287
Transport Hazard class: 6.1
Packing group II



Danger
H301-H311-H373-H411
P260-P264-P280h-P301+P310a-P330-P361+P364

Nessler's reagent solution A > RPE - For analysis

RPE

Description ..... Yellow clear liquid Identification ..... Positive Density at 20° C ..... ≥1.2

Code	Size	Packaging	Notes
464422	500 ml	Plastic bottle	

For the determination of nitrogen

Product specifications are subject to changes. Please visit our website for updates.

Nessler's reagent solution B

• Nessler reattivo soluzione B • Réactif de Nessler solution B • Nessler reactivo solución B • Natriumhydroxid

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Nessler's reagent solution B > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Assay .....20.50 ÷ 21.50 % (NaOH)

Code	Size	Packaging	Notes
464432	500 ml	Plastic bottle	

For the determination of nitrogen

Neutral red

• Rosso neutro • Rouge neutre • Rojo neutro • Neutralrot

Synonym:
3-Amino-7-dimethylamino-2-methylphenazine hydrochloride

$C_{15}H_{17}ClN_4$
Molecular Weight: 288,78
CAS: 553-24-2
EEC-N: 209-035-8

Classification transport
ONU: 3143
Transport Hazard class: 6.1
Packing group III

Warning
H302
P264-P270-P301+P312a-P330-P501a

Neutral red > RPE - For analysis - C.I. 50040

RPE

Description .....Green-brown powder Loss on drying .....≤ 5 % pH range ..... 6.8 ÷ 8.0
Identification ..... Positive Colour change..... red - yellow

Code	Size	Packaging	Notes
476951	10 g	Glass bottle	

Dye for microscopy (histology, hematology).

Nickel, powder

• Nichel, polvere • Nickel, poudre • Níquel, polvo • Nickel, Pulver

Ni
Molecular Weight: 58,71
CAS: 7440-02-0
EEC-N: 231-111-4

Danger
H317-H351-H372
P260-P264-P280-P308+P313-P362+P364-P333+P313

Nickel, powder > RPE - For analysis

RPE

Description .....Polvere metallica Identification ..... Positive Assay .....≥ 98 %

Code	Size	Packaging	Notes
464384	100 g	Glass bottle	

Nickel standard solution

• Nichel standard soluzione • Nickel solution standard • Niquel, solución patrón • Nickel-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Warning
H290-HEU208
P234-P390-P406

Nickel standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615002001	100 ml	Plastic bottle	A 0.1 ppm solution: to dilute according to Ref Ph.Eur 5002001
615002002	100 ml	Plastic bottle	A 0.2 ppm solution: to dilute according to Ref Ph.Eur 5002002
615002009	100 ml	Plastic bottle	A 10 ppm solution: to dilute according to Ref Ph.Eur 5002000

Nickel standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505752	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505755	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505753	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Nickel standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503771	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503773	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503775	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503777	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Nickel standard solution > RS - Standard solution for AAS

RS

Description ..... Green clear liquid Identification ..... Positive Titration factor ..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
507749	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507487	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
E497575	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
E497571	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Nickel standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
464271		Plastic ampoule	conc. 1.000 ppm Matrix: Nitric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Nickel standard solution > RS - Quality control standard solution for AAS (graphite furnace)

RS

Code	Size	Packaging	Notes
504363	50 ml	Plastic bottle	conc. 10 +/- 1 µg/L Matrix: 2% Nitric acid

Nickel (II) acetate tetrahydrate

• Nichel acetato oso tetraidrato • Nickel (II) acétate tétrahydraté • Níquel (II) acetato
• Nickel (II) acetattetrahydrat

Synonym:
Acetic acid nickel(II) salt

$Ni(CH_3COO)_2.4H_2O$
Molecular Weight: 248,86
CAS: 6018-89-9

Warning
H302-H332-H317-H351
P261-P271-P280-P304+P340-P308+P313-P330

Nickel (II) acetate tetrahydrate > RPE - For analysis

RPE

Description ..... Green powder
Identification ..... Positive
Co ..... ≤ 0.2 %
Cr ..... ≤ 5 ppm
Cu ..... ≤ 10 ppm
Fe ..... ≤ 20 ppm
Mn ..... ≤ 20 ppm
Pb ..... ≤ 5 ppm
Zn ..... ≤ 10 ppm
Assay ..... 23 ÷ 24 % (Ni)

Code	Size	Packaging	Notes
464474	100 g	Plastic bottle	
464476	500 g	Plastic bottle	
464477	1 kg	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Nickel (II) ammonium sulfate hexahydrate

• Nichel ammonio solfato oso • Nickel (II) ammonium sulfate hexahydraté • Níquel (II) amonio sulfato hexahidrato • Nickel (II) ammoniumsulfathexahydrat

$Ni(NH_4)_2(SO_4)_2.6H_2O$
Molecular Weight: 395
CAS: 7785-20-8

Danger
H302-H332-H334-H317-H341-H350i-H360D-H372-H410-HA26
P271-P280-P284-P304+P340-P308+P313-P342+P311a

Nickel (II) ammonium sulfate hexahydrate > RPE - For analysis

RPE

Description .. Green-blue crystalline powder
Identification ..... Positive
pH sol. 5% at 25° C ..... 4.3 ÷ 4.7
Chloride ..... ≤ 10 ppm
Water-insoluble matter ..... ≤ 30 ppm
Nitrate ..... ≤ 100 ppm
Subst. not ppt. (NH4)2S ..... ≤ 0.2 %
Ca ..... ≤ 50 ppm
Cd ..... ≤ 20 ppm
Co ..... ≤ 10 ppm
Cu ..... ≤ 10 ppm
Fe ..... ≤ 10 ppm
K ..... ≤ 200 ppm
Na ..... ≤ 200 ppm
Pb ..... ≤ 10 ppm
Zn ..... ≤ 20 ppm
Assay (complexometric) ..... 99 ÷ 100 %

Code	Size	Packaging	Notes
464545	250 g	Plastic bottle	
464547	1 kg	Plastic bottle	

Low content in cobalt

Nickel (II) carbonate basic

• Nichel carbonato oso basico • Nickel (II) carbonate basique • Níquel (II) carbonato básico • Nickel (II) carbonat basisch

$NiCO_3.2Ni(OH)_2.4H_2O$
Molecular Weight: 376,23
CAS: 39430-27-8
EEC-N: 235-715-9

Danger
H302-H332-H315-H334-H317-H341-H350i-H360D-H372-H410-HA26
P271-P280-P284-P304+P340-P308+P313-P342+P311a

Nickel (II) carbonate basic > RPE - For analysis

RPE

Description ..... Green powder
Identification ..... Positive
Chloride ..... ≤ 0.1 %
Diluted HCl-ins. matter ..... ≤ 500 ppm
Co ..... ≤ 0.1 %
Cu ..... ≤ 50 ppm
Fe ..... ≤ 100 ppm
Pb ..... ≤ 50 ppm
Zn ..... ≤ 50 ppm
Assay (complexometric) ..... ≥ 45 % (Ni)

Code	Size	Packaging	Notes
464604	100 g	Plastic bottle	
464605	1 kg	Plastic bottle	

Nickel (II) chloride hexahydrate

• Nichel cloruro oso esaidrato • Nickel (II) chlorure hexahydraté • Níquel (II) cloruro hexahidratado • Nickel (II) chloridhexahydrat

$NiCl_2.6H_2O$
Molecular Weight: 237,7
CAS: 7791-20-0
EEC-N: 231-743-0

Classification transport
ONU: 3288
Transport Hazard class: 6.1
Packing group III

Danger
H301-H334-H317-H351
P261-P280-P284-P301+P310a-P304+P340-P342+P311a

Nickel (II) chloride hexahydrate > RPE - For analysis

RPE

Description ..... Green crystals
Identification ..... Positive
Cd ..... ≤ 10 ppm
Co ..... ≤ 10 ppm
Cu ..... ≤ 10 ppm
Fe ..... ≤ 10 ppm
Pb ..... ≤ 5 ppm
Zn ..... ≤ 5 ppm
Assay (argentimetric) ..... ≥ 97 %
Ni ..... ≥ 24 %

Code	Size	Packaging	Notes
464644	100 g	Glass bottle	
464645	250 g	Plastic bottle	
464647	1 kg	Plastic bottle	

Low content in cobalt

Nickel (II) nitrate hexahydrate

• Nichel nitrato oso esaidrato • Nickel (II) nitrate hexahydraté • Níquel (II) nitrato hexahidratado • Nickel (II) nitrathexahydrat

$Ni(NO_3)_2.6H_2O$
Molecular Weight: 290,81
CAS: 13478-00-7
EEC-N: 236-068-5

Classification transport
ONU: 2725
Transport Hazard class: 5.1
Packing group III



Danger
H272-H302-H334-H317-H351
P210-P261-P280-P284-P304+P340-P342+P311a

Nickel (II) nitrate hexahydrate > RPE - For analysis

RPE

Description	Green crystals	As	≤ 5 ppm	Cu	≤ 10 ppm	S	≤ 10 ppm
Identification	Positive	Ca	≤ 10 ppm	Fe	≤ 10 ppm	Zn	≤ 10 ppm
Water-insoluble matter	≤ 100 ppm	Cd	≤ 1 ppm	Mg	≤ 10 ppm	Assay (complexometric)	≥ 98.5 %
Chloride	≤ 50 ppm	Co	≤ 0.5 %	Mn	≤ 20 ppm		
Al	≤ 10 ppm	Cr	≤ 1 ppm	Pb	≤ 1 ppm		

Code	Size	Packaging	Notes
464685	250 g	Plastic bottle	
464686	1 kg	Plastic bottle	

Nickel (II) nitrate 10g/l

• Nichel nitrato oso 10 g/l • Nickel (II) nitrate 10 g/l • Níquel (II) nitrato 10 gIL • Nickel (II)-nitrat 10g/l

Classification transport
ONU: 3098
Transport Hazard class: 5.1
Packing group II



Danger
H314-H334-H317-H341-H350i-H360D-H372-H411-HA26
P280-P284-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P362+P364-P342+P311a

Nickel (II) nitrate 10g/l > RS - Matrix modifiers for AAS-GTA

RS

Code	Size	Packaging	Notes
503197	50 ml	Plastic bottle	Matrix: 1% Nitric acid

Nickel (II) sulfate hexahydrate

• Nichel solfato oso esaidrato • Nickel (II) sulfate hexahydraté • Níquel (II) sulfato hexahidratado • Nickel (II) sulfathexahydrat

$NiSO_4.6H_2O$
Molecular Weight: 262,86
CAS: 10101-97-0



Danger
H302-H332-H315-H334-H317-H341-H350i-H360D-H372-H410-HA26
P271-P280-P284-P304+P340-P308+P313-P342+P311a

Nickel (II) sulfate hexahydrate > RPE - For analysis

RPE

Description	Green crystals	Water-insoluble matter	≤100 ppm	Cu	≤10 ppm	Pb	≤10 ppm
Identification	Positive	Subst. not ppt. (NH4)2S	≤0.1 %	Fe	≤80 ppm	Zn	≤20 ppm
pH sol. 5% at 25° C	2.5 ÷ 6.5	Ca	≤400 ppm	K	≤500 ppm	Assay (complexometric)	≥99 %
Total nitrogen	≤20 ppm	Cd	≤10 ppm	Mn	≤10 ppm		
Chloride	≤10 ppm	Co	≤5 ppm	Na	≤100 ppm		

Code	Size	Packaging	Notes
464775	250 g	Plastic bottle	
464777	1 kg	Plastic bottle	
464772	25 kg	Plastic bucket	

Low content in cobalt

Nickel (II) sulfate hexahydrate > RE - Pure

RE

Description	Green crystals	Water-insoluble matter	≤0.1 %	Assay (complexometric)	98 ÷ 100 %
Identification	Positive	Fe	≤100 ppm		

Code	Size	Packaging	Notes
355757	1 kg	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Nicotinamide

• Nicotinammide • Nicotinamide • Nicotinamida • Niacinamid

Synonym:
• *Vitamin B3*
• *Niacetamide*

$C_6H_6ON_2$
Molecular Weight: 122,13
CAS: 98-92-0
EEC-N: 202-713-4

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Nicotinamide > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP-DAB

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Related substances .....Conform Ph.Eur.
Ready carbonizable substances.....Conform USP-NF
Organic volatile impurities Conform USP-NF
Melting point..... 128 ÷ 131 ° C
pH solution 5% .....6.0 ÷ 7.5
Loss on drying ..... ≤0.5 %
Sulphated ash..... ≤0.1 %
Heavy metals (Pb).....≤30 ppm
Assay (non-aqueos medium) .99.0 ÷ 101.0 % s.s.
Assay (HPLC) .....98.5 ÷ 101.5 % s.s.

Code	Size	Packaging	Notes
392304	100 g	Plastic bottle	
392307	1 kg	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Nicotinic acid

• Acido nicotinico • Acide nicotinique • Acido nicotínico • Nikotinsäure

Synonym:
3-Picolinic acid

N:CHC(COOH):CHCH:CH
Molecular Weight: 123,11
CAS: 59-67-6
EEC-N: 200-441-0

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Nicotinic acid > RPE - For analysis

RPE

Description .....White crystalline powder
Identification ..... Positive
Melting point..... 235.7 ÷ 237.3 ° C
Chloride.....≤0.02 %
Heavy metals (Pb).....≤20 ppm
Loss on drying .....≤1.0 %
Residue on ignition.....≤0.1 %
Assay (non-aqueos medium).....99.5 ÷ 100.5 %
Sulphate.....≤ 0.1 %

Code	Size	Packaging	Notes
407914	100 g	Plastic bottle	

Nigrosine

• Nigrosina • Nigrosine • Nigrosina • Nigrosine

Synonym:
Acid black 2

CAS: 8005-03-6

Warning
H302-H312-H332
P261-P264-P271-P280h-P301+P312a-P304+P340

Nigrosine > RS - For microscopy - C.I. 50420

RS

Description ..... black granules
Identification ..... Positive

Code	Size	Packaging	Notes
464852	25 g	Glass bottle	
464853	50 g	Glass bottle	

Ninhydrin

• Ninidrina • Ninhydrine • Ninhidrina • Ninhydrin

Synonym:
1,2,3-Indantrione monohydrate

$C_6H_4COCOCO.H_2O$
Molecular Weight: 178,15
CAS: 485-47-2
EEC-N: 207-618-1

Warning
H302-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Ninhydrin > RPE - For analysis - ACS

RPE

Description ..... Yellow powder
Identification ..... Positive
Ident.and melting point ..... Conform
Aminoacids sensitivity ..... Conform
Solubility ..... Conform

Code	Size	Packaging	Notes
464928	5 g	Glass bottle	
464922	25 g	Glass bottle	

Ninhydrin solution

• Ninidrina solzione • Ninhydrine solution • Ninhidrina solución • Ninhydrinlösung

Synonym:
1,2,3-Indantrione monohydrate

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group III

Danger
H226-H318
P210-P241-P280-P303+P361+P353-P305+P351+P338-P310a

Ninhydrin solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611058304	100 ml	Glass bottle	Ninhydrin solution R1 Ref Ph.Eur 1058304
611058305	100 ml	Glass bottle	Ninhydrin solution R2 Ref Ph.Eur 1058305

Ninhydrin and Tin (II) chloride reagent

• Ninidrina e stagno cloruro reattivo • Réactif à la ninhydrine et au étain (II) chlorure • Ninhidrina y estaño (II) cloruro reactivo
• Ninhydrin und Zinn (II) chlorid Reagenz

Ninhydrin and Tin (II) chloride reagent > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611058301	100 ml	Glass bottle	Ref Ph.Eur 1058301

Niobium standard solution

• Niobio standard soluzione • Niobium solution standard • Niobio, solución patrón • Niob-Standardlösung

Classification transport
ONU: 2922
Transport Hazard class: 8
Packing group II

Danger
H290-H302-H311-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P361+P364

Niobium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505737	100 ml	Plastic bottle	conc. 10 ppm Matrix: Nitric acid and hydrofluoric acid
505738	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid and hydrofluoric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Product specifications are subject to changes. Please visit our website for updates.

Niobium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503751	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid
503753	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid
503755	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid and hydrofluoric acid
503757	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid and hydrofluoric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Niobium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507750	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid
507511	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Nital solution 4%

• Nital soluzione 4% • Nital solution 4% • Nital solución 4% • Nital-Lösung 4%

Classification transport
ONU: 3316
Transport Hazard class: 9
Packing group II



Danger
H225-H272-H290-H331-H314-HEU071
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Nital solution 4% > RS - Macrography

RS

Code	Size	Packaging	Notes
505021	1 l	Plastic bottle	

Composition: 4ml HNO3 65%; 100ml Ethanol

Nitrate standard solution

• Nitrati standard soluzione • Nitrate solution standard • Nitrato, solución patrón • Nitrat-Standardlösung

Nitrate standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615002101	100 ml	Plastic bottle	A 10 ppm solution: to dilute according to Ref Ph.Eur 5002101
615002102	100 ml	Plastic bottle	A 2 ppm solution: to dilute according to Ref Ph.Eur 5002102
615002109	100 ml	Plastic bottle	A 100 ppm solution: to dilute according to Ref Ph.Eur 5002100

Nitrate standard solution > RS - Standard solution for ion chromatography

RS

Code	Size	Packaging	Notes
503331	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503333	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Nitric acid fuming 99%

• Acido nitrico fumante 99% • Acide nitrique fumant 99% • Acido nítrico fumante 99% • Salpetersäure raucht zu 99%

HNO_3
Molecular Weight: 63,01
CAS: 7697-37-2
EEC-N: 231-714-2

Classification transport
ONU: 2031
Transport Hazard class: 8
Packing group I



Danger
H272-H290-H331-H314-HEU071
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Nitric acid fuming 99% > RPE - For analysis

RPE

Description	Clear colourless liquid	Sulfate	≤ 5 ppm	Cu	≤ 100 ppb	Na	≤ 300 ppb
Identification	Conform	Al	≤ 500 ppb	Fe	≤ 200 ppb	Sn	≤ 200 ppb
Assay (HNO3)	≥ 99 %	Ca	≤ 300 ppb	Mg	≤ 200 ppb	Zn	≤ 500 ppb
Chloride	≤ 1 ppm	Cr	≤ 200 ppb	K	≤ 200 ppb		

Code	Size	Packaging	Notes
408143	500 ml	Glass bottle PVC coated	
408142	2 l	Glass bottle PVC coated	

Nitric acid 69.5%

• Acido nitrico 69,5% • Acide nitrique 69.5% • Acido nítrico 69.5% • Salpetersäure 69.5%

HNO_3
Molecular Weight: 63,01
CAS: 7697-37-2

Classification transport
ONU: 2031
Transport Hazard class: 8
Packing group II



Danger
H272-H290-H331-H314-HEU071
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Nitric acid 69.5% > RS - VLSI - For electronic use

RS

Code	Size	Packaging	Notes
527671	1 l	Plastic bottle	
527670	2.5 l	Plastic bottle	

For specifications, contact our customer service for a certificate of analysis

Nitric acid 69.5% > RS - RSE - For electronic use

RS

Description	Clear colourless liquid	As	≤0.005 ppm	Fe	≤0.1 ppm	Pt	≤0.02 ppm
Identification	Positive	Au	≤0.05 ppm	Ga	≤0.02 ppm	Sb	≤0.01 ppm
Density at 20° C	1.408 ÷ 1.418	B	≤0.01 ppm	In	≤0.02 ppm	Si	≤0.1 ppm
Assay (acidimetric)	69.1 ÷ 69.9 %	Ba	≤0.05 ppm	K	≤0.1 ppm	Sn	≤0.02 ppm
Chloride	≤0.05 ppm	Be	≤0.02 ppm	Li	≤0.02 ppm	Sr	≤0.02 ppm
Phosphate	≤0.1 ppm	Bi	≤0.02 ppm	Mg	≤0.1 ppm	Ta	≤0.1 ppm
Heavy metals (Pb)	≤0.05 ppm	Ca	≤0.1 ppm	Mn	≤0.01 ppm	Ti	≤0.01 ppm
Residue on ignition	≤2 ppm	Cd	≤0.005 ppm	Mo	≤0.05 ppm	Tl	≤0.02 ppm
Sulphate	≤0.5 ppm	Co	≤0.01 ppm	Na	≤0.3 ppm	V	≤0.01 ppm
Ag	≤0.02 ppm	Cr	≤0.01 ppm	Ni	≤0.01 ppm	Zn	≤0.05 ppm
Al	≤0.05 ppm	Cu	≤0.005 ppm	Pb	≤0.02 ppm	Zr	≤0.01 ppm

Code	Size	Packaging	Notes
408097	1 l	Glass bottle	
408098	2.5 l	Glass bottle	

Nitric acid 69.5% > RS - MOS - For electronic use

RS

Description	Clear colourless liquid	As	≤0.005 ppm	Fe	≤0.1 ppm	Pt	≤0.05 ppm
Identification	Positive	Au	≤0.05 ppm	Ga	≤0.02 ppm	Sb	≤0.01 ppm
Density at 20° C	1.408 ÷ 1.418	B	≤0.01 ppm	In	≤0.02 ppm	Sn	≤0.02 ppm
Assay (acidimetric)	69.1 ÷ 69.9 %	Ba	≤0.05 ppm	K	≤0.1 ppm	Sr	≤0.02 ppm
Chloride	≤0.05 ppm	Be	≤0.02 ppm	Li	≤0.02 ppm	Ta	≤0.1 ppm
Phosphate	≤0.1 ppm	Bi	≤0.02 ppm	Mg	≤0.1 ppm	Ti	≤0.05 ppm
Heavy metals (Pb)	≤0.05 ppm	Ca	≤0.1 ppm	Mn	≤0.01 ppm	Tl	≤0.05 ppm
Residue on ignition	≤2 ppm	Cd	≤0.005 ppm	Mo	≤0.05 ppm	V	≤0.05 ppm
Sulphate	≤0.5 ppm	Co	≤0.01 ppm	Na	≤0.3 ppm	Zn	≤0.05 ppm
Ag	≤0.02 ppm	Cr	≤0.01 ppm	Ni	≤0.01 ppm	Zr	≤0.05 ppm
Al	≤0.05 ppm	Cu	≤0.005 ppm	Pb	≤0.02 ppm		

Code	Size	Packaging	Notes
408151	1 l	Glass bottle	
408152	2.5 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Nitric acid 69.5% > RPE - For analysis - ISO - ACS - Reag.Ph.Eur. - Reag.USP

RPE

Description	Clear liquid	Nitrite	≤35 ppm	Cd	≤0.005 ppm	Na	≤0.5 ppm
Colour (APHA)	≤10	Silicate	≤1 ppm	Co	≤0.01 ppm	Ni	≤0.05 ppm
Identification	Positive	Sulphate	≤0.5 ppm	Cr	≤0.1 ppm	Pb	≤0.02 ppm
Hg	≤10 ppb	Ag	≤0.02 ppm	Cu	≤0.01 ppm	Sr	≤0.02 ppm
Density at 20° C	1.408 ÷ 1.416	Al	≤0.05 ppm	Fe	≤0.2 ppm	Ti	≤0.1 ppm
Chloride	≤0.1 ppm	As	≤0.005 ppm	K	≤0.1 ppm	Tl	≤0.05 ppm
Residue on evaporation	≤5 ppm	Ba	≤0.1 ppm	Li	≤0.02 ppm	V	≤0.05 ppm
Phosphate	≤0.5 ppm	Be	≤0.02 ppm	Mg	≤0.1 ppm	Zn	≤0.01 ppm
Heavy metals (Pb)	≤0.2 ppm	Bi	≤0.1 ppm	Mn	≤0.01 ppm	Assay (acidimetric)	69.1 ÷ 69.9 %
Sulphated ash	≤4 ppm	Ca	≤5 ppm	Mo	≤0.05 ppm		

Code	Size	Packaging	Notes
408071	1 l	Glass bottle	
408076	1 l	Glass bottle PVC coated	
524530	1 l	Plastic bottle	
408072	2.5 l	Glass bottle	
524531	2.5 l	Plastic bottle	
408075	34 kg	Plastic drum	

Nitric acid 67-70%

• Acido nitrico 67-70% • Acide nitrique 67-70% • Acido nítrico 67-70% • Salpetersäure 67-70%

HNO_3
Molecular Weight: 63,01
CAS: 7697-37-2

Classification transport
ONU: 2031
Transport Hazard class: 8
Packing group II



Danger
H272-H290-H331-H314-HEU071
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Nitric acid 67-70% > RS - Superpure - For trace analysis at ppb level

RS

Description	Clear liquid	Li	≤ 0.1 ppb	Ce	≤ 0.1 ppb	Pt	≤ 0.5 ppb
Identification	Positive	Mg	≤ 1 ppb	Colour (APHA)	≤ 10	Pr	≤ 0.1 ppb
Chloride	≤ 0.2 ppm	Mn	≤ 0.1 ppb	Cs	≤ 0.1 ppb	Re	≤ 0.1 ppb
Total phosphorus	≤ 0.01 ppm	Hg	≤ 0.1 ppb	Dy	≤ 0.1 ppb	Rh	≤ 0.5 ppb
Total sulphur	≤ 0.3 ppm	Mo	≤ 0.1 ppb	Er	≤ 0.1 ppb	Rb	≤ 0.1 ppb
Co	≤ 0.5 ppb	Ni	≤ 0.5 ppb	Eu	≤ 0.1 ppb	Ru	≤ 0.5 ppb
Sb	≤ 0.5 ppb	Se	≤ 1 ppb	Gd	≤ 0.1 ppb	Sm	≤ 0.1 ppb
Fe	≤ 1 ppb	Ag	≤ 0.1 ppb	Ga	≤ 0.1 ppb	Sc	≤ 0.1 ppb
Ba	≤ 0.1 ppb	Na	≤ 1 ppb	Ge	≤ 0.1 ppb	Te	≤ 0.1 ppb
Be	≤ 0.1 ppb	Sr	≤ 0.1 ppb	Au	≤ 0.1 ppb	Tb	≤ 0.1 ppb
Bi	≤ 0.1 ppb	Th	≤ 0.1 ppb	Hf	≤ 0.1 ppb	Tl	≤ 0.1 ppb
B	≤ 1 ppb	Sn	≤ 0.5 ppb	Ho	≤ 0.1 ppb	Tm	≤ 0.1 ppb
Cd	≤ 0.5 ppb	Ti	≤ 0.5 ppb	In	≤ 0.1 ppb	W	≤ 0.1 ppb
Ca	≤ 1 ppb	U	≤ 0.1 ppb	La	≤ 0.1 ppb	Yb	≤ 0.1 ppb
Cr	≤ 1 ppb	V	≤ 0.5 ppb	Lu	≤ 0.1 ppb	Y	≤ 0.1 ppb
K	≤ 1 ppb	Zn	≤ 0.5 ppb	Nd	≤ 0.1 ppb		
Cu	≤ 0.5 ppb	Zr	≤ 0.1 ppb	Nb	≤ 0.1 ppb		
Pb	≤ 0.1 ppb	Assay (acidimetric)	67 ÷ 70 %	Pd	≤ 0.5 ppb		

Code	Size	Packaging	Notes
408115	500 ml	Plastic bottle	
408116	1 l	Plastic bottle	
408117	2.5 l	Plastic bottle	

Nitric acid 67-69%

• Acido nitrico 67-69% • Acide nitrique 67-69% • Acido nítrico 67-69% • Salpetersäure 67-69%

HNO_3
Molecular Weight: 63,01
CAS: 7697-37-2

Classification transport
ONU: 2031
Transport Hazard class: 8
Packing group II



Danger
H272-H290-H331-H314-HEU071
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Nitric acid 67-69% > RS - Ultrapure - For trace analysis at ppt level

RS

Description	Clear colourless liquid	Mg	≤ 10 ppt	Ce	≤ 10 ppt	Pd	≤ 20 ppt
Identification	Positive	Mn	≤ 10 ppt	Cs	≤ 10 ppt	Pt	≤ 20 ppt
Ag	≤ 10 ppt	Mo	≤ 10 ppt	Dy	≤ 1 ppt	Pr	≤ 1 ppt
Al	≤ 20 ppt	Na	≤ 10 ppt	Er	≤ 1 ppt	Re	≤ 10 ppt
As	≤ 20 ppt	Ni	≤ 20 ppt	Eu	≤ 1 ppt	Rh	≤ 10 ppt
B	≤ 10 ppt	Pb	≤ 10 ppt	Gd	≤ 1 ppt	Rb	≤ 10 ppt
Ba	≤ 10 ppt	Sn	≤ 20 ppt	Ga	≤ 10 ppt	Ru	≤ 20 ppt
Be	≤ 10 ppt	Sr	≤ 10 ppt	Ge	≤ 10 ppt	Sm	≤ 1 ppt
Bi	≤ 10 ppt	Tl	≤ 10 ppt	Au	≤ 20 ppt	Sc	≤ 10 ppt
Ca	≤ 10 ppt	Ti	≤ 10 ppt	Hf	≤ 10 ppt	Te	≤ 1 ppt
Cd	≤ 10 ppt	V	≤ 10 ppt	Ho	≤ 1 ppt	Tb	≤ 1 ppt
Co	≤ 10 ppt	Zn	≤ 10 ppt	In	≤ 1 ppt	Tm	≤ 1 ppt
Cr	≤ 10 ppt	Zr	≤ 10 ppt	La	≤ 1 ppt	W	≤ 10 ppt
Cu	≤ 10 ppt	Assay (acidimetric)	67 ÷ 69 %	Li	≤ 10 ppt	Yb	≤ 1 ppt
Fe	≤ 10 ppt	U	≤ 1 ppt	Lu	≤ 1 ppt	Y	≤ 1 ppt
Hg	≤ 50 ppt	Th	≤ 1 ppt	Nd	≤ 1 ppt		
K	≤ 10 ppt	Sb	≤ 10 ppt	Nb	≤ 1 ppt		

Code	Size	Packaging	Notes
408052	250 ml	Plastic bottle	
408051	500 ml	Plastic bottle	

Nitric acid 67.5 (42° Be)

• Acido nitrico 67.5% (42° Be) • Acide nitrique 67.5% (42° Be) • Acido nítrico 67.5% (42° Be) • Salpetersäure 67.5% (42° Be)

HNO_3
Molecular Weight: 63,01
CAS: 7697-37-2

Classification transport
ONU: 2031
Transport Hazard class: 8
Packing group II



Danger
H272-H290-H331-H314-HEU071
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Nitric acid 67.5 (42° Be) > RPE - For analysis

RPE

Assay (acidimetric)	65.5 - 69.7 %	Heavy metals (Pb)	≤ 0.2 ppm	Co	≤ 0.01 ppm	Na	≤ 0.2 ppm
Identification (I.R.)	Positive	Ag	≤ 0.02 ppm	Cr	≤ 0.1 ppm	Ni	≤ 0.02 ppm
Colour	≤ 10 APHA	Al	≤ 0.05 ppm	Cu	≤ 0.01 ppm	Pb	≤ 0.02 ppm
Density at 20°C	1.39 - 1.42	As	≤ 0.005 ppm	Fe	≤ 0.2 ppm	Sr	≤ 0.02 ppm
Residue on evaporation	≤ 4 ppm	Ba	≤ 0.1 ppm	K	≤ 0.05 ppm	Ti	≤ 0.01 ppm
Silicate	≤ 1 ppm	Be	≤ 0.02 ppm	Li	≤ 0.02 ppm	Tl	≤ 0.02 ppm
Chloride	≤ 0.1 ppm	Bi	≤ 0.1 ppm	Mg	≤ 0.05 ppm	V	≤ 0.01 ppm
Sulphate	≤ 0.5 ppm	Ca	≤ 0.5 ppm	Mn	≤ 0.02 ppm	Zn	≤ 0.05 ppm
Phosphate	≤ 0.5 ppm	Cd	≤ 0.005 ppm	Mo	≤ 0.02 ppm		

Code	Size	Packaging	Notes
528530	5 l	Tank	

Nitric acid 67.5 (42° Be) > RE - Pure

RE

Description	Clear colourless liquid	Residue on ignition	≤50 ppm	Heavy metals (Pb)	≤10 ppm
Identification	Positive	Chloride	≤30 ppm	Fe	≤10 ppm
Density at 20° C	1.395 ÷ 1.415	Sulphate	≤50 ppm	Assay (acidimetric)	65.0 ÷ 70.0 %

Code	Size	Packaging	Notes
305502	2.5 l	Glass bottle	
305501	40 kg	Plastic tank	
305505	70 kg	Plastic tank	

Product specifications are subject to changes. Please visit our website for updates.

Nitric acid 65%

• Acido nitrico 65% • Acide nitrique 65% • Acido nítrico 65% • Salpetersäure 65%

HNO_3
Molecular Weight: 63,01
CAS: 7697-37-2

Classification transport
ONU: 2031
Transport Hazard class: 8
Packing group II

Danger
H272-H290-H331-H314-HEU071
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Nitric acid 65% > RS - RSE - For electronic use

RS

Description	Clear colourless liquid	As	≤0.005 ppm	Fe	≤0.1 ppm	Pt	≤0.05 ppm
Identification	Positive	Au	≤0.05 ppm	Ga	≤0.02 ppm	Sb	≤0.01 ppm
Density at 20° C	1.390 ÷ 1.410	B	≤0.01 ppm	In	≤0.02 ppm	Sn	≤0.02 ppm
Assay (acidimetric)	64.0 ÷ 66.0 %	Ba	≤0.05 ppm	K	≤0.1 ppm	Sr	≤0.02 ppm
Chloride	≤0.05 ppm	Be	≤0.02 ppm	Li	≤0.02 ppm	Ta	≤0.1 ppm
Phosphate	≤0.1 ppm	Bi	≤0.02 ppm	Mg	≤0.1 ppm	Ti	≤0.05 ppm
Heavy metals (Pb)	≤0.05 ppm	Ca	≤0.1 ppm	Mn	≤0.01 ppm	Tl	≤0.05 ppm
Residue on ignition	≤2 ppm	Cd	≤0.005 ppm	Mo	≤0.05 ppm	V	≤0.05 ppm
Sulphate	≤0.5 ppm	Co	≤0.01 ppm	Na	≤0.3 ppm	Zn	≤0.05 ppm
Ag	≤0.02 ppm	Cr	≤0.01 ppm	Ni	≤0.01 ppm	Zr	≤0.05 ppm
Al	≤0.05 ppm	Cu	≤0.005 ppm	Pb	≤0.02 ppm		

Code	Size	Packaging	Notes
408101	1 l	Glass bottle	
408102	2.5 l	Glass bottle	

Nitric acid 65% > RS - For enviromental analysis - ISO

RS

Description	Clear colourless liquid	Ag	≤0.02 ppm	Cr	≤0.1 ppm	Na	≤0.5 ppm
Identification	Positive	Al	≤0.05 ppm	Cu	≤0.01 ppm	Ni	≤0.05 ppm
Density at 20° C	1.390 ÷ 1.410	As	≤0.005 ppm	Fe	≤0.2 ppm	Pb	≤0.02 ppm
Chloride	≤0.1 ppm	Ba	≤0.1 ppm	Hg	≤0.005 ppm	Sr	≤0.02 ppm
Phosphate	≤0.5 ppm	Be	≤0.02 ppm	K	≤0.1 ppm	Ti	≤0.1 ppm
Heavy metals (Pb)	≤0.2 ppm	Bi	≤0.1 ppm	Li	≤0.02 ppm	Tl	≤0.05 ppm
Residue on ignition	≤4 ppm	Ca	≤0.5 ppm	Mg	≤0.1 ppm	V	≤0.05 ppm
Silicate	≤1 ppm	Cd	≤0.005 ppm	Mn	≤0.01 ppm	Zn	≤0.05 ppm
Sulphate	≤0.5 ppm	Co	≤0.01 ppm	Mo	≤0.05 ppm	Assay (acidimetric)	64 ÷ 66 %

Code	Size	Packaging	Notes
407951	1 l	Glass bottle	
407952	2.5 l	Glass bottle	

Low content in Hg

Nitric acid 65% > RPE - For analysis - ISO

RPE

Description	Clear colourless liquid	Ag	≤0.02 ppm	Cr	≤0.1 ppm	Ni	≤0.05 ppm
Identification	Positive	Al	≤0.05 ppm	Cu	≤0.01 ppm	Pb	≤0.02 ppm
Density at 20° C	1.390 ÷ 1.410	As	≤0.005 ppm	Fe	≤0.2 ppm	Sr	≤0.02 ppm
Chloride	≤0.1 ppm	Ba	≤0.1 ppm	K	≤0.1 ppm	Ti	≤0.1 ppm
Phosphate	≤0.5 ppm	Be	≤0.02 ppm	Li	≤0.02 ppm	Tl	≤0.05 ppm
Heavy metals (Pb)	≤0.2 ppm	Bi	≤0.1 ppm	Mg	≤0.1 ppm	V	≤0.05 ppm
Residue on ignition	≤4 ppm	Ca	≤0.5 ppm	Mn	≤0.01 ppm	Zn	≤0.05 ppm
Silicate	≤1 ppm	Cd	≤0.005 ppm	Mo	≤0.05 ppm	Assay (acidimetric)	64 ÷ 66 %
Sulphate	≤0.5 ppm	Co	≤0.01 ppm	Na	≤0.5 ppm		

Code	Size	Packaging	Notes
408021	1 l	Glass bottle PVC coated	
408022	1 l	Glass bottle	
524535	1 l	Plastic bottle	
408025	2.5 l	Glass bottle	
524536	2.5 l	Plastic bottle	
408027	34 kg	Plastic drum	

Nitric acid 65% > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 1.390 ÷ 1.410
Chloride ..... ≤10 ppm
Heavy metals (Pb) ..... ≤10 ppm
Residue on ignition ..... ≤500 ppm
Sulphate ..... ≤100 ppm
As ..... ≤1 ppm
Fe ..... ≤50 ppm
Assay (acidimetric) ..... 64 ÷ 66 %

Code	Size	Packaging	Notes
305201	1 l	Glass bottle	
305207	2.5 l	Glass bottle	
305202	34 kg	Plastic drum	
305205	34 kg	Drum	

Nitric acid 18%

• Acido nitrico 18% • Acide nitrique 18% • Acido nítrico 18% • Salpetersäure 18%

HNO_3
Molecular Weight: 63,01
CAS: 7697-37-2

Classification transport
ONU: 2031
Transport Hazard class: 8
Packing group II



Danger
H290-H332-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Nitric acid 18% > RS - RSE - For electronic use

RS

Assay ..... 17 - 19 %

Code	Size	Packaging	Notes
408191	1 l	Plastic bottle	

Nitric acid 10%

• Acido nitrico 10% • Acide nitrique 10% • Acido nítrico 10% • Salpetersäure 10%

HNO_3
Molecular Weight: 63,01
CAS: 7697-37-2

Classification transport
ONU: 2031
Transport Hazard class: 8
Packing group II



Danger
H290-H332-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Nitric acid 10% > RS - For analysis

RS

Assay ..... 9 - 11 %
Density d20/4 ..... 1.048 - 1.06

Code	Size	Packaging	Notes
PS0568/41	10 l	Plastic tank	

Nitric acid 8 mol/l (8N)

• Acido nitrico 8 mol/l (8N) • Acide nitrique 8 mol/l (8N) • Acido nítrico 8 mol/l (8N) • Salpetersäure 8 mol/l (8N)

HNO_3
Molecular Weight: 63,01
CAS: 7697-37-2

Classification transport
ONU: 2031
Transport Hazard class: 8
Packing group II



Danger
H290-H331-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Nitric acid 8 mol/l (8N) > RPE - For analysis

RPE

Assay (potentiometry) ..... 7.984 - 8.016 N

Code	Size	Packaging	Notes
PS0311/20	2.5 l	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Nitric acid 2 mol/l (2N)

• Acido nitrico 2 mol/l (2N) • Acide nitrique 2 mol/l (2N) • Acido nítrico 2 mol/l (2N) • Salpetersäure 2 mol/l (2N)

HNO_3
Molecular Weight: 63,01
CAS: 7697-37-2

Classification transport
ONU: 2031
Transport Hazard class: 8
Packing group II



Danger
H290-H332-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Nitric acid 2 mol/l (2N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 1.998 - 2.002 N NIST 723 ..... e

Code	Size	Packaging	Notes
408185000	5 l	Plastic tank	Certified with NIST traceability

Volumetric solution ready-to-use

Nitric acid 1 mol/l (1N)

• Acido nitrico 1 mol/l (1N) • Acide nitrique 1 mol/l (1N) • Acido nítrico 1 mol/l (1N) • Salpetersäure 1 mol/l (1N)

HNO_3
Molecular Weight: 63,01
CAS: 7697-37-2

Classification transport
ONU: 2031
Transport Hazard class: 8
Packing group II



Danger
H290-H332-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Nitric acid 1 mol/l (1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.999 - 1.001 N NIST 723 ..... e

Code	Size	Packaging	Notes
408176000	500 ml	Plastic bottle	Certified with NIST traceability
408171000	1 l	Plastic bottle	Certified with NIST traceability

Volumetric solution ready-to-use

Nitric acid 0.1 mol/l (0.1N)

• Acido nitrico 0.1 mol/l (0.1N) • Acide nitrique 0.1 mol/l (0.1N) • Acido nítrico 0.1 mol/l (0.1N) • Salpetersäure 0.1 mol/l (0.1N)

HNO_3
Molecular Weight: 63,01
CAS: 7697-37-2

Nitric acid 0.1 mol/l (0.1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.0999 - 0.1001 N NIST 723 ..... e

Code	Size	Packaging	Notes
408206000	500 ml	Plastic bottle	Certified with NIST traceability

6.301 g of HNO3. Volumetric solution ready-to-use

Nitric acid 0.1 mol/l (0.1N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
408231		Glass ampoule	Volume: 55 ml

6,301 g HNO3. Volumetric concentrated solution to prepare 1 L of solution 0,1 N

Nitric acid, dilute

• Acido nitrico diluito • Acide nitrique dilué • Acido nítrico diluido • Salpetersäure verdünnt

HNO_3
Molecular Weight: 63,01
CAS: 7697-37-2

Classification transport
ONU: 2031
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Nitric acid, dilute > RS - For analysis according to Ph. Eur. Chap. 4.1.1 — RS

Code	Size	Packaging	Notes
611058402	100 ml	Plastic bottle	Ref Ph.Eur 1058402
611058409	250 ml	Plastic bottle	Ref Ph.Eur 1058402

Nitric acid cerium salt ▶ Cerium (III) nitrate hexahydrate

Nitrilotriacetic acid

• Acido nitrilotriacetico • Acide nitrilotriacétique • Acido nitrilotriacético • Nitrilotriessigsäure

Synonym:
• *Nitrilotriacetic acid*
• *NTA*

$N(CH_2COOH)_3$
Molecular Weight: 191,15
CAS: 139-13-9
EEC-N: 205-355-7

Warning
H302
P264-P270-P301+P312a-P330-P501a

Nitrilotriacetic acid > RPE - For analysis — RPE

Description ..... White powder Identification (I.R.) ..... Conform Assay (complexometric) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
408242	100 g	Glass bottle	

Nitrite standard solution

• Nitriti standard soluzione • Nitrite solution standard • Nitrito, solución patrón • Nitrit-Standardlösung

Nitrite standard solution > RS - Standard solution for ion chromatography — RS

Code	Size	Packaging	Notes
503321	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503323	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

m-Nitrobenzaldehyde

• m-Nitrobenzaldeide • m-Nitrobenzaldéhyde • m-Nitrobenzaldehyde • m-Nitrobenzaldehyd

$NO_2C_6H_4CHO$
Molecular Weight: 151,12
CAS: 99-61-6
EEC-N: 202-772-6

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

m-Nitrobenzaldehyde > RE - Pure — RE

Description ..... Yellowish crystalline powder Melting point ..... 56 ÷ 59 ° C Assay (GLC) ..... ≥ 98.5 %
Identification ..... Positive Water ..... ≤ 0.5 %

Code	Size	Packaging	Notes
465142	25 g	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Nitrobenzene

• Nitrobenzene • Nitrobenzène • Nitrobenceno • Nitrobenzol

$C_6H_5NO_2$
Molecular Weight: 123,11
CAS: 98-95-3
EEC-N: 202-716-0

Classification transport
ONU: 1662
Transport Hazard class: 6.1
Packing group II



Danger
H301-H311-H331-H351-H360F-H372-H412-HA26
P271-P280-P304+P340-P308+P313-P330-P361+P364-P403+P233

Nitrobenzene > RPE - For analysis - ACS

RPE

Description ..... Yellow clear liquid | Identification ..... Positive | Assay (GLC) ..... ≥ 98.5 % | Refractive index at 20°C. 1.5500 ÷ 1.5530

Code	Size	Packaging	Notes
465222	1 l	Glass bottle	

m-Nitrobenzoic acid

• Acido m-nitrobenzoico • Acide m-nitrobenzoique • Acido m-nitrobenzoico • m-Nitrobenzoesäure

Synonym:
3-Nitrobenzoic acid

$NO_2C_6H_4COOH$
Molecular Weight: 167,12
CAS: 121-92-6
EEC-N: 204-508-5



Warning
H302
P264-P270-P301+P312a-P330-P501a

m-Nitrobenzoic acid > RPE - For analysis

RPE

Description ..... Light yellow to light green powder or crystals | Identification (I.R.) ..... Positive | Melting point ..... 139 - 143 °C | Assay (GC) ..... ≥ 97.5 %

Code	Size	Packaging	Notes
408414	100 g	Glass bottle	

p-Nitrophenol

• p-Nitrofenolo • p-Nitrophénol • p-Nitrofenol • p-Nitrophenol

Synonym:
4-Nitrophenol

$NO_2C_6H_4OH$
Molecular Weight: 139,11
CAS: 100-02-7
EEC-N: 202-811-7

Classification transport
ONU: 1663
Transport Hazard class: 6.1
Packing group III



Warning
H302-H312-H332-H373
P260-P264-P271-P280h-P301+P312a-P304+P340

p-Nitrophenol > RPE - For analysis

RPE

Description ..... Yellow crystals | Identification ..... Positive | Sensitivity (pH 5.6-7.6) ..... Conform | Melting point ..... 109 ÷ 114 ° C | Water ..... ≤ 0.9 % | Colour change ..... Colourless-yellow

Code	Size	Packaging	Notes
465744	100 g	Glass bottle	

Acid-base indicator

p-Nitrophenol solution 0,1% in water

• p-Nitrofenolo soluzione 0.1% in acqua • p-Nitrophénol solution 0.1% dans l'eau
• p-Nitrofenol solución 0.1% en agua • p-Nitrophenol-Lösung 0.1% in Wasser

Synonym:
4-Nitrophenol

p-Nitrophenol solution 0,1% in water > RPE - For analysis

RPE

Description ..... Yellow clear liquid | Identification ..... Positive | pH range ..... 5.0 - 7.0

Code	Size	Packaging	Notes
E465776	500 ml	Bottle	

Nonylphenol ethoxylated 10 ETO

• Nonilfenolo etossilato 10 ETO • Nonylphénol éthoxylate 10 ETO • Nonilfenol etoxilato 10 ETO
• Nonylphenol ethoxylated - 10 ETO

Synonym:
4-Nonylphenyl-polyethylene glycol

CAS: 9016-45-9
EEC-N: 500-024-6



Warning
H302-H315-H319-H411
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Nonylphenol ethoxylated 10 ETO > RS - Standard for detection of surfactants

RS

Description .....Yellow liquid Identification ..... Positive Assay (at production) .....≥90 %

Code	Size	Packaging	Notes
466361	10 g	Glass bottle	

Nuclear fast red

• Rosso solido nucleare • Rouge nucléaire • Rojo nuclear sólido • Nuklear schnell rot

Synonym:
4-Amino-9,10-dihydro-1,3-dihydroxy-9,10-dioxo-2-anthracenesulfonic acid sodium salt

$C_{14}H_8NNaO_7S$
Molecular Weight: 357,28
CAS: 6409-77-4
EEC-N: 229-088-0

Nuclear fast red > RS - For microscopy - C.I. 60760

RS

Description .....Red brown powder Identification ..... Positive

Code	Size	Packaging	Notes
477011	10 g	Glass bottle	
477012	25 g	Glass bottle	

Dye for cytology

Product specifications are subject to changes. Please visit our website for updates.

n-Octane
• n-Ottano • n-Octane • n-Octano • n-Octan

$CH_3(CH_2)_6CH_3$
Molecular Weight: 114,23
CAS: 111-65-9
EEC-N: 203-892-1

Classification transport
ONU: 1262
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H336-H304-H410
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

n-Octane > RPE - For analysis - Reag. Ph. Eur.

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Ready carbonizable substances..... Conform
Density at 20° C ..... 0.697 ÷ 0.707
Refractive index at 20°C. 1.3940 ÷ 1.4010
Boiling point..... 125.0 ÷ 126.0 ° C
Water (K.F.) ..... ≤100 ppm
Residue on evaporation ..... ≤10 ppm
Acidity (caprilic acid)..... ≤17 ppm
Alcalinity (NH3)..... ≤0.2 ppm
Subst. reducing KMnO4..... ≤20 ppm (5m)
Total sulphur ..... ≤50 ppm
Assay (GLC) ..... ≥ 99 %

Code	Size	Packaging	Notes
467562	1 l	Glass bottle	

n-Octane > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Colour ..... ≤ 10 APHA
Density at 20° C ..... 0.697 ÷ 0.707
Refractive index at 20°C. 1.3925 ÷ 1.4025
Boiling point..... 124.5 ÷ 126.5 ° C
Residue on evaporation ..... ≤30 ppm
Acidity (caprilic acid)..... ≤50 ppm
Total sulphur ..... ≤50 ppm
Assay (GLC) ..... ≥95 %

Code	Size	Packaging	Notes
356661	1 l	Glass bottle	
356663	2.5 l	Glass bottle	

Octane 80 blend
• Miscela ottano 80 • Mélange octane 80 • Mezcla octano 80 • Octan 80-Mischung

$CH_3(CH_2)_6CH_3$
Molecular Weight: 114,23
CAS: 111-65-9

Classification transport
ONU: 3295
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H336-H304-H410
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Octane 80 blend > RE - ASTM

RE

Clear,colourless liquid..... Conform
Isooctane (ASTM) ..... 79.9 - 80.1 % (V/V)
n-Heptane (ASTM)..... 19.9 - 20.1 % (V/V)

Code	Size	Packaging	Notes
525992	5 l	Plastic tank	
525993	25 l	Metal drum	
525994	140 kg	Drum	

Suitable for ASTM methods D2700 and D2699

1-Octanesulphonic acid sodium salt
• Acido 1-ottansolfonico sale sodico • Acide 1-octanesulfonique sel sodique • Acido 1-octanosulfónico sal sódica • 1-Octansulfonsäure-Natriumsalz

$C_8H_{17}NaO_3S$
Molecular Weight: 216,28
CAS: 5324-84-5
EEC-N: 226-195-4

1-Octanesulphonic acid sodium salt > RS - For ion pair chromatography

RS

Description ..... White crystalline powder
Water (K.F.) ..... ≤ 2 %
Assay ..... ≥ 98 %
Absorbance (0,25M)
At 200 nm ..... ≤ 0.10 AU
At 210 nm ..... ≤ 0.05 AU
At 220 nm ..... ≤ 0.04 AU
At 230 nm ..... ≤ 0.03 AU
At 240 nm ..... ≤ 0.01 AU
At 250 nm ..... ≤ 0.01 AU
At 260 nm ..... ≤ 0.01 AU

Code	Size	Packaging	Notes
405861	25 g	Glass bottle	
405862	100 g	Plastic bottle	
405863	1 kg	Plastic bottle	

1-Octanesulfonic acid sodium salt monohydrate

• Acido 1-octansolfonico sale sodico monoidrato • Acide octanesulfonique sel sodique monohydraté
• Acido 1-octanosulfónico sal sódica monohidrato • 1-Octansulfonsäure-Natriumsalz-Monohydrat

Synonym:
Sodium 1-octanesulfonate monohydrate

$CH_3(CH_2)_7SO_3Na.H_2O$
Molecular Weight: 234,29
CAS: 207596-29-0

1-Octanesulfonic acid sodium salt monohydrate > RS - For ion pair chromatography

RS

Description ..... White crystalline powder
Loss on drying ..... 7.0 - 9.0 %
Assay ..... ≥ 99.0 %

Absorbance (0,25M)
At 200 nm ..... ≤ 0.3 AU
At 210 nm ..... ≤ 0.2 AU
At 220 nm ..... ≤ 0.1 AU
At 230 nm ..... ≤ 0.075 AU
At 240 nm ..... ≤ 0.05 AU
At 250 nm ..... ≤ 0.05 AU
At 260 nm ..... ≤ 0.04 AU

Code	Size	Packaging	Notes
405931	25 g	Glass bottle	
405932	100 g	Plastic bottle	

Octanoic acid ▶ n-Caprylic acid

Octanol-1

• Alcole n-ottilico • Octanol-1 • Octanol-1 • 1-Octanol

Synonym:
• *1-Octanol*
• *Capryl alcohol*

$CH_3(CH_2)_6CH_2OH$
Molecular Weight: 130,23
CAS: 111-87-5
EEC-N: 203-917-6

Warning
H302-H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Octanol-1 > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive

Density at 25° C ..... 0.815 ÷ 0.830
Refractive index at 20°C ..... 1.425 ÷ 1.440

Boiling point ..... 194.0 ÷ 196.0 °C
Melting point ..... -16 ÷ -14 °C

Assay (GLC) ..... ≥ 99 %

Code	Size	Packaging	Notes
415002	100 ml	Glass bottle	
415003	1 l	Glass bottle	
415004	30 l	Plastic drum	

Oil of cedar wood

• Olio di legno cedro condensato • Huile de cèdre condensée • Aceite de cedro condensado
• Öl aus Zedernholz

Synonym:
• *Cedar oil*
• *Cedarwood oil*

CAS: 8002-27-9

Warning
H317-H412
P261-P280g-P302+P352a-P362+P364-P333+P313-P501a

Oil of cedar wood > RS - For microscopy

RS

Description ..... Yellow colourless liquid
Identification ..... Positive
Density at 20°C ..... 0.990 ÷ 1.010
Refractive index at 20°C ..... 1.515 - 1.520

Code	Size	Packaging	Notes
466753	100 ml	Glass bottle	
466757	1 l	Glass bottle	

Immersion medium for microscopy

Oil Red O ▶ Red for oils O

Product specifications are subject to changes. Please visit our website for updates.

Oil refined of almonds

• Olio di mandorla raffinato • Huile d'amande raffinée • Aceite refinado de almendras • Aus Mandeln raffiniertes Öl

Synonym:
Almond oil

CAS: 8007-69-0

Oil refined of almonds > ERBApharm - According to pharmacopoeia: NF

ERBApharm

Description ..... Yellow colourless liquid
Identification ..... Positive
Relative density ..... 0.910 ÷ 0.915
Acid value ..... ≤ 0.5
Peroxide value ..... ≤ 5.0
Not saponificable matt. ..... ≤ 0.9 %
Composition of fatty acids (GC). Conform NF
Sterol composition ..... Conform NF

Code	Size	Packaging	Notes
356251	1 l	Glass bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Oleic acid

• Acido oleico • Acide oléique • Acido oleico • Olsäure

Synonym:
• *cis-9-Octadecenoic acid*
• *Elainic acid*

$CH_3(CH_2)_7CH{:}CH(CH_2)_7COOH$
Molecular Weight: 282,45
CAS: 112-80-1
EEC-N: 204-007-1

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Oleic acid > RE - Pure

RE

Description ..... Clear yellow liquid
Identification ..... Positive
Density at 20° C ..... 0.890 ÷ 0.910
Acid value ..... ≥ 195 mg KOH / g
Iodine value ..... ≥ 89 g / 100g

Code	Size	Packaging	Notes
305704	1 l	Glass bottle	
305701	24 kg	Metal drum	

Orange G

• Arancio G • Orange G • Naranja G • Orange G

Synonym:
• *1-Phenylazo-2-naphthol-6,8-disulfonic acid disodium salt*
• *Acid orange 10*

$C_{16}H_{10}N_2Na_2O_7S_2$
Molecular Weight: 452,36
CAS: 1936-15-8
EEC-N: 217-705-6

Orange G > RS - For microscopy - C.I. 16230

RS

Description ..... Red-orange powder
Identification ..... Positive
Maximum absorption ..... 470 ÷ 520 nm
Loss on drying at 110°C ..... ≤ 10.00 %

Code	Size	Packaging	Notes
423432	25 g	Glass bottle	

Dye for cytology

Orange II

• Arancio II • Orange II • Naranja II • Orange II

Synonym:
Acid Orange 7

$C_{16}H_{11}N_2NaO_4S$
Molecular Weight: 350,33
CAS: 633-96-5
EEC-N: 211-199-0

Orange II > RPE - For analysis - C.I. 15510

RPE

Description ..... Orange red powder
Identification ..... Positive

Code	Size	Packaging	Notes
423341	10 g	Glass bottle	

Dye for microscopy (histology). Indicator acid - base (pH 11.0 ÷ 13.0). For the extraction and determination of cationic surfactants

Orcein

• Orceina • Orcéine • Orceina • Orcein

Synonym:
Natural Red 28

CAS: 1400-62-0
EEC-N: 215-750-6

Warning
H302
P264-P270-P301+P312a-P330-P501a

Orcein > RS - For microscopy - C.I. Natural Red 28

RS

Description ..... Brown powder Identification ..... Positive

Code	Size	Packaging	Notes
466858	5 g	Glass bottle	
466859	25 g	Glass bottle	

Dye for botanical and histology

Orcinol monohydrate

• Orcina monoidrato • Orcinol • Orcina monohidrato • Orcinol monohydrate

Synonym:
3,5-Dihydroxytoluene|5-Methylresorcinol

$C_7H_8O_2.H_2O$
Molecular Weight: 142,15
CAS: 6153-39-5

Warning
H302-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Orcinol monohydrate > RPE - For analysis

RPE

Description ..... White powder or pinkish
Identification ..... Positive
Melting point ..... 56 ÷ 61 °C
Water (K.F.) ..... 10 ÷ 14 %
Residue on ignition ..... ≤ 0.1 %
Assay (GLC) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
466908	5 g	Glass bottle	

Organic standard PCB

• Standard organico PCB • Standard organique PCB • Patrones orgánicos PCB • Organisch standard: PCB

Organic standard PCB > RS - For environmental analysis

RS

Code	Size	Packaging	Notes
507127	1 ml	Glass ampoule	PCB 29 [15862-07-4] 10µg/ml in isooctane
507128	1 ml	Glass ampoule	PCB 30 [35693-92-6] 10µg/ml in isooctane
507129	1 ml	Glass ampoule	PCB 73 [74338-23-1] 10µg/ml in isooctane
507131	1 ml	Glass ampoule	PCB 89 [73575-57-2] 10µg/ml in isooctane
507132	1 ml	Glass ampoule	PCB 90 [68194-07-0] 10µg/ml in isooctane
507133	1 ml	Glass ampoule	PCB 106 [70424-69-0] 10µg/ml in isooctane
507134	1 ml	Glass ampoule	PCB 164 [74472-45-0] 10µg/ml in isooctane
507135	1 ml	Glass ampoule	PCB 143 [68194-15-0] 10µg/ml in isooctane
507136	1 ml	Glass ampoule	PCB 155[33979-03-2] 10µg/ml in isooctane
507137	1 ml	Glass ampoule	PCB 198 [68194-17-2] 10µg/ml in isooctane
507138	1 ml	Glass ampoule	PCB 207 [52663-79-3] 10µg/ml in isooctane
507139	1 ml	Glass ampoule	PCB 209 [2051-24-3] 10µg/ml in isooctane
507154	1 ml	Glass ampoule	PCB 209 [2051-24-3] 100µg/ml in isooctane

Custom formulations of organic substances are available. Contact us for more details.

Product specifications are subject to changes. Please visit our website for updates.

Organic Standard: PAH multielement mixture

• Standard organico: PAH mix • Standard organique: mélange multiéléments HAP • Patrones orgánicos : mezcla multi-HAP
• Organisch standard: PAH Multielement-Mischung

Organic Standard: PAH multielement mixture > RS - For environmental analysis

RS

Code	Size	Packaging	Notes
507063	5 x 1 ml	Glass ampoule	15 components 50ug/ml each in acetonitrile: Acenaphthene [CAS:83-32-9]; Anthracene [CAS:120-12-7]; Benzo(a)anthracene [CAS:56-55-3]; Benzo(a)pyrene [CAS:50-32-8]; Benzo(b) fluoranthene [CAS:205-99-2]; Benzo(g,h,i)perylene [CAS:191-24-2]; Benzo(k)fluoranthene [CAS:207-08-9]; Chrysene [CAS:218-01-9]; Dibenzo(a,h)anthracene [CAS:53-70-3]; Fluoranthene [CAS:206-44-0]; Fluorene [CAS:86-73-7]; Indeno(1,2,3-c,d) pyrene [CAS:193-39-5]; Naphthalene [CAS:91-20-3]; Phenanthrene [CAS:85-01-8]; Pyrene [CAS:129-00-0]
506835	10 x 1 ml	Glass ampoule	13 components in acetonitrile/ acetone 95/5: Phenanthrene [CAS:85-01-8] 600µg/ml; Anthracene [CAS:120-12-7] 40µg/ml; Fluoranthene [CAS:206-44-0] 160µg/ml; Pyrene [CAS:129-00-0] 160µg/ml; Benzo(a)anthracene [CAS:56-55-3] 20µg/ml; Chrysene [CAS:218-01-9] 80µg/ml; (95/5) Benzo(b)fluoranthene [CAS:205-99-2] 20µg/ml; Benzo(k)fluoranthene [CAS:207-08-9] 10µg/ml; Benzo(j)fluoranthene [CAS:205-82-3] 20µg/ml; Benzo(a)pyrene [CAS:50-32-8] 20µg/ml; Dibenzo(a,h)anthracene [CAS:53-70-3] 10µg/ml; Indeno(1,2,3-c,d)pyrene [CAS:193-39-5] 12µg/ml; Benzo(g,h,i)perylene [CAS:191-24-2] 20µg/ml
506878	1 ml	Glass ampoule	19 components 100 µg/ml each in acetonitrile: Naphthalene [CAS:91-20-3]; Acenaphthylene [CAS:208-96-8]; 1-Methylnaphthalene [CAS:90-12-0]; 2-Methylnaphthalene [CAS:91-57-6]; Acenaphthene [CAS:83-32-9]; Fluorene [CAS:86-73-7]; Phenanthrene [CAS:85-01-8]; Anthracene [CAS:120-12-7]; Fluoranthene [CAS:206-44-0]; Pyrene [CAS:129-00-0]; 2-Methyl-Fluoranthene [CAS:33543-31-6]; Benzo(a)anthracene [CAS:56-55-3;] Chrysene [CAS:218-01-9]; Benzo(b)fluoranthene [CAS:205-99-2]; Benzo(k)fluoranthene [CAS:207-08-9]; Benzo(a) pyrene [CAS:50-32-8]; Benzo(g,h,i)perylene [CAS:191-24-2]; Indeno(1,2,3-c,d)pyrene [CAS:193-39-5]; Dibenzo(a,h)anthracene [CAS:53-70-3]
506938	1 ml	Glass ampoule	23 components 1 µg/ml each in methanol: Anthracene [CAS:120-12-7]; Fluoranthene [CAS:206-44-0]; Naphthalene [CAS:91-20-3]; Benzo(a)pyrene [CAS:50-32-8]; Benzo(b)fluoranthene [CAS:205-99-2]; Benzo(g,h,i)perylene [CAS:191-24-2]; Benzo(k) fluoranthene [CAS:207-08-9]; Indeno(1,2,3-c,d)pyrene [CAS:193-39-5]; Biphenyl [CAS:92-52-4]; Acenaphthene [CAS:83-32-9]; Acenaphthylene [CAS:208-96-8]; Ben zo(a)anthracene [CAS:56-55-3]; Chrysene [CAS:218-01-9]; Dibenzo(a,h)anthracene [CAS:53-70-3]; Fluorene [CAS:86-73-7]; 2-Methylnaphthalene [CAS:91-57-6]; 2-Methyl-Fluoranthene [CAS:33543-31-6]; Phenanthrene [CAS:85-01-8]; Pyrene [CAS:129-00-0]; 1-Benzothiophene [CAS:95-15-8]; Dibenzothiophene [CAS:132-65-0]; Benzo(e)pyrene [CAS:192-97-2]; Perylene [CAS:198-55-0]
506979	1 ml	Glass ampoule	6 components in acetonitrile: Fluoranthene [CAS:206-44-0] 2mg/l; Benzo(b)fluoranthene [CAS:205-99-2] 2mg/l; Benzo(k) fluoranthene [C AS:207-08-9] 2mg/l; Benzo(a)pyrene [CAS:50-32-8] 2mg/l; Benzo(g,h,i)perylene [CAS:191-24-2] 10mg/l; Indeno(1,2,3-c,d)pyrene [CAS:193-39-5] 10mg/l
506980	1 ml	Glass ampoule	15 components in acetonitrile: Acenaphthene [CAS:83-32-9] 5mg/l; Fluorene [CAS:86-73-7] 5mg/l; Fluoranthene [CAS:206-44-0] 5mg/l; Benzo(a)anthracene [CAS:56-55-3] 5mg/l; Chrysene [CAS:218-01-9] 5mg/l; Benzo(b)fluoranthene [CAS:205-99-2] 5mg/l; Benzo(a)pyrene [CAS:50-32-8] 5mg/l; Dibenzo(a,h) anthracene [CAS:53-70-3] 5mg/l; Benzo(g,h,i)perylene [CAS:191-24-2] 10mg/l; Indeno(1,2,3-c,d)pyrene [CAS:193-39-5]10mg/l; 2-Methylnaphthalene [CAS:91-57-6] 10mg/l; 2-Methyl-Fluoranthene [CAS:33543-31-6] 10mg/l; Anthracene [CAS:120-12-7] 2mg/l; Benzo(k)fluoranthene [CAS:207-08-9] 2mg/l; Pyrene [CAS:129-00-0] 20mg/l

Code	Size	Packaging	Notes
507094	1 ml	Glass ampoule	19 components in Methanol: 2-Methylnaphthalene [CAS:91-57-6] 40mg/l; Anthracene [CAS:120-12-7] 20mg/l; Fluoranthene [CAS:206-44-0] 20mg/l; 2-Methyl-Fluoranthene [CAS:33543316] 20mg/l Benzo(a)anthracene [CAS:56-55-3] 20mg/l; Benzo(b) fluoranthene [CAS:205-99-2] 20mg/l Benzo(k)fluoranthene [CAS:207-08-9] 20mg/l; Benzo(a)pyrene [CAS:50-32-8] 20mg/l; Dibenzo(a,h)anthracene [CAS:53-70-3] 20mg/l; Benzo(g,h,i) perylene [CAS:191-24-2] 20mg/l; Indeno(1,2,3-c,d)pyrene [CAS:193-39-5] 20mg/l; Benzo(b)chrysene [CAS:214-17-5] 2mg/l; Naphthalene [CAS:91-20-3] 40mg/l; Acenaphthene [CAS:83-32-9] 40mg/l; Fluorene [CAS:86-73-7] 20mg/l; Phenanthrene [CAS:85-01-8] 20mg/l; Pyrene [CAS:129-00-0] 20mg/l; Chrysene [CAS:218-01-9] 20mg/l; Acenaphthylene [CAS:208-96-8] 400mg/l
507859	1.5 ml	Glass ampoule	15 components 1 µg/ml each in dichloromethane: Naphthalene [CAS:91-20-3]; Acenaphthene [CAS:83-32-9]; Fluorene [CAS:86-73-7]; Phenanthrene [CAS:85-01-8]; Anthracene [CAS:120-12-7]; Fluoranthene [CAS:206-44-0]; Pyrene [CAS:129-00-0]; Benzo(a)anthracene [CAS:56-55-3]; Chrysene [CAS:218-01-9] Benzo(b)fluoranthene [CAS:205-99-2]; Benzo(k)fluoranthene [CAS:207-08-9]; Benzo(a)pyrene [CAS:50-32-8]; Dibenzo(a,h) anthracene [CAS:53-70-3]; Benzo(g,h,i)perylene [CAS:191-24-2]; Indeno(1,2,3-c,d)pyrene [CAS:193-39-5]
507899	1.5 ml	Glass ampoule	15 components 1 µg/ml each in dichloromethane: Naphthalene [CAS:91-20-3]; Acenaphthene [CAS:83-32-9]; Fluorene [CAS:86-73-7]; Phenanthrene [CAS:85-01-8]; Anthracene [CAS:120-12-7]; Fluoranthene [CAS:206-44-0]; Pyrene [CAS:129-00-0]; Benzo(a)anthracene [CAS:56-55-3]; Chrysene [CAS:218-01-9] Benzo(b)fluoranthene [CAS:205-99-2]; Benzo(k)fluoranthene [CAS:207-08-9]; Benzo(a)pyrene [CAS:50-32-8]; Dibenzo(a,h) anthracene [CAS:53-70-3]; Benzo(g,h,i)perylene [CAS:191-24-2]; Indeno(1,2,3-c,d)pyrene [CAS:193-39-5]
506821	10 ml	Glass ampoule	19 components 10mg/l each in methanol: Naphthalene [CAS:91-20-3]; Acenaphthylene [CAS:208-96-8]; Acenaphthene [CAS:83-32-9]; Fluorene [CAS:86-73-7]; Phenanthrene [CAS:85-01-8]; Anthracene [CAS:120-12-7]; Fluoranthene [CAS:206-44-0]; Pyrene [CAS:129-00-0]; Benzo(a)anthracene [CAS:56-55-3]; Chrysene [CAS:218-01-9]; Benzo(b)fluoranthene [CAS:205-99-2]; Benzo(k)fluoranthene [CAS:207-08-9]; Benzo(a)pyrene [CAS:50-32-8]; Dibenzo(a,h)anthracene [CAS:53-70-3]; Benzo(g,h,i) perylene [CAS:191-24-2]; Indeno(1,2,3-c,d)pyrene [CAS:193-39-5]; Benzo(e)pyrene [CAS:192-97-2]; 2-Methylnaphthalene [CAS:91-57-6]; 2-Methyl-Fluoranthene [C AS:33543-31-6]

Custom formulations of organic substances are available. Contact us for more details.

Organic Standard: PCB multielement mixture

• Standard organico: PCB mix • Standard organique: mélange multiéléments PCB • Patrones orgánicos : mezcla multi-PCB
• Organisch standard: PCB Multielement-Mischung

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H336-H304-H410
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Organic Standard: PCB multielement mixture > RS - For environmental analysis

RS

Code	Size	Packaging	Notes
507103	5 x 1 ml	Glass ampoule	7 components 10 µg/ml each in isooctane: PCB 28 [CAS:7012-37-5]; PCB 52 [CAS:35693-99-3]; PCB 101 [CAS:37680-73-2]; PCB 118 [CAS:31508-00-6]; PCB 138 [CAS:35065-28-2]; PCB 153 [CAS:35065-27-1]; PCB 180 [CAS:35065-29-3]
507609	1 ml	Glass ampoule	2 components 100 µg/ml each in isooctane: PCB 30 [CAS:35693-92-6]; PCB 155 [CAS:33979-03-2]
507679	1 ml	Glass ampoule	8 components 100 µg/ml each in hexane: PCB 28 [CAS:7012-37-5]; PCB 52 [CAS:35693-99-3]; PCB 101 [CAS:37680-73-2]; PCB 118 [CAS:31508-00-6]; PCB 138 [CAS:35065-28-2]; PCB 153 [CAS:35065-27-1]; PCB 180 [CAS:35065-29-3]; PCB 194 [CAS:35694-08-7]
507062	1.2 ml	Glass ampoule	14 components 10 mg/l each in hexane: PCB 28 [CAS:7012-37-5]; PCB 31 [CAS:16606-02-3]; PCB 52 [CAS:35693-99-3]; PCB 101 [CAS:37680-73-2]; PCB 105 [CAS:32598-14-4]; PCB 118 [CAS:31508-00-6]; PCB 132 [CAS:38380-05-1]; PCB 138 [CAS:35065-28-2]; PCB 149 [CAS:38380-04-0]; PCB 153 [CAS:35065-27-1]; PCB 160 [CAS:41411-62-5]; PCB 163 [CAS:74472-44-9]; PCB 180 [CAS:35065-29-3]; PCB 193 [CAS:69782-91-8]
506732	5 ml	Glass ampoule	19 components in ethyle acetate: PCB 18 [CAS:37680-65-2] 0.34mg/l; PCB 28 [CAS:7012-37-5] 0.6mg/l; PCB 52 [CAS:35693-99-3] 0.9mg/l; PCB 77 [CAS:32598-13-3] 2.2mg/l; PCB 81 [CAS:70362-50-4] 3.3mg/l; PCB 101 [CAS:37680-73-2] 2.28mg/l; PCB 105 [CAS:32598-14-4] 2.6mg/l; PCB 114 [CAS:74472-37-0] 9.6mg/l; PCB 118 [CAS:31508-00-6] 2.6mg/l; PCB 123 [CAS:65510-44-3] 2.7mg/l; PCB 126 [CAS:57465-28-8] 3mg/l; PCB 138 [CAS:35065-28-2] 6mg/l; PCB 153 [CAS:35065-27-1] 5mg/l; PCB 156 [CAS:38380-08-4] 5mg/l; PCB 157 [CAS:69782-90-7] 7mg/l; PCB 167 [CAS:52663-72-6] 8mg/l; PCB 169 [CAS:32774-16-6] 10mg/l; PCB 180 [CAS:35065-29-3] 10mg/l; PCB 189 [CAS:39635-31-9] 7mg/l
507889	5 ml	Glass ampoule	14 components 1 mg/l each in hexane: PCB 28 [CAS:7012-37-5]; PCB 52 [CAS:35693-99-3]; PCB 101 [CAS:37680-73-2]; PCB 153 [CAS:35065-27-1]; PCB 138 [CAS:35065-28-2]; PCB 180 [CAS:35065-29-3] PCB 194 [CAS:35694-08-7]; 1,2,4-Trichlorobenzene [CAS:120-82-1]; Hexachloro-1, 3-butadiene [CAS:87-68-3]; Hexachlorobenzene [CAS:118-74-1]; Alpha-HCH [CAS:319-84-6]; Beta-HCH [CAS:319-85-7]; Gamma-HCH (Lindane) [CAS:58-89-9]; Delta-HCH [CAS:319-86-8]
507115	10 ml	Glass ampoule	14 components 10 µg/ml each in isooctane according to EN 61619: PCB 18 [CAS:37680-65-2]; PCB 28 [CAS:7012-37-5]; PCB 31 [CAS:16606-02-3]; PCB 44 [CAS:41464-39-5]; PCB 52 [CAS:35693-99-3]; PCB 101 [CAS:37680-73-2]; PCB 118 [CAS:31508-00-6]; PCB 138 [CAS:35065-28-2]; PCB 149 [CAS:38380-04-0]; PCB 153 [CAS:35065-27-1]; PCB 170 [CAS:35065-30-6]; PCB 180 [CAS:35065-29-3]; PCB 194 [CAS:35694-08-7]; PCB 209 [CAS:2051-24-3]
509144	10 ml	Glass ampoule	PCB 138 [35065-28-2] 50µg/ml in isooctane
509145	10 ml	Glass ampoule	PCB 153 [35065-27-1] 50µg/ml in isooctane
509146	10 ml	Glass ampoule	PCB 18 [37680-65-2] 50µg/ml in isooctane
509147	10 ml	Glass ampoule	PCB 28 [7012-37-5] 50µg/ml in isooctane
509148	10 ml	Glass ampoule	PCB 52 [35693-99-3] 50µg/ml in isooctane

Custom formulations of organic substances are available. Contact us for more details.

Organic standard: Pesticide mixture

• Standard organico: Pesticidi mix • Standard organique: mélange de pesticides • Patrones orgánicos : mezcla de pesticidas
• Organisch standard: Pestizidmischung

Organic standard: Pesticide mixture > RS - For environmental analysis

RS

Code	Size	Packaging	Notes
506897	1 ml	Glass ampoule	45 components 10mg/ml each in cyclohexane/acetone: Acetochlor [CAS:34256-82-1]10mg/l; Aclonifen [CAS:74070-46-5] 10mg/l; Alachlor [CAS:15972-60-8] 10mg/l; Bifenthrin [CAS:82657-04-3] 10mg/l;Cadusafos [CAS:95465-99-9] 10mg/l; Captan [CAS:133-06-2] 10mg/l; Carbofuran [CAS:1563-66-2] 10mg/l; Chlorfenvinphos [CAS:470-90-6]10mg/l; Chlormephos [CAS:24934-91-6] 10mg/l; Chlorothalonil [CAS:1897-45-6] 10mg/l Chlorpyrifos [CAS:2921-88-2] 10mg/l; Chlorpyrifos methyl [CAS:5598-13-0] 10mg/l; lambda-Cyhalothrin [CAS:91465-08-6] 10mg/l; Cypermethrin [CAS:52315-07-8]10mg/l; Delta-HCH [CAS:319-86-8] 10mg/l; Diazinon [CAS:333-41-5] 10mg/l; Dichlobenil [CAS:1194-65-6] 10mg/l; Dinoterb [CAS:1420-07-1] 10mg/l; Endosulfan-total (sulfate) [CAS:1031-07-8] 10mg/l; Fipronil [CAS:120068-37-3] 10mg/l; Folpet [CAS:133-07-3] 10mg/l; Heptachlor-endo-epoxide [CAS:28044-83-9] 10mg/l; Hexachloro-1,3-butadiene [CAS:87-68-3]10mg/l; Iprodione [CAS:36734-19-7] 10mg/l; Isofenphos [CAS:25311-71-1]10mg/l; Malathion [CAS:121-75-5] 10mg/l; Metazachlor [CAS:67129-08-2] 10mg/l; Oxadiazon [CAS:19666-30-9] 10mg/l; Oxyfluorfen [CAS:42874-03-3] 10mg/l; Parathion (Parathion-ethyl) [CAS:56-38-2] 10mg/l; Parathion-methyl [CAS:298-00-0] 10mg/l; Pendimethalin [CAS:40487-42-1] 10mg/l; Pentachlorobenzene [CAS:608-93-5] 10mg/l; Procymidone [CAS:32809-16-8] 10mg/l; Propachlor [CAS:1918-16-7] 10mg/l; Tebutam [CAS:35256-85-0] 10mg/l; Tefluthrin [CAS:79538-32-2] 10mg/l; Terbufos [CAS:13071-79-9] 10mg/l; Tolylfluanid [CAS:731-27-1] 10mg/l; Triazophos [CAS:24017-47-8] 10mg/l; Trifluralin [CAS:1582-09-8] 10mg/l; Vinclozolin [CAS:50471-44-8] 10mg/l; Piperonyl butoxide [CAS:51-03-6] 10mg/l; Metolachlor [CAS:51218-45-2] 5mg/l; S-Metolachlor [CAS:87392-12-9] 5mg/l

Product specifications are subject to changes. Please visit our website for updates.

Code	Size	Packaging	Notes
506905	1 ml	Glass ampoule	79 components in acetone: Bifenthrin [CAS:82657-04-3] 120µg/ml; lambda-Cyhalothrin [CAS:91465-08-6] 100µg/ml; Cypermethrin [CAS:5231 5-07-8] 130µg/ml; Deltamethrin [CAS:52918-63-5] 130µg/ml; Fenvalerate [CAS:51630-58-1] 105µg/ml; Permethrin [CAS:52645-53-1] 100µg/ml; tau-Fluvalinate [CAS:102851-06-9] 100µg/ml; Tetramethrin [CAS:7696-12-0] 100µg/ml; Aldrin [CAS:309-00-2] 20µg/ml; cis-Chlordane [CAS:5103-71-9] 20µg/ml; trans-Chlordane [CAS:5103-74-2] 20µg/ml; 2,4'-DDD [CAS:53-19-0] 20µg/ml; 4,4'-DDD (TDE) [CAS:72-54-8] 20µg/ml; 2,4'-DDE [CAS:3424-82-6] 20µg/ml; 4,4'-DDE [CAS:72-55-9] 20µg/ml; 2,4'-DDT [CAS:789-02-6] 20µg/ml; 4,4'-DDT [CAS:50-29-3] 20µg/ml; Dieldrin [CAS:60-57-1] 20µg/ml; Endosulfan-alpha [CAS:959-98-8] 20µg/ml; Endosulfan-beta [CAS:33213-65-9] 20µg/ml; Endosulfan-total (sulfate) [CAS:1031-07-8] 20µg/ml; Endrin [CAS:72-20-8] 20µg/ml; Endrin aldehyde [CAS:7421-93-4] 20µg/ml; Alp ha-HCH [CAS:319-84-6] 20µg/ml; Beta-HCH [CAS:319-85-7] 20µg/ml; Delta-HCH [CAS:319-86-8] 20µg/ml; Gamma-HCH (Lindane) [CAS:58-89-9] 20µg/ml; Heptachlor [CAS:76-44-8] 20µg/ml; Heptachlor-endo-epoxide [CAS:28044-83-9] 20µg/ml; Heptachlor-exo-epoxide [CAS:1024-57-3] 20µg/ml;Hexachlorobenzene [CAS:118-74-1] 20µg/ml; PCB 209 [CAS:2051-24-3] 20µg/ml; PCB 29 [CAS:15862-07-4] 20µg/ml; Vinclozolin [CAS:50471-44-8] 20µg/ml; Alachlor [CAS:15972-60-8] 100µg/ml; Bromopropylate [CAS:18181-80-1] 50µg/ml; Chlorothalonil [CAS:1897-45-6] 25µg/ml; Dicofol [CAS:115-32-2] 75µg/ml; Iprodione [CAS:36734-19-7] 200µg/ml; Nitrofen [CAS:1836-75-5] 20µg/ml; oxy-Chlordane [CAS:27304-13-8] 20µg/ml; Phosalone [CAS:2310-17-0] 20µg/ml; Procymidone [CAS:32809-16-8] 150µg/ml; Tetradifon [CAS:116-29-0] 20µg/ml; Bromophos-ethyl [CAS:4824-78-6] 100µg/ml; Bromophos-methyl [CAS:2104-96-3] 100µg/ml; Chlorfenvinphos [CAS:470-90-6] 100µg/ml; Chlorpyrifos (Chlorpyrifos-ethyl) [CAS:2921-88-2] 100µg/ml; Chlorpyrifos methyl [CAS:5598-13-0] 100µg/ml; Diazinon [CAS:333-41-5] 100µg/ml; Dichlorvos [CAS:62-73-7] 100µg/ml; Dimethoate [CAS:60-51-5] 100µg/ml; Disulfoton [CAS:298-04-4] 50µg/ml; Fenchlorphos [CAS:299-84-3] 100µg/ml; Fenthion [CAS:55-38-9] 100µg/ml; Malathion [CAS:121-75-5] 100µg/ml; Parathion (Parathion-ethyl) [CAS:56-38-2] 100µg/ml; Parathion-methyl [CAS:298-00-0] 100µg/ml; Pirimiphos-methyl [CAS:29232-93-7] 100µg/ml; Terbufos [CAS:13071-79-9] 100µg/ml; Acephate [CAS:30560-19-1] 100µg/ml; Azinphos-ethyl [CAS:2642-71-9] 400µg/ml; Azinphos-methyl [CAS:86-50-0] 400µg/ml; Demeton-S-methyl [CAS:919-86-8] 100µg/ml; Ethion [CAS:563-12-2] 20µg/ml; Fenamiphos [CAS:22224-92-6] 50µg/ml; Fenitrothion [CAS:122-14-5] 50µg/ml; Fonofos [CAS:944-22-9] 40µg/ml; Metalaxyl [CAS:57837-19-1] 600µg/ml; Methamidophos [CAS:10265-92-6] 100µg/ml; Methidathion [CAS:950-37-8] 100µg/ml; Mevinphos [CAS:7786-34-7] 100µg/ml; Monocrotophos [CAS:6923-22-4] 100µg/ml; Oxa dixyl [CAS:77732-09-3] 400µg/ml; Phorate [CAS:298-02-2] 50µg/ml; Phosphamidon [CAS:13171-21-6] 100µg/ml; Pirimiphos-ethyl [CAS:23505-41-1] 50µg/ml; Triazophos [CAS:24017-47-8] 100µg/ml; Tefluthrin [CAS:79538-32-2] 10µg/ml
506948	1 ml	Glass ampoule	29 component 20µg/ml each in toluene/acetone: Aldrin [CAS:309-00-2]; Alpha-HCH [CAS:319-84-6]; Beta-HCH [CAS:319-85-7]; Delta-HCH [CAS:319-86-8]; cis-Chlordane (alpha-Chlordane) [CAS:5103-71-9]; Dieldrin [CAS:60-57-1]; Endosulfan-alpha [CAS:959-98-8]; Endos ulfan-beta [CAS:33213-65-9]; Endosulfan-total (sulfate) [CAS:1031-07-8]; Endrin [CAS:72-20-8]; Endrin aldehyde [CAS:7421-93-4]; Endrin ketone [CAS:53494-70-5]; Gamma-HCH (Lindane) [CAS:58-89-9]; trans-Chlordane (Gamma-Chlordane) [CAS:5103-74-2]; Heptachlor [CAS:76-44-8]; Heptachlor-exo-epoxide (cis-Heptachlorepoxide (cis-, exo-,) [CAS:1024-57-3]; Methoxychlor (DMTD) [CAS:72-43-5]; 4,4'-DDD (TDE) [CAS:72-54-8]; 4,4'-DDE [CAS:72-55-9]; 4,4'-DDT [CAS:50-29-3]; Dicofol [CAS:115-32-2]; Nitrofen [CAS:1836-75-5]; Isodrin [CAS:465-73-6]; Alachlor [CAS:15972-60-8]; Hexachlorobenzene (HCB) [CAS:118-74-1]; 2,4'-DDE [CAS:3424-82-6]; 2,4'-DDD [CAS:53-19-0]; 2,4'-DDT [CAS:789-02-6]; oxy-Chlordane [CAS:27304-13-8]; trans-Nonachlor [CAS:39765-80-5]

Code	Size	Packaging	Notes
506950	1 ml	Glass ampoule	17 components 20µg/ml each in toluene/acetone: Cyfluthrin [CAS:68359-37-5]; Cypermethrin [CAS:52315-07-8]; Fenvalerate [CAS:51630-58-1]; Permethrin [CAS:52645-53-1]; Phenothrin [CAS:26002-80-2]; Tetramethrin [CAS:7696-12-0]; lambda-Cyhalothrin [CAS:91465-08-6]; Piperonyl butoxide [CAS:51-03-6]; Bifenthrin [CAS:82657-04-3]; Chlorothalonil [CAS:1897-45-6]; Quintozene [CAS:82-68-8]; Tecnazene [CAS:117-18-0]; Chlorobenzilate [CAS:510-15-6]; Vinclozolin [CAS:50471-44-8]; Chlordecone hydrate [CAS:143-50-0]; Captan [CAS:133-06-2]
506803	10 ml	Bottle	12 components 10µg/ml each in acetonitrile: Azoxystrobin [CAS:131860-33-8]; Boscalid [CAS:188425-85-6]; Carbendazim [CAS:10605-21-7]; Chlorpyrifos [CAS:2921-88-2]; Cyprodinil [CAS:121552-61-2]; Linuron [CAS:330-55-2]; Metalaxyl [CAS:57837-19-1]; Methomyl [CAS:16752-77-5]; Myclobutanil [CAS:88671-89-0]; Pyrimethanil [CAS:53112-28-0]; Pirimicarb [CAS:23103-98-2]; Thiabendazole [CAS:148-79-8]

Custom formulations of organic substances are available. Contact us for more details.

Organic standard: Mixture for hydrocarbon analysis

• Soluzione standard per analisi degli idrocarburi • Mélange standard pour analyse des hydrocarbures • Mezcla estándar para análisis de hidrocarburos • Standardgemisch für die Kohlenwasserstoffanalyse

Organic standard: Mixture for hydrocarbon analysis > RS - For enviromental analysis

RS

Code	Size	Packaging	Notes
506736	1 ml	Glass ampoule	5 components 5000µg/ml each in methanol: Benzene; Toluene; o-Xylene; m-Xylene; p-Xylene
506742	1 ml	Glass ampoule	14 components in methanol: 1,1-Dichloroethene 1000µg/ml; Dichloromethane 5000µg/ml; trans-1,2-Dichloroethene 5000µg/ml; 1,1-Dichloroethane 5000µg/ml; cis-1,2-Dichloroethene 5000µg/ml; 1,2-Dichloroethane 5000µg/ml; Chloroform 500µg/ml; 1,1,1-Trichloroethane 500µg/ml; Trichloroethene 500µg/ml; Bromodichloromethane 500µg/ml; Dibromochloromethane 500µg/ml;Tribromomethane 500µg/ml; Tetrachloromethane 100µg/ml; Tetrachloroethene 100µg/ml; Hydrocarbons Mixture Benzene; 1-Bromo-2-chloroethane; Chlorobenzene; 1,1-Dichloroethane; 1,2-dichloroethane
507189	1 ml	Glass ampoule	6 components 1000 µg/ml each in methanol: Benzene; Toluene; o-Xylene; m-Xylene; p-Xylene; Ethylbenzene
507190	1 ml	Glass ampoule	6 components 1000 µg/ml each in methanol: Trichloroethene; Tetrachloroethene; 1,2-Dichloroethane; Tetrachloromethane; 1,1,1-Trich loroethane; 1,1,2-Trichloroethane
507191	1 ml	Glass ampoule	4 components 1000 µg/ml each in methanol: Tribromomethane; Chloroform; Bromodichloromethane; Dibromochloromethane
507474	2 ml	Glass ampoule	6 components 1000 µg/ml each in carbon disulfide: Benzene; Ethylbenzene; Toluene; m-Xylene; o-Xylene; p-Xylene
506614	100 ml	Bottle	22 components 2500mg/Kg each 1-Bromo-2-chloroethane; Chlorobenzene; 1,1-Dichloroethane; 1,1-Dichloroethene; cis-1,2-Dichloroethene; trans-1,2-Dichloroethene; Dichloromethane; Pentachloroethane; 1,1,2,2-Tetrachloroethane; Tetrachloroethene; Tetrachloromethane; 1,1,2-Trichloroethane; Trichloroethene: Chloroprene; Chloromethane; Vinylchloride; 1,3-Butadiene; Chloroethane; 1,2-Dichlorobutane; Ethylene; Chloroform; Matrix: Benzene

Custom formulations of organic substances are available. Contact us for more details.

Organic standard: Mixture for hydrocarbon analysis > RS - For enviromental analysis according to NF EN ISO 9377-2

RS

Code	Size	Packaging	Notes
506002	1 ml	Glass ampoule	Standard quality control of two mineral oils in acetone
506010	1 ml	Glass ampoule	Mixture of mineral oil without additive 2 to 5 mg / ml each in hexane
506020	1 ml	Glass ampoule	Standard mixture of n-alkanes (C10 to C40 in pairs) of 50 mcg / ml each in hexane
506012	5 ml	Glass ampoule	2 component: 5g/l each of Mineral Oil [CAS:8042-47-5]; Diesel Oil in n-Hexane
506040	5 ml	Glass ampoule	N-tetracontane mixture (20 mg / l) and n-decane (20 mg / l) in hexane
506011	10 ml	Glass ampoule	Mixture of mineral oil without additive 2 to 1 mg / ml each in hexane
506013	10 ml	Glass ampoule	2 components: 5g/l each of Mineral Oil [CAS:8042-47-5]; Diesel Oil in n-Hexane
506021	10 ml	Glass ampoule	Standard mixture of n-alkanes (C10 to C40 in pairs) of 50 µg / ml each in hexane
506030	10 ml	Glass ampoule	Mother solution stearyl stearate 2 g / l in hexane

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Orthophosphoric acid 99%

• Acido fosforico 99% • Acide phosphorique 99% • Acido orto-fosforico 99% • ortho-Phosphorsäure 99%

Synonym: *Phosphoric acid*

H_3PO_4
Molecular Weight: 98
CAS: 7664-38-2
EEC-N: 231-633-2

Classification transport
ONU: 3453
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Orthophosphoric acid 99% > RPE - For analysis

RPE

Description ..... White deliquescent crystals	Nitrate .....≤2 ppm	Co .....≤5 ppm	Na .....≤50 ppm
Identification ..... Positive	Reducing subst as H3PO3 .....≤50 ppm	Cu .....≤5 ppm	Ni .....≤5 ppm
Chloride .....≤2 ppm	Silicate .....≤500 ppm	Fe .....≤10 ppm	Pb .....≤5 ppm
Fluoride .....≤5 ppm	Sulphate .....≤10 ppm	K .....≤20 ppm	Zn .....≤10 ppm
Ca, Mg and ppt by NH4OH .....≤50 ppm	As .....≤0.5 ppm	Mg .....≤10 ppm	Assay (acidimetric) ..... ≥99 %
Heavy metals (Pb) .....≤10 ppm	Cd .....≤5 ppm	Mn .....≤0.5 ppm	

Code	Size	Packaging	Notes
405967	1 kg	Plastic bottle	
405961	10 kg	Plastic tank	

Orthophosphoric acid 85%

• Acido fosforico 85% • Acide phosphorique 85% • Acido orto-fosforico 85% • ortho-Phosphorsäure 85%

Synonym: *Phosphoric acid*

H_3PO_4
Molecular Weight: 98
CAS: 7664-38-2

Classification transport
ONU: 1805
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Orthophosphoric acid 85% > RS - VLSI - For electronic use

RS

Code	Size	Packaging	Notes
527592	1 l	Plastic bottle	
527591	2.5 l	Plastic bottle	

For specifications, contact our customer service for a certificate of analysis

Orthophosphoric acid 85% > RS - RSE - For electronic use

RS

Description	Clear liquid	Heavy metals (Pb)	≤5 ppm	Ca	≤20 ppm	Mg	≤5 ppm
Colour (APHA)	≤10	Nitrate	≤3 ppm	Cd	≤5 ppm	Mn	≤0.5 ppm
Identification	Positive	Subst. reducing KMnO4	≤10 ppm	Co	≤1 ppm	Na	≤30 ppm
Density at 20° C	1.689 ÷ 1.701	Sulphate	≤5 ppm	Cu	≤2 ppm	Ni	≤3 ppm
Assay (acidimetric)	85.0 ÷ 87.0 %	Volatile acid	≤3 ppm	Fe	≤5 ppm	Pb	≤1 ppm
Ammonium	≤5 ppm	Al	≤0.5 ppm	Ga	≤0.1 ppm	Sr	≤5 ppm
Chloride	≤1 ppm	As + Sb (as As)	≤0.5 ppm	K	≤5 ppm	Zn	≤10 ppm
Fluoride	≤5 ppm	Bi	≤1 ppm	Li	≤1 ppm		

Code	Size	Packaging	Notes
406022	1 l	Plastic bottle	
406021	2.5 l	Plastic bottle	

Orthophosphoric acid 85% > RPE - For analysis - ISO - ACS - Reag.Ph.Eur. - Reag.USP

RPE

Description	Clear liquid	Nitrate	≤3 ppm	K	≤50 ppm	Appearance of solution	Conform Ph. Eur.
Colour (APHA)	≤10	Reducing substances	≤25 ppm	Mg	≤10 ppm	Substances precipitated with ammonia	Conform Ph. Eur.
Identification	Positive	Sulphate	≤5 ppm	Mn	≤0.5 ppm	Phosphorous and hypophosphorous acid	Conform Ph. Eur.
Density at 20° C	1.689 ÷ 1.701	As	≤0.5 ppm	Na	≤250 ppm	Alkali phosphates	Conform USP-NF
Volatile acid	≤10 ppm	Ca	≤20 ppm	Ni	≤5 ppm		
Water-insoluble matter	≤10 ppm	Cd	≤5 ppm	Pb	≤1 ppm		
Chloride	≤1 ppm	Co	≤1 ppm	Sb	≤4 ppm		
Fluoride	≤5 ppm	Cu	≤2 ppm	Zn	≤10 ppm		
Heavy metals (Pb)	≤5 ppm	Fe	≤5 ppm	Assay (acidimetric)	85 ÷ 87 %		

Code	Size	Packaging	Notes
406002	1 l	Plastic bottle	
406005	2.5 l	Plastic bottle	
406003	40 kg	Plastic drum	

Orthophosphoric acid 85% > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP-DAB

ERBApharm

Description	Clear colourless liquid	Subst. ppt by NH4OH	Conform Ph.Eur.	Chloride	≤50 ppm	Fe	≤50 ppm
Identification	Positive	Nitrate	Conform USP-NF	Sulphate	≤100 ppm	Assay (acidimetric)	85 ÷ 88 %
Appearance of solution	Conform Ph.Eur.	Alkali phosphates	Conform USP-NF	Heavy metals (Pb)	≤10 ppm	Origin (BSE/TSE)	Synthesis
Hypophos. phosphor acid	Conform Ph.Eur.	Density at 20° C	1.689 ÷ 1.701	As	≤2 ppm	Residual solvents (Current ICH)	Conform

Code	Size	Packaging	Notes
304061	1 l	Plastic bottle	
304062	2.5 l	Plastic bottle	
304063	40 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Orthophosphoric acid 85% > RE - Pure

RE

Description	Clear colourless liquid	Density at 20°C	1.69 ÷ 1.71	Iron (Fe)	≤ 20 ppm
Assay	85.0 ÷ 86.0 %	Residue on ignition	≤ 20 ppm	Chloride (Cl)	≤ 5 ppm

Code	Size	Packaging	Notes
528535	5 l	Tank	

Product specifications are subject to changes. Please visit our website for updates.

Orthophosphoric acid 75%

• Acido fosforico 75% • Acide phosphorique 75% • Acido orto-fosfórico 75% • ortho-Phosphorsäure 75%

Synonym:
Phosphoric acid

H_3PO_4
Molecular Weight: 98
CAS: 7664-38-2

Classification transport
ONU: 1805
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Orthophosphoric acid 75% > RE - Pure

RE

Description ..... Clear liquid
Identification ..... Positive
Density at 20° C ..... 1.568 - 1.589
Chloride ..... ≤5 ppm
Heavy metals (Pb) ..... ≤10 ppm
Sulphate ..... ≤120 ppm
As ..... ≤1 ppm
Fe ..... ≤10 ppm
Assay (acidimetric) ..... 74.0 - 76.0 %

Code	Size	Packaging	Notes
304051	1 l	Plastic bottle	
304054	2.5 l	Plastic bottle	
304052	85 kg	Plastic tank	

Orthophosphoric acid 10%

• Acido fosforico 10% • Acide phosphorique 10% • Acido orto-fosfórico 10% • Ortho-phosphorsäure 10%

Synonym:
Phosphoric acid

H_3PO_4
Molecular Weight: 98
CAS: 7664-38-2
EEC-N: 231-633-2

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Orthophosphoric acid 10% > RS - For analysis

RS

Density d20/4 ..... 1.048 - 1.059
Assay ..... 9 - 11 %

Code	Size	Packaging	Notes
PS0084/22	5 l	Plastic tank	
PS0084/42	20 l	Plastic tank	

Orthophosphoric acid-d3 85% in D2O

• Acido fosforico-d3 85% in D20 • Acide phosphorique-d3 85% dans D20
• Acido orto-fosforico-d3 85% in D20 • ortho-Phosphorsäure-d3 85% in D20

Synonym:
• *Phosphoric acid-d3 solution*
• *Trideuterophosphoric acid*

D_3O_4P
Molecular Weight: 101,02
CAS: 14335-33-2

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group II

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Orthophosphoric acid-d3 85% in D20 > RS - For NMR - min 99%

RS

Code	Size	Packaging	Notes
P5055	25 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Osmium standard solution

• Osmio standard soluzione • Osmium solution standard • Osmio, solución patrón • Osmium-Standardlösung

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Osmium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505758	100 ml	Plastic bottle	conc. 100 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Osmolality Standards

• Osmolalità Standards • Etalons d'osmolalité • Patrones de osmolalidad • Osmolalitätsstandards

Osmolality Standards > RS - For calibration

RS

Code	Size	Packaging	Notes
540301	12 x 5 ml	Glass ampoule	100mOsm/Kg H2O
540302	12 x 5 ml	Glass ampoule	1500mOsm/Kg H2O
540303	12 x 5 ml	Glass ampoule	200mOsm/Kg H2O
540304	12 x 5 ml	Glass ampoule	2000mOsm/Kg H2O
540305	12 x 5 ml	Glass ampoule	290mOsm/Kg H2O
540306	12 x 5 ml	Glass ampoule	300mOsm/Kg H2O
540307	12 x 5 ml	Glass ampoule	400mOsm/Kg H2O
540308	12 x 5 ml	Glass ampoule	500mOsm/Kg H2O
540309	12 x 5 ml	Glass ampoule	850mOsm/Kg H2O
540310	12 x 5 ml	Glass ampoule	900mOsm/Kg H2O

Osmolality Standards Protein Based

• Osmolalità standard a base di proteine • Etalons d'osmolalité à base de protéines • Patrones de osmolalidad a base de proteínas
• Osmolalitätsstandards basierend auf Proteinen

Osmolality Standards Protein Based > RS - For calibration

RS

Code	Size	Packaging	Notes
540351	12 x 5 ml	Glass ampoule	240mOsm/Kg H2O
540352	12 x 5 ml	Glass ampoule	280mOsm/Kg H2O
540353	12 x 5 ml	Glass ampoule	320mOsm/Kg H2O

Osmolality Standards Urine Based

• Osmolalità Standards a base di urina • Etalons d'osmolalité à base d'urine • Patrones de osmolalidad a base de orina
• Osmolalitätsstandards basierend auf Urin

Osmolality Standards Urine Based > RS - For calibration

RS

Code	Size	Packaging	Notes
540354	12 x 5 ml	Glass ampoule	300mOsm/Kg H2O
540355	12 x 5 ml	Glass ampoule	800mOsm/Kg H2O

Oxalic acid dihydrate

• Acido ossalico diidrato • Acide oxalique dihydraté • Acido oxálico dihidrato • Oxalsäuredihydrat

Synonym:
Ethanedioic acid dihydrate

$(COOH)_2.2H_2O$
Molecular Weight: 126,07
CAS: 6153-56-6
EEC-N: 205-634-3



Danger
H302-H318
P264-P280i-P301+P312a-P305+P351+P338-P310a-P501a

Oxalic acid dihydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White crystals
Identification ..... Positive
Substances darkened by sulphuric acid ..... Conform
Total nitrogen ..... ≤ 10 ppm
Chloride ..... ≤ 20 ppm
Water-insoluble matter ..... ≤ 50 ppm
Heavy metals (Pb) ..... ≤ 5 ppm
Residue on ignition ..... ≤ 100 ppm
Ca ..... ≤ 10 ppm
Fe ..... ≤ 2 ppm
Assay (oxidimetric) ..... 99.5 ÷ 102.5 %
Sulphate ..... ≤ 50 ppm

Code	Size	Packaging	Notes
408736	500 g	Plastic bottle	
408737	1 kg	Plastic bottle	
408731	5 kg	Plastic tank	
408733	25 kg	Plastic bucket	

Product specifications are subject to changes. Please visit our website for updates.

Oxalic acid dihydrate > RE - Pure

RE

Description ..... White crystals Identification ..... Positive Fe ..... ≤ 15 ppm Assay (acidimetric) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
305757	1 kg	Plastic bottle	
305758	5 kg	Plastic tank	

Oxalic acid 0.5 mol/l (1N)

• Acido ossalico 0.5 mol/l (1N) • Acide oxalique 0.5 mol/l (1N) • Acido oxálico 0.5 mol/l (1N) • Oxalsäure 0.5 mol/l (1N)

$(COOH)_2$
CAS: 144-62-7

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Warning
H290-H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Oxalic acid 0.5 mol/l (1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.998 - 1.002 N

Code	Size	Packaging	Notes
408826	500 ml	Plastic bottle	

45,02 g of C2H2O4. Volumetric solution ready-to-use. For oxydometry. Stabilized with sulfuric acid

Oxalic acid 0.05 mol/l (0.1N)

• Acido ossalico 0.05 mol/l (0.1N) • Acide oxalique 0.05 mol/l (0.1N) • Acido oxálico 0.05 mol/l (0.1N) • Oxalsäure 0.05 mol/l (0.1N)

$(COOH)_2$
Molecular Weight: 90,03
CAS: 144-62-7

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Warning
H290-H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Oxalic acid 0.05 mol/l (0.1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.0998 - 0.1002 N

Code	Size	Packaging	Notes
408856	500 ml	Plastic bottle	

4,502 g of C2H2O4. Volumetric solution ready-to-use. For oxydometry. Stabilized with sulfuric acid

Oxalic acid 0.05 mol/l (0.1N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
408871		Plastic ampoule	Volume: 165 ml

Volumetric concentrated solution to prepare 1 L of solution 0,1 N

Oxalic acid 0.005 mol/l (0.01N)

• Acido ossalico 0.005 mol/l (0.01N) • Acide oxalique 0.005 mol/l (0.01N) • Acido oxálico 0.005 mol/l (0.01N) • Oxalsäure 0.005 mol/l (0.01N)

$(COOH)_2$
CAS: 144-62-7

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Warning
H290-H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Oxalic acid 0.005 mol/l (0.01N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
408901		Plastic ampoule	Volume: 55 ml

Volumetric concentrated solution to prepare 1 L of solution 0,01 N

Oxalic acid diammonium salt ▶ Ammonium oxalate monohydrate

2,2'-Oxydiethanol ▶ Diethylene glycol

Product specifications are subject to changes. Please visit our website for updates.

Palladium standard solution

• Palladio standard soluzione • Palladium solution standard • Paladio, solución patrón • Palladium-Standardlösung

Palladium standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2 — RS

Code	Size	Packaging	Notes
615003600	100 ml	Plastic bottle	A 500 ppm solution Ref Ph.Eur 5003600

Palladium standard solution > RS - Standard solution for ICP-MS — RS

Code	Size	Packaging	Notes
505772	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505775	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Palladium standard solution > RS - Standard solution for ICP — RS

Code	Size	Packaging	Notes
503811	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503813	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503815	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid
503817	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Palladium standard solution > RS - Standard solution for AAS — RS

Code	Size	Packaging	Notes
507751	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507512	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Palladium (II) chloride

• Palladio cloruro oso • Palladium (II) chlorure • Paladio (II) cloruro • Palladium (II) chlorid

$PdCl_2$
Molecular Weight: 177,31
CAS: 7647-10-1
EEC-N: 231-596-2

Classification transport
ONU: 3260
Transport Hazard class: 8
Packing group III



Danger
H290-H318-H317
P261-P280a-P305+P351+P338-P310a-P362+P364-P333+P313

Palladium (II) chloride > RPE - For analysis — RPE

Description ..... Red brown powder Identification ..... Positive Assay (gravimetric) ..... ≥59.5 % Pd

Code	Size	Packaging	Notes
467737	1 g	Glass bottle	
467731	10 g	Glass bottle	

Palladium nitrate 2 g/l solution

• Palladio nitrato 2 g/L soluzione • Palladium nitrate 2 g/l • Paladio nitrato solución 2 g/l • Palladiumnitrat 2 g/l Lösung

$Pd(NO_3)_2$
Molecular Weight: 230,43 (an,)
CAS: 10102-05-3

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group II



Danger
H290-H332-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Palladium nitrate 2 g/l solution > RS - Matrix modifiers for AAS-GTA

RS

Code	Size	Packaging	Notes
503198	50 ml	Plastic bottle	Matrix: 1% Nitric acid
503202	50 ml	Bottle	Matrix: 15% Nitric acid

Papanicolaou Haematoxylin solution according to Harris

• Papanicolaou Ematossilina soluzione secondo Harris • Papanicolaou Hématoxyline selon Harris • Papanicolaou Hematoxilina en solución según Harris • Papanicolaou Haematoxylin-Lösung nach Harris



Danger
H318
P280i-P305+P351+P338-P310a

Papanicolaou Haematoxylin solution according to Harris > RS - For histology

RS

Description .....Dark red liquid Identification ..... Positive Max absorbance wave-length....555 ÷ 565 nm Absorbance(lamba max) ..... ≥0.52 uA

Code	Size	Packaging	Notes
446462	500 ml	Plastic bottle	In Vitro Diagnostic Medical Device
446464	6 x 500 ml	Plastic bottle	In Vitro Diagnostic Medical Device
446461	1 l	Plastic bottle	In Vitro Diagnostic Medical Device
446465	6 x 1 l	Plastic bottle	In Vitro Diagnostic Medical Device
446463	2.5 l	Plastic bottle	In Vitro Diagnostic Medical Device
446466	4 x 2.5 l	Plastic bottle	In Vitro Diagnostic Medical Device

Papanicolaou solution EA 50

• Papanicolaou soluzione EA 50 • Papanicolaou solution EA 50 • Papanicolaou solución EA 50 • Papanicolaou-Lösung EA 50

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II



Danger
H225-H319-H370-H373
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

Papanicolaou solution EA 50 > RS - For histology

RS

Description .....Green clear liquid Maximum abs.Lambda max1(water,dil.1:200) .630 - 634 nm Maximum abs.Lambda max2(water,dil.1:200) .515 - 518 nm
Identification ..... Positive Absorbance (Lambda max1,water,dil.1:200) ..... 0.24 - 0.30 Absorbance (Lambda max2,water,dil.1:200) ..... 1.25 - 1.30

Code	Size	Packaging	Notes
467782	500 ml	Plastic bottle	In Vitro Diagnostic Medical Device
E467784	6 x 500 ml	Plastic bottle	In Vitro Diagnostic Medical Device
467781	1 l	Plastic bottle	In Vitro Diagnostic Medical Device
E467785	6 x 1 l	Plastic bottle	In Vitro Diagnostic Medical Device
467783	2.5 l	Plastic bottle	In Vitro Diagnostic Medical Device
E467786	4 x 2.5 l	Plastic bottle	In Vitro Diagnostic Medical Device

Dye for cytological diagnosis and oncology. Contains ethanol and methanol

Product specifications are subject to changes. Please visit our website for updates.

Papanicolaou solution OG 6

• Papanicolaou soluzione OG 6 • Papanicolaou solution OG 6 • Papanicolaou solución OG 6 • Papanicolaou Lösung OG 6

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II

Danger
H225-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

Papanicolaou solution OG 6 > RS - For histology

RS

Description ..... Orange clear liquid
Identification ..... Positive
Density at 20°C ..... ~ 0.83
Empirical test ..... Positive
Absorbance at 480 nm ..... 0.6 ÷ 0.8

Code	Size	Packaging	Notes
467792	500 ml	Plastic bottle	In Vitro Diagnostic Medical Device
E467794	6 x 500 ml	Plastic bottle	In Vitro Diagnostic Medical Device
467791	1 l	Plastic bottle	In Vitro Diagnostic Medical Device
E467795	6 x 1 l	Plastic bottle	In Vitro Diagnostic Medical Device
467793	2.5 l	Plastic bottle	In Vitro Diagnostic Medical Device
E467796	4 x 2.5 l	Plastic bottle	In Vitro Diagnostic Medical Device

Dye for cytological diagnosis and oncology. Contains ethanol and methanol

Paraffin 56-58°C - Erbaplast (without DMSO)

• Paraffina 56-58°C Erbaplast (senza DMSO) • Paraffine 56-58°C Erbaplast (sans DMSO) • Parafina 56-58°C Erbaplast (sin DMSO)
• Paraffin 56-58°C Erbaplast (ohne DMSO)

CAS: 92045-76-6
EEC-N: 295-458-3

Paraffin 56-58°C - Erbaplast (without DMSO) > RS - For histology - CE - IVD

RS

Description ..... White pellets
Melting point ..... 56 ÷ 58 °C

Code	Size	Packaging	Notes
467958	4 x 2 kg	Bag	In Vitro Diagnostic Medical Device

Paraffin oil

• Olio di vaselina • Huile de vaseline • Aceite de vaselina • Paraffinöl

Synonym:
Mineral oil

CAS: 8012-95-1
EEC-N: 232-455-8

Paraffin oil > RS - For optical spectroscopy

RS

Description ..... Colourless oily liquid
Identification (I.R.) ..... Conform
Density at 20°C ..... ~ 0.880

Code	Size	Packaging	Notes
466792	100 ml	Glass bottle	

Paraffin oil > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP

ERBApharm

Description ..... Colourless oily liquid
Identification (I.R.) ..... Conform Ph.Eur.
Acidity or alkalinity ..... Conform Ph.Eur.
Polycyclic Aromatic Hydrocarbons . Conform Ph.Eur.
Solid paraffins ..... Conform Ph.Eur.
Ready carbonizable substances ..... Conform Ph.Eur.
Density at 20°C ..... 0.827 ÷ 0.890
Density at 25°C ..... 0.845 ÷ 0.905
Viscosity at 20°C ..... 110 ÷ 230 mPa x s
Viscosity at 40°C ..... 34.5 ÷ 150 mm 2/s

Code	Size	Packaging	Notes
356601	1 l	Glass bottle	
356608	5 l	Aluminium can	
356603	23 kg	Metal drum	
356607	185 kg	Metal drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Paraffin white soft

• Vaselina bianca • Vaseline blanche • Vaselina blanca • Weiße Vaseline

Synonym:
Vaseline

CAS: 8009-03-8
EEC-N: 232-373-2

Paraffin white soft > ERBApharm - According to pharmacopoeia: BP-NF

ERBApharm

Description	White mass	Aspetto sostanza fusa	Conform BP	USP-NF		Melting point	47.0 ÷ 65.0 ° C
Identification	Positive	Consistenza	60 ÷ 300 BP	Acidity or alkalinity	Conform BP	Sulphated ash	≤0.05 %
Reaction	Conform USP-NF	Ready carbonizable substances	Conform	Polynuclear hydrocarbon	≤ 300 ppm		

Code	Size	Packaging	Notes
388407	1 kg	Metallic can	
388409	5 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Paraffin white soft > RE - Pure

RE

Appearance	Soft whitish mass	Identification	Conform	Drop melting point	50 - 60 °C

Code	Size	Packaging	Notes
388607	1 kg	Plastic bottle	
388609	5 kg	Plastic tank	

Paraformaldehyde

• Paraformaldeide • Paraformaldéhyde • Paraformaldehído • Paraformaldehyd

Synonym:
Polyoxymethylene

$(CH_2O)n$
CAS: 30525-89-4

Classification transport
ONU: 2213
Transport Hazard class: 4.1
Packing group III



Danger
H228-H302-H332-H315-H318-H317-H351-H412
P210-P280-P304+P340-P310a-P305+P351+P338-P330-P362+P364

Paraformaldehyde > RE - Pure

RE

Description	White powder	Sulphated ash	≤ 0.1 %	Acidity or alcalinity	Passes test	Insoluble in NH4OH	Passes test
Identification	Positive	Assay (oxidimetric)	≥ 95 %	Heavy metals (Pb)	≤ 0.001 %		

Code	Size	Packaging	Notes
387507	1 kg	Plastic bottle	
387503	25 kg	Fibre drum	

Pararosaniline solution, decolorised

• Pararosanilina soluzione, decolorata • Pararosaniline décolorée en solution • Pararrosanilina solución, decolorado • Pararosanilin in Lösung verfärbt



Danger
H350-HA26
P201-P202-P280-P308+P313-P405-P501a

Pararosaniline solution, decolorised > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611062201	100 ml	Glass bottle	Ref Ph.Eur 1062201

Storage: protected from light

Product specifications are subject to changes. Please visit our website for updates.

n-Pentane 99%

• n-Pentano 99% • n-Pentane 99% • n-Pentano 99% • n-Pentan 99%

$CH_3(CH_2)_3CH_3$
Molecular Weight: 72,15
CAS: 109-66-0
EEC-N: 203-692-4

Classification transport
ONU: 1265
Transport Hazard class: 3
Packing group II

Danger
H225-H336-H304-H411-HEU066
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

n-Pentane 99% > RS - For GC-MS

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C ..... 1.355 - 1.359
Water (K.F.) ..... ≤ 50 ppm
Residue on evaporation ..... ≤ 2 ppm
Colour ..... ≤ 5 APHA
Acidity (acetic acid) ..... ≤ 10 ppm
Assay (GC) ..... ≥ 99.0 %
GC-MS.Individual peak (n-hexadecane). ≤ 2 µg/L
Ret.range n-undecane to n-tetracontane (scanning area 30-600amu)

Code	Size	Packaging	Notes
468172	1 l	Glass bottle	

n-Pentane 99% > RS - ATRASOL - For traces analysis

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C ..... 1.355 - 1.359
Water content (K.F.) ..... ≤ 50 mg/Kg
Colour ..... ≤ 5 Hazen
Assay (GC) ..... ≥ 99 %
Non volatile residue ..... ≤ 2 mg/Kg
Free acid (as CH3COOH) ..... ≤ 10 mg/Kg
Total sulphur (S) ..... ≤ 10 ppm
GC (FID) - NC Atrasol ..... Conform
GC-ECD.Individual peak (CCl4) ..... ≤ 1 µg/l
Ret.range dichloromethane to 1,2,4-trichlorobenzene
GC-ECD.Individual peak (Lindane) .≤ 2 ng/L
Ret.range 1,2,4-trichlorobenzene to decachlorobiphenyle
GC-FID.Individual peak (n-hexadecane). ≤ 2 µg/L
Ret.range n-undecane to n-tetracontane

Code	Size	Packaging	Notes
P064323016	1 l	Glass bottle	
P064323021	2.5 l	Glass bottle	

n-Pentane 99% > RS - SPECTROSOL - For optical spectroscopy

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Colour ..... ≤10 APHA
Density at 20° C ..... 0.623 ÷ 0.629
Refractive index at 20°C. 1.3552 ÷ 1.3606
Boiling point ..... 35.8 ÷ 36.3 ° C
Acidity or alkalinity ..... ≤0.0002 meq/g
Water (K.F.) ..... ≤50 ppm
Residue on evaporation ..... ≤2 ppm
Total sulphur ..... ≤10 ppm
Aromatic compounds ..... ≤5 ppm
Assay (GLC) ..... ≥99 %
U.V. Transmittance
at 210 nm ..... ≥45 %
at 220 nm ..... ≥89 %
at 230 nm ..... ≥95 %
Trasmittance from 240 nm ..... ≥98 %

Code	Size	Packaging	Notes
468142	1 l	Glass bottle	
468141	2.5 l	Glass bottle	

n-Pentane 99% > RS - For enviromental analysis

RS

Description ..... Clear liquid
Identification ..... Positive
Colour (APHA) ..... ≤ 10
Water ..... ≤ 50 ppm
Not volatile residue ..... ≤ 5 ppm
Free acids (CH3COOH) ..... ≤ 10 ppm
Total sulphur ..... ≤ 10 ppm
GC-ECD (Carbonio tetracloruro) ..... ≤ 1 µg/l
GC-ECD (Lindano) ..... ≤ 2 ng/l
GC-FID (Esadecano) ..... ≤ 5 µg/l
Assay (GLC) ..... ≥ 99 %

Code	Size	Packaging	Notes
468204	1 l	Glass bottle	

n-Pentane 99% > RPE - For analysis - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Colour (APHA) ..... ≤10
Density at 15° C ..... 0.629 ÷ 0.633
Identification (I.R.) ..... Positive
Aromatics ..... ≤0.001 %
Alcohol miscibility ..... Complete
Residue on evaporation ..... ≤0.001 % p/v
Diethyl ether miscib. ..... Complete
Assay (GLC) ..... ≥99 %
Chloroform miscibility ..... Complete
Density at 20° C ..... 0.623 ÷ 0.629
Refractive index at 20°C. 1.3552 ÷ 1.3606
Boiling point ..... 35.8 ÷ 36.3 ° C
Acidity or alkalinity ..... ≤0.0001 meq/g
Water (K.F.) ..... ≤100 ppm
Total phosphorus ..... ≤0.5 ppm
Total silicon ..... ≤0.02 ppm
Total sulphur ..... ≤5 ppm
Ca ..... ≤0.5 ppm
Cu ..... ≤0.05 ppm
Fe ..... ≤0.2 ppm
K ..... ≤0.2 ppm
Mg ..... ≤0.1 ppm
Na ..... ≤1 ppm
Pb ..... ≤0.05 ppm
Zn ..... ≤0.1 ppm

Code	Size	Packaging	Notes
468151000	1 l	Glass bottle	

n-Pentane 99% > RE - Pure

RE

Description ..... Clear liquid	Density at 20°C ..... 0.621 ÷ 0631	Residue on evaporation ..... ≤ 20 ppm	Assay (GLC) ..... ≥ 99 %
Identification ..... Positive	Refractive index at 20°C. 1.3529 ÷ 1.3629	Water (K.F.) ..... ≤ 150 ppm	
Colour ..... ≤ 10 APHA	Boiling point..... 34.75 ÷ 37.25 °C	Aromatic compounds..... ≤ 10 ppm	

Code	Size	Packaging	Notes
528994	1 l	Glass bottle	
528993	2.5 l	Glass bottle	
528995	5 l	Aluminium can	
528996	25 l	Drum	
528997	200 l	Metal drum	

n-Pentane

• n-Pentano • n-Pentane • n-Pentano • n-Pentan

$CH_3(CH_2)_3CH_3$
Molecular Weight: 72,15
CAS: 109-66-0
EEC-N: 203-692-4

Classification transport
ONU: 1265
Transport Hazard class: 3
Packing group II



Danger
H225-H336-H304-H411-HEU066
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

n-Pentane > RS - For HPLC - Isocratic Grade

RS

Refractive index at 20°C..... 1.355 - 1.359	UV transmittance at 210 nm ..... ≥ 5 %	UV transmittance at 300 nm ..... ≥ 98 %	Assay (GC) ..... ≥ 95 %
Water content (K.F.) ..... ≤ 100 mg/Kg	UV transmittance at 230 nm ..... ≥ 80 %	Aromatic compounds..... ≤ 5 mg/Kg	Total sulphur (S) ..... ≤ 2 ppm
Colour ..... ≤ 10 Hazen	UV transmittance at 290 nm ..... ≥ 85 %	Non volatile residue ..... ≤ 5 mg/Kg	Free acid (as CH3COOH)..... ≤ 10 mg/Kg

Code	Size	Packaging	Notes
P0643716	1 l	Glass bottle	
P0643721	2.5 l	Glass bottle	

n-Pentane > RS - For GC-MS

RS

Appearance ..... Clear colourless liquid	Water (K.F.) ..... ≤ 50 ppm	Assay (GC) ..... ≥ 96.5 %	Ret.range n-undecane to n-tetracontane (scanning area 30-600amu)
Refractive index at 20°C..... 1.355 - 1.359	Residue on evaporation ..... ≤ 2 ppm	GC-MS.Individual peak (n-hexadecane). ≤ 2 µg/L	
Density d20/4 ..... 0.621 - 0.631	Colour ..... ≤ 5 APHA		

Code	Size	Packaging	Notes
468182	1 l	Glass bottle	

n-Pentane > RS - ATRASOL - For trace analysis, Suitable for Hydrocarbon index determination

RS

Appearance ..... Clear colourless liquid	Non volatile residue ..... ≤ 2 mg/Kg	GC-FID.Hydrocarbon oil index...≤ 0.05 mg/l	Ret.range 1,2,4-trichlorobenzene to decachlorobiphenyle
Refractive index at 20°C..... 1.355 - 1.359	Acidity (acetic acid)..... ≤ 10 ppm	Retention time n-decane - n-tetracontane	
Water content (K.F.) ..... ≤ 50 mg/Kg	Assay (GC) ..... ≥ 96.5 %	GC-FID.Individual peak (C10-C40)..≤ 2 µg/l	
Colour ..... ≤ 5 Hazen	Density d20/4 ..... 0.621 - 0.631	GC-ECD.Individual peak (Lindane) .≤ 2 ng/L	

Code	Size	Packaging	Notes
P0643216	1 l	Glass bottle	
P0643221	2.5 l	Glass bottle	

According to NF-EN-ISO9377-2 for hydrocarbon index determination

n-Pentane > RS - PESTIPUR - For pesticide analysis

RS

Description ..... Clear liquid	Water ..... ≤ 50 ppm	Total sulphur ..... ≤ 2 ppm	Assay (GLC) ..... ≥ 95 %
Identification ..... Positive	Acidity (acetic acid)..... ≤ 10 ppm	GC-ECD (Lindano) ..... ≤ 3 ng/l	
Colour ..... ≤ 10 hazen	Not volatile residue..... ≤ 5 ppm	GC-NPD (Ethylparation)..... ≤ 3 ng/l	

Code	Size	Packaging	Notes
468161	1 l	Glass bottle	
468162	2.5 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

n-Pentane > RS - Anhydrous - For analysis

RS

Refractive index at 20°C ..... 1.355 - 1.359	Colour ..... ≤ 10 Hazen	n-hexane ..... ≤ 0.4 %	Total sulphur (S) ..... ≤ 2 ppm
Water content (K.F.) ..... ≤ 50 mg/Kg	Assay (GC) ..... ≥ 95 %	Cyclopentane ..... ≤ 1 %	Free acid (as CH3COOH) ..... ≤ 10 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg	Aromatic compounds ..... ≤ 20 mg/Kg	2,2-dimethylbutane ..... ≤ 1 %	

Code	Size	Packaging	Notes
P0641016	1 l	Glass bottle	

n-Pentane > RPE - For analysis

RPE

Refractive index at 20°C ..... 1.355 - 1.359	Colour ..... ≤ 10 Hazen	n-hexane ..... ≤ 0.4 %	Total sulphur (S) ..... ≤ 2 ppm
Water content (K.F.) ..... ≤ 150 mg/Kg	Assay (GC) ..... ≥ 95 %	Cyclopentane ..... ≤ 1 %	Free acid (as CH3COOH) ..... ≤ 10 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg	Aromatic compounds ..... ≤ 20 mg/Kg	2,2-dimethylbutane ..... ≤ 1 %	

Code	Size	Packaging	Notes
468121	1 l	Glass bottle	
468122	2.5 l	Glass bottle	
468123	200 l	Metal drum	

n-Pentane > RE - Pure

RE

Description ..... Clear colourless liquid	Refractive index at 20°C ..... 1.3529 ÷ 1.3629	Assay (GLC) ..... ≥95 %
Identification ..... Positive	Water (K.F.) ..... ≤200 ppm	Colour ..... ≤ 10 APHA
Density at 20° C ..... 0.621 ÷ 0.631	Residue on evaporation ..... ≤ 25 ppm	n-Hexane ..... ≤ 0.4 %

Code	Size	Packaging	Notes
356951	1 l	Glass bottle	
356954	5 l	Aluminium can	
356952	16 kg	Drum	
356953	200 l	Metal drum	

2,4-Pentanedione ▶ Acetylacetone

1-Pentanesulphonic acid sodium salt

• Acido 1-pentansolfonico sale sodico • Acide 1-pentanesulfonique sel sodique
• Acido 1-pentanosulfónico sal sódica • 1-Pentansulfonsäure-Natriumsalz

Synonym:
Sodium pentanesulfonate

$CH_3(CH_2)_4SO_3Na$
Molecular Weight: 174,19
CAS: 22767-49-3
EEC-N: 245-208-4

1-Pentanesulphonic acid sodium salt > RS - For ion pair chromatography

RS

Description ..... White crystalline powder	Absorbance (0,25M)	At 220 nm ..... ≤ 0.04 AU	At 250 nm ..... ≤ 0.01 AU
Water (K.F.) ..... ≤ 2 %	At 200 nm ..... ≤ 0.10 AU	At 230 nm ..... ≤ 0.03 AU	At 260 nm ..... ≤ 0.01 AU
Assay ..... ≥ 98 %	At 210 nm ..... ≤ 0.05 AU	At 240 nm ..... ≤ 0.01 AU	

Code	Size	Packaging	Notes
405841	25 g	Glass bottle	
405842	100 g	Plastic bottle	

1-Pentanesulphonic acid sodium salt > RPE - For analysis

RPE

Description ..... White powder	Identification ..... Positive	Assay ..... ≥95 %

Code	Size	Packaging	Notes
409064	25 g	Glass bottle	

1-Pentanesulphonic acid sodium salt monohydrate

• Acido 1-pentansolfonico sale sodico monoidrato • Acide 1-pentanesulfonique sel sodique monohydrate
• Acido 1-pentanosulfónico sal sódica monohidrato • 1-Pentansulfonsäure-Natriumsalz-Monohydrat

Synonym:
Sodium 1-pentanesulfonate monohydrate

$CH_3(CH_2)_4SO_3Na.H_2O$
Molecular Weight: 192,19
CAS: 207605-40-1

1-Pentanesulphonic acid sodium salt monohydrate > RS - For ion pair chromatography

RS

Description ..... White crystalline powder	Absorbance (0,25M)	At 220 nm ..... ≤ 0.03 AU	At 250 nm ..... ≤ 0.01 AU
Loss on drying ..... 7.0 - 9.0 %	At 200 nm ..... ≤ 0.1 AU	At 230 nm ..... ≤ 0.02 AU	
Assay ..... ≥ 99.0 %	At 210 nm ..... ≤ 0.05 AU	At 240 nm ..... ≤ 0.01 AU	

Code	Size	Packaging	Notes
405891	25 g	Glass bottle	
405892	100 g	Plastic bottle	

Pentanol-1 ▶ n-Amyl alcohol

Pepsin HCl

• Pepsina HCl • Pepsine • Pepsina HCl • Pepsin HCl



Danger
H290-H314-HEU208
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Pepsin HCl > RPE - For analysis

RPE

Code	Size	Packaging	Notes
468245	250 ml	Bottle	

Cleaning solution proteins

Perchloric acid 65-71%

• Acido perclorico 65-71% • Acide perchlorique 65-71% • Acido perclórico 65-71% • Perchlorsäure 65-71%

$HClO_4$
Molecular Weight: 100,47
CAS: 7601-90-3

Classification transport
ONU: 1873
Transport Hazard class: 5.1
Packing group I



Danger
H271-H302-H314
P210-P280-P283-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Perchloric acid 65-71% > RS - Superpure - For trace analysis at ppb level

RS

Description ..... Clear liquid	Dy ..... ≤ 0.5 ppb	Mo ..... ≤ 0.5 ppb	Tb ..... ≤ 0.5 ppb
Identification ..... Positive	Er ..... ≤ 0.5 ppb	Nd ..... ≤ 0.5 ppb	Tl ..... ≤ 0.5 ppb
Colour ..... ≤ 10 APHA	Eu ..... ≤ 0.5 ppb	Ni ..... ≤ 1 ppb	Th ..... ≤ 1 ppb
Al ..... ≤ 1 ppb	Gd ..... ≤ 0.5 ppb	Pd ..... ≤ 0.5 ppb	Tm ..... ≤ 0.5 ppb
Sb ..... ≤ 0.5 ppb	Ga ..... ≤ 0.5 ppb	Pt ..... ≤ 0.5 ppb	Sn ..... ≤ 1 ppb
As ..... ≤ 0.5 ppb	Au ..... ≤ 0.5 ppb	K ..... ≤ 1 ppb	Ti ..... ≤ 1 ppb
Ba ..... ≤ 1 ppb	Ho ..... ≤ 0.5 ppb	Pr ..... ≤ 0.5 ppb	U ..... ≤ 0.5 ppb
Be ..... ≤ 0.5 ppb	In ..... ≤ 0.5 ppb	Rh ..... ≤ 0.5 ppb	V ..... ≤ 0.5 ppb
Bi ..... ≤ 0.5 ppb	Fe ..... ≤ 1 ppb	Rb ..... ≤ 0.5 ppb	Yb ..... ≤ 0.5 ppb
Cd ..... ≤ 1 ppb	La ..... ≤ 0.5 ppb	Sm ..... ≤ 0.5 ppb	Y ..... ≤ 0.5 ppb
Ca ..... ≤ 1 ppb	Pb ..... ≤ 1 ppb	Sc ..... ≤ 0.5 ppb	Zn ..... ≤ 1 ppb
Ce ..... ≤ 0.5 ppb	Li ..... ≤ 0.5 ppb	Ag ..... ≤ 1 ppb	Zr ..... ≤ 0.5 ppb
Cs ..... ≤ 0.5 ppb	Lu ..... ≤ 0.5 ppb	Na ..... ≤ 1 ppb	Assay (acidimetric) ..... 65 ÷ 71 %
Co ..... ≤ 0.5 ppb	Mg ..... ≤ 1 ppb	Sr ..... ≤ 0.5 ppb	
Cu ..... ≤ 0.5 ppb	Mn ..... ≤ 1 ppb	Te ..... ≤ 0.5 ppb	

Code	Size	Packaging	Notes
409193	1 l	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Perchloric acid 65%

• Acido perclorico 65% • Acide perchlorique 65% • Acido perclórico 65% • Perchlorsäure 65%

$HClO_4$
Molecular Weight: 100,47
CAS: 7601-90-3

Classification transport
ONU: 1873
Transport Hazard class: 5.1
Packing group I



Danger
H271-H302-H314
P210-P280-P283-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Perchloric acid 65% > RS - For enviromental analysis - ISO

RS

Description	Clear liquid	Residue on ignition	≤30 ppm	Cr	≤0.2 ppm	Pb	≤0.05 ppm
Colour (APHA)	≤10	Sulphate	≤4 ppm	Cu	≤0.1 ppm	Sr	≤0.02 ppm
Identification	Positive	Ag	≤0.1 ppm	Fe	≤0.5 ppm	Ti	≤0.1 ppm
Density at 20° C	1.587 ÷ 1.607	Al	≤0.05 ppm	Hg	≤0.02 ppm	Tl	≤0.05 ppm
Total nitrogen	≤10 ppm	As	≤0.05 ppm	K	≤0.1 ppm	V	≤0.05 ppm
Chlorate	≤10 ppm	Ba	≤0.02 ppm	Li	≤0.02 ppm	Zn	≤0.1 ppm
Free chlorine	≤0.5 ppm	Bi	≤0.1 ppm	Mg	≤0.5 ppm	Assay (acidimetric)	64 ÷ 66 %
Chloride	≤1 ppm	Ca	≤0.5 ppm	Mn	≤0.02 ppm		
Fluoride	≤1 ppm	Cd	≤0.005 ppm	Mo	≤0.05 ppm		
Phospate-silicate(SiO2)	≤1 ppm	Co	≤0.05 ppm	Ni	≤0.1 ppm		

Code	Size	Packaging	Notes
409121	1 l	Glass bottle	

Low content in Hg

Perchloric acid 65% > RPE - For analysis - ISO

RPE

Description	Clear liquid	Heavy metals (Pb)	≤1 ppm	Cd	≤0.05 ppm	Ni	≤0.1 ppm
Colour (APHA)	≤10	Residue on ignition	≤30 ppm	Co	≤0.05 ppm	Pb	≤0.05 ppm
Identification	Positive	Sulphate	≤5 ppm	Cu	≤0.1 ppm	Sr	≤0.02 ppm
Density at 20° C	1.587 ÷ 1.607	Ag	≤0.1 ppm	Fe	≤0.5 ppm	Ti	≤0.1 ppm
Total nitrogen	≤10 ppm	Al	≤0.05 ppm	K	≤0.1 ppm	Tl	≤0.05 ppm
Chlorate	≤10 ppm	As	≤0.05 ppm	Li	≤0.02 ppm	V	≤0.05 ppm
Free chlorine	≤0.5 ppm	Ba	≤0.02 ppm	Mg	≤0.5 ppm	Zn	≤0.1 ppm
Chloride	≤1 ppm	Bi	≤0.1 ppm	Mn	≤0.02 ppm	Assay (acidimetric)	64 ÷ 66 %
Phospate-silicate(SiO2)	≤5 ppm	Ca	≤0.5 ppm	Mo	≤0.05 ppm		

Code	Size	Packaging	Notes
409111	1 l	Glass bottle	
409113	35 kg	Drum	

Perchloric acid 65% > RE - Pure

RE

Description	Clear colourless liquid	Total nitrogen	≤50 ppm	Heavy metals (Pb)	≤1 ppm	Fe	≤5 ppm
Identification	Positive	Chloride	≤10 ppm	Residue on ignition	≤50 ppm	Assay (acidimetric)	64 ÷ 66 %
Density at 20° C	1.587 ÷ 1.607	Phospate-silicate(SiO2)	≤50 ppm	Sulphate	≤50 ppm		

Code	Size	Packaging	Notes
306091	1 l	Glass bottle	

Perchloric acid 0.1 mol/l (0.1N)

• Acido perclorico 0.1 mol/l (0.1N) • Acide perchlorique 0.1 mol/l (0.1N) • Acido perclórico 0.1 mol/l (0.1N) • Perchlorsäure 0.1 mol/l (0.1N)

Classification transport
ONU: 2789
Transport Hazard class: 8
Packing group II



Danger
H226-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Perchloric acid 0.1 mol/l (0.1N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613003900	1 l	Glass bottle	Ref Ph.Eur 3003900

Perchloric acid 0.1 mol/l (0.1N) in acetic acid

• Acido perclorico 0.1 mol/l (0.1N) in acido acetico • Acide perchlorique 0.1 mol/l (0.1N) dans l'acide acétique
• Acido perclórico 0.1 mol/l (0.1N) en acido acético • Perchlorsäure 0.1 mol/l (0.1N) essigsäure

Classification transport
ONU: 2789
Transport Hazard class: 8
Packing group II



Danger
H226-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Perchloric acid 0.1 mol/l (0.1N) in acetic acid > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.0998 - 0.1002 N NIST 84 I

Code	Size	Packaging	Notes
409136	500 ml	Glass bottle	Certified with NIST traceability
409131	1 l	Glass bottle	Certified with NIST traceability

10,046 g of HClO4. Volumetric solution ready-to-use. Content is guaranteed for standardized volumes at 20 °C

Perchloric acid 0.05 mol/l (0.05N)

• Acido perclorico 0.05 mol/l (0.05N) • Acide perchlorique 0.05 mol/l (0.05N) • Acido perclórico 0.05 mol/l (0.05N) • Perchlorsäure 0.05 mol/l (0.05N)

Classification transport
Transport Hazard class: 8
Packing group II



Perchloric acid 0.05 mol/l (0.05N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613004000	1 l	Glass bottle	Ref Ph.Eur 3004000

Perchloric acid 0.01 mol/l (0.01N)

• Acido perclorico 0.01 mol/l (0.01N) • Acide perchlorique 0.01 mol/l (0.01N) • Acido perclórico 0.01 mol/l (0.01N) • Perchlorsäure 0.01 mol/l (0.01N)

Classification transport
ONU: 2789
Transport Hazard class: LQ



Danger
H226-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Perchloric acid 0.01 mol/l (0.01N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
E409141	500 ml	Glass bottle	

1.0046 g of HClO4. Volumetric solution ready-to-use. Solution in acetic anhydride

Perchloric acid solution

• Acido perclorico soluzione • Acide perchlorique solution • Acido perclórico solución • Perchlorsäurelösung

Perchloric acid solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611062901	100 ml	Plastic bottle	Ref Ph.Eur 1062901

Product specifications are subject to changes. Please visit our website for updates.

Periodic acid

• Acido periodico • Acide périodique • Acido periódico • Periodsäure

$HIO_4.2H_2O$
Molecular Weight: 227,94
CAS: 10450-60-9
EEC-N: 233-937-0

Classification transport
ONU: 3084
Transport Hazard class: 8
Packing group II



Danger
H272-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Periodic acid > RPE - For analysis

RPE

Description ..... Colorless or White crystals
Assay ..... ≥ 99 %
Iodide ..... ≤ 0.005 %
Water insolubles ..... ≤ 0.02 %
Chloride ..... ≤ 0.02 %
Residue on ignition (SO4) ..... ≤ 0.1 %

Code	Size	Packaging	Notes
409182	25 g	Glass bottle	
409184	100 g	Glass bottle	
409185	250 g	Glass bottle	

Petroleum

• Petrolio • Pétrole • Petróleo • Petroleum

CAS: 64771-72-8
EEC-N: 929-018-5



Danger
H304-HEU066
P301+P310a-P331-P405-P501a

Petroleum > RE - Pure

RE

Description ..... Clear colourless liquid
Identification (I.R.) ..... Positive
Density at 15°C ..... 0.746 ÷ 0.752
Boiling point ..... 185 ÷ 245 °C
Residue on ignition ..... ≤ 100 ppm

Code	Size	Packaging	Notes
357151	1 l	Glass bottle	
357155	21 kg	Metal drum	

Petroleum benzin E

• Benzina E • Essence E • Benzina de Petroleo E • Benzin E

CAS: 64742-49-0

Classification transport
ONU: 3295
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H336-H304-H411
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Petroleum benzin E > RE - Pure

RE

Density d15/4 ..... 0.720 - 0.745
Colour ..... ≤ 10 Hazen
Aromatic compounds ..... ≤ 100 mg/Kg
Boiling point ..... 100 - 140 °C
Total sulphur (S) ..... ≤ 5 ppm

Code	Size	Packaging	Notes
P0370048	25 l	Metal drum	

Petroleum ether 100 - 140°C

• Etere di petrolio (Ligroina) 100 - 140°C • Ether de pétrole 100 - 140°C • Eter de petróleo 100 - 140°C • Petrolether 100 - 140°C

CAS: 64742-49-0
EEC-N: 265-151-9

Classification transport
ONU: 3295
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H336-H304-H411-HEU066
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Petroleum ether 100 - 140°C > RE - Pure

RE

Description ..... Clear colourless liquid
Water (K.F.) ..... ≤ 200 ppm
Identification ..... Positive
Residue on evaporation ..... ≤ 50 ppm
Density at 15° C ..... 0.725 ÷ 0.740
Total sulphur ..... ≤ 10 ppm
Aromatics (Thiophene) ..... ≤ 0.06 %
Boling point min. ..... ≥100 °C
Boiling point max. ..... ≤140 °C

Code	Size	Packaging	Notes
348913	1 l	Glass bottle	
348912	2.5 l	Glass bottle	
508230	5 l	Plastic tank	
348914	20 kg	Metal drum	
508232	25 l	Metal drum	

Petroleum ether 80 - 120°C

• Etere di petrolio (Ligroina) 80 - 120°C • Ether de pétrole 80 - 120°C • Eter de petróleo 80 - 120°C • Petrolether 80 - 120°C

CAS: 64742-49-0

Classification transport
ONU: 1268
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H336-H411
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Petroleum ether 80 - 120°C > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 15° C ..... ≥0.723
Residue on evaporation ..... ≤100 ppm
Boling point min. ..... ≥80 °C
Boiling point max. ..... ≤120 °C

Code	Size	Packaging	Notes
348901	1 l	Glass bottle	
348905	2.5 l	Glass bottle	

Petroleum ether 80 - 100°C

• Etere di petrolio (Benzina) 80 - 100°C • Ether de pétrole 80 - 100°C • Eter de petróleo 80 - 100°C • Petrolether 80 - 100°C

CAS: 64742-49-0

Classification transport
ONU: 1268
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H336-H304-H411
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Petroleum ether 80 - 100°C > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Boiling point ..... 80 ÷ 100 ° C
Carb.sulf. miscibility ..... Complete
CCl4 miscibility ..... Complete
Anhyd.Ethyl alc.miscib. ..... Complete
Benzene miscibility ..... Complete
Diethyl ether miscib. ..... Complete
Ready carbonizable substances ..... Conform
Density at 15° C ..... ~ 0.709
Water (K.F.) ..... ≤0.01 %
Residue on evaporation ..... ≤10 ppm
Acidity (acetic acid) ..... ≤5 ppm
Aromatics ..... ≤100 ppm
Subst. reducing KMnO4 ..... ≤20 ppm (5m)
Total sulphur ..... ≤50 ppm

Code	Size	Packaging	Notes
427031	1 l	Glass bottle	
427036	18 kg	Metal drum	

Product specifications are subject to changes. Please visit our website for updates.

Petroleum ether 80 - 100°C > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 15° C ..... ~ 0.708
Residue on evaporation ..... ≤100 ppm
Boling point min. ..... ≥80 °C
Boiling point max ..... ≤100 °C

Code	Size	Packaging	Notes
323501	1 l	Glass bottle	
323503	2.5 l	Glass bottle	
323502	19 kg	Metal drum	

Petroleum ether 75 - 120°C

• Etere di petrolio (Ligroina) 75 - 120°C • Ether de pétrole 75 - 120°C • Eter de petróleo 75 - 120°C • Petrolether 75 - 120°C

CAS: 64742-49-0

Classification transport
ONU: 1268
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H336-H304-H411
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Petroleum ether 75 - 120°C > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Diethyl ether miscib. ..... Complete
Misc.with Abs.Ethanol ..... Complete
Boling point min. ..... ≥75 °C
Boiling point max ..... ≤120 °C
Ready carbonizable substances ..... Conform
Density at 15° C ..... ≥0.715
Water (K.F.) ..... ≤50 ppm
Residue on evaporation ..... ≤10 ppm
Acidity (acetic acid) ..... ≤0.7 ppm
Alcalinity (NH3) ..... ≤0.2 ppm
Subst. reducing KMnO4 ..... ≤20 ppm (5m)
Total sulphur ..... ≤50 ppm

Code	Size	Packaging	Notes
458001	1 l	Glass bottle	
458003	2.5 l	Glass bottle	

Petroleum ether 60 - 80°C

• Etere di petrolio (Benzina) 60 - 80°C • Ether de pétrole 60 - 80°C • Eter de petróleo 60 - 80°C • Petrolether 60 - 80°C

CAS: 64742-49-0

Classification transport
ONU: 1268
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H361f-H336-H373-H304-H411
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Petroleum ether 60 - 80°C > RPE - For analysis

RPE

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Boiling point ..... 60 ÷ 80 ° C
Ready carbonizable substances ..... Conform
Density at 20° C ..... 0.660 ÷ 0.690
Water (K.F.) ..... ≤100 ppm
Residue on evaporation ..... ≤10 ppm
Acidity (acetic acid) ..... ≤5 ppm
Subst. reducing KMnO4 ..... ≤20 ppm (5m)
Total sulphur ..... ≤50 ppm

Code	Size	Packaging	Notes
427001	1 l	Glass bottle	
427003	2.5 l	Glass bottle	
427007	18 kg	Metal drum	

Petroleum ether 55 - 85°C

• Etere di petrolio (Benzina) 55 - 85°C • Ether de pétrole 55 - 85°C • Eter de petróleo 55 - 85°C • Petrolether 55 - 85°C

CAS: 64742-49-0

Classification transport
ONU: 1268
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H361f-H336-H373-H411
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Petroleum ether 55 - 85°C > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20°C ..... 0.660 ÷ 0.690
Residue on evaporation ..... ≤100 ppm
Boling point min. ..... ≥55 °C
Boiling point max ..... ≤85 °C

Code	Size	Packaging	Notes
323401	1 l	Glass bottle	
323403	2.5 l	Glass bottle	
323402	18 kg	Metal drum	

Petroleum ether 40 - 70°C

• Etere di petrolio 40 - 70°C • Ether de pétrole 40 - 70°C • Eter de petróleo 40 - 70°C • Petrolether 40 - 70°C

CAS: 64742-49-0

Classification transport
ONU: 1268
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H336-H304-H411
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Petroleum ether 40 - 70°C > RPE - For analysis

RPE

Description	Clear liquid	Residue on evaporation	≤20 ppm	Total sulphur	≤50 ppm	Fe	≤0.1 ppm
Colour (APHA)	≤10	Water (K.F.)	≤100 ppm	Al	≤0.5 ppm	Mg	≤0.1 ppm
Identification	Positive	Boling point min.	≥40 °C	Ba	≤0.1 ppm	Mn	≤0.02 ppm
Alcohol miscibility	Complete	Boiling point max	≤70 °C	Ca	≤0.5 ppm	Ni	≤0.02 ppm
Diethyl ether miscib.	Complete	Acidity (acetic acid)	≤0.7 ppm	Cd	≤0.05 ppm	Pb	≤0.1 ppm
Fat oils	Conform	Alcalinity (NH3)	≤0.2 ppm	Co	≤0.02 ppm	Sn	≤0.1 ppm
Ready carbonizable substances	Conform	Benzene	≤100 ppm	Cr	≤0.02 ppm	Zn	≤0.1 ppm
Density at 15° C	0.645 ÷ 0.670	Subst. reducing KMnO4	≤20 ppm (5m)	Cu	≤0.02 ppm		

Code	Size	Packaging	Notes
447821	1 l	Glass bottle	
447824	5 l	Aluminium can	
447822	19 kg	Aluminium can	

Petroleum ether 40 - 70°C > RE - Pure

RE

Description	Clear liquid	Density at 15° C	0.645 ÷ 0.670	Benzene	≤200 ppm
Colour (APHA)	≤ 10	Residue on evaporation	≤20 ppm	Boling point min.	≥40 °C
Identification	Positive	Water (K.F.)	≤100 ppm	Boiling point max	≤70 °C

Code	Size	Packaging	Notes
341024	1 l	Glass bottle	
341022	19 kg	Aluminium can	

Petroleum ether 40 - 65°C

• Etere di petrolio 40 - 65°C • Ether de pétrole 40 - 65°C • Eter de petróleo 40 - 65°C • Petrolether 40 - 65°C

CAS: 64742-49-0

Classification transport
ONU: 1268
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H336-H304-H411
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Petroleum ether 40 - 65°C > RS - PESTIPUR - For pesticide analysis

RS

Appearance	Clear colourless liquid	Water content (K.F.)	≤ 100 mg/Kg	Ret.range 1,2,4-trichlorobenzene to decachlorobiphenyle		Retention time Atrazin to Coumaphos	
Boiling point	40 - 65 °C	Non volatile residue	≤ 2 mg/Kg	GC-NPD.Individual peak (Ethylparathion)	≤ 3 ng/l		
Density d15/4	0.640 - 0.655	GC chromatogram	Conform				
Colour	≤ 10 Hazen	GC-ECD.Individual peak (Lindane) .	≤ 3 ng/L				

Code	Size	Packaging	Notes
447851	1 l	Glass bottle	
447852	2.5 l	Glass bottle	

For chlorinated compounds analysis

Petroleum ether 40 - 65°C > RPE - For analysis

RPE

Description	Clear colourless liquid	Refractive index at 20°C	1.366 ÷ 1.376	Water (K.F.)	≤ 100 ppm	Aromatic compounds	≤ 100 ppm
Colour	≤ 10 APHA	Boiling point	40 ÷ 65 °C	Assay (CPG)	Conform		
Density at 15°C	0.645 ÷ 0.660	Residue on evaporation	≤ 10 ppm	n-Hexane	≤ 2 %		

Code	Size	Packaging	Notes
447811	1 l	Glass bottle	
447812	2.5 l	Glass bottle	
447813	5 l	Plastic tank	
447814	10 l	Metal tank	
447815	25 l	Metal drum	
447816	200 l	Metal drum	

Product specifications are subject to changes. Please visit our website for updates.

Petroleum ether 40 - 60°C

• Etere di petrolio 40 - 60°C • Ether de pétrole 40 - 60°C • Eter de petróleo 40 - 60°C • Petrolether 40 - 60°C

CAS: 64742-49-0

Classification transport
ONU: 1268
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H336-H304-H411
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Petroleum ether 40 - 60°C > RPE - For analysis

RPE

Description	Clear liquid	Ready carbonizable substances	Conform	Total sulphur	≤50 ppm	Fe	≤0.1 ppm
Colour (APHA)	≤10	Density at 15° C	0.647 ÷ 0.654	Al	≤0.5 ppm	Mg	≤0.1 ppm
Identification	Positive	Residue on evaporation	≤10 ppm	Ba	≤0.1 ppm	Mn	≤0.02 ppm
Alcohol miscibility	Complete	Water (K.F.)	≤100 ppm	Ca	≤0.5 ppm	Ni	≤0.02 ppm
Diethyl ether miscib.	Complete	Acidity (acetic acid)	≤0.7 ppm	Cd	≤0.05 ppm	Pb	≤0.1 ppm
Boling point min.	≥40 °C	Alcalinity (NH3)	≤0.2 ppm	Co	≤0.02 ppm	Sn	≤0.1 ppm
Boiling point max	≤60 °C	Benzene	≤100 ppm	Cr	≤0.02 ppm	Zn	≤0.1 ppm
Fat oils	Conform	Subst. reducing KMnO4	≤20 ppm (5m)	Cu	≤0.02 ppm	Bromine rating	≤ 1

Code	Size	Packaging	Notes
447833	1 l	Glass bottle	
447831	2.5 l	Glass bottle	
447832	5 l	Aluminium can	
447836	5 l	Plastic tank	
447834	19 kg	Aluminium can	

Petroleum ether 40 - 60°C > RE - Pure

RE

Description	Clear liquid	Density at 20°C	0.643 ÷ 0.673	Water (K.F.)	≤ 200 ppm
Colour	≤ 10 APHA	Refractive index at 20°C	1.368 ÷ 1.378	Aromatics	≤100 ppm
Identification	Positive	Residue on evaporation	≤ 20 ppm	Distillation intervalle	40 ÷ 60 °C

Code	Size	Packaging	Notes
528283	5 l	Plastic tank	

Petroleum ether 35 - 60°C

• Etere di petrolio 35 - 60°C • Ether de pétrole 35 - 60°C • Eter de petróleo 35 - 60°C • Petrolether 35 - 60°C

CAS: 109-66-0

Classification transport
ONU: 1265
Transport Hazard class: 3
Packing group II



Danger
H225-H336-H304-H411-HEU066
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Petroleum ether 35 - 60°C > RS - ATRASOL - For trace analysis, Suitable for Hydrocarbon index determination

RS

Appearance	Clear colourless liquid	Colour	≤ 5 Hazen	Ret.range 1,2,4-trichlorobenzene to decachlorobiphenyle	Retention time n-decane - n-tetracontane
Refractive index at 20°C	1.355 - 1.359	Non volatile residue	≤ 2 mg/Kg	GC-FID.Hydrocarbon oil index	≤ 0.05 mg/l	GC-FID.Individual peak (C10-C40)	≤ 5 µg/l
Water content (K.F.)	≤ 50 mg/Kg	GC-ECD.Individual peak (Lindane)	≤ 3 ng/L		

Code	Size	Packaging	Notes
P0883216	1 l	Glass bottle	
P0883221	2.5 l	Glass bottle	

According to NF-EN-ISO9377-2 for hydrocarbon index determination. Contains n-Pentane, iso-Pentane, 1-pentene, 2,2-Dimethylbutane and cyclopentane

Petroleum ether 35 - 60°C > RS - PESTIPUR - For pesticide analysis

RS

Clear,colourless liq.appearance	Conform	Refractive index at 20°C	1.355 - 1.359	Retention time trichlorobenzene to mirex	GC-NPD.Individual peak (Ethylparathion) ≤ 3 ng/l
Identification	Conform	Water content (K.F.)	≤ 100 mg/Kg	Non volatile residue	≤ 2 mg/Kg	Retention time Atrazin to Coumaphos
Colour	≤ 10 Apha	GC-ECD.Individual peak (Lindane)	≤ 3 ng/l	Total sulphur (S)	≤ 10 ppm	

Code	Size	Packaging	Notes
447862	1 l	Glass bottle	
447861	2.5 l	Glass bottle	

For chlorinated and nitrogenous compounds analysis. Contains n-Pentane, iso-Pentane, 1-pentene, 2,2-Dimethylbutane and cyclopentane

Petroleum ether 35 - 60°C > RPE - For analysis

RPE

Description	Clear liquid	Density at 20°C	0.643 ÷ 0.673	Water (K.F.)	≤ 150 ppm	Bromine rating	≤ 1
Colour	≤ 10 APHA	Refractive index at 20°C	1.368 ÷ 1.378	Aromatics	≤ 20 ppm	Distillation intervalle	35 ÷ 60 °C
Identification	Positive	Residue on evaporation	≤ 10 ppm	Total sulphur	≤ 10 ppm	Assay (CPG)	Conform

Code	Size	Packaging	Notes
528070	1 l	Glass bottle	
528071	2.5 l	Glass bottle	
528280	5 l	Plastic tank	
528281	25 l	Metal drum	
528282	200 l	Metal drum	

Petroleum ether 30 - 50°C

• Etere di petrolio 30 - 50°C • Ether de pétrole 30 - 50°C • Eter de petróleo 30 - 50°C • Petrolether 30 - 50°C

CAS: 109-66-0

Classification transport
ONU: 1268
Transport Hazard class: 3
Packing group II



Danger
H225-H336-H304-H411
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Petroleum ether 30 - 50°C > RPE - For analysis

RPE

Description	Clear liquid	Residue on evaporation	≤10 ppm	Subst. reducing KMnO4	≤20 ppm (5m)	Cu	≤0.02 ppm
Colour (APHA)	≤10	Water (K.F.)	≤100 ppm	Total sulphur	≤50 ppm	Fe	≤0.1 ppm
Identification	Positive	Boling point min.	≥30 °C	Al	≤0.5 ppm	Mg	≤0.1 ppm
Alcohol miscibility	Complete	Boiling point max	≤50 °C	Ba	≤0.1 ppm	Mn	≤0.02 ppm
Diethyl ether miscib.	Complete	Acidity (acetic acid)	≤0.7 ppm	Ca	≤0.5 ppm	Ni	≤0.02 ppm
Fat oils	Conform	Alcalinity (NH3)	≤0.2 ppm	Cd	≤0.05 ppm	Pb	≤0.1 ppm
Ready carbonizable substances	Conform	Benzene	≤100 ppm	Co	≤0.02 ppm	Sn	≤0.1 ppm
Density at 15° C	≥0.633	Unsaturated hydrocarbon	≤0.4 %	Cr	≤0.02 ppm	Zn	≤0.1 ppm

Code	Size	Packaging	Notes
447801	1 l	Glass bottle	
447804	5 l	Aluminium can	
447802	18 kg	Aluminium can	

Petroleum ether 30 - 50°C > RE - Pure

RE

Description	Clear colourless liquid	Density at 15° C	≥0.633	Boling point min.	≥30 °C
Identification	Positive	Residue on evaporation	≤100 ppm	Boiling point max	≤50 °C

Code	Size	Packaging	Notes
341034	1 l	Glass bottle	
341032	18 kg	Aluminium can	

Product specifications are subject to changes. Please visit our website for updates.

Petroleum ether 30 - 40°C

• Etere di petrolio 30 - 40°C • Ether de pétrole 30 - 40°C • Eter de petróleo 30 - 40°C • Petrolether 30 - 40°C

CAS: 109-66-0

Classification transport
ONU: 1265
Transport Hazard class: 3
Packing group II



Danger
H225-H336-H304-H411-HEU066
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Petroleum ether 30 - 40°C > RPE - For analysis

RPE

Description	Clear liquid	Ready carbonizable substances	Conform	Subst. reducing KMnO4	≤20 ppm (5m)	Cu	≤0.02 ppm
Colour (APHA)	≤10	Density at 15° C	≥0.630	Total sulphur	≤50 ppm	Fe	≤0.1 ppm
Identification	Positive	Water (K.F.)	≤100 ppm	Al	≤0.5 ppm	Mg	≤0.1 ppm
Alcohol miscibility	Complete	Residue on evaporation	≤10 ppm	Ba	≤0.1 ppm	Mn	≤0.02 ppm
Diethyl ether miscib.	Complete	Acidity (acetic acid)	≤0.7 ppm	Ca	≤0.5 ppm	Ni	≤0.02 ppm
Boling point min.	≥30 °C	Alcalinity (NH3)	≤0.2 ppm	Cd	≤0.05 ppm	Pb	≤0.1 ppm
Boiling point max	≤40 °C	Benzene	≤200 ppm	Co	≤0.02 ppm	Sn	≤0.1 ppm
Fat oils	Conform	Unsaturated hydrocarbon	≤0.4 %	Cr	≤0.02 ppm	Zn	≤0.1 ppm

Code	Size	Packaging	Notes
447793	1 l	Glass bottle	
447795	5 l	Aluminium can	
447792	18 kg	Aluminium can	

Phenacetin ▶ p-Acetylphenetidine

o-Phenantroline monohydrate

• o-Fenantrolina monoidrata • o-Phénanthroline monohydraté • o-Fenantrolina monohidrato
• o-Phenanthrolinmonohydrat

Synonym:
1,10-Phenanthroline monohydrate

$C_{12}H_8N_2.H_2O$
Molecular Weight: 198,21
CAS: 5144-89-8
EEC-N: 200-629-2

Classification transport
ONU: 3143
Transport Hazard class: 6.1
Packing group III



Danger
H301-H410
P264-P270-P301+P310a-P330-P405-P501a

o-Phenantroline monohydrate > RPE - For analysis - ACS

RPE

Description	White - pink powder	Suitab. as indicator	Conform	Assay	≥ 99.0 %
Identification	Positive	Iron sensitivity	Conform		

Code	Size	Packaging	Notes
450038	5 g	Glass bottle	
450039	25 g	Glass bottle	

Redox indicator

1,10-Phenanthroline iron(II) sulfate complex ▶ Ferroin 0.025 mol/l solution

o-Phenantroline-Iron (II) sulphate solution in sulphuric acid

• o-Fenantrolina - Ferro solfato oso soluzione in acido solforico • o-Phénanthroline-Fer (II) sulfate
• o-Fenantrolina-Hierro (II) sulfato solución en acido sulfúrico • o-Phenanthrolin-Eisen (II) sulfat

Synonym:
• *1,10-Phenanthroline iron(II) sulfate complex*
• *Ferroin*

$[Fe(C_{12}H_8N_2)_3]SO_4$
Molecular Weight: 692,52
CAS: 14634-91-4

H412
P273-P501a

o-Phenantroline-Iron (II) sulphate solution in sulphuric acid > RPE - For analysis

RPE

Description ..... Red clear liquid
Identification ..... Positive

Code	Size	Packaging	Notes
E450043	100 ml	Glass bottle	

Phenol

• Fenolo • Phénol • Fenol • Phenol

Synonym:
Hydroxybenzene

C_6H_5OH
Molecular Weight: 94,11
CAS: 108-95-2
EEC-N: 203-632-7

Classification transport
ONU: 1671
Transport Hazard class: 6.1
Packing group II



Danger
H301-H311-H331-H314-H341-H373
P280-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338-P308+P313-P361+P364-P403+P233

Phenol > RPE - For analysis - ACS

RPE

Code	Size	Packaging	Notes
451271	1 kg	Plastic bottle	

For specifications, contact our customer service for a certificate of analysis

Phenol > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP-DAB

ERBApharm

Description ..... Crystalline mass
Identification ..... Positive
Acidity ..... Conform Ph.Eur.
Appearance of solution ..... Conform Ph.Eur.
Reaction, solution app. ..... Conform USP-NF
Freezing point ..... ≥ 39.5 °C
Water (K.F.) ..... ≤0.5 %
Not volatile residue ..... ≤0.05 %
Assay (iodometric) ..... 99.0 ÷ 100.5 %

Code	Size	Packaging	Notes
343407	1 kg	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Phenol liquified 85%

• Fenolo liquido 85% • Phénol liquide 85% • Fenol liquido 85% • Phenol flüssiggas 85%

Synonym:
Hydroxybenzene

C_6H_5OH
Molecular Weight: 94,11
CAS: 108-95-2

Classification transport
ONU: 2821
Transport Hazard class: 6.1
Packing group II



Danger
H301-H311-H330-H314-H341-H373
P280-P284-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Phenol liquified 85% > RE - Pure

RE

Description ..... Clear liquid
Identification ..... Positive
Reaction ..... Conform
Not volatile residue ..... ≤0.05 %
Assay ..... 82.0 ÷ 86.5 %

Code	Size	Packaging	Notes
343411	1 l	Glass bottle	

Phenol red

• Rosso fenolo • Rouge de phénol • Rojo de fenol • Phenolrot

Synonym:
Phenolsulfonphthalein

$C_{19}H_{14}O_5S$
Molecular Weight: 354,38
CAS: 143-74-8
EEC-N: 205-609-7



Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Phenol red > RPE - For analysis - ACS

RPE

Description ..... Red crystalline powder
Identification ..... Positive
Appearance of solution ..... Conform ACS
Colour change ..... yellow red
pH range ..... 6.8 - 8.2

Code	Size	Packaging	Notes
476838	5 g	Glass bottle	
476839	25 g	Glass bottle	

Clark indicator series. Dye for microscopy

Product specifications are subject to changes. Please visit our website for updates.

Phenol Red solution 0.2% in ethanol

• Rosso fenolo soluzione 0,2% in alcole etilico • Rouge de phénol solution 0.2% dans l'éthanol • Rojo de fenol solución 0.2% en alcohol etilico • Phenolrote Lösung 0.2% in Ethanol

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Phenol Red solution 0.2% in ethanol > RPE - For analysis

RPE

Description .....Red clear liquid Identification ..... Positive Sensitivity(6.8-8.4) .....Conform Colour change..... yellow red

Code	Size	Packaging	Notes
E476845	250 ml	Glass bottle	

Indicator series Clark indicator acid-base

Phenol red solution

• Rosso fenolo solzione • Rouge de phénol solution • Rojo de fenol solución • Phenolrote Lösung

HEU210

Phenol red solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611063601	100 ml	Plastic bottle	Ref Ph.Eur 1063601 / Colour change: pH 6.8 (yellow) to pH 8.4 (reddish-violet)
611063603	500 ml	Plastic bottle	Phenol red solution R2 Ref Ph.Eur 1063603

Phenolphthalein

• Fenolftaleina • Phénolphtaléine • Fenolftaleína • Phenolphthalein

Synonym:
3,3-bis(4-Hydroxyphenyl)-1(3H)isobenzofuranone

$C_{20}H_{14}O_4$
Molecular Weight: 318,33
CAS: 77-09-8
EEC-N: 201-004-7

Danger
H341-H350-H361f-HA26
P201-P202-P280-P308+P313-P405-P501a

Phenolphthalein > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White powder Appear. of alcohol sol..... Conform pH range ..... 8.0 ÷ 10.0
Identification ..... Positive Colour change..... Colourless-red

Code	Size	Packaging	Notes
451154	100 g	Plastic bottle	
451156	500 g	Plastic bottle	

Phenolphthalein solution 1% in ethanol

• Fenolftaleina soluzione 1% in alcole etilico • Phénolphthaléine solution à 1% dans l'éthanol • Fenolftaleína solución 1% en alcohol etilico • Phenolphthaleinlösung 1% in Ethanol

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group III

Danger
H226-H319-H341-H350-HA26
P210-P241-P280-P303+P361+P353-P305+P351+P338-P308+P313

Phenolphthalein solution 1% in ethanol > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611063703	100 ml	Plastic bottle	Phenolphtalein solution R1 Ref Ph.Eur 1063703

Phenolphthalein solution 1% in ethanol > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
617000131	100 ml	Plastic bottle	Phenolphtalein TS

Phenolphthalein solution 1% in ethanol > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Sensit.(pH 8.3-10.0)..... Conform Colour change..... colurless-purple

Code	Size	Packaging	Notes
451191	250 ml	Glass bottle	
E451191	250 ml	Glass bottle	Only for italian market
451192	1 l	Glass bottle	
E451192	1 l	Glass bottle	Only for italian market

Acid-base indicator

Phenolphthalein solution 0.1%

• Fenolftaleina soluzione 0.1% • Phénolphtaléine solution 0.1% • Fenolftaleína solución 0.1% • Phenolphthaleinlösung 0.1%

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II

Danger
H225
P210-P240-P241-P280-P303+P361+P353-P501a

Phenolphthalein solution 0.1% > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611063702	100 ml	Plastic bottle	Ref Ph.Eur 1063702

Phenolphtalein solution 0.1% in ethanol

• Fenolftaleina soluzione 0.1% in alcole etilico • Phénolphthaléine solution 0.1% dans l'éthanol • Fenolftaleína solución 0.1% en etanol • Phenolphthaleinlösung 0.1% in Ethanol

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II

Danger
H225
P210-P240-P241-P280-P303+P361+P353-P501a

Phenolphtalein solution 0.1% in ethanol > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611063709	1 l	Bottle	Ref Ph.Eur 1063702

L-Phenylalanine

• L-Fenilalanina • L-Phénylalanine • L-Fenilalanina • L-Phenylalanin

Synonym:
(S)-2-Amino-3-phenylpropionic acid

$C_6H_5CH_2CH(NH_2)COOH$
Molecular Weight: 165,19
CAS: 63-91-2
EEC-N: 200-568-1

L-Phenylalanine > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Potere rotator. specif.(C=2 in Acqua) .....-33.0 ÷ -35.2 ° s.s.
Pb ..... ≤ 10 ppm
pH solution 1% ..... 5.4 ÷ 6.0
Loss on drying ..... ≤ 0.3 %
Heavy metals (Pb)..... ≤ 20 ppm
Residue on ignition..... ≤ 0.1 %
Assay (non-aqueos medium)..... 98.5 ÷ 102.0 % (s.s.)

Code	Size	Packaging	Notes
450328	5 g	Glass bottle	
450329	100 g	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

2-Phenylethanol

• Alcole 2-feniletilico • 2-Phényléthanol • Alcohol 2-feniletílico • 2-Phenylethanol

Synonym:
2-Phenylethyl alcohol

$C_6H_5CH_2CH_2OH$
Molecular Weight: 122,17
CAS: 60-12-8
EEC-N: 200-456-2

Warning
H302-H319
P264-P280i-P301+P312a-P305+P351+P338-P337+P313-P501a

2-Phenylethanol > ERBApharm - According to pharmacopoeia: USP

ERBApharm

Description ..... Clear colourless liquid
Identification ..... Positive
Chlorinated compounds ..... Conform USP-NF
Aldehyde ..... Conform USP-NF
Density at 25° C ..... 1.017 ÷ 1.020
Refractive index at 20°C ..... 1.531 ÷ 1.534
Sulphated ash ..... ≤50 ppm
Residue solvents ..... Conform USP-NF
Origin (BSE/TSE) ..... Vegetable

Code	Size	Packaging	Notes
529022	1 l	Glass bottle	
529021	2.4 l	Glass bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

2-Phenylethanol > RE - Pure

RE

Description ..... Liquido limpido incolore
Identification ..... Positive
Density at 20° C ..... 1.017 ÷ 1.023
Refractive index at 20° C ..... 1.529 ÷ 1.535
Assay (CPG) ..... ≥ 98 %

Code	Size	Packaging	Notes
308731	1 l	Glass bottle	

Phenylhydrazine hydrochloride

• Fenilidrazina cloridrato • Phénylhydrazine chlorhydrate • Fenilhidracina clorhidrato • Phenylhydrazinhydrochlorid

$C_6H_5NHNH_2.HCl$
Molecular Weight: 144,61
CAS: 59-88-1
EEC-N: 200-444-7

Classification transport
ONU: 2811
Transport Hazard class: 6.1
Packing group III

Danger
H301-H311-H331-H315-H319-H317-H341-H350-H372-H400-HA26
P280-P304+P340-P305+P351+P338-P308+P313-P330-P361+P364-P403+P233

Phenylhydrazine hydrochloride > RPE - For analysis

RPE

Description ..... White to light yellow to pink-beige
Identification ..... Positive
Sulphated ash ..... ≤ 0.2 %
Assay (acidimetric) ..... ≥ 99.0 %

Code	Size	Packaging	Notes
450843	50 g	Plastic bottle	

Phenylhydrazine hydrochloride solution

• Fenilidrazina cloridrato soluzione • Phénylhydrazine chlorhydrate solution • Fenilhidracina clorhidrata solución • Phenylhydrazinhydrochloridlösung

$C_6H_8N_2.HCl$
Molecular Weight: 144,6
CAS: 59-88-1

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Danger
H290-H350-HEU208-HA26
P234-P280-P308+P313-P390-P406-P501a

Phenylhydrazine hydrochloride solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611064501	100 ml	Glass bottle	Ref Ph.Eur 1064501

Phloroglucinol
• Floroglucina • Phloroglucinol • Floroglucina • Phloroglucin

Synonym:
1,3,5-Trihydroxybenzene

$1,3,5\text{-}(OH)_3C_6H_3$
Molecular Weight: 126,11
CAS: 108-73-6
EEC-N: 203-611-2

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Phloroglucinol > RE - Pure

RE

Description ..... Polvere o scaglie bianche
Identification ..... Positive
Melting point ..... 213 - 218 ° C
Water (K.F.) ..... ≤ 2 %
Assay (HPLC) ..... ≥99 %

Code	Size	Packaging	Notes
452031	10 g	Glass bottle	
452033	50 g	Glass bottle	

Phloxin B
• Floxina B • Phloxine B • Floxina B • Phloxin B

Synonym:
• *Acid Red 92*
• *2',4',5',7'-Tetrabromo-4,5,6,7-tetrafluoresceine disodium salt*

$C_{20}H_2Br_4Cl_4Na_2O_5$
Molecular Weight: 829,64
CAS: 18472-87-2
EEC-N: 242-355-6

Phloxin B > RS - For microscopy - C.I. 45410

RS

Description ..... Red brown powder
Identification ..... Positive

Code	Size	Packaging	Notes
452051	10 g	Glass bottle	
452052	25 g	Glass bottle	

Dye for botanical, cytology and histology

Phosphate standard solution
• Fosfati standard soluzione • Phosphate solution standard • Fosfato, solución patrón • Phosphat-Standardlösung

Phosphate standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615002200	100 ml	Plastic bottle	A 5 ppm solution: to dilute according to Ref Ph.Eur 5002200
615004200	1 l	Plastic bottle	A 200 ppm solution Ref Ph.Eur 5004200

Phosphate standard solution > RS - Standard solution for ion chromatography

RS

Code	Size	Packaging	Notes
503341	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503343	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Phosphate buffer pH 9.0
• Tampone fosfato pH 9.0 • Tampon phosphate pH 9.0 • Tampon fosfato pH 9.0 • Pufferlosung phosphat pH 9.0

Phosphate buffer pH 9.0 > RS - For analysis according to Ph. Eur. Chap. 4.1.3

RS

Code	Size	Packaging	Notes
614008300	1 l	Plastic bottle	Ref Ph.Eur 4008300

DASIT GROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

Phosphate buffer pH 7.4
• Tampone fosfato pH 7.4 • Tampon phosphate pH 7.4 • Tampon fosfato pH 7.4 • Pufferlosung phosphat pH 7.4

Phosphate buffer pH 7.4 > RS - For analysis
RS

Temperature of measurement ..... 19 - 21 °C pH..... 7.35 - 7.45 unité pH

Code	Size	Packaging	Notes
524965	5 l	Kubidos	
PS0740/95	5 l	Kubidos	

Composition: Potassium dihydrogen phosphate 1.9g/l, disodium hydrogen phosphate 19.3g/l, deionized water 992.5 g/l. Traceable to NI ST

Phosphate buffer pH 7.4 > RS - For analysis according to Ph. Eur. Chap. 4.1.3
RS

Code	Size	Packaging	Notes
614004800	1 l	Plastic bottle	Ref Ph.Eur 4004800

Phosphate buffer pH 7.2
• Tampone fosfato pH 7.2 • Tampon phosphate pH 7.2 • Tampón fosfato pH 7.2 • Pufferlosung phosphat pH 7.2

Phosphate buffer pH 7.2 > RS - For analysis
RS

pH..... 7.16 ÷ 7.24 Temperature ..... 19 ÷ 21 °C

Code	Size	Packaging	Notes
525925	2.5 l	Glass bottle	
525921	25 l	Plastic tank	

Phosphate buffer pH 6.8
• Tampone fosfato pH 6.8 • Tampon phosphate pH 6.8 • Tampon fosfato pH 6.8 • Pufferlosung phosphat pH 6.8

Phosphate buffer pH 6.8 > RS - For analysis
RS

pH..... 6.75 - 6.85 Temperature ..... 19 ÷ 21 °C

Code	Size	Packaging	Notes
524952	10 l	Plastic tank	

Composition: Potassium dihydrogen phosphate 6.8g/l, sodium hydroxide 0.9g/l, deionized water 992.3 g/l. Traceable to NIST

Phosphate buffer pH 6.8 > RS - For analysis according to Ph. Eur. Chap. 4.1.3
RS

Code	Size	Packaging	Notes
614003400	1 l	Plastic bottle	Ref Ph.Eur 4003400

Phosphate buffer pH 6.0
• Tampone fosfato pH 6.0 • Tampon phosphate pH 6.0 • Tampon fosfato pH 6.0 • Pufferlosung phosphat pH 6.0

Phosphate buffer pH 6.0 > RS - For analysis according to Ph. Eur. Chap. 4.1.3
RS

Code	Size	Packaging	Notes
614002400	1 l	Plastic bottle	Ref Ph.Eur 4002400

Phosphate buffer pH 5.5
• Tampone fosfato pH 5.5 • Tampon phosphate pH 5.5 • Tampon fosfato pH 5.5 • Pufferlosung phosphat pH 5.5

Phosphate buffer pH 5.5 > RS - For analysis according to Ph. Eur. Chap. 4.1.3
RS

Code	Size	Packaging	Notes
614002000	1 l	Plastic bottle	Ref Ph.Eur 4002000

Phosphate buffer pH 3.0

• Tampone fosfato pH 3.0 • Tampon phosphate pH 3.0 • Tampon fosfato pH 3.0 • Pufferlosung phosphat pH 3.0

Phosphate buffer pH 3.0 > RS - For analysis according to Ph. Eur. Chap. 4.1.3

RS

Code	Size	Packaging	Notes
614000501	100 ml	Plastic bottle	Ref Ph.Eur 4000500
614000500	1 l	Plastic bottle	Ref Ph.Eur 4000500

Phosphate buffer pH 2.0

• Tampone fosfato pH 2.0 • Tampon phosphate pH 2.0 • Tampon fosfato pH 2.0 • Pufferlosung phosphat pH 2.0

Phosphate buffer pH 2.0 > RS - For analysis according to Ph. Eur. Chap. 4.1.3

RS

Code	Size	Packaging	Notes
614007900	1 l	Plastic bottle	Ref Ph.Eur 4007900

Phosphomolybdic acid

• Acido fosfomolibdico • Acide phosphomolybdique • Acido fosfomolíbdico • Phosphormolybdänsäure

Synonym:
Molybdophosphoric acid

$2H_3PO_4.20MoO_3.48H2O$
Molecular Weight: 3939,5
CAS: 51429-74-4

Classification transport
ONU: 3084
Transport Hazard class: 8
Packing group II



Danger
H272-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Phosphomolybdic acid > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Yellow crystals	Chloride ..... ≤ 200 ppm	Sulphate ..... ≤ 250 ppm
Identification ..... Positive	Water-insoluble matter ..... ≤ 100 ppm	Ca ..... ≤ 200 ppm
Ammonium ..... ≤ 100 ppm	Heavy metals (Pb) ..... ≤ 50 ppm	Fe ..... ≤ 50 ppm

Code	Size	Packaging	Notes
405913	50 g	Glass bottle	
405915	250 g	Glass bottle	

Phosphomolybdotungstic reagent

• Reattivo fosfomolibdotungstico • Réactif phosphomolybdotungstique • Fosfomolibdotungstenico reactivo • Phosphomolybdotungstica-Reagenz

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III



Warning
H290-H302-H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Phosphomolybdotungstic reagent > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611065000	100 ml	Plastic bottle	Ref Ph.Eur 1065000

Storage: at 2 °C to 8 °C

Product specifications are subject to changes. Please visit our website for updates.

Phosphonic acid

• Acido fosfonico • Acide phosphonique • Acido fosfonico • Phosphorige Säure

Synonym:
Phosphorous acid

H_3PO_3
Molecular Weight: 82
CAS: 13598-36-2
EEC-N: 237-066-7

Classification transport
ONU: 2834
Transport Hazard class: 8
Packing group III

Danger
H302-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Phosphonic acid > RPE - For analysis

RPE

Description ..... White crystals
Identification ..... Positive
Chloride ..... ≤ 100 ppm
Sulphate ..... ≤ 80 ppm
Fe ..... ≤ 20 ppm
Assay (acidimetric) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
406056	500 g	Glass bottle	
406053	20 kg	Plastic bucket	

Phosphoric acid ▶ Orthophosphoric acid 99%

Phosphoric anhydride ▶ Phosphorus pentoxide

Phosphorous acid ▶ Phosphonic acid

Phosphorus standard solution

• Fosforo standard soluzione • Phosphore solution standard • Fósforo, solución patrón • Phosphor-Standardlösung

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Phosphorus standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505762	100 ml	Plastic bottle	conc. 10 ppm Matrix: Water
505765	100 ml	Plastic bottle	conc. 100 ppm Matrix: Water
505763	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Phosphorus standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503791	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503793	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503795	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Water
503797	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Phosphorus pentoxide

• Anidride fosforica • Anhydride phosphorique • Anhídrido fosfórico • Phosphorpentoxid

Synonym:
Phosphoric anhydride

P_2O_5
Molecular Weight: 141,94
CAS: 1314-56-3
EEC-N: 215-236-1

Classification transport
ONU: 1807
Transport Hazard class: 8
Packing group II

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Phosphorus pentoxide > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White powder
Identification ..... Positive
Water-insoluble matter ..... ≤200 ppm
Ammonium ..... ≤100 ppm
Heavy metals (Pb) ..... ≤100 ppm
Assay (acidimetric) ..... ≥98.0 %

Code	Size	Packaging	Notes
421808	100 g	Glass bottle	
421802	250 g	Glass bottle	

Phosphorus pentoxide > RE - Pure

RE

Description ..... White powder
Identification ..... Positive
Heavy metals (Pb) ..... ≤500 ppm
Assay (acidimetric) ..... ≥98 %

Code	Size	Packaging	Notes
317753	250 g	Glass bottle	

Phosphosulfuric acid

• Acido fosfosolforico • Acide phosphosulfurique • Acido fosfosulfúrico • Phosphosulfonsäure

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Phosphosulfuric acid > RS - For nitrogen detection according to Kjeldahl

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Ammonium ..... ≤1 ppm
Nitrate ..... ≤0.2 ppm

Code	Size	Packaging	Notes
E406101	1 l	Glass bottle	

Phosphotungstic acid

• Acido fosfotungstico • Acide phosphotungstique • Acido fosfotúngstico • Phosphorwolframsäure

Synonym:
Tungstophosphoric acid

$H_3PO_4.12WO_3.H_2O$
Molecular Weight: 2880,05
CAS: 12501-23-4

Classification transport
ONU: 3260
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Phosphotungstic acid > RPE - For analysis

RPE

Description ..... White powder
Identification ..... Positive
Chloride ..... ≤ 20 ppm
Sulphate ..... ≤ 50 ppm
Ammonium ..... ≤ 50 ppm
Residue on calcination ..... ≤ 17 %
Heavy metals (Pb) ..... ≤ 40 ppm
Fe ..... ≤ 30 ppm
Na ..... ≤ 100 ppm
Assay ..... ≥ 82 %

Code	Size	Packaging	Notes
406154	100 g	Glass bottle	

Phosphotungstic acid solution

• Acido fosfotungstico soluzione • Acide phosphotungstique solution • Acido fosfotúngstico solución
• Phosphorwolframsäure-Lösung

Synonym:
• *Phosphotungstic acid hydrate*
• *Tungstophosphoric acid*

Phosphotungstic acid solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611065200	100 ml	Plastic bottle	Ref Ph.Eur 1065200

Product specifications are subject to changes. Please visit our website for updates.

o-Phthalaldehyde

• o-Ftalaldeide • o-Phtalaldéhyde • o-Ftaldialdehído • o-Phthalaldehyd

Synonym: *Benzene-1,2-dicarboxaldehyde*

$C_6H_4(CHO)_2$
Molecular Weight: 134,13
CAS: 643-79-8
EEC-N: 211-402-2

Danger
H301
P264-P270-P301+P310a-P330-P405-P501a

o-Phthalaldehyde > RPE - For analysis

RPE

Description ..... Yellow-orange crystals
Identification ..... Positive
Melting point ..... 55 - 58 ° C
Acidity (Phtalic acid) ..... ≤ 0.3 %
Assay (GLC) ..... ≥ 98.0 %

Code	Size	Packaging	Notes
452751	10 g	Glass bottle	

Phthalic acid

• Acido ftalico • Acide phtalique • Acido ftálico • Phthalsäure

Synonym: *1,2-Benzenedicarboxylic acid*

$1,2\text{-}(COOH)_2C_6H_4$
Molecular Weight: 166,13
CAS: 88-99-3
EEC-N: 201-873-2

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Phthalic acid > RPE - For analysis - Reag. Ph. Eur.

RPE

Description ..... White crystalline powder
Identification ..... Positive
Melting point ..... ~ 210 °C
Assay (acidimetric) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
406205	250 g	Plastic bottle	

Phthalic anhydride

• Anidride ftalica • Anhydride phtalique • Anhídrido ftálico • Phthalsäureanhydrid

$C_6H_4(CO)_2O$
Molecular Weight: 148,12
CAS: 85-44-9
EEC-N: 201-607-5

Classification transport
ONU: 2214
Transport Hazard class: 8
Packing group III

Danger
H302-H315-H318-H334-H317-H335
P284-P304+P340-P310a-P305+P351+P338-P330-P362+P364-P342+P311a-P403+P233

Phthalic anhydride > RE - Pure

RE

Description ..... White flakes
Identification ..... Positive
Melting point ..... 129 ÷ 132 °C
Assay ..... ≥ 98.5 %

Code	Size	Packaging	Notes
318007	1 kg	Plastic bottle	

3-Picolinic acid ▶ Nicotinic acid

Picric acid solution

• Acido picrico soluzione • Acide picrique solution • Acido pícrico solución • Pikrinsäure

Synonym: *2,4,6-Trinitrophenol*

$C_6H_3N_3O_7$
Molecular Weight: 229,11
CAS: 88-89-1

HEU210

Picric acid solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611065801	100 ml	Plastic bottle	Ref Ph.Eur 1065801
611065802	100 ml	Plastic bottle	Picric acid solution R1 Ref Ph.Eur 1065802

Picric acid solution > RPE - For analysis

RPE

Description ..... Yellow clear liquid Identification ..... Positive Density at 20° C ..... 1.00 ÷ 1.02 Assay ..... 1.1 ÷ 1.3 %

Code	Size	Packaging	Notes
409302	500 ml	Plastic bottle	
409305	2.5 l	Plastic bottle	

Saturated aqueous solution ~ 1.2%

Piperidine

• Piperidina • Pipéridine • Piperidina • Piperidin

Synonym:
Hexahydropyridine

$NH(CH_2)_4CH_2$
Molecular Weight: 85,15
CAS: 110-89-4
EEC-N: 203-813-0

Classification transport
ONU: 2401
Transport Hazard class: 8
Packing group I

Danger
H225-H302-H311-H331-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P361+P364-P403+P233

Piperidine > RS - For peptide synthesis

RS

Refractive index at 20°C ..... 1.45 - 1.454 Water content (K.F.) ..... ≤ 3000 mg/Kg Colour ..... ≤ 10 Hazen Assay (GC) ..... ≥ 99 %

Code	Size	Packaging	Notes
P0663518	500 ml	Glass bottle	
P0663516	1 l	Glass bottle	
P0663521	2.5 l	Glass bottle	

Piperidine > RPE - For analysis

RPE

Description ..... Clear liquid Density at 20° C ..... 0.856 ÷ 0.866 Water (K.F.) ..... <0.3 %
Colour ..... < 50 APHA Assay (GLC) ..... ≥99.0 %

Code	Size	Packaging	Notes
469551	100 ml	Glass bottle	
469552	500 ml	Glass bottle	

Piperidine > RE - Pure

RE

Refractive index at 20°C ..... 1.450 - 1.454 Density d20/4 ..... 0.856 - 0.866 Colour ..... ≤ 20 Hazen
Identification (IR) ..... Conform Water content (K.F.) ..... ≤ 5000 mg/Kg Assay (GC) ..... ≥ 99 %

Code	Size	Packaging	Notes
P0660216	1 l	Glass bottle	
P0660221	2.5 l	Glass bottle	
P0660229	5 l	Plastic tank	
P0660248	25 l	Metal drum	

Platinum standard solution

• Platino standard soluzione • Platine solution standard • Platino, solución patrón • Platin-Standardlösung

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group III

Platinum standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505787	100 ml	Plastic bottle	conc. 10 ppm Matrix: Hydrochloric acid
505788	100 ml	Plastic bottle	conc. 100 ppm Matrix: Hydrochloric acid
505789	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Product specifications are subject to changes. Please visit our website for updates.

Platinum standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503831	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503833	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503835	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid
503837	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Polyvinylpyrrolidone

• Polivinilpirrolidone • Polyvinylpyrrolidone • Polivinilpirrolidona • Polyvinylpyrrolidon

Synonym:
• *PVP*
• *Polyvidone*

$(C_6H_9NO)n$
Molecular Weight: 25000-30000
CAS: 9003-39-8

Polyvinylpyrrolidone > RPE - For analysis

RPE

Description ..... White powder Identification ..... Positive Loss on drying .....≤5 %

Code	Size	Packaging	Notes
470071	500 g	Plastic bottle	
470072	1.5 kg	Plastic bottle	

Ponceau red BS

• Rosso Ponceau BS • Rouge Ponceau BS • Rojo Ponceau BS • Ponceau rot BS

Synonym:
• *Ponceau BS*
• *Acid Red 66*

$C_{22}H_{14}N_4Na_2O_7S_2$
Molecular Weight: 556,48
CAS: 4196-99-0
EEC-N: 224-084-5

Ponceau red BS > RS - For microscopy - C.I. 26905

RS

Description ..... Red brick powder Identification ..... Positive

Code	Size	Packaging	Notes
476941	10 g	Glass bottle	

Dye for histology

Ponceau red S

• Rosso Ponceau S • Rouge Ponceau S • Rojo Ponceau S • Ponceau rot S

Synonym:
3-Hydroxy-4-(2-sulfo-4-[4-sulfophenylazo] phenylazo)-2,7-naphthalenedisulfonic acid sodium salt

$C_{22}H_{12}N_4Na_4O_{13}S_4$
Molecular Weight: 760,56
CAS: 6226-79-5
EEC-N: 228-319-2

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Ponceau red S > RS - For microscopy - C.I. 27195

RS

Description ..... Brown powder Identification ..... Positive

Code	Size	Packaging	Notes
476981	5 g	Glass bottle	
476982	25 g	Glass bottle	

Dye for histochemistry

Potassium standard solution

• Potassio standard soluzione • Potassium solution standard • Potasio, solución patrón • Kalium-Standardlösung

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Potassium standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615002401	100 ml	Plastic bottle	A 20 ppm solution: to dilute according to Ref Ph.Eur 5002401
615002402	100 ml	Plastic bottle	A 0.1 % solution Ref Ph.Eur 5002402
615002409	100 ml	Plastic bottle	A 100 ppm solution: to dilute according to Ref Ph.Eur 5002400
615005100	100 ml	Plastic bottle	A 600 ppm solution: to dilute according to Ref Ph.Eur 5005100

Potassium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505682	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505685	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505683	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Potassium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503671	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503673	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503675	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503677	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Potassium standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
507753	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
506960	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
E497605	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
E497601	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Potassium standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
470081		Plastic ampoule	conc. 1.000 ppm Matrix: Nitric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Product specifications are subject to changes. Please visit our website for updates.

Potassium standard solution > RS - Standard solution for ion chromatography

RS

Code	Size	Packaging	Notes
503271	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503273	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Potassium acetate

• Potassio acetato • Potassium acétate • Potasio acetato • Kaliumacetat

CH_3COOK
Molecular Weight: 98,15
CAS: 127-08-2
EEC-N: 204-822-2

Potassium acetate > RPE - For analysis

RPE

Description	White granular powder	Sulphate	≤ 50 ppm	Assay (non-aqueos medium)	≥ 99 %	Mg	≤ 100 ppm
Identification	Positive	Heavy metals (Pb)	≤ 10 ppm	Water	≤ 1.0 %	Hg	≤ 1 ppm
pH sol. 5% at 20°C	7.5 ÷ 8.5	Fe	≤ 20 ppm	Ca	≤ 100 ppm		
Chloride	≤ 50 ppm	Zn	≤ 20 ppm	Na	≤ 0.4 %		

Code	Size	Packaging	Notes
470145	100 g	Plastic bottle	
470146	500 g	Plastic bottle	
470147	1 kg	Plastic bottle	
470143	25 kg	Plastic bucket	

Potassium acetate > ERBApharm - According to pharmacopoeia: BP-Ph.Eur.

ERBApharm

Description	White crystalline powder	pH solution 5%	7.5 ÷ 9.0	Sulphate	≤200 ppm	Assay (non-aqueos medium)	99.0 ÷ 101.0 % s.s.
Identification	Positive	Loss on drying	≤3.0 %	Al	≤1 ppm	Origin (BSE/TSE)	Synthesis
Appearance of solution	Conform Ph.Eur.	Chloride	≤200 ppm	Fe	≤20 ppm		
Reducing substances	Conform Ph.Eur.	Heavy metals (Pb)	≤4 ppm	Na	≤0.5 %		

Code	Size	Packaging	Notes
358907	1 kg	Plastic bottle	
358908	5 kg	Plastic tank	
358903	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Potassium aluminum sulfate dodecahydrate ▶ Aluminum potassium sulfate dodecahydrate

Potassium antimonyl tartrate trihydrate ▶ Antimony potassium tartrate

Potassium bicarbonate

• Potassio bicarbonato • Potassium bicarbonate • Potasio bicarbonato • Kaliumbicarbonat

Synonym:
Potassium hydrogen carbonate

$KHCO_3$
Molecular Weight: 100,12
CAS: 298-14-6
EEC-N: 206-059-0

Potassium bicarbonate > RPE - For analysis - ACS

RPE

Description	White crystalline powder	Chloride	≤10 ppm	Ca	≤20 ppm	Assay (alkalimetric)	99.7 ÷ 100.5 % s.s.
Identification	Positive	Phosphate	≤5 ppm	Fe	≤5 ppm		
Water-insoluble matter	≤100 ppm	Total sulphur	≤30 ppm	Mg	≤10 ppm		
Ammonium	≤5 ppm	Heavy metals (Pb)	≤5 ppm	Na	≤300 ppm		

Code	Size	Packaging	Notes
470285	100 g	Plastic bottle	
470286	500 g	Plastic bottle	
470287	1 kg	Plastic bottle	
470289	5 kg	Plastic jar	

Potassium bisulfate

• Potassio bisolfato • Potassium bisulfate • Potasio bisulfato • Kaliumbisulfat

Synonym:
Potassium hydrogen sulfate

$KHSO_4$
Molecular Weight: 136,17
CAS: 7646-93-7
EEC-N: 231-594-1

Classification transport
ONU: 2509
Transport Hazard class: 8
Packing group II



Danger
H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Potassium bisulfate > RPE - For analysis

RPE

Description	White crystals	Ammonium	≤ 20 ppm	Heavy metals (Pb)	≤ 20 ppm	Fe	≤ 20 ppm
Identification	Positive	Chloride	≤ 20 ppm	Ca	≤ 200 ppm	Assay (acidimetric)	98 ÷ 102 %

Code	Size	Packaging	Notes
470556	100 g	Plastic bottle	
470557	1 kg	Plastic bottle	
470552	25 kg	Plastic bucket	

Potassium bitartrate ▶ Potassium L-tartrate monobasic

Potassium bromate

• Potassio bromato • Potassium bromate • Potasio bromato • Kaliumbromat

$KBrO_3$
Molecular Weight: 167,01
CAS: 7758-01-2
EEC-N: 231-829-8

Classification transport
ONU: 1484
Transport Hazard class: 5.1
Packing group II



Danger
H271-H301-H350-HA26
P210-P280-P283-P301+P310a-P306+P360-P308+P313

Potassium bromate > RS - For analysis according to Ph. Eur. Chap. 4.2.1

RS

Code	Size	Packaging	Notes
612000300	50 g	Plastic bottle	Ref Ph.Eur 2000300

Potassium bromate > RPE - For analysis - ISO - ACS - Reag.Ph.Eur. - Reag.USP

RPE

Description	White crystals	Bromide	Conform ACS	Sulphate	≤ 50 ppm	Assay (oxidimetric)	≥ 99.8 %
Identification	Positive	Water-insoluble matter	≤ 50 ppm	Fe	≤ 20 ppm		
pH sol. 5% at 25° C	5.0 ÷ 9.0	Heavy metals (Pb)	≤ 5 ppm	Na	≤ 100 ppm		

Code	Size	Packaging	Notes
470655	250 g	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Potassium bromate 0.033 mol/l (0.198N)

• Potassio bromato 0.033 mol/l (0.198N) • Potassium bromate 0.033 mol/l (0.198N) • Potasio bromato 0.033 mol/l (0.198N)
• Kaliumbromat 0.033 mol / l (0.198 N)



Danger
H350-HA26
P201-P202-P280-P308+P313-P405-P501a

Potassium bromate 0.033 mol/l (0.198N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2 — RS

Code	Size	Packaging	Notes
613004200	1 l	Plastic bottle	Ref Ph.Eur 3004200

Potassium bromate 0.02 mol/l (0.12N)

• Potassio bromato 0.02 mol/l (0.12N) • Potassium bromate 0.02 mol/l (0.12N) • Potasio bromato 0.02 mol/l (0.12N) • Kaliumbromat 0.02 mol/l (0.12N)

Potassium bromate 0.02 mol/l (0.12N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2 — RS

Code	Size	Packaging	Notes
613004300	1 l	Plastic bottle	Ref Ph.Eur 3004300

Potassium bromate 0.0167 mol/l (0.1N)

• Potassio bromato 0.0167 mol/l (0.1N) • Potassium bromate 0.0167 mol/l (0.1N) • Potasio bromato 0.0167 mol/l (0.1N) • Kaliumbromat 0.0167 mol/l (0.1N)

$KBrO_3$
CAS: 7758-01-2



Danger
H350-HA26
P201-P202-P280-P308+P313-P405-P501a

Potassium bromate 0.0167 mol/l (0.1N) > RPE - NORMEX - For analysis — RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
470681		Glass ampoule	Volume: 60 ml

Volumetric concentrated solution to prepare 1 L of solution 0,1 N

Potassium bromide

• Potassio bromuro • Potassium bromure • Potasio bromuro • Kaliumbromid

KBr
Molecular Weight: 119,01
CAS: 7758-02-3
EEC-N: 231-830-3



Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Potassium bromide > RS - For optical spectroscopy — RS

Description ..... White cryst. powder Identification (I.R.)..... Conform

Code	Size	Packaging	Notes
470701	100 g	Glass bottle	

Potassium bromide > RPE - For analysis - ACS — RPE

Description ..... White crystals	Bromate ..... ≤ 10 ppm	Sulphate ..... ≤ 50 ppm	Fe ..... ≤ 5 ppm
Identification ..... Positive	Chloride ..... ≤ 0.2 %	Heavy metals (Pb)..... ≤ 5 ppm	Mg ..... ≤ 10 ppm
pH sol. 5% at 25° C ..... 5.0 ÷ 8.8	Iodate ..... ≤ 10 ppm	Ba ..... ≤ 20 ppm	Na ..... ≤ 0.02 %
Water-insoluble matter ..... ≤ 50 ppm	Iodide ..... ≤ 10 ppm	Ca ..... ≤ 20 ppm	Assay (argentimetric) ..... ≥ 99.0 %

Code	Size	Packaging	Notes
470734	100 g	Plastic bottle	
470735	250 g	Plastic bottle	
470737	1 kg	Plastic bottle	
470733	25 kg	Plastic bucket	

Potassium bromide > ERBApharm - According to pharmacopoeia: BP-DAB-Ph.Eur.-Ph.Franc.

ERBApharm

Description ..... White crystalline powder	Bromate .....Conform Ph.Eur.	Heavy metals (Pb).....≤ 10 ppm	Assay (argentimetric) ..98.0 ÷ 100.5 % s.s.
Identification ..... Positive	Iodide .....Conform Ph.Eur.	Mg,alkal.earth met.(Ca) .....≤ 200 ppm	
Appearance of solution .....Conform Ph.Eur.	Loss on drying ..... ≤ 1.0 %	Sulphate .....≤ 100 ppm	
Acidity or alkalinity.....Conform Ph.Eur.	Chloride..... ≤ 0.6 %	Fe .....≤ 20 ppm	

Code	Size	Packaging	Notes
359707	1 kg	Plastic bottle	
359702	5 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Potassium carbonate

• Potassio carbonato • Potassium carbonate • Potasio carbonato • Kaliumcarbonat

K_2CO_3
Molecular Weight: 138,21
CAS: 584-08-7

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Potassium carbonate > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611068900	100 g	Glass bottle	Ref Ph.Eur 1068900

Potassium carbonate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White powder	Phosphate .....≤ 10 ppm	Ca.....≤ 50 ppm	Assay (alkalimetric)..... ≥ 99.0 %
Identification ..... Positive	Silicate .....≤ 50 ppm	Fe .....≤ 5 ppm	
Water-insoluble matter .....≤ 100 ppm	Total sulphur .....≤ 40 ppm	Mg .....≤ 20 ppm	
Chloride.....≤ 30 ppm	Heavy metals (Pb).....≤ 5 ppm	Na.....≤ 200 ppm	

Code	Size	Packaging	Notes
470805	250 g	Plastic bottle	
470807	1 kg	Plastic bottle	
470801	5 kg	Plastic tank	

Potassium carbonate > RE - Pure

RE

Description ..... White crystalline powder	Sulphate .....≤ 50 ppm	Na..... ≤ 0.25 %	KOH ..... ≤ 0.15 %
Identification ..... Positive	Heavy metals (Pb).....≤ 1 ppm	Assay (alkalimetric).....99.0 ÷ 100.0 %	
Chloride.....≤ 20 ppm	Fe .....≤ 5 ppm	Loss on drying ..... ≤ 0.8 %	

Code	Size	Packaging	Notes
359808	1 kg	Plastic bottle	
359809	5 kg	Plastic tank	
359803	25 kg	Plastic bucket	

Potassium chloride

• Potassio cloruro • Potassium chlorure • Potasio cloruro • Kaliumchlorid

KCl
Molecular Weight: 74,55
CAS: 7447-40-7
EEC-N: 231-211-8

Potassium chloride > RS - For soils analysis

RS

Assay (argentimetric)..... ≥ 99.0 %	Water insoluble substances..... ≤ 0.005 %	Ammonium (NH4)..... ≤ 0.00007 %	Fe ..... ≤ 0.0002 %
pH sol. 5% at 25°C .....5.4 ÷ 8.6	Sulphate..... ≤ 0.002 %	Phosphate ..... ≤ 0.001 %	Heavy metals (Pb)..... ≤ 0.0002 %

Code	Size	Packaging	Notes
471181	5 kg	Plastic jar	

DASIT GROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

Potassium chloride > RPE - For analysis - ISO - ACS - Reag.Ph.Eur. - Reag.USP

RPE

Description .....White crystals
Identification ..... Positive
Acidity or alkalinity.....Conform
Appearance of solution .....Conform
pH sol. 5% at 25° C .....5.4 ÷ 8.6
Loss on drying ..... ≤ 1.0 %
Water-insoluble matter .....≤ 50 ppm
Mg and alkaline-earth metals (Ca).... ≤ 200 ppm
Bromide.....≤ 100 ppm
Nitrate,Chlorate (NO3) .....≤ 30 ppm
Phosphate .....≤ 5 ppm
Iodide .....≤ 20 ppm
Sulphate.....≤ 10 ppm
Heavy metals (Pb).....≤ 5 ppm
Al .....≤ 1 ppm
Ba .....≤ 10 ppm
Ca .....≤ 20 ppm
Fe .....≤ 3 ppm
Mg .....≤ 10 ppm
Na.....≤ 50 ppm
Assay (argentimetric) ..99.0 ÷ 100.5 % s.s.

Code	Size	Packaging	Notes
471175	100 g	Plastic bottle	
471176	500 g	Plastic bottle	
471177	1 kg	Plastic bottle	
471171	5 kg	Plastic jar	
471173	25 kg	Plastic bucket	

Potassium chloride > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP-DAB

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution ..... Conform Ph.Eur
Acidity or alkalinity.....Conform Ph.Eur.
Iodide .....≤ 50 ppm
Barium .....Conform Ph.Eur.
Na..... Conform USP-NF
Calcium + Magnesium .... Conform USP-NF
Residue solvents ..... Conform USP-NF
Loss on drying ..... ≤ 1.0 %
Bromide..... ≤ 0.1 %
Heavy metals (Pb).....≤ 10 ppm
Mg,alkal.earth met.(Ca) .....≤ 200 ppm
Sulphate .....≤ 300 ppm
Fe .....≤ 20 ppm
Assay (argentimetric).....99.0 ÷ 100.5 %
Origin (BSE/TSE)..... Synthesis
Residual solvents (Current ICH).....Conform

Code	Size	Packaging	Notes
360107	1 kg	Plastic bottle	
360109	5 kg	Plastic tank	
360106	25 kg	Plastic bucket	
360104	50 kg	Fibre drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Potassium chloride 3.5 mol/l (3.5N)

• Potassio cloruro 3.5 mol/l (3.5N) • Potassium chlorure 3.5 mol/l (3.5N) • Potasio cloruro 3.5 mol/l (3.5N) • Kaliumchlorid 3.5 mol/l (3.5 N)

KCl
Molecular Weight: 74,55
CAS: 7447-40-7

Potassium chloride 3.5 mol/l (3.5N) > RPE - For analysis

RPE

Code	Size	Packaging	Notes
471225	250 ml	Plastic bottle	

Content is guaranteed for standardized volumes at 20 °C

Potassium chloride 3.5 mol/l (3.5N) + silver chloride

• Potassio cloruro 3.5 mol/l (3.5N) e argento cloruro • Potassium chlorure 3.5 mol/l (3.5N) avec argent chlorure
• Potasio cloruro 3.5 mol/l (3.5N) con plata cloruro • Kaliumchlorid 3.5 mol/l (3.5N) mit Silberchlorid

KCl
Molecular Weight: 74,55
CAS: 7447-40-7

Potassium chloride 3.5 mol/l (3.5N) + silver chloride > RPE - For analysis

RPE

Code	Size	Packaging	Notes
471245	250 ml	Plastic bottle	

Electrolytic solution filling. Content is guaranteed for standardized volumes at 20 °C

Potassium chloride 3 mol/l (3N)

• Potassio cloruro 3 mol/l (3N) • Potassium chlorure 3 mol/l (3N) • Potasio cloruro 3 mol/l (3N) • Kaliumchlorid 3 mol/l (3N)

KCl
Molecular Weight: 74,55
CAS: 7447-40-7

Potassium chloride 3 mol/l (3N) > RPE - For analysis

RPE

Code	Size	Packaging	Notes
471215	250 ml	Plastic bottle	

Electrolyte for the reference electrode.
Content is guaranteed for standardized volumes at 20 °C.

Potassium chloride 3 mol/l (3N) water-glycerol solution

• Potassio cloruro 3 mol/l (3N) soluzione acqua-glicerina • Potassium chlorure 3 mol/l (3N) solution eau/glycérine
• Potasio cloruro 3 mol/l (3N) solución agua-glicerina • Kaliumchlorid 3 mol/l (3N) Lösung aus Wasser und Glycerin

KCl
Molecular Weight: 74,55
CAS: 7447-40-7

Potassium chloride 3 mol/l (3N) water-glycerol solution > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titolo (KCl) ..... 1.485 ÷ 1.515 Density at 20°C ..... 1.190 ÷ 1.210

Code	Size	Packaging	Notes
471275	250 ml	Plastic bottle	

Electrolytic solution filling. Content is guaranteed for standardized volumes at 20 °C

Potassium chloride 3 mol/l (3N) + silver chloride

• Potassio cloruro 3 mol/l (3N) e argento cloruro • Potassium chlorure 3 mol/l (3N) avec argent chlorure • Potasio cloruro 3 mol/l (3N) con plata cloruro
• Kaliumchlorid 3 mol/l (3N) mit Silberchlorid

KCl
Molecular Weight: 74,55
CAS: 7447-40-7

Potassium chloride 3 mol/l (3N) + silver chloride > RPE - For analysis

RPE

Code	Size	Packaging	Notes
471235	250 ml	Plastic bottle	

Electrolyte for the reference electrode.
Content is guaranteed for standardized volumes at 20 °C

Potassium chloride 1 mol/l (1N)

• Potassio cloruro 1 mol/l (1N) • Potassium chlorure 1 mol/l (1N) • Potasio cloruro 1 mol/l (1N) • Kaliumchlorid 1 mol/l (1N)

KCl
Molecular Weight: 74,55
CAS: 7447-40-7

Potassium chloride 1 mol/l (1N) > RS - For soils analysis

RS

Assay (potentiometry) ..... 0.98 - 1.02 N

Code	Size	Packaging	Notes
PS0772/79	100 l	Plastic drum	

Product specifications are subject to changes. Please visit our website for updates.

Potassium chloride 0.1 mol/l (0.1N)

• Potassio cloruro 0.1 mol/l (0.1N) • Potassium chlorure 0.1 mol/l (0.1N) • Potasio cloruro 0.1 mol/l (0.1N) • Kaliumchlorid 0.1 mol/l (0.1N)

KCl
Molecular Weight: 74,55
CAS: 7447-40-7

Potassium chloride 0.1 mol/l (0.1N) > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611069101	1 l	Plastic bottle	Ref Ph.Eur 1069101

Potassium chloride 0.01 mol/l (0.01N)

• Potassio cloruro 0.01 mol/l (0.01N) • Potassium chlorure 0.01 mol/l (0.01N) • Potasio cloruro 0.01 mol/l (0.01N) • Kaliumchlorid 0.01 mol/l (0.01N)

KCl
Molecular Weight: 74,55
CAS: 7447-40-7

Potassium chloride 0.01 mol/l (0.01N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
471255	1 l	Plastic bottle	Potassium chloride 0.01 M

Electrolytic solution filling

Potassium chloride 25g/l in HCl

• Potassio cloruro soluzione 25g/l in acido cloridrico • Potassium chlorure 25g/l dans HCl • Potasio cloruro solución 25g/l en HCl • Kaliumchlorid 25g/l (HCl)

KCl
Molecular Weight: 74,55
CAS: 7447-40-7

Potassium chloride 25g/l in HCl > RS - Ionisation standard solution for AAS

RS

Code	Size	Packaging	Notes
504538	500 ml	Plastic bottle	Matrix: 2% Hydrochloric acid

Potassium chloride 12g/l

• Potassio cloruro soluzione 12g/l • Potassium chlorure 12g/l • Potasio cloruro solución 12g/l • Kaliumchlorid 12g/l

KCl
Molecular Weight: 74,55
CAS: 7447-40-7

Potassium chloride 12g/l > RS - For analysis according to Ph. Eur.Chap. 2.2.25

RS

Code	Size	Packaging	Notes
506432	10 ml	Sealed cuvette	
506433	100 ml	Glass bottle	

Potassium chloride saturated solution

• Potassio cloruro soluzione satura • Potassium chlorure solution saturée • Potasio cloruro solución saturada • Kaliumchlorid-Lösung gesättigt

KCl
Molecular Weight: 74,55
CAS: 7447-40-7

Potassium chloride saturated solution > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Density at 20° C ..... 1.169 ÷ 1.173

Code	Size	Packaging	Notes
471265	250 ml	Plastic bottle	

Electrolytic solution filling

Potassium chloride solution

• Potassio cloruro soluzione • Potassium chlorure solution • Potasio cloruro solución • Kaliumchlorid-Lösung

KCl
Molecular Weight: 74,55
CAS: 7447-40-7

Potassium chloride solution > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Density at 20° C ..... 1.12 ÷ 1.14 pH at 20° C ..... 6.6 ÷ 6.8

Code	Size	Packaging	Notes
471285	250 ml	Plastic bottle	

Conservation solution for electrodes

Potassium chromate

• Potassio cromato • Potassium chromate • Potasio cromato • Kaliumchromat

K_2CrO_4
Molecular Weight: 194,2
CAS: 7789-00-6
EEC-N: 232-140-5

Classification transport
ONU: 3288
Transport Hazard class: 6.1
Packing group III



Danger
H315-H319-H317-H340-H350i-H335-H410-HA26
P261-P271-P280-P304+P340-P305+P351+P338-P403+P233

Potassium chromate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Yellow crystals
Identification ..... Positive
pH sol. 5% at 25° C ..... 8.6 ÷ 9.8
Water-insoluble matter ..... ≤50 ppm
Chloride ..... ≤50 ppm
Sulphate ..... ≤300 ppm
Ca ..... ≤50 ppm
Na ..... ≤200 ppm
Assay (oxidimetric) ..... ≥99.0 %

Code	Size	Packaging	Notes
471295	250 g	Plastic bottle	
471297	1 kg	Plastic bottle	

Potassium chromate solution 10%

• Potassio cromato soluzione 10% • Potassium chromate 10% solution • Potasio cromato solución 10% • Kaliumchromatlösung bei 10%

K_2CrO_4
Molecular Weight: 194,2
CAS: 7789-00-6

Classification transport
ONU: 3287
Transport Hazard class: 6.1
Packing group II



Danger
H315-H319-H317-H340-H350i-H411-HA26
P261-P280-P305+P351+P338-P308+P313-P362+P364-P337+P313

Potassium chromate solution 10% > RPE - For analysis

RPE

Code	Size	Packaging	Notes
505032	1 l	Bottle	

Product specifications are subject to changes. Please visit our website for updates.

Potassium chromate 5% solution

• Potassio cromato soluzione 5% • Potassium chromate solution 5% • Potasio cromato solución 5% • Kaliumchromatlösung bei 5%

K_2CrO_4
Molecular Weight: 194,2
CAS: 7789-00-6



Danger
H317-H340-H350i-H411-HA26
P261-P280-P308+P313-P362+P364-P333+P313-P501a

Potassium chromate 5% solution > RS - For agroalimentary analysis

RS

Appearance ..... Conform Assay ..... 4.75 ÷ 5.25 %

Code	Size	Packaging	Notes
502681	1 l	Plastic bottle	

Composition according to NF V04-314: K2CrO4 50g water QSP 1 L

Potassium chromate 5% solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611069201	1 l	Plastic bottle	Ref Ph.Eur 1069201

Potassium chromium (III)sulfate dodecahydrate ▶ Chromium (III) potassium sulfate dodecahydrate

Potassium citrate tribasic monohydrate

• Potassio citrato tribasico monoidrato • Potassium citrate tribasique monohydraté
• Potasio citrato tribásico monohidrato • Tribasisches Kaliumcitrat-Monohydrat

Synonym:
Tripotassium citrate

$K_3C_6H_5O_7.H_2O$
Molecular Weight: 324,42
CAS: 6100-05-6
EEC-N: 212-755-5

Potassium citrate tribasic monohydrate > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Reducing substances ..... Conform
Ready carbonizable substances..... Conform
pH sol. 5% at 25° C ..... 7.5 ÷ 9.5
Chloride..... ≤10 ppm
Water-insoluble matter ..... ≤30 ppm
Heavy metals (Pb)..... ≤20 ppm
Oxalate..... ≤100 ppm
Sulfate..... ≤ 50 ppm
As ..... ≤0.4 ppm
Ca..... ≤50 ppm
Cu..... ≤5 ppm
Fe ..... ≤5 ppm
Na..... ≤1500 ppm
Ni..... ≤5 ppm
Zn ..... ≤2 ppm
Assay (non-aqueos medium) ..... ≥99.5 %

Code	Size	Packaging	Notes
471025	250 g	Plastic bottle	
471027	1 kg	Plastic bottle	

Potassium citrate tribasic monohydrate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-Ph. Franc.-BP

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution ..... Conform Ph.Eur.
Acidity or alkalinity..... Conform Ph.Eur.
Ready carbonizable substances..... Conform
Ph.Eur.
Tartrate..... Conform USP-NF
Organic volatile impurities Conform USP-NF
Water (K.F.) ..... 4.0 ÷ 7.0 %
Loss on drying ..... 3.0 ÷ 6.0 %
Chloride..... ≤50 ppm
Heavy metals (Pb)..... ≤10 ppm
Oxalate..... ≤300 ppm
Sulphate..... ≤150 ppm
Na..... ≤0.3 %
Assay (non-aqueos medium) . 99.0 ÷ 101.0 % s.s.

Code	Size	Packaging	Notes
359956	500 g	Plastic bottle	
359957	1 kg	Plastic bottle	
359958	2.5 kg	Plastic bottle	
359959	5 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Potassium dichromate

• Potassio bicromato • Potassium dichromate • Potasio dicromato • Kaliumdichromat

$K_2Cr_2O_7$
Molecular Weight: 294,19
CAS: 7778-50-9
EEC-N: 231-906-6

Classification transport
ONU: 3086
Transport Hazard class: 6.1
Packing group I



Danger
H272-H301-H312-H330-H314-H334-H317-H340-H350-H360FD-H372-H410-HA26
P210-P280-P284-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338-P342+P311a-P403+P233

Potassium dichromate > RPE - For analysis — RPE

Description ..... Orange crystals	Chloride ..... ≤ 10 ppm	Fe ..... ≤ 10 ppm
Identification ..... Positive	Sulphate ..... ≤ 50 ppm	Na ..... ≤ 200 ppm
Loss on drying at 105°C ..... ≤ 0.05 %	Ca ..... ≤ 30 ppm	Assay (oxydimetric) ..... ≥ 99.0 %

Code	Size	Packaging	Notes
470336	500 g	Plastic bottle	
470337	1 kg	Plastic bottle	

Potassium dichromate - Sulfuric acid solution

• Potassio bicromato - Soluzione di acido solforico • Potassium dichromate - Solution dans l'acide sulfurique • Potasio dicromato - Solución de ácido sulfúrico • Kaliumdichromat - Lösung in Schwefelsäure

HEU203

Potassium dichromate - Sulfuric acid solution > RS - For analysis according to Ph. Eur.Chap. 2.2.25 — RS

Code	Size	Packaging	Notes
506442	2 x 10 ml	Sealed cuvette	conc. 60 mg/l
506452	2 x 10 ml	Sealed cuvette	conc. 600 mg/l
506443	100 ml	Glass bottle	conc. 60 mg/l
506453	100 ml	Glass bottle	conc. 600 mg/l

Potassium dichromate solution 0.5%

• Potassio bicromato soluzione 0.5% • Potassium dichromate solution 0.5% • Potasio dicromato solución 0.5% • Kaliumdichromatlösung 0.5%

$K_2Cr_2O_7$
Molecular Weight: 294,19
CAS: 7778-50-9
EEC-N: 231-906-6



Danger
H340-H350-H360FD-H412-HEU203-HEU208-HA26
P201-P273-P280-P308+P313-P405-P501a

Potassium dichromate solution 0.5% > RS - For analysis according to Ph. Eur. Chap. 4.1.1 — RS

Code	Size	Packaging	Notes
611069509	100 ml	Plastic bottle	Potassium dichromate solution R1 Ref Ph.Eur 1069502
611069502	1 l	Plastic bottle	Potassium dichromate solution R1 Ref Ph.Eur 1069502

Potassium dichromate 0.167 mol/l (1N)

• Potassio bicromato 0.167 mol/l (1 N) • Potassium dichromate 0.167 mol/l (1 N) • Potasio dicromato 0.167 mol/l (1 N) • Kaliumdichromat 0.167 mol/l (1N)

$K_2Cr_2O_7$
Molecular Weight: 294,19
CAS: 7778-50-9

Classification transport
ONU: 3082
Transport Hazard class: 9
Packing group III



Danger
H302-H332-H315-H319-H334-H317-H340-H350-H360FD-H373-H411-HA26
P271-P280-P284-P304+P340-P305+P351+P338-P342+P311a

Potassium dichromate 0.167 mol/l (1N) > RPE - For analysis — RPE

Code	Size	Packaging	Notes
507536	1 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Potassium dichromate 0.0417 mol/l (0.25N)

• Potassio bicromato 0.0417 mol/l (0.25N) • Potassium dichromate 0.0417 mol/l (0.25N) • Potasio dicromato 0.0417 mol/l (0.25N)
• Kaliumdichromat 0.0417 mol/l (0.25N)

$K_2Cr_2O_7$
Molecular Weight: 294,19
CAS: 7778-50-9

Classification transport
ONU: 3082
Transport Hazard class: 9
Packing group III



Danger
H302-H332-H315-H319-H334-H317-H340-H350-H360FD-H373-H411-HA26
P271-P280-P284-P304+P340-P305+P351+P338-P342+P311a

Potassium dichromate 0.0417 mol/l (0.25N) > RS - For environmental analysis (COD determination) — RS

Description ..... Orange clear liquid Assay (potentiometry) ..... 0.2495 - 0.2505 N

Code	Size	Packaging	Notes
470451	1 l	Glass bottle	

Content is guaranteed for standardized volumes at 20 °C

Potassium dichromate 0.04 mol/l (0.24N) in 80 g/l HgSO4

• Potassio bicromato 0.04 mol/l (0.24N) in 80 g/l HgSO4 • Potassium dichromate 0.04 mol/l (0.24N) dans 80 g/l HgSO4
• Potasio dicromato 0.04 mol/l (0.24N) en 80 g/l HgSO4 • Kaliumdichromat 0.04 mol/l(0.24N) in 80 g/l HgSO4

$K_2Cr_2O_7$
Molecular Weight: 294,18
CAS: 7778-50-9

Classification transport
ONU: 2922
Transport Hazard class: 8
Packing group II



Danger
H290-H302-H331-H314-H334-H317-H340-H350-H360FD-H373-H411-HA26
P280-P284-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P362+P364-P342+P311a-P403+P233

Potassium dichromate 0.04 mol/l (0.24N) in 80 g/l HgSO4 > RS - For environmental analysis (COD determination) — RS

Assay ..... 0,0398 ÷ 0,0402 mol/L

Code	Size	Packaging	Notes
526711	1 l	Glass bottle	
526712	2.5 l	Glass bottle	

Potassium dichromate 0.0167 mol/l (0.1N)

• Potassio bicromato 0.0167 mol/l (0.1 N) • Potassium dichromate 0.0167 mol/l (0.1 N) • Potasio dicromato 0.0167 mol/l (0.1 N)
• Kaliumdichromat 0.0167 mol/l (0.1N)

$K_2Cr_2O_7$
Molecular Weight: 294,19
CAS: 7778-50-9



Danger
H340-H350-H360FD-H412-HEU203-HEU208-HA26
P201-P273-P280-P308+P313-P405-P501a

Potassium dichromate 0.0167 mol/l (0.1N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2 — RS

Code	Size	Packaging	Notes
613004600	1 l	Plastic bottle	Ref Ph.Eur 3004600

Potassium dichromate 0.0167 mol/l (0.1N) > RPE - NORMEX - For analysis — RPE

Description ..... Orange clear liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
470501		Glass ampoule	Volume: 60 ml

Volumetric concentrated solution to prepare 1 L of solution 0,1 N

Potassium dichromate solution 106 g/l

• Potassio bicromato soluzione 106 g/l • Potassium dichromate solution 106g/l • Potasio dicromato solución 106 g/l • Kaliumdichromatlösung 106 g/l

$K_2Cr_2O_7$
Molecular Weight: 294,18
CAS: 7778-50-9

Classification transport
ONU: 3287
Transport Hazard class: 6.1
Packing group II



Danger
H301-H330-H314-H334-H317-H340-H350-H360FD-H335-H372-H411-HA26
P280-P284-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338-P342+P311a-P403+P233

Potassium dichromate solution 106 g/l > RS - For analysis according to Ph. Eur. Chap. 4.1.1 — RS

Code	Size	Packaging	Notes
611069501	1 l	Plastic bottle	Ref Ph.Eur 1069501

Potassium dichromate 0.1414 g/l

• Potassio bicromato 0.1414 g/l • Potassium dichromate 0.1414 g/l • Potasio dicromato 0.1414 g/l • Kaliumdichromat 0.1414 g/l

$K_2Cr_2O_7$
Molecular Weight: 294,18
CAS: 7778-50-9

HEU203

Potassium dichromate 0.1414 g/l > RS - For analysis — RS

Code	Size	Packaging	Notes
504594	1 l	Plastic bottle	

Potassium dihydrogen phosphate ▶ Potassium phosphate monobasic

Potassium disulfate ▶ Potassium pyrosulphate

Potassium ferricyanide

• Potassio ferricianuro • Potassium ferricyanure • Potasio ferricianuro • Kaliumferricyanid

Synonym:
• *Potassium hexacyanoferrate(III)*
• *Red prussiate*

$K_3[Fe(CN)_6]$
Molecular Weight: 329,24
CAS: 13746-66-2
EEC-N: 237-323-3

HEU032

Potassium ferricyanide > RS - For analysis according to Ph. Eur. Chap. 4.1.1 — RS

Code	Size	Packaging	Notes
611069700	500 g	Glass bottle	Ref Ph.Eur 1069800

Potassium ferricyanide > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP — RPE

Description ..... Red orange crystals
Identification ..... Positive
Water-insoluble matter ..... ≤50 ppm
Chloride ..... ≤100 ppm
Ferrocyanide ..... ≤500 ppm
Sulphate ..... ≤100 ppm
Assay (oxidimetric) ..... ≥99.0 %

Code	Size	Packaging	Notes
471364	100 g	Plastic bottle	
471365	250 g	Plastic bottle	
471367	1 kg	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Potassium ferricyanide > RE - Pure

RE

Description ..... Red orange crystals
Moisture ..... ≤ 0.1 %
Water insoluble substances ..... ≤ 0.1 %
K4Fe(CN)6.3H2O ..... ≤ 0.35 %
KCl ..... ≤ 0.3 %
Assay (K3Fe(CN)6) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
360257	1 kg	Plastic bottle	
360258	5 kg	Plastic tank	
360252	25 kg	Plastic bucket	

Potassium ferrocyanide trihydrate

• Potassio ferrocianuro triidrato • Potassium ferrocyanure trihydraté • Potasio ferrocianuro trihidrato
• Kaliumferrocyanid-Trihydrat

Synonym:
• *Potassium hexacyanoferrate(II) trihydrate*
• *Yellow prussiate*

$K_4[Fe(CN)_6],3H_2O$
Molecular Weight: 368,34
CAS: 14459-95-1
EEC-N: 237-722-2

HEU032

Potassium ferrocyanide trihydrate > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611069800	500 g	Glass bottle	Ref Ph.Eur 1069800

Potassium ferrocyanide trihydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Yellow crystals
Identification ..... Positive
Chloride ..... ≤100 ppm
Water-insoluble matter ..... ≤50 ppm
Sulphate ..... Conform ACS
Assay (oxidimetric) ..... 98.5 ÷ 102.0 %

Code	Size	Packaging	Notes
471484	100 g	Plastic bottle	
471485	250 g	Plastic bottle	
471487	1 kg	Plastic bottle	
471488	2.5 kg	Plastic bottle	
471483	25 kg	Drum	

Potassium ferrocyanide trihydrate > RE - Pure

RE

Description ..... Yellow crystals
Identification ..... Positive
Water-insoluble matter ..... ≤ 0.1 %
Assay (oxidimetric) ..... ≥ 98 %

Code	Size	Packaging	Notes
360557	1 kg	Plastic bottle	
360558	5 kg	Plastic tank	
360552	25 kg	Plastic bucket	

Potassium ferrocyanide solution 10%

• Potassio ferrocianuro soluzione 10% • Potassium ferrocyanure solution 10% • Potasio ferrocianuro solución 10% • Kaliumferrocyanidlösung 10%

$K_4Fe(CN)_6.3H_2O$
Molecular Weight: 422,41
CAS: 14459-95-1

HEU032

Potassium ferrocyanide solution 10% > RPE - For analysis

RPE

Description ..... Yellow clear liquid
Identification ..... Positive
Density at 20° C ..... ≥1.06
Assay (oxidimetric) ..... 9.5 ÷ 10.5 %

Code	Size	Packaging	Notes
E471501	1 l	Glass bottle	

Potassium ferrocyanide solution 53 g/l

• Potassio ferricianuro soluzione 53 g/l • Potassium ferrocyanure solution 53 g/l • Potasio ferrocianuro solución 53 g/l • Kaliumferrocyanidlösung 53 g / l

$K_4Fe(CN)_6.3H_2O$
Molecular Weight: 422,41
CAS: 14459-95-1

HEU032

Potassium ferrocyanide solution 53 g/l > RS - For analysis according to Ph. Eur. Chap. 4.1.1 — RS

Code	Size	Packaging	Notes
611069801	100 ml	Plastic bottle	Ref Ph.Eur 1069801

Potassium fluoride

• Potassio fluoruro • Potassium fluorure • Potasio fluoruro • Kaliumfluorid

KF
Molecular Weight: 58,1
CAS: 7789-23-3
EEC-N: 232-151-5

Classification transport
ONU: 1812
Transport Hazard class: 6.1
Packing group III

Danger
H301-H311-H331
P261-P304+P340-P311a-P330-P361+P364-P403+P233

Potassium fluoride > RPE - For analysis — RPE

Description ..... White powder
Identification ..... Positive
Acidity (Hydrofluor ac) ..... ≤0.1 %
Alcalinity (as K2CO3) ..... ≤ 0.1 %
Water not sol. matter ..... ≤0.01 %
Chloride ..... ≤0.05 %
Fluosilicates ..... ≤0.05 %
Sulphate ..... ≤100 ppm
Sulphite ..... ≤100 ppm
Heavy metals (Pb) ..... ≤20 ppm
Cu ..... ≤20 ppm
Fe ..... ≤10 ppm
Ni ..... ≤20 ppm
Pb ..... ≤20 ppm
Zn ..... ≤20 ppm
Assay (non-aqueos medium) ..... ≥98.0 %

Code	Size	Packaging	Notes
471564	100 g	Plastic bottle	
471561	250 g	Plastic bottle	
471562	1 kg	Plastic bottle	
471563	10 kg	Plastic bucket	

Potassium fluoride dihydrate

• Potassio fluoruro biidrato • Potassium fluorure dihydraté • Potasio fluoruro dihidrato • Kaliumfluorid-Dihydrat

$KF.2H_2O$
Molecular Weight: 94,13
CAS: 13455-21-5
EEC-N: 232-151-5

Classification transport
ONU: 1812
Transport Hazard class: 6.1
Packing group III

Danger
H301-H311-H331
P261-P304+P340-P311a-P330-P361+P364-P403+P233

Potassium fluoride dihydrate > RPE - For analysis — RPE

Description ..... White crystals
Identification ..... Positive
Assay ..... ≥ 98.5 %

Code	Size	Packaging	Notes
471555	250 g	Plastic bottle	

Potassium guaiacolsulfonate

• Potassio solfoguaiacolato • Potassium sulfogaïacolate • Potasio sulfoguayacolato • Kaliumguajacolsulfonat

$C_7H_7KO_5S$
Molecular Weight: 242,29
CAS: 1321-14-8
EEC-N: 215-314-5

Potassium guaiacolsulfonate > RE - Pure — RE

Description ..... White crystalline powder
Identification ..... Positive
Water (K.F.) ..... ≤6.0 %
Guaiacol (TLC) ..... ≤0.5 %
Heavy metals (Pb) ..... ≤20 ppm
Sulphate ..... ≤200 ppm
Assay (non-aqueos medium) ..... ≥95.5 %

Code	Size	Packaging	Notes
363807	1 kg	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

di-Potassium hexachloroplatinate

• di-Potassio esacloroplatinato • Potassium chloroplatinate • Di-Potasio hexacloroplatinato
• Dikaliumhexachloroplatinat

Synonym:
Potassium hexachloroplatinate(IV)

K_2PtCl_6
Molecular Weight: 486,01
CAS: 16921-30-5
EEC-N: 240-979-3

Classification transport
ONU: 3290
Transport Hazard class: 6.1
Packing group II

Danger
H301-H318-H334-H317
P261-P284-P301+P310a-P304+P340-P305+P351+P338-P342+P311a

di-Potassium hexachloroplatinate > RPE - For analysis

RPE

Description	Yellow powder	Ca	≤ 20 ppm	Mg	≤ 20 ppm	Ru	≤ 20 ppm
Identification	Positive	Cu	≤ 10 ppm	Pd	≤ 20 ppm	Si	≤ 20 ppm
Ag	≤ 20 ppm	Fe	≤ 20 ppm	Pb	≤ 20 ppm	Assay (gravimetric)	40 ÷ 40.2 % Pt
Au	≤ 20 ppm	Ir	≤ 20 ppm	Rh	≤ 20 ppm		

Code	Size	Packaging	Notes
471127	1 g	Glass ampoule	

Potassium hexacyanoferrate(III) ▶ Potassium ferricyanide

Potassium hexacyanoferrate(II) trihydrate ▶ Potassium ferrocyanide trihydrate

Potassium hydrogen carbonate ▶ Potassium bicarbonate

Potassium hydrogen iodate

• Potassio iodato acido • Potassium iodate acide • Potasio yodato acido • Kaliumjodatsäure

Synonym:
Potassium biiodate

$KIO_3.HIO_3$
Molecular Weight: 389,92
CAS: 13455-24-8
EEC-N: 236-650-9

Classification transport
ONU: 3085
Transport Hazard class: 5.1
Packing group II

Danger
H272-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Potassium hydrogen iodate > RPE - For analysis

RPE

Description	White crystalline powder	Br-BrO3-Cl-ClO3 (Cl)	≤500 ppm	Fe	≤20 ppm	Assay (oxidimetric)	≥99 %
Identification	Positive	Sulphate	≤100 ppm	Ni	≤20 ppm		
Water-insoluble matter	≤500 ppm	Cu	≤20 ppm	Pb	≤20 ppm		

Code	Size	Packaging	Notes
472641	50 g	Glass bottle	

di-Potassium hydrogen phosphate ▶ Potassium phosphate dibasic anhydrous

di-Potassium hydrogen phosphate trihydrate ▶ Potassium phosphate dibasic trihydrate

Potassium hydrogen phthalate

• Potassio ftalato acido • Potassium phtalate acide • Potasio ftalato acido • Kaliumphthalsäure

Synonym:
• *Phtalic acid monopotassium salt*
• *Potassium biphtalate*

$HOOCC_6H_4COOK$
Molecular Weight: 204,23
CAS: 877-24-7
EEC-N: 212-889-4

Potassium hydrogen phthalate > RS - For analysis according to Ph. Eur. Chap. 4.2.1

RS

Code	Size	Packaging	Notes
612000400	50 g	Plastic bottle	Ref Ph.Eur 2000400

Potassium hydrogen phthalate > RS - For volumetry and pHmetry

RS

Description .....White crystals Identification .....Positive pH sol. M/20 at 25° C .....4.001 ÷ 4.011 Assay .....≥ 99.5 %

Code	Size	Packaging	Notes
471913	25 g	Glass bottle	

Potassium hydrogen phthalate > RPE - For analysis

RPE

Description .....White crystals
Identification .....Positive
Loss on drying (110°C).....≤500 ppm
Total nitrogen.....≤10 ppm
Chloride.....≤20 ppm
Water-insoluble matter .....≤30 ppm
Heavy metals (Pb).....≤5 ppm
Total sulphur .....≤20 ppm
Ca.....≤10 ppm
Cd.....≤5 ppm
Co.....≤5 ppm
Cr.....≤10 ppm
Cu.....≤5 ppm
Fe .....≤5 ppm
Mg .....≤10 ppm
Mn .....≤5 ppm
Na.....≤100 ppm
Ni.....≤5 ppm
Pb .....≤5 ppm
Zn .....≤5 ppm
Assay (acidimetric) .....≥99.5 %

Code	Size	Packaging	Notes
471865	250 g	Plastic bottle	
471866	1 kg	Plastic bottle	
471867	2.5 kg	Plastic bottle	

Potassium hydrogen phthalate 0.2 mol/l (0.2N)

• Potassio ftalato acido 0.2 mol/l (0.2N) • Potassium phtalate acide 0.2 mol/l (0.2N)
• Potasio ftalato acido 0.2 mol/l (0.2N) • Kaliumphthalsäure 0.2 mol/l (0.2N)

Synonym:
• *Phtalic acid monopotassium salt*
• *Potassium biphtalate*

$HOOCC_6H_4COOK$
Molecular Weight: 204,23
CAS: 877-24-7

Potassium hydrogen phthalate 0.2 mol/l (0.2N) > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611070001	1 l	Plastic bottle	Ref Ph.Eur 1070001

Potassium hydrogen phthalate 0.1 mol/l (0.1N)

• Potassio ftalato acido 0.1 mol/l (0.1N) • Potassium phtalate acide 0.1 mol/l (0.1N)
• Potasio ftalato acido 0.1 mol/l (0.1N) • Kaliumphthalsäure 0.1 mol/l (0.1N)

Synonym:
• *Phtalic acid monopotassium salt*
• *Potassium biphtalate*

$HOOCC_6H_4COOK$
Molecular Weight: 204,23
CAS: 877-24-7

Classification transport
ONU: 2789
Transport Hazard class: 8
Packing group II



Danger
H226-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Potassium hydrogen phthalate 0.1 mol/l (0.1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor.....0.998 ÷ 1.002

Code	Size	Packaging	Notes
E471926	500 ml	Bottle	

Ready-to-use solution in acetic anhydride

Potassium hydrogen sulfate ▶ Potassium bisulfate

Product specifications are subject to changes. Please visit our website for updates.

Potassium hydroxide, flakes

• Potassio idrossido, scaglie • Potassium hydroxyde, écailles • Potasio hidróxido, escamas
• Kaliumhydroxid, Flocken

Synonym:
Caustic potash

KOH
Molecular Weight: 56,1
CAS: 1310-58-3
EEC-N: 215-181-3

Classification transport
ONU: 1813
Transport Hazard class: 8
Packing group II

Danger
H290-H302-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Potassium hydroxide, flakes > ERBApharm - According to pharmacopoeia: BP-Ph.Eur. — ERBApharm

Description ..... White flakes
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Carbonate..... ≤ 2.0 %
Chloride..... ≤ 200 ppm
Phosphate ..... ≤ 100 ppm
Heavy metals (Pb)..... ≤ 10 ppm
Sulphate..... ≤ 200 ppm
Fe ..... ≤ 10 ppm
Na..... ≤ 1.0 %
Assay (total alkalin.).....85.0 ÷ 100.5 %

Code	Size	Packaging	Notes
362201	25 kg	Sack	
362202	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Potassium hydroxide, flakes > RE - Pure — RE

Description ..... White flakes
Identification ..... Positive
Chloride..... ≤ 80 ppm
Potassio carbonato..... ≤ 1 %
Na (NaOH)..... ≤ 0.9 %
Sulphate..... ≤ 20 ppm
Assay (acidimetric) ..... ≥ 85 %

Code	Size	Packaging	Notes
362257	1 kg	Plastic bottle	
362258	5 kg	Plastic tank	
362251	25 kg	Sack	

Potassium hydroxide, pellets

• Potassio idrossido, gocce • Potassium hydroxyde, pastilles • Potasio hidróxido, gotas
• Kaliumhydroxid, Pellets

Synonym:
Caustic potash

KOH
Molecular Weight: 56,1
CAS: 1310-58-3
EEC-N: 215-181-3

Classification transport
ONU: 1813
Transport Hazard class: 8
Packing group II

Danger
H290-H302-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Potassium hydroxide, pellets > RS - RSE - For electronic use — RS

Description ..... White pellets
Identification ..... Positive
Total nitrogen..... ≤3 ppm
Carbonate..... ≤5000 ppm
Chloride..... ≤10 ppm
Phosphate ..... ≤5 ppm
Heavy metals (Pb)..... ≤2 ppm
Silicate ..... ≤15 ppm
Sulphate..... ≤5 ppm
As ..... ≤0.5 ppm
Ca ..... ≤5 ppm
Cd..... ≤0.1 ppm
Cu..... ≤0.5 ppm
Fe ..... ≤5 ppm
Hg..... ≤0.1 ppm
Mg ..... ≤5 ppm
Mn ..... ≤0.1 ppm
Na..... ≤300 ppm
Ni..... ≤1 ppm
Pb..... ≤1 ppm
Zn..... ≤1 ppm
Assay (alkalimetric)..... ≥86 %

Code	Size	Packaging	Notes
472097	1 kg	Plastic bottle	
472092	5 kg	Plastic jar	

Potassium hydroxide, pellets > RS - For microanalysis — RS

Description ..... White pellets Identification ..... Positive

Code	Size	Packaging	Notes
472086	500 g	Plastic bottle	

Potassium hydroxide, pellets > RPE - For analysis - ACS - ISO

RPE

Description ..... White pellets	Carbonate..... ≤ 2.0 %	Heavy metals (Ag) .....≤ 10 ppm	Assay (alkalimetric)..... ≥ 85 %
Identification ..... Positive	Chloride.....≤ 100 ppm	Fe .....≤ 10 ppm	Mg .....≤ 20 ppm
Ca.....≤ 50 ppm	Phosphate .....≤ 5 ppm	Na.....≤ 500 ppm	
Total nitrogen.....≤ 10 ppm	Sulphate .....≤ 30 ppm	Ni.....≤ 10 ppm	

Code	Size	Packaging	Notes
472171	100 g	Plastic bottle	
472172	500 g	Plastic bottle	
472173	1 kg	Plastic bottle	
472175	5 kg	Plastic jar	

Low content in sodium

Potassium hydroxide, pellets > RPE - For analysis

RPE

Description ..... White pellets	Phosphate .....≤4 ppm	Ca.....≤10 ppm	Cu.....≤ 5 ppm
Identification ..... Positive	Silicate .....≤5 ppm	Insoluble in water .....≤ 50 ppm	Zn.....≤ 5 ppm
Total nitrogen.....≤5 ppm	Sulphate.....≤10 ppm	Heavy metals (as Pb) .....≤ 5 ppm	
Carbonate..... ≤0.6 %	Ag.....≤0.5 ppm	Precipitable with NH4OH.....≤ 100 ppm	
Chloride.....≤10 ppm	Al.....≤10 ppm	As .....≤ 1 ppm	

Code	Size	Packaging	Notes
472057	1 kg	Plastic bottle	
472059	5 kg	Plastic jar	
472056	25 kg	Plastic bucket	

Potassium hydroxide, pellets > ERBApharm - According to pharmacopoeia: FU-Ph.Eur.

ERBApharm

Description ..... White pellets	Na.....≤1.0 %	Phosphate .....≤100 ppm	Fe .....≤10 ppm
Identification ..... Positive	Carbonate.....≤2.0 %	Heavy metals (Pb).....≤10 ppm	Titolo (alcalinità totale) .....85.0 ÷ 100.5 %
Appearance of solution .....Conform Ph.Eur.	Chloride.....≤200 ppm	Sulphate.....≤200 ppm	

Code	Size	Packaging	Notes
362237	1 kg	Plastic bottle	
362239	5 kg	Plastic tank	
362235	25 kg	Sack	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Potassium hydroxide solution 45%

• Potassio idrossido soluzione 45% • Potassium hydroxyde solution 45% • Potasio hidróxido solución 45% • Kaliumhydroxid 45%

Synonym: *Caustic potash*

KOH
Molecular Weight: 56,1
CAS: 1310-58-3

Classification transport
ONU: 1814
Transport Hazard class: 8
Packing group II

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Potassium hydroxide solution 45% > RS - RSE - For electronic use

RS

Description .....Clear liquid	Ni.....≤ 1 ppm	Sulphate.....≤ 3 ppm	Mg .....≤ 3 ppm
Colour (APHA) ..... ≤ 20	Assay (acidimetric) .....45.0 ÷ 46.0 %	Silicate .....≤ 10 ppm	Mn .....≤ 0.5 ppm
Carbonate..... ≤ 1.0 %	Subst. ppt by NH4OH.....≤ 50 ppm	Heavy metals (Pb).....≤ 3 ppm	Hg.....≤ 0.05 ppm
Cr.....≤ 1 ppm	Total nitrogen.....≤ 0.5 ppm	Al.....≤ 1 ppm	Pb.....≤ 1 ppm
Cu.....≤ 0.5 ppm	Chloride.....≤ 5 ppm	As .....≤ 0.5 ppm	Zn.....≤ 1 ppm
Fe .....≤ 1 ppm	Phosphate .....≤ 3 ppm	Ca.....≤ 3 ppm	

Code	Size	Packaging	Notes
472103	5 l	Plastic tank	

Product specifications are subject to changes. Please visit our website for updates.

Potassium hydroxide solution 38% (40° Bé) in water

• Potassio idrossido soluzione 38% (40° Bé) in acqua • Potassium hydroxyde 38% (40°Bé)
• Potasio hidróxido solución 38% (40° Bé) en agua • Kaliumhydroxid 38% (40 ° Bé) in Wasser

Synonym:
Caustic potash

KOH
Molecular Weight: 56,1
CAS: 1310-58-3

Classification transport
ONU: 1814
Transport Hazard class: 8
Packing group II

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Potassium hydroxide solution 38% (40° Bé) in water > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Carbonate ..... ≤1.5 %
Chloride ..... ≤40 ppm
Nitrogen compounds (N) ..... ≤10 ppm
Phosphate ..... ≤10 ppm
Heavy metals (Pb) ..... ≤20 ppm
Subst. ppt by NH4OH ..... ≤200 ppm
Silicate ..... ≤100 ppm
Sulphate ..... ≤20 ppm
Fe ..... ≤5 ppm
Assay ..... 38 - 40 %

Code	Size	Packaging	Notes
E472151	1 l	Plastic bottle	
E472152	35 kg	Plastic drum	

Potassium hydroxide solution 33%

• Potassio idrossido soluzione 33% • Potassium hydroxyde 33% • Potasio hidróxido solución 33% • Kaliumhydroxid 33%

KOH
Molecular Weight: 56,1
CAS: 1310-58-3

Classification transport
ONU: 1814
Transport Hazard class: 8
Packing group II

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Potassium hydroxide solution 33% > RS - For analysis

RS

Code	Size	Packaging	Notes
PS0766/22	5 l	Plastic tank	

Potassium hydroxide solution 28%

• Potassio idrossido soluzione 28% • Potassium hydroxyde 28% • Potasio hidróxido solución 28%
• Kaliumhydroxid 28%

Synonym:
Caustic potash

KOH
Molecular Weight: 56,1
CAS: 1310-58-3

Classification transport
ONU: 1814
Transport Hazard class: 8
Packing group II

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Potassium hydroxide solution 28% > RS - For gas analysis according to Orsat

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... ~ 1.27
Assay (alkalimetric) ..... 27 ÷ 29 %

Code	Size	Packaging	Notes
E472221	1 l	Plastic bottle	

Potassium hydroxide solution 3% in ethanol

• Potassio idrossido soluzione 3% in etanolo • Potassium hydroxyde solution 3% dans l'ethanol
• Potasio hidróxido solución 3% en etanol • Kaliumhydroxidlösung 3% in Ethanol

Synonym:
Caustic potash

KOH
Molecular Weight: 56,1
CAS: 1310-58-3

Classification transport
ONU: 2920

Potassium hydroxide solution 3% in ethanol > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611070303	100 ml	Glass bottle	Ref Ph.Eur 1070303

Potassium hydroxide 2 mol/l (2N) in ethanol

• Potassio idrossido 2 mol/l (2N) in etanolo • Potassium hydroxyde 2 mol/l (2N) dans l'ethanol
• Potasio hidróxido 2 mol/l (2N) en etanol • Kaliumhydroxid 2 mol / l (2N) in Ethanol

Synonym:
Caustic potash

KOH
Molecular Weight: 56,1
CAS: 1310-58-3

Classification transport
ONU: 2733
Transport Hazard class: 3
Packing group II



Danger
H225-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Potassium hydroxide 2 mol/l (2N) in ethanol > RS - For analysis according to Ph. Eur. Chap. 4.1.1 — RS

Code	Size	Packaging	Notes
611070301	100 ml	Plastic bottle	Ref Ph.Eur 1070301

Potassium hydroxide 1 mol/l (1N)

• Potassio idrossido 1 mol/l (1N) • Potassium hydroxyde 1 mol/l (1N) • Potasio hidróxido 1 mol/l (1N)
• Kaliumhydroxid 1 mol/l (1N)

Synonym:
Caustic potash

KOH
Molecular Weight: 56,1
CAS: 1310-58-3

Classification transport
ONU: 1814
Transport Hazard class: 8
Packing group II



Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Potassium hydroxide 1 mol/l (1N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2 — RS

Code	Size	Packaging	Notes
613009100	1 l	Plastic bottle	Ref Ph.Eur 3009100

Potassium hydroxide 1 mol/l (1N) > RPE - For analysis — RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.999 - 1.001 N NIST 84 I

Code	Size	Packaging	Notes
472287000	1 l	Plastic bottle	Certified with NIST traceability
472282000	5 l	Kubidos	Certified with NIST traceability
472281000	10 l	Kubidos	Certified with NIST traceability

56.11 g of KOH. Volumetric solution ready-to-use. Content is guaranteed for standardized volumes at 20°C.

Potassium hydroxide 1 mol/l (1N) > RPE - NORMEX - For analysis — RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
472311		Plastic ampoule	Volume: 165 ml

56,11 g of KOH. Volumetric concentrated solution to prepare 1 L of solution 1 N

Potassium hydroxide 0.5 mol/l (0.5N)

• Potassio idrossido 0.5 mol/l (0.5N) • Potassium hydroxyde 0.5 mol/l (0.5N)
• Potasio hidróxido 0.5 mol/l (0.5N) • Kaliumhydroxid 0.5 mol/l (0.5N)

Synonym:
Caustic potash

KOH
Molecular Weight: 56,1
CAS: 1310-58-3

Classification transport
ONU: 1814
Transport Hazard class: 8
Packing group II



Danger
H302-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Potassium hydroxide 0.5 mol/l (0.5N) > RPE - For analysis — RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.4995 - 0.5005 N NIST 84 I

Code	Size	Packaging	Notes
472337000	1 l	Plastic bottle	Certified with NIST traceability
472332000	5 l	Kubidos	Certified with NIST traceability
472331000	10 l	Kubidos	Certified with NIST traceability

28.055 g of KOH. Volumetric solution ready-to-use.

Product specifications are subject to changes. Please visit our website for updates.

Potassium hydroxide 0.5 mol/l (0.5N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
472391		Plastic ampoule	Volume: 55 ml

28,055 g of KOH. Volumetric concentrated solution to prepare 1 L of solution 0,5 N

Potassium hydroxide 0.5 mol/l (0.5N) in ethanol

• Potassio idrossido 0.5 mol/l (0.5N) in alcole etilico • Potassium hydroxyde 0.5 mol/l (0.5N) dans l'éthanol
• Potassium hydroxide 0.5 mol/l (0.5N) in ethanol • Kaliumhydroxid 0.5 mol/l (0.5N) in ethanol

Synonym:
Caustic potash

KOH
Molecular Weight: 56,1
CAS: 1310-58-3

Classification transport
ONU: 2924
Transport Hazard class: 3
Packing group II



Danger
H225-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Potassium hydroxide 0.5 mol/l (0.5N) in ethanol > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613004900	1 l	Glass bottle	Ref Ph.Eur 3004900

Potassium hydroxide 0.5 mol/l (0.5N) in ethanol > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611070302	1 l	Glass bottle	Ref Ph.Eur 1070302

Potassium hydroxide 0.5 mol/l (0.5N) in ethanol > RPE - For analysis

RPE

Description ..... Colourless to orange liquid Assay (potentiometry) ..... 0.499 - 0.501 N NIST 84 I

Code	Size	Packaging	Notes
472021000	1 l	Plastic bottle	Certified with NIST traceability
472022000	1 l	Glass bottle	Certified with NIST traceability

28.055 g of KOH. Volumetric solution ready-to-use

Potassium hydroxide 0.5 mol/l (0.5N) in methanol

• Potassio idrossido 0.5 mol/l (0.5N) in alcole metilico
• Potassium hydroxyde 0.5 mol/l (0.5N) dans le méthanol • Potasio hidróxido 0.5 mol/l (0.5N) en metanol
• Kaliumhydroxid 0.5 mol/l (0.5N) in methanol

Synonym:
Caustic potash

KOH
Molecular Weight: 56,1
CAS: 1310-58-3

Classification transport
ONU: 3286
Transport Hazard class: 3
Packing group II



Danger
H225-H301-H314-H370
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338

Potassium hydroxide 0.5 mol/l (0.5N) in methanol > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.499 - 0.501 N NIST 84 I

Code	Size	Packaging	Notes
472366000	500 ml	Plastic bottle	Certified with NIST traceability
472364000	1 l	Glass bottle	Certified with NIST traceability
472367000	5 l	Plastic tank	Certified with NIST traceability

28.055 g of KOH. Volumetric solution ready-to-use

Potassium hydroxide 0.46 mol/l (0.46N)

• Potassio idrossido 0.46 mol/l (0.46N) • Potassium hydroxyde 0.46 mol/l (0.46N)
• Potasio hidróxido 0.46 mol/l (0.46N) • Kaliumhydroxid 0.46 mol/l (0.46N)

Synonym:
Caustic potash

KOH
Molecular Weight: 56,1
CAS: 1310-58-3

Classification transport
ONU: 1814
Transport Hazard class: 8
Packing group II



Danger
H302-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Potassium hydroxide 0.46 mol/l (0.46N) > RS - For agroalimentary analysis

RS

Description ..... Clear colourless liquid Colour ..... ≤ 10 APHA Assay ..... 0.455 ÷ 0.465 N

Code	Size	Packaging	Notes
502212	5 l	Plastic tank	

Potassium hydroxide 0.25 mol/l (0.25N)

• Potassio idrossido 0.25 mol/l (0.25N) • Potassium hydroxyde 0.25 mol/l (0.25N)
• Potasio hidróxido 0.25 mol/l (0.25N) • Kaliumhydroxid 0.25 mol/l (0.25N)

Synonym:
Caustic potash

KOH
Molecular Weight: 56,1
CAS: 1310-58-3

Classification transport
ONU: 1814
Transport Hazard class: 8
Packing group III



Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Potassium hydroxide 0.25 mol/l (0.25N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.2498 - 0.2503 N NIST 84 I

Code	Size	Packaging	Notes
472427000	1 l	Plastic bottle	Certified with NIST traceability
472422000	5 l	Kubidos	Certified with NIST traceability
472421000	10 l	Kubidos	Certified with NIST traceability

14.027 g of KOH. Volumetric solution ready-to-use

Potassium hydroxide 0.23 mol/l (0.23N)

• Potassio idrossido 0.23 mol/l (0.23N) • Potassium hydroxyde 0.23 mol/l (0.23N)
• Potasio hidróxido 0.23 mol/l (0.23N) • Kaliumhydroxid 0.23 mol/l (0.23N)

Synonym:
Caustic potash

KOH
Molecular Weight: 56,1
CAS: 1310-58-3

Classification transport
ONU: 1814
Transport Hazard class: 8
Packing group III



Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Potassium hydroxide 0.23 mol/l (0.23N) > RS - For agroalimentary analysis

RS

Description ..... Clear colourless liquid Assay ..... 0.225 ÷ 0.235 N Colour ..... ≤ 10 APHA

Code	Size	Packaging	Notes
502092	5 l	Plastic tank	

Potassium hydroxide 0.1 mol/l (0.1N)

• Potassio idrossido 0.1 mol/l (0.1N) • Potassium hydroxyde 0.1 mol/l (0.1N)
• Potasio hidróxido 0.1 mol/l (0.1N) • Kaliumhydroxid 0.1 mol/l (0.1N)

Synonym:
Caustic potash

KOH
Molecular Weight: 56,1
CAS: 1310-58-3

Classification transport
ONU: 1814
Transport Hazard class: 8
Packing group III



Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Potassium hydroxide 0.1 mol/l (0.1N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613004800	1 l	Plastic bottle	Ref Ph.Eur 3004800

Product specifications are subject to changes. Please visit our website for updates.

Potassium hydroxide 0.1 mol/l (0.1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.0999 - 0.1001 N NIST 84 I

Code	Size	Packaging	Notes
472457000	1 l	Plastic bottle	Certified with NIST traceability
472452000	5 l	Kubidos	Certified with NIST traceability
472451000	10 l	Kubidos	Certified with NIST traceability

5.61 g of KOH. Volumetric solution ready-to-use

Potassium hydroxide 0.1 mol/l (0.1N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
472511		Plastic ampoule	Volume: 55 ml

5,611 g KOH. Volumetric concentrated solution to prepare 1 L of solution 0,1 N

Potassium hydroxide 0.1 mol/l (0.1N) in ethanol

• Potassio idrossido 0.1 mol/l (0.1N) in alcole etilico • Potassium hydroxyde 0.1 mol/l (0.1N) dans l'ethanol
• Potasio hidróxido 0.1 mol/l (0.1N) en etanol • Kaliumhydroxid 0.1 mol/l (0.1N) in ethanol

Synonym:
Caustic potash

KOH
Molecular Weight: 56,1
CAS: 1310-58-3

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P332+P313

Potassium hydroxide 0.1 mol/l (0.1N) in ethanol > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613005100	1 l	Glass bottle	Ref Ph.Eur 3005100

Potassium hydroxide 0.1 mol/l (0.1N) in ethanol > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.0998 - 0.1002 N NIST 84 I

Code	Size	Packaging	Notes
472041000	1 l	Glass bottle	Certified with NIST traceability
472042000	1 l	Plastic bottle	Certified with NIST traceability

5.61 g of KOH. Volumetric solution ready-to-use

Potassium hydroxide 0.1 mol/l (0.1N) in methanol

• Potassio idrossido 0.1 mol/l (0.1N) in alcole metilico
• Potassium hydroxyde 0.1 mol/l (0.1N) dans le méthanol • Potasio hidróxido 0.1 mol/l (0.1N) en metanol
• Kaliumhydroxid 0.1 mol/l (0.1N) in methanol

Synonym:
Caustic potash

KOH
Molecular Weight: 56,1
CAS: 1310-58-3

Classification transport
ONU: 1230
Transport Hazard class: 3
Packing group II



Danger
H225-H301-H311-H331-H315-H319-H370
P210-P280-P303+P361+P353-P304+P340-
P305+P351+P338-P403+P233

Potassium hydroxide 0.1 mol/l (0.1N) in methanol > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.0998 - 0.1002 N NIST 84 I

Code	Size	Packaging	Notes
472486000	500 ml	Plastic bottle	Certified with NIST traceability
472484000	1 l	Glass bottle	Certified with NIST traceability

5.61 g of KOH. Volumetric solution ready-to-use

Potassium hydroxide in solution

• Potassio idrossido in soluzione • Potassium hydroxyde en solution • Potasio hidróxido en solución • Kaliumhydroxid in Lösung

KOH
Molecular Weight: 56,11
CAS: 1310-58-3

Classification transport
ONU: 2924
Transport Hazard class: 3
Packing group II



Danger
H225-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Potassium hydroxide in solution > RS - For analysis according to Ph. Eur. Chap. 4.2.2 — RS

Code	Size	Packaging	Notes
613005001	500 ml	Glass bottle	Ref Ph.Eur 3005001

Potassium Iodate

• Potassio iodato • Potassium iodate • Potasio yodato • Kaliumiodat

KIO_3
Molecular Weight: 214
CAS: 7758-05-6
EEC-N: 231-831-9

Classification transport
ONU: 1479
Transport Hazard class: 5.1
Packing group II



Danger
H272
P210-P220-P280-P370+P378a-P501a

Potassium Iodate > RPE - For analysis — RPE

Description ..... white crystalline powder
Identification (I.R.)..... Positive
Acidity or alcalinity..... Passes test
Loss on drying ..... ≤ 0.5 %
Iodide ..... ≤ 20 ppm
Heavy metals (Pb)..... ≤ 20 ppm
As ..... ≤ 3 ppm
Assay (dried base)..... 99.0 ÷ 101.0 %
SO4..... < 50 ppm

Code	Size	Packaging	Notes
472563	50 g	Glass bottle	
472565	250 g	Glass bottle	

Potassium iodate 0.05 mol/l (0.3N)

• Potassio iodato 0.05 mol/l (0.3N) • Potassium iodate 0.05 mol/l (0.3N) • Potasio yodato 0.05 mol/l (0.3N) • Kaliumjodat 0.05 mol/l (0.3N)

KIO_3
Molecular Weight: 214
CAS: 7758-05-6

Potassium iodate 0.05 mol/l (0.3N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2 — RS

Code	Size	Packaging	Notes
613005200	1 l	Glass bottle	Ref Ph.Eur 3005200

Potassium iodate 0.0167 mol/l (0.1N)

• Potassio iodato 0.0167 mol/l (0.1N) • Potassium iodate 0.0167 mol/l (0.1N) • Potasio yodato 0.0167 mol/l (0.1N) • Kaliumjodat 0.0167 mol/l (0.1N)

KIO_3
Molecular Weight: 214
CAS: 7758-05-6

Potassium iodate 0.0167 mol/l (0.1N) > RPE - NORMEX - For analysis — RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
472601		Glass ampoule	Volume: 60 ml

Volumetric concentrated solution to prepare 1 L of solution 0.1 N

Product specifications are subject to changes. Please visit our website for updates.

Potassium iodate 0.00167 mol/l (0.01N)

• Potassio iodato 0.00167 mol/l (0.01N) • Potassium iodate 0.00167 mol/l (0.01N) • Potasio yodato 0.00167 mol/l (0.01N)
• Kaliumjodat 0.00167 mol/l (0.01N)

KIO_3
Molecular Weight: 214
CAS: 7758-05-6

Potassium iodate 0.00167 mol/l (0.01N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
472631		Glass ampoule	Volume: 60 ml

0,3567 g KIO3. Volumetric concentrated solution to prepare 1 L of solution 0,01 N

Potassium iodide

• Potassio ioduro • Potassium iodure • Potasio yoduro • Kaliumiodid

KI
Molecular Weight: 166,01
CAS: 7681-11-0
EEC-N: 231-659-4

Danger
H372
P260-P264-P270-P314-P501a

Potassium iodide > RS - For microanalysis

RS

Description ..... White cryst. powder Identification ..... Positive

Code	Size	Packaging	Notes
472821	100 g	Glass bottle	

Potassium iodide > RPE - For analysis - ACS

RPE

Description ..... White or colourless solid
Identification ..... Positive
pH sol. 5% at 25° C ..... 6.0 ÷ 9.2
Loss on drying at 150°C..... ≤ 0.2 %
Water-insoluble matter ..... ≤ 50 ppm
Chloride + bromide (Cl) ..... ≤ 100 ppm
Iodate..... ≤ 3 ppm
Phosphate ..... ≤ 10 ppm
Sulphate..... ≤ 50 ppm
Heavy metals (Pb)..... ≤ 5 ppm
Ba..... ≤ 20 ppm
Ca..... ≤ 20 ppm
Fe ..... ≤ 3 ppm
Mg ..... ≤ 10 ppm
Na..... ≤ 50 ppm
Assay (oxidimetric) ..... ≥ 99.0 %

Code	Size	Packaging	Notes
472735	250 g	Plastic bottle	
472737	1 kg	Plastic bottle	
472736	25 kg	Plastic bucket	

Potassium iodide > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP-DAB

ERBApharm

DescriptionPolvere bianca o quasi, o cristalli incolori
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Alcalinity.....Conform Ph.Eur.
Thiosulphate .....Conform Ph.Eur.
Nitrat.nitrit.and NH4OH .... Conform USP-NF
Thiosulfates and barium .. Conform USP-NF
Loss on drying ..... ≤1.0 %
Iodate.....≤4 ppm
Heavy metals (Pb).....≤10 ppm
Sulphate.....≤150 ppm
Fe .....≤20 ppm
Assay (oxidimetric) .....99.0 ÷ 100.5 % s.s.

Code	Size	Packaging	Notes
362405	250 g	Plastic bottle	
362407	1 kg	Plastic bottle	
362409	5 kg	Plastic tank	
362403	10 kg	Plastic tank	
362402	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Potassium iodide solution 10%

• Potassio ioduro soluzione 10% • Potassium iodure solution 10% • Potasio yoduro solución 10% • Kaliumiodid 10%

KI
Molecular Weight: 166,01
CAS: 7681-11-0

Danger
H372
P260-P264-P270-P314-P501a

Potassium iodide solution 10% > RPE - For analysis

RPE

Description ..... Clear colourless liquid Density at 20° C ..... 1.072 - 1.080 Assay ..... 9.5 ÷ 10.5 % p/p

Code	Size	Packaging	Notes
472831	500 ml	Glass bottle	

Potassium iodide solution 3.9%

• Potassio ioduro soluzione 3.9% • Potassium iodure solution 3.9% • Potasio yoduro solución 3.9% • Kaliumiodid 3.9%

KI
Molecular Weight: 166,01
CAS: 7681-11-0

Warning
H373
P260-P314-P501a

Potassium iodide solution 3.9% > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay .....3.7 - 4.1 % p/v

Code	Size	Packaging	Notes
472815000	250 ml	Glass bottle	

Content is guaranteed for standardized volumes at 20°C. Keep tightly closed

Potassium iodide solution

• Potassio ioduro soluzione • Potassium iodure solution • Potasio yoduro solución • Kaliumjodidlösung

KI
Molecular Weight: 166,01
CAS: 7681-11-0

Potassium iodide solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611070504	100 ml	Glass bottle	Solution saturated Ref Ph.Eur 1070504
611070505	100 ml	Glass bottle	Solution iodinated R1 Ref Ph.Eur 1070505
611070502	1 l	Plastic bottle	Potassium iodide solution 166 g/l Ref Ph.Eur 1070502

Storage: protected from light

Potassium iodide starch paper

• Cartina di amido di potassio • Papier potassium iodure d' amidon • Almidón de yoduro de potasio papel • Kaliumiodid-Stärke-Papier

Potassium iodide starch paper > RS - For pHmetry

RS

Code	Size	Packaging	Notes
434980000	1 roll	Dispenser	Paper starch iodide, Color change: White --> Blue-purple

Product specifications are subject to changes. Please visit our website for updates.

Potassium iodobismuthate solution

• potassio iodobismutato soluzione • Potassium iodobismuthate solution • Potasio yodobismutato solución • Kaliumjodobismutatlösung



Danger
H318-H372
P260-P264-P280i-P305+P351+P338-P310a-P501a

Potassium iodobismuthate solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611070600	100 ml	Glass bottle	Ref Ph.Eur 1070600
611070602	100 ml	Glass bottle	Potassium iodobismuthate solution R2 Ref Ph.Eur 1070602

Potassium metabisulfite

• Potassio metabisolfito • Potassium métabisulfite • Potasio metabisulfito • Kaliumdisulfit

Synonym:
Potassium disulfite

$K_2S_2O_5$
Molecular Weight: 222,33
CAS: 16731-55-8
EEC-N: 240-795-3



Warning
H319-H335-HEU031
P261-P271-P304+P340-P305+P351+P338-P337+P313-P403+P233

Potassium metabisulfite > ERBApharm - According to pharmacopoeia: NF

ERBApharm

Description ..... White powder or pieces
Identification ..... Positive
Heavy metals (Pb) ..... ≤10 ppm
Fe ..... ≤10 ppm
Assay (SO2) ..... 51.8 ÷ 57.6 %

Code	Size	Packaging	Notes
362627	1 kg	Plastic bottle	
362629	5 kg	Plastic tank	
362622	10 kg	Plastic tank	
362623	25 kg	Plastic bucket	

This product should be used in compliance with the current legislation.

Potassium nitrate

• Potassio nitrato • Potassium nitrate • Potasio nitrato • Kaliumnitrat

KNO_3
Molecular Weight: 101,1
CAS: 7757-79-1
EEC-N: 231-818-8

Classification transport
ONU: 1486
Transport Hazard class: 5.1
Packing group III



Danger
H272
P210-P220-P280-P370+P378a-P501a

Potassium nitrate > RPE - For analysis - ISO - ACS - Reag.Ph.Eur. - Reag.USP

RPE

Description ..... White crystalline powder
Identification ..... Positive
pH sol. 5% in H20 ..... 4.5 ÷ 8.5
Water-insoluble matter ..... ≤50 ppm
Chloride ..... ≤20 ppm
Phosphate ..... ≤5 ppm
Iodate ..... ≤5 ppm
Nitrite ..... ≤10 ppm
Sulphate ..... ≤30 ppm
Heavy metals (Pb) ..... ≤5 ppm
Ca ..... ≤50 ppm
Fe ..... ≤3 ppm
Mg ..... ≤20 ppm
Na ..... ≤50 ppm
Assay (acidimetric) ..... ≥99.0 %

Code	Size	Packaging	Notes
473006	100 g	Plastic bottle	
473007	1 kg	Plastic bottle	
473009	5 kg	Plastic tank	
473001	25 kg	Drum	

Potassium nitrate > ERBApharm - According to pharmacopoeia: BP-Ph.Eur.

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Acidity or alkalinity.....Conform Ph.Eur.
Reducing substances .....Conform Ph.Eur.
Ammonium .....≤ 100 ppm
Heavy metals (Pb).....≤ 10 ppm
Sulphate .....≤ 150 ppm
Ca.....≤ 100 ppm
Fe .....≤ 20 ppm
Na.....≤ 0.10 %
Loss on drying ..... ≤ 0.5 %
Assay (acidimetric) .....99.0 ÷ 101.0 % s.s.

Code	Size	Packaging	Notes
363007	1 kg	Plastic bottle	
363009	5 kg	Plastic tank	
363002	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Potassium nitrate 1 mol/l (1N)

• Potassio nitrato 1 mol/l (1N) • Potassium nitrate 1 mol/l (1N) • Potasio nitrato 1 mol/l (1N) • Kaliumnitrat 1 mol/l (1N)

KNO_3
Molecular Weight: 101,1
CAS: 7757-79-1

Potassium nitrate 1 mol/l (1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20°C ..... 1.057 ÷ 1.061

Code	Size	Packaging	Notes
473045	250 ml	Plastic bottle	

Potassium nitrite

• Potassio nitrito • Potassium nitrite • Potasio nitrito • Nitrit-Kalium

KNO_2
Molecular Weight: 85,1
CAS: 7758-09-0
EEC-N: 231-832-4

Classification transport
ONU: 1488
Transport Hazard class: 5.1
Packing group II



Danger
H272-H301-H400
P210-P220-P264-P280-P301+P310a-P330

Potassium nitrite > RPE - For analysis

RPE

Description .....Yellowish crystals
Identification ..... Positive
Chloride.....≤ 200 ppm
Sulphate .....≤ 200 ppm
Fe .....≤ 10 ppm
Pb .....≤ 5 ppm
Assay ..... ≥ 97 %

Code	Size	Packaging	Notes
473084	100 g	Glass bottle	

Potassium oxalate monohydrate

• Potassio ossalato monoidrato • Potassium oxalate monohydraté • Potasio oxalato monohidrato
• Kaliumoxalat-Monohydrat

Synonym:
• *Ethanedioic acid*
• *Oxalic acid potassium salt*

$(COOK)_2.H_2O$
Molecular Weight: 184,23
CAS: 6487-48-5
EEC-N: 209-506-8



Warning
H302-H312
P264-P270-P280h-P301+P312a-P330-P501a

Potassium oxalate monohydrate > RPE - For analysis - ACS

RPE

Description ..... White crystalline powder
Identification ..... Positive
Neutrality .....Conform
Ready carbonizable substances.....Conform
Water-insoluble matter .....≤100 ppm
Ammonium .....≤20 ppm
Chloride.....≤20 ppm
Sulphate .....≤100 ppm
Heavy metals (Pb).....≤20 ppm
Fe .....≤10 ppm
Na.....≤200 ppm
Assay (oxidimetric) .....98.5 ÷ 101.0 %

Code	Size	Packaging	Notes
473135	250 g	Plastic bottle	
473137	1 kg	Plastic bottle	
473133	25 kg	Plastic bucket	

Product specifications are subject to changes. Please visit our website for updates.

Potassium m-periodate

• Potassio (meta)-periodato • Potassium m-périodate • Potasio (meta)-peryodato • Kaliumperiodat

Synonym:
Periodic acid potassium salt

KIO_4
Molecular Weight: 230
CAS: 7790-21-8
EEC-N: 232-196-0

Classification transport
ONU: 1479
Transport Hazard class: 5.1
Packing group II

Danger
H272
P210-P220-P280-P370+P378a-P501a

Potassium m-periodate > RPE - For analysis - Reag. Ph. Eur.

RPE

Description ..... White crystalline powder
Identification ..... Positive
Mn ..... ≤ 1 ppm
Assay (iodometric) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
473332	25 g	Glass bottle	
473334	100 g	Glass bottle	

Potassium permanganate

• Potassio permanganato • Potassium permanganate • Potasio permanganato • Kaliumpermanganat

$KMnO_4$
Molecular Weight: 158,04
CAS: 7722-64-7
EEC-N: 231-760-3

Classification transport
ONU: 1490
Transport Hazard class: 5.1
Packing group II

Danger
H272-H302-H410
P210-P220-P264-P280-P301+P312a-P501a

Potassium permanganate > RS - For enviromental analysis - ACS

RS

Description ..... Dark violet crystals
Identification ..... Positive
Chloride & Chlorate(Cl) ..... ≤ 50 ppm
Water-insoluble matter ..... ≤ 0.2 %
Sulphate ..... ≤ 200 ppm
Hg ..... ≤ 0.05 ppm
Assay (oxidimetric) ..... ≥ 99.0 %

Code	Size	Packaging	Notes
476671	100 g	Glass bottle	

Low content in Hg

Potassium permanganate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Dark violet crystals
Identification ..... Positive
Chloride & Chlorate(Cl) ..... ≤ 50 ppm
Water-insoluble matter ..... ≤ 0.2 %
Sulphate ..... ≤ 200 ppm
Assay (oxidimetric) ..... 99.0 ÷ 100.5 %
Appearance of solution ..... Conform Ph.Eur.

Code	Size	Packaging	Notes
473384	100 g	Plastic bottle	
473385	250 g	Plastic bottle	
473387	1 kg	Plastic bottle	
473381	25 kg	Metal drum	

Potassium permanganate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP-DAB

ERBApharm

Description ..... Dark-violet crystals
Identification ..... Positive
Appearance of solution ..... Conform Ph.Eur.
Loss on drying ..... ≤ 0.5 %
Water not sol. matter ..... ≤ 0.2 %
Chloride ..... ≤ 200 ppm
Sulphate ..... ≤ 500 ppm
Assay (oxidimetric) ..... 99.0 ÷ 100.5 %

Code	Size	Packaging	Notes
363107	1 kg	Plastic bottle	
363109	5 kg	Plastic tank	
363101	25 kg	Metal drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Potassium permanganate solution 3%

• Potassio permanganato soluzione 3% • Potassium permanganate solution 3% • Potasio permanganato solución 3% • Kaliumpermanganatlösung 3%

$KMnO_4$
Molecular Weight: 158,04
CAS: 7722-64-7

H411
P273-P391-P501a

Potassium permanganate solution 3% > RS - For analysis according to Ph. Eur. Chap. 4.1.1 — RS

Code	Size	Packaging	Notes
611070902	1 l	Glass bottle	Ref Ph.Eur 1070902

Potassium permanganate 0.2 mol/l (1N)

• Potassio permanganato 0.2 mol/l (1N) • Potassium permanganate 0.2 mol/l (1N) • Potasio permanganato 0.2 mol/l (1N) • Kaliumpermanganat 0.2 mol/l (1N)

$KMnO_4$
Molecular Weight: 158,04
CAS: 7722-64-7

H411
P273-P391-P501a

Potassium permanganate 0.2 mol/l (1N) > RPE - For analysis — RPE

Description ..... Clear purple liquid Assay (potentiometry) ..... 0.99 - 1.01 N NIST 136 ..... e

Code	Size	Packaging	Notes
473514000	1 l	Glass bottle	Certified with NIST traceability

3.1606 g of KMnO4. Volumetric solution ready-to-use

Potassium permanganate 0.02 mol/l (0.1N)

• Potassio permanganato 0.02 mol/l (0.1N) • Potassium permanganate 0.02 mol/l (0.1N) • Potasio permanganato 0.02 mol/l (0.1N) • Kaliumpermanganat 0.02 mol/l (0.1N)

$KMnO_4$
Molecular Weight: 158,04
CAS: 7722-64-7

Potassium permanganate 0.02 mol/l (0.1N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2 — RS

Code	Size	Packaging	Notes
613005301	100 ml	Glass bottle	Ref Ph.Eur 3005300
613005309	250 ml	Glass bottle	Ref Ph.Eur 3005300
613005300	1 l	Glass bottle	Ref Ph.Eur 3005300

Potassium permanganate 0.02 mol/l (0.1N) > RPE - For analysis — RPE

Description ..... Purple clear liquid Assay (potentiometry) ..... 0.0998 - 0.1002 N NIST 136 ..... e

Code	Size	Packaging	Notes
473567000	1 l	Glass bottle	Certified with NIST traceability
473565000	5 l	Kubidos	Certified with NIST traceability

3.1606 g of KMnO4. Volumetric solution ready-to-use

Potassium permanganate 0.02 mol/l (0.1N) > RPE - NORMEX - For analysis — RPE

Description ..... Clear purple liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
473591		Glass ampoule	Volume: 65 ml

3,161 g KMnO4. Volumetric concentrated solution to prepare 1 L of solution 0,1 N

Product specifications are subject to changes. Please visit our website for updates.

Potassium permanganate 0.002 mol/l (0.01N)

• Potassio permanganato 0.002 mol/l (0.01N) • Potassium permanganate 0.002 mol/l (0.01N) • Potasio permanganato 0.002 mol/l (0.01N)
• Kaliumpermanganat 0.002 mol/l (0.01N)

$KMnO_4$
Molecular Weight: 158,04
CAS: 7722-64-7

H412
P273-P501a

Potassium permanganate 0.002 mol/l (0.01N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear purple liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
473661		Glass ampoule	Volume: 60 ml

0,3161 g KMnO4. Volumetric concentrated solution to prepare 1 L of solution 0,01 N

Potassium permanganate and phosphoric acid solution

• Potassio permanganato e acido fosforico soluzione • Potassium permanganate - Solution phosphorique • Potasio permanganato y ácido fosfórico solución
• Kaliumpermanganat - Phosphorlösung

$KMnO_4$
Molecular Weight: 158,04
CAS: 7722-64-7

Potassium permanganate and phosphoric acid solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611070901	100 ml	Glass bottle	Ref Ph.Eur 1070901

Potassium persulfate

• Potassio persolfato • Potassium persulfate • Potasio persulfato • Kaliumpersulfat

Synonym:
Potassium peroxodisulfate

$K_2S_2O_8$
Molecular Weight: 270,31
CAS: 7727-21-1
EEC-N: 231-781-8

Classification transport
ONU: 1492
Transport Hazard class: 5.1
Packing group III



Danger
H272-H302-H315-H319-H334-H317-H335
P210-P280-P284-P304+P340-P305+P351+P338-P342+P311a-P403+P233

Potassium persulfate > RE - Pure

RE

Description ..... White crystalline powder Identification ..... Positive Acidity(Sulphuric acid) ..... ≤ 0.15 % Fe ..... ≤ 5 ppm Assay (oxidimetric) ..... ≥ 99 %

Code	Size	Packaging	Notes
473701	1 kg	Plastic bottle	

Potassium phosphate dibasic anhydrous

• Potassio fosfato bibasico anidro • Potassium phosphate dibasique anhydre
• Potasio fosfato dibásico anhidro • Dikaliumhydrogenphosphat wasserfrei

Synonym:
Dipotassium hydrogenphosphate

K_2HPO_4
Molecular Weight: 174,18
CAS: 7758-11-4
EEC-N: 231-834-5

Potassium phosphate dibasic anhydrous > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White crystalline powder	Loss on drying (105°C)..... ≤1.0 %	Chloride..... ≤30 ppm	Fe ..... ≤10 ppm
Identification ..... Positive	Water-insoluble matter ..... ≤100 ppm	Sulphate..... ≤50 ppm	Na..... ≤500 ppm
pH sol. 5% in H2O..... 8.5 ÷ 9.6	Total nitrogen..... ≤10 ppm	Heavy metals (Pb)..... ≤5 ppm	Assay (potetiometric) ..... ≥98.0 %

Code	Size	Packaging	Notes
471785	100 g	Plastic bottle	
471786	500 g	Plastic bottle	
471787	1 kg	Plastic bottle	
471782	5 kg	Plastic tank	
471781	25 kg	Plastic drum	
471788	25 kg	Plastic bucket	

Potassium phosphate dibasic anhydrous > RE - Pure

RE

Description ..... white crystalline powder	Chloride..... ≤100 ppm	Sulphate..... ≤200 ppm	Assay (acidimetric) ..... ≥98 %
Identification ..... Positive	Heavy metals (Pb)..... ≤30 ppm	Fe ..... ≤50 ppm	

Code	Size	Packaging	Notes
361757	1 kg	Plastic bottle	
361752	5 kg	Plastic tank	
361751	25 kg	Plastic bucket	

Potassium phosphate dibasic trihydrate

• Potassio fosfato bibasico triidrato • Potassium phosphate dibasique trihydraté
• Potasio fosfato dibasico trihidrato • Zweiwertiges Kaliumphosphat-Trihydrat

Synonym:
Dipotassium hydrogen phosphate trihydrate

$K_2HPO_4.3H_2O$
Molecular Weight: 228,23
CAS: 16788-57-1
EEC-N: 231-834-5

Potassium phosphate dibasic trihydrate > RPE - For analysis

RPE

Description ..... White crystals	Total nitrogen..... ≤50 ppm	Cu..... ≤30 ppm	Pb..... ≤30 ppm
Identification ..... Positive	Chloride..... ≤50 ppm	Fe ..... ≤30 ppm	
pH sol. 5% in H2O..... 8.5 ÷ 9.6	Sulphate..... ≤100 ppm	Assay (non-aqueos medium)..... ≥ 99.0 %	
Water-insoluble matter ..... ≤100 ppm	As ..... ≤1 ppm	Ni..... ≤30 ppm	

Code	Size	Packaging	Notes
471764	100 g	Plastic bottle	
471766	500 g	Plastic bottle	
471767	1 kg	Plastic bottle	
471761	25 kg	Fibre drum	

Product specifications are subject to changes. Please visit our website for updates.

Potassium phosphate monobasic

• Potassio fosfato monobasico • Potassium phosphate monobasique • Potasio fosfato monobásico
• Kaliumphosphat monobasisch

Synonym:
• *Monopotassium phosphate*
• *Potassium dihydrogen phosphate*

KH_2PO_4
Molecular Weight: 136,09
CAS: 7778-77-0
EEC-N: 231-913-4

Potassium phosphate monobasic > RPE - For analysis - ISO - ACS - Reag.Ph.Eur. - Reag.USP

RPE

Description ..... White crystals
Identification ..... Positive
pH sol. 5% in H2O ..... 4.2 ÷ 4.5
Loss on drying (105°C) ..... ≤0.2 %
Water-insoluble matter ..... ≤100 ppm
Chloride ..... ≤10 ppm
Sulphate ..... ≤30 ppm
Heavy metals (Pb) ..... ≤10 ppm
Fe ..... ≤10 ppm
Na ..... ≤50 ppm
Assay (potetiometric) 99.0 ÷ 100.5 (s-s) %
Appearance of solution ..... Clear colourless
liquid Ph. Eur.
Reducing substances ..... Pass test Ph.Eur.
As ..... ≤ 2 ppm
Loss on drying 130°C ..... ≤ 2.0 %

Code	Size	Packaging	Notes
471685	100 g	Plastic bottle	
471686	500 g	Plastic bottle	
471687	1 kg	Plastic bottle	
471682	5 kg	Plastic jar	
471681	25 kg	Plastic bucket	

Potassium phosphate monobasic > ERBApharm - According to pharmacopoeia: NF

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Organic volatile impurities Conform USP-NF
Loss on drying ..... ≤1.0 %
Fluoride ..... ≤10 ppm
Heavy metals (Pb) ..... ≤20 ppm
Not soluble matter ..... ≤0.2 %
As ..... ≤3 ppm
Pb ..... ≤5 ppm
Assay ..... 98.0 ÷ 100.5 % s.s.
Origin (BSE/TSE) ..... Synthesis
Residual solvents (Current ICH) ..... Conform

Code	Size	Packaging	Notes
361507	1 kg	Plastic bottle	
361509	5 kg	Plastic tank	
361503	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Potassium phosphate monobasic 0.2 mol/l (0.2N)

• Potassio fosfato monobasico 0.2 mol/l (0.2N) • Potassium phosphate monobasique 0.2 mol/l (0.2N)
• Potasio fosfato monobásico 0.2 mol/l (0.2N) • Monobasisches Kaliumphosphat 0.2 mol/l (0.2N)

Synonym:
• *Monopotassium phosphate*
• *Potassium dihydrogen phosphate*

KH_2PO_4
Molecular Weight: 136,09
CAS: 7778-77-0

Potassium phosphate monobasic 0.2 mol/l (0.2N) > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611069601	1 l	Plastic bottle	Ref Ph.Eur 1069601

Potassium pyroantimonate acid

• Potassio piroantimoniato acido • Potassium pyroantimoniate acide • Potasio piroantimoniato acido
• Kaliumsäurepyroantimoniat

Synonym:
• *Potassium hexahydroxoantimonate(V)*
• *Potassium antimonate, hydrated*

$KSb(OH)_6$
Molecular Weight: 262,9
CAS: 12208-13-8
EEC-N: 235-387-7

Classification transport
ONU: 1549
Transport Hazard class: 6.1
Packing group III



Warning
H302-H332-H411
P261-P264-P271-P301+P312a-P304+P340-P501a

Potassium pyroantimonate acid > RPE - For analysis

RPE

Description ..... White to off white crystals powder
Identification ..... Positive
Assay (as K3BO3 on calc.mat.) ..... ≥ 99 %
Loss on ignition ..... ≤ 22 %

Code	Size	Packaging	Notes
473835	250 g	Plastic bottle	

Potassium pyroantimonate solution

• Potassio piroantimoniato soluzione • Potassium pyroantimoine solution • Potasio piroantimoniato solución • Kaliumpyroantimonlösung

Synonym:
• *Potassium hexahydroxoantimonate(V)*
• *Potassium antimonate, hydrated*

$KSb(OH)_6$
Molecular Weight: 262,9
CAS: 12208-13-8

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group II

Potassium pyroantimonate solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611071303	100 ml	Plastic bottle	Ref Ph.Eur 1071300
611071309	250 ml	Plastic bottle	Ref Ph.Eur 1071300
611071302	1 l	Plastic bottle	Ref Ph.Eur 1071300

Potassium pyroantimonate solution > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
617000311	100 ml	Bottle	

tetra-Potassium pyrophosphate

• Potassio pirofosfato • Potassium pyrophosphate • Potasio pirofosfato • Kaliumpyrophosphat

$K_4P_2O_7$
Molecular Weight: 330,35
CAS: 7320-34-5
EEC-N: 230-785-7

tetra-Potassium pyrophosphate > RPE - For analysis

RPE

Description	White powder	Fluoride	≤ 5 ppm	As	≤ 1 ppm	Assay (acidimetric)	≥ 95.0 %
Identification	Positive	Heavy metals (Pb)	≤ 20 ppm	Fe	≤ 30 ppm		
Water-insoluble matter	≤ 0.1 %	pH sol. 1% at 25° C	10.0 ÷ 10.5	Pb	≤ 1 ppm		

Code	Size	Packaging	Notes
473915	250 g	Plastic bottle	

Potassium pyrosulphate

• Potassio pirosolfato • Potassium pyrosulfate • Potasio pirosulfato • Kaliumpyrosulfat

Synonym:
Potassium disulfate

$K_2O_7S_2$
Molecular Weight: 254,33
CAS: 7790-62-7
EEC-N: 232-216-8

Classification transport
ONU: 3260
Transport Hazard class: 8
Packing group II

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Potassium pyrosulphate > RPE - For analysis - ACS

RPE

Description	White granules	Chloride	≤ 20 ppm	Fe	≤ 20 ppm	Ca	≤ 20 ppm
Identification	Positive	Phosphate	≤ 10 ppm	Na	≤ 100 ppm	Mg	≤ 10 ppm
Water (K.F.)	≤ 2.5 %	Heavy metals (Pb)	≤ 10 ppm	Assay (acidimetric)	37.5 ÷ 38.6 %(H2SO4)	Water-insoluble matter	≤ 100 ppm

Code	Size	Packaging	Notes
474016	100 g	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Potassium sodium tartrate tetrahydrate

• Potassio sodio tartrato tetraidrato • Potassium sodium tartrate tétrahydraté
• Potasio y sodio tartrato tetrahidrato • Kalium natrium tartrat tetrahydrate

Synonym:
Rochelle salt

$C_4H_4O_6KNa.4H_2O$
Molecular Weight: 282,23
CAS: 6381-59-5
EEC-N: 205-698-2

Potassium sodium tartrate tetrahydrate > RPE - For analysis - ISO - ACS - Reag.Ph.Eur. - Reag.USP

RPE

Description	White semitransparent crystals	Chloride	≤10 ppm	Sulphate	≤50 ppm	%	
Identification	Positive	Phosphate	≤20 ppm	Ca	≤50 ppm		
pH sol. 5% at 25° C	6.0 ÷ 8.5	Water-insoluble matter	≤50 ppm	Fe	≤10 ppm		
Ammonium	≤20 ppm	Heavy metals (Pb)	≤5 ppm	Assay (non-aqueos medium)	99.0 ÷ 102.0		

Code	Size	Packaging	Notes
474115	100 g	Plastic bottle	
474116	500 g	Plastic bottle	
474117	1 kg	Plastic bottle	
474119	5 kg	Plastic jar	
474112	25 kg	Plastic bucket	
474114	50 kg	Plastic bucket	

Potassium sodium tartrate tetrahydrate > ERBApharm - According to pharmacopoeia: USP

ERBApharm

Description	White crystalline powder	Ammonia	Conform USP-NF	Assay (alkalimetric)	99.0 ÷ 102.0 % s.s.
Identification	Positive	Water (K.F.)	21.0 ÷ 27.0 %	Origin (BSE/TSE)	Conform
Alcalinity	Conform USP-NF	Heavy metals (Pb)	≤10 ppm	Residual solvents (Current ICH)	Conform

Code	Size	Packaging	Notes
363457	1 kg	Plastic bottle	
363459	5 kg	Plastic tank	
363455	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Potassium sorbate

• Potassio sorbato • Potassium sorbate • Potasio sorbato • Kaliumsorbat

Synonym:
• *Potassium 2,4-hexadienoate*
• *Sorbic acid potassium salt*

$CH_3(CH:CH)_2COOK$
Molecular Weight: 150,22
CAS: 24634-61-5
EEC-N: 246-376-1

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Potassium sorbate > RE - Pure

RE

Description	White rode-like granules	Loss on drying in vacuo	≤1 %	As	≤3 ppm	Assay (non-aqueos medium)	≥99 % s s
Identification	Positive	Aldehydes(Formaldehyde)	≤0.1 %	Pb	≤10 ppm		
M.p. extr. Sorbic acid	133 ÷ 135 ° C	Cu + Zn	≤50 ppm	Zn	≤25 ppm		

Code	Size	Packaging	Notes
363884	1 kg	Plastic bottle	

Potassium sulfate

• Potassio solfato • Potassium sulfate • Potasio sulfato • Kaliumsulfat

K_2SO_4
Molecular Weight: 174,27
CAS: 7778-80-5
EEC-N: 231-915-5

Potassium sulfate > RS - For microanalysis

RS

Description ..... White crystals Identification ..... Positive

Code	Size	Packaging	Notes
474205	250 g	Plastic bottle	

Potassium sulfate > RPE - For analysis - ISO - ACS - Reag.Ph.Eur. - Reag.USP

RPE

Description ..... White crystals
Identification ..... Positive
pH sol. 5% in H20 ..... 5.5 ÷ 8.5
Total nitrogen ..... ≤5 ppm
Chloride ..... ≤10 ppm
Water-insoluble matter ..... ≤100 ppm
Heavy metals (Pb) ..... ≤5 ppm
Fe ..... ≤5 ppm
Na ..... ≤200 ppm
Assay (acidimetric) ..... ≥99.0 %
Ca ..... ≤100 ppm
Mg ..... ≤50 ppm

Code	Size	Packaging	Notes
474166	100 g	Plastic bottle	
474167	1 kg	Plastic bottle	
474169	5 kg	Plastic jar	

Potassium sulfate > RE - Pure

RE

Description ..... White crystalline powder
Identification ..... Positive
Loss on drying ..... ≤ 0.1 %
Assay (gravimetric) ..... ≥99 %
Na ..... ≤ 300 ppm
Chloride ..... ≤ 40 ppm
Mg ..... ≤ 300 ppm

Code	Size	Packaging	Notes
363607	1 kg	Plastic bottle	
363608	5 kg	Plastic tank	
363602	25 kg	Plastic bucket	

Potassium tartrate

• Potassio tartrato • Potassium tartrate • Potasio tartrato • Kaliumtartrat

Synonym:
Dipotassium tartrate hemihydrate

$C_4H_4K_2O_6 \cdot 1/2H_2O$
Molecular Weight: 235,28
CAS: 6100-19-2

Potassium tartrate > RPE - For analysis

RPE

Description ..... White crystals
Identification ..... Positive
pH sol. 5% at 25° C ..... 7.0 ÷ 9.0
Ammonium ..... ≤ 10 ppm
Chloride ..... ≤ 10 ppm
Phosphate ..... ≤ 5 ppm
Water-insoluble matter ..... ≤ 50 ppm
Heavy metals (Pb) ..... ≤ 5 ppm
Sulphate ..... ≤ 50 ppm
As ..... ≤ 0.4 ppm
Ca ..... ≤ 50 ppm
Cu ..... ≤ 2 ppm
Fe ..... ≤ 5 ppm
Na ..... ≤ 200 ppm
Ni ..... ≤ 2 ppm
Pb ..... ≤ 2 ppm
Zn ..... ≤ 2 ppm
Assay (non-aqueos medium) ..... ≥ 99 %

Code	Size	Packaging	Notes
474465	250 g	Plastic bottle	
474467	1 kg	Plastic bottle	

DASIT GROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

Potassium L-tartrate monobasic

• Potassio tartrato acido • Potassium L-tartrate monobasique • Potasio tartrato acido
• Kalium L-Tartrat monobasisch

Synonym:
Potassium bitartrate

$COOK(CHOH)_2COOH$
Molecular Weight: 188,18
CAS: 868-14-4
EEC-N: 212-769-1

Potassium L-tartrate monobasic > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Loss on drying ..... ≤ 0.2 %
Pot. rotat. spec. a 20°C (c=10; NaOH 1N) ..... +32 ÷ +33 °
Chloride ..... ≤ 200 ppm
Heavy metals (Pb) ..... ≤ 10 ppm
Sulphate ..... ≤ 100 ppm
As ..... ≤ 1 ppm
Assay (non-aqueos medium) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
474515	250 g	Plastic bottle	
474517	1 kg	Plastic bottle	

Potassium L-tartrate monobasic > RE - Pure

RE

Description ..... White crystalline powder
Identification ..... Positive
Loss on drying ..... ≤0.5 %
Ac. tartarico libero ..... ≤0.2 %
Chloride ..... ≤350 ppm
Sulphate ..... ≤0.12 %
As ..... ≤1 ppm
Pb ..... ≤10 ppm
Assay (acidimetric) ..... ≥99.5 % s.s.

Code	Size	Packaging	Notes
363907	1 kg	Plastic bottle	

Potassium tellurite

• Potassio tellurito • Potassium tellurite • Potasio telurito • Kaliumtellurit

$K_2TeO_3.nH_2O$
Molecular Weight: 253,8
CAS: 7790-58-1
EEC-N: 232-213-1



Danger
H301
P264-P270-P301+P310a-P330-P405-P501a

Potassium tellurite > RE - Pure

RE

Description ..... Yellowish crystals
Identification ..... Positive
Assay (oxidimetric) ..... ≥95 %

Code	Size	Packaging	Notes
474652	25 g	Glass bottle	

Potassium tetraiodomercurate solution, alkaline

• Potassio tetraiodomercurato soluzione alcalina • Potassium tétraiodomercurate - Solution alcaline • Potasio tetraiodomercurato solución alcalina
• Kaliumtetraiodomercurat - alkalische Lösung

Classification transport
ONU: 3316
Transport Hazard class: 9
Packing group II



Danger
H290-H300-H310-H330-H314-H373-H410
P264-P273-P301+P310a-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338-P320-P361+P364

Potassium tetraiodomercurate solution, alkaline > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611071600	200 ml	Bottle	Ref Ph.Eur 1071600

Potassium thiocyanate

• Potassio solfocianuro • Potassium sulfocyanure • Potasio sulfocianuro • Kaliumthiocyanat

Synonym:
Potassium rhodanide

KSCN
Molecular Weight: 97,18
CAS: 333-20-0
EEC-N: 206-370-1

Warning
H302-H312-H332-H412-HEU032
P261-P264-P271-P280h-P301+P312a-P304+P340

Potassium thiocyanate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description	White crystals	Water-insoluble matter	≤50 ppm	Chloride	≤50 ppm	Fe	≤2 ppm
Identification	Positive	Reducing iodine	Conform ACS	Sulphate	≤50 ppm	Na	≤50 ppm
pH sol. 5% at 25° C	5.3 ÷ 8.7	Ammonium	≤30 ppm	Heavy metals (Pb)	≤5 ppm	Assay (argentimetric)	≥99.0 %

Code	Size	Packaging	Notes
474355	250 g	Plastic bottle	
474357	1 kg	Plastic bottle	

Potassium thiocyanate > RE - Pure

RE

Description	White crystalline powder	pH sol. 5% at 20°C	5.0 ÷ 8.7	Sulphate	≤ 0.1 %	Assay (argentimetric)	≥ 98 %
Identification	Positive	Chloride	≤ 500 ppm	Fe	≤ 20 ppm		

Code	Size	Packaging	Notes
363756	500 g	Plastic bottle	
363752	25 kg	Plastic bucket	

Potassium thiocyanate solution 5%

• Potassio solfocianuro soluzione 5% • Potassium sulfocyanure solution 5%
• Potasio sulfocianuro solución 5% • Kaliumthiocyanat 5%

Synonym:
Potassium rhodanide

KSCN
Molecular Weight: 97,18
CAS: 333-20-0

HEU031-HEU210

Potassium thiocyanate solution 5% > RPE - For analysis

RPE

Description	Clear colourless liquid	Identification	Positive	Density at 20° C	1.020 ÷ 1.030

Code	Size	Packaging	Notes
E474381	1 l	Plastic bottle	

Stabilized with methyle p-hydroxybenzoate and n-Propyle p-hydroxybenzoate

Potassium thiocyanate 0.1 mol/l (0.1N)

• Potassio solfocianuro 0.1 mol/l (0.1N) • Potassium sulfocyanure 0.1 mol/l (0.1N)
• Potasio sulfocianuro 0.1 mol/l (0.1N) • Kaliumthiocyanat 0.1 mol/l (0.1N)

Synonym:
Potassium rhodanide

KSCN
Molecular Weight: 97,18
CAS: 333-20-0

HEU031-HEU210

Potassium thiocyanate 0.1 mol/l (0.1N) > RPE - For analysis

RPE

Description	Clear colourless liquid	Identification	Positive	Titration factor	0.998 ÷ 1.002

Code	Size	Packaging	Notes
E474417	1 l	Plastic bottle	

9.718 g of KSCN. Volumetric solution ready-to-use. Stabilized with p-oxybenzoate

Product specifications are subject to changes. Please visit our website for updates.

Potassium thiocyanate solution

• Potassio solfocianuro soluzione • Potassium sulfocyanure solution • Potasio sulfocianuro solución • Kaliumsulfocyanidlösung

Synonym:
Potassium rhodanide

KSCN
Molecular Weight: 97,18
CAS: 333-20-0

Potassium thiocyanate solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611071801	1 l	Plastic bottle	A 97 g/l solution Ref Ph.Eur 1071801

Praseodymium standard solution

• Praseodimio standard soluzione • Praséodyme solution standard • Praseodimio, solución patrón • Praseodym-Standardlösung

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Praseodymium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505782	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505785	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Praseodymium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503821	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503823	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503825	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503827	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Primary opalescent suspension

• Sospensione primaria di opalescenza • Suspension mère d'opalescence • Suspensión opalescente primaria • Suspensionsmutter der Opaleszenz

Danger
H317-H350-HA26
P261-P280-P308+P313-P362+P364-P333+P313-P501a

Primary opalescent suspension > RS - For analysis according to USP - Ph. Eur. Chap. 2.2.1

RS

Description ..... Opalescent liquid Assay (15ml in 1l of purif.water)..... 57 - 63 NTU

Code	Size	Packaging	Notes
612201100	100 ml	Glass bottle	Formazin suspension
612201101	1 l	Glass bottle	Formazin suspension

L(-)Proline

• L(-)Prolina • L(-)Proline • L(-)Prolina • L-Prolin

Synonym:
(S)-Pyrrolidine-2-carboxylic acid

$NH(CH_2)_3CHCOOH$
Molecular Weight: 115,13
CAS: 147-85-3
EEC-N: 205-702-2

L(-)Proline > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Potere rotat. spec. a 20°C (C=4; H2O) ..... -86.8 ÷ -84.5 °
pH 2,5% at 25° C ..... 5.5 ÷ 7
Loss on drying ..... ≤ 0.3 %
Assay (non-aqueos medium) ..... 99.0 ÷ 101.0 %

Code	Size	Packaging	Notes
474708	5 g	Glass bottle	

Propan-1-ol

• Propan-1-olo • Propane-1-ol • Propan-1-ol • 1-Propanol

Synonym:
Propyl alcohol

$CH_3CH_2CH_2OH$
Molecular Weight: 60,1
CAS: 71-23-8
EEC-N: 200-746-9

Classification transport
ONU: 1274
Transport Hazard class: 3
Packing group III



Danger
H225-H302-H318-H336
P210-P280-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Propan-1-ol > RS - For HPLC - Isocratic Grade

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.803 ÷ 0.805
Refractive index at 20°C ..... 1.3840 ÷ 1.3860
Boiling point ..... 96.9 ÷ 97.4 ° C
Acidity or alkalinity ..... ≤0.00015 meq/g
Water (K.F.) ..... ≤500 ppm
Residue on evaporation ..... ≤10 ppm
Assay (GLC) ..... ≥99.5 %
U.V. Transmittance
at 220 nm ..... ≥20 %
at 230 nm ..... ≥56 %
at 240 nm ..... ≥79 %
at 250 nm ..... ≥89 %
at 270 nm ..... ≥96 %
at 290 nm ..... ≥98 %

Code	Size	Packaging	Notes
412541000	1 l	Glass bottle	
412542000	2.5 l	Glass bottle	

Propan-1-ol > RS - Anhydrous - For analysis

RS

Refractive index at 20°C ..... 1.384 - 1.386
Water content (K.F.) ..... ≤ 300 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 99.5 %
Free acid (as CH3COOH) ..... ≤ 0.03 % m/m

Code	Size	Packaging	Notes
P0941016	1 l	Glass bottle	
P0941021	2.5 l	Glass bottle	

Propan-1-ol > RPE - For analysis - Reag. Ph. Eur.

RPE

Description ..... Clear liquid
Identification (I.R.) ..... Conform
Colour ..... ≤ 10 APHA
Density at 20° C ..... 0.803 ÷ 0.805
Refractive index at 20°C ..... 1.3840 ÷ 1.3860
Boiling point ..... 96 ÷ 98 ° C
Distillation range ..... 96 ÷ 99°C
Water (K.F.) ..... ≤ 1000 ppm
Residue on evaporation ..... ≤ 5 ppm
Acidity (acetic acid) ..... ≤ 0.03 %
Assay (GLC) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
415104	1 l	Glass bottle	
415102	2.5 l	Glass bottle	
415108	10 l	Plastic tank	
415106	25 l	Aluminium can	

Product specifications are subject to changes. Please visit our website for updates.

Propan-1-ol > RE - Pure

RE

Description ..... Clear liquid	Density at 20° C ..... 0.802 ÷ 0.806	Water (K.F.) ..... ≤0.5 %	Assay (GLC) ..... ≥99 %
Identification ..... Positive	Refractive index at 20°C 1.3830 ÷ 1.3870	Residue on evaporation ..... ≤ 50 ppm	Colour ..... ≤ 10 APHA

Code	Size	Packaging	Notes
309351	1 l	Glass bottle	
309352	2.5 l	Glass bottle	
309354	5 l	Plastic tank	
309353	25 l	Plastic tank	
309358	165 kg	Metal drum	

Propan-2-ol

• Propan-2-olo • Propane-2-ol • Propan-2-ol • 2-Propanol

Synonym:
• *Isopropanol*
• *IPA*

$CH_3CHOHCH_3$
Molecular Weight: 60,1
CAS: 67-63-0
EEC-N: 200-661-7

Classification transport
ONU: 1219
Transport Hazard class: 3
Packing group II

Danger
H225-H319-H336
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Propan-2-ol > RS - For LC/MS

RS

Description ..... Clear colourless liquid	Acidity (acetic acid) ..... ≤ 0.0010 %	At 260 nm ..... ≥ 98.5 %	Al ..... ≤ 50 ppb
Colour ..... ≤ 10 APHA	Alcalinity (NH3) ..... ≤ 0.0005 %	HPLC Gradient	Fe ..... ≤ 50 ppb
Identification (I.R.) ..... Positive	Assay (CPG) ..... ≥ 99.95 %	At 254 nm ..... ≤ 2 mAU	Na ..... ≤ 50 ppb
Refractive index at 20°C ..... 1.375 - 1.379	Transmittance	Test LC-MS TIC (50-2000m/z) ES I(+)	Ca ..... ≤ 50 ppb
Water (K.F.) ..... ≤ 200 ppm	At 220 nm ..... ≥ 64 %	Sensitive Impurities (reserpine) ... ≤ 100 ppb	Mg ..... ≤ 50 ppb
Residue on evaporation ..... ≤ 2 ppm	At 230 nm ..... ≥ 80 %	Metals compounds	K ..... ≤ 50 ppb

Code	Size	Packaging	Notes
415183	1 l	Glass bottle	
415184	2.5 l	Glass bottle	

Propan-2-ol > RS - For HPLC PLUS Gradient grade

RS

Description ..... Clear colourless liquid	Acidity or alkalinity ..... ≤0.0001 meq/g	at 254 nm ..... ≤2 ppb	at 240 nm ..... ≥89 %
Identification ..... Positive	Water (K.F.) ..... ≤0.1 %	at 365 nm ..... ≤2 ppb	at 250 nm ..... ≥96 %
Density at 20° C ..... 0.784 ÷ 0.786	Residue on evaporation ..... ≤2 ppm	U.V. Transmittance	at 260 nm ..... ≥98 %
Refractive index at 20°C 1.3766 ÷ 1.3786	Assay (GLC) ..... ≥99.9 %	at 220 nm ..... ≥63 %	at 270 nm ..... ≥99 %
Boiling point ..... 82.1 ÷ 82.6 ° C	Fluorescence	at 230 nm ..... ≥79 %	

Code	Size	Packaging	Notes
412711000	1 l	Glass bottle	
412712000	2.5 l	Glass bottle	

Propan-2-ol > RS - For HPLC Isocratic grade - ACS - Reag.Ph.Eur - Reag.USP

RS

Appearance ..... Clear liquid	Density at 20°C ..... 0.784 ÷ 0.786	Assay (GLC) ..... ≥ 99.9 %	at 275 nm ..... ≤ 0.03 AU
Identification ..... Positive	Refractive index at 20°C 1.3766 ÷ 1.3786	Absorbance UV (ACS - USP)	at 300 nm ..... ≤ 0.02 AU
Color ..... ≤ 10 APHA	Carbonyl comp. (propionald. and acetone) ≤ 20 ppm	at 210 nm ..... ≤ 1.00 AU	from 400 to 330 nm ..... ≤ 0.01 AU
Solubility in water ..... Passes test	Residue after evaporation ..... ≤ 5 ppm	at 220 nm ..... ≤ 0.40 AU	Transmittance UV (RS CLHP - Reag Ph Eur)
Miscibility in alcohol ..... Passes test	Water (H2O) ..... ≤ 500 ppm	at 230 nm ..... ≤ 0.20 AU	at 250 nm ..... ≥ 95 %
Miscibility in water ..... Passes test	Titrable acid or base ..... ≤ 0.0001 meq/g	at 245 nm ..... ≤ 0.08 AU	
Boiling point ..... 82.05 ÷ 82.55 °C		at 260 nm ..... ≤ 0.04 AU	

Code	Size	Packaging	Notes
412821	1 l	Glass bottle	
525161	2.5 l	Glass bottle	

Propan-2-ol > RS - For HPLC - Isocratic Grade

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.784 ÷ 0.786
Refractive index at 20°C . 1.3766 ÷ 1.3786
Boiling point ..... 82.1 ÷ 82.6 ° C
Acidity or alkalinity ..... ≤0.0001 meq/g
Water (K.F.) ..... ≤0.1 %
Residue on evaporation ..... ≤2 ppm
Assay (GLC) ..... ≥99.9 %
U.V. Transmittance
At 210 nm ..... ≥ 20 %
at 220 nm ..... ≥63 %
at 230 nm ..... ≥79 %
at 240 nm ..... ≥89 %
at 250 nm ..... ≥96 %
at 260 nm ..... ≥98 %
at 270 nm ..... ≥99 %

Code	Size	Packaging	Notes
412421000	1 l	Glass bottle	
412422000	2.5 l	Glass bottle	

Propan-2-ol > RS - For preparative HPLC

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.784 ÷ 0.786
Refractive index at 20°C . 1.3766 ÷ 1.3786
Boiling point ..... 82.1 ÷ 82.6 °C
Water (K.F.) ..... ≤500 ppm
Residue on evaporation ..... ≤5 ppm
Alcalinity (NH3) ..... ≤0.0002 meq/g
Assay (GLC) ..... ≥99.5 %
U.V. Transmittance
at 220 nm ..... ≥50 %
at 255 nm ..... ≥97 %

Code	Size	Packaging	Notes
415112	2.5 l	Glass bottle	

Propan-2-ol > RS - PESTIPUR - For pesticide analysis

RS

Clear,colourless liq.appearance ..... Conform
Identification ..... Conform
Refractive index at 20°C ..... 1.375 - 1.379
Colour ..... ≤ 10 Apha
Water content (K.F.) ..... ≤ 500 mg/Kg
Non volatile residue ..... ≤ 2 mg/Kg
Free acid (as CH3COOH) ..... ≤ 20 mg/Kg
Assay (GC) ..... ≥ 99.9 %
GC-ECD.Individual peak (Lindane) .. ≤ 3 ng/l
Retention time trichlorobenzene to mirex
GC-NPD.Individual peak (Ethylparathion) ≤ 3 ng/l
Retention time Atrazin to Coumaphos

Code	Size	Packaging	Notes
415281	1 l	Glass bottle	

For chlorinated and nitrogenous compounds analysis.

Propan-2-ol > RS - SPECTROSOL - For optical spectroscopy

RS

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Density at 20° C ..... 0.784 ÷ 0.786
Refractive index at 20°C . 1.3766 ÷ 1.3786
Boiling point ..... 82.1 ÷ 82.6 ° C
Water (K.F.) ..... ≤500 ppm
Residue on evaporation ..... ≤5 ppm
Acidity ..... ≤0.0005 meq/g
Alcalinity ..... ≤0.0002 meq/g
Assay (GLC) ..... ≥99.8 %
Fluorescence
at 254 nm ..... ≤2 ppb
at 365 nm ..... ≤2 ppb
U.V. Transmittance
at 205 nm ..... ≥10 %
at 215 nm ..... ≥50 %
at 230 nm ..... ≥80 %
at 250 nm ..... ≥95 %
at 260 nm ..... ≥98 %

Code	Size	Packaging	Notes
415213	1 l	Glass bottle	
P0952721	2.5 l	Glass bottle	

Propan-2-ol > RS - Anhydrous - For analysis

RS

Refractive index at 20°C ..... 1.375 - 1.379
Water content (K.F.) ..... ≤ 300 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 99.8 %
Free acid (as CH3COOH) ..... ≤ 10 mg/Kg

Code	Size	Packaging	Notes
P0951010	200 ml	Bottle with septum	
P0951016	1 l	Glass bottle	

Propan-2-ol > RS - VLSI - For electronic use

RS

Code	Size	Packaging	Notes
527696	1 l	Plastic bottle	
527690	2.5 l	Plastic bottle	
527691	30 l	Plastic tank	

For specifications, contact our customer service for a certificate of analysis

Product specifications are subject to changes. Please visit our website for updates.

Propan-2-ol > RS - RSE - For electronic use

RS

Description	White flakes	Chloride	≤0.2 ppm	Ca	≤0.2 ppm	Ni	≤0.01 ppm
Identification	Positive	Total phosphorus	≤0.1 ppm	Cd	≤0.01 ppm	Pb	≤0.01 ppm
Colour (APHA)	≤10	Heavy metals (Pb)	≤0.1 ppm	Co	≤0.01 ppm	Pt	≤0.05 ppm
Assay (HRGC)	≥99 %	Subst. reducing KMnO4	≤2.5 ppm	Cr	≤0.01 ppm	Sb	≤0.01 ppm
Water miscibility	Conform	Phosphate	≤0.5 ppm	Cu	≤0.01 ppm	Sn	≤0.02 ppm
Assay (GLC)	≥99.7 %	Total sulphur	≤1 ppm	Fe	≤0.1 ppm	Sr	≤0.02 ppm
Resistivity	≥10 Mohm.cm	Ag	≤0.02 ppm	Ga	≤0.02 ppm	Ti	≤0.05 ppm
Density at 20° C	0.784 ÷ 0.786	Al	≤ 50 ppb	In	≤0.02 ppm	Tl	≤0.05 ppm
Boiling point	82.1 ÷ 82.6 ° C	As	≤0.01 ppm	K	≤0.1 ppm	V	≤0.05 ppm
Water (K.F.)	≤500 ppm	Au	≤0.05 ppm	Li	≤0.02 ppm	Zn	≤ 10 ppb
Residue on evaporation	≤5 ppm	B	≤0.01 ppm	Mg	≤0.1 ppm	Zr	≤0.05 ppm
Acidity (propionic ac.)	≤10 ppm	Ba	≤0.1 ppm	Mn	≤0.01 ppm		
Alcalinity (NH3)	≤2 ppm	Be	≤0.02 ppm	Mo	≤0.05 ppm		
Aldeides - ketones	≤50 ppm	Bi	≤0.02 ppm	Na	≤0.5 ppm		

Code	Size	Packaging	Notes
415237	1 l	Glass bottle	
415235	2.5 l	Glass bottle	
415231	5 l	Metal tank	
415238	5 l	Plastic bottle	
415236	27 l	Metal drum	
415233	200 l	Metal drum	

Propan-2-ol > RS - MOS - For electronic use

RS

Description	Clear liquid	Chloride	≤0.2 ppm	Ca	≤0.2 ppm	Na	≤0.5 ppm
Colour (APHA)	≤10	Phosphate	≤0.5 ppm	Cd	≤0.01 ppm	Ni	≤0.01 ppm
Identification	Positive	Heavy metals (Pb)	≤0.1 ppm	Co	≤0.01 ppm	Pb	≤0.01 ppm
Water miscibility	Conform	Subst. reducing KMnO4	≤2.5 ppm	Cr	≤0.01 ppm	Pt	≤0.05 ppm
Resistivity	≥10 Mohm.cm	Total sulphur	≤1 ppm	Cu	≤0.01 ppm	Sb	≤0.01 ppm
Assay (GLC)	≥99.7 %	Ag	≤0.02 ppm	Fe	≤0.1 ppm	Sn	≤0.02 ppm
Density at 20° C	0.784 ÷ 0.786	Al	≤0.05 ppm	Ga	≤0.02 ppm	Sr	≤0.02 ppm
Boiling point	82.1 ÷ 82.6 ° C	As	≤0.01 ppm	In	≤0.02 ppm	Ti	≤0.05 ppm
Water (K.F.)	≤500 ppm	Au	≤0.05 ppm	K	≤0.1 ppm	Tl	≤0.05 ppm
Residue on evaporation	≤5 ppm	B	≤0.01 ppm	Li	≤0.02 ppm	V	≤0.05 ppm
Acidity (propionic ac.)	≤10 ppm	Ba	≤0.1 ppm	Mg	≤0.1 ppm	Zn	≤0.01 ppm
Alcalinity (NH3)	≤2 ppm	Be	≤0.02 ppm	Mn	≤0.01 ppm	Zr	≤0.05 ppm
Aldeides - ketones	≤50 ppm	Bi	≤0.02 ppm	Mo	≤0.05 ppm		

Code	Size	Packaging	Notes
415162	1 l	Glass bottle	
415161	2.5 l	Glass bottle	

Propan-2-ol > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description	Clear liquid	Residue on evaporation	≤ 10 ppm	B	≤ 0.02 ppm	Mg	≤ 0.1 ppm
Colour (APHA)	≤ 10	Acidity or alkalinity	≤ 0.0001 meq/g	Ba	≤ 0.5 ppm	Mn	≤ 0.01 ppm
Identification (I.R.)	Conform	Carbonyl compounds (acetone)	≤ 20 ppm	Ca	≤ 0.5 ppm	Ni	≤ 0.01 ppm
Water miscibility	Conform	Carbonyl compounds (propionaldehyde)	≤ 20 ppm	Cd	≤ 0.01 ppm	Pb	≤ 0.01 ppm
Density at 20° C	0.785 ÷ 0.789	Subst. reducing KMnO4	≤ 5 ppm	Co	≤ 0.01 ppm	Zn	≤ 0.01 ppm
Refractive index at 20°C	1.3766 ÷ 1.3786	Heavy metals (Pb)	≤ 1 ppm	Cr	≤ 0.01 ppm	Assay (GLC)	≥ 99.9 %
Boiling point	82.1 ÷ 82.6 ° C	Al	≤ 0.5 ppm	Cu	≤ 0.01 ppm		
Water (K.F.)	≤ 0.1 %			Fe	≤ 0.1 ppm		

Code	Size	Packaging	Notes
415154	1 l	Glass bottle	
524171	1 l	Spray bottle	6 units / box
415156	2.5 l	Glass bottle	
415158	2.5 l	Plastic bottle	
415173	5 l	Plastic bottle	
529174	5 l	Plastic tank	
415153	10 l	Plastic tank	
415157	25 l	Aluminium can	
524170	25 l	Plastic tank	
415152	200 l	Metal drum	

Propan-2-ol > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-Ph.Franc.-BP

ERBApharm

Description ..... Clear colourless liquid
Identification ..... Conform
Appearance of solution ..... Conform Ph. Eur.
Acidity or alcalinity (ml NaOH 0,01M). ≤ 0.6 ml
Acidity (ml NaOH 0,020N) ..... ≤ 0.70 ml
Peroxide ..... Conform Ph.Eur.
Density at 20°C ..... 0.785 - 0.789 Ph.Eur.
Density at 25°C ..... 0.783 - 0.787
Refractive index at 20°C ..... 1.376 - 1.378
Water (K.F.) ..... ≤ 0.1 %
Benzene ..... ≤ 2 ppm
Non volatile substances ..... ≤ 20 ppm
Limit of volatiles impurities Conform USP-NF
Assay (CPG) ..... ≥ 99.9 %
Related compounds ..... ≤ 0.3 %
Any single impurity ..... ≤ 0.1 %
Total impurities (GC) ..... ≤ 1.0 %
UV Absorbance
at 230 nm ..... ≤ 0.30 AU
at 250 nm ..... ≤ 0.10 AU
at 270 nm ..... ≤ 0.03 AU
at 290 nm ..... ≤ 0.02 AU
at 310 nm ..... ≤ 0.01 AU
UV Absorbance curve from 230 to 310 nm ..... Curve with no observable peaks or shoulders Ph.Eur.
Origin (BSE/TSE) ..... Synthesis
Residual solvents (Current ICH) ..... Conform

Code	Size	Packaging	Notes
309501	1 l	Glass bottle	
309505	2.5 l	Glass bottle	
529165	5 l	Plastic tank	
309506	10 l	Plastic tank	
309504	25 l	Metal drum	
309507	25 l	Plastic drum	
309500	200 l	Metal drum	
309509	200 l	Plastic drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Propan-2-ol > RE - Pure

RE

Description ..... Clear colourless liquid
Water (K.F.) ..... ≤ 2000 ppm
Residue on evaporation ..... ≤ 20 ppm
Refractive index at 20°C ..... 1.375 ÷ 1.379
Acidity (acetic acid) ..... ≤ 20 ppm
Colour ..... ≤ 10 APHA

Code	Size	Packaging	Notes
529093	5 l	Plastic tank	
529092	25 l	Plastic tank	
529091	200 l	Metal drum	

Propan-2-ol 70%

• Propan-2-olo 70% • Propane-2-ol 70% • Propan-2-ol 70% • 2-Propanol 70%

Synonym:
• *2-Propanol*
• *Isopropanol*

$CH_3CHOHCH_3$
Molecular Weight: 60,1
CAS: 67-63-0

Classification transport
ONU: 1219
Transport Hazard class: 3
Packing group II

Danger
H225-H319-H336
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Propan-2-ol 70% > RPE - For analysis - ACS

RPE

Description ..... Clear colourless liquid
Density at 20°C ..... 0.856 ÷ 0.862

Code	Size	Packaging	Notes
524161	25 l	Plastic tank	

Propan-2-ol 70% > ERBApharm - According to pharmacopoeia: Ph.Eur.- Microbiological tested

ERBApharm

Description ..... Clear colourless liquid
Density at 20°C ..... 0.856 - 0.862
Assay (alcohometric) at 20°C ..... 69 - 71 %(m/m)
Related substances (CPG) ..... ≤ 0.3 %
Absorbance ..... Conform Ph.Eur.
Benzene ..... ≤ 2 ppm
Origin (BSE/TSE) ..... Synthesis
Residual solvents (Current ICH) ..... Conform
Total aerobic microbial count (TAMC) .... ≤ 5 CFU/100ml
Total yeasts/mould count (TYMC) ..... ≤ 5 CFU/100ml
Test of specified micro-organisms
Enterobacteriaceae ..... Absent/100 ml
Staphylococcus aureus ..... Absent/100 ml
Pseudomonas aeruginosa ... Absent/100 ml

Code	Size	Packaging	Notes
524195	5 l	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Product specifications are subject to changes. Please visit our website for updates.

Propan-2-ol 70% > ERBApharm - Prepared from raw material according Ph.Eur

ERBApharm

Description ..... Clear colourless liquid Density at 20°C ..... 0.856 ÷ 0.862 Assay (alcohometric) at 20°C ..... 69 ÷ 71 %(m/m)

Code	Size	Packaging	Notes
524182	1 l	Plastic bottle	
524183	1 l	Spray bottle	6 units / box
524184	2.5 l	Plastic bottle	
524181	5 l	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Propanedioic acid ▶ Malonic acid

1,2-Propanediol ▶ Propylene glycol

1,3-Propanediol

• 1,3-Propandiolo • 1,3-Propanediol • 1,3-Propanodiol • 1,3-Propandiol

Synonym:
• *1,3-Dihydroxypropane*
• *Trimethylene glycol*

$CH_2OHCH_2CH_2OH$
Molecular Weight: 76,1
CAS: 504-63-2
EEC-N: 207-997-3

1,3-Propanediol > RE - Pure

RE

Appearance ..... Clear and viscous liquid Refractive index at 20°C ..... 1.438 - 1.442 Colour ..... ≤ 15 Hazen
Identification (IR) ..... Conform Water content (K.F.) ..... ≤ 1000 mg/Kg Assay (GC) ..... ≥ 99.7 %

Code	Size	Packaging	Notes
P8040216	1 l	Glass bottle	
P8040222	5 l	Plastic tank	
P8040268	190 l	Metal drum	

1,2,3-Propanetriol ▶ Glycerol (30°Bé)

1-Propanesulfonic acid sodium salt

• Acido 1-propansolfonico sale sodico • Acide 1-propanesulfonique sel sodique • Acido 1-propanosulfónico sal sódica • 1-Propansulfonsäure-Natriumsalz

$CH_3CH_2CH_2SO_3Na$
Molecular Weight: 146,14
CAS: 14533-63-2

1-Propanesulfonic acid sodium salt > RS - For ion pair chromatography

RS

Description ..... White crystalline powder Absorbance At 220 nm ..... ≤ 0.02 AU At 260 nm ..... ≤ 0.01 AU
Water (K.F.) ..... ≤ 2.0 % At 200 nm ..... ≤ 0.05 AU At 230 nm ..... ≤ 0.02 AU
Assay ..... ≥ 98.0 % At 210 nm ..... ≤ 0.03 AU At 250 nm ..... ≤ 0.01 AU

Code	Size	Packaging	Notes
405901	25 g	Glass bottle	
405902	100 g	Plastic bottle	

Propanoic acid ▶ Propionic acid

2-Propanone ▶ Acetone

Propionaldehyde

• Aldeide propionica • Aldéhyde propionique • Aldehído propiónico • Propionaldehyd

Synonym:
Propanal

CH_3CH_2CHO
Molecular Weight: 58,08
CAS: 123-38-6
EEC-N: 204-623-0

Classification transport
ONU: 1275
Transport Hazard class: 3
Packing group II



Danger
H225-H302-H315-H319-H335
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Propionaldehyde > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Refractive index at 20°C. 1.3610 ÷ 1.3660
Assay (GLC) ..... ≥ 96.0 %

Code	Size	Packaging	Notes
310504	100 ml	Glass bottle	

Propionic acid

• Acido propionico • Acide propionique • Acido propiónico • Propionsäure

Synonym:
• *Propanoic acid*
• *Propanyl acid*

CH_3CH_2COOH
Molecular Weight: 74,08
CAS: 79-09-4
EEC-N: 201-176-3

Classification transport
ONU: 3463
Transport Hazard class: 8
Packing group II



Danger
H226-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Propionic acid > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Water miscibility ..... Conform
Alcohol miscibility ..... Conform
Density at 20° C ..... 0.992 ÷ 0.994
Refractive index at 20°C. 1.3864 ÷ 1.3884
Boiling point ..... 140 ÷ 142 ° C
Water (K.F.) ..... ≤0.25 %
Chloride ..... ≤5 ppm
Heavy metals (Pb) ..... ≤5 ppm
Residue on evaporation ..... ≤50 ppm
Sulphate ..... ≤10 ppm
Fe ..... ≤2 ppm
Assay (acidimetric) ..... ≥99.5 %

Code	Size	Packaging	Notes
409551	250 ml	Glass bottle	
409553	1 l	Glass bottle	

Propionic acid > ERBAPharm - According to pharmacopoeia: USP-NF

ERBApharm

Description ..... Clear colourless liquid
Identification (I.R.) ..... Positive
Density at 25°C ..... 0.988 - 0.993
Distillation range ..... 138.5 - 142.5 °C
Nonvolatil residue ..... ≤ 0.01 %
Heavy metals (Pb) ..... ≤ 0.001 %
Readily oxidizable substances ..... Conform
Aldehyde ..... Conform
Assay (acidimetric) ..... 99.5 - 100.5 %
Origin (BSE/TSE) ..... Synthesis
Residual solvents (Current ICH) ..... Conform

Code	Size	Packaging	Notes
529050	20 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Propionic acid > RE - Pure

RE

Description ..... Yellow clear liquid
Identification ..... Positive
Density at 20° C ..... 0.988 ÷ 0.998
Refractive index at 20°C ..... 1.3854 ÷ 1.3894
Residue on evaporation ..... ≤ 100 ppm
Assay (acidimetric) ..... ≥98 %

Code	Size	Packaging	Notes
306254	1 l	Glass bottle	

n-Propyl acetate

• n-Propile acetato • n-Propyle acétate • n-Propil acetato • n-Propylacetat

$CH_3COO(CH_2)_2CH_3$
Molecular Weight: 102,13
CAS: 109-60-4
EEC-N: 203-686-1

Classification transport
ONU: 1276
Transport Hazard class: 3
Packing group II



Danger
H225-H319-H336-HEU066
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

n-Propyl acetate > RPE - For analysis

RPE

Description ..... Clear liquid
Identification ..... Positive
Colour ..... ≤ 15 Apha
Refractive index at 20°C ..... 1.3812 ÷ 1.3882
Water (K.F.) ..... ≤ 1000 ppm
Assay (GLC) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
474807	1 l	Glass bottle	

DASITGROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

n-Propyl alcohol ▶ Propan-1-ol

Propyl p-hydroxybenzoate

• Propile p-ossibenzoato • Propyle p-oxybenzoate • Propil p-hidroxibenzoato • Propyl-p-oxybenzoat

Synonym: *Propylparaben*

$HOC_6H_4COO(CH_2)_2CH_3$
Molecular Weight: 180,21
CAS: 94-13-3
EEC-N: 202-307-7

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Propyl p-hydroxybenzoate > ERBApharm - According to pharmacopoeia: DAB-BP-FU-NF-Ph.Eur.-Ph.Franc. ERBApharm

Description ..... White crystalline powder	Acidity .....Conform Ph.Eur.	Melting point.....96 ÷ 99 °C	Origin (BSE/TSE)..... Synthesis
Identification .....Positive	Related compounds.....Conform Ph.Eur.	Sulphated ash..... ≤ 0.1 %	
Appearance of solution .....Conform Ph.Eur.	Organic volatile impurities Conform USP-NF	Assay (saponification) .....98.0 ÷ 102.0 %	

Code	Size	Packaging	Notes
363953	50 g	Glass bottle	
363956	500 g	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Propylene carbonate

• Propilene carbonato • Propylène carbonate • Propileno carbonato • Propylencarbonat

Synonym: *1,2-Propanediol cyclic carbonate*

$OCH(CH_3)CH_2OCO$
Molecular Weight: 102,09
CAS: 108-32-7
EEC-N: 203-572-1

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Propylene carbonate > RPE - For analysis RPE

Description .....Clear colourless liquid	Density at 20° C .....1.200 ÷ 1.206	Boiling point.....241.0 ÷ 242.5 °C	Assay (acidimetric) ..... ≥99 %
Identification .....Positive	Refractive index at 20°C.1.4199 ÷ 1.4219	Residue on ignition.....≤100 ppm	

Code	Size	Packaging	Notes
474871	1 l	Glass bottle	

Propylene glycol

• Glicol propilenico • Propylène glycol • Propilenglicol • Propylenglycol

Synonym: *1,2-Propanediol*

$CH_2OHCHOHCH_3$
Molecular Weight: 76,09
CAS: 57-55-6
EEC-N: 200-338-0

Propylene glycol > RPE - For analysis RPE

Description .....Clear colourless liquid	Density at 20° C .....1.034 ÷ 1.038	Alkalinity (NH4OH) .....≤0.85 ppm	Residue on ignition.....≤70 ppm
Identification .....Positive	Refractive index at 20°C.1.4309 ÷ 1.4339	Chloride.....≤20 ppm	Sulphate.....≤20 ppm
Water miscibility.....Conform	Boiling point.....188.0 ÷ 190.0 ° C	Carbonyl Compounds (CO).....≤100 ppm	As .....≤2 ppm
Miscb. with Acetone ..... Complete	Water (K.F.) ..... ≤0.1 %	Heavy metals (Pb).....≤2 ppm	Fe .....≤2 ppm
Alcohol miscibility..... Complete	Acidity (acetic acid).....≤ 20 ppm	Peroxides (H2O2) .....≤5 ppm	Assay (GLC) ..... ≥99.5 %

Code	Size	Packaging	Notes
454054	1 l	Glass bottle	
454053	2.5 l	Glass bottle	
454052	30 kg	Plastic drum	

Propylene glycol > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP

ERBApharm

Description ..... Clear colourless liquid
Identification ..... Positive Ph. Eur.
Identification A (IR)..... Conform USP
Identification B:
Diethylene glycol ..... ≤ 0.10 %
Ethylene glycol ..... ≤ 0.10 %
Identification C .Same RT to standard by GC USP
Acidity ..... Conform Ph.Eur.
Reducing substances ..... Conform Ph.Eur.
Oxidizing substances ..... Conform Ph.Eur.
Density at 20° C .....1.035 - 1.040 Ph.Eur.
Density at 25°C .....1.035 - 1.037 USP
Boiling point..... 184 ÷ 189 °C
Refractive index at 20°C.....1.431 ÷ 1.433
Water (K.F.) ..... ≤0.2 %
Sulphated ash..... ≤70 ppm
Chloride..... ≤70 ppm
Heavy metals (Pb)..... ≤ 5 ppm (m/V)
Sulphate ..... ≤60 ppm
Assay (GLC) ..... ≥99.5 %

Code	Size	Packaging	Notes
346701	1 l	Glass bottle	
346703	2.5 l	Glass bottle	
346705	60 kg	Plastic tank	
346708	200 kg	Plastic drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Pumice stone

• Pomice • Pierre ponce • Piedra pómez • Bimsstein

CAS: 1332-09-8

Pumice stone > RPE - For analysis

RPE

Description ..... Granuli grigiastri
Identification ..... Positive

Code	Size	Packaging	Notes
469971	250 g	Plastic bottle	

Pyridine

• Piridina • Pyridine • Piridina • Pyridin

N:CHCH:CHCH:CH
CAS: 110-86-1
EEC-N: 203-809-9

Classification transport
ONU: 1282
Transport Hazard class: 3
Packing group II



Danger
H225-H302-H312-H332
P210-P241-P261-P280-P303+P361+P353-P304+P340

Pyridine > RS - Anhydrous - For analysis

RS

Water content (K.F.) ..... ≤ 200 mg/Kg
Colour ..... ≤ 20 Hazen
Assay (GC) ..... ≥ 99.8 %
Non volatile residue ..... ≤ 10 mg/Kg
Refractive index at 20°C ..... 1.508 - 1.512

Code	Size	Packaging	Notes
P0671010	200 ml	Bottle with septum	
P0671016	1 l	Glass bottle	
P0671046	1 l	Glass bottle PVC coated	
P0671021	2.5 l	Glass bottle	
P0671068	200 l	Metal drum	

Pyridine > RS - For peptide synthesis

RS

Water content (K.F.) ..... ≤ 100 mg/Kg
Colour ..... ≤ 10 Hazen
Kaiser Test ..... Conform
Assay (GC) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
P0673516	1 l	Glass bottle	
P0673521	2.5 l	Glass bottle	

Pyridine > RS - For potentiometry

RS

Water content (K.F.) ..... ≤ 500 mg/Kg
Colour ..... ≤ 10 Hazen
Tetrabutylamm.hydroxyde test ..... Conform
Assay (GC) ..... ≥ 99.5 %
Non volatile residue ..... ≤ 20 mg/Kg

Code	Size	Packaging	Notes
P06725P16	1 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Pyridine > RS - For titration according to Karl Fischer

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.979 ÷ 0.985
Refractive index at 20°C. 1.5050 ÷ 1.5140
Boiling point..... 114.2 ÷ 116.2 ° C
Water (K.F.) ..... ≤500 ppm
Assay (GLC) ..... ≥99.6 %

Code	Size	Packaging	Notes
469651	250 ml	Glass bottle	
469652	1 l	Glass bottle	

Pyridine > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Water solubility..... Conform
Boiling point..... 114.2 - 116.2 ° C
Subst. reducing KMnO4..... Conform
Water (K.F.) ..... ≤0.1 %
Residue on evaporation ..... ≤20 ppm
Ammonia..... ≤20 ppm
Chloride..... ≤10 ppm
Sulphate..... ≤10 ppm
Cu..... ≤5 ppm
Assay (GLC) ..... ≥99.0 %

Code	Size	Packaging	Notes
469622	500 ml	Glass bottle	
469629	1 l	Glass bottle	
469624	2.5 l	Glass bottle	
469626	20 kg	Plastic drum	
469621	25 kg	Metal drum	
469623	200 l	Metal drum	

Pyridine > RE - Pure

RE

Description ..... Clear colourless liquid or yellowish
Colour ..... ≤ 20 APHA
Identification ..... Positive
Density at 20° C ..... 0.979 ÷ 0.985
Refractive index at 20°C..... 1.5055 ÷ 1.5135
Boiling point..... 113.7 ÷ 116.7 °C
Water (K.F.) ..... ≤0.1 %
Residue on evaporation ..... ≤50 ppm
Assay (GLC) ..... ≥99.8 %

Code	Size	Packaging	Notes
358752	1 l	Glass bottle	
358754	25 kg	Metal drum	
528257	200 l	Metal drum	

Pyridine-d5

• Piridina-d5 • Pyridine-d5 • Piridina-d5 • Pyridin-d5

Synonym:
Pentadeuteropyridine

C_5D_5N
Molecular Weight: 84,14
CAS: 7291-22-7
EEC-N: 230-720-2

Classification transport
ONU: 1282
Transport Hazard class: 3
Packing group II

Danger
H225-H302-H312-H332
P210-P241-P261-P280-P303+P361+P353-P304+P340

Pyridine-d5 > RS - For NMR - min 99.95%

RS

Code	Size	Packaging	Notes
P5370	2 x 0.6 ml	Glass ampoule	

For specifications, contact our customer service for a certificate of analysis

Pyridine-d5 > RS - For NMR - min 99.5%

RS

Code	Size	Packaging	Notes
P5369A	2 x 0.75 ml	Glass ampoule	
P5364A	10 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

1-(2-Pyridylazo)-2-naphthol

• 1-(2-Piridile-azo)-2-naftolo • 1-(2-Pyridyl-azo)-2-napthol • 1-(2-Piridilazo)-2-naftol
• 1-(2-pyridyl-azo)-2-naphthol

Synonym:
PAN

N:CHCH:CHCH:CN:$NC_{10}H_6OH$
Molecular Weight: 249,27
CAS: 85-85-8
EEC-N: 201-637-9

1-(2-Pyridylazo)-2-naphthol > RPE - For analysis

RPE

Description ..... Orange red powder Identification ..... Positive Assay ..... ≥ 97.5 %

Code	Size	Packaging	Notes
469592	5 g	Glass bottle	

Complexometric indicator. For extraction and spectrophotometric determination of the transition metals

Pyrocatechol

• Pirocatecolo • Pyrocatéchol • Pirocatequina • Brenzkatechin

Synonym:
• 1,2-Benzenediol
• Catechol

$1,2-(OH)_2C_6H_4$
Molecular Weight: 110,11
CAS: 120-80-9
EEC-N: 204-427-5

Classification transport
ONU: 2811
Transport Hazard class: 6.1
Packing group III

Warning
H302-H312-H315-H319
P264-P280-P305+P351+P338-P332+P313-P362+P364-P337+P313

Pyrocatechol > RPE - For analysis

RPE

Description ..... Grey-brown flakes Identification ..... Positive Melting point ..... 103 ÷ 105 ° C Assay (GLC) ..... ≥ 99.0 %

Code	Size	Packaging	Notes
469753	50 g	Glass bottle	

Pyrocatechol violet

• Violetto pirocatechina • Violet de pyrocatéchol • Violeta de pirocatequina • Brenzcatechinviolett

Synonym:
Pyrocatecholsulfonphthalein

$C_{19}H_{14}O_7S$
Molecular Weight: 386,39
CAS: 115-41-3
EEC-N: 204-088-3

Pyrocatechol violet > RPE - For analysis

RPE

Description ..... Brown crystalline powder Identification ..... Positive Sensitivity as indicat. ..... Conform

Code	Size	Packaging	Notes
491871	1 g	Glass bottle	
491872	25 g	Glass bottle	

Complexometric indicator

Pyrogallol

• Pirogallolo • Acide pyrogallique • Pirogalol • Pyrogallol

Synonym:
1,2,3-Trihydroxybenzene

$1,2,3-(OH)_3C_6H_3$
Molecular Weight: 126,11
CAS: 87-66-1
EEC-N: 201-762-9

Warning
H302-H312-H332-H341-H412
P261-P271-P280-P304+P340-P308+P313-P330

Pyrogallol > RPE - For analysis - ACS - Reag. Ph.Eur.

RPE

Description ..... White crystalline powder Melting point ..... 131 ÷ 135 °C Fe ..... ≤10 ppm Sulphate ..... ≤50 ppm
Identification (I.R.) ..... Positive Sulphated ash ..... ≤ 0.005 % Heavy metals (Pb) ..... ≤5 ppm Chloride ..... ≤10 ppm

Code	Size	Packaging	Notes
409435	250 g	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Pyrrolidine dithiocarbamic acid ammonium salt

• Acido pirrolidinditiocarbammico sale ammonico • Acide pyrrolidinedithiocarboxylique-1,sel ammoniacal
• Acido pirrolidinditiocarbámico sal de amonio • Pyrrolidindithiocarbaminsäure-Ammoniumsalz

Synonym:
Ammonium pyrrolidinedithiocarbamate

$CH_2CH_2CH_2CH_2NCSSNH_4$
Molecular Weight: 164,29
CAS: 5108-96-3
EEC-N: 225-834-4

Pyrrolidine dithiocarbamic acid ammonium salt > RPE - For analysis - Reag. Ph. Eur.

RPE

Description ..... White yellowish powder Identification ..... Positive Assay ..... ≥ 99 % (NH3)

Code	Size	Packaging	Notes
409471	10 g	Glass bottle	

Quinaldine red

• Rosso chinaldina • Rouge de quinaldine • Rojo de quinaldina • Chinaldinrot

Synonym:
2-(4-Dimethylaminostyryl)-1-ethylquinolinium iodide

$C_{21}H_{23}IN_2$
Molecular Weight: 430,33
CAS: 117-92-0
EEC-N: 204-221-5

Quinaldine red > RPE - For analysis

RPE

Description .....Dark green powder
Identification ..... Positive
Non acq.media ind.sens. ..... Conform
Loss on drying .....≤5 %
Colour change.....Incolore - rosso

Code	Size	Packaging	Notes
476687	1 g	Glass bottle	
476688	25 g	Glass bottle	

Acid-base indicator (pH 1.4 ÷ 3.2)

Quinidine sulfate

• Chinidina solfato • Quinidine sulfate • Quinidina sulfato • Chinidinsulfat

$(C_{20}H_{24}N_2O_2)_2.H_2SO_4.2H_2O$
Molecular Weight: 782,96
CAS: 6591-63-5
EEC-N: 200-046-3



Warning
H302
P264-P270-P301+P312a-P330-P501a

Quinidine sulfate > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Ready carbonizable substances.....Conform
pH sol. 1% at 25° C .....6.0 ÷ 6.8
Specific optical rotation+275.5 ÷ +280.5 °
Loss on drying .....4.2 ÷ 5.0 %
Chloride.....≤50 ppm
Alchool/Chlorof.isolub. .....≤500 ppm
Heavy metals (Pb).....≤10 ppm
Residue on ignition.....≤500 ppm
Fe .....≤10 ppm
Assay (non-aqueos medium) ..... ≥99 %

Code	Size	Packaging	Notes
436701	10 g	Glass bottle	

Quinoline

• Chinolina • Quinoléine • Quinoleína • Chinolin

Synonym:
1-Benzazine

C_9H_7N
Molecular Weight: 129,16
CAS: 91-22-5
EEC-N: 202-051-6

Classification transport
ONU: 2656
Transport Hazard class: 6.1
Packing group III



Danger
H302-H312-H315-H319-H341-H350-H411-HA26
P280-P305+P351+P338-P308+P313-P330-P332+P313-P337+P313

Quinoline > RE - Pure

RE

Description .....Yellow-brown clear liquid
Identification ..... Positive
Density at 20° C .....1.088 ÷ 1.100
Refractive index at 20°C.1.6218 ÷ 1.6318
Boiling point.....235.6 ÷ 238.6 °C
Residue on ignition..... ≤0.1 %
Assay (GLC) ..... ≥95 %

Code	Size	Packaging	Notes
333701	100 ml	Glass bottle	
333707	1 l	Glass bottle	

Quinone ▶ p-Benzoquinone

Product specifications are subject to changes. Please visit our website for updates.

Raffinose
• Raffinosio • Raffinose • Rafinosa • Raffinose

$C_{18}H_{32}O_{16}.5H_2O$
Molecular Weight: 594,51
CAS: 17629-30-0
EEC-N: 208-146-9

Raffinose > RPE - For analysis — RPE

Description ..... White crystalline powder
Identification ..... Positive
Melting point..... 78 - 80 °C
Specific optical rotation..... +103 - +107 °
Water (K.F.) ..... 13 - 17 %
Residue on ignition..... ≤ 0.1 %
Total nitrogen ..... ≤100 ppm
Chloride..... ≤50 ppm
Heavy metals (Pb)..... ≤10 ppm
Sulphate..... ≤50 ppm
As ..... ≤2 ppm
Fe ..... ≤10 ppm

Code	Size	Packaging	Notes
475132	25 g	Glass bottle	

Raney's alloy
• Lega Raney • Alliage de Raney • Aleación según Raney • Raney-Legierung

CAS: 12003-78-0

Classification transport
ONU: 3089
Transport Hazard class: 4.1
Packing group III

Danger
H228-H260-H317-H351-H372
P210-P223-P231a+P232-P241-P280-P308+P313

Raney's alloy > RPE - For analysis — RPE

Description ..... Greyish metallic powder
Identification ..... Positive
Al..... ~50 %
Ni..... ~50 %

Code	Size	Packaging	Notes
457675	250 g	Plastic bottle	

Reagent for lipolysis
• Reattivo unico per lipolisi • Réactif unique pour lipolyse • Reactivo unico para lipolisi • Reagenz für die Lipolyse

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group III

Warning
H290
P234-P390-P406

Reagent for lipolysis > RPE - For analysis — RPE

Density at 20°C ..... 1.151 ÷ 1.161
pH at 20°C ..... 7.9 ÷ 8.3

Code	Size	Packaging	Notes
524910	2.5 l	Glass bottle	

Composition: Reagent for copper: 90.9% HCl 0,7N: 4.55% EDTA(8% w/v): 4.55%

Reagent TAN
• Reattivo TAN • Réactif TAN • Reactivo TAN • TAN-Reagenz

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II

Danger
H225-H315-H319-H361d-H336-H373
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Reagent TAN > RS - For analysis — RS

Water content (K.F.) ..... 4500 - 5500 mg/Kg
Free acid (as CH3COOH)..... ≤ 5 mg/Kg
Refractive index at 20°C..... 1.433 - 1.437

Code	Size	Packaging	Notes
PS0327/21	2.5 l	Glass bottle	
PS0327/29	5 l	Plastic tank	
PS0327/39	10 l	Plastic tank	

Composition: 495 ml propanol-2, 500ml Toluene, 5 ml water

Reagent TBN ASTM D2896

• Reattivo TBN ASTM D2896 • Réactif TBN ASTM D2896 • Reactivo TBN ASTM D2896 • TBN-Reagenz D2896

Classification transport
ONU: 2924
Transport Hazard class: 3
Packing group III



Danger
H226-H314-H411
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Reagent TBN ASTM D2896 > RS - For analysis

RS

Refractive index at 20°C ..... 1.464 - 1.469

Code	Size	Packaging	Notes
PS0423/21	2.5 l	Glass bottle	
PS0423/29	5 l	Plastic tank	
PS0423/39	10 l	Plastic tank	

Composition: 333 ml acetic acid, 667 ml chlorobenzene

Reagent TBN ASTM D4739

• Reattivo TBN ASTM D4739 • Réactif TBN ASTM D4739 • Reactivo TBN ASTM D4739 • TBN-Reagenz D4739

Classification transport
ONU: 1992
Transport Hazard class: 3
Packing group II



Danger
H225-H302-H315-H319-H351-H361d-H336-H372-HEU301
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Reagent TBN ASTM D4739 > RS - For analysis

RS

Aspect ..... Conform Density at 20°C ..... 1.036 ÷ 1.044 Water (K.F.) ..... ~ 1 % (m/m)

Code	Size	Packaging	Notes
526615	5 l	Plastic tank	

Composition: Bidistillated water: 0.05 L, Isopropyl alcohol: 1.650 L, Toluene: 1.650 L, Chloroform: 1.650 L

Red for oils O

• Rosso per olio O • Rouge pour l'huile O • Rojo para aceites O • Oelrot o

Synonym:
- *Oil Red O*
- *1-[2,5-Dimethyl-4-(2,5-dimethylphenylazo) phenylazo]-2-naphthol*

$C_{26}H_{24}N_4O$
Molecular Weight: 408,5
CAS: 1320-06-5
EEC-N: 215-295-3

Red for oils O > RPE - For analysis - C.I. 26125

RPE

Description ..... Red brick powder Identification ..... Positive

Code	Size	Packaging	Notes
476961	25 g	Glass bottle	

Dye for lipoprotein in acetate cellulose capsule

Redox solution 220 mV at 25°C

• Soluzione Redox 220 mV a 25°C • Solution réductice 220 mV à 25°C • Solución Redox 220 mV a 20°C • Reduktionslösung 220 mV bei 25 ° C

Redox solution 220 mV at 25°C > RS - For electrochemistry

RS

Description ..... Yellow clear liquid Identification ..... Positive pH at 25°C ..... 6.95 ÷ 7.05 Redox potential at 25°C ..... 215 ÷ 225 mV

Code	Size	Packaging	Notes
478032	500 ml	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Redox solution 468 mV at 25°C

• Soluzione Redox 468 mV a 25°C • Solution réductice 468 mV à 25°C • Solución Redox 468 mV a 25°C • Reduktionslösung 468 mV bei 25 ° C

Classification transport
ONU: 3264
Transport Hazard class: LQ

Warning
H290-H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Redox solution 468 mV at 25°C > RS - For electrochemistry

RS

Description .....Liquido limpido verde-marrone chiaro Identification ..... Positive Redox potential at 25°C..... 458 ÷ 478 mV

Code	Size	Packaging	Notes
478052	500 ml	Plastic bottle	

Refractive Index standards

• Indice di rifrazione standard • Etalons d'indice de réfraction • Patrones de índice de refracción • Brechungsindex-Standards

Refractive Index standards > RS - For calibration

RS

Code	Size	Packaging	Notes
540101	15 ml	Bottle	1.34325 at 20°C
540102	15 ml	Bottle	1.34782 at 20°C
540103	15 ml	Bottle	1.35171 at 20°C
540104	15 ml	Bottle	1.37233 at 20°C
540105	15 ml	Bottle	1.38115 at 20°C
540106	15 ml	Bottle	1.40978 at 20°C
540107	15 ml	Bottle	1.42009 at 20°C
540108	15 ml	Bottle	1.44193 at 20°C

Store between 2 - 8 °C

Reinecke salt

• Sale di Reinecke • Sel de Reinecke • Sal de Reinecke • Reinecke salt

Synonym:
Ammonium tetrarhodanatodiamminechromate(III)

$NH_4[Cr(NH_3)_2(SCN)_4].H_2O$
Molecular Weight: 354,44
CAS: 13573-16-5
EEC-N: 237-003-3

Warning
H302-H312-H332-HEU032
P261-P264-P271-P280h-P301+P312a-P304+P340

Reinecke salt > RPE - For analysis - ACS

RPE

Description ..... Dark red crystalline powder Diluted HCl-ins. matter .....≤500 ppm Assay (gravimetric) .....≥93.0 %
Identification ..... Positive Sens.(Choline chloride) .....≥0.5 mg/ml

Code	Size	Packaging	Notes
420672	25 g	Glass bottle	

Resorcinol

• Resorcina • Résorcine • Resorcina • Resorcin

Synonym:
1,3-Benzenediol

$1,3\text{-}(OH)_2C_6H_4$
Molecular Weight: 110,11
CAS: 108-46-3
EEC-N: 203-585-2

Classification transport
ONU: 2876
Transport Hazard class: 6.1
Packing group III

Warning
H302-H315-H319-H400
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Resorcinol > RPE - For analysis

RPE

Description .....White flakes Melting point..... 109 ÷ 111 ° C Assay (GLC) .....≥ 97.5 %
Identification ..... Positive Residue on ignition.....≤ 0.1 %

Code	Size	Packaging	Notes
476565	250 g	Plastic bottle	

L(+)Rhamnose

• L(+)Ramnosio • L(+)Rhamnose • L(+)Ramnosa • L-(+)-Rhamnose

Synonym:
6-Deoxy-L-mannose

$C_6H_{12}O_5.H_2O$
Molecular Weight: 182,17
CAS: 10030-85-0
EEC-N: 222-793-4

L(+)Rhamnose > RPE - For analysis

RPE

Appearance ..... White to pale cream crystalline powder
Assay (HPLC) ..... ≥ 98.5 %
Identification (FTIR) ..... Conforms
Water (K.F.) ..... 9 - 11 %
Melting point ..... 87 - 94 °C
Optical Rotation (C=10 in water, 20h) ..... +7 - +9 °

Code	Size	Packaging	Notes
476312	25 g	Glass bottle	

Rhenium standard solution

• Renio standard soluzione • Rhénium solution standard • Renio, solución patrón • Rhenium-Standardlösung

Classification transport
ONU: 3264

Rhenium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505802	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505805	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Rhenium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507754	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507513	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Rhodium standard solution

• Rodio standard soluzione • Rhodium solution standard • Rodio, solución patrón • Rhodium-Standardlösung

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group III

Warning
H290
P234-P390-P406

Rhodium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505807	100 ml	Plastic bottle	conc. 10 ppm Matrix: Hydrochloric acid
505808	100 ml	Plastic bottle	conc. 100 ppm Matrix: Hydrochloric acid
505809	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Product specifications are subject to changes. Please visit our website for updates.

Rhodium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503861	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503863	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503865	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Rhodol solution

• Rhodol soluzione • Solution Rhodol • Rhodol solución • Rhodol-Lösung

H412-HEU208
P273-P501a

Rhodol solution > RS - For analysis

RS

Refractive index at 20°C..... 1.342 - 1.346

Code	Size	Packaging	Notes
PS0158/15	1 l	Plastic bottle	
PS0158/29	5 l	Plastic tank	

Riboflavine

• Riboflavina • Riboflavine • Riboflavina • Riboflavin

Synonym:
Vitamin B2

$C_{17}H_{20}O_6N_4$
Molecular Weight: 376,37
CAS: 83-88-5
EEC-N: 201-507-1

Riboflavine > RE - Pure

RE

Description ..... Yellow-orange powder
Identification ..... Positive
Absorbance ..... Conform
Pot. rot. spec. a20°C (c=0.5;NaOH 0.05N) .....-115 ÷-135 °
s.s.
Loss on drying .....≤ 1.5 %
Sulphated ash.....≤ 0.1 %
Lumiflavine (TLC).....≤ 0.025 %
Lumiflavine (spectr.) .....≤ 0.025 A
Titolo (fluorimetrico).....98.0 ÷ 102.0 % s.s.

Code	Size	Packaging	Notes
389511	10 g	Glass bottle	

D(-)Ribose

• D(-)Ribosio • D(-)Ribose • D(-)Ribosa • D(-)-Ribose

$C_5H_{10}O_5$
Molecular Weight: 150,13
CAS: 50-69-1
EEC-N: 200-059-4

D(-)Ribose > RPE - For analysis

RPE

Description ..... Yellowish powder
Identification ..... Positive
Loss on drying ..... ≤ 1 %
Potere rotat. spec. a 20°C (C=4; H2O)-19.2
÷-20.8 °
Separazione (TLC) ..... ≥ 99.50 %
Heavy metals (Pb).....≤ 10 ppm
Residue on ignition..... ≤ 0.1 %

Code	Size	Packaging	Notes
476608	5 g	Glass bottle	

Rice starch
• Amido di riso • Amidon de riz • Almidón de arroz • Reisstärke

$(C_6H_{10}O_5)n$
CAS: 9005-25-8
EEC-N: 232-679-6

Rice starch > ERBApharm - According to pharmacopoeia: BP-FU-Ph.Eur.-Ph.Franc.

ERBApharm

Description	White powder	Loss on drying	≤ 15.0 %	Zolfo biossido	≤ 50 ppm	TYMC	≤ 100 CFU/g
Identification	Positive	Sulphated ash	≤ 0.6 %	Oxidizing substances	≤ 0.002 %	Escherichia coli	Absent Ph.Eur.
Microscopic test	Conform Ph.Eur.	pH (sosp. 20%)	5.0 ÷ 8.0	Microbial tests		Salmonella	Absent Ph. Eur.
Foreingn cellular elem.	Conform Ph.Eur.	Fe	≤ 10 ppm	TAMC	≤ 1000 CFU/g		

Code	Size	Packaging	Notes
313107	1 kg	Plastic bottle	
313108	2.5 kg	Plastic bottle	
313109	5 kg	Plastic tank	
313102	25 kg	Fibre drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Rosolic acid
• Acido p-rosolico • Acide p-rosolique • Acido p-rosólico • Rosolsäure

Synonym:
• *4-[Bis(4-hydroxyphenyl)methylene]-2,5-cyclohexadienone*
• *Aurin*

$C_{19}H_{14}O_3$
Molecular Weight: 290,32
CAS: 603-45-2
EEC-N: 210-041-8

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Rosolic acid > RPE - For analysis - C.I. 43800

RPE

Description	Red - brown crystalline powder	Loss on drying	≤ 10 %	pH range	6.2 ÷ 8.2
Identification	Positive	Colour change	yellow red		

Code	Size	Packaging	Notes
409702	25 g	Glass bottle	

Dye for microscopy (bacteriology). Indicator acid - base

Rubidium standard solution
• Rubidio standard soluzione • Rubidium solution standard • Rubidio, solución patrón • Rubidium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Rubidium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505792	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505795	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Rubidium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503841	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503843	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503845	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503847	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

DASIT GROUP | CARLO ERBA REAGENTS
Product specifications are subject to changes. Please visit our website for updates.

Rubidium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507755	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507514	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Ruthenium standard solution

• Rutenio standard soluzione • Ruthénium solution standard • Rutenio, solución patrón • Ruthenium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Ruthenium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505812	100 ml	Plastic bottle	conc. 10 ppm Matrix: Hydrochloric acid
505815	100 ml	Plastic bottle	conc. 100 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Ruthenium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503871	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503873	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503875	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid
503877	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Saccharose ▶ D(+)-Sucrose

Safranine T
• Safranina T • Safranine T • Safranina T • Safranin T

Synonym:
Basic Red 2

$C_{20}H_{19}ClN_4$
Molecular Weight: 350,85
CAS: 477-73-6
EEC-N: 207-518-8

Warning
H302-H319
P264-P280i-P301+P312a-P305+P351+P338-P337+P313-P501a

Safranine T > RS - For microscopy - C.I. 50420 — RS

Description ..... Red brown powder Identification ..... Positive

Code	Size	Packaging	Notes
477232	25 g	Glass bottle	

Dye for bacteriology, cytology

Safranine T hydroalcoholic solution for Gram-Hucker Kit
• Safranina T soluzione idroalcolica per kit Gram-Hucker
• Safranine T solution hydroalcoolique pour kit de Gram-Hucker
• Safranina T solución hidroalcohólica para kit Gram-Hucker
• Safranin T hydroalkoholische Lösung für Gram-Hucker Kit

Synonym:
Basic Red 2

$C_{20}H_{19}ClN_4$
Molecular Weight: 350,85
CAS: 477-73-6

HEU210

Safranine T hydroalcoholic solution for Gram-Hucker Kit > RS - For bacteriology — RS

Description ..... Red clear liquid Identification ..... Positive

Code	Size	Packaging	Notes
477241	250 ml	Glass bottle	In Vitro Diagnostic Medical Device

Salicylaldehyde azine
• Salicilaldeide azina • Salicylaldéhyde-azine • Salicilaldehido azina • Salicylaldehyd Azin

Salicylaldehyde azine > RS - For analysis according to Ph. Eur. Chap. 4.1.1 — RS

Code	Size	Packaging	Notes
611075500	100 ml	Glass bottle	Ref Ph.Eur 1075500

Salicylic acid
• Acido salicilico • Acide salicylique • Acido salicílico • Salicylsäure

Synonym:
2-Hydroxybenzoic acid

2-HOC_6H_4COOH
Molecular Weight: 138,12
CAS: 69-72-7
EEC-N: 200-712-3

Danger
H302-H315-H318-H335
P304+P340-P310a-P305+P351+P338-P330-P362+P364-P403+P233

Salicylic acid > RPE - For analysis - ACS — RPE

Description ..... White crystalline powder
Identification ..... Positive
Ready carbonizable substances..... Conform
Melting point..... 158.0 ÷ 161.0 °C
Chloride..... ≤10 ppm
Heavy metals (Pb)..... ≤5 ppm
Residue on ignition..... ≤100 ppm
Sulphate..... ≤30 ppm
Fe ..... ≤2 ppm
Assay (HPLC) ..... ≥99.0 %

Code	Size	Packaging	Notes
409773	100 g	Plastic bottle	
409777	1 kg	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Salicylic acid > ERBApharm - According to pharmacopoeia: FU

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .....Conform F.U.
Related substances (HPLC) .....Conform
Melting point..... 158.0 ÷ 161.0 °C
Loss on drying ..... ≤0.5 %
Sulphated ash..... ≤0.1 %
Chloride.....≤100 ppm
Sulphate.....≤200 ppm
Heavy metals (Pb).....≤20 ppm
Assay (acidimetric) .....99.0 ÷ 100.5 % s.s.
Origin (BSE/TSE)..... Synthesis
Residual solvents (Current ICH).....Conform

Code	Size	Packaging	Notes
306381	1 kg	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Salicylic acid > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Related compounds..... Conform PhEur
Melting point..... 158.0 ÷ 161.0 ° C
Loss on drying ..... ≤0.5 %
Sulphated ash.....≤500 ppm
Chloride.....≤100 ppm
Sulphate.....≤200 ppm
Heavy metals (Pb).....≤20 ppm
Assay (acidimetric) .....99.5 ÷ 100.5 % s.s.
Origin (BSE/TSE)..... Synthesis
Residual solvents (Current ICH).....Conform

Code	Size	Packaging	Notes
306377	1 kg	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Salmiac ▶ Ammonium chloride

Samarium standard solution

• Samario standard soluzione • Samarium solution standard • Samario, solución patrón • Samarium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Samarium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505852	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505855	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Samarium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503931	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503933	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503935	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503937	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Samarium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507756	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507515	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Sand of Fontainebleau

• Sabbia di Fontainebleau • Sable de Fontainebleau • Arena de Fontainebleau • Sand von Fontainebleau

Synonym:
Silica

SiO_2
Molecular Weight: 60,09
CAS: 14808-60-7
EEC-N: 238-878-4

Sand of Fontainebleau > RS - For agroalimentary analysis

RS

Density at 20°C ..... 2 ÷ 3 Granulometry ..... 180 ÷ 500 micrometer

Code	Size	Packaging	Notes
502064	1 kg	Plastic bottle	
502063	5 kg	Plastic bucket	
502062	25 kg	Plastic bucket	

Sand purified

• Sabbia purificata • Sable purifié • Arena purificada • Sand gereinigt

Synonym:
Silica

SiO_2
Molecular Weight: 60,09
CAS: 14808-60-7
EEC-N: 238-878-4

Sand purified > RS - For flash chromatography

RS

Description ..... Hazel granules Identification ..... Positive Particle size (40ö100) ..... Conform mesh

Code	Size	Packaging	Notes
477153	1 kg	Plastic bottle	

Saponin

• Saponina • Saponine • Saponina • Saponin

CAS: 8047-15-2
EEC-N: 232-462-6

Warning
H335
P261-P271-P304+P340-P312a-P403+P233-P501a

Saponin > RE - Pure

RE

Description ..... Yellow powder Identification ..... Positive pH 1% at 25°C ..... 5.0 - 6.5 Water ..... ≤ 5 %

Code	Size	Packaging	Notes
365755	250 g	Plastic bottle	
365757	1 kg	Plastic bottle	
365758	5 kg	Plastic tank	

Saybolt Colour Standards

• Standard del colore Saybolt • Etalons couleurs Saybolt • Patrones de color Saybolt • Saybolt Farbstandards

Saybolt Colour Standards > RS - For calibration

RS

Code	Size	Packaging	Notes
540709	100 ml	Glass bottle	-15
540710	100 ml	Glass bottle	+0
540711	100 ml	Glass bottle	+12
540712	100 ml	Glass bottle	+15
540713	100 ml	Glass bottle	+19
540714	100 ml	Glass bottle	+25
540715	100 ml	Glass bottle	+30

Product specifications are subject to changes. Please visit our website for updates.

Scandium standard solution

• Scandio standard soluzione • Scandium solution standard • Escandio, solución patrón • Scandium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III



Scandium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505837	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505838	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505839	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Scandium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503901	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503903	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503905	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503907	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Scandium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507757	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507516	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Schiff's reagent for Aldehydes

• Schiff reattivo soluzione per Aldedi • Réactif de Schiff pour Aldéhydes • Schiff reactivo solución para Aldehídos • Schiff-Reagens für Aldehyde



Danger
H314
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Schiff's reagent for Aldehydes > RS - For analysis

RS

Description ..... Clear colouless or light yellow liquid Identification ..... Positive

Code	Size	Packaging	Notes
477601	500 ml	Glass bottle	

Schiff's reagent for PAS coloration

• Schiff reattivo per colorazione PAS • Réactif de Schiff pour PAS coloration • Schiff reactivo solución para PAS coloratión • Schiff-Reagenz für NICHT Färbung



Warning
H290-H319
P234-P264-P280i-P305+P351+P338-P337+P313-P406

Schiff's reagent for PAS coloration > RS - For histology

RS

Description ..... Yellow clear liquid
Identification ..... Positive

Code	Size	Packaging	Notes
477591	500 ml	Glass bottle	In Vitro Diagnostic Medical Device
477592	6 x 500 ml	Glass bottle	In Vitro Diagnostic Medical Device

Sebacic acid

• Acido sebacico • Acide sébacique • Acido sebácico • Sebacinsäure

Synonym:
Decanedioic acid

$HOOC(CH_2)_8COOH$
Molecular Weight: 202,25
CAS: 111-20-6
EEC-N: 203-845-5

Sebacic acid > RE - Pure

RE

Description ..... White granular powder
Identification ..... Positive
Melting point ..... 132.5 ÷ 136.5 °C
Residue on ignition ..... ≤0.1 %
Assay (GLC) ..... ≥94 %

Code	Size	Packaging	Notes
409875	250 g	Plastic bottle	

Selenic mixture

• Miscela selenica • Mélange sélénique • Mezcla selenica • Selenmischung



Warning
H319-H411
P264-P280i-P305+P351+P338-P337+P313-P391-P501a

Selenic mixture > RS - For nitrogen detection according to Wieninger

RS

Description ..... Pads or dark gray powder
Identification ..... Positive

Code	Size	Packaging	Notes
463421	250 g	Plastic bottle	
463422	1 kg	Plastic bottle	

Selenium, powder

• Selenio, polvere • Sélénium, poudre • Selenio, polvo • Selen

Se
Molecular Weight: 78,96
CAS: 7782-49-2
EEC-N: 231-957-4

Classification transport
ONU: 3288
Transport Hazard class: 6.1
Packing group II



Danger
H301-H331-H373-H413
P271-P301+P310a-P304+P340-P311a-P330-P403+P233

Selenium, powder > RPE - For analysis

RPE

Description ..... Blackish powder
Identification ..... Positive
As ..... ≤ 100 ppm
Cu ..... ≤ 100 ppm
Fe ..... ≤ 100 ppm
Hg ..... ≤ 100 ppm
Pb ..... ≤ 500 ppm
Te ..... ≤ 500 ppm
Assay ..... ≥ 99.50 % (Se)

Code	Size	Packaging	Notes
477702	25 g	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Selenium standard solution

• Selenio standard soluzione • Sélénium solution standard • Selenio, solución patrón • Selen-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group II



Danger
H290-H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Selenium standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615002501	100 ml	Plastic bottle	A 1 ppm solution Ref Ph.Eur 5002501
615002500	1 l	Plastic bottle	A 100 ppm solution Ref Ph.Eur 5002500

Selenium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505842	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505845	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505843	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Selenium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503911	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503913	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503915	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503917	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Selenium standard solution > RS - Standard solution for AAS

RS

Description ..... Clear pinky liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
507758	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507491	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
E497625	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
E497621	500 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Selenium standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear pinkish liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
477691		Plastic ampoule	conc. 1.000 ppm Matrix: Water - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Selenium dioxide
• Selenio biossido • Sélenium dioxyde • Selenio dióxido • Selendioxid

SeO_2
Molecular Weight: 110,96
CAS: 7446-08-4
EEC-N: 231-194-7

Classification transport
ONU: 3283
Transport Hazard class: 6.1
Packing group II



Danger
H301-H331-H373-H410
P271-P301+P310a-P304+P340-P311a-P330-P403+P233

Selenium dioxide > RPE - For analysis

RPE

Description ..... White-pink crystals
Identification ..... Positive
Chloride ..... ≤ 500 ppm
Fe ..... ≤ 50 ppm
Pb ..... ≤ 50 ppm
Assay (oxidimetric) ..... ≥ 98 %

Code	Size	Packaging	Notes
477762	25 g	Glass bottle	

Selenous acid
• Acido selenioso • Acide sélénieux • Acido selenioso • Selenige Säure

H_2SeO_3
Molecular Weight: 128,98
CAS: 7783-00-8
EEC-N: 231-974-7

Classification transport
ONU: 3283
Transport Hazard class: 6.1
Packing group II



Danger
H301-H331-H373-H410
P271-P301+P310a-P304+P340-P311a-P330-P403+P233

Selenous acid > RPE - For analysis

RPE

Description ..... Whitish powder
Identification ..... Positive
Assay (iodometric) ..... ≥ 97.5 %

Code	Size	Packaging	Notes
409964	100 g	Glass bottle	

Silica gel 60A 6 - 35µ
• Gel di silice 60A 6 - 35µ • Gel de silice 60A 6 - 35µ • Gel de sílice 60A 6 - 35µ • Kieselgel 60A 6 - 35µ

SiO_2
Molecular Weight: 60,09
CAS: 7631-86-9

Silica gel 60A 6 - 35µ > RS - For chromatography

RS

Particle size ..... 6 - 35 µm
Loss on drying ..... ≤ 6 %
Particle size > 10.1 µm ..... ≥ 90 %
Particle size > 40.3 µm ..... ≤ 10 %
Particle size > 6.4 µm ..... ≥ 98 %

Code	Size	Packaging	Notes
P2010017	1 kg	Plastic bottle	
P2010027	5 kg	Plastic bucket	
P2010044	25 kg	Plastic drum	

400-2500 mesh. Stir before use

Silica gel 60A 20 - 45µ
• Gel di silice 60A 20 - 45µ • Gel de silice 60A 20 - 45µ • Gel de sílice 60A 20 - 45µ • Kieselgel 60A 20 - 45µ

SiO_2
Molecular Weight: 60,09
CAS: 7631-86-9

Silica gel 60A 20 - 45µ > RS - For chromatography

RS

Particle size ..... 20 - 45 µm
Loss on drying ..... ≤ 6 %
Particle size > 16 µm ..... ≥ 98 %
Particle size > 20.2 µm ..... ≥ 90 %
Particle size > 64.0 µm ..... ≤ 4 %

Code	Size	Packaging	Notes
P2200017	1 kg	Plastic bottle	
P2200027	5 kg	Plastic bucket	

325 - 625 mesh. Stir before use

Product specifications are subject to changes. Please visit our website for updates.

Silica gel 60A 35 - 70µ

• Gel di silice 60A 35 - 70µ • Gel de silice 60A 35 - 70µ • Gel de sílice 60A 35 - 70µ • Kieselgel 60A 35 - 70µ

SiO_2
Molecular Weight: 60,09
CAS: 7631-86-9

Silica gel 60A 35 - 70µ > RS - For chromatography

RS

Particle size ..... 35 ÷ 70 µm
Particle size <40 µm ..... ≤ 15 %
Particle size >63 µm ..... ≤ 10 %
Loss on drying ..... ≤ 8.0 %
pH (5%m/m solution) ..... 6.5 ÷ 7.5

Code	Size	Packaging	Notes
P2000017	1 kg	Plastic bottle	
P2000026	2 kg	Plastic bottle	
P2000027	5 kg	Plastic bucket	
P2000044	25 kg	Plastic drum	

220 - 400 mesh. Stir before use

Silica gel 60A 35 - 70µ > RS - For Flash chromatography

RS

Description ..... White powder
Identification ..... Positive
pH sospension 10% H20 ..... 6.2 ÷ 7.2
Chloride ..... ≤100 ppm
Apparent density ..... 380 ÷ 420 g/l
Nitrate ..... ≤150 ppm
Sulphate ..... ≤700 ppm
Cd ..... ≤5 ppm
Cu ..... ≤5 ppm
Fe ..... ≤50 ppm
Ni ..... ≤5 ppm
Pb ..... ≤5 ppm
Zn ..... ≤5 ppm
Granulometry
> 20.2 µm ..... ≥ 99.0 %
> 32.0 µm ..... ≥ 87.0 %
> 80.7 µm ..... ≤ 3.0 %

Code	Size	Packaging	Notes
453351	100 g	Plastic bottle	
453352	500 g	Plastic bottle	
453353	1 kg	Plastic bottle	
453355	5 kg	Plastic tank	

220 - 400 mesh

Silica gel 60A 40 - 63µ

• Gel di silice 60A 40 - 63µ • Gel de silice 60A 40 - 63µ • Gel de sílice 60A 40 - 63µ • Kieselgel 60A 40 - 63µ

SiO_2
Molecular Weight: 60,09
CAS: 7631-86-9

Silica gel 60A 40 - 63µ > RS - For chromatography

RS

Particle size ..... 40 ÷ 63 µm
Particle size <40 µm ..... ≤ 15 %
Particle size >63 µm ..... ≤ 10 %
Loss on drying ..... ≤ 6 %
pH (5%m/m solution) ..... 6.5 ÷ 7.5

Code	Size	Packaging	Notes
P2050017	1 kg	Plastic bottle	
P2050027	5 kg	Plastic bucket	
P2050044	25 kg	Metal bucket	

230 - 400 mesh. Stir before use

Silica gel 60A 70 - 200µ

• Gel di silice 60A 70 - 200µ • Gel de silice 60A 70 - 200µ • Gel de sílice 60A 70 - 200µ • Kieselgel 60A 70 - 200µ

SiO_2
Molecular Weight: 60,09
CAS: 7631-86-9

Silica gel 60A 70 - 200µ > RS - For chromatography

RS

Particle size ..... 70 ÷ 200 µm
Particle size <60 µm ..... ≤ 10 %
Particle size >200 µm ..... ≤ 10 %
Loss on drying ..... ≤ 6.0 %
pH (5%m/m solution) ..... 6.5 - 7.5

Code	Size	Packaging	Notes
P2100017	1 kg	Plastic bottle	
P2100026	2 kg	Plastic bottle	
P2100027	5 kg	Plastic bucket	
P2100044	25 kg	Plastic drum	

70 -220 mesh. Stir before use

Silica gel 60A 0,06÷0,20 mm

• Gel di silice 60A 0,06÷0,20 mm • Gel de silice 60A 0,06÷0,20 mm • Gel de sílice 60A 0,06÷0,20 mm • Kieselgel 60A 0.06 ÷ 0.20 mm

Synonym:
• *Silica*
• *Silicon dioxide*

SiO_2
Molecular Weight: 60,09
CAS: 7631-86-9

Silica gel 60A 0,06÷0,20 mm > RS - For chromatography

RS

Description ..... White powder
Identification ..... Positive
pH sospension 10% H2O ..... 6.5 ÷ 7.5
< 0.06 mm ..... ≤5 %
> 0.20 mm ..... ≤5 %

Code	Size	Packaging	Notes
453336	500 g	Plastic bottle	
453337	1 kg	Plastic bottle	
453332	5 kg	Plastic bucket	
453331	20 kg	Plastic bucket	

70 - 230 mesh

Silica gel granular

• Gel di silice granulare • Gel de silice granulés • Gel de sílice granulado • Kieselgel granuliert

Synonym:
• *Silica*
• *Silicon dioxide*

SiO_2
Molecular Weight: 60,09
CAS: 7631-86-9

Silica gel granular > RE - Pure

RE

Description ..... White granules
Identification ..... Positive
Functionality ..... Conform

Code	Size	Packaging	Notes
453278	10 x 50 g	Carton box	
453272	10 x 100 g	Carton box	
453273	10 x 250 g	Carton box	
453275	10 x 500 g	Carton box	
453277	1 kg	Plastic bottle	
453279	5 kg	Plastic tank	

Product specifications are subject to changes. Please visit our website for updates.

Silica gel granular with indicator cobalt free

• Gel di silice granulare con indicatore esente da cobalto
• Gel de silice granulés avec indicateur exempt de cobalt • Gel de sílice granulado sin cobalto
• Kieselgelgranulat mit kobaltfreiem Indikator

Synonym:
• *Silica*
• *Silicon dioxide*

SiO_2
Molecular Weight: 60,09
CAS: 7631-86-9

Silica gel granular with indicator cobalt free > RE - Pure

RE

Description ..... Small bags containing yellowish granules Identification ..... Positive Functionality ..... Conform

Code	Size	Packaging	Notes
453317	1 kg	Plastic bottle	
453319	5 kg	Plastic tank	
453315	25 kg	Plastic bucket	

Silicon standard solution

• Silicio standard soluzione • Silicium solution standard • Silicio, solución patrón • Silizium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Silicon standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505847	100 ml	Plastic bottle	conc. 10 ppm Matrix: Water
505848	100 ml	Plastic bottle	conc. 100 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Silicon standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503921	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503923	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503925	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Water
503927	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Water
504271	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid
504273	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid
504275	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid and hydrofluoric acid
504277	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid and hydrofluoric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Silicon standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
E497635	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Water
E497631	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Silicon standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
477961		Plastic ampoule	conc. 1.000 ppm Matrix: Water - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Silicon carbide ▶ Carborundum, granules

Silicon dioxide

• Silicio biossido • Silicium dioxyde • Silicio dióxido • Siliciumdioxid

Synonym:
Silica

SiO_2
Molecular Weight: 60,09
CAS: 14808-60-7
EEC-N: 238-878-4

Warning
H373
P260-P314-P501a

Silicon dioxide > RPE - For analysis

RPE

Description ..... White powder Identification ..... Positive Loss on ignition..... ≤ 0.5 % Assay ..... ≥ 99.5 %

Code	Size	Packaging	Notes
422104	100 g	Plastic bottle	
422106	500 g	Plastic bottle	

Silicotungstic acid

• Acido silicotungstico • Acide silicotungstique • Acido silicotúngstico • Kieselwolframsäure

Synonym:
Tungstosilicic acid hydrate

$SiO_2.12WO_3.26H_2O$
Molecular Weight: 3310,66
CAS: 12027-43-9

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Silicotungstic acid > RPE - For analysis

RPE

Description ..... White powder Identification ..... Positive Solubility..... Clear colourless liquid Loss on ignition..... 14.0 - 16.2 %

Code	Size	Packaging	Notes
410051	10 g	Glass bottle	

Silver, sheet

• Argento, lamina • Argent, lames • Plata, hojas • Silber, Blatt

Ag
Molecular Weight: 107,87
CAS: 7440-22-4
EEC-N: 231-131-3

Silver, sheet > RPE - For analysis

RPE

Description ..... Lamina Identification ..... Positive Assay ..... ≥99.9 %

Code	Size	Packaging	Notes
423752	25 g	Bag	0.5 mm

Product specifications are subject to changes. Please visit our website for updates.

Silver, wool

• Argento, lana • Argent, laine • Plata, lana • Silber, Wolle

Ag
Molecular Weight: 107,87
CAS: 7440-22-4
EEC-N: 231-131-3

Silver, wool > RS - For microanalysis

RS

Description ..... Lana Identification ..... Positive Assay ..... ≥99.9 %

Code	Size	Packaging	Notes
423791	5 g	Glass bottle	

Silver standard solution

• Argento standard soluzione • Argent solution standard • Plata, solución patrón • Silberstandardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Silver standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615002609	100 ml	Glass bottle	A 5 ppm solution: to dilute according to Ref Ph.Eur 5002600

Silver standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505302	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505305	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505303	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Silver standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503401	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503403	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503405	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503407	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Silver standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
507526	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507480	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Silver standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
423611		Plastic ampoule	conc. 1.000 ppm Matrix: Nitric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Silver acetate

• Argento acetato • Argent acétate • Plata acetato • Silberacetat

Synonym:
Acetic acid silver salt

CH_3COOAg
Molecular Weight: 166,92
CAS: 563-63-3
EEC-N: 209-254-9

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Silver acetate > RE - Pure

RE

Description ..... White greyish powder Identification ..... Positive Assay (argentimetric)..... ≥98 %

Code	Size	Packaging	Notes
319502	25 g	Glass bottle	
319507	250 g	Plastic bottle	

Silver carbonate

• Argento carbonato • Argent carbonate • Plata carbonato • Silbercarbonat

Ag_2CO_3
Molecular Weight: 275,75
CAS: 534-16-7
EEC-N: 208-590-3

Danger
H318-H410
P273-P280i-P305+P351+P338-P310a-P391-P501a

Silver carbonate > RE - Pure

RE

Description ..... Yellow-green powder Identification ..... Positive Substances not ppt HCl ..... ≤1 % Assay (argentimetric)..... 99.5 ÷ 101.5 %

Code	Size	Packaging	Notes
320002	25 g	Glass bottle	
320007	250 g	Glass bottle	

Silver chloride

• Argento cloruro • Argent chlorure • Plata cloruro • Silberchlorid

AgCl
Molecular Weight: 143,32
CAS: 7783-90-6
EEC-N: 232-033-3

Classification transport
ONU: 3077
Transport Hazard class: 9
Packing group III

Warning
H410
P273-P391-P501a

Silver chloride > RE - Pure

RE

Description ..... Whitish irregular granules
Identification ..... Positive
Al ..... ≤ 60 ppm
Ca ..... ≤ 60 ppm
Cu ..... ≤ 100 ppm
Fe ..... ≤ 100 ppm
Mg ..... ≤ 60 ppm
Mn ..... ≤ 60 ppm
Ni ..... ≤ 60 ppm
Pb ..... ≤ 60 ppm

Code	Size	Packaging	Notes
320502	25 g	Glass bottle	
320504	100 g	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Silver diethyldithiocarbamate

• Argento dietilditiocarbammato • Argent diéthyldithiocarbamate • Plata dietilditiocarbamato
• Silberdiethyldithiocarbaminat

Synonym:
• *DETC*
• *Diethyldiocarbamic acid silver salt*

$(C_2H_5)_2NCSSAg$
Molecular Weight: 256,14
CAS: 1470-61-7
EEC-N: 216-003-7

Silver diethyldithiocarbamate > RPE - For analysis

RPE

Appearance Light yellow crystalline powder — Assay (C5H10AgNS2) ..... ≥ 99 % — Insoluble matter in water ..... ≤ 0.1 % — Loss on drying ..... ≤ 0.1 %

Code	Size	Packaging	Notes
423913	10 g	Glass bottle	

Silver manganese paper

• Carta di manganese d'argento • Papier argent manganèse • Papel de plata de manganeso • Silber-Mangan-Papier

H412
P273-P501a

Silver manganese paper > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611078200	50 stripes	Tube	Ref Ph.Eur 1078200

Silver nitrate

• Argento nitrato • Argent nitrate • Plata nitrato • Silbernitrat

Synonym:
Nitric acid silver(I) salt

$AgNO_3$
Molecular Weight: 169,87
CAS: 7761-88-8
EEC-N: 231-853-9

Classification transport
ONU: 1493
Transport Hazard class: 5.1
Packing group II



Danger
H272-H290-H314-H410
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Silver nitrate > RPE - For analysis - ACS

RPE

Description ..... White crystals
Identification ..... Positive
Appearance of solution ..... Conform ACS
Fee acidity ..... Conform ACS
Substances not ppt HCl ..... ≤100 ppm
Chloride ..... ≤5 ppm
Sulphate ..... ≤20 ppm
Cu ..... ≤2 ppm
Fe ..... ≤2 ppm
Pb ..... ≤10 ppm
Assay (argentimetric) ..... ≥99.0 %

Code	Size	Packaging	Notes
423952	25 g	Glass bottle	
423954	100 g	Glass bottle	
423955	250 g	Glass bottle	
423957	1 kg	Plastic bottle	

Silver nitrate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP-DAB

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .... Conform USP-NF
Acidity or alkalinity ..... Conform Ph.Eur.
Al,Pb,Cu,Bi ..... Conform Ph.Eur.
Foreign salts ..... ≤ 0.3 %
Assay (argentimetric) ..... 99.8 ÷ 100.5 %

Code	Size	Packaging	Notes
320904	100 g	Glass bottle	
320907	1 kg	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Silver nitrate solution 5%

• Argento nitrato soluzione 5% • Argent nitrate solution 5% • Plata nitrato solución 5% • Silbernitrat 5%

Synonym:
Nitric acid silver(I) salt

$AgNO_3$
Molecular Weight: 169,87
CAS: 7761-88-8

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group II



Danger
H314-H411
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Silver nitrate solution 5% > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Density at 15° C ..... 1.036 ÷ 1.044 Assay ..... 4.0 ÷ 6.0 %

Code	Size	Packaging	Notes
E423982	1 l	Bottle	

Silver nitrate solution 2.9075%

• Argento nitrato soluzione 2.9075% • Argent nitrate solution 2.90756% • Plata nitrato solución 2.90756% • Silbernitrat 2.90756%

Synonym:
Nitric acid silver(I) salt

$AgNO_3$
Molecular Weight: 169,87
CAS: 7761-88-8

H412
P273-P501a

Silver nitrate solution 2.9075% > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
E424001	1 l	Bottle	

Silver nitrate 1 mol/l (1N)

• Argento nitrato 1 mol/l (1N) • Argent nitrate 1 mol/l (1N) • Plata nitrato 1 mol/l (1N) • Silbernitrat 1 mol/l (1N)

Synonym:
Nitric acid silver(I) salt

$AgNO_3$
Molecular Weight: 169,87
CAS: 7761-88-8

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group II



Danger
H314-H411
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Silver nitrate 1 mol/l (1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.999 - 1.001 N NIST 399 c

Code	Size	Packaging	Notes
424036000	500 ml	Glass bottle	Certified with NIST traceability
424035000	1 l	Glass bottle	Certified with NIST traceability

169.87 g of AgNO3. Volumetric solution ready-to-use

Silver nitrate 0.5 mol/l (0.5N)

• Argento nitrato 0.5 mol/l (0.5N) • Argent nitrate 0.5 mol/l (0.5N) • Plata nitrato 0.5 mol/l (0.5N) • Silbernitrat 0.5 mol/l (0.5N)

Synonym:
Nitric acid silver(I) salt

$AgNO_3$
Molecular Weight: 169,87
CAS: 7761-88-8

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group II



Danger
H314-H411
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Silver nitrate 0.5 mol/l (0.5N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.4995 - 0.5005 N NIST 399 c

Code	Size	Packaging	Notes
424051000	1 l	Glass bottle	Certified with NIST traceability

Volumetric solution ready-to-use

Product specifications are subject to changes. Please visit our website for updates.

Silver nitrate 0.1 mol/l (0.1N)

• Argento nitrato 0.1 mol/l (0.1N) • Argent nitrate 0.1 mol/l (0.1N) • Plata nitrato 0.1 mol/l (0.1N)
• Silbernitrat 0.1 mol/l (0.1N)

Synonym:
Nitric acid silver(I) salt

$AgNO_3$
Molecular Weight: 169,87
CAS: 7761-88-8

Warning
H315-H319-H412
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Silver nitrate 0.1 mol/l (0.1N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613005600	1 l	Glass bottle	Ref Ph.Eur 3005600

Storage: protected from light

Silver nitrate 0.1 mol/l (0.1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.0999 - 0.1001 N NIST 399 c

Code	Size	Packaging	Notes
424067000	1 l	Plastic bottle	Certified with NIST traceability
424062000	5 l	Kubidos	Certified with NIST traceability
424063000	5 l	Plastic tank	Certified with NIST traceability
424061000	10 l	Kubidos	Certified with NIST traceability

16.987 g of AgNO3. Volumetric solution ready-to-use

Silver nitrate 0.1 mol/l (0.1N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
424081		Glass ampoule	Volume: 60 ml

16,987 g of AgNO3. Volumetric concentrated solution to prepare 1 L of solution 0,1 N

Silver nitrate 0.1 mol/l (0.1N) in 2-propanol

• Argento nitrato 0.1 mol/l (0.1N) in isopropanolo • Argent nitrate 0.1 mol/l (0.1N) dans le propanol-2 • Plata nitrato 0.1 mol/l (0.1N) en propanol-2
• Silbernitrat 0.1 mol/l (0.1N) in 2-Propanol

Classification transport
ONU: 1219
Transport Hazard class: 3
Packing group II

Danger
H225-H315-H319-H336-H412
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Silver nitrate 0.1 mol/l (0.1N) in 2-propanol > RS - For analysis

RS

Assay (potentiometry) ..... 0.0998 - 0.1002 N

Code	Size	Packaging	Notes
PS0250/16	1 l	Glass bottle	

Silver nitrate 0.05 mol/l (0.05N)

• Argento nitrato 0.05 mol/l (0.05N) • Argent nitrate 0.05 mol/l (0.05N) • Plata nitrato 0.05 mol/l (0.05N)
• Silbernitrat 0.05 mol/l (0.05N)

Synonym:
Nitric acid silver(I) salt

$AgNO_3$
Molecular Weight: 169,87
CAS: 7761-88-8

H412
P273-P501a

Silver nitrate 0.05 mol/l (0.05N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.04995 - 0.05005 N NIST 399 c

Code	Size	Packaging	Notes
424101000	1 l	Plastic bottle	Certified with NIST traceability

Volumetric solution ready-to-use

Silver nitrate 0.01 mol/l (0.01N)

• Argento nitrato 0.01 mol/l (0.01N) • Argent nitrate 0.01 mol/l (0.01N) • Plata nitrato 0.01 mol/l (0.01N)
• Silbernitrat 0.01 mol/l (0.01N)

Synonym:
Nitric acid silver(I) salt

$AgNO_3$
Molecular Weight: 169,87
CAS: 7761-88-8

H412
P273-P501a

Silver nitrate 0.01 mol/l (0.01N) > RPE - For analysis

RPE

Assay (potentiometry) ..... 0.00998 - 0.01002 N

Code	Size	Packaging	Notes
PS0030/15	1 l	Plastic bottle	

Silver nitrate 0.01 mol/l (0.01N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
424161		Glass ampoule	Volume: 60 ml

1,6987 g of AgNO3. Volumetric concentrated solution to prepare 1 L of solution 0,01 N

Silver nitrate 0.01 mol/l (0.01N) in propanol-2

• Argento nitrato 0.01 mol/l (0.01N) in isopropanolo • Argent nitrate 0.01 mol/l (0.01N) dans le propanol-2 • Plata nitrato 0.01 mol/l (0.01N) en propanol-2
• Silbernitrat 0.01 mol/l (N/10) in Propanol-2

$AgNO_3$
Molecular Weight: 169,87
CAS: 7761-88-8

Classification transport
ONU: 1219
Transport Hazard class: 3
Packing group II



Danger
H225-H319-H336
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Silver nitrate 0.01 mol/l (0.01N) in propanol-2 > RS - For analysis

RS

Assay (potentiometry) ..... 0.00998 - 0.01002 N

Code	Size	Packaging	Notes
PS0252/16	1 l	Glass bottle	

Silver nitrate solution

• Argento nitrato soluzione • Argent nitrate solution • Plata nitrato solución • Silbernitratlösung

Synonym:
Nitric acid silver(I) salt

$AgNO_3$
Molecular Weight: 169,87
CAS: 7761-88-8



Danger
H315-H318-H411
P264-P280a-P305+P351+P338-P310a-P362+P364-P332+P313

Silver nitrate solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611078307	100 ml	Glass bottle	Silver nitrate R1 Ref Ph.Eur 1078301
611078306	100 ml	Glass bottle	Silver nitrate R2 Ref Ph.Eur 1078302
611078301	1 l	Glass bottle	Silver nitrate R1 Ref Ph.Eur 1078301
611078302	1 l	Glass bottle	Silver nitrate R2 Ref Ph.Eur 1078302

Storage: protected from light

Silver nitrate solution > RS - For analysis according to JP

RS

Code	Size	Packaging	Notes
616001001	1 l	Plastic bottle	Silver nitrate TS

Storage: protected from light

Product specifications are subject to changes. Please visit our website for updates.

Silver nitrate solution > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
617000201	1 l	Glass bottle	Silver nitrate TS

Silver oxide

• Argento ossido • Argent oxyde • Plata óxido • Silberoxid

Synonym:
Silver(I) oxide

Ag_2O
Molecular Weight: 231,74
CAS: 20667-12-3
EEC-N: 243-957-1

Classification transport
ONU: 1479
Transport Hazard class: 5.1
Packing group II



Danger
H271-H318-H410
P210-P280-P283-P305+P351+P338-P310a-P306+P360

Silver oxide > RPE - For analysis

RPE

Description ..... Grey powder
Identification ..... Positive
Loss on drying ..... ≤ 0.5 %
Ca ..... ≤ 50 ppm
Cd ..... ≤ 50 ppm
Cu ..... ≤ 0.08 %
Fe ..... ≤ 150 ppm
Mg ..... ≤ 50 ppm
Ni ..... ≤ 50 ppm
Pb ..... ≤ 50 ppm
Zn ..... ≤ 50 ppm
Assay (argentimetric) ..... ≥ 99 %

Code	Size	Packaging	Notes
424181	25 g	Glass bottle	
424182	250 g	Plastic bottle	

Silver sulfate

• Argento solfato • Argent sulfate • Plata sulfato • Silbersulfat

Synonym:
Sulfuric acid disilver(I) salt

Ag_2SO_4
Molecular Weight: 311,79
CAS: 10294-26-5
EEC-N: 233-653-7



Danger
H318-H410
P273-P280i-P305+P351+P338-P310a-P391-P501a

Silver sulfate > RPE - For analysis

RPE

Description ..... Whitish powder
Identification ..... Positive
Umidità (H2O) ..... ≤ 0.5 %
Cu ..... ≤ 500 ppm
Fe ..... ≤ 150 ppm
Assay (argentimetric) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
424201	25 g	Glass bottle	
424203	100 g	Glass bottle	

Silver sulfate solution 0.7% in sulfuric acid

• Argento solfato soluzione 0.7% in acido solforico • Argent sulfate solution 0.7% dans l'acide sulfurique
• Plata sulfato solución 0.7% en acido sulfúrico • Silbersulfat 0.7% in Schwefelsäure

Synonym:
Sulfuric acid disilver(I) salt

Ag_2SO_4
Molecular Weight: 311,79
CAS: 10294-26-5

Classification transport
ONU: 1830
Transport Hazard class: 8
Packing group II



Danger
H290-H314-H412
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Silver sulfate solution 0.7% in sulfuric acid > RS - For environmental analysis (COD determination)

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Assay ..... 0.685 ÷ 0.715 %

Code	Size	Packaging	Notes
424191	1 l	Glass bottle	
424192	2.5 l	Glass bottle	

Soap solution in ethanol

• Sapone soluzione in alcole etilico • Savon en solution dans l'éthanol • Jabón solución en alcohol etilico • Seife in Ethanollösung

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II



Warning
H226-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

Soap solution in ethanol > RS - For hydrotimetry according to Boutron-Boudet

RS

Description ..... Yellowish clear liquid Identification ..... Positive

Code	Size	Packaging	Notes
E477507	1 l	Glass bottle	

Soda lime

• Calce sodata • Chaux sodée • Cal sodada • Natronkalk

CAS: 8006-28-8



Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Soda lime > RS - For anesthesia and basal metabolism test

RS

Description ..... White granules
Identification ..... Positive
Water ..... 12.0 ÷ 19.0 %
Alcalinity (NaOH) ..... ≤ 3.5 %
Activity ..... ≥ 19.0 %
Hardness ..... ≥ 75.0 %
Umidità assorbita ..... ≤ 7.5 %
Diameter ..... 2.5 ÷ 5.0 mm

Code	Size	Packaging	Notes
432873	1 kg	Plastic bottle	
432874	4,5 kg	Plastic tank	

With ethyl violet indicator

Soda lime > RS - For CO2 absorption

RS

Description ..... Granuli white
Identification ..... Positive
CO2 Absorption ..... > 19 %
Diameter: ..... Conform
> 2.80 mm ..... < 1.0 %
> 2.00 mm ..... < 30.0 %
> 1.40 mm ..... Balance
> 0.60 mm ..... < 20.0 %
< 0.60 mm ..... < 1.0 %

Code	Size	Packaging	Notes
432861	1 kg	Plastic bottle	
432862	5 kg	Plastic tank	

With indicator manganese salt

Soda lime > RS - For microanalysis

RS

Description ..... White granules
Identification ..... Positive
Water ..... 16 ÷ 19 %
Hardness ..... ≥80 %
> 4.75 mm ..... ≤7.0 %
< 0.6 mm ..... ≤1.0 %

Code	Size	Packaging	Notes
432851	500 g	Plastic bottle	

With indicator ethyl violet

Soda lime > RPE - For analysis

RPE

Description ..... White granules
Identification ..... Positive
Water ..... 12 ÷ 19 %
Hardness ..... ≥ 75 %
Activity ..... ≥ 19.0 %
Alcalinity (NaOH) ..... ≤ 3.5 %
> 8.0 mm ..... Nil
> 4.75 mm ..... ≤ 7.0 %
< 0.425 mm ..... ≤ 2.0 %

Code	Size	Packaging	Notes
432801	1 kg	Plastic bottle	
432802	2.5 kg	Plastic bottle	
432803	25 kg	Plastic bucket	

With ethyl purple indicator. Diameter 2,5 - 6 mm

Product specifications are subject to changes. Please visit our website for updates.

Sodium standard solution

• Sodio standard soluzione • Sodium solution standard • Sodio solución patrón • Natriumstandardlösung

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sodium standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615002701	100 ml	Plastic bottle	A 50 ppm solution: to dilute according to Ref Ph.Eur 5002701
615002709	100 ml	Plastic bottle	A 200 ppm solution: to dilute according to Ref Ph.Eur 5002700
615005700	1 l	Plastic bottle	A 1000 ppm solution Ref Ph.Eur 5005700

Sodium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505732	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505735	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505733	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Sodium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503741	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503743	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503745	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503747	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Sodium standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
507759	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503749	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
E497645	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
E497641	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Sodium standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
478101		Plastic ampoule	conc. 1.000 ppm Matrix: Nitric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Sodium standard solution > RS - Standard solution for ion chromatography

RS

Code	Size	Packaging	Notes
503301	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503303	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Sodium acetate anhydrous

• Sodio acetato anidro • Sodium acétate anhydre • Sodio acetato anhidro • Natriumacetat wasserfrei

Synonym: *Acetic acid sodium salt anhydrous*

CH_3COONa
Molecular Weight: 82,035
CAS: 127-09-3
EEC-N: 204-823-8

Sodium acetate anhydrous > RPE - For analysis

RPE

Description	White hygroscopic powder	Heavy metals (Pb)	≤10 ppm	Ca	≤100 ppm	Ni	≤2 ppm
Identification	Positive	Nitrate	≤10 ppm	Cu	≤2 ppm	Zn	≤2 ppm
pH sol. 5% at 25° C	7.5 ÷ 9.2	Subst. reducing KMnO4	≤50 ppm (1h)	Fe	≤3 ppm	Assay (non-aqueos medium)	99.0 ÷ 101.0 % (s.s.)
Loss on drying	≤1.0 %	Sulphate	≤100 ppm	K	≤0.1 %		
Ca, Mg and ppt by NH4OH	≤50 ppm	Al	≤5 ppm	Mg	≤5 ppm		
Chloride	≤350 ppm	As	≤1 ppm	Mn	≤5 ppm		

Code	Size	Packaging	Notes
478165	100 g	Plastic bottle	
478166	500 g	Plastic bottle	
478167	1 kg	Plastic bottle	
478163	25 kg	Fibre drum	

Sodium acetate anhydrous > ERBApharm - According to pharmacopoeia: USP

ERBApharm

Description	White crystalline powder	pH	7.5 ÷ 9.2 USP-NF	Heavy metals (Pb)	≤10 ppm	% s.s.	
Identification	Positive	Loss on drying	≤1.0 %	Sulphate	≤50 ppm	Origin (BSE/TSE)	Synthesis
Calcium + Magnesium	Conform USP-NF	Water not sol. matter	≤500 ppm	Al	≤0.2 ppm		
K	Conform USP-NF	Chloride	≤350 ppm	Assay (non-aqueos medium)	99.0 ÷ 101.0		

Code	Size	Packaging	Notes
366377	1 kg	Plastic bottle	
366372	5 kg	Plastic tank	
366371	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium acetate trihydrate

• Sodio acetato triidrato • Sodium acétate trihydraté • Sodio acetato trihidrato • Natriumacetat-trihydrat

Synonym: *Acetic acid sodium salt trihydrate*

$CH_3COONa.3H_2O$
Molecular Weight: 136,08
CAS: 6131-90-4
EEC-N: 612-115-9

Sodium acetate trihydrate > RPE - For analysis - ISO - ACS - Reag.USP

RPE

Description	White crystals	Chloride	≤10 ppm	Sulphate	≤20 ppm	%	
Identification	Positive	Phosphate	≤5 ppm	Fe	≤5 ppm	Ca	≤50 ppm
pH sol. 5% at 25° C	7.5 ÷ 9.2	Water-insoluble matter	≤50 ppm	K	≤50 ppm	Mg	≤20 ppm
Subst. reducing KMnO4	Conform	Heavy metals (Pb)	≤5 ppm	Assay (non-aqueos medium)	99.0 ÷ 101.0		

Code	Size	Packaging	Notes
478135	100 g	Plastic bottle	
478136	500 g	Plastic bottle	
478137	1 kg	Plastic bottle	
478139	5 kg	Plastic tank	
478132	25 kg	Plastic bucket	

Product specifications are subject to changes. Please visit our website for updates.

Sodium acetate trihydrate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-BP

ERBApharm

Description ..... White or almost white,cryst. powder or colourless cryst.
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Reducing substances .....Conform Ph.Eur.
K..... Conform USP-NF
pH sol. 5% at 25° C .....7.5 ÷ 9.0
Loss on drying 130° C.....39.0 ÷ 40.5 %
Chloride.....≤200 ppm
Heavy metals (Pb).....≤10 ppm
Sulphate.....≤50 ppm
Not soluble matter.....≤500 ppm
As .....≤2 ppm
Calcium + Magnesium .....≤50 ppm
Fe .....≤10 ppm
Al .....≤0.2 ppm
Assay (non-aqueos medium) .99.0 ÷ 101.0 % (s.s.)

Code	Size	Packaging	Notes
366207	1 kg	Plastic bottle	
366209	5 kg	Plastic tank	
366205	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium alginate

• Sodio alginato • Sodium alginate • Sodio alginato • Natriumalginat

Synonym:
• *Algin*
• *Alginic acid*

$(C_6H_7O_6Na)n$
Molecular Weight: >200000
CAS: 9005-38-3

Sodium alginate > ERBApharm - According to pharmacopoeia: FU-Ph.Eur.

ERBApharm

Description ..... Beige powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Loss on drying ..... ≤15.0 %
Sulphated ash.....30.0 ÷ 36.0 %
Chloride.....≤1.0 %
Heavy metals (Pb).....≤20 ppm
Ca..... ≤1.50 %
Microbial tests
TAMC ..... ≤1000 CFU/g
TYMC ..... ≤100 CFU/g
Escherichia coli ..... Absent
Salmonella..... Absent

Code	Size	Packaging	Notes
366551	100 g	Plastic bottle	
366552	1 kg	Plastic bottle	
366553	5 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium aluminate

• Sodio alluminato • Sodium aluminate • Sodio aluminato • Natriumaluminat

$NaAlO_2$
Molecular Weight: 81,97
CAS: 11138-49-1
EEC-N: 234-391-6

Classification transport
ONU: 2812
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sodium aluminate > RE - Pure

RE

Description .....White-grey powder
Identification ..... Positive
Iron (as Fe203) .....≤ 0.05 %
Aluminium (as Al203)..... 50 - 60 %
Sodium (as Na20) ..... 37 - 45 %

Code	Size	Packaging	Notes
478237	1 kg	Plastic bottle	
478232	25 kg	Plastic bucket	

Sodium ammonium hydrogen phosphate

• Sodio ammonio idrogeno fosfato • Sodium ammonium hydrogénophosphate
• Sodio y amonio hidrógeno fosfato • Natriumammoniumhydrogenphosphat

Synonym:
• *Ammonium sodium phosphate dibasic tetrahydrate*
• *Ammonium sodium hydrogen phosphate*

$NaNH_4HPO_4.4H_2O$
Molecular Weight: 209,07
CAS: 13011-54-6
EEC-N: 235-860-8

Sodium ammonium hydrogen phosphate > RPE - For analysis

RPE

Description	White crystalline powder	Fluoride	≤10 ppm	As	≤0.5 ppm	Mg	≤20 ppm
Identification	Positive	Water-insoluble matter	≤50 ppm	Ca	≤50 ppm	Ni	≤25 ppm
pH sol. 5% at 25° C	7.5 ÷ 8.5	Heavy metals (Pb)	≤5 ppm	Cu	≤25 ppm	Zn	≤25 ppm
Carbonate	≤10 ppm	Nitrate	≤10 ppm	Fe	≤5 ppm	Assay (acidimetric)	≥99.5 %
Chloride	≤10 ppm	Sulphate	≤50 ppm	K	≤300 ppm		

Code	Size	Packaging	Notes
478357	1 kg	Plastic bottle	

Sodium L-ascorbate

• Sodio L-ascorbato • Sodium L-ascorbate • Sodio L-ascorbato • Natrium L-ascorbat

Synonym:
• *L(+)Ascorbic acid sodium salt*
• *Vitamin C sodium salt*

$C_6H_7O_6Na$
Molecular Weight: 198,11
CAS: 134-03-2
EEC-N: 205-126-1

Sodium L-ascorbate > RE - Pure

RE

Description ..... white to yellow crystalline powder Identification ..... Positive Assay (HClO4) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
366681	100 g	Plastic bottle	
366684	1 kg	Plastic bottle	

Sodium arsenite 0.1 mol/l (0.2N)

• Sodio arsenito 0.1 mol/l (0.2N) • Sodium arsenite 0.1 mol/l (0.2N) • Sodio arsenito 0.1 mol/l (0.2N)
• Natriumarsenit 0.1 mol/l (0.2N)

Synonym:
• *Sodium (meta)arsenite*
• *Sodium dioxoarsenate*

$AsNaO_2$
Molecular Weight: 129,91
CAS: 7784-46-5

Classification transport
ONU: 1686
Transport Hazard class: 6.1
Packing group III



Danger
H302-H350-H411-HA26
P264-P280-P301+P312a-P330-P308+P313-P501a

Sodium arsenite 0.1 mol/l (0.2N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613005800	1 l	Plastic bottle	Ref Ph.Eur 3005800

Sodium arsenite 0.05 mol/l (0.1N)

• Sodio arsenito 0.05 mol/l (0.1N) • Sodium arsenite 0.05 mol/l (0.1N) • Sodio arsenito 0.05 mol/l (0.1N)
• Natriumarsenit 0.05 mol/l (0.1N)

Synonym:
• *Sodium (meta)arsenite*
• *Sodium dioxoarsenate*

$AsNaO_2$
Molecular Weight: 129,91
CAS: 7784-46-5

Classification transport
ONU: 1686
Transport Hazard class: 6.1
Packing group II



Danger
H302-H319-H350-H411-HA26
P264-P280-P301+P312a-P305+P351+P338-P308+P313-P337+P313

Sodium arsenite 0.05 mol/l (0.1N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
402381		Glass ampoule	Volume: 60 ml

6,494 g NaAsO2. Volumetric concentrated solution to prepare 1 L of solution 0,1 N

Product specifications are subject to changes. Please visit our website for updates.

Sodium azide

• Sodio azide • Sodium azide • Sodio azida • Natriumazid

NaN_3
Molecular Weight: 65,01
CAS: 26628-22-8
EEC-N: 247-852-1

Classification transport
ONU: 1687
Transport Hazard class: 6.1
Packing group II

Danger
H300-H310-H373-H410-HEU032
P260-P264-P280h-P301+P310a-P330-P361+P364

Sodium azide > RE - Pure

RE

Description ..... White crystalline powder
Identification ..... Positive
Loss on drying ..... ≤ 0.5 %
Water-insoluble matter ..... ≤ 500 ppm
pH solution 5% ..... ≥ 9
Heavy metals (Pb) ..... ≤ 20 ppm
Sodium carbonate ..... ≤ 0.15 %
Assay (oxidimetric) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
478484	25 g	Glass bottle	
478482	250 g	Glass bottle	
478481	2.5 kg	Plastic bottle	

Sodium benzoate

• Sodio benzoato • Sodium benzoate • Sodio benzoato • Natriumbenzoat

Synonym:
Benzoic acid sodium salt

C_6H_5COONa
Molecular Weight: 144,11
CAS: 532-32-1
EEC-N: 208-534-8

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Sodium benzoate > ERBApharm - According to pharmacopoeia: BP-FU-NF-Ph.Eur.

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution ..... Conform Ph.Eur.
Acidity or alkalinity ..... Conform Ph.Eur.
Water (K.F.) ..... ≤1.5 %
Loss on drying ..... ≤2.0 %
Halog. comp. ionized Cl ..... ≤200 ppm
Halog. comp. total Cl ..... ≤300 ppm
Heavy metals (Pb) ..... ≤10 ppm
Assay (non-aqueos medium) .99.0 ÷ 100.5 % s.s.

Code	Size	Packaging	Notes
366757	1 kg	Plastic bottle	
366759	5 kg	Plastic bucket	
366754	25 kg	Fibre drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium bicarbonate

• Sodio bicarbonato • Sodium bicarbonate • Sodio bicarbonato • Natriumbicarbonat

Synonym:
Sodium hydrogen carbonate

$NaHCO_3$
Molecular Weight: 84,01
CAS: 144-55-8
EEC-N: 205-633-8

Sodium bicarbonate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White crystalline powder
Identification ..... Positive
Water-insoluble matter ..... ≤150 ppm
Ammonium ..... ≤5 ppm
Chloride ..... ≤30 ppm
Phosphate ..... ≤10 ppm
Total sulphur ..... ≤30 ppm
Heavy metals (Pb) ..... ≤5 ppm
Fe ..... ≤10 ppm
K ..... ≤50 ppm
Assay (alkalimetric) ..... 99.7 ÷ 100.3 % s.s.
Ca ..... ≤100 ppm
Mg ..... ≤50 ppm
As ..... ≤2 ppm
Appearance of solution ..... Conform
Carbonate ..... Conform
Loss on drying ..... ≤0.25 %
Sulphate ..... ≤150 ppm

Code	Size	Packaging	Notes
478535	100 g	Plastic bottle	
478536	500 g	Plastic bottle	
478537	1 kg	Plastic bottle	
478531	5 kg	Plastic jar	
478532	25 kg	Plastic bucket	

Sodium bicarbonate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-BP-DAB

ERBApharm

Description ..... White crystalline powder	Not soluble matter ..... Conform USP-NF	Heavy metals (Pb) ..... ≤5 ppm	Ca ..... ≤100 ppm
Identification ..... Positive	Loss (silica gel) ..... ≤0.25 %	Sulphate ..... ≤150 ppm	Fe ..... ≤20 ppm
Appearance of solution ..... Conform Ph.Eur.	Ammonium ..... ≤20 ppm	Sulfur compounds ..... ≤150 ppm	Assay (alkalimetric) ..... 99.0 ÷ 100.5 % s.s.
Carbonate ..... Conform Ph.Eur.	Chloride ..... ≤150 ppm	As ..... ≤2 ppm	Normal carbonate ..... Conform

Code	Size	Packaging	Notes
366908	1 kg	Plastic bottle	
366909	5 kg	Plastic tank	
366902	25 kg	Plastic bucket	
366904	50 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium bisulfate monohydrate

• Sodio bisolfato monoidrato • Sodium bisulfate monohydraté • Sodio bisulfato monohidrato • Natriumbisulfat-Monohydrat

Synonym: *Sodium hydrogen sulfate monohydrate*

$NaHSO_4.H_2O$
Molecular Weight: 138,07
CAS: 10034-88-5
EEC-N: 231-665-7

Classification transport
ONU: 3260
Transport Hazard class: 8
Packing group III



Danger
H318
P280i-P305+P351+P338-P310a

Sodium bisulfate monohydrate > RPE - For analysis

RPE

Description ..... White crystals	Chloride ..... ≤ 20 ppm	Fe ..... ≤ 50 ppm
Identification ..... Positive	Heavy metals (Pb) ..... ≤ 50 ppm	Assay ..... ≥ 96.0 %

Code	Size	Packaging	Notes
478675	100 g	Plastic bottle	
478676	500 g	Plastic bottle	
478677	1 kg	Plastic bottle	
478673	25 kg	Drum	

Sodium borohydride

• Sodio boroidruro • Sodium borohydrure • Sodio boro hidruro • Natriumtetrahydroborat

Synonym: *Sodium tetrahydroborate*

$NaBH_4$
Molecular Weight: 37,83
CAS: 16940-66-2
EEC-N: 241-004-4

Classification transport
ONU: 1426
Transport Hazard class: 4.3
Packing group I



Danger
H260-H301-H314-H360FD-HEU014-HA26
P223-P280-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338

Sodium borohydride > RE - Pure - Powder

RE

Description ..... White powder | Identification ..... Positive | Assay (oxidimetric) ..... ≥ 95 %

Code	Size	Packaging	Notes
478953	50 g	Metallic can	
478955	250 g	Metallic can	
478957	1 kg	Metallic can	

Sodium borohydride > RE - Pure - Pearls

RE

Description ..... White pearls | Identification ..... Positive | Assay ..... ≥ 97.5 %

Code	Size	Packaging	Notes
478964	100 g	Metallic can	

Product specifications are subject to changes. Please visit our website for updates.

Sodium bromide

• Sodio bromuro • Sodium bromure • Sodio bromuro • Bromnatrium

NaBr
Molecular Weight: 102,9
CAS: 7647-15-6
EEC-N: 231-599-9

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Sodium bromide > RPE - For analysis - ACS

RPE

Description	White crystalline powder	Chloride	≤ 0.2 %	Ba	≤ 20 ppm	Mg	≤ 10 ppm
Identification	Positive	Water-insoluble matter	≤ 50 ppm	Ca	≤ 20 ppm	Assay (argentimetric)	≥ 99.0 %
pH sol. 5% at 25° C	5.0 ÷ 8.8	Heavy metals (Pb)	≤ 5 ppm	Fe	≤ 5 ppm		
Bromate	≤ 10 ppm	Sulphate	≤ 20 ppm	K	≤ 0.1 %		

Code	Size	Packaging	Notes
479055	250 g	Plastic bottle	
479057	1 kg	Plastic bottle	

Sodium bromide > ERBApharm - According to pharmacopoeia: Ph.Eur.-Ph.Franc.

ERBApharm

Description	White crystalline powder	Bromate	Conform Ph.Eur.	Mg,alkal.earth met.(Ca)	≤200 ppm	Assay (argentimetric)	98.0 ÷ 100.5 % s.s.
Identification	Positive	Iodide	Conform Ph.Eur.	Heavy metals (Pb)	≤10 ppm		
Appearance of solution	Conform Ph.Eur.	Loss on drying	≤3.0 %	Sulphate	≤100 ppm		
Acidity or alkalinity	Conform Ph.Eur.	Chloride	≤0.6 %	Fe	≤20 ppm		

Code	Size	Packaging	Notes
367357	1 kg	Plastic bottle	
367359	5 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium 1-butanesulfonate ▶ 1-Butanesulfonic acid sodium salt

Sodium carbonate anhydrous

• Sodio carbonato anidro • Sodium carbonate anhydre • Sodio carbonato anhidro
• Natriumcarbonat wasserfrei

Synonym:
• *Calcined soda*
• *Carbonic acid disodium salt*

Na_2CO_3
Molecular Weight: 105,99
CAS: 497-19-8
EEC-N: 207-838-8

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Sodium carbonate anhydrous > RS - Standard for volumetry

RS

Description	White crystals	Identification	Positive	Assay	≥99.7 %

Code	Size	Packaging	Notes
479331	50 g	Glass bottle	

Sodium carbonate anhydrous > RPE - For analysis - ISO - ACS - Reag.Ph.Eur. - Reag.USP

RPE

Description	White crystalline powder	Phosphate	≤ 10 ppm	Total sulphur	≤ 30 ppm	Ca	≤ 300 ppm
Identification	Positive	Water-insoluble matter	≤ 100 ppm	Fe	≤ 5 ppm	Mg	≤ 50 ppm
Loss on drying (285°C)	≤ 1.0 %	Heavy metals (Pb)	≤ 5 ppm	K	≤ 50 ppm		
Chloride	≤ 10 ppm	Silicate	≤ 50 ppm	Assay (alkalimetric)	≥ 99.5 % s.s.		

Code	Size	Packaging	Notes
479305	100 g	Plastic bottle	
479306	500 g	Plastic bottle	
479307	1 kg	Plastic bottle	
479301	5 kg	Plastic tank	
479302	25 kg	Drum	

Sodium carbonate anhydrous > ERBApharm - According to pharmacopoeia: Ph.Eur.-NF

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Loss on drying ..... ≤0.5 %
Appearance of solution .....Conform Ph.Eur.
Alkali hydroxides and bicarbonates Conform Ph.Eur.
Chloride.....≤ 125 ppm
Sulphate.....≤ 250 ppm
As .....≤ 5 ppm
Fe .....≤ 50 ppm
Heavy metals (Pb).....≤ 10 ppm
Residue solvents ..... Conform USP-NF
Assay (acidimetric) .....99.5 ÷ 100.5 %s.s.

Code	Size	Packaging	Notes
367707	1 kg	Plastic bottle	
367703	5 kg	Plastic tank	
367705	25 kg	Plastic bucket	
367704	50 kg	Fibre drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium carbonate decahydrate

• Sodio carbonato decaidrato • Sodium carbonate décahydraté • Sodio carbonato decahidrato • Sodium carbonate decahydrate

$Na_2CO_3.10H_2O$
Molecular Weight: 286,14
CAS: 6132-02-1
EEC-N: 207-838-8

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Sodium carbonate decahydrate > RPE - For analysis - ISO

RPE

Description .....White crystals
Identification ..... Positive
Free alkalies (NaOH) .....≤400 ppm
Total nitrogen.....≤5 ppm
Bicarbonate ..... ≤0.2 %
Chloride.....≤5 ppm
Phosphate .....≤5 ppm
Water-insoluble matter .....≤25 ppm
Heavy metals (Pb).....≤3 ppm
Subst. ppt by NH4OH.....≤100 ppm
Reducing iodine .....≤50 ppm
Silicate .....≤20 ppm
Total sulphur .....≤10 ppm
Al .....≤5 ppm
As .....≤0.1 ppm
Ca .....≤20 ppm
Cu .....≤2 ppm
Fe .....≤2 ppm
K.....≤50 ppm
Mg .....≤2 ppm
Ni.....≤2 ppm
Pb .....≤2 ppm
Zn .....≤2 ppm
Assay (alkalimetric)..... ≥99.5 %

Code	Size	Packaging	Notes
479125	100 g	Plastic bottle	
479126	500 g	Plastic bottle	
479127	1 kg	Plastic bottle	
479121	5 kg	Plastic tank	
479122	25 kg	Plastic bucket	

Sodium carbonate decahydrate > ERBApharm - According to pharmacopoeia: BP-FU-Ph.Eur.-Ph.Franc.

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Alc.hydroxides + bicar. .....Conform Ph.Eur.
Chloride.....≤50 ppm
Heavy metals (Pb).....≤20 ppm
Sulphate.....≤100 ppm
As .....≤2 ppm
Fe .....≤20 ppm
Assay (alkalimetric)36.7 ÷ 40.0 % Na2C03

Code	Size	Packaging	Notes
367608	1 kg	Plastic bottle	
367609	5 kg	Plastic tank	
367601	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium carbonate monohydrate

• Sodio carbonato monoidrato • Sodium carbonate monohydraté • Sodio carbonato monohidrato
• Natriumcarbonat-Monohydrat

Synonym:
Carbonic acid disodium salt

$Na_2CO_3.H_2O$
Molecular Weight: 105,99
CAS: 5968-11-6
EEC-N: 207-838-8

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Sodium carbonate monohydrate > RS - For analysis according to Ph. Eur. Chap. 4.2.1

RS

Code	Size	Packaging	Notes
612000500	50 g	Plastic bottle	Ref Ph.Eur 2000500

Product specifications are subject to changes. Please visit our website for updates.

Sodium carbonate monohydrate > RPE - For analysis - ACS

RPE

Description .....White, granular crystals
Identification ..... Positive
Loss on drying .....13.0 ÷ 15.0 %
Ca..... ≤ 0.03 %
Chloride.....≤ 10 ppm
Phosphate .....≤ 5 ppm
Water-insoluble matter .....≤ 100 ppm
Heavy metals (Pb).....≤ 5 ppm
Mg .....≤ 50 ppm
Silicate .....≤ 50 ppm
Total sulphur .....≤ 40 ppm
Fe .....≤ 5 ppm
K.....≤ 50 ppm
Assay (alkalimetric)..... ≥ 99.5 %

Code	Size	Packaging	Notes
479255	100 g	Plastic bottle	
479256	500 g	Plastic bottle	
479257	1 kg	Plastic bottle	

Sodium carbonate monohydrate > ERBApharm - According to pharmacopoeia: FU-Ph.Eur.-Ph.Franc.

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Alkali hydroxides and bicarbonates Conform Ph.Eur.
Chloride.....≤ 125 ppm
Sulphate.....≤ 250 ppm
Heavy metals (Pb).....≤ 50 ppm
As .....≤ 5 ppm
Fe .....≤ 50 ppm
Assay (alkalimetric).....83.0 ÷ 87.5 %

Code	Size	Packaging	Notes
367691	1 kg	Plastic bottle	
367692	5 kg	Plastic tank	
367693	25 kg	Plastic bucket	
367694	50 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium carbonate 0.5 mol/l (1N)

• Sodio carbonato 0.5 mol/l (1N) • Sodium carbonate 0.5 mol/l (1N) • Sodio carbonato 0.5 mol/l (1N)
• Natriumcarbonat 0.5 mol/l (1N)

Synonym:
• *Calcined soda*
• *Carbonic acid disodium salt*

Molecular Weight: 105,99

HEU210

Sodium carbonate 0.5 mol/l (1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Assay (potentiometry) ..... 0.998 - 1.002 N

Code	Size	Packaging	Notes
479186	500 ml	Plastic bottle	

52,995 g of Na2CO3. Volumetric solution ready-to-use

Sodium carbonate 0.05 mol/l (0.1N)

• Sodio carbonato 0.05 mol/l (0.1N) • Sodium carbonate 0.05 mol/l (0.1N) • Sodio carbonato 0.05 mol/l (0.1N)
• Natriumcarbonat 0.05 mol/l (0.1N)

Synonym:
• *Calcined soda*
• *Carbonic acid disodium salt*

HEU210

Sodium carbonate 0.05 mol/l (0.1N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
479211		Plastic ampoule	Volume: 55 ml

5,299 g of Na2CO3. Volumetric concentrated solution to prepare 1 L of solution 0,1 N

Sodium carbonate solution 20%

• Sodio carbonato soluzione 20% • Sodium carbonate solution 20% • Sodio carbonato solución 20%
• Natriumcarbonat 20%

Synonym:
• *Calcined soda*
• *Carbonic acid disodium salt*

HEU210

Sodium carbonate solution 20% > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Density at 20° C ..... 1.072 ÷ 1.078

Code	Size	Packaging	Notes
479151	1 l	Plastic bottle	

Mass percentage based on Na2CO3.10H2O content

Sodium carbonate solution

• Sodio carbonato soluzione • Sodium carbonate anhydre solution • Sodio carbonato solución
• Natriumcarbonat wasserfreie Lösung

Synonym:
• *Calcined soda*
• *Carbonic acid disodium salt*

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Sodium carbonate solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611079301	1 l	Plastic bottle	A 106 g/l solution ref Ph.Eur 1079301

Sodium chloride

• Sodio cloruro • Sodium chlorure • Sodio cloruro • Natriumchlorid

Synonym:
Halite

NaCl
Molecular Weight: 58,44
CAS: 7647-14-5
EEC-N: 231-598-3

Sodium chloride > RS - For enviromental analysis

RS

Description ..... White crystals Identification ..... Positive Hg ..... ≤0.005 ppm Assay (argentimetric) ..... ≥99.5 %

Code	Size	Packaging	Notes
479671	100 g	Glass bottle	

Low content in Hg

Sodium chloride > RS - For analysis according to Ph. Eur. Chap. 4.2.1

RS

Code	Size	Packaging	Notes
612000600	250 g	Plastic bottle	Ref Ph.Eur 2000600

Sodium chloride > RS - Standard for volumetry

RS

Description ..... White crystals Identification ..... Positive Assay ..... ≥99.5 %

Code	Size	Packaging	Notes
479652	50 g	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Sodium chloride > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description .....White crystals
Identification ..... Positive
pH sol. 5% at 25° C .....5.0 ÷ 9.0
Loss on drying ..... ≤ 0.5 %
Appearance of solution .....Conform
Bromide.....≤ 100 ppm
Nitrate,Chlorate (NO3) .....≤ 30 ppm
Phosphate .....≤ 5 ppm
Water-insoluble matter .....≤ 50 ppm
Iodide .....≤ 20 ppm
Heavy metals (Pb).....≤ 5 ppm
Sulphate .....≤ 40 ppm
Al.....≤ 0.2 ppm
As.....≤ 1 ppm
Ca.....≤ 20 ppm
Ba.....Conform
Fe .....≤ 2 ppm
K.....≤ 50 ppm
Mg .....≤ 10 ppm
Assay (argentimetric)..99.0 ÷ 100.5 % s.s.
Acidity or alkalinity.....Conform
Nitrite .....Conform
Ferrocyanide.....≤ 1 ppm
Magnesium and alkali metals....≤ 100 ppm

Code	Size	Packaging	Notes
479685	100 g	Plastic bottle	
479686	500 g	Plastic bottle	
479687	1 kg	Plastic bottle	
479689	5 kg	Plastic jar	
479681	25 kg	Plastic bucket	

Sodium chloride > RPE - For analysis - According to ASTM B117 ISO 9227/2006

RPE

Description ..... White crystalline powder
Identification ..... Positive
Heavy metals (Pb).....≤ 5 ppm
Loss on drying ..... ≤ 0.5 %
Cu.....≤ 0.3 ppm
Ni.....≤ 10 ppm
Halogen (Iodide+Bromide+fluoride)≤ 0.1 %
Impurezze totali..... ≤ 0.3 % s.s.
Assay (argentimetric) ..... ≥ 99.8 % s.s.
NaI ..... ≤ 0.1 %

Code	Size	Packaging	Notes
479663	1 kg	Plastic bottle	
479662	5 kg	Plastic bucket	
479661	25 kg	Plastic bucket	

For salt spray tests

Sodium chloride > ERBApharm - According to pharmacopoeia: Ph.Eur.- Microbiological tested

ERBApharm

Description ..... White crystalline powder
Identification (I.R.).....Positive
Appearance of solution .....Pass test
Acidity or alcalinity.....Pass test
Barium .....Pass test
Iodide .....Pass test
Ferrocyanide.....Pass test
Nitrite .....Pass test
Loss on drying ..... ≤ 0.5 %
Mg,alkal.earth met.(Ca) .....≤ 100 ppm
Bromide.....≤ 100 ppm
Phosphate .....≤ 25 ppm
Heavy metals (Pb).....≤ 5 ppm
Sulphate .....≤ 200 ppm
Al.....≤ 0.2 ppm
As .....≤ 1 ppm
Fe .....≤ 2 ppm
K.....≤ 500 ppm
Assay (argentimetric)...99.0 ÷ 100.5 % s.s
Total aerobic microbial count (TAMC) ≤ 100 CFU/g
Total yeasts/mould count (TYMC)≤ 10 CFU/g

Code	Size	Packaging	Notes
368281	10 kg	Drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium chloride > ERBApharm-According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP-DAB

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Acidity or alkalinity.....Conform Ph.Eur.
Barium .....Conform Ph.Eur.
Iodide .....Conform Ph.Eur.
Ferrocyanide.....Conform Ph.Eur.
Nitrite .....Conform Ph.Eur.
Loss on drying ..... ≤ 0.5 %
Mg,alkal.earth met.(Ca) .....≤ 100 ppm
Bromide.....≤ 100 ppm
Phosphate .....≤ 25 ppm
Heavy metals (Pb).....≤ 5 ppm
Sulphate .....≤ 200 ppm
Al.....≤ 0.2 ppm
As.....≤ 1 ppm
Fe .....≤ 2 ppm
Assay (argentimetric)..99.0 ÷ 100.5 % s.s.

Code	Size	Packaging	Notes
368257000	1 kg	Plastic bottle	
368253000	25 kg	Sack	
368256000	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium chloride > ERBApharm - According to pharmacopoeia: Ph.Eur.-FU-Ph.Franc.-BP-DAB-USP-JP

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Acidity or alkalinity.....Conform Ph.Eur.
Residue solvents .....Conform USP
Barium .....Conform Ph.Eur.
Iodide .....Conform Ph.Eur.
Ferrocyanide.....Conform Ph.Eur.
Nitrite .....Conform Ph.Eur.
Loss on drying ..... ≤ 0.5 %
Mg,alkal.earth met.(Ca) .....≤ 100 ppm
Bromide.....≤ 100 ppm
Phosphate .....≤ 25 ppm
Heavy metals (Pb).....≤ 3 ppm
Sulphate .....≤ 200 ppm
Al.....≤ 0.2 ppm
As.....≤ 1 ppm
Fe .....≤ 2 ppm
K.....≤ 500 ppm
Assay (argentimetric)..99.0 ÷ 100.5 % s.s.

Code	Size	Packaging	Notes
368257	1 kg	Plastic bottle	
368259	5 kg	Plastic tank	
368253	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium chloride 5 mol/l (5N)

• Sodio cloruro 5 mol/l (5N) • Sodium chlorure 5 mol/l (5N) • Sodio cloruro 5 mol/l (5N)
• Natriumchlorid 5 mol/l (5N)

Synonym:
Halite

NaCl
Molecular Weight: 58,44
CAS: 7647-14-5

Sodium chloride 5 mol/l (5N) > RPE - For analysis

RPE

Clear,colourless liquid ..... Conform Assay (potentiometry) ..... 4.990 - 5.010 N

Code	Size	Packaging	Notes
502131	1 l	Plastic bottle	

Content is guaranteed for standardized volumes at 20 °C

Sodium chloride 0.1 mol/l (0,1N)

• Sodio cloruro 0.1 mol/l (0,1N) • Sodium chlorure 0.1 mol/l (0,1N) • Sodio cloruro 0.1 mol/l (0,1N)
• Natriumchlorid 0.1 mol/l (0.1N)

Synonym:
Halite

NaCl
Molecular Weight: 58,44
CAS: 7647-14-5

Sodium chloride 0.1 mol/l (0,1N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
479781		Plastic ampoule	Volume: 55 ml

5,844 g NaCl. Volumetric concentrated solution to prepare 1 L of solution 0,1 N

Sodium citrate dibasic sesquihydrate

• Sodio citrato bibasico • Sodium citrate dibasique • Sodio citrato dibásico
• Dibasisches Natriumcitrat-Sesquihydrat

Synonym:
Citric acid disodium salt

$C_6H_6O_7Na_2.1.5H_2O$
Molecular Weight: 263,1
CAS: 144-33-2
EEC-N: 205-623-3

Sodium citrate dibasic sesquihydrate > ERBApharm - According to pharmacopoeia: BP

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Ready carbonizable substances Conform BP
pH solution 3% ..... 4.9 ÷ 5.2
Chloride ..... ≤330 ppm
Heavy metals (Pb) ..... ≤20 ppm
Oxalate ..... ≤150 ppm
Sulphate ..... ≤0.12 %
As ..... ≤2 ppm
Assay (acidimetric) ..... 98.0 ÷ 104.0 %
Origin (BSE/TSE) ..... Vegetable
Residual solvents (Current ICH) ..... Conform

Code	Size	Packaging	Notes
367951	1 kg	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium citrate tribasic anhydrous

• Sodio citrato tribasico anidro • Sodium citrate tribasique anhydre • Sodio citrato tribásico anhidro • Natriumcitrat tribasisch wasserfrei

$Na_3C_6H_5O_7$
Molecular Weight: 258,12
CAS: 68-04-2
EEC-N: 200-675-3

Sodium citrate tribasic anhydrous > ERBApharm - According to pharmacopoeia: USP

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Alcalinity ..... Conform USP-NF
Tartrate ..... Conform USP-NF
Loss at 180°C ..... ≤1.0 %
Heavy metals (Pb) ..... ≤10 ppm
Assay (protonometric) 99.0 ÷ 100.5 % s.s.

Code	Size	Packaging	Notes
368107	1 kg	Plastic bottle	
368102	20 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Product specifications are subject to changes. Please visit our website for updates.

Sodium citrate tribasic dihydrate

• Sodio citrato tribasico biidrato • Sodium citrate tribasique bihydraté • Sodio citrato tribásico dihidrato
• Natrium-basisches Citrat-Bihydrat

Synonym:
• *Trisodium citrate dihydrate*
• *Citric acid trisodium salt dihydrate*

$Na_3C_6H_5O_7.2H_2O$
Molecular Weight: 294,1
CAS: 6132-04-3
EEC-N: 200-675-3

Sodium citrate tribasic dihydrate > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Reducing iodine ..... Conform
Ready carbonizable substances..... Conform
pH sol. 5% at 25° C ..... 7.5 ÷ 8.7
Ammonium ..... ≤10 ppm
Chloride ..... ≤10 ppm
Total phosphorus ..... ≤10 ppm
Water-insoluble matter ..... ≤30 ppm
Heavy metals (Pb) ..... ≤5 ppm
Oxalate ..... ≤100 ppm
Total sulphur ..... ≤20 ppm
As ..... ≤0.2 ppm
Ca ..... ≤20 ppm
Cu ..... ≤2 ppm
Fe ..... ≤5 ppm
K ..... ≤250 ppm
Ni ..... ≤2 ppm
Pb ..... ≤2 ppm
Zn ..... ≤2 ppm
Assay (non-aqueos medium) ..... ≥99 %

Code	Size	Packaging	Notes
479485	250 g	Plastic bottle	
479487	1 kg	Plastic bottle	
479488	2.5 kg	Plastic bottle	
479486	25 kg	Plastic bucket	
479484	50 kg	Plastic bucket	

Sodium citrate tribasic dihydrate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-BP-DAB

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution ..... Conform Ph.Eur.
Acidity or alkalinity ..... Conform Ph.Eur.
Ready carbonizable substances..... Conform Ph.Eur.
Tartrate ..... Conform USP-NF
Water (K.F.) ..... 11.0 ÷ 13.0 %
Chloride ..... ≤50 ppm
Heavy metals (Pb) ..... ≤10 ppm
Oxalate ..... ≤300 ppm
Sulphate ..... ≤150 ppm
Assay (protonometric) 99.0 ÷ 100.5 % s.s.

Code	Size	Packaging	Notes
368057	1 kg	Plastic bottle	
368058	5 kg	Plastic tank	
368052	10 kg	Plastic tank	
368051	25 kg	Plastic bucket	
368054	50 kg	Fibre drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium cobalt nitrite

• Sodio cobalto nitrito • Sodium cobalt nitrite • Sodio cobalto nitrito • Natriumkobaltnitrit

Synonym:
Sodium hexanitrocobaltate(III)

$Na_3Co(NO_2)_6$
Molecular Weight: 403,94
CAS: 13600-98-1
EEC-N: 237-077-7

Classification transport
ONU: 1479
Transport Hazard class: 5.1
Packing group II

Danger
H272-H315-H319-H334-H317-H351-H335
P210-P280-P284-P304+P340-P305+P351+P338-P342+P311a-P403+P233

Sodium cobalt nitrite > RPE - For analysis

RPE

Description ..... Dark orange powder
Identification ..... Positive
Diluted acetic acid insoluble matter . ≤ 0.02 %
Suitabily for K determ. ..... Conform

Code	Size	Packaging	Notes
479833	50 g	Glass bottle	

Sodium cyanoborohydride

• Sodio cianoboroidruro • Sodium cyanoborohydrure • Sodio cianoborohidruro • Natrium-cyanoborhydrid

Synonym:
Sodium cyanotrihydridoborate

$Na(H_3BCN)$
Molecular Weight: 62,84
CAS: 25895-60-7
EEC-N: 247-317-2

Classification transport
ONU: 1409
Transport Hazard class: 4.3
Packing group I



Danger
H260-H314-HEU032
P223-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sodium cyanoborohydride > RPE - For analysis

RPE

Description ..... White powder | Identification ..... Positive | Assay (iodometric) ..... ≥90 % s s

Code	Size	Packaging	Notes
479371	25 g	Glass bottle	

Sodium 1-decanesulfonate ▶ 1-Decanesulfonic acid sodium salt

Sodium dichloroisocyanurate dihydrate

• Sodio dicloroisocianurato biidrato • Sodium dichloroisocyanurate dihydraté
• Sodio dicloroisocianurato dihidrato • Natriumdichlorisocyanatdihydrat

Synonym:
Dichloroisocyanuric acid sodium salt dihydrate

$C_3Cl_2N_3NaO_3.2H_2O$
Molecular Weight: 255,98
CAS: 51580-86-0
EEC-N: 220-767-7

Classification transport
ONU: 3077
Transport Hazard class: 9
Packing group III



Warning
H302-H319-H335-H410-HEU031
P261-P271-P304+P340-P305+P351+P338-P337+P313-P403+P233

Sodium dichloroisocyanurate dihydrate > RPE - For analysis

RPE

Description ..... White granular powder | Identification ..... Positive | Assay ..... ≥ 98.0 %

Code	Size	Packaging	Notes
479921	10 g	Glass bottle	

Sodium diethyldithiocarbamate trihydrate

• Sodio dietilditiocarbammato triidrato • Sodium diéthyldithiocarbamate trihydraté
• Sodio dietilditiocarbamato trihidrato • Natriumdiethyldithiocarbaminat-Trihydrat

Synonym:
Diethyldithiocarbamic acid sodium salt

$(C_2H_5)_2NCSSNa.3H_2O$
Molecular Weight: 225,23
CAS: 20624-25-3

Classification transport
ONU: 3077
Transport Hazard class: 9
Packing group III



Warning
H302-H315-H319-H335-H400
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Sodium diethyldithiocarbamate trihydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Whyte-yellow crystals | Copper sensitivity ..... Conform | Sodium (Na2SO4) ..... 30.5 ÷ 32.5 %
Identification (I.R.) ..... Conform | Water solubility ..... Conform

Code	Size	Packaging	Notes
405144	100 g	Plastic bottle	

Indicator for the determination of heavy metals

Sodium dihydrogen phosphate dihydrate ▶ Sodium phosphate monobasic dihydrate

Sodium dihydrogen phosphate monohydrate ▶ Sodium phosphate monobasic monohydrate

Sodium dithionite ▶ Sodium hydrosulfite

Sodium 1-dodecanesulfonate ▶ 1-Dodecanesulfonic acid sodium salt

Sodium dodecylbenzenesulfonate ▶ Dodecylbenzenesulphonic acid sodium salt

Product specifications are subject to changes. Please visit our website for updates.

Sodium dodecyl sulfate ▶ Sodium laurylsulfate

Sodium fluoride

• Sodio fluoruro • Sodium fluorure • Sodio fluoruro • Natriumfluorid

NaF
Molecular Weight: 41,99
CAS: 7681-49-4
EEC-N: 231-667-8

Classification transport
ONU: 1690
Transport Hazard class: 6.1
Packing group III

Danger
H301-H315-H319-HEU032
P301+P310a-P330-P305+P351+P338-P362+P364-P332+P313-P337+P313

Sodium fluoride > RPE - For analysis - ACS - ISO

RPE

Description ..... White crystalline powder	Alcalinity ..... ≤0.01 meq/g	Heavy metals (Pb) ..... ≤30 ppm	K ..... ≤200 ppm
Identification ..... Positive	Chloride ..... ≤50 ppm	Sulphate ..... ≤300 ppm	Assay (acidimetric) ..... ≥99 %
Loss on drying ..... ≤0.3 %	Fluosilicates ..... ≤0.1 %	Sulphite ..... ≤50 ppm	
Acidity ..... ≤0.03 meq/g	Water-insoluble matter ..... ≤200 ppm	Fe ..... ≤30 ppm	

Code	Size	Packaging	Notes
479955	250 g	Plastic bottle	
479957	1 kg	Plastic bottle	
479954	25 kg	Plastic bucket	

Sodium fluoride > RE - Pure

RE

Description ..... White crystalline powder	Chloride ..... ≤ 100 ppm	Sulphate ..... ≤ 5000 ppm	Assay ..... ≥ 97 %
Identification ..... Positive	Heavy metals (Pb) ..... ≤ 100 ppm	Fe ..... ≤ 200 ppm	

Code	Size	Packaging	Notes
368457	1 kg	Plastic bottle	
368458	5 kg	Plastic tank	
368451	25 kg	Plastic bucket	

Sodium formate

• Sodio formiato • Sodium formiate • Sodio formiato • Natriumformiat

Synonym:
Formic acid sodium salt

HCOONa
Molecular Weight: 68,01
CAS: 141-53-7
EEC-N: 205-488-0

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Sodium formate > RPE - For analysis - ACS

RPE

Description ..... White crystalline powder	Chloride ..... ≤ 10 ppm	Fe ..... ≤ 5 ppm
Identification ..... Positive	Sulphate ..... ≤ 10 ppm	Heavy metals (Pb) ..... ≤ 5 ppm
Not soluble matter ..... ≤ 50 ppm	Ca ..... ≤ 50 ppm	Assay (oxidimetric) ..... ≥ 99.0 %

Code	Size	Packaging	Notes
480045	100 g	Plastic bottle	
480046	500 g	Plastic bottle	

Sodium gluconate

• Sodio gluconato • Sodium gluconate • Sodio gluconato • Natriumgluconat

Synonym:
• *2,3,4,5,6-Pentahydroxycaproic acid sodium salt*
• *D-gluconate sodium salt*

$C_6H_{11}NaO_7$
Molecular Weight: 218,13
CAS: 527-07-1
EEC-N: 208-407-7

Sodium gluconate > RE - Pure

RE

Description ..... White powder
Identification ..... Positive
Chloride ..... ≤500 ppm
Red.ing sugars(Glucose) ..... ≤1 %
Assay (non-aqueos medium) ..... ≥98 %

Code	Size	Packaging	Notes
369582	1 kg	Plastic bottle	
369581	10 kg	Carton box	

Sodium glutamate acid

• Sodio glutammato acido • Sodium glutamate • Sodio glutamato acido • Natriumglutamatsäure

Synonym:
• *L-Glutamic acid monosodium salt hydrate*
• *L-2-Aminopentanedioic acid*

$C_5H_8NO_4Na.H_2O$
Molecular Weight: 187
CAS: 142-47-2
EEC-N: 205-538-1

Sodium glutamate acid > RE - Pure

RE

Description ..... White cryst. need.sha
Identification ..... Positive
Potere rotator. specif.(C=1 HCl 6N).. +23 ÷ +25.3 °
Assay (ex nitrogen) ..... ≥ 98 %

Code	Size	Packaging	Notes
369667	1 kg	Plastic bottle	
369663	25 kg	Plastic bucket	

Sodium glycerophosphate pentahydrate

• Sodio glicerofosfato pentaidrato • Sodium glycérophosphate pentahydraté
• Sodio glicerofosfato pentahidrato • Natriumglycerophosphatpentahydrat

Synonym:
• *Glycerol-2-phosphate disodium salt hydrate*
• *BGP*

$C_3H_7O_6PNa_2.5H_2O$
Molecular Weight: 306
CAS: 13408-09-8

Sodium glycerophosphate pentahydrate > ERBApharm - According to pharmacopoeia: Ph.Eur.

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Phosphate ..... ≤ 0.1 %
Heavy metals (Pb) ..... ≤ 20 ppm
Fe ..... ≤ 20 ppm
Appearance of solution ..... Conform Ph.Eur.
Alcalinity ..... Conform Ph.Eur.
Water (K.F.) ..... 25.0 ÷ 35.0 %
Chloride ..... ≤ 0.02 %
Glycerol and alcohol-soluble substances .. ≤ 1.0 %
Sulphate ..... ≤ 0.05 %
Assay ..... 98.0 ÷ 105.0 % anidro

Code	Size	Packaging	Notes
369447	1 kg	Plastic bottle	
369449	5 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium 1-heptanesulfonate ▶ 1-Heptanesulphonic acid sodium salt

Product specifications are subject to changes. Please visit our website for updates.

Sodium hexafluorosilicate

• Sodio esafluorosilicato • Sodium hexafluorosilicate • Sodio hexafluorosilicato • Natriumhexafluorosilikat

Na_2SiF_6
Molecular Weight: 188,06
CAS: 16893-85-9
EEC-N: 240-934-8

Classification transport
ONU: 2674
Transport Hazard class: 6.1
Packing group III

Danger
H301-H311-H331
P261-P304+P340-P311a-P330-P361+P364-P403+P233

Sodium hexafluorosilicate > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Chloride ..... ≤200 ppm
Sulphate ..... ≤200 ppm
Heavy metals (Pb) ..... ≤50 ppm
Fe ..... ≤50 ppm
Assay (acidimetric) ..... ≥98.5 %

Code	Size	Packaging	Notes
480005	250 g	Plastic bottle	

Sodium hexametaphosphate

• Sodio esametafosfato • Sodium hexamétaphosphate • Sodio hexametafosfato • Natriumhexametaphosphat

Synonym:
Sodium metaphosphate

$(NaPO_3)_6$
Molecular Weight: 611,76
CAS: 10124-56-8
EEC-N: 233-343-1

Sodium hexametaphosphate > RE - Pure

RE

Description ..... white hygroscopic powder
Identification ..... Positive
pH (1% solution) ..... 5.8 ÷ 6.5
Water insoluble substances ..... ≤ 0.06 %
Fe ..... ≤500 ppm
Assay (as P203) ..... ≥ 68 %

Code	Size	Packaging	Notes
368357	1 kg	Plastic bottle	
368351	10 kg	Plastic tank	
368352	25 kg	Plastic bucket	

Sodium 1-hexanesulfonate ▶ 1-Hexanesulphonic acid sodium salt

Sodium 1-hexanesulfonate monohydrate ▶ 1-Hexanesulphonic acid sodium salt monohydrate

Sodium hexanitrocobaltate(III) ▶ Sodium cobalt nitrite

Sodium hydrogen carbonate ▶ Sodium bicarbonate

di-Sodium hydrogen phosphate ▶ Sodium phosphate dibasic anhydrous

di-Sodium hydrogen phosphate dihydrate ▶ Sodium phosphate dibasic dihydrate

di-Sodium hydrogen phosphate dodecahydrate ▶ Sodium phosphate dibasic dodecahydrate

Sodium hydrogen sulfate monohydrate ▶ Sodium bisulfate monohydrate

Sodium hydrogen tartrate monohydrate

• Sodio tartrato acido • Sodium tartrate acide • Sodio tartrato acido • Natriumhydrogentartratmonohydrat

Synonym: *Sodium bitartrate monohydrate*

$NaOOC(CHOH)_2COOH.H_2O$
Molecular Weight: 190,09
CAS: 526-94-3
EEC-N: 208-400-9

Sodium hydrogen tartrate monohydrate > RPE - For analysis

RPE

Description	White crystals	Chloride	≤ 10 ppm	As	≤ 0.4 ppm	Ni	≤ 2 ppm
Identification	Positive	Total phosphorus	≤ 10 ppm	Ca	≤ 100 ppm	Pb	≤ 2 ppm
pH sol. 5% at 25° C	3.30 ÷ 3.60	Water-insoluble matter	≤ 50 ppm	Cu	≤ 2 ppm	Zn	≤ 2 ppm
Loss on drying	9 ÷ 10 %	Heavy metals (Pb)	≤ 10 ppm	Fe	≤ 10 ppm	Assay (acidimetric)	≥ 99.5 %
Ammonium	≤ 50 ppm	Total sulphur	≤ 50 ppm	K	≤ 100 ppm		

Code	Size	Packaging	Notes
483706	500 g	Plastic bottle	
483703	25 kg	Fibre drum	

Sodium hydrosulfite

• Sodio idrosolfito • Sodium hydrosulfite • Sodio hidrosulfito • Natriumdithionit

Synonym: *Sodium dithionite*

$Na_2S_2O_4$
Molecular Weight: 174,11
CAS: 7775-14-6
EEC-N: 231-890-0

Classification transport
ONU: 1384
Transport Hazard class: 4.2
Packing group II



Danger
H251-H319-HEU031
P235-P264-P280-P305+P351+P338-P337+P313-P420

Sodium hydrosulfite > RE - Pure

RE

Description	White crystalline powder	Identification	Positive	Assay (oxidimetric)	≥ 80 %

Code	Size	Packaging	Notes
370011	1 kg	Metallic can	
370014	2.5 kg	Metallic can	

Sodium hydroxide, pearls

• Sodio idrossido, perline • Sodium hydroxyde, perles • Sodio hidróxido, perlas • Natriumhydroxid, Perlen

Synonym: *Caustic soda*

NaOH
Molecular Weight: 40
CAS: 1310-73-2
EEC-N: 215-185-5

Classification transport
ONU: 1823
Transport Hazard class: 8
Packing group II



Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sodium hydroxide, pearls > ERBApharm - According to pharmacopoeia: Ph.Eur.-NF

ERBApharm

Description	Pearls white	pH	≥ 11.0	Sulphate	≤ 200 ppm	Origin (BSE/TSE)	Synthesis
Identification	Positive	Carbonate	≤ 2.0 %	Fe	≤ 10 ppm	Residual solvents (Current ICH)	Conform
Appearance of solution	Conform Ph. Eur.	Chloride	≤ 200 ppm	Assay (total alkalinity)	97.0 ÷ 100.5 %		
Not sol.matter,org.mat	Conform NF	Heavy metals (Pb)	≤ 20 ppm	Content of sodium	54.0 ÷ 59.8 %		

Code	Size	Packaging	Notes
369743	1 kg	Plastic bottle	
369741	5 kg	Plastic tank	
369742	25 kg	Sack	
369744	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Product specifications are subject to changes. Please visit our website for updates.

Sodium hydroxide, pellets

• Sodio idrossido, gocce • Sodium hydroxyde, pastilles • Sodio hidróxido, lentejas • Natriumhydroxid, Pellets

Synonym: *Caustic soda*

NaOH
Molecular Weight: 40
CAS: 1310-73-2
EEC-N: 215-185-5

Classification transport
ONU: 1823
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sodium hydroxide, pellets > RS - RSE - For electronic use

RS

Description	White pellets	Heavy metals (Pb)	≤2 ppm	Cd	≤0.1 ppm	Mn	≤0.5 ppm
Identification	Positive	Silicate	≤20 ppm	Cu	≤0.5 ppm	Ni	≤2 ppm
Total nitrogen	≤3 ppm	Sulphate	≤5 ppm	Fe	≤3 ppm	Pb	≤1 ppm
Carbonate	≤5000 ppm	Al	≤2 ppm	Hg	≤0.1 ppm	Zn	≤1 ppm
Chloride	≤10 ppm	As	≤0.5 ppm	K	≤100 ppm	Assay (alkalimetric)	≥98.5 %
Phosphate	≤5 ppm	Ca	≤5 ppm	Mg	≤5 ppm		

Code	Size	Packaging	Notes
480527	1 kg	Plastic bottle	
480522	5 kg	Plastic jar	
480525	25 kg	Plastic bucket	

Sodium hydroxide, pellets > RPE - For analysis - ACS - ISO

RPE

Description	White pellets	Chloride	≤ 50 ppm	Ca	≤ 50 ppm	Mg	≤ 20 ppm
Identification	Positive	Phosphate	≤ 10 ppm	Fe	≤ 10 ppm	Ni	≤ 10 ppm
Total nitrogen	≤ 10 ppm	Sulphate	≤ 30 ppm	Hg	≤ 0.1 ppm	Assay (alkalimetric)	≥ 97.0 %
Carbonate	≤ 1.0 %	Heavy metals (Ag)	≤ 20 ppm	K	≤ 0.02 %		

Code	Size	Packaging	Notes
480505	100 g	Plastic bottle	
480501	500 g	Plastic bottle	
480507	1 kg	Plastic bottle	
480509	5 kg	Plastic jar	
480502	10 kg	Plastic tank	
480508	25 kg	Plastic bucket	

Sodium hydroxide, pellets > ERBApharm - According to pharmacopoeia: BP-NF-Ph.Eur.-JP

ERBApharm

Description	white or almost white pellets	pH	≥11.0	Fe	≤10 ppm	Assay (total alkalinity)	97.0 ÷ 100.5 %
Identification	Positive	Carbonate	≤2.0 %	K	≤ 0.5 %	Content of sodium	54.0 - 59.8 %
Appearance of solution	Conform Ph.Eur.	Chloride	≤200 ppm	Sulphate	≤200 ppm	Origin (BSE/TSE)	Synthesis
Not sol.matter,org.mat.	Conform USP-NF	Heavy metals (Pb)	≤ 30 ppm (JP)	Hg	≤ 0.1 ppm	Residual solvents (Current ICH)	Conform

Code	Size	Packaging	Notes
369777	1 kg	Plastic bottle	
369772	5 kg	Plastic tank	
369771	20 kg	Plastic bucket	
369774	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium hydroxide on silica

• Sodio idrossido su silice • Sodium hydroxyde sur silice • Sodio hidróxido sobre sílice
• Natriumhydroxid auf Kieselgel

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2

Classification transport
ONU: 1823
Transport Hazard class: 8
Packing group II



Danger
H314-H335
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338-P403+P233

Sodium hydroxide on silica > RS - For microanalysis

RS

Description .....Dark grey granules Identification .....Positive Average grain diameter.....1.6 ÷ 3 mm ca. CO2 absorption .....≥30 %

Code	Size	Packaging	Notes
424494	100 g	Glass bottle	
424497	1 kg	Plastic bottle	

Sodium hydroxide solution 50%

• Sodio idrossido soluzione 50% • Sodium hydroxyde solution 50% • Sodio hidróxido solución 50%
• Natronlauge 50 %

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2
EEC-N: 215-185-5

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II



Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Sodium hydroxide solution 50% > RPE - For analysis

RPE

Assay .....48 - 51.5 % Density d20/4 .....1.5 - 1.54 NaCl .....≤ 0.02 %
Iron (Fe).....≤ 15 mg/Kg Na2CO3 .....≤ 0.7 %

Code	Size	Packaging	Notes
P4540041	10 l	Plastic tank	
P4540049	25 l	Plastic tank	

Sodium hydroxide solution 40%

• Sodio idrossido soluzione 40% • Sodium hydroxyde solution 40% • Sodio hidróxido solución 40%
• Natronlauge 40%

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2
EEC-N: 215-185-5

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II



Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Sodium hydroxide solution 40% > RS - For agroalimentary analysis

RS

Colour .....≤ 10 APHA Density at 20°C .....1.420 ÷ 1.440 Total nitrogen .....≤ 10 ppm
Description .....Clear colourless liquid Assay .....39.0 ÷ 41.0 %

Code	Size	Packaging	Notes
502721	5 l	Plastic tank	
502722	10 l	Plastic tank	

572 g of NaOH for 1 L

Product specifications are subject to changes. Please visit our website for updates.

Sodium hydroxide solution 35-37%

• Sodio idrossido soluzione 35-37% • Sodium hydroxyde solution 35-37% • Sodio hidróxido solución 35-37%
• Natronlauge 35-37%

Synonym: *Caustic soda*

NaOH
Molecular Weight: 40
CAS: 1310-73-2
EEC-N: 215-185-5

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Sodium hydroxide solution 35-37% > RS - For agroalimentary analysis

RS

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Density at 20°C ..... 1.380 ÷ 1.400
Nitrogen compounds ..... ≤ 10 ppm

Code	Size	Packaging	Notes
502112	5 l	Plastic tank	

500g of NaOH for 1L . According to normative T90-110

Sodium hydroxide solution 35%

• Sodio idrossido soluzione 35% • Sodium hydroxyde solution 35% • Sodio hidróxido solución 35%
• Natronlauge 35 %

Synonym: *Caustic soda*

NaOH
Molecular Weight: 40
CAS: 1310-73-2
EEC-N: 215-185-5

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Sodium hydroxide solution 35% > RPE - For analysis

RPE

Density d20/4 ..... 1.38 - 1.39
Description ..... Clear colourless liquid
Assay ..... 35 - 36 %
Density at 20° C ..... 1.38 - 1.39
Assay (alkalimetric) ..... 35 - 36 %

Code	Size	Packaging	Notes
480591	1 l	Plastic bottle	
480593	25 kg	Plastic tank	

Sodium hydroxide solution 32%

• Sodio idrossido soluzione 32% • Sodium hydroxyde solution 32% • Sodio hidróxido solución 32%
• Natronlauge 32%

Synonym: *Caustic soda*

NaOH
Molecular Weight: 40
CAS: 1310-73-2
EEC-N: 215-185-5

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Sodium hydroxide solution 32% > RS - For agroalimentary analysis

RS

Assay ..... 30 - 34 %
Density at 20°C ..... 1.322 - 1.374
Fe ..... ≤ 10 ppm

Code	Size	Packaging	Notes
524510	25 kg	Plastic tank	

Sodium hydroxide solution 32% > RS - For Kjeldahl

RS

Description ..... Slightly opalescent liquid
Colour ..... ≤ 10 APHA
Assay ..... 31.0 - 33.0 %
Density at 20°C ..... 1.339 - 1.359
Nitrogen compounds ..... ≤ 1 ppm

Code	Size	Packaging	Notes
480561	1 l	Plastic bottle	
480566	2.5 l	Plastic bottle	
526521	5 l	Plastic tank	
480564	10 l	Plastic tank	
480562	25 kg	Plastic tank	
480563	30 kg	Plastic tank	

Sodium hydroxide solution 30%

• Sodio idrossido soluzione 30% • Sodium hydroxyde solution 30% • Sodio hidróxido solución 30%
• Natronlauge 30%

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2
EEC-N: 215-185-5

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Sodium hydroxide solution 30% > RS - For agroalimentary analysis

RS

Description ..... Slightly opalescent liquid
Colour ..... ≤ 10 APHA
Density at 20°C ..... ≥ 1.323

Code	Size	Packaging	Notes
502741	5 l	Plastic tank	

Sodium hydroxide solution 30% > RPE - For analysis

RPE

Na2CO3 ..... ≤ 0.5 %
Silicate (as SiO2) ..... ≤ 30 mg/Kg
Chloride (Cl-) ..... ≤ 10 mg/Kg
Rapprochement BA/FT ..... Conform
Sodium hydroxyde content ..... 30 - 31 %

Code	Size	Packaging	Notes
P1050552	30 l	Plastic tank	

Sodium hydroxide solution 30% > RPE - For nitrogen dosing

RPE

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Density at 20°C ..... 1.323 ÷ 1.333
Total nitrogen ..... ≤ 1 ppm
Assay ..... 29.5 ÷ 30.5 %

Code	Size	Packaging	Notes
502731	1 l	Plastic bottle	

Sodium hydroxide solution 30% > ERBApharm - Prepared from raw material according Ph.Eur

ERBApharm

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Identification ..... Positive
Density at 20°C ..... 1.311 ÷ 1.344
Chloride ..... ≤ 200 ppm
Sulphate ..... ≤ 200 ppm
Iron ..... ≤10 ppm
Heavy metals (Pb) ..... ≤ 20 ppm
Carbonate ..... ≤ 0.6 %
Assay ..... 29.5 ÷ 30.5 %

Code	Size	Packaging	Notes
369704	1 l	Plastic bottle	
369702	20 l	Plastic tank	
369701000	10 kg	Plastic tank	
369706	200 l	Plastic drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium hydroxide solution 30% > RE - Pure

RE

Description ..... Opalescent liquid
Identification ..... Positive
Density at 20° C ..... 1.306 - 1.349
Alcalinity (NaOH) ..... 28.0 ÷ 32.0 %
Carbonate ..... ≤1.2 %
Chloride ..... ≤300 ppm
Sulphate ..... ≤150 ppm
Fe ..... ≤30 ppm

Code	Size	Packaging	Notes
369762	10 kg	Plastic tank	
369761	30 kg	Plastic tank	
369766	50 kg	Plastic tank	

Product specifications are subject to changes. Please visit our website for updates.

Sodium hydroxide solution 20% w/v

• Sodio idrossido soluzione 20% p/v • Sodium hydroxyde solution 20% m/v
• Sodio hidróxido solución 20% p/v • Natronlauge 20% w/v

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2
EEC-N: 215-185-5

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sodium hydroxide solution 20% w/v > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611081401	1 l	Plastic bottle	Ref Ph.Eur 1081401

Sodium hydroxide solution 20% w/v > RPE - For analysis

RPE

Description ..... Slightly opalescent liquid Colour ..... ≤ 10 APHA Assay (NaOH) ..... 19.9 ÷ 20.1 %m/v

Code	Size	Packaging	Notes
524505	10 l	Plastic tank	

Sodium hydroxide solution 20% w/w

• Sodio idrossido soluzione 20% p/p • Sodium hydroxyde solution à 20% m/m • Sodio hidróxido solución 20% p/p • Natronlauge 20% w/w

NaOH
Molecular Weight: 40
CAS: 1310-73-2
EEC-N: 215-185-5

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sodium hydroxide solution 20% w/w > RPE - For analysis

RPE

Description ..... Clear colourless liquid Density at 20° C ..... 1.213 - 1.225 Assay (alkalimetric) ..... 19.5 - 20.5 % Carbonate ..... ≤1 %

Code	Size	Packaging	Notes
480621	1 l	Plastic bottle	
480622	30 kg	Plastic tank	

Sodium hydroxide solution 20% w/w > ERBAPharm - Prepared from raw material according Ph.Eur

ERBApharm

Description ..... Clear colourless liquid Assay (alkalimetric) ..... 19.5 - 20.5 % Origine (BSE-TSE) ..... Conform
Density at 20° C ..... 1.213 - 1.225 Carbonate ..... ≤1 % Residual solvents (Current ICH) ..... Conform

Code	Size	Packaging	Notes
480631	1 l	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium hydroxide solution 10% w/v

• Sodio idrossido soluzione 10% p/v • Sodium hydroxyde solution 10% m/v
• Sodio hidróxido solución 10% p/v • Natronlauge 10% w/v

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2
EEC-N: 215-185-5

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sodium hydroxide solution 10% w/v > RS - For agroalimentary analysis

RS

Description ..... Slightly opalescent liquid Identification ..... Positive Assay (NaOH) ..... 9 - 11 %m/v

Code	Size	Packaging	Notes
508615	5 l	Plastic tank	

Sodium hydroxide solution 10% w/v > RPE - For analysis

RPE

Description ..... Slightly opalescent liquid Colour ..... ≤ 10 APHA Assay (NaOH) ..... 9.9 ÷ 10.1 %m/v

Code	Size	Packaging	Notes
524506	5 l	Plastic tank	
524507	10 l	Plastic tank	

Sodium hydroxide solution 10% w/v > RE - Pure

RE

Description ..... Slightly opalescent liquid Colour ..... ≤ 10 APHA Assay (NaOH) ..... 9.5 - 10.5 %m/v

Code	Size	Packaging	Notes
526642	5 l	Plastic tank	
526641	10 l	Plastic tank	
526644	50 l	Plastic tank	

Sodium hydroxide solution 5% w/v

• Sodio idrossido soluzione 5% p/v • Sodium hydroxyde solution 5% m/v • Sodio hidróxido solución 5% w/v • Natronlauge 5% w/v

Synonym: *Caustic soda*

NaOH
Molecular Weight: 40
CAS: 1310-73-2
EEC-N: 215-185-5

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sodium hydroxide solution 5% w/v > RPE - For analysis

RPE

Description ..... Sightly opalescent liquid Colour ..... ≤ 10 APHA Assay (NaOH) ..... 4.9 ÷ 5.1 %m/v

Code	Size	Packaging	Notes
524502	5 l	Plastic tank	
524501	10 l	Plastic tank	

Sodium hydroxide solution 5% w/v > RE - Pure

RE

Description ..... Sightly opalescent liquid Colour ..... ≤ 10 APHA Assay (NaOH) ..... 4.5 - 5.5 %m/v

Code	Size	Packaging	Notes
526632	5 l	Plastic tank	
526634	10 l	Plastic tank	

Sodium hydroxide 6 mol/l (6N)

• Sodio idrossido 6 mol/l (6N) • Sodium hydroxyde 6 mol/l (6N) • Sodio hidróxido 6 mol/l (6N) • Natronlauge 6 mol/l (6N)

NaOH
Molecular Weight: 40
CAS: 1310-73-2
EEC-N: 215-185-5

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sodium hydroxide 6 mol/l (6N) > ERBApharm - Prepared from raw material according Ph.Eur

ERBApharm

Assay (potentiometry) ..... 5.982 - 6.018 N Colour ..... ≤ 100 Hazen

Code	Size	Packaging	Notes
524651	2 l	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Product specifications are subject to changes. Please visit our website for updates.

Sodium hydroxide 5 mol/l (5N)

• Sodio idrossido 5 mol/l (5N) • Sodium hydroxyde 5 mol/l (5N) • Sodio hidróxido 5 mol/l (5N)
• Natronlauge 5 mol/l (5N)

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2
EEC-N: 215-185-5

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Sodium hydroxide 5 mol/l (5N) > RS - For agroalimentary analysis

RS

Description ..... Clear colourless liquid Assay ..... 4.995 ÷ 5.005 N

Code	Size	Packaging	Notes
526513	1 l	Plastic bottle	
526512	5 l	Plastic tank	

Content is guaranteed for standardized volumes at 20°C

Sodium hydroxide 4 mol/l (4N)

• Sodio idrossido 4 mol/l (4N) • Sodium hydroxyde 4 mol/l (4N) • Sodio hidróxido 4 mol/l (4N) • Natronlauge 4 mol/l (4N)

NaOH
Molecular Weight: 40
CAS: 1310-73-2
EEC-N: 215-185-5

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Sodium hydroxide 4 mol/l (4N) > RPE - For agroalimentary analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 3.992 - 4.008 N

Code	Size	Packaging	Notes
502662	2.5 l	Plastic bottle	
502664	10 l	Plastic tank	

Sodium hydroxide 3 mol/l (3N)

• Sodio idrossido 3 mol/l (3N) • Sodium hydroxyde 3 mol/l (3N) • Sodio hidróxido 3 mol/l (3N)
• Natronlauge 3 mol/l (3N)

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2
EEC-N: 215-185-5

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Sodium hydroxide 3 mol/l (3N) > ERBApharm - Prepared from raw material according Ph.Eur

ERBApharm

Identification (Ph.Eur)..... Conform Assay (Ph.Eur)..... 2.85 - 3.15 N Origine (BSE-TSE)..... Conform Residual solvents (Current ICH)..... Conform

Code	Size	Packaging	Notes
524732	500 ml	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium hydroxide 2 mol/l (2N)

• Sodio idrossido 2 mol/l (2N) • Sodium hydroxyde 2 mol/l (2N) • Sodio hidróxido 2 mol/l (2N)
• Natronlauge 2 mol/l (2N)

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2
EEC-N: 215-185-5

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Sodium hydroxide 2 mol/l (2N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 1.998 - 2.002 N NIST 84 I

Code	Size	Packaging	Notes
480686000	500 ml	Plastic bottle	Certified with NIST traceability
480687000	1 l	Plastic bottle	Certified with NIST traceability
480682000	5 l	Plastic tank	Certified with NIST traceability
480681000	10 l	Kubidos	Certified with NIST traceability
480684000	20 l	Plastic tank	Certified with NIST traceability

80 g de NaOH. Volumetric solution ready-to-use

Sodium hydroxide 2 mol/l (2N) > ERBApharm - Prepared from raw material according Ph.Eur

ERBApharm

Identification (Ph.Eur)..... Conform Assay (Ph.Eur)..... 1.9 -2.1 N Origine (BSE-TSE)..... Conform Residual solvents (Current ICH)..... Conform

Code	Size	Packaging	Notes
524671	1 l	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium hydroxide 1.2 mol/l (1.2N)

• Sodio idrossido 1.2 mol/l (1.2N) • Sodium hydroxyde 1.2 mol/l (1.2N) • Sodio hidróxido 1.2 mol/l (1.2N)
• Natronlauge 1.2 mol/l (1.2N)

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H314
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Sodium hydroxide 1.2 mol/l (1.2N) > RS - For analysis

RS

Assay (potentiomety) ..... 1.1976 - 1.2024 N

Code	Size	Packaging	Notes
PS0736/41	10 l	Plastic tank	
PS0736/42	20 l	Plastic tank	
PS0736/49	25 l	Plastic tank	

Sodium hydroxide 1 mol/l (1N)

• Sodio idrossido 1 mol/l (1N) • Sodium hydroxyde 1 mol/l (1N) • Sodio hidróxido 1 mol/l (1N)
• Natronlauge 1 mol/l (1N)

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H314
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Sodium hydroxide 1 mol/l (1N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613006301	500 ml	Plastic bottle	Ref Ph.Eur 3006300
613006300	1 l	Plastic bottle	Ref Ph.Eur 3006300

Product specifications are subject to changes. Please visit our website for updates.

Sodium hydroxide 1 mol/l (1N) > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
617000121	500 ml	Plastic bottle	

Sodium hydroxide 1 mol/l (1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Assay (potentiometry) ..... 0.999 - 1.001 N NIST 84 I

Code	Size	Packaging	Notes
480717000	1 l	Plastic bottle	Certified with NIST traceability
480711000	5 l	Kubidos	Certified with NIST traceability
480714000	5 l	Plastic tank	Certified with NIST traceability
480713000	10 l	Kubidos	Certified with NIST traceability

40 g de NaOH. Volumetric solution ready-to-use

Sodium hydroxide 1 mol/l (1N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
480741		Plastic ampoule	Volume: 165 ml

40 g NaOH.Volumetric concentrated solution to prepare 1 L of solution 1 N

Sodium hydroxide 1 mol/l (1N) > ERBApharm - Prepared from raw material according USP - NF

ERBApharm

Identification (USP) ..... Conform Origine (BSE-TSE) ..... Conform
Assay (USP) ..... 0.95 - 1.05 N Residual solvents (Current ICH) ..... Conform

Code	Size	Packaging	Notes
524761	1 l	Plastic bottle	

Sodium hydroxide 1 mol/l (1N) > ERBApharm - Prepared from raw material according Ph.Eur

ERBApharm

Identification (Ph.Eur) ..... Conform Origine (BSE-TSE) ..... Conform
Assay (Ph.Eur) ..... 0.95 - 1.05 N Residual solvents (Current ICH) ..... Conform

Code	Size	Packaging	Notes
524621	1 l	Plastic bottle	

Sodium hydroxide 0.7 mol/l (N/1.4)

• Sodio idrossido 0.7 mol/l (N/1.4) • Sodium hydroxyde 0.7 mol/l (N/1.4) • Sodio hidróxido 0.7 mol/l (N/1.4)
• Natronlauge 0.7 mol/l (N/1.4)

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sodium hydroxide 0.7 mol/l (N/1.4) > RS - For agroalimentary analysis

RS

Description ..... Clear colourless liquid Colour ..... ≤ 10 APHA Assay ..... 0.710 ÷ 0.718 N

Code	Size	Packaging	Notes
526511	10 l	Kubidos	

Sodium hydroxide 0.5 mol/l (0.5N)

• Sodio idrossido 0.5 mol/l (0.5N) • Sodium hydroxyde 0.5 mol/l (0.5N) • Sodio hidróxido 0.5 mol/l (0.5N)
• Natronlauge 0.5 mol/l (0.5N)

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sodium hydroxide 0.5 mol/l (0.5N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.4995 - 0.5005 N NIST 84 I

Code	Size	Packaging	Notes
480777000	1 l	Plastic bottle	Certified with NIST traceability
480771000	5 l	Kubidos	Certified with NIST traceability
480772000	10 l	Kubidos	Certified with NIST traceability
480773000	10 l	Plastic tank	Certified with NIST traceability

20 g of NaOH. Volumetric solution ready-to-use

Sodium hydroxide 0.5 mol/l (0.5N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
480801		Plastic ampoule	Volume: 55 ml

20 g of NaOH. Volumetric concentrated solution to prepare 1 L of solution 0,5 N

Sodium hydroxide 0.357 mol/l (0.357N)

• Sodio idrossido 0.357 mol/l (0.357N) • Sodium hydroxyde 0.357 mol/l (0.357N)
• Sodio hidróxido 0.357 mol/l (0.357N) • Natronlauge 0.357 mol/l (0.357N)

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group III

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Sodium hydroxide 0.357 mol/l (0.357N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.3566 - 0.3574 N NIST 84 I

Code	Size	Packaging	Notes
480837000	1 l	Plastic bottle	Certified with NIST traceability

14.28 g of NaOH. Volumetric solution ready-to-use: 1/2.82N

Sodium hydroxide 0.25 mol/l (0.25N)

• Sodio idrossido 0.25 mol/l (0.25N) • Sodium hydroxyde 0.25 mol/l (0.25N)
• Sodio hidróxido 0.25 mol/l (0.25N) • Natronlauge 0.25 mol/l (0.25N)

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group III

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Sodium hydroxide 0.25 mol/l (0.25N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.2498 - 0.2503 N NIST 84 I

Code	Size	Packaging	Notes
480867000	1 l	Plastic bottle	Certified with NIST traceability
480861000	5 l	Kubidos	Certified with NIST traceability
480862000	10 l	Kubidos	Certified with NIST traceability
480863000	25 l	Plastic tank	Certified with NIST traceability

10 g of NaOH. Volumetric solution ready-to-use

Product specifications are subject to changes. Please visit our website for updates.

Sodium hydroxide 0.25 mol/l (0.25N) > ERBApharm - Prepared from raw material according Ph.Eur

ERBApharm

Colour ..... ≤ 100 APHA Assay ..... 0.2495 ÷ 0.2505 N

Code	Size	Packaging	Notes
369812	5 l	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade. Content is guaranteed for standardized volumes at 20°C. Keep tightly closed

Sodium hydroxide 0.2 mol/l (0.2N)

• Sodio idrossido 0.2 mol/l (0.2N) • Sodium hydroxyde 0.2 mol/l (0.2N) • Sodio hidróxido 0.2 mol/l (0.2N)
• Natronlauge 0.2 mol/l (0.2N)

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group III

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Sodium hydroxide 0.2 mol/l (0.2N) > RS - For agroalimentary analysis

RS

Description ..... Clear colourless liquid Colour ..... ≤ 10 APHA Assay ..... 0.19 ÷ 0.21 M

Code	Size	Packaging	Notes
502782	500 ml	Plastic bottle	
502781000	10 l	Plastic tank	

Sodium hydroxide 0.2 mol/l (0.2N) > RPE - For analysis

RPE

Assay (potentiometry) ..... 0.1998 - 0.2002 N

Code	Size	Packaging	Notes
P3440015	1 l	Plastic bottle	

Sodium hydroxide 1/9 mol/l (N/9)

• Sodio idrossido 1/9 mol/l (N/9) • Sodium hydroxyde 1/9 mol/l (N/9) • Sodio hidróxido 1/9 mol/l (N/9)
• Natronlauge 1/9 mol/l (N/9)

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group III

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Sodium hydroxide 1/9 mol/l (N/9) > RPE - For analysis

RPE

Assay (potentiometry) ..... 0.1109 - 0.1113 N

Code	Size	Packaging	Notes
P4500022	5 l	Plastic tank	

Sodium hydroxide 0.1 mol/l (0.1N)

• Sodio idrossido 0.1 mol/l (0.1N) • Sodium hydroxyde 0.1 mol/l (0.1N) • Sodio hidróxido 0.1 mol/l (0.1N)
• Natronlauge 0.1 mol/l (0.1N)

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2

Sodium hydroxide 0.1 mol/l (0.1N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613006601	500 ml	Plastic bottle	Ref Ph.Eur 3006600
613006600	1 l	Plastic bottle	Ref Ph.Eur 3006600

Sodium hydroxide 0.1 mol/l (0.1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.0999 - 0.1001 N NIST 84 I

Code	Size	Packaging	Notes
480897000	1 l	Plastic bottle	Certified with NIST traceability
480891000	5 l	Kubidos	Certified with NIST traceability
480892000	10 l	Kubidos	Certified with NIST traceability
480893000	10 l	Plastic tank	Certified with NIST traceability
480895000	50 l	Plastic drum	Certified with NIST traceability

4 g of NaOH. Volumetric solution ready-to-use

Sodium hydroxide 0.1 mol/l (0.1N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
480921		Plastic ampoule	Volume: 55 ml

4 g of NaOH. Volumetric concentrated solution to prepare 1 L of solution 0,1 N

Sodium hydroxide 0.1 mol/l (0.1N) > ERBApharm - Prepared from raw material according Ph.Eur

ERBApharm

Identification (Ph.Eur) ..... Conform Assay (Ph.Eur) ..... 0.095 - 0.105 N Origine (BSE-TSE) ..... Conform Residual solvents (Current ICH) ..... Conform

Code	Size	Packaging	Notes
524631	1 l	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium hydroxide 0.1 mol/l (0.1N) in ethanol

• Sodio idrossido 0.1 mol/l (N/10) in etanolo • Sodium hydroxyde 0.1 mol/l (0.1N) dans l'éthanol
• Sodio hidróxido 0.1 mol/l (N/10) en etanol • Natronlauge 0.1 mol/l (0.1N) in Ethanol

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II



Danger
H225-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

Sodium hydroxide 0.1 mol/l (0.1N) in ethanol > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613007001	100 ml	Glass bottle	Ref Ph.Eur 3007000
613007000	1 l	Glass bottle	Ref Ph.Eur 3007000

Sodium hydroxide 0.01 mol/l (0.01N)

• Sodio idrossido 0.01 mol/l (0.01N) • Sodium hydroxyde 0.01 mol/l (0.01N)
• Sodio hidróxido 0.01 mol/l (0.01N) • Natronlauge 0.01 mol/l (0.01N)

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2

HEU210

Sodium hydroxide 0.01 mol/l (0.01N) > RPE - For analysis

RPE

Assay (potentiometry) ..... 0.00998 - 0.01002 N

Code	Size	Packaging	Notes
PS0215/15	1 l	Plastic bottle	

Sodium hydroxide 0.01 mol/l (0.01N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
481001		Plastic ampoule	Volume: 55 ml

0,4 g of NaOH. Volumetric concentrated solution to prepare 1 L of solution 0,01 N

Product specifications are subject to changes. Please visit our website for updates.

Sodium hydroxide solution

• Sodio idrossido soluzione • Sodium hydroxyde solution • Sodio hidróxido solución • Natronlauge

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2
EEC-N: 215-185-5

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sodium hydroxide solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1 — RS

Code	Size	Packaging	Notes
611081402	1 l	Plastic bottle	Ref Ph.Eur 1081402
611081404	1 l	Plastic bottle	Sodium hydroxide solution, strong Ref Ph.Eur 1081404

Sodium hydroxide solution, methanolic

• Sodio idrossido soluzione, metanolica • Sodium hydroxyde solution dans le méthanol
• Sodio hidróxido solución, metanólico • Natronlauge in Methanol

Synonym:
Caustic soda

NaOH
Molecular Weight: 40
CAS: 1310-73-2

Classification transport
ONU: 2810

Sodium hydroxide solution, methanolic > RS - For analysis according to Ph. Eur. Chap. 4.1.1 — RS

Code	Size	Packaging	Notes
611081405	100 ml	Glass bottle	Sodium hydroxide solution, methanolic R1 Ref Ph.Eur 1081405

Sodium hydroxide-d 1 30%

• Sodio idrossido-d 30% • Sodium hydroxyde-d 30% • Sodio hidróxido-d 30% • Natronlauge-d 30%

Synonym:
Sodium deuteroxide solution

NaOD
Molecular Weight: 41
CAS: 14014-06-3
EEC-N: 237-825-2

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sodium hydroxide-d 1 30% > RS - For NMR - min 99.5% — RS

Code	Size	Packaging	Notes
P5675	25 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Sodium hydroxide-d 1 mol/l

• Sodio idrossido-d 1 mol/l • Sodium hydroxyde-d 1N • Sodio hidróxido-d 1 mol/l • Natronlauge-d 1N

Synonym:
Sodium deuteroxide solution

NaOD
Molecular Weight: 41
CAS: 14014-06-3

Classification transport
ONU: 1824
Transport Hazard class: 8
Packing group II

Sodium hydroxide-d 1 mol/l > RS - For NMR - min 99.8% — RS

Code	Size	Packaging	Notes
P5665	25 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Sodium hypochlorite solution 12.5%

• Sodio ipoclorito soluzione 12.5% • Sodium hypochlorite solution 12.5% • Sodio hipoclorito solución 12.5% • Natriumhypochloritlösung 12.5%

NaClO
Molecular Weight: 74,44
CAS: 7681-52-9

Classification transport
ONU: 1791
Transport Hazard class: 8
Packing group III



Danger
H290-H314-H400-H412-HEU031
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sodium hypochlorite solution 12.5% > RPE - For analysis — RPE

Appearance ..... Slight yellow liquid French chlorometric degree ..... ≥ 44 ° Active chlorine percentage ..... ≥ 11.6 % (m/m)

Code	Size	Packaging	Notes
P9350015	1 l	Plastic bottle	
P9350046	1 l	Glass bottle PVC coated	
P9350049	25 l	Plastic tank	

Store at ambient temperature

Sodium hypochlorite solution in water

• Sodio ipoclorito soluzione in acqua • Sodium hypochlorite solution aqueuse • Sodio hipoclorito solución en agua • Natriumhypochloritösung

NaClO
Molecular Weight: 74,44
CAS: 7681-52-9

Classification transport
ONU: 1791
Transport Hazard class: 8
Packing group III



Danger
H314-H412-HEU031
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sodium hypochlorite solution in water > RS - For analysis according to Ph. Eur. Chap. 4.1.1 — RS

Code	Size	Packaging	Notes
611081609	250 ml	Glass bottle	Sodium hypochlorite solution, strong Ref Ph.Eur 1081600
611081600	1 l	Glass bottle	Sodium hypochlorite solution, strong Ref Ph.Eur 1081600

Sodium hypochlorite solution in water > RS - For analysis according to JP — RS

Code	Size	Packaging	Notes
616001008	100 ml	Glass bottle	Sodium hypochlorite TS

Sodium hypochlorite solution in water > RS - For analysis according to USP — RS

Code	Size	Packaging	Notes
617000181	1 l	Glass bottle	

Sodium hypochlorite solution in water > RPE - For analysis — RPE

Description ..... Yellow clear liquid
Density at 20°C ..... ~ 1.1
Alcalinity (NaOH)..... ≤1.8 % m/m
Cd..... ≤10 ppm
Cr..... ≤10 ppm
Cu..... ≤10 ppm
Hg..... ≤10 ppm
Mn..... ≤10 ppm
Ni..... ≤10 ppm
Zn..... ≤10 ppm
Assay (iodometric)..... 5 ÷ 9 %(Cl)m/m

Code	Size	Packaging	Notes
481181	1 l	Plastic bottle	
481185	30 kg	Plastic drum	

Store at ambient temperature

Sodium hypochlorite solution in water > RE - Pure — RE

Description ..... Yellow clear liquid Alcalinity (NaOH)..... ≤ 1.8 % m/m Assay (iodometric) ..... 5 - 9 % Cl Density at 20°C ..... 1.1 - 1.2

Code	Size	Packaging	Notes
370321	1 l	Plastic bottle	
370323	5 l	Plastic bottle	
370322	30 kg	Plastic drum	

Store at ambient temperature

Product specifications are subject to changes. Please visit our website for updates.

Sodium hypophosphite

• Sodio ipofosfito • Sodium hypophosphite • Sodio hipofosfito • Natriumhypophosphit

$NaH_2PO_2.H_2O$
Molecular Weight: 106,06
CAS: 10039-56-2
EEC-N: 231-669-9

Sodium hypophosphite > RE - Pure

RE

Description .....Semitransparent crystals
Identification .....Positive
Chloride.....≤100 ppm
Heavy metals (Pb).....≤10 ppm
Sulphate.....≤500 ppm
As .....≤2 ppm
Fe .....≤50 ppm
Assay (oxidimetric) .....≥101 % t.q.

Code	Size	Packaging	Notes
481201	1 kg	Plastic bottle	
481202	5 kg	Plastic tank	

Sodium iodide

• Sodio ioduro • Sodium iodure • Sodio yoduro • Natriumiodid

NaI
Molecular Weight: 149,89
CAS: 7681-82-5
EEC-N: 231-679-3



Danger
H315-H319-H334-H317-H335
P261-P284-P304+P340-P305+P351+P338-P342+P311a-P403+P233

Sodium iodide > RPE - For analysis

RPE

Description .....White crystals
Identification .....Positive
Loss on drying ..... ≤ 2 %
Iodate.....≤ 4 ppm
Sulphate.....≤ 150 ppm
Fe .....≤ 20 ppm
Assay (oxidimetric) .....99 ÷ 101.5 %

Code	Size	Packaging	Notes
481163	50 g	Glass bottle	
481164	500 g	Plastic bottle	
481162	25 kg	Plastic bucket	

Sodium iodide > ERBApharm - According to pharmacopoeia: BP-FU-Ph.Eur.-Ph.Franc.

ERBApharm

Description ..... White crystalline powder
Identification .....Positive
Alcalinity.....Conform Ph.Eur.
Appearance of solution .....Conform Ph.Eur.
Thiosulphate .....Conform Ph.Eur.
Iodate.....Conform Ph.Eur.
Loss on drying ..... ≤ 3.0 %
Heavy metals (Pb).....≤ 10 ppm
Sulphate.....≤ 150 ppm
Fe .....≤ 20 ppm
Assay (oxidimetric) .....99.0 ÷ 100.5 % s.s.

Code	Size	Packaging	Notes
370305	250 g	Plastic bottle	
370307	1 kg	Plastic bottle	
370309	5 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium laurylsulfate

• Sodio laurilsolfato • Sodium laurylsulfate • Sodio laurilsulfato • Natriumdodecylsulfat

Synonym:
• *Sodium dodecyl sulfate*
• *Dodecyl sodium sulfate*

$CH_3(CH_2)_{11}OSO_3Na$
Molecular Weight: 288,38
CAS: 151-21-3
EEC-N: 273-257-1



Danger
H315-H318-H335-H412
P261-P304+P340-P310a-P305+P351+P338-P362+P364-P403+P233

Sodium laurylsulfate > RS - For surfactants detection

RS

Description .....White crystalline powder
Identification (I.R.)..... Positive
pH sol. 1% at 20°C ..... 8.5 ÷ 10.5
Free sulphate (Na2SO4).....≤ 2.5 %
Assay (acidimetric) .....≥ 92.0 %

Code	Size	Packaging	Notes
481231	250 g	Plastic bottle	
481233	10 kg	Plastic bucket	
481235	25 kg	Drum	

Sodium laurylsulfate > RPE - For analysis

RPE

Anionic surfactant ..... 93 - 98 %
Unsulphated matter ..... ≤ 1 %
Loss on drying ..... ≤ 1 %
Na2SO4 ..... ≤ 3 %
NaCl ..... ≤ 1.5 %
pH (1% solution) ..... 8.5 - 10.5

Code	Size	Packaging	Notes
P7600513	100 g	Plastic bottle	
P7600514	500 g	Plastic bottle	
P7600517	1 kg	Plastic bottle	

Sodium metabisulfite

• Sodio metabisolfito • Sodium métabisulfite • Sodio metabisulfito • Dinatriumdisulfit

Synonym:
• *Sodium bisulfite*
• *Sodium hydrogensulfite*

$Na_2O_5S_2$
Molecular Weight: 190,1
CAS: 7681-57-4
EEC-N: 231-673-0



Danger
H302-H318-HEU031
P264-P280i-P301+P312a-P305+P351+P338-P310a-P501a

Sodium metabisulfite > RPE - For analysis - ACS

RPE

Description ..... White crystalline powder
Identification ..... Positive
Chloride ..... ≤ 0.05 %
Water-insoluble matter ..... ≤ 50 ppm
Heavy metals (Pb) ..... ≤ 10 ppm
Thiosulphate ..... ≤ 0.05 %
Fe ..... ≤ 20 ppm
Assay (oxidimetric) ..... ≥ 97.0 %

Code	Size	Packaging	Notes
481286	100 g	Plastic bottle	
481287	1 kg	Plastic bottle	
481288	2.5 kg	Plastic bottle	
481283	25 kg	Plastic bucket	

Sodium metabisulfite > ERBApharm - According to pharmacopoeia: BP-NF-Ph.Eur.-FU

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution ..... Conform Ph.Eur.
pH solution 5% ..... 3.5 ÷ 5.0
Chloride ..... ≤500 ppm
Thiosulphate ..... ≤500 ppm
Heavy metals (Pb) ..... ≤20 ppm
As ..... ≤5 ppm
Fe ..... ≤20 ppm
Assay (oxidimetric) ..... 95.0 ÷ 100.5 %
Assay (SO2) ..... 65.0 ÷ 67.4 %
Residual solvents (Current ICH) ..... Conform
Origin (BSE/TSE) ..... Synthesis

Code	Size	Packaging	Notes
370751	1 kg	Plastic bottle	
370752	2.5 kg	Plastic bottle	
370753	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium metaperiodate

• Sodio (meta)-periodato • Sodium m-périodate • Sodio metaperyodato • Natriummetaperiodat

Synonym:
Sodium periodate

$NaIO_4$
Molecular Weight: 213,89
CAS: 7790-28-5
EEC-N: 232-197-6

Classification transport
ONU: 3085
Transport Hazard class: 5.1
Packing group I



Danger
H271-H314-H372-H400
P210-P280-P283-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sodium metaperiodate > RPE - For analysis - ACS

RPE

Description ..... White crystalline powder
Identification ..... Positive
Other halogens (Cl) ..... ≤ 0.02 %
Mn ..... ≤ 3 ppm
Assay (iodometric) ..... 99.8 ÷ 100.3 % (s.s.)

Code	Size	Packaging	Notes
482234	100 g	Glass bottle	
482236	1 kg	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Sodium metaphosphate

• Sodio metafosfato • Sodium métaphosphate • Sodio metafosfato • Natriummetaphosphat

$NaPO_3$
Molecular Weight: 102,2
CAS: 10361-03-2
EEC-N: 233-782-9

Sodium metaphosphate > RE - Pure

RE

Description ..... White crystalline powder
Identification ..... Positive
Chloride ..... ≤500 ppm
Heavy metals (Pb) ..... ≤20 ppm
Sulphate ..... ≤0.1 %
Fe ..... ≤500 ppm
Assay ..... ≥68.0 % P2O5

Code	Size	Packaging	Notes
481557	1 kg	Plastic bottle	
481552	25 kg	Plastic bucket	

Sodium methoxide 0.1 mol/l

• Sodio metossido 0.1 mol/l • Sodium méthanolate 0.1 mol/l • Sodio metóxido 0.1 mol/l
• Natriummethanolat 0.1 mol/l

Synonym:
Sodium methylate

CH_3NaO
Molecular Weight: 54,02
CAS: 124-41-4

Classification transport
ONU: 1992
Transport Hazard class: 3
Packing group II



Danger
H225-H301-H370
P210-P241-P280-P301+P310a-P303+P361+P353-P403+P235

Sodium methoxide 0.1 mol/l > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613007101	100 ml	Glass bottle	Ref Ph.Eur 3007100
613007100	1 l	Glass bottle	Ref Ph.Eur 3007100

Sodium molybdate dihydrate

• Sodio molibdato diidrato • Sodium molybdate dihydraté • Sodio molibdato dihidrato • Natriummolybdat

Synonym:
Molybdic acid sodium salt dihydrate

$Na_2MoO_4.2H_2O$
Molecular Weight: 241,95
CAS: 10102-40-6
EEC-N: 231-551-7



Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Sodium molybdate dihydrate > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Water insoluble substances ..... ≤ 0.05 %
Pb ..... ≤ 50 ppm
Assay (Mo) ..... ≥ 39.5 %

Code	Size	Packaging	Notes
481684	100 g	Glass bottle	
481685	250 g	Plastic bottle	
481687	1 kg	Plastic bottle	

Sodium nitrate

• Sodio nitrato • Sodium nitrate • Sodio nitrato • Natriumnitrat

$NaNO_3$
Molecular Weight: 84,99
CAS: 7631-99-4
EEC-N: 231-554-3

Classification transport
ONU: 1498
Transport Hazard class: 5.1
Packing group III



Danger
H272-H319
P210-P220-P264-P280-P305+P351+P338-P337+P313

Sodium nitrate > RPE - For analysis - ISO - ACS - Reag.Ph.Eur. - Reag.USP

RPE

Description	White crystals	Chloride	≤ 10 ppm	Sulphate	≤ 30 ppm	Ca	≤ 50 ppm
Identification	Positive	Phosphate	≤ 5 ppm	Heavy metals (Pb)	≤ 5 ppm	Mg	≤ 20 ppm
pH sol. 5% in H20	5.5 ÷ 8.3	Iodate	≤ 5 ppm	Fe	≤ 3 ppm		
Water-insoluble matter	≤ 50 ppm	Nitrite	≤ 10 ppm	Assay (acidimetric)	≥ 99.0 %		

Code	Size	Packaging	Notes
481755	100 g	Plastic bottle	
481756	500 g	Plastic bottle	
481757	1 kg	Plastic bottle	
481759	5 kg	Plastic jar	
481751	25 kg	Drum	

Sodium nitrate > RE - Pure

RE

Description	Yellowish crystals	Chloride	≤0.5 %	Heavy metals (Pb)	≤50 ppm	Fe	≤50 ppm
Identification	Positive	Water-insoluble matter	≤500 ppm	Sulphate	≤0.5 %	Assay (non-aqueos medium)	≥96 %

Code	Size	Packaging	Notes
371809	5 kg	Plastic tank	
371802	25 kg	Plastic bucket	
371804	50 kg	Fibre drum	

Sodium nitrite

• Sodio nitrito • Sodium nitrite • Sodio nitrito • Natriumnitrit

$NaNO_2$
Molecular Weight: 68,99
CAS: 7632-00-0
EEC-N: 231-555-9

Classification transport
ONU: 1500
Transport Hazard class: 5.1
Packing group III



Danger
H272-H301-H400
P210-P220-P264-P280-P301+P310a-P330

Sodium nitrite > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description	Yellow crystals	Chloride	≤50 ppm	Ca	≤100 ppm	Assay (oxidimetric)	≥97.0 %
Identification	Positive	Sulphate	≤100 ppm	Fe	≤10 ppm		
Water-insoluble matter	≤100 ppm	Heavy metals (Pb)	≤10 ppm	K	≤50 ppm		

Code	Size	Packaging	Notes
481825	100 g	Plastic bottle	
481826	500 g	Plastic bottle	
481827	1 kg	Plastic bottle	
481829	5 kg	Plastic jar	

Sodium nitrite > ERBAPharm - According to pharmacopoeia: Ph.Eur.-USP-BP

ERBApharm

Description	Yellow crystals	Heavy metals (Pb)	≤20 ppm	Appearance of solution	Conform Ph. Eur.	Sulfate	≤ 200 ppm
Identification	Positive	Assay (oxidimetric)	97.0 ÷ 101.0 % s.s.	Acidity or alcalinity	Conform Ph.Eur.		
Loss on drying	≤0.25 %	Assay (Ph. Eur.)	98.5 - 100.5 % (ds)	Chloride	≤ 50 ppm		

Code	Size	Packaging	Notes
371901	1 kg	Plastic bottle	
371902	5 kg	Plastic tank	
371903	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Product specifications are subject to changes. Please visit our website for updates.

Sodium nitrite > RE - Pure

RE

Description ..... Yellow crystals
Identification ..... Positive
Chloride ..... ≤0.1 %
Heavy metals (Pb) ..... ≤50 ppm
Sulphate ..... ≤0.1 %
Fe ..... ≤50 ppm
Assay (oxidimetric) ..... ≥95 %

Code	Size	Packaging	Notes
372109	5 kg	Plastic tank	

Sodium nitrite 0.1 mol/l (0.1N)

• Sodio nitrito 0.1 mol/l (0.1N) • Sodium nitrite 0.1 mol/l (0.1N) • Sodio nitrito 0.1 mol/l (0.1N) • Natriumnitrit 0.1 mol/l (0.1N)

$NaNO_2$
Molecular Weight: 68,99
CAS: 7632-00-0

Sodium nitrite 0.1 mol/l (0.1N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613007200	1 l	Plastic bottle	Ref Ph.Eur 3007200

Sodium nitrite solution 500 g/l

• Sodio nitrito 500 g/l soluzione • Sodium nitrite 500g/l • Sodio nitrito solución 500 g/L • Natriumnitrit 500 g/l

$NaNO_2$
Molecular Weight: 68,99
CAS: 7632-00-0

Classification transport
ONU: 3219
Transport Hazard class: 5.1
Packing group III



Danger
H272-H301-H400
P210-P220-P264-P280-P301+P310a-P330

Sodium nitrite solution 500 g/l > RPE - For analysis

RPE

Assay ..... 490 - 510 g/L

Code	Size	Packaging	Notes
524725	2 l	Plastic bottle	

Sodium nitroprusside dihydrate

• Sodio nitroprussiato diidrato • Sodium nitroprussiate dihydraté • Sodio nitroprusiato dihidrato • Natriumnitroprussid-Dihydrat

$Na_2Fe(CN)_5NO.2H_2O$
Molecular Weight: 297,95
CAS: 13755-38-9
EEC-N: 238-373-9

Classification transport
ONU: 3288
Transport Hazard class: 6.1
Packing group III



Danger
H301
P264-P270-P301+P310a-P330-P405-P501a

Sodium nitroprusside dihydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Red brick crystals
Identification ..... Positive
Chloride ..... ≤ 200 ppm
Sulphate ..... Conform
Not soluble matter ..... ≤ 0.01 %

Code	Size	Packaging	Notes
481932	50 g	Glass bottle	
481934	100 g	Glass bottle	

Sodium 1-octanesulfonate monohydrate ▶ 1-Octanesulfonic acid sodium salt monohydrate

Sodium oxalate

• Sodio ossalato • Sodium oxalate • Sodio oxalato • Natriumoxalat

Synonym:
• *Ethandioic acid sodium salt*
• *Oxalic acid disodium salt*

$(COONa)_2$
Molecular Weight: 134
CAS: 62-76-0
EEC-N: 200-550-3

Warning
H302-H312
P264-P270-P280h-P301+P312a-P330-P501a

Sodium oxalate > RS - Standard for volumetry

RS

Description ..... White crystalline powder Identification ..... Positive Assay ..... ≥99.8 %

Code	Size	Packaging	Notes
482101	50 g	Glass bottle	

Sodium oxalate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White crystalline powder
Identification ..... Positive
Neutrality ..... Conform
Ready carbonizable substances ..... Conform
Loss on drying ..... ≤100 ppm
Water-insoluble matter ..... ≤50 ppm
Ammonium ..... ≤20 ppm
Chloride ..... ≤20 ppm
Sulphate ..... ≤20 ppm
Heavy metals (Pb) ..... ≤20 ppm
Fe ..... ≤10 ppm
K ..... ≤50 ppm
Assay (oxidimetric) ..... ≥99.5 %

Code	Size	Packaging	Notes
482064	100 g	Plastic bottle	
482065	250 g	Plastic bottle	
482067	1 kg	Plastic bottle	

Sodium oxalate > RE - Pure

RE

Description ..... White crystalline powder
Identification ..... Positive
Chloride ..... ≤ 300 ppm
Sulphate ..... ≤ 100 ppm
Fe ..... ≤ 100 ppm
Assay (oxidimetric) ..... 94 ÷ 96 %

Code	Size	Packaging	Notes
372201	1 kg	Plastic bottle	
372203	5 kg	Plastic tank	

Sodium 1-pentanesulfonate ▶ 1-Pentanesulphonic acid sodium salt

Sodium 1-pentanesulfonate monohydrate ▶ 1-Pentanesulphonic acid sodium salt monohydrate

Sodium perborate tetrahydrate

• Sodio perborato tetraidrato • Sodium perborate tétrahydraté • Sodio perborato tetrahidratado • Natriumperborat-Tetrahydrat

$NaBO_3.4H_2O$
Molecular Weight: 153,88
CAS: 10486-00-7
EEC-N: 231-556-4

Danger
H302-H332-H318-H360Df-H335-HA26
P261-P280-P304+P340-P310a-P305+P351+P338-P403+P233

Sodium perborate tetrahydrate > RE - Pure

RE

Description ..... White crystalline powder
Identification ..... Positive
Chloride ..... ≤500 ppm
Heavy metals (Pb) ..... ≤100 ppm
Sulphate ..... ≤0.1 %
As ..... ≤20 ppm
Fe ..... ≤500 ppm
Assay (oxidimetric) ..... ≥96 %

Code	Size	Packaging	Notes
482183	1 kg	Plastic bottle	
482185	5 kg	Plastic tank	
482187	25 kg	Plastic bucket	

Product specifications are subject to changes. Please visit our website for updates.

Sodium perchlorate monohydrate

• Sodio perclorato monoidrato • Sodium perchlorate monohydraté • Sodio perclorato • Natriumperchlorat-Monohydrat

Synonym: *Hyperchloric acid sodium salt*

$NaClO_4.H_2O$
Molecular Weight: 140,46
CAS: 7791-07-3
EEC-N: 231-511-9

Classification transport
ONU: 1502
Transport Hazard class: 5.1
Packing group II



Danger
H272-H302
P210-P220-P264-P280-P301+P312a-P501a

Sodium perchlorate monohydrate > RPE - For analysis

RPE

Description ..... White crystal	Chloride ..... ≤ 50 ppm	Heavy metals (Pb) ..... ≤ 10 ppm	Fe ..... ≤ 50 ppm
Identification (I.R.) ..... Positive	Chlorates ..... ≤ 0.10 %	Sulfate ..... ≤ 50 ppm	Assay ..... ≥ 99.0 %

Code	Size	Packaging	Notes
482204	250 g	Glass bottle	

Sodium peroxide

• Sodio perossido • Sodium péroxyde • Sodio peróxido • Natriumperoxid

Synonym: *Sodium superoxide*

Na_2O_2
Molecular Weight: 77,98
CAS: 1313-60-6
EEC-N: 215-209-4

Classification transport
ONU: 1504
Transport Hazard class: 5.1
Packing group I



Danger
H271-H314
P210-P280-P283-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sodium peroxide > RPE - For analysis

RPE

Description ..... Pale yellow powder	Chloride ..... ≤40 ppm	Ca ..... ≤500 ppm	Pb ..... ≤20 ppm
Identification ..... Positive	Phosphate ..... ≤20 ppm	Fe ..... ≤20 ppm	Assay (oxidimetric) ..... ≥97 %
Total nitrogen ..... ≤20 ppm	Total sulphur ..... ≤200 ppm	K ..... ≤200 ppm	

Code	Size	Packaging	Notes
482252	1 kg	Metallic can	

Sodium persulfate

• Sodio persolfato • Sodium persulfate • Sodio persulfato • Natriumpersulfat

$Na_2S_2O_8$
Molecular Weight: 238,1
CAS: 7775-27-1
EEC-N: 231-892-1

Classification transport
ONU: 1505
Transport Hazard class: 5.1
Packing group III



Danger
H272-H302-H315-H334-H317-H335
P210-P280-P284-P304+P340-P342+P311a-P403+P233

Sodium persulfate > RPE - For analysis

RPE

Assay ..... ≥ 99 %	Iron (Fe) ..... ≤ 5 mg/Kg	Active oxygen ..... ≥ 6.65 %	Sulfuric acid ..... ≤ 0.1 %

Code	Size	Packaging	Notes
P1650517	1 kg	Plastic bottle	

Sodium persulfate > RE - Pure

RE

Description ..... White crystalline powder	Identification ..... Positive	Fe ..... ≤10 ppm	Assay (oxidimetric) ..... ≥98 %

Code	Size	Packaging	Notes
482365	250 g	Plastic bottle	
482367	2.5 kg	Plastic bottle	
482363	25 kg	Plastic bucket	

Sodium persulfate 1 mol/l

• Sodio persolfato 1 mol/l • Sodium persulfate 1 mol/l • Sodio persulfato 1 mol/l • Natriumpersulfat 1 mol/l

$Na_2S_2O_8$
Molecular Weight: 238.10
CAS: 7775-27-1

Classification transport
ONU: 3216
Transport Hazard class: 5.1
Packing group III



Danger
H272-H315-H334-H317
P210-P261-P280-P284-P304+P340-P342+P311a

Sodium persulfate 1 mol/l > RPE - For analysis

RPE

Refractive index at 20°C..... 1.35 - 1.354

Code	Size	Packaging	Notes
PS0083/42	20 l	Plastic tank	

Sodium phosphate dibasic anhydrous

• Sodio fosfato bibasico anidro • Sodium phosphate dibasique anhydre • Sodio fosfato dibásico anhidro • Natriumhydrogenphosphat wasserfrei

Synonym:
Disodium hydrogen phosphate

Na_2HPO_4
Molecular Weight: 141,96
CAS: 7558-79-4
EEC-N: 231-448-7

Sodium phosphate dibasic anhydrous > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White powder
Identification ..... Positive
pH sol. 5% at 25° C .....8.7 ÷ 9.3
Loss on drying ..... ≤ 0.2 %
Chloride.....≤ 20 ppm
Water-insoluble matter .....≤ 100 ppm
Heavy metals (Pb).....≤ 10 ppm
Sulphate.....≤ 50 ppm
Fe .....≤ 20 ppm
Assay (alkalimetric)..... ≥ 99.0 %

Code	Size	Packaging	Notes
480143	100 g	Plastic bottle	
480144	500 g	Plastic bottle	
480141	1 kg	Plastic bottle	
480142	5 kg	Plastic jar	

Sodium phosphate dibasic anhydrous > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP

ERBApharm

Description ..... White powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Reducing substances .....Conform Ph.Eur.
Fosfato monobasico .....Conform Ph.Eur.
Origin (BSE/TSE)..... Synthesis
Residual solvents (Current ICH).....Conform
Not soluble matter..... ≤ 0.4 %
Chloride.....≤ 200 ppm
Sulphate.....≤ 500 ppm
As .....≤ 2 ppm
Fe .....≤ 20 ppm
Heavy metals (Pb).....≤ 10 ppm
Loss on drying at 105°C..... ≤ 1.0 %
Assay (potentiometric) 98.0 ÷ 100.5 % s.s.
Loss on drying at 130°C..... ≤ 5.0 %

Code	Size	Packaging	Notes
369212	1 kg	Plastic bottle	
369213	5 kg	Plastic bucket	
369211	25 kg	Sack	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium phosphate dibasic anhydrous > ERBApharm - According to pharmacopoeia: USP

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Loss on drying ..... ≤ 5.0 %
Loss on drying 130° C..... ≤5.0 %
Water-insoluble matter ..... ≤0.4 %
Chloride..... ≤0.06 %
Not soluble matter..... ≤0.4 %
Sulphate..... ≤0.2 %
As .....≤16 ppm
Heavy metals (Pb).....≤20 ppm
Assay .....98.0 ÷ 100.5 % (s.s.)
Assay (potetiometric) 98.0 ÷ 100.5 % (s.s.)

Code	Size	Packaging	Notes
369275	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Product specifications are subject to changes. Please visit our website for updates.

Sodium phosphate dibasic anhydrous > RE - Pure

RE

Description ..... White powder
Identification ..... Positive
Loss on drying ..... ≤8 %
Heavy metals (Pb)..... ≤50 ppm
As ..... ≤10 ppm
Assay (acidimetric) ..... ≥92.0 %

Code	Size	Packaging	Notes
369257	1 kg	Plastic bottle	
369258	5 kg	Plastic tank	
369252	25 kg	Plastic bucket	

Sodium phosphate dibasic dihydrate

• Sodio fosfato bibasico biidrato • Sodium phosphate dibasique dihydraté • Sodio fosfato dibásico dihidrato • Natriumhydrogenphosphat Dihydrat

Synonym:
Disodium hydrogen phosphate dihydrate

$Na_2HPO_4.2H_2O$
Molecular Weight: 177,99
CAS: 10028-24-7
EEC-N: 231-448-7

Sodium phosphate dibasic dihydrate > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Loss on drying ..... 18.5 ÷ 21.5 %
Chloride..... ≤ 10 ppm
Fluoride ..... ≤ 3 ppm
Water-insoluble matter ..... ≤ 50 ppm
Heavy metals (Pb)..... ≤ 10 ppm
Sulphate..... ≤ 50 ppm
As ..... ≤ 1 ppm
Fe ..... ≤ 5 ppm
Pb ..... ≤ 1 ppm
Assay (potetiometric) ..... 98.0 ÷ 100.5 %
Loss on ignition..... 25.1 ÷ 25.5 %
pH solution 1% ..... 9.0 ÷ 9.3
Hg..... ≤ 1 ppm
Ammonium..... ≤ 10 ppm
Reducing substances ..... Conform
Sodium dihydrogen phosphate ..... ≤ 2.5 %

Code	Size	Packaging	Notes
480225	100 g	Plastic bottle	
480226	500 g	Plastic bottle	
480227	1 kg	Plastic bottle	
480222	5 kg	Plastic jar	

Sodium phosphate dibasic dihydrate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution ..... Conform Ph.Eur.
Reducing substances ..... Conform Ph.Eur.
Fosfato monosodico ..... ≤ 2.5 %
Water not sol. matter..... ≤ 0.4 %
Chloride..... ≤ 400 ppm
Sulphate..... ≤ 0.1 %
Heavy metals (Pb)..... ≤ 20 ppm
As ..... ≤ 4 ppm
Fe ..... ≤ 40 ppm
Loss on drying 130° C..... 19.5 ÷ 21.0 %
Titolo (alcalimetrico)..98.0 ÷ 100.5 % s.s.

Code	Size	Packaging	Notes
369185	5 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium phosphate dibasic dodecahydrate

• Sodio fosfato bibasico dodecaidrato • Sodium phosphate dibasique dodécahydraté
• Sodio fosfato dibasico dodecahidrato • Dibasisches Natriumphosphat-Dodecahydrat

Synonym:
Disodium hydrogen phosphate dodecahydrate

$Na_2HPO_4.12H_2O$
Molecular Weight: 358,14
CAS: 10039-32-4
EEC-N: 231-448-7

Sodium phosphate dibasic dodecahydrate > RPE - For analysis

RPE

Description ..White or gray crystalline mass	Heavy metals (Pb).....≤ 5 ppm	Fe .....≤ 5 ppm	Assay (potetiometric) .....98 ÷ 102 %
Identification .....Positive	Sulphate.....≤ 50 ppm	K.....≤ 100 ppm	Mono or Tribasic salt .....Conform
pH sol. 5% at 25° C .....9.0 ÷ 9.4	As .....≤ 0.5 ppm	Mg .....≤ 10 ppm	Co.....≤ 5 ppm
Total nitrogen.....≤ 10 ppm	Ca.....≤ 10 ppm	Ni.....≤ 5 ppm	Cr.....≤ 5 ppm
Chloride.....≤ 5 ppm	Cd.....≤ 5 ppm	Pb .....≤ 5 ppm	Mn .....≤ 5 ppm
Water-insoluble matter .....≤ 50 ppm	Cu.....≤ 5 ppm	Zn .....≤ 5 ppm	

Code	Size	Packaging	Notes
480133	100 g	Plastic bottle	
480136	500 g	Plastic bottle	
480137	1 kg	Plastic bottle	
480131	5 kg	Plastic tank	
480132	10 kg	Plastic bucket	
480135	25 kg	Plastic bucket	

Sodium phosphate dibasic dodecahydrate > ERBApharm - According to pharmacopoeia: BP-DAB-FU-Ph. Eur.-Ph.Franc.-USP

ERBApharm

Description .....White crystals	Monobasic phosphate.....Conform Ph.Eur.	As .....≤2 ppm	Loss on drying 130° C.....57.0 ÷ 61.0 %
Identification .....Positive	Chloride.....≤200 ppm	Fe .....≤20 ppm	
Appearance of solution .....Conform Ph.Eur.	Heavy metals (Pb).....≤10 ppm	Assay (alkalimetric).....98.5 ÷ 100.5 % s.s.	
Reducing substances .....Conform Ph.Eur.	Sulphate.....≤500 ppm	Not soluble matter..... ≤ 0.4 %	

Code	Size	Packaging	Notes
369158	1 kg	Plastic bottle	
369159	5 kg	Plastic tank	
369152	25 kg	Plastic bucket	
369154	50 Kg	Fibre drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium phosphate dibasic dodecahydrate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.

ERBApharm

Description .....White crystals	Monobasic phosphate.....Conform Ph.Eur.	As .....≤2 ppm	Loss on drying 130° C.....57.0 ÷ 61.0 %
Identification .....Positive	Chloride.....≤200 ppm	Fe .....≤20 ppm	
Appearance of solution .....Conform Ph.Eur.	Heavy metals (Pb).....≤10 ppm	Assay (alkalimetric).....98.5 ÷ 100.5 % s.s.	
Reducing substances .....Conform Ph.Eur.	Sulphate.....≤500 ppm	Not soluble matter..... ≤ 0.4 %	

Code	Size	Packaging	Notes
529510	1 kg	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Sodium phosphate monobasic dihydrate

• Sodio fosfato monobasico biidrato • Sodium phosphate monobasique dihydraté
• Sodio fosfato monobásico dihidrato • Natriumphosphat-einbasiges Dihydrat

Synonym:
Sodium dihydrogen phosphate dihydrate

$NaH_2PO_4.2H_2O$
Molecular Weight: 156,01
CAS: 13472-35-0
EEC-N: 231-449-1

Sodium phosphate monobasic dihydrate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Reducing substances .....Conform Ph.Eur.
Al,Ca and related sub. ..... Conform USP-NF
pH solution 5% .....4.2 ÷ 4.5
Water (K.F.) .....18.0 ÷ 26.5 %
Loss on drying 130° C.....21.5 ÷ 24.0 %
Chloride.....≤140 ppm
Heavy metals (Pb).....≤10 ppm
Sulphate.....≤300 ppm
Not soluble matter..... ≤0.2 %
As .....≤2 ppm
Fe .....≤10 ppm
Assay (alkalimetric).....98.0 ÷ 100.5 % s.s.

Code	Size	Packaging	Notes
369138	1 kg	Plastic bottle	
369139	5 kg	Plastic tank	
369132	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium phosphate monobasic monohydrate

• Sodio fosfato monobasico monoidrato • Sodium phosphate monobasique monohydraté
• Sodio fosfato monobásico monohidrato • Natriumphosphat monobasisch Monohydrat

Synonym:
• *Monosodium phosphate*
• *Sodium dihydrogen phosphate monohydrate*

$NaH_2PO_4.H_2O$
Molecular Weight: 137,99
CAS: 10049-21-5
EEC-N: 231-449-2

Sodium phosphate monobasic monohydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White crystalline powder
Identification ..... Positive
pH sol. 5% in H2O.....4.1 ÷ 4.5
Water-insoluble matter .....≤ 100 ppm
Chloride.....≤ 5 ppm
Heavy metals (Pb).....≤ 10 ppm
Sulphate.....≤ 30 ppm
Ca.....≤ 50 ppm
Fe .....≤ 10 ppm
K.....≤ 100 ppm
Assay .....98.0 ÷ 102.0 %

Code	Size	Packaging	Notes
480085	100 g	Plastic bottle	
480086	500 g	Plastic bottle	
480087	1 kg	Plastic bottle	
480082	5 kg	Plastic jar	
480081	25 kg	Plastic bucket	

Sodium phosphate monobasic monohydrate > ERBApharm - According to pharmacopoeia: USP

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
Al,Ca and related sub. ..... Conform USP-NF
pH (1:20) .....4.1 ÷ 4.5
Water (K.F.) .....10.0 ÷ 15.0 %
Chloride.....≤140 ppm
Sulphate..... ≤0.15 %
Not soluble matter..... ≤0.2 %
As .....≤8 ppm
Assay (alkalimetric).....98.0 ÷ 103.0 % s.s.
Origin (BSE/TSE)..... Synthesis
Residual solvents (Current ICH).....Conform
Heavy metals (Pb).....≤ 20 ppm

Code	Size	Packaging	Notes
369143	1 kg	Plastic bottle	
369141	5 kg	Plastic tank	
369142	25 kg	Sack	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium phosphate tribasic dodecahydrate

• Sodio fosfato tribasico dodecaidrato • Sodium phosphate tribasique dodécahydraté
• Sodio fosfato tribásico dodecahidratado • Natriumphosphat tribasisch Dodecahydrat

Synonym:
• *Trisodium phosphate dodecahydrate*
• *TSP*

$Na_3PO_4.12H_2O$
Molecular Weight: 380,12
CAS: 10101-89-0
EEC-N: 231-509-8

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Sodium phosphate tribasic dodecahydrate > RS - Nuclear

RS

Code	Size	Packaging	Notes
526001	25 kg	Drum	

For specifications, contact our customer service for a certificate of analysis

Sodium phosphate tribasic dodecahydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White crystals
Identification ..... Positive
Free alkalies (NaOH) ..... ≤ 2.5 %
Water-insoluble matter ..... ≤ 100 ppm
Chloride ..... ≤ 10 ppm
Sulphate ..... ≤ 100 ppm
Heavy metals (Pb) ..... ≤ 10 ppm
Fe ..... ≤ 10 ppm
Assay (alkalimetric) ..... 98.0 ÷ 102.0 %

Code	Size	Packaging	Notes
480275	100 g	Plastic bottle	
480276	500 g	Plastic bottle	
480277	1 kg	Plastic bottle	
480272	5 kg	Plastic tank	
480271	25 kg	Fibre drum	

Sodium phosphate tribasic dodecahydrate > RE - Pure

RE

Description ..... White powder
Identification ..... Positive
pH sol. 1% ..... 11.8 ÷ 12.5
Water-insoluble matter ..... ≤ 0.2 %
Assay (alkalimetric) ..... ≥ 95.0 %

Code	Size	Packaging	Notes
369309	5 kg	Plastic tank	
369301	25 kg	Plastic bucket	

Sodium o-Phosphite pentahydrate

• Sodio o-fosfito pentaidrato • Sodium o-Phosphite pentahydraté • Sodio o-Fosfito pentahidrato
• Natrium-o-phosphitpentahydrat

Synonym:
• *Sodium phosphite dibasic pentahydrate*
• *di-Sodium hydrogen phosphite*

$Na_2HPO_3.5H_2O$
Molecular Weight: 216
CAS: 13517-23-2

Sodium o-Phosphite pentahydrate > RPE - For analysis

RPE

Description .Colorless transparent granule crystal, weak odor
Identification ..... Positive
Water ..... 38.0 ÷ 45.0 %
Heavy metals (Pb) ..... ≤25 ppm
Assay (oxidimetric) ..... 97.0 ÷ 102.0 %

Code	Size	Packaging	Notes
482042	25 kg	Plastic bucket	
482041	50 kg	Fibre drum	

Product specifications are subject to changes. Please visit our website for updates.

Sodium pyrophosphate decahydrate

• Sodio pirofosfato decaidrato • Sodium pyrophosphate décahydraté • Sodio pirofosfato decahidrato • Natriumpyrophosphatdecahydrat

Synonym:
Sodium pyrophosphate tetrabasic decahydrate

$Na_4P_2O_7.10H_2O$
Molecular Weight: 446,06
CAS: 13472-36-1
EEC-N: 231-767-1

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Sodium pyrophosphate decahydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description .....White crystals
Identification .....Positive
pH sol. 5% at 25° C .....9.5 ÷ 10.5
Total nitrogen .....≤10 ppm
Chloride .....≤20 ppm
Water-insoluble matter .....≤100 ppm
Heavy metals (Pb) .....≤10 ppm
Sulphate .....≤50 ppm
Fe .....≤10 ppm
Assay (acidimetric) .....99.0 ÷ 103.0 %

Code	Size	Packaging	Notes
482426	100 g	Plastic bottle	
482427	1 kg	Plastic bottle	
482422	10 kg	Box	
482421	25 kg	Drum	

Sodium salicylate

• Sodio salicilato • Sodium salicylate • Sodio salicilato • Natriumsalicylat

Synonym:
• *Salicylic acid sodium salt*
• *Sodium-2-hydroxybenzoate*

HOC_6H_4COONa
Molecular Weight: 160,11
CAS: 54-21-7
EEC-N: 200-198-0

Warning
H302
P264-P270-P301+P312a-P330-P501a

Sodium salicylate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP-DAB

ERBApharm

Description ..... White crystalline powder
Identification .....Positive
Appearance of solution .....Conform Ph.Eur.
Acidity .....Conform Ph.Eur.
Sulphite+thiosul. (SO2) .... Conform USP-NF
Organic volatile impurities Conform USP-NF
Loss on drying ..... ≤ 0.5 %
Chloride..... ≤ 0.02 %
Heavy metals (Pb)..... ≤ 20 ppm
Sulphate..... ≤ 0.06 %
Assay (non-aqueos medium) .99.0 ÷ 100.5 % s.s.
Water (K.F.) ..... ≤ 0.5 %

Code	Size	Packaging	Notes
373607	1 kg	Plastic bottle	
373608	5 kg	Plastic tank	
373603	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium silicate

• Sodio silicato • Sodium silicate • Sodio silicato • Natriumsilicat

Synonym:
Water glass

$Na_2O.xSiO_2$
Molecular Weight: 182,13 (an,)
CAS: 1344-09-8
EEC-N: 215-687-4

Danger
H315-H318-H335
P261-P304+P340-P310a-P305+P351+P338-P362+P364-P403+P233

Sodium silicate > RE - Pure

RE

Description ..... White powder
Identification ..... Positive
Loss on drying ..... ~ 16 %
Assay(Na20 2Si02 anydr)..... ≥72 %

Code	Size	Packaging	Notes
373908	2.5 kg	Plastic bucket	
373909	5 kg	Plastic bucket	
373902	25 kg	Fibre drum	

Sodium succinate hexahydrate

• Sodio succinato esaidrato • Sodium succinate hexahydraté • Sodio succinato hexahidratado
• Natriumsuccinat-Hexahydrat

Synonym:
• *Butanedioic acid disodium salt*
• *Succinic acid disodium salt*

$(CH_2COONa)_2.6H_2O$
Molecular Weight: 270,15
CAS: 6106-21-4
EEC-N: 205-778-7

Sodium succinate hexahydrate > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
pH sol. 5% at 25° C ..... 8.4 ÷ 9.2
Ammonium ..... ≤10 ppm
Chloride ..... ≤10 ppm
Phosphate ..... ≤20 ppm
Water-insoluble matter ..... ≤50 ppm
Heavy metals (Pb) ..... ≤5 ppm
Sulphate ..... ≤50 ppm
As ..... ≤1 ppm
Ca ..... ≤50 ppm
Cu ..... ≤2 ppm
Fe ..... ≤5 ppm
K ..... ≤50 ppm
Ni ..... ≤2 ppm
Pb ..... ≤2 ppm
Zn ..... ≤2 ppm
Assay (non-aqueos medium) ..... ≥99 %

Code	Size	Packaging	Notes
483554	100 g	Plastic bottle	
483555	250 g	Plastic bottle	
483557	2.5 kg	Plastic bottle	
483551	25 kg	Fibre drum	

Sodium sulfate anhydrous

• Sodio solfato anidro • Sodium sulfate anhydre • Sodio sulfato anhidro • Natriumsulfat wasserfrei

Na_2SO_4
Molecular Weight: 142,04
CAS: 7757-82-6
EEC-N: 231-820-9

Sodium sulfate anhydrous > RS - For anhydrification

RS

Description ..... Crystals or white crystalline powder
Identification ..... Positive
Loss on drying 130°C ..... ≤ 1 %
pH (10% solution) ..... 6 - 10
NaCl ..... ≤ 0.25 %
Fe ..... ≤ 5 ppm
Assay (acidimetric) ..... ≥ 99 %

Code	Size	Packaging	Notes
P1320017	1 kg	Plastic bottle	
P1320027	5 kg	Plastic tank	
P1320044	25 kg	Plastic bucket	

Sodium sulfate anhydrous > RS - For residual pesticides analysis

RS

Description ..... Crystals or white crystalline powder
Identification ..... Positive
pH sol. 5% in H2O ..... 5.2 ÷ 9.2
Loss on ignition ..... ≤0.5 %
Total nitrogen ..... ≤5 ppm
Calcium + Magnesium ..... ≤ 150 ppm
Chloride ..... ≤10 ppm
Water-insoluble matter ..... ≤50 ppm
Heavy metals (Pb) ..... ≤5 ppm
As ..... ≤1 ppm
Fe ..... ≤10 ppm
Assay (acidimetric) ..... ≥99.0 %

Code	Size	Packaging	Notes
483025	500 g	Glass bottle	

Sodium sulfate anhydrous > RPE - Crystalline powder - For analysis - ACS - ISO

RPE

Description ..... Crystals or white crystalline powder
Identification ..... Positive
pH sol. 5% in H2O ..... 5.2 ÷ 9.2
Loss on ignition ..... ≤0.5 %
Water-insoluble matter ..... ≤100 ppm
Total nitrogen ..... ≤5 ppm
Chloride ..... ≤10 ppm
Phosphate ..... ≤10 ppm
Heavy metals (Pb) ..... ≤5 ppm
Ca ..... ≤100 ppm
Fe ..... ≤10 ppm
K ..... ≤100 ppm
Mg ..... ≤50 ppm
Assay (acidimetric) ..... ≥99.0 %

Code	Size	Packaging	Notes
483006	100 g	Plastic bottle	
483005	500 g	Plastic bottle	
483007	1 kg	Plastic bottle	
483009	5 kg	Plastic jar	
483001	25 kg	Plastic bucket	

DASIT GROUP | CARLO ERBA REAGENTS
Product specifications are subject to changes. Please visit our website for updates.

Sodium sulfate anhydrous > ERBApharm - According to pharmacopoeia: Ph.Eur.

ERBApharm

Description ..... White powder	Chloride ..... ≤ 450 ppm	Mg ..... ≤ 200 ppm	Residual solvents (Current ICH) ..... Conform
Identification ..... Positive	Heavy metals (Pb) ..... ≤ 45 ppm	Loss on drying 130°C ..... ≤ 0.5 %	
Appearance of solution ..... Conform Ph. Eur.	Ca ..... ≤450 ppm	Assay ..... 98.5 ÷ 101.0 %	
Acidity or alcalinity ..... Conform Ph.Eur.	Fe ..... ≤ 90 ppm	Origin (BSE/TSE) ..... Synthesis	

Code	Size	Packaging	Notes
375713	1 kg	Plastic bottle	
375716	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium sulfate decahydrate

• Sodio solfato decaidrato • Sodium sulfate décahydraté • Sodio sulfato decahidrato • Natriumsulfat-Decahydrat

Synonym: *Glauber's salt*

$Na_2SO_4.10H_2O$
Molecular Weight: 322,19
CAS: 7727-73-3
EEC-N: 231-820-9

Sodium sulfate decahydrate > RPE - For analysis - ACS

RPE

Description ..... White crystals	Chloride ..... ≤5 ppm	Ca ..... ≤50 ppm	Assay (acidimetric) ..... ≥99 %
Identification ..... Positive	Phosphate ..... ≤5 ppm	Fe ..... ≤5 ppm	
pH solution 5% ..... 5.2 ÷ 9.2	Water-insoluble matter ..... ≤100 ppm	K ..... ≤50 ppm	
Total nitrogen ..... ≤3 ppm	Heavy metals (Pb) ..... ≤3 ppm	Mg ..... ≤30 ppm	

Code	Size	Packaging	Notes
482957	1 kg	Plastic bottle	
482959	5 kg	Plastic tank	

Sodium sulfide nonahydrate

• Sodio solfuro nonaidrato • Sodium sulfure nonahydraté • Sodio sulfuro nonahidrato • Natriumsulfid-Nonahydrat

$Na_2S.9H_2O$
Molecular Weight: 240,18
CAS: 1313-84-4
EEC-N: 215-211-5

Classification transport
ONU: 1849
Transport Hazard class: 8
Packing group II



Danger
H302-H311-H314-H400
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P361+P364

Sodium sulfide nonahydrate > RPE - For analysis - ACS

RPE

Description ..... Crystals or Chunks	Ammonium ..... ≤ 50 ppm	Sulphite and thiosulphate (SO4) ..... ≤ 0.1 %
Identification ..... Positive	Fe ..... Conform	Assay (oxidimetric) ..... ≥ 98.0 %

Code	Size	Packaging	Notes
483484	100 g	Plastic bottle	
483485	250 g	Plastic bottle	
483487	1 kg	Plastic bottle	
483489	5 kg	Plastic tank	

Sodium sulfide nonahydrate solution

• Sodio solfuro nonaidrato soluzione • Sodium sulfure nonahydraté solution • Sodio sulfuro nonahidrato solución • Natriumsulfid-Nonahydrat-Lösung

$Na_2S.9H_2O$
Molecular Weight: 240,18
CAS: 1313-84-4

Classification transport
ONU: 2922
Transport Hazard class: 8
Packing group III



Danger
H302-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sodium sulfide nonahydrate solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611083901	100 ml	Plastic bottle	Ref Ph.Eur 1083901
611083902	100 ml	Plastic bottle	Sodium sulfide solution R1 Ref Ph.Eur 1083902

Sodium sulfide trihydrate

• Sodio solfuro triidrato • Sodium sulfure trihydraté • Sodio sulfuro trihidrato • Natriumsulfid-Trihydrat

$Na_2S.3H_2O$
Molecular Weight: 132
CAS: 42607-30-7
EEC-N: 215-211-5

Classification transport
ONU: 1849
Transport Hazard class: 8
Packing group II



Danger
H301-H314-HEU031
P280-P301+P310a-P301+P330+P331-
P303+P361+P353-P304+P340-P305+P351+P338

Sodium sulfide trihydrate > RE - Pure

RE

Description ..... Scaglie gialle
Identification ..... Positive
Assay (oxidimetric) ..... ≥97 %

Code	Size	Packaging	Notes
376403	25 kg	Fibre drum	

Sodium sulfite anhydrous

• Sodio solfito anidro • Sodium sulfite anhydre • Sodio sulfito anhidro • Natriumsulfit wasserfrei

Na_2SO_3
Molecular Weight: 126,04
CAS: 7757-83-7
EEC-N: 231-821-4

HEU031

Sodium sulfite anhydrous > RPE - For analysis - ACS

RPE

Description ..... White powder
Identification ..... Positive
Acidity ..... Conform
Free alkalinity ..... ≤0.03 meq/g
Chloride ..... ≤200 ppm
Water-insoluble matter ..... ≤50 ppm
Heavy metals (Pb) ..... ≤10 ppm
Fe ..... ≤10 ppm
Assay (oxidimetric) ..... ≥98.0 %

Code	Size	Packaging	Notes
483256	100 g	Plastic bottle	
483257	1 kg	Plastic bottle	
483258	2.5 kg	Plastic bottle	
483252	25 kg	Plastic bucket	

Sodium sulfite anhydrous > ERBApharm - According to pharmacopoeia: BP-Ph.Eur.

ERBApharm

Description ..... White powder
Identification ..... Positive
Appearance of solution ..... Conform Ph.Eur.
Heavy metals (Pb) ..... ≤10 ppm
Thiosulphate ..... ≤0.1 %
Fe ..... ≤10 ppm
Se ..... ≤10 ppm
Zn ..... ≤25 ppm
Assay (oxidimetric) ..... 95.0 ÷ 100.5 %

Code	Size	Packaging	Notes
376006	1 kg	Plastic bottle	
376008	2.5 kg	Plastic bottle	
376009	5 kg	Plastic tank	
376002	10 kg	Plastic tank	
376003	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium sulfocyanate

• Sodio solfocianuro • Sodium sulfocyanure • Sodio sulfocianuro • Natriumsulfocyanid

Synonym:
• *Sodium thiocyanate*
• *Sodium rhodanate*

NaSCN
Molecular Weight: 81,07
CAS: 540-72-7
EEC-N: 208-754-4

Warning
H302-HEU032
P264-P270-P301+P312a-P330-P501a

Sodium sulfocyanate > RPE - For analysis - ACS

RPE

Description .....White crystals
Identification .....Positive
Ammonium .....≤20 ppm
Carbonate .....≤0.2 %
Chloride .....≤100 ppm
Water-insoluble matter .....≤50 ppm
Heavy metals (Pb) .....≤5 ppm
Sulphate .....≤100 ppm
Sulphide .....≤10 ppm
Fe .....≤2 ppm
Assay (argentimetric) .....≥98.0 %

Code	Size	Packaging	Notes
483354	100 g	Plastic bottle	
483356	500 g	Plastic bottle	

Sodium tartrate dihydrate

• Sodio tartrato diidrato • Sodium tartrate dihydraté • Sodio tartrato dihidrato • Natriumtartratdihydrat

Synonym:
• *L-(+)tartaric acid disodium salt*
• *Disodium tartrate dihydrate*

$(CHOHCOONa)_2.2H_2O$
Molecular Weight: 230,08
CAS: 6106-24-7
EEC-N: 212-773-3

Sodium tartrate dihydrate > RS - For Karl Fischer's reagent standardization - ACS

RS

Description ..... White crystalline powder
Identification .....Positive
pH sol. 5% at 25° C .....7.0 ÷ 9.0
Loss on drying 150° C....15.61 ÷ 15.71 %
Ammonium .....≤30 ppm
Chloride .....≤5 ppm
Phosphate .....≤5 ppm
Water-insoluble matter .....≤50 ppm
Heavy metals (Pb) .....≤5 ppm
Sulphate .....≤50 ppm
Ca .....≤100 ppm
Fe .....≤10 ppm
Assay (non-aqueos medium) .99.0 ÷ 101.0 %

Code	Size	Packaging	Notes
483561	100 g	Plastic bottle	

Sodium tartrate dihydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White crystalline powder
Identification .....Positive
pH sol. 5% at 25° C .....7.0 ÷ 9.0
Loss on drying at 150°C.15.61 ÷ 15.71 %
Ammonium .....≤ 30 ppm
Chloride .....≤ 5 ppm
Phosphate .....≤ 5 ppm
Water-insoluble matter .....≤ 50 ppm
Heavy metals (Pb) .....≤ 5 ppm
Sulphate .....≤ 50 ppm
Ca .....≤ 100 ppm
Fe .....≤ 10 ppm
Assay (non-aqueos medium) .99.0 ÷ 101.0 %

Code	Size	Packaging	Notes
483635	100 g	Plastic bottle	
483636	500 g	Plastic bottle	
483637	1 kg	Plastic bottle	

Sodium tetraborate anhydrous

• Sodio tetraborato anidro • Sodium tétraborate anhydre • Sodio tetraborato anhidro
• Natriumtetraborat wasserfrei

Synonym:
Borax

$Na_2B_4O_7$
Molecular Weight: 201,22
CAS: 1330-43-4
EEC-N: 215-540-4

Danger
H360FD-HA26
P201-P202-P280-P308+P313-P405-P501a

Sodium tetraborate anhydrous > RPE - For analysis

RPE

Description .....White crystals
Identification .....Positive
Chloride .....≤ 500 ppm
Heavy metals (Pb) .....≤ 50 ppm
Ca .....≤ 500 ppm
Fe .....≤ 100 ppm
Assay (acidimetric) .....≥ 98.0 %

Code	Size	Packaging	Notes
483735	100 g	Plastic bottle	
483736	1 kg	Plastic bottle	

Sodium tetraborate decahydrate

• Sodio tetraborato decaidrato • Sodium tétraborate décahydraté • Sodio tetraborato decahidrato • Natriumtetraboratdecahydrat

Synonym:
Borax decahydrate

$Na_2B_4O_7.10H_2O$
Molecular Weight: 381,37
CAS: 1303-96-4
EEC-N: 215-540-4

Danger
H319-H360FD-HA26
P264-P280-P305+P351+P338-P308+P313-P337+P313-P501a

Sodium tetraborate decahydrate > RPE - For analysis - ACS - ISO

RPE

Description	White crystals	Phosphate	≤5 ppm	Ca	≤50 ppm	Ni	≤2 ppm
Identification	Positive	Water-insoluble matter	≤30 ppm	Cu	≤2 ppm	Pb	≤2 ppm
pH sol. M/100 at 25° C	9.00 ÷ 9.50	Heavy metals (Pb)	≤5 ppm	Fe	≤5 ppm	Zn	≤2 ppm
Carbonate	≤50 ppm	Sulphate	≤10 ppm	K	≤50 ppm	Assay (acidimetric)	≥99.5 %
Chloride	≤5 ppm	As	≤1 ppm	Mg	≤10 ppm		

Code	Size	Packaging	Notes
478815	100 g	Plastic bottle	
478816	500 g	Plastic bottle	
478817	1 kg	Plastic bottle	
478819	5 kg	Plastic tank	
478812	25 kg	Plastic bucket	

Sodium tetraborate decahydrate > ERBApharm - According to pharmacopoeia: BP-FU-NF-Ph.Eur.-Ph.Franc.

ERBApharm

Description	White crystalline powder	Carbonates and bicarb.	Conform USP-NF	Heavy metals (Pb)	≤ 20 ppm	Ca	≤ 100 ppm
Identification	Positive	pH sol. 4% at 25 °C	9.0 ÷ 9.6	Sulphate	≤ 50 ppm	Fe	≤ 4 ppm
Appearance of solution	Conform Ph.Eur	Ammonium	≤ 10 ppm	As	≤ 5 ppm	Assay (alkalimetric)	99.0 ÷ 103.0 %

Code	Size	Packaging	Notes
367207	1 kg	Plastic bottle	
367209	5 kg	Plastic tank	
367201	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium tetraphenylborate

• Sodio tetrafenilborato • Sodium tétraphénylborate • Sodio tetrafenilborato • Natriumtetraphenylborat

Synonym:
Tetraphenylboron sodium salt

$Na[B(C_6H_5)_4]$
Molecular Weight: 342,23
CAS: 143-66-8
EEC-N: 205-605-5

Classification transport
ONU: 2811
Transport Hazard class: 6.1
Packing group III

Danger
H301
P264-P270-P301+P310a-P330-P405-P501a

Sodium tetraphenylborate > RPE - For analysis - ACS

RPE

Description	White crystalline powder	Appearance of solution	Conform	Assay (gravimetric)	≥99.5 %
Identification	Positive	Loss on drying	≤0.5 %		

Code	Size	Packaging	Notes
483758	5 g	Glass bottle	
483751	25 g	Glass bottle	

Reagent for the precipitation titration and potassium

Sodium thiocyanate ▶ Sodium sulfocyanate

Product specifications are subject to changes. Please visit our website for updates.

Sodium thiosulfate anhydrous

• Sodio tiosolfato anidro • Sodium thiosulfate anhydre • Sodio tiosulfato anhidro • Natriumthiosulfat wasserfrei

$Na_2S_2O_3$
Molecular Weight: 158,11
CAS: 7772-98-7
EEC-N: 231-867-5

Sodium thiosulfate anhydrous > RE - Pure

RE

Description ..... White crystalline powder
Identification ..... Positive
pH sol 10%.....6.5 ÷ 9.5
Sulphat + sulphit (SO4) ..... ≤ 1 %
Loss on drying ..... ≤ 2 %
Heavy metals (Pb).....≤ 50 ppm
Assay (oxidimetric) ..... ≥98 %
Fe .....≤ 50 ppm

Code	Size	Packaging	Notes
378377	1 kg	Plastic bottle	
378378	5 kg	Plastic tank	
378372	10 kg	Plastic tank	

Sodium thiosulfate pentahydrate

• Sodio tiosolfato pentaidrato • Sodium thiosulfate pentahydraté • Sodio tiosulfato pentahidrato • Natriumthiosulfat pentahydrat

$Na_2S_2O_3.5H_2O$
Molecular Weight: 248,18
CAS: 10102-17-7
EEC-N: 231-867-5

Sodium thiosulfate pentahydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description .....Semitransparent crystals
Identification ..... Positive
pH sol. 5% at 25° C .....6.0 ÷ 8.4
Total nitrogen .....≤20 ppm
Water-insoluble matter .....≤50 ppm
Sulphat + sulphit (SO4) ..... ≤0.1 %
Sulphide .....≤1 ppm
Assay (iodometric).....99.5 ÷ 101.0 %
Appearance of solution .....Conform
pH sol. 10% at 20°C .....6.0 ÷ 8.4
Heavy metals (Pb).....≤10 ppm

Code	Size	Packaging	Notes
483825	100 g	Plastic bottle	
483826	500 g	Plastic bottle	
483827	1 kg	Plastic bottle	
483829	5 kg	Plastic tank	
483821	25 kg	Plastic bucket	

Sodium thiosulfate pentahydrate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP

ERBApharm

Description .....Colourless crystals
Identification ..... Positive
Appearance of solution .....Conform Ph.Eur.
Sulphide .....Conform Ph.Eur.
Ca..... Conform USP-NF
pH sol 10%.....6.0 ÷ 8.4
Loss at 45°C .....32.0 ÷ 37.0 %
Sulphat + sulphit (SO4) ..... ≤0.2 %
Heavy metals (Pb).....≤10 ppm
Assay (iodometric).....99.0 ÷ 100.5 % s.s.

Code	Size	Packaging	Notes
377907	1 kg	Plastic bottle	
377909	5 kg	Plastic tank	
377901	25 kg	Drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sodium thiosulfate pentahydrate > RE - Pure

RE

Description ..... Colourless crystals
Identification ..... Positive
Heavy metals (Pb)..... ≤ 10 ppm
S..... ≤ 20 ppm
Fe ..... ≤ 5 ppm
Assay (oxidimetric) .....≥ 99 %

Code	Size	Packaging	Notes
378207	1 kg	Plastic bottle	
378209	5 kg	Plastic tank	
378202	25 kg	Plastic bucket	

Sodium thiosulfate 1 mol/l (1N)

• Sodio tiosolfato 1 mol/l (1N) • Sodium thiosulfate 1 mol/l (1N) • Sodio tiosulfato 1 mol/l (1N) • Natriumthiosulfat 1 mol/l (1N)

$Na_2S_2O_3.5H_2O$
Molecular Weight: 248,18
CAS: 10102-17-7

Sodium thiosulfate 1 mol/l (1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.99 - 1.01 N

Code	Size	Packaging	Notes
484026	500 ml	Plastic bottle	

158,11 g of Na2S2O3. Volumetric solution ready-to-use

Sodium thiosulfate 0.5 mol/l (0.5N)

• Sodio tiosolfato 0.5 mol/l (0.5N) • Sodium thiosulfate 0.5 mol/l (0.5N) • Sodio tiosulfato 0.5 mol/l (0.5N) • Natriumthiosulfat 0.5 mol/l (0.5N)

$Na_2S_2O_3.5H_2O$
Molecular Weight: 248,18
CAS: 10102-17-7

Sodium thiosulfate 0.5 mol/l (0.5N) > RPE - For analysis

RPE

Assay (potentiometry) ..... 0.499 - 0.501 N

Code	Size	Packaging	Notes
P3530015	1 l	Plastic bottle	

Sodium thiosulfate 0.2 mol/l (0.2N)

• Sodio tiosolfato 0.2 mol/l (0.2N) • Sodium thiosulfate 0.2 mol/l (0.2N) • Sodio tiosulfato 0.2 mol/l (0.2N) • Natriumthiosulfat 0.2 mol/l (0.2N)

$Na_2S_2O_3.5H_2O$
Molecular Weight: 248,18
CAS: 10102-17-7

Sodium thiosulfate 0.2 mol/l (0.2N) > RPE - For analysis

RPE

Assay (potentiometry) ..... 0.1998 - 0.2002 N

Code	Size	Packaging	Notes
P3520022	5 l	Plastic tank	

Sodium thiosulfate 0.1 mol/l (0.1N)

• Sodio tiosolfato 0.1 mol/l (0.1N) • Sodium thiosulfate 0.1 mol/l (0.1N) • Sodio tiosulfato 0.1 mol/l (0.1N) • Natriumthiosulfat 0.1 mol/l (0.1N)

$Na_2S_2O_3.5H_2O$
Molecular Weight: 248,18
CAS: 10102-17-7

Sodium thiosulfate 0.1 mol/l (0.1N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613007301	500 ml	Glass bottle	Ref Ph.Eur 3007300
613007300	1 l	Glass bottle	Ref Ph.Eur 3007300

Sodium thiosulfate 0.1 mol/l (0.1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.0999 - 0.1001 N NIST 136 ..... e

Code	Size	Packaging	Notes
484077000	1 l	Glass bottle	Certified with NIST traceability
484072000	5 l	Kubidos	Certified with NIST traceability
484071000	10 l	Kubidos	Certified with NIST traceability

15.811 g of Na2S2O3. Volumetric solution ready-to-use

Product specifications are subject to changes. Please visit our website for updates.

Sodium thiosulfate 0.1 mol/l (0.1N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
484121		Plastic ampoule	Volume: 55 ml

Volumetric concentrated solution to prepare 1 L of solution 0,1 M

Sodium thiosulfate 0.0394 mol/l (0.0394N)

• Sodio tiosolfato 0.0394 mol/l (0.0394N) • Sodium thiosulfate 0.0394 mol/l (0.0394N) • Sodio tiosulfato 0.0394 mol/l (0.0394N)
• Natriumthiosulfat 0.0394 mol/l (0.0394N)

$Na_2S_2O_3.5H_2O$
Molecular Weight: 248,18
CAS: 10102-17-7

Sodium thiosulfate 0.0394 mol/l (0.0394N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.0390 - 0.0398 N

Code	Size	Packaging	Notes
484141	2.5 l	Glass bottle	

Sodium thiosulfate 0.0197 mol/l (0.0197N)

• Sodio tiosolfato 0.0197 mol/l (0.0197N) • Sodium thiosulfate 0.0197 mol/l (0.0197N) • Sodio tiosulfato 0.0197 mol/l (0.0197N)
• Natriumthiosulfat 0.0197 mol/l (0.0197N)

$Na_2S_2O_3.5H_2O$
CAS: 10102-17-7

Sodium thiosulfate 0.0197 mol/l (0.0197N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
484155	2.5 l	Glass bottle	

Sodium thiosulfate 0.01 mol/l (0.01N)

• Sodio tiosolfato 0.01 mol/l (0.01N) • Sodium thiosulfate 0.01 mol/l (0.01N) • Sodio tiosulfato 0.01 mol/l (0.01N) • Natriumthiosulfat 0.01 mol/l (0.01N)

$Na_2S_2O_3.5H_2O$
Molecular Weight: 248,18
CAS: 10102-17-7

Sodium thiosulfate 0.01 mol/l (0.01N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
484161		Plastic ampoule	Volume: 55 ml

1,581 g Na2S2O3. Volumetric concentrated solution to prepare 1 L of solution 0,01 N

Sodium p-toluenesulfonate ▶ p-Toluenesulfonic acid sodium salt

Sodium tungstate dihydrate

• Sodio tungstato diidrato • Sodium tungstate dihydraté • Sodio tungstato dihidrato • Natriumwolframatdihydrat

Synonym:
Tungstic acid sodium salt dihydrate

$Na_2WO_4.2H_2O$
Molecular Weight: 329,86
CAS: 10213-10-2
EEC-N: 236-743-4

Warning
H302-H319
P264-P280i-P301+P312a-P305+P351+P338-P337+P313-P501a

Sodium tungstate dihydrate > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description	White crystals	Chloride	≤ 50 ppm	Sulphate	≤ 100 ppm
Identification	Positive	Water-insoluble matter	≤ 100 ppm	Mo	≤ 10 ppm
Free alkalinity	≤ 0.02 meq/g	Heavy metals and Fe(Pb)	≤ 10 ppm	Assay (gravimetric)	99.0 ÷ 101.0 %

Code	Size	Packaging	Notes
484233	50 g	Glass bottle	
484236	500 g	Plastic bottle	

Solvent Plus

• Solvente Plus • Solvant Plus • Disolvente Plus • Lösungsmittel Plus

Molecular Weight: 176
CAS: 68551-19-9
EEC-N: 271-369-5

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group III

Danger
H226-H304-HEU066
P210-P241-P280-P301+P310a-P303+P361+P353-P403+P235

Solvent Plus > RS - For histology

RS

Description	Clear liquid	Identification	Positive	Aromatic compounds	≤ 0.05 %	Distillation range	175 ÷ 213 °C

Code	Size	Packaging	Notes
446187	2.5 l	Glass bottle	
446181	5 l	Plastic tank	

Isoparaffins based histological clearing agent

Sorbitol

• Sorbitolo • Sorbitol • Sorbitol • Sorbit

Synonym:
D-Glucitol

$C_6H_{14}O_6$
Molecular Weight: 182,17
CAS: 50-70-4
EEC-N: 200-061-5

Sorbitol > ERBApharm - According to pharmacopoeia: FU-Ph.Eur.

ERBApharm

Description	White crystalline powder	Reducing sugar	Conform Ph.Eur.	Sostanze analoghe (HPLC)	Conform Ph.Eur.	Escherichia coli	Absent Ph. Eur.
Identification	Positive	Assay (HPLC)	97.0 ÷ 102.0 % anidro	Conductivity	≤ 20 µS.cm-1	Salmonella	Absent Ph. Eur.
Water (K.F.)	≤ 1.5 %	Appearance of solution	Conform Ph.Eur.	Microbial tests			
Ni	≤ 1 ppm	Specific optical rotation (anhydrous)	+4.0 ÷ +7.0 °	TAMC	≤ 1000 CFU/g		
Pb	≤ 0.5 ppm			TYMC	≤ 100 CFU/g		

Code	Size	Packaging	Notes
379013	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Product specifications are subject to changes. Please visit our website for updates.

Sorbitol (no crystallizable) solution 70%

• Sorbitolo (non cristallizzabile) soluzione 70% • Sorbitol solution 70%
• Sorbitol (no cristalizable) solución 70% • Sorbitollösung 70%

Synonym:
D-Glucitol

$C_6H_{14}O_6$
Molecular Weight: 182,17
CAS: 50-70-4
EEC-N: 200-061-5

Sorbitol (no crystallizable) solution 70% > ERBApharm - According to pharmacopoeia: Ph.Eur.-NF

ERBApharm

Description ..... Clear colourless liquid
Identification ..... Positive
Appearance of solution ..... Conform Ph.Eur.
Water (K.F.) ..... 28.5 ÷ 31.5 %
Conductivity ..... ≤ 10 µS.cm-1
Ni ..... ≤ 1 ppm
Reducing sugar ..... ≤ 0.2 %
Zuccheri riducenti dopo idrolisi ..... ≤ 9.3 %
Specific optical rotation ..... +1.5 - +3.5 °
Residue on ignition ..... ≤ 0.1 %
pH sol. 14% ..... 5.0 - 7.5
Assay (HPLC)
Sostanza anidra ..... 68.0 ÷ 72.0 %
D- Sorbitolo ..... 72.0 ÷ 92.0 %

Code	Size	Packaging	Notes
379021	1 l	Plastic bottle	
379022	5 l	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

D-Sorbitol

• D-Sorbitolo • D-Sorbitol • D-Sorbitol • Sorbit D

Synonym:
D-Glucitol

$C_6H_{14}O_6$
Molecular Weight: 182,17
CAS: 50-70-4
EEC-N: 200-061-5

D-Sorbitol > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Loss on drying ..... ≤1 %
Chloride ..... ≤50 ppm
Heavy metals (Pb) ..... ≤10 ppm
Residue on ignition ..... ≤0.1 %
Sulphate ..... ≤100 ppm
Red.ing sugars(Glucose) ..... ≤0.1 %
Totol sugars(Glucose) ..... ≤0.3 %
As ..... ≤2 ppm
Ca ..... ≤50 ppm
Fe ..... ≤10 ppm
Assay (oxidimetric) ..... ≥98 %

Code	Size	Packaging	Notes
484704	100 g	Plastic bottle	
484705	250 g	Plastic bottle	
484701	1 kg	Plastic bottle	

Standard solution 1.30 µS/cm

• Soluzione standard da 1.30 µS/cm • Etalon de conductivité 1.30 µS/cm • Solución patrón de 1.30 µS/cm • Leitfähigkeitsstandard 1.30 µS/cm

Classification transport
ONU: 1274
Transport Hazard class: 3
Packing group III



Danger
H226-H302-H318-H336
P210-P280-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Standard solution 1.30 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Conductivity at 25°C ..... 1.25 ÷ 1.35 µS/cm

Code	Size	Packaging	Notes
575231	250 ml	Glass bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 5 µS/cm

• Soluzione standard da 5 µS/cm • Etalon de conductivité 5 µS/cm • Solución patrón de 5 µS/cm • Leitfähigkeitsstandard 5 µS/cm

Classification transport
ONU: 1274
Transport Hazard class: LQ



Danger
H226-H302-H318-H336
P210-P280-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Standard solution 5 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....4.95 ÷ 5.05 µS/cm

Code	Size	Packaging	Notes
575001	250 ml	Glass bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 10 µS/cm

• Soluzione standard da 10 µS/cm • Etalon de conductivité 10µS/cm • Solución patrón de 10 µS/cm • Leitfähigkeitsstandard 10 µS/cm

Classification transport
ONU: 1274
Transport Hazard class: 3
Packing group III



Danger
H226-H302-H318-H336
P210-P280-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Standard solution 10 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....9.80 ÷ 10.20 µS/cm

Code	Size	Packaging	Notes
575011	250 ml	Glass bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 20 µS/cm

• Soluzione standard da 20 µS/cm • Etalon de conductivité 20 µS/cm • Solución patrón de 20 µS/cm • Leitfähigkeitsstandard 20 µS/cm

Classification transport
ONU: 1274
Transport Hazard class: LQ



Danger
H226-H302-H318-H336
P210-P280-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Standard solution 20 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....19.80 ÷ 20.20 µS/cm

Code	Size	Packaging	Notes
575021	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 50 µS/cm

• Soluzione standard da 50 µS/cm • Etalon de conductivité 50 µS/cm • Solución patrón de 50 µS/cm • Leitfähigkeitsstandard 50 µS/cm

Classification transport
ONU: 1274
Transport Hazard class: 3
Packing group III



Danger
H226-H318-H336
P210-P280-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Standard solution 50 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....49.0 - 51.0 µS/cm

Code	Size	Packaging	Notes
575031	500 ml	Glass bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Product specifications are subject to changes. Please visit our website for updates.

Standard solution 84 µS/cm
• Soluzione standard da 84 µS/cm • Etalon de conductivité 84 µS/cm • Solución patrón de 84 µS/cm • Leitfähigkeitsstandard 84 µS/cm

Standard solution 84 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....83.16 ÷ 84.84 µS/cm

Code	Size	Packaging	Notes
575041	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 100 µS/cm
• Soluzione standard da 100 µS/cm • Etalon de conductivité 100 µS/cm • Solución patrón de 100 µS/cm • Leitfähigkeitsstandard 100 µS/cm

Standard solution 100 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....99.0 ÷ 101.0 µS/cm

Code	Size	Packaging	Notes
575051	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 147 µS/cm
• Soluzione standard da 147 µS/cm • Etalon de conductivité 147 µS/cm • Solución patrón de 147 µS/cm • Leitfähigkeitsstandard 147 µS/cm

Standard solution 147 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....145.5 ÷ 148.5 µS/cm

Code	Size	Packaging	Notes
575061	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 200 µS/cm
• Soluzione standard da 200 µS/cm • Etalon de conductivité 200 µS/cm • Solución patrón de 200 µS/cm • Leitfähigkeitsstandard 200 µS/cm

Standard solution 200 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....198.0 ÷ 202.0 µS/cm

Code	Size	Packaging	Notes
575071	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 500 µS/cm
• Soluzione standard da 500 µS/cm • Etalon de conductivité 500 µS/cm • Solución patrón de 500 µS/cm • Leitfähigkeitsstandard 500 µS/cm

Standard solution 500 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....495.0 ÷ 505.0 µS/cm

Code	Size	Packaging	Notes
575081	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 1000 µS/cm

• Soluzione standard da 1000 µS/cm • Etalon de conductivité 1000 µS/cm • Solución patrón de 1000 µS/cm • Leitfähigkeitsstandard 1000 µS/cm

Standard solution 1000 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....990.0 ÷ 1010.0 µS/cm

Code	Size	Packaging	Notes
575091	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 1413 µS/cm

• Soluzione standard da 1413 µS/cm • Etalon de conductivité 1413 µS/cm • Solución patrón de 1413 µS/cm • Leitfähigkeitsstandard 1413 µS/cm

Standard solution 1413 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....1399 ÷ 1427 µS/cm

Code	Size	Packaging	Notes
575101	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 5000 µS/cm

• Soluzione standard da 5000 µS/cm • Etalon de conductivité 5000 µS/cm • Solución patrón de 5000 µS/cm • Leitfähigkeitsstandard 5000 µS/cm

Warning
H302
P264-P270-P301+P312a-P330-P501a

Standard solution 5000 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....4950 ÷ 5050 µS/cm

Code	Size	Packaging	Notes
575111	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 10000 µS/cm

• Soluzione standard da 10000 µS/cm • Etalon de conductivité 10 000 µS/cm • Solución patrón de 10000 µS/cm • Leitfähigkeitsstandard 10.000 µS/cm

Warning
H302
P264-P270-P301+P312a-P330-P501a

Standard solution 10000 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....9900 ÷ 10100 µS/cm

Code	Size	Packaging	Notes
575121	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 12880 µS/cm

• Soluzione standard da 12880 µS/cm • Etalon de conductivité 12 880 µS/cm • Solución patrón de 12880 µS/cm • Leitfähigkeitsstandard 12 880 µS/cm

Standard solution 12880 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....12751 ÷ 13009 µS/cm

Code	Size	Packaging	Notes
575131	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Product specifications are subject to changes. Please visit our website for updates.

Standard solution 20000 µS/cm

• Soluzione standard da 20000 µS/cm • Etalon de conductivité 20 000 µS/cm • Solución patrón de 20000 µS/cm • Leitfähigkeitsstandard 20.000 µS/cm

Warning
H302
P264-P270-P301+P312a-P330-P501a

Standard solution 20000 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....19800 ÷ 20200 µS/cm

Code	Size	Packaging	Notes
575141	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 50000 µS/cm

• Soluzione standard da 50000 µS/cm • Etalon de conductivité 50 000 µS/cm • Solución patrón de 50000 µS/cm • Leitfähigkeitsstandard 50.000 µS/cm

Warning
H302
P264-P270-P301+P312a-P330-P501a

Standard solution 50000 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....49500 ÷ 50500 µS/cm

Code	Size	Packaging	Notes
575151	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 100000 µS/cm

• Soluzione standard da 100000 µS/cm • Etalon de conductivité 100 000 µS/cm • Solución patrón de 100000 µS/cm • Leitfähigkeitsstandard 100.000 µS/cm

Warning
H302
P264-P270-P301+P312a-P330-P501a

Standard solution 100000 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....99000 ÷ 101000 µS/cm

Code	Size	Packaging	Notes
575161	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 150000 µS/cm

• Soluzione standard da 150000 µS/cm • Etalon de conductivité 150 000 µS/cm • Solución patrón de 150000 µS/cm • Leitfähigkeitsstandard 150.000 µS/cm

Warning
H302
P264-P270-P301+P312a-P330-P501a

Standard solution 150000 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....148500 ÷ 151500 µS/cm

Code	Size	Packaging	Notes
575171	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 200000 µS/cm

• Soluzione standard da 200000 µS/cm • Etalon de conductivité 200 000 µS/cm • Solución patrón de 200000 µS/cm • Leitfähigkeitsstandard 200.000 µS/cm



Warning
H302
P264-P270-P301+P312a-P330-P501a

Standard solution 200000 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....198000 ÷ 202000 µS/cm

Code	Size	Packaging	Notes
575181	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 300000 µS/cm

• Soluzione standard da 300000 µS/cm • Etalon de conductivité 300 000 µS/cm • Solución patrón de 300000 µS/cm • Leitfähigkeitsstandard 300 000 µS/cm



Warning
H302
P264-P270-P301+P312a-P330-P501a

Standard solution 300000 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....297000 ÷ 303000 µS/cm

Code	Size	Packaging	Notes
575191	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 350000 µS/cm

• Soluzione standard da 350000 µS/cm • Etalon de conductivité 350 000 µS/cm • Solución patrón de 350000 µS/cm • Leitfähigkeitsstandard 350.000 µS/cm

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group II



Warning
H290-H302-H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Standard solution 350000 µS/cm > RS - For conductivity

RS

Code	Size	Packaging	Notes
575201	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Standard solution 450000 µS/cm

• Soluzione standard da 450000 µS/cm • Etalon de conductivité 450 000 µS/cm • Solución patrón de 450000 µS/cm • Leitfähigkeitsstandard 450.000 µS/cm

Classification transport
ONU: 1789
Transport Hazard class: 8
Packing group III



Warning
H290
P234-P390-P406

Standard solution 450000 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....445500 ÷ 454500 µS/cm

Code	Size	Packaging	Notes
575211	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Product specifications are subject to changes. Please visit our website for updates.

Standard solution 500000 µS/cm

• Soluzione standard da 500000 µS/cm • Etalon de conductivité 500 000 µS/cm • Solución patrón de 500000 µS/cm • Leitfähigkeitsstandard 500.000 µS/cm

Classification transport
ONU: 2796
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Standard solution 500000 µS/cm > RS - For conductivity

RS

Description ..... Clear colourless liquid Identification ..... Positive Conductivity at 25°C .....495000 ÷ 505000 µS/cm

Code	Size	Packaging	Notes
575221	500 ml	Plastic bottle	

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Starch paste solution 1%

• Salda d'amido soluzione 1% • Empois d'amidon solution 1% • Almidón en pasta solución 1% • Stärkelösung 1%

$(C_6H_{10}O_5)n$
CAS: 9005-84-9

Starch paste solution 1% > RPE - For analysis

RPE

Description ..... Colourless opaline liquid Identification ..... Positive

Code	Size	Packaging	Notes
E477301	250 ml	Glass bottle	
E477302	1 l	Bottle	

Stabilized

Starch soluble

• Amido solubile • Amidon soluble • Almidón soluble • Stärke

$(C_6H_{10}O_5)n$
CAS: 9005-84-9
EEC-N: 232-686-4

Starch soluble > RPE - For analysis - Reag. Ph. Eur.

RPE

Description ..... White powder pH solution 2% ..... 5.0 ÷ 7.0 Sulphated ash .....≤ 1.5 %
Identification ..... Positive Loss on drying 100° C.....≤ 20 %

Code	Size	Packaging	Notes
417585	250 g	Plastic bottle	
417587	1 kg	Plastic bottle	

Starch soluble solution

• Amido solubile soluzione • Amidon soluble solution • Almidón soluble solución • Lösliche Stärkelösung

Starch soluble solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611085103	100 ml	Plastic bottle	Ref Ph.Eur 1085103
611085104	1 l	Bottle	Ref Ph.Eur 1085103

Stearic acid

• Acido stearico • Acide stéarique • Acido estearico • Stearinsäure

Synonym:
1-Heptadecanecarboxylic acid

$CH_3(CH_2)_{16}COOH$
Molecular Weight: 284,48
CAS: 57-11-4
EEC-N: 200-313-4

Classification transport
ONU: 1325
Transport Hazard class: 4.1
Packing group III

Warning
H228-H315-H319-H335
P210-P241-P280-P304+P340-P305+P351+P338-P403+P233

Stearic acid > ERBApharm - Vegetal origin - According to pharmacopoeia: Ph.Eur.-NF

ERBApharm

Description ..... White flakes
Identification ..... Positive
Appearance ..... Conform Ph.Eur.
Color of solution ..... Pass test USP-NF

Acidity ..... Conform Ph.Eur.
Acidity index ..... 194 ÷ 212
Iodine value ..... ≤ 4,0
Sulphated ash ..... ≤0,1 %

Heavy metals (Pb) ..... ≤10 ppm
Ni ..... ≤1 ppm
Melting point ..... 53 ÷ 59 ° C
Assay (Stearic Acid) ..... 40.0 ÷ 60.0 %

Assay (Stearic Acid+Palmitic Acid) ≥90.0 %

Code	Size	Packaging	Notes
307112	2.5 kg	Plastic bottle	
307115	25 kg	Fibre drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Stearic acid calcium salt ▶ Calcium stearate

Stearic acid magnesium salt ▶ Magnesium stearate

Strontium standard solution

• Stronzio standard soluzione • Strontium solution standard • Estroncio, solución patrón • Strontium-Standardlösung

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Strontium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505867	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505868	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505869	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Strontium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503951	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503953	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503955	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503957	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Product specifications are subject to changes. Please visit our website for updates.

Strontium standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
507760	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507493	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
E497665	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
E497661	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Strontium standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
485391		Plastic ampoule	conc. 1.000 ppm Matrix: Nitric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Strontium standard solution > RS - Standard solution for ion chromatography

RS

Code	Size	Packaging	Notes
503361	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Strontium acetate

• Stronzio acetato • Strontium acétate • Estroncio acetato • Strontium acetate

$Sr(CH_3COO)_2$
Molecular Weight: 205,71
CAS: 543-94-2
EEC-N: 208-854-8

Strontium acetate > RPE - For analysis

RPE

Description ..... White cryst. powder
Identification ..... Positive
pH sol. 5% at 25° C ..... 6.5 ÷ 8.5
Chloride..... ≤50 ppm
Insol.in dil.acetic ac. ..... ≤50 ppm
Heavy metals (Pb)..... ≤5 ppm
Nitrate ..... ≤30 ppm
Substanc. not ppt H_2SO_4..... ≤0.15 %
Sulphate ..... ≤30 ppm
Ba ..... ≤50 ppm
Ca ..... ≤0.1 %
Cu ..... ≤2.5 ppm
Fe ..... ≤5 ppm
K..... ≤500 ppm
Na..... ≤0.15 %
Ni..... ≤2.5 ppm
Pb ..... ≤2.5 ppm
Zn ..... ≤2.5 ppm
Assay (complexometric) ..... ≥99 %

Code	Size	Packaging	Notes
485304	100 g	Glass bottle	

Strontium bromide monohydrate

• Stronzio bromuro monoidrato • Strontium bromure monohydraté • Estroncio bromuro monohidrato • Strontiumbromidmonohydrat

$SrBr_2.H_2O$
Molecular Weight: 265,43
CAS: 14519-13-2
EEC-N: 233-969-5

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Strontium bromide monohydrate > RPE - For analysis

RPE

Description ..... Polvere granular bianca
Identification ..... Positive
pH sol. 5% at 25° C ..... 4.8 ÷ 8.0
Bromate ..... ≤10 ppm
Chloride..... ≤0.2 %
Water-insoluble matter ..... ≤100 ppm
Iodide ..... ≤50 ppm
Heavy metals (Pb)..... ≤5 ppm
Substanc. not ppt H_2SO_4..... ≤0.2 %
Sulphate ..... ≤50 ppm
Sulphide ..... ≤10 ppm
Ba ..... ≤50 ppm
Ca ..... ≤0.1 %
Cu ..... ≤25 ppm
Fe ..... ≤5 ppm
Ni..... ≤25 ppm
Pb ..... ≤25 ppm
Zn ..... ≤25 ppm
Assay (complexometric) ..... 99 ÷ 100 %

Code	Size	Packaging	Notes
485354	100 g	Glass bottle	

Strontium carbonate

• Stronzio carbonato • Strontium carbonate • Estroncio carbonato • Strontiumcarbonat

$SrCO_3$
Molecular Weight: 147,63
CAS: 1633-05-2
EEC-N: 216-643-7

Strontium carbonate > RPE - For analysis

RPE

Description	White powder	HCl-insoluble matter	≤100 ppm	Ca	≤0.2 %	Zn	≤5 ppm
Identification	Positive	Heavy metals (Pb)	≤60 ppm	Cu	≤5 ppm	Assay (complexometric)	99 ÷ 100 %
Alkalinity(SrOH)	≤300 ppm	Substanc. not ppt H2SO4	≤0.3 %	Fe	≤10 ppm		
Chloride	≤10 ppm	Sulphate	≤50 ppm	Ni	≤5 ppm		
Phosphate	≤10 ppm	Ba	≤200 ppm	Pb	≤5 ppm		

Code	Size	Packaging	Notes
485404	100 g	Glass bottle	
485407	1 kg	Plastic bottle	

Strontium chloride hexahydrate

• Stronzio cloruro esaidrato • Strontium chlorure hexahydraté • Estroncio cloruro hexahidratado • Strontiumchlorid-Hexahydrat

$SrCl_2.6H_2O$
Molecular Weight: 266,62
CAS: 10025-70-4
EEC-N: 233-971-6

Danger
H318
P280i-P305+P351+P338-P310a

Strontium chloride hexahydrate > RPE - For analysis - ACS

RPE

Description	crist. bianchi	Water-insoluble matter	≤50 ppm	Ba	≤500 ppm	Mg	≤2 ppm
Identification	Positive	Heavy metals (Pb)	≤5 ppm	Ca	≤500 ppm	Assay (complexometric)	99.0 ÷ 103.0 %
pH sol. 5% at 25° C	5.0 ÷ 7.0	Sulphate	≤10 ppm	Fe	≤5 ppm		

Code	Size	Packaging	Notes
485455	250 g	Plastic bottle	
485457	1 kg	Plastic bottle	

Strontium nitrate

• Stronzio nitrato • Strontium nitrate • Estroncio nitrato • Strontiumnitrat

$Sr(NO_3)_2$
Molecular Weight: 211,63
CAS: 10042-76-9
EEC-N: 233-131-9

Classification transport
ONU: 1507
Transport Hazard class: 5.1
Packing group III

Danger
H272-H302
P210-P220-P264-P280-P301+P312a-P501a

Strontium nitrate > RPE - For analysis - ACS

RPE

Description	White crystalline powder	Chloride	≤20 ppm	Sulphate	≤50 ppm	Assay (complexometric)	≥99.0 %
Identification	Positive	Water-insoluble matter	≤100 ppm	Ba	≤500 ppm		
pH sol. 5% at 25° C	5.0 ÷ 7.0	Mg and Alkali salts	≤0.15 %	Ca	≤500 ppm		
Loss on drying	≤0.1 %	Heavy metals (Pb)	≤5 ppm	Fe	≤5 ppm		

Code	Size	Packaging	Notes
485605	250 g	Plastic bottle	
485607	1 kg	Plastic bottle	

Strontium nitrate > RE - Pure

RE

Description	White cryst. powder	Chloride	≤500 ppm	Heavy metals (Pb)	≤50 ppm	Fe	≤50 ppm
Identification	Positive	Water-insoluble matter	≤0.1 %	Sulphate	≤500 ppm	Assay (complexometric)	≥98 %

Code	Size	Packaging	Notes
379707	1 kg	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Strontium sulfate

• Stronzio solfato • Strontium sulfate • Estroncio sulfato • Strontiumsulfat

Synonym: *Celestite*

$SrSO_4$
Molecular Weight: 183,68
CAS: 7759-02-6
EEC-N: 231-850-2

Strontium sulfate > RPE - For analysis

RPE

Description	White powder	Chloride	≤20 ppm	Ca	≤0.2 %	Zn	≤5 ppm
Identification	Positive	Heavy metals (Pb)	≤10 ppm	Cu	≤5 ppm	Assay (complexometric)	≥99 %
Loss on ignition	≤1 %	Nitrate	≤0.1 %	Fe	≤10 ppm		
Acidity(Sulphuric acid)	≤200 ppm	Soluble salts	≤0.3 %	Ni	≤5 ppm		
Alkalinity(SrOH)	≤30 ppm	Ba	≤100 ppm	Pb	≤5 ppm		

Code	Size	Packaging	Notes
485705	250 g	Plastic bottle	

Succinic acid

• Acido succinico • Acide succinique • Ácido succínico • Bernsteinsäure

Synonym: *Butanedioic acid*

$HOOC(CH_2)_2COOH$
Molecular Weight: 118,09
CAS: 110-15-6
EEC-N: 203-740-4

Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Succinic acid > RPE - For analysis

RPE

Description	white crystalline powder	Water (K.F.)	≤ 0.5 %	Fe	≤ 5 ppm
Melting point	185 ÷ 190 °C	Assay (complexometry)	≥ 99.5 % s.s	Other organic acid	≤ 0.5 %

Code	Size	Packaging	Notes
411025	250 g	Plastic bottle	
411027	1 kg	Plastic bottle	
411023	25 kg	Plastic bucket	

Succinic anhydride

• Anidride succinica • Anhydride succinique • Anhídrido succínico • Bernsteinsäureanhydrid

Synonym: *Dihydro-2,5-furandione*

$(CH_2CO)_2O$
Molecular Weight: 100,07
CAS: 108-30-5
EEC-N: 203-570-0

Warning
H302-H319-H335
P261-P271-P304+P340-P305+P351+P338-P337+P313-P403+P233

Succinic anhydride > RPE - For analysis

RPE

Description	White crystal powder	Melting point	118 ÷ 120 °C	Heavy metals (Pb)	≤ 20 ppm	Assay (dried base)	≥ 99.0 %
Identification (I.R.)	Positive	Residue on ignition	≤ 500 ppm	Fe	≤ 20 ppm		

Code	Size	Packaging	Notes
422204	100 g	Glass bottle	

Succinic anhydride > RE - Pure

RE

Description	White crystal powder	Melting point	118 ÷ 120 °C	Heavy metals (Pb)	≤ 20 ppm
Identification (I.R.)	Positive	Residue on ignition	≤ 0.1 %	Assay (dried base)	≥ 99 %

Code	Size	Packaging	Notes
318507	1 kg	Plastic bottle	

D(+)-Sucrose

• D(+)-Saccarosio • D(+)-Saccharose • D(+)-Sacarosa • D(+)-Saccharose

Synonym:
α-D-Glucopyranosyl β-D-fructofuranoside

$C_{12}H_{22}O_{11}$
Molecular Weight: 342,3
CAS: 57-50-1
EEC-N: 200-334-9

D(+)-Sucrose > RPE - For analysis - ACS

RPE

Description ..... White crystals	Acidity ..... ≤0.0008 meq/g	Water-insoluble matter ..... ≤50 ppm	Sulphat + sulphit (SO4) ..... ≤50 ppm
Identification ..... Positive	Loss on drying ..... ≤300 ppm	Heavy metals (Pb) ..... ≤5 ppm	Inver.sugar (Glucose) ..... ≤500 ppm
Specific optical rotation... +66.3 ÷ +66.8 °	Chloride ..... ≤50 ppm	Residue on ignition ..... ≤100 ppm	Fe ..... ≤5 ppm

Code	Size	Packaging	Notes
477186	100 g	Plastic bottle	
477187	1 kg	Plastic bottle	
477182	5 kg	Plastic tank	
477183	25 kg	Plastic bucket	

D(+)-Sucrose > ERBApharm - According to pharmacopoeia: BP-FU-NF-Ph.Eur.-Ph.Franc.

ERBApharm

Description ..... White crystalline powder	Dextrine ..... Conform Ph.Eur.	Colore (A) ..... ≤ 45 Ph.Eur.	Sulphite ..... ≤ 10 ppm
Identification ..... Positive	Reducing sugar ..... Conform Ph.Eur.	Specific optical rotation... +66.3 ÷ +67.0 °	Pb ..... ≤ 0.5 ppm
Appearance of solution ..... Conform Ph.Eur.	Conductivity ..... Conform Ph.Eur.	Loss on drying ..... ≤ 0.1 %	

Code	Size	Packaging	Notes
365157	1 kg	Plastic bottle	
365158	5 kg	Plastic tank	
365152	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sudan black B

• Nero Sudan B • Noir Soudan B • Negro Sudan B • Schwarzer Sudan B

Synonym:
Ceres black BN

$C_{29}H_{24}N_6$
Molecular Weight: 456,55
CAS: 4197-25-5
EEC-N: 224-087-1

Sudan black B > RS - For microscopy - C.I. 26150

RS

Description ..... Black powder Identification ..... Positive

Code	Size	Packaging	Notes
464241	25 g	Glass bottle	

Dye for histology

Sudan III

• Sudan III • Soudan III • Sudan III • Sudan III

Synonym:
1-[4-(Phenylazo)phenylazo]-2-naphthol

$C_{22}H_{16}N_4O$
Molecular Weight: 352,4
CAS: 85-86-9
EEC-N: 201-638-4

Sudan III > RS - For microscopy - C.I. 26100

RS

Description ..... Red brick powder Identification ..... Positive

Code	Size	Packaging	Notes
485902	25 g	Glass bottle	

Dye for histology

Product specifications are subject to changes. Please visit our website for updates.

Sudan III hydroalcoholic saturated solution

• Sudan III soluzione idroalcolica satura • Soudan III solution hydroalcoolique
• Sudan III solución hidroalcohólica • Sudan III hydroalkoholische Lösung

Synonym:
1-[4-(Phenylazo)phenylazo]-2-naphthol

$C_{22}H_{16}N_4O$
Molecular Weight: 352,4
CAS: 85-86-9

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II

Danger
H225-H319-H336
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Sudan III hydroalcoholic saturated solution > RS - For the colouring of neutral fats — RS

Description ..... Red clear liquid Identification ..... Positive Density at 20° C ..... 0.855 ÷ 0.861

Code	Size	Packaging	Notes
E485952	250 ml	Glass bottle	

Sudan yellow

• Giallo Sudan • Jaune Soudan • Amarillo Sudán • Sudangelb

Synonym:
• *Sudan I*
• *1-Phenylazo-2-naphthol*

$C_{16}H_{12}N_2O$
Molecular Weight: 248,28
CAS: 842-07-9
EEC-N: 212-668-2

Warning
H317-H341-H351-H413
P261-P280-P308+P313-P362+P364-P333+P313-P501a

Sudan yellow > RPE - For analysis - C.I. 12055 — RPE

Description ..... Red-orange crystalline powder Identification ..... Positive

Code	Size	Packaging	Notes
453581	10 g	Glass bottle	
453582	25 g	Glass bottle	

Sulfamic acid

• Acido solfammico • Acide sulfamique • Acido sulfamico • Sulfaminsäure

Synonym:
Amidosulfonic acid

NH_2SO_3H
Molecular Weight: 97,09
CAS: 5329-14-6
EEC-N: 226-218-8

Classification transport
ONU: 2967
Transport Hazard class: 8
Packing group III

Warning
H315-H319-H412
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Sulfamic acid > RPE - For analysis - ACS — RPE

Description ..... White crystals
Identification ..... Positive
Chloride ..... ≤ 10 ppm
Water-insoluble matter ..... ≤ 100 ppm
Heavy metals (Pb) ..... ≤ 10 ppm
Residue on ignition ..... ≤ 100 ppm
Sulphate ..... ≤ 0.05 %
Fe ..... ≤ 5 ppm
Assay (acidimetric) ..... 99.3 ÷ 100.3 %

Code	Size	Packaging	Notes
410105	250 g	Plastic bottle	
410106	500 g	Plastic bottle	
410104	25 kg	Plastic bucket	

Sulfamic acid > RE - Pure — RE

Description ..... White crystals
Identification ..... Positive
Heavy metals (Pb) ..... ≤20 ppm
Sulphate ..... ≤0.3 %
Fe ..... ≤20 ppm
Assay ..... ≥99.5 %

Code	Size	Packaging	Notes
306507	1 kg	Plastic bottle	
306508	5 kg	Plastic tank	
306503	25 kg	Plastic bucket	

Sulfanilamide
• Solfanilammide • Sulfanilamide • Sulfanilamida • Sulfanilsäureamid

Synonym:
p-Aminobenzenesulfonamide

$C_6H_8O_2N_2S$
Molecular Weight: 172,2
CAS: 63-74-1
EEC-N: 200-563-4

Sulfanilamide > RS - For microanalysis

RS

Description ..... White crystalline powder
Identification ..... Positive
Melting point ..... 164 ÷ 167 ° C
Loss on drying ..... ≤ 0.5 %
Assay ..... ≥ 97.5 % (s.s.)

Code	Size	Packaging	Notes
485961	2 g	Glass bottle	

Sulfanilamide > RPE - For analysis

RPE

Description ..... White to off white crystals or powder
Assay ..... ≥ 99.0 % (ds)
IR ..... Conforms to structure
Acidity ..... Passes test
Loss on drying 105°C ..... ≤ 0.5 %
Sulphated ash ..... ≤ 0.1 %

Code	Size	Packaging	Notes
485971	100 g	Plastic bottle	

Sulfanilic acid
• Acido solfanilico • Acide sulfanilique • Acido sulfanílico • Sulfanilsäure

Synonym:
• *4-Aminobenzenesulfonic acid*
• *Aniline-4-sulfonic acid*

$4\text{-}NH_2C_6H_4SO_3H$
Molecular Weight: 173,19
CAS: 121-57-3
EEC-N: 204-482-5

Warning
H315-H319-H317
P261-P264-P280a-P305+P351+P338-P362+P364-P337+P313

Sulfanilic acid > RS - For analysis according to Ph. Eur. Chap. 4.2.1

RS

Code	Size	Packaging	Notes
612000700	50 g	Plastic bottle	Ref Ph.Eur 2000700

Sulfanilic acid > RPE - For analysis - ACS - Reag. Ph.Eur.

RPE

Description ..... White powder
Identification ..... Positive
Sodium carbonate 5% ins ..... ≤200 ppm
Residue on ignition ..... ≤100 ppm
Chloride ..... ≤20 ppm
Nitrite ..... ≤0.5 ppm
Sulphate ..... ≤100 ppm
Assay (acidimetric) ..... 98.0 ÷ 102.0 %

Code	Size	Packaging	Notes
410154	100 g	Plastic bottle	
410156	500 g	Plastic bottle	

Sulfate standard solution
• Solfati standard soluzione • Sulfate solution standard • Solfato, solución patrón • Sulfat-Standardlösung

Sulfate standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615002801	100 ml	Plastic bottle	A 10 ppm solution R1: to dilute according to Ref Ph.Eur 5002801
615002802	100 ml	Plastic bottle	A 100 ppm solution: to dilute according to Ref Ph.Eur 5002802
615002809	100 ml	Plastic bottle	A 10 ppm solution: to dilute according to Ref Ph.Eur 5002800

Product specifications are subject to changes. Please visit our website for updates.

Sulfate standard solution > RS - Standard solution for ion chromatography

RS

Code	Size	Packaging	Notes
503351	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
503353	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Sulfite standard solution

• Solfito standard soluzione • Sulfite standard solution • Sulfito, solución patrón • Sulfit-Standardlösung

Sulfite standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615002900	100 ml	Plastic bottle	A 1,5 ppm solution Ref Ph.Eur 5002900

Sulfolane

• Sulfolano • Sulfolane • Sulfolano • Sulfolano

Synonym:
Tetrahydrothiophene 1,1-dioxide

$CH_2(CH_2)_3.SO_2$
Molecular Weight: 120,17
CAS: 126-33-0
EEC-N: 204-783-1

Warning
H302
P264-P270-P301+P312a-P330-P501a

Sulfolane > RS - Anhydrous - For analysis

RS

Water content (K.F.) ..... ≤ 200 mg/Kg Assay (GC) ..... ≥ 98.5 % Colourless to light yellow appearance ..... Conform

Code	Size	Packaging	Notes
P932SP16	1 l	Glass bottle	

Sulfomolybdic reagent

• Reattivo solfomolibdico • Réactif sulfomolybdique • Reactivo sulfomolibdico • Sulfomolybdat-Reagens

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfomolybdic reagent > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611086500	100 ml	Plastic bottle	Sulfomolybdic reagent R3 Ref Ph.Eur 1086500

Sulfosalicylic acid

• Acido solfosalicilico • Acide sulfosalycilique • Acido sulfosalicílico • Sulfosalicylsäure

Synonym:
2-Hydroxy-5-sulfobenzoic acid

$HO.C_6H_3(COOH)SO_3H.2H_2O$
Molecular Weight: 254,2
CAS: 5965-83-3
EEC-N: 202-555-6

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Sulfosalicylic acid > RPE - For analysis - ACS

RPE

Description ..... Solido bianco
Identification ..... Positive
Salicilic acid ..... ≤ 0.04 %
Chloride ..... ≤ 10 ppm
Water-insoluble matter ..... ≤ 0.02 %
Heavy metals (Pb) ..... ≤ 20 ppm
Residue on ignition ..... ≤ 0.1 %
Sulphate ..... ≤ 0.02 %
Fe ..... ≤ 10 ppm
Assay (acidimetric) ..... 99.0 ÷ 101.0 %

Code	Size	Packaging	Notes
410894	100 g	Glass bottle	
410896	500 g	Plastic bottle	

Sulfur sublimed and washed

• Zolfo sublimato lavato ventilato • Soufre sublimé lavé ventilé • Azufre sublimado lavado • Schwefel sublimiert und gewaschen

S
Molecular Weight: 32,06
CAS: 7704-34-9
EEC-N: 231-722-6

Classification transport
ONU: 1350
Transport Hazard class: 4.1
Packing group III

Warning
H315
P264-P280g-P302+P352a-P332+P313-P362+P364

Sulfur sublimed and washed > RE - Pure

RE

Description ..... Yellow powder
Identification ..... Positive
Residue on ignition ..... ≤ 0.1 %
Acidity (H_2SO_4) ..... ≤ 0.1 %
Assay (gravimetric) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
378807	1 kg	Plastic bottle	
378809	5 kg	Plastic tank	
378802	25 kg	Plastic bucket	

Sulfur standard solution

• Zolfo standard soluzione • Soufre solution standard • Azufre, solución patrón • Schwefel-Standardlösung

Sulfur standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505822	100 ml	Plastic bottle	conc. 10 ppm Matrix: Water
505825	100 ml	Plastic bottle	conc. 100 ppm Matrix: Water
505823	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Sulfur standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
504291	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
504293	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
504295	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Water
504297	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Water

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Sulfuric acid 98%

• Acido solforico 98% • Acide sulfurique 98% • Acido sulfúrico 98% • Schwefelsäure 98%

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9
EEC-N: 231-639-5

Classification transport
ONU: 1830
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid 98% > RS - For microanalysis

RS

Description ..... Clear liquid
Identification ..... Positive
Density at 20° C ..... ~ 1.835
Assay ..... 95 ÷ 98 %

Code	Size	Packaging	Notes
410421	1 l	Glass bottle	

Sulfuric acid 98% > RPE - For nitrogen dosing

RPE

Description ..... Clear colourless liquid
Density at 20°C ..... ~ 1.84
Total nitrogen (N) ..... ≤ 2 ppm
Assay (acidimetric) ..... 95 - 98 %

Code	Size	Packaging	Notes
502641	2.5 l	Glass bottle	

Content is guaranteed for standardized volumes at 20 °C

Product specifications are subject to changes. Please visit our website for updates.

Sulfuric acid 96% (66°Be)

• Acido solforico 96% (66°Bé) • Acide sulfurique 96% (66°Bé) • Acido sulfúrico 96% (66°Bé) • Schwefelsäure 96% (66° Be)

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9
EEC-N: 231-639-5

Classification transport
ONU: 1830
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid 96% (66°Be) > RE - Pure

RE

Description Clear or opaline colourless liquid
Identification ..... Positive
Density at 20° C ..... 1.830 ÷ 1.836
Chloride ..... ≤500 ppm
Cu ..... ≤ 0.02 ppm
Assay (acidimetric) ..... ≥96 %
Pb ..... ≤ 0.03 ppm
Zn ..... ≤ 0.25 ppm

Code	Size	Packaging	Notes
306751	2.5 l	Glass bottle	
306755	25 kg	Plastic bucket	
306752	50 kg	Plastic drum	

Sulfuric acid 96%

• Acido solforico 96% • Acide sulfurique 96% • Acido sulfurico 96% • Schwefelsäure 96%

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9
EEC-N: 231-639-5

Classification transport
ONU: 1830
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid 96% > RS - VLSI - For electronic use

RS

Code	Size	Packaging	Notes
527631	1 l	Plastic bottle	
527630	2.5 l	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Sulfuric acid 96% > RS - RSE - For electronic use

RS

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Density at 20° C ..... 1.834 ÷ 1.836
Assay (acidimetric) ..... 95.0 ÷ 97.0 %
Ammonium ..... ≤0.5 ppm
Chloride ..... ≤0.1 ppm
Heavy metals (Pb) ..... ≤0.4 ppm
Nitrate ..... ≤0.1 ppm
Phosphate ..... ≤0.5 ppm
Residue on ignition ..... ≤3 ppm
Subst. reducing KMn04 ..... ≤2 ppm

Ag ..... ≤0.02 ppm
Al ..... ≤0.05 ppm
As ..... ≤0.005 ppm
Au ..... ≤0.05 ppm
B ..... ≤0.01 ppm
Ba ..... ≤0.05 ppm
Be ..... ≤0.02 ppm
Bi ..... ≤0.02 ppm
Ca ..... ≤0.2 ppm
Cd ..... ≤0.005 ppm
Co ..... ≤0.01 ppm
Cr ..... ≤0.01 ppm

Cu ..... ≤0.01 ppm
Fe ..... ≤0.1 ppm
Ga ..... ≤0.02 ppm
In ..... ≤0.02 ppm
K ..... ≤0.1 ppm
Li ..... ≤0.02 ppm
Mg ..... ≤0.1 ppm
Mn ..... ≤0.01 ppm
Mo ..... ≤0.01 ppm
Na ..... ≤0.5 ppm
Ni ..... ≤0.01 ppm
Pb ..... ≤0.02 ppm

Pt ..... ≤0.05 ppm
Sb ..... ≤0.01 ppm
Se ..... ≤0.5 ppm
Sn ..... ≤0.02 ppm
Sr ..... ≤0.02 ppm
Ti ..... ≤0.05 ppm
Tl ..... ≤0.05 ppm
V ..... ≤0.05 ppm
Zn ..... ≤0.02 ppm
Zr ..... ≤0.05 ppm

Code	Size	Packaging	Notes
410374	1 l	Glass bottle	
410371	2.5 l	Glass bottle	

Sulfuric acid 96% > RS - MOS - For electronic use

RS

Description	Clear liquid	Ag	≤0.02 ppm	Cu	≤0.01 ppm	Pt	≤0.05 ppm
Colour (APHA)	≤10	Al	≤0.05 ppm	Fe	≤0.1 ppm	Sb	≤0.01 ppm
Identification	Positive	As	≤0.005 ppm	Ga	≤0.02 ppm	Se	≤0.5 ppm
Density at 20° C	1.834 ÷ 1.836	Au	≤0.05 ppm	In	≤0.02 ppm	Sn	≤0.02 ppm
Assay (acidimetric)	95.0 ÷ 97.0 %	B	≤0.01 ppm	K	≤0.1 ppm	Sr	≤0.02 ppm
Ammonium	≤0.5 ppm	Ba	≤0.05 ppm	Li	≤0.02 ppm	Ti	≤0.05 ppm
Chloride	≤0.1 ppm	Be	≤0.02 ppm	Mg	≤0.1 ppm	Tl	≤0.05 ppm
Heavy metals (Pb)	≤0.4 ppm	Bi	≤0.02 ppm	Mn	≤0.01 ppm	V	≤0.05 ppm
Nitrate	≤0.1 ppm	Ca	≤0.2 ppm	Mo	≤0.01 ppm	Zn	≤0.02 ppm
Phosphate	≤0.5 ppm	Cd	≤0.005 ppm	Na	≤0.5 ppm	Zr	≤0.05 ppm
Residue on ignition	≤3 ppm	Co	≤0.01 ppm	Ni	≤0.01 ppm		
Subst. reducing KMn04	≤2 ppm	Cr	≤0.01 ppm	Pb	≤0.02 ppm		

Code	Size	Packaging	Notes
410382	1 l	Glass bottle	
410381	2.5 l	Glass bottle	

Sulfuric acid 96% > RS - For enviromental analysis - ACS - Reag. Ph.Eur. - Reag. USP

RS

Description	Clear oily liquid	Subst. reducing KMn04	≤2 ppm	Cu	≤0.01 ppm	Se	≤3 ppm
Colour (APHA)	≤10	Ag	≤0.02 ppm	Fe	≤0.1 ppm	Sr	≤0.02 ppm
Colour of 2N solution(APHA)	≤10	Al	≤0.05 ppm	Hg	≤0.005 ppm	Ti	≤0.05 ppm
Identification	Positive	As	≤0.005 ppm	K	≤0.1 ppm	Tl	≤0.05 ppm
Density at 20° C	1.834 ÷ 1.836	Ba	≤0.1 ppm	Li	≤0.02 ppm	V	≤0.05 ppm
Ammonium	≤1 ppm	Be	≤0.02 ppm	Mg	≤0.1 ppm	Zn	≤0.05 ppm
Chloride	≤0.1 ppm	Bi	≤0.02 ppm	Mn	≤0.01 ppm	Zr	≤0.05 ppm
Phosphate	≤0.5 ppm	Ca	≤0.2 ppm	Mo	≤0.01 ppm	Assay (acidimetric)	95.0 ÷ 97.0 %
Heavy metals (Pb)	≤0.8 ppm	Cd	≤0.005 ppm	Na	≤0.5 ppm		
Nitrate	≤0.2 ppm	Co	≤0.01 ppm	Ni	≤0.02 ppm		
Residue on ignition	≤4 ppm	Cr	≤0.05 ppm	Pb	≤0.02 ppm		

Code	Size	Packaging	Notes
410261	1 l	Glass bottle	

Low content in Hg

Sulfuric acid 96% > RPE - For analysis - ISO

RPE

Description	Clear oily liquid	Heavy metals (Pb)	≤0.8 ppm	Cu	≤0.01 ppm	Pb	≤0.02 ppm
Colour (APHA)	≤10	Nitrate	≤0.2 ppm	Fe	≤0.1 ppm	Se	≤3 ppm
Colour of 2N solution(APHA)	≤10	Residue on ignition	≤4 ppm	K	≤0.1 ppm	Sr	≤0.02 ppm
Identification	Positive	Subst. reducing KMn04	≤2 ppm	Li	≤0.02 ppm	Zn	≤0.05 ppm
Density at 20° C	1.834 ÷ 1.836	As	≤0.01 ppm	Mg	≤0.2 ppm	Assay (acidimetric)	95 ÷ 97 %
Ammonium	≤1 ppm	Ca	≤0.2 ppm	Na	≤0.5 ppm		
Chloride	≤0.1 ppm	Cd	≤0.005 ppm	Ni	≤0.05 ppm		

Code	Size	Packaging	Notes
410301	1 l	Glass bottle	
410303	1 l	Glass bottle PVC coated	
524540	1 l	Plastic bottle	
410306	2.5 l	Glass bottle	
524541	2.5 l	Plastic bottle	
410304	5 l	Plastic tank	
524543	25 l	Plastic tank	
410307	30 kg	Plastic drum	
410302	50 kg	Plastic drum	

Content is guaranteed for standardized volumes at 20 °C

Product specifications are subject to changes. Please visit our website for updates.

Sulfuric acid 96% > ERBApharm - According to pharmacopoeia: BP-NF-Ph.Eur.

ERBApharm

Description ..... Clear colourless liquid	Nitrate ..... Conform Ph.Eur.	Chloride ..... ≤ 50 ppm	Fe ..... ≤ 25 ppm
Identification ..... Positive	Density at 20° C ..... ~ 1.84	Heavy metals (Pb) ..... ≤ 5 ppm	Assay (acidimetric) ..... 95.0 ÷ 98.0 %
Appearance of solution ..... Conform Ph.Eur.	Sulphated ash ..... ≤ 50 ppm	As ..... ≤ 1 ppm	Subst. reducing KMnO4 ... Conform USP-NF

Code	Size	Packaging	Notes
306651	1 l	Glass bottle	
306657	2.5 l	Glass bottle	
306653	50 kg	Plastic drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Sulfuric acid 95-97%

• Acido solforico 95-97% • Acide sulfurique 95-97% • Acido sulfúrico 95-97% • Schwefelsäure 95-97%

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9
EEC-N: 231-639-5

Classification transport
ONU: 1830
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid 95 - 97 % > RPE - For nitrogen dosing - Reag.Ph.Eur.

RPE

Description ..... Clear oily liquid	Ammonium ..... ≤ 1 ppm	Residue on ignition ..... ≤ 4 ppm	Subs. reducing KMnO4 ..... Pass test
Identification ..... Positive	Chloride ..... ≤ 0.1 ppm	As ..... ≤ 0.01 ppm	Test nitrate (Ph.Eur.) ..... Pass test
Colour ..... ≤ 10 APHA	Heavy metals (Pb) ..... ≤ 0.8 ppm	Fe ..... ≤ 0.1 ppm	Assay (acidimetric) ..... 95.0 ÷ 97.0 %
Density at 20°C ..... 1.834 ÷ 1.836	Nitrate ..... ≤ 0.2 ppm	Total nitrogen ..... ≤ 2 ppm	

Code	Size	Packaging	Notes
502302	1 l	Glass bottle	

Sulfuric acid 93-98%

• Acido solforico 93-98% • Acide sulfurique 93-98% • Acido sulfúrico 93-98% • Schwefelsäure 93-98%

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9
EEC-N: 231-639-5

Classification transport
ONU: 1830
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid 93-98% > RS - Ultrapure - For trace analysis at ppt level

RS

Description ..... Clear colourless liquid	Mg ..... ≤ 50 ppt	Dy ..... ≤ 10 ppt	Rb ..... ≤ 10 ppt
Identification ..... Positive	Mn ..... ≤ 10 ppt	Er ..... ≤ 10 ppt	Sm ..... ≤ 10 ppt
Ag ..... ≤ 50 ppt	Mo ..... ≤ 10 ppt	Eu ..... ≤ 10 ppt	Sc ..... ≤ 10 ppt
Al ..... ≤ 50 ppt	Na ..... ≤ 50 ppt	Gd ..... ≤ 10 ppt	Se ..... ≤ 500 ppt
As ..... ≤ 500 ppt	Ni ..... ≤ 50 ppt	Ga ..... ≤ 10 ppt	Te ..... ≤ 50 ppt
Ba ..... ≤ 10 ppt	Pb ..... ≤ 10 ppt	Ge ..... ≤ 100 ppt	Tb ..... ≤ 10 ppt
Be ..... ≤ 10 ppt	Sn ..... ≤ 50 ppt	Hf ..... ≤ 10 ppt	Tl ..... ≤ 10 ppt
Bi ..... ≤ 10 ppt	Sr ..... ≤ 10 ppt	Ho ..... ≤ 10 ppt	Tm ..... ≤ 10 ppt
Ca ..... ≤ 50 ppt	Ti ..... ≤ 50 ppt	In ..... ≤ 10 ppt	W ..... ≤ 10 ppt
Cd ..... ≤ 10 ppt	V ..... ≤ 10 ppt	La ..... ≤ 10 ppt	U ..... ≤ 10 ppt
Co ..... ≤ 10 ppt	Zn ..... ≤ 50 ppt	Li ..... ≤ 10 ppt	Yb ..... ≤ 10 ppt
Cr ..... ≤ 10 ppt	Assay (acidimetric) ..... 93 ÷ 98 %	Lu ..... ≤ 10 ppt	Y ..... ≤ 10 ppt
Cu ..... ≤ 10 ppt	Th ..... ≤ 10 ppt	Nd ..... ≤ 10 ppt	Zr ..... ≤ 10 ppt
Fe ..... ≤ 50 ppt	Sb ..... ≤ 50 ppt	Nb ..... ≤ 10 ppt	
Hg ..... ≤ 100 ppt	Ce ..... ≤ 10 ppt	Pr ..... ≤ 10 ppt	
K ..... ≤ 50 ppt	Cs ..... ≤ 10 ppt	Rh ..... ≤ 50 ppt	

Code	Size	Packaging	Notes
410351	500 ml	Plastic bottle	

Sulfuric acid 93-98% > RS - Superpure - For trace analysis at ppb level

RS

Description	Clear liquid	Mn	≤ 0.5 ppb	Cs	≤ 0.1 ppb	Sm	≤ 0.1 ppb
Identification	Positive	Mo	≤ 0.5 ppb	Dy	≤ 0.1 ppb	Sc	≤ 0.1 ppb
Ag	≤ 1 ppb	Na	≤ 1 ppb	Er	≤ 0.1 ppb	Te	≤ 0.1 ppb
Al	≤ 1 ppb	Ni	≤ 0.5 ppb	Eu	≤ 0.1 ppb	Tl	≤ 0.1 ppb
As	≤ 0.5 ppb	Pb	≤ 0.1 ppb	Gd	≤ 0.1 ppb	Tm	≤ 0.1 ppb
Ba	≤ 0.1 ppb	Sb	≤ 1 ppb	Ga	≤ 0.1 ppb	W	≤ 0.5 ppb
Be	≤ 0.1 ppb	Se	≤ 10 ppb	Ge	≤ 1 ppb	U	≤ 0.1 ppb
Bi	≤ 0.1 ppb	Sn	≤ 1 ppb	Hf	≤ 0.1 ppb	Yb	≤ 0.1 ppb
Ca	≤ 1 ppb	Sr	≤ 0.5 ppb	Ho	≤ 0.1 ppb	Y	≤ 0.1 ppb
Cd	≤ 0.5 ppb	Ti	≤ 1 ppb	In	≤ 0.1 ppb	Tb	≤ 0.1 ppb
Co	≤ 0.5 ppb	V	≤ 0.5 ppb	La	≤ 0.1 ppb	Chloride	≤ 0.7 ppm
Cr	≤ 0.5 ppb	Zn	≤ 1 ppb	Lu	≤ 0.1 ppb	Nitrate	≤ 0.2 ppm
Cu	≤ 0.5 ppb	Zr	≤ 0.5 ppb	Nd	≤ 0.1 ppb	Total phosphorus	≤ 0.05 ppm
Hg	≤ 0.1 ppb	Assay (acidimetric)	93 ÷ 98 %	Nb	≤ 0.1 ppb	Reducing substances KMn04	≤ 20 ppm
K	≤ 1 ppb	Th	≤ 0.1 ppb	Pr	≤ 0.1 ppb		
Li	≤ 0.5 ppb	Colour (APHA)	≤ 10	Rh	≤ 0.5 ppb		
Mg	≤ 1 ppb	Ce	≤ 0.1 ppb	Rb	≤ 0.5 ppb		

Code	Size	Packaging	Notes
410405	500 ml	Plastic bottle	
410406	1 l	Plastic bottle	
410407	2.5 l	Plastic bottle	

Sulfuric acid 90%

• Acido solforico 90% • Acide sulfurique 90% • Acido sulfúrico 90% • Schwefelsäure 90%

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9
EEC-N: 231-639-5

Classification transport
ONU: 1830
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid 90% > RS - For analysis according to Gerber

RS

Description ..... Clear colourless liquid Density at 20° C ..... 1.815 - 1.825 Assay (acidimetric) ..... 90 - 92 %

Code	Size	Packaging	Notes
410391	1 l	Plastic bottle	
410394	2.5 l	Glass bottle	

Sulfuric acid 85%

• Acido solforico 85% • Acide sulfurique 85% • Acido sulfúrico 85% • Schwefelsäure 85%

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

Classification transport
ONU: 1830
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid 85% > RS - For analysis

RS

Density d20/4 ..... 1.774 - 1.783 Sulfuric acid content ..... 84.5 - 85.5 %

Code	Size	Packaging	Notes
PS0433/15	1 l	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Sulfuric acid 72%

• Acido solforico 72% • Acide sulfurique 72% • Acido sulfúrico 72% • Schwefelsäure 72%

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

Classification transport
ONU: 1830
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid 72% > RS - For agroalimentary analysis

RS

Description ..... Clear colourless liquid Colour ..... ≤ 10 APHA Density at 20°C ..... 1.629 ÷ 1.639 Assay ..... 71.50 ÷ 72.50 %

Code	Size	Packaging	Notes
502771	2.5 l	Glass bottle	

Sulfuric acid 69%

• Acido solforico 69% • Acide sulfurique 69% • Acido sulfúrico 69% • Schwefelsäure 69%

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

Classification transport
ONU: 1830
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid 69% > RS - For milk analysis

RS

Clear,colourless solution ..... Conform Sulfuric acid content ..... 68.0 - 70.0 % Density d20/4 ..... 1.587- 1.611

Code	Size	Packaging	Notes
PS0893/21	2.5 l	Glass bottle	

Sulfuric acid 62%

• Acido solforico 62% • Acide sulfurique 62% • Acido sulfúrico 62% • Schwefelsäure 62%

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

Classification transport
ONU: 1830
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid 62% > RS - For milk analysis

RS

Clear,colourless solution ..... Conform Sulfuric acid content ..... 61.0 - 63.0 % Density d20/4 ..... 1.509 - 1.531

Code	Size	Packaging	Notes
PS0894/21	2.5 l	Glass bottle	

Sulfuric acid 50%

• Acido solforico 50% • Acide sulfurique 50% • Acido sulfúrico 50% • Schwefelsäure 50%

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

Classification transport
ONU: 2796
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid 50% > RE - Pure

RE

Description ..... Clear colourless liquid Density at 20° C ..... 1.385 ÷ 1.405 Assay (acidimetric) ..... 49.0 ÷ 51.0 %
Identification ..... Positive Density at 20°C ..... 1.385 ÷ 1.405

Code	Size	Packaging	Notes
E306702	1 l	Bottle	
528541	5 l	Plastic tank	
E306704	35 kg	Plastic tank	

Sulfuric acid 35% (30°Be)

• Acido solforico 35% (30°Bé) • Acide sulfurique 35% (30°Be) • Acido sulfúrico 35% (30°Bé) • Schwefelsäure 35% (30°Be)

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

Classification transport
ONU: 2796
Transport Hazard class: 8
Packing group II



Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid 35% (30°Be) > RE - Pure — RE

Description ..... Clear or opaline colourless liquid Density at 20° C ..... 1.252 - 1.260 Assay ..... 34.0 - 35.0 %

Code	Size	Packaging	Notes
307001000	30 kg	Plastic tank	

Sulfuric acid 30%

• Acido solforico 30% • Acide sulfurique 30% • Acido sulfúrico 30% • Schwefelsäure 30%

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

Classification transport
ONU: 2796
Transport Hazard class: 8
Packing group II



Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid 30% > RS - For analysis — RS

Density d20/4 ..... 1.206 - 1.230 H2SO4 content ..... 29 - 31 %

Code	Size	Packaging	Notes
PS0009/15	1 l	Plastic bottle	

Sulfuric acid 25%

• Acido solforico 25% • Acide sulfurique 25% • Acido sulfúrico 25% • Schwefelsäure 25%

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

Classification transport
ONU: 2796
Transport Hazard class: 8
Packing group II



Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid 25% > RE - Pure — RE

Density d20/4 ..... 1.174 - 1.182 Sulfuric acid content ..... 24.5 - 25.5 %

Code	Size	Packaging	Notes
504562	1 l	Glass bottle	
PS0212/21	2.5 l	Glass bottle	

Sulfuric acid 20%

• Acido solforico 20% • Acide sulfurique 20% • Acido sulfúrico 20% • Schwefelsäure 20%

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

Classification transport
ONU: 2796
Transport Hazard class: 8
Packing group II



Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid 20% > RPE - For analysis — RPE

Description ..... Clear colourless liquid Density at 20° C ..... 1.136 - 1.158 Assay ..... 20 - 22 %

Code	Size	Packaging	Notes
410511000	1 l	Plastic bottle	
410516	20 kg	Plastic drum	

Content is guaranteed for standardized volumes at 20°C. Keep tightly

Product specifications are subject to changes. Please visit our website for updates.

Sulfuric acid 10% v/v

• Acido solforico 10% v/v • Acide sulfurique 10% v/v • Acido sulfúrico 10% V/V • Schwefelsäure 10% v/v

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

Classification transport
ONU: 2796
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid 10% v/v > RPE - For analysis

RPE

Density at 20°C ..... 1.12 - 1.14

Code	Size	Packaging	Notes
502591	1 l	Bottle	

Sulfuric acid 4 mol/l (8N)

• Acido solforico 4 mol/l (8N) • Acide sulfurique 4 mol/l (8N) • Acido sulfúrico 4 mol/l (8N) • Schwefelsäure 4 mol/l (8N)

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

Classification transport
ONU: 2796
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid 4 mol/l (8N) > RS - For environmental analysis (COD determination)

RS

Description ..... Clear colourless liquid Assay ..... 7.984 - 8.016 N

Code	Size	Packaging	Notes
526741	1 l	Bottle	

Sulfuric acid 2.5 mol/l (5N)

• Acido solforico 2.5 mol/l (5N) • Acide sulfurique 2.5 mol/l (5N) • Acido sulfúrico 2.5 mol/l (5N) • Schwefelsäure 2.5 mol/l (5N)

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

Classification transport
ONU: 2796
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid 2.5 mol/l (5N) > RPE - For analysis

RPE

Assay (potentiomety) ..... 4.995 - 5.005 N

Code	Size	Packaging	Notes
P3240015	1 l	Plastic bottle	

Sulfuric acid 1 mol/l (2N)

• Acido solforico 1 mol/l (2N) • Acide sulfurique 1 mol/l (2N) • Acido sulfúrico 1 mol/l (2N) • Schwefelsäure 1 mol/l (2N)

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

Classification transport
ONU: 2796
Transport Hazard class: 8
Packing group II

Warning
H290-H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Sulfuric acid 1 mol/l (2N) > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
617000281	1 l	Bottle	

Sulfuric acid 1 mol/l (2N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 1.998 - 2.002 N NIST 723 ..... e

Code	Size	Packaging	Notes
410547000	1 l	Plastic bottle	Certified with NIST traceability
410548000	10 l	Plastic tank	Certified with NIST traceability

98.06 g of H2SO4. Volumetric solution ready-to-use

Sulfuric acid 0.5 mol/l (1N)

• Acido solforico 0.5 mol/l (1N) • Acide sulfurique 0,5 mol/l (1N) • Acido sulfúrico 0.5 mol/l (1N) • Schwefelsäure 0.5 mol/l (1N)

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

HEU210

Sulfuric acid 0.5 mol/l (1N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613007800	1 l	Plastic bottle	Ref Ph.Eur 3007800

Sulfuric acid 0.5 mol/l (1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.999 - 1.001 N NIST 723 ..... e

Code	Size	Packaging	Notes
410577000	1 l	Plastic bottle	Certified with NIST traceability
410572000	5 l	Kubidos	Certified with NIST traceability
410575000	5 l	Plastic tank	Certified with NIST traceability
410571000	10 l	Kubidos	Certified with NIST traceability

49.03 g of H2SO4. Volumetric solution ready-to-use

Sulfuric acid 0.5 mol/l (1N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
410591		Plastic ampoule	Volume: 165 ml

49,03 g of H2SO4 . Volumetric concentrated solution to prepare 1 L of solution 1 N

Sulfuric acid 0.33 mol/l (2N/3)

• Acido solforico 0.33 mol/l (2N/3) • Acide sulfurique 0.33 mol/l (2N/3) • Acido sulfúrico 0.33 mol/l (2N/3) • Schwefelsäure 0.33 mol/l (2N/3)

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

HEU210

Sulfuric acid 0.33 mol/l (2N/3) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.6653 - 0.6680 N

Code	Size	Packaging	Notes
410634	1 l	Plastic bottle	

32,363 g of H2SO4. Volumetric solution ready-to-use. Content is guaranteed for standardized volumes at 20 °C

Product specifications are subject to changes. Please visit our website for updates.

Sulfuric acid 0.26 mol/l (0.52N)

• Acido solforico 0.26 mol/l (0.52N) • Acide sulfurique 0.26 mol/l (0.52N) • Acido sulfúrico 0.26 mol/l (0.52N) • Schwefelsäure 0.26 mol/l (0.52N)

H_2SO_4 HEU210
Molecular Weight: 98,08
CAS: 7664-93-9

Sulfuric acid 0.26 mol/l (0.52N) > RS - For agroalimentary analysis

RS

Description ..... Clear colourless liquid Aspect ..... Conform Assay ..... 0.515 ÷ 0.525 N

Code	Size	Packaging	Notes
502202	5 l	Plastic tank	

Sulfuric acid 0.25 mol/l (0.5N)

• Acido solforico 0.25 mol/l (0.5N) • Acide sulfurique 0.25 mol/l (0.5N) • Acido sulfúrico 0.25 mol/l (0.5N) • Schwefelsäure 0.25 mol/l (0.5N)

H_2SO_4 HEU210
Molecular Weight: 98,08
CAS: 7664-93-9

Sulfuric acid 0.25 mol/l (0.5N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.4995 - 0.5005 N NIST 723 ..... e
Identification ..... Positive Titration factor ..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
410667000	1 l	Plastic bottle	Certified with NIST traceability
410663000	5 l	Kubidos	Certified with NIST traceability
410662000	10 l	Kubidos	Certified with NIST traceability

24.52 g of H2SO4. Volumetric solution ready-to-use

Sulfuric acid 0.25 mol/l (0.5N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
410681		Plastic ampoule	Volume: 55 ml

24,52 g H2SO4. Volumetric concentrated solution to prepare 1 L of solution 0,5 N

Sulfuric acid 0.166 mol/l (0.333N)

• Acido solforico 0.166 mol/l (0.333N) • Acide sulfurique 0.166 mol/l (0.333N) • Acido sulfúrico 0.166 mol/l (0.333N) • Schwefelsäure 0.166 mol/l (0.333N)

H_2SO_4 HEU210
Molecular Weight: 98,08
CAS: 7664-93-9

Sulfuric acid 0.166 mol/l (0.333N) > RPE - For analysis

RPE

Assay (potentiometry) ..... 0.3331 - 0.3337 N

Code	Size	Packaging	Notes
PS0217/15	1 l	Plastic bottle	

Sulfuric acid 0.13 mol/l (0.26N)

• Acido solforico 0.13 mol/l (0.26N) • Acide sulfurique 0.13 mol/l (0.26N) • Acido sulfúrico 0.13 mol/l (0.26N) • Schwefelsäure 0.13 mol/l (0.26N)

H_2SO_4 HEU210
Molecular Weight: 98,08
CAS: 7664-93-9

Sulfuric acid 0.13 mol/l (0.26N) > RS - For agroalimentary analysis

RS

Description ..... Clear colourless liquid Assay ..... 0.255 ÷ 0.265 N

Code	Size	Packaging	Notes
502651	5 l	Plastic tank	

Content is guaranteed for standardized volumes at 20 °C. Keep tightly closed

Sulfuric acid 0.125 mol/l (0.25N)

• Acido solforico 0.125 mol/l (0.25N) • Acide sulfurique 0.125 mol/l (0.25N) • Acido sulfúrico 0.125 mol/l (0.25N) • Schwefelsäure 0.125 mol/l (0.25N)

H_2SO_4 HEU210
Molecular Weight: 98,08
CAS: 7664-93-9

Sulfuric acid 0.125 mol/l (0.25N) > RPE - For analysis

RPE

Assay (potentiometry) ..... 0.2495 - 0.2505 N

Code	Size	Packaging	Notes
PS0445/22	5 l	Plastic tank	

Sulfuric acid 0.1 mol/l (0.2N)

• Acido solforico 0.1 mol/l (0.2N) • Acide sulfurique 0.1 mol/l (0.2N) • Acido sulfúrico 0.1 mol/l (0.2N) • Schwefelsäure 0.1 mol/l (0.2N)

H_2SO_4 HEU210
Molecular Weight: 98,08
CAS: 7664-93-9

Sulfuric acid 0.1 mol/l (0.2N) > RS - For agroalimentary analysis

RS

Description ..... Clear liquid Aspect ..... Conform Assay ..... 0.195 ÷ 0.205 N

Code	Size	Packaging	Notes
502100000	1 l	Plastic bottle	

Sulfuric acid 0.05 mol/l (0.1N)

• Acido solforico 0.05 mol/l (0.1N) • Acide sulfurique 0.05 mol/l (0.1N) • Acido sulfúrico 0.05 mol/l (0.1N) • Schwefelsäure 0.05 mol/l (0.1N)

H_2SO_4 HEU210
Molecular Weight: 98,08
CAS: 7664-93-9

Sulfuric acid 0.05 mol/l (0.1N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613008001	500 ml	Plastic bottle	Ref Ph.Eur 3008000
613008000	1 l	Plastic bottle	Ref Ph.Eur 3008000

Product specifications are subject to changes. Please visit our website for updates.

Sulfuric acid 0.05 mol/l (0.1N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.0999 - 0.1001 N NIST 723 ..... e

Code	Size	Packaging	Notes
410717000	1 l	Plastic bottle	Certified with NIST traceability
410712000	5 l	Kubidos	Certified with NIST traceability
410711000	10 l	Kubidos	Certified with NIST traceability
410715000	10 l	Plastic tank	Certified with NIST traceability
410714000	20 l	Plastic tank	Certified with NIST traceability

4.904 g of H2SO4. Volumetric solution ready-to-use

Sulfuric acid 0.05 mol/l (0.1N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor ..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
410731		Plastic ampoule	Volume: 55 ml

4,904 g of H2SO4. Volumetric concentrated solution to prepare 1 L of solution 0,1 N

Sulfuric acid 0.025 mol/l (0.05N)

• Acido solforico 0.025 mol/l (0.05N) • Acide sulfurique 0.025 mol/l (0.05N) • Acido sulfúrico 0.025 mol/l (0.05N) • Schwefelsäure 0.025 mol/l (0.05N)

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

HEU210

Sulfuric acid 0.025 mol/l (0.05N) > RPE - For analysis

RPE

Assay (potentiometry) ..... 0.0499 - 0.0501 N

Code	Size	Packaging	Notes
PS0016/96	10 l	Kubidos	

Sulfuric acid 0.02 mol/l (0.04N)

• Acido solforico 0.02 mol/l (0.04N) • Acide sulfurique 0.02 mol/l (0.04N) • Acido sulfúrico 0.02 mol/l (0.04N) • Schwefelsäure 0.02 mol/l (0.04N)

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

HEU210

Sulfuric acid 0.02 mol/l (0.04N) > RPE - For analysis

RPE

Assay (potentiometry) ..... 0.03992 - 0.04008 N

Code	Size	Packaging	Notes
PS0219/15	1 l	Plastic bottle	
PS0219/95	5 l	Kubidos	
PS0219/96	10 l	Kubidos	

Sulfuric acid 0.01 mol/l (0.02N)

• Acido solforico 0.01 mol/l (0.02N) • Acide sulfurique 0.01 mol/l (0.02N) • Acido sulfúrico 0.01 mol/l (0.02N) • Schwefelsäure 0.01 mol/l (0.02N)

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

HEU210

Sulfuric acid 0.01 mol/l (0.02N) > RPE - For analysis

RPE

Assay (potentiometry) ..... 0.01996 - 0.02004 N

Code	Size	Packaging	Notes
PS0047/15	1 l	Plastic bottle	

Sulfuric acid 0.005 mol/l (0.01N)

• Acido solforico 0.005 mol/l (0.01N) • Acide sulfurique 0.005 mol/l (0.01N) • Acido sulfúrico 0.005 mol/l (0.01N) • Schwefelsäure 0.005 mol/l (0.01N)

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

HEU210

Sulfuric acid 0.005 mol/l (0.01N) > RPE - For analysis

RPE

Assay (potentiomety) ..... 0.00998 - 0.01002 N

Code	Size	Packaging	Notes
PS0026/95	5 l	Kubidos	

Sulfuric acid 0.005 mol/l (0.01N) > RPE - NORMEX - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
410791		Plastic ampoule	Volume: 55 ml

0,490 g of H2SO4. Volumetric concentrated solution to prepare 1 L of solution 0,01 N

Sulfuric acid 0.0025 mol/l (0.005N)

• Acido solforico 0.0025 mol/l (0.005N) • Acide sulfurique 0.0025 mol/l (0.005N) • Acido sulfúrico 0.0025 mol/l (0.005N)
• Schwefelsäure 0.0025 mol/l (0.005N)

H_2SO_4
Molecular Weight: 98.08
CAS: 7664-93-9

Sulfuric acid 0.0025 mol/l (0.005N) > RS - For analysis

RS

Assay (potentiometry) ..... 0.00495 - 0.00505 N

Code	Size	Packaging	Notes
424111	10 l	Kubidos	

Content is guaranteed for standardized volumes at 20 °C

Sulfuric acid with 10 g/l Silver sulfate

• Argento solfato 10 g/l in acido solforico • Acide sulfurique à 10 g/l d'argent sulfate • Acido sulfúrico con 10 g/l de plata sulfato
• Schwefelsäure 10g/l Silbersulfat

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

Classification transport
ONU: 1830
Transport Hazard class: 8
Packing group II

Danger
H290-H314-H412
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid with 10 g/l Silver sulfate > RS - For environmental analysis (COD determination)

RS

Assay ..... 9.0 ÷ 11.0 g/l

Code	Size	Packaging	Notes
526605	1 l	Glass bottle	
526606	2.5 l	Glass bottle	

According to NF T90101 of 02/2001

Product specifications are subject to changes. Please visit our website for updates.

Sulfuric acid with 6.6 g/l Silver sulfate

• Argento solfato 6.6 g/l in acido solforico • Acide sulfurique à 6.6 g/l d'argent sulfate • Acido sulfúrico con 6.6 g/l de plata sulfato • Schwefelsäure 6.6 g/l Silbersulfat

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

Classification transport
ONU: 1830
Transport Hazard class: 8
Packing group II

Danger
H290-H314-H412
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid with 6.6 g/l Silver sulfate > RS - For environmental analysis (COD determination)

RS

Assay ..... 6.0 ÷ 7.2 g/l

Code	Size	Packaging	Notes
526602	2.5 l	Glass bottle	

According to NF T90101 of 02/2001

Sulfuric acid, dilute

• Acido solforico, diluito • Acide sulfurique dilué • Acido sulfúrico, diluido • Verdünnte Schwefelsäure

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9

Classification transport
ONU: 2796
Transport Hazard class: 8
Packing group II

Warning
H290-H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Sulfuric acid, dilute > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611086804	1 l	Plastic bottle	A 98 g/l solution Ref Ph.Eur 1086804

Sulfuric acid d=1.820

• Acido solforico d=1,820 • Acide sulfurique d=1,820 • Acido sulfúrico d=1,820 • Schwefelsäure d=1.820

H_2SO_4
Molecular Weight: 98,08
CAS: 7664-93-9
EEC-N: 231-639-5

Classification transport
ONU: 1830
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Sulfuric acid d=1.820 > RS - For agroalimentary analysis

RS

Description ..... Clear colourless liquid Colour (APHA) ..... ≤ 10 APHA Density at 20°C ..... 1.815 ÷ 1.825

Code	Size	Packaging	Notes
502020	5 l	Plastic tank	

According to NF V04-263 and V04-210

Talc
• Talco • Talc • Talco • Talkum

Synonym:
Hydrous magnesium silicate

$3MgO.4SiO_2.H_2O$
Molecular Weight: 379,29
CAS: 14807-96-6
EEC-N: 238-877-9

Talc > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP

ERBApharm

Description	White powder	Al	≤ 2.0 %	Pb	≤ 10 ppm	TAMC	≤ 100 CFU/g
Identification	Positive	Ca	≤ 0.90 %	Loss on ignition	≤ 7.0 %	TYMC	≤ 50 CFU/g
Acidity or alkalinity	Conform	Fe	≤ 0.25 %	Asbestos	Absent		
Sostanze idrosolubili	≤ 0.1 %	Mg	17.0 ÷ 19.5 %	Microbial tests			

Code	Size	Packaging	Notes
382107	1 kg	Plastic bottle	
382109	5 kg	Plastic tank	
382105	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Tannic acid
• Acido tannico • Acide tannique • Acido tánico • Gerbsäure

Synonym:
Gallotannin

$C_{76}H_{52}O_{46}$
Molecular Weight: 1701,23
CAS: 1401-55-4
EEC-N: 215-753-2

Warning
H302-H412
P264-P270-P273-P301+P312a-P330-P501a

Tannic acid > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU

ERBApharm

Description	Yellowish powder	Loss on drying	≤ 12.0 %	Resins	Conform Ph.Eur.
Identification	Positive	Sulphated ash	≤ 0.1 %	Heavy metals (as Pb)	≤ 40 ppm
Appearance of solution	Conform Ph.Eur.	Dextrins, gum, salts, sugars	Conform Ph.Eur.	As	≤ 3 ppm

Code	Size	Packaging	Notes
307157	1 kg	Bag	
307152	5 kg	Plastic bucket	
307153	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Tannic acid > RE - Pure

RE

Description	Brown powder	Water solubility	Conform	Water	< 7 %	Assay (gravimetric)	> 94 % s.s.
Identification	Positive	Alcohol solubility	Conform	Sulphated ash	< 0.3 %		

Code	Size	Packaging	Notes
411074	100 g	Plastic bottle	
411076	500 g	Plastic bottle	

Tantalum standard solution
• Tantalio standard soluzione • Tantale solution standard • Tántalo, solución patrón • Tantal-Standardlösung

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group III



Tantalum standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505872	100 ml	Plastic bottle	conc. 10 ppm Matrix: Nitric acid and hydrofluoric acid
505875	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid and hydrofluoric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Product specifications are subject to changes. Please visit our website for updates.

Tantalum standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503961	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid
503963	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid
503965	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrofluoric acid and nitric acid
503967	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrofluoric acid and nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Tantalum standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507761	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507517	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

L(+) Tartaric Acid

• Acido L(+)tartarico • Acide L(+)tartrique • Acido L(+)tartárico • L(+)-Weinsäure

$HOOC(CHOH)_2COOH$
Molecular Weight: 150,09
CAS: 87-69-4
EEC-N: 201-766-0

Danger
H318
P280i-P305+P351+P338-P310a

L(+) Tartaric Acid > RPE - For analysis - ACS - ISO

RPE

Description ..... White crystalline powder	Residue on ignition .....≤200 ppm	Oxalate .....Conform	Fe .....≤5 ppm
Identification .....Positive	Chloride.....≤10 ppm	Sulfur compounds (as SO4).....≤20 ppm	Assay (acidimetric) ..... ≥99.0 %
Water-insoluble matter .....≤50 ppm	Phosphate .....≤10 ppm	Heavy metals (Pb).....≤5 ppm	

Code	Size	Packaging	Notes
411125	250 g	Plastic bottle	
411127	1 kg	Plastic bottle	
411121	25 kg	Plastic bucket	

L(+) Tartaric Acid > ERBApharm - According to pharmacopoeia: DAB-BP-FU-NF-Ph.Eur.-Ph.Franc.

ERBApharm

Description ..... White crystalline powder	Loss on drying ..... ≤0.2 %	Sulphate .....≤150 ppm	Assay (acidimetric) .....99.7 ÷ 100.5 % s.s.
Identification .....Positive	Sulphated ash..... ≤0.1 %	Heavy metals (Pb).....≤10 ppm	Residual solvents (Current ICH).....Conform
Appearance of solution .....Conform Ph.Eur.	Oxalic acid.....≤ 360 ppm	Calcium .....≤200 ppm	
Specific optical rotation... +12.0 ÷ +12.8 °	Chloride.....≤100 ppm	Origin (BSE/TSE).....Vegetable	

Code	Size	Packaging	Notes
307357	1 kg	Plastic bottle	
307359	5 kg	Plastic tank	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

L(+) Tartaric Acid > ERBApharm - Crystals - According to pharmacopoeia: Ph.Eur.-NF-FU-Ph.Franc.-BP-

ERBApharm

Description .....Colourless crystals	Specific optical rotation... +12.0 ÷ +12.8 °	Chloride.....≤ 100 ppm	Assay (acidimetric) .....99.7 ÷ 100.5 % s.s.
Identification .....Positive	Loss on drying ..... ≤ 0.2 %	Sulphate.....≤ 150 ppm	Origin (BSE/TSE).....Vegetable
Appearance of solution .....Conform Ph.Eur.	Sulphated ash..... ≤ 0.1 %	Heavy metals (Pb).....≤ 10 ppm	
Residual solvents (Current ICH).....Conform	Oxalic acid.....≤ 360 ppm	Calcium .....≤ 200 ppm	

Code	Size	Packaging	Notes
307307	1 kg	Plastic bottle	
307309	5 kg	Plastic tank	
307301	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

L(+) Tartaric acid solution 20% in water

• Acido L(+)tartarico soluzione 20% in acqua • Acide L(+)tartrique 20% • Acido L(+)tartárico solución 20% en agua • L (+) Weinsäure 20%

$HOOC(CHOH)_2COOH$
Molecular Weight: 150,09
CAS: 87-69-4



Danger
H318
P280i-P305+P351+P338-P310a

L(+) Tartaric acid solution 20% in water > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Density at 15° C ..... ~ 1.10 Assay (acidimetric) ..... 20 ÷ 21 %

Code	Size	Packaging	Notes
E411131	1 l	Bottle	

L-(+)-Tartaric acid diammonium salt ▶ Ammonium L(+)tartrate

Tartrazine

• Tartrazina • Tartrazine • Tartracina • Tartrazin

Synonym:
Acid Yellow 23

$C_{16}H_9N_4Na_3O_9S_2$
Molecular Weight: 534,39
CAS: 1934-21-0
EEC-N: 217-699-5

Tartrazine > RS - For microscopy - C.I. 19140

RS

Description ..... Orange powder Identification ..... Positive Water insoluble substances..... ≤ 0.2 % Assay (oxidimetric) ..... ≥85 %

Code	Size	Packaging	Notes
486903	50 g	Glass bottle	

Dye for histology

Tauber reagent

• Tauber reattivo • Réactif de Tauber • Tauber reactivo • Tauber-Reagenz

HEU210

Tauber reagent > RS - For microscopy

RS

Description ..... Yellow clear liquid Identification ..... Positive Sensibilità acido ascorbico ..... Conform

Code	Size	Packaging	Notes
490422	500 ml	Glass bottle	

TCA ▶ Trichloroacetic acid

Tellurium lumps

• Tellurio pezzi • Tellure en morceaux • Teluro trozos • Tellur in Stücken

Te
Molecular Weight: 127,61
CAS: 13494-80-9
EEC-N: 236-813-4



Danger
H332-H319-H360-H335-HA26
P261-P271-P280-P304+P340-P305+P351+P338-P403+P233

Tellurium lumps > RPE - For analysis

RPE

Description ..... Pezzi lucenti Identification ..... Positive Assay (oxidimetric) ..... 99 ÷ 100 %

Code	Size	Packaging	Notes
487002	25 g	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Tellurium powder
• Tellurio polvere • Tellure en poudre • Teluro polvo • Tellurpulver

Te
Molecular Weight: 127,6
CAS: 13494-80-9
EEC-N: 236-813-4

Danger
H332-H319-H360-H335-HA26
P261-P271-P280-P304+P340-P305+P351+P338-P403+P233

Tellurium powder > RPE - For analysis — RPE

Description ..... Black powder Identification ..... Positive Assay (oxidimetric) ..... ≥99 %

Code	Size	Packaging	Notes
487023	50 g	Glass bottle	

Tellurium standard solution
• Tellurio standard soluzione • Tellurium solution standard • Teluro, solución patrón • Tellur-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group II

Danger
H290-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Tellurium standard solution > RS - Standard solution for ICP-MS — RS

Code	Size	Packaging	Notes
505887	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505888	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Tellurium standard solution > RS - Standard solution for ICP — RS

Code	Size	Packaging	Notes
503981	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503983	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503985	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid
503987	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Tellurium standard solution > RS - Standard solution for AAS — RS

Code	Size	Packaging	Notes
507762	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507518	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Terbium standard solution

• Terbio standard soluzione • Terbium solution standard • Terbio, solución patrón • Terbium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Terbium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505882	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505885	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505883	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Terbium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503971	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503973	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
503975	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
503977	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Tetrabutylammonium bisulfate

• Tetrabutilammonio bisolfato • Tetrabutylammonium bisulfate • Tetrabutilamonio bisulfato • Tetrabutylammonium bisulfate

$C_{16}H_{37}NO_4S$
Molecular Weight: 339,53
CAS: 32503-27-8
EEC-N: 251-068-5

Warning
H302-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Tetrabutylammonium bisulfate > RS - For ion pair chromatography

RS

Absorbance UV curve (10%)
A210nm (1M) ..... ≤ 0.06 AU
A220nm (1M) ..... ≤ 0.05 AU
A230nm (1M) ..... ≤ 0.03 AU
A260nm (1M) ..... ≤ 0.02 AU
A500nm (1M) ..... ≤ 0.02 AU
Assay ..... ≥ 99.0 %

Code	Size	Packaging	Notes
405971	25 g	Plastic bottle	
405972	100 g	Plastic bottle	

Tetrabutylammonium bisulfate > RPE - For analysis

RPE

Description ..... White crystals
Identification ..... Positive
Assay (acidimetric) ..... ≥ 97.5 %
Melting point ..... 168 ÷ 172 °C
Water (K.F.) ..... ≤ 0.15 %

Code	Size	Packaging	Notes
487101	250 g	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Tetrabutylammonium bromide

• Tetrabutilammonio bromuro • Tétrabutylammonium bromure • Tetrabutilamonio bromuro • Tetrabutylammoniumbromid

$(C_4H_9)_4NBr$
Molecular Weight: 322,37
CAS: 1643-19-2
EEC-N: 216-699-2

Warning
H302
P264-P270-P301+P312a-P330-P501a

Tetrabutylammonium bromide > RS - For polarography

RS

Description ..... White crystals Identification ..... Positive Melting point..... 100÷104 °C

Code	Size	Packaging	Notes
487051	10 g	Glass bottle	

Tetrabutylammonium hydroxide solution 0.1 mol/l (0.1N)

• Tetrabutilammonio idrossido soluzione 0.1 mol/l (0.1N) • Tétrabutylammonium hydroxyde 0.1 mol/l (0.1N)
• Tetrabutilamonio hidróxido solución 0.1 mol/l (0.1N) • Tetrabutylammoniumhydroxid 0.1 mol/l (0.1N)

$(C_4H_9)_4NOH$
Molecular Weight: 259,48
CAS: 2052-49-5

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II

Danger
H225-H302-H315-H319-H361d-H371-H336-H373
P210-P241-P280-P303+P361+P353-P305+P351+P338-P308+P313

Tetrabutylammonium hydroxide solution 0.1 mol/l (0.1N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613008300	1 l	Glass bottle	Ref Ph.Eur 3008300

Tetrabutylammonium hydroxide 0.1 mol/l (0.1N) in 2-propanol

• Tetrabutilammonio idrossido 0.1 mol/l (0.1N) in isopropanolo • Tétrabutylammonium hydroxyde 0.1 mol/l (0.1N) dans propanol-2
• Tetrabutilamonio hidróxido 0.1 mol/l (0.1N) en 2-propanol • Tetrabutylammoniumhydroxid 0.1 mol/l (0.1N) in 2-Propanol

$(C_4H_9)_4NOH$
Molecular Weight: 259,48
CAS: 2052-49-5

Classification transport
ONU: 1992
Transport Hazard class: 3
Packing group II

Danger
H225-H301-H314-H370-H336
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338

Tetrabutylammonium hydroxide 0.1 mol/l (0.1N) in 2-propanol > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613008400	1 l	Glass bottle	Ref Ph.Eur 3008400

Tetrabutylammonium hydroxide 0.1 mol/l (0.1N) in 2-propanol > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive Density at 20° C ..... ~ 0.80 Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
E487031	500 ml	Glass bottle	

Tetrabutylammonium hydroxide 0.1 mol/l (0.1N) in methanol / propanol-2 (50/50)

- Tetrabutilammonio idrossido 0.1 mol/l (0.1N) in metanolo/propan-2-olo (50/50)
- Tetrabutylammonium hydroxyde 0.1 mol/l (0.1N) dans méthanol / propanol-2 (50/50)
- Tetrabutilamonio hidróxido 0.1 mol/l (0.1N) en methanol / propanol-2 (50/50)
- Tetrabutylammoniumhydroxid 0.1 mol/l (0.1 N) in Methanol / Propanol-2 (50/50)

$(C_4H_9)_4NOH$
Molecular Weight: 259,48
CAS: 2052-49-5

Classification transport
ONU: 1992
Transport Hazard class: 3
Packing group II

Danger
H225-H301-H315-H318-H370-H336
P210-P241-P280-P301+P310a-P303+P361+P353-P305+P351+P338

Tetrabutylammonium hydroxide 0.1 mol/l (0.1N) in methanol / propanol-2 (50/50) > RPE - For analysis — RPE

Assay (potentiometry) ..... 0.0998 - 0.1002 N

Code	Size	Packaging	Notes
P3840016	1 l	Glass bottle	

Tetrachloroauric(III) acid ▶ Gold(III) chloride trihydrate

Tetrachloro-1,4-benzoquinone ▶ Chloranil

Tetrachloroethylene

- Tetracloroetilene • Tétrachloroéthylène • Tetracloroetileno • Tetrachlorethylen

Synonym:
Perchloroethylene

C_2Cl_4
Molecular Weight: 165,83
CAS: 127-18-4
EEC-N: 204-825-9

Classification transport
ONU: 1897
Transport Hazard class: 6.1
Packing group III

Warning
H315-H319-H317-H351-H336-H411
P261-P271-P280-P304+P340-P305+P351+P338-P403+P233

Tetrachloroethylene > RS - SPECTROSOL - For optical spectroscopy — RS

Clear,colourless liq.appearance ..... Conform
Identification ..... Conform
Refractive index at 20°C .. 1.5014 - 1.5074
Density at 20°C ..... 1.618 - 1.628
Boiling point ..... 120.3 - 121.8 °C
Water content (K.F.) ..... ≤ 100 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
Assay (GC) ..... ≥ 99.9 %
Free acid ..... ≤ 0.0003 meq/g
Free alkali ..... ≤ 0.00006 meq/g
Fluorescence (quinine) at 365 nm .. ≤ 2 ppb
UV transmittance at 300 nm ..... ≥ 75 %
UV transmittance at 310 nm ..... ≥ 83 %
UV transmittance at 320 nm ..... ≥ 88 %
UV transmittance at 340 nm ..... ≥ 88 %
UV transmittance at 360 nm ..... ≥ 88 %
UV transmittance at 400 nm ..... ≥ 97 %

Code	Size	Packaging	Notes
P0682716	1 l	Glass bottle	
P0682721	2.5 l	Glass bottle	

Tetrachloroethylene > RPE - For analysis - Stabilized — RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Alcohol miscibility ..... Complete
Benzene miscibility ..... Complete
Diethyl ether miscib. ..... Complete
Phosgene ..... Conform
Ready carbonizable substances ..... Conform
Density at 20° C ..... 1.618 ÷ 1.628
Refractive index at 20°C . 1.5014 ÷ 1.5074
Boiling point ..... 120.3 ÷ 121.8 °C
Water (K.F.) ..... ≤ 200 ppm
Residue on evaporation ..... ≤ 10 ppm
Acidity (HCl) ..... ≤ 5 ppm
Alcalinity (NH3) ..... ≤ 0.5 ppm
Free chlorine ..... ≤ 0.1 ppm
Chloride ..... ≤ 1 ppm
Assay (GLC) ..... ≥ 99.5 %

Code	Size	Packaging	Notes
449671	1 l	Glass bottle	
449672	2.5 l	Glass bottle	
449673	35 kg	Drum	

Tetrachloroethylene > RE - Pure - Stabilized — RE

Refractive index at 20°C ..... 1.503 - 1.507
Description ..... Clear colourless liquid
Water content (K.F.) ..... ≤ 100 mg/Kg
Identification ..... Positive
Non volatile residue ..... ≤ 50 mg/Kg
Density at 20° C ..... 1.618 ÷ 1.628
Colour ..... ≤ 10 Hazen
Residue on evaporation ..... ≤20 ppm
Assay (GC) ..... ≥ 99 %
Assay (GLC) ..... ≥98 %

Code	Size	Packaging	Notes
343001	1 l	Glass bottle	
P0680228	5 l	Plastic tank	
343003	40 kg	Metal drum	

Product specifications are subject to changes. Please visit our website for updates.

Tetrachloroethane-d2

• Tetracloroetano-d2 • Tétrachloroéthane-d2 • Tetracloroetano-d2 • Tetrachlorethan-d2

Synonym:
• *1,1,2,2-Tetrachloroethane-d2*
• *1,2-Dideutero-1,1,2,2-tetrachloroethane*

$C_2D_2Cl_4$
Molecular Weight: 169,86
CAS: 33685-54-0
EEC-N: 251-634-1

Classification transport
ONU: 1702
Transport Hazard class: 6.1
Packing group II

Danger
H310-H330-H411
P284-P304+P340-P310a-P320-P361+P364-P403+P233

Tetrachloroethane-d2 > RS - For NMR - min 99.5%

RS

Code	Size	Packaging	Notes
P5435	25 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Tetraethylammonium bromide

• Tetraetilammonio bromuro • Tétraéthylammonium bromure • Tetraetilamonio bromuro
• Tetraethylammoniumbromid

Synonym:
TEA Bromide

$(C_2H_5)_4NBr$
Molecular Weight: 210,17
CAS: 71-91-0
EEC-N: 200-769-4

Tetraethylammonium bromide > RS - For polarography

RS

Description ..... White powder Identification ..... Positive

Code	Size	Packaging	Notes
487152	25 g	Glass bottle	

Tetrahydrofuran

• Tetraidrofurano • Tétrahydrofuranne • Tetrahidrofurano • Tetrahydrofuran

$OCH_2CH_2CH_2CH_2$
Molecular Weight: 72,11
CAS: 109-99-9
EEC-N: 203-726-8

Classification transport
ONU: 2056
Transport Hazard class: 3
Packing group II

Danger
H225-H302-H319-H351-H335-HEU019
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Tetrahydrofuran > RS - For HPLC - Isocratic grade - Not stabilized

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.885 ÷ 0.893
Refractive index at 20°C ..... 1.4040 ÷ 1.4100
Boiling point ..... 64.0 ÷ 65.0 ° C
Water (K.F.) ..... ≤200 ppm
Residue on evaporation ..... ≤2 ppm
Peroxide ..... ≤200 ppm
Acidity ..... ≤0.00035 meq/g
Alcalinity ..... ≤0.0006 meq/g
Assay (GLC) ..... ≥99.9 %
U.V. Transmittance
at 230 nm ..... ≥15 %
at 240 nm ..... ≥40 %
at 250 nm ..... ≥55 %
at 260 nm ..... ≥70 %
at 270 nm ..... ≥85 %
At 275 nm ..... ≥ 90 %
at 280 nm ..... ≥92 %
At 310 nm ..... ≥ 98 %
At 315 nm ..... ≥ 99 %

Code	Size	Packaging	Notes
412451000	1 l	Glass bottle	
412453000	1 l	Glass bottle PVC coated	
412452000	2.5 l	Glass bottle	

Tetrahydrofuran > RS - For HPLC - Isocratic grade - Stabilized with BHT

RS

Clear,colourless liq.appearance ..... Conform
Identification ..... Conform
Colour ..... ≤ 10 Apha
Refractive index at 20°C ..... 1.405 - 1.409
Water content (K.F.) ..... ≤ 200 mg/Kg
Peroxides (as H2O2) ..... ≤ 50 mg/Kg
Stabilizer (ionol) ..... 40 - 60 mg/Kg
UV transmittance at 240 nm ..... ≥ 10 %
UV transmittance at 250 nm ..... ≥ 40 %
UV transmittance at 280 nm ..... ≥ 30 %
UV transmittance at 300 nm ..... ≥ 90 %
UV transmittance at 320 nm ..... ≥ 95 %
Free acid (as CH3COOH) ..... ≤ 20 mg/Kg
Assay (GC) ..... ≥ 99.8 %
Non volatile residue (without stab.) ≤ 5 mg/Kg

Code	Size	Packaging	Notes
412471	1 l	Glass bottle	
412472	2.5 l	Glass bottle	

Tetrahydrofuran > RS - For preparative HPLC - Stabilized with BHT

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.885 ÷ 0.893
Refractive index at 20°C. 1.4040 ÷ 1.4100
Boiling point..... 64.0 ÷ 65.0 ° C
Water (K.F.) ..... ≤200 ppm
Residue on evaporation ..... ≤5 ppm
Alcalinity..... ≤0.0002 meq/g
Peroxide ..... ≤50 ppm
Assay (GLC) ..... ≥99.8 %
U.V. Transmittance
at 320 nm ..... ≥90 %
Stabilizer (Ionol)..... 40 ÷ 60 ppm

Code	Size	Packaging	Notes
487352	2.5 l	Glass bottle	

Tetrahydrofuran > RS - SPECTROSOL - For optical spectroscopy - Not stabilized - Reag.Ph.Eur.

RS

Description ..... Clear colourless liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Density at 20° C ..... 0.885 ÷ 0.893
Refractive index at 20°C. 1.4040 ÷ 1.4100
Boiling point..... 64.0 ÷ 65.0 ° C
Water (K.F.) ..... ≤200 ppm
Residue on evaporation ..... ≤5 ppm
Acidity ..... ≤0.0005 meq/g
Alcalinity..... ≤0.0002 meq/g
Peroxides (H2O2) ..... ≤300 ppm
Assay (GLC) ..... ≥99.9 %
U.V. Transmittance
at 240 nm ..... ≥30 %
at 250 nm ..... ≥50 %
at 300 nm ..... ≥95 %
at 320 nm ..... ≥98 %
UV Absorbance at 255 nm ..... ≤ 0.70 AU
UV Absorbance at 270 nm ..... ≤ 0.10 AU
UV Absorbance at 310 nm ..... ≤ 0.01 AU

Code	Size	Packaging	Notes
487345	1 l	Glass bottle	
487346	2.5 l	Glass bottle	

Tetrahydrofuran > RS - Anhydrous - For analysis - Stabilized with BHT

RS

Refractive index at 20°C..... 1.405 - 1.409
Water content (K.F.) ..... ≤ 100 mg/Kg
Colour ..... ≤ 10 Hazen
Peroxides (as H2O2)..... ≤ 20 mg/Kg
Stabilizer (ionol)..... 200 - 400 mg/Kg
Assay (GC) ..... ≥ 99.9 %
Free acid (as CH3COOH)..... ≤ 20 mg/Kg
Clear,colourless liq.appearance ..... Conform
Identification (IR)..... Conform
Density d20/4 ..... 0.884 - 0.894
Non volatile residue (without stab.)≤ 10 mg/Kg

Code	Size	Packaging	Notes
P0701010	200 ml	Bottle with septum	
P07010T10	200 ml	Bottle with septum	On molecular sieves 4A
P0701016	1 l	Glass bottle	
P07010T16	1 l	Glass bottle	On molecular sieves 4A
P0701021	2.5 l	Glass bottle	

Tetrahydrofuran > RPE - For analysis - ACS - Reag.Ph.Eur. - Reag.USP - Stabilized with BHT

RPE

Description ..... Clear colourless liquid
Colour (APHA) ..... ≤ 10
Identification (I.R.)..... Conform
Density at 20° C ..... 0.885 ÷ 0.893
Refractive index at 20°C. 1.4040 ÷ 1.4100
Boiling point..... 64.0 ÷ 65.0 ° C
Water (K.F.) ..... ≤150 ppm
Residue on evaporation ..... ≤300 ppm
Acidity (acetic acid)..... ≤20 ppm
Alcalinity (NH3)..... ≤7 ppm
Peroxides (H2O2) ..... ≤20 ppm
Al ..... ≤0.5 ppm
Ba ..... ≤0.1 ppm
Ca ..... ≤0.5 ppm
Cd ..... ≤0.05 ppm
Co ..... ≤0.02 ppm
Cr..... ≤0.02 ppm
Cu..... ≤0.02 ppm
Fe ..... ≤0.1 ppm
Mg ..... ≤0.1 ppm
Mn ..... ≤0.02 ppm
Ni..... ≤0.02 ppm
Pb ..... ≤0.1 ppm
Sn ..... ≤0.1 ppm
Zn ..... ≤0.1 ppm
Assay (GLC) ..... ≥99.9 %
Stabilized with BHT..... 200 ÷ 350 ppm
Residue on evaporation (without stab). ≤ 10 ppm

Code	Size	Packaging	Notes
487308	1 l	Glass bottle	
487303	2.5 l	Glass bottle	
487305	5 l	Plastic tank	
487307	5 l	Aluminium can	
487301	23 kg	Metal drum	
487309	200 l	Metal drum	

Product specifications are subject to changes. Please visit our website for updates.

Tetrahydrofuran > RE - Pure - Stabilized with BHT

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.884 ÷ 0.894
Refractive index at 20°C. 1.4020 ÷ 1.4120
Boiling point..... 63.8 ÷ 65.3 ° C
Water (K.F.) ..... ≤200 ppm
Acidity ..... ≤50 ppm
Peroxides (H2O2) ..... ≤100 ppm
Assay (GLC) ..... ≥99.5 %
Residue on evaporation (without stab). ≤ 50 ppm
Stabilized with BHT..... 200 ÷ 350 ppm

Code	Size	Packaging	Notes
382981	1 l	Glass bottle	
382985	2.5 l	Glass bottle	
382982	5 l	Aluminium can	
528481	5 l	Plastic tank	
382986	23 kg	Metal drum	
382983	200 l	Metal drum	

Tetrahydrofuran-d8

• Tetraidrofurano-d8 • Tétrahydrofuranne-d8 • Tetrahidrofurano-d8 • Tetrahydrofuran-d8

Synonym:
Octadeuterotetrahydrofuran

C_4D_8O
Molecular Weight: 80,16
CAS: 1693-74-9
EEC-N: 216-898-4

Classification transport
ONU: 2056
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H319-H351-H335-HEU019
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Tetrahydrofuran-d8 > RS - For NMR - min 99.5%

RS

Code	Size	Packaging	Notes
P5380	2 x 0.6 ml	Glass ampoule	
P5385	25 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Tetramethylammonium hydroxide 10%

• Tetrametilammonio idrossido 10% • Tétraméthylammonium hydroxyde 10% • Tetrametilamonio hidróxido 10% • Tetramethylammoniumhydroxid 10%

$C_4H_{13}NO$
Molecular Weight: 91,15
CAS: 75-59-2

Classification transport
ONU: 1835
Transport Hazard class: 8
Packing group II



Danger
H301-H311-H314-H370-H372
P280-P301+P310a-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338

Tetramethylammonium hydroxide 10% > RS - For polarography and steroids detection

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Assay (acidimetric) ..... 9.5 ÷ 10.5 %

Code	Size	Packaging	Notes
487491	100 ml	Glass bottle	
487492	250 ml	Glass bottle	

N,N,N',N'-Tetramethyl-p-phenylenediamine dihydrochloride

• N,N,N',N'-Tetrametil-p-fenilendiammina bicloridrato
• N,N,N',N'-Tétraméthyl-p-phénylènediamine dichlorhydraté
• N,N,N',N'-Tetrametil -p-fenilendiamina diclorhidrato
• N,N,N',N'-Tetramethyl-p-phenylendiamin-Dihydrochlorid

Synonym:
• *TMPPD*
• *Wurster's reagent*

$C_6H_4[N(CH_3)_2]_2.2HCl$
Molecular Weight: 237,17
CAS: 637-01-4
EEC-N: 211-274-8



Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

N,N,N',N'-Tetramethyl-p-phenylenediamine dihydrochloride > RPE - For analysis

RPE

Description . White-hazel crystalline powder
Identification ..... Positive
Melting point..... 219 ÷ 222 °C
Assay (non-aqueos medium)..... ≥ 98.5 %

Code	Size	Packaging	Notes
487601	5 g	Glass bottle	

TFA ▶ Trifluoroacetic acid

TFAA ▶ Trifluoroacetic anhydride

Thallium standard solution

• Tallio standard soluzione • Thallium solution standard • Talio, solución patrón • Thallium-Standardlösung

Thallium standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615003000	100 ml	Plastic bottle	A 10 ppm solution Ref Ph.Eur 5003000

Thallium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505912	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505915	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505913	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Thallium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
504011	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504013	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504015	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
504017	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Thiazole Yellow G ▶ Clayton's yellow

Thioacetamide

• Tioacetammide • Thioacétamide • Tioacetamida • Thioacetamid

Synonym:
Ethanethioamide

C_2H_5NS
Molecular Weight: 75,13
CAS: 62-55-5



Thioacetamide > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
617000211	100 ml	Plastic bottle	Thioacetamide TS

Thioacetamide > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White crystalline powder	Melting point ..... 111 ÷ 114 ° C	Assay (argentimetric) ..... ≥ 99.0 %
Identification ..... Positive	Residue on ignition ..... ≤ 500 ppm	Appearance of solution 2% ..... Conform

Code	Size	Packaging	Notes
487803	50 g	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Thioacetamide solution 40 g/l

• Tioacetammide soluzione 40 g/l • Thioacétamide solution 40 g/l • Tioacetamida solución 40 g/l
• Thioacetamidlösung 40 g/l

Synonym:
Ethanethioamide

CH_3CSNH_2
Molecular Weight: 75,13
CAS: 62-55-5

Danger
H350-HA26
P201-P202-P280-P308+P313-P405-P501a

Thioacetamide solution 40 g/l > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611089603	100 ml	Plastic bottle	Ref Ph.Eur 1089602
611089602	1 l	Plastic bottle	Ref Ph.Eur 1089602

2-Thiobarbituric acid

• Acido 2-tiobarbiturico • Acide 2-thiobarbiturique • Acido 2-tiobarbitúrico • 2-Thiobarbitursäure

Synonym:
• *4,6-Dihydroxy-2-mercaptopyrimidine*
• *4,6-Dihydroxypyrimidine-2-thiol*

$NHCOCH_2CONHCS$
Molecular Weight: 144,15
CAS: 504-17-6
EEC-N: 207-985-8

2-Thiobarbituric acid > RPE - For analysis

RPE

Description ..... Yellowish crystalline powder
Identification ..... Positive
Loss on drying ..... ≤ 2 %
Assay (acidimetric) ..... ≥ 97.5 % (s.s.)

Code	Size	Packaging	Notes
411271	5 g	Glass bottle	
411272	25 g	Glass bottle	

Thiocarbamide ▶ Thiourea

Thioethylene glycol ▶ 2-Mercaptoethanol

Thioglycolic acid 80%

• Acido tioglicolico 80% • Acide thioglycolique 80% • Acido tioglicólico 80% • Thioglykolsäure 80%

Synonym:
Mercaptoacetic acid

$C_2H_4O_2S$
Molecular Weight: 92,12
CAS: 68-11-1

Classification transport
ONU: 1940
Transport Hazard class: 8
Packing group II

Danger
H301-H311-H331-H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P311a-P305+P351+P338-P361+P364-P403+P233

Thioglycolic acid 80% > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611089700	10 ml	Glass bottle	Ref Ph.Eur 1089700

Thioglycolic acid 80% > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... ≥ 1.260
Heavy metals (Pb) ..... ≤10 ppm
Iron sensitivity ..... ≥0.1 µg/ml
Sulphate ..... ≤50 ppm
Fe ..... ≤50 ppm
Assay (oxidimetric) ..... ≥ 78 %

Code	Size	Packaging	Notes
411385	500 ml	Glass bottle	

Thiourea

• Tiourea • Thiourée • Tiourea • Thioharnstoff

Synonym:
• *Thiocarbamide*
• *Sulfourea*

NH_2CSNH_2
Molecular Weight: 76,12
CAS: 62-56-6
EEC-N: 200-543-5

Warning
H302-H351-H361d-H411
P264-P280-P301+P312a-P330-P308+P313-P501a

Thiourea > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... White crystalline powder
Identification ..... Positive
Water solubility..... Conform
Melting point..... 174 ÷ 177 °C
Loss on drying ..... ≤0.5 %
Residue on ignition..... ≤0.1 %
Assay (argentimetric)..... ≥99.0 % s.s.

Code	Size	Packaging	Notes
488104	100 g	Plastic bottle	
488105	250 g	Plastic bottle	
488107	1 kg	Plastic bottle	
488102	25 kg	Plastic bucket	
488101	50 kg	Fibre drum	

Thiourea > RE - Pure

RE

Description ..... White crystalline powder
Identification ..... Positive
Melting point..... 173 ÷ 178 ° C
Loss on drying ..... ≤2 %
Residue on ignition..... ≤0.1 %
Fe ..... ≤10 ppm
Assay (argentimetric)..... ≥98 %

Code	Size	Packaging	Notes
385407	1 kg	Plastic bottle	
385409	5 kg	Plastic tank	
385403	25 kg	Plastic bucket	

Thorium standard solution

• Torio standard soluzione • Thorium solution standard • Torio, solución patrón • Thorium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Thorium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
504281	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504283	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Thulium standard solution

• Tulio standard soluzione • Thullium solution standard • Tulio, solución patrón • Thullium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Thulium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505917	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505918	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Product specifications are subject to changes. Please visit our website for updates.

Thulium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507763	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507519	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Thymol

• Timolo • Thymol • Timol • Thymol

Synonym:
2-Isopropyl-5-methylphenol

$C_{10}H_{14}O$
Molecular Weight: 150,22
CAS: 89-83-8
EEC-N: 201-944-8

Classification transport
ONU: 2430
Transport Hazard class: 8
Packing group II



Danger
H302-H314-H411
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Thymol > ERBApharm - According to pharmacopoeia: BP-DAB-NF-Ph.Eur.-FU

ERBApharm

Description ..... Colourless crystals
Identification ..... Positive
Appearance of solution ..... Conform Ph.Eur.
Acidity ..... Conform Ph.Eur.
Related compounds ..... Conform Ph.Eur.
Organic volatile impurities ..... Conform NF
Melting point ..... 48 ÷ 51 °C
Not volatile residue ..... ≤500 ppm
Assay ..... 99.0 ÷ 101.0 %

Code	Size	Packaging	Notes
384205	250 g	Plastic bottle	
384201	1 kg	Plastic bottle	
384202	2.5 kg	Plastic bottle	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Thymol blue

• Blu timolo • Bleu de thymol • Azul de timol • Thymolblau

Synonym:
Thymolsulfonphthalein

$C_{27}H_{30}O_5S$
Molecular Weight: 466,59
CAS: 76-61-9
EEC-N: 200-973-3

Thymol blue > RPE - For analysis

RPE

Description ..... Polvere verde bruna
Identification ..... Positive
pH range ..... 1.2 ÷ 2.8
Colour change ..... rosso-giallo

Code	Size	Packaging	Notes
429228	5 g	Glass bottle	
429222	25 g	Glass bottle	
429223	50 g	Plastic bottle	

Thymol blue 0.4% in ethanol

• Blu timolo soluzione 0,4% in alcole etilico • Bleu de thymol solution 0.4% dans l'éthanol
• Azul de timol solucion 0.4% en alcohol etilico • Thymolblau lösung 0.4% in Ethanol

Synonym:
Thymolsulfonphthalein

$C_{27}H_{30}O_5S$
Molecular Weight: 466,59
CAS: 76-61-9

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group III



Warning
H226-H319
P210-P241-P280-P303+P361+P353-P305+P351+P338-P337+P313

Thymol blue 0.4% in ethanol > RPE - For analysis

RPE

Description ..... Blue-brown liquid
Identification ..... Positive
Sensitivity(pH 1.2-2.8) ..... Conform
Sensitivity(pH 8.2-9.6) ..... Conform
Colour change ..... red-yellow

Code	Size	Packaging	Notes
E429235	250 ml	Glass bottle	

Acid-basis indicator

Thymol blue indicator
• Blu timolo indicatore • Indicateur bleu de thymol • Indicador azul de timol • Thymol blauer Indikator

Thymol blue indicator > RS - For analysis — RS

pH..... 6.8 - 7 unite pH Temperature of measurement ..... 15 - 25 °C

Code	Size	Packaging	Notes
PS0270/15	1 l	Plastic bottle	

Thymol blue solution
• Blu timolo soluzione • Bleu de thymol solution • Azul de timol solución • Thymolblau lösung

Classification transport
ONU: 2924
Transport Hazard class: 3
Packing group III



Warning
H226-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

Thymol blue solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1 — RS

Code	Size	Packaging	Notes
611090601	100 ml	Plastic bottle	Ref Ph.Eur 1090601

Colour change: pH 1.2 (red) to pH 2.8 (yellow); pH 8.0 (olive-green) to pH 9.6 (blue)

Thymol blue TA indicator
• Blu timolo TA indicatore • Indicateur TA au bleu de thymol • TA indicador azul de timol • Thymolblau TA-Indikator

Thymol blue TA indicator > RS - For analysis — RS

pH..... 6.5 - 7.5 unite pH Temperature of measurement ..... 15 - 25 °C Coca-cola colour ..... Conform

Code	Size	Packaging	Notes
PS0187/15	1 l	Plastic bottle	
PS0187/16	1 l	Glass bottle	

Thymolphthalein
• Timolftaleina • Thymolphtaléine • Timolftaleína • Thymolphthalein

Synonym:
5',5''-Diisopropyl-2',2''-dimethylphenolphthalein

$C_{28}H_{30}O_4$
Molecular Weight: 430,55
CAS: 125-20-2
EEC-N: 204-729-7

Thymolphthalein > RPE - For analysis - ACS — RPE

Description ..... White crystalline powder Identification ..... Positive Colour change..... incolore ÷ blu pH range ..... 8.8 ÷ 10.5

Code	Size	Packaging	Notes
487728	5 g	Glass bottle	
487729	25 g	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Thymolphthalein 0.1% hydroalcoholic solution

• Timolftaleina 0.1% soluzione idroalcolica • Thymolphtaléine 0.1% solution hydroalcoolique
• Timolftaleína 0.1% solución hidroalcohólica • Thymolphthalein 0.1% ige hydroalkoholische Lösung

Synonym:
5',5''-Diisopropyl-2',2''-dimethylphenolphthalein

$C_{28}H_{30}O_4$
Molecular Weight: 430,55
CAS: 125-20-2

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II

Danger
H225-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

Thymolphthalein 0.1% hydroalcoholic solution > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.877 ÷ 0.883
Colour change..... incolore blu
pH range ..... 9.3 - 10.5

Code	Size	Packaging	Notes
E487755	250 ml	Glass bottle	

Michaelis indicator series

Thymolphthalein solution 0.1% in ethanol

• Timolftaleina soluzione 0.1% in etanolo • Thymolphtaléine solution 0.1% dans l'éthanol
• Timolftaleína solución 0.1% en etanol • Thymolphthaleinlösung 0.1% in Ethanol

Synonym:
5',5''-Diisopropyl-2',2''-dimethylphenolphthalein

$C_{28}H_{30}O_4$
Molecular Weight: 430,55
CAS: 125-20-2

Classification transport
ONU: 1170
Transport Hazard class: 3
Packing group II

Danger
H225-H319
P210-P241-P280-P303+P361+P353-
P305+P351+P338-P337+P313

Thymolphthalein solution 0.1% in ethanol > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611090701	100 ml	Plastic bottle	Ref Ph.Eur 1090701

Colour change: pH 9.3 (colourless) to pH 10.5 (blue)

Tin, powder

• Stagno, polvere • Etain, poudre • Estaño, polvo • Zinnpulver

Sn
Molecular Weight: 118,69
CAS: 7440-31-5
EEC-N: 231-141-8

Tin, powder > RPE - For analysis

RPE

Description ..... Grey powder
Identification ..... Positive
Assay (gravimetric) ..... ≥99 %
Particle size >106 µm ..... ≥ 0 %
Particle size >75 µm ..... ≤ 0.5 %
Particle size >45 µm ..... 3 - 15 %

Code	Size	Packaging	Notes
484914	100 g	Glass bottle	
484917	1 kg	Plastic bottle	

Tin foil

• Stagno lastra • Etain en feuilles • Estaño hojas • Zinn geht

Sn
Molecular Weight: 118,69
CAS: 7440-31-5
EEC-N: 231-141-8

Tin foil > RPE - For analysis

RPE

Description ..... Metallic foil
Identification ..... Positive
Assay (gravimetric) ..... ≥99 %

Code	Size	Packaging	Notes
484887	1 kg	Box	

Tin standard solution

• Stagno standard soluzione • Etain solution standard • Estaño, solución patrón • Zinn-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Warning
H290
P234-P390-P406

Tin standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615003101	100 ml	Plastic bottle	A 0.1 ppm solution: to dilute according to Ref Ph.Eur 5003101
615003109	100 ml	Plastic bottle	A 5 ppm solution: to dilute according to Ref Ph.Eur 5003100

Tin standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505862	100 ml	Plastic bottle	conc. 10 ppm Matrix: Nitric acid and hydrofluoric acid
505865	100 ml	Plastic bottle	conc. 100 ppm Matrix: Hydrofluoric acid and nitric acid
505863	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Tin standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
503941	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503943	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
503945	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid
503947	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Tin standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
503949	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrofluoric acid and nitric acid
507492	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid
E497655	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
E497651	500 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Tin standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
484861		Plastic ampoule	conc. 1.000 ppm Matrix: Hydrochloric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Product specifications are subject to changes. Please visit our website for updates.

Tin (II) chloride dihydrate

• Stagno cloruro oso diidrato • Chlorure d'étain (II) dihydrate • Estaño (II) cloruro dihidrato
• Zinn (II) chloriddihydrat

Synonym:
Stannous chloride dihydrate

$SnCl_2.2H_2O$
Molecular Weight: 225,63
CAS: 10025-69-1
EEC-N: 231-868-0

Classification transport
ONU: 3260
Transport Hazard class: 8
Packing group III



Danger
H314-H335-H373
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338-P403+P233

Tin (II) chloride dihydrate > RPE - For analysis - ACS

RPE

Description .....White crystals
Identification .....Positive
HCl solubility .....Conform
Sulphate .....Conform
Ca .....≤ 50 ppm
Fe .....≤ 30 ppm
K .....≤ 50 ppm
Na .....≤ 100 ppm
Pb .....≤ 100 ppm
Assay (oxidimetric) .....98.0 ÷ 103.0 %

Code	Size	Packaging	Notes
485004	100 g	Plastic bottle	
485005	250 g	Plastic bottle	
485007	1 kg	Plastic bottle	
485002	5 kg	Plastic jar	

Tin (II) chloride dihydrate > RE - Pure

RE

Description .....White crystals
Identification .....Positive
Assay (oxidimetric) .....≥97.0 %

Code	Size	Packaging	Notes
379406	500 g	Plastic bottle	
379407	5 kg	Plastic tank	
379403	25 kg	Plastic bucket	

Tin (II) chloride solution

• Stagno cloruro oso soluzione • Etain (II) chlorure solution • Estaño (II) cloruro solución
• Zinn (II) chloridlösung

Synonym:
Stannous chloride

$SnCl_2.2H_2O$
Molecular Weight: 189,62
CAS: 10025-69-1



Danger
H314-H373
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Tin (II) chloride solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611085001	100 ml	Glass bottle	Ref Ph.Eur 1085001

Tin (II) chloride solution 10%

• Stagno cloruro oso soluzione 10% • Chlorure d'étain (II) dihydrate solution 10%
• Estaño (II) cloruro solución 10% • Zinn (II) -chloridlösung 10%

Synonym:
Stannous chloride

$SnCl_2.2H_2O$
Molecular Weight: 189,62
CAS: 10025-69-1



Danger
H314-H373
P280-P301+P330+P331-P303+P361+P353-
P304+P340-P310a-P305+P351+P338

Tin (II) chloride solution 10% > RPE - For analysis

RPE

Description .....Clear colourless liquid
Identification .....Positive
Density at 20° C .....1.075 ÷ 1.085
Assay .....9.5 ÷ 10.5 %

Code	Size	Packaging	Notes
E485041	1 l	Bottle	

Tin (II) sulfate

• Stagno solfato oso • Etain (II) sulfate • Estaño (II) sulfato • Zinn (II) sulfat

Synonym:
Stannous sulfate

$SnSO_4$
Molecular Weight: 214,75
CAS: 7488-55-3
EEC-N: 231-302-2

Tin (II) sulfate > RE - Pure

RE

Description ..... Yellowish crystalline powder
Identification ..... Positive
Co ..... ≤50 ppm
Cu ..... ≤20 ppm
Fe ..... ≤100 ppm
Ni ..... ≤20 ppm
Pb ..... ≤200 ppm
Assay (oxidimetric) ..... ≥95 %

Code	Size	Packaging	Notes
379601	1 kg	Plastic bottle	

Tin (IV) chloride pentahydrate

• Stagno cloruro ico pentaidrato • Etain (IV) chlorure pentahydrate • Estaño (IV) cloruro pentahidratado
• Zinn (IV) chloridpentahydrat

Synonym:
Tin tetrachloride

$SnCl_4.5H_2O$
Molecular Weight: 350,58
CAS: 10026-06-9
EEC-N: 231-588-8

Classification transport
ONU: 2440
Transport Hazard class: 8
Packing group III

Danger
H314-H412
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Tin (IV) chloride pentahydrate > RPE - For analysis

RPE

Description .... White to off-white crystalline powder or lumps
Identification ..... Positive
Fe ..... ≤ 50 ppm
Assay (argentimetric) ..... ≥ 98.0 %
Pb ..... ≤ 50 ppm

Code	Size	Packaging	Notes
485074	100 g	Glass bottle	
485076	500 g	Plastic bottle	

Tin (IV) oxide

• Stagno ossido ico • Etain (IV) oxyde • Estaño (IV) óxido • Zinn (IV) oxid

SnO_2
Molecular Weight: 150,69
CAS: 18282-10-5
EEC-N: 242-159-0

Tin (IV) oxide > RPE - For analysis

RPE

Description ..... White to light grey powder
Identification ..... Positive
Assay ..... ≥ 99.9 %

Code	Size	Packaging	Notes
485154	100 g	Glass bottle	

Tisab

• Tisab • Tisab • Tisab • TISAB

Danger
H318
P280i-P305+P351+P338-P310a

Tisab > RS - For analysis

RS

pH ..... 4.7 - 4.9 unite pH
Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
PS0363/21	2.5 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Tisab > RS - For fluorides analysis

RS

pH..... 5.00 - 6.00 unité pH Temperature of measurement ..... 19 - 21 °C

Code	Size	Packaging	Notes
PS0809/22	5 l	Plastic tank	

Tisab III solution

• Tisab III soluzione • Tisab III solution • Tisab III solución • Tisab III-Lösung

Danger
H318
P280i-P305+P351+P338-P310a

Tisab III solution > RS - For fluorides analysis

RS

Description ..... Clear colourless liquid Identification ..... Positive pH at 20° C ..... 5.0 ÷ 5.5

Code	Size	Packaging	Notes
488162	500 ml	Plastic bottle	

Titanium standard solution

• Titanio standard soluzione • Titane solution standard • Titanio, solución patrón • Titan-Standardlösung

Titanium standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615003200	1 l	Plastic bottle	A 100 ppm solution Ref Ph.Eur 5003200

Titanium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505907	100 ml	Plastic bottle	conc. 10 ppm Matrix: Nitric acid and hydrofluoric acid
505908	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid and hydrofluoric acid
505909	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrofluoric acid and nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Titanium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
504001	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid and hydrofluoric acid
504003	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid and hydrofluoric acid
504005	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid and hydrofluoric acid
504007	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid and hydrofluoric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Titanium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507764	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid
507520	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Titanium dioxide

• Titanio biossido • Titane dioxyde • Titanio dióxido • Titandioxid

Synonym:
- *Titanium(IV) oxide*
- *Titania*

TiO_2
Molecular Weight: 79,9
CAS: 13463-67-7
EEC-N: 236-675-5

Titanium dioxide > RPE - For analysis

RPE

Description	White powder	Phosphate	≤0.1 %	Water solubility	≤0.4 %	Zn	≤50 ppm
Identification	Positive	H2SO4-insoluble matter	≤0.1 %	As	≤2 ppm	Assay (oxidimetric)	≥98.5 %
Loss on drying	≤0.5 %	Heavy metals (Pb)	≤10 ppm	Cu	≤5 ppm		
Loss on ignition	≤1.0 %	Sulphate	≤0.1 %	Fe	≤50 ppm		
Chloride	≤200 ppm	Solubility in dil. HCl	≤0.5 %	Pb	≤10 ppm		

Code	Size	Packaging	Notes
488256	100 g	Plastic bottle	
488257	1 kg	Plastic bottle	
488251	10 kg	Carton box	

Titanium dioxide > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-BP

ERBApharm

Description	White powder	Barium	Conform Ph.Eur.	Loss on ignition	≤0.5 %	Fe	≤200 ppm
Identification	Positive	Water-soluble subst.	≤0.25 %	Heavy metals (Pb)	≤20 ppm	Assay (oxidimetric)	99.0 ÷ 100.5 %s.s.
Appearance of solution	Conform Ph.Eur.	Acid soluble matter	≤0.5 %	Sb	≤100 ppm		
Acidity or alkalinity	Conform Ph.Eur.	Loss on drying	≤0.5 %	As	≤1 ppm		

Code	Size	Packaging	Notes
385751	1 kg	Plastic bottle	
385752	5 kg	Plastic tank	
385753	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Titanium dioxide > RE - Pure

RE

Description	Greyish powder	Loss on drying	≤1 %	Assay (oxidimetric)	≥98 %
Identification	Positive	Solubility in dil. HCl	≤1 %		

Code	Size	Packaging	Notes
385707	1 kg	Plastic bottle	
385709	5 kg	Plastic tank	
385702	25 kg	Plastic bucket	

Titanium isopropylate

• Titanio isopropilato • Titane isopropylate • Titanio tetrapropilato • Titanisopropylat

Synonym:
- *Titanium(IV) isopropoxide*
- *Tetraisopropyl orthotitanate*

$Ti[OCH(CH_3)_2]_4$
Molecular Weight: 284,26
CAS: 546-68-9
EEC-N: 208-909-6

Classification transport
ONU: 2413
Transport Hazard class: 3
Packing group III



Danger
H226-H331-H319-H336
P210-P280-P303+P361+P353-P304+P340-P311a-P305+P351+P338-P403+P233

Titanium isopropylate > RE - Pure

RE

Description	Clear slightly yellow liq.	Density at 20° C	~ 0.965	Assay (gravimetric)	16.6 - 17.3 % Ti
Identification	Positive	Melting point	≥ 15 ° C		

Code	Size	Packaging	Notes
488421	100 ml	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Titanium trichloride-sulfuric acid reagent

• Reattivo titanio tricloruro-acido solforico • Titane trichlorure-acide sulfurique • Titanio tricloruro-ácido sulfúrico • Titantrichlorid-Schwefelsäure

$TiCl_3$
Molecular Weight: 154,22
CAS: 7705-07-9

Titanium trichloride-sulfuric acid reagent > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611091202	100 ml	Glass bottle	Ref Ph.Eur 1091202

Titanium trichloride-sulfuric acid reagent > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
617000221	100 ml	Plastic bottle	Titanium Trichloride-Sulfuric Acid TS

o-Tolidine solution 0.1%

• o-Tolidina soluzione 0.1% • o-Tolidine en solution à 0,1% • o-Tolidina solución 0.1% • o-Tolidinlösung 0.1%

Synonym:
3,3'-Dimethylbenzidine

$C_{14}H_{16}N_2$
Molecular Weight: 212,28
CAS: 119-93-7

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group III



Danger
H290-H350-HA26
P234-P280-P308+P313-P390-P406-P501a

o-Tolidine solution 0.1% > RPE - For analysis

RPE

Description ..... Clear colourless liquid Identification ..... Positive

Code	Size	Packaging	Notes
488461	1 l	Glass bottle	

For the determination of Au, Ce, Cl, halogen free, Mn

o-Tolidine solution

• o-Tolidina soluzione • o-Tolidine en solution • o-Tolidina solución • o-Tolidin in Lösung

$C_{14}H_{16}N_2$
Molecular Weight: 212,28
CAS: 119-93-7

o-Tolidine solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611123001	500 ml	Glass bottle	Ref Ph.Eur 1123001

Toluene

• Toluene • Toluène • Tolueno • Toluol

$C_6H_5CH_3$
Molecular Weight: 92,14
CAS: 108-88-3
EEC-N: 203-625-9

Classification transport
ONU: 1294
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H361d-H336-H373-H304
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Toluene > RS - For HPLC - Isocratic Grade

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.865 ÷ 0.869
Refractive index at 20°C. 1.4931 ÷ 1.4991
Boiling point ..... 110.1 ÷ 111.1 ° C

Acidity or alkalinity ..... ≤0.0015 meq/g
Water (K.F.) ..... ≤100 ppm
Residue on evaporation ..... ≤2 ppm
Assay (GLC) ..... ≥99.8 %
U.V. Transmittance

At 285 nm ..... ≥ 10 %
at 290 nm ..... ≥55 %
at 300 nm ..... ≥80 %
at 310 nm ..... ≥90 %
at 320 nm ..... ≥94 %

at 330 nm ..... ≥96 %
at 350 nm ..... ≥99 %

Code	Size	Packaging	Notes
412641000	1 l	Glass bottle	
412642000	2.5 l	Glass bottle	

Toluene > RS - For preparative HPLC

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.865 ÷ 0.869
Refractive index at 20°C. 1.4931 ÷ 1.4991
Boiling point..... 110.1 ÷ 111.1 °C
Water (K.F.) ..... ≤200 ppm
Residue on evaporation ..... ≤5 ppm
Alcalinity..... ≤0.0002 meq/g
Assay (GLC) ..... ≥99.8 %
U.V. Transmittance
at 300 nm ..... ≥75 %
at 350 nm ..... ≥98 %

Code	Size	Packaging	Notes
488531	2.5 l	Glass bottle	

Toluene > RS - ATRASOL - For traces analysis

RS

Appearance ..... Clear colourless liquid
Refractive index at 20°C..... 1.494 - 1.498
Water content (K.F.) ..... ≤ 50 mg/Kg
Colour ..... ≤ 10 Hazen
Free acid (as HCl) ..... ≤ 10 mg/Kg
Non volatile residue ..... ≤ 5 mg/Kg
Assay (GC) ..... ≥ 99.9 %
GC (FID) - NC Atrasol ..... Conform
GC-ECD.Individual peak (Lindane) .≤ 2 ng/L
Ret.range 1,2,4-trichlorobenzene
to decachlorobiphenyle
GC-FID.Individual peak (n-hexadecane). ≤ 5
µg/L
Ret.range n-dodecane to n-tetracontane

Code	Size	Packaging	Notes
P0713216	1 l	Glass bottle	
P0713221	2.5 l	Glass bottle	
P0713282	4 l	Glass bottle	

Toluene > RS - PESTIPUR - For pesticide analysis

RS

Description ..... Clear liquid
Colour ..... ≤ 10 hazen
Identification ..... Positive
Water ..... ≤ 100 ppm
Assay (GC) ..... ≥ 99.8 %
Acidity (HCl) ..... ≤ 10 ppm
Not volatile residue ..... ≤ 5 ppm
GC-ECD (Lindano) ..... ≤ 3 ng/l
Non volatile residue ..... ≤ 5 mg/Kg
GC-NPD (Ethylparation)..... ≤ 3 ng/l
GC-ECD.Individual peak (Lindane) .≤ 3 ng/L
Assay (GLC) ..... ≥ 99.8 %

Code	Size	Packaging	Notes
488591	1 l	Glass bottle	
488592	2.5 l	Glass bottle	
488594	4 l	Glass bottle	

Toluene > RS - SPECTROSOL - For optical spectroscopy

RS

Description ..... Clear liquid
Colour (APHA) ..... ≤10
Identification ..... Positive
Density at 20° C ..... 0.865 ÷ 0.869
Refractive index at 20°C..... 1.4931 ÷ 1.4991
Boiling point..... 110.1 ÷ 111.1 ° C
Acidity or alkalinity..... ≤0.0015 meq/g
Water (K.F.) ..... ≤100 ppm
Residue on evaporation ..... ≤5 ppm
Assay (GLC) ..... ≥99.8 %
Free acid (as HCl) ..... ≤ 10 mg/Kg
U.V. Transmittance at 285 nm ..... ≥10 %
U.V. Transmittance at 290 nm ..... ≥55 %
U.V. Transmittance at 300 nm ..... ≥80 %
U.V. Transmittance at 310 nm ..... ≥90 %
U.V. Transmittance at 320 nm ..... ≥93 %
U.V. Transmittance at 350 nm ..... ≥98 %

Code	Size	Packaging	Notes
488601	1 l	Glass bottle	
488602	2.5 l	Glass bottle	

Toluene > RS - Anhydrous - For analysis

RS

Refractive index at 20°C..... 1.494 - 1.498
Water content (K.F.) ..... ≤ 50 mg/Kg
Non volatile residue ..... ≤ 10 mg/Kg
Colour ..... ≤ 10 Hazen
Assay (GC) ..... ≥ 99.8 %
Free acid (as HCl) ..... ≤ 10 mg/Kg
Benzene ..... ≤ 200 mg/Kg
Styrene..... ≤ 10 mg/Kg
Ethylbenzene + xylene..... ≤ 500 mg/Kg

Code	Size	Packaging	Notes
P0711010	200 ml	Bottle with septum	
P07110T10	200 ml	Bottle with septum	On molecular sieves 4A
P0711016	1 l	Glass bottle	
P0711021	2.5 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Toluene > RPE - For analysis - ISO - ACS - Reag.Ph.Eur. - Reag.USP

RPE

Description ..... Clear colourless liquid
Identification (I.R.)..... Positive
Colour (APHA) ..... ≤10
Density at 20° C .....0.865 ÷ 0.869
Refractive index at 20°C.1.4931 ÷ 1.4991
Boiling point.....110.1 ÷ 111.1 ° C
Alcohol miscibility..... Complete
Chloroform miscibility ..... Complete
Diethyl ether miscib..... Complete

Acidity (benzoic acid) .....≤14 ppm
Alcalinity (NH3).....≤2 ppm
Water (K.F.) .....≤150 ppm
Residue on evaporation .....≤10 ppm
Ready carbonizable substances.....Conform
Benzene ..... ≤ 0.02 %
Tiophene .....≤1 ppm
Total sulphur .....≤3 ppm
Assay (GLC) ..... ≥99.8 %

Al.....≤0.5 ppm
B.....≤0.01 ppm
Ba.....≤0.1 ppm
Ca.....≤0.5 ppm
Cd.....≤0.01 ppm
Co.....≤0.01 ppm
Cr.....≤0.01 ppm
Cu.....≤0.01 ppm
Fe .....≤0.05 ppm

Mg .....≤0.1 ppm
Mn .....≤0.01 ppm
Ni.....≤0.01 ppm
Pb .....≤0.01 ppm
Sn .....≤0.1 ppm
Zn .....≤0.01 ppm
Styrene.....≤ 10 ppm
Ethyl benzene + xylene.....≤ 500 ppm

Code	Size	Packaging	Notes
488551	1 l	Glass bottle	
488555	2.5 l	Glass bottle	
488552	5 l	Aluminium can	
488557	24 kg	Metal drum	
488556	170 kg	Metal drum	

Toluene > RE - Pure

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C .....0.862 ÷ 0.872

Refractive index at 20°C.....1.491 ÷ 1.501
Boiling point.....109.9 ÷ 111.4 °C
Residue on evaporation .....≤50 ppm

Water (K.F.) .....≤300 ppm
Benzene ..... ≤ 0.02 %
Total sulphur .....≤100 ppm

Assay (GLC) ..... ≥99.5 %
Ethyl benzene + xylene.....≤ 750 ppm

Code	Size	Packaging	Notes
386002	1 l	Glass bottle	
386001	2.5 l	Glass bottle	
386003	23 kg	Metal drum	
386009	170 kg	Metal drum	

Toluene > RE - Pure - Low content in benzene

RE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C .....0.862 ÷ 0.872

Refractive index at 20°C.....1.491 ÷ 1.501
Boiling point.....109.9 ÷ 111.4 °C
Residue on evaporation .....≤50 ppm

Water (K.F.) .....≤300 ppm
Benzene ..... ≤ 0.02 %
Total sulphur .....≤100 ppm

Assay (GLC) ..... ≥99.5 %
Ethyl benzene + xylene.....≤ 750 ppm

Code	Size	Packaging	Notes
528231	5 l	Plastic tank	
528233	25 l	Metal drum	
528232	200 l	Metal drum	

Toluene > RE - ASTM

RE

Appearance ..... Clear liquid
Colour ..... ≤ 10 APHA

Peroxide .....≤ 5 ppm
Density at 15°C

Water (K.F.) .....≤ 200 ppm
Toluene (ASTM)..... ≥ 99.5 % (V/V)

Antioxidant..... 1500 - 2500 mg/L

Code	Size	Packaging	Notes
386102	5 l	Plastic tank	
386104	25 l	Metal drum	
386106	180 kg	Drum	

Suitable for ASTM methods D2700 and D2699

Toluene in solution in hexane

• Toluene soluzione in esano • Toluène en solution dans l'hexane • Tolueno en solución de hexano • Toluol in Lösung in Hexan

$C_6H_5CH_3$
Molecular Weight: 92,14
CAS: 108-88-3

Classification transport
ONU: 1993
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H361f-H336-H373-H411
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Toluene in solution in hexane > RS - For analysis according to Ph. Eur.Chap. 2.2.25 — RS

Code	Size	Packaging	Notes
506462	10 ml	Sealed cuvette	
506463	100 ml	Glass bottle	

Toluene-d8

• Toluene-d8 • Toluène-d8 • Tolueno-d8 • Toluol-d8

$C_6D_5CD_3$
Molecular Weight: 100,19
CAS: 2037-26-5
EEC-N: 218-009-5

Classification transport
ONU: 1294
Transport Hazard class: 3
Packing group II



Danger
H225-H315-H361d-H336-H373-H304
P210-P241-P280-P303+P361+P353-P304+P340-P403+P233

Toluene-d8 > RS - For NMR - min 99.5% — RS

Code	Size	Packaging	Notes
P5399A	2 x 0.75 ml	Glass ampoule	
P5393A	5 ml	Glass ampoule	
P5395	25 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

p-Toluene sulfonamide

• p-Toluensolfonammide • p-Toluènesulfonamide • p-Toluenosulfonamida • p-Toluolsulfonamid

$CH_3C_6H_4SO_2NH_2$
Molecular Weight: 171,22
CAS: 70-55-3
EEC-N: 200-741-1

p-Toluene sulfonamide > RPE - For analysis — RPE

Description ..... White crystalline powder
Identification ..... Positive
Melting point ..... ~ 136 ° C
Assay (ex nitrogen) ..... ≥98 %

Code	Size	Packaging	Notes
488661	100 g	Plastic bottle	

p-Toluenesulfonic acid

• Acido p-toluensolfonico monoidrato • Acide p-toluenesulfonique monohydrate • Acido p-toluenosulfónico • Toluol-4-sulfonsäure

$CH_3C_6H_4SO_3H.H_2O$
Molecular Weight: 190,22
CAS: 6192-52-5
EEC-N: 203-180-0

Classification transport
ONU: 2585
Transport Hazard class: 8
Packing group III



Danger
H314-H335
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

p-Toluenesulfonic acid > RPE - For analysis — RPE

Description ..... White crystals
Identification ..... Positive
Melting point ..... 99 ÷ 103 °C
Water (K.F.) ..... ≤13 %
Acidity(Sulphuric acid) ..... ≤ 1 %
Fe ..... ≤50 ppm
Assay (oxidimetric) ..... ≥ 97 %

Code	Size	Packaging	Notes
411436	500 g	Plastic bottle	
411432	20 kg	Plastic bucket	

DASIT GROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

p-Toluenesulfonic acid > RE - Pure

RE

Description ..... White crystals Identification ..... Positive Fe ..... ≤100 ppm Assay (acidimetric) ..... ≥97 %

Code	Size	Packaging	Notes
307508	10 kg	Carton box	

p-Toluenesulfonic acid sodium salt

• Acido p-toluensolfonico sale sodico • Acide p-toluenesulfonique sel sodique
• Acido p-toluenosulfónico sal sódica • p-Toluolsulfonsäure-Natriumsalz

Synonym:
Sodium p-toluenesulfonate

$CH_3C_6H_4SO_3Na$
Molecular Weight: 194,19
CAS: 657-84-1
EEC-N: 211-522-5

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

p-Toluenesulfonic acid sodium salt > RPE - For analysis

RPE

Description ..... White crystals Identification ..... Positive Assay ..... ≥ 97.5 %

Code	Size	Packaging	Notes
411504	100 g	Glass bottle	

Toluidine blue

• Blu toluidina • Bleu de toluidine • Azul de toluidina • Toluidinblau

Synonym:
Basic Blue 17

$C_{15}H_{16}ClN_3S$
Molecular Weight: 305,83
CAS: 92-31-9
EEC-N: 202-146-2

Toluidine blue > RS - For microscopy - C.I. 52040

RS

Description ..... Black powder Identification ..... Positive Spettro (UV) ..... Conform

Code	Size	Packaging	Notes
429282	25 g	Glass bottle	

Dye for cytology-histochemistry

Total-ionic-strength-adjustment buffer

• Tampone forza ionica totale • Tampon pour ajustement de la force ionique totale • Tampón fuerza iónica total
• Puffer zur Einstellung der Gesamtionenstärke

Total-ionic-strength-adjustment buffer > RS - For analysis according to Ph. Eur. Chap. 4.1.3

RS

Code	Size	Packaging	Notes
614007700	1 l	Plastic bottle	Ref Ph.Eur 4007700
614008800	1 l	Plastic bottle	Ref Ph.Eur 4008800

Triacetin

• Triacetina • Triacétine • Triacetina • Triacetin

Synonym:
Glyceryl triacetate

$C_9H_{14}O_6$
Molecular Weight: 218,21
CAS: 102-76-1
EEC-N: 203-051-9

Triacetin > RPE - For analysis

RPE

Description ..... Clear colourless liquid Acidity (acetic acid) ..... ≤100 ppm Assay (GLC) ..... ≥99.0 %
Identification ..... Positive Water ..... ≤0.1 %

Code	Size	Packaging	Notes
489152	1 l	Glass bottle	

2,2,2 - Trichlorethanol

• 2,2,2 - Tricloroetanolo • 2,2,2 - Trichloroéthanol • 2,2,2 - Tricloroetanol • 2,2,2 - Trichlorethanol

Cl_3CCH_2OH
Molecular Weight: 149,4
CAS: 115-20-8
EEC-N: 204-071-0



Danger
H302-H318-H336
P261-P271-P304+P340-P310a-P305+P351+P338-P403+P233

2,2,2 - Trichlorethanol > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Refractive index at 20°C. 1.4880 ÷ 1.4910
Assay (GLC) ..... ≥98.5 %

Code	Size	Packaging	Notes
415271	100 ml	Glass bottle	

Trichloroacetic acid

• Acido tricloroacetico • Acide trichloroacétique • Acido tricloroacetico • Trichloressigsaure

Synonym:
TCA

CCl_3COOH
Molecular Weight: 163,39
CAS: 76-03-9
EEC-N: 200-927-2

Classification transport
ONU: 1839
Transport Hazard class: 8
Packing group II



Danger
H314-H335-H410
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Trichloroacetic acid > RPE - For analysis - ACS - Reag. Ph.Eur. - Reag. USP

RPE

Description ..... Colourless crystals
Identification ..... Positive
Ready carbonizable substances ..... Conform
Residue on ignition ..... ≤ 0.03 %
Chloride ..... ≤ 20 ppm
Phosphate ..... ≤ 5 ppm
Nitrate ..... ≤ 20 ppm
Sulphate ..... ≤ 0.02 %
Heavy metals (Pb) ..... ≤ 20 ppm
Fe ..... ≤ 10 ppm
Assay (acidimetric) ..... ≥ 99.0 %
Water ..... ≤ 0.5 %

Code	Size	Packaging	Notes
411524	100 g	Glass bottle	
411525	250 g	Plastic bottle	
411527	1 kg	Plastic bottle	

Trichloroacetic acid > RE - Pure

RE

Description ..... Colourless crystals
Identification ..... Positive
Water ..... ≤ 0.5 %
Fe ..... ≤ 10 ppm
Assay (acidimetric) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
307557	1 kg	Plastic bottle	
307558	5 kg	Plastic tank	

Trichloroacetic acid solution 20%

• Acido tricloroacetico soluzione 20% • Acide trichloroacétique 20% • Acido tricloroacético solución 20% • Trichloressigsäure 20%

Synonym:
TCA

CCl_3COOH
Molecular Weight: 163,39
CAS: 76-03-9

Classification transport
ONU: 2564
Transport Hazard class: 8
Packing group II



Danger
H314-H335-H400-H411
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Trichloroacetic acid solution 20% > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Assay (acidimetric) ..... 19.5 - 20.5 %

Code	Size	Packaging	Notes
502073	100 ml	Bottle	
411554000	1 l	Glass bottle	

For the determination of iron in the blood according Heimayer. For protein precipitation. Content is guaranteed for standardized volumes at 20°C. Keep tightly closed

Product specifications are subject to changes. Please visit our website for updates.

1,2,4-Trichlorobenzene

• 1,2,4-Triclorobenzene • 1,2,4-Trichlorobenzène • 1,2,4-Triclorobenceno • 1,2,4-Trichlorbenzol

$C_6H_3Cl_3$
Molecular Weight: 181,45
CAS: 120-82-1
EEC-N: 204-428-0

Classification transport
ONU: 2321
Transport Hazard class: 6.1
Packing group III



Warning
H302-H315-H410
P264-P280g-P301+P312a-P332+P313-P362+P364-P501a

1,2,4-Trichlorobenzene > RS - SPECTROSOL - For optical spectroscopy

RS

Refractive index at 20°C..1.5697 - 1.5737
Water content (K.F.).....≤ 100 mg/Kg
Colour ..... ≤ 10 Hazen
UV transmittance at 310 nm ..... ≥ 40 %
UV transmittance at 315 nm ..... ≥ 80 %
UV transmittance at 385 nm ..... ≥ 98 %
Assay (GC)..... ≥ 99 %
Non volatile residue .....≤ 10 mg/Kg

Code	Size	Packaging	Notes
P0722721	2.5 l	Glass bottle	

1,2,4-Trichlorobenzene > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 1.451 ÷ 1.457
Refractive index at 20°C..... 1.5687 ÷ 1.5747
Boiling point..... 212.5 ÷ 213.5 ° C
Melting point..... 16.0 ÷ 18.0 ° C
Water (K.F.) .....≤0.1 %
Residue on ignition .....≤10 ppm
Assay (GLC) .....≥98.5 %

Code	Size	Packaging	Notes
489382	1 l	Glass bottle	

Should be stored at not less than 20 °C

Trichloromethane ▶ Chloroform

Tricresol

• Tricresolo • Tricrésol • Tricresol • Trikresol

$HOC_6H_4CH_3$
Molecular Weight: 108,14
CAS: 1319-77-3
EEC-N: 215-293-2

Classification transport
ONU: 2022
Transport Hazard class: 6.1
Packing group II



Danger
H301-H311-H314
P280-P301+P310a-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338

Tricresol > RE - Pure

RE

Description ..... Brown clear liquid
Identification ..... Positive
Density at 25° C .....≥1.035

Code	Size	Packaging	Notes
386202	1 l	Glass bottle	

Mixture of isomers

Triethanolamine

• Trietanolammina • Triéthanolamine • Trietanolamina • Triethanolamin

Synonym:
• *2,2',2''-Nitrilotriethanol*
• *Tris(2-hydroxyethyl)amine*

$(CH_2OHCH_2)_3N$
Molecular Weight: 149,19
CAS: 102-71-6
EEC-N: 203-049-8

Triethanolamine > RPE - For analysis

RPE

Description ..... Yellowish liquid
Identification ..... Positive
Water miscibility.....Conform
Alcohol miscibility..... Complete
Density at 20° C .....1.120 ÷ 1.128
Refractive index at 20°C..1.4797 ÷ 1.4907
Melting point.....20.0 ÷ 22.0 ° C
Water (K.F.) ..... ≤0.3 %
Chloride.....≤10 ppm
Diethanolamine ..... ≤1.5 %
Monoethanolamine..... ≤0.5 %
Heavy metals (Pb).....≤2 ppm
Residue on ignition.....≤50 ppm
Sulphate .....≤20 ppm
Fe .....≤2 ppm
Assay (non-aqueos medium)..... ≥98 %

Code	Size	Packaging	Notes
489504	1 l	Glass bottle	
489501	30 kg	Metal drum	

Keep in a dark place

Triethanolamine > ERBApharm - According to pharmacopoeia: FU-Ph.Eur.

ERBApharm

Description ..... Clear colourless liquid or yellowish
Identification ..... Positive
Identification B ..... pass test
Identification C ..... Pass test
Density at 20° C ..... 1.120 - 1.128
Refractive index at 20°C ..... 1.482 ÷ 1.485
Water (K.F.) ..... ≤ 0.5 %
Sulphated ash ..... ≤ 0.05 %
Diethanolamine ..... ≤ 0.5 %
Monoethanolamine ..... ≤ 0.1 %
Related substances ..... ≤ 1.0 %
Appearance of solution ..... Conform Ph.Eur.
Total basis ..... 99.0 - 103.0 % anidro
Heavy metals (Pb) ..... ≤ 10 ppm
N-Nitrosodiethanolamine ..... ≤ 24 ppb

Code	Size	Packaging	Notes
386301	1 l	Glass bottle	
386303	2.5 l	Glass bottle	
386304	30 kg	Metal drum	
386305	220 kg	Metal drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Triethylamine

• Trietilammina • Triéthylamine • Trietilamina • Triethylamin

Synonym:
N,N-Diethylethanamine

$(C_2H_5)_3N$
Molecular Weight: 101,19
CAS: 121-44-8
EEC-N: 204-469-4

Classification transport
ONU: 1296
Transport Hazard class: 3
Packing group II



Danger
H225-H302-H311-H331-H314-H335
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P361+P364-P403+P233

Triethylamine > RS - For HPLC - Isocratic Grade

RS

Clear,colourless liquid ..... Conform
Water content (K.F.) ..... ≤ 0.05 % (m/m)
Residue on evaporation ..... ≤ 0.001 % (m/m)
Assay (GC) ..... ≥ 99.7 %
UV transmittance at 250nm (0.1M) ..... ≥ 10 %
UV transmittance at 254nm (0.1M) ..... ≥ 75 %

Code	Size	Packaging	Notes
489631	1 l	Glass bottle	
489633	2.5 l	Glass bottle	

Triethylamine > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 0.725 ÷ 0.729
Refractive index at 20°C. 1.3983 ÷ 1.4023
Boiling point ..... 89.0 ÷ 90.0 °C
Residue on evaporation ..... ≤100 ppm
Assay (GLC) ..... ≥99.5 %

Code	Size	Packaging	Notes
489556	1 l	Glass bottle	

Triethylamine > RE - Pure

RE

Description ..... Clear colourless liquid
Water content (K.F.) ..... ≤ 1000 mg/Kg
Identification ..... Positive
Density at 20° C ..... 0.724 ÷ 0.730
Refractive index at 20°C. 1.3953 ÷ 1.4053
Boiling point ..... 88.5 ÷ 90.5 °C
Residue on evaporation ..... ≤0.02 %
Assay (GLC) ..... ≥99.5 %
Diethylamine ..... ≤ 0.1 %

Code	Size	Packaging	Notes
386601	1 l	Glass bottle	
386603	5 l	Plastic tank	
386602	20 kg	Plastic tank	
P0790266	200 l	Combined drum	

Product specifications are subject to changes. Please visit our website for updates.

Triethylene glycol

• Glicol trietilenico • Glycol triéthylénique • Trietilenglicol • Triethylenglycol

Synonym: *Triglycol*

$(CH_2OHCH_2OCH_2)_2$
Molecular Weight: 150,18
CAS: 112-27-6
EEC-N: 203-953-2

Triethylene glycol > RPE - For analysis

RPE

Description	Clear colourless liquid	Density at 20° C	1.123 ÷ 1.131	Chloride	≤ 2 ppm	Sulphate	≤ 20 ppm
Identification	Positive	Refractive index at 20°C	1.4553 ÷ 1.4603	Heavy metals (Pb)	≤ 2 ppm	Fe	≤ 2 ppm
Water miscibility	Conform	Water (K.F.)	≤ 0.2 %	Peroxides (H2O2)	≤ 50 ppm	Assay (GLC)	≥ 98 %
Alcohol miscibility	Complete	Acidity (acetic acid)	≤ 150 ppm	Residue on ignition	≤ 30 ppm		

Code	Size	Packaging	Notes
454111	1 l	Glass bottle	
454112	30 kg	Plastic drum	

Trifluoroacetic acid

• Acido trifluoroacetico • Acide trifluoroacétique • Acido trifluoroacético • Trifluoressigsäure

Synonym: *TFA*

CF_3COOH
Molecular Weight: 114,02
CAS: 76-05-1
EEC-N: 200-929-3

Classification transport
ONU: 2699
Transport Hazard class: 8
Packing group I

Danger
H332-H314-H412
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Trifluoroacetic acid > RS - For LC/MS

RS

Description	Clear colourless liquid	Water (K.F.)	≤ 0.05 %	Fluoride	≤ 50 ppm	Assay (acidimetric)	≥ 99.9 %
Identification	Positive	Chloride	≤ 10 ppm	Sulphate	≤ 10 ppm	Suitability	LC-MS Tested

Code	Size	Packaging	Notes
411541	10 x 1 ml	Glass ampoule	
411542	10 x 2.5 ml	Glass ampoule	
411543	50 ml	Plastic bottle	

Eluent phase additive

Trifluoroacetic acid > RS - SPECTROSOL - For optical spectroscopy

RS

Appearance	Clear colourless liquid	Water content (K.F.)	≤ 0.05 % m/m	UV Absorbance at 280 nm	≤ 0.10 AU	UV Absorbance at 320 nm	≤ 0.05 AU
Identification (IR)	Conform	UV Absorbance at 260 nm	≤ 1.2 AU	UV Absorbance at 290 nm	≤ 0.09 AU	Content (Acidimetry)	99.9 - 101.0 % m/m
Density d20/4	1.480 - 1.500	UV Absorbance at 270 nm	≤ 0.15 AU	UV Absorbance at 300 nm	≤ 0.08 AU	Residue on evaporation	≤ 0.002 % m/m

Code	Size	Packaging	Notes
P0082746	1 l	Glass bottle PVC coated	
P0082747	2.5 l	Glass bottle PVC coated	

Trifluoroacetic acid > RS - For peptide synthesis

RS

Identification (IR)	Conform	Colour	≤ 10 Hazen	Content (Acidimetry)	≥ 99.9 % m/m	Sulphate (SO4--)	≤ 20 mg/Kg
Density d20/4	1.480 - 1.500	Water content (K.F.)	≤ 500 mg/Kg	Chloride (Cl-)	≤ 20 mg/Kg	Fluoride	≤ 50 mg/Kg

Code	Size	Packaging	Notes
P0082103	100 ml	Glass bottle	
P0082112	1 l	Glass bottle PVC coated	
P0082147	2.5 l	Glass bottle PVC coated	

Trifluoroacetic acid > RPE - For analysis

RPE

Description	Clear colourless liquid	Water (K.F.)	≤ 0.05 %	Sulphate	≤ 10 ppm
Identification	Positive	Assay (acidimetric)	≥ 99.9 %		

Code	Size	Packaging	Notes
411561	100 ml	Glass bottle	
411564	250 ml	Glass bottle	

Trifluoroacetic acid > RE - Pure

RE

Identification (IR)..... Conform
Density d20/4 ..... 1.480 - 1.500
Colour ..... ≤ 10 Hazen
Water content (K.F.) ..... ≤ 1000 mg/Kg
Content (Acidimetry) ..... ≥ 99 % m/m
Chloride (Cl-)..... ≤ 20 mg/Kg
Sulphate (S04--) ..... ≤ 20 mg/Kg
Fluoride ..... ≤ 50 mg/Kg

Code	Size	Packaging	Notes
P0080247	2.5 l	Glass bottle PVC coated	
P0080212	1 kg	Glass bottle PVC coated	
P0080297	30 kg	Plastic bucket	

Trifluoroacetic acid-d

• Acido trifluoroacetico-d • Acide trifluoroacétique-d • Acido trifluoroacético-d • Trifluoressigsäure-d

CF_3COOD
Molecular Weight: 115,03
CAS: 599-00-8
EEC-N: 209-961-2

Classification transport
ONU: 2699
Transport Hazard class: 8
Packing group I



Danger
H332-H314-H412
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Trifluoroacetic acid-d > RS - For NMR - min 99.5%

RS

Code	Size	Packaging	Notes
P5419A	2 x 0.75 ml	Glass ampoule	
P5413A	5 ml	Glass bottle	

For specifications, contact our customer service for a certificate of analysis

Trifluoroacetic anhydride

• Anidride trifluoroacetica • Anhydride trifluoroacétique • Anhídrido trifluoroacético
• Trifluoressigsäureanhydrid

Synonym:
TFAA

$(CF_3CO)_2O$
Molecular Weight: 210,04
CAS: 407-25-0
EEC-N: 206-982-9

Classification transport
ONU: 3265
Transport Hazard class: 8
Packing group I



Danger
H301-H314
P280-P301+P310a-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338

Trifluoroacetic anhydride > RPE - For analysis

RPE

Description ..... Clear colourless liquid
Identification ..... Positive
Density at 20° C ..... 1.511 ÷ 1.515
Residue on evaporation ..... ≤10 ppm
Assay (as anhydride)..... ≥98 %

Code	Size	Packaging	Notes
422225	500 ml	Glass bottle	

For derivatization

1,2,3-Trihydroxybenzene ▶ Pyrogallol

1,3,5-Trihydroxybenzene ▶ Phloroglucinol

3,4,5-Trihydroxybenzoic acid monohydrate ▶ Gallic acid monohydrate

2,4,6-Trihydroxypyrimidine ▶ Barbituric acid

Triiodomethane ▶ Iodoform

Product specifications are subject to changes. Please visit our website for updates.

Trimethylcetylammonium bromide

• Trimetilcetilammonio bromuro • Triméthylcétylammonium bromure • Trimetilcetilamonio bromuro
• Trimethylcetylammoniumbromid

Synonym:
• *Hexadecyltrimethylammonium bromide*
• *CTAB*

$CH_3(CH_2)_{15}N(CH_3)_3Br$
Molecular Weight: 364,46
CAS: 57-09-0
EEC-N: 200-311-3

Classification transport
ONU: 3077
Transport Hazard class: 9
Packing group III

Danger
H302-H315-H318-H335-H410
P304+P340-P310a-P305+P351+P338-P330-P362+P364-P403+P233

Trimethylcetylammonium bromide > RPE - For analysis

RPE

Description ..... White crystalline powder Identification ..... Positive Assay (non-aqueos medium) ..... ≥ 99.0 %

Code	Size	Packaging	Notes
489833	50 g	Plastic bottle	
489831	500 g	Plastic bottle	

Trimethylcetylammonium p-toluenesulfonate

• Trimetilcetilammonio p-toluensolfonato • Triméthylcétylammonium p-toluènesulfonate
• Trimetilcetilamonio p-toluensulfonato • Trimethylcethylammonium p-toluenesulphonate

Synonym:
• *Hexadecyltrimethylammonium p-toluenesulfonate*
• *Cetyltrimethylammonium p-toluenesulfonate*

$CH_3C_6H_4SO_3N(CH_3)_3(CH_2)_{15}CH_3$
Molecular Weight: 455,74
CAS: 138-32-9
EEC-N: 205-324-8

Warning
H302-H312-H332-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Trimethylcetylammonium p-toluenesulfonate > RE - Pure

RE

Description ..... White powder Soluz. in acqua (1:200) ..... Complete Assay ..... ≥ 99.0 %
Identification ..... Positive pH sol. 1% ..... 5.0 ÷ 8.0

Code	Size	Packaging	Notes
387203	25 kg	Fibre drum	

2,2,4-Trimethylpentane ▶ Isooctane

1,3,7-Trimethylxanthine ▶ Caffeine anhydrous

2,4,6-Trinitrophenol ▶ Picric acid solution

Trioctylphosphine oxide

• Triottilfosfina ossido • Trioctylphosphine oxyde • Trioctilfosfina óxido • Trioctylphosphinoxid

Synonym:
TOPO

$(C_8H_{17})_3PO$
Molecular Weight: 386,65
CAS: 78-50-2
EEC-N: 201-121-3

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Trioctylphosphine oxide > RPE - For analysis

RPE

Description ..... White crystals Identification ..... Positive Melting point ..... 52 ÷ 57 °C Assay (GLC) ..... ≥ 98.5 %

Code	Size	Packaging	Notes
489581	50 g	Glass bottle	

Triphenylphosphine

• Trifenilfosfina • Triphénylphosphine • Trifenilfosfina • Triphenylphosphin

Synonym:
Phosphorustriphenyl

$(C_6H_5)_3P$
Molecular Weight: 262,29
CAS: 603-35-0
EEC-N: 210-036-0



Warning
H302-H317-H410
P261-P264-P280g-P301+P312a-P333+P313-P501a

Triphenylphosphine > RPE - For analysis

RPE

Description ..... White powder
Identification ..... Positive
Melting point..... 80.0 ÷ 82.0 ° C
Assay (GLC) ..... ≥98 %

Code	Size	Packaging	Notes
489591	100 g	Glass bottle	

2,3,5-Triphenyltetrazolium chloride

• 2,3,5-Trifeniltetrazolo cloruro • 2,3,5-Triphényltétrazolium chlorure • 2,3,5-Trifeniltetrazolio cloruro
• 2,3,5-Triphenyltetrazolium-chlorid

Synonym:
• *TPTZ*
• *Tetrazolium Red*

$(C_6H_5)_3CN_4Cl$
Molecular Weight: 334,81
CAS: 298-96-4
EEC-N: 206-071-6



Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

2,3,5-Triphenyltetrazolium chloride > RPE - For analysis

RPE

Description ..... White yellowish powder
Identification ..... Positive
Water (K.F.) ..... ≤ 1.0 %
Assay (non-aqueos medium) 99 ÷ 100 % s
s

Code	Size	Packaging	Notes
489651	10 g	Glass bottle	

Tris (hydroxymethyl)-aminomethane

• Tris idrossimetilamminometano • Tris(hydroxyméthyl)aminométhane • Tris (hidroximetil)aminometano
• Tris(hydroxymethyl)-aminomethan

Synonym:
Tris base

$NH_2C(CH_2OH)_3$
Molecular Weight: 121,14
CAS: 77-86-1
EEC-N: 201-064-4



Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Tris (hydroxymethyl)-aminomethane > RS - For pHmetry

RS

Description ..... White crystalline powder
Identification ..... Positive
pH solution 5% ..... 10.5 ÷ 11.5
Solution colour ..... ≤25 APHA
Melting point..... ≥168.5 °C
Water (K.F.) ..... ≤0.5 %
Residue on ignition..... ≤0.1 %
Heavy metals (Pb)..... ≤5 ppm
Absorbance (1M aq, 290 nm)..... ≤ 0.050
Assay (alkalimetric)..... ≥99.3 % s.s.

Code	Size	Packaging	Notes
489973	1 kg	Plastic bottle	
489971	25 kg	Plastic bucket	

Tris (hydroxymethyl)-aminomethane > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Melting point..... 168 ÷ 172 ° C
pH (1M a 25 °C)..... 10.5 ÷ 11.5
Loss on drying ..... ≤0.5 %
Water-insoluble matter ..... ≤100 ppm
Heavy metals (Pb)..... ≤2 ppm
Residue on ignition..... ≤ 0.1 %
Fe ..... ≤1 ppm
Assay (alkalimetric)..... ≥99.5 %

Code	Size	Packaging	Notes
489981	100 g	Plastic bottle	
489984	500 g	Plastic bottle	
489983	1 kg	Plastic bottle	
489985	25 kg	Plastic bucket	

Product specifications are subject to changes. Please visit our website for updates.

Tris (hydroxymethyl)-aminomethane > ERBAPharm - According to pharmacopoeia: USP

ERBApharm

Description ..... White crystalline powder
Identification ..... Positive
pH solution 5% ..... 10.0 ÷ 11.5
Melting point..... 168 ÷ 172 °C
Loss on drying ..... ≤ 1.0 %
Residue on ignition..... ≤ 0.1 %
Heavy metals (Pb)..... ≤ 10 ppm
Titolo (alcalimetrico)99.0 ÷ 101.0 % (s.s.)

Code	Size	Packaging	Notes
313441	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Tris (hydroxymethyl)-aminomethane > RE - Pure

RE

Description ..... White crystalline powder
Identification ..... Positive
Melting point..... 168 ÷ 172 ° C
pH (1M a 25 °C)..... 10.5 ÷ 11.5
Loss on drying ..... ≤1 %
Water-insoluble matter ..... ≤100 ppm
Heavy metals (Pb)..... ≤2 ppm
Residue on ignition..... ≤ 0.1 %
Fe ..... ≤2 ppm
Assay (alkalimetric)..... ≥99 %

Code	Size	Packaging	Notes
313432	1 kg	Plastic bottle	
313431	25 kg	Plastic bucket	

Tris(hydroxymethyl)aminomethane solution

• Tris idrossimetilamminometano soluzione • Tris(hydroxyméthyl)aminométhane solution • Tris (hidroximetil)aminometano solución
• Tris (hydroxymethyl) aminomethanlösung

Tris(hydroxymethyl)aminomethane solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611094201	100 ml	Plastic bottle	Ref Ph.Eur 1094201

Tris(hydroxymethyl)aminomethane buffer solution pH 8.1

• Tampone tris(idrossimetil)aminometano tampn soluzione pH 8.1 • Tampon tris(hydroxyméthyl)aminométhane pH 8.1
• Tris (hidroximetil)aminometano solución tampón pH 8.1 • Tris (hydroxymethyl) aminomethan-Pufferlösung pH 8.1

Tris(hydroxymethyl)aminomethane buffer solution pH 8.1 > RS - For analysis according to Ph. Eur. Chap. 4.1.3

RS

Code	Size	Packaging	Notes
614006200	1 l	Plastic bottle	Ref Ph.Eur 4006200

Tris(hydroxymethyl)aminomethane-EDTA buffer solution pH 8.4

• Tampone tris(idrossimetil)aminometano - EDTA pH 7.4 • Tampon tris(hydroxyméthyl)aminométhane-EDTA pH 8.4
• Tampón tris(hidroximetil)aminometano - EDTA pH 7.4 • Tris (hydroxymethyl) aminomethan-EDTA-Pufferlösung pH 8.4

Tris(hydroxymethyl)aminomethane-EDTA buffer solution pH 8.4 > RS - For analysis according to Ph. Eur. Chap. 4.1.3

RS

Code	Size	Packaging	Notes
614006600	1 l	Plastic bottle	Ref Ph.Eur 4006600

Tris (hydroxymethyl)-aminomethane hydrochloride

• Tris idrossimetilamminometano cloridrato • Tris(hydroxyméthyl)aminométhane chlorhydraté
• Tris (hidroximetil)aminometano hidrocloruro • Tris(hydroxymethyl)-aminomethan HCl

Synonym:
2-Amino-2(hydroxymethyl)1,3-propanediol hydrochloride

$NH_2C(CH_2OH)_3.HCl$
Molecular Weight: 157,6
CAS: 1185-53-1
EEC-N: 214-684-5

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Tris (hydroxymethyl)-aminomethane hydrochloride > RPE - For analysis

RPE

Description ..... White crystals
Identification ..... Positive
Pb ..... ≤ 5 ppm
As ..... ≤ 10 ppm
Ba ..... ≤ 1 ppm
Ca ..... ≤ 10 ppm
Cu ..... ≤ 1 ppm
Fe ..... ≤ 2 ppm
Mg ..... ≤ 10 ppm
Mn ..... ≤ 1 ppm
Assay ..... ≥ 99 %
Absorbance (1M acq. 280 nm) ..... ≤ 0.05
Melting point ..... 150 ÷ 154 °C

Code	Size	Packaging	Notes
479911	100 g	Plastic bottle	
479912	500 g	Plastic bottle	
479913	1 kg	Plastic bottle	

Tris(hydroxymethyl)aminomethane sodium chloride buffer solution pH 7.4

• Tampone tris(idrossimetil)aminometano sodio cloruro pH 7.4 • Tampon tris(hydroxyméthyl)aminométhane-chlorure de sodium pH 7.4
• Tampón tris(hidroximetil)aminometano sodio cloruro pH 7.4 • Tris (hydroxymethyl) aminomethan-Natriumchlorid-Pufferlösung pH 7.4

Tris(hydroxymethyl)aminomethane sodium chloride buffer solution pH 7.4 > RS - For analysis according to Ph. Eur. Chap. 4.1.3

RS

Code	Size	Packaging	Notes
614004900	1 l	Plastic bottle	Ref Ph.Eur 4004900

Triton© X100 solution

• Triton© X100 soluzione • Triton© X100 solution • Triton© X100 solución • Triton© X100-Lösung

Synonym:
4-(1,1,3,3-Tetramethylbutyl)phenyl-polyethylene glycol

$4\text{-}(C_8H_{17})C_6H_4(OCH_2CH_2)nOH$
CAS: 9002-93-1

Danger
H302-H318-H412
P264-P280i-P301+P312a-P305+P351+P338-P310a-P501a

Triton© X100 solution > RE - Pure

RE

Clear to light unclear liquid ..... Conform
Appearance ..... Without particles in suspension
Density d20/4 ..... 1.055 - 1.075
Water content (K.F.) ..... ≤ 2000 mg/Kg
Colour ..... ≤ 60 Hazen
Cloud point (1% in water) ..... 63 - 69 °C

Code	Size	Packaging	Notes
P0120041	10 l	Plastic tank	

Tropaeolin O

• Tropeolina O • Tropéoline O • Tropeolina O • Tropaeolin O

Synonym:
• *Tropaeolin O sodium salt*
• *Acid Orange 6*

$C_{12}H_9N_2NaO_5S$
Molecular Weight: 316,27
CAS: 547-57-9
EEC-N: 208-924-8

Warning
H315-H319
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Tropaeolin O > RPE - For analysis - C.I. 14270

RPE

Description ..... Polvere marrone arancio
Identification ..... Positive
Loss on drying ..... ≤15 %
Colour change ..... giallo arancio
pH range ..... 11.1 - 12.7

Code	Size	Packaging	Notes
490001	10 g	Glass bottle	
490002	25 g	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Tropaeolin O solution 0.1%

• Tropeolina O soluzione 0,1% • Tropéoline O solution 0.1% • Tropeolina O solución 0.1%
• Tropaeolin 0 0.1 %

Synonym:
• *Tropaeolin O sodium salt*
• *Acid Orange 6*

$C_{12}H_9N_2NaO_5S$
Molecular Weight: 316,27
CAS: 547-57-9

Tropaeolin O solution 0.1% > RPE - For analysis

RPE

Description ..... Orange clear liquid
Identification ..... Positive
Sensitivity (pH 11-13) ..... Conform
Colour change ..... giallo arancio
pH range ..... 11.1 - 12.7

Code	Size	Packaging	Notes
E490056	500 ml	Bottle	

Tungsten standard solution

• Tungsteno standard soluzione • Tungstène solution standard • Tungsteno, solución patrón • Wolfram-Standardlösung

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group III



Tungsten standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505932	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Ammonium hydroxyde
505935	100 ml	Plastic bottle	conc. 100 ppm. Matrix: Ammonium hydroxyde

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Tungsten standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
504058	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Water
504051	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: 4% ammonia
504053	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: 4% ammonia
504055	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: 10% ammonia
504057	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: 10% ammonia

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Tungsten standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507765	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid
507521	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Tungstic acid
• Acido tungstico • Acide tungstique • Acido túngstico • Wolframsäure

H_2WO_4
Molecular Weight: 249,86
CAS: 7783-03-1
EEC-N: 231-975-2

Tungstic acid > RPE - For analysis

RPE

Description ..... Yellow-greenish powder
Identification ..... Positive
Loss on ignition ..... 6 ÷ 8 %
Fe ..... ≤ 20 ppm
Assay (WO3) ..... ≥ 92.0 % t.q.
Al ..... ≤ 10 ppm
As ..... ≤ 10 ppm
Cu ..... ≤ 5 ppm
Mo ..... ≤ 50 ppm
Na ..... ≤ 30 ppm
Si ..... ≤ 30 ppm
Ti ..... ≤ 10 ppm

Code	Size	Packaging	Notes
411628	250 g	Plastic bottle	

Tungsten (VI) oxide
• Anidride tungstica • Anhydride tungstique • Anhídrido túngstico • Wolfram (VI) oxid

WO_3
Molecular Weight: 231,85
CAS: 1314-35-8
EEC-N: 215-231-4

Warning
H302
P264-P270-P301+P312a-P330-P501a

Tungsten (VI) oxide > RS - For microanalysis

RS

Description ..... Yellow powder
Identification ..... Positive
Assay (gravimetric) ..... ≥99 %

Code	Size	Packaging	Notes
422241	100 g	Glass bottle	

Turbidity std 4000NTU formazine
• Torbidità std 4000NTU formazina • Standard turbidité 4000NTU Formazine • Turbiedad std 4000NTU formazina • Standardtrübung 4000NTU Formazin

Danger
H332-H317-H350-HA26
P261-P271-P280-P304+P340-P308+P313-P362+P364

Turbidity std 4000NTU formazine > RS - For turbidimetry

RS

Description ..... Liquido opalescente
Identification ..... Positive

Code	Size	Packaging	Notes
489162	500 ml	Glass bottle	

Turk's reagent
• Turk reattivo • Réactif de Turk • Türk reactivo • Türkisches Reagenz

HEU210

Turk's reagent > RS - For microscopy

RS

Description ..... Purple liquid
Identification ..... Positive

Code	Size	Packaging	Notes
E490451	500 ml	Plastic bottle	In Vitro Diagnostic Medical Device

Dye for hematology

Product specifications are subject to changes. Please visit our website for updates.

L(-)Tyrosine

• L(-)Tirosina • L(-)Tyrosine • L(-)Tirosina • L-Tyrosin

Synonym:
(S)-2-Amino-3-(4-hydroxyphenyl)propionic acid

$4\text{-}HOC_6H_4CH_2CHNH_2COOH$
Molecular Weight: 181,19
CAS: 60-18-4
EEC-N: 200-460-4

L(-)Tyrosine > RPE - For analysis

RPE

Description ..... White crystalline powder	Loss on drying ..... ≤ 0.5 %	Chloride.....≤ 400 ppm	Assay ($HClO_4$).....≥ 99.0 % (d.s.)
Identification ..... Positive	Residue on ignition..... ≤ 0.1 %	Heavy metals (Pb).....≤10 ppm	
Spec. opt. rot. (C=5 in HCl 1N)-12.3÷-9.8°C	Ammonium .....≤ 200 ppm	Fe .....≤10 ppm	

Code	Size	Packaging	Notes
488152	25 g	Glass bottle	

Uranine ▶ Fluorescein sodium salt

Uranium standard solution

• Uranio standard soluzione • Uranium solution standard • Uranio, solución patrón • Uran-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Danger
H314
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Uranium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505922	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505923	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Uranium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
504031	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504033	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504035	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
504037	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Urea

• Urea • Urée • Urea • Urea

Synonym:
• *Carbamide*
• *Carbonyldiamide*

$(NH_2)_2CO$
Molecular Weight: 60,06
CAS: 57-13-6
EEC-N: 200-315-5

Urea > RPE - For analysis - ACS

RPE

Description ..... White crystalline powder
Identification ..... Positive
Melting point..... 132 ÷ 135 °C
Chloride..... ≤ 5 ppm
Sulphate..... ≤ 10 ppm
Water-insoluble matter ..... ≤ 100 ppm
Heavy metals (Pb)..... ≤ 10 ppm
Residue on ignition..... ≤ 100 ppm
Fe ..... ≤ 10 ppm
Assay (non-aqueos medium) .99.0 ÷ 100.5 %

Code	Size	Packaging	Notes
490757	100 g	Plastic bottle	
490758	500 g	Plastic bottle	
490759	1 kg	Plastic bottle	
490751	25 kg	Sack	

Urea > RE - Pure

RE

Description ..... White granules
Identification ..... Positive
Melting point..... 130 ÷ 135 °C
Ammoniacal salts..... ≤0.2 %
Water-insoluble matter ..... ≤200 ppm
Fe ..... ≤20 ppm

Code	Size	Packaging	Notes
387807	1 kg	Plastic bottle	
387809	5 kg	Plastic tank	
387805	25 kg	Plastic bucket	
387801	50 kg	Sack	

Product specifications are subject to changes. Please visit our website for updates.

Vanadium standard solution

• Vanadio standard soluzione • Vanadium solution standard • Vanadio, solución patrón • Vanadium-Standardlösung

Classification transport	H314
ONU: 3264	P280-P301+P330+P331-P303+P361+P353-
Transport Hazard class: LQ	P304+P340-P310a-P305+P351+P338

Vanadium standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615003300	100 ml	Plastic bottle	A 1 g/l solution Ref Ph.Eur 5003300

Vanadium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505927	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505928	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505929	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Vanadium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
504041	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504043	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504045	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
504047	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Vanadium standard solution > RS - Standard solution for AAS

RS

Description ..... Yellow clear liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
507766	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504187	500 ml	Plastic bottle	conc. 1000 ppm Matrix: Nitric acid
E497675	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Sulfuric acid
E497671	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Sulfuric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Vanadium standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Blue clear liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
491091		Plastic ampoule	conc. 1.000 ppm Matrix: Nitric acid - Volume: 50 ml

Vanadium (V) oxide

• Vanadio pentossido • Vanadium (V) pentaoxyde • Vanadio (V) pentóxido • Vanadium (V) oxid

Synonym:
Vandia

V_2O_5
Molecular Weight: 181,88
CAS: 1314-62-1
EEC-N: 215-239-8

Classification transport
ONU: 2862
Transport Hazard class: 6.1
Packing group III



Danger
H302-H332-H341-H361d-H335-H372-H411
P260-P271-P280-P304+P340-P308+P313-P403+P233

Vanadium (V) oxide > RPE - For analysis

RPE

Description ..... Ochre powder
Identification ..... Positive
Cr ..... ≤ 100 ppm
Fe ..... ≤ 300 ppm
Si ..... ≤ 200 ppm
Assay (oxidimetric) ..... ≥ 99.60 %

Code	Size	Packaging	Notes
491103	50 g	Glass bottle	

Vanillin

• Vanillina • Vanilline • Vanillina • Vanillin

Synonym:
4-Hydroxy-3-methoxybenzaldehyde

$C_8H_8O_3$
Molecular Weight: 152,15
CAS: 121-33-5
EEC-N: 204-465-2



Warning
H302
P264-P270-P301+P312a-P330-P501a

Vanillin > ERBApharm - According to pharmacopoeia: BP-DAB-NF-Ph.Eur.-FU

ERBApharm

Description ..... Yellowish crystals
Identification ..... Positive
Appearance of solution ..... Conform Ph.Eur.
Related compounds ..... Conform Ph.Eur.
React. w. sulphuric ac. ..... Conform Ph.Eur.
Melting point ..... 81 ÷ 83 ° C
Loss (silica gel) ..... ≤1.0 % (4h)
Sulphated ash ..... ≤500 ppm
Assay (alkalimetric) ..... 99.0 ÷ 101.0 % s.s.

Code	Size	Packaging	Notes
388104	100 g	Plastic bottle	
388107	1 kg	Plastic bottle	
388108	5 kg	Plastic tank	
388102	10 kg	Carton box	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Vanillin solution, phosphoric

• Vanillina soluzione, fosforico • Vaniline phosphorique solution • Vanillina solución fosfórica • Vanillinlösung, Phosphorsäure

Classification transport
ONU: 2733
Transport Hazard class: 3
Packing group II



Danger
H225-H314
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Vanillin solution, phosphoric > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611095302	100 ml	Plastic bottle	Ref Ph.Eur 1095302

Vaseline ▶ Paraffin white soft

Product specifications are subject to changes. Please visit our website for updates.

Victoria blue

• Blu vittoria • Bleu Victoria • Azul Victoria • Viktoriablau

$C_{34}H_{34}N_3Cl$
Molecular Weight: 506,08
CAS: 1325-85-5
EEC-N: 215-408-6

Classification transport
ONU: 3143
Transport Hazard class: 6.1
Packing group III

Warning
H302
P264-P270-P301+P312a-P330-P501a

Victoria blue > RS - For microscopy - C.I. 42595

RS

Description ..... Green granules Identification ..... Positive

Code	Size	Packaging	Notes
429381	10 g	Glass bottle	
429382	25 g	Glass bottle	

Dye for microscopy (Anderson Method)

Viscosity standards

• Standard di viscosità • Standards de viscosité • Patrones de viscosidad • Viskositätsstandards

Viscosity standards > RS - For calibration

RS

Code	Size	Packaging	Notes
540801	500 ml	Glass bottle	6.7cSt@20°C, 5.8cSt@25°C, 4.2cSt@37.78°C, 4cSt@40°C, 3.2cSt@50°C
540802	500 ml	Glass bottle	14cSt@20°C, 12cSt@25°C, 8cSt@37.78°C, 7.5cSt@40°C, 5.8cSt@50°C
540803	500 ml	Glass bottle	20cSt@20°C, 16cSt@25°C, 11cSt@37.78°C, 10cSt@40°C, 7.5cSt@50°C
540804	500 ml	Glass bottle	30cSt@20°C, 24cSt@25°C, 15cSt@37.78°C, 14cSt@40°C, 10cSt@50°C
540805	500 ml	Glass bottle	88cSt@20°C, 66cSt@25°C, 35cSt@37.78°C, 32cSt@40°C, 21cSt@50°C
540806	500 ml	Glass bottle	160cSt@20°C, 120cSt@25°C, 60cSt@37.78°C, 54cSt@40°C, 35cSt@50°C
540807	500 ml	Glass bottle	790cSt@20°C, 580cSt@25°C, 280cSt@37.78°C, 250cSt@40°C, 160cSt@50°C
540808	500 ml	Glass bottle	3300cSt@20°C, 2300cSt@25°C, 1100cSt@37.78°C, 940cSt@40°C, 560cSt@50°C
540809	500 ml	Glass bottle	19000cSt@20°C, 12000cSt@25°C, 4000cSt@37.78°C, 3400cSt@40°C, 1700cSt@50°C
540810	500 ml	Glass bottle	28000cSt@20°C, 17000cSt@25°C, 6000cSt@37.78°C, 5100cSt@40°C, 2500cSt@50°C
540811	500 ml	Glass bottle	41000cSt@20°C, 25000cSt@25°C, 8000cSt@37.78°C, 6700cSt@40°C, 3200cSt@50°C

Water
• Acqua • Eau • Agua • Wasser

H_2O
Molecular Weight: 18,02
CAS: 7732-18-5
EEC-N: 231-791-2

Water > RS - For UHPLC-MS

RS

Description ..... Clear colourless liquid
Colour ..... ≤ 5 APHA
Identification (I.R.) ..... Positive
Conductivity ..... ≤ 0.09 µS/cm
Residue on evaporation ..... ≤ 0.4 ppm
Total organic carbon ..... ≤ 10 ppb
Acidity ..... ≤ 0.0002 %
Alkalinity ..... ≤ 0.00005 %
Transmittance
At 200 nm ..... ≥ 95 %
At 230 nm ..... ≥ 99 %
Fluorescence (quinine)
At 254 nm ..... ≤ 0.3 ppb
At 365 nm ..... ≤ 0.3 ppb
UHPLC gradient peak
At 210 nm ..... ≤ 2 mAU
Drift at 210 nm ..... ≤ 8 mAU
Drift at 254 nm ..... ≤ 3 mAU
Sensitive Impurities (reserpine) ..... ≤ 30 ppb
Al ..... ≤ 20 ppb
Fe ..... ≤ 30 ppb
Ca ..... ≤ 50 ppb
Mg ..... ≤ 20 ppb
Na ..... ≤ 50 ppb
K ..... ≤ 50 ppb

Code	Size	Packaging	Notes
412091	1 l	Glass bottle	
412092	2.5 l	Glass bottle	

Water > RS - For LC/MS

RS

Description ..... Clear colourless liquid
Colour ..... ≤ 5 APHA
Identification (I.R.) ..... Positive
Conductivity ..... ≤ 0.09 µS/cm
Residue on evaporation ..... ≤ 0.5 ppm
Total organic carbon ..... ≤ 10 ppb
Gradiente HPLC (Test) ..... Conform
Acidity ..... ≤ 0.0002 %
Alkalinity ..... ≤ 0.00005 %
HPLC Gradient
at 210 nm ..... ≤ 2 mAU
at 254 nm ..... ≤ 0.5 mAU
Drift HPLC
Idon. test grad. LC-MS (TIC,100-2000m/z)
Impurezze sensibili (reserpina) ..... ≤ 50 ppb
Contenuto metalli
Al ..... ≤ 20 ppb
Fe ..... ≤ 30 ppb
Ca ..... ≤ 50 ppb
Mg ..... ≤ 20 ppb
Na ..... ≤ 100 ppb
K ..... ≤ 50 ppb

Code	Size	Packaging	Notes
412111	1 l	Glass bottle	
412112	2.5 l	Glass bottle	

Water > RS - For HPLC PLUS

RS

Description ..... Clear colourless liquid
Identification ..... Positive
Residue on evaporation ..... ≤0.5 ppm
Total organic carbon ..... ≤0.1 ppm
Fluorescence
at 254 nm ..... ≤1 ppb
at 365 nm ..... ≤0.5 ppb
UV Abs.max elut.peak
at 210 nm ..... ≤5 mAu
at 220 nm ..... ≤3 mAu
Conductivity during production ≤ 0.1 µS/cm
Heavy metals (Pb) ..... ≤ 0.1 ppm
Nitrate ..... ≤ 0.1 ppm
CO2 ..... Not detectable

Code	Size	Packaging	Notes
412141	1 l	Glass bottle	
412142	2.5 l	Glass bottle	

Water > RS - For Headspace chromatography

RS

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Residue on evaporation ..... ≤ 2 ppm
GC/HS
Residual solvent of class 1(acc. to ICH) ≤ 1 µg/g
Residual solvent of class 2(acc. to ICH) ≤ 10 µg/g
Residual solvent of class 3(acc. to ICH) ≤ 50 µg/g

Code	Size	Packaging	Notes
412011	1 l	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Water > RS - Ultrapure - For trace analysis at ppt level

RS

Description	Clear liquid	Co	≤ 10 ppt	Mn	≤ 10 ppt	Na	≤ 10 ppt
Colour (APHA)	≤ 10	Cu	≤ 10 ppt	Hg	≤ 20 ppt	Sr	≤ 10 ppt
Identification	Positive	Dy	≤ 1 ppt	Mo	≤ 10 ppt	Ta	≤ 10 ppt
Chloride	≤ 1 ppb	Er	≤ 1 ppt	Nd	≤ 1 ppt	Te	≤ 1 ppt
Phosphate	≤ 1 ppb	Eu	≤ 1 ppt	Ni	≤ 10 ppt	Tb	≤ 10 ppt
Sulphate	≤ 1 ppb	Gd	≤ 1 ppt	Nb	≤ 10 ppt	Tl	≤ 10 ppt
Al	≤ 20 ppt	Ga	≤ 10 ppt	Pd	≤ 10 ppt	Th	≤ 1 ppt
Sb	≤ 10 ppt	Ge	≤ 10 ppt	Pt	≤ 10 ppt	Tm	≤ 10 ppt
As	≤ 10 ppt	Au	≤ 10 ppt	K	≤ 10 ppt	Sn	≤ 10 ppt
Ba	≤ 10 ppt	Hf	≤ 1 ppt	Pr	≤ 10 ppt	Ti	≤ 10 ppt
Be	≤ 10 ppt	Ho	≤ 1 ppt	Re	≤ 10 ppt	W	≤ 10 ppt
Bi	≤ 10 ppt	In	≤ 1 ppt	Rh	≤ 10 ppt	U	≤ 1 ppt
B	≤ 20 ppt	Fe	≤ 10 ppt	Rb	≤ 10 ppt	V	≤ 10 ppt
Cd	≤ 10 ppt	La	≤ 1 ppt	Ru	≤ 10 ppt	Yb	≤ 10 ppt
Ca	≤ 10 ppt	Pb	≤ 10 ppt	Sm	≤ 10 ppt	Y	≤ 1 ppt
Ce	≤ 10 ppt	Li	≤ 10 ppt	Sc	≤ 10 ppt	Zn	≤ 10 ppt
Cs	≤ 10 ppt	Lu	≤ 1 ppt	Se	≤ 50 ppt	Zr	≤ 10 ppt
Cr	≤ 10 ppt	Mg	≤ 10 ppt	Ag	≤ 10 ppt		

Code	Size	Packaging	Notes
412185	500 ml	Plastic bottle	

Water > RS - Superpure - For trace analysis at ppb level

RS

Description	Clear liquid	Resistivity	≥ 18 Mohm	Mg	≤ 50 ppb
Identification	Positive	Total organic carbon	≤ 20 ppb	K	≤ 50 ppb
Colour	≤ 5 APHA	Al	≤ 50 ppb	Assorbimento massimo (a 254 nm)	0.002
Residue on evaporation	≤ 1 ppm	Ca	≤ 50 ppb	Absorbance	
Fluorescenza (Chinina)	≤ 100 ppt	Fe	≤ 50 ppb	at 200 nm	≤ 0.01
				at 254 nm	≤ 0.005
				at 300 nm	≤ 0.005
				at 400 nm	≤ 0.005
				LC Gradiente (Suitability Test)	Conform

Code	Size	Packaging	Notes
412151	1 l	Glass bottle	

Water > RS - For analysis according to Ph. Eur.Chapter 2.2.25

RS

Code	Size	Packaging	Notes
506411	100 ml	Glass bottle	Spectrophotometry Stray Light Blank

Water > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611095501	1 l	Plastic bottle	Water ammonium-free Ref Ph.Eur 1095501
611095506	1 l	Plastic bottle	Water nitrate-free Ref Ph.Eur 1095506

Water > RPE - For analysis

RPE

Description	Clear colourless liquid	As	≤ 0.01 ppm	Cu	≤ 0.01 ppm	Pb	≤ 0.01 ppm
pH at 25°C	5 ÷ 7	Au	≤ 0.01 ppm	Fe	≤ 0.01 ppm	Si	≤ 0.01 ppm
Conductivity	≤ 0.1 µS/cm	B	≤ 0.01 ppm	In	≤ 0.01 ppm	Sn	≤ 0.01 ppm
Residue on evaporation	≤ 1 ppm	Ba	≤ 0.01 ppm	K	≤ 0.01 ppm	Sr	≤ 0.01 ppm
Residue on calcination	≤ 1 ppm	Be	≤ 0.01 ppm	Li	≤ 0.01 ppm	Tl	≤ 0.01 ppm
Chloride	≤ 0.1 ppm	Bi	≤ 0.01 ppm	Mg	≤ 0.01 ppm	V	≤ 0.01 ppm
Phosphate	≤ 0.1 ppm	Ca	≤ 0.01 ppm	Mn	≤ 0.01 ppm	Zn	≤ 0.01 ppm
Nitrate	≤ 0.1 ppm	Cd	≤ 0.01 ppm	Mo	≤ 0.01 ppm	Zr	≤ 0.01 ppm
Sulphate	≤ 0.1 ppm	Co	≤ 0.01 ppm	Na	≤ 0.1 ppm	Oxidizing substances (O)	≤ 0.4 mg/l
Ag	≤ 0.01 ppm	Cr	≤ 0.01 ppm	Ni	≤ 0.01 ppm		

Code	Size	Packaging	Notes
307592	1 l	Plastic bottle	
307593	2.5 l	Plastic bottle	
307582	5 l	Plastic bottle	
307591	5 l	Plastic tank	
307581	10 l	Plastic tank	
307586	10 l	Kubidos	
307584	20 l	Kubidos	
307583	25 kg	Plastic tank	
307587	25 kg	Plastic tank	With tap
307585	50 kg	Plastic drum	
307589	1000 kg	Combined drum	

Conform to EN-ISO 3696 grade 3

Water purified

• Acqua purificata • Eau purifiee • Agua purificada • Gereinigtes Wasser

H_2O
Molecular Weight: 18,02
CAS: 7732-18-5
EEC-N: 231-791-2

Water purified > ERBApharm® - According to pharmacopoeia: Ph.Eur.-FU-Ph.Franc.-BP-DAB-USP-JP

ERBApharm

Description	Clear colourless liquid	Chloride	Pass test	Residue on evaporation	≤10 ppm	Origin (BSE/TSE)	Mineral
Identification	Positive	Sulphate	Pass test	Ammonium	≤0.2 ppm	Total Organic Carbon	≤ 0.05 ppm
pH	5.0 ÷ 7.0	Calcium	Conform	Nitrate	≤0.2 ppm	Conductivity at 20°C	≤ 1.1 µS/cm
Acidity or alkalinity	Conform Ph.Eu.	Calcium + Magnesium	Conform Ph.Eu.	Heavy metals (Pb)	≤0.1 ppm		
Oxidizable substances	Conform Ph.Eu.	Carbon dioxide	Conform	Total aerobic bacteria	Conform Ph.Eu.		

Code	Size	Packaging	Notes
307606	1 l	Plastic bottle	
307601	5 l	Plastic tank	
307602	10 l	Kubidos	
307603	25 kg	Plastic tank	
307604	200 l	Plastic drum	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Water + 0.1% v/v formic acid

• Acqua + 0.1%(v/v) acido formico • Eau + 0.1%(v/v) acide formique • Agua + 0.1%(v/v) acido fórmico • Wasser + 0.1% (v / v) Ameisensäure

H_2O
Molecular Weight: 18,02
CAS: 7732-18-5

Water + 0.1% v/v formic acid > RS - For LC/MS

RS

Description ..... Clear colourless liquid
Colour ..... ≤ 10 APHA
Acidity (formic acid) ..... 0.095 ÷ 0.105 %
HPLC Gradient
At 210 nm ..... ≤ 50 mAU
At 254 nm ..... ≤ 10 mAU

Transmittance
At 230 nm ..... ≥ 45 %
pH at 20°C ..... 2.6 - 2.8
Test LC-MS TIC (100-2000m/z)
Sensitive Impurities (reserpine) ..... ≤ 50 ppb
Metals content

Al ..... ≤ 20 ppb
Fe ..... ≤ 30 ppb
Ca ..... ≤ 50 ppb
Mg ..... ≤ 20 ppb
Na ..... ≤ 100 ppb
K ..... ≤ 50 ppb

Raw material used
Water (code 412110) ..... Batch number
Formic acid 98-99% (code 405820) . Batch number

Code	Size	Packaging	Notes
412121	1 l	Glass bottle	
412122	2.5 l	Glass bottle	

Water + 0.1% v/v trifluoroacetic acid

• Acqua+ 0.1% v/v acido trifluoroacetico • Eau + 0.1% v/v d'acide trifluoroacétique • Agua + 0,1% v/v acido trifluoroacético • Wasser + 0.1% v/v trifluoressigsäure

H_2O
CAS: 7732-18-5

Water + 0.1% v/v trifluoroacetic acid > RS - For HPLC

RS

Clear,colourless solution ..... Conform
TFA content ..... 0.095 - 0.105 %

Code	Size	Packaging	Notes
412031	1 l	Glass bottle	

Water chlorine

• Acqua di cloro • Eau de chlore • Agua de cloro • Chlorwasser

Cl_2
Molecular Weight: 70,91
CAS: 7782-50-5



Warning
H400
P273-P391-P501a

Water chlorine > RPE - For analysis

RPE

Description ..... Yellow clear liquid
Identification ..... Positive
Assay (ex chloryne) ..... 0.4 ÷ 0.7 % (p/p)

Code	Size	Packaging	Notes
411981	1 l	Glass bottle	

Water deionized and acidified

• Acqua deionizzata acidificata • Eau désionisée acidifiée • Agua desionizada acidificada • Gesäuertes entionisiertes Wasser

Molecular Weight: 63,01

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III



Warning
H290-H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Water deionized and acidified > RS - Blanks for AAS, ICP, ICP-MS

RS

Code	Size	Packaging	Notes
504550	1 l	Bottle	Matrix: 2 % Nitric acid
504551	1 l	Bottle	Matrix: 5 % Nitric acid
504552	1 l	Bottle	Matrix: 10 % Nitric acid
504553	1 l	Bottle	Matrix: 2 % Hydrochloric acid
504554	1 l	Bottle	Matrix: 5 % Hydrochloric acid
504557	1 l	Bottle	Matrix: 10 % Hydrochloric acid

Wijs' reagent

• Wijs reattivo • Réactif de Wijs • Wijs reactivo • Wijs 'Reagenz

Classification transport
ONU: 2920
Transport Hazard class: 8
Packing group II



Danger
H226-H314-H373
P210-P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Wijs' reagent > RS - For detection of iodine index

RS

Description ..... Brown clear liquid Identification ..... Positive Iodine value .....≥90

Code	Size	Packaging	Notes
E491901	250 ml	Glass bottle	
E491902	1 l	Glass bottle	

Wool fat ▶ Lanolin anhydrous

Wright's stain solution in methanol

• Wright colorante soluzione in alcole metilico • Colorant de Wright solution methanolique • Colorante de Wright solución in alcohol metilico • Wright's Färbelösung in Methanol

Classification transport
ONU: 1992
Transport Hazard class: 3
Packing group II



Danger
H225-H301-H370
P210-P241-P280-P301+P310a-P303+P361+P353-P403+P235

Wright's stain solution in methanol > RS - For microscopy

RS

Description ..... Blue clear liquid Identification ..... Positive Assorbanza a 518 nm ..... ≥ 0.17 Assorbanza a 660 nm ..... ≥ 0.15

Code	Size	Packaging	Notes
492011	100 ml	Glass bottle	In Vitro Diagnostic Medical Device

Dye for hematology

Product specifications are subject to changes. Please visit our website for updates.

Xylene, mix of isomers

• Xilene, isomeri misti • Xylène, mélange d'isomères • Xileno, mezcla de isómeros • Xylol, isomerengemisch

C_8H_{10}
Molecular Weight: 106,17
EEC-N: 905-588-0

Classification transport
ONU: 1307
Transport Hazard class: 3
Packing group III

Danger
H226-H312-H332-H315-H319-H335-H373-H304
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

Xylene, mix of isomers > RS - RSE - For electronic use

RS

Description	Clear liquid	Phosphate	≤1 ppm	Cd	≤0.005 ppm	Ni	≤0.01 ppm
Colour (APHA)	≤10	Heavy metals (Pb)	≤0.1 ppm	Co	≤0.01 ppm	Pb	≤0.01 ppm
Identification	Positive	Toluene	≤5000 ppm	Cr	≤0.01 ppm	Pt	≤0.05 ppm
Ready carbonizable substances	Conform	Total sulphur	≤3 ppm	Cu	≤0.01 ppm	Sb	≤0.01 ppm
Density at 20° C	0.864 ÷ 0.870	Ag	≤0.02 ppm	Fe	≤0.05 ppm	Sn	≤0.02 ppm
Boiling point	137.0 ÷ 140.0 °	Al	≤0.05 ppm	Ga	≤0.02 ppm	Sr	≤0.02 ppm
Resistivity	≥1 Mohm.cm	As	≤0.01 ppm	In	≤0.02 ppm	Ti	≤0.05 ppm
Water (K.F.)	≤100 ppm	Au	≤0.05 ppm	K	≤0.1 ppm	Tl	≤0.05 ppm
Residue on evaporation	≤5 ppm	B	≤0.01 ppm	Li	≤0.02 ppm	V	≤0.05 ppm
Acidity (HCl)	≤5 ppm	Ba	≤0.1 ppm	Mg	≤0.1 ppm	Zn	≤0.01 ppm
Alcalinity (NH3)	≤1 ppm	Be	≤0.02 ppm	Mn	≤0.01 ppm	Zr	≤0.05 ppm
Benzene	≤100 ppm	Bi	≤0.02 ppm	Mo	≤0.05 ppm		
Chloride	≤3 ppm	Ca	≤0.2 ppm	Na	≤0.1 ppm		

Code	Size	Packaging	Notes
492358	1 l	Glass bottle	
492359	2.5 l	Glass bottle	

Xylene, mix of isomers > RPE - For analysis - ISO - ACS - Reag.Ph.Eur. - Reag.USP

RPE

Description	Clear liquid	Residue on evaporation	≤10 ppm	Al	≤0.5 ppm	Mg	≤0.1 ppm
Colour (APHA)	≤10	Acidity (benzoic acid)	≤14 ppm	Ba	≤0.1 ppm	Mn	≤0.02 ppm
Identification (I.R.)	Positive	Alcalinity (NH3)	≤2 ppm	Ca	≤0.5 ppm	Ni	≤0.02 ppm
Boiling point	137 - 140 °C	Benzene	≤ 50 ppm	Cd	≤0.05 ppm	Pb	≤0.05 ppm
Ready carbonizable substances	Conform	Ethylbenzene	≤25 %	Co	≤0.02 ppm	Sn	≤0.1 ppm
Density at 20° C	0.864 ÷ 0.870	Tiophene	≤1 ppm	Cr	≤0.02 ppm	Zn	≤0.05 ppm
Refractive index at 20°C	1.4947 ÷ 1.4987	Toluene	≤0.5 %	Cu	≤0.02 ppm	Assay(isomeric mixture)	≥99.0 %
Water (K.F.)	≤200 ppm	Total sulphur	≤3 ppm	Fe	≤0.1 ppm	Any single impurity	≤ 0.5 %

Code	Size	Packaging	Notes
492301	1 l	Glass bottle	
492306	2.5 l	Glass bottle	
492305	5 l	Plastic tank	
492303	23 kg	Metal drum	
492304	170 kg	Metal drum	

Xylene, mix of isomers > RE - Pure - Low content in benzene

RE

Description	Clear liquid	Refractive index at 20°C	1.4917 ÷ 1.5017	Water (K.F.)	≤ 200 ppm	Assay (isomeric mix)	≥ 98.5 %
Identity (IR)	Positive	Boiling point	137.5 ÷ 139.5 °C	Benzene	≤ 50 ppm		
Density at 20°C	0.862 ÷ 0.872	Residue on evaporation	≤ 100 ppm	Total sulphur	≤ 100 ppm		

Code	Size	Packaging	Notes
392602	1 l	Glass bottle	
392603	2.5 l	Glass bottle	
528251	5 l	Plastic tank	
392605	23 kg	Metal drum	
392608	170 kg	Metal drum	
528252	200 l	Metal drum	

o-Xylene

• o-Xilene • o-xylène • o-Xileno • o-Xylol

Synonym:
1,2-Dimethylbenzene

$C_6H_4(CH_3)_2$
Molecular Weight: 106,17
CAS: 95-47-6
EEC-N: 202-422-2

Classification transport
ONU: 1307
Transport Hazard class: 3
Packing group II



Danger
H226-H312-H332-H315-H319-H335-H304
P210-P280-P303+P361+P353-P304+P340-P305+P351+P338-P403+P233

o-Xylene > RPE - For analysis

RPE

Description	Clear colourless liquid	Ready carbonizable substances	Conform	Water (K.F.)	≤300 ppm	Tiophene	≤1 ppm
Identification	Positive	Density at 20° C	0.875 ÷ 0.885	Acidity (benzoic acid)	≤14 ppm	Toluene	≤0.15 %
Alcohol miscibility	Complete	Refractive index at 20°C	1.5028 ÷ 1.5088	Alcalinity (NH3)	≤2 ppm	Total sulphur	≤5 ppm
Chloroform miscibility	Complete	Boiling point	139.0 ÷ 149.0 ° C	Benzene	≤ 500 ppm	Assay (GLC)	≥99 %
Diethyl ether miscib.	Complete	Freezing point	-23.5 ÷ -26.5 °C	Ethylbenzene	≤0.15 %		

Code	Size	Packaging	Notes
492403	1 l	Glass bottle	
492404	2.5 l	Glass bottle	
492401	24 kg	Metal drum	

Xylenecyanol

• Xilencianolo FF • Xylènecyanol FF • Xilenocianol FF • Xylolcyanol

Synonym:
• *Xylene Cyanol FF*
• *Acid blue 147*

$C_{25}H_{27}N_2NaO_6S_2$
Molecular Weight: 538,61
CAS: 2650-17-1
EEC-N: 220-167-5



Warning
H319
P264-P280i-P305+P351+P338-P337+P313

Xylenecyanol > RPE - For analysis

RPE

Description	Cristalli verde intenso	Identification	Positive	Dye content	≥ 70 %	UV Lambda max	612 to 616 nm

Code	Size	Packaging	Notes
492211	1 g	Glass bottle	
492212	25 g	Glass bottle	

Dye for microscopy (histology). Acid base indicator. C.I. 42135

2,4-Xylenol

• 2,4-Xilenolo • 2,4-Xylenol • 2,4-Xilenol • 2,4-Xylenol

Synonym:
2,4-Dimethylphenol | 4-Hydroxy-m-xylene

$(CH_3)_2C_6H_3OH$
Molecular Weight: 122,17
CAS: 105-67-9
EEC-N: 203-321-6

Classification transport
ONU: 3430
Transport Hazard class: 6.1
Packing group II



Danger
H301-H311-H314-H411
P280-P301+P310a-P301+P330+P331-P303+P361+P353-P304+P340-P305+P351+P338

2,4-Xylenol > RE - Pure

RE

Description	Yellow clear liquid	Density at 20° C	1.015 ÷ 1.021	Residue on ignition	≤100 ppm
Identification	Positive	Boiling point	210.5 ÷ 212.5 °C	Assay (GC)	≥ 98.0 %

Code	Size	Packaging	Notes
492661	25 ml	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Xylenol orange

• Arancio xilenolo • Orange de xylénol • Naranja de xilenol • Xylenolorange

Synonym:
• *Xylenol Orange disodium salt*
• *Cresolsulfonphtalein disodium salt*

$C_{31}H_{30}O_{13}N_2SNa_2$
Molecular Weight: 716,63
CAS: 1611-35-4
EEC-N: 216-553-8

Warning
H315-H319-H335
P261-P271-P304+P340-P305+P351+P338-P332+P313-P403+P233

Xylenol orange > RPE - For analysis

RPE

Description ..... Orange crystalline powder Identification ..... Positive Sens.as complex.indicat ..... Conform

Code	Size	Packaging	Notes
423597	1 g	Glass bottle	
423598	5 g	Glass bottle	
423599	25 g	Glass bottle	

Complexometric indicator

D(+)Xylose

• D(+)Xilosio • D(+)Xylose • D(+)Xilosa • D(+)-Xylose

$C_5H_{10}O_5$
Molecular Weight: 150,13
CAS: 58-86-6
EEC-N: 200-400-7

D(+)Xylose > RPE - For analysis

RPE

Description ..... White crystalline powder
Identification ..... Positive
Specific optical rotation..... +19 ÷ +21 °
Water (K.F.) ..... ≤0.2 %
Chloride..... ≤50 ppm
Heavy metals (Pb)..... ≤10 ppm
Residue on ignition..... ≤0.1 %
Sulphate..... ≤50 ppm

Code	Size	Packaging	Notes
492803	50 g	Glass bottle	
492804	100 g	Glass bottle	

D(+)Xylose > RE - Pure

RE

Description .....Whyte crystalline powder
Identification ..... Positive
Melting point..... 144 ÷ 148 ° C
Specific optical rotation... +19.0 ÷ +20.0 °
Loss on drying ..... ≤0.3 %
Acidity (acetic acid)..... ≤300 ppm
Chloride..... ≤50 ppm
Heavy metals (Pb)..... ≤10 ppm
Residue on ignition..... ≤0.1 %
Sulphate..... ≤50 ppm
As..... ≤1 ppm

Code	Size	Packaging	Notes
392631	500 g	Plastic bottle	
392635	25 kg	Plastic bucket	

Yeast dried
• Lievito di birra secco • Levure de bière sèche • Levdura de cerveza seco • Bierhefe trocken

Yeast dried > RE - Pure

RE

Description ..... Pv-gran nocciola odore ca Identification ..... Positive

Code	Size	Packaging	Notes
348854	100 g	Plastic bottle	

Ytterbium standard solution
• Itterbio standard soluzione • Ytterbium solution standard • Iterbio, solución patrón • Ytterbium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Ytterbium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505947	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505948	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Ytterbium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
504071	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504073	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504075	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
504077	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Ytterbium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507768	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507523	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Yttrium standard solution
• Ittrio standard soluzione • Yttrium solution standard • Itrio, solución patrón • Yttrium-Standardlösung

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Yttrium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505942	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505945	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505943	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Product specifications are subject to changes. Please visit our website for updates.

Yttrium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
504061	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504063	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
504065	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid
504067	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Yttrium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507767	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
507522	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Ziehl-Neelsen's reagent

• Ziehl-Neelsen reattivo soluzione • Réactif de Ziehl-Neelsen en solution • Ziehl-Neelsen reactivo solución • Ziehl-Neelsens Reagenz

Classification transport
ONU: 2810
Transport Hazard class: 6.1
Packing group II



Danger
H331-H314-H341
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P308+P313-P403+P233

Ziehl-Neelsen's reagent > RS - For bacteriology

RS

Description ..... Red liquid Identification ..... Positive Assorbanza a 554 nm ..... ≥ 0.14

Code	Size	Packaging	Notes
493101	250 ml	Plastic bottle	In Vitro Diagnostic Medical Device
493102	1 l	Plastic bottle	In Vitro Diagnostic Medical Device

Dye for bacteriology. Contains carbol fuchsin.

Zinc, foil

• Zinco, lastra • Zinc, lames • Zinc, hojas • Zink, Folie

Zn
Molecular Weight: 65,38
CAS: 7440-66-6
EEC-N: 231-175-3

Zinc, foil > RPE - For analysis

RPE

Description ..... Grey foil Identification ..... Positive Assay (oxidimetric) ..... ≥98 %

Code	Size	Packaging	Notes
493507	1 kg	Bag	

Zinc, granular

• Zinco, granuli • Zinc, granule • Zinc, granulada • Zink, körnig

Zn
Molecular Weight: 65,38
CAS: 7440-66-6
EEC-N: 231-175-3

Classification transport
ONU: 1436
Transport Hazard class: 4.3
Packing group II



Danger
H250-H260-H410
P210-P222-P223-P231a+P232-P280-P402+P404

Zinc, granular > RPE - For analysis

RPE

Description ..... Granuli grigi As ..... ≤0.15 ppm Pb ..... ≤500 ppm
Identification ..... Positive Fe ..... ≤200 ppm

Code	Size	Packaging	Notes
493451	500 g	Plastic bottle	0.5 - 1 mm
493307	1 kg	Plastic bottle	1 - 7 mm
493309	5 kg	Plastic bottle	1 - 7 mm
493303	25 kg	Glass bottle	1 - 7 mm

Product specifications are subject to changes. Please visit our website for updates.

Zinc, powder
• Zinco, polvere • Zinc, poudre • Zinc, polvo • Zink, Pulver

Zn
Molecular Weight: 65,38
CAS: 7440-66-6

Classification transport
ONU: 3077
Transport Hazard class: 9
Packing group III

Warning
H410
P273-P391-P501a

Zinc, powder > RE - Pure — RE

Description ..... Grey powder Identification ..... Positive Assay (oxidimetric) ..... ≥85 %

Code	Size	Packaging	Notes
493705	250 g	Glass bottle	
493707	1 kg	Plastic bottle	
493702	25 kg	Metal drum	

Zinc, activated
• Zinco, attivato • Zinc activé • Zinc, activado • Zink aktiviert

Zn
CAS: 7440-66-6

Classification transport
ONU: 3077
Transport Hazard class: 9
Packing group III

Warning
H410
P273-P391-P501a

Zinc, activated > RS - For analysis according to Ph. Eur. Chap. 4.1.1 — RS

Code	Size	Packaging	Notes
611096501	100 g	Plastic bottle	Ref Ph.Eur 1096501

Zinc standard solution
• Zinco standard soluzione • Zinc solution standard • Zinc, solución patrón • Zink-Standardlösung

Zn
Molecular Weight: 65,38
CAS: 7440-66-6

Classification transport
ONU: 3264
Transport Hazard class: 8
Packing group III

Warning
H290
P234-P390-P406

Zinc standard solution > RS - For analysis according to Ph. Eur. Chap. 4.2.1 — RS

Code	Size	Packaging	Notes
612000800	100 g	Plastic bottle	Ref Ph.Eur 2000800

Zinc standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2 — RS

Code	Size	Packaging	Notes
615003402	100 ml	Plastic bottle	A 10 ppm solution: to dilute according to Ref Ph.Eur 5003402
615003403	100 ml	Plastic bottle	A 5 ppm solution: to dilute according to Ref Ph.Eur 5003403
615003409	100 ml	Plastic bottle	A 100 ppm solution: to dilute according to Ref Ph.Eur 5003401

Zinc standard solution > RS - Standard solution for ICP-MS — RS

Code	Size	Packaging	Notes
505952	100 ml	Plastic bottle	conc. 10 ppm. Matrix: Nitric acid
505955	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid
505953	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Zinc standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
504081	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
504085	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid
504083	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
504087	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Zinc standard solution > RS - Standard solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
507769	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
E497685	100 ml	Glass bottle	conc. 1.000 ppm Matrix: Hydrochloric acid
507477	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid
E497681	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Zinc standard solution > RS - NORMEX - Concentrated solution for AAS

RS

Description ..... Clear colourless liquid Identification ..... Positive Titration factor..... 0.995 ÷ 1.005

Code	Size	Packaging	Notes
493151		Plastic ampoule	conc. 1.000 ppm Matrix: Nitric acid - Volume: 50 ml

Volumetric concentrated solution to be diluted accordingly with the instructions for use printed on the package

Zinc acetate dihydrate

• Zinco acetato diidrato • Zinc acétate dihydraté • Zinc acetato dihidrato • Zinkacetatdihydrat

$Zn(CH_3COO)_2.2H_2O$
Molecular Weight: 219,49
CAS: 5970-45-6

Warning
H319-H410
P264-P280i-P305+P351+P338-P337+P313-P391-P501a

Zinc acetate dihydrate > RPE - For analysis

RPE

Description ..... White shining crystals
Identification ..... Positive
pH sol. 5% at 25° C ..... 6.2 ÷ 6.6
Chloride ..... ≤5 ppm
Phosphate ..... ≤20 ppm
Water-insoluble matter ..... ≤30 ppm
Sulphate ..... ≤10 ppm
As ..... ≤0.4 ppm
Cd ..... ≤5 ppm
Cu ..... ≤5 ppm
Fe ..... ≤ 5 ppm
K ..... ≤100 ppm
Mn ..... ≤5 ppm
Ni ..... ≤5 ppm
Pb ..... ≤5 ppm
Assay (complexometric) ..... ≥99.5 %

Code	Size	Packaging	Notes
493806	100 g	Plastic bottle	
493807	1 kg	Plastic bottle	
493803	25 kg	Plastic bucket	

DASIT GROUP | CARLO ERBA REAGENTS
Product specifications are subject to changes. Please visit our website for updates.

Zinc carbonate basic

• Zinco carbonato basico • Zinc carbonate basique • Zinc carbonato básico • Zinkcarbonat basisch

$Zn_5(CO_3)_2(OH)_6$
Molecular Weight: 548,97
CAS: 5263-02-5
EEC-N: 226-076-7



Warning
H315-H319-H410
P264-P280a-P305+P351+P338-P332+P313-P362+P364-P337+P313

Zinc carbonate basic > RPE - For analysis

RPE

Description	White powder	Cd	≤ 10 ppm	Mn	≤ 20 ppm	Loss on drying	≤ 3 %
Identification	Positive	Cu	≤ 5 ppm	Ni	≤ 5 ppm	Loss on ignition	25.5 ± 0.5 %
Sulphate	≤ 0.5 %	Fe	≤ 20 ppm	Pb	≤ 5 ppm	Assay (alkalimetric)	≥ 69 % (ZnO)

Code	Size	Packaging	Notes
494006	500 g	Plastic bottle	

Zinc chloride anhydrous

• Zinco cloruro anidro • Zinc chlorure anhydre • Zinc cloruro anhidro • Zinkchlorid wasserfrei

Synonym:
Dichlorozinc

$ZnCl_2$
Molecular Weight: 136,28
CAS: 7646-85-7
EEC-N: 231-592-0

Classification transport
ONU: 2331
Transport Hazard class: 8
Packing group III



Danger
H302-H314-H410
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338

Zinc chloride anhydrous > RPE - For analysis - ACS

RPE

Description	White crystals	Nitrate	≤30 ppm	Fe	≤10 ppm	Pb	≤50 ppm
Identification	Positive	Oxichloride(ZnO)	Conform	K	≤200 ppm	Assay (argentimetric)	≥97.0 %
Ammonium	≤50 ppm	Sulphate	≤100 ppm	Mg	≤100 ppm		
Diluted HCl-ins. matter	≤50 ppm	Ca	≤600 ppm	Na	≤500 ppm		

Code	Size	Packaging	Notes
494104	100 g	Plastic bottle	
494105	250 g	Plastic bottle	
494107	1 kg	Plastic bottle	
494106	10 kg	Plastic tank	

Zinc chloride anhydrous > RE - Pure

RE

Description	White crystalline powder	Sulphate	≤0.05 %	Assay (complexometric)	≥97 %
Identification	Positive	Fe	≤10 ppm		

Code	Size	Packaging	Notes
393007	1 kg	Plastic bottle	
393009	5 kg	Plastic tank	

Zinc chloride solution 60%

• Zinco cloruro soluzione 60% • Zinc chlorure solution à 60% • Zinc cloruro solución 60%
• Zinkchloridlösung 60%

Synonym:
Dichlorozinc

$ZnCl_2$
Molecular Weight: 136,28
CAS: 7646-85-7

Classification transport
ONU: 1840
Transport Hazard class: 8
Packing group III



Danger
H302-H314-H335-H410
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Zinc chloride solution 60% > RPE - For analysis

RPE

Description	Liquido incolore	Identification	Positive	Density at 20° C	≥1.750

Code	Size	Packaging	Notes
E494301	1 l	Bottle	

Zinc chloride solution, iodinated

• Zinco cloruro soluzione, iodata • Chlorure de zinc - Solution iodée • Zinc cloruro solución, yodada
• Zinkchlorid - Jodlösung

Synonym:
Dichlorozinc

Molecular Weight: 136,28

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group III



Danger
H302-H314-H335-H372-H410
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Zinc chloride solution, iodinated > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611096603	500 ml	Glass bottle	Ref Ph.Eur 1096602

Storage: protected from light

Zinc chloride-formic acid solution

• Zinco cloruro - acido formico soluzione • Chlorure de zinc - acide formique - Solution • Zinc cloruro - ácido fórmico solución
• Zinkchlorid-Ameisensäure - Lösung

Classification transport
ONU: 1760
Transport Hazard class: 8
Packing group II



Danger
H302-H332-H314-H335-H411
P280-P301+P330+P331-P303+P361+P353-P304+P340-P310a-P305+P351+P338-P403+P233

Zinc chloride-formic acid solution > RS - For analysis according to Ph. Eur. Chap. 4.1.1

RS

Code	Size	Packaging	Notes
611096601	1 l	Glass bottle	Ref Ph.Eur 1096601

Zinc dibenzyldithiocarbamate

• Zinco dibenzilditiocarbammato • Zinc dibenzyldithiocarbamate • Zinc dibencilditiocarbamato • n,n-Zinkdibenzyldithiocarbamat

$[(C_6H_5CH_2)_2NCSS]_2Zn$
Molecular Weight: 610,2
CAS: 14726-36-4
EEC-N: 238-778-0

Zinc dibenzyldithiocarbamate > RPE - For analysis

RPE

Description ..... White powder
Identification ..... Positive
Melting point ..... 183 ÷ 185 ° C
Assay (complexometric) ..... ≥94 %

Code	Size	Packaging	Notes
494311	10 g	Glass bottle	

Zinc nitrate hexahydrate

• Zinco nitrato esaidrato • Zinc nitrate hexahydraté • Zinc nitrato hexahidratado • Zinknitrathexahydrat

$Zn(NO_3)_2.6H_2O$
Molecular Weight: 297,47
CAS: 10196-18-6
EEC-N: 231-943-8

Classification transport
ONU: 1514
Transport Hazard class: 5.1
Packing group II



Danger
H272-H302-H315-H319-H335
P210-P261-P280-P304+P340-P305+P351+P338-P403+P233

Zinc nitrate hexahydrate > RPE - For analysis

RPE

Description ..... White semitransparent crystals
Identification ..... Positive
Chloride ..... ≤ 100 ppm
Fe ..... ≤ 20 ppm
Pb ..... ≤ 100 ppm
Assay (complexometric) ..... ≥ 97.5 %

Code	Size	Packaging	Notes
494506	100 g	Plastic bottle	
494507	1 kg	Plastic bottle	

Product specifications are subject to changes. Please visit our website for updates.

Zinc oxide

• Zinco ossido • Zinc oxyde • Zinc óxido • Zinkoxid

ZnO
Molecular Weight: 81,37
CAS: 1314-13-2
EEC-N: 215-222-5

Warning
H410
P273-P391-P501a

Zinc oxide > RPE - For analysis

RPE

Description	White powder	Phosphate	≤5 ppm	As	≤0.5 ppm	Ni	≤10 ppm
Identification	Positive	Dil. H2SO4-ins. matter	≤100 ppm	Ca	≤50 ppm	Pb	≤50 ppm
Alcalinity	Conform	Nitrate	≤20 ppm	Cd	≤10 ppm	Zn	≤20 ppm
Loss on ignition	≤0.5 %	Subst. not ppt. (NH4)2S	≤0.1 %	Cu	≤5 ppm	Assay (alkalimetric)	≥99.0 %
Carbonate	≤0.2 %	Subst. reducing KMnO4	≤10 ppm(15m)	Fe	≤10 ppm		
Chloride	≤5 ppm	Total sulphur	≤50 ppm	Mn	≤5 ppm		

Code	Size	Packaging	Notes
494606	100 g	Plastic bottle	
494607	1 kg	Plastic bottle	
494602	25 kg	Plastic bucket	

Zinc oxide > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP

ERBApharm

Description	Yellowish powder	Fe - other heavy metals	Conform USP-NF	Fe	≤200 ppm	Origin (BSE/TSE)	Synthesis
Identification	Positive	Loss on calcin. 500°C	≤1.0 %	Pb	≤50 ppm	Residual solvents (Current ICH)	Conform
Alcalinity	Conform Ph.Eur.	As	≤5 ppm	Assay (complexometric) s.s.c.	99.0 ÷ 100.5 %		
Coal and acid ins.matt.	Conform Ph.Eur.	Cd	≤10 ppm				

Code	Size	Packaging	Notes
393507	1 kg	Plastic bottle	
393509	5 kg	Plastic tank	
393503	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Zinc stearate

• Zinco stearato • Zinc stéarate • Zinc estearato • Zinkstearat

Synonym:
Stearic acid zinc salt

$C_{36}H_{70}O_4Zn$
Molecular Weight: 632,33
CAS: 557-05-1
EEC-N: 209-151-9

Warning
H335
P261-P271-P304+P340-P312a-P403+P233-P501a

Zinc stearate > ERBApharm - Vegetal origin - According to pharmacopoeia: Ph.Eur.-USP-FU

ERBApharm

Description	White powder	Acidity ind. fat acids	195 ÷ 210	As	≤ 1.5 ppm	Assay as ZnO (complexometric)	12.5 ÷ 14.0 %
Identification	Positive	Freezing point	≥54 °C	Cd	≤ 5.0 ppm		
Appearance of solution	Conform Ph.Eur.	Alkaly-alkaline earth	≤ 1.0 %	Pb	≤ 10 ppm		
Solution appea. fat ac.	Conform Ph.Eur.	Chloride	≤ 250 ppm	Assay as Zn (complexometric)	10.0 ÷ 12.0 %		
Acidity or alkalinity	Conform Ph.Eur.	Sulphate	≤ 0.6 %				

Code	Size	Packaging	Notes
395451	1 kg	Plastic bottle	
395452	10 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Zinc sulfate heptahydrate

• Zinco solfato eptaidrato • Zinc sulfate heptahydraté • Zinc sulfato heptahidratado • Zink-Sulfat heptahydrat

$ZnSO_4.7H_2O$
Molecular Weight: 287,54
CAS: 7446-20-0



Danger
H302-H318-H410
P264-P280i-P301+P312a-P305+P351+P338-P310a-P501a

Zinc sulfate heptahydrate > RPE - For analysis - ACS

RPE

Description	White crystals	Chloride	≤5 ppm	Fe	≤10 ppm	Na	≤500 ppm
Identification	Positive	Water-insoluble matter	≤100 ppm	K	≤100 ppm	Pb	≤30 ppm
pH sol. 5% at 25° C	4.4 ÷ 6.0	Nitrate	≤20 ppm	Mg	≤50 ppm	Assay (complexometric)	99.0 ÷ 103.0 %
Ammonium	≤10 ppm	Ca	≤50 ppm	Mn	≤3 ppm		

Code	Size	Packaging	Notes
494905	100 g	Plastic bottle	
494906	500 g	Plastic bottle	
494907	1 kg	Plastic bottle	
494909	5 kg	Plastic jar	
494901	25 kg	Drum	

Zinc sulfate heptahydrate > ERBApharm - According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP

ERBApharm

Description	White crystalline powder	Acidity	Conform USP-NF	Alkaline,alk.earth met.	≤0.9 %	Pb	≤20 ppm
Identification	Positive	pH sol. 5% at 25° C	4.4 ÷ 5.6	As	≤14 ppm	Assay (complexometric)	99.0 ÷ 104.0 %
Appearance of solution	Conform Ph.Eur.	Chloride	≤300 ppm	Fe	≤100 ppm		

Code	Size	Packaging	Notes
394007	1 kg	Plastic bottle	
394009	5 kg	Plastic tank	
394001	25 kg	Plastic bucket	

In case of Excipient use, please contact us to inquire about the availability of this product in XcipharmTM grade

Zinc sulfate monohydrate

• Zinco solfato monoidrato • Zinc sulfate monohydrate • Zinc sulfato monohidrato • Zinksulfat-Monohydrat

$ZnSO_4.H_2O$
Molecular Weight: 179,45
CAS: 7446-19-7
EEC-N: 231-793-3

Classification transport
ONU: 3077
Transport Hazard class: 9
Packing group III



Danger
H302-H318-H410
P264-P280i-P301+P312a-P305+P351+P338-P310a-P501a

Zinc sulfate monohydrate > RPE - For analysis

RPE

Description	White fine powder	pH sol. 5% at 25°C	4 ÷ 6	Pb	≤ 10 ppm	As	≤ 1 ppm
Identification	Positive	Chloride	≤ 1 %	Cu	≤ 5 ppm	Hg	≤ 0.1 ppm
Assay (Zn)	≥ 35 %	Mn	≤ 30 ppm	Cd	≤ 10 ppm		
Water (K.F.)	≤ 1 %	Fe	≤ 30 ppm	Ni	≤ 15 ppm		

Code	Size	Packaging	Notes
495005	250 g	Plastic bottle	
495007	1 kg	Plastic bottle	

Zinc sulfate 0.1 mol/l (0.2N)

• Zinco solfato 0.1 mol/l (0.2N) • Zinc sulfate 0.1 mol/l (0.2N) • Zinc sulfato 0.1 mol/l (0.2N) • Zink-Sulfat 0.1 mol/l (0.2N)

$ZnSO_4.7H_2O$
Molecular Weight: 287,54
CAS: 7446-20-0

Zinc sulfate 0.1 mol/l (0.2N) > RS - For analysis according to Ph. Eur. Chap. 4.2.2

RS

Code	Size	Packaging	Notes
613008601	500 ml	Plastic bottle	Ref Ph.Eur 3008600
613008600	1 l	Plastic bottle	Ref Ph.Eur 3008600

Product specifications are subject to changes. Please visit our website for updates.

Zinc sulfate 0.1 mol/l (0.2N) > RS - For analysis according to JP

RS

Code	Size	Packaging	Notes
616000050	1 l	Plastic bottle	Zinc (standard reagent)/ hydrochloric acid / bromine TS /water

Zinc sulfate 0.1 mol/l (0.2N) > RS - For analysis according to USP

RS

Code	Size	Packaging	Notes
617000291	1 l	Bottle	

Zinc sulfate 0.1 mol/l (0.2N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Assay (potentiometry) ..... 0.1996 - 0.2004 N

Code	Size	Packaging	Notes
494921	1 l	Plastic bottle	

Zinc sulfate 0.05 mol/l (0.05N)

• Zinco solfato 0.05 mol/l (0.05N) • Zinc sulfate 0.05 mol/l (0.05N) • Zinc sulfato 0.05 mol/l (0.05N) • Zink-Sulfat 0.05 mol/l

$ZnSO_4.7H_2O$
Molecular Weight: 287,54
CAS: 7446-20-0

Warning
H319-H412
P264-P273-P280i-P305+P351+P338-P337+P313-P501a

Zinc sulfate 0.05 mol/l (0.05N) > RPE - For analysis

RPE

Description ..... Clear colourless liquid Titration factor ..... 0.998 ÷ 1.002

Code	Size	Packaging	Notes
494931	1 l	Plastic bottle	

Zinc sulfide

• Zinco solfuro • Zinc sulfure • Zinc sulfuro • Zinksulfid

Synonym:
Zinc sulphide

ZnS
Molecular Weight: 97,43
CAS: 1314-98-3
EEC-N: 215-251-3

HEU031

Zinc sulfide > RPE - For analysis

RPE

Description ..... White-green powder
Identification ..... Positive
Ammonium ..... ≤500 ppm
Chloride ..... ≤500 ppm
Heavy metals (Pb) ..... ≤20 ppm
Fe ..... ≤50 ppm
Assay (complexometric) ..... ≥98 %

Code	Size	Packaging	Notes
495105	250 g	Plastic bottle	
495107	1 kg	Plastic bottle	

Zirconium powder

• Zirconio polvere • Zirconium en poudre • Zirconio polvo • Zirkoniumpulver

Zr
Molecular Weight: 91,22
CAS: 7440-67-7
EEC-N: 231-176-9

Classification transport
ONU: 2008
Transport Hazard class: 4.2
Packing group I

Danger
H250-H260
P210-P222-P223-P231a+P232-P280-P402+P404

Zirconium powder > RPE - For analysis

RPE

Description ..... Dark grey powder Identification ..... Positive Assay (gravimetric) ..... ≥97 %

Code	Size	Packaging	Notes
495202	25 g	Glass bottle	

Zirconium standard solution

• Zirconio standard soluzione • Zirconium solution standard • Zirconio, solución patrón • Zirkonium-Standardlösung

Zirconium standard solution > RS - For analysis according to Ph. Eur. Chap. 4.1.2

RS

Code	Size	Packaging	Notes
615003500	100 ml	Plastic bottle	A 1 g/l solution Ref Ph.Eur 5003500

Zirconium standard solution > RS - Standard solution for ICP-MS

RS

Code	Size	Packaging	Notes
505957	100 ml	Plastic bottle	conc. 10 ppm Matrix: Nitric acid and hydrofluoric acid
505958	100 ml	Plastic bottle	conc. 100 ppm Matrix: Nitric acid and hydrofluoric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Zirconium standard solution > RS - Standard solution for ICP

RS

Code	Size	Packaging	Notes
504091	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid and hydrofluoric acid
504095	100 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid and hydrofluoric acid
504093	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Hydrochloric acid and hydrofluoric acid
504097	500 ml	Plastic bottle	conc. 10.000 ppm Matrix: Hydrochloric acid and hydrofluoric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Zirconium standard solution > RS - Standard solution for AAS

RS

Code	Size	Packaging	Notes
507770	100 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid
507524	500 ml	Plastic bottle	conc. 1.000 ppm Matrix: Nitric acid and hydrofluoric acid

Provided with a Certificate of analysis with references on the analytical method, the N.I.S.T. Standard Reference Materials and the confidence interval

Zirconium (IV) oxide

• Zirconio ossido • Zirconium (IV) oxyde • Zirconio (IV) óxido • Zirkonium (IV) oxid

Synonym:
Zirconia

ZrO_2
Molecular Weight: 123,22
CAS: 1314-23-4
EEC-N: 215-227-2

Zirconium (IV) oxide > RPE - For analysis

RPE

Description ..... Polvere bianca
Identification ..... Positive
Chloride ..... ≤100 ppm
Fe ..... ≤200 ppm
Assay (complexometric) ..... ≥99 %

Code	Size	Packaging	Notes
495305	250 g	Glass bottle	

Product specifications are subject to changes. Please visit our website for updates.

Zirconyl nitrate

• Zirconile nitrato • Zirconyle nitrate • Zirconio oxinitrato • Zirconyl nitrat

Synonym:
Zirconium(IV) oxynitrate hydrate

$ZrO(NO_3)_2.nH_2O$
Molecular Weight: 231,23 (an,)
CAS: 14985-18-3
EEC-N: 237-529-3

Classification transport
ONU: 1477
Transport Hazard class: 5.1
Packing group II

Danger
H272-H315-H319-H335
P210-P261-P280-P304+P340-P305+P351+P338-P403+P233

Zirconyl nitrate > RE - Pure

RE

Description ..... White powder Identification ..... Positive Hf..... ≤ 4 % Assay (gravimetric) ..... ≥ 99.4 %

Code	Size	Packaging	Notes
396105	50 g	Glass bottle	

Zirconyl nitrate

• Zirconile nitrato • Zirconyle nitrate • Zirconio oxinitrato • Zirconyl nitrat

Synonym:
Zirconium(IV) oxynitrate hydrate

$ZrO(NO_3)_2.nH_2O$
Molecular Weight: 231,23 (an,)
CAS: 14985-18-3
EEC-N: 237-529-3

Classification transport
ONU: 1477
Transport Hazard class: 5.1
Packing group II

Danger
H272-H315-H319-H335
P210-P261-P280-P304+P340-P305+P351+P338-P403+P233

Zirconyl nitrate > RE - Pure

RE

Description ..... White powder Identification ..... Positive Hf..... ≤ 4 % Assay (gravimetric) ..... ≥ 99.4 %

Code	Size	Packaging	Notes
396105	50 g	Glass bottle	

BEST CHOICE

Xcipharm™

Your requirements are our mission

The new product range for excipient use:



QUALITY ASSURANCE

- Flow chart
- Change control

PROCESS

- Validated cleaning procedures and dedicated equipment
- Labels consolidation
- Certificate of analysis with the name of the producer and the date of production of the raw material
- Stability studies

DOCUMENTATION

- BSE/TSE statement
- OGM statement
- Residual solvents statement
- ICH Q3D
- Risk assessment (2015/C95/02)

www.carloerbareagents.com

DASIT GROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

CHEMICAL TABLES




BEST CHOICE


Shuttle Service

Stainless steel returnable containers to optimize solvent quality and the management of packaging waste

- Improved safety
- Environmental friendly
- Quality guaranteed
- Efficient Logistic



www.carloerbareagents.com


PERIODIC TABLE OF THE ELEMENTS
Legend
Categories
Common Constants
Phase at STP
Gas Liquid Solid Synthetic
Period
Lanthanides Series
Actinides Series
Lanthanides
Actinides
Chemicals, Labware and Life Sciences
for Laboratories and Industry
DASITGROUP
CARLO ERBA
REAGENTS
www.carloerbareagents.com

MEASUREMENTS AND SYMBOLS

Unit Name	Symbols	Dimension
Ångström	Å	Length
Becquerel	Bq	Activity
Calorie	cal	Heat quantity
Candela	cd	Luminous intensity
Kilogram	kg	Mass
Coulomb	C	Eletric carge
Second	s	Time
Dyne	dyn	Force
Dyne per centimetre	dyn/cm	Surface tension
Erg	erg	Work, energy
Farad	F	Eletric capacitance
Ampere	A	Electric current
Herz	Hz	Frequency
Joule	J	Energy, work
Metre	m	Length
Micron	µ o µm	Length
Mole	mol	Amount of substance
Newton	N	Force
Pascal	Pa	Pressure
Poise	P o Po	Dynamic viscosity
Liter	L	Volume
Stokes	St	Kinematic viscosity
Kelvin	K	Temperature
Volt	V	Eletric potential
Watt	W	Power

DECIMAL UNIT MULTIPLES AND SUBMULTIPLES

Multiples

Factor	Name	Decimal number	Prefix	Symbol
10^{18}	trillion	1 000 000 000 000 000 000	exa	E
10^{15}	–	1 000 000 000 000 000	penta	P
10^{12}	one thousand bilion	1 000 000 000 000	tera	T
10^{9}	bilion	1 000 000 000	giga	G
10^{6}	milion	1 000 000	mega	M
10^{3}	one thousand	1 000	kilo	k
10^{2}	one hundred	100	hecto	h
10^{1}	ten	10	deca	da
10^{0}	one	1	–	–

Submultiples

Factor	Name	Decimal number	Prefix	Symbol
10^{0}	one	1	–	–
10^{-1}	one tenth	0,1	deci	d
10^{-2}	one hundredth	0,01	centi	c
10^{-3}	one thousandth	0,001	milli	m
10^{-6}	one milionth	0,000 001	micro	µ
10^{-9}	one bilionth	0,000 000 001	nano	n
10^{-12}	one thousand bilionth	0,000 000 000 001	piko	p
10^{-15}	–	0,000 000 000 000 001	femto	f
10^{-18}	one trilionth	0,000 000 000 000 000 001	atto	a

APHA COLOUR

APHA	Pt-Co ml	Water ml
10	1.0	49.0
20	2.0	48.0
30	3.0	47.0
40	4.0	46.0
60	6.0	44.0
80	8.0	42.0
100	10.0	40.0
150	15.0	35.0
200	20.0	30.0
300	30.0	20.0
400	40.0	10.0
500	50.0	–

The colour of the liquid (pure substance or solutions) might be conventionally expressed in APHA units (American Public Health Association).

The determination is carried out comparing the colour of the liquid with that of reference solution prepared under specific conditions. For the comparison, two identical 50 ml Nessler cylinders of transparent glass, containing equal volumes of the liquid and reference solution were used. This was prepared by diluting a certain amount of Platinum-Cobalt, so as to obtain the given APHA value according to the ratios given in the table below.

Platinum-Cobalt (500 APHA)

Dissolve 1,246 g of Potassium chloroplatinate RPE and 1,000 g of Cobaltous chloride hexahydrate RPE in 200 ml of distilled water. Add 100 ml of hydrochloric acid 37% and dilute to 1000 ml with distilled water. This solution has a conventional colorimetric value of 500 APHA units.

VISCOSITY - UNIT OF MEASURE

International System (SI)

Shear stress τ ..... Pascal (Pa)=Newton/m² (N/m²)
Velocity gradient g..... m/s
Dynamic viscosity h ..... Pascal x second (Pa·s)
..... Millipascal x second (mPa·s)
Kinematic viscosity n..... m²/s =104 Stokes
..... mm²/s

CGS System

Shear stress τ ..... dine/cm²
Velocity gradient g..... cm/s
Dynamic viscosity h ..... Poise (P) = dine·s/cm² = 1 Pa
..... Centipoise (cP)
Kinematic viscosity n..... Stokes (St) = 0,1 Pa·s
..... Centistokes (cSt) = 1 mPa·s

Tab. 1 - Distilled water – specific viscosity at different temperatures (1)

Temperature	cP (centipoise)	Viscosity
0°	0,0179	1,000
5°	0,0151	0,843
10°	0,0130	0,730
15°	0,0114	0,637
17,5°	0,0107	0,599
20°	0,0100	0,561
30°	0,0080	0,446
50°	0,0054	0,307
70°	0,0040	0,226
100°	0,0028	0,158

(1) From: Küster F.W. - Thiel A., Tabelle Logaritmiche, ed. Hoepli, 1965

Tab. 2 - Table of viscosity in increasing order (cP at 20° C)

Solvent	Viscosity (cP)
Pentane	0.23
Diethyl ether	0.23
Methyl-tert-butyl ether	0.27
Petroleum ether	0.30
Hexane	0.31
Acetone	0.32
Acetonitrile	0.36
Heptane	0.41
Dichlorometane	0.43
1-Chlorobutane	0.45
Ethyl acetate	0.45
2,2,4 Trimethylpentane	0.51

Solvent	Viscosity (cP)
Methanol	0.55
Tetrahydrofuran	0.55
Chloroform	0.57
Toluene	0.59
Benzene	0.65
1,1,2-Tricloro 1,2,2-trifluoroethane	0.71
1,2 Dichloroethane	0.79
Dimethylformamide	0.85
Tetrachloroethylene	0.93
Pyridine	0.95
Carbon tetrachloride	0.97
Cyclohexane	1.00

Solvent	Viscosity (cP)
Water	1.00
Hethanol absolute	1.20
Acetic acid glacial	1.29
1,4 Dioxane	1.54
2-Methoxyethanol	1.72
Dimethylsulfoxide	2.24
Propan-1-ol	2.26
Propan-2-ol	2.30
Octan-1-ol	approx. 10.64

DENSITY

Ammonium hydroxide

d 15°C 4°C	°Bé	% NH3
1.000	10	–
0.992	11	1.61
0.986	12	3.30
0.979	13	4.80
0.972	14	6.55
0.966	15	8.33
0.959	16	9.91
0.953	17	11.60
0.947	18	13.31
0.941	19	15.04
0.935	20	17.12
0.929	21	18.64
0.923	22	20.08
0.917	23	22.39
0.912	24	24.34
0.906	25	26.31
0.900	26	27.99
0.895	27	29.69
0.889	28	31.75

Hydrochloric acid

15°C d 4°C	°Bé	% m HCl
1.0069	1	1.56
1.014	2	2.99
1.021	3	4.55
1.028	4	5.99
1.036	5	7.56
1.043	6	9.14
1.050	7	10.59
1.058	8	12.17
1.066	9	13.61
1.074	10	15.16
1.082	11	16.70
1.090	12	18.30
1.098	13	20.00
1.106	14	21.60
1.115	15	23.05
1.124	16	24.79
1.133	17	26.55
1.142	18	28.15
1.151	19	29.95
1.160	20	32.10
1.169	21	33.65
1.179	22	35.40
1.189	23	37.25
1.199	24	39.10

Nitric acid

d 15°C 4°C	°Bé	% m HNO_3
1.0069	1	1,39
1.014	2	2,69
1.021	3	4,08
1.028	4	5,37
1.036	5	6.76
1.043	6	8,13
1.050	7	9,35
1.058	8	10,68
1.066	9	11,88
1.074	10	13,15
1.082	11	14,47
1.090	12	15,70
1.098	13	17,11
1.106	14	18,46
1.115	15	19,61
1.124	16	21,00
1.133	17	22,40
1.142	18	23,70
1.151	19	25,15
1.160	20	26,65
1.169	21	28,03
1.179	22	29,38
1.189	23	30,88
1.199	24	32,36
1.209	25	33,80
1.219	26	35,28
1.229	27	36,96
1.240	28	38,44
1.250	29	40,12
1.261	30	41,81
1.273	31	43,49
1.284	32	45,18
1.295	33	46,98
1.307	34	48,72
1.319	35	50,71
1.331	36	52,80
1.344	37	54,93
1.356	38	57,13
1.369	39	59,39
1.382	40	61,92
1.396	41	64,71
1.409	42	67,50
1.423	43	70,80
1.437	44	74,32
1.452	45	78,18
1.467	46	82,48
1.482	47	87,23
1.498	48	93,45
1.513	49	99,07

Phosphoric acid

d 15°C 4°C	°Bé	% m H_3PO_4
1.0069	1	1.38
1.014	2	2.76
1.021	3	4.13
1.028	4	5.51
1.036	5	6.90
1.043	6	8.26
1.050	7	9.64
1.058	8	11.02
1.066	9	12.40
1.074	10	13.77
1.082	11	15.15
1.090	12	16.53
1.098	13	17.91
1.106	14	19.28
1.115	15	20.66
1.124	16	22.04
1.133	17	23.42
1.142	18	24.80
1.151	19	26.17
1.160	20	27.55
1.169	21	28.93
1.179	22	30.31
1.189	23	31.68
1.199	24	33.06
1.209	25	34.44
1.219	26	35.82
1.229	27	37.19
1.240	28	38.57
1.250	29	39.95
1.261	30	41.33
1.273	31	42.70
1.284	32	44.08
1.295	33	45.46
1.307	34	46.84
1.319	35	48.21
1.331	36	49.59
1.344	37	50.97
1.356	38	52.04
1.369	39	53.72
1.382	40	55.10
1.396	41	56.48
1.409	42	57.86
1.423	43	59.23
1.437	44	60.61
1.452	45	61.99
1.467	46	63.37
1.482	47	64.75
1.498	48	66.12
1.513	49	67.50
1.529	50	68.88
1.545	51	70.26
1.562	52	71.63
1.579	53	73.01
1.597	54	74.39
1.615	55	75.77
1.633	56	77.14
1.652	57	78.52
1.671	58	79.90
1.690	59	81.28
1.710	60	82.65
1.731	61	83.03
1.752	62	85.41
1.773	63	86.80
1.795	64	88.16
1.818	65	89.55
1.841	66	90.92

Potassium hydroxide

d 15°C 4°C	°Bé	% m KOH
1.0069	1	0.9
1.014	2	1.7
1.021	3	2.6
1.028	4	3.5
1.036	5	4.5
1.043	6	5.5
1.050	7	6.4
1.058	8	7.4
1.066	9	8.3
1.074	10	9.2
1.082	11	10.1
1.090	12	11.0
1.098	13	12.0
1.106	14	12.9
1.115	15	13.8
1.124	16	14.8
1.133	17	15.7
1.142	18	16.6
1.151	19	17.6
1.160	20	18.6
1.169	21	19.5
1.179	22	20.5
1.189	23	21.4
1.199	24	22.4
1.209	25	23.3
1.219	26	24.2
1.229	27	25.1
1.240	28	26.1
1.250	29	27.0
1.261	30	28.0
1.273	31	28.9
1.284	32	29.8
1.295	33	30.7
1.307	34	31.7
1.319	35	32.7
1.331	36	33.7
1.344	37	34.9
1.356	38	35.9

1.369	39	36.9
1.382	40	37.9
1.396	41	38.9
1.409	42	39.9
1.423	43	40.9
1.437	44	42.1
1.452	45	43.4
1.467	46	44.6
1.482	47	45.8
1.498	48	47.1
1.513	49	48.3
1.529	50	49.4
1.545	51	50.6

Sodium hydroxide

d 15°C 4°C	°Bé	% m NaOH
1.0069	1	0.59
1.014	2	1.20
1.021	3	1.85
1.028	4	2.50
1.036	5	3.15
1.043	6	3.79
1.050	7	4.50
1.058	8	5.20
1.066	9	5.86
1.074	10	6.58
1.082	11	7.30
1.090	12	8.07
1.098	13	8.78
1.106	14	9.50
1.115	15	10.30
1.124	16	11.06
1.133	17	11.90
1.142	18	12.69
1.151	19	13.50
1.160	20	14.35
1.169	21	15.15
1.179	22	16.00
1.189	23	16.90
1.199	24	17.81
1.209	25	18.71
1.219	26	19.65
1.229	27	20.60
1.240	28	21.55
1.250	29	22.50
1.261	30	23.50
1.273	31	24.48
1.284	32	25.50
1.295	33	26.58
1.307	34	27.65
1.319	35	28.83
1.331	36	30.00
1.344	37	31.20
1.356	38	32.50
1.369	39	33.73
1.382	40	35.00
1.396	41	36.36
1.409	42	37.65
1.423	43	39.06
1.437	44	40.47
1.452	45	42.02
1.467	46	43.58
1.482	47	45.16
1.498	48	46.73
1.513	49	48.41
1.529	50	50.10

Sulphuric acid

d 15°C 4°C	°Bé	% m H_2SO_4
1.0069	1	1.20
1.014	2	2.20
1.021	3	3.35
1.028	4	4.40
1.036	5	5.54
1.043	6	6.67
1.050	7	7.67
1.058	8	8.77
1.066	9	9.78
1.074	10	10.90
1.082	11	12.06
1.090	12	13.13
1.098	13	14.35
1.106	14	15.48
1.115	15	16.49
1.124	16	17.66
1.133	17	18.85
1.142	18	19.93
1.151	19	21.17
1.160	20	22.45
1.169	21	23.60
1.179	22	24.76
1.189	23	26.04
1.199	24	27.32
1.209	25	28.58
1.219	26	29.84
1.229	27	31.23
1.240	28	32.40
1.250	29	33.66
1.261	30	34.90
1.273	31	36.17
1.284	32	37.45
1.295	33	38.84
1.307	34	40.12
1.319	35	41.50
1.331	36	42.98
1.344	37	44.28
1.356	38	45.62
1.369	39	46.94
1.382	40	48.35
1.396	41	49.85
1.409	42	51.15
1.423	43	52.51
1.437	44	53.91
1.452	45	55.34
1.467	46	56.74
1.482	47	58.13
1.498	48	59.54
1.513	49	61.12
1.529	50	62.53
1.545	51	64.05
1.562	52	65.50
1.579	53	66.95
1.597	54	68.41
1.615	55	70.00
1.633	56	71.70
1.652	57	73.18
1.671	58	74.80
1.690	59	76.50
1.710	60	78.04
1.731	61	80.02
1.752	62	81.83
1.773	63	84.00
1.795	64	86.30
1.818	65	90.05
1.841	66	95.69

POLARITY

Table of polarity in increasing order (p)

Solvent	Polarity (p)
Heptane	0.1
Hexane	0.1
Petroleum ether	0.1
2,2,4 Trimethylpentane	0.1
Cyclohexane	0.2
1-Chlorobutane	1.0
Carbon tetrachloride	1.6
Toluene	2.4
Metyl-tert butyle ether	2.5
Benzene	2.7
Diethyl ether	2.8
Dichloromethane	3.1
Octan-1-ol	3.4
1,2 Dichloroethane	3.5
Propan-1-ol	3.9
Propan-2-ol	4.0
Tetrahydrofuran	4.0
Chloroform	4.1
Ethanol absolute	4.3
Ethyl acetate	4.4
1,4 Dioxane	4.8
Acetone	5.1
Methanol	5.1
Pyridine	5.3
2-Methoxyethanol	5.5
Acetonitrile	5.8
Acetic acid glacial	6.0
Dimethylformamide	6.4
Dimethysulfoxide	7.2
Water	10.2

INDICATORS

Table of pH range and colour shades

Indicator	pH range	Acid	Basic
Malachite green	0.0-2.0	yellow	green-blue
Brilliant green	0.0-2.6	yellow	green
Eosin Y	0.0-3.0	yellow	green
Erythrosin B	0.0-3.6	orange	red
Methyl green	0.1-2.3	yellow	blue
Methyl violet	0.1-2.7	yellow	violet
Picric acid	0.2-1.0	colourless	yellow
Cresol red	0.2-1.8	red	yellow
Crystal violet	0.8-2.6	yellow	blue-violet
Thymol blue	1.2-2.8	red	yellow
Tropaeolin OO	1.3-3.2	red	yellow
Eosin B	1.4-2.4	colourless	rose
Quinaldine red	1.4-3.2	colourless	rose
2,4-Dinitrophenol	2.4-4.0	colourless	yellow
Methyl yellow	2.9-4.0	red	yellow
Bromophenol blue	3.0-4.6	yellow	blue-violet
Congo red	3.0-5.2	blue	yellow-orange
Methyl orange	3.1-4.4	red	orange
Alizarine sodium sulphonate	3.7-5.2	yellow	violet
a-Naphtil red	3.7-5.0	red	yellow
Bromocresol green	4.0-5.6	yellow	blue
2,5-Dinitrophenol	4.0-5.8	colourless	yellow
Alizarine red	4.3-6.3	yellow	violet
Methyl red	4.4-6.2	red	yellow
Chlorophenol red	4.8-6.4	yellow	red
Bromocresol purple	5.2-6.8	yellow	purple
p-Nitrophenol	5.4-7.5	colourless	yellow
Bromoxylenol blue	5.7-7.5	yellow	blue
Alizarine	5.8-7.2	yellow	red
Bromothymol blue	6.0-7.6	yellow	blue
Bromophenol blue	6.2-7.6	yellow	blue
Phenol red	6.4-8.2	yellow	red
3-Nitrophenol	6.6-8.6	colourless	yellow-orange
Neutral red	6.8-8.0	red	yellow
Rosolic acid	6.8-8.0	yellow	red
Cresol red	7.2-8.8	yellow	red
a-Naphtolphtalein	7.3-8.7	rose	green
Cresol purple	7.4-9.0	yellow	purple
Tropaeolin OOO	7.6-8.9	yellow	rose-red
Thymol blue	8.0-9.6	yellow	blue
Phenolphtalein	8.0-10.0	colourless	red
a-Naphtolbenzein	9.0-11.0	yellow	blue
Thymolphtalein	9.4-10.6	colourless	blue
Alkali blue 6B	9.4-14.0	violet	rose
Alizarin	10.0-12.0	yellow	purple
Nilo blue	10.1-11.1	blue	red
Diazoviolet	10.1-12.0	yellow	violet
Tropaeolin O	11.0-13.0	yellow	orange-brown
Nitramine	11.0-13.0	colourless	orange-brown
Poirrier blue	11.0-13.0	blue	violet-rose
Clayton's yellow O	12.0-13.0	yellow	red
Trinitrobenzoic acid	12.0-13.4	colourless	orange-red
Indigo carmine dried	11.5-13.0	blue	yellow
Epsilon blue	11.6-13.0	orange	violet

COLOUR INDEX

Colour Index	Colour Index Name	Commercial Name	Synonyms Index
10316	Acid Yellow 1	Naphthol yellow S	Acid yelow S
11020	Solvent Yellow 1,2	Methyl yellow	
11270	Basic Orange 2	Chrysoidin Y	Brown salt R
12055	Solvent Yellow 14	Sudan yellow	Sudan I Sudan yellow R
12140	Solvent Orange 7	Sudan II	Sudan red Sudan Orange RR
13020	Acid Red 2	Methyl red	
13025	Acid Orange 52	Methyil orange	Orange III Helianthin
13065	Acid Yellow 36	Methanyl yellow	Tropaeolin G
13080	Acid Orange 5	Tropaeolin 00	Orange IV
14030	Mordant orange 1	Alizarin yellow R	Alizarin yellow G Orange R
14270	Acid Orange 6	Tropaeolin O	Tropeolina Y
14645	Mordant black 11	Heriochrome black T	Superchrome black T
15510	Acid Orange 7	Orange II	Tropaeolin 000
15705	Mordant black 17	Calcon	Palatine chrome balck Eriochrome blue black B
16150	Acid red 26	Ponceau de Xilidine	Ponceau 2 R Brilliant Ponceau
16185	Acid Red 27	Amaranth	Naphthol red S, C o O Solid red O
16230	Acid Orange 10	Orange G	Orange GG
16570	Acid red 29	Chromotrope 2R	Acid phloxin GR
19140	Acid yellow 23	Tartrazine	Acid yellow T
19540	Direct yellow 9	Titan yellow	Thiazole yellow G Clayton yellow
20470	Acid Black 1	Naphthalene black 12 B	Naftol blue black Amido black 10B Pontacyl black blue SX
21010	Basic brown 4	Bismark Brown R	Vesuvine BL
22120	Direct red 28	Congo red	Cotton redB
23850	Direct blue 14	Trypan blue	Congo blue 3B
23860	Direct blue 53	Evans blue	Diazol pure blue Geigy blue 536 med
24890	Direct yellow 4	Brilliant yellow	Yellow paper
26050	Solvent red 19	Sudan red 7B	Fast red 7B
26100	Solvent red 23	Sudan III	
26105	Solvent red 24	Sudan IV	Scarlett R (Michaelis) Fat ponceau
26125	Solvent red 27	Oil red 0	Sudan red 5B
26150	Solvent black 3	Sudan black B	Ceres black BN
26905	Acid red 66	Scarlett Biebrich	Imperial scarlett Brilliant ponceau S Ponceau red BS
27195	Acid red 112	Ponceau red S	Java scarlet
37025	Azoic Diazo No. 6	o-Nitroaniline	Orange GRS
37030		m-Nitroaniline	
37035	Azoic Diazo No. 37	p-Nitroaniline	Nitrazol CF Nitrosamine red
37235	Azoic diazo N. 48	Fast blue B salt	Dianisidine blue Diazo blue B salt Blue salt BNS
41000	Basic yellow 2	Auramine 0	Pyoctanine yellow
42000	Basic green 4	Malachite green	Vittoria green B China green
42040	Basic green 1	Brilliant green	Aniline green Diamone green Emerald green
42045	Acid blu 1	Eriogalucine	Disulphine blue V Sulphon blue
42053	Food green 3	Fast green FCF	
42090	Acid Blue 9	Erioglaucine	Alphazurine FG
42095	Acid green 5	Light green SF	Acid green F Acid green G Lissamine green SF
42135	Acid blue 147	Xilencyanol FF	Cyanol FF
42510	Basic violet 14	Rosanilin	Fuchsin brilliant
		Basic fuchsin	Rosaniline hydrocholride Magenta I
42535	Basic violet 1	Gentian violet	Methyl violet 2R

Colour Index	Colour Index Name	Commercial Name	Synonyms Index
42535	Basic violet 1	Violetto Metile 2 B	
42555	Basic violet 3	Crystal violet	Methyl violet 6B
42556	Basic green	Iodine green	
42563	Basic blue 8	Vittoria blue 4R	Fast Blue 4R
42585	Basic blue 20	Methyl green	
42600	Basic violet 4	Ethyl violet	Ethyl purple 6B
42655	Acid blue 90	Brilliant Indocyanin G	Coomassie brilliant blue G250 Eriosin Cyanin brilliant G
42660	Acid blue 83	Brilliant Indocyanin 6 B	Coomassie brilliant blue R Brilliant acid cyanine 6B
42685	Acid violet 19	Acid fuchsin	Fuchsin S Rubin S Acid Magenta
42755	Acid blue 22	Aniline blue (water soluble)	China blue Cotton blue Blu di Hofman Opal blue Water blue I
42765	Acid blue 119	Alkali blue 6B	Reflex blue AG
42775	Solvent blue 3	Aniline blue (alcohol soluble)	Light blue Lyon's blue Paris blue Gentian blue
42780	Acid blue 93	Methyl blue	Helvetia blue Soluble blue 8B Poirier's blue C4B
43800		Rosolic acid sodiu salt	Aurine (water soluble) Corollin (water soluble)
43820	Mordant blue 3	Chromoxane canine R	Cyanin R Solochrome Eriochrome canine R
43825	Mordant blue 29	Cromoxane pure blue BLD	Cromeazurol S
45005	Basic dye	Pyronine G	Pyronine Y
45170	Basic violet 10	Rhodamine B	Rhodamine O Brilliant rhodamine B
45350	Acid yellow 73	Fluorescein sodium salt	Uranin
45380	Acid red 87	Eosin Y (yellowish)	Tetrabromofluoresceina sodica
45386	Solvent red 45	Ethyl Eosin (alcohol soluble)	Eosin S
45400	Acid red 91	Eosina B (blue shade)	Eosin scarlet
45410	Acid red 92 (soluble in acqua)	Phloxin B	Cyanosin Magdala red Tetrabromotetrachlorofluoroscein Sodium salt
45430	Acid red 51	Erythrosin B	Erythrosin J
45440	Acid red 94	Rose Bengal	
46005	Basic orange 14	Acridine orange	Euchrysine
49700		Indophenol	Indophenol blue
50040	Basic red 5	Neutral red	Toluylene red Neophospine
50240	Basic red 2	Safranine O	Cotton red
50420	Acid black 2	Nigrosine (water soluble)	Aniline blue black
51010	Basic dye	Brilliant cresyl blue	Cresyl blue BBS
51050	Mordant Blue 14	Celestine blue B	Coerin 2R
51180	Basic blue 12	Nilo blue A	Nilo blue BX
52000	Basic violet	Thionine acetate	Lauth's violet
52015	Basic blue 9	Methylene blue	
52040	Basic blue 17	Toluidine blue	
56085	Mordant dye	Murexide	
58000	Mordant red 11	Alizarin	
58005	Mordant red 3	Alizarin red S	
58500	Mordant violet-26	Quinizarin	Alizarin orange A Alizarina cianina 3R
60760	Pigment dyes	Nuclear fast red	Calcium red Kerneckrot Helio fast rubin BBL
61515	Solvent blue 19	Blu Oracet B	
73000	Vat blue 1	Indigo	Indigo blue
73015	Acid blue 74	Indigo carmine	Sodium indigo disulphonate
74240	Ingrain blue 1	Alcian blue 8GX	Alcian blue
75290	Natural black 1	Hematoxylin	Hematein
75300	Natural yellow 3	Curcumin	Curcuma
75470	Natural red 4	Acido carminico	Carminio Cocciniglia
75660	Natural Yellow 11	Morin	Fustic

SOLUTIONS CHEMISTRY


Freezing mixtures

Mixture	Solution concentration	Temperature °C
Ammonium chloride	solution 23 %	- 3° C
Potassium chloride	solution 20 %	- 12° C
Ammonium nitrate	solution 50 %	- 15° C
Sodium chloride	solution 25 %	- 21° C
Sodium nitrate	solution 33 %	- 24° C
Calcium chloride 6 H2O	solution 62 % - with ice	- 39° C
Calcium chloride 6 H2O	solution 59 % - with ice	- 55° C
Methanol or Acetone with dry ice		- 77° C

Molarity and normality chart for common acid and base solutions

			Volume required for a liter	
Acid	**Molarity**	**Normality**	*1 M solution*	*1 N solution*
Acetic acid 99.5%	17,4 M	17,4 N	57,5 ml	57,5 ml
Ammonia sol. 25%	13,2 M	13,2 N	75,6 ml	75,6 ml
Ammonia sol. 35%	18,5 M	18,5 N	54,0 ml	54,0 ml
Hydrochloric acid 37%	11,6 M	11,6 N	85,8 ml	85,8 ml
Hydrochloric acid 32%	10,2 M	10,2 N	98,2 ml	98,2 ml
Hydrofluoric acid 40%	22,6 M	22,6 N	44,2 ml	44,2 ml
Nitric acid 70%	15,7 M	15,7 N	63,7 ml	63,7 ml
Phosphoric acid 85%	14,7 M	44,2 N	67,8 ml	22,6 ml
Sulphuric acid 96%	18,0 M	36,0 N	55,5 ml	27,8 ml

Miscibility table



Solvent	b.p. (°C)	UV (nm) cut off 1AU	d (g/mL) at 20°C	Refractive index at 20°C	Viscosity (cP) at 20°C
Acetone	56	330	0,786	1,359	0,32
Acetonitrile	82	190	0,786	1,344	0,37
Water	100	190	0,998	1,333	1,00
Ethanol	78	210	0,789	1,360	1,20
Isopropanol	82	205	0,785	1,377	2,30
Methanol	65	205	0,791	1,329	0,60
n-butanol	125	215	0,81	1,394	0,73
n-propanol	97	210	0,803	1,384	2,27
Benzene	80	280	0,879	1,501	0,65
Butyl Acetate	125	254	0,882	1,399	2,98
Cyclohexane	81	200	0,779	1,426	1,00
Chloroform	61	245	1,498	1,946	0,57
1,2-Dichloroethane	84	225	1,257	1,444	0,79
Dichloromethane	41	233	1,326	1,424	0,44
Diethylether	35	215	0,713	1,353	0,23
Dimethylformamide	155	268	0,944	1,431	0,85
Dimethysulfoxide	189	268	1,092	1,478	2,24
Dioxane	101	215	1,033	1,422	1,54
Heptane	98	200	0,684	1,387	0,41
Hexane	69	195	0,655	1,375	0,31
Ethyl Acetate	77	256	0,894	1,372	0,45
Isoctane	99	215	0,692	1,392	0,51
Methyl-ethyl ketone	80	329	0,806	1,379	0,45
Methyl-tert-butyl ether	55	210	0,741	1,369	0,27
Pentane	36	190	0,626	1,358	0,23
Tetrahydrofurane	65	212	0,886	1,407	0,55
Toluene	111	284	0,867	1,496	0,59
Trichloroethylene	87	273	1,462	1,477	0,57
Xylene	139	288	0,861	1,500	0,61

full square means immiscible

CONVERSION TABLES

US and British measuring units and conversion factors

Lenght		
1 mil	=	25,4 µm
1 inch	=	2,54000 centimetres
1 foot	=	30,48006 centimetres
1 yard	=	0,91440 metres
1 mile	=	1609,34 metres
1 mile (nautical)	=	1853,25 metres
1 millimetre	=	0,03937 inches
1 centimetre	=	0,3937 inches
1 metre	=	39,37 inches
1 metre	=	3,2028 foots
1 metre	=	1,09361 yards
1 kilometer	=	0,62136 miles
1 kilometer	=	0,53959 miles (nautical)

Volume		
1 cubic inch	=	16,38716 cubic centimetres
1 cubic foot	=	28,31625 cubic decimetres
1 cubic yard	=	0,76456 cubic metres
1 cubic centimetre	=	0,06102 cubic inches
1 cubic decimetre	=	0,03531 cubic foots
1 cubic metre	=	1,30794 cubic yards

Capacity		
1 quart (USA liquid)	=	0,94633 liters
1 gallon (USA)	=	3,78533 liters
1 barrel (USA)	=	0,11562 cubic metres
1 quart (UK)	=	1,13650 litres
1 gallon (UK)	=	4,5596 litres
1 barrel (UK)	=	0,16366 cubic metres
1 litre	=	1,056681 quarts (USA)
1 litre	=	0,264177 gallons (USA)
1 litre	=	0,87990 quarts (UK)
1 litre	=	0,219976 gallons (UK)

Weight (Mass)		
1 grain	=	64,80 milligrams
1 dramma	=	1,772 grams
1 ounce (US)	=	28,3495 grams
1 pound	=	453,5924 grams
1 ton (short) (US)	=	907,18486 kilograms
1 ton (long) (UK)	=	1016,0470 kilograms
1 grams	=	15,4324 grains
1 gram	=	0,03527 ounces
1 kilogram	=	2,20462 pounds
1 metric ton.	=	1,10231 short tons (US)
1 metric ton	=	0,98420 long tons (UK)

Concentrations

		ppm	**ppb**	**ppt**	**Potency**	**Proportion**
%	**g/kg** **mg/g** **µg/mg**	**mg/kg** **µg/g** **ng/mg**	**µg/kg** **ng/g** **pg/mg**	**ng/kg** **pg/g** **fg/mg**		
1	10	10.000			1×10^{-2}	1: 100
0,5	5	5.000			5×10^{-3}	1: 500
0,1	1	1.000			1×10^{-3}	1: 1.000
0,05	0,5	500			5×10^{-4}	1: 5.000
0,01	0,1	100			1×10^{-4}	1: 10.000
0,005	0,05	50			5×10^{-5}	1: 50.000
0,001	0,01	10	10.000		1×10^{-5}	1: 100.000
0,000.5	0,005	5	5.000		5×10^{-6}	1: 500.000
0,000.1	0,001	1	1.000		1×10^{-6}	1: 1.000.000
0,000.05	0,0005	0,5	500		5×10^{-7}	1: 5.000.000
0,000.01	0,000.1	0,1	100		1×10^{-7}	1: 10.000.000
0,000.001	0,000.01	0,01	10	10.000	1×10^{-8}	1: 100.000.000
0,000.0001	0,000.001	0,001	1	1.000	1×10^{-9}	1: 1.000.000.000
			0,1	100	1×10^{-10}	1: 10.000.000.000
			0,01	10	1×10^{-11}	1: 100.000.000.000
			0,001	1	1×10^{-12}	1: 1.000.000.000.000

Transmittance vs Absorbance unit

% T	A.U.	% T	A.U.	% T	A.U.	% T	A.U.	% T	A.U.	% T	A.U.
0.5	2.301	17.5	0.757	34.5	0.462	51.5	0.288	68.5	0.164	85.5	0.068
1.0	2.000	18.0	0.745	35.0	0.456	52.0	0.284	69.0	0.161	86.0	0.066
1.5	1.824	18.5	0.733	35.5	0.450	52.5	0.280	69.5	0.158	86.5	0.063
2.0	1.699	19.0	0.721	36.0	0.444	53.0	0.276	70.0	0.155	87.0	0.060
2.5	1.620	19.5	0.710	36.5	0.438	53.5	0.271	70.5	0.152	87.5	0.058
3.0	1.523	20.0	0.699	37.0	0.432	54.0	0.268	71.0	0.149	88.0	0.056
3.5	1.469	20.5	0.688	37.5	0.426	54.5	0.263	71.5	0.146	88.5	0.053
4.0	1.398	21.0	0.678	38.0	0.420	55.0	0.260	72.0	0.143	89.0	0.051
4.5	1.347	21.5	0.667	38.5	0.414	55.5	0.256	72.5	0.140	89.5	0.048
5.0	1.301	22.0	0.658	39.0	0.409	56.0	0.252	73.0	0.137	90.0	0.046
5.5	1.260	22.5	0.647	39.5	0.403	56.5	0.248	73.5	0.134	90.5	0.043
6.0	1.222	23.0	0.638	40.0	0.398	57.0	0.244	74.0	0.131	91.0	0.041
6.5	1.187	23.5	0.628	40.5	0.392	57.5	0.240	74.5	0.128	91.5	0.039
7.0	1.155	24.0	0.620	41.0	0.387	58.0	0.237	75.0	0.125	92.0	0.036
7.5	1.125	24.5	0.611	41.5	0.382	58.5	0.233	75.5	0.122	92.5	0.034
8.0	1.097	25.0	0.602	42.0	0.377	59.0	0.229	76.0	0.119	93.0	0.032
8.5	1.071	25.5	0.593	42.5	0.372	59.5	0.225	76.5	0.116	93.5	0.030
9.0	1.046	26.0	0.585	43.0	0.367	60.0	0.222	77.0	0.114	94.0	0.027
9.5	1.022	26.5	0.577	43.5	0.361	60.5	0.218	77.5	0.111	94.5	0.025
10.0	1.000	27.0	0.569	44.0	0.357	61.0	0.215	78.0	0.108	95.0	0.022
10.5	0.979	27.5	0.561	44.5	0.351	61.5	0.211	78.5	0.105	95.5	0.020
11.0	0.959	28.0	0.553	45.0	0.347	62.0	0.208	79.0	0.102	96.0	0.018
11.5	0.943	28.5	0.545	45.5	0.342	62.5	0.204	79.5	0.099	96.5	0.015
12.0	0.921	29.0	0.538	46.0	0.337	63.0	0.201	80.0	0.097	97.0	0.013
12.5	0.903	29.5	0.530	46.5	0.332	63.5	0.197	80.5	0.094	97.5	0.011
13.0	0.886	30.0	0.523	47.0	0.327	64.0	0.194	81.0	0.092	98.0	0.009
13.5	0.870	30.5	0.516	47.5	0.323	64.5	0.190	81.5	0.089	98.5	0.006
14.0	0.854	31.0	0.509	48.0	0.319	65.0	0.187	82.0	0.086	99.0	0.004
14.5	0.838	31.5	0.502	48.5	0.314	65.5	0.184	82.5	0.083	99.5	0.002
15.0	0.824	32.0	0.495	49.0	0.310	66.0	0.180	83.0	0.081	100.0	0.000
15.5	0.810	32.5	0.488	49.5	0.305	66.5	0.177	83.5	0.078		
16.0	0.796	33.0	0.482	50.0	0.301	67.0	0.174	84.0	0.076		
16.5	0.782	33.5	0.475	50.5	0.297	67.5	0.171	84.5	0.073		
17.0	0.770	34.0	0.469	51.0	0.292	68.0	0.168	85.0	0.071		

Baumé vs specific gravity

Conversion rules at a temperature of 60°F:

For liquids more dense than water:

$$s.g.= \frac{145}{145- \text{degrees Baumé}}$$

For liquids less dense than water:

$$s.g.= \frac{140}{130 + \text{degrees Baumé}}$$

°Be	Specific gravity	°Be	Specific gravity	°Be	Specific gravity	°Be	Specific gravity	°Be	Specific gravity	°Be	Specific gravity
103,33	0,60	36,67	0,84	10,74	1,08	34,73	1,32	51,45	1,55	63,77	1,79
101,40	0,61	35,68	0,85	11,36	1,09	35,15	1,32	51,75	1,56	63,99	1,79
99,51	0,61	34,71	0,85	11,97	1,09	35,57	1,33	52,05	1,56	64,22	1,80
97,64	0,62	33,74	0,86	12,58	1,10	35,98	1,33	52,35	1,57	64,44	1,80
95,81	0,62	32,79	0,86	13,18	1,10	36,39	1,34	52,64	1,57	64,67	1,81
94,00	0,63	31,85	0,87	13,78	1,11	36,79	1,34	52,94	1,58	64,89	1,81
92,22	0,63	30,92	0,87	14,37	1,11	37,19	1,35	53,23	1,58	65,11	1,82
90,47	0,64	30,00	0,88	14,96	1,12	37,59	1,35	53,52	1,59	65,33	1,82
88,75	0,64	29,09	0,88	15,54	1,12	37,99	1,36	53,81	1,59	65,55	1,83
87,05	0,65	28,19	0,89	16,11	1,13	38,38	1,36	54,09	1,60	65,77	1,83
85,38	0,65	27,30	0,89	16,68	1,13	38,77	1,37	54,38	1,60	65,98	1,84
83,74	0,66	26,42	0,90	17,25	1,14	39,16	1,37	54,66	1,61	66,20	1,84
82,12	0,66	25,56	0,90	17,81	1,14	39,55	1,38	54,94	1,61	66,41	1,85
80,53	0,67	24,70	0,91	18,36	1,15	39,93	1,38	55,22	1,62	66,62	1,85
78,96	0,67	23,85	0,91	18,91	1,15	40,31	1,39	55,49	1,62	66,83	1,86
77,41	0,68	23,01	0,92	19,46	1,16	40,68	1,39	55,77	1,63	67,04	1,86
75,88	0,68	22,17	0,92	20,00	1,16	41,06	1,40	56,04	1,63	67,25	1,87
74,38	0,69	21,35	0,93	20,54	1,17	41,43	1,40	56,31	1,64	67,46	1,87
72,90	0,69	20,54	0,93	21,07	1,17	41,80	1,41	56,59	1,64	67,67	1,88
71,44	0,70	19,73	0,94	21,60	1,18	42,16	1,41	56,85	1,65	67,87	1,88
70,00	0,70	18,94	0,94	22,12	1,18	42,53	1,42	57,12	1,65	68,08	1,89
68,58	0,71	18,15	0,95	22,64	1,19	42,89	1,42	57,39	1,66	68,28	1,89
67,18	0,71	17,37	0,95	23,15	1,19	43,25	1,43	57,65	1,66	68,48	1,90
65,80	0,72	16,60	0,96	23,66	1,20	43,60	1,43	57,91	1,67	68,68	1,90
64,44	0,72	15,83	0,96	24,17	1,20	43,95	1,44	58,17	1,67	68,88	1,91
63,10	0,73	15,08	0,97	24,67	1,21	44,31	1,44	58,43	1,68	69,08	1,91
61,78	0,73	14,33	0,97	25,17	1,21	44,65	1,45	58,69	1,68	69,28	1,92
60,48	0,74	13,59	0,98	25,66	1,22	45,00	1,45	58,95	1,69	69,48	1,92
59,19	0,74	12,86	0,98	26,15	1,22	45,34	1,46	59,20	1,69	69,68	1,93
57,92	0,75	12,13	0,99	26,63	1,23	45,68	1,46	59,45	1,70	69,87	1,93
56,67	0,75	11,41	0,99	27,11	1,23	46,02	1,47	59,71	1,70	70,06	1,94
55,43	0,76	10,70	1,00	27,59	1,24	46,36	1,47	59,96	1,71	70,26	1,94
54,21	0,76	0,72	1,01	28,06	1,24	46,69	1,48	60,20	1,71	70,45	1,95
53,01	0,77	1,44	1,01	28,53	1,25	47,03	1,48	60,45	1,72	70,64	1,95
51,82	0,77	2,14	1,02	29,00	1,25	47,36	1,49	60,70	1,72	70,83	1,96
50,65	0,78	2,84	1,02	29,46	1,26	47,68	1,49	60,94	1,73	71,02	1,96
49,49	0,78	3,54	1,03	29,92	1,26	48,01	1,50	61,18	1,73	71,21	1,97
48,34	0,79	4,22	1,03	30,38	1,27	48,33	1,50	61,43	1,74	71,40	1,97
47,22	0,79	4,90	1,04	30,83	1,27	48,65	1,51	61,67	1,74	71,58	1,98
46,10	0,80	5,58	1,04	31,27	1,28	48,97	1,51	61,91	1,75	71,77	1,98
43,91	0,81	6,24	1,05	31,72	1,28	49,29	1,52	62,14	1,75	71,95	1,99
42,84	0,81	6,90	1,05	32,16	1,29	49,61	1,52	62,38	1,76	72,14	1,99
41,78	0,82	7,56	1,06	32,60	1,29	49,92	1,53	62,61	1,76	72,32	2,00
40,73	0,82	8,21	1,06	33,03	1,30	50,23	1,53	62,85	1,77	72,50	2,00
39,70	0,83	8,85	1,07	33,46	1,30	50,54	1,54	63,08	1,77		
38,67	0,83	9,49	1,07	33,89	1,31	50,84	1,54	63,31	1,78		
37,66	0,84	10,12	1,08	34,31	1,31	51,15	1,55	63,54	1,78		

Normality and Molarity chart for common volumetric solutions

Volumetric solution	Normality	Molarity
Acetic acid	0.01 N	0.01 M
Acetic acid	0.1 N	0.1 M
Ammonium thiocyanate	0.01 N	0.01 M
Ammonium thiocyanate	0.1 N	0.1 M
Bromine	0.1 N	0.05 M
Cerium (IV) sulphate	0.1 N	0.1 M
EDTA sodium salt	0.2 N	0.1 M
EDTA sodium salt	0.1 N	0.05 M
EDTA sodium salt	0.02 N	0.01 M
Hydrochloric acid	0.01 N	0.01 M
Hydrochloric acid	0.1 N	0.1 M
Hydrochloric acid	0.5 N	0.5 M
Hydrochloric acid	1 N	1 M
Hydrochloric acid	2 N	2 M
Iodine	0.01 N	0.005 M
Iodine	0.1 N	0.05 M
Iodine	1 N	0.5 M
Litium methoxide	0.1 N	0.1 M
Mercuric perchlorate	0.01 N	0.01 N
Nitric acid	0.1 N	0.1 M
Nitric acid	1 N	1 M
Oxalic acid	0.01 N	0.005 M
Oxalic acid	0.1 N	0.05 M
Oxalic acid	1 N	0.5 M
Perchloric acid	0.01 N	0.01 M
Perchloric acid	0.1 N	0.1 M
Potassium bromate	0.1 N	0.0167 M
Potassium dichromate	0.1 N	0.0167 M
Potassium hydr. phtalate	0.1 N	0.1 M
Potassium hydroxide	0.1 N	0.1 M
Potassium hydroxide	0.25 N	0.25 M
Potassium hydroxide	0.5 N	0.5 M
Potassium hydroxide	1 N	1 M
Potassium iodate	0.01 N	0.00167 M
Potassium iodate	0.1 N	0.0167 M
Potassium permanganate	0.01 N	0.002 M
Potassium permanganate	0.1 N	0.02 M
Potassium permanganate	1 N	0.2 M
Potassium thiocyanate	0.1 N	0.1 M
Silver nitrate	0.01 N	0.01 M
Silver nitrate	0.1 N	0.1 M
Silver nitrate	1 N	1 M
Sodium arsenite	0.1 N	0.05 M
Sodium carbonate	0.1 N	0.05 M
Sodium carbonate	1 N	0.5 M
Sodium chloride	0.1 N	0.1 M
Sodium hydroxide	0.01 N	0.01 M
Sodium hydroxide	0.1 N	0.1 M
Sodium hydroxide	0.25 N	0.25 M
Sodium hydroxide	1/2.82 N	1/2.82 M (0.357 mol/l)
Sodium hydroxide	0.5 N	0.5 M
Sodium hydroxide	1 N	1 M
Sodium hydroxide	2 N	2 M
Sodium thiosulphate	0.01 N	0.01 M
Sodium thiosulphate	0.1 N	0.1 M
Sodium thiosulphate	1 N	1 M
Sulphuric acid	0.01 N	0.005 M
Sulphuric acid	0.1 N	0.05 M
Sulphuric acid	0.5 N	0.25 M
Sulphuric acid	1 N	0.5 M
Sulphuric acid	2 N	1 M

Mesh size conversion table

Mesh size	Micron size approximate	Millimeters approximate	Inches
4	4760	4,760	0,185
6	3360	3,360	0,131
8	2380	2,380	0,093
12	1680	1,680	0,065
16	1190	1,190	0,046
20	840	0,840	0,0328
30	590	0,590	0,0232
40	420	0,420	0,0164
50	297	0,297	0,0116
60	250	0,250	0,0097
70	210	0,210	0,0082
80	177	0,177	0,0069
100	149	0,149	0,0058
140	105	0,105	0,0041
200	74	0,074	0,0029
230	62	0,062	0,0024
270	53	0,053	0,0021
325	44	0,044	0,0017
400	37	0,037	0,0015
625	20	0,020	0,0008
1250	10	0,010	0,0004
2500	5	0,005	0,0002

SAFETY IN THE LABORATORY

Alterable chemical products

The expiration date of all of our reagents is printed on both the label and the certificate of analysis. It applies to products stored in their original and intact packaging and away from heat and light as specified in the safety data sheet.

Once that an alterable product has been opened, the final user should determine the expiry date of the product basis on a risk analysis that includes the following parameters:

- Systematic risk
- Chemical risk
- Utilisation risk

Our range of products includes some alterable chemical products that, due to their own chemical properties, may be subject to an alteration during the time.

We indicate here below with some abbreviations the most common types of alteration found on our alterable products.

A: Alterable molecule
C: Colour change
D: Deliquescent
F: Interaction with the content
ID: Hydrolisable
IG: Hygroscope
M: Formation of an insoluble precipation
O: Oxidation
P: Polymerisation

Products	Alterability
Diethylene Glycol dimethylether	O
Diethylene Glycol monobutylether	O
Diethylene Glycol monoethylether	O
Diethylether not stabilized	O
Diethylsulfate	ID-C
Dihydroxyacetone	ID-D
Dimethylaminonaphtalene-5 -sulfonyl chloride	ID
Dimethylsulfate	ID-C
Diphenylamine	C
Diphenyldithiocarbazide	O-C
Ergometrine maleate	A
Ergometrine tartrate	A
Ethyl formate	ID
Ethyl-5-methylpyridine-2	C
Ethylaniline	C
Ethylchloroacetate	O
Ethylene Glycol dimethylether	O
Ethylene Glycol monoethylether	O
Ethylene Glycol monoethylether acetate	O
Ethylene Glycol monomethylether	O
Formaldehyde 40% w/v	O-P
Formic acid 99%	O
Furan	C-O
Furfural	C
Furfuryl alcohol	C
Gaiacol	C
Hanus's reagent	A
Heptanal	O
Hydrazine hydroxyde	A
Hydrazine solution	A
Hydrogen peroxide	A
Hydrogen sulfide saturated solution	O-M
Hydroquinone monomethylether	O

Products	Alterability
Iodine trichloride	F
Karl Fisher reagent	IG
Lutidine (2,4)	C
Magnesium perchlorate	IG
Magnesium peroxide	A
Methyl isobutylacetone	C
Methyldichloroacetate	ID
mono-Ethanolamine	C
n,n-Diethylaniline	C
n,n-Dimethylphenylenediamine	
n-Ethyl piperidine	C
Nitric acid fuming 90%	F
n-Methylaniline	C
Orthophosphoric acid 99%	O
Orthophosphoric acid 99%	IG
p-Dimethylaminobenzaldehyde	C
Phenol	C
Phosphomolybdic acid	C
Phosphorous pentachloride	F
Phosphorus pentoxide	IG
Picoline	C
Piperidine	C
Piruvic acid	P
Potassium ethyl xantogenate	ID
Potassium metabisulfite	O
p-Oxalate	C
p-phenetidin	C
Propionaldehyde	O
Propionitrile	C
Protonaldehyde	C-O
p-Toluidine	C
Pyridine hydrochloride	IG
Pyrrole	C
Pyrrolidine	C

Products	Alterability
Salicylaldehyde	O-C
Silver diethylcarbamate	C
Soda lime	A
Sodium acetate anhydrous	IG
Sodium citrate tribasic anhydrous	IG
Sodium cyanide solution	A
Sodium hydrate and hypochloride solution	A-O
Sodium hydrosulfite	O
Sodium hypochloride solution	A-O
Sodium metabisulfite	O
Sodium methylate	ID
Sodium sulfide nonahydrate	C-IG
Starch solution 1%	M
Strontium iodide	C
Styrene	P
Succinic anhydride	ID
Sulfuric acid 96%	C
Sulphurous acid saturated solution	O
sym-Diphenylcarbazide	C
Sym-Tetrabromoethane	A
Sym-Tetrachloroethane	A
Tetrahydrofurfuryl alcohol	C-D
Tin chloride anhydrous	A
Titanium trichloride 15%	M
Trichloroacetic acid solution 20%	A
Triphenylchloromethane	F
Vitamin A acetate	A
Vitamin A palmitate	A
Water chlorine saturated solution	A
Zinc oxide	IG

Incompatible chemical products

ACETIC ACID	Nitric acid Perchloric acid Alcohols Chromium oxide Ethylen glycol Permanganates Peroxides
HYDROFLUORIC ACID	Ammonia
NITRIC ACID CONCENTRATED	Acetic acid Hydrocyanic acid Aniline Chromium oxide Hydrogen sulphide Flammable liquids and gasses
OXALIC ACID	Silver Mercury
PERCHLORIC ACID	Alcohols Acetic anhydride Bismuth and its alloys Paper Wood
SULPHURIC ACID	Potassium chlorate Potassium perchlorate Potassium permanganate
AMMONIA	Hydrofluoric acid Bromine Chlorine Iodine Calcium hypoclorite Mercury
AMMONIUM NITRATE	Acids Chlorates Flammable liquids Nitrates Metallic powders Flammable organic substances Sulphur
ANILINE	Nitric acid Hydrogen peroxide
SILVER	Acetiyene Oxalic acid Tartaric acid Ammonium salts
BROMINE AND CHLORINE	Acetylene Ammonia Benzene Petroleum ether Butadiene Butane Hydrogen Methane Propane Metallic powders
CYANIDES	Acids
CHLORATES	Acidi Metallic powders Ammonium salts Flammable organic substances Sulphur
CHROMIUM OXIDE	Acetic acid Alcohols Petroleum ether Canphor Glycerol Flammable liquids Naphthalin
PHOSPHORUS	Oxigenated combinations Sulpur
HYDROCARBONS	Bromine Chlorine Chromium oxide Fluorine Sodium peroxide
IODINE	Acetylene Ammonia
FLAMMABLE LIQUIDS	Nitric acid Halogens Ammonium nitrate Chromium oxide Hydrogen peroxide Sodium peroxide
MERCURY	Acetylene Ammonia
ALKALY METALS	Water Halogens Carbon dioxide Carbon tetrachloride other Halogenated alkanes
HYDROGEN PEROXIDE	Acetone Alcohols Aniline Chromium Copper Iridium Metals and metallic salts Nitromethane Organic substances Flammable substances
POTASSIUM PERMANGANATE	Sulphuric acid Benzaldehyde Ethylen glycol Glycerol
COPPER	Acetilene Hydrogen peroxide
SODIUM PEROXIDE	Acetic acid Ethyl alcohol Methyl alcohol Acetic anhydride Benzaldehyde Ethyil acetate Ethylen glycol Furfurol Carbon disulfide

HAZARD STATEMENTS

H200	Unstable explosives.
H201	Explosive; mass explosion hazard.
H202	Explosive, severe projection hazard.
H203	Explosive; fire, blast or projection hazard.
H204	Fire or projection hazard.
H205	May mass explode in fire.
H220	Extremely flammable gas.
H221	Flammable gas.
H222	Extremely flammable aerosol.
H223	Flammable aerosol.
H224	Extremely flammable liquid and vapour.
H225	Highly flammable liquid and vapour.
H226	Flammable liquid and vapour.
H228	Flammable solid.
H240	Heating may cause an explosion.
H241	Heating may cause a fire or explosion.
H242	Heating may cause a fire.
H250	Catches fire spontaneously if exposed to air.
H251	Self-heating: may catch fire.
H252	Self-heating in large quantities; may catch fire.
H260	In contact with water releases flammable gases which may ignite spontaneously.
H261	In contact with water releases flammable gases.
H270	May cause or intensify fire; oxidiser.
H271	May cause fire or explosion; strong oxidiser.
H272	May intensify fire; oxidiser.
H280	Contains gas under pressure; may explode if heated.
H281	Contains refrigerated gas; may cause cryogenic burns or injury.
H290	May be corrosive to metals.
H300	Fatal if swallowed.
H301	Toxic if swallowed.
H302	Harmful if swallowed.
H304	May be fatal if swallowed and enters airways.
H310	Fatal in contact with skin.
H311	Toxic in contact with skin
H312	Harmful in contact with skin.
H314	Causes severe skin burns and eye damage.
H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H318	Causes serious eye damage.
H319	Causes serious eye irritation.
H330	Fatal if inhaled.
H331	Toxic if inhaled.
H332	Harmful if inhaled.
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H335	May cause respiratory irritation.
H336	May cause drowsiness or dizziness.
H340	May cause genetic defects.
H341	Suspected of causing genetic defects.
H350	May cause cancer.
H351	Suspected of causing cancer.
H360	May damage fertility or the unborn child.
H361	Suspected of damaging fertility or the unborn child.
H362	May cause harm to breast-fed children.
H370	Causes damage to organs.
H371	May cause damage to organs.
H372	Causes damage to organs through prolonged or repeated exposure exposure cause the hazard.
H373	May cause damage to organs through prolonged or repeated exposure exposure cause the hazard.
H400	Very toxic to aquatic life.
H410	Very toxic to aquatic life with long lasting effects.
H411	Toxic to aquatic life with long lasting effects.
H412	Harmful to aquatic life with long lasting effects.
H413	May cause long lasting harmful effects to aquatic life.

PRECAUTIONARY STATEMENTS

P101	If medical advice is needed, have product container or label at hand.
P102	Keep out of reach of children.
P103	Read label before use.
P201	Obtain special instructions before use.
P202	Do not handle until all safety precautions have been read and understood.
P210	Keep away from heat/sparks/open flames/hot surfaces. - No smoking.
P211	Do not spray on an open flame or other ignition source.
P220	Keep/Store away from clothing/... /combustible materials.
P221	Take any precaution to avoid mixing with combustibles...
P222	Do not allow contact with air.
P223	Keep away from any possible contact with water, because of violent reaction and possible flash fire.
P230	Keep wetted with...
P231	Handle under inert gas.
P232	Protect from moisture.
P231+P232	Handle under inert gas. Protect from moisture.
P233	Keep container tightly closed.
P234	Keep only in original container.
P235	Keep cool.
P235+P410	Keep cool. Protect from sunlight.
P240	Ground/bond container and receiving equipment.
P241	Use explosion-proof electrical/ventilating/lighting/.../equipment.
P242	Use only non-sparking tools.
P243	Take precautionary measures against static discharge.
P244	Keep reduction valves free from grease and oil.
P250	Do not subject to grinding/shock/.../friction.
P251	Pressurized container: Do not pierce or burn, even after use.
P260	Do not breathe dust/fume/gas/mist/vapours/spray.
P261	Avoid breathing dust/fume/gas/mist/vapours/spray.
P262	Do not get in eyes, on skin, or on clothing.
P263	Avoid contact during pregnancy/while nursing.
P264	Wash ... thoroughly after handling.
P270	Do no eat, drink or smoke when using this product.
P271	Use only outdoors or in a well-ventilated area.
P272	Contaminated work clothing should not be allowed out of the workplace.
P273	Avoid release to the environment.
P280	Wear protective gloves/protective clothing/eye protection/face protection.
P281	Use personal protective equipment as required.
P282	Wear cold insulating gloves/face shield/eye protection.
P283	Wear fire/flame resistant/retardant clothing.
P284	Wear respiratory protection.
P285	In case of inadequate ventilation wear respiratory protection.
P301	IF SWALLOWED.
P301+P310	IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician.
P301+P312	IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.
P301+P330+P331	IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
P302	IF ON SKIN.
P302+P334	IF ON SKIN: Gently wash with plenty of soap and water.
P302+P350	IF ON SKIN: Wash with plenty of soap and water.
P302+P352	IF ON SKIN: Immerse in cool water/wrap in wet bandages.
P303	IF ON SKIN (or hair).
P303+P361+P353	IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.
P304	IF INHALED.
P304+P340	IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
P304+P341	IF INHALED: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing.
P305	IF IN EYES.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P306	IF ON CLOTHING.
P306+P360	IF ON CLOTHING: rinse immediately contaminated clothing and skin with plenty of water before removing clothes.
P307	IF exposed
P307+P311	IF exposed: Call a POISON CENTER or doctor/physician.


P308	IF exposed or concerned.
P308+P313	IF exposed or concerned: Get medical advice/attention.
P309	IF exposed or if you feel unwell.
P309+P311	IF exposed or if you feel unwell: Call a POISON CENTER or doctor/physician.
P310	Immediately call a POISON CENTER or doctor/physician.
P311	Call a POISON CENTER or doctor/physician.
P312	Call a POISON CENTER or doctor/physician if you feel unwell.
P313	Get medical advice/attention.
P314	Get medical advice/attention if you feel unwell.
P315	Get immediate medical advice/attention.
P320	Specific treatment is urgent (see... on this label).
P321	Specific treatment (se ... on this label).
P322	Specific measures (see... on this label).
P330	Rinse mouth.
P331	Do NOT induce vomiting.
P332	If skin irritation occurs.
P332+P313	If skin irritation occurs: Get medical advice/attention.
P333	If skin irritation or rash occurs.
P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
P334	Immerse in cool water/wrap in wet bandages.
P335	Brush off loose particles from skin.
P335+P334	Brush off loose particles from skin. Immerse in cool water/wrap in wet bandages.
P336	Thaw frosted parts with lukewarm water. Do no rub affected area.
P337	If eye irritation persists.
P337+P313	If eye irritation persists: Get medical advice/attention.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.
P340	Remove victim to fresh air and keep at rest in a position comfortable for breathing.
P341	If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing.
P342	If experiencing respiratory symptoms.
P342+P311	If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician.
P350	Gently wash with plenty of soap and water.
P351	Rinse cautiously with water for several minutes.
P352	Wash with plenty of soap and water.
P353	Rinse skin with water/shower.
P360	Rinse immediately contaminated clothing and skin with plenty of water before removing clothes.
P361	Remove/Take off immediately all contaminated clothing.
P362	Take off contaminated clothing and wash before reuse.
P363	Wash contaminated clothing before reuse.
P370	In case of fire.
P370+P376	In case of fire: Stop leak if safe to do so.
P370+P378	In case of fire: Use ... for extinction.
P370+P380	In case of fire: Evacuate area.
P370+P380+P375	In case of fire: Evacuate area. Fight fire remotely due to the risk of explosion.
P371	In case of major fire and large quantities.
P371+P380+P375	In case of major fire and large quantities: Evacuate area. Fight fire remotely due to the risk of explosion.
P372	Explosion risk in case of fire.
P373	DO NOT fight fire when fire reaches explosives.
P374	Fight fire with normal precautions from a reasonable distance.
P375	Fight fire remotely due to the risk of explosion.
P376	Stop leak if safe to do so.
P377	Leaking gas fire: Do not extinguish, unless leak can be stopped safely.
P378	Use ... for extinction.
P380	Evacuate area.
P381	Eliminate all ignition sources if safe to do so.
P390	Absorb spillage to prevent material damage.
P391	Collect spillage.
P401	Store ...
P402	Store in a dry place.
P402+P404	Store in a dry place. Store in a closed container.
P403	Store in a well-ventilated place.
P403+P233	Store in a well-ventilated place. Keep container tightly closed.
P403+P235	Store in a well-ventilated place. Keep cool.
P404	Store in a closed container.
P405	Store locked up.

P406	Store in corrosive resistant/... container with a resistant inner liner.
P407	Maintain air gap between stacks/pallets.
P410	Protect from sunlight.
P410+P412	Protect from sunlight. Do no expose to temperatures exceeding 50 °C/122 °F.
P411	Store at temperatures not exceeding ... °C.
P411+P235	Store at temperatures not exceeding ... °C.
P412	Do not expose to temperatures exceeding 50 °C/122 °F.
P413	Store bulk masses greater than ... kg.
P420	Store away from other materials.
P422	Store contents under ...
P501	Dispose of contents/container to ...

ADDITIONAL STATEMENTS

EUH 001	Explosive when dry.
EUH 006	Explosive with or without contact with air.
EUH 014	Reacts violently with water.
EUH 018	In use, may form flammable/explosive vapour-air mixture.
EUH 019	May form explosive peroxides.
EUH 029	Contact with water liberates toxic gas .
EUH 031	Contact with acids liberates toxic gas.
EUH 032	Contact with acids liberates very toxic gas.
EUH 044	Risk of explosion if heated under confinement.
EUH 059	Hazardous to the ozone layer.
EUH 066	Repeated exposure may cause skin dryness or cracking.
EUH 070	Toxic by eye contact.
EUH 071	Corrosive to the respiratory tract.
EUH 201	Contains lead. Should not be used on surfaces liable to be chewed or sucked by children.
EUH 201A	Warning! Contains lead.
EUH 202	Cyanoacrylate. Danger. Bonds skin and eyes in seconds. Keep out of the reach of children.
EUH 203	Contains chromium (VI). May produce an allergic reaction.
EUH 204	Contains isocyanates. May produce an allergic reaction.
EUH 205	Contains epoxy constituents. May produce an allergic reaction.
EUH 206	Warning! Do not use together with other products. May release dangerous gases (chlorine).
EUH 207	Warning! Contains cadmium. Dangerous fumes are formed during use. See information supplied by the manufacturer. Comply with the safety instructions.
EUH 208	Contains (name of sensitising substance). May produce an allergic reaction.
EUH 209	Can become highly flammable in use.
EUH 209A	Can become flammable in use.
EUH 201	Safety data sheet available on request.
EUH 401	To avoid risks to human health and the environment, comply with the instructions for use.

INDEXES



INDEX BY PRODUCT CODE

300101	Absorbent for split liquids	129
300102	Absorbent for split liquids	129
300857	p-Acetylphenetidine	146
301251	Acetophenone	146
301356	Chlorobutanol hemihydrate	259
301357	Chlorobutanol hemihydrate	259
301501	Acetone	140
301502	Acetone	140
301503	Acetone	140
301504	Acetone	140
301505	Acetone	140
301506	Acetone	140
301507	Acetone	140
301852	Acetic acid 80%	134
301853	Acetic acid 80%	134
301855	Acetic acid 80%	134
302002	Acetic acid 96%	133
302003	Acetic acid 96%	133
302005	Acetic acid 96%	133
302007	Acetic acid 96%	133
302011	Acetic acid glacial	133
302013	Acetic acid glacial	133
302014	Acetic acid glacial	133
302015	Acetic acid glacial	133
302016	Acetic acid glacial	133
302021	Acetic acid glacial	133
302031	Acetic acid glacial	133
302032	Acetic acid glacial	133
302033	Acetic acid glacial	133
302034	Acetic acid glacial	133
302082	Benzoic acid	199
302087	Benzoic acid	199
302089	Benzoic acid	199
302177	Boric acid	206
302178	Boric acid	206
302179	Boric acid	206
302185	Boric acid	205
302484	Citric acid anhydrous	269
302485	Citric acid anhydrous	269
302486	Citric acid anhydrous	269
302487	Citric acid anhydrous	269
302488	Citric acid anhydrous	269
302501	Citric acid monohydrate	270
302504	Citric acid monohydrate	270
302507	Citric acid monohydrate	270
302509	Citric acid monohydrate	270
302551	Citric acid monohydrate	269
302554	Citric acid monohydrate	269
302557	Citric acid monohydrate	269
302559	Citric acid monohydrate	269
302575	Kerosene	449
302591	Hydrochloric acid 10%	400
302601	Hydrochloric acid 32%	397
302602	Hydrochloric acid 32%	397
302604	Hydrochloric acid 32%	397
302621	Hydrochloric acid 37%	395
302622	Hydrochloric acid 37%	395
302623	Hydrochloric acid 37%	395
302624	Hydrochloric acid 37%	395
302625	Hydrochloric acid 37%	395
302626	Hydrochloric acid 37%	395
302627	Hydrochloric acid 37%	395
302643	Hydrochloric acid 37%	395
302664	Hydrochloric acid 32% (20°Bé)	397
303151	Dichloroacetic acid	297
303201	Ethylenediaminetetraacetic acid disodium salt	345
303202	Ethylenediaminetetraacetic acid disodium salt	345
303203	Ethylenediaminetetraacetic acid disodium salt	345
303225	Ethylenediaminetetraacetic acid disodium salt	346
303226	Ethylenediaminetetraacetic acid disodium salt	346
303227	Ethylenediaminetetraacetic acid disodium salt	346
303251	Ethylenediaminetetraacetic acid	344
303252	Ethylenediaminetetraacetic acid	344
303731	Hydrofluoric acid 39.5%	412
303871	Hydrochloric acid 37%	395
303901	Formic acid 85%	363
303905	Formic acid 85%	363
303911	Formic acid 99%	363
303912	Formic acid 99%	363
303913	Formic acid 99%	363
304051	Orthophosphoric acid 75%	549
304052	Orthophosphoric acid 75%	549
304054	Orthophosphoric acid 75%	549
304061	Orthophosphoric acid 85%	548
304062	Orthophosphoric acid 85%	548
304063	Orthophosphoric acid 85%	548
304201	Gallic acid monohydrate	367
304505	L(+)Glutamic acid	371
304507	L(+)Glutamic acid	371
304651	L(+)-Lactic acid	452
304652	L(+)-Lactic acid	452
304653	L(+)-Lactic acid	452
305201	Nitric acid 65%	532
305202	Nitric acid 65%	532
305205	Nitric acid 65%	532
305207	Nitric acid 65%	532
305501	Nitric acid 67.5 (42° Be)	530
305502	Nitric acid 67.5 (42° Be)	530
305505	Nitric acid 67.5 (42° Be)	530
305701	Oleic acid	539
305704	Oleic acid	539
305757	Oxalic acid dihydrate	551
305758	Oxalic acid dihydrate	551
306091	Perchloric acid 65%	561
306254	Propionic acid	628
306377	Salicylic acid	643
306381	Salicylic acid	643
306503	Sulfamic acid	727
306507	Sulfamic acid	727
306508	Sulfamic acid	727
306651	Sulfuric acid 96%	733
306653	Sulfuric acid 96%	733
306657	Sulfuric acid 96%	733
306751	Sulfuric acid 96% (66°Be)	731
306752	Sulfuric acid 96% (66°Be)	731
306755	Sulfuric acid 96% (66°Be)	731
307001000	Sulfuric acid 35% (30°Be)	736
307112	Stearic acid	722
307115	Stearic acid	722
307152	Tannic acid	744
307153	Tannic acid	744
307157	Tannic acid	744
307301	L(+) Tartaric Acid	745
307307	L(+) Tartaric Acid	745
307309	L(+) Tartaric Acid	745
307357	L(+) Tartaric Acid	745
307359	L(+) Tartaric Acid	745
307508	p-Toluenesulfonic acid	769
307557	Trichloroacetic acid	770
307558	Trichloroacetic acid	770
307581	Water	788
307582	Water	788
307583	Water	788
307584	Water	788
307585	Water	788
307586	Water	788
307587	Water	788
307589	Water	788
307591	Water	788
307592	Water	788
307593	Water	788
307601	Water purified	788
307602	Water purified	788
307603	Water purified	788
307604	Water purified	788
307606	Water purified	788
307641	Hydrogen peroxide solution 6%	415
307671	Hydrogen peroxide solution 3%	416
307678	Hydrogen peroxide solution 3%	416
307685	Hydrogen peroxide solution 30%	415
307701	Hydrogen peroxide solution 40% w/v	413
307708	Hydrogen peroxide solution 40% w/v	413
307709	Hydrogen peroxide solution 40% w/v	413
307742	Hydrogen peroxide solution 35%	413
307901	n-Amyl alcohol	180
308001	Isoamyl alcohol	437
308003	Isoamyl alcohol	437
308131	Benzyl alcohol	200
308132	Benzyl alcohol	200
308137	Benzyl alcohol	200
308138	Benzyl alcohol	200
308251	Butanol-1	226
308257	Butanol-1	226
308259	Butanol-1	226
308301	Isobutanol	437
308303	Isobutanol	437
308357	Cetyl alcohol	255
308358	Cetyl alcohol	255
308359	Cetyl alcohol	255
308600	Ethanol absolute anhydrous	333
308601	Ethanol absolute anhydrous	333
308602	Ethanol absolute anhydrous	333
3086022	Ethanol absolute anhydrous	333
308603	Ethanol absolute anhydrous	333
3086032	Ethanol absolute anhydrous	333
308604	Ethanol absolute anhydrous	333
308605	Ethanol absolute anhydrous	333
3086052	Ethanol absolute anhydrous	333
308607	Ethanol absolute anhydrous	333
3086072	Ethanol absolute anhydrous	333
308608	Ethanol absolute anhydrous	333
308609	Ethanol absolute anhydrous	333
3086092	Ethanol absolute anhydrous	333
308621	Ethanol 94° denaturated	339
308623	Ethanol 94° denaturated	339
308624	Ethanol 94° denaturated	339
308625	Ethanol 94° denaturated	339
308641	Ethanol 96°	336
3086412	Ethanol 96°	336
308643	Ethanol 96°	336
308644	Ethanol 96°	336
3086442	Ethanol 96°	336
308645	Ethanol 96°	336
3086452	Ethanol 96°	336
308646	Ethanol 96°	336
3086462	Ethanol 96°	336
308647	Ethanol 96°	336
3086472	Ethanol 96°	336
308648	Ethanol 96°	336
3086482	Ethanol 96°	336
308649	Ethanol 96°	336
3086492	Ethanol 96°	336
308651	Ethanol absolute denaturated	338
308652	Ethanol absolute denaturated	338
308653	Ethanol absolute denaturated	338
308655	Ethanol absolute denaturated	338
308656	Ethanol absolute denaturated	338
308661	Ethanol absolute anhydrous	333

Product specifications are subject to changes. Please visit our website for updates.

INDEX BY PRODUCT CODE

3086612	Ethanol absolute anhydrous	333
308662	Ethanol absolute anhydrous	333
3086622	Ethanol absolute anhydrous	333
308663	Ethanol absolute anhydrous	333
308664	Ethanol absolute anhydrous	333
3086642	Ethanol absolute anhydrous	333
308667	Ethanol absolute anhydrous	333
308681	Ethanol 90° denaturated	339
308682	Ethanol 90° denaturated	339
308683	Ethanol 90° denaturated	339
308687	Ethanol 90° denaturated	339
308731	2-Phenylethanol	573
308771	Ethanol 70% v/v	337
308775	Ethanol 70% v/v	337
3087752	Ethanol 70% v/v	337
309000	Methanol	492
309001	Methanol	492
309002	Methanol	492
309004	Methanol	492
309008	Methanol	492
309009	Methanol	492
309201	Methanol	492
309203	Methanol	492
309204	Methanol	492
309351	Propan-1-ol	623
309352	Propan-1-ol	623
309353	Propan-1-ol	623
309354	Propan-1-ol	623
309358	Propan-1-ol	623
309500	Propan-2-ol	626
309501	Propan-2-ol	626
309504	Propan-2-ol	626
309505	Propan-2-ol	626
309506	Propan-2-ol	626
309507	Propan-2-ol	626
309509	Propan-2-ol	626
310348	Formaldehyde 35% w/w	359
310349	Formaldehyde 35% w/w	359
310351	Formaldehyde 35% w/w	359
310355	Formaldehyde 35% w/w	359
310356	Formaldehyde 35% w/w	359
310358	Formaldehyde 35% w/w	359
310504	Propionaldehyde	628
311002	Aluminum ammonium sulfate dodecahydrate	152
311009	Aluminum ammonium sulfate dodecahydrate	152
311252	Aluminum chloride hexahydrate	153
311254	Aluminum chloride hexahydrate	153
311256	Aluminum chloride hexahydrate	153
311257	Aluminum chloride hexahydrate	153
311734	Aluminum hydroxide	154
312001	Aluminum nitrate nonahydrate	154
312007	Aluminum nitrate nonahydrate	154
312008	Aluminum nitrate nonahydrate	154
312252	Aluminum oxide	155
312258	Aluminum oxide	155
312259	Aluminum oxide	155
312261	Aluminum oxide activated	156
312401	Aluminum potassium sulfate dodecahydrate	156
312402	Aluminum potassium sulfate dodecahydrate	156
312404	Aluminum potassium sulfate dodecahydrate	156
312508	Aluminum potassium sulfate dodecahydrate	157
312751	Aluminum sulfate	157
312752	Aluminum sulfate	157
312753	Aluminum sulfate	157
313071	Maize starch	476
313072	Maize starch	476
313073	Maize starch	476
313102	Rice starch	640
313107	Rice starch	640
313108	Rice starch	640
313109	Rice starch	640
313251	Isoamyl acetate	436
313252	Isoamyl acetate	436
313431	Tris (hydroxymethyl)-aminomethane	777
313432	Tris (hydroxymethyl)-aminomethane	777
313441	Tris (hydroxymethyl)-aminomethane	777
313502	Ammonium acetate	166
313504	Ammonium acetate	166
313507	Ammonium acetate	166
313508	Ammonium acetate	166
313601	Ammonium bicarbonate	167
313895	tri-Ammonium citrate	170
313951	Ammonium chloride	169
313952	Ammonium chloride	169
313954	Ammonium chloride	169
313956	Ammonium chloride	169
313957	Ammonium chloride	169
314001	Ammonium chloride	169
314002	Ammonium chloride	169
314261	Ammonium hydrogen difluoride	172
314504	Ammonium phosphate monobasic	176
314505	Ammonium phosphate monobasic	176
314506	Ammonium phosphate monobasic	176
314507	Ammonium phosphate monobasic	176
314757	Ammonium phosphate dibasic	176
314758	Ammonium phosphate dibasic	176
314861	Ammonia solution 28%	161
314863	Ammonia solution 28%	161
314866	Ammonia solution 28%	161
314871	Ammonia solution 30%	161
314873	Ammonia solution 30%	161
315502	Ammonium nitrate	174
315509	Ammonium nitrate	174
316002	Ammonium persulfate	175
316008	Ammonium persulfate	175
316251	Ammonium sulfate	177
316252	Ammonium sulfate	177
316257	Ammonium sulfate	177
316303	Ammonium thiocyanate	178
316307	Ammonium thiocyanate	178
316308	Ammonium thiocyanate	178
316501	Acetic anhydride	137
316502	Acetic anhydride	137
316503	Acetic anhydride	137
317511	Chromium (VI) oxide solution	267
317753	Phosphorus pentoxide	578
318007	Phthalic anhydride	579
318507	Succinic anhydride	725
319502	Silver acetate	654
319507	Silver acetate	654
320002	Silver carbonate	654
320007	Silver carbonate	654
320502	Silver chloride	654
320504	Silver chloride	654
320904	Silver nitrate	655
320907	Silver nitrate	655
321502	Barium carbonate	192
321507	Barium carbonate	192
321553	Canada balsam	247
321554	Canada balsam	247
321752	Barium chloride dihydrate	193
321757	Barium chloride dihydrate	193
321758	Barium chloride dihydrate	193
322001	Barium hydroxide octahydrate	194
322004	Barium hydroxide octahydrate	194
322007	Barium hydroxide octahydrate	194
322009	Barium hydroxide octahydrate	194
322607	Barium sulfate	196
322721	Benzalkonium chloride	197
322737	Benzalkonium chloride	197
322738	Benzalkonium chloride	197
323101	Benzyl benzoate	201
323102	Benzyl benzoate	201
323401	Petroleum ether 55 - 85°C	565
323402	Petroleum ether 55 - 85°C	565
323403	Petroleum ether 55 - 85°C	565
323501	Petroleum ether 80 - 100°C	565
323502	Petroleum ether 80 - 100°C	565
323503	Petroleum ether 80 - 100°C	565
324185	Bismuth(III) nitrate pentahydrate	204
325601	n-Butyl acetate	229
325602	n-Butyl acetate	229
325603	n-Butyl acetate	229
325604	n-Butyl acetate	229
325631	Isobutyl acetate	438
325633	Isobutyl acetate	438
325701	Di-n-butylphthalate	297
325741	Cadmium chloride monohydrate	233
326356	Caffeine anhydrous	234
326357	Caffeine anhydrous	234
326358	Caffeine anhydrous	234
326454	Calcium hydroxide	241
326458	Calcium hydroxide	241
326503	Calcium acetate anhydrous	236
326507	Calcium acetate anhydrous	236
326511	Calcium acetate anhydrous	236
326512	Calcium acetate anhydrous	236
326513	Calcium acetate anhydrous	236
327002	Calcium carbonate	237
327003	Calcium carbonate	237
327059	Calcium carbonate	237
327101	Calcium carbonate	237
327105	Calcium carbonate	237
327507	Calcium chloride hexahydrate	239
327509	Calcium chloride hexahydrate	239
327603	Calcium chloride dihydrate	239
327607	Calcium chloride dihydrate	239
327609	Calcium chloride dihydrate	239
328252	Calcium chloride anhydrous	238
328257	Calcium chloride anhydrous	238
328258	Calcium chloride anhydrous	238
328757	Calcium chloride anhydrous	238
328759	Calcium chloride anhydrous	238
328807	Calcium chloride anhydrous	238
328809	Calcium chloride anhydrous	238
330303	Calcium phosphate dibasic dihydrate	243
330307	Calcium phosphate dibasic dihydrate	243
330403	Calcium phosphate tribasic	244
330407	Calcium phosphate tribasic	244
330409	Calcium phosphate tribasic	244
330601	Calcium gluconate	240
330608	Calcium gluconate	240
330609	Calcium gluconate	240
331003	Calcium hydroxide	241
331007	Calcium hydroxide	241
331008	Calcium hydroxide	241
331407	Calcium lactate	241
331408	Calcium lactate	241
331501	Calcium nitrate tetrahydrate	242
331509	Calcium nitrate tetrahydrate	242
331555	Calcium oxide, lumps	242
331557	Calcium oxide, lumps	242
331564	Calcium oxide, powder	242
331567	Calcium oxide, powder	242
331602	Calcium pantothenate	242

INDEX BY PRODUCT CODE

331751	Calcium sulfate dihydrate	244
331752	Calcium sulfate dihydrate	244
331761	Calcium sulfate hemihydrate	245
331762	Calcium sulfate hemihydrate	245
331763	Calcium sulfate hemihydrate	245
332261	Calcium stearate	244
332262	Calcium stearate	244
332265	Calcium stearate	244
332356	Camphor natural	246
332401	Camphor synthetic	246
332402	Camphor synthetic	246
332406	Camphor synthetic	246
332573	Kaolin washed	443
332574	Kaolin washed	443
332658	Charcoal vegetable	256
332659	Charcoal vegetable	256
333701	Quinoline	634
333707	Quinoline	634
333751	Cyclohexane	287
333752	Cyclohexane	287
333753	Cyclohexane	287
333801	Cyclohexanol	287
333901	Cyclohexanone	288
333902	Cyclohexanone	288
333905	Cyclohexanone	288
334251	Chlorobenzene	259
334254	Chlorobenzene	259
334255	Chlorobenzene	259
334351	Chloroform	262
334353	Chloroform	262
334354	Chloroform	262
334356	Chloroform	262
337001	Diacetone alcohol	295
337002	Diacetone alcohol	295
337251	o-Dichlorobenzene	298
337303	p-Dichlorobenzene	298
337307	p-Dichlorobenzene	298
337315	Dichloromethane	304
337331	Dichloromethane	304
337333	Dichloromethane	304
337335	Dichloromethane	304
337337	Dichloromethane	304
337501	Diethylamine	306
337801	Diethanolamine	306
338001	1,4-Dioxane	321
338002	1,4-Dioxane	321
338003	1,4-Dioxane	321
338004	1,4-Dioxane	321
338005	1,4-Dioxane	321
338112	Diethyl phthalate	310
338113	Diethyl phthalate	310
338114	Diethyl phthalate	310
338115	Diethyl phthalate	310
338121	n,n-Diethyl-p-phenylenediamine sulfate	310
338124	n,n-Diethyl-p-phenylenediamine sulfate	310
339381	n-Heptane 99%	381
339382	n-Heptane 99%	381
339385	n-Heptane 99%	381
339386	n-Heptane 99%	381
339751	n-Hexane	387
339752	n-Hexane	387
339755	n-Hexane	387
339756	n-Hexane	387
339757	n-Hexane	387
339758	n-Hexane	387
339759	n-Hexane	387
339851	Hexane mixture of isomers	388
339852	Hexane mixture of isomers	388
339856	Hexane mixture of isomers	388
340151	1,2-Dichloroethane	300
340155	1,2-Dichloroethane	300
340731	Diethyl ether	309
340751	Diethyl ether	309
340752	Diethyl ether	309
340759	Diethyl ether	309
340762	Diethyl ether	310
340765	Diethyl ether	310
341022	Petroleum ether 40 - 70°C	566
341024	Petroleum ether 40 - 70°C	566
341032	Petroleum ether 30 - 50°C	568
341034	Petroleum ether 30 - 50°C	568
341502	Ethyl acetate	342
341503	Ethyl acetate	342
341506	Ethyl acetate	342
341511	Ethyl acetate	342
341751	Ethyl acetoacetate	343
342101	Ethyl formate	350
343001	Tetrachloroethylene	750
343003	Tetrachloroethylene	750
343407	Phenol	570
343411	Phenol liquified 85%	570
343441	Iron (III) ammonium citrate red	431
343442	Iron (III) ammonium citrate red	431
343443	Iron (III) ammonium citrate red	431
343605	Iron (III) ammonium citrate green	431
343606	Iron (III) ammonium citrate green	431
343607	Iron (III) ammonium citrate green	431
343757	Iron (III) ammonium oxalate	432
344003	Iron (II) ammonium sulfate hexahydrate	429
344007	Iron (II) ammonium sulfate hexahydrate	429
344008	Iron (II) ammonium sulfate hexahydrate	429
344107	Iron (III) ammonium sulfate dodecahydrate	432
344108	Iron (III) ammonium sulfate dodecahydrate	432
344201	Iron (III) citrate	434
344202	Iron (III) citrate	434
344504	Iron (III) chloride hexahydrate	434
344507	Iron (III) chloride hexahydrate	434
344508	Iron (III) chloride hexahydrate	434
344509	Iron (III) chloride hexahydrate	434
344957	Iron (II) sulfate heptahydrate	431
344959	Iron (II) sulfate heptahydrate	431
345355	Fluorescein sodium salt	357
345356	Fluorescein sodium salt	357
345357	Fluorescein sodium salt	357
346102	Glycerol (30°Bé)	373
346106	Glycerol (30°Bé)	373
346131	Glycerol 90% (28° Bé)	373
346132	Glycerol 90% (28° Bé)	373
346134	Glycerol 90% (28° Bé)	373
346161	Glycerol (30°Bé)	373
346162	Glycerol (30°Bé)	373
346164	Glycerol (30°Bé)	373
346165	Glycerol (30°Bé)	373
346167	Glycerol (30°Bé)	373
346205	Glycine	374
346207	Glycine	374
346208	Glycine	374
346301	Diethylene glycol	307
346303	Diethylene glycol	307
346304	Diethylene glycol	307
346501	Ethylene glycol	349
346502	Ethylene glycol	349
346503	Ethylene glycol	349
346504	Ethylene glycol	349
346508	Ethylene glycol	349
346509	Ethylene glycol	349
346701	Propylene glycol	630
346703	Propylene glycol	630
346705	Propylene glycol	630
346708	Propylene glycol	630
346971	D(+)Glucose monohydrate	371
346972	D(+)Glucose monohydrate	371
346973	D(+)Glucose monohydrate	371
346983	D(+)Glucose anhydrous	370
346987	D(+)Glucose anhydrous	370
346989	D(+)Glucose anhydrous	370
347107	Gum arabic	377
347355	Lanolin anhydrous	453
347357	Lanolin anhydrous	453
347359	Lanolin anhydrous	453
348124	Hydroquinone	416
348126	Hydroquinone	416
348128	Hydroquinone	416
348129	Hydroquinone	416
348354	Inositol	423
348451	Iodine	424
348452	Iodine	424
348454	Iodine	424
348455	Iodine	424
348457	Iodine	424
348554	Iodoform	427
348558	Iodoform	427
348702	Lactose monohydrate	453
348703	Lactose monohydrate	453
348707	Lactose monohydrate	453
348708	Lactose monohydrate	453
348754	Lecithin of soya	461
348854	Yeast dried	794
348901	Petroleum ether 80 - 120°C	564
348905	Petroleum ether 80 - 120°C	564
348912	Petroleum ether 100 - 140°C	564
348913	Petroleum ether 100 - 140°C	564
348914	Petroleum ether 100 - 140°C	564
348955	Lithium carbonate	463
348957	Lithium carbonate	463
349257	Magnesium carbonate basic	470
349272	Magnesium carbonate basic	470
349279	Magnesium carbonate basic	470
349355	Magnesium chloride hexahydrate	471
349357	Magnesium chloride hexahydrate	471
349359	Magnesium chloride hexahydrate	471
349372	Magnesium chloride hexahydrate	471
349377	Magnesium chloride hexahydrate	471
349407	Magnesium glycerophosphate	471
349455	Magnesium hydroxide	472
349557	Magnesium nitrate hexahydrate	472
349653	Magnesium oxide heavy	473
349655	Magnesium oxide heavy	473
349656	Magnesium oxide heavy	473
349753	Magnesium peroxide	474
349757	Magnesium peroxide	474
349851	Magnesium sulfate heptahydrate	475
349852	Magnesium sulfate heptahydrate	475
349859	Magnesium sulfate heptahydrate	475
350032	Magnesium stearate	474
350033	Magnesium stearate	474
350035	Magnesium stearate	474
350401	Maltose monohydrate	478
351502	Manganese (II) chloride tetrahydrate	480
351507	Manganese (II) chloride tetrahydrate	480
351508	Manganese (II) chloride tetrahydrate	480

DASIT GROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

INDEX BY PRODUCT CODE

352002	Manganese (II) sulfate monohydrate	481
352007	Manganese (II) sulfate monohydrate	481
352008	Manganese (II) sulfate monohydrate	481
352051	D-Mannitol	482
352052	D-Mannitol	482
352053	D-Mannitol	482
352103	L-Menthol	483
352106	L-Menthol	483
352654	Mercury (I) chloride	485
352657	Mercury (I) chloride	485
354007	Methyl 4-hydroxybenzoate	499
354008	Methyl 4-hydroxybenzoate	499
354152	Methyl salicylate	505
354155	Methyl salicylate	505
354251	Ethyl methyl ketone	350
354253	Ethyl methyl ketone	350
354254	Ethyl methyl ketone	350
354501	Dichloromethane	303
354751	1-Naphthol	517
355757	Nickel (II) sulfate hexahydrate	524
356251	Oil refined of almonds	539
356351	Castor oil	249
356352	Castor oil	249
356353	Castor oil	249
356601	Paraffin oil	555
356603	Paraffin oil	555
356607	Paraffin oil	555
356608	Paraffin oil	555
356661	n-Octane	537
356663	n-Octane	537
356951	n-Pentane	559
356952	n-Pentane	559
356953	n-Pentane	559
356954	n-Pentane	559
357151	Petroleum	563
357155	Petroleum	563
357253	Lead (II) acetate trihydrate	458
358007	Lead (II) nitrate	459
358008	Lead (II) nitrate	459
358252	Lead (II) oxide	460
358257	Lead (II) oxide	460
358259	Lead (II) oxide	460
358752	Pyridine	631
358754	Pyridine	631
358903	Potassium acetate	583
358907	Potassium acetate	583
358908	Potassium acetate	583
359702	Potassium bromide	586
359707	Potassium bromide	586
359803	Potassium carbonate	586
359808	Potassium carbonate	586
359809	Potassium carbonate	586
359956	Potassium citrate tribasic monohydrate	591
359957	Potassium citrate tribasic monohydrate	591
359958	Potassium citrate tribasic monohydrate	591
359959	Potassium citrate tribasic monohydrate	591
360104	Potassium chloride	587
360106	Potassium chloride	587
360107	Potassium chloride	587
360109	Potassium chloride	587
360252	Potassium ferricyanide	595
360257	Potassium ferricyanide	595
360258	Potassium ferricyanide	595
360552	Potassium ferrocyanide trihydrate	595
360557	Potassium ferrocyanide trihydrate	595
360558	Potassium ferrocyanide trihydrate	595
361503	Potassium phosphate monobasic	615
361507	Potassium phosphate monobasic	615
361509	Potassium phosphate monobasic	615
361751	Potassium phosphate dibasic anhydrous	614
361752	Potassium phosphate dibasic anhydrous	614
361757	Potassium phosphate dibasic anhydrous	614
362201	Potassium hydroxide, flakes	599
362202	Potassium hydroxide, flakes	599
362235	Potassium hydroxide, pellets	600
362237	Potassium hydroxide, pellets	600
362239	Potassium hydroxide, pellets	600
362251	Potassium hydroxide, flakes	599
362257	Potassium hydroxide, flakes	599
362258	Potassium hydroxide, flakes	599
362402	Potassium iodide	607
362403	Potassium iodide	607
362405	Potassium iodide	607
362407	Potassium iodide	607
362409	Potassium iodide	607
362622	Potassium metabisulfite	609
362623	Potassium metabisulfite	609
362627	Potassium metabisulfite	609
362629	Potassium metabisulfite	609
363002	Potassium nitrate	610
363007	Potassium nitrate	610
363009	Potassium nitrate	610
363101	Potassium permanganate	611
363107	Potassium permanganate	611
363109	Potassium permanganate	611
363455	Potassium sodium tartrate tetrahydrate	617
363457	Potassium sodium tartrate tetrahydrate	617
363459	Potassium sodium tartrate tetrahydrate	617
363602	Potassium sulfate	618
363607	Potassium sulfate	618
363608	Potassium sulfate	618
363752	Potassium thiocyanate	620
363756	Potassium thiocyanate	620
363807	Potassium guaiacolsulfonate	596
363884	Potassium sorbate	617
363907	Potassium L-tartrate monobasic	619
363953	Propyl p-hydroxybenzoate	629
363956	Propyl p-hydroxybenzoate	629
364007	Copper (II) acetate hydrate	278
364008	Copper (II) acetate hydrate	278
364507	Copper (II) chloride dihydrate	278
364508	Copper (II) chloride dihydrate	278
364611	Copper (I) iodide	277
364631	Copper (I) iodide	277
364637	Copper (I) iodide	277
364752	Copper (II) sulfate pentahydrate	281
364757	Copper (II) sulfate pentahydrate	281
364759	Copper (II) sulfate pentahydrate	281
365002	Copper (II) sulfate anhydrous	280
365006	Copper (II) sulfate anhydrous	280
365007	Copper (II) sulfate anhydrous	280
365152	D(+)-Sucrose	726
365157	D(+)-Sucrose	726
365158	D(+)-Sucrose	726
365755	Saponin	644
365757	Saponin	644
365758	Saponin	644
366205	Sodium acetate trihydrate	663
366207	Sodium acetate trihydrate	663
366209	Sodium acetate trihydrate	663
366371	Sodium acetate anhydrous	662
366372	Sodium acetate anhydrous	662
366377	Sodium acetate anhydrous	662
366551	Sodium alginate	663
366552	Sodium alginate	663
366553	Sodium alginate	663
366681	Sodium L-ascorbate	664
366684	Sodium L-ascorbate	664
366754	Sodium benzoate	665
366757	Sodium benzoate	665
366759	Sodium benzoate	665
366902	Sodium bicarbonate	666
366904	Sodium bicarbonate	666
366908	Sodium bicarbonate	666
366909	Sodium bicarbonate	666
367201	Sodium tetraborate decahydrate	710
367207	Sodium tetraborate decahydrate	710
367209	Sodium tetraborate decahydrate	710
367357	Sodium bromide	667
367359	Sodium bromide	667
367601	Sodium carbonate decahydrate	668
367608	Sodium carbonate decahydrate	668
367609	Sodium carbonate decahydrate	668
367691	Sodium carbonate monohydrate	669
367692	Sodium carbonate monohydrate	669
367693	Sodium carbonate monohydrate	669
367694	Sodium carbonate monohydrate	669
367703	Sodium carbonate anhydrous	668
367704	Sodium carbonate anhydrous	668
367705	Sodium carbonate anhydrous	668
367707	Sodium carbonate anhydrous	668
367951	Sodium citrate dibasic sesquihydrate	672
368051	Sodium citrate tribasic dihydrate	673
368052	Sodium citrate tribasic dihydrate	673
368054	Sodium citrate tribasic dihydrate	673
368057	Sodium citrate tribasic dihydrate	673
368058	Sodium citrate tribasic dihydrate	673
368102	Sodium citrate tribasic anhydrous	672
368107	Sodium citrate tribasic anhydrous	672
368253	Sodium chloride	671
368253000	Sodium chloride	671
368256000	Sodium chloride	671
368257	Sodium chloride	671
368257000	Sodium chloride	671
368259	Sodium chloride	671
368281	Sodium chloride	671
368351	Sodium hexametaphosphate	677
368352	Sodium hexametaphosphate	677
368357	Sodium hexametaphosphate	677
368451	Sodium fluoride	675
368457	Sodium fluoride	675
368458	Sodium fluoride	675
369132	Sodium phosphate monobasic dihydrate	703
369138	Sodium phosphate monobasic dihydrate	703

INDEX BY PRODUCT CODE

369139	Sodium phosphate monobasic dihydrate	703
369141	Sodium phosphate monobasic monohydrate	703
369142	Sodium phosphate monobasic monohydrate	703
369143	Sodium phosphate monobasic monohydrate	703
369152	Sodium phosphate dibasic dodecahydrate	702
369154	Sodium phosphate dibasic dodecahydrate	702
369158	Sodium phosphate dibasic dodecahydrate	702
369159	Sodium phosphate dibasic dodecahydrate	702
369185	Sodium phosphate dibasic dihydrate	701
369211	Sodium phosphate dibasic anhydrous	700
369212	Sodium phosphate dibasic anhydrous	700
369213	Sodium phosphate dibasic anhydrous	700
369252	Sodium phosphate dibasic anhydrous	701
369257	Sodium phosphate dibasic anhydrous	701
369258	Sodium phosphate dibasic anhydrous	701
369275	Sodium phosphate dibasic anhydrous	700
369301	Sodium phosphate tribasic dodecahydrate	704
369309	Sodium phosphate tribasic dodecahydrate	704
369447	Sodium glycerophosphate pentahydrate	676
369449	Sodium glycerophosphate pentahydrate	676
369581	Sodium gluconate	676
369582	Sodium gluconate	676
369663	Sodium glutamate acid	676
369667	Sodium glutamate acid	676
369701000	Sodium hydroxide solution 30%	682
369702	Sodium hydroxide solution 30%	682
369704	Sodium hydroxide solution 30%	682
369706	Sodium hydroxide solution 30%	682
369741	Sodium hydroxide, pearls	678
369742	Sodium hydroxide, pearls	678
369743	Sodium hydroxide, pearls	678
369744	Sodium hydroxide, pearls	678
369761	Sodium hydroxide solution 30%	682
369762	Sodium hydroxide solution 30%	682
369766	Sodium hydroxide solution 30%	682
369771	Sodium hydroxide, pellets	679
369772	Sodium hydroxide, pellets	679
369774	Sodium hydroxide, pellets	679
369777	Sodium hydroxide, pellets	679
369812	Sodium hydroxide 0.25 mol/l (0.25N)	689
370011	Sodium hydrosulfite	678
370014	Sodium hydrosulfite	678
370305	Sodium iodide	693
370307	Sodium iodide	693
370309	Sodium iodide	693
370321	Sodium hypochlorite solution in water	692
370322	Sodium hypochlorite solution in water	692
370323	Sodium hypochlorite solution in water	692
370751	Sodium metabisulfite	694
370752	Sodium metabisulfite	694
370753	Sodium metabisulfite	694
371802	Sodium nitrate	696
371804	Sodium nitrate	696
371809	Sodium nitrate	696
371901	Sodium nitrite	696
371902	Sodium nitrite	696
371903	Sodium nitrite	696
372109	Sodium nitrite	697
372201	Sodium oxalate	698
372203	Sodium oxalate	698
373603	Sodium salicylate	705
373607	Sodium salicylate	705
373608	Sodium salicylate	705
373902	Sodium silicate	705
373908	Sodium silicate	705
373909	Sodium silicate	705
375713	Sodium sulfate anhydrous	707
375716	Sodium sulfate anhydrous	707
376002	Sodium sulfite anhydrous	708
376003	Sodium sulfite anhydrous	708
376006	Sodium sulfite anhydrous	708
376008	Sodium sulfite anhydrous	708
376009	Sodium sulfite anhydrous	708
376403	Sodium sulfide trihydrate	708
377901	Sodium thiosulfate pentahydrate	711
377907	Sodium thiosulfate pentahydrate	711
377909	Sodium thiosulfate pentahydrate	711
378202	Sodium thiosulfate pentahydrate	711
378207	Sodium thiosulfate pentahydrate	711
378209	Sodium thiosulfate pentahydrate	711
378372	Sodium thiosulfate anhydrous	711
378377	Sodium thiosulfate anhydrous	711
378378	Sodium thiosulfate anhydrous	711
378802	Sulfur sublimed and washed	730
378807	Sulfur sublimed and washed	730
378809	Sulfur sublimed and washed	730
379013	Sorbitol	714
379021	Sorbitol (no crystallizable) solution 70%	715
379022	Sorbitol (no crystallizable) solution 70%	715
379403	Tin (II) chloride dihydrate	761
379406	Tin (II) chloride dihydrate	761
379407	Tin (II) chloride dihydrate	761
379601	Tin (II) sulfate	762
379707	Strontium nitrate	724
382105	Talc	744
382107	Talc	744
382109	Talc	744
382981	Tetrahydrofuran	753
382982	Tetrahydrofuran	753
382983	Tetrahydrofuran	753
382985	Tetrahydrofuran	753
382986	Tetrahydrofuran	753
384201	Thymol	757
384202	Thymol	757
384205	Thymol	757
385403	Thiourea	756
385407	Thiourea	756
385409	Thiourea	756
385702	Titanium dioxide	764
385707	Titanium dioxide	764
385709	Titanium dioxide	764
385751	Titanium dioxide	764
385752	Titanium dioxide	764
385753	Titanium dioxide	764
386001	Toluene	767
386002	Toluene	767
386003	Toluene	767
386009	Toluene	767
386102	Toluene	767
386104	Toluene	767
386106	Toluene	767
386202	Tricresol	771
386301	Triethanolamine	772
386303	Triethanolamine	772
386304	Triethanolamine	772
386305	Triethanolamine	772
386601	Triethylamine	772
386602	Triethylamine	772
386603	Triethylamine	772
387203	Trimethylcetylammonium p-toluenesulfonate	775
387503	Paraformaldehyde	556
387507	Paraformaldehyde	556
387801	Urea	782
387805	Urea	782
387807	Urea	782
387809	Urea	782
388102	Vanillin	784
388104	Vanillin	784
388107	Vanillin	784
388108	Vanillin	784
388407	Paraffin white soft	556
388409	Paraffin white soft	556
388607	Paraffin white soft	556
388609	Paraffin white soft	556
388701	Gentian violet	368
388703	Gentian violet	368
389511	Riboflavine	639
391804	p-Aminobenzoic acid	158
391805	p-Aminobenzoic acid	158
392304	Nicotinamide	525
392307	Nicotinamide	525
392602	Xylene, mix of isomers	791
392603	Xylene, mix of isomers	791
392605	Xylene, mix of isomers	791
392608	Xylene, mix of isomers	791
392631	D(+)Xylose	793
392635	D(+)Xylose	793
393007	Zinc chloride anhydrous	799
393009	Zinc chloride anhydrous	799
393503	Zinc oxide	801
393507	Zinc oxide	801
393509	Zinc oxide	801
394001	Zinc sulfate heptahydrate	802
394007	Zinc sulfate heptahydrate	802
394009	Zinc sulfate heptahydrate	802
395451	Zinc stearate	801
395452	Zinc stearate	801
396105	Zirconyl nitrate	805
400155	Acetal	129
400204	Acetamide	129
400255	Acetanilide	129
400305	Acetylacetone	146
400307	Acetylacetone	146
400522	n-Acetyl-L-cysteine	146
400932	Acetone	138
400952	Acetone	138
400961	Acetone	140
400962	Acetone	140
400963	Acetone	140
400965	Acetone	140
400971	Acetone	140
400972	Acetone	140
400974	Acetone	140
400978	Acetone	140
400979	Acetone	140
400991	Acetone	138
400992000	Acetone	138
400994	Acetone	138

INDEX BY PRODUCT CODE

401032	Acetone	139
401034	Acetone	139
401041	Acetone	140
401042	Acetone	140
401051	Acetone	139
401052	Acetone	139
401054	Acetone	139
401055	Acetone	139
401058	Acetone	139
401183000	Acetonitrile	144
401185000	Acetonitrile	144
401212	Acetonitrile	143
401216	Acetonitrile	143
401241	Acetonitrile	143
401242	Acetonitrile	143
401243	Acetonitrile	143
401361	Acetic acid glacial	131
401391	Acetic acid glacial	132
401392	Acetic acid glacial	132
401396	Acetic acid glacial	132
401397	Acetic acid glacial	132
401405	Acetic acid glacial	131
401406	Acetic acid glacial	131
401407	Acetic acid glacial	131
401411	Acetic acid glacial	130
401412	Acetic acid glacial	130
401413	Acetic acid glacial	130
401414	Acetic acid glacial	130
401421	Acetic acid glacial	132
401422	Acetic acid glacial	132
401424	Acetic acid glacial	132
401425	Acetic acid glacial	132
401431	Acetic acid glacial	130
401432	Acetic acid glacial	130
401453	Acetic acid glacial	132
401455	Acetic acid glacial	132
401462	Acetic acid glacial	131
401463	Acetic acid glacial	131
401531	Acetic acid 12%	135
401561	Acetic acid 0.1 mol/l (0.1N)	136
401785	Adipic acid	147
402032	1-Amino-2-naphthol-4-sulfonic acid	158
402133	Anisic acid	181
402381	Sodium arsenite 0.05 mol/l (0.1N)	664
402404	L(+)Ascorbic acid	186
402406	L(+)Ascorbic acid	186
402407	L(+)Ascorbic acid	186
402442	L(+)Aspartic acid	186
402532	Barbituric acid	191
402535	Barbituric acid	191
402621	Benzoic acid	199
402635	Benzoic acid	199
402636	Benzoic acid	199
402637	Benzoic acid	199
402762	Boric acid	205
402763	Boric acid	205
402764	Boric acid	205
402765	Boric acid	205
402766	Boric acid	205
402767	Boric acid	205
402922	Hydrobromic acid 48%	393
402925	Hydrobromic acid 48%	393
403236	n-Butyric acid	231
403308	Calconcarbonic acid	245
403421	n-Caprylic acid	247
403473	n-Caproic acid	247
403721	Citric acid monohydrate	269
403722	Citric acid monohydrate	269
403724	Citric acid monohydrate	269
403725	Citric acid monohydrate	269
403727	Citric acid monohydrate	269
403821	Chloranilic acid	257
403822	Chloranilic acid	257
403871	Hydrochloric acid 37%	395
403872	Hydrochloric acid 37%	395
403874	Hydrochloric acid 37%	395
403876	Hydrochloric acid 37%	395
403878	Hydrochloric acid 37%	395
403891	Hydrochloric acid 32-35%	396
403901	Hydrochloric acid 23%	399
403905	Hydrochloric acid 23%	399
403909	Hydrochloric acid 23%	399
403915	Hydrochloric acid 34-37%	396
403916	Hydrochloric acid 34-37%	396
403917	Hydrochloric acid 34-37%	396
403921	Hydrochloric acid 29-31%	398
403941	Hydrochloric acid 37%	394
403942	Hydrochloric acid 37%	394
403971	Hydrochloric acid 37%	394
403974	Hydrochloric acid 37%	394
403977	Hydrochloric acid 37%	394
403981	Hydrochloric acid 32%	397
403982	Hydrochloric acid 32%	397
403984	Hydrochloric acid 32%	397
403986	Hydrochloric acid 32%	397
403988	Hydrochloric acid 32%	397
404033	Hydrochloric acid 8%	400
404036	Hydrochloric acid 8%	400
404061000	Hydrochloric acid 2 mol/l (2N)	403
404062000	Hydrochloric acid 2 mol/l (2N)	403
404067000	Hydrochloric acid 2 mol/l (2N)	403
404091000	Hydrochloric acid 1 mol/l (1N)	404
404092000	Hydrochloric acid 1 mol/l (1N)	404
404094000	Hydrochloric acid 1 mol/l (1N)	404
404097000	Hydrochloric acid 1 mol/l (1N)	404
404111	Hydrochloric acid 1 mol/l (1N)	404
404141000	Hydrochloric acid 0.5 mol/l (0.5N)	405
404142000	Hydrochloric acid 0.5 mol/l (0.5N)	405
404147000	Hydrochloric acid 0.5 mol/l (0.5N)	405
404161	Hydrochloric acid 0.5 mol/l (0.5N)	405
404191000	Hydrochloric acid 0.1 mol/l (0.1N)	406
404192000	Hydrochloric acid 0.1 mol/l (0.1N)	406
404195000	Hydrochloric acid 0.1 mol/l (0.1N)	406
404197000	Hydrochloric acid 0.1 mol/l (0.1N)	406
404211	Hydrochloric acid 0.1 mol/l (0.1N)	406
404251	Hydrochloric acid 0.01 mol/l (0.01N)	408
404267	Hydrochloric acid 0.01 mol/l (0.01N)	408
404308	Chloroacetic acid	258
404872	Chromotropic acid disodium salt	268
405101	Dichloroacetic acid	297
405103	Dichloroacetic acid	297
405144	Sodium diethyldithiocarbamate trihydrate	674
405192	Diethylenetriaminepentacetic acid	308
405351	Dodecylbenzenesulphonic acid sodium salt	325
405352	Dodecylbenzenesulphonic acid sodium salt	325
405421	Ethylenediaminetetraacetic acid disodium salt 0.1 mol/l (0.2N)	346
405424	Ethylenediaminetetraacetic acid tripotassic salt	348
405431	Ethylenediaminetetraacetic acid disodium salt 0.01 mol/l (0.02N)	347
405442000	Ethylenediaminetetraacetic acid disodium salt 0.01 mol/l (0.02N)	347
405443000	Ethylenediaminetetraacetic acid disodium salt 0.01 mol/l (0.02N)	347
405451	Heptafluorobutyric acid	380
405463	Ethylenediaminetetraacetic acid	344
405465	Ethylenediaminetetraacetic acid	344
405482	Ethylenediaminetetraacetic acid tetrasodium salt tetrahydrate	348
405486	Ethylenediaminetetraacetic acid tetrasodium salt tetrahydrate	348
405491	Ethylenediaminetetraacetic acid disodium salt	345
405492	Ethylenediaminetetraacetic acid disodium salt	345
405494	Ethylenediaminetetraacetic acid disodium salt	345
405497	Ethylenediaminetetraacetic acid disodium salt	345
405501000	Ethylenediaminetetraacetic acid disodium salt 0.05 mol/l (0.1N)	347
405502000	Ethylenediaminetetraacetic acid disodium salt 0.05 mol/l (0.1N)	347
405511000	Ethylenediaminetetraacetic acid disodium salt 0.1 mol/l (0.2N)	346
405512000	Ethylenediaminetetraacetic acid disodium salt 0.1 mol/l (0.2N)	346
405513000	Ethylenediaminetetraacetic acid disodium salt 0.1 mol/l (0.2N)	346
405514000	Ethylenediaminetetraacetic acid disodium salt 0.1 mol/l (0.2N)	346
405521	Ethylene glycol bis(2-aminoethyl ether)-n,n,n',n'-tetraacetic acid	349
405522	Ethylene glycol bis(2-aminoethyl ether)-n,n,n',n'-tetraacetic acid	349
405531	Ethylenediaminetetraacetic acid dipotassium salt dihydrate	345
405541	Ethylenediaminetetraacetic acid potassium and magnesium salt dihydrate	347
405582	Ethylenediaminetetraacetic acid dipotassium salt dihydrate	345
405611	Hydrofluoric acid 47-51%	411
405621	1-Hexanesulphonic acid sodium salt	389
405622	1-Hexanesulphonic acid sodium salt	389
405631	1-Butanesulfonic acid sodium salt	225
405632	1-Butanesulfonic acid sodium salt	225
405653	Hydrofluoric acid 50%	410
405716	Hydrofluoric acid 47-51%	411
405722	Hydrofluoric acid 50%	410
405737	Hydrofluoric acid 50%	410
405761	Hydrofluoric acid 39.5%	412

INDEX BY PRODUCT CODE

405765	Hydrofluoric acid 39.5%	412
405775	Hydrofluoric acid diluted	412
405792	Formic acid 99%	363
405793	Formic acid 99%	363
405794	Formic acid 99%	363
405821	Formic acid 99%	362
405822	Formic acid 99%	362
405823	Formic acid 99%	362
405824	Formic acid 99%	362
405832	Formic acid 85%	363
405833	Formic acid 85%	363
405835	Formic acid 85%	363
405841	1-Pentanesulphonic acid sodium salt	559
405842	1-Pentanesulphonic acid sodium salt	559
405851	1-Heptanesulphonic acid sodium salt	383
405852	1-Heptanesulphonic acid sodium salt	383
405861	1-Octanesulphonic acid sodium salt	537
405862	1-Octanesulphonic acid sodium salt	537
405863	1-Octanesulphonic acid sodium salt	537
405871	1-Decanesulfonic acid sodium salt	290
405872	1-Decanesulfonic acid sodium salt	290
405881	1-Dodecanesulfonic acid sodium salt	324
405882	1-Dodecanesulfonic acid sodium salt	324
405891	1-Pentanesulphonic acid sodium salt monohydrate	560
405892	1-Pentanesulphonic acid sodium salt monohydrate	560
405901	1-Propanesulfonic acid sodium salt	627
405902	1-Propanesulfonic acid sodium salt	627
405913	Phosphomolybdic acid	576
405915	Phosphomolybdic acid	576
405921	1-Hexanesulphonic acid sodium salt monohydrate	389
405922	1-Hexanesulphonic acid sodium salt monohydrate	389
405931	1-Octanesulfonic acid sodium salt monohydrate	538
405932	1-Octanesulfonic acid sodium salt monohydrate	538
405941	Dodecyltrimethylammonium bromide	325
405961	Orthophosphoric acid 99%	547
405967	Orthophosphoric acid 99%	547
405971	Tetrabutylammonium bisulfate	748
405972	Tetrabutylammonium bisulfate	748
406002	Orthophosphoric acid 85%	548
406003	Orthophosphoric acid 85%	548
406005	Orthophosphoric acid 85%	548
406021	Orthophosphoric acid 85%	548
406022	Orthophosphoric acid 85%	548
406053	Phosphonic acid	577
406056	Phosphonic acid	577
406154	Phosphotungstic acid	578
406205	Phthalic acid	579
406284	Fumaric acid	365
406287	Fumaric acid	365
406335	Gallic acid monohydrate	367
406336	Gallic acid monohydrate	367
406434	Glycolic acid	374
406485	L(+)Glutamic acid	371
406831	Hydriodic acid 57%	393
406918	Ferron	354
406961	Hypophosphorous acid 50%	418
406962	Hypophosphorous acid 50%	418
407261	Maleic acid	477
407263	Maleic acid	477
407266	Maleic acid	477
407314	DL-Malic acid	477
407316	DL-Malic acid	477
407363	Malonic acid	478
407421	D(-)Mandelic acid	478
407441	alpha-Methoxyphenylacetic acid	495
407465	Metaphosphoric acid	487
407467	Metaphosphoric acid	487
407481	Methanesulfonic acid	487
407483	Methanesulfonic acid	487
407914	Nicotinic acid	525
407951	Nitric acid 65%	531
407952	Nitric acid 65%	531
408021	Nitric acid 65%	531
408022	Nitric acid 65%	531
408025	Nitric acid 65%	531
408027	Nitric acid 65%	531
408051	Nitric acid 67-69%	530
408052	Nitric acid 67-69%	530
408071	Nitric acid 69.5%	529
408072	Nitric acid 69.5%	529
408075	Nitric acid 69.5%	529
408076	Nitric acid 69.5%	529
408097	Nitric acid 69.5%	528
408098	Nitric acid 69.5%	528
408101	Nitric acid 65%	531
408102	Nitric acid 65%	531
408115	Nitric acid 67-70%	529
408116	Nitric acid 67-70%	529
408117	Nitric acid 67-70%	529
408142	Nitric acid fuming 99%	528
408143	Nitric acid fuming 99%	528
408151	Nitric acid 69.5%	528
408152	Nitric acid 69.5%	528
408171000	Nitric acid 1 mol/l (1N)	533
408176000	Nitric acid 1 mol/l (1N)	533
408185000	Nitric acid 2 mol/l (2N)	533
408191	Nitric acid 18%	532
408206000	Nitric acid 0.1 mol/l (0.1N)	533
408231	Nitric acid 0.1 mol/l (0.1N)	533
408242	Nitrilotriacetic acid	534
408414	m-Nitrobenzoic acid	535
408731	Oxalic acid dihydrate	550
408733	Oxalic acid dihydrate	550
408736	Oxalic acid dihydrate	550
408737	Oxalic acid dihydrate	550
408826	Oxalic acid 0.5 mol/l (1N)	551
408856	Oxalic acid 0.05 mol/l (0.1N)	551
408871	Oxalic acid 0.05 mol/l (0.1N)	551
408901	Oxalic acid 0.005 mol/l (0.01N)	551
409064	1-Pentanesulphonic acid sodium salt	559
409111	Perchloric acid 65%	561
409113	Perchloric acid 65%	561
409121	Perchloric acid 65%	561
409131	Perchloric acid 0.1 mol/l (0.1N) in acetic acid	562
409136	Perchloric acid 0.1 mol/l (0.1N) in acetic acid	562
409182	Periodic acid	563
409184	Periodic acid	563
409185	Periodic acid	563
409193	Perchloric acid 65-71%	560
409302	Picric acid solution	580
409305	Picric acid solution	580
409435	Pyrogallol	632
409471	Pyrrolidine dithiocarbamic acid ammonium salt	633
409551	Propionic acid	628
409553	Propionic acid	628
409702	Rosolic acid	640
409773	Salicylic acid	642
409777	Salicylic acid	642
409875	Sebacic acid	646
409964	Selenous acid	648
410051	Silicotungstic acid	652
410104	Sulfamic acid	727
410105	Sulfamic acid	727
410106	Sulfamic acid	727
410154	Sulfanilic acid	728
410156	Sulfanilic acid	728
410261	Sulfuric acid 96%	732
410301	Sulfuric acid 96%	732
410302	Sulfuric acid 96%	732
410303	Sulfuric acid 96%	732
410304	Sulfuric acid 96%	732
410306	Sulfuric acid 96%	732
410307	Sulfuric acid 96%	732
410351	Sulfuric acid 93-98%	733
410371	Sulfuric acid 96%	731
410374	Sulfuric acid 96%	731
410381	Sulfuric acid 96%	732
410382	Sulfuric acid 96%	732
410391	Sulfuric acid 90%	734
410394	Sulfuric acid 90%	734
410405	Sulfuric acid 93-98%	734
410406	Sulfuric acid 93-98%	734
410407	Sulfuric acid 93-98%	734
410421	Sulfuric acid 98%	730
410511000	Sulfuric acid 20%	736
410516	Sulfuric acid 20%	736
410547000	Sulfuric acid 1 mol/l (2N)	738
410548000	Sulfuric acid 1 mol/l (2N)	738
410571000	Sulfuric acid 0.5 mol/l (1N)	738
410572000	Sulfuric acid 0.5 mol/l (1N)	738
410575000	Sulfuric acid 0.5 mol/l (1N)	738
410577000	Sulfuric acid 0.5 mol/l (1N)	738
410591	Sulfuric acid 0.5 mol/l (1N)	738
410634	Sulfuric acid 0.33 mol/l (2N/3)	738
410662000	Sulfuric acid 0.25 mol/l (0.5N)	739
410663000	Sulfuric acid 0.25 mol/l (0.5N)	739
410667000	Sulfuric acid 0.25 mol/l (0.5N)	739
410681	Sulfuric acid 0.25 mol/l (0.5N)	739
410711000	Sulfuric acid 0.05 mol/l (0.1N)	741
410712000	Sulfuric acid 0.05 mol/l (0.1N)	741
410714000	Sulfuric acid 0.05 mol/l (0.1N)	741
410715000	Sulfuric acid 0.05 mol/l (0.1N)	741
410717000	Sulfuric acid 0.05 mol/l (0.1N)	741
410731	Sulfuric acid 0.05 mol/l (0.1N)	741
410791	Sulfuric acid 0.005 mol/l (0.01N)	742
410894	Sulfosalicylic acid	729
410896	Sulfosalicylic acid	729
411023	Succinic acid	725
411025	Succinic acid	725
411027	Succinic acid	725
411074	Tannic acid	744
411076	Tannic acid	744
411121	L(+) Tartaric Acid	745
411125	L(+) Tartaric Acid	745
411127	L(+) Tartaric Acid	745
411271	2-Thiobarbituric acid	755
411272	2-Thiobarbituric acid	755
411385	Thioglycolic acid 80%	755
411432	p-Toluenesulfonic acid	768
411436	p-Toluenesulfonic acid	768
411504	p-Toluenesulfonic acid sodium salt	769
411524	Trichloroacetic acid	770
411525	Trichloroacetic acid	770

Product specifications are subject to changes. Please visit our website for updates.

INDEX BY PRODUCT CODE

411527	Trichloroacetic acid	770
411541	Trifluoroacetic acid	773
411542	Trifluoroacetic acid	773
411543	Trifluoroacetic acid	773
411554000	Trichloroacetic acid solution 20%	770
411561	Trifluoroacetic acid	773
411564	Trifluoroacetic acid	773
411628	Tungstic acid	780
411981	Water chlorine	789
412011	Water	786
412031	Water + 0.1% v/v trifluoroacetic acid	789
412041	Acetonitrile	142
412042	Acetonitrile	142
412051	Hydrogen peroxide solution 30-32%	413
412071	Hydrogen peroxide solution 30%	415
412072	Hydrogen peroxide solution 30%	415
412074	Hydrogen peroxide solution 30%	415
412076	Hydrogen peroxide solution 30%	415
412077	Hydrogen peroxide solution 30%	415
412081	Hydrogen peroxide solution 30%	414
412091	Water	786
412092	Water	786
412102	Hydrogen peroxide solution 30%	414
412111	Water	786
412112	Water	786
412121	Water + 0.1% v/v formic acid	789
412122	Water + 0.1% v/v formic acid	789
412141	Water	786
412142	Water	786
412151	Water	787
412161	Hydrogen peroxide solution 30%	414
412162	Hydrogen peroxide solution 30%	414
412163	Hydrogen peroxide solution 30%	414
412185	Water	787
412321	Acetonitrile + 0.1% v/v trifluoroacetic acid	145
412322	Acetonitrile + 0.1% v/v trifluoroacetic acid	145
412331	Acetonitrile + 0.1% v/v formic acid	145
412332	Acetonitrile + 0.1% v/v formic acid	145
412341	Acetonitrile	142
412342	Acetonitrile	142
412371000	Acetonitrile	142
412372000	Acetonitrile	142
412374	Acetonitrile	142
412381	Methanol	488
412383	Methanol	488
412391000	Acetonitrile	142
412392000	Acetonitrile	142
412393	Acetonitrile	142
412395	Acetonitrile	142
412407	Acetonitrile	143
412409	Acetonitrile	143
412411000	Acetonitrile	143
412412000	Acetonitrile	143
412413000	Acetonitrile	143
412421000	Propan-2-ol	624
412422000	Propan-2-ol	624
412431000	Cyclohexane	286
412432000	Cyclohexane	286
412441000	Isooctane	439
412442000	Isooctane	439
412451000	Tetrahydrofuran	751
412452000	Tetrahydrofuran	751
412453000	Tetrahydrofuran	751
412471	Tetrahydrofuran	751
412472	Tetrahydrofuran	751
412501	Acetone	138
412502	Acetone	138
412511000	Butanol-1	226
412512000	Butanol-1	226
412521	Ethanol absolute anhydrous	331
4125212	Ethanol absolute anhydrous	331
412522	Ethanol absolute anhydrous	331
4125222	Ethanol absolute anhydrous	331
412531	Methanol	489
412532	Methanol	489
412533	Methanol	489
412535	Methanol	489
412541000	Propan-1-ol	622
412542000	Propan-1-ol	622
412571	Chloroform	259
412572	Chloroform	259
412591000	n-Heptane 99%	380
412592000	n-Heptane 99%	380
412601000	n-Hexane	386
412602000	n-Hexane	386
412611000	Ethyl acetate	340
412612000	Ethyl acetate	340
412621000	Dichloromethane	300
412622000	Dichloromethane	300
412631	Hexane mixture of isomers	388
412632	Hexane mixture of isomers	388
412641000	Toluene	765
412642000	Toluene	765
412652	Chloroform	259
412653	Chloroform	259
412661	Dichloromethane	300
412662	Dichloromethane	300
412671	Diethyl ether	308
412672	Diethyl ether	308
412674	Diethyl ether	308
412681	2-Methyltetrahydrofuran	505
412682	2-Methyltetrahydrofuran	505
412691	n-Hexane 99%	385
412692	n-Hexane 99%	385
412701	Ethanol absolute anhydrous	331
4127012	Ethanol absolute anhydrous	331
412702	Ethanol absolute anhydrous	331
4127022	Ethanol absolute anhydrous	331
412703	Ethanol absolute anhydrous	331
4127032	Ethanol absolute anhydrous	331
412704	Ethanol absolute anhydrous	331
412711000	Propan-2-ol	623
412712000	Propan-2-ol	623
412721	Methanol	488
412722	Methanol	488
412724	Methanol	488
412725	Methanol	488
412821	Propan-2-ol	623
413603	ß-Alanine	148
413654	Albumin from eggs, dried	148
413656	Albumin from eggs, dried	148
413671	Albumin from eggs powder	148
413672	Albumin from eggs powder	148
413783	n-Amyl alcohol	180
413801	Isoamyl alcohol	436
413832	Isoamyl alcohol	437
413833	Isoamyl alcohol	437
413836	Isoamyl alcohol	437
413892	Isoamyl alcohol	436
413941	tert-Amyl alcohol	180
413944	tert-Amyl alcohol	180
413945	tert-Amyl alcohol	180
414022	Benzyl alcohol	200
414024	Benzyl alcohol	200
414052	Benzyl alcohol	200
414131	Butanol-1	226
414132	Butanol-1	226
414133	Butanol-1	226
414211	Isobutanol	437
414213	Isobutanol	437
414214	Isobutanol	437
414261	Butanol-2	227
414264	Butanol-2	227
414266	Butanol-2	227
414341	tert-Butanol	227
414342	tert-Butanol	227
414343	tert-Butanol	227
414346	tert-Butanol	227
414427	Cetyl alcohol	254
414511	Mixture Ethanol 99° / Isopropanol	508
414541	Ethanol 96°	334
4145412	Ethanol 96°	334
414542	Ethanol 96°	334
4145422	Ethanol 96°	334
414551	Mixture Ethanol 95° / Isopropanol	508
414583	Ethanol absolute anhydrous	332
4145832	Ethanol absolute anhydrous	332
414587	Ethanol absolute anhydrous	332
4145872	Ethanol absolute anhydrous	332
414601	Ethanol absolute anhydrous	332
4146012	Ethanol absolute anhydrous	332
414602	Ethanol absolute anhydrous	332
4146022	Ethanol absolute anhydrous	332
414603	Ethanol absolute anhydrous	332
4146032	Ethanol absolute anhydrous	332
414604	Ethanol absolute anhydrous	332
414605	Ethanol absolute anhydrous	332
4146052	Ethanol absolute anhydrous	332
414606	Ethanol absolute anhydrous	332
4146062	Ethanol absolute anhydrous	332
414607	Ethanol absolute anhydrous	332
4146072	Ethanol absolute anhydrous	332
414608	Ethanol absolute anhydrous	332
4146082	Ethanol absolute anhydrous	332
414609	Ethanol absolute anhydrous	332
4146092	Ethanol absolute anhydrous	332
414622	Ethanol 96°	335
414631	Ethanol 96°	335
4146312	Ethanol 96°	335
414632	Ethanol 96°	335
4146322	Ethanol 96°	335
414633	Ethanol 96°	335
414634	Ethanol 96°	335
4146342	Ethanol 96°	335
414635	Ethanol 96°	335
4146352	Ethanol 96°	335
414637	Ethanol 96°	335
4146372	Ethanol 96°	335
414638	Ethanol 96°	335
414639	Ethanol 96°	335
414667	Ethanol 96°	334
4146672	Ethanol 96°	334
414677	Ethanol absolute anhydrous	331
4146772	Ethanol absolute anhydrous	331
414813	Methanol	491
414814	Methanol	491
414815	Methanol	491
414816	Methanol	491
414817	Methanol	491
414818	Methanol	491
414819	Methanol	491
414821	Methanol	491
414822	Methanol	491

INDEX BY PRODUCT CODE

414831	Methanol	488
414832	Methanol	488
414854	Methanol	491
414855	Methanol	491
414861	Methanol + 0.1% v/v formic acid	492
414862	Methanol + 0.1% v/v formic acid	492
414871	Methanol + 0.1% v/v trifluoroacetic acid	493
414872	Methanol + 0.1% v/v trifluoroacetic acid	493
414881	Methanol	491
414883	Methanol	491
414902	Methanol	490
414903	Methanol	490
414914	Methanol	490
414917	Methanol	490
414930	Methanol	489
414932	Methanol	489
414941	Methanol	488
414942	Methanol	488
414952	Methanol	489
414981	Methanol	490
415002	Octanol-1	538
415003	Octanol-1	538
415004	Octanol-1	538
415102	Propan-1-ol	622
415104	Propan-1-ol	622
415106	Propan-1-ol	622
415108	Propan-1-ol	622
415112	Propan-2-ol	624
415152	Propan-2-ol	625
415153	Propan-2-ol	625
415154	Propan-2-ol	625
415156	Propan-2-ol	625
415157	Propan-2-ol	625
415158	Propan-2-ol	625
415161	Propan-2-ol	625
415162	Propan-2-ol	625
415173	Propan-2-ol	625
415183	Propan-2-ol	623
415184	Propan-2-ol	623
415213	Propan-2-ol	624
415231	Propan-2-ol	625
415233	Propan-2-ol	625
415235	Propan-2-ol	625
415236	Propan-2-ol	625
415237	Propan-2-ol	625
415238	Propan-2-ol	625
415271	2,2,2 - Trichlorethanol	770
415281	Propan-2-ol	624
415312	Anisaldehyde	181
415362	Benzaldehyde	197
415631	Formaldehyde 4% w/v buffered at pH 6.9	361
415633	Formaldehyde 4% w/v buffered at pH 6.9	361
415634	Formaldehyde 4% w/v buffered at pH 6.9	361
415636	Formaldehyde 4% w/v buffered at pH 6.9	361
415661	Formaldehyde 37% w/v	358
415662	Formaldehyde 37% w/v	358
415666	Formaldehyde 37% w/v	358
415667	Formaldehyde 37% w/v	358
415669	Formaldehyde 37% w/v	358
415672	Formaldehyde 5% w/v buffered at pH 6.9	360
415674	Formaldehyde 5% w/v buffered at pH 6.9	360
415682	Formaldehyde 37% w/v neutralized	359
415683	Formaldehyde 37% w/v neutralized	359
415684	Formaldehyde 37% w/v neutralized	359
415685	Formaldehyde 37% w/v neutralized	359
415686	Formaldehyde 37% w/v neutralized	359
415691	Formaldehyde 4% w/v buffered at pH 6.9	361
415692	Formaldehyde 4% w/v buffered at pH 6.9	361
415693	Formaldehyde 4% w/v buffered at pH 6.9	361
415694	Formaldehyde 4% w/v buffered at pH 6.9	361
415695	Formaldehyde 4% w/v buffered at pH 6.9	361
415696	Formaldehyde 4% w/v buffered at pH 6.9	361
415721	Glutardialdehyde solution 50%	372
415892	Alizarin	149
416002	Alizarin red	149
416281	Allylthiourea	150
416283	Allylthiourea	150
416531	Alumina white	151
416581	Aluminum standard solution	152
416815	Aluminum, powder	151
416817	Aluminum, powder	151
416892	Aluminum ammonium sulfate dodecahydrate	152
416895	Aluminum ammonium sulfate dodecahydrate	152
416897	Aluminum ammonium sulfate dodecahydrate	152
416942	Aluminum chloride hexahydrate	153
416943	Aluminum chloride hexahydrate	153
416945	Aluminum chloride hexahydrate	153
416947	Aluminum chloride hexahydrate	153
416949	Aluminum chloride hexahydrate	153
416996	Aluminum chloride anhydrous	153
417046	Aluminum hydroxide	154
417047	Aluminum hydroxide	154
417095	Aluminum nitrate nonahydrate	154
417096	Aluminum nitrate nonahydrate	154
417097	Aluminum nitrate nonahydrate	154
417144	Aluminum oxide	155
417145	Aluminum oxide	155
417146	Aluminum oxide	155
417147	Aluminum oxide	155
417182	Aluminum oxide (acid)	155
417185	Aluminum oxide (acid)	155
417214	Aluminum oxide (basic)	155
417217	Aluminum oxide (basic)	155
417241	Aluminum oxide (neutral)	156
417245	Aluminum oxide (neutral)	156
417295	Aluminum potassium sulfate dodecahydrate	156
417296	Aluminum potassium sulfate dodecahydrate	156
417297	Aluminum potassium sulfate dodecahydrate	156
417424	Aluminum sulfate	157
417425	Aluminum sulfate	157
417427	Aluminum sulfate	157
417585	Starch soluble	721
417587	Starch soluble	721
417781	Isoamyl acetate	436
417782	Isoamyl acetate	436
418381	4-Aminophenazone	159
418564	m-Aminophenol	159
418594	p-Aminophenol	159
418771	Ammonium acetate	166
418772	Ammonium acetate	166
418773	Ammonium acetate	166
418775	Ammonium acetate	166
418776	Ammonium acetate	166
418777	Ammonium acetate	166
418781	Ammonium acetate	166
418925	Ammonium bicarbonate	167
418926	Ammonium bicarbonate	167
418927	Ammonium bicarbonate	167
418929	Ammonium bicarbonate	167
419174	Ammonium bromide	167
419175	Ammonium bromide	167
419176	Ammonium bromide	167
419177	Ammonium bromide	167
419201	Ammonium carbamate	168
419202	Ammonium carbamate	168
419204	Ammonium carbamate	168
419232	Ammonium carbonate	168
419235	Ammonium carbonate	168
419236	Ammonium carbonate	168
419237	Ammonium carbonate	168
419239	Ammonium carbonate	168
419312	Ammonium citrate dibasic	170
419313	Ammonium citrate dibasic	170
419315	Ammonium citrate dibasic	170
419316	Ammonium citrate dibasic	170
419317	Ammonium citrate dibasic	170
419412	Ammonium chloride	169
419415	Ammonium chloride	169
419416	Ammonium chloride	169
419417	Ammonium chloride	169
419419	Ammonium chloride	169
419634	Ammonium fluoride	171
419635	Ammonium fluoride	171
419637	Ammonium fluoride	171
419638	Ammonium fluoride	171
419733	Ammonium formate	172
419734	Ammonium formate	172
419735	Ammonium formate	172
419736	Ammonium formate	172
419737	Ammonium formate	172
419741	Ammonium formate	171
419785	Ammonium phosphate monobasic	176
419786	Ammonium phosphate monobasic	176
419787	Ammonium phosphate monobasic	176
419831	Ammonium phosphate dibasic	176
419832	Ammonium phosphate dibasic	176
419834	Ammonium phosphate dibasic	176
419835	Ammonium phosphate dibasic	176
419836	Ammonium phosphate dibasic	176
419837	Ammonium phosphate dibasic	176
419941	Ammonia solution 30%	161
419943	Ammonia solution 30%	161
419945	Ammonia solution 30%	161
419946	Ammonia solution 30%	161
419948	Ammonia solution 30%	161
419981	Ammonia solution 20 - 22%	163
419983	Ammonia solution 20 - 22%	163
419984	Ammonia solution 20 - 22%	163
419993	Ammonia solution 25%	162
420031	Ammonia solution 25%	162
420051	Ammonia solution 25%	162
420052	Ammonia solution 25%	162
420071	Ammonia solution 30%	160
420073	Ammonia solution 30%	160
420075	Ammonia solution 30%	160
420077	Ammonia solution 30%	160
420084	Ammonia solution 25%	162
420085	Ammonia solution 25%	162
420133	Ammonium iodide	173
420135	Ammonium iodide	173
420161	Ammonia solution 20 - 22%	163

INDEX BY PRODUCT CODE

420175	Ammonia solution 20 - 22%	163
420234	Ammonium molybdate tetrahydrate	173
420236	Ammonium molybdate tetrahydrate	173
420238	Ammonium molybdate tetrahydrate	173
420391	Ammonium molybdate tetrahydrate	173
420422	Ammonium nitrate	174
420424	Ammonium nitrate	174
420425	Ammonium nitrate	174
420426	Ammonium nitrate	174
420427	Ammonium nitrate	174
420429	Ammonium nitrate	174
420473	Ammonium oxalate monohydrate	175
420475	Ammonium oxalate monohydrate	175
420476	Ammonium oxalate monohydrate	175
420477	Ammonium oxalate monohydrate	175
420478	Ammonium oxalate monohydrate	175
420623	Ammonium persulfate	175
420625	Ammonium persulfate	175
420626	Ammonium persulfate	175
420627	Ammonium persulfate	175
420629	Ammonium persulfate	175
420672	Reinecke salt	637
420724	Ammonium sulfamate	177
420725	Ammonium sulfamate	177
420771	Ammonium sulfate	177
420772	Ammonium sulfate	177
420774	Ammonium sulfate	177
420775	Ammonium sulfate	177
420776	Ammonium sulfate	177
420777	Ammonium sulfate	177
420885	Ammonium thiocyanate	178
420886	Ammonium thiocyanate	178
420887	Ammonium thiocyanate	178
420946	Ammonium thiocyanate 1 mol/l (1N)	179
420977	Ammonium thiocyanate 0.1 mol/l (0.1N)	179
421001	Ammonium thiocyanate 0.1 mol/l (0.1N)	179
421061	Ammonium thiocyanate 0.01 mol/l (0.01N)	179
421101	Ammonium sulfide solution 20%	178
421206	Ammonium L(+)tartrate	178
421491	Acetic anhydride	137
421493	Acetic anhydride	137
421496	Acetic anhydride	137
421735	Chromium (VI) oxide	267
421802	Phosphorus pentoxide	578
421808	Phosphorus pentoxide	578
421955	Maleic anhydride	477
422004	Molybdenum (VI) oxide	511
422005	Molybdenum (VI) oxide	511
422104	Silicon dioxide	652
422106	Silicon dioxide	652
422204	Succinic anhydride	725
422225	Trifluoroacetic anhydride	774
422241	Tungsten (VI) oxide	780
422251	Magnesium perchlorate	474
422252	Magnesium perchlorate	474
422254	Magnesium perchlorate	474
422376	Aniline hydrochloride	181
422731	Antimony standard solution	183
422834	Antimony trichloride	183
422835	Antimony trichloride	183
423035	Antimony potassium tartrate	183
423036	Antimony potassium tartrate	183
423037	Antimony potassium tartrate	183
423281	Anthrone	182
423282	Anthrone	182
423341	Orange II	539
423432	Orange G	539
423461	Acridine orange	147
423501	Methyl orange	501
423502	Methyl orange	501
423503	Methyl orange	501
423504	Methyl orange	501
423505	Methyl orange	501
423597	Xylenol orange	793
423598	Xylenol orange	793
423599	Xylenol orange	793
423611	Silver standard solution	654
423752	Silver, sheet	652
423791	Silver, wool	653
423913	Silver diethyldithiocarbamate	655
423952	Silver nitrate	655
423954	Silver nitrate	655
423955	Silver nitrate	655
423957	Silver nitrate	655
424035000	Silver nitrate 1 mol/l (1N)	656
424036000	Silver nitrate 1 mol/l (1N)	656
424051000	Silver nitrate 0.5 mol/l (0.5N)	656
424061000	Silver nitrate 0.1 mol/l (0.1N)	657
424062000	Silver nitrate 0.1 mol/l (0.1N)	657
424063000	Silver nitrate 0.1 mol/l (0.1N)	657
424067000	Silver nitrate 0.1 mol/l (0.1N)	657
424081	Silver nitrate 0.1 mol/l (0.1N)	657
424101000	Silver nitrate 0.05 mol/l (0.05N)	657
424111	Sulfuric acid 0.0025 mol/l (0.005N)	742
424161	Silver nitrate 0.01 mol/l (0.01N)	658
424181	Silver oxide	659
424182	Silver oxide	659
424191	Silver sulfate solution 0.7% in sulfuric acid	659
424192	Silver sulfate solution 0.7% in sulfuric acid	659
424201	Silver sulfate	659
424203	Silver sulfate	659
424268	L(+)Arginine monohydrochloride	184
424271	L(+)Arginine	184
424281	Arsenazo III	184
424282	Arsenazo III	184
424494	Sodium hydroxide on silica	680
424497	Sodium hydroxide on silica	680
424544	L(+)Asparagine	186
424547	L(+)Asparagine	186
424691	Azomethine H	189
424692	Azomethine H	189
424721	Azure II	189
424731	Azure II eosin	190
424861	Barium standard solution	192
424895	Barium acetate	192
424896	Barium acetate	192
424897	Barium acetate	192
424943	Barium carbonate	192
424945	Barium carbonate	192
425022	Barium chloride dihydrate	193
425025	Barium chloride dihydrate	193
425026	Barium chloride dihydrate	193
425027	Barium chloride dihydrate	193
425029	Barium chloride dihydrate	193
425245	Barium chromate	194
425292	Barium hydroxide octahydrate	194
425295	Barium hydroxide octahydrate	194
425296	Barium hydroxide octahydrate	194
425297	Barium hydroxide octahydrate	194
425341	Barium nitrate	195
425342	Barium nitrate	195
425347	Barium nitrate	195
425411	Barium perchlorate trihydrate	195
425497	Barium sulfate	196
426113	Benzene	198
426231	Benzenesulfonyl chloride	198
426761	Benzyl benzoate	200
426763	Benzyl benzoate	200
427001	Petroleum ether 60 - 80°C	565
427003	Petroleum ether 60 - 80°C	565
427007	Petroleum ether 60 - 80°C	565
427031	Petroleum ether 80 - 100°C	564
427036	Petroleum ether 80 - 100°C	564
427345	Dibenzoyl peroxide	296
427347	Dibenzoyl peroxide	296
428064	Bismuth, granules	202
428071	Bismuth standard solution	203
428103	Bismuth(III) carbonate basic	203
428105	Bismuth(III) carbonate basic	203
428234	Bismuth(III) nitrate pentahydrate	204
428236	Bismuth(III) nitrate pentahydrate	204
428284	Bismuth(III) nitrate basic	203
428286	Bismuth(III) nitrate basic	203
428294	Bismuth(III) nitrate basic	203
428432	Biuret 97%	205
428532	Alkali blue 6B	150
428551	Alcian Blue 8GS 1%	148
428561	Alcian blue 8GX	149
428582	Aniline blue soluble in water	180
428642	Coomassie brilliant blue R 250	274
428653	Bromophenol blue	213
428655	Bromophenol blue	213
428658	Bromophenol blue	213
428659	Bromophenol blue	213
428691	Bromophenol blue solution 0.02%	214
428702	Bromothymol blue	215
428703	Bromothymol blue	215
428708	Bromothymol blue	215
428731	Bromothymol blue 0.02%	215
428811	Brillant cresyl blue	208
428812	Brillant cresyl blue	208
428901	Lactophenol blue solution	452
428932	Methyl blue	497
428984	Methylene blue	498
428991	Carbolated Methylene Blue hydroalcoholic solution	248
429021	Methylthymol blue sodium salt	506
429022	Methylthymol blue sodium salt	506
429222	Thymol blue	757
429223	Thymol blue	757
429228	Thymol blue	757
429282	Toluidine blue	769
429291	Carbolated Toluidine Blue hydroalcoholic solution	248
429381	Victoria blue	785
429382	Victoria blue	785
429641	Boron standard solution	208
429751	Fixative Bouin liquid	355
429981	Methylene blue	498
429982	Methylene blue	498
430652	alpha-Bromonaphthalene	213
431252	Bismarck brown R	201
431601000	n-Butyl acetate	229
431602000	n-Butyl acetate	229
431604	n-Butyl acetate	229
431721	Isobutyl acetate	438
431811	n-Butyl chloride	229
431817	n-Butyl chloride	229
431821	n-Butyl chloride	229
432001	tert-Butylmethylether	230
432002	tert-Butylmethylether	230

INDEX BY PRODUCT CODE

432011	tert-Butylmethylether	231
432013	tert-Butylmethylether	231
432015	tert-Butylmethylether	231
432022000	tert-Butylmethylether	230
432031	tert-Butylmethylether	230
432032	tert-Butylmethylether	230
432034	tert-Butylmethylether	230
432061	tert-Butylmethylether	230
432062	tert-Butylmethylether	230
432121	Butylhydroxytoluene	230
432311	Cadmium standard solution	232
432344	Cadmium acetate dihydrate	233
432345	Cadmium acetate dihydrate	233
432347	Cadmium acetate dihydrate	233
432444	Cadmium carbonate	233
432446	Cadmium carbonate	233
432644	Cadmium nitrate tetrahydrate	234
432645	Cadmium nitrate tetrahydrate	234
432801	Soda lime	660
432802	Soda lime	660
432803	Soda lime	660
432851	Soda lime	660
432861	Soda lime	660
432862	Soda lime	660
432873	Soda lime	660
432874	Soda lime	660
432941	Calcium standard solution	235
432982	Calcium acetate monohydrate	236
432985	Calcium acetate monohydrate	236
432987	Calcium acetate monohydrate	236
433183	Calcium carbonate	237
433185	Calcium carbonate	237
433187	Calcium carbonate	237
433216	Calcium carbonate	236
433245	Calcium carbonate	237
433371	Calcium chloride hexahydrate	239
433375	Calcium chloride hexahydrate	239
433377	Calcium chloride hexahydrate	239
433381	Calcium chloride dihydrate	238
433382	Calcium chloride dihydrate	238
433384	Calcium chloride dihydrate	238
433386	Calcium chloride dihydrate	238
433387	Calcium chloride dihydrate	238
433403	Calcium chloride anhydrous	238
433405	Calcium chloride anhydrous	238
433406	Calcium chloride anhydrous	238
433407	Calcium chloride anhydrous	238
433535	Calcium chloride anhydrous	237
433585	Calcium fluoride	240
433587	Calcium fluoride	240
433685	Calcium phosphate monobasic monohydrate	243
433774	Calcium phosphate tribasic	243
433776	Calcium phosphate tribasic	243
433873	Calcium hydroxide	240
433874	Calcium hydroxide	240
433875	Calcium hydroxide	240
433877	Calcium hydroxide	240
433951	Calcium nitrate tetrahydrate	241
433955	Calcium nitrate tetrahydrate	241
433956	Calcium nitrate tetrahydrate	241
433957	Calcium nitrate tetrahydrate	241
434151	Calcium sulfate dihydrate	244
434155	Calcium sulfate dihydrate	244
434156	Calcium sulfate dihydrate	244
434171	Calcon	245
434181	Calmagite	246
434351	Haemalum solution according to Carazzi	378
434454	Charcoal activated	255
434455	Charcoal activated	255
434462	Charcoal activated	255
434501	Charcoal decolorizing	255
434507	Charcoal decolorizing	255
434766	Carborundum, granules	248
434786	Carborundum, powder	248
434932	Indigo carmine dried	422
434980000	Potassium iodide starch paper	608
435121000	Indicator papers	421
435131000	Indicator papers	420
435140000	Indicator papers	420
435150000	Indicator papers	420
435161000	Indicator papers	420
435220000	Congo red paper	273
435260000	Litmus paper	466
435300000	Litmus paper	466
435340000	Litmus paper	466
435421000	Indicator papers	420
435431000	Indicator papers	420
435441000	Indicator papers	420
435451000	Indicator papers	420
435493000	Indicator papers	421
435494000	Indicator papers	421
435496000	Indicator papers	421
435498000	Indicator papers	421
435502000	Indicator papers	421
435511000	Indicator papers	420
435642000	Indicator papers	421
435643000	Indicator papers	421
435644000	Indicator papers	421
435645000	Indicator papers	421
435922	Kovac reagent	451
435963	Casein	249
436061	Cellulose, powder	250
436081	Cerium (IV) ammonium nitrate	251
436082	Cerium (IV) ammonium nitrate	251
436091	Cerium (IV) ammonium sulfate dihydrate	252
436092	Cerium (IV) ammonium sulfate dihydrate	252
436203	Cerium (III) nitrate hexahydrate	251
436402	Cerium (IV) sulfate tetrahydrate	252
436404	Cerium (IV) sulfate tetrahydrate	252
436426	Cerium (IV) sulfate 0.1 mol/l	253
436501	Cesium chloride	254
436502	Cesium chloride	254
436534	Cesium sulfate	254
436701	Quinidine sulfate	634
436853	p-Benzoquinone	199
436854	p-Benzoquinone	199
436901	Cyclohexane	287
436902	Cyclohexane	287
436903	Cyclohexane	287
436905	Cyclohexane	287
436906	Cyclohexane	287
436908	Cyclohexane	287
436931	Cyclohexane	286
436932	Cyclohexane	286
436963	Cyclohexane	286
436967	Cyclohexane	286
437052	Cyclohexanone	288
437053	Cyclohexanone	288
437055	Cyclohexanone	288
437104	Cyclohexylamine	288
437251	Cinchonine	268
437308	L-Cysteine	289
437309	L-Cysteine	289
437351	L-Cystine	289
437355	L-Cystine	289
437401	Citral	268
437551	Chloramine T sodium salt	256
437554	Chloramine T sodium salt	256
437555	Chloramine T sodium salt	256
437557	Chloramine T sodium salt	256
437601	Chloranil	256
437704	Chloroacetamide	258
438203	p-Chlorobenzaldehyde	258
438251	Chlorobenzene	258
438253	Chlorobenzene	258
438255	Chlorobenzene	258
438504	4-Chlorophenol	264
438581	Chloroform	262
438582	Chloroform	262
438591	Chloroform	260
438592	Chloroform	260
438601	Chloroform	261
438603	Chloroform	261
438606	Chloroform	261
438607	Chloroform	261
438612	Chloroform	261
438613	Chloroform	261
438614	Chloroform	261
438641	Chloroform	260
438651	Chloroform	260
438652	Chloroform	260
438662	Chloroform	261
438664	Chloroform	261
438681	Chloroform	260
438682	Chloroform	260
438732	Chloroform	260
439131	Cobalt standard solution	271
439154	Cobalt (II) acetate tetrahydrate	271
439155	Cobalt (II) acetate tetrahydrate	271
439156	Cobalt (II) acetate tetrahydrate	271
439204	Cobalt (II) ammonium sulfate hexahydrate	271
439207	Cobalt (II) ammonium sulfate hexahydrate	271
439353	Cobalt (II) chloride hexahydrate	272
439355	Cobalt (II) chloride hexahydrate	272
439357	Cobalt (II) chloride hexahydrate	272
439455	Cobalt (II) nitrate hexahydrate	273
439457	Cobalt (II) nitrate hexahydrate	273
439504	Cobalt (II) nitrate hexahydrate	272
439705	Cobalt (II) sulfate heptahydrate	273
440381	o-Cresol	282
440385	o-Cresol	282
440387	o-Cresol	282
440435	m-Cresol	282
440437	m-Cresol	282
440572	Chrysoidine Y	268
440641	Chromium standard solution	265
440724	Chromium (III) chloride hexahydrate	265
440727	Chromium (III) chloride hexahydrate	265
440775	Chromium (III) nitrate nonahydrate	266
440776	Chromium (III) nitrate nonahydrate	266
440823	Chromium (III) oxide	266
440825	Chromium (III) oxide	266
440827	Chromium (III) oxide	266
440876	Chromium (III) potassium sulfate dodecahydrate	266
440877	Chromium (III) potassium sulfate dodecahydrate	266
440955	Chromium (III) sulfate	267
440957	Chromium (III) sulfate	267
441221	Decalcifying agent	290
441553	Diacetyldioxime	295
441623	Diacetyldioxime sodium salt	295
441625	Diacetyldioxime sodium salt	295
441771	Diacetone alcohol	295
441774	Diacetone alcohol	295
441873	Dibenzoylmethane	296
442261	Dichloromethane	302
442262	Dichloromethane	302
442291	Dichloromethane	301
442292000	Dichloromethane	301
442294	Dichloromethane	301
442301	o-Dichlorobenzene	298

Product specifications are subject to changes. Please visit our website for updates.

442353	m-Dichlorobenzene	298
442371	Dichloromethane	302
442406	p-Dichlorobenzene	298
442407	p-Dichlorobenzene	298
442458	2,6-Dichloroquinone-4-chlorimide	305
442508	2,6-Dichlorophenolindophenol sodium salt	305
442554	Diethanolamine	306
442557	Diethanolamine	306
442558	Diethanolamine	306
442756	Diethylamine	306
443056	Diethyl carbonate	307
443201000	1,4-Dioxane	321
443202000	1,4-Dioxane	321
443204000	1,4-Dioxane	321
443205	1,4-Dioxane	321
443206000	1,4-Dioxane	321
443231	1,4-Dioxane	321
443252	Diethylene glycol	307
443253	Diethylene glycol	307
443255	Diethylene glycol	307
443341	n,n-Diethyl-p-phenylenediamine sulfate	310
443404	Diethyl phthalate	310
443671	4-Diphenylaminesulfonic acid sodium salt	322
443752	sym-Diphenylcarbazide	323
443754	sym-Diphenylcarbazide	323
443801	sym-Diphenylcarbazone	323
444053	Diphenylthiocarbazone	323
444131	Differentiator for kit Gram-Hucker	311
444207	Digitonin	311
444252	Dimedone	313
444307	N,N-Dimethylacetamide	314
444308	N,N-Dimethylacetamide	314
444309	N,N-Dimethylacetamide	314
444311	N,N-Dimethylacetamide	313
444552	Methyl yellow	506
444603	p-Dimethylaminobenzaldehyde	314
444604	p-Dimethylaminobenzaldehyde	314
444678	p-Dimethylaminobenzylide-nerhodanine	314
444731	Neocuproine hydrochloride	519
444901	N,N-Dimethylformamide dimethylacetal	317
444923	N,N-Dimethylformamide	316
444925	N,N-Dimethylformamide	316
444926	N,N-Dimethylformamide	316
444928	N,N-Dimethylformamide	316
444941	N,N-Dimethylformamide	315
444942	N,N-Dimethylformamide	315
444956	N,N-Dimethylformamide	315
444957	N,N-Dimethylformamide	315
444981	N,N-Dimethylformamide	314
444982	N,N-Dimethylformamide	314
444991	N,N-Dimethylformamide	315
445101	Dimethylsulphoxide	318
445103	Dimethylsulphoxide	318
445106	Dimethylsulphoxide	318
445107	Dimethylsulphoxide	318
445111	Dimethylsulphoxide	318
445112	Dimethylsulphoxide	318
445121	Dimethylsulphoxide	318
445141	Dimethylsulphoxide	317
445142	Dimethylsulphoxide	317
445151	Isohexane	438
445152	Isohexane	438
445231	Dimidium bromide	320
445232	Dimidium bromide	320
445233	Dimidium bromide	320
445421	2,4-Dinitrochlorobenzene	320
445524	2,4-Dinitrophenylhydrazine (with 30% of water)	320
445551	Imidazole	419
445552	Imidazole	419
445958	2,2'-Dipyridyl	323
445981	Diisopropylamine	311
446008	Dithiooxamide	324
446181	Solvent Plus	714
446187	Solvent Plus	714
446372	Haemalum solution according to Mayer	378
446377	Haemalum solution according to Mayer	378
446461	Papanicolaou Haematoxylin solution according to Harris	554
446462	Papanicolaou Haematoxylin solution according to Harris	554
446463	Papanicolaou Haematoxylin solution according to Harris	554
446464	Papanicolaou Haematoxylin solution according to Harris	554
446465	Papanicolaou Haematoxylin solution according to Harris	554
446466	Papanicolaou Haematoxylin solution according to Harris	554
446472	Hematoxylin	379
446473	Hematoxylin	379
446475	Hematoxylin	379
446602	Eosin B	328
446632	Eosin Y	328
446634	Eosin Y	328
446644	Eosin Y 1% solution acqueous	328
446664	Eosin Y 0.5% solution alcoholic	329
446731	Hexamethyldisilazane	384
446781	n-Heptane 99%	381
446782	n-Heptane 99%	381
446783	n-Heptane 99%	381
446785	n-Heptane 99%	381
446787	n-Heptane 99%	381
446788	n-Heptane 99%	381
446789	n-Heptane 99%	381
446811	Eriochromocyanine R	330
446812	Eriochromocyanine R	330
446824	n-Heptane 99%	380
446831	n-Heptane	382
446832	n-Heptane	382
446841	Heptane mixture of isomers	382
446842	Heptane mixture of isomers	382
446875	Hexamethylenetetramine	384
446876	Hexamethylenetetramine	384
446891	Hexane mixture of isomers	388
446892	Hexane mixture of isomers	388
446893	Hexane mixture of isomers	388
446901	n-Hexane	387
446902	n-Hexane	387
446903	n-Hexane	387
446904	n-Hexane	387
446905	n-Hexane	387
446907	n-Hexane	387
446932	n-Hexane	387
446934	n-Hexane	387
446951	n-Heptane 99%	380
446952	n-Heptane 99%	380
446971	Erythrosin extra B	330
446972	Erythrosin extra B	330
446981	Esbach's reagent	330
446991	n-Hexane	387
447011	n-Hexane	386
447012	n-Hexane	386
447013	n-Hexane	386
447041	n-Hexane 99%	386
447042	n-Hexane 99%	386
447051	n-Hexane 99%	385
447052	n-Hexane 99%	385
447111	n-Hexane 99%	385
447112000	n-Hexane 99%	385
447121000	1,2-Dichloroethane	299
447131	Isohexane	439
447132	Isohexane	439
447181	Hexane mixture of isomers	388
447182	Hexane mixture of isomers	388
447191	1,2-Dichloroethane	299
447192	1,2-Dichloroethane	299
447212	n-Hexane 99%	385
447311	Isohexane	439
447312	Isohexane	439
447351	Ethanolamine	340
447352	Ethanolamine	340
447521	Diethyl ether	309
447522	Diethyl ether	309
447523	Diethyl ether	309
447525	Diethyl ether	309
447531	Diethyl ether	309
447532	Diethyl ether	309
447534	Diethyl ether	309
447539	Diethyl ether	309
447593	Diethyl ether	308
447651	Diethyl ether	308
447652	Diethyl ether	308
447792	Petroleum ether 30 - 40°C	569
447793	Petroleum ether 30 - 40°C	569
447795	Petroleum ether 30 - 40°C	569
447801	Petroleum ether 30 - 50°C	568
447802	Petroleum ether 30 - 50°C	568
447804	Petroleum ether 30 - 50°C	568
447811	Petroleum ether 40 - 65°C	566
447812	Petroleum ether 40 - 65°C	566
447813	Petroleum ether 40 - 65°C	566
447814	Petroleum ether 40 - 65°C	566
447815	Petroleum ether 40 - 65°C	566
447816	Petroleum ether 40 - 65°C	566
447821	Petroleum ether 40 - 70°C	566
447822	Petroleum ether 40 - 70°C	566
447824	Petroleum ether 40 - 70°C	566
447831	Petroleum ether 40 - 60°C	567
447832	Petroleum ether 40 - 60°C	567
447833	Petroleum ether 40 - 60°C	567
447834	Petroleum ether 40 - 60°C	567
447836	Petroleum ether 40 - 60°C	567
447851	Petroleum ether 40 - 65°C	566
447852	Petroleum ether 40 - 65°C	566
447861	Petroleum ether 35 - 60°C	567
447862	Petroleum ether 35 - 60°C	567
447931	Diisopropylether	312
447932	Diisopropylether	312
447933	Diisopropylether	312
447935	Diisopropylether	312
448211	Ethyl acetate	341
448251	Ethyl acetate	342
448252	Ethyl acetate	342
448253	Ethyl acetate	342
448254	Ethyl acetate	342
448255	Ethyl acetate	342
448256	Ethyl acetate	342
448258	Ethyl acetate	342
448271	Ethyl acetate	341
448272	Ethyl acetate	341
448306	Ethyl acetate	342
448307	Ethyl acetate	342
448308	Ethyl acetate	342
448342	Ethyl acetate	341
448351	Ethyl acetate	341
448352000	Ethyl acetate	341
448383	Ethyl acetate	340
448384	Ethyl acetate	340
449424	Ethylenediamine	344
449425	Ethylenediamine	344
449426	Ethylenediamine	344

INDEX BY PRODUCT CODE

449451	Ethylenediamine	344
449671	Tetrachloroethylene	750
449672	Tetrachloroethylene	750
449673	Tetrachloroethylene	750
449773	Eugenol	351
449895	Kieselguhr composed	449
449897	Kieselguhr composed	449
449926	Fehling's A reagent	353
449927	Fehling's A reagent	353
450038	o-Phenantroline monohydrate	569
450039	o-Phenantroline monohydrate	569
450328	L-Phenylalanine	572
450329	L-Phenylalanine	572
450843	Phenylhydrazine hydrochloride	573
451154	Phenolphthalein	571
451156	Phenolphthalein	571
451191	Phenolphthalein solution 1% in ethanol	572
451192	Phenolphthalein solution 1% in ethanol	572
451271	Phenol	570
451311	Iron standard solution	429
451373	Iron, powder	427
451377	Iron, powder	427
451395	Iron, reduced by hydrogen	428
451397	Iron, reduced by hydrogen	428
451451	Iron (II) ammonium sulfate hexahydrate	429
451452	Iron (II) ammonium sulfate hexahydrate	429
451453	Iron (II) ammonium sulfate hexahydrate	429
451457	Iron (II) ammonium sulfate hexahydrate	429
451502	Iron (III) ammonium sulfate dodecahydrate	432
451503	Iron (III) ammonium sulfate dodecahydrate	432
451504	Iron (III) ammonium sulfate dodecahydrate	432
451505	Iron (III) ammonium sulfate dodecahydrate	432
451507	Iron (III) ammonium sulfate dodecahydrate	432
451573	Iron (II) chloride tetrahydrate	430
451574	Iron (II) chloride tetrahydrate	430
451575	Iron (II) chloride tetrahydrate	430
451576	Iron (II) chloride tetrahydrate	430
451625	Iron (III) chloride hexahydrate	433
451626	Iron (III) chloride hexahydrate	433
451627	Iron (III) chloride hexahydrate	433
451692	Iron (III) chloride anhydrous sublimed	433
451695	Iron (III) chloride anhydrous sublimed	433
451696	Iron (III) chloride anhydrous sublimed	433
451722	Iron (III) nitrate nonahydrate	435
451723	Iron (III) nitrate nonahydrate	435
451725	Iron (III) nitrate nonahydrate	435
451727	Iron (III) nitrate nonahydrate	435
451824	Iron (III) oxide	435
451826	Iron (III) oxide	435
451877	Iron (II) sulfate heptahydrate	430
451878	Iron (II) sulfate heptahydrate	430
451879	Iron (II) sulfate heptahydrate	430
451925	Iron (III) sulfate	435
451926	Iron (III) sulfate	435
451927	Iron (III) sulfate	435
452031	Phloroglucinol	574
452033	Phloroglucinol	574
452051	Phloxin B	574
452052	Phloxin B	574
452083	Fluorescein	357
452086	Fluorescein	357
452087	Fluorescein	357
452112	Fluorescein sodium salt	357
452113	Fluorescein sodium salt	357
452117	Fluorescein sodium salt	357
452271	Florisil 60-100 mesh	356
452273	Florisil 60-100 mesh	356
452282	Formamide	362
452286	Formamide	362
452331	Florisil 60-100 mesh	356
452332	Florisil 60-100 mesh	356
452333	Florisil 60-100 mesh	356
452351	Florisil 100-200 mesh	357
452353	Florisil 100-200 mesh	357
452665	D(-)Fructose	364
452666	D(-)Fructose	364
452751	o-Phthalaldehyde	579
452812	Fuchsin acid	364
452814	Fuchsin acid	364
452842	Fuchsin basic	365
452844	Fuchsin basic	365
453125	D(+)Galactose	366
453126	D(+)Galactose	366
453226	Gelatine	368
453272	Silica gel granular	650
453273	Silica gel granular	650
453275	Silica gel granular	650
453277	Silica gel granular	650
453278	Silica gel granular	650
453279	Silica gel granular	650
453315	Silica gel granular with indicator cobalt free	651
453317	Silica gel granular with indicator cobalt free	651
453319	Silica gel granular with indicator cobalt free	651
453331	Silica gel 60A 0,06÷0,20 mm	650
453332	Silica gel 60A 0,06÷0,20 mm	650
453336	Silica gel 60A 0,06÷0,20 mm	650
453337	Silica gel 60A 0,06÷0,20 mm	650
453351	Silica gel 60A 35 - 70µ	649
453352	Silica gel 60A 35 - 70µ	649
453353	Silica gel 60A 35 - 70µ	649
453355	Silica gel 60A 35 - 70µ	649
453451	Alizarin yellow R	150
453518	Clayton's yellow	270
453519	Clayton's yellow	270
453542	Metanil yellow	487
453562	Naphthol yellow S	518
453581	Sudan yellow	727
453582	Sudan yellow	727
453611	Giemsa's reagent	369
453614	Giemsa's reagent	369
453616	Giemsa's reagent	369
453751	Glycerol (30°Bé)	372
453752	Glycerol (30°Bé)	372
453755	Glycerol (30°Bé)	372
453758	Glycerol (30°Bé)	372
453759	Glycerol (30°Bé)	372
453771	Glycerol (30°Bé)	372
453772	Glycerol (30°Bé)	372
453804	Glycine	373
453807	Glycine	373
453881	2-(2-Butoxyethoxy)ethanol	228
453883	2-(2-Butoxyethoxy)ethanol	228
453902	Ethylene glycol	348
453904	Ethylene glycol	348
453905	Ethylene glycol	348
453906	Ethylene glycol	348
453941	2-Butoxy ethanol	228
454021	2-Methoxy ethanol	495
454023	2-Methoxy ethanol	495
454024	2-Methoxy ethanol	495
454028	2-Methoxy ethanol	495
454052	Propylene glycol	629
454053	Propylene glycol	629
454054	Propylene glycol	629
454111	Triethylene glycol	773
454112	Triethylene glycol	773
454131	Glyoxal-bis-(2-hydroxyanil)	374
454132	Glyoxal-bis-(2-hydroxyanil)	374
454333	D(+)Glucose anhydrous	370
454335	D(+)Glucose anhydrous	370
454336	D(+)Glucose anhydrous	370
454337	D(+)Glucose anhydrous	370
454338	D(+)Glucose anhydrous	370
454441	Gram - Hucker Kit	376
454452	Griess' reagent A	376
454462	Griess' reagent B	376
454481	Griess' reagent	376
454911	Histolemon	390
454912	Histolemon	390
454915	Histolemon	390
454921	Hyamine 1622	391
455054	Hydrazine dihydrochloride	392
455056	Hydrazine dihydrochloride	392
455057	Hydrazine dihydrochloride	392
455291	Hydrindantin	392
455325	Hydroquinone	416
455523	Hydroxylamine sulphate	417
455525	Hydroxylamine sulphate	417
455527	Hydroxylamine sulphate	417
455621	Indicator for iodometry	420
455622	Indicator for iodometry	420
455801	Indole	423
455853	Inositol	423
455901	Inulin	424
455902	Inulin	424
455903	Inulin	424
455954	Iodine resublimed	425
455955	Iodine resublimed	425
455957	Iodine resublimed	425
455959	Iodine resublimed	425
456036000	Iodine 0.05 mol/l (0.1N)	425
456037000	Iodine 0.05 mol/l (0.1N)	425
456051	Iodine 0.05 mol/l (0.1N)	426
456121	Iodine 0.005 mol/l (0.01N)	426
456135000	Iodine 0.5 mol/l (1N)	425
456137000	Iodine 0.5 mol/l (1N)	425
456641	Isooctane	440
456732	Isooctane	440
456734	Isooctane	440
456753	Isooctane	439
456754	Isooctane	439
456791	Isooctane	439
456792	Isooctane	439
456851	Histamine dihydrochloride	389
456852	Histamine dihydrochloride	389
456951	L-Histidine hydrochloride monohydrate	390
456952	L-Histidine hydrochloride monohydrate	390
457502	Lanthanum nitrate hexahydrate	454
457506	Lanthanum nitrate hexahydrate	454
457511	Lanthanum oxide	455
457521	Glass wool	370
457531	Amman's lactophenol solution	160
457551	Lactose monohydrate	452
457552	Lactose monohydrate	452
457553	Lactose monohydrate	452
457557	Lactose monohydrate	452
457625	Devarda's alloy	294
457627	Devarda's alloy	294
457675	Raney's alloy	635
457928	L(+)-Leucine	461
458001	Petroleum ether 75 - 120°C	565
458003	Petroleum ether 75 - 120°C	565
458121	L(+)Lysine monohydrochloride	467

Product specifications are subject to changes. Please visit our website for updates.

INDEX BY PRODUCT CODE

458122	L(+)Lysine monohydrochloride	467
458124	L(+)Lysine monohydrochloride	467
458163	Lithium tetraborate anhydrous	465
458164	Lithium tetraborate anhydrous	465
458204	Lithium carbonate	463
458207	Lithium carbonate	463
458211	Lithium standard solution	462
458254	Lithium chloride	463
458256	Lithium chloride	463
458257	Lithium chloride	463
458271	Lithium chloride	463
458272	Lithium chloride	463
458273	Lithium chloride	463
458275	Lithium chloride	463
458281	Lithium hydroxide anhydrous	464
458282	Lithium hydroxide anhydrous	464
458292	Lithium hydroxide monohydrate	464
458303	Lithium hydride	464
458355	Lithium nitrate	465
458356	Lithium nitrate	465
458404	Lithium sulfate monohydrate	465
458405	Lithium sulfate monohydrate	465
458741	Lugol concentrated solution	466
458751	Lugol solution for Gram-Hucker kit	466
458762	Lugol's Reagent Iodine-Iodide Solution	466
458763	Lugol's Reagent Iodine-Iodide Solution	466
458772	Luminol	467
458891	Magnesium standard solution	469
459044	Magnesium, ribbon	468
459066	Magnesium, powder	468
459085	Magnesium, turnings	468
459131	Magnesium acetate tetrahydrate	470
459135	Magnesium acetate tetrahydrate	470
459137	Magnesium acetate tetrahydrate	470
459285	Magnesium carbonate basic	470
459287	Magnesium carbonate basic	470
459331	Magnesium chloride hexahydrate	471
459332	Magnesium chloride hexahydrate	471
459334	Magnesium chloride hexahydrate	471
459336	Magnesium chloride hexahydrate	471
459337	Magnesium chloride hexahydrate	471
459435	Magnesium hydrogen phosphate trihydrate	472
459437	Magnesium hydrogen phosphate trihydrate	472
459535	Magnesium nitrate hexahydrate	472
459536	Magnesium nitrate hexahydrate	472
459537	Magnesium nitrate hexahydrate	472
459584	Magnesium oxide	473
459586	Magnesium oxide	473
459587	Magnesium oxide	473
459617	Magnesium oxide	473
459662	Magnesium sulfate heptahydrate	475
459665	Magnesium sulfate heptahydrate	475
459666	Magnesium sulfate heptahydrate	475
459667	Magnesium sulfate heptahydrate	475
459669	Magnesium sulfate heptahydrate	475
459863	Maltose monohydrate	478
459865	Maltose monohydrate	478
459911	Manganese standard solution	480
459965	Manganese electrolytic	479
460001	Manganese (II) acetate tetrahydrate	480
460005	Manganese (II) acetate tetrahydrate	480
460007	Manganese (II) acetate tetrahydrate	480
460052	Manganese (IV) oxide	481
460055	Manganese (IV) oxide	481
460056	Manganese (IV) oxide	481
460156	Manganese (II) chloride tetrahydrate	480
460158	Manganese (II) chloride tetrahydrate	480
460159	Manganese (II) chloride tetrahydrate	480
460305	Manganese (II) sulfate monohydrate	481
460307	Manganese (II) sulfate monohydrate	481
460352	D-Mannitol	481
460353	D-Mannitol	481
460355	D-Mannitol	481
460357	D-Mannitol	481
460502	Mayer's reagent	482
460511	Hematoxylin solution according to Mayer	379
460512	Hematoxylin solution according to Mayer	379
460513	Hematoxylin solution according to Mayer	379
460515	Hematoxylin solution according to Mayer	379
460581	May Grünwald reagent	482
460584	May Grünwald reagent	482
460586	May Grünwald reagent	482
460691	2-Mercaptoethanol	483
460741	Mercury standard solution	485
461003	Mercury (II) chloride	485
461105	Mercury (II) iodide	486
461325	Mercury (II) oxide red	486
461405	Mercury (II) sulfate	486
461805	4-Methylaminophenol sulfate	496
461943	Methyl isobutyl ketone	500
461945	Methyl isobutyl ketone	500
462017	Methyl acetate	496
462207	Methyl benzoate	496
462238	3-Methyl-2-benzothiazolinone hydrazone hydrochloride	497
462601	Methyl iodide	500
462604	Methyl iodide	500
462701	Ethyl methyl ketone	350
462702	Ethyl methyl ketone	350
462703	Ethyl methyl ketone	350
462704	Ethyl methyl ketone	350
462872	N-Methyl-2-pyrrolidone	503
462873	N-Methyl-2-pyrrolidone	503
462874	N-Methyl-2-pyrrolidone	503
462875	N-Methyl-2-pyrrolidone	503
462881	N-Methyl-2-pyrrolidone	502
463001	Dichloromethane	303
463002	Dichloromethane	303
463003	Dichloromethane	303
463008	Dichloromethane	303
463025	Dichloromethane	302
463126	DL-Methionine	494
463161	Dichloromethane	303
463162	Dichloromethane	303
463251	Mix Diethyl ether/Ethanol 70/30 w/w	507
463255	Mix Diethyl ether/Ethanol 70/30 w/w	507
463281	Dichloromethane	300
463291	Dichloromethane	300
463311	Dichloromethane	303
463314	Dichloromethane	303
463318	Dichloromethane	303
463332	Dichloromethane	301
463342	Dichloromethane	301
463421	Selenic mixture	646
463422	Selenic mixture	646
463431	Molybdenum standard solution	511
463453	Morpholine	512
463531	Mucicarmine hydroalcoholic solution	512
463608	Murexide	516
463609	Murexide	516
463651	Naphthalene	517
463654	Naphthalene	517
463655	Naphthalene	517
463831	N-(1-Naphtyl)ethylenediamine dihydrochloride	518
463891	alpha-Naphtholbenzein	518
463935	1-Naphthol	517
463984	2-Naphthol	517
463986	2-Naphthol	517
464221	Eriochrome black T	330
464222	Eriochrome black T	330
464231	Nessler's reagent single solution	520
464232	Nessler's reagent single solution	520
464241	Sudan black B	726
464271	Nickel standard solution	522
464384	Nickel, powder	521
464422	Nessler's reagent solution A	520
464432	Nessler's reagent solution B	521
464474	Nickel (II) acetate tetrahydrate	522
464476	Nickel (II) acetate tetrahydrate	522
464477	Nickel (II) acetate tetrahydrate	522
464545	Nickel (II) ammonium sulfate hexahydrate	523
464547	Nickel (II) ammonium sulfate hexahydrate	523
464604	Nickel (II) carbonate basic	523
464605	Nickel (II) carbonate basic	523
464644	Nickel (II) chloride hexahydrate	523
464645	Nickel (II) chloride hexahydrate	523
464647	Nickel (II) chloride hexahydrate	523
464685	Nickel (II) nitrate hexahydrate	524
464686	Nickel (II) nitrate hexahydrate	524
464772	Nickel (II) sulfate hexahydrate	524
464775	Nickel (II) sulfate hexahydrate	524
464777	Nickel (II) sulfate hexahydrate	524
464852	Nigrosine	525
464853	Nigrosine	525
464922	Ninhydrin	526
464928	Ninhydrin	526
465142	m-Nitrobenzaldehyde	534
465222	Nitrobenzene	535
465744	p-Nitrophenol	535
466361	Nonylphenol ethoxylated 10 ETO	536
466753	Oil of cedar wood	538
466757	Oil of cedar wood	538
466782	Immersion oil	420
466783	Immersion oil	420
466792	Paraffin oil	555
466831	Holmium oxide	391
466858	Orcein	540
466859	Orcein	540
466908	Orcinol monohydrate	540
466961	Gold standard solution	375
467007	Gold(III) chloride trihydrate	375
467254	p-Hydroxybenzaldehyde	416
467353	8-Hydroxyquinoline	417

INDEX BY PRODUCT CODE

467355	8-Hydroxyquinoline	417
467356	8-Hydroxyquinoline	417
467562	n-Octane	537
467731	Palladium (II) chloride	553
467737	Palladium (II) chloride	553
467781	Papanicolaou solution EA 50	554
467782	Papanicolaou solution EA 50	554
467783	Papanicolaou solution EA 50	554
467791	Papanicolaou solution OG 6	555
467792	Papanicolaou solution OG 6	555
467793	Papanicolaou solution OG 6	555
467958	Paraffin 56-58°C - Erbaplast (without DMSO)	555
468121	n-Pentane	559
468122	n-Pentane	559
468123	n-Pentane	559
468141	n-Pentane 99%	557
468142	n-Pentane 99%	557
468151000	n-Pentane 99%	557
468161	n-Pentane	558
468162	n-Pentane	558
468172	n-Pentane 99%	557
468182	n-Pentane	558
468204	n-Pentane 99%	557
468245	Pepsin HCl	560
468791	Lead standard solution	456
468866	Lead	455
468932	Lead (II) acetate trihydrate	458
468934	Lead (II) acetate trihydrate	458
468935	Lead (II) acetate trihydrate	458
468937	Lead (II) acetate trihydrate	458
468984	Lead (II) acetate basic	457
468985	Lead (II) acetate basic	457
468987	Lead (II) acetate basic	457
469055	Lead (IV) oxide	460
469057	Lead (IV) oxide	460
469353	Lead (II) nitrate	458
469355	Lead (II) nitrate	458
469356	Lead (II) nitrate	458
469357	Lead (II) nitrate	458
469404	Lead (II) oxide	460
469505	Lead (II) sulfate	460
469506	Lead (II) sulfate	460
469551	Piperidine	580
469552	Piperidine	580
469592	1-(2-Pyridylazo)-2-naphthol	632
469621	Pyridine	631
469622	Pyridine	631
469623	Pyridine	631
469624	Pyridine	631
469626	Pyridine	631
469629	Pyridine	631
469651	Pyridine	631
469652	Pyridine	631
469753	Pyrocatechol	632
469971	Pumice stone	630
470017	Hexachloroplatinic acid hexahydrate	383
470038	Bromocresol purple	212
470039	Bromocresol purple	212
470067	m-Cresol purple	282
470068	m-Cresol purple	282
470071	Polyvinylpyrrolidone	581
470072	Polyvinylpyrrolidone	581
470081	Potassium standard solution	582
470143	Potassium acetate	583
470145	Potassium acetate	583
470146	Potassium acetate	583
470147	Potassium acetate	583
470285	Potassium bicarbonate	584
470286	Potassium bicarbonate	584
470287	Potassium bicarbonate	584
470289	Potassium bicarbonate	584
470336	Potassium dichromate	592
470337	Potassium dichromate	592
470451	Potassium dichromate 0.0417 mol/l (0.25N)	593
470501	Potassium dichromate 0.0167 mol/l (0.1N)	593
470552	Potassium bisulfate	584
470556	Potassium bisulfate	584
470557	Potassium bisulfate	584
470655	Potassium bromate	584
470681	Potassium bromate 0.0167 mol/l (0.1N)	585
470701	Potassium bromide	585
470733	Potassium bromide	585
470734	Potassium bromide	585
470735	Potassium bromide	585
470737	Potassium bromide	585
470801	Potassium carbonate	586
470805	Potassium carbonate	586
470807	Potassium carbonate	586
471025	Potassium citrate tribasic monohydrate	591
471027	Potassium citrate tribasic monohydrate	591
471127	di-Potassium hexachloroplatinate	597
471171	Potassium chloride	587
471173	Potassium chloride	587
471175	Potassium chloride	587
471176	Potassium chloride	587
471177	Potassium chloride	587
471181	Potassium chloride	586
471215	Potassium chloride 3 mol/l (3N)	588
471225	Potassium chloride 3.5 mol/l (3.5N)	587
471235	Potassium chloride 3 mol/l (3N) + silver chloride	588
471245	Potassium chloride 3.5 mol/l (3.5N) + silver chloride	587
471255	Potassium chloride 0.01 mol/l (0.01N)	589
471265	Potassium chloride saturated solution	590
471275	Potassium chloride 3 mol/l (3N) water-glycerol solution	588
471285	Potassium chloride solution	590
471295	Potassium chromate	590
471297	Potassium chromate	590
471364	Potassium ferricyanide	594
471365	Potassium ferricyanide	594
471367	Potassium ferricyanide	594
471483	Potassium ferrocyanide trihydrate	595
471484	Potassium ferrocyanide trihydrate	595
471485	Potassium ferrocyanide trihydrate	595
471487	Potassium ferrocyanide trihydrate	595
471488	Potassium ferrocyanide trihydrate	595
471555	Potassium fluoride dihydrate	596
471561	Potassium fluoride	596
471562	Potassium fluoride	596
471563	Potassium fluoride	596
471564	Potassium fluoride	596
471681	Potassium phosphate monobasic	615
471682	Potassium phosphate monobasic	615
471685	Potassium phosphate monobasic	615
471686	Potassium phosphate monobasic	615
471687	Potassium phosphate monobasic	615
471761	Potassium phosphate dibasic trihydrate	614
471764	Potassium phosphate dibasic trihydrate	614
471766	Potassium phosphate dibasic trihydrate	614
471767	Potassium phosphate dibasic trihydrate	614
471781	Potassium phosphate dibasic anhydrous	614
471782	Potassium phosphate dibasic anhydrous	614
471785	Potassium phosphate dibasic anhydrous	614
471786	Potassium phosphate dibasic anhydrous	614
471787	Potassium phosphate dibasic anhydrous	614
471788	Potassium phosphate dibasic anhydrous	614
471865	Potassium hydrogen phthalate	598
471866	Potassium hydrogen phthalate	598
471867	Potassium hydrogen phthalate	598
471913	Potassium hydrogen phthalate	598
472021000	Potassium hydroxide 0.5 mol/l (0.5N) in ethanol	603
472022000	Potassium hydroxide 0.5 mol/l (0.5N) in ethanol	603
472041000	Potassium hydroxide 0.1 mol/l (0.1N) in ethanol	605
472042000	Potassium hydroxide 0.1 mol/l (0.1N) in ethanol	605
472056	Potassium hydroxide, pellets	600
472057	Potassium hydroxide, pellets	600
472059	Potassium hydroxide, pellets	600
472086	Potassium hydroxide, pellets	599
472092	Potassium hydroxide, pellets	599
472097	Potassium hydroxide, pellets	599
472103	Potassium hydroxide solution 45%	600
472171	Potassium hydroxide, pellets	600
472172	Potassium hydroxide, pellets	600
472173	Potassium hydroxide, pellets	600
472175	Potassium hydroxide, pellets	600
472281000	Potassium hydroxide 1 mol/l (1N)	602
472282000	Potassium hydroxide 1 mol/l (1N)	602
472287000	Potassium hydroxide 1 mol/l (1N)	602
472311	Potassium hydroxide 1 mol/l (1N)	602
472331000	Potassium hydroxide 0.5 mol/l (0.5N)	602
472332000	Potassium hydroxide 0.5 mol/l (0.5N)	602
472337000	Potassium hydroxide 0.5 mol/l (0.5N)	602
472364000	Potassium hydroxide 0.5 mol/l (0.5N) in methanol	603
472366000	Potassium hydroxide 0.5 mol/l (0.5N) in methanol	603
472367000	Potassium hydroxide 0.5 mol/l (0.5N) in methanol	603
472391	Potassium hydroxide 0.5 mol/l (0.5N)	603
472421000	Potassium hydroxide 0.25 mol/l (0.25N)	604
472422000	Potassium hydroxide 0.25 mol/l (0.25N)	604

Product specifications are subject to changes. Please visit our website for updates.

472427000	Potassium hydroxide 0.25 mol/l (0.25N)	604
472451000	Potassium hydroxide 0.1 mol/l (0.1N)	605
472452000	Potassium hydroxide 0.1 mol/l (0.1N)	605
472457000	Potassium hydroxide 0.1 mol/l (0.1N)	605
472484000	Potassium hydroxide 0.1 mol/l (0.1N) in methanol	605
472486000	Potassium hydroxide 0.1 mol/l (0.1N) in methanol	605
472511	Potassium hydroxide 0.1 mol/l (0.1N)	605
472563	Potassium Iodate	606
472565	Potassium Iodate	606
472601	Potassium iodate 0.0167 mol/l (0.1N)	606
472631	Potassium iodate 0.00167 mol/l (0.01N)	607
472641	Potassium hydrogen iodate	597
472735	Potassium iodide	607
472736	Potassium iodide	607
472737	Potassium iodide	607
472815000	Potassium iodide solution 3.9%	608
472821	Potassium iodide	607
472831	Potassium iodide solution 10%	608
473001	Potassium nitrate	609
473006	Potassium nitrate	609
473007	Potassium nitrate	609
473009	Potassium nitrate	609
473045	Potassium nitrate 1 mol/l (1N)	610
473084	Potassium nitrite	610
473133	Potassium oxalate monohydrate	610
473135	Potassium oxalate monohydrate	610
473137	Potassium oxalate monohydrate	610
473332	Potassium m-periodate	611
473334	Potassium m-periodate	611
473381	Potassium permanganate	611
473384	Potassium permanganate	611
473385	Potassium permanganate	611
473387	Potassium permanganate	611
473514000	Potassium permanganate 0.2 mol/l (1N)	612
473565000	Potassium permanganate 0.02 mol/l (0.1N)	612
473567000	Potassium permanganate 0.02 mol/l (0.1N)	612
473591	Potassium permanganate 0.02 mol/l (0.1N)	612
473661	Potassium permanganate 0.002 mol/l (0.01N)	613
473701	Potassium persulfate	613
473835	Potassium pyroantimonate acid	615
473915	tetra-Potassium pyrophosphate	616
474016	Potassium pyrosulphate	616
474112	Potassium sodium tartrate tetrahydrate	617
474114	Potassium sodium tartrate tetrahydrate	617
474115	Potassium sodium tartrate tetrahydrate	617
474116	Potassium sodium tartrate tetrahydrate	617
474117	Potassium sodium tartrate tetrahydrate	617
474119	Potassium sodium tartrate tetrahydrate	617
474166	Potassium sulfate	618
474167	Potassium sulfate	618
474169	Potassium sulfate	618
474205	Potassium sulfate	618
474355	Potassium thiocyanate	620
474357	Potassium thiocyanate	620
474465	Potassium tartrate	618
474467	Potassium tartrate	618
474515	Potassium L-tartrate monobasic	619
474517	Potassium L-tartrate monobasic	619
474652	Potassium tellurite	619
474708	L(-)Proline	622
474756	Isopropylamine	442
474807	n-Propyl acetate	628
474821	Isopropyl acetate	441
474871	Propylene carbonate	629
475132	Raffinose	635
475151	Copper standard solution	276
475185	Copper electrolytic, wire	275
475187	Copper electrolytic, wire	275
475215	Copper electrolytic rebaked, sheet	274
475305	Copper electrolytic, turnings	274
475307	Copper electrolytic, turnings	274
475334	Copper reduced, powder	275
475337	Copper reduced, powder	275
475405	Copper (II) acetate hydrate	277
475407	Copper (II) acetate hydrate	277
475553	Copper (II) carbonate (basic)	278
475555	Copper (II) carbonate (basic)	278
475557	Copper (II) carbonate (basic)	278
475605	Copper (I) chloride	276
475607	Copper (I) chloride	276
475685	Copper (II) chloride dihydrate	278
475687	Copper (II) chloride dihydrate	278
475782	Copper (II) nitrate trihydrate	279
475783	Copper (II) nitrate trihydrate	279
475784	Copper (II) nitrate trihydrate	279
475786	Copper (II) nitrate trihydrate	279
475966	Copper (II) oxide wire	279
475994	Copper (II) oxide	279
475997	Copper (II) oxide	279
476092	Copper (II) sulfate pentahydrate	280
476096	Copper (II) sulfate pentahydrate	280
476097	Copper (II) sulfate pentahydrate	280
476099	Copper (II) sulfate pentahydrate	280
476154	Copper (II) sulfate pentahydrate	280
476243	Copper (II) sulfate anhydrous	280
476245	Copper (II) sulfate anhydrous	280
476247	Copper (II) sulfate anhydrous	280
476312	L(+)Rhamnose	638
476565	Resorcinol	637
476608	D(-)Ribose	639
476671	Potassium permanganate	611
476687	Quinaldine red	634
476688	Quinaldine red	634
476762	Congo red	273
476764	Congo red	273
476778	o-Cresol red	283
476838	Phenol red	570
476839	Phenol red	570
476881	Methyl red	503
476882	Methyl red	503
476883	Methyl red	503
476941	Ponceau red BS	581
476951	Neutral red	521
476961	Red for oils O	636
476981	Ponceau red S	581
476982	Ponceau red S	581
477011	Nuclear fast red	536
477012	Nuclear fast red	536
477153	Sand purified	644
477182	D(+)-Sucrose	726
477183	D(+)-Sucrose	726
477186	D(+)-Sucrose	726
477187	D(+)-Sucrose	726
477232	Safranine T	642
477241	Safranine T hydroalcoholic solution for Gram-Hucker Kit	642
477591	Schiff's reagent for PAS coloration	646
477592	Schiff's reagent for PAS coloration	646
477601	Schiff's reagent for Aldehydes	645
477691	Selenium standard solution	647
477702	Selenium, powder	646
477721	Molecular sieves 3 A	510
477731	Molecular sieves 3 A	510
477762	Selenium dioxide	648
477961	Silicon standard solution	652
478032	Redox solution 220 mV at 25°C	636
478052	Redox solution 468 mV at 25°C	637
478101	Sodium standard solution	661
478132	Sodium acetate trihydrate	662
478135	Sodium acetate trihydrate	662
478136	Sodium acetate trihydrate	662
478137	Sodium acetate trihydrate	662
478139	Sodium acetate trihydrate	662
478163	Sodium acetate anhydrous	662
478165	Sodium acetate anhydrous	662
478166	Sodium acetate anhydrous	662
478167	Sodium acetate anhydrous	662
478232	Sodium aluminate	663
478237	Sodium aluminate	663
478357	Sodium ammonium hydrogen phosphate	664
478481	Sodium azide	665
478482	Sodium azide	665
478484	Sodium azide	665
478531	Sodium bicarbonate	665
478532	Sodium bicarbonate	665
478535	Sodium bicarbonate	665
478536	Sodium bicarbonate	665
478537	Sodium bicarbonate	665
478673	Sodium bisulfate monohydrate	666
478675	Sodium bisulfate monohydrate	666
478676	Sodium bisulfate monohydrate	666
478677	Sodium bisulfate monohydrate	666
478812	Sodium tetraborate decahydrate	710
478815	Sodium tetraborate decahydrate	710
478816	Sodium tetraborate decahydrate	710
478817	Sodium tetraborate decahydrate	710
478819	Sodium tetraborate decahydrate	710
478953	Sodium borohydride	666
478955	Sodium borohydride	666
478957	Sodium borohydride	666
478964	Sodium borohydride	666
479055	Sodium bromide	667
479057	Sodium bromide	667
479121	Sodium carbonate decahydrate	668
479122	Sodium carbonate decahydrate	668
479125	Sodium carbonate decahydrate	668
479126	Sodium carbonate decahydrate	668
479127	Sodium carbonate decahydrate	668
479151	Sodium carbonate solution 20%	670
479186	Sodium carbonate 0.5 mol/l (1N)	669
479211	Sodium carbonate 0.05 mol/l (0.1N)	669
479255	Sodium carbonate monohydrate	669

INDEX BY PRODUCT CODE

479256	Sodium carbonate monohydrate	669
479257	Sodium carbonate monohydrate	669
479301	Sodium carbonate anhydrous	667
479302	Sodium carbonate anhydrous	667
479305	Sodium carbonate anhydrous	667
479306	Sodium carbonate anhydrous	667
479307	Sodium carbonate anhydrous	667
479331	Sodium carbonate anhydrous	667
479371	Sodium cyanoborohydride	674
479484	Sodium citrate tribasic dihydrate	673
479485	Sodium citrate tribasic dihydrate	673
479486	Sodium citrate tribasic dihydrate	673
479487	Sodium citrate tribasic dihydrate	673
479488	Sodium citrate tribasic dihydrate	673
479652	Sodium chloride	670
479661	Sodium chloride	671
479662	Sodium chloride	671
479663	Sodium chloride	671
479671	Sodium chloride	670
479681	Sodium chloride	671
479685	Sodium chloride	671
479686	Sodium chloride	671
479687	Sodium chloride	671
479689	Sodium chloride	671
479781	Sodium chloride 0.1 mol/l (0,1N)	672
479833	Sodium cobalt nitrite	673
479911	Tris (hydroxymethyl)-aminomethane hydrochloride	778
479912	Tris (hydroxymethyl)-aminomethane hydrochloride	778
479913	Tris (hydroxymethyl)-aminomethane hydrochloride	778
479921	Sodium dichloroisocyanurate dihydrate	674
479954	Sodium fluoride	675
479955	Sodium fluoride	675
479957	Sodium fluoride	675
480005	Sodium hexafluorosilicate	677
480045	Sodium formate	675
480046	Sodium formate	675
480081	Sodium phosphate monobasic monohydrate	703
480082	Sodium phosphate monobasic monohydrate	703
480085	Sodium phosphate monobasic monohydrate	703
480086	Sodium phosphate monobasic monohydrate	703
480087	Sodium phosphate monobasic monohydrate	703
480131	Sodium phosphate dibasic dodecahydrate	702
480132	Sodium phosphate dibasic dodecahydrate	702
480133	Sodium phosphate dibasic dodecahydrate	702
480135	Sodium phosphate dibasic dodecahydrate	702
480136	Sodium phosphate dibasic dodecahydrate	702
480137	Sodium phosphate dibasic dodecahydrate	702
480141	Sodium phosphate dibasic anhydrous	700
480142	Sodium phosphate dibasic anhydrous	700
480143	Sodium phosphate dibasic anhydrous	700
480144	Sodium phosphate dibasic anhydrous	700
480222	Sodium phosphate dibasic dihydrate	701
480225	Sodium phosphate dibasic dihydrate	701
480226	Sodium phosphate dibasic dihydrate	701
480227	Sodium phosphate dibasic dihydrate	701
480271	Sodium phosphate tribasic dodecahydrate	704
480272	Sodium phosphate tribasic dodecahydrate	704
480275	Sodium phosphate tribasic dodecahydrate	704
480276	Sodium phosphate tribasic dodecahydrate	704
480277	Sodium phosphate tribasic dodecahydrate	704
480501	Sodium hydroxide, pellets	679
480502	Sodium hydroxide, pellets	679
480505	Sodium hydroxide, pellets	679
480507	Sodium hydroxide, pellets	679
480508	Sodium hydroxide, pellets	679
480509	Sodium hydroxide, pellets	679
480522	Sodium hydroxide, pellets	679
480525	Sodium hydroxide, pellets	679
480527	Sodium hydroxide, pellets	679
480561	Sodium hydroxide solution 32%	681
480562	Sodium hydroxide solution 32%	681
480563	Sodium hydroxide solution 32%	681
480564	Sodium hydroxide solution 32%	681
480566	Sodium hydroxide solution 32%	681
480591	Sodium hydroxide solution 35%	681
480593	Sodium hydroxide solution 35%	681
480621	Sodium hydroxide solution 20% w/w	683
480622	Sodium hydroxide solution 20% w/w	683
480631	Sodium hydroxide solution 20% w/w	683
480681000	Sodium hydroxide 2 mol/l (2N)	686
480682000	Sodium hydroxide 2 mol/l (2N)	686
480684000	Sodium hydroxide 2 mol/l (2N)	686
480686000	Sodium hydroxide 2 mol/l (2N)	686
480687000	Sodium hydroxide 2 mol/l (2N)	686
480711000	Sodium hydroxide 1 mol/l (1N)	687
480713000	Sodium hydroxide 1 mol/l (1N)	687
480714000	Sodium hydroxide 1 mol/l (1N)	687
480717000	Sodium hydroxide 1 mol/l (1N)	687
480741	Sodium hydroxide 1 mol/l (1N)	687
480771000	Sodium hydroxide 0.5 mol/l (0.5N)	688
480772000	Sodium hydroxide 0.5 mol/l (0.5N)	688
480773000	Sodium hydroxide 0.5 mol/l (0.5N)	688
480777000	Sodium hydroxide 0.5 mol/l (0.5N)	688
480801	Sodium hydroxide 0.5 mol/l (0.5N)	688
480837000	Sodium hydroxide 0.357 mol/l (0.357N)	688
480861000	Sodium hydroxide 0.25 mol/l (0.25N)	688
480862000	Sodium hydroxide 0.25 mol/l (0.25N)	688
480863000	Sodium hydroxide 0.25 mol/l (0.25N)	688
480867000	Sodium hydroxide 0.25 mol/l (0.25N)	688
480891000	Sodium hydroxide 0.1 mol/l (0.1N)	690
480892000	Sodium hydroxide 0.1 mol/l (0.1N)	690
480893000	Sodium hydroxide 0.1 mol/l (0.1N)	690
480895000	Sodium hydroxide 0.1 mol/l (0.1N)	690
480897000	Sodium hydroxide 0.1 mol/l (0.1N)	690
480921	Sodium hydroxide 0.1 mol/l (0.1N)	690
481001	Sodium hydroxide 0.01 mol/l (0.01N)	690
481162	Sodium iodide	693
481163	Sodium iodide	693
481164	Sodium iodide	693
481181	Sodium hypochlorite solution in water	692
481185	Sodium hypochlorite solution in water	692
481201	Sodium hypophosphite	693
481202	Sodium hypophosphite	693
481231	Sodium laurylsulfate	693
481233	Sodium laurylsulfate	693
481235	Sodium laurylsulfate	693
481283	Sodium metabisulfite	694
481286	Sodium metabisulfite	694
481287	Sodium metabisulfite	694
481288	Sodium metabisulfite	694
481552	Sodium metaphosphate	695
481557	Sodium metaphosphate	695
481684	Sodium molybdate dihydrate	695
481685	Sodium molybdate dihydrate	695
481687	Sodium molybdate dihydrate	695
481751	Sodium nitrate	696
481755	Sodium nitrate	696
481756	Sodium nitrate	696
481757	Sodium nitrate	696
481759	Sodium nitrate	696
481825	Sodium nitrite	696
481826	Sodium nitrite	696
481827	Sodium nitrite	696
481829	Sodium nitrite	696
481932	Sodium nitroprusside dihydrate	697
481934	Sodium nitroprusside dihydrate	697
482041	Sodium o-Phosphite pentahydrate	704
482042	Sodium o-Phosphite pentahydrate	704
482064	Sodium oxalate	698
482065	Sodium oxalate	698
482067	Sodium oxalate	698
482101	Sodium oxalate	698
482183	Sodium perborate tetrahydrate	698
482185	Sodium perborate tetrahydrate	698
482187	Sodium perborate tetrahydrate	698
482204	Sodium perchlorate monohydrate	699
482234	Sodium metaperiodate	694
482236	Sodium metaperiodate	694
482252	Sodium peroxide	699
482363	Sodium persulfate	699
482365	Sodium persulfate	699
482367	Sodium persulfate	699
482421	Sodium pyrophosphate decahydrate	705
482422	Sodium pyrophosphate decahydrate	705

Product specifications are subject to changes. Please visit our website for updates.

482426	Sodium pyrophosphate decahydrate	705
482427	Sodium pyrophosphate decahydrate	705
482957	Sodium sulfate decahydrate	707
482959	Sodium sulfate decahydrate	707
483001	Sodium sulfate anhydrous	706
483005	Sodium sulfate anhydrous	706
483006	Sodium sulfate anhydrous	706
483007	Sodium sulfate anhydrous	706
483009	Sodium sulfate anhydrous	706
483025	Sodium sulfate anhydrous	706
483252	Sodium sulfite anhydrous	708
483256	Sodium sulfite anhydrous	708
483257	Sodium sulfite anhydrous	708
483258	Sodium sulfite anhydrous	708
483354	Sodium sulfocyanate	709
483356	Sodium sulfocyanate	709
483484	Sodium sulfide nonahydrate	707
483485	Sodium sulfide nonahydrate	707
483487	Sodium sulfide nonahydrate	707
483489	Sodium sulfide nonahydrate	707
483551	Sodium succinate hexahydrate	706
483554	Sodium succinate hexahydrate	706
483555	Sodium succinate hexahydrate	706
483557	Sodium succinate hexahydrate	706
483561	Sodium tartrate dihydrate	709
483635	Sodium tartrate dihydrate	709
483636	Sodium tartrate dihydrate	709
483637	Sodium tartrate dihydrate	709
483703	Sodium hydrogen tartrate monohydrate	678
483706	Sodium hydrogen tartrate monohydrate	678
483735	Sodium tetraborate anhydrous	709
483736	Sodium tetraborate anhydrous	709
483751	Sodium tetraphenylborate	710
483758	Sodium tetraphenylborate	710
483821	Sodium thiosulfate pentahydrate	711
483825	Sodium thiosulfate pentahydrate	711
483826	Sodium thiosulfate pentahydrate	711
483827	Sodium thiosulfate pentahydrate	711
483829	Sodium thiosulfate pentahydrate	711
484026	Sodium thiosulfate 1 mol/l (1N)	712
484071000	Sodium thiosulfate 0.1 mol/l (0.1N)	712
484072000	Sodium thiosulfate 0.1 mol/l (0.1N)	712
484077000	Sodium thiosulfate 0.1 mol/l (0.1N)	712
484121	Sodium thiosulfate 0.1 mol/l (0.1N)	713
484141	Sodium thiosulfate 0.0394 mol/l (0.0394N)	713
484155	Sodium thiosulfate 0.0197 mol/l (0.0197N)	713
484161	Sodium thiosulfate 0.01 mol/l (0.01N)	713
484233	Sodium tungstate dihydrate	714
484236	Sodium tungstate dihydrate	714
484701	D-Sorbitol	715
484704	D-Sorbitol	715
484705	D-Sorbitol	715
484861	Tin standard solution	760
484887	Tin foil	759
484914	Tin, powder	759
484917	Tin, powder	759
485002	Tin (II) chloride dihydrate	761
485004	Tin (II) chloride dihydrate	761
485005	Tin (II) chloride dihydrate	761
485007	Tin (II) chloride dihydrate	761
485074	Tin (IV) chloride pentahydrate	762
485076	Tin (IV) chloride pentahydrate	762
485154	Tin (IV) oxide	762
485304	Strontium acetate	723
485354	Strontium bromide monohydrate	723
485391	Strontium standard solution	723
485404	Strontium carbonate	724
485407	Strontium carbonate	724
485455	Strontium chloride hexahydrate	724
485457	Strontium chloride hexahydrate	724
485605	Strontium nitrate	724
485607	Strontium nitrate	724
485705	Strontium sulfate	725
485902	Sudan III	726
485961	Sulfanilamide	728
485971	Sulfanilamide	728
486211	Buffer pH 1	216
486221	Buffer pH 1	216
486231	Buffer pH 2	217
486241	Buffer pH 2	217
486251	Buffer pH 3	217
486252	Buffer pH 3	217
486261	Buffer pH 3	217
486271	Buffer pH 4	218
486273	Buffer pH 4	218
486274	Buffer pH 4	218
486276	Buffer pH 4	218
486281	Buffer pH 4	218
486291	Buffer pH 4	219
486301	Buffer pH 5	219
486311	Buffer pH 5	219
486321	Buffer pH 6	220
486331	Buffer pH 6	220
486401	Buffer pH 6.8	220
486411	Buffer pH 7.20 Weise	221
486421	Buffer pH 7	221
486431	Buffer pH 7	221
486441	Buffer pH 7.2	221
486451	Buffer pH 7	221
486453	Buffer pH 7	221
486454	Buffer pH 7	221
486455	Buffer pH 7	221
486456	Buffer pH 7	221
486461	Buffer pH 7.4	222
486531	Buffer pH 8	222
486541	Buffer pH 8	222
486542	Buffer pH 8	222
486571	Buffer pH 9	223
486581	Buffer pH 10.06	224
486591	Buffer pH 9	223
486593	Buffer pH 9	223
486594	Buffer pH 9	223
486601	Buffer pH 10	224
486611	Buffer pH 10	224
486613	Buffer pH 10	224
486614	Buffer pH 10	224
486615	Buffer pH 10	224
486621	Buffer pH 12	225
486631	Buffer pH 11	224
486641	Buffer pH 13	225
486691	Buffer pH 12	225
486701	Buffer pH 13	225
486702	Buffer pH 13	225
486741	Buffer pH 3.56	218
486751	Buffer pH 1.68	216
486761	Buffer pH 4	218
486762	Buffer pH 4	218
486771	Buffer pH 11	224
486772	Buffer pH 11	224
486791	Buffer pH 7	221
486792	Buffer pH 7	221
486811	Buffer pH 10.06	224
486812	Buffer pH 10.06	224
486841	Buffer pH 4.62	219
486871	Buffer pH 6.88	220
486881	Buffer pH 9.22	223
486903	Tartrazine	746
487002	Tellurium lumps	746
487023	Tellurium powder	747
487051	Tetrabutylammonium bromide	749
487101	Tetrabutylammonium bisulfate	748
487152	Tetraethylammonium bromide	751
487301	Tetrahydrofuran	752
487303	Tetrahydrofuran	752
487305	Tetrahydrofuran	752
487307	Tetrahydrofuran	752
487308	Tetrahydrofuran	752
487309	Tetrahydrofuran	752
487345	Tetrahydrofuran	752
487346	Tetrahydrofuran	752
487352	Tetrahydrofuran	752
487491	Tetramethylammonium hydroxide 10%	753
487492	Tetramethylammonium hydroxide 10%	753
487601	N,N,N',N'-Tetramethyl-p-phenylenediamine dihydrochloride	753
487728	Thymolphthalein	758
487729	Thymolphthalein	758
487803	Thioacetamide	754
488101	Thiourea	756
488102	Thiourea	756
488104	Thiourea	756
488105	Thiourea	756
488107	Thiourea	756
488152	L(-)Tyrosine	781
488162	Tisab III solution	763
488251	Titanium dioxide	764
488256	Titanium dioxide	764
488257	Titanium dioxide	764
488421	Titanium isopropylate	764
488461	o-Tolidine solution 0.1%	765
488531	Toluene	766
488551	Toluene	767
488552	Toluene	767
488555	Toluene	767
488556	Toluene	767
488557	Toluene	767
488591	Toluene	766
488592	Toluene	766
488594	Toluene	766
488601	Toluene	766
488602	Toluene	766
488661	p-Toluene sulfonamide	768
489054	Litmus	465
489152	Triacetin	769
489162	Turbidity std 4000NTU formazine	780
489382	1,2,4-Trichlorobenzene	771
489501	Triethanolamine	771
489504	Triethanolamine	771
489556	Triethylamine	772
489561	N,O-Bis(trimethylsilyl)-trifluoroacetamide	204
489581	Trioctylphosphine oxide	775
489591	Triphenylphosphine	776
489631	Triethylamine	772
489633	Triethylamine	772
489651	2,3,5-Triphenyltetrazolium chloride	776
489831	Trimethylcetylammonium bromide	775
489833	Trimethylcetylammonium bromide	775

INDEX BY PRODUCT CODE

489934	N,O-Bis(trimethylsilyl) acetamide	204
489971	Tris (hydroxymethyl)-aminomethane	776
489973	Tris (hydroxymethyl)-aminomethane	776
489981	Tris (hydroxymethyl)-aminomethane	776
489983	Tris (hydroxymethyl)-aminomethane	776
489984	Tris (hydroxymethyl)-aminomethane	776
489985	Tris (hydroxymethyl)-aminomethane	776
490001	Tropaeolin O	778
490002	Tropaeolin O	778
490422	Tauber reagent	746
490751	Urea	782
490757	Urea	782
490758	Urea	782
490759	Urea	782
491091	Vanadium standard solution	783
491103	Vanadium (V) oxide	784
491152	Brilliant green	209
491207	Bromocresol green	211
491208	Bromocresol green	211
491303	Malachite green	476
491304	Malachite green	476
491351	Methyl green	499
491352	Methyl green	499
491371	Light green	461
491372	Light green	461
491391	Fast green FCF	353
491502	Crystal violet	284
491561	Crystal violet oxalate for Gram-Hucker Kit	284
491871	Pyrocatechol violet	632
491872	Pyrocatechol violet	632
492011	Wright's stain solution in methanol	790
492211	Xylenecyanol	792
492212	Xylenecyanol	792
492301	Xylene, mix of isomers	791
492303	Xylene, mix of isomers	791
492304	Xylene, mix of isomers	791
492305	Xylene, mix of isomers	791
492306	Xylene, mix of isomers	791
492358	Xylene, mix of isomers	791
492359	Xylene, mix of isomers	791
492401	o-Xylene	792
492403	o-Xylene	792
492404	o-Xylene	792
492661	2,4-Xylenol	792
492803	D(+)Xylose	793
492804	D(+)Xylose	793
493101	Ziehl-Neelsen's reagent	796
493102	Ziehl-Neelsen's reagent	796
493151	Zinc standard solution	798
493303	Zinc, granular	796
493307	Zinc, granular	796
493309	Zinc, granular	796
493451	Zinc, granular	796
493507	Zinc, foil	796
493702	Zinc, powder	797
493705	Zinc, powder	797
493707	Zinc, powder	797
493803	Zinc acetate dihydrate	798
493806	Zinc acetate dihydrate	798
493807	Zinc acetate dihydrate	798
494006	Zinc carbonate basic	799
494104	Zinc chloride anhydrous	799
494105	Zinc chloride anhydrous	799
494106	Zinc chloride anhydrous	799
494107	Zinc chloride anhydrous	799
494311	Zinc dibenzyldithiocarbamate	800
494506	Zinc nitrate hexahydrate	800
494507	Zinc nitrate hexahydrate	800
494602	Zinc oxide	801
494606	Zinc oxide	801
494607	Zinc oxide	801
494901	Zinc sulfate heptahydrate	802
494905	Zinc sulfate heptahydrate	802
494906	Zinc sulfate heptahydrate	802
494907	Zinc sulfate heptahydrate	802
494909	Zinc sulfate heptahydrate	802
494921	Zinc sulfate 0.1 mol/l (0.2N)	803
494931	Zinc sulfate 0.05 mol/l (0.05N)	803
495005	Zinc sulfate monohydrate	802
495007	Zinc sulfate monohydrate	802
495105	Zinc sulfide	803
495107	Zinc sulfide	803
495202	Zirconium powder	803
495305	Zirconium (IV) oxide	804
497551	Mercury standard solution	484
497555	Mercury standard solution	484
502000	Acetic acid 1 mol/l (1N)	136
502002	Boric acid 4%	206
502010	Hydrochloric acid 4 mol/l (4N)	402
502011	Hydrochloric acid 3 mol/l (3N)	403
502020	Sulfuric acid d=1.820	743
502044	Hydrogen peroxide solution 30%	414
502050	Amidoschwarz 10B solution	157
502051	Amidoschwarz 10B solution	157
502062	Sand of Fontainebleau	644
502063	Sand of Fontainebleau	644
502064	Sand of Fontainebleau	644
502073	Trichloroacetic acid solution 20%	770
502092	Potassium hydroxide 0.23 mol/l (0.23N)	604
502100000	Sulfuric acid 0.1 mol/l (0.2N)	740
502112	Sodium hydroxide solution 35-37%	681
502120	Kjeldahl selenium catalyst	450
502121	Kjeldahl catalyst for water analysis	450
502122	Kjeldahl catalyst for water analysis	450
502123	Kjeldahl titanium catalyst	451
502131	Sodium chloride 5 mol/l (5N)	672
502202	Sulfuric acid 0.26 mol/l (0.52N)	739
502212	Potassium hydroxide 0.46 mol/l (0.46N)	604
502302	Sulfuric acid 95 - 97 %	733
502591	Sulfuric acid 10% v/v	737
502601	Boric acid 4% with indicator	206
502611	Boric acid 1% with indicator	207
502612	Boric acid 1% with indicator	207
502621	Hydrochloric acid 3 mol/l (3N)	403
502622	Hydrochloric acid 3 mol/l (3N)	403
502631	Hydrochloric acid 0.2 mol/l (0.2N)	405
502641	Sulfuric acid 98%	730
502651	Sulfuric acid 0.13 mol/l (0.26N)	740
502662	Sodium hydroxide 4 mol/l (4N)	685
502664	Sodium hydroxide 4 mol/l (4N)	685
502671	Mixture for checking solderings	509
502681	Potassium chromate 5% solution	591
502701	Carrez reagent zinc salt	249
502711	Carrez reagent potassium salt	249
502721	Sodium hydroxide solution 40%	680
502722	Sodium hydroxide solution 40%	680
502731	Sodium hydroxide solution 30%	682
502741	Sodium hydroxide solution 30%	682
502761	Hydrochloric acid 1.128% m/v	401
502771	Sulfuric acid 72%	735
502781000	Sodium hydroxide 0.2 mol/l (0.2N)	689
502782	Sodium hydroxide 0.2 mol/l (0.2N)	689
502791	Kjeldahl catalyst without selenium and titanium	450
502792	Kjeldahl catalyst without selenium and titanium	450
502802	Kjeldahl titanium catalyst	451
502811	Kjeldahl antifoam	449
502821	Kjeldahl catalyst according to Wieninger	450
502831	Hydrochloric acid 6 mol/l (6N)	401
502832	Hydrochloric acid 6 mol/l (6N)	401
503171	Bromate standard solution	210
503173	Bromate standard solution	210
503181	Chlorate standard solution	257
503183	Chlorate standard solution	257
503191	Chlorite standard solution	257
503193	Chlorite standard solution	257
503194	Ammonium di-hydrogen phosphate 25 mg/L solution	172
503195	Ammonium nitrate 200 mg/l solution	174
503196	Magnesium nitrate 10 g/l solution	473
503197	Nickel (II) nitrate 10g/l	524
503198	Palladium nitrate 2 g/l solution	554
503202	Palladium nitrate 2 g/l solution	554
503211	Bromide standard solution	210
503213	Bromide standard solution	210
503221	Calcium standard solution	235
503223	Calcium standard solution	235
503231	Chloride standard solution	257
503233	Chloride standard solution	257
503241	Chromate standard solution	264
503243	Chromate standard solution	264
503251	Fluoride standard solution	358
503253	Fluoride standard solution	358
503261	Iodide standard solution	424
503263	Iodide standard solution	424
503271	Potassium standard solution	583
503273	Potassium standard solution	583
503281	Lithium standard solution	462
503283	Lithium standard solution	462
503291	Magnesium standard solution	469
503293	Magnesium standard solution	469
503301	Sodium standard solution	661
503303	Sodium standard solution	661
503311	Ammonium standard solution	165
503313	Ammonium standard solution	165
503321	Nitrite standard solution	534
503323	Nitrite standard solution	534
503331	Nitrate standard solution	527
503333	Nitrate standard solution	527
503341	Phosphate standard solution	574
503343	Phosphate standard solution	574
503351	Sulfate standard solution	729
503353	Sulfate standard solution	729
503358	Cyanide standard solution	285
503361	Strontium standard solution	723
503389	Calcium standard solution	235
503390	Magnesium standard solution	469
503401	Silver standard solution	653
503403	Silver standard solution	653
503405	Silver standard solution	653
503407	Silver standard solution	653
503411	Aluminum standard solution	152
503413	Aluminum standard solution	152
503415	Aluminum standard solution	152
503417	Aluminum standard solution	152
503421	Arsenic standard solution	185
503423	Arsenic standard solution	185
503425	Arsenic standard solution	185

Product specifications are subject to changes. Please visit our website for updates.

INDEX BY PRODUCT CODE

503427	Arsenic standard solution	185
503431	Gold standard solution	375
503433	Gold standard solution	375
503435	Gold standard solution	375
503437	Gold standard solution	375
503441	Boron standard solution	208
503443	Boron standard solution	208
503445	Boron standard solution	208
503447	Boron standard solution	208
503451	Barium standard solution	191
503453	Barium standard solution	191
503455	Barium standard solution	191
503457	Barium standard solution	191
503461	Beryllium standard solution	201
503463	Beryllium standard solution	201
503471	Bismuth standard solution	202
503473	Bismuth standard solution	202
503475	Bismuth standard solution	202
503477	Bismuth standard solution	202
503481	Calcium standard solution	235
503483	Calcium standard solution	235
503485	Calcium standard solution	235
503487	Calcium standard solution	235
503491	Cadmium standard solution	232
503493	Cadmium standard solution	232
503495	Cadmium standard solution	232
503497	Cadmium standard solution	232
503501	Cerium standard solution	250
503503	Cerium standard solution	250
503505	Cerium standard solution	250
503507	Cerium standard solution	250
503511	Cobalt standard solution	271
503513	Cobalt standard solution	271
503515	Cobalt standard solution	271
503517	Cobalt standard solution	271
503521	Chromium standard solution	265
503523	Chromium standard solution	265
503525	Chromium standard solution	265
503527	Chromium standard solution	265
503531	Cesium standard solution	253
503533	Cesium standard solution	253
503535	Cesium standard solution	253
503537	Cesium standard solution	253
503541	Copper standard solution	276
503543	Copper standard solution	276
503545	Copper standard solution	276
503547	Copper standard solution	276
503571	Europium standard solution	351
503573	Europium standard solution	351
503575	Europium standard solution	351
503577	Europium standard solution	351
503581	Iron standard solution	428
503583	Iron standard solution	428
503585	Iron standard solution	428
503587	Iron standard solution	428
503601	Gadolinium standard solution	366
503603	Gadolinium standard solution	366
503605	Gadolinium standard solution	366
503607	Gadolinium standard solution	366
503611	Lutetium standard solution	467
503613	Lutetium standard solution	467
503615	Lutetium standard solution	467
503617	Lutetium standard solution	467
503631	Mercury standard solution	484
503633	Mercury standard solution	484
503635	Mercury standard solution	484
503637	Mercury standard solution	484
503640	Mercury standard solution	484
503651	Indium standard solution	423
503653	Indium standard solution	423
503655	Indium standard solution	423
503657	Indium standard solution	423
503671	Potassium standard solution	582
503673	Potassium standard solution	582
503675	Potassium standard solution	582
503677	Potassium standard solution	582
503681	Lanthanum standard solution	453
503683	Lanthanum standard solution	453
503685	Lanthanum standard solution	453
503687	Lanthanum standard solution	453
503691	Lithium standard solution	462
503693	Lithium standard solution	462
503695	Lithium standard solution	462
503697	Lithium standard solution	462
503711	Magnesium standard solution	469
503713	Magnesium standard solution	469
503715	Magnesium standard solution	469
503717	Magnesium standard solution	469
503718	Magnesium standard solution	469
503719	Magnesium standard solution	469
503721	Manganese standard solution	479
503723	Manganese standard solution	479
503725	Manganese standard solution	479
503727	Manganese standard solution	479
503731	Molybdenum standard solution	510
503733	Molybdenum standard solution	510
503735	Molybdenum standard solution	510
503737	Molybdenum standard solution	510
503741	Sodium standard solution	661
503743	Sodium standard solution	661
503745	Sodium standard solution	661
503747	Sodium standard solution	661
503749	Sodium standard solution	661
503751	Niobium standard solution	527
503753	Niobium standard solution	527
503755	Niobium standard solution	527
503757	Niobium standard solution	527
503761	Neodymium standard solution	520
503763	Neodymium standard solution	520
503765	Neodymium standard solution	520
503767	Neodymium standard solution	520
503771	Nickel standard solution	522
503773	Nickel standard solution	522
503775	Nickel standard solution	522
503777	Nickel standard solution	522
503791	Phosphorus standard solution	577
503793	Phosphorus standard solution	577
503795	Phosphorus standard solution	577
503797	Phosphorus standard solution	577
503801	Lead standard solution	456
503803	Lead standard solution	456
503805	Lead standard solution	456
503807	Lead standard solution	456
503811	Palladium standard solution	553
503813	Palladium standard solution	553
503815	Palladium standard solution	553
503817	Palladium standard solution	553
503821	Praseodymium standard solution	621
503823	Praseodymium standard solution	621
503825	Praseodymium standard solution	621
503827	Praseodymium standard solution	621
503831	Platinum standard solution	581
503833	Platinum standard solution	581
503835	Platinum standard solution	581
503837	Platinum standard solution	581
503841	Rubidium standard solution	640
503843	Rubidium standard solution	640
503845	Rubidium standard solution	640
503847	Rubidium standard solution	640
503861	Rhodium standard solution	639
503863	Rhodium standard solution	639
503865	Rhodium standard solution	639
503871	Ruthenium standard solution	641
503873	Ruthenium standard solution	641
503875	Ruthenium standard solution	641
503877	Ruthenium standard solution	641
503891	Antimony standard solution	182
503893	Antimony standard solution	182
503895	Antimony standard solution	182
503897	Antimony standard solution	182
503898	Antimony standard solution	182
503899	Antimony standard solution	182
503901	Scandium standard solution	645
503903	Scandium standard solution	645
503905	Scandium standard solution	645
503907	Scandium standard solution	645
503911	Selenium standard solution	647
503913	Selenium standard solution	647
503915	Selenium standard solution	647
503917	Selenium standard solution	647
503921	Silicon standard solution	651
503923	Silicon standard solution	651
503925	Silicon standard solution	651
503927	Silicon standard solution	651
503931	Samarium standard solution	643
503933	Samarium standard solution	643
503935	Samarium standard solution	643
503937	Samarium standard solution	643
503941	Tin standard solution	760
503943	Tin standard solution	760
503945	Tin standard solution	760
503947	Tin standard solution	760
503949	Tin standard solution	760
503951	Strontium standard solution	722
503953	Strontium standard solution	722
503955	Strontium standard solution	722
503957	Strontium standard solution	722
503961	Tantalum standard solution	745
503963	Tantalum standard solution	745
503965	Tantalum standard solution	745
503967	Tantalum standard solution	745
503971	Terbium standard solution	748
503973	Terbium standard solution	748
503975	Terbium standard solution	748
503977	Terbium standard solution	748
503981	Tellurium standard solution	747
503983	Tellurium standard solution	747
503985	Tellurium standard solution	747
503987	Tellurium standard solution	747
504001	Titanium standard solution	763
504003	Titanium standard solution	763
504005	Titanium standard solution	763
504007	Titanium standard solution	763
504011	Thallium standard solution	754
504013	Thallium standard solution	754
504015	Thallium standard solution	754
504017	Thallium standard solution	754
504031	Uranium standard solution	782
504033	Uranium standard solution	782
504035	Uranium standard solution	782
504037	Uranium standard solution	782
504041	Vanadium standard solution	783
504043	Vanadium standard solution	783
504045	Vanadium standard solution	783
504047	Vanadium standard solution	783
504051	Tungsten standard solution	779
504053	Tungsten standard solution	779
504055	Tungsten standard solution	779
504057	Tungsten standard solution	779
504058	Tungsten standard solution	779
504061	Yttrium standard solution	795
504063	Yttrium standard solution	795
504065	Yttrium standard solution	795
504067	Yttrium standard solution	795
504071	Ytterbium standard solution	794
504073	Ytterbium standard solution	794
504075	Ytterbium standard solution	794
504077	Ytterbium standard solution	794

INDEX BY PRODUCT CODE

504081	Zinc standard solution	798
504083	Zinc standard solution	798
504085	Zinc standard solution	798
504087	Zinc standard solution	798
504091	Zirconium standard solution	804
504093	Zirconium standard solution	804
504095	Zirconium standard solution	804
504097	Zirconium standard solution	804
504186	Aluminum standard solution	152
504187	Vanadium standard solution	783
504190	Aluminum standard solution	152
504194	Iron standard solution	429
504195	Chromium standard solution	265
504221	Hafnium standard solution	378
504223	Hafnium standard solution	378
504225	Hafnium standard solution	378
504227	Hafnium standard solution	378
504231	Dysprosium standard solution	325
504233	Dysprosium standard solution	325
504235	Dysprosium standard solution	325
504237	Dysprosium standard solution	325
504241	Erbium standard solution	329
504243	Erbium standard solution	329
504245	Erbium standard solution	329
504247	Erbium standard solution	329
504251	Germanium standard solution	369
504253	Germanium standard solution	369
504255	Germanium standard solution	369
504257	Germanium standard solution	369
504261	Holmium standard solution	390
504263	Holmium standard solution	390
504265	Holmium standard solution	390
504267	Holmium standard solution	390
504271	Silicon standard solution	651
504273	Silicon standard solution	651
504275	Silicon standard solution	651
504277	Silicon standard solution	651
504281	Thorium standard solution	756
504283	Thorium standard solution	756
504291	Sulfur standard solution	730
504293	Sulfur standard solution	730
504295	Sulfur standard solution	730
504297	Sulfur standard solution	730
504301	Multielement standard for ICP	514
504303	Multielement standard for ICP	514
504305	Multielement standard for ICP	514
504306	Multielement standard for ICP	515
504307	Multielement standard for ICP	515
504308	Multielement standard for ICP	515
504309	Multielement standard for ICP	515
504310	Multielement standard for ICP	515
504311	Multielement standard for ICP	515
504312	Multielement standard for ICP	515
504313	Multielement standard for ICP	515
504350	Multielement standard for ICP	514
504351	Multielement standard for ICP	514
504352	Multielement standard for ICP and ICP-MS	515
504353	Multielement standard for ICP	514
504354	Multielement standard for ICP	514
504355	Multielement standard for ICP	514
504356	Multielement standard for ICP	514
504357	Multielement standard for ICP	514
504360	Cadmium standard solution	232
504361	Copper standard solution	276
504362	Manganese standard solution	480
504363	Nickel standard solution	522
504364	Lead standard solution	456
504370	Mercury standard solution	485
504392	Multielement standard for ICP and ICP-MS	515
504393	Multielement standard for ICP and ICP-MS	515
504396	Multielement standard for ICP	513
504439	Arsenic standard solution	185
504480	Multielement standard for ICP	514
504526	Multianions standard for ion chromatography	512
504527	Multianions standard for ion chromatography	512
504530	Eluent sodium carbonate/ sodium bicarbonate	328
504531	Eluent sodium carbonate/ sodium bicarbonate	328
504532	Eluent sodium carbonate/ sodium bicarbonate	328
504533	Eluent sodium carbonate	327
504534	Eluent sodium bicarbonate	327
504536	Cesium chloride 25 g/l solution	254
504537	Lanthanum chloride 25 g/l solution	454
504538	Potassium chloride 25g/l in HCl	589
504545	Copper standard solution	276
504550	Water deionized and acidified	789
504551	Water deionized and acidified	789
504552	Water deionized and acidified	789
504553	Water deionized and acidified	789
504554	Water deionized and acidified	789
504557	Water deionized and acidified	789
504562	Sulfuric acid 25%	736
504571	Hydrochloric acid 50% v/v	393
504594	Potassium dichromate 0.1414 g/l	594
504677	Multianions standard for ion chromatography	512
504901	Lead standard solution	455
504902	Lead standard solution	455
505008	Calcium carbonate	237
505009	Calcium carbonate	237
505021	Nital solution 4%	527
505032	Potassium chromate solution 10%	590
505302	Silver standard solution	653
505303	Silver standard solution	653
505305	Silver standard solution	653
505307	Aluminum standard solution	151
505308	Aluminum standard solution	151
505309	Aluminum standard solution	151
505312	Arsenic standard solution	185
505313	Arsenic standard solution	185
505315	Arsenic standard solution	185
505317	Gold standard solution	375
505318	Gold standard solution	375
505319	Gold standard solution	375
505322	Boron standard solution	208
505323	Boron standard solution	208
505325	Boron standard solution	208
505327	Barium standard solution	191
505328	Barium standard solution	191
505329	Barium standard solution	191
505332	Beryllium standard solution	201
505333	Beryllium standard solution	201
505335	Beryllium standard solution	201
505337	Bismuth standard solution	202
505338	Bismuth standard solution	202
505339	Bismuth standard solution	202
505542	Calcium standard solution	235
505543	Calcium standard solution	235
505545	Calcium standard solution	235
505547	Cadmium standard solution	232
505548	Cadmium standard solution	232
505549	Cadmium standard solution	232
505552	Cerium standard solution	250
505555	Cerium standard solution	250
505562	Cobalt standard solution	270
505563	Cobalt standard solution	270
505565	Cobalt standard solution	270
505567	Chromium standard solution	264
505568	Chromium standard solution	264
505569	Chromium standard solution	264
505572	Cesium standard solution	253
505575	Cesium standard solution	253
505577	Copper standard solution	275
505578	Copper standard solution	275
505579	Copper standard solution	275
505582	Dysprosium standard solution	325
505585	Dysprosium standard solution	325
505592	Erbium standard solution	329
505595	Erbium standard solution	329
505602	Europium standard solution	351
505605	Europium standard solution	351
505612	Iron standard solution	428
505613	Iron standard solution	428
505615	Iron standard solution	428
505617	Gallium standard solution	367
505618	Gallium standard solution	367
505619	Gallium standard solution	367
505622	Gadolinium standard solution	366
505625	Gadolinium standard solution	366
505632	Germanium standard solution	369
505633	Germanium standard solution	369
505635	Germanium standard solution	369
505642	Hafnium standard solution	378
505645	Hafnium standard solution	378
505652	Mercury standard solution	484
505654	Mercury standard solution	484
505655	Mercury standard solution	484
505657	Holmium standard solution	390
505658	Holmium standard solution	390
505662	Indium standard solution	422
505663	Indium standard solution	422
505665	Indium standard solution	422
505675	Iridium standard solution	427
505682	Potassium standard solution	582
505683	Potassium standard solution	582
505685	Potassium standard solution	582
505692	Lanthanum standard solution	453
505693	Lanthanum standard solution	453
505695	Lanthanum standard solution	453
505702	Lithium standard solution	462
505703	Lithium standard solution	462
505705	Lithium standard solution	462
505707	Lutetium standard solution	467
505708	Lutetium standard solution	467
505709	Lutetium standard solution	467
505712	Magnesium standard solution	469
505713	Magnesium standard solution	469
505715	Magnesium standard solution	469
505717	Manganese standard solution	479
505718	Manganese standard solution	479
505719	Manganese standard solution	479
505722	Molybdenum standard solution	510
505723	Molybdenum standard solution	510
505725	Molybdenum standard solution	510
505732	Sodium standard solution	661
505733	Sodium standard solution	661
505735	Sodium standard solution	661
505737	Niobium standard solution	526
505738	Niobium standard solution	526
505742	Neodymium standard solution	519
505745	Neodymium standard solution	519
505752	Nickel standard solution	522
505753	Nickel standard solution	522
505755	Nickel standard solution	522
505758	Osmium standard solution	549
505762	Phosphorus standard solution	577
505763	Phosphorus standard solution	577
505765	Phosphorus standard solution	577
505767	Lead standard solution	456
505768	Lead standard solution	456
505769	Lead standard solution	456
505772	Palladium standard solution	553

505775	Palladium standard solution	553
505782	Praseodymium standard solution	621
505785	Praseodymium standard solution	621
505787	Platinum standard solution	580
505788	Platinum standard solution	580
505789	Platinum standard solution	580
505792	Rubidium standard solution	640
505795	Rubidium standard solution	640
505802	Rhenium standard solution	638
505805	Rhenium standard solution	638
505807	Rhodium standard solution	638
505808	Rhodium standard solution	638
505809	Rhodium standard solution	638
505812	Ruthenium standard solution	641
505815	Ruthenium standard solution	641
505822	Sulfur standard solution	730
505823	Sulfur standard solution	730
505825	Sulfur standard solution	730
505832	Antimony standard solution	182
505833	Antimony standard solution	182
505835	Antimony standard solution	182
505837	Scandium standard solution	645
505838	Scandium standard solution	645
505839	Scandium standard solution	645
505842	Selenium standard solution	647
505843	Selenium standard solution	647
505845	Selenium standard solution	647
505847	Silicon standard solution	651
505848	Silicon standard solution	651
505852	Samarium standard solution	643
505855	Samarium standard solution	643
505862	Tin standard solution	760
505863	Tin standard solution	760
505865	Tin standard solution	760
505867	Strontium standard solution	722
505868	Strontium standard solution	722
505869	Strontium standard solution	722
505872	Tantalum standard solution	744
505875	Tantalum standard solution	744
505882	Terbium standard solution	748
505883	Terbium standard solution	748
505885	Terbium standard solution	748
505887	Tellurium standard solution	747
505888	Tellurium standard solution	747
505907	Titanium standard solution	763
505908	Titanium standard solution	763
505909	Titanium standard solution	763
505912	Thallium standard solution	754
505913	Thallium standard solution	754
505915	Thallium standard solution	754
505917	Thulium standard solution	756
505918	Thulium standard solution	756
505922	Uranium standard solution	782
505923	Uranium standard solution	782
505927	Vanadium standard solution	783
505928	Vanadium standard solution	783
505929	Vanadium standard solution	783
505932	Tungsten standard solution	779
505935	Tungsten standard solution	779
505942	Yttrium standard solution	794
505943	Yttrium standard solution	794
505945	Yttrium standard solution	794
505947	Ytterbium standard solution	794
505948	Ytterbium standard solution	794
505952	Zinc standard solution	797
505953	Zinc standard solution	797
505955	Zinc standard solution	797
505957	Zirconium standard solution	804
505958	Zirconium standard solution	804
506002	Organic standard: Mixture for hydrocarbon analysis	547
506010	Organic standard: Mixture for hydrocarbon analysis	547
506011	Organic standard: Mixture for hydrocarbon analysis	547
506012	Organic standard: Mixture for hydrocarbon analysis	547
506013	Organic standard: Mixture for hydrocarbon analysis	547
506020	Organic standard: Mixture for hydrocarbon analysis	547
506021	Organic standard: Mixture for hydrocarbon analysis	547
506030	Organic standard: Mixture for hydrocarbon analysis	547
506040	Organic standard: Mixture for hydrocarbon analysis	547
506110	Multielement standard for ICP	513
506120	Multielement standard for ICP	513
506130	Multielement standard for ICP	513
506140	Multielement standard for ICP	513
506150	Multielement standard for ICP	513
506411	Water	787
506432	Potassium chloride 12g/l	589
506433	Potassium chloride 12g/l	589
506442	Potassium dichromate - Sulfuric acid solution	592
506443	Potassium dichromate - Sulfuric acid solution	592
506452	Potassium dichromate - Sulfuric acid solution	592
506453	Potassium dichromate - Sulfuric acid solution	592
506462	Toluene in solution in hexane	768
506463	Toluene in solution in hexane	768
506473	Holmium perchlorate in solution	391
506614	Organic standard: Mixture for hydrocarbon analysis	546
506732	Organic Standard: PCB multielement mixture	543
506736	Organic standard: Mixture for hydrocarbon analysis	546
506742	Organic standard: Mixture for hydrocarbon analysis	546
506803	Organic standard: Pesticide mixture	546
506821	Organic Standard: PAH multielement mixture	542
506835	Organic Standard: PAH multielement mixture	542
506878	Organic Standard: PAH multielement mixture	542
506897	Organic standard: Pesticide mixture	546
506905	Organic standard: Pesticide mixture	546
506918	Mercury standard solution	484
506938	Organic Standard: PAH multielement mixture	542
506941	Beryllium standard solution	201
506948	Organic standard: Pesticide mixture	546
506950	Organic standard: Pesticide mixture	546
506960	Potassium standard solution	582
506979	Organic Standard: PAH multielement mixture	542
506980	Organic Standard: PAH multielement mixture	542
507039	Magnesium standard solution	469
507062	Organic Standard: PCB multielement mixture	543
507063	Organic Standard: PAH multielement mixture	542
507094	Organic Standard: PAH multielement mixture	542
507103	Organic Standard: PCB multielement mixture	543
507115	Organic Standard: PCB multielement mixture	543
507127	Organic standard PCB	540
507128	Organic standard PCB	540
507129	Organic standard PCB	540
507131	Organic standard PCB	540
507132	Organic standard PCB	540
507133	Organic standard PCB	540
507134	Organic standard PCB	540
507135	Organic standard PCB	540
507136	Organic standard PCB	540
507137	Organic standard PCB	540
507138	Organic standard PCB	540
507139	Organic standard PCB	540
507154	Organic standard PCB	540
507189	Organic standard: Mixture for hydrocarbon analysis	546
507190	Organic standard: Mixture for hydrocarbon analysis	546
507191	Organic standard: Mixture for hydrocarbon analysis	546
507393	Iron standard solution	429
507410	Hydrofluoric acid 0.1 mol/l	412
507474	Organic standard: Mixture for hydrocarbon analysis	546
507476	Calcium standard solution	235
507477	Zinc standard solution	798
507478	Copper standard solution	276
507479	Antimony standard solution	183
507480	Silver standard solution	653
507481	Barium standard solution	191
507482	Bismuth standard solution	202
507483	Cadmium standard solution	232
507484	Cobalt standard solution	271
507485	Chromium standard solution	265
507486	Lithium standard solution	462
507487	Nickel standard solution	522
507488	Manganese standard solution	479
507489	Mercury standard solution	484
507490	Lead standard solution	456
507491	Selenium standard solution	647
507492	Tin standard solution	760
507493	Strontium standard solution	723
507494	Molybdenum standard solution	511
507496	Arsenic standard solution	185
507497	Beryllium standard solution	201
507498	Cerium standard solution	250
507499	Cesium standard solution	253
507500	Dysprosium standard solution	326
507501	Erbium standard solution	329
507502	Europium standard solution	352
507503	Gallium standard solution	367
507504	Gadolinium standard solution	366
507505	Germanium standard solution	369
507506	Hafnium standard solution	378
507507	Holmium standard solution	391
507508	Indium standard solution	423
507509	Lanthanum standard solution	454
507510	Neodymium standard solution	520
507511	Niobium standard solution	527
507512	Palladium standard solution	553
507513	Rhenium standard solution	638
507514	Rubidium standard solution	641
507515	Samarium standard solution	643
507516	Scandium standard solution	645
507517	Tantalum standard solution	745
507518	Tellurium standard solution	747
507519	Thulium standard solution	757

INDEX BY PRODUCT CODE

507520	Titanium standard solution	763
507521	Tungsten standard solution	779
507522	Yttrium standard solution	795
507523	Ytterbium standard solution	794
507524	Zirconium standard solution	804
507525	Antimony standard solution	183
507526	Silver standard solution	653
507527	Barium standard solution	191
507528	Bismuth standard solution	202
507529	Cadmium standard solution	232
507530	Calcium standard solution	235
507531	Cerium standard solution	250
507532	Cesium standard solution	253
507533	Cobalt standard solution	271
507536	Potassium dichromate 0.167 mol/l (1N)	592
507577	Eluent sodium carbonate	327
507578	Eluent sodium bicarbonate	327
507609	Organic Standard: PCB multielement mixture	543
507679	Organic Standard: PCB multielement mixture	543
507688	Mixtures for residual solvents analysis	509
507689	Mixtures for residual solvents analysis	509
507690	Mixtures for residual solvents analysis	509
507691	Mixtures for residual solvents analysis	509
507692	Mixtures for residual solvents analysis	509
507693	Mixtures for residual solvents analysis	509
507694	Mixtures for residual solvents analysis	509
507734	Dysprosium standard solution	326
507735	Erbium standard solution	329
507736	Europium standard solution	352
507737	Gadolinium standard solution	366
507739	Gallium standard solution	367
507740	Germanium standard solution	369
507741	Hafnium standard solution	378
507742	Holmium standard solution	391
507743	Indium standard solution	423
507744	Lanthanum standard solution	454
507745	Lithium standard solution	462
507746	Manganese standard solution	479
507747	Molybdenum standard solution	511
507748	Neodymium standard solution	520
507749	Nickel standard solution	522
507750	Niobium standard solution	527
507751	Palladium standard solution	553
507752	Lead standard solution	456
507753	Potassium standard solution	582
507754	Rhenium standard solution	638
507755	Rubidium standard solution	641
507756	Samarium standard solution	643
507757	Scandium standard solution	645
507758	Selenium standard solution	647
507759	Sodium standard solution	661
507760	Strontium standard solution	723
507761	Tantalum standard solution	745
507762	Tellurium standard solution	747
507763	Thulium standard solution	757
507764	Titanium standard solution	763
507765	Tungsten standard solution	779
507766	Vanadium standard solution	783
507767	Yttrium standard solution	795
507768	Ytterbium standard solution	794
507769	Zinc standard solution	798
507770	Zirconium standard solution	804
507859	Organic Standard: PAH multielement mixture	542
507889	Organic Standard: PCB multielement mixture	543
507899	Organic Standard: PAH multielement mixture	542
508001	Dimethylsulphoxide	319
508002	Dimethylsulphoxide	319
508111	Dimethylsulphoxide	318
508200	Acetone	141
508201	Acetone	141
508212	n-Heptane	382
508215	n-Heptane	382
508216	n-Heptane	382
508217	n-Heptane	382
508221	Ethyl acetate	343
508222	Ethyl acetate	343
508230	Petroleum ether 100 - 140°C	564
508232	Petroleum ether 100 - 140°C	564
508235	Cyclohexane	287
508320	Chloroform	262
508321	Chloroform	262
508322	Chloroform	262
508370	Dichloromethane	304
508374	Dichloromethane	304
508375	Dichloromethane	304
508615	Sodium hydroxide solution 10% w/v	683
508645	Acetic acid 27%	135
508801	N,N-Dimethylformamide	316
508802	N,N-Dimethylformamide	316
508803	N,N-Dimethylformamide	316
508804	N,N-Dimethylformamide	316
508840	Fixative AFA liquid	354
508861	Formaldehyde 4% w/v buffered at pH 6.9	361
508863	Formaldehyde 4% w/v buffered at pH 6.9	361
508871	Formaldehyde acetic	362
508881	Fixative Davidson liquid	355
509144	Organic Standard: PCB multielement mixture	543
509145	Organic Standard: PCB multielement mixture	543
509146	Organic Standard: PCB multielement mixture	543
509147	Organic Standard: PCB multielement mixture	543
509148	Organic Standard: PCB multielement mixture	543
524102	Methanol	491
524103	Methanol	491
524125	Ethanol absolute anhydrous	332
5241252	Ethanol absolute anhydrous	332
524132	Ethanol 96°	335
524135	Ethanol 96°	335
524161	Propan-2-ol 70%	626
524170	Propan-2-ol	625
524171	Propan-2-ol	625
524181	Propan-2-ol 70%	627
524182	Propan-2-ol 70%	627
524183	Propan-2-ol 70%	627
524184	Propan-2-ol 70%	627
524195	Propan-2-ol 70%	626
524263	n-Heptane 99%	381
524265	n-Heptane 99%	381
524267	n-Heptane 99%	381
524314	Dichloromethane	303
524319	Dichloromethane	303
524381	Heptane mixture of isomers	382
524391	Isopentane	441
524411	Mixture C.H.M.	507
524412	Mixture C.H.M.	507
524501	Sodium hydroxide solution 5% w/v	684
524502	Sodium hydroxide solution 5% w/v	684
524505	Sodium hydroxide solution 20% w/v	683
524506	Sodium hydroxide solution 10% w/v	684
524507	Sodium hydroxide solution 10% w/v	684
524510	Sodium hydroxide solution 32%	681
524520	Acetic acid glacial	132
524521	Acetic acid glacial	132
524525	Hydrochloric acid 37%	395
524526	Hydrochloric acid 37%	395
524530	Nitric acid 69.5%	529
524531	Nitric acid 69.5%	529
524535	Nitric acid 65%	531
524536	Nitric acid 65%	531
524540	Sulfuric acid 96%	732
524541	Sulfuric acid 96%	732
524543	Sulfuric acid 96%	732
524561	Hydrochloric acid 20%	399
524605	Acetic acid 1 mol/l (1N)	136
524611	Acetic acid 0.03 mol/l (0.03N)	136
524621	Sodium hydroxide 1 mol/l (1N)	687
524631	Sodium hydroxide 0.1 mol/l (0.1N)	690
524651	Sodium hydroxide 6 mol/l (6N)	684
524671	Sodium hydroxide 2 mol/l (2N)	686
524725	Sodium nitrite solution 500 g/l	697
524732	Sodium hydroxide 3 mol/l (3N)	685
524761	Sodium hydroxide 1 mol/l (1N)	687
524910	Reagent for lipolysis	635
524920	Formaldehyde 4% w/v buffered at pH 6.9	361
524930	Formaldehyde 30% w/v	360
524952	Phosphate buffer pH 6.8	575
524965	Phosphate buffer pH 7.4	575
524971	Buffer pH 4.5	219
525101	Methanol	489
525102	Methanol	489
525161	Propan-2-ol	623
525320	Dichloromethane	304
525321	Dichloromethane	304
525921	Phosphate buffer pH 7.2	575
525925	Phosphate buffer pH 7.2	575
525992	Octane 80 blend	537
525993	Octane 80 blend	537
525994	Octane 80 blend	537
526001	Sodium phosphate tribasic dodecahydrate	704
526151	Isopar G	440
526231	Formaldehyde acetic	362
526261	Fixative Bouin liquid	355
526263001	Fixative AFA liquid	354
526264	Fixative liquid without acetic acid	356
526267	Fixative AFA liquid	354
526268	Fixative Bouin Hollande liquid	355
526270	Fixative Bouin liquid	355
526273	Formaldehyde acetic	362
526274	Fixative FIXALL-HIS liquid	356
526277	Fixative Davidson liquid	355
526311	Fixative Bouin liquid	355
526511	Sodium hydroxide 0.7 mol/l (N/1.4)	687
526512	Sodium hydroxide 5 mol/l (5N)	685
526513	Sodium hydroxide 5 mol/l (5N)	685
526521	Sodium hydroxide solution 32%	681
526531	Hydrochloric acid 0.714 mol/l (N/1.4)	404
526533	Hydrochloric acid 0.0714 mol/l (N/14)	407
526535	Hydrochloric acid 0.2 mol/l (0.2N) in propanol-2	405

Product specifications are subject to changes. Please visit our website for updates.

526536	Hydrochloric acid 0.1 mol/l (0.1N) in propanol-2	407
526537	Hydrochloric acid 0.02 mol/l (0.02N)	408
526545	Acetic acid 45%	134
526546	Acetic acid 45%	134
526602	Sulfuric acid with 6.6 g/l Silver sulfate	743
526605	Sulfuric acid with 10 g/l Silver sulfate	742
526606	Sulfuric acid with 10 g/l Silver sulfate	742
526615	Reagent TBN ASTM D4739	636
526623	ADF Solution	147
526625	ADF Solution	147
526632	Sodium hydroxide solution 5% w/v	684
526634	Sodium hydroxide solution 5% w/v	684
526641	Sodium hydroxide solution 10% w/v	684
526642	Sodium hydroxide solution 10% w/v	684
526644	Sodium hydroxide solution 10% w/v	684
526711	Potassium dichromate 0.04 mol/l (0.24N) in 80 g/l HgSO4	593
526712	Potassium dichromate 0.04 mol/l (0.24N) in 80 g/l HgSO4	593
526741	Sulfuric acid 4 mol/l (8N)	737
526751	Ferroin 0.025 mol/l solution	354
526761	Iron (II) ammonium sulfate 0.12N	429
526911	Formaldehyde 10% v/v according to Lillie	360
526912	Formaldehyde 10% v/v according to Lillie	360
526920	NDF Solution	519
526921	NDF Solution	519
526931	Formaldehyde 4% w/v buffered at pH 6.9	361
526933	Formaldehyde 4% w/v buffered at pH 6.9	361
526934	Formaldehyde 4% w/v with sodium chloride	361
526936	Formaldehyde 4% w/v buffered at pH 6.9	361
526937	Formaldehyde 4% w/v buffered at pH 6.9	361
526941	NDF Plus solution	519
527591	Orthophosphoric acid 85%	547
527592	Orthophosphoric acid 85%	547
527600	Hydrochloric acid 37%	394
527601	Hydrochloric acid 37%	394
527620	Hydrogen peroxide solution 30%	414
527621	Hydrogen peroxide solution 30%	414
527630	Sulfuric acid 96%	731
527631	Sulfuric acid 96%	731
527640	Methanol	490
527641	Methanol	490
527650	Acetone	139
527651	Acetone	139
527655	Acetone	139
527670	Nitric acid 69.5%	528
527671	Nitric acid 69.5%	528
527680	Ethanol absolute anhydrous	332
527681	Ethanol absolute anhydrous	332
527690	Propan-2-ol	624
527691	Propan-2-ol	624
527696	Propan-2-ol	624
528070	Petroleum ether 35 - 60°C	568
528071	Petroleum ether 35 - 60°C	568
528101	Methanol	492
528131	Ethanol absolute anhydrous	333
528151	Ethanol 96°	336
528152	Ethanol 96°	336
528153	Ethanol 96°	336
528154	Ethanol 96°	336
528170	Ethanol 70% v/v	337
528191	Ethanol 70°modified	339
528192	Ethanol 70°modified	339
528201	Acetone	141
528203	Acetone	141
528204	Acetone	141
528206	Acetone	141
528215	Cyclohexane	287
528216	Cyclohexane	287
528217	Cyclohexane	287
528220	N,N-Dimethylformamide	316
528221	N,N-Dimethylformamide	316
528224	n-Heptane 99%	381
528225	n-Heptane 99%	381
528226	n-Heptane 99%	381
528228	n-Heptane 99%	381
528231	Toluene	767
528232	Toluene	767
528233	Toluene	767
528245	Heptane mixture of isomers	383
528246	Heptane mixture of isomers	383
528247	Heptane mixture of isomers	383
528251	Xylene, mix of isomers	791
528252	Xylene, mix of isomers	791
528257	Pyridine	631
528260	Methylcyclohexane	498
528261	Methylcyclohexane	498
528262	Methylcyclohexane	498
528264	Methylcyclohexane	498
528275	Diethyl ether	310
528276	Diethyl ether	310
528280	Petroleum ether 35 - 60°C	568
528281	Petroleum ether 35 - 60°C	568
528282	Petroleum ether 35 - 60°C	568
528283	Petroleum ether 40 - 60°C	567
528294	Ethyl acetate	343
528295	Ethyl acetate	343
528296	Ethyl acetate	343
528297	Ethyl acetate	343
528299	Ethyl acetate	343
528300	Butanol-1	226
528301	Butanol-1	226
528325	Chloroform	262
528326	Chloroform	262
528327	Chloroform	262
528328	Chloroform	262
528329	Chloroform	262
528332	Cyclohexanone	288
528335	Dimethylsulphoxide	319
528340	N-Methyl-2-pyrrolidone	503
528341	N-Methyl-2-pyrrolidone	503
528343	N-Methyl-2-pyrrolidone	503
528346	N-Methyl-2-pyrrolidone	503
528370	Dichloromethane	304
528371	Dichloromethane	304
528372	Dichloromethane	304
528377	Dichloromethane	304
528379	Dichloromethane	304
528461	Dichloromethane	304
528462	Dichloromethane	304
528463	Dichloromethane	304
528464	Dichloromethane	304
528481	Tetrahydrofuran	753
528491	Isopentane	441
528492	Isopentane	441
528501	Ammonia solution 32%	160
528503	Ammonia solution 32%	160
528525	Hydrochloric acid 37%	395
528530	Nitric acid 67.5 (42° Be)	530
528535	Orthophosphoric acid 85%	548
528541	Sulfuric acid 50%	735
528550000	Hydrochloric acid 6 mol/l (6N)	401
528571	Hydrochloric acid 0.1 mol/l (0.1N)	406
528573	Hydrochloric acid 0.1 mol/l (0.1N)	406
528583	Hydrochloric acid 1 mol/l (1N)	404
528584	Hydrochloric acid 1 mol/l (1N)	404
528651	Hydrochloric acid 6 mol/l (6N)	401
528661	Hydrochloric acid 0.1 mol/l (0.1N)	406
528662	Hydrochloric acid 0.1 mol/l (0.1N)	406
528671	Hydrochloric acid 1 mol/l (1N)	404
528672	Hydrochloric acid 1 mol/l (1N)	404
528673	Hydrochloric acid 1 mol/l (1N)	404
528681	Hydrochloric acid 4 mol/l (4N)	402
528691	Hydrochloric acid 2 mol/l (2N)	403
528731	Hydrochloric acid 5 mol/l (5N)	402
528761	Ethanol absolute denaturated	338
528762	Ethanol absolute denaturated	338
528763	Ethanol absolute denaturated	338
528764	Ethanol absolute denaturated	338
528765	Ethanol absolute denaturated	338
528766	Ethanol absolute denaturated	338
528771	Ethanol 95° denaturated	338
528772	Ethanol 95° denaturated	338
528773	Ethanol 95° denaturated	338
528774	Ethanol 95° denaturated	338
528775	Ethanol 95° denaturated	338
528940	Hexane mixture of isomers	388
528941	Hexane mixture of isomers	388
528942	Hexane mixture of isomers	388
528950	n-Hexane 99%	386
528951	n-Hexane 99%	386
528952	n-Hexane 99%	386
528960	Isooctane	440
528961	Isooctane	440
528962	Isooctane	440
528963	Isooctane	440
528970	tert-Butylmethylether	231
528971	tert-Butylmethylether	231
528974	tert-Butylmethylether	231
528975	Ethyl methyl ketone	350
528976	Ethyl methyl ketone	350
528977	Ethyl methyl ketone	350
528979	tert-Butylmethylether	231
528980	Methyl isobutyl ketone	501
528981	Methyl isobutyl ketone	501
528993	n-Pentane 99%	558
528994	n-Pentane 99%	558
528995	n-Pentane 99%	558
528996	n-Pentane 99%	558
528997	n-Pentane 99%	558
529021	2-Phenylethanol	573
529022	2-Phenylethanol	573
529050	Propionic acid	628
529091	Propan-2-ol	626
529092	Propan-2-ol	626
529093	Propan-2-ol	626
529100	Methanol	492
529121	Ethanol absolute anhydrous	333
529122	Ethanol absolute anhydrous	333
5291222	Ethanol absolute anhydrous	333
529124	Ethanol absolute anhydrous	333
529125	Ethanol absolute anhydrous	333
529141	Ethanol 96°	336
5291412	Ethanol 96°	336
529152	Ethanol 96°	336
529165	Propan-2-ol	626
529174	Propan-2-ol	625

INDEX BY PRODUCT CODE

529180	Ethanol 60% v/v	337
529181	Ethanol 60% v/v	337
529183	Ethanol 70% v/v	336
529184	Ethanol 70% v/v	337
529184000	Ethanol 70% v/v	337
529186	Ethanol 70% v/v	337
529187000	Ethanol 70% v/v	336
529189	Ethanol 70% v/v	336
529221	Ethyl acetate	342
529222	Ethyl acetate	342
529261	Ethanol 50% v/v	337
529301	Chloroform	262
529311	Mix Ethanol absolute/Diethyl ether 50/50 (w/w)	507
529371	Mix Diethyl ether / Ethanol 95% 50/50(w/w) - 30 mg/l phenolphtaleine	507
529381	Mix Ethanol absolute/Diethyl ether 50/50 (v/v)	508
529510	Sodium phosphate dibasic dodecahydrate	702
540001	Melting point standards	483
540002	Melting point standards	483
540003	Melting point standards	483
540004	Melting point standards	483
540005	Melting point standards	483
540006	Melting point standards	483
540007	Melting point standards	483
540008	Melting point standards	483
540009	Melting point standards	483
540010	Melting point standards	483
540011	Melting point standards	483
540012	Melting point standards	483
540013	Melting point standards	483
540014	Melting point standards	483
540101	Refractive Index standards	637
540102	Refractive Index standards	637
540103	Refractive Index standards	637
540104	Refractive Index standards	637
540105	Refractive Index standards	637
540106	Refractive Index standards	637
540107	Refractive Index standards	637
540108	Refractive Index standards	637
540201	Brix Standards	209
540202	Brix Standards	209
540203	Brix Standards	209
540204	Brix Standards	209
540205	Brix Standards	209
540206	Brix Standards	209
540207	Brix Standards	209
540208	Brix Standards	209
540209	Brix Standards	209
540210	Brix Standards	209
540220	Brix Standards	209
540221	Brix Standards	209
540222	Brix Standards	209
540223	Brix Standards	209
540224	Brix Standards	209
540225	Brix Standards	209
540301	Osmolality Standards	550
540302	Osmolality Standards	550
540303	Osmolality Standards	550
540304	Osmolality Standards	550
540305	Osmolality Standards	550
540306	Osmolality Standards	550
540307	Osmolality Standards	550
540308	Osmolality Standards	550
540309	Osmolality Standards	550
540310	Osmolality Standards	550
540351	Osmolality Standards Protein Based	550
540352	Osmolality Standards Protein Based	550
540353	Osmolality Standards Protein Based	550
540354	Osmolality Standards Urine Based	550
540355	Osmolality Standards Urine Based	550
540401	Density Standard	292
540402	Density Standard	292
540403	Density Standard	292
540404	Density Standard	292
540405	Density Standard	292
540406	Density Standard	292
540407	Density Standard	292
540408	Density Standard	292
540409	Density Standard	292
540410	Density Standard	292
540411	Density Standard	292
540412	Density Standard	292
540413	Density Standard	292
540414	Density Standard	292
540415	Density Standard	292
540416	Density Standard	292
540417	Density Standard	292
540418	Density Standard	292
540420	Density Standard	292
540421	Density Standard	292
540422	Density Standard	292
540451	Density Standard	293
540452	Density Standard	293
540453	Density Standard	293
540454	Density Standard	293
540455	Density Standard	293
540456	Density Standard	293
540457	Density Standard	293
540458	Density Standard	293
540459	Density Standard	293
540460	Density Standard	293
540461	Density Standard	293
540462	Density Standard	293
540463	Density Standard	293
540464	Density Standard	293
540465	Density Standard	293
540470	Density Standard	293
540471	Density Standard	293
540472	Density Standard	293
540473	Density Standard	293
540566	Density Standard	293
540567	Density Standard	293
540568	Density Standard	293
540569	Density Standard	293
540601	ASTM colour standards	187
540602	ASTM colour standards	187
540603	ASTM colour standards	187
540604	ASTM colour standards	187
540701	Gardner Colour Standards	368
540702	Gardner Colour Standards	368
540703	Gardner Colour Standards	368
540704	Gardner Colour Standards	368
540705	Gardner Colour Standards	368
540706	Gardner Colour Standards	368
540707	Gardner Colour Standards	368
540708	Gardner Colour Standards	368
540709	Saybolt Colour Standards	644
540710	Saybolt Colour Standards	644
540711	Saybolt Colour Standards	644
540712	Saybolt Colour Standards	644
540713	Saybolt Colour Standards	644
540714	Saybolt Colour Standards	644
540715	Saybolt Colour Standards	644
540801	Viscosity standards	785
540802	Viscosity standards	785
540803	Viscosity standards	785
540804	Viscosity standards	785
540805	Viscosity standards	785
540806	Viscosity standards	785
540807	Viscosity standards	785
540808	Viscosity standards	785
540809	Viscosity standards	785
540810	Viscosity standards	785
540811	Viscosity standards	785
554192	Eukitt	351
554193	Eukitt	351
554194	Eukitt	351
570011	Karl Fischer reagent 1 component 5 mg H2O/ml	443
570021	Karl Fischer reagent 1 component 2 mg H2O/ml	443
570031	Karl Fischer solvent for oils one component	444
570041	Karl Fischer solvent for aldehydes and ketones one component	444
570051	Karl Fischer titrant 2 component 5 mg H2O/ml	444
570061	Karl Fischer titrant 2 component 2 mg H2O/ml	444
570071	Karl Fischer solvent 2 component	445
570081	Karl Fischer titrant 1 component 5 mg H2O/ml for aldehydes and ketones	443
570091	Karl Fischer solvent 2 component for aldehydes and ketones - Methanol free	445
570101	Karl Fischer solvent for oils 2 component	445
570111	Karl Fischer 2 component buffered solvent	445
570121	Karl Fischer anolyte solution	446
570131	Karl Fischer anolyte solution for cells with and without diagrams	447
570141	Karl Fischer anolyte solution - CFC free	446
570151	Karl Fischer anolyte solution, oven	446
570161	Karl Fischer anolyte solution for ketones and aldehydes - Methanol free	446
570171	Karl Fischer anolyte solution for oils	447
570181	Karl Fischer catholyte solution	447
570191	Karl Fischer catholyte solution for aldehydes and ketones	447
570201	Karl Fischer water standard 0.10 mg/g	449
570211	Karl Fischer water standard 1.0 mg/g	448
570221	Karl Fischer water standard 10.0 mg/g	448
570231	Karl Fischer water standard 5.0 mg/ml	448
570241	Karl Fischer reagent 1 component 5 mg H2O/ml	443
570251	Karl Fischer reagent 1 component 2 mg H2O/ml	443
575001	Standard solution 5 µS/cm	716
575011	Standard solution 10 µS/cm	716
575021	Standard solution 20 µS/cm	716
575031	Standard solution 50 µS/cm	716
575041	Standard solution 84 µS/cm	717
575051	Standard solution 100 µS/cm	717
575061	Standard solution 147 µS/cm	717
575071	Standard solution 200 µS/cm	717
575081	Standard solution 500 µS/cm	717
575091	Standard solution 1000 µS/cm	718
575101	Standard solution 1413 µS/cm	718

Product specifications are subject to changes. Please visit our website for updates.

INDEX BY PRODUCT CODE

575111	Standard solution 5000 µS/cm	718
575121	Standard solution 10000 µS/cm	718
575131	Standard solution 12880 µS/cm	718
575141	Standard solution 20000 µS/cm	719
575151	Standard solution 50000 µS/cm	719
575161	Standard solution 100000 µS/cm	719
575171	Standard solution 150000 µS/cm	719
575181	Standard solution 200000 µS/cm	720
575191	Standard solution 300000 µS/cm	720
575201	Standard solution 350000 µS/cm	720
575211	Standard solution 450000 µS/cm	720
575221	Standard solution 500000 µS/cm	721
575231	Standard solution 1.30 µS/cm	715
611000401	Acetic acid 30%	134
611000402	Acetic acid 12%	135
611000501	Acetic anhydride	137
611002501	Ethanol 96°	334
611003701	Aminohippuric acid reagent	158
611004601	Aminopyrazolone solution	159
611004701	Ammonia solution 17%	164
611004702	Ammonia solution diluted	164
611004703	Ammonia solution diluted	164
611005201	Ammonium carbonate solution 158 g/l	168
611007301	Anisaldehyde solution	181
611007302	Anisaldehyde solution	181
611007701	Antimony trichloride	183
611009300	Barium chloride dihydrate	192
611009301	Barium chloride solution 61 g/l	193
611009303	Barium chloride solution 61 g/l	193
611009309	Barium chloride solution 61 g/l	193
611009401	Barium hydroxide solution 47.3 g/l	195
611009409	Barium hydroxide solution 47.3 g/l	195
611011601	Biuret reagent	205
611012401	Bromine solution	210
611012402	Bromine water	211
611012409	Bromine water	211
611012601	Bromocresol green solution	211
611012602	Bromocresol green - Methyl red solution	212
611012701	Bromocresol purple solution	213
611012801	Bromophenol blue solution	214
611012802	Bromophenol blue solution	214
611012803	Bromophenol blue solution	214
611012901	Bromothymol blue solution	215
611012903	Bromothymol blue solution	215
611015201	Calcium sulfate hemihydrate solution	245
611015202	Calcium sulfate hemihydrate solution	245
611017901	Chloral hydrate	256
611021000	Citric acid monohydrate	269
611022001	Congo red solution	274
611022501	Copper (II) sulfate solution 12.5%	281
611022600	Ammoniacal solution of copper tetrammine	165
611022801	Cresol red R solution in Sodium hydroxide 0.1N/ Ethanol 96% / Water	283
611022901	Crystal violet solution 0.5% in anhydrous acetic acid	284
611023100	Cupri-citric solution	284
611023300	Cupri-tartaric solution	285
611032001	1,4-Dioxane	321
611032002	1,4-Dioxane	321
611032003	1,4-Dioxane	321
611032101	Diphenylamine solution 1% in sulfuric acid	322
611032102	Diphenylamine solution 1% in sulfuric acid	322
611032109	Diphenylamine solution 1% in sulfuric acid	322
611036401	Ethylene oxide solution	349
611036408	Ethylene oxide solution	349
611037702	Iron (III) ammonium sulfate solution 100 g/l	433
611037703	Iron (III) ammonium sulfate solution 100 g/l	433
611038100	Ferroin 0.025 mol/l solution	354
611039101	Formaldehyde 35% w/w	359
611039401	Fuchsin solution decolorised	365
611039402	Fuchsin solution decolorised	365
611043101	Holmium perchlorate in solution	391
611043501	Hydrochloric acid 25% w/v	398
611043503	Hydrochloric acid, dilute	409
611043504	Hydrochloric acid, dilute	409
611043507	Hydrochloric acid, brominated	409
611044301	Hydroxylamine solution, alcoholic	417
611045601	Indigo carmine solution	422
611045901	Iodine bromide solution	426
611046300	Iodoplatinate reagent	427
611046309	Iodoplatinate reagent	427
611048001	Lanthanum nitrate solution 50 g/l	454
611048101	Lead (II) acetate cotton	457
611048102	Lead (II) acetate paper	457
611048103	Lead (II) acetate solution 95 g/l	458
611048301	Lead (II) nitrate solution 33 g/l	459
611048400	Lead (II) acetate basic solution	457
611049901	Magnesium oxide	473
611050501	Malachite green solution 0.5% in acetic acid anhydrous	476
611052101	Mercuric bromide paper	484
611052201	Mercury chloride solution 54 g/l	485
611052600	Mercury (II) sulfate solution	486
611053203	Methanol, hydrochloric	493
611054801	Methyl orange mixed solution	502
611054802	Methyl Orange solution 0.1%	501
611055101	Methyl red mixed solution	504
611055102	Methyl red solution	504
611055901	Dichloromethane acidified with 1% hydrochloric acid	305
611056700	Molybdovanadic reagent	511
611056801	Eriochrome black T	329
611057601	alpha-Naphtholbenzein solution 0.2% in acetic acid	518
611058301	Ninhydrin and Tin (II) chloride reagent	526
611058304	Ninhydrin solution	526
611058305	Ninhydrin solution	526
611058402	Nitric acid, dilute	534
611058409	Nitric acid, dilute	534
611062201	Pararosaniline solution, decolorised	556
611062901	Perchloric acid solution	562
611063601	Phenol red solution	571
611063603	Phenol red solution	571
611063702	Phenolphthalein solution 0.1%	572
611063703	Phenolphthalein solution 1% in ethanol	571
611063709	Phenolphtalein solution 0.1% in ethanol	572
611064501	Phenylhydrazine hydrochloride solution	573
611065000	Phosphomolybdotungstic reagent	576
611065200	Phosphotungstic acid solution	578
611065801	Picric acid solution	579
611065802	Picric acid solution	579
611068900	Potassium carbonate	586
611069101	Potassium chloride 0.1 mol/l (0.1N)	589
611069201	Potassium chromate 5% solution	591
611069501	Potassium dichromate solution 106 g/l	594
611069502	Potassium dichromate solution 0.5%	592
611069509	Potassium dichromate solution 0.5%	592
611069601	Potassium phosphate monobasic 0.2 mol/l (0.2N)	615
611069700	Potassium ferricyanide	594
611069800	Potassium ferrocyanide trihydrate	595
611069801	Potassium ferrocyanide solution 53 g/l	596
611070001	Potassium hydrogen phthalate 0.2 mol/l (0.2N)	598
611070301	Potassium hydroxide 2 mol/l (2N) in ethanol	602
611070302	Potassium hydroxide 0.5 mol/l (0.5N) in ethanol	603
611070303	Potassium hydroxide solution 3% in ethanol	601
611070502	Potassium iodide solution	608
611070504	Potassium iodide solution	608
611070505	Potassium iodide solution	608
611070600	Potassium iodobismuthate solution	609
611070602	Potassium iodobismuthate solution	609
611070901	Potassium permanganate and phosphoric acid solution	613
611070902	Potassium permanganate solution 3%	612
611071302	Potassium pyroantimonate solution	616
611071303	Potassium pyroantimonate solution	616
611071309	Potassium pyroantimonate solution	616
611071600	Potassium tetraiodomercurate solution, alkaline	619
611071801	Potassium thiocyanate solution	621
611075500	Salicylaldehyde azine	642
611078200	Silver manganese paper	655
611078301	Silver nitrate solution	658
611078302	Silver nitrate solution	658
611078306	Silver nitrate solution	658
611078307	Silver nitrate solution	658
611079301	Sodium carbonate solution	670
611081401	Sodium hydroxide solution 20% w/v	683
611081402	Sodium hydroxide solution	691
611081404	Sodium hydroxide solution	691
611081405	Sodium hydroxide solution, methanolic	691
611081600	Sodium hypochlorite solution in water	692
611081609	Sodium hypochlorite solution in water	692

INDEX BY PRODUCT CODE

611083901	Sodium sulfide nonahydrate solution	707
611083902	Sodium sulfide nonahydrate solution	707
611085001	Tin (II) chloride solution	761
611085103	Starch soluble solution	721
611085104	Starch soluble solution	721
611086500	Sulfomolybdic reagent	729
611086804	Sulfuric acid, dilute	743
611089602	Thioacetamide solution 40 g/l	755
611089603	Thioacetamide solution 40 g/l	755
611089700	Thioglycolic acid 80%	755
611090601	Thymol blue solution	758
611090701	Thymolphthalein solution 0.1% in ethanol	759
611091202	Titanium trichloride-sulfuric acid reagent	765
611094201	Tris(hydroxymethyl) aminomethane solution	777
611095302	Vanillin solution, phosphoric	784
611095501	Water	787
611095506	Water	787
611096501	Zinc, activated	797
611096601	Zinc chloride-formic acid solution	800
611096603	Zinc chloride solution, iodinated	800
611123001	o-Tolidine solution	765
612000200	Benzoic acid	199
612000300	Potassium bromate	584
612000400	Potassium hydrogen phthalate	597
612000500	Sodium carbonate monohydrate	668
612000600	Sodium chloride	670
612000700	Sulfanilic acid	728
612000800	Zinc standard solution	797
612002100	Degree of coloration of liquids: primary solutions	291
612002200	Degree of coloration of liquids: primary solutions	291
612002300	Degree of coloration of liquids: primary solutions	291
612201100	Primary opalescent suspension	621
612201101	Primary opalescent suspension	621
612202100	Degree of coloration of liquids: primary solutions	290
612202200	Degree of coloration of liquids: primary solutions	290
612202300	Degree of coloration of liquids: primary solutions	290
612202400	Hydrochloric acid, dilute	409
612202510	Degree of coloration of liquids: standard solutions	291
612202520	Degree of coloration of liquids: standard solutions	291
612202530	Degree of coloration of liquids: standard solutions	291
612202540	Degree of coloration of liquids: standard solutions	291
612202550	Degree of coloration of liquids: standard solutions	291
612203168	Buffer pH 1.68	216
612203401	Buffer pH 4	218
612203687	Buffer pH 6.88	220
612203741	Buffer pH 7.4	222
612203918	Buffer pH 9.22	223
613000100	Cerium (IV) ammonium nitrate 0.1 mol/l	251
613000200	Cerium (IV) ammonium nitrate 0.01 mol/l	251
613000300	Cerium (IV) ammonium sulfate 0.1 mol/l	252
613000301	Cerium (IV) ammonium sulfate 0.1 mol/l	252
613000400	Cerium (IV) ammonium sulfate 0.01 mol/l	252
613000500	Ammonium thiocyanate 0.1 mol/l (0.1N)	179
613000501	Ammonium thiocyanate 0.1 mol/l (0.1N)	179
613000600	Barium chloride 0.1 mol/l (0.2N)	193
613000700	Barium perchlorate 0.05 mol/l	196
613000900	Hyamine 1622 solution 0.004M	392
613000901	Hyamine 1622 solution 0.004M	392
613001000	Bromide - bromate 0.0167 mol/l	210
613001100	Cerium (IV) sulfate 0.1 mol/l	253
613001101	Cerium (IV) sulfate 0.1 mol/l	253
613001300	Iron (III) ammonium sulfate 0.1 mol/l	433
613001400	Iron (II) sulfate 0.1mol/l	431
613001500	Hydrochloric acid 6 mol/l (6N)	401
613001700	Hydrochloric acid 2 mol/l (2N)	403
613001800	Hydrochloric acid 1 mol/l (1N)	403
613001801	Hydrochloric acid 1 mol/l (1N)	403
613002100	Hydrochloric acid 0.1 mol/l (0.1N)	406
613002101	Hydrochloric acid 0.1 mol/l (0.1N)	406
613002700	Iodine 0.05 mol/l (0.1N)	425
613002900	Iodine 0.01 mol/l (0.02N)	426
613003100	Lead (II) nitrate 0.1 mol/l	459
613003101	Lead (II) nitrate 0.1 mol/l	459
613003900	Perchloric acid 0.1 mol/l (0.1N)	561
613004000	Perchloric acid 0.05 mol/l (0.05N)	562
613004200	Potassium bromate 0.033 mol/l (0.198N)	585
613004300	Potassium bromate 0.02 mol/l (0.12N)	585
613004600	Potassium dichromate 0.0167 mol/l (0.1N)	593
613004800	Potassium hydroxide 0.1 mol/l (0.1N)	604
613004900	Potassium hydroxide 0.5 mol/l (0.5N) in ethanol	603
613005001	Potassium hydroxide in solution	606
613005100	Potassium hydroxide 0.1 mol/l (0.1N) in ethanol	605
613005200	Potassium iodate 0.05 mol/l (0.3N)	606
613005300	Potassium permanganate 0.02 mol/l (0.1N)	612
613005301	Potassium permanganate 0.02 mol/l (0.1N)	612
613005309	Potassium permanganate 0.02 mol/l (0.1N)	612
613005600	Silver nitrate 0.1 mol/l (0.1N)	657
613005800	Sodium arsenite 0.1 mol/l (0.2N)	664
613005900	Ethylenediaminetetraacetic acid disodium salt 0.1 mol/l (0.2N)	346
613005901	Ethylenediaminetetraacetic acid disodium salt 0.1 mol/l (0.2N)	346
613006300	Sodium hydroxide 1 mol/l (1N)	686
613006301	Sodium hydroxide 1 mol/l (1N)	686
613006600	Sodium hydroxide 0.1 mol/l (0.1N)	689
613006601	Sodium hydroxide 0.1 mol/l (0.1N)	689
613007000	Sodium hydroxide 0.1 mol/l (0.1N) in ethanol	690
613007001	Sodium hydroxide 0.1 mol/l (0.1N) in ethanol	690
613007100	Sodium methoxide 0.1 mol/l	695
613007101	Sodium methoxide 0.1 mol/l	695
613007200	Sodium nitrite 0.1 mol/l (0.1N)	697
613007300	Sodium thiosulfate 0.1 mol/l (0.1N)	712
613007301	Sodium thiosulfate 0.1 mol/l (0.1N)	712
613007800	Sulfuric acid 0.5 mol/l (1N)	738
613008000	Sulfuric acid 0.05 mol/l (0.1N)	740
613008001	Sulfuric acid 0.05 mol/l (0.1N)	740
613008300	Tetrabutylammonium hydroxide solution 0.1 mol/l (0.1N)	749
613008400	Tetrabutylammonium hydroxide 0.1 mol/l (0.1N) in 2-propanol	749
613008600	Zinc sulfate 0.1 mol/l (0.2N)	802
613008601	Zinc sulfate 0.1 mol/l (0.2N)	802
613008700	Cupriethylenediamine hydroxide 1 mol/l	285
613008800	Hydrochloric acid 0.1 mol/l (0.1N) in ethanol	407
613009100	Potassium hydroxide 1 mol/l (1N)	602
613009400	Iodine 0.5 mol/l (1N)	425
613009601	Barium perchlorate 0.025 mol/l	196
613009700	Lead (II) nitrate 0.05 mol/l	459
613010100	Lanthanum nitrate 0.1 mol/l	454
614000100	Acetone	140
614000200	Buffer pH 2	217
614000500	Phosphate buffer pH 3.0	576
614000501	Phosphate buffer pH 3.0	576
614000600	Buffer pH 3.5	217
614000601	Buffer pH 3.5	217
614000900	Buffer pH 3.7	218
614001400	Acetate buffer pH 4.6	130
614001700	Buffer pH 5.2	219
614002000	Phosphate buffer pH 5.5	575
614002200	Acetate buffer pH 6.0	130
614002400	Phosphate buffer pH 6.0	575
614003400	Phosphate buffer pH 6.8	575
614003500	Buffer pH 7	220
614004600	Buffer pH 7.4	222
614004800	Phosphate buffer pH 7.4	575
614004900	Tris(hydroxymethyl) aminomethane sodium chloride buffer solution pH 7.4	778
614006200	Tris(hydroxymethyl) aminomethane buffer solution pH 8.1	777
614006600	Tris(hydroxymethyl) aminomethane-EDTA buffer solution pH 8.4	777
614007000	Buffer pH 9	222
614007200	Ammonium chloride buffer solution pH 9.5	170
614007300	Ammonium chloride buffer solution pH 10.0	170
614007301	Ammonium chloride buffer solution pH 10.0	170
614007700	Total-ionic-strength-adjustment buffer	769
614007900	Phosphate buffer pH 2.0	576
614008300	Phosphate buffer pH 9.0	574
614008800	Total-ionic-strength-adjustment buffer	769
614013400	Ammonium chloride buffer solution pH 10.7	169
615000200	Aluminum standard solution	151
615000201	Aluminum standard solution	151
615000202	Aluminum standard solution	151
615000203	Aluminum standard solution	151
615000301	Ammonium standard solution	165
615000302	Ammonium standard solution	165
615000309	Ammonium standard solution	165
615000400	Antimony standard solution	182

Product specifications are subject to changes. Please visit our website for updates.

INDEX BY PRODUCT CODE

615000501	Arsenic standard solution	185
615000509	Arsenic standard solution	185
615000601	Barium standard solution	191
615000609	Barium standard solution	191
615000700	Cadmium standard solution	232
615000709	Cadmium standard solution	232
615000801	Calcium standard solution	234
615000802	Calcium standard solution	234
615000803	Calcium standard solution	234
615000804	Calcium standard solution	234
615000809	Calcium standard solution	234
615000901	Chloride standard solution	257
615000909	Chloride standard solution	257
615001000	Chromium standard solution	264
615001002	Chromium standard solution	264
615001100	Copper standard solution	275
615001209	Ferrocyanide standard solution	354
615001300	Ferricyanide standard solution	353
615001401	Fluoride standard solution	358
615001409	Fluoride standard solution	358
615001601	Iron standard solution	428
615001602	Iron standard solution	428
615001603	Iron standard solution	428
615001605	Iron standard solution	428
615001606	Iron standard solution	428
615001609	Iron standard solution	428
615001700	Lead standard solution	455
615001701	Lead standard solution	455
615001702	Lead standard solution	455
615001703	Lead standard solution	455
615001704	Lead standard solution	455
615001705	Lead standard solution	455
615001706	Lead standard solution	455
615001709	Lead standard solution	455
615001801	Magnesium standard solution	468
615001802	Magnesium standard solution	468
615001803	Magnesium standard solution	468
615001809	Magnesium standard solution	468
615001900	Mercury standard solution	484
615001901	Mercury standard solution	484
615002001	Nickel standard solution	521
615002002	Nickel standard solution	521
615002009	Nickel standard solution	521
615002101	Nitrate standard solution	527
615002102	Nitrate standard solution	527
615002109	Nitrate standard solution	527
615002200	Phosphate standard solution	574
615002401	Potassium standard solution	582
615002402	Potassium standard solution	582
615002409	Potassium standard solution	582
615002500	Selenium standard solution	647
615002501	Selenium standard solution	647
615002609	Silver standard solution	653
615002701	Sodium standard solution	661
615002709	Sodium standard solution	661
615002801	Sulfate standard solution	728
615002802	Sulfate standard solution	728
615002809	Sulfate standard solution	728
615002900	Sulfite standard solution	729
615003000	Thallium standard solution	754
615003101	Tin standard solution	760
615003109	Tin standard solution	760
615003200	Titanium standard solution	763
615003300	Vanadium standard solution	783
615003402	Zinc standard solution	797
615003403	Zinc standard solution	797
615003409	Zinc standard solution	797
615003500	Zirconium standard solution	804
615003600	Palladium standard solution	553
615003700	Glyoxal standard solution	374
615003809	Iodine 10 ppm	426
615004100	Chloride standard solution	257
615004200	Phosphate standard solution	574
615004300	Cobalt standard solution	270
615004400	Germanium standard solution	369
615004500	Manganese standard solution	479
615005100	Potassium standard solution	582
615005300	Bismuth standard solution	202
615005700	Sodium standard solution	661
615005800	Manganese standard solution	479
616000001	Arsenic trioxide solution	185
616000050	Zinc sulfate 0.1 mol/l (0.2N)	803
616001001	Silver nitrate solution	658
616001008	Sodium hypochlorite solution in water	692
616001018	Barium chloride 30 g/l	194
616001028	Cobalt (II) chloride in solution	272
616001037	Copper (II) sulphate solution	281
616001038	Copper (II) sulphate solution	281
616001047	Iron chloride in solution	434
616001048	Iron chloride in solution	434
616001057	Cobalt (II) chloride in solution	272
617000001	Arsenic trioxide solution	185
617000101	Bromocresol green solution	212
617000111	Methyl red solution	504
617000121	Sodium hydroxide 1 mol/l (1N)	687
617000131	Phenolphthalein solution 1% in ethanol	571
617000141	Bromine solution	210
617000151	Ammonia solution 6N	164
617000161	Barium chloride 30 g/l	194
617000181	Sodium hypochlorite solution in water	692
617000191	Hydrochloric acid 1 mol/l (1N)	404
617000201	Silver nitrate solution	659
617000211	Thioacetamide	754
617000221	Titanium trichloride-sulfuric acid reagent	765
617000281	Sulfuric acid 1 mol/l (2N)	737
617000291	Zinc sulfate 0.1 mol/l (0.2N)	803
617000301	Lead (II) acetate cotton	457
617000311	Potassium pyroantimonate solution	616
617000321	Lead (II) nitrate	458
618000001	Degree of coloration of liquids: primary solutions	291
618000002	Degree of coloration of liquids: primary solutions	291
618000003	Degree of coloration of liquids: primary solutions	291
E306702	Sulfuric acid 50%	735
E306704	Sulfuric acid 50%	735
E307661	Hydrogen peroxide solution 3.5% w/v	415
E406101	Phosphosulfuric acid	578
E409141	Perchloric acid 0.01 mol/l (0.01N)	562
E411131	L(+) Tartaric acid solution 20% in water	746
E411921	Aqueous calcium hydroxide	184
E413903	Isoamyl alcohol	436
E415932	Alizarin saturated solution in ethanol	149
E419361	Ammonium citrate solution 20%	171
E420001	Ammonia solution 10%	164
E420002	Ammonia solution 10%	164
E420371	Ammonium molybdate solution 2.5% in nitric acid	174
E420521	Ammonium oxalate solution 4%	175
E423562	Methyl Orange solution 0.1%	501
E423982	Silver nitrate solution 5%	656
E424001	Silver nitrate solution 2.9075%	656
E425101	Barium chloride solution 10%	193
E425301	Barium hydroxide solution 5%	195
E425742	Benedict's reagent	197
E428541	Alkali Blue 6B solution 2% in ethanol	150
E428665	Bromophenol blue solution 0.4% in ethanol	213
E428715	Bromothymol blue 0.4% in ethanol	215
E429011	Methylene blue solution 1%	499
E429031	Methylene blue saturated solution	498
E429235	Thymol blue 0.4% in ethanol	757
E433427	Calcium chloride solution 0.025%	239
E441071	Cupriethylenediamine solution	285
E441581	Diacetyldioxime solution 1% in ethanol	296
E446302	Ehrlich's reagent	327
E449936	Fehling's B reagent	353
E449937	Fehling's B reagent	353
E450043	o-Phenantroline-Iron (II) sulphate solution in sulphuric acid	569
E451191	Phenolphthalein solution 1% in ethanol	572
E451192	Phenolphthalein solution 1% in ethanol	572
E451521	Iron (III) ammonium sulfate solution 33% in nitric acid	432
E451653	Iron (III) chloride solution 4.5%	434
E453612	Giemsa's reagent	369
E453613	Giemsa's reagent	369
E453615	Giemsa's reagent	369
E454872	Hanus's reagent	379
E454972	Hyamine 1622 solution 0.004M	392
E455256	Idrimeter Erba Solution A	419
E455257	Idrimeter Erba Solution A	419
E455266	Idrimeter Erba Solution B	419
E455267	Idrimeter Erba Solution B	419
E455271	Idrimeter Erba Indicator C	419
E455274	Idrimeter Erba Indicator C	419
E455651	Indicator for ammoniacal nitrogen solution	420
E455661	Indicator universal pH 0-5 hydroalcoholic solution	421
E455662	Indicator universal pH 0-5 hydroalcoholic solution	421
E455702	Indicator universal pH 1-11 hydroalcoholic solution	421
E455706	Indicator universal pH 1-11 hydroalcoholic solution	421
E455711	Indicator universal pH 1-11 water solution	422
E455712	Indicator universal pH 1-11 water solution	422
E458321	Lithium methoxide 0.1 mol/l (0.1N)	464
E458761	Lugol's Reagent Iodine-Iodide Solution	466
E460582	May Grünwald reagent	482
E460583	May Grünwald reagent	482
E460585	May Grünwald reagent	482
E463562	Folin-Ciocalteu's reagent	358
E465776	p-Nitrophenol solution 0,1% in water	535
E467784	Papanicolaou solution EA 50	554
E467785	Papanicolaou solution EA 50	554
E467786	Papanicolaou solution EA 50	554
E467794	Papanicolaou solution OG 6	555
E467795	Papanicolaou solution OG 6	555
E467796	Papanicolaou solution OG 6	555
E470045	Bromocresol purple solution 0.4% in ethanol	212
E471501	Potassium ferrocyanide solution 10%	595

INDEX BY PRODUCT CODE

E471926	Potassium hydrogen phthalate 0.1 mol/l (0.1N)	598
E472151	Potassium hydroxide solution 38% (40° Bé) in water	601
E472152	Potassium hydroxide solution 38% (40° Bé) in water	601
E472221	Potassium hydroxide solution 28%	601
E474381	Potassium thiocyanate solution 5%	620
E474417	Potassium thiocyanate 0.1 mol/l (0.1N)	620
E475632	Copper (I) chloride solution 7% in ammonia	277
E476211	Copper (II) sulfate solution 12.5%	281
E476805	o-Cresol Red solution 0.2% in ethanol	283
E476845	Phenol Red solution 0.2% in ethanol	571
E476915	Methyl red solution water/ ethanol 0.2%	503
E476921	Methyl red solution 0.1% in ethanol	504
E477301	Starch paste solution 1%	721
E477302	Starch paste solution 1%	721
E477507	Soap solution in ethanol	660
E485041	Tin (II) chloride solution 10%	761
E485952	Sudan III hydroalcoholic saturated solution	727
E487031	Tetrabutylammonium hydroxide 0.1 mol/l (0.1N) in 2-propanol	749
E487755	Thymolphthalein 0.1% hydroalcoholic solution	759
E490056	Tropaeolin O solution 0.1%	779
E490451	Turk's reagent	780
E491255	Bromocresol green 0.04% hydroalcoholic solution	211
E491551	Crystal violet solution 0.5% in anhydrous acetic acid	284
E491651	Gentian violet carbolated solution	368
E491661	Gentian violet carbolated solution	368
E491901	Wijs' reagent	790
E491902	Wijs' reagent	790
E494301	Zinc chloride solution 60%	799
E497401	Aluminum standard solution	152
E497405	Aluminum standard solution	152
E497411	Antimony standard solution	183
E497415	Antimony standard solution	183
E497441	Barium standard solution	191
E497445	Barium standard solution	191
E497451	Bismuth standard solution	202
E497455	Bismuth standard solution	202
E497461	Boron standard solution	208
E497465	Boron standard solution	208
E497481	Calcium standard solution	235
E497485	Calcium standard solution	235
E497491	Cobalt standard solution	271
E497495	Cobalt standard solution	271
E497501	Chromium standard solution	265
E497511	Iron standard solution	429
E497515	Iron standard solution	429
E497521	Lithium standard solution	462
E497525	Lithium standard solution	462
E497531	Magnesium standard solution	469
E497535	Magnesium standard solution	469
E497541	Manganese standard solution	479
E497545	Manganese standard solution	479
E497561	Molybdenum standard solution	511
E497565	Molybdenum standard solution	511
E497571	Nickel standard solution	522
E497575	Nickel standard solution	522
E497581	Gold standard solution	375
E497585	Gold standard solution	375
E497591	Lead standard solution	456
E497595	Lead standard solution	456
E497601	Potassium standard solution	582
E497605	Potassium standard solution	582
E497611	Copper standard solution	276
E497615	Copper standard solution	276
E497621	Selenium standard solution	647
E497625	Selenium standard solution	647
E497631	Silicon standard solution	651
E497635	Silicon standard solution	651
E497641	Sodium standard solution	661
E497645	Sodium standard solution	661
E497651	Tin standard solution	760
E497655	Tin standard solution	760
E497661	Strontium standard solution	723
E497665	Strontium standard solution	723
E497671	Vanadium standard solution	783
E497675	Vanadium standard solution	783
E497681	Zinc standard solution	798
E497685	Zinc standard solution	798
FG201138	Ausilab 210	187
FG201139	Ausilab 250	188
FG201140	Ausilab 500	189
FG201146	Ausilab 400	189
FG201149	Ausilab 110	187
FG201151	Ausilab 280	188
FG201200	Ausilab 290	188
FG201220	Ausilab 140	187
FG20C1156C5	Ausilab 260	188
P0011016	n-Butyl acetate	228
P0011021	n-Butyl acetate	228
P0021010	Ethyl acetate	341
P0021016	Ethyl acetate	341
P0021021	Ethyl acetate	341
P00210T21	Ethyl acetate	341
P0023216	Ethyl acetate	341
P0023221	Ethyl acetate	341
P0040228	Methyl acetate	496
P0040240	Methyl acetate	496
P0040248	Methyl acetate	496
P0040268	Methyl acetate	496
P0043721	Methyl acetate	495
P0051010	Acetone	139
P0051016	Acetone	139
P0053216	Acetone	138
P0053221	Acetone	138
P0053282	Acetone	138
P0060228	Acetonitrile	144
P0060248	Acetonitrile	144
P0060268	Acetonitrile	144
P0061010	Acetonitrile	144
P0061016	Acetonitrile	144
P0061021	Acetonitrile	144
P00610S10	Acetonitrile	144
P00610S16	Acetonitrile	144
P00610S21	Acetonitrile	144
P00610T10	Acetonitrile	144
P0063510	Acetonitrile	144
P0063516	Acetonitrile	144
P0063521	Acetonitrile	144
P00637S16	Acetonitrile	143
P00637S21	Acetonitrile	143
P00725P15	Acetic acid glacial	132
P00725P21	Acetic acid glacial	132
P0080212	Trifluoroacetic acid	774
P0080247	Trifluoroacetic acid	774
P0080297	Trifluoroacetic acid	774
P0082103	Trifluoroacetic acid	773
P0082112	Trifluoroacetic acid	773
P0082147	Trifluoroacetic acid	773
P0082746	Trifluoroacetic acid	773
P0082747	Trifluoroacetic acid	773
P0093016	Decane	290
P0120041	Triton© X100 solution	778
P013A1010	Ethanol absolute anhydrous	331
P013A1016	Ethanol absolute anhydrous	331
P013A1021	Ethanol absolute anhydrous	331
P0171010	Butanol-1	226
P0171016	Butanol-1	226
P0180241	Butanol-2	227
P0190222	tert-Butanol	227
P0190268	tert-Butanol	227
P0191016	tert-Butanol	227
P0201016	Ethyl methyl ketone	350
P02402E40	Chloroform	262
P02402E48	Chloroform	262
P02405A28	Chloroform	262
P02410A10	Chloroform	261
P02410A16	Chloroform	261
P02410A21	Chloroform	261
P02410AT21	Chloroform	261
P02410E10	Chloroform	261
P02410E16	Chloroform	261
P02410E21	Chloroform	261
P02432E16	Chloroform	260
P02432E21	Chloroform	260
P0251010	Cyclohexane	286
P0251016	Cyclohexane	286
P0260248	Cyclohexanol	288
P0260268	Cyclohexanol	288
P0280228	1,2-Dichloroethane	300
P0280268	1,2-Dichloroethane	300
P0281010	1,2-Dichloroethane	299
P0281016	1,2-Dichloroethane	299
P0281021	1,2-Dichloroethane	299
P0282716	1,2-Dichloroethane	299
P02910A10	Dichloromethane	302
P02910A16	Dichloromethane	302
P02910A21	Dichloromethane	302
P02910A48	Dichloromethane	302
P02910AT10	Dichloromethane	302
P02910AT16	Dichloromethane	302
P02910E10	Dichloromethane	302
P02910E16	Dichloromethane	302
P02910E21	Dichloromethane	302
P02927A21	Dichloromethane	302
P02932A16	Dichloromethane	301
P02932A21	Dichloromethane	301
P02932A82	Dichloromethane	301
P02932E16	Dichloromethane	301
P02932E21	Dichloromethane	301
P0300221	1,2-Dimethoxyethane	313
P0300268	1,2-Dimethoxyethane	313
P0301010	1,2-Dimethoxyethane	313
P03010T16	1,2-Dimethoxyethane	313
P0341010	N,N-Dimethylformamide	315
P0341016	N,N-Dimethylformamide	315
P0341021	N,N-Dimethylformamide	315
P03410T10	N,N-Dimethylformamide	315
P03410T16	N,N-Dimethylformamide	315
P0343216	N,N-Dimethylformamide	315
P0343221	N,N-Dimethylformamide	315
P0343516	N,N-Dimethylformamide	316
P0343521	N,N-Dimethylformamide	316
P0343522	N,N-Dimethylformamide	316
P0343541	N,N-Dimethylformamide	316
P0343549	N,N-Dimethylformamide	316
P0343550	N,N-Dimethylformamide	316
P0343567	N,N-Dimethylformamide	316
P0350239	Dimethylsulphoxide	319
P0350242	Dimethylsulphoxide	319
P0350266	Dimethylsulphoxide	319
P03502T10	Dimethylsulphoxide	319
P03502T16	Dimethylsulphoxide	319
P03502T21	Dimethylsulphoxide	319
P0353216	Dimethylsulphoxide	318

Product specifications are subject to changes. Please visit our website for updates.

INDEX BY PRODUCT CODE

P0353221	Dimethylsulphoxide	318
P0355016	Dimethylsulphoxide	318
P0355041	Dimethylsulphoxide	318
P0355049	Dimethylsulphoxide	318
P0361010	1,4-Dioxane	321
P0361016	1,4-Dioxane	321
P0361021	1,4-Dioxane	321
P0370048	Petroleum benzin E	563
P0400252	N,N-Diisopropylethylamine	312
P0400272	N,N-Diisopropylethylamine	312
P0410228	Diethylene glycol dimethyl ether	307
P0410248	Diethylene glycol dimethyl ether	307
P0430228	Diisopropylether	312
P0430240	Diisopropylether	312
P0430248	Diisopropylether	312
P0430268	Diisopropylether	312
P0431016	Diisopropylether	312
P0440508	Diethyl ether	309
P0441008	Diethyl ether	309
P0441010	Diethyl ether	309
P0441016	Diethyl ether	309
P0441021	Diethyl ether	309
P04410T10	Diethyl ether	309
P04410T16	Diethyl ether	309
P0501016	n-Heptane 99%	381
P0501021	n-Heptane 99%	381
P0502721	n-Heptane 99%	380
P0520040	Hexane mixture of isomers	388
P052053068	n-Hexane 99%	386
P0521016	n-Hexane	387
P05230S01/16	n-Hexane 99%	385
P05230S01/21	n-Hexane 99%	385
P0523216	n-Hexane	386
P0523221	n-Hexane	386
P052323016	n-Hexane 99%	385
P052323021	n-Hexane 99%	385
P0531016	Isobutanol	437
P0580516	Methylcyclohexane	498
P0581016	Methylcyclohexane	497
P0582716	Methylcyclohexane	497
P0600268	Methyl isobutyl ketone	501
P0601016	Methyl isobutyl ketone	500
P0630240	Isooctane	440
P0630268	Isooctane	440
P0641016	n-Pentane	559
P0643216	n-Pentane	558
P0643221	n-Pentane	558
P064323016	n-Pentane 99%	557
P064323021	n-Pentane 99%	557
P0643716	n-Pentane	558
P0643721	n-Pentane	558
P0651016	Isopentane	441
P0660216	Piperidine	580
P0660221	Piperidine	580
P0660229	Piperidine	580
P0660248	Piperidine	580
P0663516	Piperidine	580
P0663518	Piperidine	580
P0663521	Piperidine	580
P0671010	Pyridine	630
P0671016	Pyridine	630
P0671021	Pyridine	630
P0671046	Pyridine	630
P0671068	Pyridine	630
P06725P16	Pyridine	630
P0673516	Pyridine	630
P0673521	Pyridine	630
P0680228	Tetrachloroethylene	750
P0682716	Tetrachloroethylene	750
P0682721	Tetrachloroethylene	750
P0701010	Tetrahydrofuran	752
P0701016	Tetrahydrofuran	752
P0701021	Tetrahydrofuran	752
P07010T10	Tetrahydrofuran	752
P07010T16	Tetrahydrofuran	752
P0711010	Toluene	766
P0711016	Toluene	766
P0711021	Toluene	766
P07110T10	Toluene	766
P0713216	Toluene	766
P0713221	Toluene	766
P0713282	Toluene	766
P0722721	1,2,4-Trichlorobenzene	771
P0790266	Triethylamine	772
P0853016	Hexadecane	383
P0871010	N-Methyl-2-pyrrolidone	502
P0873516	N-Methyl-2-pyrrolidone	502
P0873521	N-Methyl-2-pyrrolidone	502
P0873541	N-Methyl-2-pyrrolidone	502
P0873549	N-Methyl-2-pyrrolidone	502
P0873566	N-Methyl-2-pyrrolidone	502
P0883216	Petroleum ether 35 - 60°C	567
P0883221	Petroleum ether 35 - 60°C	567
P0890228	Isopropyl acetate	442
P0890240	Isopropyl acetate	442
P0890268	Isopropyl acetate	442
P0900221	Methylisoamyl ketone	500
P0921016	tert-Butylmethylether	231
P0931010	Methanol	490
P0931016	Methanol	490
P0931021	Methanol	490
P0933216	Methanol	489
P0933221	Methanol	489
P0941016	Propan-1-ol	622
P0941021	Propan-1-ol	622
P0951010	Propan-2-ol	624
P0951016	Propan-2-ol	624
P0952721	Propan-2-ol	624
P1050552	Sodium hydroxide solution 30%	682
P1320017	Sodium sulfate anhydrous	706
P1320027	Sodium sulfate anhydrous	706
P1320044	Sodium sulfate anhydrous	706
P1460012	Magnesium sulfate anhydrous	475
P1460027	Magnesium sulfate anhydrous	475
P1460044	Magnesium sulfate anhydrous	475
P1460057	Magnesium sulfate anhydrous	475
P1650517	Sodium persulfate	699
P1810017	Molecular sieves 3 A	510
P1820017	Molecular sieves 4 A	510
P1820027	Molecular sieves 4 A	510
P1820047	Molecular sieves 4 A	510
P2000017	Silica gel 60A 35 - 70µ	649
P2000026	Silica gel 60A 35 - 70µ	649
P2000027	Silica gel 60A 35 - 70µ	649
P2000044	Silica gel 60A 35 - 70µ	649
P2010017	Silica gel 60A 6 - 35µ	648
P2010027	Silica gel 60A 6 - 35µ	648
P2010044	Silica gel 60A 6 - 35µ	648
P2050017	Silica gel 60A 40 - 63µ	649
P2050027	Silica gel 60A 40 - 63µ	649
P2050044	Silica gel 60A 40 - 63µ	649
P2100017	Silica gel 60A 70 - 200µ	650
P2100026	Silica gel 60A 70 - 200µ	650
P2100027	Silica gel 60A 70 - 200µ	650
P2100044	Silica gel 60A 70 - 200µ	650
P2200017	Silica gel 60A 20 - 45µ	648
P2200027	Silica gel 60A 20 - 45µ	648
P3100015	Acetic acid 0.1 mol/l (0.1N)	136
P3160015	Hydrochloric acid 5 mol/l (5N)	402
P3160095	Hydrochloric acid 5 mol/l (5N)	402
P3240015	Sulfuric acid 2.5 mol/l (5N)	737
P3250016	Iron (II) ammonium sulfate 0.1N	430
P3440015	Sodium hydroxide 0.2 mol/l (0.2N)	689
P3520022	Sodium thiosulfate 0.2 mol/l (0.2N)	712
P3530015	Sodium thiosulfate 0.5 mol/l (0.5N)	712
P3840016	Tetrabutylammonium hydroxide 0.1 mol/l (0.1N) in methanol / propanol-2 (50/50)	750
P4500022	Sodium hydroxide 1/9 mol/l (N/9)	689
P4540041	Sodium hydroxide solution 50%	680
P4540049	Sodium hydroxide solution 50%	680
P4610017	Charcoal activated	255
P5006	Chloroform-d + 0.03% TMS	263
P5039	Acetic acid-d4	137
P5044A	Acetone-d6	141
P5045	Acetone-d6	141
P5046	Acetone-d6	141
P5049	Acetone-d6	141
P5055	Orthophosphoric acid-d3 85% in D2O	549
P5060	Acetone-d6	141
P5070	Acetonitrile-d3	145
P5073A	Acetonitrile-d3	145
P5079	Acetonitrile-d3	145
P5085	Benzene-d6	198
P5086	Benzene-d6	198
P5089	Benzene-d6	198
P5115	Chloroform-d	263
P5116	Chloroform-d	263
P5117	Chloroform-d	263
P5118	Chloroform-d	263
P5119	Chloroform-d	263
P5130	Chloroform-d	263
P5140	Methanol-d4 + 0.03% TMS	494
P5151A	Cyclohexane-d12	287
P5160T	Deuterium oxide-d2 + 0.03% TSP d4	294
P5161T	Deuterium oxide-d2 + 0.5% TSP d4	294
P5164	Deuterium oxide-d2	293
P5165	Deuterium oxide-d2	293
P5165S	Deuterium oxide-d2	293
P5166	Deuterium oxide-d2	293
P5168	Deuterium oxide-d2	293
P5169	Deuterium oxide-d2	293
P5170D	Deuterium oxide-d2 + 0.01% DMSO	294
P5175	Deuterium oxide-d2	293
P5179	Deuterium oxide-d2	293
P5189A	N,N-Dimethylformamide-d7	316
P5200	Dimethylsulphoxide-d6	319
P5204A	Dimethylsulphoxide-d6	319
P5204S	Dimethylsulphoxide-d6	319
P5205	Dimethylsulphoxide-d6	319
P5206	Dimethylsulphoxide-d6	319
P5209	Dimethylsulphoxide-d6	319
P5220	Dimethylsulphoxide-d6	319
P5229	Dimethylsulphoxide-d6	319
P5262A	Ethanol-d6 anhydrous	340
P5275	Methanol-d1	494
P5280	Methanol-d4	493
P5283A	Methanol-d4	493
P5284	Methanol-d4	493
P5284S	Methanol-d4	493
P5285	Methanol-d4	493
P5289	Methanol-d4	493
P5309	Methanol-d3	494
P5310	Methanol-d4	493
P5319	Methanol-d4	493
P5325	Chloroform-d	263
P5330	Dichloromethane-d2	305
P5335	Dichloromethane-d2	305
P5339	Dichloromethane-d2	305
P5364A	Pyridine-d5	631
P5369A	Pyridine-d5	631
P5370	Pyridine-d5	631

INDEX BY PRODUCT CODE

P5380	Tetrahydrofuran-d8	753
P5385	Tetrahydrofuran-d8	753
P5393A	Toluene-d8	768
P5395	Toluene-d8	768
P5399A	Toluene-d8	768
P5413A	Trifluoroacetic acid-d	774
P5419A	Trifluoroacetic acid-d	774
P5435	Tetrachloroethane-d2	751
P5505	Chloroform-d	263
P5533A	1,2-Dichlorobenzene-d4	299
P5541	Dimethylsulphoxide-d6 + 0.03% TMS	319
P5545	Dimethylsulphoxide-d6 + 0.03% TMS	319
P5602	Dimethylsulphoxide-d6 + 0.03% TMS	320
P5605	Dimethylsulphoxide-d6 + 0.03% TMS	320
P5665	Sodium hydroxide-d 1 mol/l	691
P5675	Sodium hydroxide-d 1 30%	691
P5685	Hydrochloric acid-d 20%	409
P5695	Hydrochloric acid-d 1 mol/l	409
P5733	Formic acid-d	364
P6070268	Isopropyl myristate	442
P6080503	Hexafluoro-2-propanol	384
P6080518	Hexafluoro-2-propanol	384
P6120028	Isoamyl acetate	436
P6151010	Formamide	362
P6170541	Glycerol (30°Bé)	372
P6170567	Glycerol (30°Bé)	372
P6263216	Isohexane	438
P6263221	Isohexane	438
P7600513	Sodium laurylsulfate	694
P7600514	Sodium laurylsulfate	694
P7600517	Sodium laurylsulfate	694
P8010216	Cyclopentyl methyl ether	289
P8010229	Cyclopentyl methyl ether	289
P8010248	Cyclopentyl methyl ether	289
P8010268	Cyclopentyl methyl ether	289
P8020216	N,N'-Dimethylpropylene uree	317
P8020218	N,N'-Dimethylpropylene uree	317
P8020229	N,N'-Dimethylpropylene uree	317
P8020248	N,N'-Dimethylpropylene uree	317
P8020268	N,N'-Dimethylpropylene uree	317
P8030216	1,3-Dioxolane	322
P8030222	1,3-Dioxolane	322
P8030249	1,3-Dioxolane	322
P8030268	1,3-Dioxolane	322
P8040216	1,3-Propanediol	627
P8040222	1,3-Propanediol	627
P8040268	1,3-Propanediol	627
P8880014	Dicalite 4158	297
P8880017	Dicalite 4158	297
P8880027	Dicalite 4158	297
P932SP16	Sulfolane	729
P9350015	Sodium hypochlorite solution 12.5%	692
P9350046	Sodium hypochlorite solution 12.5%	692
P9350049	Sodium hypochlorite solution 12.5%	692
P9960216	2-Methyltetrahydrofuran	505
P9960221	2-Methyltetrahydrofuran	505
P9960229	2-Methyltetrahydrofuran	505
P9960248	2-Methyltetrahydrofuran	505
P9960268	2-Methyltetrahydrofuran	505
P9990216	4-Methyltetrahydropyran	506
P9990218	4-Methyltetrahydropyran	506
P9990221	4-Methyltetrahydropyran	506
PA160021	m-Cresol	282
PS0009/15	Sulfuric acid 30%	736
PS0016/96	Sulfuric acid 0.025 mol/l (0.05N)	741
PS0026/95	Sulfuric acid 0.005 mol/l (0.01N)	742
PS0030/15	Silver nitrate 0.01 mol/l (0.01N)	658
PS0047/15	Sulfuric acid 0.01 mol/l (0.02N)	741
PS0083/42	Sodium persulfate 1 mol/l	700
PS0084/22	Orthophosphoric acid 10%	549
PS0084/42	Orthophosphoric acid 10%	549
PS0158/15	Rhodol solution	639
PS0158/29	Rhodol solution	639
PS0187/15	Thymol blue TA indicator	758
PS0187/16	Thymol blue TA indicator	758
PS0189/15	Bromophenol blue TAC indicator	214
PS0189/16	Bromophenol blue TAC indicator	214
PS0194/22	Ammonia buffer solution pH 10	165
PS0194/95	Ammonia buffer solution pH 10	165
PS0200/15	Buffer pH 10	223
PS0206/41	Hydrochloric acid 0.04 mol/l (0.04N)	408
PS0212/21	Sulfuric acid 25%	736
PS0215/15	Sodium hydroxide 0.01 mol/l (0.01N)	690
PS0217/15	Sulfuric acid 0.166 mol/l (0.333N)	739
PS0219/15	Sulfuric acid 0.02 mol/l (0.04N)	741
PS0219/95	Sulfuric acid 0.02 mol/l (0.04N)	741
PS0219/96	Sulfuric acid 0.02 mol/l (0.04N)	741
PS0221/29	Acetic acid 12%	136
PS0222/52	Acetic acid 25%	135
PS0226/41	Boric buffer solution	207
PS0237/41	Acetic acid 20%	135
PS0250/16	Silver nitrate 0.1 mol/l (0.1N) in 2-propanol	657
PS0252/16	Silver nitrate 0.01 mol/l (0.01N) in propanol-2	658
PS0269/15	Bromophenol blue indicator	214
PS0270/15	Thymol blue indicator	758
PS0311/20	Nitric acid 8 mol/l (8N)	532
PS0313/20	Hydrochloric Acid 9 mol/l (9N)	401
PS0327/21	Reagent TAN	635
PS0327/29	Reagent TAN	635
PS0327/39	Reagent TAN	635
PS0342/15	Hydrochloric acid 0.02 mol/l (0.02N)	408
PS0347/22	Hydrochloric acid 12%	399
PS0347/49	Hydrochloric acid 12%	399
PS0347/66	Hydrochloric acid 12%	399
PS0363/21	Tisab	762
PS0423/21	Reagent TBN ASTM D2896	636
PS0423/29	Reagent TBN ASTM D2896	636
PS0423/39	Reagent TBN ASTM D2896	636
PS0427/19	Buffer pH 9	223
PS0433/15	Sulfuric acid 85%	734
PS0445/22	Sulfuric acid 0.125 mol/l (0.25N)	740
PS0492/41	Fehling's A reagent	353
PS0493/41	Fehling's B reagent	353
PS0562/22	Boric acid 20 g/l with indicator	207
PS0563/21	Boric acid 3%	206
PS0568/41	Nitric acid 10%	532
PS0573/21	Mixture for bromine index determination	508
PS0587/15	Hydrochloric acid 0.05 mol/l (0.05N)	407
PS0589/15	Hydrochloric acid 4 mol/l (4N)	402
PS0589/22	Hydrochloric acid 4 mol/l (4N)	402
PS0589/49	Hydrochloric acid 4 mol/l (4N)	402
PS0703/22	Boric acid 20g/l	207
PS0736/41	Sodium hydroxide 1.2 mol/l (1.2N)	686
PS0736/42	Sodium hydroxide 1.2 mol/l (1.2N)	686
PS0736/49	Sodium hydroxide 1.2 mol/l (1.2N)	686
PS0740/95	Phosphate buffer pH 7.4	575
PS0751/29	Hydrochloric acid 20%	399
PS0766/22	Potassium hydroxide solution 33%	601
PS0768/41	Hydrochloric acid 10%	400
PS0769/20	Hydrochloric acid 26%	398
PS0772/79	Potassium chloride 1 mol/l (1N)	588
PS0784/95	Buffer acetate pH 4.5	216
PS0809/22	Tisab	763
PS0852/29	Acetone / water 98/2 (v/v) with bromophenol blue 0.02 g/l	141
PS0864/41	Hydrochloric acid 5%	400
PS0893/21	Sulfuric acid 69%	735
PS0894/21	Sulfuric acid 62%	735

Product specifications are subject to changes. Please visit our website for updates.

INDEX BY CAS NUMBER

50-00-0	Formaldehyde 37% w/v	358
50-00-0	Formaldehyde 37% w/v neutralized	359
50-00-0	Formaldehyde 35% w/w	359
50-00-0	Formaldehyde 30% w/v	360
50-00-0	Formaldehyde 10% v/v according to Lillie	360
50-00-0	Formaldehyde 5% w/v buffered at pH 6.9	360
50-00-0	Formaldehyde 4% w/v buffered at pH 6.9	361
50-00-0	Formaldehyde 4% w/v with sodium chloride	361
50-00-0	Formaldehyde acetic	362
50-69-1	D(-)Ribose	639
50-70-4	Sorbitol	714
50-70-4	Sorbitol (no crystallizable) solution 70%	715
50-70-4	D-Sorbitol	715
50-81-7	L(+)Ascorbic acid	186
50-99-7	D(+)Glucose anhydrous	370
52-90-4	L-Cysteine	289
54-21-7	Sodium salicylate	705
55-55-0	4-Methylaminophenol sulfate	496
56-40-6	Glycine	373
56-81-5	Glycerol (30°Bé)	372
56-81-5	Glycerol 90% (28° Bé)	373
56-84-8	L(+)Aspartic acid	186
56-86-0	L(+)Glutamic acid	371
56-89-3	L-Cystine	289
56-92-8	Histamine dihydrochloride	389
57-09-0	Trimethylcetylammonium bromide	775
57-11-4	Stearic acid	722
57-13-6	Urea	782
57-48-7	D(-)Fructose	364
57-50-1	D(+)-Sucrose	726
57-55-6	Propylene glycol	629
58-08-2	Caffeine anhydrous	234
58-15-1	Aminopyrazolone solution	159
58-86-6	D(+)Xylose	793
59-23-4	D(+)Galactose	366
59-51-8	DL-Methionine	494
59-67-6	Nicotinic acid	525
59-88-1	Phenylhydrazine hydrochloride	573
59-88-1	Phenylhydrazine hydrochloride solution	573
60-00-4	Ethylenediaminetetraacetic acid	344
60-10-6	Diphenylthiocarbazone	323
60-11-7	Methyl yellow	506
60-12-8	2-Phenylethanol	573
60-18-4	L(-)Tyrosine	781
60-24-2	2-Mercaptoethanol	483
60-29-7	Diethyl ether	308
60-35-5	Acetamide	129
61-73-4	Methylene blue saturated solution	498
61-73-4	Methylene blue solution 1%	499
61-90-5	L(+)-Leucine	461
62-44-2	p-Acetylphenetidine	146
62-54-4	Calcium acetate anhydrous	236
62-55-5	Thioacetamide	754
62-55-5	Thioacetamide solution 40 g/l	755
62-56-6	Thiourea	756
62-76-0	Sodium oxalate	698
63-74-1	Sulfanilamide	728
63-91-2	L-Phenylalanine	572
64-17-5	Ethanol absolute anhydrous	331
64-17-5	Ethanol 96°	334
64-17-5	Ethanol 70% v/v	336
64-17-5	Ethanol 60% v/v	337
64-17-5	Ethanol 50% v/v	337
64-17-5	Ethanol absolute denaturated	338
64-17-5	Ethanol 95° denaturated	338
64-17-5	Ethanol 94° denaturated	339
64-17-5	Ethanol 90° denaturated	339
64-17-5	Ethanol 70°modified	339
64-18-6	Formic acid 99%	362
64-18-6	Formic acid 85%	363
64-19-7	Acetic acid glacial	130
64-19-7	Acetic acid 96%	133
64-19-7	Acetic acid 80%	134
64-19-7	Acetic acid 45%	134
64-19-7	Acetic acid 30%	134
64-19-7	Acetic acid 27%	135
64-19-7	Acetic acid 25%	135
64-19-7	Acetic acid 20%	135
64-19-7	Acetic acid 12%	135
64-19-7	Acetic acid 1 mol/l (1N)	136
64-19-7	Acetic acid 0.1 mol/l (0.1N)	136
64-19-7	Acetic acid 0.03 mol/l (0.03N)	136
65-61-2	Acridine orange	147
65-85-0	Benzoic acid	199
67-42-5	Ethylene glycol bis(2-aminoethyl ether)-n,n,n',n'-tetraacetic acid	349
67-43-6	Diethylenetriaminepentacetic acid	308
67-52-7	Barbituric acid	191
67-56-1	Methanol	488
67-56-1	Methanol + 0.1% v/v formic acid	492
67-56-1	Methanol + 0.1% v/v trifluoroacetic acid	493
67-56-1	Methanol, hydrochloric	493
67-63-0	Propan-2-ol	623
67-63-0	Propan-2-ol 70%	626
67-64-1	Acetone	138
67-66-3	Chloroform	259
67-68-5	Dimethylsulphoxide	317
68-04-2	Sodium citrate tribasic anhydrous	672
68-11-1	Thioglycolic acid 80%	755
68-12-2	N,N-Dimethylformamide	314
69-65-8	D-Mannitol	481
69-72-7	Salicylic acid	642
70-55-3	p-Toluene sulfonamide	768
71-23-8	Propan-1-ol	622
71-36-3	Butanol-1	226
71-41-0	n-Amyl alcohol	180
71-43-2	Benzene	198
71-91-0	Tetraethylammonium bromide	751
72-48-0	Alizarin	149
72-48-0	Alizarin saturated solution in ethanol	149
74-79-3	L(+)Arginine	184
74-88-4	Methyl iodide	500
75-05-8	Acetonitrile	142
75-05-8	Acetonitrile + 0.1% v/v formic acid	145
75-05-8	Acetonitrile + 0.1% v/v trifluoroacetic acid	145
75-09-2	Dichloromethane	300
75-09-2	Dichloromethane acidified with 1% hydrochloric acid	305
75-12-7	Formamide	362
75-31-0	Isopropylamine	442
75-47-8	Iodoform	427
75-59-2	Tetramethylammonium hydroxide 10%	753
75-65-0	tert-Butanol	227
75-75-2	Methanesulfonic acid	487
75-85-4	tert-Amyl alcohol	180
76-03-9	Trichloroacetic acid	770
76-03-9	Trichloroacetic acid solution 20%	770
76-05-1	Trifluoroacetic acid	773
76-22-2	Camphor synthetic	246
76-59-5	Bromothymol blue	215
76-59-5	Bromothymol blue 0.4% in ethanol	215
76-59-5	Bromothymol blue 0.02%	215
76-59-5	Bromothymol blue solution	215
76-60-8	Bromocresol green	211
76-60-8	Bromocresol green 0.04% hydroalcoholic solution	211
76-60-8	Bromocresol green solution	211
76-61-9	Thymol blue	757
76-61-9	Thymol blue 0.4% in ethanol	757
77-09-8	Phenolphthalein	571
77-86-1	Tris (hydroxymethyl)-aminomethane	776
77-92-9	Citric acid anhydrous	269
78-50-2	Trioctylphosphine oxide	775
78-78-4	Isopentane	441
78-83-1	Isobutanol	437
78-92-2	Butanol-2	227
78-93-3	Ethyl methyl ketone	350
79-07-2	Chloroacetamide	258
79-09-4	Propionic acid	628
79-11-8	Chloroacetic acid	258
79-14-1	Glycolic acid	374
79-20-9	Methyl acetate	495
79-33-4	L(+)-Lactic acid	452
79-40-3	Dithiooxamide	324
79-43-6	Dichloroacetic acid	297
83-07-8	4-Aminophenazone	159
83-88-5	Riboflavine	639
84-66-2	Diethyl phthalate	310
84-74-2	Di-n-butylphthalate	297
85-44-9	Phthalic anhydride	579
85-85-8	1-(2-Pyridylazo)-2-naphthol	632
85-86-9	Sudan III	726
85-86-9	Sudan III hydroalcoholic saturated solution	727
87-66-1	Pyrogallol	632
87-69-4	L(+) Tartaric Acid	745
87-69-4	L(+) Tartaric acid solution 20% in water	746
87-88-7	Chloranilic acid	257
87-89-8	Inositol	423
88-89-1	Picric acid solution	579
88-99-3	Phthalic acid	579
89-83-8	Thymol	757
90-11-9	alpha-Bromonaphthalene	213
90-15-3	1-Naphthol	517
90-44-8	Anthrone	182
91-20-3	Naphthalene	517
91-22-5	Quinoline	634
92-31-9	Toluidine blue	769
93-58-3	Methyl benzoate	496
94-13-3	Propyl p-hydroxybenzoate	629
94-36-0	Dibenzoyl peroxide	296
95-45-4	Diacetyldioxime	295
95-45-4	Diacetyldioxime solution 1% in ethanol	296
95-47-6	o-Xylene	792
95-48-7	o-Cresol	282
95-50-1	o-Dichlorobenzene	298
96-47-9	2-Methyltetrahydrofuran	505
97-00-7	2,4-Dinitrochlorobenzene	320
97-53-0	Eugenol	351
98-09-9	Benzenesulfonyl chloride	198
98-86-2	Acetophenone	146
98-92-0	Nicotinamide	525
98-95-3	Nitrobenzene	535
99-61-6	m-Nitrobenzaldehyde	534
99-76-3	Methyl 4-hydroxybenzoate	499
100-02-7	p-Nitrophenol	535
100-09-4	Anisic acid	181
100-10-7	p-Dimethylaminobenzaldehyde	314

INDEX BY CAS NUMBER

100-51-6	Benzyl alcohol	200
100-52-7	Benzaldehyde	197
100-97-0	Hexamethylenetetramine	384
101-38-2	2,6-Dichloroquinone-4-chlorimide	305
102-71-6	Triethanolamine	771
102-76-1	Triacetin	769
103-84-4	Acetanilide	129
103-84-4	Melting point standards	483
104-88-1	p-Chlorobenzaldehyde	258
105-57-7	Acetal	129
105-58-8	Diethyl carbonate	307
105-67-9	2,4-Xylenol	792
106-46-7	p-Dichlorobenzene	298
106-48-9	4-Chlorophenol	264
106-51-4	p-Benzoquinone	199
107-06-2	1,2-Dichloroethane	299
107-15-3	Ethylenediamine	344
107-21-1	Ethylene glycol	348
107-83-5	Isohexane	438
107-92-6	n-Butyric acid	231
107-95-9	ß-Alanine	148
108-10-1	Methyl isobutyl ketone	500
108-18-9	Diisopropylamine	311
108-19-0	Biuret 97%	205
108-19-0	Biuret reagent	205
108-20-3	Diisopropylether	312
108-21-4	Isopropyl acetate	441
108-24-7	Acetic anhydride	137
108-30-5	Succinic anhydride	725
108-31-6	Maleic anhydride	477
108-32-7	Propylene carbonate	629
108-39-4	m-Cresol	282
108-46-3	Resorcinol	637
108-73-6	Phloroglucinol	574
108-87-2	Methylcyclohexane	497
108-88-3	Toluene	765
108-88-3	Toluene in solution in hexane	768
108-90-7	Chlorobenzene	258
108-91-8	Cyclohexylamine	288
108-93-0	Cyclohexanol	287
108-94-1	Cyclohexanone	288
108-95-2	Phenol	570
108-95-2	Phenol liquified 85%	570
109-57-9	Allylthiourea	150
109-60-4	n-Propyl acetate	628
109-66-0	n-Pentane 99%	557
109-66-0	n-Pentane	558
109-66-0	Petroleum ether 35 - 60°C	567
109-66-0	Petroleum ether 30 - 50°C	568
109-66-0	Petroleum ether 30 - 40°C	569
109-69-3	n-Butyl chloride	229
109-86-4	2-Methoxy ethanol	495
109-89-7	Diethylamine	306
109-94-4	Ethyl formate	350
109-99-9	Tetrahydrofuran	751
110-12-3	Methylisoamyl ketone	500
110-15-6	Succinic acid	725
110-16-7	Maleic acid	477
110-17-8	Fumaric acid	365
110-19-0	Isobutyl acetate	438
110-27-0	Isopropyl myristate	442
110-54-3	n-Hexane 99%	385
110-54-3	n-Hexane	386
110-54-3	Hexane mixture of isomers	388
110-71-4	1,2-Dimethoxyethane	313
110-82-7	Cyclohexane	286
110-86-1	Pyridine	630
110-89-4	Piperidine	580
110-91-8	Morpholine	512
111-20-6	Sebacic acid	646
111-30-8	Glutardialdehyde solution 50%	372
111-42-2	Diethanolamine	306
111-46-6	Diethylene glycol	307
111-65-9	n-Octane	537
111-65-9	Octane 80 blend	537
111-76-2	2-Butoxy ethanol	228
111-87-5	Octanol-1	538
111-96-6	Diethylene glycol dimethyl ether	307
112-27-6	Triethylene glycol	773
112-34-5	2-(2-Butoxyethoxy)ethanol	228
112-80-1	Oleic acid	539
115-20-8	2,2,2 - Trichlorethanol	770
115-39-9	Bromophenol blue	213
115-39-9	Bromophenol blue solution 0.4% in ethanol	213
115-39-9	Bromophenol blue solution 0.02%	214
115-39-9	Bromophenol blue solution	214
115-40-2	Bromocresol purple	212
115-40-2	Bromocresol purple solution 0.4% in ethanol	212
115-41-3	Pyrocatechol violet	632
116-63-2	1-Amino-2-naphthol-4-sulfonic acid	158
117-92-0	Quinaldine red	634
118-10-5	Cinchonine	268
118-75-2	Chloranil	256
119-26-6	2,4-Dinitrophenylhydrazine (with 30% of water)	320
119-36-8	Methyl salicylate	505
119-93-7	o-Tolidine solution 0.1%	765
119-93-7	o-Tolidine solution	765
120-46-7	Dibenzoylmethane	296
120-51-4	Benzyl benzoate	200
120-72-9	Indole	423
120-80-9	Pyrocatechol	632
120-82-1	1,2,4-Trichlorobenzene	771
121-33-5	Vanillin	784
121-44-8	Triethylamine	772
121-54-0	Hyamine 1622	391
121-54-0	Hyamine 1622 solution 0.004M	392
121-57-3	Sulfanilic acid	728
121-92-6	m-Nitrobenzoic acid	535
122-39-4	Diphenylamine solution 1% in sulfuric acid	322
123-08-0	p-Hydroxybenzaldehyde	416
123-11-5	Anisaldehyde	181
123-30-8	p-Aminophenol	159
123-31-9	Hydroquinone	416
123-38-6	Propionaldehyde	628
123-42-2	Diacetone alcohol	295
123-51-3	Isoamyl alcohol	436
123-54-6	Acetylacetone	146
123-86-4	n-Butyl acetate	228
123-91-1	1,4-Dioxane	321
123-92-2	Isoamyl acetate	436
124-04-9	Adipic acid	147
124-07-2	n-Caprylic acid	247
124-18-5	Decane	290
124-41-4	Sodium methoxide 0.1 mol/l	695
125-20-2	Thymolphthalein	758
125-20-2	Thymolphthalein 0.1% hydroalcoholic solution	759
125-20-2	Thymolphthalein solution 0.1% in ethanol	759
126-33-0	Sulfolane	729
126-81-8	Dimedone	313
127-08-2	Potassium acetate	583
127-09-3	Sodium acetate anhydrous	662
127-18-4	Tetrachloroethylene	750
127-19-5	N,N-Dimethylacetamide	313
128-37-0	Butylhydroxytoluene	230
130-22-3	Alizarin red	149
134-03-2	Sodium L-ascorbate	664
135-19-3	2-Naphthol	517
137-08-6	Calcium pantothenate	242
138-32-9	Trimethylcetylammonium p-toluenesulfonate	775
139-13-9	Nitrilotriacetic acid	534
140-22-7	sym-Diphenylcarbazide	323
141-43-5	Ethanolamine	340
141-53-7	Sodium formate	675
141-78-6	Ethyl acetate	340
141-82-2	Malonic acid	478
141-97-9	Ethyl acetoacetate	343
142-04-1	Aniline hydrochloride	181
142-47-2	Sodium glutamate acid	676
142-62-1	n-Caproic acid	247
142-82-5	n-Heptane 99%	380
142-82-5	n-Heptane	382
142-82-5	Heptane mixture of isomers	382
143-66-8	Sodium tetraphenylborate	710
143-74-8	Phenol red	570
144-33-2	Sodium citrate dibasic sesquihydrate	672
144-55-8	Eluent sodium bicarbonate	327
144-55-8	Sodium bicarbonate	665
144-62-7	Oxalic acid 0.5 mol/l (1N)	551
144-62-7	Oxalic acid 0.05 mol/l (0.1N)	551
144-62-7	Oxalic acid 0.005 mol/l (0.01N)	551
145-50-6	alpha-Naphtholbenzein	518
145-50-6	alpha-Naphtholbenzein solution 0.2% in acetic acid	518
147-85-3	L(-)Proline	622
148-24-3	8-Hydroxyquinoline	417
150-13-0	p-Aminobenzoic acid	158
151-21-3	Sodium laurylsulfate	693
288-32-4	Imidazole	419
298-14-6	Potassium bicarbonate	584
298-96-4	2,3,5-Triphenyltetrazolium chloride	776
301-04-2	Lead (II) acetate cotton	457
302-17-0	Chloral hydrate	256
333-20-0	Potassium thiocyanate	620
333-20-0	Potassium thiocyanate solution 5%	620
333-20-0	Potassium thiocyanate 0.1 mol/l (0.1N)	620
333-20-0	Potassium thiocyanate solution	621
366-18-7	2,2'-Dipyridyl	323
375-22-4	Heptafluorobutyric acid	380
407-25-0	Trifluoroacetic anhydride	774
409-21-2	Carborundum, granules	248
409-21-2	Carborundum, powder	248
464-49-3	Camphor natural	246
471-34-1	Calcium carbonate	236
477-73-6	Safranine T	642
477-73-6	Safranine T hydroalcoholic solution for Gram-Hucker Kit	642
485-47-2	Ninhydrin	526
493-52-7	Methyl red	503
493-52-7	Methyl red solution	504
497-19-8	Eluent sodium carbonate	327
497-19-8	Sodium carbonate anhydrous	667
504-17-6	2-Thiobarbituric acid	755
504-63-2	1,3-Propanediol	627
513-77-9	Barium carbonate	192
513-78-0	Cadmium carbonate	233
517-28-2	Hematoxylin	379
517-28-2	Hematoxylin solution according to Mayer	379
518-47-8	Fluorescein sodium salt	357
518-67-2	Dimidium bromide	320
521-31-3	Luminol	467
526-94-3	Sodium hydrogen tartrate monohydrate	678
527-07-1	Sodium gluconate	676
532-32-1	Sodium benzoate	665
532-82-1	Chrysoidine Y	268
534-16-7	Silver carbonate	654

Product specifications are subject to changes. Please visit our website for updates.

INDEX BY CAS NUMBER

536-17-4	p-Dimethylaminobenzylide-nerhodanine	314
538-62-5	sym-Diphenylcarbazone	323
540-69-2	Ammonium formate	171
540-72-7	Sodium sulfocyanate	709
540-84-1	Isooctane	439
541-73-1	m-Dichlorobenzene	298
543-80-6	Barium acetate	192
543-94-2	Strontium acetate	723
544-76-3	Hexadecane	383
546-68-9	Titanium isopropylate	764
547-57-9	Tropaeolin O	778
547-57-9	Tropaeolin O solution 0.1%	779
547-58-0	Methyl orange	501
547-58-0	Methyl Orange solution 0.1%	501
547-91-1	Ferron	354
548-24-3	Eosin B	328
548-62-9	Crystal violet	284
548-62-9	Crystal violet solution 0.5% in anhydrous acetic acid	284
548-62-9	Gentian violet	368
548-62-9	Gentian violet carbolated solution	368
553-24-2	Neutral red	521
554-13-2	Lithium carbonate	463
557-04-0	Magnesium stearate	474
557-05-1	Zinc stearate	801
563-63-3	Silver acetate	654
569-64-2	Malachite green	476
573-58-0	Congo red	273
573-58-0	Congo red solution	274
584-08-7	Potassium carbonate	586
587-98-4	Metanil yellow	487
591-27-5	m-Aminophenol	159
599-00-8	Trifluoroacetic acid-d	774
603-35-0	Triphenylphosphine	776
603-45-2	Rosolic acid	640
611-71-2	D(-)Mandelic acid	478
616-91-1	n-Acetyl-L-cysteine	146
617-48-1	DL-Malic acid	477
620-45-1	2,6-Dichlorophenolindophenol sodium salt	305
631-61-8	Ammonium acetate	166
632-99-5	Fuchsin basic	365
632-99-5	Fuchsin solution decolorised	365
633-03-4	Brilliant green	209
633-96-5	Orange II	539
637-01-4	N,N,N',N'-Tetramethyl-p-phenylenediamine dihydrochloride	753
643-79-8	o-Phthalaldehyde	579
646-06-0	1,3-Dioxolane	322
657-27-2	L(+)Lysine monohydrochloride	467
657-84-1	p-Toluenesulfonic acid sodium salt	769
666-52-4	Acetone-d6	141
811-98-3	Methanol-d4	493
811-98-3	Methanol-d4 + 0.03% TMS	494
842-07-9	Sudan yellow	727
846-70-8	Naphthol yellow S	518
860-22-0	Indigo carmine dried	422
860-22-0	Indigo carmine solution	422
865-34-9	Lithium methoxide 0.1 mol/l (0.1N)	464
865-49-6	Chloroform-d	263
865-49-6	Chloroform-d + 0.03% TMS	263
868-14-4	Potassium L-tartrate monobasic	619
872-50-4	N-Methyl-2-pyrrolidone	502
877-24-7	Potassium hydrogen phthalate	597
877-24-7	Potassium hydrogen phthalate 0.2 mol/l (0.2N)	598
877-24-7	Potassium hydrogen phthalate 0.1 mol/l (0.1N)	598
920-66-1	Hexafluoro-2-propanol	384
925-94-0	Formic acid-d	364
927-20-8	Magnesium glycerophosphate	471
999-97-3	Hexamethyldisilazane	384
1064-48-8	Amidoschwarz 10B solution	157
1066-33-7	Ammonium bicarbonate	167
1076-43-3	Benzene-d6	198
1111-78-0	Ammonium carbamate	168
1119-34-2	L(+)Arginine monohydrochloride	184
1119-94-4	Dodecyltrimethylammonium bromide	325
1149-16-2	Glyoxal-bis-(2-hydroxyanil)	374
1185-53-1	Tris (hydroxymethyl)-aminomethane hydrochloride	778
1185-57-5	Iron (III) ammonium citrate green	431
1185-57-5	Iron (III) ammonium citrate red	431
1186-52-3	Acetic acid-d4	137
1303-96-4	Sodium tetraborate decahydrate	710
1304-85-4	Bismuth(III) nitrate basic	203
1305-62-0	Calcium hydroxide	240
1305-78-8	Calcium oxide, lumps	242
1305-78-8	Calcium oxide, powder	242
1308-38-9	Chromium (III) oxide	266
1309-37-1	Iron (III) oxide	435
1309-42-8	Magnesium hydroxide	472
1309-48-4	Magnesium oxide	473
1309-48-4	Magnesium oxide heavy	473
1309-60-0	Lead (IV) oxide	460
1310-58-3	Potassium hydroxide, flakes	599
1310-58-3	Potassium hydroxide, pellets	599
1310-58-3	Potassium hydroxide solution 45%	600
1310-58-3	Potassium hydroxide solution 38% (40° Bé) in water	601
1310-58-3	Potassium hydroxide solution 33%	601
1310-58-3	Potassium hydroxide solution 28%	601
1310-58-3	Potassium hydroxide solution 3% in ethanol	601
1310-58-3	Potassium hydroxide 2 mol/l (2N) in ethanol	602
1310-58-3	Potassium hydroxide 1 mol/l (1N)	602
1310-58-3	Potassium hydroxide 0.5 mol/l (0.5N)	602
1310-58-3	Potassium hydroxide 0.5 mol/l (0.5N) in ethanol	603
1310-58-3	Potassium hydroxide 0.5 mol/l (0.5N) in methanol	603
1310-58-3	Potassium hydroxide 0.46 mol/l (0.46N)	604
1310-58-3	Potassium hydroxide 0.25 mol/l (0.25N)	604
1310-58-3	Potassium hydroxide 0.23 mol/l (0.23N)	604
1310-58-3	Potassium hydroxide 0.1 mol/l (0.1N)	604
1310-58-3	Potassium hydroxide 0.1 mol/l (0.1N) in ethanol	605
1310-58-3	Potassium hydroxide 0.1 mol/l (0.1N) in methanol	605
1310-58-3	Potassium hydroxide in solution	606
1310-65-2	Lithium hydroxide anhydrous	464
1310-66-3	Lithium hydroxide monohydrate	464
1310-73-2	Sodium hydroxide, pearls	678
1310-73-2	Sodium hydroxide, pellets	679
1310-73-2	Sodium hydroxide on silica	680
1310-73-2	Sodium hydroxide solution 50%	680
1310-73-2	Sodium hydroxide solution 40%	680
1310-73-2	Sodium hydroxide solution 35-37%	681
1310-73-2	Sodium hydroxide solution 35%	681
1310-73-2	Sodium hydroxide solution 32%	681
1310-73-2	Sodium hydroxide solution 30%	682
1310-73-2	Sodium hydroxide solution 20% w/v	683
1310-73-2	Sodium hydroxide solution 20% w/w	683
1310-73-2	Sodium hydroxide solution 10% w/v	683
1310-73-2	Sodium hydroxide solution 5% w/v	684
1310-73-2	Sodium hydroxide 6 mol/l (6N)	684
1310-73-2	Sodium hydroxide 5 mol/l (5N)	685
1310-73-2	Sodium hydroxide 4 mol/l (4N)	685
1310-73-2	Sodium hydroxide 3 mol/l (3N)	685
1310-73-2	Sodium hydroxide 2 mol/l (2N)	686
1310-73-2	Sodium hydroxide 1.2 mol/l (1.2N)	686
1310-73-2	Sodium hydroxide 1 mol/l (1N)	686
1310-73-2	Sodium hydroxide 0.7 mol/l (N/1.4)	687
1310-73-2	Sodium hydroxide 0.5 mol/l (0.5N)	688
1310-73-2	Sodium hydroxide 0.357 mol/l (0.357N)	688
1310-73-2	Sodium hydroxide 0.25 mol/l (0.25N)	688
1310-73-2	Sodium hydroxide 0.2 mol/l (0.2N)	689
1310-73-2	Sodium hydroxide 1/9 mol/l (N/9)	689
1310-73-2	Sodium hydroxide 0.1 mol/l (0.1N)	689
1310-73-2	Sodium hydroxide 0.1 mol/l (0.1N) in ethanol	690
1310-73-2	Sodium hydroxide 0.01 mol/l (0.01N)	690
1310-73-2	Sodium hydroxide solution	691
1310-73-2	Sodium hydroxide solution, methanolic	691
1312-81-8	Lanthanum oxide	455
1313-13-9	Manganese (IV) oxide	481
1313-27-5	Molybdenum (VI) oxide	511
1313-60-6	Sodium peroxide	699
1313-84-4	Sodium sulfide nonahydrate	707
1313-84-4	Sodium sulfide nonahydrate solution	707
1314-13-2	Zinc oxide	801
1314-23-4	Zirconium (IV) oxide	804
1314-35-8	Tungsten (VI) oxide	780
1314-56-3	Phosphorus pentoxide	578
1314-62-1	Vanadium (V) oxide	784
1314-98-3	Zinc sulfide	803
1317-36-8	Lead (II) oxide	460
1317-38-0	Copper (II) oxide	279
1317-38-0	Copper (II) oxide wire	279
1319-77-3	Tricresol	771
1320-06-5	Red for oils O	636
1321-14-8	Potassium guaiacolsulfonate	596
1324-76-1	Alkali blue 6B	150
1324-76-1	Alkali Blue 6B solution 2% in ethanol	150
1325-85-5	Victoria blue	785
1327-53-3	Arsenic trioxide solution	185
1330-43-4	Sodium tetraborate anhydrous	709
1332-09-8	Pumice stone	630
1332-58-7	Kaolin washed	443
1333-82-0	Chromium (VI) oxide	267
1333-82-0	Chromium (VI) oxide solution	267
1335-32-6	Lead (II) acetate basic	457
1336-21-6	Ammonia solution 32%	160
1336-21-6	Ammonia solution 30%	160

INDEX BY CAS NUMBER

1336-21-6	Ammonia solution 28%	161
1336-21-6	Ammonia solution 25%	162
1336-21-6	Ammonia solution 20 - 22%	163
1336-21-6	Ammonia solution 17%	164
1336-21-6	Ammonia solution 10%	164
1336-21-6	Ammonia solution 6N	164
1336-21-6	Ammonia solution diluted	164
1341-49-7	Ammonium hydrogen difluoride	172
1343-88-0	Florisil 60-100 mesh	356
1343-88-0	Florisil 100-200 mesh	357
1344-09-8	Sodium silicate	705
1344-28-1	Alumina white	151
1344-28-1	Aluminum oxide	155
1344-28-1	Aluminum oxide (acid)	155
1344-28-1	Aluminum oxide (basic)	155
1344-28-1	Aluminum oxide (neutral)	156
1344-28-1	Aluminum oxide activated	156
1393-92-6	Litmus	465
1400-62-0	Orcein	540
1401-55-4	Tannic acid	744
1455-13-6	Methanol-d1	494
1465-25-4	N-(1-Naphtyl)ethylenediamine dihydrochloride	518
1470-61-7	Silver diethyldithiocarbamate	655
1516-08-1	Ethanol-d6 anhydrous	340
1592-23-0	Calcium stearate	244
1611-35-4	Xylenol orange	793
1633-05-2	Strontium carbonate	724
1634-04-4	tert-Butylmethylether	230
1643-19-2	Tetrabutylammonium bromide	749
1665-00-5	Dichloromethane-d2	305
1668-00-4	Arsenazo III	184
1693-74-9	Tetrahydrofuran-d8	753
1733-12-6	o-Cresol red	283
1733-12-6	o-Cresol Red solution 0.2% in ethanol	283
1735-17-7	Cyclohexane-d12	287
1762-95-4	Ammonium thiocyanate	178
1762-95-4	Ammonium thiocyanate 1 mol/l (1N)	179
1762-95-4	Ammonium thiocyanate 0.1 mol/l (0.1N)	179
1762-95-4	Ammonium thiocyanate 0.01 mol/l (0.01N)	179
1787-61-7	Eriochrome black T	329
1829-00-1	Clayton's yellow	270
1849-29-2	Methanol-d3	494
1934-21-0	Tartrazine	746
1936-15-8	Orange G	539
1945-77-3	Methylthymol blue sodium salt	506
2037-26-5	Toluene-d8	768
2052-49-5	Tetrabutylammonium hydroxide solution 0.1 mol/l (0.1N)	749
2052-49-5	Tetrabutylammonium hydroxide 0.1 mol/l (0.1N) in 2-propanol	749
2052-49-5	Tetrabutylammonium hydroxide 0.1 mol/l (0.1N) in methanol / propanol-2 (50/50)	750
2199-69-1	1,2-Dichlorobenzene-d4	299
2206-26-0	Acetonitrile-d3	145
2206-27-1	Dimethylsulphoxide-d6	319
2206-27-1	Dimethylsulphoxide-d6 + 0.03% TMS	319
2216-51-5	L-Menthol	483
2243-76-7	Alizarin yellow R	150
2303-01-7	m-Cresol purple	282
2321-07-5	Fluorescein	357
2338-05-8	Iron (III) citrate	434
2353-45-9	Fast green FCF	353
2386-53-0	1-Dodecanesulfonic acid sodium salt	324
2386-54-1	1-Butanesulfonic acid sodium salt	225
2538-85-4	Calcon	245
2650-17-1	Xylenecyanol	792
2832-45-3	1-Hexanesulphonic acid sodium salt	389
3012-65-5	Ammonium citrate dibasic	170
3012-65-5	Ammonium citrate solution 20%	171
3051-09-0	Murexide	516
3147-14-6	Calmagite	246
3164-29-2	Ammonium L(+)tartrate	178
3244-88-0	Fuchsin acid	364
3458-72-8	tri-Ammonium citrate	170
3564-18-9	Eriochromocyanine R	330
3737-95-9	Calconcarbonic acid	245
4196-99-0	Ponceau red BS	581
4197-25-5	Sudan black B	726
4472-41-7	N,N-Dimethylformamide-d7	316
4637-24-5	N,N-Dimethylformamide dimethylacetal	317
4717-96-8	4-Methyltetrahydropyran	506
5108-96-3	Pyrrolidine dithiocarbamic acid ammonium salt	633
5141-20-8	Light green	461
5144-89-8	o-Phenantroline monohydrate	569
5263-02-5	Zinc carbonate basic	799
5324-84-5	1-Octanesulphonic acid sodium salt	537
5329-14-6	Sulfamic acid	727
5341-61-7	Hydrazine dihydrochloride	392
5392-40-5	Citral	268
5614-37-9	Cyclopentyl methyl ether	289
5743-04-4	Cadmium acetate dihydrate	233
5743-26-0	Calcium acetate monohydrate	236
5743-47-5	Calcium lactate	241
5794-13-8	L(+)Asparagine	186
5808-22-0	Chromotropic acid disodium salt	268
5892-10-4	Bismuth(III) carbonate basic	203
5934-29-2	L-Histidine hydrochloride monohydrate	390
5949-29-1	Citric acid monohydrate	269
5950-69-6	Hydrindantin	392
5965-83-3	Sulfosalicylic acid	729
5968-11-6	Sodium carbonate monohydrate	668
5970-45-6	Zinc acetate dihydrate	798
5995-86-8	Gallic acid monohydrate	367
5996-10-1	D(+)Glucose monohydrate	371
6001-64-5	Chlorobutanol hemihydrate	259
6009-70-7	Ammonium oxalate monohydrate	175
6009-70-7	Ammonium oxalate solution 4%	175
6018-89-9	Nickel (II) acetate tetrahydrate	522
6046-93-1	Copper (II) acetate hydrate	277
6080-56-4	Lead (II) acetate trihydrate	458
6100-05-6	Potassium citrate tribasic monohydrate	591
6100-19-2	Potassium tartrate	618
6104-59-2	Coomassie brilliant blue R 250	274
6106-21-4	Sodium succinate hexahydrate	706
6106-24-7	Sodium tartrate dihydrate	709
6131-90-4	Sodium acetate trihydrate	662
6132-02-1	Sodium carbonate decahydrate	668
6132-04-3	Sodium citrate tribasic dihydrate	673
6147-53-1	Cobalt (II) acetate tetrahydrate	271
6152-67-6	4-Diphenylaminesulfonic acid sodium salt	322
6153-39-5	Orcinol monohydrate	540
6153-56-6	Oxalic acid dihydrate	550
6156-78-1	Manganese (II) acetate tetrahydrate	480
6192-52-5	p-Toluenesulfonic acid	768
6226-79-5	Ponceau red S	581
6283-63-2	n,n-Diethyl-p-phenylenediamine sulfate	310
6303-21-5	Hypophosphorous acid 50%	418
6363-53-7	Maltose monohydrate	478
6381-59-5	Potassium sodium tartrate tetrahydrate	617
6381-92-6	Ethylenediaminetetraacetic acid disodium salt	345
6381-92-6	Ethylenediaminetetraacetic acid disodium salt 0.1 mol/l (0.2N)	346
6381-92-6	Ethylenediaminetetraacetic acid disodium salt 0.05 mol/l (0.1N)	347
6381-92-6	Ethylenediaminetetraacetic acid disodium salt 0.01 mol/l (0.02N)	347
6409-77-4	Nuclear fast red	536
6484-52-2	Ammonium nitrate	174
6484-52-2	Ammonium nitrate 200 mg/l solution	174
6487-48-5	Potassium oxalate monohydrate	610
6591-63-5	Quinidine sulfate	634
7021-09-2	alpha-Methoxyphenylacetic acid	495
7080-50-4	Chloramine T sodium salt	256
7087-68-5	N,N-Diisopropylethylamine	312
7220-79-3	Methylene blue	498
7226-23-5	N,N'-Dimethylpropylene uree	317
7291-22-7	Pyridine-d5	631
7296-20-0	Neocuproine hydrochloride	519
7320-34-5	tetra-Potassium pyrophosphate	616
7429-90-5	Aluminum, powder	151
7439-89-6	Iron, powder	427
7439-89-6	Iron, reduced by hydrogen	428
7439-92-1	Lead	455
7439-95-4	Magnesium, powder	468
7439-95-4	Magnesium, ribbon	468
7439-95-4	Magnesium, turnings	468
7439-96-5	Manganese electrolytic	479
7440-02-0	Nickel, powder	521
7440-22-4	Silver, sheet	652
7440-22-4	Silver, wool	653
7440-31-5	Tin, powder	759
7440-31-5	Tin foil	759
7440-44-0	Charcoal activated	255
7440-44-0	Charcoal decolorizing	255
7440-44-0	Charcoal vegetable	256
7440-50-8	Copper electrolytic rebaked, sheet	274
7440-50-8	Copper electrolytic, turnings	274
7440-50-8	Copper electrolytic, wire	275
7440-50-8	Copper reduced, powder	275
7440-66-6	Zinc, foil	796
7440-66-6	Zinc, granular	796
7440-66-6	Zinc, powder	797
7440-66-6	Zinc, activated	797
7440-66-6	Zinc standard solution	797
7440-67-7	Zirconium powder	803
7440-69-9	Bismuth, granules	202
7440-70-2	Calcium standard solution	234
7446-08-4	Selenium dioxide	648
7446-14-2	Lead (II) sulfate	460
7446-19-7	Zinc sulfate monohydrate	802
7446-20-0	Zinc sulfate heptahydrate	802
7446-20-0	Zinc sulfate 0.1 mol/l (0.2N)	802
7446-20-0	Zinc sulfate 0.05 mol/l (0.05N)	803
7446-70-0	Aluminum chloride anhydrous	153
7447-40-7	Potassium chloride	586
7447-40-7	Potassium chloride 3.5 mol/l (3.5N)	587
7447-40-7	Potassium chloride 3.5 mol/l (3.5N) + silver chloride	587

Product specifications are subject to changes. Please visit our website for updates.

INDEX BY CAS NUMBER

7447-40-7	Potassium chloride 3 mol/l (3N)	588
7447-40-7	Potassium chloride 3 mol/l (3N) water-glycerol solution	588
7447-40-7	Potassium chloride 3 mol/l (3N) + silver chloride	588
7447-40-7	Potassium chloride 1 mol/l (1N)	588
7447-40-7	Potassium chloride 0.1 mol/l (0.1N)	589
7447-40-7	Potassium chloride 0.01 mol/l (0.01N)	589
7447-40-7	Potassium chloride 25g/l in HCl	589
7447-40-7	Potassium chloride 12g/l	589
7447-40-7	Potassium chloride saturated solution	590
7447-40-7	Potassium chloride solution	590
7447-41-8	Lithium chloride	463
7487-88-9	Magnesium sulfate anhydrous	475
7487-94-7	Mercury (II) chloride	485
7488-55-3	Tin (II) sulfate	762
7553-56-2	Iodine	424
7553-56-2	Iodine resublimed	425
7553-56-2	Iodine 0.5 mol/l (1N)	425
7553-56-2	Iodine 0.05 mol/l (0.1N)	425
7553-56-2	Iodine 0.01 mol/l (0.02N)	426
7553-56-2	Iodine 0.005 mol/l (0.01N)	426
7558-79-4	Sodium phosphate dibasic anhydrous	700
7580-67-8	Lithium hydride	464
7601-90-3	Perchloric acid 65-71%	560
7601-90-3	Perchloric acid 65%	561
7631-86-9	Silica gel 60A 6 - 35µ	648
7631-86-9	Silica gel 60A 20 - 45µ	648
7631-86-9	Silica gel 60A 35 - 70µ	649
7631-86-9	Silica gel 60A 40 - 63µ	649
7631-86-9	Silica gel 60A 70 - 200µ	650
7631-86-9	Silica gel 60A 0,06÷0,20 mm	650
7631-86-9	Silica gel granular	650
7631-86-9	Silica gel granular with indicator cobalt free	651
7631-99-4	Sodium nitrate	696
7632-00-0	Sodium nitrite	696
7632-00-0	Sodium nitrite 0.1 mol/l (0.1N)	697
7632-00-0	Sodium nitrite solution 500 g/l	697
7646-79-9	Cobalt (II) chloride in solution	272
7646-85-7	Zinc chloride anhydrous	799
7646-85-7	Zinc chloride solution 60%	799
7646-93-7	Potassium bisulfate	584
7647-01-0	Hydrochloric acid 50% v/v	393
7647-01-0	Hydrochloric acid 37%	394
7647-01-0	Hydrochloric acid 34-37%	396
7647-01-0	Hydrochloric acid 32-35%	396
7647-01-0	Hydrochloric acid 32% (20°Bé)	397
7647-01-0	Hydrochloric acid 32%	397
7647-01-0	Hydrochloric acid 29-31%	398
7647-01-0	Hydrochloric acid 26%	398
7647-01-0	Hydrochloric acid 25% w/v	398
7647-01-0	Hydrochloric acid 23%	399
7647-01-0	Hydrochloric acid 20%	399
7647-01-0	Hydrochloric acid 12%	399
7647-01-0	Hydrochloric acid 10%	400
7647-01-0	Hydrochloric acid 8%	400
7647-01-0	Hydrochloric acid 5%	400
7647-01-0	Hydrochloric acid 1.128% m/v	401
7647-01-0	Hydrochloric acid 6 mol/l (6N)	401
7647-01-0	Hydrochloric acid 5 mol/l (5N)	402
7647-01-0	Hydrochloric acid 4 mol/l (4N)	402
7647-01-0	Hydrochloric acid 3 mol/l (3N)	403
7647-01-0	Hydrochloric acid 2 mol/l (2N)	403
7647-01-0	Hydrochloric acid 1 mol/l (1N)	403
7647-01-0	Hydrochloric acid 0.714 mol/l (N/1.4)	404
7647-01-0	Hydrochloric acid 0.5 mol/l (0.5N)	405
7647-01-0	Hydrochloric acid 0.2 mol/l (0.2N)	405
7647-01-0	Hydrochloric acid 0.2 mol/l (0.2N) in propanol-2	405
7647-01-0	Hydrochloric acid 0.1 mol/l (0.1N)	406
7647-01-0	Hydrochloric acid 0.1 mol/l (0.1N) in ethanol	407
7647-01-0	Hydrochloric acid 0.1 mol/l (0.1N) in propanol-2	407
7647-01-0	Hydrochloric acid 0.0714 mol/l (N/14)	407
7647-01-0	Hydrochloric acid 0.05 mol/l (0.05N)	407
7647-01-0	Hydrochloric acid 0.04 mol/l (0.04N)	408
7647-01-0	Hydrochloric acid 0.02 mol/l (0.02N)	408
7647-01-0	Hydrochloric acid 0.01 mol/l (0.01N)	408
7647-01-0	Hydrochloric acid, dilute	409
7647-01-0	Hydrochloric acid, brominated	409
7647-10-1	Palladium (II) chloride	553
7647-14-5	Ausilab 260	188
7647-14-5	Sodium chloride	670
7647-14-5	Sodium chloride 5 mol/l (5N)	672
7647-14-5	Sodium chloride 0.1 mol/l (0,1N)	672
7647-15-6	Sodium bromide	667
7647-17-8	Cesium chloride	254
7647-17-8	Cesium chloride 25 g/l solution	254
7664-38-2	Orthophosphoric acid 99%	547
7664-38-2	Orthophosphoric acid 85%	547
7664-38-2	Orthophosphoric acid 75%	549
7664-38-2	Orthophosphoric acid 10%	549
7664-39-3	Hydrofluoric acid 50%	410
7664-39-3	Hydrofluoric acid 47-51%	411
7664-39-3	Hydrofluoric acid 39.5%	412
7664-39-3	Hydrofluoric acid diluted	412
7664-93-9	Sulfuric acid 98%	730
7664-93-9	Sulfuric acid 96% (66°Be)	731
7664-93-9	Sulfuric acid 96%	731
7664-93-9	Sulfuric acid 95 - 97 %	733
7664-93-9	Sulfuric acid 93-98%	733
7664-93-9	Sulfuric acid 90%	734
7664-93-9	Sulfuric acid 85%	734
7664-93-9	Sulfuric acid 72%	735
7664-93-9	Sulfuric acid 69%	735
7664-93-9	Sulfuric acid 62%	735
7664-93-9	Sulfuric acid 50%	735
7664-93-9	Sulfuric acid 35% (30°Be)	736
7664-93-9	Sulfuric acid 30%	736
7664-93-9	Sulfuric acid 25%	736
7664-93-9	Sulfuric acid 20%	736
7664-93-9	Sulfuric acid 10% v/v	737
7664-93-9	Sulfuric acid 4 mol/l (8N)	737
7664-93-9	Sulfuric acid 2.5 mol/l (5N)	737
7664-93-9	Sulfuric acid 1 mol/l (2N)	737
7664-93-9	Sulfuric acid 0.5 mol/l (1N)	738
7664-93-9	Sulfuric acid 0.33 mol/l (2N/3)	738
7664-93-9	Sulfuric acid 0.26 mol/l (0.52N)	739
7664-93-9	Sulfuric acid 0.25 mol/l (0.5N)	739
7664-93-9	Sulfuric acid 0.166 mol/l (0.333N)	739
7664-93-9	Sulfuric acid 0.13 mol/l (0.26N)	740
7664-93-9	Sulfuric acid 0.125 mol/l (0.25N)	740
7664-93-9	Sulfuric acid 0.1 mol/l (0.2N)	740
7664-93-9	Sulfuric acid 0.05 mol/l (0.1N)	740
7664-93-9	Sulfuric acid 0.025 mol/l (0.05N)	741
7664-93-9	Sulfuric acid 0.02 mol/l (0.04N)	741
7664-93-9	Sulfuric acid 0.01 mol/l (0.02N)	741
7664-93-9	Sulfuric acid 0.005 mol/l (0.01N)	742
7664-93-9	Sulfuric acid 0.0025 mol/l (0.005N)	742
7664-93-9	Sulfuric acid with 10 g/l Silver sulfate	742
7664-93-9	Sulfuric acid with 6.6 g/l Silver sulfate	743
7664-93-9	Sulfuric acid, dilute	743
7664-93-9	Sulfuric acid d=1.820	743
7681-11-0	Potassium iodide	607
7681-11-0	Potassium iodide solution 10%	608
7681-11-0	Potassium iodide solution 3.9%	608
7681-11-0	Potassium iodide solution	608
7681-49-4	Sodium fluoride	675
7681-52-9	Sodium hypochlorite solution 12.5%	692
7681-52-9	Sodium hypochlorite solution in water	692
7681-57-4	Sodium metabisulfite	694
7681-65-4	Copper (I) iodide	277
7681-82-5	Sodium iodide	693
7697-37-2	Nitric acid fuming 99%	528
7697-37-2	Nitric acid 69.5%	528
7697-37-2	Nitric acid 67-70%	529
7697-37-2	Nitric acid 67-69%	530
7697-37-2	Nitric acid 67.5 (42° Be)	530
7697-37-2	Nitric acid 65%	531
7697-37-2	Nitric acid 18%	532
7697-37-2	Nitric acid 10%	532
7697-37-2	Nitric acid 8 mol/l (8N)	532
7697-37-2	Nitric acid 2 mol/l (2N)	533
7697-37-2	Nitric acid 1 mol/l (1N)	533
7697-37-2	Nitric acid 0.1 mol/l (0.1N)	533
7697-37-2	Nitric acid, dilute	534
7698-05-7	Hydrochloric acid-d 20%	409
7698-05-7	Hydrochloric acid-d 1 mol/l	409
7704-34-9	Sulfur sublimed and washed	730
7705-07-9	Titanium trichloride-sulfuric acid reagent	765
7705-08-0	Iron (III) chloride anhydrous sublimed	433
7705-08-0	Iron chloride in solution	434
7722-64-7	Potassium permanganate	611
7722-64-7	Potassium permanganate solution 3%	612
7722-64-7	Potassium permanganate 0.2 mol/l (1N)	612
7722-64-7	Potassium permanganate 0.02 mol/l (0.1N)	612
7722-64-7	Potassium permanganate 0.002 mol/l (0.01N)	613
7722-64-7	Potassium permanganate and phosphoric acid solution	613
7722-76-1	Ammonium di-hydrogen phosphate 25 mg/L solution	172
7722-76-1	Ammonium phosphate monobasic	176
7722-84-1	Hydrogen peroxide solution 40% w/v	413
7722-84-1	Hydrogen peroxide solution 35%	413
7722-84-1	Hydrogen peroxide solution 30-32%	413
7722-84-1	Hydrogen peroxide solution 30%	414
7722-84-1	Hydrogen peroxide solution 6%	415
7722-84-1	Hydrogen peroxide solution 3.5% w/v	415
7722-84-1	Hydrogen peroxide solution 3%	416
7726-95-6	Bromine solution	210
7726-95-6	Bromine water	211
7727-21-1	Potassium persulfate	613

INDEX BY CAS NUMBER

7727-43-7	Barium sulfate	196
7727-54-0	Ammonium persulfate	175
7727-73-3	Sodium sulfate decahydrate	707
7732-18-5	Water	786
7732-18-5	Water purified	788
7732-18-5	Water + 0.1% v/v formic acid	789
7732-18-5	Water + 0.1% v/v trifluoroacetic acid	789
7757-79-1	Potassium nitrate	609
7757-79-1	Potassium nitrate 1 mol/l (1N)	610
7757-82-6	Sodium sulfate anhydrous	706
7757-83-7	Sodium sulfite anhydrous	708
7758-01-2	Potassium bromate	584
7758-01-2	Potassium bromate 0.0167 mol/l (0.1N)	585
7758-02-3	Potassium bromide	585
7758-05-6	Potassium Iodate	606
7758-05-6	Potassium iodate 0.05 mol/l (0.3N)	606
7758-05-6	Potassium iodate 0.0167 mol/l (0.1N)	606
7758-05-6	Potassium iodate 0.00167 mol/l (0.01N)	607
7758-09-0	Potassium nitrite	610
7758-11-4	Potassium phosphate dibasic anhydrous	614
7758-23-8	Calcium phosphate monobasic monohydrate	243
7758-87-4	Calcium phosphate tribasic	243
7758-89-6	Copper (I) chloride	276
7758-89-6	Copper (I) chloride solution 7% in ammonia	277
7758-98-7	Copper (II) sulfate anhydrous	280
7758-98-7	Copper (II) sulfate solution 12.5%	281
7758-98-7	Copper (II) sulphate solution	281
7758-99-8	Copper (II) sulfate pentahydrate	280
7759-02-6	Strontium sulfate	725
7761-88-8	Silver nitrate	655
7761-88-8	Silver nitrate solution 5%	656
7761-88-8	Silver nitrate solution 2.9075%	656
7761-88-8	Silver nitrate 1 mol/l (1N)	656
7761-88-8	Silver nitrate 0.5 mol/l (0.5N)	656
7761-88-8	Silver nitrate 0.1 mol/l (0.1N)	657
7761-88-8	Silver nitrate 0.05 mol/l (0.05N)	657
7761-88-8	Silver nitrate 0.01 mol/l (0.01N)	658
7761-88-8	Silver nitrate 0.01 mol/l (0.01N) in propanol-2	658
7761-88-8	Silver nitrate solution	658
7772-98-7	Sodium thiosulfate anhydrous	711
7773-06-0	Ammonium sulfamate	177
7774-29-0	Mercury (II) iodide	486
7774-34-7	Calcium chloride hexahydrate	239
7775-14-6	Sodium hydrosulfite	678
7775-27-1	Sodium persulfate	699
7775-27-1	Sodium persulfate 1 mol/l	700
7778-50-9	Potassium dichromate	592
7778-50-9	Potassium dichromate solution 0.5%	592
7778-50-9	Potassium dichromate 0.167 mol/l (1N)	592
7778-50-9	Potassium dichromate 0.0417 mol/l (0.25N)	593
7778-50-9	Potassium dichromate 0.04 mol/l (0.24N) in 80 g/l HgSO4	593
7778-50-9	Potassium dichromate 0.0167 mol/l (0.1N)	593
7778-50-9	Potassium dichromate solution 106 g/l	594
7778-50-9	Potassium dichromate 0.1414 g/l	594
7778-77-0	Potassium phosphate monobasic	615
7778-77-0	Potassium phosphate monobasic 0.2 mol/l (0.2N)	615
7778-80-5	Potassium sulfate	618
7782-49-2	Selenium, powder	646
7782-50-5	Water chlorine	789
7782-61-8	Iron (III) nitrate nonahydrate	435
7782-63-0	Iron (II) sulfate heptahydrate	430
7782-75-4	Magnesium hydrogen phosphate trihydrate	472
7783-00-8	Selenous acid	648
7783-03-1	Tungstic acid	780
7783-20-2	Ammonium sulfate	177
7783-28-0	Ammonium phosphate dibasic	176
7783-35-9	Mercury (II) sulfate	486
7783-83-7	Iron (III) ammonium sulfate dodecahydrate	432
7783-85-9	Iron (II) ammonium sulfate hexahydrate	429
7783-85-9	Iron (II) ammonium sulfate 0.12N	429
7783-85-9	Iron (II) ammonium sulfate 0.1N	430
7783-90-6	Silver chloride	654
7784-13-6	Aluminum chloride hexahydrate	153
7784-24-9	Aluminum potassium sulfate dodecahydrate	156
7784-26-1	Aluminum ammonium sulfate dodecahydrate	152
7784-27-2	Aluminum nitrate nonahydrate	154
7784-31-8	Aluminum sulfate	157
7784-46-5	Sodium arsenite 0.1 mol/l (0.2N)	664
7784-46-5	Sodium arsenite 0.05 mol/l (0.1N)	664
7785-20-8	Nickel (II) ammonium sulfate hexahydrate	523
7788-99-0	Chromium (III) potassium sulfate dodecahydrate	266
7789-00-6	Potassium chromate	590
7789-00-6	Potassium chromate solution 10%	590
7789-00-6	Potassium chromate 5% solution	591
7789-02-8	Chromium (III) nitrate nonahydrate	266
7789-20-0	Deuterium oxide-d2	293
7789-20-0	Deuterium oxide-d2 + 0.01% DMSO	294
7789-20-0	Deuterium oxide-d2 + 0.5% TSP d4	294
7789-20-0	Deuterium oxide-d2 + 0.03% TSP d4	294
7789-23-3	Potassium fluoride	596
7789-75-5	Calcium fluoride	240
7789-77-7	Calcium phosphate dibasic dihydrate	243
7790-21-8	Potassium m-periodate	611
7790-28-5	Sodium metaperiodate	694
7790-58-1	Potassium tellurite	619
7790-62-7	Potassium pyrosulphate	616
7790-69-4	Lithium nitrate	465
7791-07-3	Sodium perchlorate monohydrate	699
7791-13-1	Cobalt (II) chloride hexahydrate	272
7791-18-6	Magnesium chloride hexahydrate	471
7791-20-0	Nickel (II) chloride hexahydrate	523
8001-79-4	Castor oil	249
8002-27-9	Oil of cedar wood	538
8002-43-5	Lecithin of soya	461
8005-03-6	Nigrosine	525
8005-78-5	Bismarck brown R	201
8006-28-8	Soda lime	660
8006-54-0	Lanolin anhydrous	453
8007-47-4	Canada balsam	247
8007-69-0	Oil refined of almonds	539
8009-03-8	Paraffin white soft	556
8012-95-1	Paraffin oil	555
8028-48-6	Histolemon	390
8047-15-2	Saponin	644
8049-11-4	Devarda's alloy	294
9000-01-5	Gum arabic	377
9000-70-8	Gelatine	368
9000-71-9	Casein	249
9002-93-1	Triton© X100 solution	778
9003-39-8	Polyvinylpyrrolidone	581
9004-34-6	Cellulose, powder	250
9005-25-8	Rice starch	640
9005-38-3	Sodium alginate	663
9005-80-5	Inulin	424
9005-84-9	Maize starch	476
9005-84-9	Starch paste solution 1%	721
9005-84-9	Starch soluble	721
9006-59-1	Albumin from eggs powder	148
9006-59-1	Albumin from eggs, dried	148
9016-45-9	Nonylphenol ethoxylated 10 ETO	536
10022-31-8	Barium nitrate	195
10022-68-1	Cadmium nitrate tetrahydrate	234
10025-69-1	Tin (II) chloride dihydrate	761
10025-69-1	Tin (II) chloride solution	761
10025-69-1	Tin (II) chloride solution 10%	761
10025-70-4	Strontium chloride hexahydrate	724
10025-77-1	Iron (III) chloride hexahydrate	433
10025-91-9	Antimony trichloride	183
10026-06-9	Tin (IV) chloride pentahydrate	762
10026-22-9	Cobalt (II) nitrate hexahydrate	272
10026-24-1	Cobalt (II) sulfate heptahydrate	273
10028-24-7	Sodium phosphate dibasic dihydrate	701
10030-85-0	L(+)Rhamnose	638
10031-43-3	Copper (II) nitrate trihydrate	279
10034-76-1	Calcium sulfate hemihydrate	245
10034-81-8	Magnesium perchlorate	474
10034-85-2	Hydriodic acid 57%	393
10034-88-5	Sodium bisulfate monohydrate	666
10034-96-5	Manganese (II) sulfate monohydrate	481
10034-99-8	Magnesium sulfate heptahydrate	475
10035-04-8	Calcium chloride dihydrate	238
10035-06-0	Bismuth(III) nitrate pentahydrate	204
10035-10-6	Hydrobromic acid 48%	393
10039-26-6	Lactose monohydrate	452
10039-32-4	Sodium phosphate dibasic dodecahydrate	702
10039-54-0	Hydroxylamine sulphate	417
10039-56-2	Sodium hypophosphite	693
10042-76-9	Strontium nitrate	724
10043-35-3	Boric acid	205
10043-35-3	Boric acid 4%	206
10043-35-3	Boric acid 4% with indicator	206
10043-35-3	Boric acid 3%	206
10043-35-3	Boric acid 1% with indicator	207
10043-35-3	Boric acid 20g/l	207
10043-35-3	Boric acid 20 g/l with indicator	207
10043-52-4	Calcium chloride anhydrous	237
10049-21-5	Sodium phosphate monobasic monohydrate	703
10060-12-5	Chromium (III) chloride hexahydrate	265
10099-58-8	Lanthanum chloride 25 g/l solution	454
10099-74-8	Lead (II) nitrate	458
10101-41-4	Calcium sulfate dihydrate	244
10101-89-0	Sodium phosphate tribasic dodecahydrate	704

DASIT GROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

INDEX BY CAS NUMBER

10101-97-0	Nickel (II) sulfate hexahydrate	524
10102-05-3	Palladium nitrate 2 g/l solution	554
10102-17-7	Sodium thiosulfate pentahydrate	711
10102-17-7	Sodium thiosulfate 1 mol/l (1N)	712
10102-17-7	Sodium thiosulfate 0.5 mol/l (0.5N)	712
10102-17-7	Sodium thiosulfate 0.2 mol/l (0.2N)	712
10102-17-7	Sodium thiosulfate 0.1 mol/l (0.1N)	712
10102-17-7	Sodium thiosulfate 0.0394 mol/l (0.0394N)	713
10102-17-7	Sodium thiosulfate 0.0197 mol/l (0.0197N)	713
10102-17-7	Sodium thiosulfate 0.01 mol/l (0.01N)	713
10102-25-7	Lithium sulfate monohydrate	465
10102-40-6	Sodium molybdate dihydrate	695
10112-91-1	Mercury (I) chloride	485
10124-56-8	Sodium hexametaphosphate	677
10125-13-0	Copper (II) chloride dihydrate	278
10196-18-6	Zinc nitrate hexahydrate	800
10213-10-2	Sodium tungstate dihydrate	714
10277-43-7	Lanthanum nitrate hexahydrate	454
10294-26-5	Silver sulfate	659
10294-26-5	Silver sulfate solution 0.7% in sulfuric acid	659
10294-39-0	Barium perchlorate trihydrate	195
10294-40-3	Barium chromate	194
10294-41-4	Cerium (III) nitrate hexahydrate	251
10294-42-5	Cerium (IV) sulfate tetrahydrate	252
10294-42-5	Cerium (IV) sulfate 0.1 mol/l	253
10294-54-9	Cesium sulfate	254
10326-27-9	Barium chloride solution 10%	193
10326-27-9	Barium chloride 0.1 mol/l (0.2N)	193
10326-27-9	Barium chloride 30 g/l	194
10361-03-2	Sodium metaphosphate	695
10361-29-2	Ammonium carbonate	168
10361-37-2	Barium chloride dihydrate	192
10361-37-2	Barium chloride solution 61 g/l	193
10378-47-9	Cerium (IV) ammonium sulfate dihydrate	252
10378-47-9	Cerium (IV) ammonium sulfate 0.1 mol/l	252
10378-47-9	Cerium (IV) ammonium sulfate 0.01 mol/l	252
10416-59-8	N,O-Bis(trimethylsilyl) acetamide	204
10450-60-9	Periodic acid	563
10486-00-7	Sodium perborate tetrahydrate	698
11024-24-1	Digitonin	311
11138-49-1	Sodium aluminate	663
12003-78-0	Raney's alloy	635
12007-60-2	Lithium tetraborate anhydrous	465
12027-06-4	Ammonium iodide	173
12027-43-9	Silicotungstic acid	652
12054-85-2	Ammonium molybdate tetrahydrate	173
12054-85-2	Ammonium molybdate solution 2.5% in nitric acid	174
12055-62-8	Holmium oxide	391
12069-69-1	Copper (II) carbonate (basic)	278
12124-97-9	Ammonium bromide	167
12125-01-8	Ammonium fluoride	171
12125-02-9	Ammonium chloride	169
12135-76-1	Ammonium sulfide solution 20%	178
12208-13-8	Potassium pyroantimonate acid	615
12208-13-8	Potassium pyroantimonate solution	616
12230-71-6	Barium hydroxide octahydrate	194
12230-71-6	Barium hydroxide solution 5%	195
12230-71-6	Barium hydroxide solution 47.3 g/l	195
12501-23-4	Phosphotungstic acid	578
13011-54-6	Sodium ammonium hydrogen phosphate	664
13235-36-4	Ethylenediaminetetraacetic acid tetrasodium salt tetrahydrate	348
13408-09-8	Sodium glycerophosphate pentahydrate	676
13419-61-9	1-Decanesulfonic acid sodium salt	290
13446-18-9	Magnesium nitrate hexahydrate	472
13446-34-9	Manganese (II) chloride tetrahydrate	480
13455-21-5	Potassium fluoride dihydrate	596
13455-24-8	Potassium hydrogen iodate	597
13463-67-7	Titanium dioxide	764
13465-95-7	Barium perchlorate 0.05 mol/l	196
13465-95-7	Barium perchlorate 0.025 mol/l	196
13472-35-0	Sodium phosphate monobasic dihydrate	703
13472-36-1	Sodium pyrophosphate decahydrate	705
13477-34-4	Calcium nitrate tetrahydrate	241
13478-00-7	Nickel (II) nitrate hexahydrate	524
13478-10-9	Iron (II) chloride tetrahydrate	430
13494-80-9	Tellurium lumps	746
13494-80-9	Tellurium powder	747
13517-23-2	Sodium o-Phosphite pentahydrate	704
13573-16-5	Reinecke salt	637
13596-46-8	Cobalt (II) ammonium sulfate hexahydrate	271
13598-36-2	Phosphonic acid	577
13600-98-1	Sodium cobalt nitrite	673
13746-66-2	Potassium ferricyanide	594
13755-38-9	Sodium nitroprusside dihydrate	697
14014-06-3	Sodium hydroxide-d 1 30%	691
14014-06-3	Sodium hydroxide-d 1 mol/l	691
14017-54-0	Holmium perchlorate in solution	391
14335-33-2	Orthophosphoric acid-d3 85% in D2O	549
14452-57-4	Magnesium peroxide	474
14459-95-1	Potassium ferrocyanide trihydrate	595
14459-95-1	Potassium ferrocyanide solution 10%	595
14459-95-1	Potassium ferrocyanide solution 53 g/l	596
14519-13-2	Strontium bromide monohydrate	723
14533-63-2	1-Propanesulfonic acid sodium salt	627
14552-35-3	Cupriethylenediamine solution	285
14634-91-4	Ferroin 0.025 mol/l solution	354
14634-91-4	o-Phenantroline-Iron (II) sulphate solution in sulphuric acid	569
14726-36-4	Zinc dibenzyldithiocarbamate	800
14807-96-6	Talc	744
14808-60-7	Sand of Fontainebleau	644
14808-60-7	Sand purified	644
14808-60-7	Silicon dioxide	652
14985-18-3	Zirconyl nitrate	805
15187-32-3	Iron (III) ammonium oxalate	432
15244-10-7	Iron (III) sulfate	435
15244-38-9	Chromium (III) sulfate	267
15708-48-2	Ethylenediaminetetraacetic acid potassium and magnesium salt dihydrate	347
16423-68-0	Erythrosin extra B	330
16674-78-5	Magnesium acetate tetrahydrate	470
16731-55-8	Potassium metabisulfite	609
16774-21-3	Cerium (IV) ammonium nitrate	251
16774-21-3	Cerium (IV) ammonium nitrate 0.1 mol/l	251
16774-21-3	Cerium (IV) ammonium nitrate 0.01 mol/l	251
16788-57-1	Potassium phosphate dibasic trihydrate	614
16893-85-9	Sodium hexafluorosilicate	677
16921-30-5	di-Potassium hexachloroplatinate	597
16940-66-2	Sodium borohydride	666
16961-25-4	Gold(III) chloride trihydrate	375
17372-87-1	Eosin Y	328
17372-87-1	Eosin Y 1% solution acqueous	328
17372-87-1	Eosin Y 0.5% solution alcoholic	329
17629-30-0	Raffinose	635
18016-24-5	Calcium gluconate	240
18282-10-5	Tin (IV) oxide	762
18472-87-2	Phloxin B	574
18497-13-7	Hexachloroplatinic acid hexahydrate	383
20624-25-3	Sodium diethyldithiocarbamate trihydrate	674
20667-12-3	Silver oxide	659
21645-51-2	Aluminum hydroxide	154
21908-53-2	Mercury (II) oxide red	486
22383-16-0	Methyl green	499
22767-49-3	1-Pentanesulphonic acid sodium salt	559
22767-50-6	1-Heptanesulphonic acid sodium salt	383
24634-61-5	Potassium sorbate	617
25102-12-9	Ethylenediaminetetraacetic acid dipotassium salt dihydrate	345
25155-30-0	Dodecylbenzenesulphonic acid sodium salt	325
25561-30-2	N,O-Bis(trimethylsilyl)-trifluoroacetamide	204
25895-60-7	Sodium cyanoborohydride	674
26628-22-8	Sodium azide	665
28300-74-5	Antimony potassium tartrate	183
28631-66-5	Aniline blue soluble in water	180
28983-56-4	Methyl blue	497
30525-89-4	Paraformaldehyde	556
32503-27-8	Tetrabutylammonium bisulfate	748
33685-54-0	Tetrachloroethane-d2	751
33864-99-2	Alcian blue 8GX	149
35658-65-2	Cadmium chloride monohydrate	233
36653-82-4	Cetyl alcohol	254
37247-10-2	Azure II	189
37267-86-0	Metaphosphoric acid	487
38894-11-0	3-Methyl-2-benzothiazolinone hydrazone hydrochloride	497
39409-82-0	Magnesium carbonate basic	470
39430-27-8	Nickel (II) carbonate basic	523
42607-30-7	Sodium sulfide trihydrate	708
51429-74-4	Phosphomolybdic acid	576
51580-86-0	Sodium dichloroisocyanurate dihydrate	674
53092-85-6	Azure II eosin	190
63449-41-2	Benzalkonium chloride	197
64742-47-8	Kerosene	449
64742-49-0	Petroleum benzin E	563
64742-49-0	Petroleum ether 100 - 140°C	564
64742-49-0	Petroleum ether 80 - 120°C	564
64742-49-0	Petroleum ether 80 - 100°C	564
64742-49-0	Petroleum ether 75 - 120°C	565
64742-49-0	Petroleum ether 60 - 80°C	565
64742-49-0	Petroleum ether 55 - 85°C	565
64742-49-0	Petroleum ether 40 - 70°C	566

INDEX BY CAS NUMBER

64742-49-0	Petroleum ether 40 - 65°C	566
64742-49-0	Petroleum ether 40 - 60°C	567
64771-72-8	Petroleum	563
65501-24-8	Ethylenediaminetetraacetic acid tripotassic salt	348
65997-17-3	Glass wool	370
68551-19-9	Solvent Plus	714
75006-64-3	Diacetyldioxime sodium salt	295
75881-23-1	Alcian Blue 8GS 1%	148
81029-05-2	Brillant cresyl blue	208
90622-57-4	Isopar G	440
91053-39-3	Kieselguhr composed	449
92045-76-6	Paraffin 56-58°C - Erbaplast (without DMSO)	555
93763-70-3	Dicalite 4158	297
206752-32-1	Azomethine H	189
207300-91-2	1-Hexanesulphonic acid sodium salt monohydrate	389
207596-29-0	1-Octanesulfonic acid sodium salt monohydrate	538
207605-40-1	1-Pentanesulphonic acid sodium salt monohydrate	560

Product specifications are subject to changes. Please visit our website for updates.

A

Acetale 129
Acetammide 129
Acetanilide 129
Acetilacetone 146
n-Acetil-L-cisteina 146
p-Acetilfenetidina 146
Acetofenone 146
Acetone 138
Acetone / Acqua 98/2 (v/v) + Blu di Bromofenolo 0.02 g/l 141
Acetone-d6 141
Acetonitrile 142
Acetonitrile + 0.1% v/v acido formico 145
Acetonitrile + 0.1% v/v acido trifluoroacetico 145
Acetonitrile-d3 145
Acido acetico glaciale 130
Acido acetico 96% 133
Acido acetico 80% 134
Acido acetico 45% 134
Acido acetico 30% 134
Acido acetico 27% 135
Acido acetico 25% 135
Acido acetico 20% 135
Acido acetico 12% 135
Acido acetico 1 mol/l (1N) 136
Acido acetico 0.1 mol/l (0.1N) 136
Acido acetico 0.03 mol/l (0.03N) 136
Acido acetico-d4 137
Acido adipico 147
Acido p-amminobenzoico 158
Acido amminoippurico reattivo 158
Acido 1-ammino-2-naftolo-4-solfonico 158
Acido anisico 181
Acido L(+)ascorbico 186
Acido L(+)aspartico 186
Acido barbiturico 191
Acido benzoico 199
Acido borico 205
Acido borico 4% 206
Acido borico 4% con indicatore 206
Acido borico 3% 206
Acido borico 1% con indicatore 207
Acido borico 20g/l 207
Acido borico 20 g/l con indicatore 207
Acido bromidrico 48% 393
Acido 1-butansolfonico sale sodico 225
Acido n-butirrico 231
Acido calconcarbonico 245
Acido n-caprilico 247
Acido n-caproico 247
Acido citrico anidro 269
Acido citrico monoidrato 269
Acido cloranilico 257
Acido cloridrico 50% v/v 393
Acido cloridrico 37% 394
Acido cloridrico 34-37% 396
Acido cloridrico 32-35% 396
Acido cloridrico 32% (20°Bé) 397
Acido cloridrico 32% 397
Acido cloridrico 29-31 % 398
Acido cloridrico 26% 398
Acido cloridrico 25% m/v 398
Acido cloridrico 23% 399
Acido cloridrico 20% 399
Acido cloridrico 12% 399
Acido cloridrico 10% 400
Acido cloridrico 8% 400
Acido cloridrico 5% 400
Acido cloridrico 1.128% m/v 401
Acido clorhidrico 9 mol/l (9N) 401
Acido cloridrico 6 mol/l (6N) 401
Acido cloridrico 5 mol/l (5N) 402
Acido cloridrico 4 mol/l (4N) 402
Acido cloridrico 3 mol/l (3N) 403
Acido cloridrico 2 mol/l (2N) 403
Acido cloridrico 1 mol/l (1N) 403
Acido cloridrico 0.714 mol/l (N/1.4) 404
Acido cloridrico 0.5 mol/l (0.5N) 405
Acido cloridrico 0.2 mol/l (0.2N) 405
Acido cloridrico 0.2 mol/l (0.2N) in isopropanolo 405
Acido cloridrico 0.1 mol/l (0.1N) 406
Acido cloridrico 0.1 mol/l (0.1N) in etanolo 407
Acido cloridrico 0.1 mol/l (0.1N) in isopropanolo 407
Acido cloridrico 0.0714 mol/l (N/14) 407
Acido cloridrico 0.05 mol/l (0.05N) 407
Acido cloridrico 0.04 mol/l (0.04N) 408
Acido cloridrico 0.02 mol/l (0.02N) 408
Acido cloridrico 0.01 mol/l (0.01N) 408
Acido cloridrico diluito 409
Acido cloridrico bromurato 409
Acido cloridrico-d 20% 409
Acido cloridrico-d 1 mol/l 409
Acido cloroacetico 258
Acido cromotropico sale bisodico 268
Acido 1-decansolfonico sale sodico 290
Acido dicloroacetico 297
Acido dietilentriamminopentacetico 308
Acido 1-dodecansolfonico sale sodico 324
Acido dodecilbenzensolfonico sale sodico 325
Acido eptafluorobutirrico 380
Acido 1-eptansolfonico sale sodico 383
Acido esacloroplatinico esaidrato 383
Acido 1-esansolfonico sale sodico 389
Acido 1-esansolfonico sale sodico monoidrato 389
Acido etilendiamminotetracetico 344
Acido etilendiamminotetracetico sale bipotassico diidrato 345
Acido etilendiamminotetracetico sale bisodico 345
Acido etilendiamminotetracetico sale bisodico 0.1 mol/l (0.2N) 346
Acido etilendiamminotetracetico sale bisodico 0.05 mol/l (0.1N) 347
Acido etilendiamminotetracetico sale bisodico 0.01 mol/l (0.02N) 347
Acido etilendiamminotetracetico sale di potassio e magnesio diidrato 347
Acido etilendiamminotetracetico sale tetrasodico tetraidrato 348
Acido etilendiamminotetracetico sale tripotassico 348
Acido etilenglicole bis-(2-amminoetiletere) tetracetico 349
Acido fluoridrico 50% 410
Acido fluoridrico 47-51% 411
Acido fluoridrico 39.5% 412
Acido fluoridrico 0.1M 412
Acido fluoridrico diluito 412
Acido formico 99% 362
Acido formico 85% 363
Acido formico-d 364
Acido fosfomolibdico 576
Acido fosfonico 577
Acido fosforico 99% 547
Acido fosforico 85% 547
Acido fosforico 75% 549
Acido fosforico 10% 549
Acido fosforico-d3 85% in D20 549
Acido ipofosforoso 50% 418
Acido fosfosolforico 578
Acido fosfotungstico 578
Acido fosfotungstico soluzione 578
Acido ftalico 579
Acido fumarico 365
Acido gallico monoidrato 367
Acido glicolico 374
Acido L(+)glutammico 371
Acido iodidrico 57% 393
Acido L(+)-lattico 452
Acido maleico 477
Acido DL-malico 477
Acido malonico 478
Acido D(-)mandelico 478
Acido metafosforico 487
Acido metansolfonico 487
Acido alfa-metossifenilacetico 495
Acido nicotinico 525
Acido nitrico fumante 99% 528
Acido nitrico 69,5% 528
Acido nitrico 67-70% 529
Acido nitrico 67-69% 530
Acido nitrico 67.5% (42° Be) 530
Acido nitrico 65% 531
Acido nitrico 18% 532
Acido nitrico 10% 532
Acido nitrico 8 mol/l (8N) 532
Acido nitrico 2 mol/l (2N) 533
Acido nitrico 1 mol/l (1N) 533
Acido nitrico 0.1 mol/l (0.1N) 533
Acido nitrico diluito 534
Acido m-nitrobenzoico 535
Acido nitrilotriacetico 534
Acido 1-octansolfonico sale sodico 537
Acido 1-octansolfonico sale sodico monoidrato 538
Acido oleico 539
Acido ossalico diidrato 550
Acido ossalico 0.5 mol/l (1N) 551
Acido ossalico 0.05 mol/l (0.1N) 551
Acido ossalico 0.005 mol/l (0.01N) 551
Acido 1-pentansolfonico sale sodico 559
Acido 1-pentansolfonico sale sodico monoidrato 560
Acido perclorico 65-71% 560
Acido perclorico 65% 561
Acido perclorico 0.1 mol/l (0.1N) 561
Acido perclorico 0.1 mol/l (0.1N) in acido acetico 562
Acido perclorico 0.05 mol/l (0.05N) 562
Acido perclorico 0.01 mol/l (0.01N) 562
Acido perclorico soluzione 562
Acido periodico 563
Acido picrico soluzione 579
Acido pirrolidinditiocarbammico sale ammonico 633
Acido 1-propansolfonico sale sodico 627
Acido propionico 628
Acido p-rosolico 640
Acido salicilico 642
Acido sebacico 646
Acido selenioso 648
Acido silicotungstico 652
Acido solfammico 727
Acido solfanilico 728
Acido solforico 98% 730
Acido solforico d=1,820 743
Acido solforico 96% (66°Bé) 731
Acido solforico 96% 731
Acido solforico 95 - 97 % 733
Acido solforico 93-98% 733
Acido solforico 90% 734
Acido solforico 85% 734
Acido solforico 72% 735
Acido solforico 69% 735
Acido solforico 62% 735
Acido solforico 50% 735
Acido solforico 35% (30°Bé) 736
Acido solforico 30% 736
Acido solforico 25% 736
Acido solforico 20% 736

INDICE ALFABETICO

Acido solforico 10% v/v 737
Acido solforico 4 mol/l (8N) 737
Acido solforico 2.5 mol/l (5N) 737
Acido solforico 1 mol/l (2N) 737
Acido solforico 0.5 mol/l (1N) 738
Acido solforico 0.33 mol/l (2N/3) 738
Acido solforico 0.26 mol/l (0.52N) 739
Acido solforico 0.25 mol/l (0.5N) 739
Acido solforico 0.166 mol/l (0.333N) 739
Acido solforico 0.13 mol/l (0.26N) 740
Acido solforico 0.125 mol/l (0.25N) 740
Acido solforico 0.1 mol/l (0.2N) 740
Acido solforico 0.05 mol/l (0.1N) 740
Acido solforico 0.025 mol/l (0.05N) 741
Acido solforico 0.02 mol/l (0.04N) 741
Acido solforico 0.01 mol/l (0.02N) 741
Acido solforico 0.005 mol/l (0.01N) 742
Acido solforico 0.0025 mol/l (0.005N) 742
Acido solforico, diluito 743
Acido solfosalicilico 729
Acido stearico 722
Acido succinico 725
Acido tannico 744
Acido L(+)tartarico 745
Acido L(+)tartarico soluzione 20% in acqua 746
Acido 2-tiobarbiturico 755
Acido tioglicolico 80% 755
Acido p-toluensolfonico monoidrato 768
Acido p-toluensolfonico sale sodico 769
Acido tricloroacetico 770
Acido tricloroacetico soluzione 20% 770
Acido trifluoroacetico 773
Acido trifluoroacetico-d 774
Acido tungstico 780
Acqua 786
Acqua deionizzata acidificata 789
Acqua purificata 788
Acqua + 0.1%(v/v) acido formico 789
Acqua+ 0.1% v/v acido trifluoroacetico 789
Acqua di bromo 211
Acqua di calce 184
Acqua di cloro 789
Afnio standard soluzione 378
ß-Alanina 148
Albumina d'uovo polvere 148
Albumina d'uovo secca 148
Alcole isoamilico 436
Alcole n-amilico 180
Alcole ter-amilico 180
Alcole benzilico 200
Alcole ter-butilico 227
Alcole cetilico 254
Alcole 2-feniletilico 573
Alcole metilico-d1 494
Alcole metilico-d3 494
Alcole metilico-d4 493
Alcole metilico-d4 + 0.03% TMS 494
Alcole n-ottilico 538
Aldeide anisica 181
Aldeide anisica soluzione 181
Aldeide benzoica 197
Aldeide formica 37% m/v 358
Aldeide formica 37% p/v neutralizzata 359
Aldeide formica 35% m/m 359
Aldeide formica 30 % m/v 360
Aldeide formica 10% v/v (Liquido di Lille) 360
Aldeide formica 5% m/v tamponata pH 6.9 360
Aldeide formica 4% m/v tamponata pH 6.9 361
Aldeide formica 4% m/v con sodio cloruro 361
Aldeide formica e acido acetico 362
Aldeide glutarica soluzione 50% 372
Aldeide propionica 628
Alizarina 149
Alizarina soluzione satura in alcole etilico 149
Alliltiourea 150
Allumina bianca 151
Alluminio, polvere 151
Alluminio standard soluzione 151
Alluminio ammonio solfato dodecaidrato 152
Alluminio cloruro anidro 153
Alluminio cloruro esaidrato 153
Alluminio idrossido 154
Alluminio nitrato nonaidrato 154
Alluminio ossido 155
Alluminio ossido attivo 156
Alluminio ossido (acido) 155
Alluminio ossido (basico) 155
Alluminio ossido (neutro) 156
Alluminio potassio solfato dodecaidrato 156
Alluminio solfato 157
Amido solubile 721
Amido solubile soluzione 721
Amido di mais 476
Amido di riso 640
Amido nero 10B in soluzione 157
4-Amminofenazone 159
m-Amminofenolo 159
p-Amminofenolo 159
Amminopirazolone soluzione 159
Ammoniaca soluzione 32% 160
Ammoniaca soluzione 30% 160
Ammoniaca soluzione 28% 161
Ammoniaca soluzione 25% 162
Ammoniaca soluzione 20 - 22% 163
Ammoniaca soluzione 17% 164
Ammoniaca soluzione 10% 164
Ammonio soluzione 6N 164
Ammoniaca soluzione diluita 164
Ammonio standard soluzione 165
Ammonio acetato 166
Ammonio bicarbonato 167
Ammonio bifluoruro 172
Ammonio bromuro 167
Ammonio carbammato 168
Ammonio carbonato 168
Ammonio carbonato soluzione 158 g/l 168
Ammonio citrato bibasico 170
Ammonio citrato soluzione 20% 171
Ammonio citrato tribasico 170
Ammonio cloruro 169
Ammonio diidrogeno fosfato 25 mg/l soluzione 172
Ammonio ferrosolfato (II) 0.12N 429
Ammonio fluoruro 171
Ammonio formiato 171
Ammonio fosfato bibasico 176
Ammonio fosfato monobasico 176
Ammonio ioduro 173
Ammonio molibdato tetraidrato 173
Ammonio molibdato soluzione 2.5% in acido nitrico 174
Ammonio nitrato 174
Ammonio nitrato 200 mg/l soluzione 174
Ammonio ossalato monoidrato 175
Ammonio ossalato soluzione 4% 175
Ammonio persolfato 175
Ammonio solfammato 177
Ammonio solfato 177
Ammonio solfocianuro 178
Ammonio solfocianuro 1 mol/l (1N) 179
Ammonio solfocianuro 0.1 mol/l (0.1N) 179
Ammonio solfocianuro 0.01 mol/l (0.01N) 179
Ammonio solfuro soluzione 20% 178
Ammonio L(+)tartrato 178
Anidride acetica 137
Anidride fosforica 578
Anidride ftalica 579
Anidride maleica 477
Anidride molibdica 511
Anidride succinica 725
Anidride trifluoroacetica 774
Anidride tungstica 780
Anilina cloridrato 181
Antimonio standard soluzione 182
Antimonio potassio tartrato 183
Antimonio tricloruro 183
Antrone 182
Arancio acridina 147
Arancio G 539
Arancio II 539
Arancio metile 501
Arancio metile soluzione 0.1% 501
Arancio xilenolo 793
Argento, lamina 652
Argento, lana 653
Argento standard soluzione 653
Argento acetato 654
Argento carbonato 654
Argento cloruro 654
Argento dietilditiocarbammato 655
Argento nitrato 655
Argento nitrato soluzione 5% 656
Argento nitrato soluzione 2.9075% 656
Argento nitrato soluzione 658
Argento nitrato 1 mol/l (1N) 656
Argento nitrato 0.5 mol/l (0.5N) 656
Argento nitrato 0.1 mol/l (0.1N) 657
Argento nitrato 0.1 mol/l (0.1N) in isopropanolo 657
Argento nitrato 0.05 mol/l (0.05N) 657
Argento nitrato 0.01 mol/l (0.01N) 658
Argento nitrato 0.01 mol/l (0.01N) in isopropanolo 658
Argento ossido 659
Argento solfato 659
Argento solfato soluzione 0.7% in acido solforico 659
Argento solfato 10 g/l in acido solforico 742
Argento solfato 6.6 g/l in acido solforico 743
L(+)Arginina 184
L(+)Arginina monocloridrato 184
Arsenazo III 184
Arsenico standard soluzione 185
Arsenico triossido soluzione 185
L(+)Asparagina 186
Assorbente per liquidi versati 129
Ausilab 110 187
Ausilab 140 187
Ausilab 210 187
Ausilab 250 188
Ausilab 260 188
Ausilab 280 188
Ausilab 290 188
Ausilab 400 189
Ausilab 500 189
Azometina H 189
Azzurro II 189
Azzurro II eosina 190

B

Balsamo del Canada 247
Bario standard soluzione 191
Bario acetato 192
Bario carbonato 192
Bario cloruro diidrato 192
Bario cloruro soluzione 10% 193
Bario cloruro soluzione 61 g/l 193
Bario cloruro 30 g/l 194
Bario cloruro 0.1 mol/l (0.2N) 193
Bario cromato 194
Bario idrossido ottaidrato 194
Bario idrossido soluzione 5% 195
Bario idrossido soluzione 47.3 g/l 195
Bario nitrato 195

Product specifications are subject to changes. Please visit our website for updates.

INDICE ALFABETICO

Bario perclorato triidrato 195
Bario perclorato 0.025 mol/l 196
Bario perclorato 0.05 mol/l 196
Bario solfato 196
Benedict reattivo 197
Benzalconio cloruro 197
Benzene 198
Benzene-d6 198
Benzenesolfonile cloruro 198
Benzile benzoato 200
Benzina E 563
p-Benzochinone 199
Berillio standard soluzione 201
Bis (4-idrossi-N-metilanilinio) solfato 496
Bismuto, granelli 202
Bismuto standard soluzione 202
Bismuto carbonato basico 203
Bismuto nitrato basico 203
Bismuto (III) nitrato pentaidrato 204
Biureto 97% 205
Biureto reattivo 205
Blu alcali 6B 150
Blu alcali 6B soluzione 2% in alcole etilico 150
Blu alcian 8GX 149
Blu alcian 8GS 1% 148
Blu anilina solubile in acqua 180
Blu bromofenolo 213
Blu bromofenolo indicatore 214
Blu bromofenolo indicatore TAC 214
Blu bromofenolo soluzione 214
Blu bromofenolo soluzione 0.4% in alcole etilico 213
Blu bromofenolo soluzione 0.02% 214
Blu bromotimolo 215
Blu bromotimolo soluzione 215
Blu bromotimolo soluzione 0.4% in alcole etilico 215
Blu bromotimolo 0.02% 215
Blu Coomassie R250 274
Blu cresile brillante 208
Blu lattofenolo soluzione 452
Blu metile 497
Blu metilene 498
Blu metilene fenicato soluzione idroalcolica 248
Blu metilene soluzione satura in alcole etilico 498
Blu metilene soluzione 1% 499
Blu metiltimolo sale sodico 506
Blu timolo 757
Blu timolo indicatore 758
Blu timolo soluzione 758
Blu timolo TA indicatore 758
Blu timolo soluzione 0,4% in alcole etilico 757
Blu toluidina 769
Blu toluidina fenicato soluzione idroalcolica 248
Blu vittoria 785
Boro standard soluzione 208
Brix Standards 209
Bromati standard soluzione 210
Bromo soluzione 210
alfa-Bromonaftalene 213
Bromuri - bromati 0.0167 mol/l 210
Bromuri standard soluzione 210
Bruno Bismarck R 201
1-Butanolo 226
2-Butanolo 227
n-Butile acetato 228
n-Butile cloruro 229
Butile idrossitoluene 230
ter-Butilmetiletere 230
2-Butossietanolo 228
2-(2-Butossietossi)etanolo 228

C

Cadmio standard soluzione 232
Cadmio acetato diidrato 233
Cadmio carbonato 233
Cadmio cloruro monoidrato 233
Cadmio nitrato tetraidrato 234
Caffeina anidra 234
Calce sodata 660
Calcio standard soluzione 234
Calcio acetato anidro 236
Calcio acetato monoidrato 236
Calcio carbonato 236
Calcio cloruro anidro 237
Calcio cloruro biidrato 238
Calcio cloruro esaidrato 239
Calcio cloruro soluzione 0.025% 239
Calcio fluoruro 240
Calcio fosfato bibasico diidrato 243
Calcio fosfato monobasico monoidrato 243
Calcio fosfato tribasico 243
Calcio gluconato 240
Calcio idrossido 240
Calcio lattato 241
Calcio nitrato tetraidrato 241
Calcio ossido, polvere 242
Calcio ossido, pezzi 242
Calcio pantotenato 242
Calcio solfato biidrato 244
Calcio solfato emidrato 245
Calcio solfato emidrato soluzione 245
Calcio stearato 244
Calcone 245
Calmagite 246
Canfora naturale 246
Canfora sintetica 246
Caolino lavato 443
Carbone attivo 255
Carbone decolorante 255
Carbone vegetale 256
Carborundo, polvere 248
Carborundo, granulare 248
Carminio indaco secco 422
Carminio indaco soluzione 422
Carrez reattivo sale di potassio 249
Carrez reattivo sale di zinco 249
Cartina di amido di potassio 608
Carta al bromuro mercurico 484
Cartine indicatrici di pH 420
Carta di manganese d'argento 655
Cartina piombo (II) acetato 457
Cartina rosso Congo 273
Cartine tornasole 466
Caseina 249
Cellulosa, polvere 250
Cerio standard soluzione 250
Cerio ammonio nitrato ico 251
Cerio ammonio nitrato ico 0.1 mol/l 251
Cerio ammonio nitrato ico 0.01 mol/l 251
Cerio ammonio solfato ico diidrato 252
Cerio ammonio solfato ico 0.1 mol/l 252
Cerio ammonio solfato ico 0.01 mol/l 252
Cerio nitrato oso esaidrato 251
Cerio solfato ico tetraidrato 252
Cerio solfato ico 0.1 mol/l 253
Cesio standard soluzione 253
Cesio cloruro 254
Cesio cloruro 25 g/L soluzione 254
Cesio solfato 254
Cherosene 449
Chinidina solfato 634
Chinolina 634
Cianuri standard soluzione 285
Cicloesano 286
Cicloesano-d12 287
Cicloesanolo 287
Cicloesanone 288
Cicloesilammina 288
Ciclopentil-metil-etere 289
Cinconina 268
L-Cisteina 289
L-Cistina 289
Citrale 268
Cloralio idrato 256
Cloramina T sale sodico 256
Cloranile 256
Clorati standard soluzione 257
Cloriti standard soluzione 257
2-Cloroacetammide 258
p-Clorobenzaldeide 258
Clorobenzene 258
Clorobutanolo emiidrato 259
4-Clorofenolo 264
Cloroformio 259
Cloroformio-d 263
Cloroformio-d + 0.03% TMS 263
Cloruri standard soluzione 257
Cobalto standard soluzione 270
Cobalto acetato oso tetraidrato 271
Cobalto ammonio solfato oso esaidrato 271
Cobalto cloruro oso esaidrato 272
Cobalto (II) cloruro in soluzione 272
Cobalto nitrato oso esaidrato 272
Cobalto solfato oso eptaidrato 273
Cotone piombo (II) acetato 457
m-Cresolo 282
o-Cresolo 282
Crisoidina Y 268
Cromati standard soluzione 264
Cromo standard soluzione 264
Cromo cloruro ico esaidrato 265
Cromo nitrato ico nonaidrato 266
Cromo (III) di potassio solfato dodecaidrato 266
Cromo (VI) ossido 267
Cromo (VI) ossido soluzione 267
Cromo sesquiossido 266
Cromo solfato ico 267
Cupri-citrico soluzione 284
Cuprietilendiammina idrossido 1 mol/l 285
Cuprietilendiammina soluzione 285
Cupri-tartarico soluzione 285

D

Decalcificante 290
Decano 290
Densità standard 292
Deuterio ossido-d2 293
Deuterio ossido-d2 + 0.5% TSP d4 294
Deuterio ossido-d2 + 0.03% TSP d4 294
Deuterio ossido-d2 + 0.01% DMSO 294
Diacetildiossima 295
Diacetildiossima soluzione 1% in alcole etilico 296
Diacetildiossima sale sodico 295
Diacetonalcole 295
Dibenzoile perossido 296
Dibenzoilmetano 296
n-Dibutilftalato 297
Dicalite 4158 297
1,2-Diclorobenzene-d4 299
m-Diclorobenzolo 298
o-Diclorobenzolo 298
p-Diclorobenzolo 298
2,6-Diclorochinone-4-clorimmide 305
2,6-Diclorofenoloindofenolo sale sodico 305
Diclorometano 300
Diclorometano acidificato con acido cloridrico 1% 305
Diclorometano-d2 305
Dietanolammina 306
Dietilammina 306
Dietilcarbonato 307

INDICE ALFABETICO

n,n-Dietil-p-fenilendiammina solfato 310
Dietiletere 308
Dietilftalato 310
Difenilammina soluzione 1% in acido solforico 322
4-Difenilammina solfonato sodico 322
sim-Difenilcarbazide 323
sim-Difenilcarbazone 323
Difeniltiocarbazone 323
Differenziatore per kit Gram-Hucker 311
Digitonina 311
Diisopropilammina 311
Diisopropiletere 312
N,N-Diisopropiletilammina 312
Dimedone 313
N,N-Dimetilacetammide 313
p-Dimetilamminobenzaldeide 314
p-Dimetilamminobenzalrodanina 314
N,N-Dimetilformammide 314
N,N-Dimetilformammide-d7 316
N,N-Dimetilformammide dimetilacetale 317
N,N'-Dimetilpropilene urea 317
Dimetilsolfossido 317
Dimetilsolfossido-d6 319
Dimetilsolfossido-d6 + 0.03% TMS 319
1,2-Dimetossietano 313
Dimidio bromuro 320
2,4-Dinitroclorobenzene 320
2,4-Dinitrofenilidrazina (con 30% di acqua) 320
1,4-Diossano 321
1,3-Diossolano 322
2,2'-Dipiridile 323
Disprosio standard soluzione 325
Ditioossammide 324
Dodeciltrimetilammonio bromuro 325

E

Ehrlich reattivo 327
Eluente sodio bicarbonato 327
Eluente sodio carbonato 327
Eluente sodio carbonato/sodio bicarbonato 328
Emallume reattivo secondo Mayer 378
Emallume reattivo soluzione secondo Carazzi 378
Ematossilina 379
Ematossilina soluzione secondo Mayer 379
Eosina B 328
Eosina Y 328
Eosina Y 1% soluzione acquosa 328
Eosina Y 0.5% soluzione alcolica 329
n-Eptano 99% 380
n-Eptano 382
Eptano miscela di isomeri 382
Erbio standard soluzione 329
Eriocromo cianina R 330
Eritrosina extra B 330
Esadecano 383
Esafluoro-2-propanolo 384
Esametildisilazano 384
Esametilentetrammina 384
iso-Esano 438
n-Esano 99% 385
n-Esano 386
Esano miscela di isomeri 388
Esbach reattivo 330
1,2-Etano dicloro 299
Etanolammina 340
Etanolo assoluto anidro 331
Etanolo 96° 334
Etanolo 70% v/v 336
Etanolo 60% v/v 337
Etanolo 50% v/v 337
Etanolo 95%-Etere etilico 50/50(m/m) - 30 mg/l fenolftaleina 507
Etanolo assoluto denaturato 338
Etanolo 95° denaturato 338
Etanolo 94° denaturato 339
Etanolo 90° denaturato 339
Etanolo 70° modificato 339
Etanolo anidro-d6 340
Etere di petrolio (Ligroina) 100 - 140°C 564
Etere di petrolio (Ligroina) 80 - 120°C 564
Etere di petrolio (Ligroina) 75 - 120°C 565
Etere di petrolio (Benzina) 80 - 100°C 564
Etere di petrolio (Benzina) 55 - 85°C 565
Etere di petrolio (Benzina) 60 - 80°C 565
Etere di petrolio 40 - 70°C 566
Etere di petrolio 40 - 65°C 566
Etere di petrolio 40 - 60°C 567
Etere di petrolio 35 - 60°C 567
Etere di petrolio 30 - 50°C 568
Etere di petrolio 30 - 40°C 569
Etile acetato 340
Etile acetoacetato 343
Etile formiato 350
Etilendiammina 344
Eugenolo 351
Eukitt 351
Europio standard soluzione 351

F

Farina fossile composta 449
Fehling reattivo soluzione A 353
Fehling reattivo soluzione B 353
o-Fenantrolina monoidrata 569
o-Fenantrolina - Ferro solfato oso soluzione in acido solforico 569
L-Fenilalanina 572
Fenilidrazina cloridrato 573
Fenilidrazina cloridrato soluzione 573
Fenolftaleina 571
Fenolftaleina soluzione 1% in alcole etilico 571
Fenolftaleina soluzione 0.1% in alcole etilico 572
Fenolftaleina soluzione 0.1% 572
Fenolo 570
Fenolo liquido 85% 570
Ferricianuro standard soluzione 353
Ferro, polvere 427
Ferro, ridotto dall'idrogeno 428
Ferro standard soluzione 428
Ferro ammonio ossalato ico 432
Ferro ammonio solfato ico dodecaidrato 432
Ferro ammonio solfato ico soluzione 33% in acido nitrico 432
Ferro ammonio solfato ico 0.1 mol/l 433
Ferro ammonio solfato ico soluzione 100 g/l 433
Ferro ammonio solfato oso esaidrato 429
Ferro ammonio solfato oso 0.1N 430
Ferro ammonio citrato rosso 431
Ferro ammonio citrato verde 431
Ferro citrato ico 434
Ferro cloruro ico anidro sublimato 433
Ferro cloruro ico soluzione 4.5% 434
Ferro cloruro ico esaidrato 433
Ferro cloruro oso tetraidrato 430
Ferro cloruro in soluzione 434
Ferro ossido ico 435
Ferro nitrato ico nonaidrato 435
Ferro solfato ico 435
Ferro solfato oso eptaidrato 430
Ferro solfato oso 0.1mol/l 431
Ferrocianuro standard soluzione 354
Ferroina 0.025 mol/l soluzione 354
Ferron 354
Florisil 100-200 mesh 357
Florisil 60-100 mesh 356
Floroglucina 574
Floxina B 574
Fluoresceina 357
Fluoresceina sodica 357
Fluoruri standard soluzione 358
Folin-Ciocalteu reattivo 358
Formammide 362
Fosfati standard soluzione 574
Fosforo standard soluzione 577
D(-)Fruttosio 364
o-Ftalaldeide 579
Fucsina acida 364
Fucsina basica 365
Fucsina decolorata soluzione 365

G

Gadolinio standard soluzione 366
D(+)Galattosio 366
Gallio standard soluzione 367
Gel di silice 60A 70 - 200µ 650
Gel di silice 60A 35 - 70µ 649
Gel di silice 60A 40 - 63µ 649
Gel di silice 60A 20 - 45µ 648
Gel di silice 60A 6 - 35µ 648
Gel di silice 60A 0,06÷0,20 mm 650
Gel di silice granulare 650
Gel di silice granulare con indicatore esente da cobalto 651
Gelatina 368
Germanio standard soluzione 369
Giallo alizarina R 150
Giallo Clayton 270
Giallo metanile 487
Giallo metile 506
Giallo naftolo S 518
Giallo Sudan 727
Giemsa reattivo 369
Glicerina (30°Bé) 372
Glicerina 90% (28° Bé) 373
Glicocolla 373
Glicol dietilenico 307
Glicol dietilenico dimetiletere 307
Glicol etilenico 348
Glicol propilenico 629
Glicol trietilenico 773
Gliossale standard soluzione 374
Gliossale-bis(2-idrossianile) 374
D(+)Glucosio anidro 370
D(+)Glucosio monoidrato 371
Gomma arabica 377
Kit Gram-Hucker 376
Grasso di lana anidro 453
Griess reagente A 376
Griess reagente B 376
Griess reattivo soluzione in acido acetico 376

H

Hanus reattivo 379
Histolemon 390
Hyamine 1622 391
Hyamine 1622 soluzione 0.004M 392

I

Idrazina dicloridrato 392
Idrimeter Erba Soluzione A 419
Idrimeter Erba Soluzione B 419
Idrimeter Erba Indicatore C 419
Idrindantina 392
Idrochinone 416
p-IdrossibenzAldeide 416
Idrossilammina solfato 417
Idrossilammina soluzione, alcolico 417
Imidazolo 419
Indicatore misto di metile arancione 502
Indicatore misto di rosso di metile 504
Indicatore per azoto ammoniacale soluzione 420
Indicatore per iodometria 420

DASITGROUP | CARLO ERBA REAGENTI

Product specifications are subject to changes. Please visit our website for updates.

Indicatore universale pH 1-11 soluzione idroalcolica 421
Indicatore universale pH 1-11 soluzione in acqua 422
Indicatore universale pH 0-5 soluzione idroalcolica 421
Indice di rifrazione standard 637
Indio standard soluzione 422
Indolo 423
Inositolo 423
Inulina 424
Iodio bisublimato 425
Iodio 424
Iodio 0.5 mol/l (1N) 425
Iodio 0.05 mol/l (0.1N) 425
Iodio 0.01 mol/l (0.02N) 426
Iodio 0.005 mol/l (0.01N) 426
Iodio 10 ppm 426
Iodo bromuro soluzione 426
Iodoformio 427
Ioduri standard soluzione 424
Iridio standard soluzione 427
Isoamilo acetato 436
Isobutanolo 437
Isobutile acetato 438
Isoottano 439
Isopar G 440
Isopentano 441
Isopropilammina 442
Isopropile acetato 441
Istamina bicloridrato 389
L-Istidina monocloridrato monoidrato 390
Itterbio standard soluzione 794
Ittrio standard soluzione 794

K

Karl Fischer reattivo 2 mg H2O/ml monocomponente 443
Karl Fischer reattivo 5 mg H2O/ml monocomponente 443
Karl Fischer reattivo 5 mg H2O/ml monocomponente per aldeidi e chetoni 443
Karl Fischer solvente monocomponente per aldeidi e chetoni 444
Karl Fischer solvente monocomponente per olii 444
Karl Fischer titolante bicomponente 2 mg H2O/ml 444
Karl Fischer titolante bicomponente 5 mg H2O/ml 444
Karl Fischer solvente bicomponente 445
Karl Fischer solvente bicomponente tamponato 445
Karl Fischer solvente bicomponente per aldeidi e chetoni - senza metanolo 445
Karl Fischer solvente bicomponente per olii 445
Karl Fischer soluzione anolita 446
Karl Fischer soluzione anolita 446
Karl Fischer soluzione anolita per aldeidi e chetoni 446
Karl Fischer soluzione anolita per olii 447
Karl Fischer soluzione anolita, con fornetto 446
Karl Fischer soluzione catolita 447
Karl Fischer soluzione catolita per aldeidi e chetoni 447
Karl Fischer reagente per celle con e senza diaframma 447
Karl Fischer standard 0.10 mg/g acqua 449
Karl Fischer standard 1.0 mg/g acqua 448
Karl Fischer standard 10.0 mg/g acqua 448
Karl Fischer standard 5.0 mg/ml acqua 448
Kjedahl catalizzatore per analisi dell'acqua 450
Kjeldahl antischiuma 449
Kjeldahl catalizzatore 450
Kjeldahl catalizzatore al selenio 450
Kjeldahl catalizzatore al titanio 451
Kjeldahl catalizzatore sec. Wieninger 450
Kovac reattivo per indolo 451

L

Lana di vetro 370
Lantanio standard soluzione 453
Lantanio cloruro 25 g/l soluzione 454
Lantanio nitrato esaidrato 454
Lantanio nitrato soluzione 50 g/l 454
Lantanio nitrato 0.1 mol/l 454
Lantanio ossido 455
Lattofenolo d'Amman soluzione 160
Lattosio monoidrato 452
Lecitina di soia 461
Lega Devarda 294
Lega Raney 635
L(+)Leucina 461
Lievito di birra secco 794
Liquido fissatore AFA 354
Liquido fissatore di Bouin 355
Liquido fissatore di Bouin Hollande 355
Liquido fissatore di Davidson 355
Liquido fissatore FIXALL-HIS 356
Liquido fissatore senza acido acetico 356
L(+)Lisina monocloridrato 467
Litio standard soluzione 462
Litio carbonato 463
Litio cloruro 463
Litio idrossido anidro 464
Litio idrossido monoidrato 464
Litio idruro 464
Litio metilato 0.1 mol/l (0.1N) 464
Litio nitrato 465
Litio solfato monoidrato 465
Litio tetraborato anidro 465
Lugol "forte" soluzione concentrata 466
Lugol reattivo soluzione iodo-iodurata 466
Lugol soluzione in acqua per kit Gram-Hucker 466
Luminol 467
Lutezio standard soluzione 467

M

Magnesio, nastro 468
Magnesio, polvere 468
Magnesio, tornitura 468
Magnesio standard soluzione 468
Magnesio acetato tetraidrato 470
Magnesio carbonato basico 470
Magnesio cloruro esaidrato 471
Magnesio fosfato bibasico triidrato 472
Magnesio idrossido 472
Magnesio glicerofosfato 471
Magnesio nitrato esaidrato 472
Magnesio nitrato 10 g/l soluzione 473
Magnesio ossido 473
Magnesio ossido pesante 473
Magnesio perclorato 474
Magnesio perossido 474
Magnesio solfato anidro 475
Magnesio solfato eptaidrato 475
Magnesio stearato 474
Maltosio monoidrato 478
Manganese elettrolitico 479
Manganese standard soluzione 479
Manganese acetato oso tetraidrato 480
Manganese biossido 481
Manganese cloruro oso tetraidrato 480
Manganese solfato oso monoidrato 481
D-Mannitolo 481
May Grünwald reattivo 482
Mayer reattivo soluzione in acqua 482
L-Mentolo 483
2-Mercaptoetanolo 483
Mercurio standard soluzione 484
Mercurio dicloruro 485
di-Mercurio dicloruro 485
Mercurio cloruro soluzione 54 g/l 485
Mercurio ioduro ico 486
Mercurio ossido rosso 486
Mercurio solfato ico 486
Mercurio solfato soluzione 486
Metanolo 488
Metanolo + 0.1% v/v acido formico 492
Metanolo + 0.1% v/v acido trifluoroacetico 493
Metanolo, cloridrico 493
Metilcicloesano 497
Metile acetato 495
Metile benzoato 496
Metile ioduro 500
Metile isobutilchetone 500
Metile p-ossibenzoato 499
N-Metile-2-pirrolidone 502
Metile salicilato 505
3-Metile-2-benzotiazolinone idrazone cloridrato 497
Metiletilchetone 350
Metilisoamilchetone 500
2-Metiltetraidrofurano 505
4-Metiltetraidropirano 506
DL-Metionina 494
2-Metossietanolo 495
Miristato di isopropile 442
Miscela alcole-etere 30/70 m/m 507
Miscela C.H.M. 507
Miscela Etanolo assoluto/Etere Etilico 50/50 (v/v) 508
Miscela Etanolo assoluto/Etere Etilico 50/50 (p/p) 507
Miscela Etanolo 99° / Isopropanolo 508
Miscela Etanolo 95° / Isopropanolo 508
Miscela ottano 80 537
Miscela per il controllo delle saldature degli imballaggi 509
Miscela per la determinazione dell'indice di bromo 508
Miscela selenica 646
m-Nitrobenzaldeide 534
Molibdeno standard soluzione 510
Morfolina 512
Mucicarminio soluzione idroalcolica 512
Muresside 516

N

Naftalene 517
N-(1-Naftil)etilendiammina dicloridrato 518
alfa-Naftolbenzeina 518
alfa-Naftolbenzeina soluzione 0.2% in acido acetico 518
1-Naftolo 517
2-Naftolo 517
Neocuproina cloridrato 519
Neodimio standard soluzione 519
Nero eriocromo T 329
Nero Sudan B 726
Nessler reattivo soluzione A 520
Nessler reattivo soluzione B 521
Nessler reattivo soluzione unica 520
Nichel standard soluzione 521
Nichel, polvere 521
Nichel acetato oso tetraidrato 522
Nichel ammonio solfato oso 523
Nichel carbonato oso basico 523
Nichel cloruro oso esaidrato 523
Nichel nitrato oso esaidrato 524
Nichel nitrato oso 10 g/l 524
Nichel solfato oso esaidrato 524
Nicotinammide 525
Nigrosina 525

INDICE ALFABETICO

Ninidrina 526
Ninidrina solzione 526
Ninidrina e stagno cloruro reattivo 526
Niobio standard soluzione 526
Nital soluzione 4% 527
Nitrati standard soluzione 527
Nitriti standard soluzione 534
Nitrobenzene 535
p-Nitrofenolo 535
p-Nitrofenolo soluzione 0.1% in acqua 535
Nonilfenolo etossilato 10 ETO 536

O

Olio d'immersione 420
Olio di legno cedro condensato 538
Olio di mandorla raffinato 539
Olio di ricino 249
Olio di vaselina 555
Olmio standard soluzione 390
Olmio ossido 391
Olmio perclorato soluzione 391
Orceina 540
Orcina monoidrato 540
Oro standard soluzione 375
Oro (III) cloruro triidrato 375
Osmio standard soluzione 549
Osmolalità Standards 550
Osmolalità standard a base di proteine 550
Osmolalità Standards a base di urina 550
8-Ossichinolina 417
n-Ottano 537

P

Palladio standard soluzione 553
Palladio cloruro oso 553
Palladio nitrato 2 g/L soluzione 554
Papanicolaou Ematossilina soluzione secondo Harris 554
Papanicolaou soluzione EA 50 554
Papanicolaou soluzione OG 6 555
Paraffina 56-58°C Erbaplast (senza DMSO) 555
Paraformaldeide 556
Pararosanilina soluzione, decolorata 556
n-Pentano 99% 557
n-Pentano 558
Pepsina HCl 560
Perossido di idrogeno soluzione 40% m/v 413
Perossido di idrogeno soluzione 35% 413
Perossido di idrogeno soluzione 30-32% 413
Perossido di idrogeno soluzione 30% 414
Perossido di idrogeno soluzione 6% 415
Perossido di idrogeno soluzione 3.5% m/v 415
Perossido di idrogeno soluzione 3% 416
Petrolio 563
Piombo 455
Piombo standard soluzione 455
Piombo acetato basico 457
Piombo acetato basico soluzione 457
Piombo acetato soluzione 95 g/l 458
Piombo di-acetato triidrato 458
Piombo biossido 460
Piombo nitrato 458
Piombo nitrato 0.1 mol/l 459
Piombo nitrato 0.05 mol/l 459
Piombo nitrato soluzione 33 g/l 459
Piombo ossido 460
Piombo solfato 460
Piperidina 580
1-(2-Piridile-azo)-2-naftolo 632
Piridina 630
Piridina-d5 631
Pirocatecolo 632
Pirogallolo 632
Platino standard soluzione 580
Polivinilpirrolidone 581
Pomice 630
Porpora bromocresolo 212
Porpora bromocresolo soluzione 0.4% in etanolo 212
Porpora bromocresolo soluzione 213
Porpora m-cresolo 282
Potassio standard soluzione 582
Potassio acetato 583
Potassio bicarbonato 584
Potassio bicromato 592
Potassio bicromato - Soluzione di acido solforico 592
Potassio bicromato 0.167 mol/l (1 N) 592
Potassio bicromato 0.1414 g/l 594
Potassio bicromato 0.0417 mol/l (0.25N) 593
Potassio bicromato 0.0167 mol/l (0.1 N) 593
Potassio bicromato 0.04 mol/l (0.24N) in 80 g/l HgSO4 593
Potassio bicromato soluzione 0.5% 592
Potassio bicromato soluzione 106 g/l 594
Potassio bisolfato 584
Potassio bromato 584
Potassio bromato 0.033 mol/l (0.198N) 585
Potassio bromato 0.02 mol/l (0.12N) 585
Potassio bromato 0.0167 mol/l (0.1N) 585
Potassio bromuro 585
Potassio carbonato 586
Potassio citrato tribasico monoidrato 591
Potassio cloruro 586
Potassio cloruro 3.5 mol/l (3.5N) 587
Potassio cloruro 3.5 mol/l (3.5N) e argento cloruro 587
Potassio cloruro 3 mol/l (3N) 588
Potassio cloruro 3 mol/l (3N) e argento cloruro 588
Potassio cloruro 3 mol/l (3N) soluzione acqua-glicerina 588
Potassio cloruro 1 mol/l (1N) 588
Potassio cloruro 0.1 mol/l (0.1N) 589
Potassio cloruro 0.01 mol/l (0.01N) 589
Potassio cloruro soluzione 590
Potassio cloruro soluzione satura 590
Potassio cloruro soluzione 25g/l in acido cloridrico 589
Potassio cloruro soluzione 12g/l 589
Potassio cromato 590
Potassio cromato soluzione 10% 590
Potassio cromato soluzione 5% 591
di-Potassio esacloroplatinato 597
Potassio ferricianuro 594
Potassio ferricianuro soluzione 53 g/l 596
Potassio ferrocianuro triidrato 595
Potassio ferrocianuro soluzione 10% 595
Potassio fluoruro 596
Potassio fluoruro biidrato 596
Potassio fosfato bibasico anidro 614
Potassio fosfato bibasico triidrato 614
Potassio fosfato monobasico 615
Potassio fosfato monobasico 0.2 mol/l (0.2N) 615
Potassio ftalato acido 597
Potassio ftalato acido 0.2 mol/l (0.2N) 598
Potassio ftalato acido 0.1 mol/l (0.1N) 598
Potassio idrossido, gocce 599
Potassio idrossido, scaglie 599
Potassio idrossido soluzione 45% 600
Potassio idrossido soluzione 38% (40° Bé) in acqua 601
Potassio idrossido soluzione 33% 601
Potassio idrossido soluzione 28% 601
Potassio idrossido soluzione 3% in etanolo 601
Potassio idrossido 2 mol/l (2N) in etanolo 602
Potassio idrossido 1 mol/l (1N) 602
Potassio idrossido 0.5 mol/l (0.5N) 602
Potassio idrossido 0.5 mol/l (0.5N) in alcole etilico 603
Potassio idrossido 0.5 mol/l (0.5N) in alcole metilico 603
Potassio idrossido 0.46 mol/l (0.46N) 604
Potassio idrossido 0.25 mol/l (0.25N) 604
Potassio idrossido 0.23 mol/l (0.23N) 604
Potassio idrossido 0.1 mol/l (0.1N) 604
Potassio idrossido 0.1 mol/l (0.1N) in alcole etilico 605
Potassio idrossido 0.1 mol/l (0.1N) in alcole metilico 605
Potassio idrossido in soluzione 606
Potassio iodato 606
Potassio iodato 0.05 mol/l (0.3N) 606
Potassio iodato 0.0167 mol/l (0.1N) 606
Potassio iodato 0.00167 mol/l (0.01N) 607
Potassio iodato acido 597
potassio iodobismutato soluzione 609
Potassio ioduro 607
Potassio ioduro soluzione 10% 608
Potassio ioduro soluzione 3.9% 608
Potassio ioduro soluzione 608
Potassio metabisolfito 609
Potassio nitrato 609
Potassio nitrato 1 mol/l (1N) 610
Potassio nitrito 610
Potassio ossalato monoidrato 610
Potassio (meta)-periodato 611
Potassio permanganato 611
Potassio permanganato soluzione 3% 612
Potassio permanganato 0.2 mol/l (1N) 612
Potassio permanganato 0.02 mol/l (0.1N) 612
Potassio permanganato 0.002 mol/l (0.01N) 613
Potassio permanganato e acido fosforico soluzione 613
Potassio persolfato 613
Potassio piroantimoniato acido 615
Potassio piroantimoniato soluzione 616
Potassio pirofosfato 616
Potassio pirosolfato 616
Potassio sodio tartrato tetraidrato 617
Potassio solfato 618
Potassio solfocianuro 620
Potassio solfocianuro 0.1 mol/l (0.1N) 620
Potassio solfocianuro soluzione 621
Potassio solfocianuro soluzione 5% 620
Potassio solfoguaiacolato 596
Potassio sorbato 617
Potassio tartrato 618
Potassio tartrato acido 619
Potassio tellurito 619
Potassio tetraiodomercurato soluzione alcalina 619
Praseodimio standard soluzione 621
L(-)Prolina 622
Propan-1-olo 622
Propan-2-olo 623
Propan-2-olo 70% 626
1,3-Propandiolo 627
n-Propile acetato 628
Propile p-ossibenzoato 629
Propilene carbonato 629
Punto di fusione standard 483

R

Raffinosio 635
Rame elettrolitico lastra ricotto 274
Rame elettrolitico, filo 275
Rame elettrolitico, tornitura 274
Rame ridotto, polvere 275
Rame standard soluzione 275
Rame acetato ico idrato 277
Rame carbonato basico ico 278
Rame cloruro ico diidrato 278

DASITGROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

Rame cloruro oso 276
Rame cloruro oso soluzione 7% in ammoniaca 277
Rame ioduro oso 277
Rame nitrato ico triidrato 279
Rame ossido ico 279
Rame ossido filo 279
Rame solfato ico anidro 280
Rame solfato ico pentaidrato 280
Rame solfato ico soluzione 12.5% 281
Rame (II) solfato in soluzione 281
Rame tetrammina soluzione ammoniacale 165
L(+)Ramnosio 638
Reattivo fosfomolibdotungstico 576
Reattivo iodoplatinato 427
Reattivo molibdovanadico 511
Reattivo solfomolibdico 729
Reattivo TAN 635
Reattivo TBN ASTM D2896 636
Reattivo TBN ASTM D4739 636
Reattivo titanio tricloruro-acido solforico 765
Reattivo unico per lipolisi 635
Renio standard soluzione 638
Resorcina 637
Rhodol soluzione 639
Riboflavina 639
D(-)Ribosio 639
Rodio standard soluzione 638
Rosso alizarina 149
Rosso chinaldina 634
Rosso Congo 273
Rosso Congo soluzione 274
Rosso fenolo 570
Rosso fenolo soluzione 0,2% in alcole etilico 571
Rosso fenolo solzione 571
Rosso metile 503
Rosso metile soluzione 504
Rosso metile soluzione 0.2% in alcole etilico 503
Rosso metile soluzione 0,1% in alcole etilico 504
Rosso neutro 521
Rosso o-cresolo 283
Rosso o-cresolo soluzione 0.2% in alcole etilico 283
Rosso per olio O 636
Rosso Ponceau BS 581
Rosso Ponceau S 581
Rosso solido nucleare 536
Rubidio standard soluzione 640
Rutenio standard soluzione 641

S

Sabbia di Fontainebleau 644
Sabbia purificata 644
D(+)-Saccarosio 726
Safranina T 642
Safranina T soluzione idroalcolica per kit Gram-Hucker 642
Salda d'amido soluzione 1% 721
Sale di Reinecke 637
Salicilaldeide azina 642
Samario standard soluzione 643
Sapone soluzione in alcole etilico 660
Saponina 644
Scandio standard soluzione 645
Schiff reattivo per colorazione PAS 646
Schiff reattivo soluzione per Aldedi 645
Selenio, polvere 646
Selenio standard soluzione 647
Selenio biossido 648
Setacci molecolari 4 A 510
Setacci molecolari 3 A 510
Silicio standard soluzione 651
Silicio biossido 652
Sodio standard soluzione 661
Sodio acetato anidro 662
Sodio acetato triidrato 662
Sodio alginato 663
Sodio alluminato 663
Sodio ammonio idrogeno fosfato 664
Sodio arsenito 0.1 mol/l (0.2N) 664
Sodio arsenito 0.05 mol/l (0.1N) 664
Sodio L-ascorbato 664
Sodio azide 665
Sodio benzoato 665
Sodio bicarbonato 665
Sodio bisolfato monoidrato 666
Sodio boroidruro 666
Sodio bromuro 667
Sodio carbonato anidro 667
Sodio carbonato decaidrato 668
Sodio carbonato monoidrato 668
Sodio carbonato soluzione 20% 670
Sodio carbonato 0.5 mol/l (1N) 669
Sodio carbonato 0.05 mol/l (0.1N) 669
Sodio carbonato soluzione 670
Sodio cianoboroidruro 674
Sodio citrato bibasico 672
Sodio citrato tribasico anidro 672
Sodio citrato tribasico biidrato 673
Sodio cloruro 670
Sodio cloruro 5 mol/l (5N) 672
Sodio cloruro 0.1 mol/l (0,1N) 672
Sodio cobalto nitrito 673
Sodio dicloroisocianurato biidrato 674
Sodio dietilditiocarbammato triidrato 674
Sodio esafluorosilicato 677
Sodio esametafosfato 677
Sodio fluoruro 675
Sodio formiato 675
Sodio fosfato bibasico anidro 700
Sodio fosfato bibasico biidrato 701
Sodio fosfato bibasico dodecaidrato 702
Sodio fosfato monobasico biidrato 703
Sodio fosfato monobasico monoidrato 703
Sodio fosfato tribasico dodecaidrato 704
Sodio o-fosfito pentaidrato 704
Sodio glicerofosfato pentaidrato 676
Sodio gluconato 676
Sodio glutammato acido 676
Sodio idrosolfito 678
Sodio idrossido, gocce 679
Sodio idrossido, perline 678
Sodio idrossido soluzione 50% 680
Sodio idrossido soluzione 40% 680
Sodio idrossido soluzione 35-37% 681
Sodio idrossido soluzione 35% 681
Sodio idrossido soluzione 32% 681
Sodio idrossido soluzione 30% 682
Sodio idrossido soluzione 20% p/p 683
Sodio idrossido soluzione 20% p/v 683
Sodio idrossido soluzione 10% p/v 683
Sodio idrossido soluzione 5% p/v 684
Sodio idrossido 6 mol/l (6N) 684
Sodio idrossido 5 mol/l (5N) 685
Sodio idrossido 4 mol/l (4N) 685
Sodio idrossido 3 mol/l (3N) 685
Sodio idrossido 2 mol/l (2N) 686
Sodio idrossido 1.2 mol/l (1.2N) 686
Sodio idrossido 1 mol/l (1N) 686
Sodio idrossido 0.7 mol/l (N/1.4) 687
Sodio idrossido 0.5 mol/l (0.5N) 688
Sodio idrossido 0.357 mol/l (0.357N) 688
Sodio idrossido 0.25 mol/l (0.25N) 688
Sodio idrossido 0.2 mol/l (0.2N) 689
Sodio idrossido 0.1 mol/l (0.1N) 689
Sodio idrossido 0.1 mol/l (N/10) in etanolo 690
Sodio idrossido 0.01 mol/l (0.01N) 690
Sodio idrossido 1/9 mol/l (N/9) 689
Sodio idrossido soluzione 691
Sodio idrossido soluzione, metanolica 691
Sodio idrossido su silice 680
Sodio idrossido-d 30% 691
Sodio idrossido-d 1 mol/l 691
Sodio ioduro 693
Sodio ipoclorito soluzione 12.5% 692
Sodio ipoclorito soluzione in acqua 692
Sodio ipofosfito 693
Sodio laurilsolfato 693
Sodio metabisolfito 694
Sodio metafosfato 695
Sodio metossido 0.1 mol/l 695
Sodio molibdato diidrato 695
Sodio nitrato 696
Sodio nitrito 696
Sodio nitrito 0.1 mol/l (0.1N) 697
Sodio nitrito 500 g/l soluzione 697
Sodio nitroprussiato diidrato 697
Sodio ossalato 698
Sodio perborato tetraidrato 698
Sodio perclorato monoidrato 699
Sodio (meta)-periodato 694
Sodio perossido 699
Sodio persolfato 699
Sodio persolfato 1 mol/l 700
Sodio pirofosfato decaidrato 705
Sodio salicilato 705
Sodio silicato 705
Sodio solfato anidro 706
Sodio solfato decaidrato 707
Sodio solfito anidro 708
Sodio solfocianuro 709
Sodio solfuro nonaidrato 707
Sodio solfuro nonaidrato soluzione 707
Sodio solfuro triidrato 708
Sodio succinato esaidrato 706
Sodio tartrato acido 678
Sodio tartrato diidrato 709
Sodio tetraborato anidro 709
Sodio tetraborato decaidrato 710
Sodio tetrafenilborato 710
Sodio tiosolfato anidro 711
Sodio tiosolfato pentaidrato 711
Sodio tiosolfato 1 mol/l (1N) 712
Sodio tiosolfato 0.5 mol/l (0.5N) 712
Sodio tiosolfato 0.2 mol/l (0.2N) 712
Sodio tiosolfato 0.1 mol/l (0.1N) 712
Sodio tiosolfato 0.0394 mol/l (0.0394N) 713
Sodio tiosolfato 0.0197 mol/l (0.0197N) 713
Sodio tiosolfato 0.01 mol/l (0.01N) 713
Sodio tungstato diidrato 714
Solfanilammide 728
Solfati standard soluzione 728
Solfito standard soluzione 729
Soluzione ADF 147
Soluzione NDF 519
Soluzione NDF PLUS 519
Soluzione ossido di etilene 349
Soluzione Redox 468 mV a 25°C 637
Soluzione Redox 220 mV a 25°C 636
Soluzione Rosso cresolo/Sodio idrossido 0.1N / Etanolo 96%/ Acqua 283
Soluzione standard da 1.30 µS/cm 715
Soluzione standard da 5 µS/cm 716
Soluzione standard da 10 µS/cm 716
Soluzione standard da 20 µS/cm 716
Soluzione standard da 50 µS/cm 716
Soluzione standard da 84 µS/cm 717
Soluzione standard da 100 µS/cm 717
Soluzione standard da 147 µS/cm 717
Soluzione standard da 200 µS/cm 717
Soluzione standard da 500 µS/cm 717
Soluzione standard da 1000 µS/cm 718
Soluzione standard da 1413 µS/cm 718

INDICE ALFABETICO

Soluzione standard da 5000 µS/cm 718
Soluzione standard da 10000 µS/cm 718
Soluzione standard da 12880 µS/cm 718
Soluzione standard da 20000 µS/cm 719
Soluzione standard da 50000 µS/cm 719
Soluzione standard da 100000 µS/cm 719
Soluzione standard da 150000 µS/cm 719
Soluzione standard da 200000 µS/cm 720
Soluzione standard da 300000 µS/cm 720
Soluzione standard da 350000 µS/cm 720
Soluzione standard da 450000 µS/cm 720
Soluzione standard da 500000 µS/cm 721
Soluzione standard per analisi degli idrocarburi 546
Soluzioni standard per la ricerca di solventi residui 509
Soluzioni primarie per il grado di colorazione dei liquidi 290
Soluzioni standard per il grado di colorazione dei liquidi 291
Solvente Plus 714
D-Sorbitolo 715
Sorbitolo 714
Sorbitolo (non cristallizzabile) soluzione 70% 715
Sospensione primaria di opalescenza 621
Stagno standard soluzione 760
Stagno lastra 759
Stagno, polvere 759
Stagno cloruro ico pentaidrato 762
Stagno cloruro oso diidrato 761
Stagno cloruro oso soluzione 10% 761
Stagno cloruro oso soluzione 761
Stagno ossido ico 762
Stagno solfato oso 762
Standard del colore ASTM 187
Standard del colore Saybolt 644
Standard Gardner del colore 368
Standard multianione per cromatografia ionica 512
Standard multielemento per ICP 513
Standard multielemento per ICP e ICP-MS 515
Standard organico: PAH mix 542
Standard organico PCB 540
Standard organico: PCB mix 543
Standard organico: Pesticidi mix 546
Standard di viscosità 785
Stronzio standard soluzione 722
Stronzio acetato 723
Stronzio bromuro monoidrato 723
Stronzio carbonato 724
Stronzio cloruro esaidrato 724
Stronzio nitrato 724
Stronzio solfato 725
Sudan III 726
Sudan III soluzione idroalcolica satura 727
Sulfolano 729

T

Talco 744
Tallio standard soluzione 754
Tampone acetato pH 4.5 216
Tampone acetato pH 4.6 130
Tampone acetato pH 6.0 130
Tampone ammoniacale pH 10 165
Tampone borico soluzione 207
Tampone cloruro di ammonio pH 9.5 170
Tampone cloruro di ammonio pH 10.0 170
Tampone cloruro di ammonio pH 10.7 169
Tampone forza ionica totale 769
Tampone fosfato pH 2.0 576
Tampone fosfato pH 3.0 576
Tampone fosfato pH 5.5 575
Tampone fosfato pH 6.0 575
Tampone fosfato pH 6.8 575
Tampone fosfato pH 7.2 575
Tampone fosfato pH 7.4 575
Tampone fosfato pH 9.0 574
Tampone tris(idrossimetil)aminometano - EDTA pH 7.4 777
Tampone tris(idrossimetil)aminometano tampn soluzione pH 8.1 777
Tampone tris(idrossimetil)aminometano sodio cloruro pH 7.4 778
Tampone pH 1 216
Tampone pH 1.68 216
Tampone pH 2 217
Tampone pH 3 217
Tampone pH 3.5 217
Tampone pH 3.56 218
Tampone pH 3.7 218
Tampone pH 4 218
Tampone pH 4.5 219
Tampone pH 4.62 219
Tampone pH 5 219
Tampone pH 5.2 219
Tampone pH 6 220
Tampone pH 6.8 220
Tampone pH 6.88 220
Tampone pH 7 220
Tampone pH 7.2 221
Tampone pH 7.20 secondo Weise 221
Tampone pH 7.4 222
Tampone pH 8 222
Tampone pH 9 222
Tampone pH 9.22 223
Tampone pH 10 223
Tampone pH 10.06 224
Tampone pH 11 224
Tampone pH 12 225
Tampone pH 13 225
Tantalio standard soluzione 744
Tartrazina 746
Tauber reattivo 746
Tellurio standard soluzione 747
Tellurio pezzi 746
Tellurio polvere 747
Terbio standard soluzione 748
Tetrabutilammonio bisolfato 748
Tetrabutilammonio bromuro 749
Tetrabutilammonio idrossido soluzione 0.1 mol/l (0.1N) 749
Tetrabutilammonio idrossido 0.1 mol/l (0.1N) in isopropanolo 749
Tetrabutilammonio idrossido 0.1 mol/l (0.1N) in metanolo/propan-2-olo (50/50) 750
Tetracloroetano-d2 751
Tetracloroetilene 750
Tetraetilammonio bromuro 751
Tetraidrofurano 751
Tetraidrofurano-d8 753
Tetrametilammonio idrossido 10% 753
N,N,N',N'-Tetrametil-p-fenilendiammina bicloridrato 753
Timolftaleina 758
Timolftaleina 0.1% soluzione idroalcolica 759
Timolftaleina soluzione 0.1% in etanolo 759
Timolo 757
Tioacetammide 754
Tioacetammide soluzione 40 g/l 755
Tiourea 756
L(-)Tirosina 781
Tisab 762
Tisab III soluzione 763
Titanio biossido 764
Titanio isopropilato 764
Titanio standard soluzione 763
o-Tolidina soluzione 765
o-Tolidina soluzione 0.1% 765
Toluene 765
Toluene soluzione in esano 768
Toluene-d8 768
p-Toluensolfonammide 768
Torbidità std 4000NTU formazina 780
Torio standard soluzione 756
Tornasole 465
Triacetina 769
1,2,4-Triclorobenzene 771
2,2,2 - Tricloroetanolo 770
Tricresolo 771
Trietanolammina 771
Trietilammina 772
Trifenilfosfina 776
2,3,5-Trifeniltetrazolo cloruro 776
Trimetilcetilammonio bromuro 775
Trimetilcetilammonio p-toluensolfonato 775
N,O-Bis (trimetilsilil)acetammide 204
N,O-Bis (trimetilsilil)trifluoroacetammide 204
Triottilfosfina ossido 775
Tris idrossimetilamminometano 776
Tris idrossimetilamminometano soluzione 777
Tris idrossimetilamminometano cloridrato 778
Triton© X100 soluzione 778
Tropeolina O 778
Tropeolina O soluzione 0,1% 779
Tulio standard soluzione 756
Tungsteno standard soluzione 779
Turk reattivo 780

U

Uranio standard soluzione 782
Urea 782

V

Vanadio standard soluzione 783
Vanadio pentossido 784
Vanillina 784
Vanillina soluzione, fosforico 784
Vaselina bianca 556
Verde brillante 209
Verde bromocresolo 211
Verde bromocresolo soluzione 211
Verde bromocresolo soluzione 0.04% idroalcolica 211
Verde bromocresolo - Rosso metile soluzione 212
Verde luce SF 461
Verde malachite 476
Verde malachite soluzione 0.5% in acido acetico anidro 476
Verde metile 499
Verde solido FCF 353
Violetto cristalli 284
Violetto cristalli ossalato soluzione per kit Gram-Hucker 284
Violetto cristalli soluzione 0.5% in acido acetico anidro 284
Violetto genziana 368
Violetto genziana soluzione fenicata 368
Violetto pirocatechina 632

W

Wijs reattivo 790
Wright colorante soluzione in alcole metilico 790

X

o-Xilene 792
Xilene, isomeri misti 791
Xilencianolo FF 792
2,4-Xilenolo 792
D(+)Xilosio 793

Z

Ziehl -Neelsen reattivo soluzione 796
Zinco, attivato 797
Zinco, granuli 796

Product specifications are subject to changes. Please visit our website for updates.

Zinco, lastra	796
Zinco, polvere	797
Zinco standard soluzione	797
Zinco acetato diidrato	798
Zinco carbonato basico	799
Zinco cloruro anidro	799
Zinco cloruro soluzione 60%	799
Zinco cloruro soluzione, iodata	800
Zinco cloruro - acido formico soluzione	800
Zinco dibenzilditiocarbammato	800
Zinco nitrato esaidrato	800
Zinco ossido	801
Zinco solfato eptaidrato	802
Zinco solfato monoidrato	802
Zinco solfato 0.1 mol/l (0.2N)	802
Zinco solfato 0.05 mol/l (0.05N)	803
Zinco solfuro	803
Zinco stearato	801
Zirconile nitrato	805
Zirconio standard soluzione	804
Zirconio ossido	804
Zirconio polvere	803
Zolfo standard soluzione	730
Zolfo sublimato lavato ventilato	730

BEST CHOICE

Xcipharm™

Your requirements are our mission

The new product range for excipient use:



QUALITY ASSURANCE

- Flow chart
- Change control

PROCESS

- Validated cleaning procedures and dedicated equipment
- Labels consolidation
- Certificate of analysis with the name of the producer and the date of production of the raw material
- Stability studies

DOCUMENTATION

- BSE/TSE statement
- OGM statement
- Residual solvents statement
- ICH Q3D
- Risk assessment (2015/C95/02)

www.carloerbareagents.com

DASIT GROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

INDEX ALPHABÉTIQUE

A

Absorbant pour liquides répandus 129
Acétal 129
Acétamide 129
Acétanilide 129
Acétone 138
Acétone-d6 141
Acétone / eau 98/2 (v/v) + Bleu de Bromophénol 0.02 g/l 141
Acétonitrile 142
Acétonitrile-d3 145
Acétonitrile + 0.1% v/v d'acide formique 145
Acétonitrile + 0.1% v/v d'acide trifluoroacétique 145
Acétophénone 146
Acetylacétone 146
n-Acétyl-L-cystéine 146
p-Acétylphénétidine 146
Acide acétique glacial 130
Acide acétique 96% 133
Acide acétique 80% 134
Acide acétique 45% 134
Acide acétique 30% 134
Acide acétique 27% 135
Acide acétique 25% 135
Acide acétique 20% 135
Acide acétique 12% 135
Acide acétique 1 mol/l (1N) 136
Acide acétique 0.1 mol/l (0.1N) 136
Acide acétique 0.03 mol/l (0.03N) 136
Acide acétique-d4 137
Acide adipique 147
Acide p-aminobenzoique 158
Acide 1-amino-2-naphtol-4-sulfonique 158
Acide anisique 181
Acide L(+)ascorbique 186
Acide L(+)aspartique 186
Acide barbiturique 191
Acide benzoique 199
Acide borique 205
Acide borique 1% avec indicateur 207
Acide borique 3% 206
Acide borique 4% 206
Acide borique 4% avec indicateur 206
Acide borique 20g/l 207
Acide borique 20 g/l avec indicateur 207
Acide bromhydrique 48% 393
Acide 1-butanesulfonique sel sodique 225
Acide n-butyrique 231
Acide calconecarbonique 245
Acide n-caproïque 247
Acide n-caprylique 247
Acide chloranilique 257
Acide chlorhydrique 50% v/v 393
Acide chlorhydrique 37% 394
Acide chlorhydrique 34-37% 396
Acide chlorhydrique 32-35% 396
Acide chlorhydrique 32% (20°Bé) 397
Acide chlorhydrique 32% 397
Acide chlorhydrique 29 - 31 % 398
Acide chlorhydrique 26% 398
Acide chlorhydrique 25% m/v 398
Acide chlorhydrique 23% 399
Acide chlorhydrique 20% 399
Acide chlorhydrique 12% 399
Acide chlorhydrique 10% 400
Acide chlorhydrique 8% 400
Acide chlorhydrique 5% 400
Acide chlorhydrique 1.128% m/v 401
Acide chlorhydrique 9 mol/l (9N) 401
Acide chlorhydrique 6 mol/l (6N) 401
Acide chlorhydrique 5 mol/l (5N) 402
Acide chlorhydrique 4 mol/l (4N) 402
Acide chlorhydrique 3 mol/l (3N) 403
Acide chlorhydrique 2 mol/l (2N) 403
Acide chlorhydrique 1 mol/l (1N) 403
Acide chlorhydrique 0.714 mol/l (N/1.4) 404
Acide chlorhydrique 0.5 mol/l (0.5N) 405
Acide chlorhydrique 0.2 mol/l (0.2N) 405
Acide chlorhydrique 0.2 mol/l (0.2N) dans propanol-2 405
Acide chlorhydrique 0.1 mol/l (0.1N) 406
Acide chlorhydrique 0.1 mol/l (0.1N) dans l'éthanol 407
Acide chlorhydrique 0.1 mol/l (0.1N) dans propanol-2 407
Acide chlorhydrique 0.0714 mol/l (N/14) 407
Acide chlorhydrique 0.05 mol/l (0.05N) 407
Acide chlorhydrique 0.04 mol/l (0.04N) 408
Acide chlorhydrique 0.02 mol/l (0.02N) 408
Acide chlorhydrique 0.01 mol/l (0.01N) 408
Acide chlorhydrique dilué 409
Acide chlorhydrique-d 1 mol/l 409
Acide chlorhydrique-d 20% 409
Acide chlorhydrique bromé 409
Acide chloroacétique 258
Acide chromotropique sel disodique 268
Acide citrique monohydraté 269
Acide citrique anhydre 269
Acide 1-decanesulfonique sel sodique 290
Acide dichloroacétique 297
Acide diéthylènetriaminepentacétique 308
Acide 1-dodecanesulfonique sel sodique 324
Acide dodécylbenzènesulfonique sel sodique 325
Acide éthylènediaminetétracétique 344
Acide éthylènediaminetétracétique sel dipotassique dihydraté 345
Acide éthylènediaminetétracétique sel disodique 345
Acide éthylènediaminetétracétique sel disodique 0.1 mol/l (0.2N) 346
Acide éthylènediaminetétracétique sel disodique 0.05 mol/l (0.1N) 347
Acide éthylènediaminetétracétique sel disodique 0.01 mol/l (0.02N) 347
Acide éthylènediaminetétracétique sel de potassium et de magnésium dihydraté 347
Acide éthylènediaminetétracétique sel tétracétique tétrahydraté 348
Acide éthylènediaminetétracétique sel tripotassique 348
Acide éthylèneglycol bis-(2-aminoethylether) tétracétique (EGTA) 349
Acide fluorhydrique 50% 410
Acide fluorhydrique 47-51% 411
Acide fluorhydrique 39.5% 412
Acide fluorhydrique 0.1M 412
Acide fluorhydrique dilué 412
Acide formique 99% 362
Acide formique 85% 363
Acide formique-d 364
Acide fumarique 365
Acide gallique monohydraté 367
Acide L(+)glutamique 371
Acide glycolique 374
Acide heptafluorobutyrique 380
Acide 1-heptanesulfonique sel sodique 383
Acide hexachloroplatinique hexahydraté 383
Acide 1-hexanesulfonique sel sodique 389
Acide 1-hexanesulfonique sel sodique monohydrate 389
Acide hypophosphoreux 50% 418
Acide iodhydrique 57% 393
Acide L(+)-lactique 452
Acide maléique 477
Acide DL-malique 477
Acide malonique 478
Acide D(-)mandélique 478
Acide métaphosphorique 487
Acide méthanesulfonique 487
Acide alpha-méthoxyphénylacétique 495
Acide nicotinique 525
Acide nitrilotriacétique 534
Acide nitrique fumant 99% 528
Acide nitrique 69.5% 528
Acide nitrique 67-70% 529
Acide nitrique 67-69% 530
Acide nitrique 67.5% (42° Be) 530
Acide nitrique 65% 531
Acide nitrique 18% 532
Acide nitrique 10% 532
Acide nitrique 8 mol/l (8N) 532
Acide nitrique 2 mol/l (2N) 533
Acide nitrique 1 mol/l (1N) 533
Acide nitrique 0.1 mol/l (0.1N) 533
Acide nitrique dilué 534
Acide m-nitrobenzoique 535
Acide 1-octanesulfonique sel sodique 537
Acide octanesulfonique sel sodique monohydraté 538
Acide oléique 539
Acide oxalique dihydraté 550
Acide oxalique 0.5 mol/l (1N) 551
Acide oxalique 0.05 mol/l (0.1N) 551
Acide oxalique 0.005 mol/l (0.01N) 551
Acide 1-pentanesulfonique sel sodique 559
Acide 1-pentanesulfonique sel sodique monohydrate 560
Acide perchlorique 65-71% 560
Acide perchlorique 65% 561
Acide perchlorique 0.1 mol/l (0.1N) 561
Acide perchlorique 0.1 mol/l (0.1N) dans l'acide acétique 562
Acide perchlorique 0.05 mol/l (0.05N) 562
Acide perchlorique 0.01 mol/l (0.01N) 562
Acide perchlorique solution 562
Acide périodique 563
Acide phosphomolybdique 576
Acide phosphonique 577
Acide phosphorique 99% 547
Acide phosphorique 85% 547
Acide phosphorique 75% 549
Acide phosphorique 10% 549
Acide phosphorique-d3 85% dans D2O 549
Acide phosphosulfurique 578
Acide phosphotungstique 578
Acide phosphotungstique solution 578
Acide phtalique 579
Acide picrique solution 579
Acide 1-propanesulfonique sel sodique 627
Acide propionique 628
Acide pyrogallique 632
Acide pyrrolidinedithiocarboxylique-1,sel ammoniacal 633
Acide p-rosolique 640
Acide salicylique 642
Acide sébacique 646
Acide sélénieux 648
Acide silicotungstique 652
Acide stéarique 722
Acide succinique 725
Acide sulfamique 727
Acide sulfanilique 728
Acide sulfosalycilique 729
Acide sulfurique 98% 730
Acide sulfurique 96% (66°Bé) 731
Acide sulfurique 96% 731
Acide sulfurique 95 - 97% 733
Acide sulfurique 93-98% 733
Acide sulfurique 90% 734
Acide sulfurique 85% 734

INDEX ALPHABÉTIQUE

Acide sulfurique 72% 735
Acide sulfurique 69% 735
Acide sulfurique 62% 735
Acide sulfurique 50% 735
Acide sulfurique 35% (30°Be) 736
Acide sulfurique 30% 736
Acide sulfurique 25% 736
Acide sulfurique 20% 736
Acide sulfurique 10% v/v 737
Acide sulfurique 4 mol/l (8N) 737
Acide sulfurique 2.5 mol/l (5N) 737
Acide sulfurique 1 mol/l (2N) 737
Acide sulfurique 0,5 mol/l (1N) 738
Acide sulfurique 0.33 mol/l (2N/3) 738
Acide sulfurique 0.26 mol/l (0.52N) 739
Acide sulfurique 0.25 mol/l (0.5N) 739
Acide sulfurique 0.166 mol/l (0.333N) 739
Acide sulfurique 0.13 mol/l (0.26N) 740
Acide sulfurique 0.125 mol/l (0.25N) 740
Acide sulfurique 0.1 mol/l (0.2N) 740
Acide sulfurique 0.05 mol/l (0.1N) 740
Acide sulfurique 0.025 mol/l (0.05N) 741
Acide sulfurique 0.02 mol/l (0.04N) 741
Acide sulfurique 0.01 mol/l (0.02N) 741
Acide sulfurique 0.005 mol/l (0.01N) 742
Acide sulfurique 0.0025 mol/l (0.005N) 742
Acide sulfurique d=1,820 743
Acide sulfurique à 10 g/l d'argent sulfate 742
Acide sulfurique à 6.6 g/l d'argent sulfate 743
Acide sulfurique dilué 743
Acide tannique 744
Acide L(+)tartrique 745
Acide L(+)tartrique 20% 746
Acide 2-thiobarbiturique 755
Acide thioglycolique 80% 755
Acide p-toluenesulfonique monohydrate 768
Acide p-toluenesulfonique sel sodique 769
Acide trichloroacétique 770
Acide trichloroacétique 20% 770
Acide trifluoroacétique 773
Acide trifluoroacétique-d 774
Acide tungstique 780
ß-Alanine 148
Albumine d'oeuf poudre 148
Albumine d'oeuf sèche 148
Alcool n-amylique 180
Alcool tert-amylique 180
Alcool benzylique 200
Alcool cétylique 254
Alcool isoamylique 436
Aldéhyde anisique 181
Aldéhyde anisique solution 181
Aldéhyde benzoïque 197
Aldéhyde formique 37% m/v neutre 359
Aldéhyde formique 37% m/v 358
Aldéhyde formique 35% m/m 359
Aldéhyde formique 30% m/v 360
Aldéhyde formique 10% v/v selon Lillie 360
Aldéhyde formique 5% m/v tamponné à pH 6.9 360
Aldéhyde formique 4% m/v tamponné à pH 6.9 361
Aldéhyde formique 4% m/v salé 361
Aldéhyde glutarique solution 50% 372
Aldéhyde propionique 628
Alizarine 149
Alizarine solution saturée dans l'éthanol 149
Alliage de Devarda 294
Alliage de Raney 635
Allylthiourée 150
Alumine blanche 151
Aluminium, poudre 151
Aluminium standard solution 151
Aluminum ammonium sulfate dodécahydraté 152
Aluminium chlorure anhydre 153
Aluminium chlorure hexahydraté 153
Aluminium hydroxyde 154
Aluminium nitrate nonahydraté 154
Aluminium oxyde 155
Aluminium oxyde (acide) 155
Aluminium oxyde (basique) 155
Aluminium oxyde (neutre) 156
Aluminium oxyde actif 156
Aluminium potassium sulfate dodécahydraté 156
Aluminium sulfate 157
Amidon de maïs 476
Amidon de riz 640
Amidon soluble 721
Amidon soluble solution 721
4-Aminophénazone 159
m-Aminophenol 159
p-Aminophénol 159
Aminopyrazolone en solution 159
Ammoniaque solution 32% 160
Ammoniaque solution 30% 160
Ammoniaque solution 28% 161
Ammoniaque solution 25% 162
Ammoniaque solution 20 - 22% 163
Ammoniaque solution 17% 164
Ammoniaque solution 10% 164
Ammoniaque solution 6N 164
Ammoniaque solution diluée 164
Ammonium standard solution 165
Ammonium acétate 166
Ammonium bicarbonate 167
Ammonium bifluorure 172
Ammonium bromure 167
Ammonium carbamate 168
Ammonium carbonate 168
Ammonium carbonate solution 158 g/l 168
Ammonium chlorure 169
Ammonium citrate dibasique 170
Ammonium citrate solution 20% 171
tri-Ammonium citrate 170
Ammonium dihydrogénophosphate 25mg/l 172
Ammonium fluorure 171
Ammonium formiate 171
Ammonium iodure 173
Ammonium molybdate tétrahydraté 173
Ammonium molybdate solution 2.5% dans l'acide nitrique 174
Ammonium nitrate 174
Ammonium nitrate 200 mg/l 174
Ammonium oxalate monohydraté 175
Ammonium oxalate solution 4% 175
Ammonium persulfate 175
Ammonium phosphate dibasique 176
Ammonium phosphate monobasique 176
Ammonium sulfamate 177
Ammonium sulfate 177
Ammonium sulfure en solution à 20% 178
Ammonium L(+)tartrate 178
Ammonium thiocyanate 178
Ammonium thiocyanate 1 mol/l (1N) 179
Ammonium thiocyanate 0.1 mol/l (0.1N) 179
Ammonium thiocyanate 0.01 mol/l (0.01N) 179
Anhydride acétique 137
Anhydride maléique 477
Anhydride molybdique 511
Anhydride phosphorique 578
Anhydride phtalique 579
Anhydride succinique 725
Anhydride trifluoroacétique 774
Anhydride tungstique 780
Aniline chlorhydrate 181
Anthrone 182
Antimoine standard solution 182
Antimoine potassium tartrate 183
Antimoine trichlorure 183
Argent, laine 653
Argent, lames 652
Argent solution standard 653
Argent acétate 654
Argent carbonate 654
Argent chlorure 654
Argent diéthyldithiocarbamate 655
Argent nitrate 655
Argent nitrate solution 5% 656
Argent nitrate solution 2.90756% 656
Argent nitrate solution 658
Argent nitrate 1 mol/l (1N) 656
Argent nitrate 0.5 mol/l (0.5N) 656
Argent nitrate 0.1 mol/l (0.1N) 657
Argent nitrate 0.1 mol/l (0.1N) dans le propanol-2 657
Argent nitrate 0.05 mol/l (0.05N) 657
Argent nitrate 0.01 mol/l (0.01N) 658
Argent nitrate 0.01 mol/l (0.01N) dans le propanol-2 658
Argent oxyde 659
Argent sulfate 659
Argent sulfate solution 0.7% dans l'acide sulfurique 659
L(+)Arginine 184
L(+)Arginine monochlorhydratée 184
L(+)Asparagine 186
Arsenazo III 184
Arsenic standard solution 185
Arsenic trioxyde solution 185
Ausilab 110 187
Ausilab 140 187
Ausilab 210 187
Ausilab 250 188
Ausilab 260 188
Ausilab 280 188
Ausilab 290 188
Ausilab 400 189
Ausilab 500 189
Azométhine H 189
Azur II 189
Azur II Eosine 190

B

Baryum standard solution 191
Baryum acétate 192
Baryum carbonate 192
Baryum chlorure dihydraté 192
Baryum chlorure solution 10% 193
Baryum chlorure 61g/l 193
Baryum chlorure 30 g/l 194
Baryum chlorure 0.1 mol/l (0.2N) 193
Baryum chromate 194
Baryum hydroxyde octahydraté 194
Baryum hydroxyde solution 5% 195
Baryum hydroxyde solution 47.3 g/l 195
Baryum nitrate 195
Baryum perchlorate trihydraté 195
Baryum perchlorate 0.05 mol/l 196
Baryum perchlorate 0.025 mol/l 196
Baryum sulfate 196
Baume du Canada 247
Benzalkonium chlorure 197
Benzène 198
Benzène-d6 198
Benzène sulfonyle chlorure 198
p-Benzoquinone 199
Benzoyle péroxyde 296
Benzyle benzoate 200
Béryllium standard solution 201
Bismuth, granules 202
Bismuth standard solution 202
Bismuth (III) carbonate basique 203
Bismuth (III) nitrate basique 203
Bismuth (III) nitrate pentahydraté 204

Product specifications are subject to changes. Please visit our website for updates.

N,O-Bis(triméthylsilyl)acétamide 204
N,O-Bis(triméthylsilyl)trifluoroacétamide 204
Biuret 97% 205
Bleu alcalin 6B 150
Bleu alcalin 6B solution 2% dans l'éthanol 150
Bleu alcian 8GS 1% 148
Bleu alcian 8GX 149
Bleu d'aniline soluble dans l'eau 180
Bleu de bromophénol 213
Bleu de bromophénol solution 0.4% dans l'éthanol 213
Bleu de bromophénol solution 0.02% 214
Bleu de bromophénol solution 214
Bleu de bromothymol 215
Bleu de bromothymol solution 0.4% dans l'éthanol 215
Bleu de bromothymol solution 0.02% 215
Bleu de bromothymol solution 215
Bleu de Coomassie brillant R 250 274
Bleu de crésyl brillant 208
Bleu de lactophénol solution 452
Bleu de méthyle 497
Bleu de méthylène 498
Bleu de méthylène solution 1% 499
Bleu de méthylène solution saturée 498
Bleu de méthylène phéniqué solution hydroalcoolique 248
Bleu de méthylthymol sel de sodium 506
Bleu de thymol 757
Bleu de thymol solution 0.4% dans l'éthanol 757
Bleu de thymol solution 758
Bleu de toluidine 769
Bleu de toluidine phéniqué solution hydroalcoolique 248
Bleu Victoria 785
Bore standard solution 208
Bromate standard solution 210
Brome solution 210
alpha-Bromonaphtalène 213
Bromure standard solution 210
Bromure - bromate 0.0167 mol/l 210
Brun Bismarck R 201
Butanol-1 226
Butanol-2 227
tert-Butanol 227
2-Butoxyéthanol 228
2-(2-Butoxyéthoxy)éthanol 228
n-Butyle acétate 228
n-Butyle chlorure 229
Butylhydroxytoluène 230
tert-Butylméthyléther 230

C

Cadmium standard solution 232
Cadmium acetate dihydrate 233
Cadmium carbonate 233
Cadmium chlorure monohydraté 233
Cadmium nitrate tétrahydraté 234
Caféine anhydre 234
Calcium standard solution 234
Calcium acétate anhydre 236
Calcium acétate monohydrate 236
Calcium carbonate 236
Calcium chlorure anhydre 237
Calcium chlorure dihydraté 238
Calcium chlorure hexahydraté 239
Calcium chlorure solution 0.025% 239
Calcium fluorure 240
Calcium gluconate 240
Calcium hydroxyde 240
Calcium lactate 241
Calcium nitrate tétrahydraté 241
Calcium oxyde, morceaux 242
Calcium oxyde, poudre 242
Calcium pantothénate 242
Calcium phosphate dibasique dihydraté 243
Calcium phosphate monobasique monohydraté 243
Calcium phosphate tribasique 243
Calcium stéarate 244
Calcium sulfate dihydraté 244
Calcium sulfate hemihydraté 245
Calcium sulfate hemihydraté solution 245
Calcon 245
Calmagite 246
Camphre naturel 246
Camphre synthétique 246
Carborundum, poudre 248
Carborundum, granulés 248
Carmin d'indigo sec 422
Carmin d'indigo solution 422
Caséine 249
Catalyseur Kjedahl pour analyse de l'eau 450
Catalyseur Kjeldahl au sélénium 450
Catalyseur Kjeldahl au titane 451
Catalyseur Kjeldahl exempt de sélénium et titane 450
Catalyseur Kjeldahl selon Wieninger 450
Cellulose, poudre 250
Cérium standard solution 250
Cérium (IV) ammonium nitrate 251
Cérium (IV) ammonium nitrate 0.1 mol/l 251
Cérium (IV) ammonium nitrate 0.01 mol/l 251
Cérium (IV) ammonium sulfate dihydraté 252
Cérium (IV) ammonium sulfate 0.1 mol/l 252
Cérium (IV) ammonium sulfate 0.01 mol/l 252
Cérium (IV) sulfate tétrahydraté 252
Cérium (IV) sulfate 0.1 mol/l 253
Cérium (III) nitrate hexahydraté 251
Césium standard solution 253
Césium chlorure 254
Césium chlorure 25 g/l 254
Césium sulfate 254
Charbon actif 255
Charbon décolorant 255
Charbon végétal 256
Chaux sodée 660
Chloral hydraté 256
Chloramine T sel sodique 256
Chloranile 256
Chlorate standard solution 257
Chlorite standard solution 257
Chloroacétamide 258
p-Chlorobenzaldéhyde 258
Chlorobenzène 258
Chlorobutanol hémihydrate 259
Chloroforme 259
Chloroforme-d 263
Chloroforme-d + 0.03% TMS 263
p-Chlorophénol 264
Chlorure standard solution 257
Chlorure d'étain (II) dihydrate 761
Chlorure d'étain (II) dihydrate solution 10% 761
Chlorure de zinc - Solution iodée 800
Chlorure de zinc-acide formique - Solution 800
Chromate standard solution 264
Chrome standard solution 264
Chrome (III) chlorure hexahydrate 265
Chrome (III) nitrate nonahydrate 266
Chrome (III) oxyde 266
Chrome (III) potassium sulfate dodécahydraté 266
Chrome (III) sulfate 267
Chrome (VI) oxyde 267
Chrome (VI) oxyde solution 267
Chrysoïdine Y 268
Cinchonine 268
Citral 268
Cobalt standard solution 270
Cobalt (II) acétate tétrahydraté 271
Cobalt (II) ammonium sulfate hexahydraté 271
Cobalt (II) chlorure hexahydrate 272
Cobalt (II) chlorure en solution 272
Cobalt (II) nitrate hexahydrate 272
Cobalt (II) sulfate heptahydrate 273
Colorant de Wright solution methanolique 790
Coton à l'acétate de plomb (II) 457
m-Crésol 282
o-Crésol 282
Cuivre électrolytique, fils 275
Cuivre électrolytique, tournures 274
Cuivre électrolytique recuit, feuilles 274
Cuivre réduit, poudre 275
Cuivre standard solution 275
Cuivre (I) chlorure 276
Cuivre (I) chlorure solution 7% dans l'ammoniaque 277
Cuivre (I) iodure 277
Cuivre (II) acétate hydrate 277
Cuivre (II) carbonate basique 278
Cuivre (II) chlorure dihydraté 278
Cuivre (II) éthylènediamine hydroxyde 1 mol/l 285
Cuivre (II) nitrate trihydraté 279
Cuivre (II) oxyde 279
Cuivre (II) oxyde en fils 279
Cuivre (II) sulfate anhydre 280
Cuivre (II) sulfate pentahydraté 280
Cuivre (II) sulfate solution 12.5% 281
Cuivre (II) sulfate en solution 281
Cuprié thylènediamine solution 285
Cyanure solution standard 285
Cyclohexane 286
Cyclohexane-d12 287
Cyclohexanol 287
Cyclohexanone 288
Cyclohexylamine 288
Cyclopentylméthyléther 289
L-Cystéine 289
L-Cystine 289

D

Décalcifiant 290
Décane 290
Deutérium oxyde-d2 293
Deutérium oxyde-d2 + 0.01% DMSO 294
Deutérium oxyde-d2 + 0.5% TSP d4 294
Deutérium oxyde-d2 + 0.03% TSP d4 294
Diacétone alcool 295
Diacétyldioxime 295
Diacétyldioxime solution 1% dans l'éthanol 296
Diacétyldioxime sel sodique 295
Dibenzoylméthane 296
n-Dibutyle phtalate 297
Dicalite 4158 297
m-Dichlorobenzène 298
o-Dichlorobenzène 298
p-Dichlorobenzène 298
1,2-Dichlorobenzène-d4 299
1,2-Dichloroéthane 299
Dichlorométhane 300
Dichlorométhane acidifié avec 1% d'acide chlorhydrique 305
Dichlorométhane-d2 305
2,6-Dichlorophénolindophénol sel de sodium 305
2,6-Dichloroquinone-4-chloroimide 305
Diéthanolamine 306
Diéthylamine 306
Diéthylcarbonate 307
Diéthyle phtalate 310
n,n'-Diéthyl-p-phénylènediamine sulfate 310
Différenciateur pour kit de Gram-Hucker 311
Digitonine 311
Diisopropylamine 311
N,N-Diisopropyléthylamine 312

INDEX ALPHABÉTIQUE

Dimédon 313
N,N-Diméthylacétamide 313
p-Diméthylaminobenzaldéhyde 314
p-Diméthylaminobenzalrhodanine 314
1,2-Diméthoxyéthane 313
N,N-Diméthylformamide 314
N,N-Diméthylformamide-d7 316
N,N-Diméthylformamide diméthylacétal 317
N,N'-Diméthylpropylène urée 317
Diméthylsulfoxyde 317
Diméthylsulfoxyde-d6 319
Diméthylsulfoxyde-d6 + 0.03% TMS 319
Dimidium bromure 320
2,4-Dinitrochlorobenzène 320
2,4-Dinitrophénylhydrazine (avec 30% d'eau) 320
1,4-Dioxane 321
1,3-Dioxolane 322
Diphénylamine solution 1% dans l'acide sulfurique 322
4-Diphénylamine sulfonate sodique 322
Sym-Diphénylcarbazide 323
Sym-Diphénylcarbazone 323
Diphénylthiocarbazone 323
2,2'-Dipyridyle 323
Dithiooxamide 324
Dodécyltriméthylammonium bromure 325
Dysprosium solution standard 325

E

Eau 786
Eau purifiee 788
Eau + 0.1%(v/v) acide formique 789
Eau + 0.1% v/v d'acide trifluoroacétique 789
Eau de brome 211
Eau de chaux 184
Eau de chlore 789
Eau désionisée acidifiée 789
Eau oxygénée solution 40% m/v 413
Eau oxygénée solution 35% 413
Eau oxygénée solution 30-32% 413
Eau oxygénée solution 30% 414
Eau oxygénée en solution à 6% 415
Eau oxygénée solution 3.5% m/v 415
Eau oxygénée solution 3% 416
Eluant sodium carbonate 327
Eluant sodium carbonate/sodium dicarbonate 328
Eluant sodium dicarbonate 327
Empois d'amidon solution 1% 721
Eosine B 328
Eosine Y 328
Eosine Y 1% solution aqueuse 328
Eosine Y 0.5% solution alcoolique 329
Erbium solution standard 329
Eriochrome cyanine R 330
Erythrosine extra B 330
Essence E 563
Etain en feuilles 759
Etain, poudre 759
Etain solution standard 760
Etain (II) chlorure solution 761
Etain (II) sulfate 762
Etain (IV) chlorure pentahydrate 762
Etain (IV) oxyde 762
Etalons couleurs ASTM 187
Etalons couleurs Gardner 368
Etalons couleurs Saybolt 644
Etalon de conductivité 500 000 µS/cm 721
Etalon de conductivité 450 000 µS/cm 720
Etalon de conductivité 350 000 µS/cm 720
Etalon de conductivité 300 000 µS/cm 720
Etalon de conductivité 200 000 µS/cm 720
Etalon de conductivité 150 000 µS/cm 719
Etalon de conductivité 100 000 µS/cm 719
Etalon de conductivité 50 000 µS/cm 719
Etalon de conductivité 20 000 µS/cm 719
Etalon de conductivité 12 880 µS/cm 718
Etalon de conductivité 10 000 µS/cm 718
Etalon de conductivité 5000 µS/cm 718
Etalon de conductivité 1413 µS/cm 718
Etalon de conductivité 1000 µS/cm 718
Etalon de conductivité 500 µS/cm 717
Etalon de conductivité 200 µS/cm 717
Etalon de conductivité 147 µS/cm 717
Etalon de conductivité 100 µS/cm 717
Etalon de conductivité 84 µS/cm 717
Etalon de conductivité 50 µS/cm 716
Etalon de conductivité 20 µS/cm 716
Etalon de conductivité 10µS/cm 716
Etalon de conductivité 5 µS/cm 716
Etalon de conductivité 1.30 µS/cm 715
Etalons de densité 292
Etalons d'indice de réfraction 637
Etalons d'osmolalité 550
Etalons d'osmolalité à base de protéines 550
Etalons d'osmolalité à base d'urine 550
Etalons en degré Brix 209
Etalon multiéléments pour chromatographie ionique 512
Etalon multiéléments pour ICP 513
Etalon multiéléments pour ICP et ICP-MS 515
Ethanol absolu anhydre 331
Ethanol 96° 334
Ethanol 70% v/v 336
Ethanol 60% v/v 337
Ethanol 50% v/v 337
Ethanol absolu dénaturé 338
Ethanol 95° dénaturé 338
Ethanol 94° dénaturé 339
Ethanol 90° dénaturé 339
Ethanol 70° modifié 339
Ethanol anhydre-d6 340
Ethanol 95%-Ether ethylique 50/50(m/m) - 30 mg/l phénolphtaléine 507
Ethanolamine 340
Ether de pétrole 100 - 140°C 564
Ether de pétrole 80 - 120°C 564
Ether de pétrole 80 - 100°C 564
Ether de pétrole 75 - 120°C 565
Ether de pétrole 60 -80°C 565
Ether de pétrole 55 - 85°C 565
Ether de pétrole 40 - 70°C 566
Ether de pétrole 40 - 65°C 566
Ether de pétrole 40 - 60°C 567
Ether de pétrole 35 - 60°C 567
Ether de pétrole 30 - 50°C 568
Ether de pétrole 30 - 40°C 569
Ether éthylique 308
Ether isopropylique 312
Ethyle acétate 340
Ethyle acétoacétate 343
Ethyle formiate 350
Ethylènediamine 344
Ethylène glycol 348
Eugenol 351
EUKITT 351
Europium solution standard 351

F

Fer, poudre 427
Fer, réduit par l'hydrogène 428
Fer solution standard 428
Fer(II) ammonium sulfate hexahydraté 429
Fer(II) ammonium sulfate 0.1N 430
Fer (II) ammonium sulfate 0.12N 429
Fer (II) chlorure tétrahydraté 430
Fer (II) sulfate heptahydraté 430
Fer (II) sulfate 0.1 mol/l 431
Fer (III) ammonium citrate rouge 431
Fer (III) ammonium citrate vert 431
Fer (III) ammonium oxalate 432
Fer (III) ammonium sulfate dodécahydraté 432
Fer (III) ammonium sulfate 33% dans l'acide nitrique 432
Fer (III) ammonium sulfate 0.1 mol/l 433
Fer (III) ammonium sulfate 100 g/l 433
Fer (III) chlorure anhydre sublimé 433
Fer (III) chlorure hexahydraté 433
Fer (III) chlorure solution 4.5% 434
Fer chlorure en solution 434
Fer (III) citrate 434
Fer(III) nitrate nonahydraté 435
fer (III) oxyde 435
Fer (III) sulfate 435
Ferricyanures standard solution 353
Ferrocyanures standard solution 354
Ferron 354
Florisil 100-200 mesh 357
Florisil 60-100 mesh 356
Fluorescéine 357
Fluorescéine sodique 357
Fluorure solution standard 358
Formamide 362
Formol acétique 362
D(-)Fructose 364
Fuchsine acide 364
Fuchsine basique 365
Fuschine décolorée solution 365

G

Gadolinium solution standard 366
D(+)Galactose 366
Gallium standard solution 367
Gélatine 368
Gel de silice 60A 0,06÷0,20 mm 650
Gel de silice 60A 6 - 35µ 648
Gel de silice 60A 20 - 45µ 648
Gel de silice 60A 35 - 70µ 649
Gel de silice 60A 40 - 63µ 649
Gel de silice 60A 70 - 200µ 650
Gel de silice granulés 650
Gel de silice granulés avec indicateur exempt de cobalt 651
Germanium solution standard 369
D(+)Glucose anhydre 370
D(+)Glucose monohydraté 371
Glycérine (30°Bé) 372
Glycérine 90% (28°Bé) 373
Glycocolle 373
Glycol diéthylénique 307
Glycol diéthylénique diméthyléther 307
Glycol éthylénique monométhylether 495
Glycol triéthylénique 773
Glyoxal standard solution 374
Glyoxal-bis-(2-hydroxyanile) 374
Gomme arabique 377

H

Hafnium solution standard 378
Hemalun en solution selon Carazzi 378
Hemalun en solution selon Mayer 378
Hématoxyline 379
Hématoxyline en solution selon Mayer 379
n-Heptane 99% 380
n-Heptane 382
Heptane mélange d´isomères 382
Hexadécane 383
Hexafluoro-2-Propanol 384
Hexaméthyldisilazane 384
Hexaméthylènetétramine 384
n-Hexane 99% 385
n-Hexane 386
Hexane mélange d´isomères 388

Product specifications are subject to changes. Please visit our website for updates.

Histamine dichlorhydratée 389
L-Histidine monochlorhydratée monohydraté 390
Histolemon 390
Holmium solution standard 390
Holmium oxyde 391
Holmium perchlorate solution 391
Huile d'amande raffinée 539
Huile de cèdre condensée 538
Huile de ricin 249
Huile de vaseline 555
Huile pour immersion 420
Hyamine 1622 391
Hyamine 1622 solution 0.004M 392
Hydrazine dichlorhydratée 392
Hydrindantine 392
Hydroquinone 416
p-Hydroxybenzaldéhyde 416
Hydroxylamine sulfate 417
Hydroxylamine solution alcoolique 417
8-Hydroxyquinoléine 417

I

Idrimeter Erba Indicateur C 419
Idrimeter Erba Solution A 419
Idrimeter Erba Solution B 419
Imidazole 419
Indicateur bleu de bromophénol 214
Indicateur bleu de thymol 758
Indicateur Ferroïne 0.025 mol/l 354
Indicateur mixte au méthylorange 502
Indicateur mixte au rouge de méthyle 504
Indicateur pour iodométrie 420
Indicateur pour l'azote ammoniacal 420
Indicateur TA au bleu de thymol 758
Indicateur TAC au bleu de bromophénol 214
Indicateur universel pH 1-11 solution hydroalcoolique 421
Indicateur universel pH 1-11 solution aqueuse 422
Indicateur universel pH 0-5 solution hydroalcoolique 421
Indium solution standard 422
Indole 423
Inositol 423
Inuline 424
Iode 424
Iode bisublimé 425
Iode 0.5 mol/l (1N) 425
Iode 0.05 mol/l (0.1N) 425
Iode 0.01 mol/l (0.02N) 426
Iode 0.005 mol/l (0.01N) 426
Iode 10 ppm 426
Iode bromure solution 426
Iodoforme 427
Iodure solution standard 424
Iridium solution standard 427
Isoamyle acétate 436
Isobutanol 437
Isobutyle acétate 438
Isohexane 438
Isooctane 439
Isopar G 440
Isopentane 441
Isopropylamine 442
Isopropyle acétate 441

J

Jaune Clayton 270
Jaune d'alizarine R 150
Jaune de méthanyle 487
Jaune de méthyle 506
Jaune naphtol S 518
Jaune Soudan 727

K

Kaolin lavé 443
Karl Fischer anolyte solution 446
Karl Fischer anolyte solution pour cétones et aldéhydes - Sans méthanol 446
Karl Fischer anolyte solution, pour four 446
Karl Fischer anolyte solution pour les huiles 447
Karl Fischer anolyte solution - sans CFC 446
Karl Fischer catholyte solution 447
Karl Fischer catholyte solution pour aldehydes et cétones 447
Karl Fischer réactif pour la cellule avec et sans diaphragme 447
Karl Fischer solvant 2 composant avec tampon 445
Karl Fischer solvant pour huiles 1 composant 444
Karl Fischer solvant 1 composant pour aldéhydes et cétones 444
Karl Fischer solvant 2 composants pour aldéhydes et cétones - Sans méthanol 445
Karl Fischer solvent 2 component 445
Karl Fischer solvant pour huiles 2 composant 445
Karl Fischer titrant 1 composant 5 mg H2O/ml pour aldehydes et cétones 443
Karl Fischer titrant 2 composants 5 mg H2O/ml 444
Karl Fischer titrant 2 composants 2 mg H2O/ml 444
Karl Fischer standard eau 0.10 mg/g 449
Karl Fischer standard eau 1.0 mg/g 448
Karl Fischer standard eau 5.0 mg/ml 448
Karl Fischer standard eau 10.0 mg/g 448
Kérosène 449
Kit Gram-Hucker 376
Kjeldahl Antimousse 449

L

Lactophénol d'Amman en solution 160
Lactose monohydrate 452
Laine de verre 370
Lanoline anhydre 453
Lanthane solution standard 453
Lanthane chlorure 25g/l 454
Lanthane nitrate hexahydraté 454
Lanthane nitrate 50 g/l 454
Lanthane nitrate 0.1 mol/l 454
Lanthane oxyde 455
Lécithine de soja 461
L(+)-Leucine 461
Levure de bière sèche 794
Liquide de Bouin 355
Liquide de Bouin Hollande 355
Liquide fixateur AFA 354
Liquide fixateur de Davidson 355
Liquide fixateur FIXALL-HIS 356
Liquide fixateur sans acide acétique 356
Lithium solution standard 462
Lithium carbonate 463
Lithium chlorure 463
Lithium hydroxyde anhydre 464
Lithium hydroxyde monohydraté 464
Lithium hydrure 464
Lithium méthylate 0.1 mol/l (0.1N) 464
Lithium nitrate 465
Lithium sulfate monohydraté 465
Lithium tétraborate anhydre 465
Luminol 467
Lutéthium solution standard 467
L(+)Lysine monochlorhydrate 467

M

Magnésium, poudre 468
Magnesium, rubans 468
Magnésium, tournures 468
Magnésium solution standard 468
Magnésium acétate tétrahydraté 470
Magnésium carbonate basique 470
Magnésium chlorure hexahydraté 471
Magnesium glycerophosphate 471
Magnésium hydroxyde 472
Magnésium nitrate hexahydraté 472
Magnésium nitrate 10 g/l 473
Magnésium oxyde 473
Magnésium oxyde lourd 473
Magnsium perchlorate 474
Magnésium péroxyde 474
Magnesium phosphate dibasique trihydraté 472
Magnésium stéarate 474
Magnésium sulfate anhydre 475
Magnésium sulfate heptahydraté 475
Maltose monohydraté 478
Manganèse électrolytique 479
Manganèse solution standard 479
Manganèse (II) acétate tetrahydrate 480
Manganèse (II) chlorure tetrahydrate 480
Manganèse (II) sulfate monohydraté 481
Manganèse (IV) dioxyde 481
D-Mannitol 481
Mélange alcool/ether 30/70 m/m 507
Mélange CHM 507
Mélange éthanol 99° / 2-propanol 508
Mélange éthanol 95° / 2-propanol 508
Mélange Ethanol absolu/Ether Ethylique 50/50 (m/m) 507
Mélange Ethanol absolu/Ether Ethylique 50/50 (v/v) 508
Mélange octane 80 537
Mélange pour indice de brome 508
Mélanges pour la recherche des solvants résiduels 509
Mélange pour le contrôle des soudures des emballages 509
Mélange sélénique 646
Mélange standard pour analyse des hydrocarbures 546
L-Menthol 483
2-Mercaptoéthanol 483
Mercure solution standard 484
Mercure (I) chlorure 485
Mercure (II) chlorure 485
Mercure (II) chlorure solution 485
Mercure (II) iodure 486
Mercure (II) oxyde rouge 486
Mercure (II) sulfate 486
Mercure (II) sulfate solution 486
Méthanol 488
Méthanol + 0.1% v/v acide formique 492
Méthanol + 0.1% v/v acide trifluoroacétique 493
Méthanol-d1 494
Méthanol-d3 494
Méthanol-d4 493
Méthanol-d4 + 0.03% TMS 494
Méthanol - chlorhydrique 493
DL-Méthionine 494
4-Méthylaminophénol sulfate 496
3-Méthyl-2-benzothiazoline hydrazone chlorhydrate 497
Méthylcyclohexane 497
Méthyle acétate 495
Méthyle benzoate 496
Méthyle iodure 500
Méthyle p-oxybenzoate 499
N-Méthyle-2-pyrrolidone 502
Méthyle salicylate 505
Méthyléthylcétone 350
Méthylisoamylcétone 500
Méthylisobutylcétone 500
Méthylorange 501
Méthylorange solution 0.1% 501
2-Méthyltétrahydrofurane 505

INDEX ALPHABÉTIQUE

4-Méthyltétrahydropyrane 506
Molybdène solution standard 510
Morpholine 512
Mucicarmin en solution hydroalcoolique 512
Murexide 516
Myristate d'isopropyle 442

N

Naphtalène 517
1-Naphtol 517
2-Naphtol 517
alpha-Naphtolbenzéine 518
alpha-Naphtolbenzéine solution 0.2% dans acide acétique 518
N-(1-Naphtyl)éthylènediamine dichlorhydraté 518
Néocuproïne chlorhydraté 519
Néodyme solution standard 519
Nickel, poudre 521
Nickel solution standard 521
Nickel (II) acétate tétrahydraté 522
Nickel (II) ammonium sulfate hexahydraté 523
Nickel (II) carbonate basique 523
Nickel (II) chlorure hexahydraté 523
Nickel (II) nitrate hexahydraté 524
Nickel (II) nitrate 10 g/l 524
Nickel (II) sulfate hexahydraté 524
Nicotinamide 525
Nigrosine 525
Ninhydrine 526
Ninhydrine solution 526
Niobium solution standard 526
Nital solution 4% 527
Nitrate solution standard 527
Nitrite solution standard 534
m-Nitrobenzaldéhyde 534
Nitrobenzène 535
p-Nitrophénol 535
p-Nitrophénol solution 0.1% dans l'eau 535
Noir amido 10B solution 157
Noir ériochrome T 329
Noir Soudan B 726
Nonylphénol éthoxylate 10 ETO 536

O

n-Octane 537
Octanol-1 538
Or solution standard 375
Or (III) chlorure acide trihydraté 375
Orange II 539
Orange G 539
Orange acridine 147
Orange de xylénol 793
Orcéine 540
Orcinol 540
Osmium solution standard 549
Oxyde d'éthylène solution 349

P

Palladium solution standard 553
Palladium (II) chlorure 553
Palladium nitrate 2 g/l 554
Papanicolaou Hématoxyline selon Harris 554
Papanicolaou solution EA 50 554
Papanicolaou solution OG 6 555
Papier (II) à l'acétate de plomb 457
Papier argent manganèse 655
Papier bromure mercurique 484
Papier indicateur de pH 420
Papier potassium iodure d' amidon 608
Papier rouge Congo 273
Papier tournesol 466
Paraffine 56-58°C Erbaplast (sans DMSO) 555
Paraformaldéhyde 556
Pararosaniline décolorée en solution 556
n-Pentane 99% 557
n-Pentane 558
Pepsine 560
Pétrole 563
o-Phénanthroline monohydraté 569
o-Phénanthroline-Fer (II) sulfate 569
Phénol 570
Phénol liquide 85% 570
Phénolphtaléine 571
Phénolphthaléine solution à 1% dans l'éthanol 571
Phénolphtaléine solution 0.1% 572
Phénolphthaléine solution 0.1% dans l'éthanol 572
L-Phénylalanine 572
2-Phényléthanol 573
Phénylhydrazine chlorhydrate 573
Phénylhydrazine chlorhydrate solution 573
Phloroglucinol 574
Phloxine B 574
Phosphate solution standard 574
Phosphore solution standard 577
o-Phtalaldéhyde 579
Pierre ponce 630
Pipéridine 580
Platine solution standard 580
Plomb solution standard 455
Plomb 455
Plomb (II) acétate trihydraté 458
Plomb (II) acétate solution 95g/l 458
Plomb (II) acétate basique 457
Plomb (II) acétate basique soluiton 457
Plomb (II) nitrate 458
Plomb (II) nitrate 0.1 mol/l 459
Plomb (II) nitrate 0.05 mol/l 459
Plomb (II) nitrate solution 33g/l 459
Plomb (II) oxyde 460
Plomb (II) sulfate 460
Plomb (IV) oxyde 460
Polyvinylpyrrolidone 581
Potassium solution standard 582
Potassium acétate 583
Potassium bicarbonate 584
Potassium bisulfate 584
Potassium bromate 584
Potassium bromate 0.033 mol/l (0.198N) 585
Potassium bromate 0.02 mol/l (0.12N) 585
Potassium bromate 0.0167 mol/l (0.1N) 585
Potassium bromure 585
Potassium carbonate 586
Potassium chloroplatinate 597
Potassium chlorure 586
Potassium chlorure solution saturée 590
Potassium chlorure solution 590
Potassium chlorure 3.5 mol/l (3.5N) avec argent chlorure 587
Potassium chlorure 3.5 mol/l (3.5N) 587
Potassium chlorure 3 mol/l (3N) solution eau/glycérine 588
Potassium chlorure 3 mol/l (3N) avec argent chlorure 588
Potassium chlorure 3 mol/l (3N) 588
Potassium chlorure 1 mol/l (1N) 588
Potassium chlorure 0.1 mol/l (0.1N) 589
Potassium chlorure 0.01 mol/l (0.01N) 589
Potassium chlorure 25g/l dans HCl 589
Potassium chlorure 12g/l 589
Potassium chromate 590
Potassium chromate 10% solution 590
Potassium chromate solution 5% 591
Potassium citrate tribasique monohydraté 591
Potassium dichromate 592
Potassium dichromate 0.0167 mol/l (0.1 N) 593
Potassium dichromate 0.0417 mol/l (0.25N) 593
Potassium dichromate 0.04 mol/l (0.24N) dans 80 g/l HgSO4 593
Potassium dichromate 0.167 mol/l (1 N) 592
Potassium dichromate solution 106g/l 594
Potassium dichromate 0.1414 g/l 594
Potassium dichromate solution 0.5% 592
Potassium dichromate - Solution dans l'acide sulfurique 592
Potassium ferricyanure 594
Potassium ferrocyanure trihydraté 595
Potassium ferrocyanure solution 10% 595
Potassium ferrocyanure solution 53 g/l 596
Potassium fluorure 596
Potassium fluorure dihydraté 596
Potassium hydroxyde, pastilles 599
Potassium hydroxyde, écailles 599
Potassium hydroxyde solution 45% 600
Potassium hydroxyde 38% (40°Bé) 601
Potassium hydroxyde 33% 601
Potassium hydroxyde 28% 601
Potassium hydroxyde 2 mol/l (2N) dans l'ethanol 602
Potassium hydroxyde 1 mol/l (1N) 602
Potassium hydroxyde 0.5 mol/l (0.5N) 602
Potassium hydroxyde 0.5 mol/l (0.5N) dans l'éthanol 603
Potassium hydroxyde 0.5 mol/l (0.5N) dans le méthanol 603
Potassium hydroxyde 0.46 mol/l (0.46N) 604
Potassium hydroxyde 0.25 mol/l (0.25N) 604
Potassium hydroxyde 0.23 mol/l (0.23N) 604
Potassium hydroxyde 0.1 mol/l (0.1N) 604
Potassium hydroxyde 0.1 mol/l (0.1N) dans l'ethanol 605
Potassium hydroxyde 0.1 mol/l (0.1N) dans le méthanol 605
Potassium hydroxyde en solution 606
Potassium hydroxyde solution 3% dans l'ethanol 601
Potassium iodate 606
Potassium iodate 0.05 mol/l (0.3N) 606
Potassium iodate 0.0167 mol/l (0.1N) 606
Potassium iodate 0.00167 mol/l (0.01N) 607
Potassium iodate acide 597
Potassium iodobismuthate solution 609
Potassium iodure 607
Potassium iodure solution 10% 608
Potassium iodure solution 3.9% 608
Potassium iodure solution 608
Potassium métabisulfite 609
Potassium nitrate 609
Potassium nitrate 1 mol/l (1N) 610
Potassium nitrite 610
Potassium oxalate monohydraté 610
Potassium m-périodate 611
Potassium permanganate 611
Potassium permanganate 0.2 mol/l (1N) 612
Potassium permanganate 0.02 mol/l (0.1N) 612
Potassium permanganate 0.002 mol/l (0.01N) 613
Potassium permanganate solution 3% 612
Potassium permanganate - Solution phosphorique 613
Potassium persulfate 613
Potassium phosphate dibasique anhydre 614
Potassium phosphate dibasique trihydraté 614
Potassium phosphate monobasique 615
Potassium phosphate monobasique 0.2 mol/l (0.2N) 615
Potassium phtalate acide 597
Potassium phtalate acide 0.2 mol/l (0.2N) 598
Potassium phtalate acide 0.1 mol/l (0.1N) 598
Potassium pyroantimoniate acide 615
Potassium pyroantimoine solution 616
Potassium pyrophosphate 616
Potassium pyrosulfate 616
Potassium sodium tartrate tétrahydraté 617

DASIT GROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

Potassium sorbate 617
Potassium sulfate 618
Potassium sulfocyanure 620
Potassium sulfocyanure 0.1 mol/l (0.1N) 620
Potassium sulfocyanure solution 5% 620
Potassium sulfocyanure solution 621
Potassium sulfogaïacolate 596
Potassium tartrate 618
Potassium L-tartrate monobasique 619
Potassium tellurite 619
Potassium tétraiodomercurate - Solution alcaline 619
Pourpre de bromocrésol 212
Pourpre de bromocrésol solution 0.4% dans l'éthanol 212
Pourpre de bromocrésol solution 213
Pourpre de m-cresol 282
Praséodyme solution standard 621
L(-)Proline 622
Propane-1-ol 622
Propane-2-ol 623
Propane-2-ol 70% 626
1,3-Propanediol 627
n-Propyle acétate 628
Propyle p-oxybenzoate 629
Propylène carbonate 629
Propylène glycol 629
Pyridine 630
Pyridine-d5 631
1-(2-Pyridyl-azo)-2-napthol 632

Q

Pyrocatéchol 632
Quinidine sulfate 634

R

Quinoléine 634
Raffinose 635
Réactif à la ninhydrine et au étain (II) chlorure 526
Réactif à l'acide aminohippurique 158
Réactif à l'iodoplatinate 427
Réactif au biuret 205
Réactif de Bénédict 197
Réactif de Carrez sel de zinc 249
Réactif de Carrez sel de potassium 249
Réactif d'Ehrlich 327
Réactif d'Esbach 330
Réactif de Fehling solution A 353
Réactif de Fehling solution B 353
Réactif de Folin-Ciocalteu 358
Réactif de Giemsa 369
Réactif de Griess 376
Réactif A de Griess 376
Réactif B de Griess 376
Réactif de Hanus 379
Réactif de Karl Fischer 1 composant 2 mgH2O/ml 443
Réactif de Karl Fischer 1 composant 5 mg H2O/ml 443
Réactif de Lugol solution 466
Réactif de Lugol solution concentrée 466
Réactif de Lugol solution pour kit de Gram-Hucker 466
Réactif de May Grünwald 482
Réactif de Mayer 482
Réactif molybdovanadique 511
Réactif de Nessler solution 520
Réactif de Nessler solution A 520
Réactif de Nessler solution B 521
Réactif phosphomolybdotungstique 576
Réactif pour l'indole selon Kovac 451
Réactif de Schiff pour Aldéhydes 645
Réactif de Schiff pour PAS coloration 646
Réactif sulfomolybdique 729
Réactif TAN 635
Réactif TBN ASTM D2896 636
Réactif TBN ASTM D4739 636
Réactif de Tauber 746
Réactif de Turk 780
Réactif unique pour lipolyse 635
Réactif de Wijs 790
Réactif de Ziehl-Neelsen en solution 796
Résorcine 637
L(+)Rhamnose 638
Rhénium solution standard 638
Rhodium solution standard 638
Riboflavine 639
D(-)Ribose 639
Rouge Alizarine 149
Rouge Congo 273
Rouge Congo solution 274
Rouge de crésol solution dans l'hydroxyde de sodium 0.1N / éthanol 96% / eau 283
Rouge de o-crésol 283
Rouge de o-crésol solution 0.2% dans l'éthanol 283
Rouge de phénol 570
Rouge de phénol solution 0.2% dans l'éthanol 571
Rouge de phénol solution 571
Rouge de quinaldine 634
Rouge méthyle 503
Rouge méthyle solution 0.1% dans l'éthanol 504
Rouge méthyle solution 0.2% dans l'éthanol 503
Rouge de méthyle solution 504
Rouge neutre 521
Rouge nucléaire 536
Rouge Ponceau BS 581
Rouge Ponceau S 581
Rouge pour l'huile O 636
Rubidium solution standard 640

S

Ruthénium solution standard 641
Sable de Fontainebleau 644
Sable purifié 644
D(+)-Saccharose 726
Safranine T 642
Safranine T solution hydroalcoolique pour kit de Gram-Hucker 642
Salicylaldéhyde-azine 642
Samarium solution standard 643
Saponine 644
Savon en solution dans l'éthanol 660
Scandium solution standard 645
Sel de Reinecke 637
Sélénium, poudre 646
Sélenium dioxyde 648
Sélénium solution standard 647
Silicium solution standard 651
Silicium dioxyde 652
Sodium solution standard 661
Sodium acétate anhydre 662
Sodium acétate trihydraté 662
Sodium alginate 663
Sodium aluminate 663
Sodium ammonium hydrogénophosphate 664
Sodium L-ascorbate 664
Sodium arsenite 0.1 mol/l (0.2N) 664
Sodium arsenite 0.05 mol/l (0.1N) 664
Sodium azide 665
Sodium benzoate 665
Sodium bicarbonate 665
Sodium bisulfate monohydraté 666
Sodium borohydrure 666
Sodium bromure 667
Sodium carbonate anhydre 667
Sodium carbonate anhydre solution 670
Sodium carbonate décahydraté 668
Sodium carbonate monohydraté 668
Sodium carbonate 0.5 mol/l (1N) 669
Sodium carbonate 0.05 mol/l (0.1N) 669
Sodium carbonate solution 20% 670
Sodium chlorure 670
Sodium chlorure 5 mol/l (5N) 672
Sodium chlorure 0.1 mol/l (0,1N) 672
Sodium citrate dibasique 672
Sodium citrate tribasique anhydre 672
Sodium citrate tribasique bihydraté 673
Sodium cobalt nitrite 673
Sodium cyanoborohydrure 674
Sodium dichloroisocyanurate dihydraté 674
Sodium diéthyldithiocarbamate trihydraté 674
Sodium fluorure 675
Sodium formiate 675
Sodium gluconate 676
Sodium glutamate 676
Sodium glycérophosphate pentahydraté 676
Sodium hexafluorosilicate 677
Sodium hexamétaphosphate 677
Sodium hydrosulfite 678
Sodium hydroxyde, perles 678
Sodium hydroxyde, pastilles 679
Sodium hydroxyde sur silice 680
Sodium hydroxyde solution 50% 680
Sodium hydroxyde solution 40% 680
Sodium hydroxyde solution 35-37% 681
Sodium hydroxyde solution 35% 681
Sodium hydroxyde solution 32% 681
Sodium hydroxyde solution 30% 682
Sodium hydroxyde solution 20% m/v 683
Sodium hydroxyde solution à 20% m/m 683
Sodium hydroxyde solution 10% m/v 683
Sodium hydroxyde solution 5% m/v 684
Sodium hydroxyde 6 mol/l (6N) 684
Sodium hydroxyde 5 mol/l (5N) 685
Sodium hydroxyde 4 mol/l (4N) 685
Sodium hydroxyde 3 mol/l (3N) 685
Sodium hydroxyde 2 mol/l (2N) 686
Sodium hydroxyde 1.2 mol/l (1.2N) 686
Sodium hydroxyde 1 mol/l (1N) 686
Sodium hydroxyde 0.7 mol/l (N/1.4) 687
Sodium hydroxyde 0.5 mol/l (0.5N) 688
Sodium hydroxyde 0.357 mol/l (0.357N) 688
Sodium hydroxyde 0.25 mol/l (0.25N) 688
Sodium hydroxyde 0.2 mol/l (0.2N) 689
Sodium hydroxyde 1/9 mol/l (N/9) 689
Sodium hydroxyde 0.1 mol/l (0.1N) 689
Sodium hydroxyde 0.1 mol/l (0.1N) dans l'éthanol 690
Sodium hydroxyde 0.01 mol/l (0.01N) 690
Sodium hydroxyde solution 691
Sodium hydroxyde solution dans le méthanol 691
Sodium hydroxyde-d 30% 691
Sodium hydroxyde-d 1N 691
Sodium hypochlorite solution aqueuse 692
Sodium hypochlorite solution 12.5% 692
Sodium hypophosphite 693
Sodium iodure 693
Sodium laurylsulfate 693
Sodium métabisulfite 694
Sodium métaphosphate 695
Sodium méthanolate 0.1 mol/l 695
Sodium molybdate dihydraté 695
Sodium nitrate 696
Sodium nitrite 696
Sodium nitrite 0.1 mol/l (0.1N) 697
Sodium nitrite 500g/l 697
Sodium nitroprussiate dihydraté 697
Sodium oxalate 698
Sodium perborate tétrahydraté 698
Sodium perchlorate monohydraté 699
Sodium m-périodate 694

INDEX ALPHABÉTIQUE

Sodium péroxyde 699
Sodium persulfate 699
Sodium persulfate 1 mol/l 700
Sodium phosphate dibasique anhydre 700
Sodium phosphate dibasique dihydraté 701
Sodium phosphate dibasique dodécahydraté 702
Sodium phosphate monobasique monohydraté 703
Sodium phosphate monobasique dihydraté 703
Sodium phosphate tribasique dodécahydraté 704
Sodium o-Phosphite pentahydraté 704
Sodium pyrophosphate décahydraté 705
Sodium salicylate 705
Sodium silicate 705
Sodium succinate hexahydraté 706
Sodium sulfate anhydre 706
Sodium sulfate décahydraté 707
Sodium sulfite anhydre 708
Sodium sulfocyanure 709
Sodium sulfure nonahydraté 707
Sodium sulfure nonahydraté solution 707
Sodium sulfure trihydraté 708
Sodium tartrate acide 678
Sodium tartrate dihydraté 709
Sodium tétraborate anhydre 709
Sodium tétraborate décahydraté 710
Sodium tétraphénylborate 710
Sodium thiosulfate anhydre 711
Sodium thiosulfate pentahydraté 711
Sodium thiosulfate 1 mol/l (1N) 712
Sodium thiosulfate 0.5 mol/l (0.5N) 712
Sodium thiosulfate 0.2 mol/l (0.2N) 712
Sodium thiosulfate 0.1 mol/l (0.1N) 712
Sodium thiosulfate 0.0394 mol/l (0.0394N) 713
Sodium thiosulfate 0.0197 mol/l (0.0197N) 713
Sodium thiosulfate 0.01 mol/l (0.01N) 713
Sodium tungstate dihydraté 714
Solution ADF 147
Solution amoniacale de tétramminecuivre 165
Solution cupri-citrique 284
Solution cupri-tartrique 285
Solutions étalon pour le degré de coloration des liquides 291
Solutions primaires pour le degré de coloration des liquides 290
Solution NDF 519
Solution NDF PLUS 519
Solution réductice 220 mV à 25°C 636
Solution réductice 468 mV à 25°C 637
Solution Rhodol 639
Solution tampon borique 207
Solvant Plus 714
D-Sorbitol 715
Sorbitol 714
Sorbitol solution 70% 715
Soudan III 726
Soudan III solution hydroalcoolique 727
Soufre sublimé lavé ventilé 730
Soufre solution standard 730
Standard organique: mélange multiéléments HAP 542
Standard organique: mélange multiéléments PCB 543
Standard organique: mélange de pesticides 546
Standard organique PCB 540
Standards de point de fusion 483
Standard turbidité 4000NTU Formazine 780
Standards de viscosité 785
Strontium solution standard 722
Strontium acétate 723
Strontium bromure monohydraté 723
Strontium carbonate 724
Strontium chlorure hexahydraté 724
Strontium nitrate 724
Strontium sulfate 725
Sulfanilamide 728
Sulfate solution standard 728
Sulfite standard solution 729
Sulfolane 729

T

Suspension mère d'opalescence 621
Talc 744
Tamis moléculaire 3 A 510
Tamis moléculaire 4 A 510
Tampon pH 1 216
Tampon pH 1.68 216
Tampon pH 2 217
Tampon pH 3 217
Tampon pH 3.5 217
Tampon pH 3.56 218
Tampon pH 3.7 218
Tampon pH 4 218
Tampon pH 4.5 219
Tampon pH 4.62 219
Tampon pH 5 219
Tampon pH 5.2 219
Tampon pH 6 220
Tampon pH 6.8 220
Tampon pH 6.88 220
Tampon pH 7 220
Tampon pH 7.2 221
Tampon pH 7.20 selon Weise 221
Tampon pH 7.4 222
Tampon pH 8 222
Tampon pH 9 222
Tampon pH 9.22 223
Tampon pH 10 223
Tampon pH 10.06 224
Tampon pH 11 224
Tampon pH 12 225
Tampon pH 13 225
Tampon acétate pH 6.0 130
Tampon acétate pH 4.6 130
Tampon acétate pH 4.5 216
Tampon ammoniacal pH10 165
Tampon chlorure d'ammonium pH 9.5 170
Tampon chlorure d'ammonium pH 10.0 170
Tampon chlorure d'ammonium pH 10.7 169
Tampon phosphate pH 2.0 576
Tampon phosphate pH 3.0 576
Tampon phosphate pH 5.5 575
Tampon phosphate pH 6.0 575
Tampon phosphate pH 6.8 575
Tampon phosphate pH 7.2 575
Tampon phosphate pH 7.4 575
Tampon phosphate pH 9.0 574
Tampon pour ajustement de la force ionique totale 769
Tampon tris(hydroxyméthyl)aminométhane pH 8.1 777
Tampon tris(hydroxyméthyl)aminométhane-chlorure de sodium pH 7.4 778
Tampon tris(hydroxyméthyl)aminométhane-EDTA pH 8.4 777
Tantale solution standard 744
Tartrazine 746
Tellure en poudre 747
Tellure en morceaux 746
Tellurium solution standard 747
Terbium solution standard 748
Terre de silice composée 449
Tetrabutylammonium bisulfate 748
Tétrabutylammonium bromure 749
Tétrabutylammonium hydroxyde 0.1 mol/l (0.1N) 749
Tétrabutylammonium hydroxyde 0.1 mol/l (0.1N) dans propanol-2 749
Tetrabutylammonium hydroxyde 0.1 mol/l (0.1N) dans méthanol / propanol-2 (50/50) 750
Tétrachloroéthane-d2 751
Tétrachloroéthylène 750
Tétraéthylammonium bromure 751
Tétrahydrofuranne 751
Tétrahydrofuranne-d8 753
Tétraméthylammonium hydroxyde 10% 753
N,N,N',N'-Tétraméthyl-p-phénylènediamine dichlorhydraté 753
Thallium solution standard 754
Thioacétamide 754
Thioacétamide solution 40 g/l 755
Thiourée 756
Thorium solution standard 756
Thullium solution standard 756
Thymol 757
Thymolphtaléine 758
Thymolphtaléine 0.1% solution hydroalcoolique 759
Thymolphtaléine solution 0.1% dans l'éthanol 759
Tisab 762
Tisab III solution 763
Titane solution standard 763
Titane dioxyde 764
Titane isopropylate 764
Titane trichlorure-acide sulfurique 765
o-Tolidine en solution à 0,1% 765
o-Tolidine en solution 765
Toluène 765
Toluène en solution dans l'hexane 768
Toluène-d8 768
p-Toluènesulfonamide 768
Tournesol 465
Triacétine 769
1,2,4-Trichlorobenzène 771
2,2,2 - Trichloroéthanol 770
Tricrésol 771
Triéthanolamine 771
Triéthylamine 772
Triméthylcétylammonium bromure 775
Triméthylcétylammonium p-toluènesulfonate 775
Trioctylphosphine oxyde 775
Triphénylphosphine 776
2,3,5-Triphényltétrazolium chlorure 776
Tris(hydroxyméthyl)aminométhane 776
Tris(hydroxyméthyl)aminométhane chlorhydraté 778
Tris(hydroxyméthyl)aminométhane solution 777
Triton© X100 solution 778
Tropéoline O 778
Tropéoline O solution 0.1% 779
Tungstène solution standard 779

U

L(-)Tyrosine 781
Uranium solution standard 782

V

Urée 782
Vanadium solution standard 783
Vanadium (V) pentaoxyde 784
Vanilline 784
Vaniline phosphorique solution 784
Vaseline blanche 556
Vert brillant 209
Vert de bromocrésol 211
Vert de bromocrésol 0.04% solution hydroalcoolique 211
Vert de bromocrésol solution 211
Vert de bromocrésol - Rouge de méthyle solution 212
Vert de malachite 476

DASITGROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

Vert malachite solution 0.5% dans l'acide acétique anhydre	476
Vert de méthyle	499
Vert lumière SF	461
Vert solide FCF	353
Violet cristallisé	284
Violet cristallisé solution 0.5% dans l'acide acétique anhydre	284
Violet cristallisé oxalate pour kit de Gram-Hucker	284
Violet de gentiane	368
Violet de gentiane solution phéniquée	368

X

Violet de pyrocatéchol	632
o-xylène	792
Xylène, mélange d'isomères	791
Xylènecyanol FF	792
2,4-Xylenol	792

Y

D(+)Xylose	793
Ytterbium solution standard	794

Z

Yttrium solution standard	794
Zinc, granule	796
Zinc, lames	796
Zinc, poudre	797
Zinc solution standard	797
Zinc acétate dihydraté	798
Zinc activé	797
Zinc carbonate basique	799
Zinc chlorure anhydre	799
Zinc chlorure solution à 60%	799
Zinc dibenzyldithiocarbamate	800
Zinc nitrate hexahydraté	800
Zinc oxyde	801
Zinc stéarate	801
Zinc sulfate heptahydraté	802
Zinc sulfate monohydrate	802
Zinc sulfate 0.1 mol/l (0.2N)	802
Zinc sulfate 0.05 mol/l (0.05N)	803
Zinc sulfure	803
Zirconium en poudre	803
Zirconium solution standard	804
Zirconium (IV) oxyde	804
Zirconyle nitrate	805

BEST CHOICE

Tailor made

Chemicals such as you imagine

CARLO ERBA Reagents thanks to its experience and production flexibility manufacture tailor made products for industry and laboratories.



- Custom mixture
- Purification
- Custom packaging
- Deliveries in bulk tanks and isotanks
- Shuttle service
- Analytical control
- Batch management
- Quality assurance

www.carloerbareagents.com

DASIT GROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

ÍNDICE ALFABÉTICO

A

Absorbente para líquidos derramados 129
Aceite de cedro condensado 538
Aceite de inmersión 420
Aceite de ricino 249
Aceite de vaselina 555
Aceite refinado de almendras 539
Acetal 129
Acetamida 129
Acetanilida 129
Acetilacetona 146
n-Acetil-L-cisteína 146
p-Acetilfenetidina 146
Acetofenona 146
Acetona 138
Acetona-d6 141
Acetona / Agua 98/2 (v/v) + Azul de bromofenol 0.02 g/l 141
Acetonitrilo 142
Acetonitrilo+ 0.1% v/v acido formico 145
Acetonitrilo+ 0.1% v/v acido trifluoroacético 145
Acetonitrilo-d3 145
Acido acético glacial 130
Acido acético 96% 133
Acido acético 80% 134
Acido acético 45 % 134
Acido acético 30% 134
Acido acético 27 % 135
Acido acético 25 % 135
Acido acético 20 % 135
Acido acético 12 % 135
Acido acético 1 mol/l (1N) 136
Acido acético 0.1 mol/l (0.1N) 136
Acido acético 0.03 mol/l (0.03N) 136
Acido acético-d4 137
Acido adípico 147
Acido p-aminobenzoico 158
Acido 1-amino-2-naftol-4-sulfonico 158
Acido anísico 181
Acido L(+)ascórbico 186
Acido L(+)aspártico 186
Acido barbitúrico 191
Acido benzóico 199
Acido bórico 205
Acido bórico 1% con indicador 207
Acido bórico 3% 206
Acido bórico 4% 206
Acido bórico 4% con indicador 206
Acido bórico 20g/l 207
Acido bórico 20 g/l con indicador 207
Acido bromhídrico 48% 393
Acido 1-butanosulfónico sal sódica 225
Acido n-butírico 231
Acido calconcarbónico 245
Acido n-caprílico 247
Acido n-caproico 247
Acido citrico anhidro 269
Acido cítrico monohidrato 269
Acido cloranílico 257
Acido clorhidrico 50% v/v 393
Acido clorhidrico 37% 394
Acido clorhidrico 34-37% 396
Acido clorhidrico 32-35% 396
Acido clorhidrico 32% (20°Bé) 397
Acido clorhidrico 32% 397
Acido clorhidrico 29-31 % 398
Acido clorhidrico 26% 398
Acido clorhidrico 25% p/v 398
Acido clorhidrico 23% 399
Acido clorhidrico 20% 399
Acido clorhidrico 12% 399
Acido clorhidrico 10% 400
Acido clorhidrico 8% 400
Acido cloridrico 5% 400
Acido clorhidrico 1.128% p/v 401
Acido cloridrico 9 mol/l (9N) 401
Acido clorhidrico 6 mol/l (6N) 401
Acido clorhidrico 5 mol/l (5N) 402
Acido clorhidrico 4 mol/l (4N) 402
Acido clorhidrico 3 mol/l (3N) 403
Acido clorhidrico 2 mol/l (2N) 403
Acido clorhidrico 1 mol/l (1N) 403
Acido clorhídrico 0.714 mol/l (N/1.4) 404
Acido clorhidrico 0.5 mol/l (0.5N) 405
Acido clorhidrico 0.2 mol/l (0.2N) 405
Acido clorhidrico 0.2 mol/l (0.2N) en propan-2-ol 405
Acido clorhidrico 0.1 mol/l (0.1N) 406
Acido clorhidrico 0.1 mol/l (0.1N) en etanol 407
Acido clorhidrico 0.1 mol/l (0.1N) en propan-2-ol 407
Acido clorhídrico 0.0714 mol/l (N/14) 407
Acido clorhidrico 0.05 mol/l (0.05N) 407
Acido clorhidrico 0.04 mol/l (0.04N) 408
Acido clorhidrico 0.02 mol/l (0.02N) 408
Acido clorhidrico 0.01 mol/l (0.01N) 408
Acido clorhídrico diluido 409
Acido clorhidrico-d 1 mol/l 409
Acido clorhidrico-d 20% 409
Acido clorhídrico brominado 409
Acido cloroacetico 258
Acido cromotrópico sal disódica 268
Acido 1-decanosulfónico sal sódica 290
Acido dicloroacético 297
Acido dietilentriaminopentaacético 308
Acido 1-dodecanosulfónico sal sódica 324
Acido dodecilbencensulfónico sal sódica 325
Acido estearico 722
Acido etilendiaminotetracético 344
Acido etilendiaminotetracético sal dipotasica dihidrato 345
Acido etilendiaminotetracético sal dipotasica-magnésica dihidrato 347
Acido etilendiaminotetracético sal disódica 345
Acido etilendiaminotetracético sal disódica 0.1 mol/l (0.2N) 346
Acido etilendiaminotetracético sal disódica 0.05 mol/l (0.1N) 347
Acido etilendiaminotetracético sal disódica 0.01 mol/l (0.02N) 347
Acido etilendiaminotetracético sal tetrasódica tetrahidrato 348
Acido etilenglicol-bis (2-aminoetil-éter) tetracético 349
Acido etilenodiaminatetracético sal tripotasíca 348
Acido fluorhídrico 50% 410
Acido fluorhídrico 47-51% 411
Acido fluorhídrico 39.5% 412
Acido fluorhídrico 0.1M 412
Acido fluorhídrico diluido 412
Acido fórmico 99% 362
Acido fórmico 85% 363
Acido fórmico-d 364
Acido fosfomolíbdico 576
Acido fosfonico 577
Acido fosfosulfúrico 578
Acido fosfotúngstico 578
Acido fosfotúngstico solución 578
Acido ftálico 579
Acido fumárico 365
Acido gálico monohidrato 367
Acido glicólico 374
Acido L(+)glutámico 371
Acido heptafluorobutírico 380
Acido 1-heptanosulfónico sal sódica 383
Acido hexacloroplatínico hexahidrato 383
Acido 1-hexanosulfónico sal sódica 389
Acido 1-hexanosulfónico sal sódica monohidrato 389
Acido hipofosforoso 50% 418
Acido idríodico 57% 393
Ferron 354
Acido L(+)-láctico 452
Acido maleico 477
Acido DL-málico 477
Acido malónico 478
Acido D(-)mandélico 478
Acido metafosfórico 487
Acido metanosulfónico 487
Acido alfa-metoxifenilacético 495
Acido nicotínico 525
Acido nítrico fumante 99% 528
Acido nítrico 69.5% 528
Acido nítrico 67-70% 529
Acido nítrico 67-69% 530
Acido nítrico 67.5% (42° Be) 530
Acido nítrico 65% 531
Acido nítrico 18% 532
Acido nítrico 10% 532
Acido nítrico 8 mol/l (8N) 532
Acido nítrico 2 mol/l (2N) 533
Acido nítrico 1 mol/l (1N) 533
Acido nítrico 0.1 mol/l (0.1N) 533
Acido nítrico diluido 534
Acido nitrilotriacético 534
Acido m-nitrobenzoico 535
Acido 1-octanosulfónico sal sódica 537
Acido 1-octanosulfónico sal sódica monohidrato 538
Acido oleico 539
Acido orto-fosforico 99% 547
Acido orto-fosforico 85% 547
Acido orto-fosfórico 75% 549
Acido orto-fosfórico 10% 549
Acido orto-fosforico-d3 85% in D2O 549
Acido oxálico dihidrato 550
Acido oxálico 0.5 mol/l (1N) 551
Acido oxálico 0.05 mol/l (0.1N) 551
Acido oxálico 0.005 mol/l (0.01N) 551
Acido 1-pentanosulfónico sal sódica 559
Acido 1-pentanosulfónico sal sódica monohidrato 560
Acido perclórico 65-71% 560
Acido perclórico 65% 561
Acido perclórico 0.1 mol/l (0.1N) 561
Acido perclórico 0.1 mol/l (0.1N) en acido acético 562
Acido perclórico 0.05 mol/l (0.05N) 562
Acido perclórico 0.01 mol/l (0.01N) 562
Acido perclórico solución 562
Acido periódico 563
Acido pícrico solución 579
Acido pirrolidinditiocarbámico sal de amonio 633
Acido 1-propanosulfónico sal sódica 627
Acido propiónico 628
Acido p-rosólico 640
Acido salicílico 642
Acido sebácico 646
Acido selenioso 648
Acido silicotúngstico 652
Acido succínico 725
Acido sulfamico 727
Acido sulfanílico 728
Acido sulfosalicílico 729
Acido sulfúrico 98% 730
Acido sulfúrico 96% (66°Bé) 731
Acido sulfurico 96% 731
Acido sulfúrico 95 - 97 % 733
Acido sulfúrico 93-98% 733
Acido sulfúrico 90% 734
Acido sulfúrico 85% 734

ÍNDICE ALFABÉTICO

Acido sulfúrico 72% 735
Acido sulfúrico 69% 735
Acido sulfúrico 62% 735
Acido sulfúrico 50% 735
Acido sulfúrico 35% (30°Bé) 736
Acido sulfúrico 30% 736
Acido sulfúrico 25% 736
Acido sulfúrico 20% 736
Acido sulfúrico 10% V/V 737
Acido sulfúrico 4 mol/l (8N) 737
Acido sulfúrico 2.5 mol/l (5N) 737
Acido sulfúrico 1 mol/l (2N) 737
Acido sulfúrico 0.5 mol/l (1N) 738
Acido sulfúrico 0.33 mol/l (2N/3) 738
Acido sulfúrico 0.26 mol/l (0.52N) 739
Acido sulfúrico 0.25 mol/l (0.5N) 739
Acido sulfúrico 0.166 mol/l (0.333N) 739
Acido sulfúrico 0.13 mol/l (0.26N) 740
Acido sulfúrico 0.125 mol/l (0.25N) 740
Acido sulfúrico 0.1 mol/l (0.2N) 740
Acido sulfúrico 0.05 mol/l (0.1N) 740
Acido sulfúrico 0.025 mol/l (0.05N) 741
Acido sulfúrico 0.02 mol/l (0.04N) 741
Acido sulfúrico 0.01 mol/l (0.02N) 741
Acido sulfúrico 0.005 mol/l (0.01N) 742
Acido sulfúrico 0.0025 mol/l (0.005N) 742
Acido sulfúrico d=1,820 743
Acido sulfúrico con 10 g/l de plata sulfato 742
Acido sulfúrico con 6.6 g/l de plata sulfato 743
Acido sulfúrico, diluido 743
Acido tánico 744
Acido L(+)tartárico 745
Acido L(+)tartárico solución 20% en agua 746
Acido 2-tiobarbitúrico 755
Acido tioglicólico 80% 755
Acido p-toluenosulfónico 768
Acido p-toluenosulfónico sal sódica 769
Acido tricloroacetico 770
Acido tricloroacético solución 20% 770
Acido trifluoroacético 773
Acido trifluoroacético-d 774
Acido túngstico 780
Agua 786
Agua purificada 788
Agua + 0.1%(v/v) acido fórmico 789
Agua + 0,1% v/v acido trifluoroacético 789
Agua de bromo 211
Agua de cloro 789
Agua desionizada acidificada 789
ß-Alanina 148
Albúmina de huevo polvo 148
Albúmina de huevo desecada 148
Alcanfor natural 246
Alcanfor sintetico 246
Alcohol n-amílico 180
Alcohol ter-amílico 180
Alcohol bencílico 200
Alcohol ter-butílico 227
Alcohol cetílico 254
Alcohol 2-feniletílico 573
Alcohol isoamílico 436
Aldehído anísico 181
Aldehído benzoico 197
Aldehído propiónico 628
Aldehído anísico solución 181
Aleación según Devarda 294
Aleación según Raney 635
Aliltiourea 150
Alizarina 149
Alizarina solución saturada en alcohol etilico 149
Almidón en pasta solución 1% 721
Almidón de arroz 640
Almidón de maiz 476
Almidón de yoduro de potasio papel 608
Almidón soluble 721
Almidón soluble solución 721
Alúmina blanca 151
Aluminum, solución patrón 151
Aluminio, polvo 151
Aluminio cloruro anhidro 153
Aluminio cloruro hexahidrato 153
Aluminio hidróxido 154
Aluminio nitrato nonahydrate 154
Aluminio óxido 155
Aluminio óxido (ácido) 155
Aluminio óxido (básico) 155
Aluminio óxido (neutro) 156
Aluminio óxido activo 156
Aluminio sulfato 157
Aluminio y amonio sulfato dodecahidrato 152
Aluminio y potasio sulfato dodecahidrato 156
Amarillo Clayton 270
Amarillo de alizarina R 150
Amarillo de metanilo 487
Amarillo de metilo 506
Amarillo naftol S 518
Amarillo Sudán 727
4-Aminofenazona 159
m-Aminofenol 159
p-Aminofenol 159
Aminopirazolone en solución 159
Amonio, solución patrón 165
Amonio acetato 166
Amonio bicarbonato 167
Amonio bifluoruro 172
Amonio bromuro 167
Amonio carbamato 168
Amonio carbonato 168
Amonio carbonato solución 158 g/l 168
Amonio cloruro 169
Amonio citrato dibásico 170
Amonio citrato solución 20% 171
Tri-amonio citrato 170
Amonio di-Hidrógeno Fosfato solución 25 mg/l 172
Amonio hidróxido 32% 160
Amonio hidróxido 30% 160
Amonio hidróxido 28% 161
Amonio hidróxido 25% 162
Amonio hidróxido 20 - 22% 163
Amonio hidróxido 17% 164
Amonio hidróxido 10% 164
Amonio hidróxido solución diluido 164
Amonio fluoruro 171
Amonio formiato 171
Amonio fosfato dibásico 176
Amonio fosfato monobásico 176
Amonio (II) hierro II sulfato 0.1 N 430
Amonio Hierro (II) sulfato 0.12N 429
Amonio ioduro 173
Amonio molibdato tetrahidratado 173
Amonio molibdato solución 2.5% en acido nítrico 174
Amonio nitrato 174
Amonio nitrato 200 mg/l solución 174
Amonio oxalato monohidrato 175
Amonio oxalato solución 4% 175
Amonio persulfato 175
Amonio solución 6N 164
Amonio sulfamato 177
Amonio sulfato 177
Amonio sulfocianuro 178
Amonio sulfocianuro 1 mol/l (1N) 179
Amonio sulfocianuro 0.1 mol/l (0.1N) 179
Amonio sulfocianuro 0.01 mol/l (0.01N) 179
Amonio sulfuro solución 20% 178
Amonio L(+)tartrato 178
Anaranjado de acridina 147
Anhídrido acético 137
Anhídrido fosfórico 578
Anhídrido ftálico 579
Anhídrido maleico 477
Anhídrido molibdénico 511
Anhídrido succínico 725
Anhídrido trifluoroacético 774
Anhídrido túngstico 780
Anilina clorhidrato 181
Antimonio, solución patrón 182
Antimonio tricloruro 183
Antimonio y potasio tartrato 183
Antrona 182
Arena de Fontainebleau 644
Arena purificada 644
L(+)Arginina 184
L(+)Arginina monoclorhidrato 184
L(+)Asparagina 186
Arsenazo III 184
Arsenic, solución patrón 185
Arsénico trióxido solución 185
Ausilab 110 187
Ausilab 140 187
Ausilab 210 187
Ausilab 250 188
Ausilab 260 188
Ausilab 280 188
Ausilab 290 188
Ausilab 400 189
Ausilab 500 189
Azometino H 189
Azufre, solución patrón 730
Azufre sublimado lavado 730
Azul Coomassie R 250 274
Azul de alcali 6B 150
Azul de alcali 6B solución 2% en alcohol etilico 150
Azul de alcian 8GS 1% 148
Azul de alcian 8GX 149
Azul de anilina soluble en agua 180
Azul de bromofenol 213
Azul de bromofenol solución 0.4% en alcohol etilico 213
Azul de bromofenol solución 0.02% 214
Azul de bromofenol solución 214
Azul de bromotimol 215
Azul de bromotimol solución 0.4% en alcohol etilico 215
Azul de bromotimol 0.02% 215
Azul de bromotimol solución 215
Azul de cresol brillante 208
Azul lactofenol solución 452
Azul de metileno 498
Azul de metileno solución 1% 499
Azul de metileno solución saturada en alcohol etilico 498
Azul de metileno fenicado solución hidroalcohólica 248
Azul de metilo 497
Azul de metiltimol sal sódica 506
Azul de metileno fenicado solución hidroalcohólica 248
Azul de timol 757
Azul de timol solucion 0.4% en alcohol etilico 757
Azul de timol solución 758
Azul de toluidina 769
Azul Victoria 785
Azur II 189
Azur II Eosina 190

B

Balsamo del Canada 247
Bario, solución patrón 191
Bario acetato 192
Bario carbonato 192

DASITGROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

Bario cloruro dihidratado 192
Bario cloruro 0.1 mol/l (0.2N) 193
Bario cloruro solución 10% 193
Bario cloruro solución 61 g/l 193
Bario cloruro 30 g/l 194
Bario cromato 194
Bario hidróxido octahidrato 194
Bario hidróxido solución 5% 195
Bario hidróxido solución 47.3 g/l 195
Bario nitrato 195
Bario perclorato trihidrato 195
Bario perclorato 0.05 mol/l 196
Bario perclorato 0.025 mol/l 196
Bario sulfato 196
Benceno 198
Benceno-d6 198
Bencenosulfonilo cloruro 198
Bencilo benzoato 200
Benzalconio cloruro 197
Benzina de Petroleo E 563
p-benzoquinona 199
Berilio, solución patrón 201
N,O-Bis (trimetilsilil) acetamida 204
N,O-Bis (trimetilsilil) trifluoroacetamida 204
Bismuto, gránulos 202
Bismuto, solución patrón 202
Bismuto (III) carbonato basico 203
Bismuto (III) nitrato basico 203
Bismuto (III) nitrato pentahidratado 204
Biuret 97% 205
Biuret reactivo 205
Boro, solución patrón 208
Brix patrones 209
Bromato, solución patrón 210
Bromo, solución patrón 210
alpha-Bromonaftaleno 213
Bromuro, solución patrón 210
Bromuro - bromato 0.0167 mol/l 210
Butanol-1 226
Butan-2-olo 227
n-Butilo acetate 228
n-Butilo cloruro 229
Butil hidroxitolueno 230
ter-Butilmetiléter 230
2-Butoxietanol 228
2-(2-Butoxietoxi) etanol 228

C

Cadmio, solución patrón 232
Cadmio acetato dihidrato 233
Cadmio carbonato 233
Cadmio cloruro monohidrato 233
Cadmio nitrato tetrahidratado 234
Cafeína anhidra 234
Cal sodada 660
Calcio, solución patrón 234
Calcio acetato anhidro 236
Calcio acetato monohidrato 236
Calcio carbonato 236
Calcio cloruro anhidro 237
Calcio cloruro dihidrato 238
Calcio cloruro hexahidrato 239
Calcio cloruro solución 0.025% 239
Calcio estearato 244
Calcio fluoruro 240
Calcio fosfato dibásico dihidrato 243
Calcio fosfato monobásico monohidrato 243
Calcio fosfato tribásico 243
Calcio gluconato 240
Calcio hidróxido 240
Calcio hidróxido acuoso 184
Calcio lactato 241
Calcio nitrato tetrahidrato 241
Calcio óxido, polvo 242
Calcio óxido, trozos 242
Calcio pantotenato 242
Calcio sulfato dihidrato 244
Calcio sulfato hemihidrato 245
Calcio sulfato hemihidrato solución 245
Calcón 245
Calmagita 246
Caolín lavado 443
Carbón activado 255
Carbón decolorante 255
Carbón vegetal 256
Carborundo, polvo 248
Carborundo, gránulos 248
Carrez reactivo sal de zinc 249
Carrez reactivo sal de potasio 249
Caseina 249
Catalizador Kjeldahl al selenio 450
Catalizador Keldahl al titanio 451
Catalizador Kjeldahl para análisis del agua 450
Catalizador Kjeldahl según Wieninger 450
Catalizador Kjeldahl sin selenio y potassio 450
Celulosa, polvo 250
Cerio, solución patrón 250
Cerio (III) nitrato hexahidratado 251
Cerio (IV) amonio nitrato 251
Cerio (IV) amonio nitrato 0.01 mol/l 251
Cerio (IV) amonio nitrato 0.1 mol/l 251
Cerio (IV) amonio sulfato dihidrato 252
Cerio (IV) amonio sulfato 0.01 mol/l 252
Cerio (IV) amonio sulfato 0.1 mol/l 252
Cerio (IV) sulfato tetrahidrato 252
Cerio (IV) sulfato 0.1 mol/l 253
Cesio, solución patrón 253
Cesio cloruro 254
Cesio cloruro solución 25 g/l 254
Cesio sulfato 254
Cianuro, solución patrón 285
Ciclohexano 286
Ciclohexano-d12 287
Ciclohexanol 287
Ciclohexanona 288
Ciclohexilamina 288
Ciclopentil metil eter 289
Cinconina 268
L-Cisteína 289
L-Cistina 289
Citral 268
Cloral hidrato 256
Cloramina T sódica 256
Cloranilo 256
Clorato, solución patrón 257
Clorita, solución patrón 257
Cloroacetamida 258
Clorobenceno 258
p-Clorobenzaldehído 258
Clorobutanol hemihidrato 259
p-Clorofenol 264
Cloroformo 259
Cloroformo-d 263
Cloroformo-d + 0.03% TMS 263
Cloruro, solución patrón 257
Cobalto, solución patrón 270
Cobalto (II) acetato tetrahidrato 271
Cobalto (II) y amonio sulfato hexahidrato 271
Cobalto (II) cloruro hexahidrato 272
Cobalto (II) cloruro en solución 272
Cobalto (II) nitrato hexahidrato 272
Cobalto (II) sulfato heptahidrato 273
Cobre electrolítico, hilos 275
Cobre electrolítico, virutas 274
Cobre electrolítico recocido, hojas 274
Cobre reducido, polvo 275
Cobre, solución patrón 275
Cobre (I) cloruro 276
Cobre (I) cloruro solución 7% en amonio hidróxido 277
Cobre (I) ioduro 277
Cobre (II) acetato hidrato 277
Cobre (II) carbonato básico 278
Cobre (II) cloruro dihidrato 278
Cobre (II) dietilendiaminico hidroxido 1 mol/l 285
Cobre etilendiamina solución 285
Cobre (II) nitrato trihidrato 279
Cobre (II) óxido 279
Cobre (II) óxido hilos 279
Cobre (II) sulfato anhidro 280
Cobre (II) sulfato pentahidrato 280
Cobre (II) sulfato solución 12.5% 281
Cobre (II) sulfato de en solución 281
Cobre-cítrico solución 284
Cobre-tartárico solución 285
Cobre tetrammina solución amoniacale 165
Colorante de Wright solución in alcohol metilico 790
m-Cresol 282
o-Cresol 282
Crisoidina Y 268
Cromato, solución patrón 264
Cromo, solución patrón 264
Cromo (III) cloruro hexahidrato 265
Cromo (III) nitrato nonahidrato 266
Cromo (III) de potasio sulfato dodecahidrato 266
Cromo (III) óxido 266
Cromo (III) sulfato 267
Cromo (VI) óxido 267
Cromo (VI) óxido solución 267

D

Decano 290
Descalcificante 290
Deuterio oxido-d2 293
Deuterio oxido-d2 + 0.01% DMSO 294
Deuterio oxido-d2 + 0.5% TSP d4 294
Deuterio oxido-d2 + 0.03% TSP d4 294
Diacetildioxima 295
Diacetildioxima solución 1% en alcohol etilico 296
Diacetildioxima sal sódica 295
Diacetonalcohol 295
Dibenzoilmetano 296
Dibenzoilo peróxido 296
n-Dibutiloftalato 297
Dicalite 4158 297
m-Diclorobenceno 298
o-Diclorobenceno 298
p-Diclorobenceno 298
1,2-Diclorobenceno-d4 299
1,2-Dicloroetano 299
2,6-Diclorofenol-indofenol sal sódica 305
Diclorometano 300
Diclorometano-d2 305
Diclorometano acidificado con ácido clorhídrico 1% 305
2,6-Dicloroquinona-4-clorimida 305
Dietanolamina 306
Dietilamina 306
Dietilcarbonato 307
Dietileter 308
Dietil ftalato 310
n,n-Dietil-p-fenilendiamonio sulfato 310
Difenilamina solución 1% en ácido sulfúrico 322
4-Difenilamina sulfonato sal sódica 322
sim-Difenilcarbacida 323
sim-Difenilcarbazona 323
Difeniltiocarbazona 323
Diferenciador para kit Gram-Hucker 311
Digitonina 311
Diisopropilamina 311
N,N-Diisopropiletilamina 312

ÍNDICE ALFABÉTICO

Dimedona 313
N,N-Dimetilacetamida 313
p-Dimetilaminobencilidenrodanina 314
p-Dimetilaminobenzaldehído 314
N,N-Dimetil formamida 314
N,N-Dimetilformamida-d7 316
N,N-Dimetilformamida dimetilacetal 317
N,N'-Dimetilpropileno urea 317
Dimetilsulfóxido 317
Dimetilsulfóxido-d6 319
Dimetilsulfóxido-d6 + 0.03% TMS 319
1,2-Dimetoxietano 313
Dimidio bromuro 320
2,4-Dinitro-1-clorobenceno 320
2,4-Dinitrofenilhidracina (con 30% de agua) 320
1,4-Dioxano 321
1,3-Dioxolano 322
2,2'-Dipiridilo 323
Disolvente Plus 714
Disprosio, solución patrón 325
Ditiooxamida 324
Dodeciltrimetilamonio bromuro 325

E

Ehrlich reactivo 327
Eluyente sodio carbonato 327
Eluyente sodio carbonato/sodio dicarbonato 328
Eluyente sodio dicarbonato 327
Emalume reactivo solución según Carazzi 378
Emalume reactivo solución según Mayer 378
Eosina B 328
Eosina Y 328
Eosina Y 1% solución acuosa 328
Eosina Y 0.5% solución alcohólica 329
Erbio, solución patrón 329
Eriocromocianina R 330
Eritrosina extra B 330
Esbach reactivo 330
Escandio, solución patrón 645
Estaño, solución patrón 760
Estaño, polvo 759
Estaño hojas 759
Estaño (II) cloruro dihidrato 761
Estaño (II) cloruro solución 10% 761
Estaño (II) cloruro solución 761
Estaño (II) sulfato 762
Estaño (IV) cloruro pentahidratado 762
Estaño (IV) óxido 762
Estroncio, solución patrón 722
Estroncio acetato 723
Estroncio bromuro monohidrato 723
Estroncio carbonato 724
Estroncio cloruro hexahidratado 724
Estroncio nitrato 724
Estroncio sulfato 725
Etanol absoluto anhidro 331
Etanol 96° 334
Etanol 70% v/v 336
Etanol 60% v/v 337
Etanol absoluto desnaturalizado 338
Etanol 95° desnaturalizado 338
Etanol 94° desnaturalizado 339
Etanol 90° desnaturalizado 339
Etanol 70° modificada 339
Etanol 50% v/v 337
Etanol anhidro-d6 340
Etanol 95% - Etil éter 50/50 (m/m) - 30 mg/l fenolftaleína 507
Etanolamina 340
Eter de petróleo 100 - 140°C 564
Eter de petróleo 80 - 120°C 564
Eter de petróleo 80 - 100°C 564
Eter de petróleo 75 - 120°C 565
Eter de petróleo 60 - 80°C 565
Eter de petróleo 55 - 85°C 565
Eter de petróleo 40 - 70°C 566
Eter de petróleo 40 - 65°C 566
Eter de petróleo 40 - 60°C 567
Eter de petróleo 35 - 60°C 567
Eter de petróleo 30 - 50°C 568
Eter de petróleo 30 - 40°C 569
Eter di-isopropílico 312
Etilendiamina 344
Etilenglicol 348
Etileno óxido solución 349
Etilo acetato 340
Etilo acetoacetato 343
Etilo formiato 350
Eugenol 351
Eukitt 351
Europio, solución patrón 351

F

Fehling reactivo solución A 353
Fehling reactivo solución B 353
o-Fenantrolina monohidrato 569
o-Fenantrolina-Hierro (II) sulfato solución en acido sulfúrico 569
L-Fenilalanina 572
Fenilhidracina clorhidrato 573
Fenilhidracina clorhidrata solución 573
Fenol 570
Fenol liquido 85% 570
Fenolftaleína 571
Fenolftaleína solución 0.1% 572
Fenolftaleína solución 0.1% en etanol 572
Fenolftaleína solución 1% en alcohol etilico 571
Ferricianuro, solución patrón 353
Ferrocianuro, solución patrón 354
Ferroína solución 0.025 mol/l 354
Florisil 100-200 mesh 357
Florisil 60-100 mesh 356
Floroglucina 574
Floxina B 574
Fluoresceína 357
Fluoresceína sódica 357
Fluoruro, solución patrón 358
Folin-Ciocalteu reactivo 358
Formaldehído 37% p/v 358
Formaldehído 37% p/v neutralizata 359
Formaldehyde 35% p/p 359
Formaldehído 30% p/v 360
Formaldehído 10% v/v según Lillie 360
Formaldehído 5% p/v tamponata pH 6.9 360
Formaldehído 4% p/v con sodio cloruro 361
Formaldehído 4% p/v tamponata pH 6.9 361
Formamida 362
Formol acético 362
Fosfato, solución patrón 574
Fosfomolibdotungstenico reactivo 576
Fósforo, solución patrón 577
D(-)Fructosa 364
o-Ftaldialdehído 579
Fucsina ácida 364
Fucsina básica 365
Fucsina solución decolorizada 365

G

Gadolinio, solución patrón 366
D(+)Galactosa 366
Galio solución estándar 367
Gel de sílice 60A 0,06÷0,20 mm 650
Gel de sílice 60A 6 - 35µ 648
Gel de sílice 60A 20 - 45µ 648
Gel de sílice 60A 35 - 70µ 649
Gel de sílice 60A 40 - 63µ 649
Gel de sílice 60A 70 - 200µ 650
Gel de sílice granulado 650
Gel de sílice granulado sin cobalto 651
Gelatina 368
Germanio, solución patrón 369
Giemsa reactivo 369
Glicerina (30°Bé) 372
Glicerina 90% (28° Bé) 373
Glicina 373
Glicol dietilénico 307
Glicol dietilénico dimetiléter 307
Glioxal, solución patrón 374
Glioxal-bis-(2-hidroxianilo) 374
D(+)Glucosa anhidra 370
D(+)Glucosa monohidrato 371
Glutaraldehído solución 50% 372
Goma arábica 377
Grasa de lana anhidra 453
Griess reactivo solución en acido acético 376
Griess reactivo A 376
Griess reactivo B 376

H

Hafnio, solución patrón 378
Hanus reactivo 379
Harina fósil compuesta 449
Hematoxilina 379
Hematoxilina solución según Mayer 379
n-Heptano 99% 380
n-Heptano 382
Heptano mezcla de isómero 382
Hexadecano 383
Hexafluoro-2-Propanol 384
Hexametildisilazano 384
Hexametilentetramina 384
n-Hexano 99% 385
n-Hexano 386
Hexano mezcla de isómeros 388
Hidrazina diclorhidrato 392
Hidrindantina 392
Hidrógeno peróxido solución 40% p/v 413
Hidrógeno peróxido solución 35% 413
Hidrógeno peróxido solución 30-32% 413
Hidrógeno peróxido solución 30% 414
Hidrógeno peróxido solución 6% 415
Hidrógeno peróxido solución 3.5% p/v 415
Hidrógeno peróxido solución 3% 416
Hidroquinona 416
p-Hidroxibenzaldehído 416
Hidroxilamina solución, alcohólico 417
Hidroxilamonio sulfato 417
Hierro, solución patrón 428
Hierro, polvo 427
Hierro, reducido para el hidrógeno 428
Hierro (II) amonio sulfato hexahidratado 429
Hierro (II) cloruro tetrahidratado 430
Hierro (II) sulfato heptahidratado 430
Hierro (II) sulfato 0.1mol/l 431
Hierro (III) y amonio citrato rojo 431
Hierro (III) y amonio citrato verde 431
Hierro (III) y amonio oxalato 432
Hierro (III) amonio sulfato dodecahidrato 432
Hierro (III) amonio sulfato solución 33% en acido nítrico 432
Hierro (III) amonio sulfato 0.1 mol/l 433
Hierro (III) amonio sulfato solución 100 g/l 433
Hierro (III) citrato 434
Hierro (III) cloruro anhidro sublimado 433
Hierro (III) cloruro hexahidratado 433
Hierro (III) cloruro solución 4.5% 434
Hierro Cloruro en solución 434
Hierro (III) nitrato nonahidratado 435
Hierro (III) óxido 435
Hierro (III) sulfato 435
Histamina diclorhidrato 389
L-Histidina monoclorhidrato monohidrato 390

DASITGROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

Histolemon 390
Holmio, solución patrón 390
Holmio oxido 391
Holmio perclorato solución 391
Hyamina 1622 391
Hyamina 1622 solución 0.004M 392

I

Idrimeter Erba Indicador C 419
Idrimeter Erba Solución A 419
Idrimeter Erba Solución B 419
Imidazol 419
Indicador azul de bromofenol 214
Indicador azul de timol 758
Indicador para iodometría 420
Indicador para nitrógeno amoniacal solución 420
Indicador universal pH 1-11 solución hidroalcohólica 421
Indicador universal pH 1-11 solución en agua 422
Indicador universal pH 0-5 solución hidroalcohólica 421
Indigo Carmin seco 422
Indigo carmín solución 422
Indio, solución patrón 422
Indol 423
Inositol 423
Inulina 424
Iridio, solución patrón 427
Isoamilo Acetato 436
Isobutanol 437
Isobutilo acetato 438
Isohexano 438
Isooctano 439
Isopar G 440
Isopentano 441
Isopropilamina 442
Isopropilo acetato 441
Isopropilo Miristato 442
Iterbio, solución patrón 794
Itrio, solución patrón 794

J

Jabón solución en alcohol etilico 660

K

Karl Fischer agua estándar 0.10 mg/g 449
Karl Fischer agua estándar 1.0 mg/g 448
Karl Fischer agua estándar 5.0 mg/ml 448
Karl Fischer agua estándar 10.0 mg/g 448
Karl Fischer anolito solución - sin CFC 446
Karl Fischer anolito solución para aldehídos y cetonas -libre de metanol 446
Karl Fischer anolito solución para horno 446
Karl Fischer anolito solución para los aceites 447
Karl Fischer catolito solución para aldehídos y cetonas 447
Karl Fischer catolito solución 447
Karl Fischer catolito solución 446
Karl Fischer reactivo 1 compuesto 2 mgH2O/ml 443
Karl Fischer reactivo 1 compuesto 5 mgH2O/ml 443
Karl Fischer reactivo de valoración 1 compuesto 5 mg H2O/ml para aldehídos y cet 443
Karl Fischer reactivo de valoración 2 compuestos 2 mg H2O/ml 444
Karl Fischer reactivo de valoración 2 compuestos 5 mg H2O/ml 444
Karl Fischer reactivo para la celda con y sin diafragma 447
Karl Fischer solvente 2 compuestos con tampón 445
Karl Fischer solvente 1 compuesto para aldehídos y cetonas 444
Karl Fischer solvente 2 componentes 445
Karl Fischer solvente 2 compuestos para aldehídos y cetonas - sin metanol 445
Karl Fischer solvente para los aceites 1 compuesto 444
Karl Fischer solvente para los aceites 2 compuestos 445
Kit Gram-Hucker 376
Kjeldahl antiespumante 449
Kovac reactivo 451

L

Lactofenol de Amman solución 160
Lactosa monohidrato 452
Lana de vidrio 370
Lantano, solución patrón 453
Lantano cloruro solución 25 g/l 454
Lantano nitrato hexahidratado 454
Lantano nitrato 0.1 mol/l 454
Lantano nitrato solución 50 g/l 454
Lantano oxido 455
Lecitina de soja 461
Levdura de cerveza seco 794
L(+)-Leucina 461
Líquido fijador AFA 354
Líquido fijador de Bouin 355
Líquido fijador de Bouin Holland 355
Líquido fijador de Davidson 355
Líquido fijador FIXALL-HIS 356
Liquido fijador sin acido acético 356
L(+)Lisina monoclorhidrato 467
Litio, solución patrón 462
Litio carbonato 463
Litio cloruro 463
Litio hidróxido anhidro 464
Litio hidróxido monohidrato 464
Litio hidruro 464
Litio metilato 0.1 mol/l (0.1N) 464
Litio nitrato 465
Litio sulfato monohidrato 465
Litio tetraborato anhidro 465
Lugol reactivo solución concentrada 466
Lugol reactivo solución yodo-yodurata 466
Lugol solución para kit Gram - Hucker 466
Luminol 467
Lutecio, solución patrón 467

M

Magnesio, solución patrón 468
Magnesio, polvo 468
Magnesio, tiras 468
Magnesio, virutas 468
Magnesio acetato tetrahidratado 470
Magnesio carbonato básico 470
Magnesio cloruro hexahidrato 471
Magnesio estearato 474
Magnesio fosfato dibásico trihidrato 472
Magnesio glicerofosfato 471
Magnesio hidróxido 472
Magnesio nitrato hexahidratado 472
Magnesio nitrato solución 10 g/l 473
Magnesio óxido 473
Magnesio óxido pesado 473
Magnesio perclorato 474
Magnesio peróxido 474
Magnesio sulfato anhidro 475
Magnesio sulfato heptahidratado 475
Maltosa monohidrato 478
Manganeso, solución patrón 479
Manganeso electrolítico 479
Manganeso (II) acetato tetraidrato 480
Manganeso (II) cloruro tetrahidratado 480
Manganeso (II) sulfato monohidrato 481
Manganeso (IV) dióxido 481
D-Manitol 481
May Grünwald reactivo 482
Mayer reactivo solución en agua 482
L-Mentol 483
2-Mercaptoetanol 483
Mercurio, solución patrón 484
Mercurio bromuro papel 484
Mercure (I) chlorure 485
Mercurio (II) dicloruro 485
Mercurio (II) dicloruro solución 54 g/l 485
Mercurio (II) óxido rojo 486
Mercurio (II) sulfato 486
Mercurio (II) sulfato solución 486
Mercurio (II) yoduro 486
Metanol 488
Metanol + 0.1% v/v acido formico 492
Metanol + 0.1% v/v acido trifluoroacético 493
Metanol-d1 494
Metanol-d3 494
Metanol-d4 493
Metanol-d4 + 0.03% TMS 494
Metanol, clorhídrico 493
4-Metil aminofenol sulfato 496
3-Metil-2-benzotiazolinona hidrazona clorhidrato 497
Metilciclohexano 497
Metiletilcetona 350
Metilo acetato 495
Metilo benzoato 496
Metilo isoamilcetona 500
Metilo isobutilcetona 500
Metilo p-oxibenzoato 499
Metilo salicilato 505
Metilo yoduro 500
N-Metil-2-pirrolidona 502
2-Metiltetrahidrofurano 505
4-Metiltetrahidropyrano 506
DL-Metionina 494
2-Metoxietanol 495
Mezcla alcohol-éter 30/70 p/p 507
Mezcla C.H.M. 507
Mezcla estándar para análisis de hidrocarburos 546
Mezcla etanol 95°/ 2-propanol 508
Mezcla etanol 99°/2-propanol 508
Mezcla Etanol absoluto/Eter Etilico 50/50 (v/v) 508
Mezcla Etanol absoluto/Eter Etilico 50/50 (p/p) 507
Mezcla octano 80 537
Mezclas para disolventes residuales analisis 509
Mezcla para el control de soldaduras de envases 509
Mezcla para la determinación de índice de bromo 508
Mezcla selenica 646
Molibdeno, solución patrón 510
Morfolina 512
Mucicarmin solución hidroalcohólica 512
Murexida 516

N

Naftaleno 517
N-(1-Naftil)etilendiamina diclorhidrato 518
1-Naftol 517
2-Naftol 517
alpha-Naftolbenceína 518
alpha-Naftolbenceína solución 0.2% en acido acético 518
Naranja II 539
Naranja de metilo 501
Naranja de metilo solución mixta 502
Naranja de metilo solución 0.1% 501
Naranja de xilenol 793
Naranja G 539
Negro Amido 10B solucion 157

ÍNDICE ALFABÉTICO

Negro de eriocromo T 329
Negro Sudan B 726
Neocuproina clorhidrato 519
Neodimio, solución patrón 519
Nessler reactivo solución única 520
Nessler reactivo solución A 520
Nessler reactivo solución B 521
Nicotinamida 525
Nigrosina 525
Ninhidrina 526
Ninhidrina solución 526
Ninhidrina y estaño (II) cloruro reactivo 526
Niobio, solución patrón 526
Niquel, solución patrón 521
Níquel, polvo 521
Níquel (II) acetato 522
Níquel (II) amonio sulfato hexahidrato 523
Níquel (II) carbonato básico 523
Níquel (II) cloruro hexahidratado 523
Níquel (II) nitrato 10 gIL 524
Níquel (II) nitrato hexahidratado 524
Níquel (II) sulfato hexahidratado 524
Nital solución 4% 527
Nitrato, solución patrón 527
Nitrito, solución patrón 534
m-Nitrobenzaldehyde 534
Nitrobenceno 535
p-Nitrofenol 535
p-Nitrofenol solución 0.1% en agua 535
Nonilfenol etoxilato 10 ETO 536

O

n-Octano 537
Octanol-1 538
Orceina 540
Orcina monohidrato 540
Oro, solución patrón 375
Oro tricloruro acido 375
Osmio, solución patrón 549
8-Oxiquinoleína 417

P

Paladio, solución patrón 553
Paladio (II) cloruro 553
Paladio nitrato solución 2 g/l 554
Papanicolaou Hematoxilina en solución según Harris 554
Papanicolaou solución EA 50 554
Papanicolaou solución OG 6 555
Papel de plata de manganeso 655
Papel indicador de pH 420
Papel tornasol 466
Parafina 56-58°C Erbaplast (sin DMSO) 555
Paraformaldehído 556
Pararrosanilina solución, decolorado 556
Pardo de Bismarck R 201
Patrones de color ASTM 187
Patrones de color Gardner 368
Patrones de color Saybolt 644
Patrones de densidad 292
Patrones de índice de refracción 637
Patrón multielementos para cromatografía iónica 512
Patrón multielementos para ICP 513
Patrón multielementos para ICP e ICP-MS 515
Patrones orgánicos PCB 540
Patrones orgánicos: mezcla multi-HAP 542
Patrones orgánicos: mezcla multi-PCB 543
Patrones orgánicos: mezcla de pesticidas 546
Patrones de osmolalidad 550
Patrones de osmolalidad a base de proteínas 550
Patrones de osmolalidad a base de orina 550
Patrones de puntos de fusión 483
Patrones de viscosidad 785
n-Pentano 99% 557
n-Pentano 558
Pepsina HCl 560
Petróleo 563
Piedra pómez 630
Piperidina 580
Piridina 630
Piridina-d5 631
1-(2-Piridilazo)-2-naftol 632
Pirocatequina 632
Pirogalol 632
Plata, solución patrón 653
Plata, hojas 652
Plata, lana 653
Plata acetato 654
Plata carbonato 654
Plata cloruro 654
Plata dietilditiocarbamato 655
Plata nitrato 655
Plata nitrato 1 mol/l (1N) 656
Plata nitrato 0.5 mol/l (0.5N) 656
Plata nitrato 0.1 mol/l (0.1N) 657
Plata nitrato 0.1 mol/l (0.1N) en propanol-2 657
Plata nitrato 0.05 mol/l (0.05N) 657
Plata nitrato 0.01 mol/l (0.01N) 658
Plata nitrato 0.01 mol/l (0.01N) en propanol-2 658
Plata nitrato solución 5% 656
Plata nitrato solución 2.90756% 656
Plata nitrato solución 658
Plata óxido 659
Plata sulfato 659
Plata sulfato solución 0.7% en acido sulfúrico 659
Platino, solución patrón 580
Plomo, solución patrón 455
Plomo 455
Plomo (II) acetato algodón 457
Plomo (II) acetato básico 457
Plomo (II) acetato solución 95 g/l 458
Plomo (II) acetato básico solución 457
Plomo (II) acetato 3-hidratado 458
Plomo (II) acetato papel 457
Plomo (II) nitrato 458
Plomo (II) nitrato 0.1 mol/l 459
Plomo (II) nitrato 0.05 mol/l 459
Plomo (II) nitrato solución 33 g/l 459
Plomo (II) óxido 460
Plomo (II) sulfato 460
Plomo (IV) óxido 460
Polivinilpirrolidona 581
Potasio, solución patrón 582
Potasio acetato 583
Potasio bicarbonato 584
Potasio bisulfato 584
Potasio bromato 584
Potasio bromato 0.033 mol/l (0.198N) 585
Potasio bromato 0.02 mol/l (0.12N) 585
Potasio bromato 0.0167 mol/l (0.1N) 585
Potasio bromuro 585
Potasio carbonato 586
Potasio citrato tribásico monohidrato 591
Potasio cloruro 586
Potasio cloruro 3.5 mol/l (3.5N) 587
Potasio cloruro 3.5 mol/l (3.5N) con plata cloruro 587
Potasio cloruro 3 mol/l (3N) 588
Potasio cloruro 3 mol/l (3N) solución agua-glicerina 588
Potasio cloruro 3 mol/l (3N) con plata cloruro 588
Potasio cloruro 1 mol/l (1N) 588
Potasio cloruro 0.1 mol/l (0.1N) 589
Potasio cloruro 0.01 mol/l (0.01N) 589
Potasio cloruro solución 25g/l en HCl 589
Potasio cloruro solución 12g/l 589
Potasio cromato 590
Potasio cromato solución 10% 590
Potasio cromato solución 5% 591
Potasio cloruro solución saturada 590
Potasio cloruro solución 590
Potasio dicromato 592
Potasio dicromato 0.167 mol/l (1 N) 592
Potasio dicromato 0.0167 mol/l (0.1 N) 593
Potasio dicromato 0.0417 mol/l (0.25N) 593
Potasio dicromato 0.04 mol/l (0.24N) en 80 g/l HgSO4 593
Potasio dicromato solución 106 g/l 594
Potassio dicromato 0.1414 g/l 594
Potasio dicromato solución 0.5% 592
Potasio dicromato - Solución de ácido sulfúrico 592
Potasio ferricianuro 594
Potasio ferrocianuro trihidrato 595
Potasio ferrocianuro solución 10% 595
Potasio ferrocianuro solución 53 g/l 596
Potasio fluoruro 596
Potasio fluoruro dihidrato 596
Potasio fosfato dibásico anhidro 614
Potasio fosfato dibasico trihidrato 614
Potasio fosfato monobásico 615
Potasio fosfato monobásico 0.2 mol/l (0.2N) 615
Potasio ftalato acido 597
Potasio ftalato acido 0.2 mol/l (0.2N) 598
Potasio ftalato acido 0.1 mol/l (0.1N) 598
Di-Potasio hexacloroplatinato 597
Potasio hidróxido, gotas 599
Potasio hidróxido, escamas 599
Potasio hidróxido solución 45% 600
Potasio hidróxido solución 38% (40° Bé) en agua 601
Potasio hidróxido solución 33% 601
Potasio hidróxido solución 28% 601
Potasio hidróxido 2 mol/l (2N) en etanol 602
Potasio hidróxido 1 mol/l (1N) 602
Potasio hidróxido 0.5 mol/l (0.5N) 602
Potassium hydroxide 0.5 mol/l (0.5N) in ethanol 603
Potasio hidróxido 0.5 mol/l (0.5N) en metanol 603
Potasio hidróxido 0.46 mol/l (0.46N) 604
Potasio hidróxido 0.25 mol/l (0.25N) 604
Potasio hidróxido 0.23 mol/l (0.23N) 604
Potasio hidróxido 0.1 mol/l (0.1N) 604
Potasio hidróxido 0.1 mol/l (0.1N) en etanol 605
Potasio hidróxido 0.1 mol/l (0.1N) en metanol 605
Potasio hidróxido solución 3% en etanol 601
Potasio hidróxido en solución 606
Potasio metabisulfito 609
Potasio nitrato 609
Potasio nitrato 1 mol/l (1N) 610
Potasio nitrito 610
Potasio oxalato monohidrato 610
Potasio permanganato 611
Potasio permanganato 0.2 mol/l (1N) 612
Potasio permanganato 0.02 mol/l (0.1N) 612
Potasio permanganato 0.002 mol/l (0.01N) 613
Potasio permanganato solución 3% 612
Potasio permanganato y ácido fosfórico solución 613
Potasio persulfato 613
Potasio (meta)-peryodato 611
Potasio piroantimoniato acido 615
Potasio piroantimoniato solución 616
Potasio pirofosfato 616
Potasio pirosulfato 616
Potasio y sodio tartrato tetrahidrato 617
Potasio sorbato 617
Potasio sulfato 618
Potasio sulfocianuro 620
Potasio sulfocianuro 0.1 mol/l (0.1N) 620
Potasio sulfocianuro solución 5% 620
Potasio sulfocianuro solución 621

Product specifications are subject to changes. Please visit our website for updates.

Potasio sulfoguayacolato 596
Potasio tartrato 618
Potasio tartrato acido 619
Potasio telurito 619
Potasio tetraiodomercurato solución alcalina 619
Potasio yodato 606
Potasio yodato 0.05 mol/l (0.3N) 606
Potasio yodato 0.0167 mol/l (0.1N) 606
Potasio yodato 0.00167 mol/l (0.01N) 607
Potasio yodato acido 597
Potasio yodobismutato solución 609
Potasio yoduro 607
Potasio yoduro solución 10% 608
Potasio yoduro solución 3.9% 608
Potasio yoduro solución 608
Praseodimio, solución patrón 621
L(-)Prolina 622
Propan-1-ol 622
Propan-2-ol 623
Propan-2-ol 70% 626
1,3-Propanodiol 627
n-Propil acetato 628
Propilenglicol 629
Propileno carbonato 629
Propil p-hidroxibenzoato 629
Púrpura de bromocresol 212
Púrpura de bromocresol solución 0.4% en alcohol etilico 212
Púrpura de bromocresol solución 213
Púrpura de m-cresol 282

Q

Queroseno 449
Quinidina sulfato 634
Quinoleína 634

R

Rafinosa 635
L(+)Ramnosa 638
Reactivo de ácido aminohipúrico 158
Reactivo de Benedict 197
Reactivo molibdenovanádico 511
Reactivo sulfomolibdico 729
Reactivo TAN 635
Reactivo TBN ASTM D2896 636
Reactivo TBN ASTM D4739 636
Reactivo unico para lipolisi 635
Renio, solución patrón 638
Resorcina 637
Rhodol solución 639
Riboflavina 639
D(-)Ribosa 639
Rodio, solución patrón 638
Rojo Alizarina 149
Rojo Congo 273
Rojo Congo papel 273
Rojo Congo solución 274
Rojo de o-cresol 283
Rojo de o-cresol solución 0.2% en alcohol etilico 283
Rojo de o-cresol solución en sodio hidróxido 0.1N / alcohol etilico 96°/agua 283
Rojo de fenol 570
Rojo de fenol solución 0.2% en alcohol etilico 571
Rojo de fenol solución 571
Rojo de metilo 503
Rojo de metilo solución 0.1% en alcohol etilico 504
Rojo de metilo solución 0.2% en alcohol etilico 503
Rojo de metilo solución 504
Rojo de metilo solución mixta 504
Rojo de quinaldina 634
Rojo neutro 521
Rojo nuclear sólido 536
Rojo Ponceau BS 581
Rojo Ponceau S 581
Rojo para aceites O 636
Rubidio, solución patrón 640
Rutenio, solución patrón 641

S

D(+)-Sacarosa 726
Safranina T 642
Safranina T solución hidroalcohólica para kit Gram-Hucker 642
Sal de Reinecke 637
Salicilaldehido azina 642
Samario, solución patrón 643
Saponina 644
Schiff reactivo solución para Aldehídos 645
Schiff reactivo solución para PAS coloratión 646
Selenio, solución patrón 647
Selenio, polvo 646
Selenio dióxido 648
Silicio, solución patrón 651
Silicio dióxido 652
Sodio solución patrón 661
Sodio acetato anhidro 662
Sodio acetato trihidrato 662
Sodio alginato 663
Sodio aluminato 663
Sodio L-ascorbato 664
Sodio arsenito 0.1 mol/l (0.2N) 664
Sodio arsenito 0.05 mol/l (0.1N) 664
Sodio azida 665
Sodio benzoato 665
Sodio bicarbonato 665
Sodio bisulfato monohidrato 666
Sodio boro hidruro 666
Sodio bromuro 667
Sodio carbonato anhidro 667
Sodio carbonato solución 670
Sodio carbonato monohidrato 668
Sodio carbonato decahidrato 668
Sodio carbonato 0.5 mol/l (1N) 669
Sodio carbonato 0.05 mol/l (0.1N) 669
Sodio carbonato solución 20% 670
Sodio cianoborohidruro 674
Sodio citrato dibásico 672
Sodio citrato tribásico anhidro 672
Sodio citrato tribásico dihidrato 673
Sodio cloruro 670
Sodio cloruro 5 mol/l (5N) 672
Sodio cloruro 0.1 mol/l (0,1N) 672
Sodio cobalto nitrito 673
Sodio dicloroisocianurato dihidrato 674
Sodio dietilditiocarbamato trihidrato 674
Sodio fluoruro 675
Sodio formiato 675
Sodio fosfato dibásico anhidro 700
Sodio fosfato dibásico dihidrato 701
Sodio fosfato dibasico dodecahidrato 702
Sodio fosfato monobásico monohidrato 703
Sodio fosfato monobásico dihidrato 703
Sodio fosfato tribásico dodecahidratado 704
Sodio o-Fosfito pentahidrato 704
Sodio glicerofosfato pentahidrato 676
Sodio gluconato 676
Sodio glutamato acido 676
Sodio hexafluorosilicato 677
Sodio hexametafosfato 677
Sodio hidrosulfito 678
Sodio hidróxido, perlas 678
Sodio hidróxido, lentejas 679
Sodio hidróxido sobre sílice 680
Sodio hidróxido solución 50% 680
Sodio hidróxido solución 40% 680
Sodio hidróxido solución 35-37% 681
Sodio hidróxido solución 35% 681
Sodio hidróxido solución 32% 681
Sodio hidróxido solución 30% 682
Sodio hidróxido solución 20% p/p 683
Sodio hidróxido solución 20% p/v 683
Sodio hidróxido solución 10% p/v 683
Sodio hidróxido solución 5% w/v 684
Sodio hidróxido 6 mol/l (6N) 684
Sodio hidróxido 5 mol/l (5N) 685
Sodio hidróxido 4 mol/l (4N) 685
Sodio hidróxido 3 mol/l (3N) 685
Sodio hidróxido 2 mol/l (2N) 686
Sodio hidróxido 1.2 mol/l (1.2N) 686
Sodio hidróxido 1 mol/l (1N) 686
Sodio hidróxido 0.7 mol/l (N/1.4) 687
Sodio hidróxido 0.5 mol/l (0.5N) 688
Sodio hidróxido 0.357 mol/l (0.357N) 688
Sodio hidróxido 0.25 mol/l (0.25N) 688
Sodio hidróxido 0.2 mol/l (0.2N) 689
Sodio hidróxido 1/9 mol/l (N/9) 689
Sodio hidróxido 0.1 mol/l (0.1N) 689
Sodio hidróxido 0.1 mol/l (N/10) en etanol 690
Sodio hidróxido 0.01 mol/l (0.01N) 690
Sodio hidróxido solución 691
Sodio hidróxido solución, metanólico 691
Sodio hidróxido-d 30% 691
Sodio hidróxido-d 1 mol/l 691
Sodio hipoclorito solución en agua 692
Sodio hipoclorito solución 12.5% 692
Sodio hipofosfito 693
Sodio laurilsulfato 693
Sodio metabisulfito 694
Sodio metafosfato 695
Sodio metaperyodato 694
Sodio metóxido 0.1 mol/l 695
Sodio molibdato dihidrato 695
Sodio nitrato 696
Sodio nitrito 696
Sodio nitrito 0.1 mol/l (0.1N) 697
Sodio nitrito solución 500 g/L 697
Sodio nitroprusiato dihidrato 697
Sodio oxalato 698
Sodio perborato tetrahidratado 698
Sodio perclorato 699
Sodio peróxido 699
Sodio persulfato 699
Sodio persulfato 1 mol/l 700
Sodio pirofosfato decahidrato 705
Sodio salicilato 705
Sodio silicato 705
Sodio succinato hexahidratado 706
Sodio sulfato anhidro 706
Sodio sulfato decahidrato 707
Sodio sulfito anhidro 708
Sodio sulfocianuro 709
Sodio sulfuro nonahidrato 707
Sodio sulfuro nonahidrato solución 707
Sodio sulfuro trihidrato 708
Sodio tartrato acido 678
Sodio tartrato dihidrato 709
Sodio tetraborato anhidro 709
Sodio tetraborato decahidrato 710
Sodio tetrafenilborato 710
Sodio tiosulfato anhidro 711
Sodio tiosulfato pentahidrato 711
Sodio tiosulfato 1 mol/l (1N) 712
Sodio tiosulfato 0.5 mol/l (0.5N) 712
Sodio tiosulfato 0.2 mol/l (0.2N) 712
Sodio tiosulfato 0.1 mol/l (0.1N) 712
Sodio tiosulfato 0.0394 mol/l (0.0394N) 713
Sodio tiosulfato 0.0197 mol/l (0.0197N) 713
Sodio tiosulfato 0.01 mol/l (0.01N) 713
Sodio tungstato dihidrato 714
Sodio y amonio hidrógeno fosfato 664
Sodio yoduro 693

ÍNDICE ALFABÉTICO

Solfato, solución patrón 728
Solución ADF 147
Soluciónes estándar para el grado de coloración de líquidos 291
Solución NDF 519
Solución NDF PLUS 519
Solución patrón de 500000 µS/cm 721
Solución patrón de 450000 µS/cm 720
Solución patrón de 350000 µS/cm 720
Solución patrón de 300000 µS/cm 720
Solución patrón de 200000 µS/cm 720
Solución patrón de 150000 µS/cm 719
Solución patrón de 100000 µS/cm 719
Solución patrón de 50000 µS/cm 719
Solución patrón de 20000 µS/cm 719
Solución patrón de 12880 µS/cm 718
Solución patrón de 10000 µS/cm 718
Solución patrón de 5000 µS/cm 718
Solución patrón de 1413 µS/cm 718
Solución patrón de 1000 µS/cm 718
Solución patrón de 500 µS/cm 717
Solución patrón de 200 µS/cm 717
Solución patrón de 147 µS/cm 717
Solución patrón de 100 µS/cm 717
Solución patrón de 84 µS/cm 717
Solución patrón de 50 µS/cm 716
Solución patrón de 20 µS/cm 716
Solución patrón de 10 µS/cm 716
Solución patrón de 5 µS/cm 716
Solución patrón de 1.30 µS/cm 715
Soluciónes primarias para el grado de coloración de líquidos 290
Solución Redox 220 mV a 20°C 636
Solución Redox 468 mV a 25°C 637
D-Sorbitol 715
Sorbitol 714
Sorbitol (no cristalizable) solución 70% 715
Sudan III 726
Sudan III solución hidroalcohólica 727
Sulfanilamida 728
Sulfito, solución patrón 729
Sulfolano 729
Suspensión opalescente primaria 621

T

TA indicador azul de timol 758
TAC indicador azul de bromofenol 214
Talco 744
Talio, solución patrón 754
Tamices moleculares 3 A 510
Tamices moleculares 4 A 510
Tampón pH 1 216
Tampón pH 1.68 216
Tampón pH 2 217
Tampón pH 3 217
Tampón pH 3.5 217
Tampón pH 3.56 218
Tampón pH 3.7 218
Tampón pH 4 218
Tampón pH 4.5 219
Tampón pH 4.62 219
Tampón pH 5 219
Tampón pH 5.2 219
Tampón pH 6 220
Tampón pH 6.8 220
Tampón pH 6.88 220
Tampón pH 7 220
Tampón pH 7.2 221
Tampón pH 7.20 según Weise 221
Tampón pH 7.4 222
Tampón pH 8 222
Tampón pH 9 222
Tampón pH 9.22 223
Tampón pH 10 223
Tampón pH 10.06 224
Tampón pH 11 224
Tampón pH 12 225
Tampón pH 13 225
Tampón acetato pH 6.0 130
Tampón acetato pH 4.6 130
Tampón acetato pH 4,5 216
Tampón bórico solución 207
Tampón cloruro di ammonio pH 9.5 170
Tampón cloruro di ammonio pH 10.0 170
Tampón cloruro di ammonio pH 10.7 169
Tampón de amoníaco pH 10 165
Tampon fosfato pH 9.0 574
Tampon fosfato pH 7.4 575
Tampón fosfato pH 7.2 575
Tampon fosfato pH 6.8 575
Tampon fosfato pH 6.0 575
Tampon fosfato pH 5.5 575
Tampon fosfato pH 3.0 576
Tampon fosfato pH 2.0 576
Tampón fuerza iónica total 769
Tampón tris(hidroximetil)aminometano sodio cloruro pH 7.4 778
Tampón tris(hidroximetil)aminometano - EDTA pH 7.4 777
Tántalo, solución patrón 744
Tartracina 746
Tauber reactivo 746
Teluro polvo 747
Teluro trozos 746
Teluro, solución patrón 747
Terbio, solución patrón 748
Tetrabutilamonio bisulfato 748
Tetrabutilamonio bromuro 749
Tetrabutilamonio hidróxido solución 0.1 mol/l (0.1N) 749
Tetrabutilamonio hidróxido 0.1 mol/l (0.1N) en methanol / propanol-2 (50/50) 750
Tetrabutilamonio hidróxido 0.1 mol/l (0.1N) en 2-propanol 749
Tetracloroetano-D2 751
Tetracloroetileno 750
Tetraetilamonio bromuro 751
Tetrahidrofurano 751
Tetrahidrofurano-d8 753
Tetrametilamonio hidróxido 10% 753
N,N,N',N'-Tetrametil -p-fenilendiamina diclorhidrato 753
Timol 757
Timolftaleína 758
Timolftaleína 0.1% solución hidroalcohólica 759
Timolftaleína solución 0.1% en etanol 759
Tioacetamida 754
Tioacetamida solución 40 g/l 755
Tiourea 756
L(-)Tirosina 781
Tisab 762
Tisab III solución 763
Titanio, solución patrón 763
Titanio dióxido 764
Titanio tetrapropilato 764
Titanio tricloruro-ácido sulfúrico 765
o-Tolidina solución 0.1% 765
o-Tolidina solución 765
Tolueno 765
Tolueno-d8 768
Tolueno en solución de hexano 768
p-Toluenosulfonamida 768
Torio, solución patrón 756
Tornasol 465
Triacetina 769
1,2,4-Triclorobenceno 771
2,2,2 - Tricloroetanol 770
Tricresol 771
Trietanolamina 771
Trietilamina 772
Trietilenglicol 773
Trifenilfosfina 776
2,3,5-Trifeniltetrazolio cloruro 776
Trimetilcetilamonio bromuro 775
Trimetilcetilamonio p-toluensulfonato 775
Trioctilfosfina óxido 775
Tris (hidroximetil)aminometano 776
Tris (hidroximetil)aminometano hidrocloruro 778
Tris (hidroximetil)aminometano solución 777
Tris (hidroximetil)aminometano solución tampón pH 8.1 777
Triton© X100 solución 778
Tropeolina O 778
Tropeolina O solución 0.1% 779
Tungsteno, solución patrón 779
Tulio, solución patrón 756
Turbiedad std 4000NTU formazina 780
Türk reactivo 780

U

Uranio, solución patrón 782
Urea 782

V

Vanadio, solución patrón 783
Vanadio (V) pentóxido 784
Vanillina 784
Vanillina solución fosfórica 784
Vaselina blanca 556
Verde brillante 209
Verde de bromocresol 211
Verde de bromocresol 0.04% solución hidroalcohólica 211
Verde de bromocresol solución 211
Verde de bromocresol - Rojo de metilo solución 212
Verde claro SF 461
Verde de malaquita 476
Verde de malaquita solución 0.5% en acido acético anhidro 476
Verde de metilo 499
Verde sólido FCF 353
Violeta cristal 284
Violeta cristal solución 0.5% en acido acético anhidro 284
Violeta cristal oxalato solución para kit Gram-Hucker 284
Violeta de genciana 368
Violeta de genciana solución fenicada 368
Violeta de pirocatequina 632

W

Wijs reactivo 790

X

o-Xileno 792
Xileno, mezcla de isómeros 791
Xilenocianol FF 792
2,4-Xilenol 792
D(+)Xilosa 793

Y

Yodo 424
Yodo bisublimado 425
Yodo 0.5 mol/l (1N) 425
Yodo 0.05 mol/l (0.1N) 425
Yodo 0.01 mol/l (0.02N) 426
Yodo 0.005 mol/l (0.01N) 426
Yodo 10 ppm 426
Yodo bromuro solución 426
Yodoformo 427
Yodoplatinato reactivo 427

DASITGROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

Yoduro, solución patrón 424

Z

Ziehl -Neelsen reactivo solución 796
Zinc, granulada 796
Zinc, hojas 796
Zinc, polvo 797
Zinc, activado 797
Zinc, solución patrón 797
Zinc acetato dihidrato 798
Zinc carbonato básico 799
Zinc cloruro anhidro 799
Zinc cloruro solución 60% 799
Zinc cloruro solución, yodada 800
Zinc cloruro - ácido fórmico solución 800
Zinc dibencilditiocarbamato 800
Zinc estearato 801
Zinc nitrato hexahidratado 800
Zinc óxido 801
Zinc sulfato heptahidratado 802
Zinc sulfato monohidrato 802
Zinc sulfato 0.1 mol/l (0.2N) 802
Zinc sulfato 0.05 mol/l (0.05N) 803
Zinc sulfuro 803
Zirconio polvo 803
Zirconio, solución patrón 804
Zirconio (IV) óxido 804
Zirconio oxinitrato 805

BEST CHOICE

Labware Consumable & Instruments

LLG Lab Logistics Group



- Glassware
- Plasticware
- Liquid Handling
- Chromatography
- Analitycal Instruments

www.carloerbareagents.com

DASITGROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

ALPHABETISCHER INDEX

A

Absorbens für verschüttete Flüssigkeiten 129
Acetal 129
Acetamid 129
Acetanilid 129
Aceton 138
Aceton-d6 141
Aceton / Wasser 98/2 (v/v) mit Bromphenolblau 0.02 g/l 141
Acetonitril 142
Acetonitril + 0,1% v/v Ameisensäure 145
Acetonitril + 0.1% v/v trifluoressigsäure 145
Acetonitril-d3 145
Acetophenon 146
Acetylaceton 146
n-Acetyl-L-cystein 146
Acridinorange 147
ADF Lösung 147
Adipinsäure 147
ß-Alanin 148
Albumin aus Eipulver 148
Alcianblau 8GX 149
Alcianblau 8GS 1% 148
Alizarin 149
Alizarin gesättigte Lösung in Ethanol 149
Alizaringelb R 150
Alizarinrot S 149
Alkaliblau 6B 150
Alkaliblaue 6B 2% ige Lösung in Ethanol 150
Alkohol Bouin Holland 355
Allylthioharnstoff 150
Aluminium, pulver 151
Aluminium-Standardlösung 151
Aluminiumammoniumsulfat-Dodecahydrat 152
Aluminiumchlorid. wasserfrei 153
Aluminiumchlorid-Hexahydrat 153
Aluminiumhydroxid 154
Aluminiumkaliumsulfat-Dodecahydrat 156
Aluminiumnitrat-Nonahydrat 154
Aluminiumoxid 155
Aluminiumoxid (basisch) 155
Aluminiumoxid (neutral) 156
Aluminiumoxid aktiviert 156
Aluminiumoxid (Säure) 155
Aluminiumsulfat 157
Ameisensäure 85% 363
Ameisensäure 99% 362
Ameisensäure -d 364
Amidoschwarz 10B losung 157
p-Aminobenzoesäure 158
1-Amino-2-naphthol-4-sulfonsäure 158
4-Aminophenazon 159
m-Aminophenol 159
p-Aminophenol 159
Aminopyrazol in Lösung 159
Ammoniaklösung 32% 160
Ammoniaklösung 30% 160
Ammoniaklösung 28% 161
Ammoniaklösung 25% 162
Ammoniaklösung 20 - 22 % 163
Ammoniaklösung 17% 164
Ammoniaklösung 10% 164
Ammoniaklösung 6N 164
Ammoniaklösung verdünnt 164
Ammoniaklösung von Tetramminecuivre 165
Ammoniakpufferlösung pH 10 165
Ammonium-Standardlösung 165
Ammoniumacetat 166
Ammoniumbicarbonat 167
Ammoniumbromid 167
Ammoniumcarbamat 168
Ammoniumcarbonat 168
Ammoniumcarbonatlösung 158 g/l 168
Ammoniumchlorid 169
tri-Ammoniumcitrat 170
Ammoniumcitrat zweibasisch 170
Ammoniumcitratlösung 20% 171
Ammoniumdihydrogenphosphat 25 mg/l 172
Ammoniumfluorid 171
Ammoniumformiat 171
Ammoniumhydrogendifluorid 172
Ammoniumiodid 173
Ammoniummolybdat tetrahydrat 173
Ammonium molybdat Losung 2.5 % 174
Ammoniumoxalat 4% 175
Ammoniumoxalatmonohydrat 175
Ammoniumpersulfat 175
Ammoniumphosphat dibasisch 176
Ammoniumphosphat einbasig 176
Ammoniumsulfamat 177
Ammoniumsulfat 177
Ammoniumsulfid in 20% iger Lösung 178
Ammoniumtartrat 178
Ammoniumthiocyanat 178
Ammoniumthiocyanat 0.01 mol/l (0.01N) 179
Ammoniumthiocyanat 0.1 mol/l (0.1N) 179
Ammoniumthiocyanat 1 mol/l (1N) 179
Ammonsalpeter 174
Ammonsalpeter 200 mg/l 174
n-Amylalkohol 180
Aniliniumchlorid 181
Anisaldehyd 181
Anisaldehydlösung 181
Anissäure 181
Anthron 182
Antimonstandardlösung 182
Antimon-Kalium-Tartrat 183
Antimontrichlorid 183
DL-Apfelsäure 477
L(+)Arginin 184
L-Argininmonohydrochlorid 184
Arsenazo III 184
Arsen-Standardlösung 185
Arsentrioxid-Lösung 185
L(+)-Ascorbinsäure 186
L(+)Asparagin 186
L(+)Asparaginsäure 186
ASTM-Farbstandards 187
Aus Mandeln raffiniertes Öl 539
Ausilab 110 187
Ausilab 140 187
Ausilab 210 187
Ausilab 250 188
Ausilab 260 188
Ausilab 280 188
Ausilab 290 188
Ausilab 400 189
Ausilab 500 189
Azomethin H 189
Azur II 189
Azur II Eosin 190

B

Barbitursäure 191
Barium-Standardlösung 191
Bariumacetat 192
Bariumcarbonat 192
Bariumchlorid 10% 193
Bariumchlorid 61 g/l 193
Bariumchlorid 30 g/l 194
Bariumchlorid 0.1 mol / l (0.2 N) 193
Bariumchromat 194
Bariumchlorid-Dihydrat 192
Bariumhydroxid-Octahydrat 194
Bariumhydroxidlösung 47.3 g / l 195
Bariumhydroxidlösung 5% 195
Bariumnitrat 195
Bariumperchlorat- trihydrat 195
Bariumperchlorat 0.05 mol/l 196
Bariumperchlorat 0.025 mol/l 196
Bariumsulfat 196
Benzaldehyd 197
Benzalkoniumchlorid 197
Benzin E 563
Benzoesäure 199
Benzol 198
Benzol-d6 198
Benzolsulfonylchlorid 198
Benzylalkohol 200
Benzylbenzoat 200
Bernsteinsäure 725
Bernsteinsäureanhydrid 725
Beryllium-Standardlösung 201
Bierhefe trocken 794
Bimsstein 630
Bismarck brown R 201
Biuret 97% 205
Biuret-Reagenz 205
Blaues Anilin. wasserlöslich 180
Blei-Standardlösung 455
Blei (II) acetat basisch 457
Blei (II) acetat basische Lösung 457
Blei (II) acetat Baumwolle 457
Blei (II) acetatlösung 95 g / l 458
Blei (II) acetatpapier 457
Blei (II) acetat Trihydrat 458
Blei (II) nitrat 458
Blei (II) nitrat 0.1 mol/l 459
Blei (II) nitrat 0.05 mol/l 459
Blei (II) nitratlösung 33 g/l 459
Blei (II) oxid 460
Blei (II) sulfat 460
Blei (IV) oxid 460
Bleibleche 455
Bor Standardlösung 208
Boric Pufferlösung 207
Borsäure 205
Borsäure 1% 207
Borsäure 3% 206
Borsäure 4% 206
Borsäure 4% mit Indikator 206
Borsäure 20 g/l 207
Borsäure 20 g/l mit Indikator 207
Bouins Pikroformol Wasserlisch 355
Brechungsindex-Standards 637
Brenzcatechinviolett 632
Brenzkatechin 632
Brillantes Kresylblau 208
Brillantgrün 209
Bromat-Standardlösung 210
Brome lösung 210
Bromidbromat 0.0167 mol/l 210
Bromid-Standardlösung 210
Bromkresolgrün 211
Bromkresolgrün - Methylrotlösung 212
Bromkresolgrün 0.04% ige hydroalkoholische Lösung 211
Bromkresolgrünlösung 211
Bromkresolpurpur 212
Bromkresol-Purpur-Lösung 213
Bromkresol-Purpur-Lösung 0.4% in Ethanol 212
alpha-Bromnaphthalin 213
Bromnatrium 667
Bromphenolblau 213
Bromphenolblau 0.02% 214
Bromphenolblau Lösung 0.4% in Ethanol 213
Bromphenolblau-Indikator 214
Bromphenolblau-Lösung 214
Bromphenolblau-TAC-Indikator 214
Bromthymolblau 215
Bromthymolblau 0.02% 215

ALPHABETISCHER INDEX

Bromthymolblau 0.4% in Ethanol 215
Bromthymolblau-Lösung 215
Bromwasser 211
Bromwasserstoffsäure 48% 393
Butanol-1 226
Butanol-2 227
tert-Butanol 227
Buttersäure 231
Butylhydroxytoluol 230
tert-Butylmethylether 230
p-Benzochinon 199
1-Butansulfonsäure-Natriumsalz 225
2-Butoxyethanol 228
2-(2-Butoxyethoxy)ethanol 228
n-Butylacetat 228
n-Butylchlorid 229

C

Cadmium-Standardlösung 232
Cadmiumacetat-Dihydrat 233
Cadmiumcarbonat 233
Cadmiumchlorid monohydrat 233
Cadmiumnitrat-Tetrahydrat 234
Caesium-Standardlösung 253
Caesiumchlorid 25 g/l Lösung 254
Calcium-Standardlösung 234
Calciumacetat monohydrat 236
Calciumacetat wasserfrei 236
Calciumcarbonat 236
Calciumchlorid wasserfrei 237
Calciumchloriddihydrat 238
Calciumchloridhexahydrat 239
Calciumchlorid 0.025% 239
Calciumfluorid 240
Calciumgluconat 240
Calciumhydroxid 240
Calciumlactat 241
Calcium nitrate tetrahydrat 241
Calciumoxid, Klumpen 242
Calciumoxid, Pulver 242
Calciumpantothenat 242
Calciumphosphat tribasisch 243
Calciumphosphat zweiwertiges Dihydrat 243
Calciumstearat 244
Calciumsulfat-Dihydrat 244
Calciumsulfat-Halbhydrat 245
Calciumsulfat-Halbhydratlösung 245
Calcon 245
Calconcarbonsäure 245
Calmagite 246
n-Capronsäure 247
n-Caprylsäure 247
Carborundum-Granulat 248
Carborundum-Pulver 248
Carrez Reagenz Kaliumsalz 249
Carrez Reagenz Zinksalz 249
Cäsiumchlorid 254
Cäsiumsulfat 254
Cer-Standardlösung 250
Cer(III)-nitrat-Hexahydrat 251
Cer (IV) ammoniumnitrat 251
Cer (IV) ammoniumnitrat 0.01 mol/l 251
Cer (IV) ammoniumnitrat 0.1 mol/l 251
Cer (IV) ammoniumsulfatdihydrat 252
Cer (IV) ammoniumsulfat 0.01 mol/l 252
Cer (IV) ammoniumsulfat 0.1 mol/l 252
Cer (IV) sulfattetrahydrat 252
Cer (IV) sulfat 0.1 mol/l 253
Cetylalkohol 254
Chinaldinrot 634
Chinidinsulfat 634
Chinolin 634
Chloracetamid 258
Chloramins-T-Natriumsalz 256
Chloranil 256
Chloranilsäure 257
Chlorat-Standardlösung 257
p-Chlorbenzaldehyd 258
Chlorbenzol 258
Chloressigsäure 258
Chlorid-Standardlösung 257
Chlorit-Standardlösung 257
Chlorobutanol Hemihydrat 259
Chloroform 259
Chloroform-d + 0.03% TMS 263
Chloroform-d1 263
p-Chlorophenol 264
Chlorwasser 789
Chrome-Standardlösung 264
Chrom (III) chloridhexahydrat 265
Chrom (III) kaliumsulfatdodecahydrat 266
Chrom (III) nitrat Nonahydrat 266
Chrom (III) oxid 266
Chrom (III) sulfat 267
Chrom (VI) oxid 267
Chrom (VI) oxidlösung 267
Chromatierte Standardlösung 264
Chromotropsäure Dinatriumsalz-Dihydrat 268
Chrysoidin Y 268
Cinchonin 268
Citral 268
Cupriethylendiamin-Lösung 285
Cupri-Weinsäure-Lösung 285
Cupri-Zitronensäure-Lösung 284
Cyanid-Standardlösung 285
Cyclohexan 286
Cyclohexan-d12 287
Cyclohexanol 287
Cyclohexanon 288
Cyclohexylamin 288
Cyclopentylmethyl Ether 289
L-Cystein 289
L-Cystin 289

D

Davidson Fixiererflüssigkeit 355
Decan 290
1-Decansulfonsäure-Natriumsalz 290
Deuteriumoxid-d 293
Deuteriumoxid-d + 0.01% DMSO 294
Deuteriumoxid-d + 0.03% TSP d4 294
Deuteriumoxid-d + 0.5% TSP d4 294
Deuteriumchlorid 20% 409
Deuteriumchlorid 1 mol/l 409
Devarda-Legierung 294
Diacetonalkohol 295
Diacetyldioxim 295
Diacetyldioximlösung 1% in Ethanol 296
Diacetyldioximnatriumsalz 295
Dibasisches Natriumcitrat-Sesquihydrat 672
Dibasisches Natriumphosphat-Dodecahydrat 702
Dibenzoylmethan 296
Dibenzoylperoxid 296
Dibutylphthalat 297
Dicalite 4158 297
m-Dichlorbenzol 298
o-Dichlorbenzol 298
p-Dichlorbenzol 298
1,2-Dichlorbenzol-d4 299
2,6-Dichlorchinon-4-chlorimid 305
Dichloressigsäure 297
1,2-Dichlorethan 299
Dichlormethan 300
Dichlormethan mit 1% Salzsäure angesäuert 305
Dichlormethan-d2 305
2,6-Dichlorphenolindophenol-Natriumsalz 305
Dichtestandards 292
Diethanolamin 306
N,N-Diethylacetamid 313
Diethylamin 306
Diethylcarbonat 307
Diethylenglycoldimethylether 307
Diethylenglykol 307
Diethylentriaminpentaessigsäure 308
Diethylether 308
n,n - Diethyl-p-phenylendiamin-sulfat 310
Diethylphthalat 310
Digitonin 311
Diisopropylamin 311
Di-isopropylether 312
N,N-Diisopropylethylamin 312
Dikaliumhexachloroplatinat 597
Dikaliumhydrogenphosphat wasserfrei 614
Dimedon 313
1,2-Dimethoxyethan 313
p-Dimethylaminobenzaldehyd 314
p-Diméthylaminobenzalrhodanine 314
N,N-Dimethylformamid 314
N,N-Dimethylformamid-d7 316
N,N-Dimethylformamiddimethylacetal 317
N,N'-Dimethylpropylenharnstoff 317
Dimethylsulfoxyd 317
Dimethylsulfoxid-d6 319
Dimethylsulfoxid-d6 + 0.03% TMS 319
Dimidiumbromid 320
Dinatriumdisulfit 694
2,4-Dinitrochlorbenzol 320
2,4-Dinitrophenylhydrazin (mit 30% Wasser) 320
1,4-Dioxan 321
1,3-Dioxolan 322
Diphenylaminlösung 1% in Schwefelsäure 322
4-Diphenylamin-Natriumsulfonat 322
sym-Diphenylcarbazid 323
sym-Diphenylcarbazon 323
Diphenylthiocarbazon 323
2,2'-Dipyridyl 323
Dithiooxamids 324
1-Dodecansulfonsäure-Natriumsalz 324
Dodecylbenzolsulfonsäure-Natrium-Salz 325
Dodecyltrimethylammoniumbromid 325
Dysprosium-Standardlösung 325

E

Ehrlich Reagenz 327
Einbasiges Calciumphosphat-Monohydrat 243
Eisen, reduziert durch Wasserstoff 428
Eisenpulver 427
Eisenstandardlösung 428
Eisen (II) ammoniumsulfathexahydrat 429
Eisen (II) ammoniumsulfat 0.1 N 430
Eisen (II) ammoniumsulfat 0.12 N 429
Eisen chlorid in Lösung 434
Eisen(II)-chlorid-Tetrahydrat 430
Eisen(II) sulfat heptahydrat 430
Eisen (II) sulfat 0.1 mol / l 431
Eisen (III) ammoniumcitrat rot 431
Eisen (III) -ammoniumcitratgrün 431
Eisen (III) ammoniumoxalat 432
Eisen (III) ammoniumsulfat dodecahydrat 432
Eisen (III) ammoniumsulfat 0.1 mol / l 433
Eisen (III) ammoniumsulfat 100 g / l 433
Eisen (III) ammoniumsulfatlösung 33% in Salpetersäure 432
Eisen (III) chloridhexahydrat 433
Eisen(III) chlorid wasserfreies sublimiertes 433
Eisen (III) chlorid 4.5% 434
Eisen (III) citrat 434
Eisen (III) nitratnonahydrat 435
Eisen (III) oxid 435
Eisen (III) sulfat 435
Eiweiß aus Eiern, getrocknet 148
Elektrolytisches Mangan 479

Product specifications are subject to changes. Please visit our website for updates.

ALPHABETISCHER INDEX

Elutionsmittel Natriumbicarbonat 327
Elutionsmittel Natriumcarbonat 327
Elutionsmittel Natriumcarbonat / Natriumbicarbonat 328
Entfärbung von Holzkohle 255
Entkalkung 290
Enzianviolett 368
Enzianviolettkarbollösung 368
Eosin B 328
Eosin Y 328
Eosin Y 0.5% alkoholische Lösung 329
Eosin Y 1% ige wässrige Lösung 328
Erbium-Standardlösung 329
Eriochromcyanin R 330
Erythrosin extra B 330
Esbach-Reagenz 330
Essigsäure 130
Essigsäureanhydrid 137
Essigsäure 12% 135
Essigsäure 20% 135
Essigsäure 25% 135
Essigsäure 27% 135
Essigsäure 30% 134
Essigsäure 45% 134
Essigsäure 80% 134
Essigsäure 96% 133
Essigsäure 0.03 mol/l (0.03N) 136
Essigsäure 0.1 mol/l (0.1N) 136
Essigsäure 1 mol/l (1N) 136
Essigsäure d-4 137
Ethanol absolut wasserfrei 331
Ethanol 96° 334
Ethanol 70% v/v 336
Ethanol 60° 337
Ethanol 50% v/v 337
Ethanol absolut vergällt 338
Ethanol vergällt 90° 339
Ethanol vergällt 94° 339
Ethanol vergällt 95° 338
Ethanol 70° geändert 339
Ethanol 95% -Diethylether 50/50 m/m -30 mg/l Phenolphthalin 507
Ethanolamin 340
Ethanol-d6 340
Ethylacetat 340
Ethylacetoacetat 343
Ethylendiamin 344
Ethylendiamintetraessigsäure 344
Ethylendiamintetraessigsäure-Dikaliumsalz-Dihydrat 345
Ethylendiamintetraessigsäure dinatriumsalz 345
Ethylendiamintetraessigsäure dinatriumsalz 0.01 mol/l (0.02N) 347
Ethylendiamintetraessigsäure dinatriumsalz 0.05 mol/l (0.1N) 347
Ethylendiamintetraessigsäure dinatriumsalz 0.1 mol/l (0.2N) 346
Ethylendiamintetraessigsäure-Kaliumsalz und Magnesiumdihydrat 347
Ethylendiamin-tetraessigsäure Tetranatriumsalz Tetrahydrat 348
Ethylendiamintetraessigsäure Tripotassium salt 348
Ethylenglycol 348
Ethylenglykol-bis-(2-aminoethyl)-n,n,n',n'-tetraessigsäure 349
Ethylenoxidlösung 349
Ethylformiat 350
Ethylmethylketon 350
Eugenol 351
EUKITT 351
Europium-Standardlösung 351

F

Fehling-Reagenzlösung A 353
Fehling-Reagenzlösung B 353
Ferricyanid-Standardlösung 353
Ferrocyanid-Standardlösung 354
Ferron 354
Fixierflüssigkeit FIXALL-HIS 356
Fixierflüssigkeit ohne Essigsäure 356
AFA-Fixiererflüssigkeit 354
Florisil 100-200 mesh 357
Florisil 60-100 mesh 356
Fluorescein 357
Fluorescein Natriumsalz 357
Fluorid-Standardlösung 358
Fluorwasserstoff Verdünnte 412
Fluorwasserstoff 39.5% 412
Fluorwasserstoff 47-51% 411
Fluorwasserstoff 50% 410
Fluorwasserstoff 0.1 mol/l 412
Folin-Ciocalteu-Reagenz 358
Formaldehyd 30% m/v 360
Formaldehyd 35% m/m 359
Formaldehyd 37% m / v neutral 359
Formaldehyd 37% w/v 358
Formaldehyd 10% v/v Lillie 360
Formaldehyd 4% m/v mit natriumchlorid 361
Formaldehyd 4%m/v pH 6.9 gepuffert 361
Formaldehyd 5%m/v pH 6.9 gepuffert 360
Formalin Essigsäure 362
Formamid 362
Fuchsin 365
Fumarsäure 365
Fuschine verfärbte Lösung 365

G

Gadolinium-Standardlösung 366
D(+)Galactose 366
Gallium-Standardlösung 367
Gallussäure monohydrat 367
Gardner Farbstandards 368
Gelatine 368
Gelbes Naphthol S 518
Gemische zum Nachweis von Restlösemitteln 509
Gerbsäure 744
Gereinigtes Wasser 788
Germanium-Standardlösung 369
Gesäuertes entionisiertes Wasser 789
Gewaschenes Kaolin 443
Giemsa-Reagenz 369
Glaswolle 370
D(+)Glucose wasserfrei 370
D(+)-Glucose-Monohydrat 371
L-Glutaminsäure 371
Glutaraldehydlösung 50% 372
Glycerin (30°Bé) 372
Glycerin 90% (28° Bé) 373
Glycin 373
Glycolsäure 374
Glyoxal-Standardlösung 374
Glyoxalbis(2-hydroxyanil) 374
Goldenes (III) säurechloridtrohydrat 375
Goldstandardlösung 375
Gram - Hucker Kit 376
Griess Reagenz 376
Griess Reagenz A 376
Griess Reagenz B 376
Grüne Malachitlösung 0.5% in wasserfreier Essigsäure 476
Gummi arabisch 377

H

Hafnium-Standardlösung 378
Hämatoxylin 379
Hämatoxylinlösung nach Mayer 379
Hanus-Reagenz 379
Helles Coomassie Blue R 250 274
Hellgrün 461
Hemalun in Lösung nach Carazzi 378
Hemalun in Lösung nach Mayer 378
Hengste im Brix-Grad 209
n-Heptan 382
n-Heptan 99% 380
1-Heptansulfonsäure-Natriumsalz 383
Heptan Isomerengemisch 382
Hexachloroplatinsäure-Hexahydrat 383
Hexadecan 383
Hexafluor-2-propanol 384
Hexamethyldisilazan 384
Hexamethylentetramin 384
Hexan Isomerengemisch 388
Histamindihydrochlorid 389
Histolemon 390
Holmium-Standardlösung 390
Holmiumoxid 391
Holmiumperchloratlösung 391
Holzkohle aktiviert 255
Holzkohlegemüse 256
Hyamine 1622 391
Hyamin 1622-Lösung 0.004M 392
Hydraziniumchlorid 392
Hydriertes Chloral 256
Hydrindantine 392
Hydroalkoholische Universalanzeige für pH 1-11 421
Hydrochinon 416
Hydroxylaminlösung, alkoholisch 417
Hydroxylaminsulfat 417
Hypophosphorige Säure 50% 418
n-Hexan 386
n-Hexan 99% 385
1-Hexansulfonsäure-Natriumsalz 389
1-Hexansulfonsäure-Natriumsalz-Monohydrat 389
L-Histidinmonohydrochloridmonohydrat 390
p-Hydroxybenzaldehyd 416
8-Hydroxychinolin 417

I

Idrimeter Erba Anzeige C 419
Idrimeter Erba Lösung A 419
Idrimeter Erba Lösung B 419
Imidazol 419
Immersionsöl 420
Indigo Carmine getrocknet 422
Indigo-Karmin-Lösung 422
Indikator Ferroin 0.025 mol/l 354
Indikator für ammoniakalischen Stickstoff 420
Indikator für die Jodometrie 420
Indium-Standardlösung 422
indol 423
Inosit 423
Inulin 424
Iodoform 427
Iodoplatinat-Reagenz 427
Iridium-Standardlösung 427
Isoamylacetat 436
Isoamylalkohol 436
Isobutanol 437
Isobutylacetat 438
Isobutylmethylketon 500
Isohexan 438
Isooctan 439
Isopar G 440
Isopentan 441
Isopropylacetat 441
Isopropylamin 442
Isopropylmyristat 442

J

Jod 424
Jod 0.005 mol/l (0.01N) 426

ALPHABETISCHER INDEX

Jod 0.01 mol / l (0.02 N) 426
Jod 0.05 mol/l (0.1N) 425
Jod 0.5 mol/l (1N) 425
Jod 10 ppm 426
Jod wieder sublimiert 425
Jodbromidlösung 426
Jodid-Standardlösung 424
Jodwasserstoffsäure 57% 393

K

Kalium-Standardlösung 582
Kaliumacetat 583
Kaliumbicarbonat 584
Kaliumbisulfat 584
Kaliumbromat 584
Kaliumbromat 0.0167 mol/l (0.1N) 585
Kaliumbromat 0.02 mol/l (0.12N) 585
Kaliumbromat 0.033 mol / l (0.198 N) 585
Kaliumbromid 585
Kaliumcarbonat 586
Kaliumchlorid 586
Kaliumchlorid-Lösung 590
Kaliumchlorid-Lösung gesättigt 590
Kaliumchlorid 0.01 mol/l (0.01N) 589
Kaliumchlorid 0.1 mol/l (0.1N) 589
Kaliumchlorid 1 mol/l (1N) 588
Kaliumchlorid 3 mol/l (3N) 588
Kaliumchlorid 3 mol/l (3N) Lösung aus Wasser und Glycerin 588
Kaliumchlorid 3 mol/l (3N) mit Silberchlorid 588
Kaliumchlorid 3.5 mol/l (3.5 N) 587
Kaliumchlorid 3.5 mol/l (3.5N) mit Silberchlorid 587
Kaliumchlorid 12g/l 589
Kaliumchlorid 25g/l (HCl) 589
Kaliumchromat 590
Kaliumchromatlösung bei 10% 590
Kaliumchromatlösung bei 5% 591
Kaliumdichromat 592
Kaliumdichromat 0.0167 mol/l (0.1N) 593
Kaliumdichromat 0.04 mol/l(0.24N) in 80 g/l HgSO4 593
Kaliumdichromat 0.0417 mol/l (0.25N) 593
Kaliumdichromat 0.167 mol/l (1N) 592
Kaliumdichromat 0.1414 g/l 594
Kaliumdichromatlösung 0.5% 592
Kaliumdichromatlösung 106 g/l 594
Kaliumdichromat - Lösung in Schwefelsäure 592
Kaliumferricyanid 594
Kaliumferrocyanid-Trihydrat 595
Kaliumferrocyanidlösung 10% 595
Kaliumferrocyanidlösung 53 g / l 596
Kaliumfluorid 596
Kaliumfluorid-Dihydrat 596
Kaliumguajacolsulfonat 596
Kaliumhydroxid, Flocken 599
Kaliumhydroxid, Pellets 599
Kaliumhydroxid 28% 601
Kaliumhydroxid 33% 601
Kaliumhydroxid 38% (40 ° Bé) in Wasser 601
Kaliumhydroxid 45% 600
Kaliumhydroxid 0.1 mol/l (0.1N) 604
Kaliumhydroxid 0.1 mol/l (0.1N) in ethanol 605
Kaliumhydroxid 0.1 mol/l (0.1N) in methanol 605
Kaliumhydroxid 0.23 mol/l (0.23N) 604
Kaliumhydroxid 0.25 mol/l (0.25N) 604
Kaliumhydroxid 0.46 mol/l (0.46N) 604
Kaliumhydroxid 0.5 mol/l (0.5N) 602
Kaliumhydroxid 0.5 mol/l (0.5N) in ethanol 603
Kaliumhydroxid 0.5 mol/l (0.5N) in methanol 603
Kaliumhydroxid 1 mol/l (1N) 602
Kaliumhydroxid 2 mol / l (2N) in Ethanol 602
Kaliumhydroxidlösung 3% in Ethanol 601
Kaliumhydroxid in Lösung 606
Kaliumiodat 606
Kaliumiodid 607
Kaliumiodid 10% 608
Kaliumiodid 3.9% 608
Kaliumiodid-Stärke-Papier 608
Kaliumjodatsäure 597
Kaliumjodidlösung 608
Kaliumjodat 0.05 mol/l (0.3N) 606
Kaliumjodat 0.0167 mol/l (0.1N) 606
Kaliumjodat 0.00167 mol/l (0.01N) 607
Kaliumjodobismutatlösung 609
Kaliumdisulfit 609
Kalium natrium tartrat tetrahydrate 617
Kaliumnitrat 609
Kaliumnitrat 1 mol/l (1N) 610
Kaliumoxalat-Monohydrat 610
Kaliumperiodat 611
Kaliumpermanganat 611
Kaliumpermanganat 0.002 mol/l (0.01N) 613
Kaliumpermanganat 0.02 mol/l (0.1N) 612
Kaliumpermanganat 0.2 mol/l (1N) 612
Kaliumpermanganatlösung 3% 612
Kaliumpermanganat - Phosphorlösung 613
Kaliumpersulfat 613
Kaliumphosphat monobasisch 615
Kaliumphthalsäure 597
Kaliumphthalsäure 0.1 mol/l (0.1N) 598
Kaliumphthalsäure 0.2 mol/l (0.2N) 598
Kaliumpyroantimonlösung 616
Kaliumpyrophosphat 616
Kaliumpyrosulfat 616
Kaliumsäurepyroantimoniat 615
Kaliumsorbat 617
Kaliumsulfat 618
Kaliumsulfocyanidlösung 621
Kaliumtartrat 618
Kalium L-Tartrat monobasisch 619
Kaliumtellurit 619
Kaliumtetraiodomercurat - alkalische Lösung 619
Kaliumthiocyanat 620
Kaliumthiocyanat 0.1 mol/l (0.1N) 620
Kaliumthiocyanat 5% 620
Kalkwasser 184
Kampfer, Natürlicher 246
Kampfer, Synthetischer 246
Kanadabalsam 247
Karl Fischer Anolyt-Lösung - ohne FCKW 446
Karl-Fischer-Anolytlösung - ohne Pyridin 446
Karl-Fischer-Anolytlösung für den Ofen 446
Karl-Fischer-Anolytlösung für Ketone und Aldehyde - ohne Methanol 446
Karl-Fischer-Anolytlösung für Öle 447
Karl-Fischer-Katholytlösung 447
Karl-Fischer-Katholyt-Lösung - Pyridinfrei 447
Karl Fischer Lösungsmittel 1 Komponente für Aldehyde und Ketone 444
Karl Fischer Lösungsmittel 2 Komponente 445
Karl Fischer Lösungsmittel für 1 Komponentenöle 444
Karl Fischer Lösungsmittel für 2 Komponentenöle 445
Karl Fischer Pufferlösung 445
Karl Fischer Reagenz 1 Komponente 2 mg H 2 O / ml 443
Karl-Fischer-Reagenz 1 Komponente 5 mg H2O / ml 443
Karl-Fischer-Reagenz für Zellen ohne Diaphragma 447
Karl Fischer Solvent 2 Komponenten für Aldehyde und Ketone - Methanolfrei 445
Karl Fischer Standardwasser 0.10 mg / g 449
Karl Fischer Standardwasser 1.0 mg / g 448
Karl Fischer Standardwasser 10.0 mg / g 448
Karl Fischer Standardwasser 5.0 mg / ml 448
Karl Fischer Titriermittel 2 Komponenten 5 mg H2O / ml - nicht hygroskopisch 444
Karl Fischer titriert 1 Komponente 5 mg H2O / ml für Aldehyde und Ketone 443
Karl Fischer titriert 2 Komponenten 2 mg H2O / ml - nicht hygroskopisch 444
Kasein 249
Kerosin 449
Kieselgel 60A 0.06 ÷ 0.20 mm 650
Kieselgel 60A 20 - 4µ 648
Kieselgel 60A 35 - 70µ 649
Kieselgel 60A 40 - 63µ 649
Kieselgel 60A 6 - 35µ 648
Kieselgel 60A 70 - 200µ 650
Kieselgelgranulat mit kobaltfreiem Indikator 651
Kieselgel granuliert 650
Kieselgur besteht 449
Kieselwolframsäure 652
Kjedahl-Katalysator für die Wasseranalyse 450
Kjeldahl entschäumer 449
Kjeldahl-Katalysator nach Wieninger 450
Kjeldahl-Katalysator Selen- und Titanfrei 450
Kjeldahl-Selen-Katalysator 450
Kjeldahl-Titankatalysator 451
Kobalt-Standardlösung 270
Kobalt (II) acetattetrahydrat 271
Kobalt (II) ammoniumsulfathexahydrat 271
Kobalt (II) chloridhexahydrat 272
Kobalt (II) chlorid in Lösung 272
Kobalt (II) nitrathexahydrat 272
Kobalt (II) sulfatheptahydrat 273
Kongo-Rot 273
Kongorotes Papier 273
Kovac-Reagenz 451
m-Kresol 282
o-Kresol 282
o-Kresolrot 283
Kresol-rote Lösung Natronlauge 0.1N/ ethanol 96% / wasser 283
O-Kresol-rote Lösung 0.2% in Ethanol 283
m-Kresolpurpur 282
Kristallviolett 284
Kristallviolettlösung 0.5% in wasserfreier Essigsäure 284
Kristallviolettoxalat für Gram-Hucker Kit 284
Kupfer reduziert, Pulver 275
Kupferelektrolyt überbacken, Blatt 274
Kupferelektrolyt, Draht 275
Kupferelektrolyt, Späne 274
Kupfer-Standardlösung 275
Kupfer (I) chlorid 276
Kupfer (I) chloridlösung 7% in Ammoniak 277
Kupfer (I) iodid 277
Kupfer (II) acetathydrat 277
Kupfer(II)-chlorid-Dihydrat 278
Kupfer (II) ethylendiaminhydroxid 1mol/l 285
Kupfer(II)-hydroxidcarbonat 278
Kupfer(II) nitrat-trihydrat 279
Kupfer(II) oxid 279
Kupfer (II) oxiddraht 279
Kupfer(II) sulfat wasserfrei 280
Kupfer(II)-sulfat-Pentahydrat 280
Kupfersulfat 12.5 % 281
Kupfer(II)-sulfat in Lösung 281

L

Lackmus 465
Lactophenol von Amman in Lösung 160
Lactophenolblau-Lösung 452
Lactose Monohydrat 452
Lanolin Wasserfreies 453
Lanthan-Standardlösung 453
Lanthanchlorid 25 g/l 454
Lanthannitrat 0.1 mol/l 454

Product specifications are subject to changes. Please visit our website for updates.

Lanthannitrat-Hexahydrat 454
Lanthannitratlösung 50 g/l 454
Lanthanoxid 455
Lävulose 364
Leitfähigkeitsstandard 1.30 µS/cm 715
Leitfähigkeitsstandard 5 µS/cm 716
Leitfähigkeitsstandard 10 µS/cm 716
Leitfähigkeitsstandard 20 µS/cm 716
Leitfähigkeitsstandard 50 µS/cm 716
Leitfähigkeitsstandard 84 µS/cm 717
Leitfähigkeitsstandard 100 µS/cm 717
Leitfähigkeitsstandard 147 µS/cm 717
Leitfähigkeitsstandard 200 µS/cm 717
Leitfähigkeitsstandard 500 µS/cm 717
Leitfähigkeitsstandard 1000 µS/cm 718
Leitfähigkeitsstandard 1413 µS/cm 718
Leitfähigkeitsstandard 5000 µS/cm 718
Leitfähigkeitsstandard 10.000 µS/cm 718
Leitfähigkeitsstandard 12 880 µS/cm 718
Leitfähigkeitsstandard 20.000 µS/cm 719
Leitfähigkeitsstandard 50.000 µS/cm 719
Leitfähigkeitsstandard 100.000 µS/cm 719
Leitfähigkeitsstandard 150.000 µS/cm 719
Leitfähigkeitsstandard 200.000 µS/cm 720
Leitfähigkeitsstandard 300 000 µS/cm 720
Leitfähigkeitsstandard 350.000 µS/cm 720
Leitfähigkeitsstandard 450.000 µS/cm 720
Leitfähigkeitsstandard 500.000 µS/cm 721
L (+) - Leucin 461
Lithium-Standardlösung 462
Lithiumcarbonat 463
Lithiumchlorid 463
Lithiumhydrid 464
Lithiumhydroxid wasserfrei 464
Lithiumhydroxidmonohydrat 464
Lithiummethylat 0.1 mol/l (0.1 N) 464
Lithiumnitrat 465
Lithiumsulfat-Monohydrat 465
Lithiumtetraborat wasserfrei 465
Lösliche Stärkelösung 721
Lösung NDF 519
Lösung NDF PLUS 519
Lösungsmittel Plus 714
Lugol Reagenz konzentrierte Lösung 466
Lugol-Reagenzlösung 466
Lugol-Reagenzlösung für das Gram-Hucker-Kit 466
Luminol 467
Lutethium-Standardlösung 467
L (+) Lysinmonohydrochlorid 467

M

Magnesiumpulver 468
Magnesiumstandardlösung 468
Magnesiumbänder 468
Magnesiumacetat-Tetrahydrat 470
Magnesiumcarbonat basisch 470
Magnesiumchlorid-Hexahydrat 471
Magnesium dreht sich um 468
Magnesiumglycerophosphat 471
Magnesium hydrogenphosphat Trihydrat 472
Magnesiumhydroxid 472
Magnesiumnitrat-Hexahydrat 472
Magnesiumnitrat 10 g/l 473
Magnesiumoxid 473
Magnesiumoxid schwer 473
Magnesiumperchlorat 474
Magnesiumperoxid 474
Magnesiumstearat 474
Magnesiumsulfat heptahydrat 475
Magnesiumsulfat wasserfrei 475
Maisstärke 476
Malachitgrün 476
Maleinsäure 477
Maleinsäureanhydrid 477
Malonsäure 478
Maltose-Monohydrat 478
D (-) Mandelsäure 478
Mangan-Standardlösung 479
Mangan (II) acetat Tetrahydrat 480
Mangan (II) chlorid Tetrahydrat 480
Mangan(II) sulfat Monohydrat 481
Mangan (IV) dioxid 481
D - Mannit 481
L-Menthol 483
2-Mercaptoethanol 483
Metanilgelb 487
Metaphosphorsäure 487
Methanol 488
Methanol-d1 494
Methanol-d3 494
Methanol-d4 493
Methanol-d4 + 0.03% TMS 494
Methanol - Salzsäure 493
Methanol + 0.1% v/v Ameisensäure 492
Methanol + 0.1% v/v trifluoressigsäure 493
Methansulfonsäure 487
DL-Methionin 494
2-Methoxy-ethanol 495
alpha-Methoxy(phenyl)essigsäure 495
Methylacetat 495
4-Methylaminophenolsulfat 496
Methylbenzoat 496
3-Methyl-2-benzothiazolinhydrazonhydrochlorid 497
Methylblau 497
2-Methylbutanol-2 180
Methylcyclohexan 497
Methylenblau 498
Methylenblau gesättigte Lösung 498
Methylenblaulösung 1% 499
Methylenblau-phenolische hydroalkoholische Lösung 248
Methylgelb 506
Methylgrün 499
Methyl-4-hydroxybenzoat 499
Methylisoamylketon 500
Methyljodid 500
Methylorangenlösung 0.1% 501
N-Methyl-2-pyrrolidon 502
Methylrot 503
Methylrotlösung 504
Methylrotlösung 0.1% in Ethanol 504
Methylrotlösung Wasser / Ethanol 0.2% 503
Methylsalicylat 505
2-Methyltetrahydrofuran 505
4-Methyltetrahydropyran 506
Methylthymolblau Natriumsalz 506
L-(+)-Milchsäure 452
Mische Diethylether / Ethanol 70/30 m/m 507
Mischindikator mit Methylorange 502
Mischindikator mit Methylrot 504
Mischung C.H.M. 507
Mischung Ethanol absolut/Diethylether 50/50 (v/v) 508
Mischung Ethanol absolut/Diethylether 50/50 (m/m) 507
Mischung Ethanol 95 ° / Isopropanol 508
Mischung Ethanol 99 ° / Isopropanol 508
Mischung für Bromindex 508
Mischung zur Kontrolle von Lötungen 509
Molekularsieb 3 A 510
Molekularsieb 4A 510
Molybdän-Standardlösung 510
Molybdän (VI) oxid 511
Molybdovanadisches Reagens 511
Monobasisches Kaliumphosphat 0.2 mol/l (0.2N) 615
Morpholin 512
Mucicarmin in hydroalkoholischer Lösung 512
Multianionen-Standard für die Ionenchromatographie 512
Multielement-Standard für ICP 513
Multielement-Standard für ICP und ICP-MS 515
Murexid 516

N

Naphthalin 517
1-Naphthol 517
2-Naphthol 517
alpha-Naphtolbenzol 518
alpha-Naphtolbenzollösung 0.2% in Essigsäure 518
N-(1-Naphthyl) ethylendiamin-dihydrochlorid 518
Natriumstandardlösung 661
Natriumacetat wasserfrei 662
Natriumacetat-trihydrat 662
Natriumalginat 663
Natriumaluminat 663
Natriumammoniumhydrogenphosphat 664
Natriumarsenit 0.05 mol/l (0.1N) 664
Natriumarsenit 0.1 mol/l (0.2N) 664
Natrium L-ascorbat 664
Natriumazid 665
Natrium-basisches Citrat-Bihydrat 673
Natriumbenzoat 665
Natriumbicarbonat 665
Natriumbisulfat-Monohydrat 666
Natriumcarbonat wasserfrei 667
Natriumcarbonat-Monohydrat 668
Natriumcarbonat 0.05 mol/l (0.1N) 669
Natriumcarbonat 0.5 mol/l (1N) 669
Natriumcarbonat 20% 670
Natriumcarbonat wasserfreie Lösung 670
Natriumchlorid 670
Natriumchlorid 0.1 mol/l (0.1N) 672
Natriumchlorid 5 mol/l (5N) 672
Natriumcitrat tribasisch wasserfrei 672
Natrium-cyanoborhydrid 674
Natriumdichlorisocyanatdihydrat 674
Natriumdiethyldithiocarbaminat-Trihydrat 674
Natriumdithionit 678
Natriumdodecylsulfat 693
Natriumfluorid 675
Natriumformiat 675
Natriumgluconat 676
Natriumglutamatsäure 676
Natriumglycerophosphatpentahydrat 676
Natriumhexafluorosilikat 677
Natriumhexametaphosphat 677
Natriumhydrogenphosphat Dihydrat 701
Natriumhydrogenphosphat wasserfrei 700
Natriumhydrogentartratmonohydrat 678
Natriumhydroxid, Pellets 679
Natriumhydroxid, Perlen 678
Natriumhydroxid 521
Natriumhydroxid auf Kieselgel 680
Natriumhypochloritlösung 12.5% 692
Natriumhypochloritösung 692
Natriumhypophosphit 693
Natriumiodid 693
Natriumkobaltnitrit 673
Natriummetaperiodat 694
Natriummetaphosphat 695
Natriummethanolat 0.1 mol/l 695
Natriummolybdat 695
Natriumnitrat 696
Natriumnitrit 696
Natriumnitrit 0.1 mol/l (0.1N) 697
Natriumnitrit 500 g/l 697
Natriumnitroprussid-Dihydrat 697
Natriumoxalat 698
Natriumperborat-Tetrahydrat 698

ALPHABETISCHER INDEX

Natriumperchlorat-Monohydrat 699
Natriumperoxid 699
Natriumpersulfat 699
Natriumpersulfat 1 mol/l 700
Natriumphosphat monobasisch Monohydrat 703
Natriumphosphat tribasisch Dodecahydrat 704
Natriumphosphat-einbasiges Dihydrat 703
Natrium-o-phosphitpentahydrat 704
Natriumpyrophosphatdecahydrat 705
Natriumsalicylat 705
Natriumsilicat 705
Natriumsuccinat-Hexahydrat 706
Natriumsulfat wasserfrei 706
Natriumsulfat-Decahydrat 707
Natriumsulfid-Nonahydrat 707
Natriumsulfid-Nonahydrat-Lösung 707
Natriumsulfid-Trihydrat 708
Natriumsulfit wasserfrei 708
Natriumsulfocyanid 709
Natriumtartratdihydrat 709
Natriumtetraboratdecahydrat 710
Natriumtetraborat wasserfrei 709
Natriumtetrahydroborat 666
Natriumtetraphenylborat 710
Natriumthiosulfat pentahydrat 711
Natriumthiosulfat wasserfrei 711
Natriumthiosulfat 0.01 mol/l (0.01N) 713
Natriumthiosulfat 0.0197 mol/l (0.0197N) 713
Natriumthiosulfat 0.0394 mol/l (0.0394N) 713
Natriumthiosulfat 0.1 mol/l (0.1N) 712
Natriumthiosulfat 0.2 mol/l (0.2N) 712
Natriumthiosulfat 0.5 mol/l (0.5N) 712
Natriumthiosulfat 1 mol/l (1N) 712
Natriumwolframatdihydrat 714
Natronkalk 660
Natronlauge 691
Natronlauge 5% w/v 684
Natronlauge 10% w/v 683
Natronlauge 20% w/v 683
Natronlauge 20% w/w 683
Natronlauge 30% 682
Natronlauge 32% 681
Natronlauge 35 % 681
Natronlauge 35-37% 681
Natronlauge 40% 680
Natronlauge 50 % 680
Natronlauge 0.01 mol/l (0.01N) 690
Natronlauge 0.1 mol/l (0.1N) in Ethanol 690
Natronlauge 0.1 mol/l (0.1N) 689
Natronlauge 0.2 mol/l (0.2N) 689
Natronlauge 0.25 mol/l (0.25N) 688
Natronlauge 0.357 mol/l (0.357N) 688
Natronlauge 0.5 mol/l (0.5N) 688
Natronlauge 0.7 mol/l (N/1.4) 687
Natronlauge 1 mol/l (1N) 686
Natronlauge 1.2 mol/l (1.2N) 686
Natronlauge 1/9 mol/l (N/9) 689
Natronlauge 2 mol/l (2N) 686
Natronlauge 3 mol/l (3N) 685
Natronlauge 4 mol/l (4N) 685
Natronlauge 5 mol/l (5N) 685
Natronlauge 6 mol/l (6N) 684
Natronlauge in Methanol 691
Natronlauge-d 1N 691
Natronlauge-d 30% 691
Neocupronhydrochlorid 519
Neodym-Standardlösung 519
Nessler Reagenzlösung A 520
Nesslers Reagenzlösung 520
Neutralrot 521
Niacinamid 525
Nickel, Pulver 521
Nickel-Standardlösung 521
Nickel (II) acetattetrahydrat 522
Nickel (II) ammoniumsulfathexahydrat 523
Nickel (II) carbonat basisch 523
Nickel (II) chloridhexahydrat 523
Nickel (II) nitrathexahydrat 524
Nickel (II)-nitrat 10g/l 524
Nickel (II) sulfathexahydrat 524
Nigrosine 525
Nikotinsäure 525
Ninhydrin 526
Ninhydrinlösung 526
Ninhydrin und Zinn (II) chlorid Reagenz 526
Niob-Standardlösung 526
Nital-Lösung 4% 527
Nitrat-Standardlösung 527
Nitrilotriessigsäure 534
Nitrit-Standardlösung 534
Nitrit-Kalium 610
m-Nitrobenzaldehyd 534
m-Nitrobenzoesäure 535
Nitrobenzol 535
p-Nitrophenol 535
p-Nitrophenol-Lösung 0.1% in Wasser 535
Nonylphenol ethoxylated - 10 ETO 536
Nuklear schnell rot 536

O

n-Octan 537
Octan 80-Mischung 537
1-Octanol 538
1-Octansulfonsäure-Natriumsalz 537
1-Octansulfonsäure-Natriumsalz-Monohydrat 538
Oelrot o 636

Ö

Öl aus Zedernholz 538

O

Olsäure 539
Orange G 539
Orange II 539
Orangenschnaps 501
Orcein 540
Orcinol monohydrate 540
Organisch standard: PAH Multielement-Mischung 542
Organisch standard: PCB 540
Organisch standard: PCB Multielement-Mischung 543
Organisch standard: Pestizidmischung 546
Osmium-Standardlösung 549
Osmolalitätsstandards 550
Osmolalitätsstandards basierend auf Proteinen 550
Osmolalitätsstandards basierend auf Urin 550
Oxalsäuredihydrat 550
Oxalsäure 0.005 mol/l (0.01N) 551
Oxalsäure 0.05 mol/l (0.1N) 551
Oxalsäure 0.5 mol/l (1N) 551

P

Palladium-Standardlösung 553
Palladium (II) chlorid 553
Palladiumnitrat 2 g/l Lösung 554
Papanicolaou Haematoxylin-Lösung nach Harris 554
Papanicolaou Lösung OG 6 555
Papanicolaou-Lösung EA 50 554
Paraffin 56-58 ° C Erbaplast (ohne DMSO) 555
Paraffinöl 555
Paraformaldehyd 556
Pararosanilin in Lösung verfärbt 556
n-Pentan 558
n-Pentan 99% 557
1-Pentansulfonsäure-Natriumsalz 559
1-Pentansulfonsäure-Natriumsalz-Monohydrat 560
Pepsin HCl 560
Perchlorsäure 65% 561
Perchlorsäure 65-71% 560
Perchlorsäure 0.01 mol/l (0.01N) 562
Perchlorsäure 0.05 mol/l (0.05N) 562
Perchlorsäure 0.1 mol/l (0.1N) 561
Perchlorsäure 0.1 mol/l (0.1N) essigsäure 562
Perchlorsäurelösung 562
Perfluorbuttersäure 380
Periodsäure 563
Petrolether 100 -140 ° C 564
Petrolether 30 - 50 ° C 568
Petrolether 30-40 ° C 569
Petrolether 35 - 60°C 567
Petrolether 40 - 60°C 567
Petrolether 40 - 65 ° C 566
Petrolether 40 - 70 ° C 566
Petrolether 55 - 85 ° C 565
Petrolether 60 -8 0 ° C 565
Petrolether 75 -120 ° C 565
Petrolether 80 -100 ° C 564
Petrolether 80 -120 ° C 564
Petroleum 563
Phenacetin 146
o-Phenanthrolin-Eisen (II) sulfat 569
o-Phenanthrolinmonohydrat 569
Phenol 570
Phenol flüssiggas 85% 570
Phenolische Toluidinblau-Hydroalkohol-Lösung 248
Phenolphthalein 571
Phenolphthaleinlösung 0.1% 572
Phenolphthaleinlösung 0.1% in Ethanol 572
Phenolphthaleinlösung 1% in Ethanol 571
Phenolrot 570
Phenolrote Lösung 571
Phenolrote Lösung 0.2% in Ethanol 571
L-Phenylalanin 572
2-Phenylethanol 573
Phenylhydrazinhydrochlorid 573
Phenylhydrazinhydrochloridlösung 573
pH-Indikatorpapier 420
Phloroglucin 574
Phloxin B 574
Phosphat-Standardlösung 574
Phosphomolybdotungstica-Reagenz 576
Phosphor-Standardlösung 577
Phosphorige Säure 577
Phosphormolybdänsäure 576
Phosphorpentoxid 578
Ortho-phosphorsäure 10% 549
ortho-Phosphorsäure 75% 549
ortho-Phosphorsäure 85% 547
ortho-Phosphorsäure 99% 547
ortho-Phosphorsäure-d3 85% in D2O 549
Phosphorwolframsäure 578
Phosphorwolframsäure-Lösung 578
Phosphosulfonsäure 578
o-Phthalaldehyd 579
Phthalsäure 579
Phthalsäureanhydrid 579
Pikrinsäure 579
Piperidin 580
Platin-Standardlösung 580
Polyvinylpyrrolidon 581
Ponceau rot BS 581
Ponceau rot S 581
Praseodym-Standardlösung 621
Primärlösungen für den Einfärbungsgrad von Flüssigkeiten 290
L-Prolin 622
1,3-Propandiol 627
1-Propanol 622
2-Propanol 623
2-Propanol 70% 626

DASITGROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

ALPHABETISCHER INDEX

1-Propansulfonsäure-Natriumsalz 627
Propionaldehyd 628
Propionsäure 628
n-Propylacetat 628
Propylencarbonat 629
Propylenglycol 629
Propyl-p-oxybenzoat 629
Puffer acetate pH 4.5 216
Puffer acetate pH 4.6 130
Puffer acetate pH 6.0 130
Puffer Ammoniumchlorid pH 10.0 170
Puffer Ammoniumchlorid pH 10.7 169
Puffer Ammoniumchlorid pH 9.5 170
Puffer pH 1 216
Puffer pH 1.68 216
Puffer pH 2 217
Puffer pH 3 217
Puffer pH 3.5 217
Puffer pH 3.56 218
Puffer pH 3.7 218
Puffer pH 4 218
Puffer pH 4.5 219
Puffer pH 4.62 219
Puffer pH 5 219
Puffer pH 5.2 219
Puffer pH 6 220
Puffer pH 6.8 220
Puffer pH 6.88 220
Puffer pH 7 220
Puffer pH 7.2 221
Puffer pH 7.2 Weise 221
Puffer pH 7.4 222
Puffer pH 8 222
Puffer pH 9 222
Puffer pH 9.22 223
Puffer pH 10 223
Puffer pH 10.06 224
Puffer pH 11 224
Puffer pH 12 225
Puffer pH 13 225
Puffer zur Einstellung der Gesamtionenstärke 769
Pufferlosung phosphat pH 2.0 576
Pufferlosung phosphat pH 3.0 576
Pufferlosung phosphat pH 5.5 575
Pufferlosung phosphat pH 6.0 575
Pufferlosung phosphat pH 6.8 575
Pufferlosung phosphat pH 7.2 575
Pufferlosung phosphat pH 7.4 575
Pufferlosung phosphat pH 9.0 574
Pyridin 630
Pyridin-d5 631
1- (2-pyridyl-azo) -2-naphthol 632
Pyrogallol 632
Pyrrolidindithiocarbaminsäure-Ammoniumsalz 633

Q

Quecksilberbromidpapier 484
Quecksilber-Standardlösung 484
Quecksilber (I) chlorid 485
Quecksilber (II) chlorid 485
Quecksilber (II) chloridlösung 485
Quecksilber (II) iodid 486
Quecksilber (II) rotes Oxid 486
Quecksilber (II) sulfatlösung 486
Quecksilber (II) sulfat 486

R

Raffinose 635
Raney-Legierung 635
Reagenz für die Lipolyse 635
Reagenz mit Aminohippursäure 158
Reagenz von Benedict 197
Reagenz von Mai Grünwald 482
Reagenz von Mayer 482
Reduktionslösung 220 mV bei 25 ° C 636
Reduktionslösung 468 mV bei 25 ° C 637
Reinecke salt 637
Reisstärke 640
Resorcin 637
L-(+)-Rhamnose 638
Rhenium-Standardlösung 638
Rhodium-Standardlösung 638
Rhodol-Lösung 639
Riboflavin 639
D(-)-Ribose 639
Rizinusöl 249
Rosolsäure 640
Rote Kongo-Lösung 274
Rubidium-Standardlösung 640
Ruthenium-Standardlösung 641

S

D(+)-Saccharose 726
Safranin T 642
Safranin T hydroalkoholische Lösung für Gram-Hucker Kit 642
Salicylaldehyd Azin 642
Salicylsäure 642
Salpetersäure 10% 532
Salpetersäure 18% 532
Salpetersäure 65% 531
Salpetersäure 67.5% (42° Be) 530
Salpetersäure 67-69% 530
Salpetersäure 67-70% 529
Salpetersäure 69.5% 528
Salpetersäure raucht zu 99% 528
Salpetersäure 0.1 mol/l (0.1N) 533
Salpetersäure 1 mol/l (1N) 533
Salpetersäure 2 mol/l (2N) 533
Salpetersäure 8 mol/l (8N) 532
Salpetersäure verdünnt 534
Salzsäure 50% v/v 393
Salzsäure 37% 394
Salzsäure 34-37% 396
Salzsäure 32-35% 396
Salzsäure 32% (20°Be) 397
Salzsäure 32% 397
Salzsäure 29 - 31 % 398
Salzsaure 26% 398
Salzsäure 25% m/v 398
Salzsäure 23% 399
Salzsäure 20% 399
Salzsäure 12% 399
Salzsäure 10% 400
Salzsäure 8% 400
Salzsäure 5% 400
Salzsäure 1.128% m/v 401
Salzsäure 9 mol/l (9N) 401
Salzsäure 6 mol/l (6N) 401
Salzsäure 5 mol/l (5N) 402
Salzsäure 4 mol/l (4N) 402
Salzsäure 3 mol/l (3N) 403
Salzsäure 2 mol/l (2N) 403
Salzsäure 1 mol/l (1N) 403
Salzsäure 0.714 mol/l (N/1.4) 404
Salzsäure 0.5 mol/l (0.5N) 405
Salzsäure 0.2 mol/l (0.2N) 405
Salzsäure 0.2 mol / l (0.2 N) in Propanol-2 405
Salzsäure 0.1 mol/l (0.1N) 406
Salzsäure 0.1 mol / l (0.1 N) in Ethanol 407
Salzsäure 0.1 mol / l (0.1 N) in Propanol-2 407
Salzsäure 0.0714 mol/l (N/14) 407
Salzsäure 0.05 mol / l (0.05 N) 407
Salzsäure 0.04 mol/l (0.04N) 408
Salzsäure 0.02 mol/l (0.02N) 408
Salzsäure 0.01 mol/l (0.01N) 408
Salzsäure verdünnen 409
Salzsäure, bromiert 409
Samarium-Standardlösung 643
Sand gereinigt 644
Sand von Fontainebleau 644
Saponin 644
Säurefuchsin 364
Saybolt Farbstandards 644
Scandium-Standardlösung 645
Schiff-Reagens für Aldehyde 645
Schiff-Reagenz für NICHT Färbung 646
Schmelzpunktstandards 483
Schneller grüner FCF 353
Schwarzer Sudan B 726
Schwarzes Eriochrom T 329
Schwefel sublimiert und gewaschen 730
Schwefelsäure 10% v/v 737
Schwefelsäure 20% 736
Schwefelsäure 25% 736
Schwefelsäure 30% 736
Schwefelsäure 35% (30°Be) 736
Schwefelsäure 50% 735
Schwefelsäure 62% 735
Schwefelsäure 69% 735
Schwefelsäure 72% 735
Schwefelsäure 85% 734
Schwefelsäure 90% 734
Schwefelsäure 93-98% 733
Schwefelsäure 95 - 97 % 733
Schwefelsäure 96% 731
Schwefelsäure 96% (66° Be) 731
Schwefelsäure 98% 730
Schwefelsäure 0.0025 mol/l (0.005N) 742
Schwefelsäure 0.005 mol/l (0.01N) 742
Schwefelsäure 0.01 mol/l (0.02N) 741
Schwefelsäure 0.02 mol/l (0.04N) 741
Schwefelsäure 0.025 mol/l (0.05N) 741
Schwefelsäure 0.05 mol/l (0.1N) 740
Schwefelsäure 0.1 mol/l (0.2N) 740
Schwefelsäure 0.125 mol/l (0.25N) 740
Schwefelsäure 0.13 mol/l (0.26N) 740
Schwefelsäure 0.166 mol/l (0.333N) 739
Schwefelsäure 0.25 mol/l (0.5N) 739
Schwefelsäure 0.26 mol/l (0.52N) 739
Schwefelsäure 0.33 mol/l (2N/3) 738
Schwefelsäure 0.5 mol/l (1N) 738
Schwefelsäure 1 mol/l (2N) 737
Schwefelsäure 2.5 mol/l (5N) 737
Schwefelsäure 4 mol/l (8N) 737
Schwefelsäure 10g/l Silbersulfat 742
Schwefelsäure 6.6 g/l Silbersulfat 743
Schwefelsäure d=1.820 743
Schwefel-Standardlösung 730
Sebacinsäure 646
Seife in Ethanollösung 660
Selen 646
Selen-Standardlösung 647
Selendioxid 648
Selenige Säure 648
Selenmischung 646
Silber, Blatt 652
Silber, Wolle 653
Silberstandardlösung 653
Silberacetat 654
Silbercarbonat 654
Silberchlorid 654
Silberdiethyldithiocarbaminat 655
Silber-Mangan-Papier 655
Silbernitrat 655
Silbernitrat 0.1 mol/l (0.1N) in 2-Propanol 657
Silbernitrat 0.01 mol/l (N/10) in Propanol-2 658
Silbernitrat 0.01 mol/l (0.01N) 658
Silbernitrat 0.05 mol/l (0.05N) 657
Silbernitrat 0.1 mol/l (0.1N) 657
Silbernitrat 0.5 mol/l (0.5N) 656
Silbernitrat 1 mol/l (1N) 656

ALPHABETISCHER INDEX

Silbernitrat 2.90756% 656
Silbernitrat 5% 656
Silbernitratlösung 658
Silberoxid 659
Silbersulfat 659
Silbersulfat 0.7% in Schwefelsäure 659
Silizium-Standardlösung 651
Siliciumdioxid 652
Sodium carbonate decahydrate 668
Sojalecithin 461
Sonnenblumenpapier 466
Sorbit 714
Sorbit D 715
Sorbitollösung 70% 715
Standardgemisch für die Kohlenwasserstoffanalyse 546
Standardlösungen für den Einfärbungsgrad von Flüssigkeiten 291
Standardtrübung 4000NTU Formazin 780
Stärke 721
Stärkelösung 1% 721
Stearinsäure 722
Strontium-Standardlösung 722
Strontium acetate 723
Strontiumbromidmonohydrat 723
Strontiumcarbonat 724
Strontiumchlorid-Hexahydrat 724
Strontiumnitrat 724
Strontiumsulfat 725
Sudan III 726
Sudan III hydroalkoholische Lösung 727
Sudangelb 727
Sulfaminsäure 727
Sulfanilsäure 728
Sulfanilsäureamid 728
Sulfat-Standardlösung 728
Sulfit-Standardlösung 729
Sulfolano 729
Sulfomolybdat-Reagens 729
Sulfosalicylsäure 729
Suspensionsmutter der Opaleszenz 621

T

Talkum 744
TAN-Reagenz 635
Tantal-Standardlösung 744
Tartrazin 746
Tauber-Reagenz 746
TBN-Reagenz D2896 636
TBN-Reagenz D4739 636
Tellur in Stücken 746
Tellurpulver 747
Tellur-Standardlösung 747
Terbium-Standardlösung 748
Tetrabutylammonium bisulfate 748
Tetrabutylammoniumbromid 749
Tetrabutylammoniumhydroxid 0.1 mol/l (0.1 N) in Methanol / Propanol-2 (50/50) 750
Tetrabutylammoniumhydroxid 0.1 mol/l (0.1N) 749
Tetrabutylammoniumhydroxid 0.1 mol/l (0.1N) in 2-Propanol 749
Tetrachlorethan-d2 751
Tetrachlorethylen 750
Tetraethylammoniumbromid 751
Tetrahydrofuran 751
Tetrahydrofuran-d8 753
Tetramethylammoniumhydroxid 10% 753
N,N,N',N'-Tetramethyl-p-phenylendiamin-Dihydrochlorid 753
Thallium-Standardlösung 754
Thioacetamid 754
Thioacetamidlösung 40 g/l 755
2-Thiobarbitursäure 755
Thioglykolsäure 80% 755
Thioharnstoff 756
Thorium-Standardlösung 756
Thullium-Standardlösung 756
Thymol 757
Thymolblau 757
Thymolblau lösung 758
Thymolblau lösung 0.4% in Ethanol 757
Thymol blauer Indikator 758
Thymolblau TA-Indikator 758
Thymolphthalein 758
Thymolphthalein 0.1% ige hydroalkoholische Lösung 759
Thymolphthaleinlösung 0.1% in Ethanol 759
TISAB 762
Tisab III-Lösung 763
Titan-Standardlösung 763
Titandioxid 764
Titangelb 270
Titanisopropylat 764
Titantrichlorid-Schwefelsäure 765
o-Tolidin in Lösung 765
o-Tolidinlösung 0.1% 765
Toluidinblau 769
Toluol 765
Toluol in Lösung in Hexan 768
Toluol-d8 768
p-Toluolsulfonsäure-Natriumsalz 769
Toluol-4-sulfonsäure 768
p-Toluolsulfonamid 768
Triacetin 769
Tribasisches Kaliumcitrat-Monohydrat 591
Trichloressigsaure 770
Trichloressigsäure 20% 770
1,2,4-Trichlorbenzol 771
2,2,2 - Trichlorethanol 770
Triethanolamin 771
Triethylamin 772
Triethylenglycol 773
Trifluoressigsäure 773
Trifluoressigsäure-d1 774
Trifluoressigsäureanhydrid 774
Trikresol 771
N,O-Bis (trimethylsilyl) acetamid 204
Trimethylcetylammoniumbromid 775
Trimethylcethylammonium p-toluenesulphonate 775
N,O-Bis-(trimethylsilyl) trifluoracetamid 204
Trioctylphosphinoxid 775
Triphenylphosphin 776
2,3,5-Triphenyltetrazolium-chlorid 776
Tris(hydroxymethyl)-aminomethan 776
Tris (hydroxymethyl) aminomethanlösung 777
Tris (hydroxymethyl) aminomethan-EDTA-Pufferlösung pH 8.4 777
Tris(hydroxymethyl)-aminomethan HCl 778
Tris (hydroxymethyl) aminomethan-Natriumchlorid-Pufferlösung pH 7.4 778
Tris (hydroxymethyl) aminomethan-Pufferlösung pH 8.1 777
Triton © X100-Lösung 778
Tropaeolin 0 778
Tropaeolin 0 0.1 % 779
Türkisches Reagenz 780
L-Tyrosin 781

U

Universalanzeige pH-Wert 0-5 hydroalkoholische Lösung 421
Universalindikator pH 1-11 wässrige Lösung 422
Unterscheidungsmerkmal für Gram-Hucker Kit 311
Uran-Standardlösung 782
Urea 782

V

Vanadium-Standardlösung 783
Vanadium (V) oxid 784
Vanillin 784
Vanillinlösung, Phosphorsäure 784
Verdünnte Schwefelsäure 743
Viktoriablau 785
Viskositätsstandards 785

W

Wasser 786
Wasser + 0.1% (v / v) Ameisensäure 789
Wasser + 0.1% v/v trifluoressigsäure 789
Wasserfreies Koffein 234
Wasserstoffperoxyd 3% 416
Wasserstoffperoxyd 3.5% w/v 415
Wasserstoffperoxyd 6 % 415
Wasserstoffperoxyd 30 % 414
Wasserstoffperoxyd 35% 413
Wasserstoffperoxyd 30-32% 413
Wasserstoffperoxyd 40% m/v 413
Weiße Vaseline 556
Weißes Aluminiumoxid 151
L(+)-Weinsäure 745
L (+) Weinsäure 20% 746
Wijs 'Reagenz 790
Wismut, Granulat 202
Wismut-Standardlösung 202
Wismut (III) basisches Carbonat 203
Wismut (III) nitrat basisch 203
Wismut (III) nitratpentahydrat 204
Wolfram-Standardlösung 779
Wolfram (VI) oxid 780
Wolframsäure 780
Wright's Färbelösung in Methanol 790

X

Xylenolorange 793
Xylol, isomerengemisch 791
Xylolcyanol 792
2,4-Xylenol 792
o-Xylol 792
D(+)-Xylose 793

Y

Ytterbium-Standardlösung 794
Yttrium-Standardlösung 794

Z

Zellulose, Pulver 250
Ziehl-Neelsens Reagenz 796
Zink aktiviert 797
Zink, Folie 796
Zink, körnig 796
Zink, Pulver 797
Zink-Standardlösung 797
Zinkacetatdihydrat 798
Zinkcarbonat basisch 799
Zinkchlorid - Jodlösung 800
Zinkchlorid wasserfrei 799
Zinkchlorid-Ameisensäure - Lösung 800
Zinkchloridlösung 60% 799
n,n-Zinkdibenzyldithiocarbamat 800
Zinknitrathexahydrat 800
Zinkoxid 801
Zinkstearat 801
Zink-Sulfat heptahydrat 802
Zinksulfat-Monohydrat 802
Zink-Sulfat 0.05 mol/l 803
Zink-Sulfat 0.1 mol/l (0.2N) 802
Zinksulfid 803
Zinn geht 759
Zinnpulver 759

DASITGROUP | CARLO ERBA REAGENTS

Product specifications are subject to changes. Please visit our website for updates.

Zinn-Standardlösung	760
Zinn (II) chloriddihydrat	761
Zinn (II) chloridlösung	761
Zinn (II) -chloridlösung 10%	761
Zinn (II) sulfat	762
Zinn (IV) chloridpentahydrat	762
Zinn (IV) oxid	762
Zirconyl nitrat	805
Zirkoniumpulver	803
Zirkonium-Standardlösung	804
Zirkonium (IV) oxid	804
Zitronensäure monohydrat	269
Zitronensäure wasserfrei	269
Zweiwertiges Kaliumphosphat-Trihydrat	614

NOTES

www.carloerbareagents.com

GENERAL INDEX

I / INTRODUCTION

CONTACT US 1
PRODUCTION 2
CUSTOM SOLUTIONS 3
PRODUCTS CLASSIFICATION 4
Xcipharm™ 6
QUALITY 8
CERTIFICATES OF ANALYSIS 10
LABELS 11
SAFETY - DATA SHEET 12
SAFETY - CLP CLASSIFICATION 13
HOW TO READ THE CATALOGUE 14

II / PRODUCT APPLICATIONS

LIQUID CHROMATOGRAPHY 19
GAS CHROMATOGRAPHY 25
ION PAIR CHROMATOGRAPHY 28
ION CHROMATOGRAPHY 29
TRACE ANALYSIS: METALS 31
TRACE ANALYSIS: ORGANICS 49
CHEMICAL OXYGEN DEMAND ANALYSIS 57
PHARMACEUTICAL PRODUCTION 58
PHARMACEUTICAL QUALITY CONTROL 69
PHYSICALS AND CHEMICAL CHARACTERISTICS 82
pHMETRY 86
VOLUMETRY 91
FOOD ANALYSIS 96
ELECTRONICS 100
HISTOLOGY, HEMATOLOGY AND CYTODIAGNOSTIC 103
ORGANIC CHEMISTRY 110
NMR SPECTROSCOPY 113
UV SPECTROSCOPY 115
PEPTIDE SYNTHESIS 117
GREEN CHEMISTRY 118
PETROLCHEMICAL ANALYSIS 119
KARL FISCHER TITRATION 120
INDICATORS 122
CONDUCTIMETRY 125

III / PRODUCTS 127

IV / CHEMICALS TABLES 807

V / INDEXES

INDEX BY PRODUCT CODE 830
INDEX BY CAS NUMBER 861
INDICE ALFABETICO 869
INDEX ALPHABÉTIQUE 879
ÍNDICE ALFABÉTICO 889
ALPHABETISCHER INDEX 899

DASITGROUP | CARLO ERBA REAGENTS

“partner in your choice”

Over 6000 products
featured for ***your analysis***
and synthesis needs
developed for a
sustainable future.

DASITGROUP | CARLO ERBA REAGENTS

ITALIA
CARLO ERBA Reagents S.r.l.
Via Raffaele Merendi 22
20010 Cornaredo (MI)
Servizio Clienti
servizioclienticer@dgroup.it
Tel.: +39 02 93 99 190
Fax: +39 02 93 99 10 01

FRANCE
CARLO ERBA Reagents SAS
Chaussée du Vexin,
Parc d'affaire des Portes
27106 Val de Reuil
Service Client
serviceclient@cer.dgroup.it
Tél.: +33 2 32 09 20 00
Fax: +33 2 32 59 11 89

DEUTSCHLAND
CARLO ERBA Reagents GmbH
Denzlinger Str. 27
79312 Emmendingen
Kundendienst
info.de@cer.dgroup.it
Tel.: +49 07641 46 881 90
Fax: +49 07641 46 881 919

ESPAÑA
CARLO ERBA Reagents S.A.
Calle Filadors 35, 6ª
Planta Puerta 5
08208 Sabadell (BCN)
Servicio Cliente
servicliente@cer.dgroup.it
Tel.: +34 93 693 37 35
Fax: +34 93 724 31 68

ALL OTHER COUNTRIES
Customer Service
export@cer.dgroup.it
Ph.: +33 2 32 09 20 00
Fax: +33 2 32 59 11 89